I, Paolo Bosio, Technical Director, Sentinel CH SpA, Via Robert Koch 2, Milan 20152, Italy, hereby certify that the attached documents that keep the number of pages indicated below are the Instruction for Use. These documents are in a full accordance with actual manuals of manufacturer.

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Signed: _____

Date: February 13th 2020

Dr. Paolo Bosio,
Technical Director
Sentinel CH SpA, Via Robert Koch 2,
Milan 20152, Italy

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Read Highlighted Changes: Revised May 2016.

INTENDED USE

The Ammonia Controls must only be used to verify the performance of the Ammonia Ultra assay **REF** 6K89-30.

CONTENTS / MATERIALS PROVIDED

REF 6K89-20 Ammonia Controls are prepared with ammonium sulfate in **human serum**. Preservative is also present.

Actual ammonia concentrations for this lot of Ammonia Controls are listed in the Ammonia Controls Value Sheet packaged with the controls. The color of the cap corresponding to each level of control is indicated in the table below.

CONT	Cap Color	QTY
CONTROL 1	White	1 x 5 mL
CONTROL 2	Blue	1 x 5 mL
CONTROL 3	Red	1 x 5 mL

MATERIALS REQUIRED BUT NOT PROVIDED

REF 6K89-30 Ammonia Ultra

PRECAUTIONS

Precautions for Users

- **IVD**
- For *In Vitro* Diagnostic Use.
- Do not use components beyond the expiration date.
- Do not mix materials from different kit lot numbers.

CAUTION: This product contains human sourced and/or potentially infectious components. Refer to the **CONTENTS / MATERIALS PROVIDED** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens.¹ Biosafety Level 2² or other appropriate biosafety practices^{3,4} should be used for materials that contain or are suspected of containing infectious agents.

- The human serum used in **CONTROL 1-3** is nonreactive for HBsAg, HIV-1 Ag or HIV-1 NAT, anti-HIV-1/HIV-2, and anti-HCV.

The following warnings and precautions apply to **CONTROL 1-3**:

Contains sodium azide.	
EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

- Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

FOR USE WITH

ARCHITECT

Ammonia Controls

REF 6K89-20

306958/R06

C6K890

STORAGE AND STABILITY

- Unopened controls are stable until the expiration date when stored at 2 to 8°C.
- Opened controls are stable for 1 month at 2 to 8°C, if contamination is avoided and the bottles are recapped immediately after use.

PREPARATION OF CONTROL

Ammonia Controls require no preparation prior to use.

INSTRUCTIONS FOR USE

1. Verify that the correct control values have been entered into the control file.
2. Remove controls from the refrigerator and allow to stand for 15 minutes at room temperature.
3. Mix bottles several times by gentle inversion to ensure the homogeneity of the solution. Avoid foaming.
4. Open the bottles, place an appropriate amount of each control in a separate sample cup, and place in the assigned positions.

NOTE: Immediately cap bottles tightly and return to refrigerated storage. Always return each cap to its original bottle.

5. Run controls according to the instructions in Section 5 of the **ARCHITECT System Operations Manual**.

INDICATIONS OF INSTABILITY OR DETERIORATION

Controls are clear, colorless liquid. Instability or deterioration should be suspected if there are visible signs of leakage, extreme turbidity, microbial growth, or if controls do not meet the appropriate acceptance criteria.

LIMITATIONS OF THE PROCEDURE

Accurate and reproducible results are dependent upon properly functioning instruments and reagents, storage of product as directed, and good laboratory technique.

BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. *Bloodborne Pathogens*.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections: Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

TRADEMARKS

The ARCHITECT c System family of instruments consists of c 4000, c 8000, and c 16000 instruments.

ARCHITECT, c 4000, c 8000, c 16000, c System, and MULTIGENT are trademarks of Abbott Laboratories in various jurisdictions.

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Key to Symbols	
	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
CONT	Contents of kit
CONTROL 1-3	Controls 1 through 3
CONTROL 1	Control 1
CONTROL 2	Control 2
CONTROL 3	Control 3
CONVENTIONAL UNITS	Conventional Units
CONTAINS: AZIDE	Contains sodium azide. Contact with acids liberates very toxic gas.
DISTRIBUTED IN THE USA BY	Distributed in the USA by
FOR USE WITH	Identifies products to be used together
GTIN	Global Trade Item Number
INFORMATION FOR USA ONLY	Information needed for United States of America only
IVD	In Vitro Diagnostic Medical Device
LOT	Lot Number
PRODUCT OF USA	Product of USA
QTY	Quantity
RANGE	Range
REF	List Number
SI UNITS	SI Units
VALUE	Value

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com.

SENTINEL CH. SpA
Via Robert Koch, 2
Milan 20152 Italy

DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA

May 2016

ADDENDUM

to the Instruction for Use Ammonia Controls (for Russia)

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WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. The maximum shelf life of the medical device is 24 months from date of manufacturing. Device must be discarded if stored with deviations in manufacturer's requirements.

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- US Department of Labor, Occupational Safety and Health Administration. 29 CFR Part 1910.1030. *Bloodborne Pathogens*.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office, December 2009.
- World Health Organization. *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization, 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

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www.roszdravnadzor.gov.ru

----- CERTIFIED SIGNATORY -----
I, the undersigned Valerio Tacchini, Notary Public in Milan,
Italy, -----

----- CERTIFY THAT: -----
- BOSIO PAOLO, born in Alba on April 13th 1963, domiciled
for his office in Milan, Robert Koch n. 2, acting as attor-
ney and legal representative of "SENTINEL CH. S.p.A." with
register office in Milan, Via Robert Koch n. 2, registered
at the Trade Register of Milan - fiscal code 07118040158, du-
ly empowered to do so, -----
of whom identity and powers I, as Notary, am certain, he si-
gned, before me, the above document. -----
Milan, Largo Donegani, 2, february 13th 2020 -----

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FOR USE WITH
ARCHITECT

REF 4R29-21

Read Highlighted Changes: Revised March 2018.

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

NAME

UIBC

INTENDED USE

The UIBC assay is used for the quantitative determination of unsaturated iron-binding capacity (UIBC) in serum or plasma on the ARCHITECT cSystems.

SUMMARY AND EXPLANATION OF THE TEST

Iron status is used in the diagnosis and management of some hematologic and hepatic conditions. In addition to serum iron levels, several tests may be ordered including UIBC, total iron binding capacity (TIBC) and ferritin. TIBC is a surrogate for transferrin and has a direct relationship to it. $TIBC = (UIBC + \text{serum iron})$.¹⁻³ UIBC and TIBC are increased in low iron states such as uncomplicated anemia and decreased in high iron conditions such as hemochromatosis. The lone exception to the preceding is the case of anemia of chronic disease where the patient may be anemic but has adequate iron reserves and a low UIBC.¹

PRINCIPLES OF THE PROCEDURE

Sample is added to an alkaline buffer containing a known concentration of iron to saturate the available binding sites on transferrin. The iron that remains free after transferrin saturation is reduced to a ferrous state and then complexed by Ferene-S* to form a stable complex, of which the color intensity is measured at 604 nm. The color intensity is directly proportional to the unbound excess iron concentration and indirectly proportional to the unsaturated iron-binding capacity. UIBC is therefore determined by subtracting the quantity of unbound iron from the total added quantity.

*Ferene-S = 3-(2-pyridyl)-5,6-bis-[2-(5-furylsulfonic acid)]-1,2,4-triazine

Methodology: Ferene

REAGENTS

Kit Contents

UIBC 4R29

Supplied as a liquid, ready-to use three-component reagent kit.

REF	4R29-21
	850*
R1	5 x 36 mL
R2	5 x 15 mL
CAL	3 x 2 mL

*Calculation is based on the minimum reagent fill volume per kit.

Reactive Ingredients	Concentration
R1	
Glycyl glycine buffer (pH 8.6)	255 mmol/L
Sodium bicarbonate	46 mmol/L
Iron (II) chloride	20 µmol/L
R2	
Ascorbic acid (pH 2.5)	300 mmol/L
Ferene-S	12.5 mmol/L
CAL	
Iron solution	500 µg/dL (89.5 µmol/L)

R1 and R2 contain stabilizers, preservatives, and detergents. CAL contains stabilizer.

Warnings and Precautions

- IVD
- For *In Vitro* Diagnostic Use
- Rx ONLY

Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.⁴⁻⁷

The following warnings and precautions apply to: R1 and R2	
WARNING	Contains methylisothiazolone and polyethylene glycol octylphenyl ether*.
H317	May cause an allergic skin reaction.
H402*	Harmful to aquatic life.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
P273*	Avoid release to the environment.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

* Not applicable where regulation EU 1272/2008 (CLP) has been implemented.

0.0555% of the ingredients in R2 are of unknown toxicity.

The following warnings and precautions apply to the CAL	
WARNING	Contains hydroxylamine hydrochloride.
H351	Suspected of causing cancer.
H317	May cause an allergic skin reaction.
H402*	Harmful to aquatic life.
Prevention	
P201	Obtain special instructions before use.
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P273*	Avoid release to the environment.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
P308+P313	IF exposed or concerned: Get medical advice / attention.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

* Not applicable where regulation EU 1272/2008 (CLP) has been implemented.

Safety Data Sheets are available at www.abbottiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

Reagent Handling

- Do not use reagents beyond the expiration date.
- Do not pool reagents within a kit or between kits.
- Do not use components from one lot with components from another lot.
- When either the **R1** or **R2** reagent cartridge becomes empty, replace both cartridges and validate the system by analyzing controls.
- Remove any air bubbles present in the reagents with a new applicator stick, or allow the reagents to sit at the appropriate storage temperature to allow the bubbles to dissipate. To minimize volume depletion, do not use a transfer pipette to remove bubbles.

CAUTION: Bubbles may interfere with proper detection of reagent level in the cartridge and cause insufficient reagent aspiration which could impact results.

For a detailed discussion of handling precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 7.

Reagent Storage

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened reagent and Calibrator	2-8°C	Until expiration date	May be used immediately after removal from 2-8°C storage.
On board	System temperature	28 days	After 28 days (672 hours), the reagent kit must be discarded.
Opened Calibrator	2-8°C	60 days	Recap vials immediately after use to avoid contamination.

Indications of Reagent Deterioration

Reagents and calibrator must be clear.

Instability or deterioration should be suspected if there are precipitates, visible signs of leakage or contamination, turbidity, or if calibration or controls do not meet the appropriate package insert and/or ARCHITECT System Operations Manual criteria.

INSTRUMENT PROCEDURE

The UIBC assay file must be installed on the ARCHITECT cSystem prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.

For information on printing assay parameters or for a detailed description of system procedures, refer to the ARCHITECT System Operations Manual, Section 5.

Alternate Result Units

The Conventional result unit for the UIBC assay is $\mu\text{g/dL}$. The corresponding SI result unit is $\mu\text{mol/L}$. To convert results from $\mu\text{g/dL}$ to $\mu\text{mol/L}$, multiply $\mu\text{g/dL}$ by 0.179.⁹

When converting to units other than those listed, refer to the ARCHITECT System Operations Manual, Section 2.

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Specimen Types

Verified specimen types to be used with this assay:

Specimen Types	Collection Tubes
Serum	Serum (with or without gel barrier)
Plasma	Lithium heparin (with or without gel barrier) Sodium heparin

- Ethylenediaminetetraacetic acid (EDTA) and sodium fluoride/oxalate are not acceptable anticoagulants.⁹
- Other specimen types and collection tube types have not been verified with this assay.
- The instrument does not provide the capability to verify specimen type. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

Specimen Conditions

Liquid anticoagulants may have a dilution effect resulting in lower concentrations.

- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- For accurate results, plasma specimens should be free of platelets and other particulate matter. Ensure centrifugation is adequate to remove platelets.

- For additional information on specimen conditions, refer to the Interference section of this package insert.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

Preparation for Analysis

Serum: Ensure complete clot formation has taken place prior to centrifugation. Centrifuge according to tube manufacturer's specifications to ensure proper separation of serum from blood cells.

Plasma: Centrifuge according to tube manufacturer's specifications to remove platelets and ensure proper separation of plasma from blood cells.

Frozen specimens must be completely thawed before mixing.

Mix thawed specimens thoroughly.

Visually inspect thawed specimens. If layering or stratification is observed, mix until specimens are visibly homogeneous. If specimens are not mixed thoroughly, inconsistent results may be obtained.

NOTE: Stored specimens must be inspected for particulates. If present, mix and centrifuge the specimen to remove particulates prior to testing.

For total sample volume requirements, refer to the ARCHITECT System Operations Manual, Section 5.

Specimen Storage

Specimen Type	Storage Temperature	Maximum Storage Time
Serum/Plasma	20-25°C	7 days ¹⁰
	2-8 °C	3 weeks ^{10, 11}
	-20°C	1 year ¹⁰

Each laboratory may establish a range around -20°C from either the freezer manufacturer's specifications or your laboratory standard operating procedure(s) for specimen storage.

Avoid multiple freeze/thaw cycles.

PROCEDURE

Materials Provided

4R29 UIBC Reagent Kit

Materials Required but not Provided

- Control material
- Saline (0.85% to 0.90% NaCl) for specimen dilution

Assay Procedure

For a detailed description of how to run an assay, refer to the ARCHITECT System Operations Manual, Section 5.

For sample volume requirements, refer to the ASSAY PARAMETERS section of this package insert. Minimum sample volume is calculated by the system and printed on the Order List Report. Ensure adequate sample volume is present prior to running the test.

Specimen Dilution Procedures

Specimens with UIBC values exceeding 500 µg/dL (89.5 µmol/L) are flagged and may be diluted by following the Manual Dilution Procedure.

Manual Dilution Procedure

1. Dilute the specimen with saline (0.85% to 0.90% NaCl).
2. Enter the dilution factor in the Patient or Control order screen. The system uses this dilution factor to automatically correct the concentration by multiplying the result by the entered factor. If the operator does not enter the dilution factor, the result must be manually multiplied by the appropriate dilution factor before reporting the result.

$$\text{Manual Dilution Factor} = \frac{(\text{Volume of Specimen} + \text{Volume of Dilution Reagent})}{\text{Volume of Specimen}}$$

If a diluted specimen result is flagged indicating it is less than the linear low limit, do not report the result. Rerun using an appropriate dilution.

For detailed information on ordering dilutions, refer to the ARCHITECT System Operations Manual, Section 5.

Calibration

NOTE: The calibrator value for the UIBC Calibrator **must be entered as a negative number**. Refer to the ASSAY PARAMETERS section of this package insert.

Calibration is stable for approximately 168 hours, but **is required with each change in reagent cartridge**. Verify calibration with at least two levels of controls according to the established quality control requirements for your laboratory. If control results fall outside acceptable ranges, recalibration may be necessary.

The calibrator value is verified using a National Institute of Standards and Technology (NIST), SRM 3126.

For detailed information on how to perform an assay calibration, refer to the ARCHITECT System Operations Manual, Section 6.

Quality Control Procedures

As appropriate, refer to your laboratory standard operating procedure(s) and/or quality assurance plan for additional quality control requirements and potential corrective actions.

- At least two levels of controls [normal and abnormal] are to be run every 24 hours.
- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, patient values may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

Specimen Shipping

Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.

RESULTS

Representative performance data are given in the EXPECTED VALUES and SPECIFIC PERFORMANCE CHARACTERISTICS sections of this package insert. Results obtained in individual laboratories may be different from the data presented.

LIMITATIONS OF THE PROCEDURE

Use disposable plastic pipettes or glassware soaked in 1 N hydrochloric acid solution and rinsed with distilled water.

For diagnostic purposes, the test findings should always be assessed in conjunction with the patient's medical history, clinical examinations, and other findings.

Significant interference may be observed with hemolyzed samples. See Interference section for additional information.

Refer to the SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS and SPECIFIC PERFORMANCE CHARACTERISTICS sections of this package insert.

EXPECTED VALUES

Reference Range

Serum and Plasma

	Range (µg/dL)	Range (µmol/L)
Male	69 to 240	12.4 to 43.0
Female	70 to 310	12.5 to 55.5

A confirmation study was conducted using 22 male and 22 female volunteers from a healthy donor population in Milan, Italy. Data were analyzed as described by Clinical and Laboratory Standards Institute (CLSI) protocol EP28-A3c.¹² From this study, the central 90% of specimens fell within the ranges above.

The circulated iron level can show a 30% diurnal variation, with a peak early in the morning.¹³

It is recommended that each laboratory determine its own reference range based upon its particular locale and population characteristics.

SPECIFIC PERFORMANCE CHARACTERISTICS

Representative performance data are given in this section. Results obtained in individual laboratories may vary.

Linearity

Linearity for the UIBC assay is 25 to 500 µg/dL (4.5 to 89.5 µmol/L). Linearity was verified using CLSI protocol EP06-A.¹⁴

Limit of Quantitation (LOQ)

The LOQ for UIBC is ≤ 25 µg/dL (4.5 µmol/L). The LOQ is the analyte concentration at which the %CV = 20%. Studies were conducted on the ARCHITECT System using CLSI Document EP17-A2.¹⁵

Precision

The precision of the UIBC assay is ≤ 7% Total CV at concentrations ≥ 143 µg/dL (25.6 µmol/L) and ≤ 10 µg/dL (1.8 µmol/L) SD at concentration < 143 µg/dL (25.6 µmol/L). The study was conducted on the ARCHITECT System using CLSI Document EP05-A3¹⁶ and summarized below.

Control		Level 1	Level 2	Level 3
N		132	132	132
Mean (µg/dL)		152	219	271
Mean (µmol/L)		27.1	39.1	48.5
Within Day	SD (µg/dL)	3.2	2.7	2.3
	SD (µmol/L)	0.6	0.5	0.4
	%CV	2.1	1.2	0.9
Between Day	SD (µg/dL)	5.4	7.7	10.1
	SD (µmol/L)	1.0	1.4	1.8
	%CV	3.6	3.5	3.7
Total	SD (µg/dL)	6.3	8.2	10.4
	SD (µmol/L)	1.1	1.5	1.9
	%CV	4.2	3.7	3.8

Interference

Interference studies were conducted on the ARCHITECT System using CLSI document EP07-A2.¹⁷ Representative data are summarized below.

Interference studies were conducted using an acceptance criteria of ± 10% or ± 14 µg/dL (2.5 µmol/L) deviation, whichever is greater, from the target. UIBC is not affected by the presence of the following interferents up to the concentrations indicated below.

Interfering Substance	Interferent Concentration	N	UIBC-Target (µg/dL)	Difference from Target (µg/dL)	%Difference from Target
Bilirubin (conjugated)	59 mg/dL	12	141	-10.6	—
	62 mg/dL	12	298	—	-4.0
Bilirubin (unconjugated)	53 mg/dL	12	156	—	-1.8
	56 mg/dL	12	284	—	-2.5
Copper	1000 µg/dL	3	139	-6.0	—
	2000 µg/dL	3	289	—	-7.8
Hemoglobin	62 mg/dL	12	168	—	-8.1
	125 mg/dL	12	323	—	-9.4
Intralipid	1000 mg/dL	12	142	-9.7	—
	1000 mg/dL	12	295	—	-1.6
Oxytetracycline	100 mg/L	3	136	13.2	—
	100 mg/L	3	303	—	4.4
Rheumatoid Factor	100 IU/mL	3	145	—	2.4
	100 IU/mL	3	300	—	-0.2
Triglycerides*	901 mg/dL	12	147	—	-2.8
	1157 mg/dL	12	269	—	-6.7
Total Protein	13.2 g/dL	12	152	—	5.2
	13.4 g/dL	12	252	—	10.0

Interfering Substance	Interferent Concentration	N	UIBC-Target (µmol/L)	Difference from Target (µmol/L)	%Difference from Target
Bilirubin (conjugated)	1008.9 µmol/L	12	25.2	-1.9	—
	1060.2 µmol/L	12	53.3	—	-4.0
Bilirubin (unconjugated)	906.3 µmol/L	12	27.9	—	-1.8
	957.6 µmol/L	12	50.8	—	-2.5
Copper	157.0 µmol/L	3	24.9	-1.1	—
	314.0 µmol/L	3	51.7	—	-7.8
Hemoglobin	0.62 g/L	12	30.1	—	-8.1
	1.25 g/L	12	57.8	—	-9.4
Intralipid	10 g/L	12	25.4	-1.7	—
	10 g/L	12	52.8	—	-1.6
Oxytetracycline	220 µmol/L	3	24.3	2.4	—
	220 µmol/L	3	54.2	—	4.4
Rheumatoid Factor	100 IU/mL	3	26.0	—	2.4
	100 IU/mL	3	53.7	—	-0.2
Triglycerides*	10.2 mmol/L	12	26.3	—	-2.8
	13.1 mmol/L	12	48.2	—	-6.7
Total Protein	132 g/L	12	27.2	—	5.2
	134 g/L	12	45.1	—	10.0

*Triglycerides testing was performed using commercially available, human-sourced lipemic material. The results displayed here may not adequately represent results observed on other instruments with similarly-sourced, but different lipemic material.^{18, 19}

Clinical samples with high triglycerides can be clear while samples with moderate triglycerides can be turbid. The light scatter for Intralipid does not always match with that of lipoproteins in the patient samples due to differences in the size and distribution of these lipoproteins.²⁰

The effect of interfering substances has only been evaluated for those listed above.

The following compounds may interfere with the UIBC assay:

- Gadolinium Magnetic Resonance Contrast Agents^{21, 22}
- Deferasirox²³

Interferences from medications or endogenous substances may affect results.

Method Comparison

Results from the UIBC assay on an ARCHITECT cSystem were compared with the results from a commercially available ferrozine methodology, the UIBC 4P79, and between ARCHITECT cSystem instruments. Studies were conducted on the ARCHITECT System using CLSI Document EP09-A3²⁴ and summarized below.

	4R29 vs. Comparative Method	4R29 vs. 4P79	Between ARCHITECT cSystems
N	110	109	110
Y-Intercept (µg/dL)	-4.97	8.03	-4.98
Y-Intercept (µmol/L)	-0.89	1.44	-0.89
Correlation Coefficient	0.995	0.995	0.996
Slope	0.97	0.97	0.96
Range (µg/dL)*	41.7 to 478.1	31.5 to 458.8	39.9 to 447.6
Range (µmol/L)*	7.5 to 85.6	5.6 to 82.1	7.1 to 80.1

*Indicated ranges are from the comparative method, UIBC 4P79, and the c8000, respectively.

BIBLIOGRAPHY

1. Jacobs DS, DeMott WR, Oxley DK. *Laboratory Test Handbook*, 5th ed. Hudson, OH: Lexi-Comp; 2001:203–204.

2. Murtagh LJ, Whiley M, et al. Unsaturated iron binding capacity and transferrin saturation are equally reliable in detection of HFE hemochromatosis. *Am J Gastroenterology* 2002;97:2093-2099.

3. Gambino R, Desvarieux E, Orth M, et al. The relationship between chemically measured total iron-binding capacity concentrations and immunologically measured transferrin concentrations in human serum. *Clin Chem* 1997;43(12):2408-2412.

4. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.

5. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.

6. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.

7. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

8. McPherson RA, Pincus MR. *Henry's Clinical Diagnosis and Management by Laboratory Methods*, 21st ed. Philadelphia, PA: WB Saunders; 2007:1407.

9. Siek G, Lawlor J, Pelczar D, et al. Direct serum total iron-binding capacity assay suitable for automated analyzers. *Clin Chem* 2002;48:1:161-166.

10. Guder WG, da Fonseca-Wollheim F, Heil W, et al. *The Quality of Diagnostic Samples*. Darmstadt, Germany: GIT Verlag; 2009:53.

11. US Pharmacopeial Convention, Inc. General notices. In: *US Pharmacopeia National Formulary*. 1995 ed (USP 23/NF18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.

12. Clinical and Laboratory Standards Institute (CLSI). *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition*. CLSI document EP28-A3c. Wayne, PA: CLSI; 2010.

13. Ricós C, Alvarez V, Cava F, García-Lario JV, et al. Current databases on biological variation: pros, cons and progress. *Scand J Clin Lab Invest*. 1999;59:491-500.

14. CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.

15. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.

16. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision of Quantitative Measurement Procedures: Approved Guideline - Third Edition*. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.

17. CLSI. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI document EP07-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.

18. Kroll MH. Evaluating interference caused by lipemia. *Clin Chem* 2004;50:1968-1969.

19. Bornhorst JA, Roberts RF, Roberts WL. Assay-specific differences in lipemic interference in native and intralipid-supplemented samples. *Clin Chem* 2004;50:2197-2201.

20. Nikolac N. Lipemia: causes, interference mechanisms, detection and management. *Biochemia Medica*. 2014;24(1):57-67.








21. Lippi G, Daves M, Mattiuzzi C. Interference of medical contrast media on laboratory testing. *Biochem Med* 2014;24:80-88.

22. Proctor K, Rao L, Roberts W. Gadolinium magnetic resonance contrast agents produce analytic interference in multiple serum assays. *Am J Clin Pathol* 2004;121:282-292.

23. Ikuta K, Ito S, Tanaka H, et al. Interference of deferasirox with assays for serum iron and serum unsaturated iron binding capacity during iron chelating therapy. *Clin Chim Acta* 2011;412:2261–2266.

24. Clinical and Laboratory Standards Institute (CLSI). *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

Key to Symbols

	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
	Calibrator
	Distributed in the USA by
	Identifies products to be used together
	Information needed for United States of America only
	In Vitro Diagnostic Medical Device
	Lot Number
	Product of Italy
	Reagent 1
	Reagent 2
	List Number
	For use by or on the order of a physician only (applicable to USA classification only).
	Serial number
	Warning: Suspected of causing cancer.
	Warning: May cause an allergic reaction.

ARCHITECT, c4000, c8000, c16000, cSystem, and SmartWash are trademarks of Abbott Laboratories in various jurisdictions.

The ARCHITECT cSystem family of instruments consists of c4000, c8000, and c16000 instruments.

All other trademarks are property of their respective owners.

 SENTINEL CH. S.p.A.
Via Robert Koch, 2
Milano (MI)
20152, Italy
+39 02 345 514 1

Abbott Laboratories
Abbott Park, IL 60064 USA

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Revised March 2018.

©2016, 2018 Abbott Laboratories

ASSAY PARAMETERS

UIBC Serum/Plasma—Conventional and SI Units

Configure assay parameters — General									
● General		○ Calibration		○ SmartWash		○ Results		○ Interpretation	
Assay: UIBC-L		Type: Photometric		Version: †					
Number: 2981		Run controls for onboard reagents by: Lot							
● Reaction definition		○ Reagent / Sample		○ Validity checks					
Reaction mode: End up									
Wavelength: 604 / 700		Primary Secondary		Read times Main: 31 – 33					
Last required read: 33									
Absorbance range:		Color correction: — —							
Sample blank type: Self		Blank: 13 – 15							
○ Reaction definition		● Reagent / Sample		○ Validity checks					
Reagent: UIBCL		Reagent volume: 175		R1 R2 50					
Diluent: Saline		Water volume: —							
Diluent dispense mode: Type 0		Dispense mode: Type 0		Type 0					
Dilution name	Sample	Diluted sample	Diluent	Water	Dilution factor	Default dilution			
STANDARD :	25.0	—	—	—	1:1.00	●			
—	—	—	—	—	—	○			
—	—	—	—	—	—	○			
○ Reaction definition		○ Reagent / Sample		● Validity checks					
Reaction check: None									
Maximum absorbance variation: —									

Configure assay parameters — Calibration									
○ General		● Calibration		○ SmartWash		○ Results		○ Interpretation	
Assay: UIBC-L		Calibration method: Linear							
● Calibrators		○ Volumes		○ Intervals		○ Validity checks			
Calibrator set: UIBC-L		Blank: Water		Calibrator level: 0*		Concentration: ††			
Replicates: 2 [Range 1 – 3]		Cal 1: UIBC-L1							
○ Calibrators		● Volumes		○ Intervals		○ Validity checks			
Calibrator: UIBC-L		Calibrator level		Sample		Diluted sample		Diluent Water	
Blank: Water		25.0		—		—		—	
Cal 1: UIBC-L1		25.0		—		—		—	
○ Calibrators		○ Volumes		● Intervals		○ Validity checks			
Calibration intervals:		Full interval: 168 (hours)							
Calibration type:		Adjust type: None							
○ Calibrators		○ Volumes		○ Intervals		● Validity checks			
Blank absorbance range:		Span: Blank – Blank							
Span absorbance range:		Expected cal factor: 0.00							
Expected cal factor tolerance %:		0							

Configure assay parameters — SmartWash									
○ General		○ Calibration		● SmartWash		○ Results		○ Interpretation	
Assay: UIBC-L									
COMPONENT	REAGENT / ASSAY	WASH	Volume	Replicates					
R1	AMIK9	Detergent A	345	1					
R1	DIG00	Detergent A	345	1					
R1	DGT0B	Detergent A	345	1					
R1	GENT9	Detergent A	345	1					
R1	TOBRA	Detergent A	345	1					
R1	VANCO	Detergent A	345	1					
R1	FRU0S	0.5% Acid Wash	345	1					
R1	FRUL9	Detergent A	345	1					
R2	AMIK9	Detergent A	345	1					
R2	DIG00	Detergent A	345	1					
R2	DGT0B	Detergent A	345	1					
R2	GENT9	Detergent A	345	1					
R2	TOBRA	Detergent A	345	1					
R2	VANCO	Detergent A	345	1					
R2	FRU0S	0.5% Acid Wash	345	1					
R2	FRUL9	Detergent A	345	1					

UIBC Serum/Plasma—Conventional Units

Configure assay parameters — Results									
○ General		○ Calibration		○ SmartWash		● Results		○ Interpretation	
Assay: UIBC-L		Dilution default range:		Low-Linearity: 25		High-Linearity: 500		Assay number: 2981 Result units: µg/dL	
Gender and age specific ranges:***									
GENDER	AGE (UNITS)	NORMAL				EXTREME			
Male	0 – 130 (Y)	69 – 240							
Female	0 – 130 (Y)	70 – 310							

Configure result units	
Assay: UIBC-L	
Version: †	
Result units: µg/dL	
Decimal places: 0	[Range 0 – 4]
Correlation factor: 1.0000	
Intercept: 0.0000	

UIBC Serum/Plasma—SI Units

Configure assay parameters — Results									
○ General		○ Calibration		○ SmartWash		● Results		○ Interpretation	
Assay: UIBC-L		Dilution default range:		Low-Linearity: 4.5		High-Linearity: 89.5		Assay number: 2981 Result units: µmol/L	
Gender and age specific ranges:***									
GENDER	AGE (UNITS)	NORMAL				EXTREME			
Male	0 – 130 (Y)	12.4 – 43.0							
Female	0 – 130 (Y)	12.5 – 55.5							
Configure result units									
Assay: UIBC-L									
Version: †									
Result units: µmol/L									
Decimal places: 1	[Range 0 – 4]								
Correlation factor: 1.0000									
Intercept: 0.0000									

† Due to differences in instrument systems and unit configurations, version numbers may vary.

‡ Displays the number of decimal places defined in the decimal places parameter field.

†† The value for the UIBC Liquid Calibrator must be entered as a negative number: -500 µg/dL (-89.5 µmol/L). These values are defined on the * Configure calibrator set screen.

*** Non-gender specific (Either) range may be user defined.

ADDENDUM

to the Instruction for Use UIBC (for Russia)

SHIPPING

Recommended shipping conditions for device is 2-8°C. Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. The maximum shelf life of the medical device is 14 months from date of manufacturing. Device must be discarded if stored with deviations in manufacturer's requirements.

MAINTENANCE

Not required.

REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste".

Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

98/79/EC, EN ISO 9001, ISO 13485, EN ISO 18113-1, EN ISO 18113-2

POTENTIAL USER

Professional user

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device cannot be re-applied

MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

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- Jacobs DS, DeMott WR, Oxley DK. *Laboratory Test Handbook*, 5th ed. Hudson, OH: Lexi-Comp; 2001:203–204.
- Murtagh LJ, Whitley M, et al. Unsaturated iron binding capacity and transferrin saturation are equally reliable in detection of HFE hemochromatosis. *Am J Gastroenterology* 2002;97:2093-2099.
- Gambino R, Desvarieux E, Orth M, et al. The relationship between chemically measured total iron-binding capacity concentrations and immunologically measured transferrin concentrations in human serum. *Clin Chem* 1997;43(12):2408-2412.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- McPherson RA, Pincus MR. *Henry's Clinical Diagnosis and Management by Laboratory Methods*, 21st ed. Philadelphia, PA: WB Saunders; 2007:1407.
- Siek G, Lawlor J, Pelczar D, et al. Direct serum total iron-binding capacity assay suitable for automated analyzers. *Clin Chem* 2002;48:1:161-166.
- Guder WG, da Fonseca-Wollheim F, Heff W, et al. *The Quality of Diagnostic Samples*. Darmstadt, Germany: GIT Verlag; 2009:53.
- US Pharmacopeial Convention, Inc. General notices. In: *US Pharmacopeia National Formulary*. 1995 ed (USP 23/NF18). Rockville, MD: The US Pharmacopeial Convention, Inc.; 1994:11.
- Clinical and Laboratory Standards Institute (CLSI). *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition*. CLSI document EP28-A3c. Wayne, PA: CLSI; 2010.

13. Ricós C, Alvarez V, Cava F, García-Lario JV, et al. Current databases on biological variation: pros, cons and progress. *Scand J Clin Lab Invest.* 1999;59:491-500.
14. CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline.* CLSI document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.
15. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition.* CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
16. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision of Quantitative Measurement Procedures: Approved Guideline - Third Edition.* CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.
17. CLSI. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition.* CLSI document EP07-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
18. Kroll MH. Evaluating interference caused by lipemia. *Clin Chem* 2004;50:1968-1969.
19. Bornhorst JA, Roberts RF, Roberts WL. Assay-specific differences in lipemic interference in native and intralipid-supplemented samples. *Clin Chem* 2004;50:2197-2201.
20. Nikolac N. Lipemia: causes, interference mechanisms, detection and management. *Biochemia Medica.* 2014;24(1):57-67.
21. Lippi G, Daves M, Mattiuzzi C. Interference of medical contrast media on laboratory testing. *Biochem Med* 2014;24:80-88.
22. Proctor K, Rao L, Roberts W. Gadolinium magnetic resonance contrast agents produce analytic interference in multiple serum assays. *Am J Clin Pathol* 2004;121:282-292.
23. Ikuta K, Ito S, Tanaka H, et al. Interference of deferasirox with assays for serum iron and serum unsaturated iron binding capacity during iron chelating therapy. *Clin Chim Acta* 2011;412:2261–2266.
24. Clinical and Laboratory Standards Institute (CLSI). *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition.* CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

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Configure assay parameters — General			
<input checked="" type="radio"/> General	<input type="radio"/> Calibration	<input type="radio"/> SmartWash	<input type="radio"/> Results
Assay: UIBC-L	Type: Photometric	Version: †	
Number: 2981			
Run controls for onboard reagents by: Lot			
<input checked="" type="radio"/> Reaction definition	<input type="radio"/> Reagent / Sample	<input type="radio"/> Validity checks	
Reaction mode: End up			
Primary		Secondary	Read times
Wavelength: 604 / 700		Main: 31 – 33	
Last required read: 33			
Absorbance range:		Color correction: — —	
Sample blank type: Self		Blank: 13 – 15	

O Reaction definition		● Reagent / Sample		O Validity checks	
Reagent: UIBCL		Reagent volume: 175 50		R1 R2	
Diluent: Saline		Water volume: —			
Diluent dispense mode: Type 0		Dispense mode: Type 0 Type 0			
Dilution name	Sample	Diluted sample	Diluent	Water	Dilution factor
STANDARD	25.0	—	—	—	1:1.00
—	—	—	—	—	—
—	—	—	—	—	—
—	—	—	—	—	—

O Reaction definition		O Reagent / Sample		● Validity checks	
Reaction check: None					
Maximum absorbance variation: —					

Configure assay parameters — Calibration			
<input type="radio"/> General	<input checked="" type="radio"/> Calibration	<input type="radio"/> SmartWash	<input type="radio"/> Results
Assay: UIBC-L	Calibration method: Linear		
<input checked="" type="radio"/> Calibrators	<input type="radio"/> Volumes	<input type="radio"/> Intervals	<input type="radio"/> Validity checks
Calibrator set: UIBC-L	Blank: Water	Calibrator level: 0†	Concentration: ††
Replicates: 2 [Range 1 – 3]	Cal 1: UIBC-L1		

O Calibrators		● Volumes		O Intervals		O Validity checks	
Calibrator: UIBC-L	Calibrator level	Sample	Diluted sample	Diluent	Water		
Blank: Water	25.0	—	—	—	—		
Cal 1: UIBC-L1	25.0	—	—	—	—		

O Calibrators		O Volumes		● Intervals		O Validity checks	
Calibration intervals:							
Full interval: 168 (hours)							
Calibration type:							
Adjust type: None							

O Calibrators		O Volumes		O Intervals		● Validity checks	
Blank absorbance range:							
Span: Blank – Blank							
Span absorbance range: — – —							
Expected cal factor: 0.00							
Expected cal factor tolerance %: 0							

Configure assay parameters — SmartWash				
<input type="radio"/> General	<input type="radio"/> Calibration	<input checked="" type="radio"/> SmartWash	<input type="radio"/> Results	<input type="radio"/> Interpretation
Assay: UIBC-L				
COMPONENT	REAGENT / ASSAY	WASH	Volume	Replicates
R1	AMIK9	Detergent A	345	1
R1	DIG00	Detergent A	345	1
R1	DGT0B	Detergent A	345	1
R1	GENT9	Detergent A	345	1
R1	TOBRA	Detergent A	345	1
R1	VANCO	Detergent A	345	1
R1	FRU0S	0.5% Acid Wash	345	1
R1	FRUL9	Detergent A	345	1
R2	AMIK9	Detergent A	345	1
R2	DIG00	Detergent A	345	1
R2	DGT0B	Detergent A	345	1
R2	GENT9	Detergent A	345	1
R2	TOBRA	Detergent A	345	1
R2	VANCO	Detergent A	345	1
R2	FRU0S	0.5% Acid Wash	345	1
R2	FRUL9	Detergent A	345	1

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Configure assay parameters — Results			
<input type="radio"/> General	<input type="radio"/> Calibration	<input type="radio"/> SmartWash	<input checked="" type="radio"/> Results
Assay: UIBC-L		Assay number: 2981	
Dilution default range:		Result units: µg/dL	
Low-Linearity: 25			
High-Linearity: 500			
Gender and age specific ranges:***			
GENDER	AGE (UNITS)	NORMAL	EXTREME
Male	0 – 130 (Y)	69 – 240	
Female	0 – 130 (Y)	70 – 310	

Configure result units	
Assay: UIBC-L	
Version: †	
Result units: µg/dL	
Decimal places: 0 [Range 0 – 4]	
Correlation factor: 1.0000	
Intercept: 0.0000	

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Configure assay parameters — Results				
<input type="radio"/> General	<input type="radio"/> Calibration	<input type="radio"/> SmartWash	<input checked="" type="radio"/> Results	<input type="radio"/> Interpretation
Assay: UIBC-L		Assay number: 2981		
Dilution default range:		Result units: µmol/L		
Low-Linearity: 4.5				
High-Linearity: 89.5				
Gender and age specific ranges:***				
GENDER	AGE (UNITS)	NORMAL	EXTREME	
Male	0 – 130 (Y)	12.4 – 43.0		
Female	0 – 130 (Y)	12.5 – 55.5		

Configure result units	
Assay: UIBC-L	
Version: †	
Result units: µmol/L	
Decimal places: 1 [Range 0 – 4]	
Correlation factor: 1.0000	
Intercept: 0.0000	

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www.goszdravnadzor.gov.ru

CERTIFIED SIGNATORY

I, the undersigned Valerio Tacchini, Notary Public in Milan, Italy, -----

CERTIFY THAT: -----

- BOSIO PAOLO, born in Alba on April 13th 1963, domiciled for his office in Milan, Robert Koch n. 2, acting as attorney and legal representative of "SENTINEL CH. S.P.A." with register office in Milan, Via Robert Koch n. 2, registered at the Trade Register of Milan - fiscal code 07118040158, duly empowered to do so, -----

of whom identity and powers I, as Notary, am certain, he signed, before me, the above document. -----

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