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OPERATOR MANUAL

for Medical Device:

Arctic Sun 5000 system and accessories for management of normothermia and hypothermia

Manufacturer: Medivance Inc., USA

State of Colorado
County of Boulder
Signed before me on Magh 14, 2022
by Scott Light
tricine Moan
(Notary's official signature)
(Notary Sofficial signature)
Notary Public
(Title of office)
(Title of office)
Ctober 04, 3022
(Commission Expiration)
(Commission Expiration)
LEANNE SLOAN
NOTARY PUBLIC
STATE OF COLORADO
NOTARY ID 20184039166
MY COMMISSION EXPIRES OCTOBER 04, 2022
0

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М.П. / Stamp

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Medivance

a wholly owned subsidiary of Bard Medical Division C.R. Bard, Inc.

321 South Taylor Ave., Suite 200 Louisville, CO 80027 (303) 926-1917

1. NAME OF THE MEDICAL DEVICE

Arctic Sun 5000 system for management of normothermia and hypothermia with accessories:

I. Basic kit:

- 1. Control module 1 ea.
- 2. Filling tube for tank 1ea.
- 3. Fluid delivery line with straps -1 ea.
- 4. Cable to connect system to the sensor YSI 400- 1
- 5. Operator manual on the CD or/and paper 1ea.
- 6. Power Cord 1 ea.
- 7. Hydrogel pads ArcticGel of different configurations (if necessary):
- 7.1. ArcticGel Universal hydrogel pad- 1 ea.;
- 7.2. ArcticGel XX Small hydrogel pad- 1 kit pack:

Pad back, right - 1 ea.,

Pad back, left - 1 ea.,

Pad tight, right -1 ea.,

Pad tight, left -1 ea.

7.3. ArcticGel Extra Small hydrogel pad- 1 kit packs

Pad back, right - 1 ea.,

Pad back, left - 1 ea.,

Pad tight, right -1 ea.,

Pad tight, left - 1 ea.

7.4. ArcticGel Extra Small hydrogel pad– 2 kits pack:

Pad back, right - 2 ea.,

Pad back, left - 2 ea.,

Pad tight, right -2 ea.,

Pad tight, left - 2 ea.

7.5. ArcticGel Small hydrogel pad- 1 kit pack:

Pad back, right - 1 ea.,

Pad back, left - 1 ea.,

Pad tight, right -1 ea.,

Pad tight, left - 1 ea.

7.6. ArcticGel Small hydrogel pad-2 kit pack:

Pad back, right - 2 ea.,

Pad back, left - 2 ea.,

Pad tight, right -2 ea.,

Pad tight, left - 2 ea.

7.7. ArcticGel Medium hydrogel pad – 1 kit pack:

Pad back, right - 1 ea.,

Pad back, left - 1 ea.,

Pad tight, right -1 ea.,

Pad tight, left – 1 ea.

7.8. ArcticGel Medium hydrogel pad- 2 kit pack:

Pad back, right - 2 ea.,

Pad back, left - 2 ea.,

Pad tight, right -2 ea.,

Pad tight, left -2 ea.

7.9. ArcticGel Large hydrogel pack – 1 kit pack:

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Pad back, left - 1 ea., Pad tight, right - 1 ea.,

Pad tight, left - 1 ea.

7.10. ArcticGel Large hydrogel pack- 2 kit pack:

Pad back, right - 1 ea., Pad back, left - 1 ea., Pad tight, right - 1 ea., Pad tight, left - 1 ea.

7.11. ArcticGel Small Universal hydrogel pad – 1 ea.

II. Accessories:

- 1. Calibration test unit − 1 ea.;
- 2. Screen protection kit 1 ea.;
- 3. Drainage tube 1 ea.;
- 4. Shunt tube -1 ea.;
- 5. Cleaning solution 1 ea.;
- 6. Temperature simulator 1 ea.;
- 7. Service kit 1 kit pack:
 - drainage tube 1 ea.;
 - shunt tube -1 ea.;
 - filling tube for tank 1 ea.;

(Further in the text - Arctic Sun 5000 system)

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2. INFORMATION ABOUT MANUFACTURER OF MEDICAL DEVICE

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Authorized representative of the manufacture in the Russian Federation:

Limited Liability Company Bard Rus Russian Federation, 121596, Moscow, Gorbunovastreet, 2, bld. 3, floor 5, premises II, room 1, office A-500.1 cabinet 5.

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1. Dymax Corporation

141 Zehner School Road Zelienople, PA 16063, CIIIA

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Тел: +7 495 775 85 82 complaints CIS@bd.com

3.THE MAIN INFORMATION ABOUT MEDICAL DEVICE

The Arctic Sun 5000 system for managing normothermia and hypothermia is a non-invasive, thermal regulating system, indicated for monitoring and controlling patient temperature between 32°C to 38,5°C.

The Arctic Sun system is composed of the Control Module and disposable non-sterile ArcticGel Pads, which are applied to areas of the patient's skin (body, thigh).

The control module circulates temperature-controlled water to the ArcticGel Pads. A commercially available medical temperature probe (not supplied by the manufacturer with the the system) connected to the control module senses the patient's core temperature. The atomatic mode increases or decreases the circulating water temperature to achieve a pre-set patient target temperature as determined by the clinician. The clinician can also adjust the water temperature (manual mode) to obtain the desired patient temperature.

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3.1.Intended use

Arctic Sun temperature management system allows to conduct artificial hypothermia to sever patients category in order to decrease metabolic process and reduce demand of the tissues in oxygen.

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The Arctic Sun 5000 Temperature Management System is a thermal regulating system, indicated for monitoring and controlling temperature patients of all ages.

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3.3. Classification of the device

Class of potential risk of use of the medical device	26
Product classification by type and duration of contact	System and accessories –no direct contact with patient skin, pads -long-term up to 30 days.
Frequency of use	Multiple use medical device, pads – single use
Software classification (ISO 62304)	Class B
Protection aginst electric shock	System – class 1, accessories – type BF
protection against ingress of water and particulate matter IPX0	IPXO
Mode of operation	Continiuous
ME EQUIPMENT and ME SYSTEMS intended for use in an OXYGEN RICH ENVIRONMENT	Not intended for use in an- oxygen rich environment
Software	version 3.0.2, March 2018

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3.4 WARNINGS

- Do not use the Arctic Sun system in the presence of lammable agents because an explosion and/or fire may esult.
- Do not use high frequency surgical instruments or ndocardial catheters while the Arctic Sun 5000 remperature Management System is in use.
- There is a risk of electrical shock and hazardous moving parts. There are no user serviceable parts inside. Do not emove covers. Refer servicing to qualified personnel.

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Power cord has a hospital grade plug. Grounding cliability can only be achieved when connected to an quivalent receptacle marked "hospital use" or "hospital rade".

When using the Arctic Sun 5000 system, note that all ther thermal conductive systems, such as water blankets and water gels, in use while warming or cooling with the arctic Sun 5000 Temperature Management System may ctually alter or interfere with patient temperature control. Do not place ArcticGel Pads over transdermal medication atches as warming can increase drug delivery, resulting in ossible harm to the patient.

The Arctic Sun 5000 system is not intended for use in the perating room environment.

5 CAUTIONS

This product is to be used by or under the supervision of ained, qualified medical personnel.

Use only sterile water. The use of other fluids will damage the Arctic Sun 5000 Temperature Management System.

When moving the Arctic Sun 5000 Temperature fanagement System always use the handle to lift the ontroller over an obstacle to avoid over balancing.

The patient's bed surface should be located between 75 n and 150 cm above the floor to ensure proper flow and hinimize risk of leaks.

The clinician is responsible to determine the ppropriateness of custom parameters. When the system is owered off, all changes to parameters will revert to the efault unless the new settings have been saved as new efaults in the Advanced Setup screen. For small patients 30 kg) it is recommended to use the following settings: Water Temperature High Limit \leq 40°C;

Vater Temperature Low Limit ≥10°C; Control Strategy = 1. It is recommended to use the Patient Temperature High and Patient Temperature Low alert settings.

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- The operator must continuously monitor patient temperature when using Manual Control and adjust the temperature of the water flowing through the pads accordingly. Patient temperature will not be controlled by the Arctic Sun 5000 Temperature Management System in Manual Control.
- Due to the system's high efficiency, Manual Control is not recommended for long duration use. The operator is advised to use the automatic therapy modes (e.g. Control Patient, Cool Patient, Rewarm Patient) for automatic patient temperature monitoring and control.
- The Arctic Sun 5000 Temperature Management System will monitor and control patient core temperature based on the temperature probe attached to the system. The clinician responsible for correctly placing the temperature probe and verifying the accuracy and placement of the patient probe at the start of the procedure.
- Manufacturer supplies temperature simulators (fixed value resistors) for testing, training and demonstration purposes. Never use this device, or other method, to circumvent the normal patient temperature feedback control when the system is connected to the patient. Doing so exposes the patient to the hazards associated with severe hypo- or hyper-thermia.
- Manufacturer recommends measuring patient emperature from a second site to verify patient emperature. Medivance recommends the use of a second patient temperature probe connected to the Arctic Sun 5000 Temperature Management System Temperature 2 input as a provides continuous monitoring and safety alarm patures. Alternatively, patient temperature may be verified beriodically with independent instrumentation.
- The displayed temperature graph is for general nformation purposes only and is not intended to replace tandard medical record documentation for use in therapy lecisions.
- Patient temperature will not be controlled, and alarms are not enabled in Stop Mode. Patient temperature may increase or decrease with the Arctic Sun 5000 Temperature Management System in Stop Mode.
- Carefully observe the system for air leaks before and luring use. If the pads fail to prime or a significant ontinuous air leak is observed in the pad return line, check onnections. If needed, replace the leaking pad.

eakage may result in lower flow rates and potentially ecrease the performance of the system.

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The ArcticGel Pads are non-sterile for single patient use. Do not reprocess or sterilize. If used in a sterile invironment, pads should be placed according to the physician's request, either prior to the sterile preparation or sterile draping. ArcticGel Pads should not be placed on a sterile field.

Use pads immediately after opening. Do not store pads once the kit has been opened.

Do not place ArcticGel Pads on skin that has signs of alceration, burns, hives, or rash.

While there are no known allergies to hydrogel materials, aution should be exercised with any patient who has a history of skin allergies or sensitivities. Examine the tackaging before use. Do not use in case of damage.

- Do not allow circulating water to contaminate the sterile ield when patient lines are disconnected.
- The water content of the hydrogel affects the pad's adhesion to the skin and conductivity, and therefore, the efficiency of controlling patient temperature. Periodically check that pads remain moist and adherent. Replace pads when the hydrogel no longer uniformly adheres to the skin. Replacing pads at least every 5 days is recommended.
- Do not puncture the ArcticGel Pads with sharp objects. unctures will result in air entering the fluid pathway and may reduce performance.
- If accessible, examine the patient's skin under the ArcticGel Pads often, especially those at higher risk of skin njury. Skin injury may occur as a cumulative result of pressure, time and temperature. Possible skin injuries include bruising, tearing, skin ulcerations, blistering, and necrosis. Do not place bean bag or other firm positioning devices under the ArcticGel Pads.

Do not place positioning devices under the pad manifolds or patient lines.

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Treatment application, monitoring and results are the responsibility of the attending physician.

f the patient does not reach target temperature in a easonable time or the patient is not able to be maintained at the target temperature, the skin may be exposed to low or high water temperatures for an extended period of time which may increase the risk for skin injury. Ensure that pad sizing/coverage and custom parameter settings are correct for the patient and treatment goals, water flow is greater han or equal to 2.3 liters per minute and the patient emperature probe is in the correct place. For patient cooling, ensure environmental factors such as excessively not rooms, heat lamps, and heated nebulizers are eliminated and patient shivering is controlled. Otherwise, consider ncreasing minimum water temperature, modifying target emperature to an attainable setting or discontinuing eatment. For patient warming, consider decreasing water temperature, modifying emperature to an attainable setting or discontinuing reatment.

Due to underlying medical or physiological conditions, some patients are more susceptible to skin damage from pressure and heat or cold. Patients at risk include those with poor tissue perfusion or poor skin integrity due to diabetes, peripheral vascular disease, poor nutritional status, steroid use or high dose vasopressor therapy. If warranted, use pressure relieving or pressure reducing devices under the patient to protect from skin injury.

Do not allow urine, antibacterial solutions or other agents to pool underneath the ArcticGel Pads. Urine and antibacterial agents can absorb into the pad hydrogel and cause chemical injury and loss of pad adhesion. Replace pads immediately if these fluids come into contact with the hydrogel.

Do not place ArcticGel Pads over an electrosurgical rounding pad. The combination of heat sources may result in skin burns.

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If needed, place defibrillation pads between the ArcticGel Pads and the patient's skin.

• Carefully remove ArcticGel Pads from the patient's skin at the completion of use. Discard used ArcticGel Pads in accordance with hospital procedures for medical waste.

The USB data port is to be used only with a standalone USB flash drive. Do not connect to another mains powered device during patient treatment.

Users should not use cleaning or decontamination methods different from those recommended by the manufacturer without first checking with the manufacturer hat the proposed methods will not damage the equipment. Do not use bleach (sodium hypochlorite) as it may damage the system.

Medivance will not be responsible for patient safety or auipment performance if the procedures to operate, maintain, modify or service the Medivance Arctic Sun 5000 Temperature Management System are other than hose specified by Medivance. Anyone performing the procedures must be appropriately trained and qualified.

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GETTING STARTED. SYSTEM SETUP

Unpack

Unpack the Arctic Sun 5000 Temperature Management System Control Module and accessories.

2) Allow the control module to remain upright for at least 2 hours prior to completing the installation and setup procedure in order to allow the chiller oil to settle. Damage to the chiller compressor may result otherwise.

Connections

Use only Medivance approved cables and accessories with the Arctic Sun 5000 Temperature Management System Control Module. Connect the Fluid Delivery Line, Patient Temp 1 cable, Patient Temp 2 cable (optional) and Fill Tube to the back of the control module.

2) Plug the Power Cord into the wall outlet. Position Arctic Sun 5000 Temperature Management System so that access to the power cord is not restricted.

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
Power On

- 1) Turn the power ON by activating the Power Switch.
- 2) The control module will automatically go through a brief self-test of the independent safety alarm.
- 3) A New User Training module option is available from the start up screen.
- 4) When the self-test is complete, the **Patient Therapy** Selection screen will appear on the control panel.

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Fill Reservoir

- 1) Fill the reservoir with sterile water only.
- 2) Four liters of water will be required to fill the reservoir at initial installation.
- 3) Add one vial of Arctic Sun 5000 Temperature Management System Cleaning Solution to the sterile water.
- 4) From the **Patient Therapy Selection** screen, press either the **Normothermia** or **Hypothermia** button, under the New Patient heading.

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- 5) From the **Hypothermia** or **Normothermia** therapy screen, press the **Fill Reservoir** button.
- 6) The **Fill Reservoir** screen will appear. Follow the directions on the screen.

Functional Verification

Perform the following functional verification procedure after initial setup and installation of the control module.

- 1) Power On the control module
- 2) From the **Patient Therapy Selection** screen, press the **Hypothermia** button to display the **Hypothermia** therapy screen.
- From the Hypothermia therapy screen, press the lanual Control button to open the Manual Control window.
- 4) Use the Up and Down arrows to set the **Manual Control** water target temperature to 40°C and the duration to 30 minutes.
- 5) Press the **Start** button to initiate **Manual Control**. Allow at least 3 minutes for the system to stabilize.
- 6) Monitor the flow rate and water temperature in the System status area on the **Hypothermia** therapy screen.
- 7) Verify that the flow rate reaches at least 1.5 liters/minute.
- 8) Verify that the water temperature increases to 30°C.
- 9) Press the Stop button.
- 10) Set the Manual Control water target temperature to 4°C and the duration to 30 minutes.
- 11) Press the Start button to initiate Manual Control.

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- 12) Monitor the flow rate and water temperature in the **System** status area of the **Hypothermia** therapy screen. Verify that the water temperature drops to 6°C.
- 13) Press the **Stop** button to stop **Manual Control**
- 14) Press the **Cancel** button to close the **Manual Control** window
- (5) Power Off the control module.

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Place ArcticGel Pads

4.1. Patient Therapy

Examine each pad for damage prior to placement.

Connect ArcticGel Pads

Orient the blue and white colors on the pad line connector and Fluid Delivery Line. While holding the pad line tubing, insert the clear pad line connector into the Fluid Delivery Line manifold. Do not press or squeeze the wings when connecting. The connector will click into place. Check the surface of the device for mechanical damage prior to use.

Directions for use

l ArcticGel Pads are only for use with an Arctic Sun 5000 Temperature Management System Control Module.

2. Select the proper number, size and style pad for the patient size and clinical indication. However, the rate of temperature change and potentially the final achievable temperature is affected by pad surface area, patient size, ad placement and water temperature range. Best system performance will be achieved by using the entire pad set (4). If the entire set of pads is not used, the minimum flow rate may not be achieved.

3. For patient comfort, the pads may be prewarmed using Water Temperature Control Mode (manual) prior to application.

4. Place the pads on healthy, clean skin only. Remove any creams or lotions from patient's skin before pad application. Remove the release liner from each pad and apply to the appropriate area. The pads may be overlapped or folded adhesive-to-adhesive to achieve proper placement. The pads may be removed and reapplied if necessary. The pad surface must be contacting the skin for optimal energy transfer efficiency. Place pads to allow for full respiratory excursion.

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5. Attach the pad's line connectors to the patient line manifolds. Begin circulating water through the pads using either Patient Temperature Control Mode (automatic) or Water Temperature Control Mode (manual). If the pads fail to prime or a significant continuous air leak is observed in the pad return line, check connections, then if needed replace the leaking pad.

6. Once the pads are primed, assure the flow rate displayed on the control panel is greater than 1.7 liters per minute, which is the minimum flow rate for a full pad kit (4).

7. When finished, empty water from pads. Cold temperature increases the adhesiveness of the hydrogel. For ease of removal, leave pads on the patient for approximately 15 minutes to allow the hydrogel to warm. Slowly remove pads from the patient and discard.

Directions for ArcticGel Small Universal Pad

• The ArcticGel Small Universal Pads are only for use with the Arctic Sun Temperature Management System. See Operators Manual for detailed instructions on system use.

Select the proper number of pads for the patient size and clinical indication. However, the rate of temperature change and potentially the final achievable temperature is affected by pad surface area, patient size, pad placement and water temperature range. Best system performance will be achieved by using the maximum number of pads.

The following sizing chart is provided as guidance. Modify the number of pads and locations as necessary to meet specific clinical needs.

Patient Weight	Number of Pads
2.5 - 5 kg (5.5 - 10 lbs)	1
5 - 10 kg (10 - 22 lbs)	2
10 - 16 kg (22 - 35 lbs)	3

Place the pads on healthy, clean skin only. Locate pad edges away from articulating areas of the body to avoid irritation. Place pads to allow for full respiratory excursion.

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- The pad surface must be contacting the skin for optimal energy transfer efficiency. The Small Universal Pads are provided with a cloth liner over the hydrogel. The pads may be used with the cloth liner in place or removed to expose the hydrogel adhesive. If the cloth liner is used, secure the pads with the Velcro straps to ensure good contact with the skin. The pads may be removed and reapplied if necessary.
- Due to the small patient size and the potential for rapid patient temperature change
- It is recommended to use the following settings: Water Temperature High Limit: ≤40°C (104°F) Water Temperature Low Limit: ≥10°C (50°F) Control Strategy: 2
- It is recommended to use the Patient Temperature High and Patient Temperature Low alert settings.
- Place a core patient temperature probe and connect to the Arctic Sun Temperature Management System patient temperature 1 connector for continuous patient temperature feedback. A rectal or esophageal temperature probe is recommended.
- Attach the pad's line connectors to the fluid delivery line manifolds. Begin treating the patient.
- If the pads fail to prime or a significant continuous air leak is observed in the pad return line, check connections, then if needed replace the leaking pad.
- Once all the pads are primed, assure the steady state flow rate displayed on the control panel is appropriate for the number of pads connected.

Number Pads 1 2 3 Minimum Flow Rate 1/m 0.9 1.7 2.4

•When finished, purge water from pads. Slowly remove pads from the patient and discard.

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Temperature Probe Placement

Patient temperature control with the Arctic Sun 5000 Temperature Management System requires patient temperature feedback provided by an indwelling patient temperature probe connected to the Patient Temperature 1 connector on the back of the Control Module.

Patient Therapy Selection

Use the Patient Therapy Selection screen to initiate a New Patient, Continue a Current Patient, or access the Advanced Setup screen.

New Patient - Normothermia

Select **Normothermia** if the therapy goal is to maintain a patient temperature at a pre-defined target temperature for n indefinite period of time. Press the **Normothermia** button to display the **Normothermia** therapy screen.

New Patient - Hypothermia

Select **Hypothermia** to reduce and maintain a patient temperature at a set target temperature for a defined period of time, then slowly re-warm the patient at a controlled rewarming rate. Press the **Hypothermia** button to display the **Hypothermia** therapy screen.

Additional Protocol Option

Two additional protocols (**Hypothermia** or **Normothermia**) may be visible on the **Patient Therapy Selection** screen.

urrent Patient

The Continue Current Patient button and the date and time that the current therapy was paused will display on the Patient Therapy Selection screen if a patient therapy was paused within the past 6 hours.

Press the **Continue Current Patient** button to resume a paused patient therapy.

Initiate Normothermia (Control Patient & Rewarm

Normothermia therapy is initiated and managed, and patient temperature is automatically controlled to a set target temperature from the Control Patient window in the Normothermia therapy screen. The Control Patient window displays the patient target temperature and the duration since the

initiation of normothermia therapy.

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To initiate Normothermia therapy:

- l) From the **Patient Therapy Selection** screen, press the **Normothermia** button to display the **Normothermia** therapy screen.
- 2) The default patient target temperature will display in the Control Patient window.
- 3) To modify the patient target temperature, press the Adjust button to display the Control Patient-Adjust window.
- 4) Control Patient to: Use the Up and Down arrows to set the desired patient target temperature to control the patient.
- 5) Rewarm at a Rate of: Use the Up and Down arrows on the right of the screen to set the rewarming rate.
- (i) Press the **Save** button to save the new settings and close the **Control Patient-Adjust** window
- 7) Press Start, in the Control Patient window to initiate herapy. You will hear a tone and then a voice stating "Therapy Started". Additionally, the Control Patient window and the Arctic Sun 5000 Temperature Management System icon will blink, indicating that therapy is in progress.

Initiate Hypothermia (Cool Patient and Rewarm Patient)

Hypothermia therapy is initiated and managed, and patient imperature is automatically controlled to a set target imperature from the Cool Patient and Rewarm Patient mindows in the Hypothermia therapy screen.

The Cool Patient window displays the cooling phase patient target temperature and the length of time remaining in the cooling phase of the

Hypothermia therapy.

the **Rewarm Patient** window displays the rewarming phase patient target temperature and the length of time remaining in the rewarming phase of the **Hypothermia** therapy.

To initiate hypothermia therapy:

From the **Patient Therapy Selection** screen, press the **Hypothermia** button to display the **Hypothermia** therapy screen.

1. Cool Patient Settings

- The default patient target temperature and duration will display in the **Cool Patient** window.
- To modify the patient target temperature and duration, press the **Adjust** button to display the **Cool Patient-Adjust** window.

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- Cool Patient To: Use the Up and Down arrows on the left side to set the desired patient target temperature to cool the patient
- •Cool Patient For: Use the Up and Down arrows on the right side to set the cooling duration to cool the patient before rewarming begins.
- Press the **Save** button to save the new settings and close the **Cool Patient-Adjust** window

1. Rewarm Patient Settings

- The default patient target temperature and duration will display in the **Rewarm Patient** window.
- 'To change the rewarming phase patient target temperature and rewarming rate, press the **Adjust** button in the **Rewarm Patient** window
- nd Down arrows on the left side to set the desired final patient target temperature.
- Rewarm Patient To: Use the Up and Down arrows on the right side to set the desired final patient target temperature.
- Rewarm at a Rate of: Use the Up and Down arrows in the center of the screen to set the rewarming rate.
- Rewarm Patient From: When cooling a patient, adjustment of the Rewarm Patient From setting on the left side of the screen is disabled and defaults to the Cool Patient target temperature.
- When rewarming a patient, the **Rewarm Patient From** sujustment is enabled and the value can be modified. The **Rewarm Patient From** setting is the temperature to which the system is currently controlling the patient. The **Rewarm Patient From** temperature will automatically increase as the rewarming process continues. This feature allows the rewarming procedure to be optimized by clowing complete control of the rewarming ramp.
- Using the Rewarm Patient From temperature, the Rewarm Patient To temperature and the rewarming rate settings, the system will calculate and display the rewarming duration and the date/time at which the patient will reach the final rewarming target temperature.
- Press the **Save** button to save the new settings and close the **Rewarm Patient-Adjust** window.

3. Initiate Patient Cooling

Press Start, in the Cool Patient window to initiate therapy. You will hear a tone and then a voice stating Therapy Started". Additionally, the Cool Patient window and the Arctic Sun 5000 Temperature Management System ion will blink, indicating that therapy is in progress.

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Initiate Patient Rewarming

Upon completion of the cooling phase, there are two options for initiation of patient rewarming, either Automatically or Manually, depending on the **Rewarming** Begins setting in **Hypothermia Settings**.

If Rewarming Begins is set to Automatically, the ewarming process starts automatically when the Cool Patient therapy is complete and the duration reaches zero. If Rewarming Begins is set to Manually, the rewarming process starts when the Start button is pressed in the Rewarm Patient window. The

woling process will continue until the Rewarm Patient Start button is pressed. An Alert will occur when the Cool Patient duration reaches zero. When the Rewarm Patient duration timer reaches zero, the system will continue to control the patient to the target temperature until the Stop button is pressed. Once in Normothermia, the timer will reset and begin tracking the duration of Normothermia therapy.

Ind Therapy

From the **Normothermia** therapy or **Hypothermia** herapy screen, press the **Stop** button to terminate water inculation to the pads.

Press the **Empty Pads** button, and follow the instructions purge the pads of water.

Disconnect the pads from the Fluid Delivery Line.

Slowly and carefully remove pads from the patient skin.

Discard the used pads in accordance with hospital medical waste.

Press the power switch Off.

Ipower is lost while the power switch is in the On position, naudible alert will be issued until it is switched Off. This alerts the user that the treatment may have been acidentally stopped.

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
Normothermia Settings

Use the **Normothermia Settings** screen to view the current settings and modify the settings for the following parameters. To modify any parameter setting, press the adjust button to the right of the parameter.

Normothermia Settings screen parameters:

Water Temperature Settings

- Pre-Condition Water
- Manual Control
- High Water Limit
- ·Low Water Limit

Patient Temperature Settings

- High Patient Alert
- ·Low Patient Alert
- ·Control Strategy

Display Settings

- ·Temperature Units
- ·Temperature Units Adjust
- Patient Temp 2

To access the Normothermia Settings screen:

- 1) Press Adjust on the Control Patient window.
- 2) Press the More button on the Control Patient Adjust window.
- The Normothermia Settings screen will be displayed.
- 4) To save the new settings as the current patient therapy settings, press the Close button. For instructions on saving the settings as the system defaults, see Advanced Setup.

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NOTE: Verify that the Patient Temperature 1 probe is correctly positioned and correctly connected to the system. If the patient temperature change is less than 0.15°C after the first hour of therapy, the system will generate Alert 116 – Patient Temperature 1 Change Not Detected. If Alert 116 is not acknowledged after 5 minutes, the system will generate Alarm 117 – Patient Temperature 1 Change Not Detected. Alarm 117 will stop therapy and an audible alarm will sound. After Alarm 117 is acknowledged, therapy will require a restart.

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Hypothermia Settings

Use the **Hypothermia Settings** screen to view the current settings and modify the settings for the following parameters. To modify any parameter setting, press the **Adjust** button to the right of the parameter.

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Temperature Units Adjust

Patient Temp 2

To access the Hypothermia Settings screen:

Press Adjust on the Cool Patient window or the Rewarm Patient window.

?) Press the More button on the Cool Patient Adjust window or Rewarm Patient Adjust window.

The **Hypothermia Settings** screen will be displayed. To save the new settings as the current patient therapy witings, press the **Close** button. For instructions on saving the settings as the system defaults, see **Advanced Setup**.

NOTE: Verify that the Patient Temperature 1 probe is wrectly positioned and correctly connected to the system. If the patient temperature change is less than 0.15°C after the first hour of therapy, the system will generate

llert 116 – Patient Temperature 1 Change Not Detected. If Alert 116 is not acknowledged after 5 minutes, the system will generate Alarm 117 – Patient Temperature 1 Change Not Detected. Alarm 117 will stop therapy and an udible alarm will sound. After Alarm 117 is knowledged, therapy will require a restart.

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4.4. Advanced Setup

Use the **Advanced Setup** screen to view the current settings and modify the settings for the following parameters. To modify any parameter setting, press the **Adjust** button to the right of the parameter.

Location / Time Settings

Language

Number Format

Current Time

Date Format

Current Date

the following functions can be initiated from the dvanced Setup screen.

Download Patient Data: The Patient Data for the last 10 km) cases are stored on the Arctic Sun 5000 Temperature Management System hard drive. This data is maintained then the Arctic Sun 5000 Temperature Management System is powered down, or in the event of a total loss of lower.

Calibration

Total Drain

Save All Settings As Default

Upload custom file

dditionally, the following information can be viewed in Advanced Setup screen.

Software Versions

Last Calibration Date

Next Calibration Due

access the Advanced Setup screen:

Press Advanced Setup button on the Patient Therapy Section screen.

The Advanced Setup screen will be displayed.

To access the Additional Protocol Selection screen:

lefer to the Arctic Sun 5000 Temperature Management system Help screens for

formation regarding additional protocol setup.

5. Alarms and Alerts

he Arctic Sun 5000 Temperature Management System afety system continually monitors the state of the device and the patient, and issues alarms or alerts to notify the user of conditions that may interfere with patient safety or system performance. There are two types of conditions:

An Alarm notifies the user that a condition that may potentially pose an unsafe situation with respect to the patient or the device. An Alarm is a High Priority condition that requires immediate operator response.

An Alert informs the user about patient and device status without interrupting the procedure. An Alert is a Medium Priority condition that requires prompt operator response.

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Alarms

An Alarm is denoted by an audio signal that repeats every seconds until the Alarm is cleared. The Alarm screen will appear that displays the alarm number, alarm title, a description of the problem or conditions that triggered the slarm, and solutions and instructions for troubleshooting and resolving

the alarm condition. If certain Alarm conditions are not acknowledged by the operator within 2 minutes, a deminder tone will sound. All Alarm settings are maintained in the event of a power interruption.

Main Safety Alarms

While there are multiple alarms and safety features in the Arctic Sun 5000 Temperature Management System, there are five main safety alarms that will place the device into Mop mode until the condition is addressed.

Alarm Specification

ligh Patient Temperature 39.5°C

Low Patient Temperature 31.0°C

High Water Temperature 42.5°C / 44°C

Low Water Temperature 3.0°C / 3.5°C

System Self-Test Failure At device power ON

System is powered On, a system self test for the independent safety alarm is automatically run.

This test simulates a "water high temperature" fault intuation on both the primary and secondary water emperature sensors. Both the primary and secondary afety systems must respond to the fault and be verified by the opposing safety system. If either safety systems do not appropriately either an alarm 80 or 81 will be issued. Contact Customer Support.

Non-Recoverable Alarms

If an Alarm condition occurs that prevents proper use of the device or proper patient treatment (such as the five main afety alarms discussed above), the system is placed into stop mode and will not allow therapy to continue.

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This type of Alarm is known as Non-Recoverable. If this situation occurs, cycle the device power (turn device Off then On). If the alarm recurs contact Customer Support.

Recoverable Alarms

Other Alarms that temporarily Stop the device until the user is able to correct the cause and clear the Alarm are classified as Recoverable. If the condition that initiated the ularm is not addressed and problem persists, the Alarm will recur.

fa Recoverable Alarm occurs:

-) When an alarm is issued the device is placed into **Stop**
- (1) Read the displayed instructions.
- Note the Alarm number.
- Press the Close button to clear the alarm.
- f) Follow the instructions to correct the alarm condition. Perform the actions in the order listed until the alarm modition is resolved.
- note you have cleared the alarm, press the Start button in the therapy window to restart therapy. You will hear a une and a voice stating "Therapy Started". Additionally, the active therapy window and the Arctic Sun 5000 Temperature Management System icon will blink.
- The condition does not resolve, contact Customer Support.

Merts

Alerts are denoted by an audio signal that repeats every 25 monds. The Alert screen will appear that displays the alert number, alert title, a description of the problem that tiggered the alert, and solutions and instructions for mubleshooting and resolving the alert condition.

If an Alert occurs:

- 1) Read the displayed instructions.
- 2) Note the Alert number.
- 3) Press the Close button to clear the alert.
- 4) Follow the instructions to correct the alert condition. Perform the actions in the order listed until the alarm condition is resolved. If the condition does not resolve, contact Customer Support.

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5) Refer to the Arctic Sun 5000 Temperature Management System Help screens for additional information regarding plarms and alerts.

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5. MAINTENANCE AND SERVICE

Maintenance Schedule

Procedure	Interval	Who performs
		Medical
Clean	As required	
external		personnel
surfaces		
Inspect	6 month	Authorized
connectors		service support
and cables		
Clean the	6 month	Authorized
condenser		service support
Replenish	6 month	Authorized
Cleaning		service support
solution		
Inspect	6 month	Authorized
screen		service support
protector		
Calibration	Every 2000	Authorized
	hours or 250	service support
	uses, whichever	
	occurs first, as	
	indicated by	O
	system display	
Inspect fluid	6 month	Authorized
delivery line		service support

Cleaning and Maintenance

loutine cleaning and preventive maintenance should be reformed on the Arctic Sun 5000 Temperature lanagement System control module every 6 months at a minimum. This consists of cleaning the external surfaces, accessories and chiller condenser, inspecting the device, and replenishing the internal cleaning solution that appresses microorganism growth in the water reservoir and hydraulic circuit.

Ixternal Surfaces

Clean the exterior body of the control module, fluid elivery lines, power cords and temperature cables using a

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of cloth and mild detergent or disinfectant according to uspital protocol.

Medivance has qualified and approves the use of the blowing types of disinfectants for exterior surfaces: wium hypochlorite, isopropyl alcohol, and quaternary mmonium

Condenser

A dirty chiller condenser will significantly reduce the poling capacity of the control module.

To clean the condenser, wipe the dust from the exterior pill using a soft cloth. Depending on the quality of your institution's air, periodically remove the back cover and acuum or brush the condenser fins. At a minimum the andenser fins should be cleaned annually. Maintenance civities should be performed by qualified personnel.

Device Inspection

Periodically inspect the external areas of the device for amaged, loose or missing parts, and frayed or twisted over cords and cables.

Discontinue using the device displaying one or more of the above conditions until the problem is corrected and has then verified to be operating correctly.

Replenish Internal Cleaning Solution

Contact Customer Service to order internal cleaning

To replenish the internal cleaning solution:

Drain the reservoir.

Turn control module power Off.

Attach the drain line to the two drain ports on the back of the control module. Place the end of the drain line into a wntainer. The water will passively drain into the container. The reservoir.

From the Hypothermia therapy screen or the Normothermia therapy screen, press the Fill Reservoir button.

• The Fill Reservoir screen will appear. Follow the directions on the screen.

Add one vial of Arctic Sun 5000 Temperature Management System cleaning solution to the first bottle of sterile water.

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The filling process will automatically stop when the servoir is full. Continue to replace the bottles of sterile later until the filling process stops.

When the Fill Reservoir process is complete, the screen sill close.

spect Screen Protector

he Control Panel's touchscreen is supplied with a isposable screen protector. If it becomes damaged, it can be removed by lifting the edge and carefully peeling it from the screen. To ensure dust and particulates are removed, the the touchscreen using isopropyl alcohol. Remove the liner from the screen protector. Then carefully apply the protector to the screen with the liner side down against the screen.

spect Fluid Delivery Line

Power On the system

From the Patient Therapy Selection screen press the hypothermia button to display the Hypothermia therapy creen.

From the Hypothermia therapy screen, press the Manual ontrol button to open the Manual Control window.

Set the Manual Control water target temperature to 28°C and the duration to 30 minutes.

Connect a shunt to a set of fluid delivery line ports.

Press the Help button and then press the Help Index auton. Select the topic Maintenance and Service and subspic System Diagnostics then press the Display button. Werify that inlet pressure is -7 ± 0.2 .

Repeat on all valves. If inlet pressure is out of range, place the two valves that the shunt is connected to.

Ensure that the shunt is removed before device is put ack in service.

Preventative Maintenance

Use of the Arctic Sun 5000 Temperature Management System in excess of 2,000 hours, without conducting preventative maintenance, may result in failure of certain system components and failure of the system to function as intended. To maintain system performance, the Arctic Sun 5000 Temperature Management System requires periodic envice of key components.

Peventive maintenance should be performed by Authorized representatives of the service support. Service includes

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regulating of electronic components, cleanin soluton replenish, change of old components, software update (if necessary).

Contact AR "Bard Rus" LLc to confirm contacts of Authorized service provider.

Software Update

Software updates will be provided on a flash drive. Installation of software updates will be performed via the USB port on the front of the control module. The software update feature will automatically initiate if the control module detects the proper files on a flash drive inserted in the USB port at power On.

To install software update:

)Insert flash drive into the USB port.

An image of a timer will display while the software pdate in being installed, and will disappear when the oftware installation process is completed.

After installation, the new software version will display the Software Version field in Advanced Setup.

Service

ontact AR "Bard Rus" Llc for technical support and ustomer service instructions to enable appropriately palified technical personnel to repair those parts of the puipment that are considered repairable.

I. Calibration

11.1. Calibration Test Unit

Aseparate device, the Calibration Test Unit (CTU) is equired to perform periodic calibration on the Arctic Sun 5000 Temperature Management System.

When to Perform a Calibration or Calibration Check

l. Calibration is recommended after 2000 hours of peration or 250 uses, whichever occurs first. The status of Calibration is available in the Advanced Settings screen.

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2. In addition, Calibration may be required after replacing certain components (Manifold Harness, Tank Temperature Sensor Harness, Processor Circuit Card, Isolation Circuit Card, Input/Output Circuit Card).

3. A Calibration Check confirms the device flow, ability to heat and cool, and the temperature sensing systems are all within specification. During the calibration check errors may be displayed with diagnostic information assisting with performance or calibration issues. After successful completion of a calibration check, a report is displayed showing a pass or fail status of all parameters checked.

51.2. Performing a Calibration

To perform a calibration on the Arctic Sun 5000 Temperature Management System, press the Advanced Setup button on the Therapy

Selection screen. Press the Start button next to Calibration and follow the on-screen instructions.

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ELECTOMAGNETIC COMPATABILITY

Medical electrical equipment needs special precautions regarding electromagnetic compatibility. Ensure that the Arctic Sun 5000 Temperature Management System is installed and used according to the electromagnetic management in provided. The following are

mmpatibility information provided. The following are uidance and manufacturer's declarations regarding dectromagnetic compatibility for the Arctic Sun 5000 imperature Management System.

The use of accessories or cables other than those specified a sold by Medivance (shown below) is not recommended. Use of unapproved accessories or cables may result in accessed emissions or in decreased immunity of the Arctic sun 5000 Temperature Management System.

If the Arctic Sun 5000 Temperature Management System sused directly adjacent to or stacked with other quipment, the user should periodically observe the Arctic Sun 5000 Temperature Management System device to crify it operates normally in that environment.

Portable and mobile RF communications equipment can fleet Medical Electrical Equipment.

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ccessories and cables approved by Medivance for use with the Arctic Sun® Temperature Management System	Part#
emperature In-Cable - Nellcor	735-02
emperature In-Cable - BARD®	735-03
emperature In-Cable - Rusch	735-04
emperature In-Cable - GE	735-05
emperature In-Cable - Philips	735-06
emperature Out-Cable - Nellcor	735-52
emperature Out-Cable - BARD®	735-53
emperature Out-Cable - Rusch	735-54
Temperature Out-Cable - GE	735-55
Temperature Out-Cable - Philips	735-56
Power Cord, US, Canada, Mexico	733-00
Power Cord, Continental Europe	733-01
Power Cord, UK, Ireland	733-02
Power Cord, Australia, New Zealand	733-03
Power Cord, Mainland China	733-04
Power Cord, Brazil	733-05
Power Cord, Switzerland	733-07
Power Cord, South Africa	733-08
Transmission Interface Module (TIM) Kit	760-00
Transmission Interface Module (TIM)	761-00
R\$232 Cord	762-00
1.1 EN/IEC 60601-1-2 Table 1	
Guidance and Manufacturer's Declaration - Electromagnetic Emissi	ons

electromagnetic environ	nment specified bel	ow. The customer or the end user of the ystem should assure that it is used in such
Emissions test	Compliance	Electromagnetic environment - guidance

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The ARCTIC SUN® Temperature Management System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

	_	
1.1 EN/IEC 60601-1-2 Table	1	(cont.)

Quidance and Manufacturer's Declaration	n - Electromagnetic Emissions
--	-------------------------------

Emissions Test	Compliance	Electromagnetic environment guidance
RF emissions CISPR 11	Class A	The ARCTIC SUN® Temperature Management System unit is suitable for use in all establishments other than domestic, establishments and those directly connected to the public
Harmonic emissions IEC 61000-3-2		
Voltage fluctuations/ Ficker emissions IEC #1000-3-3	Complies	low-voltage power supply network that supplies buildings for domestic purposes.

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
Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The ARCTIC SUN® Temperature Management System unit is intended for use in the electromagnetic environment specified below. The customer or the end user of the ARCTIC SUN® Temperature Management System unit should assure it is used only in such an environment.

such an environment.				
Immunity Test	IEC60601 test level	Compliance Level	Intended Electromagnetic Environment	
Electromagnetic discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.	
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/ output lines	±2kV for power supply lines ±1kV for input/ output lines	Mains power quality should be that of a typical commercial or hospital environment.	
Surge IEC 61000-4-5	±1kV differential mode (line-line) ±2kV common mode (line-earth)	±1kV differential mode (line-line) ±2kV common mode (line-earth)	Mains power quality should be that of a typical commercial or hospital environment.	
Voltage dps, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 seconds	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the ARCTIC SUN® Temperature Management System unit requires continued operation during power mains interruptions, it is recommended that the ARCTIC SUN® Temperature Management System unit be powered from an uninterruptible power supply with sufficient capacity to run the unit for the maximum required time of interruption.	
Power frequency (50/60Hz) magnetic field EC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or baselist environment	

Note: UT is the a.c. mains voltage prior to application of the test level.

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13 EN/IEC 60601-1-2:2007 Sub-clause 5.2.2.2 Table 4:

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

Immunity	IEC60601	Compliance	Intended Electromagnetic Environment
Test	test level	Level	
Conducted	3Vrms	3Vms	Portable and mobile RF communications equipment should be used no closer to any part of the Arcnic Sun® Temperature Management System unit, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance d = 1.2\p 80MHz to 800 MHz d = 2.3\p 800MHz to 8.06Hz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended minimum separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey*, should be less than the compliance level in each frequency range.* Interference may occur in the vicinity of equipment marked with the following symbol:
RF	150kHz to	150kHz to	
EC 61000-46	80MHz	80MHz	
REadiated	3V/m	3V/m	
RF	80MHz to	80MHz to	
EC 61000-43	2.5GHz	2.5GHz	

IOTE 1 At 80MHz and 800MHz, the higher frequency range applies
IOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation
infected by absorption and reflection from objects, structures and people.

Field strengths from fixed transmitters, such as base stations for radio (cellular) subsets telephones and land mobile radios, amateur radio, AM and FM radio substants and TV broadcast cannot be predicted theoretically with accuracy. To assess telectromagnetic environment due to fixed RF transmitters, an electromagnetic site say should be considered. If the measured field strength in the location in which the kinc Sun® Temperature Management System unit is used exceeds the applicable Frompliance level above, the Arcmc Sun® Temperature Management System unit study be observed to verify normal operation. If abnormal performance is observed, utilized measures may be necessary, such as re-orienting or relocating the licen Sun® Temperature Management System unit.

1. Over the frequency range 150kHz to 80MHz, field strengths should be less than 3V/m.

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Temperature Range

Operating:.....10°C to 27°C
Storage:.....-30°C to 50°C

At operating temperatures higher than 27°C, the efrigeration system's cooling capacity and therefore its bility to cool a patient is compromised.

3.DISPOSAL

Ipon end of life, dispose of in accordance with local WEEE regulations.

lystem, pads and accessories in accordance with SanPiN 1.3684-21 medical waste – class A.

4.WARRANTIES

Manufacturer warrants to the original purchaser that Arctic Sun 5000 system for managable normo and hypothermia with accessories ("Equipment") and ArcticGel Pad "Disposables") will be free of defects in workmanship and materials for the period set forth in the labeling and if no such period is set forth, then one year from the date of installation. If the Equipment or a Disposable proves to be o defective, such Equipment or Disposables may be epaired, replaced, refunded or credited, at manufacturer's ption. An extended warranty for Equipment is available ir purchase. The warranty covers all parts and labor associated with defects in material and workmanship of the quipment. Authorized service will, at its discretion, letermine if the Equipment is to be repaired on site, or at he manufacturers' service center. If Equipment is to be sturned for service. Any unauthorized equipment repair erformed during the warranty period will void the varranty. All returns must be authorized in advance by authorized service. The liability of Medivance Inc. under is product warranty does not extend to any abuse, ecidental damage, misuse, improper storage, alteration,

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15.CLAIM

For questions about the circulation of the medical device in the Russian Federation, please contact Authorized Representative of the Manufacturer in the Russian Federation:

Limited Liability Company « Bard Rus»

Russian Federation, 121596, Moscow, Gorbunova street, 2, bld. 3, floor 5, premises II, room 1, office A-500.1 cabinet

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Parameter	Specification
Therapy Modes	Normothermia: Control Patient, Rewarm Patient Hypothermia: Cool Patient, Rewarm Patient
Heater Capacity	2500 BTU/hr / 750 Watts
Circulating Fluid	Sterile Water
Reservoir Capacity	3.5 liters
Water Flow Rate	5 liters per minute
Patient Probe Type	YSI 400 Series compatible
Patient Temperature Inputs	Patient Temp 1: control, monitor, alarm Patient Temp 2: monitor, alarm
Patient Temperature Display Range	10°C to 44°C 50°F to 111.2°F 0.1°C /°F increments
Patient Temperature Measurement Accuracy	±0.4°C (10°C to 32°C) ±0.2°C (32°C to 38°C) ±0.4°C (38°C to 44°C) Includes ±0.1°C external probe
Responses of the PCLCS (Physiologic Closed-Loop Control System)	Settling Time: ~4.5 hrs Relative Overshoot: <0.5°C Command Overshoot: <0.5°C Response Time: Warming (max) 33°C to 37°C: ~6 hrs Cooling 37°C to 33°C: ~2 hrs Steady State Deviation: 0 Tracking Error: 0 Note: All values derived from testing in simulated use.
Patient Temperature Control Range	32°C to 38.5°C 89.6°F to 101.3°F 0.1°C /°F increments
Water Temperature Display Range	3°C to 45°C / 37.4°F to 113.0°F 0.1°C /°F Increments
Water Temperature Control Range (Manual)	4°C to 42°C / 39.2°F to 107.6°F 1°C /°F increments
High Water Temperature Limit	36°C to 42°C / 96.8°F to 107.6°F 1°C /°F increments
Low Water Temperature Limit	4°C to 25°C / 39.2°F to 77°F 1°C /°F increments
Time to heat water from 20°C to 37°C	8 minutes (approximate)
Sound Pressure	Alarm Tone: 70dB to 80dB at 1 meter, repeats every 10 seconds Alert Tone: 63dB to 71dB at 1 meter, repeats every 25 seconds Reminder Tone: 65dB at 3 meters, 0.5 seconds on/20 seconds off
Mains Input	100-120VAC, 50-50Hz, 11A 220-240VAC, 50-60Hz, 5.5A
Leakage Current	<300 µA
Operating Relative Humidity Range	5% to 70% non-condensing
Storage Relative Humidity Range	5% to 95% non-condensing
Operating Temperature Range	10°C to 27°C / 50°F to 80°F
Storage Temperature Range	-30°C to 50°C / -20°F to 120°F
Atmospheric Pressure Range	60 kPa to 110 kPa
Dimensions	Height: 35 inches (89 cm) Width: 14 inches (36 cm) Depth: 18.5 inches (47 cm)
Weight	Empty: 43 kg / 95 lbs ; Filled: 47 kg / 103 lbs

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Acrtic Sun 5000 temperature management system is indendent for use with ArcticGel pads only.

• Arctic Sun 5000 system can be used with temperature probes YSI 400 Series compatible and probes of other manufacturers in accordance with the technical specifications and instruction for use. Temperaure probes are not supplied with the system.

Table to of temperature and resistance correlation for comparatable sensors:

Temperature, °C	Resistance, Ω
32	1,667
33	1,599
34 35	1,533
35	1,471
36	1,412
37	1,355
38	1,301
39	1,249

Cable for sensor YSI 400 series 1 pc. is included into the kit of the device and is supplied with the ArcticSun 5000 system. For connection of the sensors cabels with the connectors GE, Philips, Rusch are used and ordered separately and are not supplied with the system.

List of the temperature sensors compatable with the ArcticSun 5000 system consist of, but not limited to the following manufacturers:

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
- Temperature probe sensor for general indication 9F/12F (not more that 10 pcs) accessory to the medical device "Patient monitor Carescape B450 with accessories". Registration certificate RZN 2016/3956 date 15.04.2016, p/78, Manufacturer: GE Healthcare Finland Oy, Kuortaneenkatu 2, 00510 Helsinki, Finland
- Sensor for general indication, series 400, pediatric/adult (not more than 10 pc.) accessory to the medical device "Patient monitor Carescape B450 with accessories". Registration certificate RZN 2016/3956 date 15.04.2016, p/78, Manufacturer: GE Healthcare Finland Oy, Kuortaneenkatu 2, 00510 Helsinki, Finland
- Temperature sensors multiple use (not more than 5 pc.), Accessory to the medical device "Module Patient monitor IntelliVue Mx with accessories", p.46, Registration certificate FSZ 2012/13060 date 13.05.2021, manufacturer: Philips Medizin systeme Boblingen Gmbh, Hewlett-Packard-Strasse 2, 71034 Boblingen, Germany
- Temperature sensors single use (not more than 5 pc.), Accessory to the medical device "Module Patient monitor IntelliVue Mx with accessories", p.47, Registration certificate FSZ 2012/13060 date 13.05.2021, manufacturer: Philips Medizin systeme Boblingen Gmbh, Hewlett-Packard-Strasse 2, 71034 Boblingen, Germany
- Temperature sensors YSI 400 series, 5 types, manufacturer "YSI Incorporated", USA, accessory to the medical device "Resuscitation and anesthology Monitor for control of the of the physiology parameters MITAR-01-"R-D" per TU 9441-002-24149103-2002", registration certificate FSR 2012/13343 date 19.02.2021, Manufacturer "NPP Monitor" LLC, 344068, Russia, Rostov-on-Don, Krasnokursantskaya street, 104a

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
1.1 Symbols present on the system, pads and accessories labeling:

③	For the safe and effective use of this device, the operator must consult the accompanying documents prior to use.
EC REP	Indicates the Authorized Representative in the European Community.
1 1	This symbol adjacent to the patient connections means that the thermal probe connection is a "Defibrillator-Proof, Type BF Applied Part", per standard IEC 60601-1 and affords the degree of patient protection defined in that standard for this type of applied part.
c Us Intertek	Per ETL Intertek, models of the Arctic Sun® Temperature Management System that bear the ETL Monogram conform to AAMI ES60601-1, IEC 60601-1-8, IEC 60601-10, IEC 80601-2-35 and are certified to CSA C22.2 No. 60601-1
<u> </u>	Indicates high temperature part or component. The maximum temperature of this internal component allowed by the protective system is also listed.
STERILE STERILE	Indicates that only sterile water should be used when filling the Arctic Sun® Temperature Management System Control Module.
•	Identifies Patient Temperature 1, the patient temperature probe input for monitoring and control.
2	Identifies Patient Temperature 2, the patient temperature probe input for monitoring
	Identifies Patient Temperature Out, the patient temperature output to an external hospital monitor.
₩	Identifies the drain port.
A	Indicates electrical hazard
	Identifies the storage temperature range.
<u></u> %	Identifies the storage relative humidity range.
***	Manufacturer
	Date of Manufacture
2	Do not re-use.
3	Risk of overbalance due to pushing, leaning, resting, etc.
X	ARCTIC SUN® Temperature Management System must be disposed of properly. DO NOT dispose of into the garbage.
	Identifies mechanical hazard

P _{XOnly}	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
á	Unit quantity
REF	Catalog number
SN	Serial number
CE	CE mark
T	Fragile
<u>11</u>	Upward
†	Keep dry
\subseteq	Use by date
(F)	Time of application
2	Do not reuse
8	Do not puncture or cut
MR S	MR-safe
LOT	Lot
NON STERRE	Non-sterile C
	Not made with natural latex
(II) O	Consult instruction for use
MD	Medical device

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