INSTRUCTION FOR USE

1. Remove the protective cap.
2. Open the capsule housing by rotating the mouthpiece in the direction of the arrow.
3. Remove one capsule from the blister right before using it. Do not store capsules in the device.
4. Insert the capsule into the capsule housing. Do not place a capsule directly into the mouthpiece. Do not swallow the capsule. Do not use the device in combination with capsules that are not explicitly intended to be used in combination with this device.
5. Close the mouthpiece by rotating it in the opposite direction, maintaining the device in upright position (mouthpiece pointing up).
6. Maintaining the device upright, firmly press both the side buttons at the same time for capsule piercing.
7. Release the side buttons before inhalation.
8. Exhale before placing the mouthpiece in your mouth.
9. Hold the device without pressing the push buttons and without obstructing the air intakes. Take the mouthpiece between lips and inhale deeply. The capsule moves from the capsule seat to the whirling chamber and starts to spin around. Depending on the delivered drug you may require one or two inhalations in order to empty the capsule. Do not blow into the mouthpiece.
10. Remove the mouthpiece from the lips and hold breath as long as you comfortably can. Then exhale. Do not blow into the mouthpiece.
11. Rotate the mouthpiece in the direction of the arrow.
12. Remove the capsule and visually check if the capsule is punctured and empty. If the capsule is not punctured and still full, start the procedure again from point 4. If the capsule is punctured but not empty, re-insert the capsule into the chamber and only repeat the inhalation step from point 8. Do not press the push buttons again. If the capsule is punctured and empty, dispose it.
13. Close the mouthpiece and replace the protective cap.

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Plastiape S.p.A. con socio unico cap. soc. € 1.000.000 int. vers.