Thopaz⁺™

Cardiothoracic drainage

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19 Technical documentation

WARNINGS

Do not use other accessories than those specified or sold by the manufacturer as replacement parts for internal components as it may result in increased emissions or decreased immunity of the Thopaz+ device.

EMC

Thopaz+ is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 and can be used in the vicinity of other EMC-tested devices that fulfill the requirements of the relevant IEC 60601-1-2 standard. Untested HF (high-frequency) sources, radio networks or the like can influence the operation of the device and may not be operated in combination with the system. Thopaz+ is a medical device that requires special safety precautions and must be installed and placed in operation in accordance with the attached EMC information. Portable and mobile RF communication devices (mobile telephones) can affect Thopaz+.

Electromagnetic emissions		
Thopaz+ is intended for use in the electromagnetic environment specified below. The customer or the user of Thopaz+ should assure that they are used in such an environment.		
Emission Tests	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	Thopaz+ uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference with nearby electronic equipment.
RF emissions CISPR 11	Class B	Thopaz+ is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

WARNINGS

Thopaz+ should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, Thopaz+ should be observed to verify normal operation in the configuration in which it will be used.

Electromagnetic compatibility (EMC, IEC 60601-1-2:2007, Table 2)

Electromagnetic immunity

Thopaz+ is intended for use in the electromagnetic environment specified below. The customer or the user of Thopaz+ should assure that they are used in such an environment.

Immunity Tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 4 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV line(s) to earth	± 1 kV differential mode ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of Thopaz+ requires continued operation during power mains interruptions, it is recommended that Thopaz+ be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels of a typical commercial or hospital environment.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Electromagnetic immunity

Thopaz⁺ is intended for use in the electromagnetic environment specified below. The customer or the user of Thopaz⁺ should assure that Thopaz⁺ is used in such an environment.

Immunity Tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3Vrms 150 kHz to 80 MHz	10Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of Thopaz⁺, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 0.35 \sqrt{P}$ $d = 0.35 \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = 0.7 \sqrt{P}$ <p>800 MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m)</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3V/m 80 MHz to 2.5 GHz	10V/m	

NOTE 1 At 80 MHz and 800 MHz, the high frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Thopaz⁺ are used exceeds the applicable RF compliance level above, Thopaz⁺ should be observed to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as reorienting or relocating Thopaz⁺.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

Recommended separation distance between portable and mobile RF communications equipment and Thopaz+

Thopaz+ is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Thopaz+ can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Thopaz+ as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter	Separation distance according to frequency of transmitter m (meters)		
	150kHz to 80 MHz $d = 0.35 \sqrt{P}$	80 MHz to 800 MHz $d = 0.35 \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$
W			
0.01	0.04	0.04	0.07
0.1	0.11	0.11	0.22
1	0.35	0.35	0.70
10	1.11	1.11	2.21
100	3.50	3.50	7.00

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Appendix B

EN / IEC 62353: 2008

Hospital	Thopaz ⁺ serial number:
	Firmware version:
Location:	Run time (Settings 2/2): h

Visual checks	NO	YES
Cracks, brittle areas, discolouration on housing and display		
Covers for ports are present and can be clicked in		
Standard rail and bed holders are present and can be opened		
Seal is present, is clean and not brittle		
Release button can be pressed, canister can be clicked in		
Thopaz ⁺ is clean with no dirty areas		

Tests	OK	Not OK
FPT test Carry out test as described in Instructions for use. Thopaz ⁺ must be disconnected from the mains (battery operation). Required test equipment: vacuum gauge -0-10 kPa Class 2.5 % or better (E.G.077.1456)	☐	☐
Acoustic signal (beep) Audible?	☐	☐
Battery Connect Thopaz ⁺ to mains adapter or place it in docking station.	☐ Battery is being charged.	☐ Battery is not being charged.

Overall evaluation	Errors	Actions
No safety or functional technical errors were detected.	☐	Thopaz ⁺ is ready for use.
Errors that are not safety relevant were detected. No direct risk.	☐ 1 _____ 2 _____ 3 _____	Thopaz ⁺ is ready for use, errors must be corrected by ____/____/____
Safety relevant errors were detected.	☐ 1 _____ 2 _____ 3 _____	Do not use Thopaz ⁺ until the errors are corrected. Send Thopaz ⁺ to Medela Service Centre.

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

NPWT University

Medela.com website provides all the educational tools you need to set-up and operate the Invia Liberty NPWT system. The Medela University page includes instructions for use, quick cards and training videos to support you in becoming a confident and knowledgeable NPWT user.

www.medela.com

NPWT University can be found at www.medela.com via the "University" tab. It is available on all formats including desktop, smartphone and tablet. We recommend saving NPWT University as a shortcut on desktop, smartphone and tablets for an even quicker access.

Table of contents

Introduction	6
To consider before use	6
Intended use	6
Indications for use	7
Contraindications	7
Warnings, cautions and safety instructions	7
WARNINGS	8
CAUTIONS	9
Safety instructions	10
Wound assessment	10
Safety-related checks	11
Dressing technique	11
Invia Liberty pump	12
Display	13
Invia Liberty pump disposables	14
Invia Liberty canister with solidifier 300ml and 800ml	14
Invia Liberty canister tubing/Y-connector/Drain adapter	15
Power supply	16
USB port	16
Preparation for use	17
Invia Liberty pump user modes	18
Administrative mode	18
Patient mode	18
Administrative Mode	19
Device setup	19
Change pressure level	20
Change pressure set-up mode	21
Select air leakage volume	21
Change therapy mode	23
Change settings	24
Therapy Log file	25
To open the Therapy Log file	25

Patient mode	25
Turn ON	25
Check pressure	25
Air leakage indicator	25
Standby	26
Turn OFF	26
Set up carrying case	26
Change Invia Liberty canister and Invia Liberty tubing	27
Battery charging	29
Alarms	30
“Warning”	31
“Alarm”	31
“Internal fault”	31
Alarm table	32
Accessories overview	36
Wound dressings	37
Sterility and requirements for usage	37
Cleaning and disinfection	38
Disposal	39
Maintenance/Safety-related check	39
Guarantee	39
Service life	39
Signs and symbols	40
Technical specifications	41
Electromagnetic compatibility (EMC)	42
Electromagnetic emissions	42
Electromagnetic immunity	43
Electromagnetic immunity	44
Recommended separation distance	45

Introduction

With the Invia Liberty you have selected a system for use in Negative Pressure Wound Therapy (NPWT). The lightweight Invia Liberty pump provides an adjustable negative pressure range and two therapy modes along with an electronic measuring and monitoring system. The pump is quiet during operation and has optical and acoustic status alarms for patient safety.

Invia Liberty is portable and can be operated independent of the electrical outlet thanks to a rechargeable battery. Acoustic and optical signals are triggered for variances from the set values as well as for faults.

To consider before use

These instructions for use are a general guide for the use of the Invia Liberty pump with associated products. Medical matters must be addressed by a physician.

In order to ensure safe and proper operation of Medela products, a quality management system is used. Please meet the conditions below (failing to do so will void the warranty). Invia Liberty NPWT system is to be used exclusively as described in these instructions for use.

- Before initiating NPWT treatment, read the instructions for use, indications, contraindications, warnings, cautions and safety instructions. Nonobservance and incorrect use can lead to considerable dangers and cause pain and injury to the patient.
- Safe and effective operation of this device requires specific instructions from a physician.
- For use only by persons who have been adequately trained in wound care and negative pressure wound therapy.
- Therapy changes (pressure level, constant or intermittent mode) should only be done as prescribed by a physician.
- In these instructions for use “pressure” in general implies “negative pressure”.
- Please keep in mind that each wound is unique and must be assessed by a qualified medical professional who must use his / her best clinical judgment when applying this therapy. The pressure level and therapy mode must be adapted to each individual patient according to his / her medical knowledge and according to the wound healing phase.

Intended use

Intended user

The Invia Liberty NPWT system is intended to be used by healthcare professionals or adequately trained lay users.

Healthcare professionals are responsible to train lay users according to the patient instructions for use and explain all related safety information.

Intended patient population

The Invia Liberty NPWT system is intended to be used on patients only exhibiting conditions as described in the indications for use. The device has not been studied in pediatric patients.

Intended environment

The Invia Liberty NPWT system is intended for use in acute, extended and home care settings.

Indications for use

The Invia Liberty Negative Pressure Wound Therapy (NPWT) system is indicated for patients who would benefit from a suction device (NPWT) as it creates an environment that promotes wound healing by secondary or tertiary (delayed primary) intention by preparing the wound bed for closure, reducing edema, promoting granulation tissue formation and perfusion, and by removing exudate and infectious material. When used on closed surgical incisions, the Invia Liberty NPWT system is also intended to manage the environment of surgical incisions that continue to drain following sutured or stapled closure by maintaining a closed environment and removing exudate via the application of Negative Pressure Wound Therapy.

The Invia Liberty NPWT system is appropriate for the following indications:

- Acute or subacute wounds
- Chronic wounds
- Dehiscent wounds
- Pressure ulcers
- Diabetic/Neuropathic ulcers
- Venous insufficiency ulcers
- Traumatic wounds
- Partial thickness burns
- Flaps and grafts
- Closed surgical incisions

Contraindications

The Invia Liberty NPWT system is contraindicated in the presence of:

- Necrotic tissue with eschar present
- Untreated osteomyelitis
- Non-enteric and unexplored fistulas
- Malignancy in the wound
- Exposed vasculature
- Exposed nerves
- Exposed anastomotic site of blood vessels or bypasses
- Exposed organs

Warnings, cautions and safety instructions

WARNINGS

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTIONS

Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury.

Safety related tip

Indicating useful information about the safe use of the device.

The Invia Liberty Negative Pressure Wound Therapy system (Invia Liberty NPWT system) is intended for the use described in these instructions for use.

Medela is only responsible for the effect on BASIC SAFETY, reliability and performance of the Invia Liberty NPWT system if it is used in accordance with the instructions for use.

Please read and observe these warnings and safety instructions before operation. These instructions for use must be kept with the device.

Please note that these instructions for use are a general guide for the use of the product. Medical situations must be addressed by a physician.

WARNINGS

- Do not modify this equipment without authorization from the manufacturer.
- This manual provides general guidelines for the use of the Invia NPWT system.
- The safe and effective operation of this device requires specific instruction from a physician.
- No modification of this equipment is allowed.
- For use only by healthcare professionals who have been adequately trained in suction procedures, wound care, negative pressure wound therapy and in the use of aspirators or adequately trained lay users. Healthcare professionals are responsible to train lay users according to the patient instructions for use and explain all related safety information.
Caution: Incorrect use can cause pain and injury to the patient.
- Consult the indications for use, cautions and contraindications when using the Invia Liberty as a vacuum source with the Invia NPWT system. Failure to read and follow all instructions in this manual prior to use may result in death or injury of the patient.
- Failure to obtain consent and any additional instructions from the treating physician prior to use, may lead to death or injury of the patient.
- Before you plug in the device, please check that your local power supply is the same as the voltage given on the specification plate.
- The device is not for use while bathing or showering
- Do not dry Invia Liberty in a microwave.
- Data transfer via USB is not possible in the running mode.
- A patient undergoing NPWT requires frequent supervision. Objective indications or signs of possible infection or complication must be addressed immediately (e.g. fever, pain, redness, increased warmth, swelling or purulent drainage). Monitor the device, wound, surrounding skin and patient status and comfort level frequently to ensure efficient, safe treatment and patient comfort.
- Do not place the foam/gauze dressing directly on exposed blood vessels, organs, nerves, tendons, bones or ligaments. When using the Invia Liberty NPWT system in close proximity to these structures a protective barrier, such as a non-adherent wound contact layer, must be used.
- Serious or fatal injury can result from bone fragments or sharp edges (e.g. staples or hardware) that could puncture protective barriers, vessels or organs.
- Patient must be closely monitored for bleeding. If sudden or increased bleeding is observed, immediately stop use of the pump, apply pressure on wound dressing and seek immediate Emergency Medical Attention.
- Should a spinal cord injury patient experience autonomic hyperreflexia, discontinue treatment with the Invia Liberty NPWT system and consult a physician immediately.
- Never place the Invia Liberty pump in water or liquids. Clamp the drain and disconnect from the dressing prior to bathing or showering.
- Consider the use of a protective barrier on skin that may come in contact with the tubing, especially in patients with fragile skin.
- Invia NPWT instructions advise 24 hours therapy without interruption. If therapy is discontinued for more than 2 hours using foam or gauze, the dressing should be replaced and therapy restarted by a healthcare professional.
- This device has not been studied in pediatric patients.
- Clamp the drain and disconnect the Invia Liberty pump prior to patient entering hyperbaric oxygen chamber (HBO) or Positron Emission Tomography (PET).
- The Invia Liberty NPWT system is not for use in the Magnetic Resonance (MR) environment, so do not take the Invia Liberty NPWT system into this environment.

- Explosion hazard – the Invia Liberty pump is not for use in potentially explosive environments including oxygen enriched environments and in areas of flammable anesthetics.
- In the event that defibrillation is required, disconnect the pump from the wound dressing before the patient is defibrillated.
- Special care is advised for dressing placement and removal in order to avoid situations such as unintentional gauze or foam retention.
- This product contains lithium-ion batteries which bear risk of fire, explosion and burns. Do not disassemble, crash, heat above 100°C, incinerate or dispose of in fire.

Contact your local Medela customer service representative for assistance with product operations.

CAUTIONS

The following statements describe medical conditions that may require special care for the safe and effective use of the Invia NPWT system.

- Patients at high risk for bleeding and hemorrhage.
- Patients taking anticoagulants or platelet aggregation inhibitors, or with patients experiencing active bleeding or difficult wound hemostasis.
- Wounds that involve a fistula.
- Using Invia NPWT system in close proximity to blood vessels and organs or exposed organs, vessels, nerves, tendons, or ligaments. Provide necessary protection of all vessels and organs using a protective barrier.
- Patients with a history of vascular anastomosis or friable, irradiated, sutured or infected blood vessels.
- Use near vagus nerve (bradycardia) or use on patients with a history of spinal cord injury (stimulation of sympathetic nervous system).
- Circumferential dressing application.

Incorrect use may cause pain and injury to the patient. Excessive negative pressure, a too tight adhesive cover dressing, or an infection of the wound may cause pain to the patient. In either case, the dressing must be removed and the wound assessed.

The patient should be monitored regularly according to physician instructions and facility guidelines to monitor patient comfort, therapy compliance and signs of wound infection.

Do not use an Invia Liberty canister or tubing if the sterile packaging is damaged.

The Negative Pressure Wound Therapy must be used 24 hours per day without interruption. If the pump is stopped for more than two hours, the dressing must be changed and therapy restarted.

Consider the patient's size and weight when prescribing this device.

Consider mode of therapy – intermittent versus continuous.

Safety instructions

- Before you plug in the device, please check that your local power supply is the same as the voltage given on the specification plate.
- The Invia Liberty pump is verified within the scope of conformity evaluation and is only to be used with products included in the Invia Liberty NPWT system and distributed by Medela. Medela can only guarantee the effective performance of the system with these products.
- Wireless communication equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, and walkie-talkies can interfere with the Invia Liberty pump and should be kept at a minimum distance of 30 cm away from the Invia Liberty pump.
- The patient should be regularly monitored according to facility or institution guidelines.
- Invia Liberty pump must remain in an upright position during use.
- Supervision is necessary when the Liberty pump is used in the vicinity of children.
- Do not use Invia Liberty pump if:
 - The power cord or plug are damaged
 - The device is not functioning properly
 - The device is damaged
 - The device has apparent safety defects
- Never pull the plug out of the main socket by pulling on the connecting cable.
- Keep the Invia Liberty pump with associated products away from hot surfaces.
- Never place the Invia Liberty pump, charger or docking station device in water or other liquids and keep the charger connector away from moisture or immersion in water.
- The Invia Liberty pump must not be used for suctioning explosive, easily flammable or corrosive liquids.
- The tubing connected to the canister must never come in direct contact with the wound area.

Wound assessment

CAUTION

Patient monitoring: The patient should be monitored regularly according to the physician's instructions and facility guidelines to check for patient comfort, therapy compliance and signs of infection.

WARNING

Objective indications or signs of a possible infection or complication must be addressed immediately (e.g. fever, pain, redness, increased warmth, swelling or purulent discharge). Non-observance can lead to considerable danger to the patient.

Observe wound/periwound tissue and exudate for signs of infection or other complications. Most common signs of infection include redness, tenderness, fever, swelling, itching, increased warmth in the wound area, strong odor or purulent discharge. Additional symptoms include nausea, vomiting, diarrhea, headache, dizziness, fainting, sore throat with swelling of the mucous membranes, disorientation, high fever (>38.8° C), refractory hypotension, orthostatic hypotension, or erythroedema (a sunburn-like rash). More serious complications of infection include pain, discomfort, fever, gangrene, toxic or septic shock. If more serious complications of infection occur, discontinue therapy and consult a physician immediately.

Safety-related checks

For the safety-related checks, the device should be maintained and repaired throughout its service life in compliance with the service procedures.

The Invia Liberty pump is a device in protection class II (EN IEC 60601-1), the safety-related checks are confined to visual inspection of the housing and charger for damage. This check must be performed prior to each use.

Devices of protection class II do not have a protective earth conductor; there is therefore no need to check the earth leakage current.

The Invia Liberty pump enclosure is made entirely of insulated material. Tests of the enclosure leakage current using common measuring instruments will therefore not reveal measurable values.

Even when suctioning a conductive fluid until the overflow protection device activates, measurements of the patient leakage current using common measuring instruments will not reveal measurable values.

Invia Liberty pump does not have patient circuits or functional earth connections.

Dressing technique

Consult the appropriate Invia wound dressing instructions for use for information regarding dressing applications. Perform a thorough wound cleansing per physician orders or facility protocol prior to dressing applications.

Routine dressing changes should occur every 48–72 hours. Dressing changes for infected wounds should be considered more frequently.

Dressing change



 Dressings (wound contact layer, wound filler, wound cover, internal drain, drain adapter, FitPad) should be changed every 48 to 72 hours, but no less than 3 times a week, or as instructed by the healthcare professional.

When treating infected wounds or wounds more susceptible to tissue in-growth into wound filler material, more frequent dressing changes may be needed. The frequency of dressing changes should be based on an evaluation of the wound characteristics rather than standard recommendations.

 Important things to remember:


- Routinely check that the negative pressure level and the therapy mode is at the prescribed setting.
- Monitor the integrity of the dressing and that the pump is on and running.

Invia Liberty pump


The max flow of Invia Liberty pump is 5 liters/minute with an adjustable pressure range of -40 to -200 mmHg (-5 to -27 kPa).

Display

Invia Liberty pump disposables

Invia Liberty canister with solidifier 300 ml and 800 ml

Material:	Polypropylene
Accuracy of graduation:	+/- 2.5 % (in the upright position)
Composition of solidifier:	Cross-linked sodium polyacrylate


CAUTIONS

For appropriate, safe operation the Invia Liberty pump must remain in an upright position during use.

If the pump tips over, set it upright again. The special construction of the safety chamber in the upper region of the canister protects the overflow protection / bacteria filter from clogging immediately, if tipped over.

When the canister is full and the pump tips over, this function is made inoperative, since the secretions will flow into the safety chamber and clog up the hydrophilic filter. In this case an alarm will sound and the canister will need to be replaced.

Invia Liberty canister tubing/Y-connector/Drain adapter

Material: TPE
 Length: 1.5 m
 Diameter: 2.0/3.0 mm

The canister tubing is comprised of two lumens:

The smaller lumen (measuring tubing) regulates the pressure and the larger lumen (suction tubing) removes the fluid from the wound into the canister. A hydrophilic overflow protection/bacteria filter in the tubing base helps to prevent contamination of Invia Liberty pump.

The Invia Y-connector with Quick-connector

For warnings and cautions of the Invia Y-connector with Quick-connector consult the Invia Y-connector with Quick-connector instructions for use.

Power supply

CAUTIONS

Before you charge the device, please check that your local power supply is the same as the voltage given on the specification plate.

WARNINGS

This product contains lithium-ion batteries which bear risk of fire, explosion and burns. Do not disassemble, crash, heat above 100 °C (212 °F), incinerate or dispose of in fire.

Invia Liberty pump is operable while connected to the electrical power supply or by an internal rechargeable lithium-ion battery. While in use and connected to the electrical power supply, the battery is re-charged.

The charge on the battery is dependent upon the run-time of the pump. The run time refers to the effective operation of the motor. Invia Liberty pump motor only turns on when the measured pressure is lower than the set pressure. If the Invia Liberty pump is running continuously, Medela suggests a maximum of 4 hours of battery usage before it is necessary to recharge.

The average run time of the battery exceeds 10 hours. This is influenced by the size of the wound, the air leakage in the system and the set pressure. If there is a leakage in the system, the pump motor will run more often which will reduce the battery run-time.

USB port

Additional equipment connected to medical electrical equipment must comply with the respective IEC or ISO standards (e.g. IEC 60950 for data processing equipment).

Furthermore all configurations shall comply with the requirements for medical electrical systems (see IEC 60601-1-1 or clause 16 of the 3Ed. of IEC 60601-1, respectively). Anyone connecting additional equipment to medical electrical equipment configures a medical system and is therefore responsible that the system complies with the requirements for medical electrical systems. Attention is drawn to the fact that local laws take priority over the above mentioned requirements. If in doubt, consult your local representative or the Medela customer service department.


Preparation for use

WARNING

Use only after instruction by trained personnel. Wear gloves for all operations and utilize universal precautions.

1. Check necessary parts

2. Connect Invia Liberty canister tubing

A. Open the tubing packaging


B. Insert the tubing base into the pump as shown (straight push)

3. Click in Invia Liberty canister

A. Unpack new canister and remove the safety sticker. Slip the pegs at the bottom of the canister into the slots at the base of the pump. Lift the canister towards the pump until it clicks into place.



B.

B. "Click" into place.


Canister secures tubing

4. Connect the dressing tubing to the canister tubing.

5. Switch Invia Liberty pump ON by pressing []

A.

A. The pump will do a self-test and start running immediately. You should hear the motor running for a short time until the pressure is built up. If the motor keeps running for more than 30sec. – check system and retry.

If the above check was successful, switch OFF Invia Liberty pump by pressing [ > 3 seconds]. If the self-test is not successful, follow instructions shown on the display or see the Alarms chapter.

For dressing applications, please refer to the instructions for use provided with the Invia dressings. Turn the Invia Liberty pump on as described in the user modes section.

Invia Liberty pump user modes

Administrative mode

Used by healthcare professional to either set up new patient therapy or to change pump settings such as pressure, air leakage volume, therapy Constant and Intermittent modes, intermittent times, language and time zone. You can enter the Administrative mode when first turning the pump on or while therapy is running.

Patient mode

In Patient mode the pump can be turned on and off, placed in standby mode and an alarm can be muted. In this mode the canister and tubing can be changed and the pump can be charged. When the pump is turned off in this mode, last settings are used by default.

Administrative Mode

 The factory default settings of the Invia Liberty pump are -125 mmHg and constant mode.


CAUTION

The pressure level should always be set according to prescribing healthcare professional instruction.

When an Invia FitPad is used, the set pressure is controlled at the wound site.

Device setup


New patient therapy

1. **Switch Invia Liberty pump ON** in Administrative mode, press and hold [], and then press [].

Self-test starts

If the self-test is not successful, follow instructions shown on the display or see the Alarm chapter.

2. **Acknowledge disclaimer.** Press “OK” [] to confirm.

“New patient? Yes/No”

“Yes” means that Invia Liberty pump will issue a new therapy ID number (= new patient).

This number can be noted in the patient's file.

The therapy ID number is displayed in the information field.

If yes


Press “Yes” [] to confirm.

Press “OK” [] to enter the main display.

“No” means that the therapy ID number and settings remains unchanged (= same patient).


If no

Press “No” [] to confirm and enter the main display.


3. Press the Selection buttons [] or [] to set pressure level.
4. Press “OK” [] to confirm and enter the main display.

Administrative mode

1. To enter the Administrative mode, unlock the display, press and hold [], and then press [].

Change pressure level

CAUTIONS

Set the pressure level and therapy mode according to prescribing healthcare professional instruction.

Default pressure mode is Standard with preset pressure levels:

mmHg	-40	-60	-80	-100	-125	-150	-175	-200
(kPa)	-5	-8	-11	-13	-17	-20	-23	-27

1. To enter the Administrative mode, unlock the display, press and hold [], and then press [].
2. Press the Selection buttons [] or [] to set pressure level.
3. Press “OK” [] to confirm.
4. To return to the main display, press “Back” [].

CAUTIONS

If the pressure level is not confirmed, the pump will switch back to previous settings and return automatically to the main display after 5 seconds.

Change pressure set-up mode

Detailed pressure mode allows to select pressure level from -40 to -200 mmHg in increments of 5mmHg (5kPa to 27kPa in increments of 1kPa)


1. To enter the Administrative mode, unlock the display, press and hold [▶], and then press [⏻].
2. Press the soft buttons [▶] + [◀] at the same time to enter the Settings menu.
3. To switch pressure modes, scroll down with [▼] and select Pressure set-up.
4. Press "OK" [◀] to select.
5. To change between Standard and Detailed pressure modes, press Selection buttons [▲] or [▼] to choose Standard or Detailed mode.
6. Press "OK" [▶] to confirm.
7. Refer to Change pressure level instructions to select pressure values.

Select air leakage volume

High air leakage mode is recommended when X-large dressing is used.




Default air leakage volume is Standard. If the vacuum on the wound is < 75% of the set pressure level or the air flow is > 1500ml/min, leakage indicator is "full" and flashes, indicating there is a big air leak in the system. The air leakage alarm will go off within 2 minutes if the air flow is above the defined level or after 1 minute if the set pressure is not maintained. Follow the instructions shown on the display, or see the Alarms chapter.

For wounds with air leakage volume between 1500ml/min and 2200ml/min, High air leakage volume mode is available. If the vacuum on the wound is < 75% of the set pressure level or the air flow is > 2200ml/min value, leakage indicator is "full" and flashes, indicating there is a big air leak in the system. The air leakage alarm will go off within 2 minutes if the air flow is above the defined level (2200ml/min) or after 1 minute if the set pressure is not maintained. There will be no leakage alarm during the first 5 minutes which allows establishing the set pressure on a wound with high volume without an alarm.


Selecting between Standard and High air leakage volume:

1. To enter the Administrative mode, unlock the display, press and hold [], and press [].
2. Press the soft buttons [ + ] at the same time to enter the Settings menu.
3. To change air leakage volume, scroll down with [] and select Air Leakage.
4. Press “OK” [] to select.
5. To change between Standard or High air leakage, choose the mode with the Selection buttons [] or [].

6. Press “OK” [] to confirm.
7. To return to the main display, press “Back” [].


Change therapy mode

1. To enter the Administrative mode, unlock the display, press and hold [▶], and then press [⏻].
2. To change the therapy mode, press "Standby" [◀ 3 seconds] to put the Invia Liberty pump into the Standby mode.

Constant mode – C


The default therapy mode is Constant.


3. To change from Constant mode to Intermittent mode, press [▶] "Change to Intermittent" and then press "On" [▶].

Intermittent mode – I

The default times for Intermittent mode are 5 minutes on and 2 minutes off.

4. To change from Intermittent mode to Constant mode, press [▶] "Change to Constant" and press "On" [▶].


- ☞ 1 minute after the last button has been pressed, the Invia Liberty pump switches into Patient mode automatically and the display is locked.
- ☞ 1 minute after the last button has been pressed, the backlight turns off. When pump alarms or any button is pressed the backlight will illuminate.

☞ The settings can only be changed when the pump is in the Administrative mode.

The settings possible to change are, Unit pressure, Air leakage volume, Constant and Intermittent modes, intermittent times, language and time zone. Pump number, pump run-time and pump version can only be viewed but not changed.

Change settings


- 1. To enter the Administrative mode, unlock the display, press and hold [▶], and then press [⏻].
- 2. Press the soft buttons [▶+◀] at the same time to enter the Settings menu.
- 3. To change a setting, choose with the Selection buttons [^] or [v].

Settings 1/3		Selection
Pressure	-125	-40 to -200mmHg (-5 to -27 kPa)
Unit Pressure	mmHg	kPa / mmHg
Pressure set-up	Standard	Standard / Detailed
Air Leakage	Standard	Standard / High
Language	English	Languages according to list
Back	OK	

Settings 2/3		Selection
On time	5 min	1 – 8 min
Off time	2 min	1 – 8 min
Time zone	+1 h	+12 until -12 hours GMT
Pump number	152887	as displayed
Pump run-time	27 h	as displayed
Back	OK	

Settings 3/3		Selection
Version	1.22	as displayed
Back	OK	


- 4. Press “OK” [▶] to select.
- 5. To change value, choose with the Selection buttons [^] or [v].
- 6. Press “OK” [▶] to confirm.

☞ If the changed setting is not confirmed with “OK”, the pump will switch back to the previous setting.

- 7. Press “Back” [▶] to exit settings.

☞ 1 minute after the last button has been pressed, the Invia Liberty pump switches into Patient mode automatically and the display is locked.

Therapy Log file


In the Therapy Log file information regarding runtimes (on/off), pressure settings, therapy modes, alarms, warnings and errors is listed. The last 51 events are logged and displayed.

To open the Therapy Log file

Date	Time	M	Duration
2011-06-02	06:41	Standby	0:01
2011-06-02	06:15	I (-400)	0:26
2011-06-02	06:15	Standby	0:00
2011-06-01	15:00	C (-130)	15:15
2011-06-01	14:48	Off	0:12
2011-06-01	13:40	C (-130)	1:00

Back -1- Graph

List

Graph

1. To enter the Administrative mode, unlock the display, press and hold [▶], and then press [⏻].
2. Press Selection buttons [▲] + [▼] simultaneously.
3. To view additional pages, scroll with the Selection buttons [▲] or [▼].
4. To view the Therapy Log file as a graph, press “Graph” [◀].
5. To exit the Therapy Log file, press “Back” [▶].

☞ The Invia Liberty pump switches automatically into main display (in Administrative mode) 30 seconds after the last button has been pressed. After additional 30 seconds, the pump switches into Patient mode and the display is locked.

Patient mode

Turn ON

Switch Invia Liberty pump ON by pressing [⏻]


Self-test starts

If the self-test is not successful, follow instructions shown on the display or see the Alarms chapter.

☞ When the pump is turned on in this mode, last settings (therapy mode and pressure level) are used by default.

Check pressure


The set pressure will be shown on the display. The motor will run for a few seconds to build up the pressure. If it runs continuously for more than 30 seconds, check system for leaks.


Air leakage indicator

An air leakage indicator is shown on the display to visualize if there is an air leakage in the system.

a) When the indicator is “empty”, the system is air tight.

b) When the indicator is "half full" there is an air leak in the system, but the pressure and therapy is maintained in accordance with the set pressure.

c) When the indicator is "full" and flashes, there is a big air leak in the system. When pump is first turned on and a big air leak is present the air leakage alarm will go off within 2 minutes in the Standard air leakage mode or within 5 minutes in the High air leakage mode if the set pressure is not maintained. When therapy is running and a big air leak is detected, alarm will go off within 2 minutes. Follow the instructions shown on the display, or see the Alarms chapter.

Standby

Change the Invia Liberty pump into the Standby mode.

Press "Standby" [ > 3 seconds] and the pump motor will stop running.

If the pump is in Standby mode for more than 5 minutes, a warning will go off, follow instructions shown on the display or see the Alarms chapter.

Turn OFF

Press [ > 3 seconds] and the pump will turn off.

Set up carrying case

The Invia Liberty is intended to be worn in a carrying case during operation.

1. Open the carrying case.
2. Slide the Invia Liberty into the designated pocket.
3. Close the velcro over the pump.

The shoulder strap can also be used as a belt strap.

To accomplish this, release the strap at the sides of the carrying case and pull it through the two loops located on the back of the carrying case.

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


Change Invia Liberty canister and Invia Liberty tubing minimum once a week or, when the canister is full, by visual inspection or when alarm sounds, in accordance with instructions on the pump display.


1. Provide sterile canister and sterile tubing.
2. Press "Standby" [ > 3 seconds] for the pump motor to stop running. Clamp the canister tubing.
- 3.1 Position the clamp next to the quick-connector and close the clamp on the pump tubing and close the clamp on the dressing tubing.

- 3.2 Disconnect the canister tubing from dressing tubing.
- 3.3 Protect the end of dressing tubing from contamination.


4. Release and remove canister.

5. Seal used canister with a cap located on the side of the canister.

6. Remove canister tubing in direction of the arrow.


7. Unpack new canister tubing and connect to Invia Liberty pump. Insert the tubing base into the pump (straight push).

- 8.1 Unpack new canister and remove the safety sticker.
 - 8.2 Slip the pegs at the bottom of the canister into the slots at the base of the pump. Lift the canister towards pump until it clicks into place.
 - 8.3 Connect the canister tubing to dressing tubing. **Ensure that the tubing is un-clamped.**
 - 8.4 Press "ON" with []. Pressure will begin to built up.
9. Dispose canister and canister tubing in accordance with local procedures. In the home care setting, return disposables to care giver for correct disposal.

Battery charging

The battery must be fully charged before first usage. To charge the battery, plug the charger into the electrical outlet port in the Invia Liberty pump. You can continue to use the pump while it is charging.

The capacity of the battery is not negatively affected if it is charged when partially discharged and the pump does not need to be fully charged at each occasion.

The battery can be charged when appropriate, or in accordance with alarm signal for “Battery low” or “Battery fully discharged (empty)” and instructions in the pump display.

There is approximately 30 minutes of charge remaining on the battery at the onset of the “Battery low” warning. If the battery is completely discharged it will take 3 to 4 hours to charge to full status.

If the Invia Liberty pump is fully charged AND the pump is still connected to an electrical source,  will appear in the battery icon.

If the pump is disconnected from an electrical source, 4 bars are visible in the battery icon which indicates that the battery is fully charged.

If Invia Liberty pump has not been in use, the battery must be charged approximately once every 6 months to ensure optimum function.

To disconnect Invia Liberty pump from the electrical source, remove the charger connector and close the safety cover.


WARNING

Do not pull on the cable or the anti-bend protection.

Alarms

The pump distinguishes between “Warning”, “Alarm” and “Internal fault”. If the Invia Liberty pump detects any situations where the therapy cannot be maintained, an acoustic alarm sounds, a fault number and a description of the problem appears on the display. For explanation of the fault number, see the Alarm Table in this chapter.

Example:

CAUTIONS

Invia Liberty instructions advise 24 hours therapy without interruption. If therapy is discontinued for more than 2 hours using foam or gauze, the dressing should be replaced and therapy restarted by a healthcare professional.

“Warning”

An acoustic alarm sounds and the fault number is shown on the display.

**“Alarm”**

An acoustic alarm sounds and the fault number is shown on the display.

When a Warning/Alarm goes off, an acoustic alarm sounds. A description of the “Warning” or “Alarm” will be shown on the display.

1. Press “Mute” [] to mute and acknowledge the alarm.
The acoustic alarm will resume in 1 minute if the problem is not solved.
2. Follow the instruction shown on the display or see the Alarm table.
3. If the problem cannot be solved, turn off the Invia Liberty pump [> 3 seconds] and contact the Medela customer service for further instructions.

“Internal fault”

Pump operation stops and an acoustic alarm sounds, “Internal fault” is shown on the display.











1. Press [> 3 seconds] and the pump will be turned off.
2. Restart the pump by pressing [] and the pump will be turned on.
3. If the Internal fault alarm remains, turn off [> 3 seconds] Invia Liberty pump and contact the Medela customer service.

Alarm table

	Fault number	Problem description on the display	Troubleshooting on the display	Pressure	
Warning	401			✓	
	402				
	405			✗	
	406			✓	
Alarm	301			✓	
	302			✓	
	305			✗	
	306			✓	

	Remarks/potential cause of fault
	<p>Recharge the battery either by placing the Invia Liberty pump in the docking station or plug in the charger to the electrical outlet port on the pump.</p> <p>Remaining time of battery is approximately 30 minutes.</p>
	Unplug USB cable
	<p>If the pump is in Standby mode for more than 5 minutes, an alarm will go off.</p> <p>To continue therapy press "On" [◀] or switch off the pump by pressing [⏻ > 3 seconds].</p>
	Cool Invia Liberty pump down as per instructions on the display.
	<p>Dressing: - Check dressing for air leakage. Press firmly around the edges of the dressing, around the drain tubing or on the Invia FitPad.</p> <p>- Apply some additional film dressing to seal the leaking area.</p> <p>Connectors: - Ensure that the dressing tubing is connected properly to the canister tubing.</p> <p>- Ensure that the canister tubing is inserted straight into the pump.</p> <p>Canister: - Ensure that the canister is properly inserted, release the canister and reposition.</p> <p>- Ensure that the orange O-ring/gasket, placed beside the canister tubing on the pump is not missing. Additional O-ring is available via Medela customer service.</p>
	<p>Tubing: - Ensure that the tubing is not twisted, kinked or clamped.</p> <p>- If the canister tubing is clogged, change the tubing.</p> <p>Canister: - If canister is full or filter clogged, replace canister.</p> <p>- If dressing tubing is clogged, change the dressing.</p> <p>Then press "On" to restart the therapy.</p>
	<p>Recharge the battery either by placing the Invia Liberty pump in the docking station or plug in the charger to the electrical outlet port on the pump.</p> <p>Remaining time of battery is 15 minutes.</p>
	Change the canister, see chapter "Change Invia Liberty canister and Invia Liberty tubing".

Alarm table

	Fault number	Problem description on the display	Troubleshooting on the display	Pressure	
Alarm	311			X	
	312			X	
	313			X	
	315			X	
Internal fault				X	

CAUTION
If fault repeats, note the fault number, switch off the pump and contact the Medela customer service.

CAUTION
Invia Liberty instructions advise 24 hours therapy without interruption. If therapy is discontinued for more than 2 hours using foam or gauze, the dressing should be replaced and therapy restarted by a healthcare professional.

	Remarks/potential cause of fault
	<p>Canister: - Release the canister and reposition.</p> <p>Tubing: - Ensure that the tubing is not twisted, kinked or clamped.</p>
	<p>Occurs when alarm 301 was pending for 5 minutes. For troubleshooting, refer to instructions for Alarm 301. Restart the pump by pressing “On” [].</p>
	<p>To restart the therapy after canister has been changed, press right selection button “On” [].</p>
	<p>Occurs when Warning 406 was pending for 30 minutes. Cool Invia Liberty pump down.</p>
	<p>Restart the pump. If internal fault remains, turn off by pressing [ > 3 seconds] and contact the Medela customer service.</p>

OPERATION HALTS


OPERATION CONTINUES

Accessories overview

CAUTION

The Invia Liberty pump is verified within the scope of conformity evaluation and is only to be used with products included in the Invia Liberty NPWT system and distributed by Medela. Medela can only guarantee performance of the system with these products.

* Product not available in all countries

Wound dressings

 Wound dressings to be applied and changed by healthcare professionals only.

The Invia Liberty NPWT system is intended to be used in conjunction with the Invia dressings only. For specific dressing indications, contraindications, warnings and cautions, consult the appropriate Invia dressing instruction for use.

Sterility and requirements for usage

The Invia Liberty disposables products are sterile and single use devices.

Canister tubing Canister Y-connector Drain adapter and drain Invia FitPad (tubing) Wound dressing		These are single use devices that should be disposed of after use. If reused performance of the product may deteriorate, cross-contamination may occur.
--	---	---

 **CAUTION**
 Do not re-sterilise. Do not use if inner package is damaged or opened prior to use.
 Do not reuse. If reused performance of the product may deteriorate, cross-contamination may occur.

Cleaning and disinfection

CAUTION

Invia Liberty pump with associated products (docking station, rail holder and charger) should be cleaned/disinfected after every use. Before cleaning the device, unplug the pump from the wall outlet.

	Cleaning	Disinfection	Sterilization	Dishwasher
Pump	✓	✓	✗	✗
Docking station				
Rail holder				
Charger	Wipe off with a damp cloth.	Wipe off with a disinfecting agent	Sterilization / cleaning in a dishwasher is NOT allowed.	
Carrying case*				

* If needed, the carrying case can be washed in washing machine

CAUTION

Invia Liberty pump, docking station, rail holder, charger and carrying case cannot be sterilized. Immersion disinfection, thermal disinfection and ultrasound cleaning are not permitted.

Disinfection (Pump housing, docking station and charger)

CAUTION

Invia Liberty pump can be disinfected with “alcohol”.

Disinfection

Invia Liberty can be disinfected with the disinfecting agent group “alcohol”.
Do not use other cleaning agents (e.g. Terralin) as they can damage the plastic housing.
Immersion disinfection, thermal disinfection and ultrasound cleaning are not permitted.

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CAUTION

Do not use other cleaning agents (e.g. Terralin) as they can damage the plastic material.

Cleaning procedure for Invia Liberty pump, docking station and charger

1. Wear suitable protection (clothing, gloves, face mask and goggles) according to local guidelines.
2. Apply the disinfectant agent in accordance with instructions from the manufacturer¹⁾.
Pay attention particularly to edges, narrow corners and bottom side.
3. Leave the disinfectant on. Follow the recommended residence time for the disinfectant as instructed by the manufacturer¹⁾.
4. Thoroughly clean the surface, all edges, housing niches, corners, bed holder, brackets, port covers and bottom side.

5. Wipe dry or air dry as instructed from the manufacturer¹⁾.
6. If needed, repeat step 2-5 to ensure proper cleaning.
7. Dispose contaminated material in accordance with local environmental guidelines.

CAUTION

Do not spray disinfectants directly into openings as this may harm electronic components.

For a detailed cleaning instruction, contact your Medela representative.

¹⁾ Manufacturer of the disinfectant agent.

Disposal

Invia Liberty pump is made from various metals and plastics. Before disposal the rechargeable battery and electronics must be removed according to instructions. Then Invia Liberty pump is no longer operational. Disposal of electronics and plastic components should be handled in accordance with local environmental guidelines.

Invia Liberty disposables and dressings should be handled and disposed in accordance with local environmental guidelines.

CAUTION

Must not be disposed of together with household refuse.

Maintenance / Safety-related check

Service work may only be carried out by authorized personnel. A safety-related check is confined to visual inspection of the housing and charger for damage and must be performed prior to each use. If Invia Liberty pump has not been in use, the battery must be charged approximately once every 6 months to ensure optimum function.

Guarantee

Guarantee for 2 years after date of delivery in used in accordance with these instructions. The manufacturer is not liable for any damage or consequential damage caused by incorrect operation, inappropriate usage as well as use by unauthorized persons.

Service life

The service life of the device is five years; the internal batteries life included.

Signs and symbols

This symbol indicates a safety related tip.

This symbol indicates a CAUTION or WARNING associated with the device.

This symbol indicates a class II device.

This symbol indicates to not dispose the device together with household refuse (for EU only).

This symbol indicates the date of manufacture (four digits for the year and two digits for the month).

This symbol indicates the name and the address of the manufacturer.

This symbol indicates the device is sterilized using ethylene oxide.

This symbol indicates a prescription device. CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician. (for US only).

This symbol indicates a type BF applied part.

This symbol indicates manufacturer's catalog number.

This symbol indicates manufacturer's serial number.

This symbol indicates manufacturer's batch code.

This symbol indicates the protection against the ingress of solid foreign objects and against harmful effects due to the ingress of water.

This symbol indicates that the device should not be used after the end of the year and month shown.

This symbol indicates to follow the instruction for use.

This symbol indicates the compliance with additional USA and Canada safety requirements for medical electrical equipment.

This symbol indicates to not use the device if package is damaged.

This symbol indicates the number of items.

This symbol indicates the direct current socket.

This symbol indicates a single use device. Do not reuse the device.

This symbol indicates MR unsafe.

This symbol indicates the temperature limitation for operation, transport and storage.

This symbol indicates the humidity limitation for operation and storage.

This symbol indicates the atmospheric pressure limitation for operation and storage.

This symbol indicates to keep the device dry.

This symbol indicates to handle the fragile device with care.

This symbol indicates to keep the device away from sunlight.

This symbol indicates that the device is in conformance with the Medical Device Directive 93/42/EEC.

This symbol indicates that the marked item or its material is part of a recovery or recycling process.

This symbol indicates the correct upright position of the transport package.

Technical specifications

vacuum range
- 40 to -200 mmHg
- 5 to -27 kPa

low flow
5 L/min

without canister
1000 g
2.2 lbs

IP33

H x W x D
150 x 170 x 95 mm

ISO 9001
ISO 13485
CE (93/42/EEC), IIa

max. noise level
- 42.5 dB(A) 1 l

alarm noise level
78 dB(A)

Operating
conditions

15% +5

°C

70 kPa

Transport /
Storage conditions

15% -20

°C

70 kPa

Switching adapter AC

Model: MSA-C2500IS12.0-30C-ZZ

IEC: 60601-1

Input: 100-240VAC, max. 0.8A,
50/60Hz

Output: 12.0VDC, 2.5A

Electromagnetic compatibility (EMC)

Invia Liberty is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 and IEC 60601-1-2:2014 4th Edition according clause 7 and 8.9.

Electromagnetic compatibility (EMC, IEC 60601-1-2:2007 3rd Edition and IEC 60601-1-2:2014 4th Edition, Table 1)

Electromagnetic emissions		
The Invia Liberty pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Invia Liberty pump should assure that it is used in such environment.		
Emission tests	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Invia Liberty pump uses RF energy only for it's internal function. Therefore, its RF emissions are very low and are not likely to cause any interference with nearby electronic equipment.
RF emissions CISPR 11	Class B	The Invia Liberty pump is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

**WARNING**

The Invia Liberty pump should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Invia Liberty pump should be observed to verify normal operation in the configuration in which it will be used.

**WARNING**

Invia Liberty is a medical device that requires special safety precautions and must be installed and placed in operation in accordance with the attached EMC information.

**WARNING**

Do not use other accessories than those specified or sold by the manufacturer as replacement parts for internal components as it may result in increased emissions or decreased immunity of the Invia Liberty pump.

**WARNING**

EMC

HF (high-frequency) surgical equipment, radio networks or the like can influence the operation of the device and may not be operated in combination with the Invia Liberty pump.

Electromagnetic compatibility (EMC, IEC 60601-1-2:2007 3rd Edition, Table 2)

Electromagnetic immunity			
The Invia Liberty pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Invia Liberty pump should assure that it is used in such environment.			
Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV line-to-earth	± 1 kV differential mode ± 2 kV line-to-earth not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Invia Liberty pump requires continued operation during power mains interruptions, it is recommended that the Invia Liberty pump be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels of a typical commercial or hospital environment.
NOTE: U_T is the a.c. mains voltage prior to application of the test level.			

Electromagnetic compatibility (EMC, IEC 60601-1-2:2007 3rd Edition, Table 4)

Electromagnetic immunity			
The Invia Liberty pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Invia Liberty pump should assure that it is used in such environment.			
Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Invia Liberty pump including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 0.35 \sqrt{P}$ $d = 0.35 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 0.7 \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	<p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a, should be less than the compliance level in each frequency range. ^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<p>^a Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Invia Liberty pump is used exceeds the applicable RF compliance level above, the Invia Liberty pump should be observed to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as reorienting or relocating the Invia Liberty pump.</p> <p>^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.</p>			

Electromagnetic compatibility (EMC, IEC 60601-1-2:2007 3rd Edition, Table 6)

Recommended separation distance between portable and mobile RF communications equipment and the Invia Liberty pump			
The Invia Liberty pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Invia Liberty pump can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Invia Liberty pump as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter	Separation distance according to frequency of transmitter		
	m		
W	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 0.35 \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$
0.01	0.035	0.035	0.07
0.1	0.11	0.11	0.22
1	0.35	0.35	0.7
10	1.11	1.11	2.21
100	3.5	3.5	7.0
For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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Table of contents

Introduction	4
Overview.....	4
Intended use	6
Indications for use	6
Contraindications.....	6
Warnings and safety instructions	7
Warnings	7
Cautions	8
Physician's orders.....	9
Wound assessment	9
Dressing technique	10
Initial delivery	11
Invia Motion NPWT pump accessories and associated products.....	11
Invia Motion disposables	12
Wound dressings.....	12
Charge battery for the first time	13
Disconnect from power source.....	13
Positioning	13
Preparation for use	14
Operating instructions.....	17
Power supply and battery management	19
Notifications and troubleshooting	21
Overview notifications	24
Cleaning and disinfection.....	25
Disposal	26
Warranty and maintenance	27
Signs and symbols.....	28
Technical specifications.....	30
Electromagnetic compatibility (EMC)	31

Introduction

The Invia Motion Negative Pressure Wound Therapy system, subsequently named “Invia Motion NPWT system”, is approved exclusively for the use described in these Instructions for use. Please read the information contained within these Instructions for use before operation of the Invia Motion NPWT system. These instructions for use must be kept with the device.

These instructions for use are a general guide for the use of the product. Medical situations must be addressed by a physician.

Indications, Contraindications, Warnings, Cautions and other information to support the safe use of the Invia Motion NPWT system are included in the respective chapters in these instructions for use.

The Invia Motion NPWT pump is a suction pump for Negative Pressure Wound Therapy that provides treatment status through a display and acoustic signals. The Invia Motion NPWT pump is a single patient use pump which provides continuous or intermittent operation and multiple negative pressure selection options. The Invia Motion NPWT pump is portable and can be operated using a rechargeable battery. Acoustic and optical signals are triggered for variances from the set values as well as for faults. The Invia Motion NPWT system is intended for use in all healthcare settings, including home care, by healthcare personnel and caregivers trained on the instructions for use. The patient or caregiver must have the visual and hearing acuity necessary to appropriately respond to notifications with the sensory and cognitive ability to understand directions from a healthcare professional pertaining to the proper use of the pump.


Invia Motion Negative Pressure Wound Therapy system is available in two versions:

- 1) **REF** 087.4022: Invia Motion – 15 Days Negative Pressure Wound Therapy system
- 2) **REF** 087.4024: Invia Motion – 60 Days Negative Pressure Wound Therapy system

Overview

The Invia Motion NPWT system and the Invia Motion carrying case are for a single patient’s entire therapy.

Buttons to operate the Invia Motion NPWT pump

Backlight

The backlight will be activated by pressing any button. The backlight will turn off after 10 seconds if no further interaction results.

Symbols on display

Notifications	
	Leakage
	Blockage
	Canister full
	Battery low (symbol is blinking)
	Battery fully discharged (symbol is blinking)
	Battery missing (symbol is blinking)
	Internal fault
	Pump time expired (symbol is blinking) (Example for Invia Motion – 15 Days)
Information	
	Shows runtime already used in days (📅) and hours (🕒) (Example for 6 days and 23 hours)
	Pump is running / On-state in intermittent mode (symbol is blinking)
	Off-state in intermittent mode (symbol is blinking)
	Constant mode
	Intermittent mode
	Negative pressure setting in mmHg
	Battery fully charged
	Battery charge level
	Battery charging (bars moving)

Intended use

The Invia Motion NPWT system is intended to be used by healthcare professionals or adequately trained lay users.

Healthcare professionals are responsible to train lay users according to the patient instructions for use and explain all related safety information.

Intended patient population

The Invia Motion NPWT system is intended to be used on patients only exhibiting conditions as described in the indications for use. The device has not been studied in pediatric patients.

Intended environment

The Invia Motion NPWT system is intended for use in acute, extended and home care settings.

Indications for use

The Invia Motion Negative Pressure Wound Therapy (NPWT) system is indicated for patients who would benefit from a suction device (NPWT) as when used on open wounds it creates an environment that promotes wound healing by secondary or tertiary (delayed primary) intention by preparing the wound bed for closure, reducing edema, promoting granulation tissue formation and perfusion, and by removing exudate and infectious material. When used on closed surgical incisions, the Invia Motion NPWT system is also intended to manage the environment of surgical incisions that continue to drain following sutured or stapled closure by maintaining a closed environment and removing exudate via the application of Negative Pressure Wound Therapy.

The Invia Motion NPWT system is appropriate for the following indications:

- Acute or subacute wounds
- Chronic wounds
- Dehisced wounds
- Pressure ulcers
- Diabetic/Neuropathic ulcers
- Venous insufficiency ulcers
- Traumatic wounds
- Partial thickness burns
- Flaps and grafts
- Closed surgical incisions

Contraindications

The Invia Motion NPWT system is contraindicated in the presence of:

- Necrotic tissue with eschar present
- Untreated osteomyelitis
- Non-enteric and unexplored fistulas
- Malignancy in the wound
- Exposed vasculature
- Exposed nerves
- Exposed anastomotic site of blood vessels or bypasses
- Exposed organs

Warnings and safety instructions

WARNINGS

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTIONS

Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury.

Safety related tip

Indicating useful information about the safe use of the device.

The Invia Motion Negative Pressure Wound Therapy system is intended for use as described in these instructions for use.

Medela is only responsible for the effect on BASIC SAFETY, reliability and performance of the Invia Motion NPWT system if it is used in accordance with the instructions for use.

**Please read and observe these warnings and safety instructions before operation.
These instructions for use must be kept with the device.**

Please note that these instructions for use are a general guide for the use of the product. Medical situations must be addressed by a physician.

Warnings

WARNINGS

- Do not modify this equipment without authorization from the manufacturer.
- Consult the indications for use, contraindications, warnings and cautions before using the Invia Motion NPWT pump as a negative pressure source with the Invia foam and gauze dressings. Failure to read and follow all instructions in this manual prior to use may result in serious or fatal injury of the patient.
- Invia motion therapy Instructions advise 24 hours therapy per day without interruption. If therapy is discontinued for more than 2 hours, the dressing should be replaced and therapy restarted by a healthcare professional.
- Unplug the pump before cleaning.
- Do not use accessories other than those specified or sold by the manufacturer as replacement parts for internal components as it may result in increased electromagnetic emissions or decreased electromagnetic immunity of the Invia Motion NPWT pump.
- The Invia Motion NPWT system should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Invia Motion NPWT system should be observed to verify normal operation in the configuration in which it will be used.
- The Invia Motion NPWT pump is intended for a single patient's entire therapy. Reuse on another patient might cause cross contamination.
- Patients must be closely monitored for bleeding. If sudden or increased bleeding is observed, immediately stop use of the pump, apply pressure and seek immediate emergency medical attention.
- Objective indications or signs of a possible infection or complication must be addressed immediately (e.g. fever, pain, redness, increased warmth, swelling or purulent discharge). Non-observance can lead to considerable danger to the patient.
- When using the Invia Motion NPWT system in close proximity to blood vessels, organs, nerves, tendons or ligaments a protective barrier, such as a non-adherent contact layer, must be used.
- Serious or fatal injury can result from bone fragments or sharp edges (e.g. staples or hardware) that could puncture protective barriers, vessels or organs.

- Patients with spinal cord injuries (stimulation of sympathetic nervous system) who experience Autonomic Hyperreflexia during therapy should have NPWT discontinued.
- In the event that defibrillation is required, disconnect the pump from the wound dressing before the patient is defibrillated.
- Special care is advised for dressing placement and removal in order to avoid situations such as unintentional gauze or foam retention.
- This device has not been studied in pediatric patients.
- Do not bring the device into the following areas: Positron Emission Tomography (PET) or Hyperbaric Oxygen Chamber (HBO).
- The Invia Motion NPWT system is Magnetic Resonance (MR) Unsafe. Do not take the Invia Motion NPWT system into the MR environment.
- HF (high-frequency) surgical equipment, radio networks or the like can influence the operation of the device and may not be operated in combination with the Invia Motion pump.
- Explosion hazard – the Invia Motion NPWT pump is not for use in potentially explosive environments including oxygen enriched environments and in areas of flammable anesthetics.
- Portable Radio Frequency (RF) communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 0.3 m to any part of the Invia Motion NPWT system, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- The Invia Motion NPWT system is intended to be used in conjunction with the Invia foam and gauze dressings.
- Invia Foam Dressing Kit with FitPad, X-Large should not be used in conjunction with the Invia Motion NPWT pump.

Cautions

CAUTIONS

- Incorrect use may cause pain and injury to the patient. Excessive negative pressure or an infection of the wound may cause pain to the patient. In either case, the dressing must be changed and the wound inspected in detail.
- The patient should be monitored regularly according to physician instructions and facility guidelines to monitor patient comfort, therapy compliance and signs of wound infection.
- Do not use an Invia Motion canister/tubing set if the sterile packaging is damaged.
- Wireless communication equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, and walkie-talkies can interfere with the Invia Motion pump and should be kept at a minimum distance of 30 cm away from the Invia Motion pump.
- The Invia Motion NPWT pump must remain in an upright (carrying case) or lying (display side up) position during use.
- The total wound volume connected to the Invia Motion NPWT pump shall be 1l or less.

The following statements describe medical conditions that may require special care for the safe and effective use of the Invia Motion NPWT system.

- Patients at high risk for bleeding and hemorrhage, including patients experiencing active bleeding or difficult wound hemostasis.
- Patients taking anticoagulants or platelet aggregation inhibitors.
- Patients with a history of vascular anastomosis or friable, irradiated, sutured or infected blood vessels.
- Patients with spinal cord injury (sympathetic nervous system stimulation): in the event a patient experiences autonomic dysreflexia (sudden changes in blood pressure or heart rate in response to stimulation of the sympathetic nervous system), discontinue negative pressure wound therapy to help minimize sensory stimulation and seek immediate medical assistance.

- Patients with infected wound or osteomyelitis.
- Wounds that involve an enteric fistula.
- To minimize the risk of bradycardia, NPWT must not be placed in the proximity of the vagus nerve.
- Avoid circumferential dressing applications.
- The Negative Pressure Wound Therapy must be used 24 hours per day without interruption. If the pump is stopped for more than 2 hours, the dressing must be changed and therapy restarted.
- Consider the patient's size and weight when prescribing this device.
- Consider mode of therapy – intermittent versus continuous.
- For maximum benefit on closed surgical incision, the Invia negative pressure therapy should be applied immediately post-surgery to clean surgically closed wounds. It is to be continuously applied for a minimum of two days up to a maximum of seven days, with regular dressing changes every 48 to 72 hours. All dressing changes should be applied under direct medical supervision. The Invia therapy system will not be effective in addressing complications associated with the following:
 - Ischemia to the incision or the incision area
 - Untreated or inadequately treated infection
 - Inadequate hemostasis of the incision
 - Cellulitis of the incision area

Physician's orders

A medical professional must assess each wound when applying this therapy and determining negative pressure intensity according to the wound characteristics. Therapy or setting changes may only be done under the order of a physician.

Wound assessment

CAUTION

Patient monitoring: The patient should be monitored regularly according to the physician's instructions and facility guidelines to check for patient comfort, therapy compliance and signs of infection.

WARNING

Objective indications or signs of a possible infection or complication must be addressed immediately (e.g. fever, pain, redness, increased warmth, swelling or purulent discharge). Non-observance can lead to considerable danger to the patient.

Observe the wound/periwound tissue and exudate for signs of infection or other complications. The most common signs of infection include redness, tenderness, fever, swelling, itching, increased warmth in the wound area, strong odor or purulent discharge. Additional symptoms include nausea, vomiting, diarrhea, headache, dizziness, fainting, sore throat with swelling of the mucous membranes, disorientation, high fever ($>38.8^{\circ}\text{C}$), refractory hypotension, orthostatic hypotension, or erythroedema (a sunburn-like rash). More serious complications of infection include pain, discomfort, fever, gangrene, toxic or septic shock. If more serious complications of infection occur, discontinue therapy and consult a healthcare professional immediately.

Dressing technique

 Wound dressings to be applied and changed by healthcare professionals only.

Consult the appropriate Invia dressing instructions for use for information regarding dressing applications included in every Invia dressing carton. Perform a thorough wound cleansing per physician orders prior to dressing applications.

 For closed surgical incision the recommended pressure level is -125 mmHg at constant mode.

Routine dressing changes should occur every 48–72 hours. Dressing changes for infected wounds should be considered more frequently.

Initial delivery

Upon delivery, check the package containing the Invia Motion NPWT system for completeness and general condition.

1 Invia Motion NPWT pump

1 Invia Motion carrying case

1 Invia Motion power supply

1 Invia Motion Clinician Instructions for Use

1 Invia Motion Patient Instructions for Use

1 Invia Motion Quick Card

☞ When first turned on, the Invia Motion NPWT pump needs to be attached to a wall outlet in order to fully charge the battery.

Invia Motion NPWT pump accessories and associated products

The Invia Motion NPWT system was verified in combination with the accessories listed below. For correct and safe operation use the Invia Motion NPWT system with these accessories only:

Y-connector
REF087.0023

Invia Motion NPWT pump

Power supply
REF200.8604

Carrying case
REF087.4003

Canister/Tubing Set
REF087.4016

Wound dressings

 Wound dressings to be applied and changed by healthcare professionals only.

The Invia Motion NPWT system is intended to be used in conjunction with the Invia dressings only. For specific dressing indications, contraindications, warnings and cautions, consult the appropriate Invia dressing instructions for use.

WARNING

Invia foam dressing kit with FitPad, X-Large should not be used in conjunction with the Invia Motion NPWT pump.

CAUTION

The total wound volume connected to the Invia Motion NPWT pump shall be 1l or less.

Invia foam and gauze dressing kits provide their own instructions for use.


Invia Motion disposables

Invia Motion canister/tubing set 0.15l (150 ml), with solidifier (sterile, single use)

Capacity: 0.15l (150ml)
Material: Polypropylene, Silicon, ABS
Sterilization Method: Ethylene Oxide


Intended use

The Medela Invia Motion canister/tubing set 0.15l (150 ml) is an accessory of the Medela Invia Motion NPWT system and it is intended to collect removed wound exudate and infectious material.

Charge battery for the first time

Select the plug. Attach the plug to the power supply.

Plug the Invia Motion power supply into a wall outlet and open the small protective cover of the pump for the Direct Current (DC) socket and attach the power supply. Moving bars on the battery display will indicate charging is occurring. [].

1. Battery is fully discharged

2. Battery is being charged (moving bars)

3. Battery is fully charged

Disconnect from power source

 Check the battery charge level periodically. Make sure the battery compartment on the back side of the pump is closed. Do not remove battery from device as long as therapy is ongoing.

Disconnecting the power source

When the battery is fully charged (after approximately 4 hours) the following non-blinking symbol will appear on the display . Charging of the battery can be slowed when exposed to high temperatures.

To disconnect the Invia Motion NPWT pump from the electrical source, hold the end of the plug and remove it from the pump.

Close the DC socket by pressing the small protective cover over the opening.

Positioning

**CAUTION**

The Invia Motion NPWT pump must remain in an upright (carrying case) or lying (display side up) position during use.

The Invia Motion NPWT pump is protected against the penetration from solid/fluid substances by a hydrophilic bacteria filter.

If the canister is full or the Invia Motion NPWT pump is placed in an inappropriate position, fluids will flow into the safety chamber and clog up the filter. In this situation, the pressure is interrupted and canister full notification will appear (acoustic signal and symbol on display). The Invia Motion canister/tubing set must be replaced.

Preparation for use

To be used by healthcare professionals or adequately trained caregivers only.

Checklist

Prior to using the Invia Motion NPWT system, it is important to check for the following:

- Completeness and general condition of the packaging
- Damage to the power cord and plug
- Damage to the pump
- Proper functioning of the device
- Battery status

Connect Invia Motion canister/tubing set 0.15 l (150 ml)

Make sure the correct canister is used: **REF** 087.4016.

1. Open sterile packaging.

**CAUTION**

Do not use the Invia Motion canister/tubing set if the sterile packaging is damaged or has expired.

2. Insert the Invia Motion canister/tubing set into the Invia Motion NPWT pump connecting port (1). Rotate the canister/tubing set to the right to secure (2).

3. Make sure the canister is aligned with the label to the front.

4. Connect the dressing tubing to the pump tubing.

 Make sure both parts are aligned correctly as shown in the picture below. Push together until you hear a click.


For dressing applications, please refer to the instructions for use provided with the Invia dressing kits.

Switch Invia Motion NPWT pump on

1. Press [] to switch the Invia Motion NPWT pump on.

2. When turned on for the first time, the Invia Motion NPWT pump will start up in the default factory setting (constant mode with pressure level of -125mmHg).

If the Invia Motion NPWT pump is turned off, the setting will be as previously set when the pump is turned back on.


Run-time, timer and safety time

REF 087.4022 Invia Motion – 15 Days Negative Pressure Wound Therapy system

REF 087.4024 Invia Motion – 60 Days Negative Pressure Wound Therapy system

These two Invia Motion NPWT pump versions have a runtime of 15 and 60 days. After one hour of use, the timer is activated and will continue even when the pump is turned off. During the first hour of use, the Invia Motion NPWT pump can be turned on and off without starting the timer. If the timer function is used to readout the therapy time, one hour needs to be added to the displayed runtime. Press [] for information about the used operation time.

After 15 or 60 days, a blinking symbol on the display alike [] will indicate the expiration of the Invia Motion NPWT pump.

Example for 15 days version (blinking)

☞ The Invia Motion NPWT pump runtime versions of 15 or 60 days provide an additional safety usage time of 72 hours. During this time the device can be operated as usual. An acoustic signal will occur after 72 hours and the negative pressure will not be maintained.

Setup carrying case

The Invia Motion NPWT pump is intended to be worn in a carrying case during operation. An Invia Motion carrying case is provided with the Invia Motion NPWT pump.

Attach shoulder/belt strap onto the carrying case as a first step.

1. Open the carrying case.
2. Slide the Invia Motion NPWT pump into the designated upper pocket and close the Velcro strap.
3. Place the excess tubing into the designated lower pocket and close the Velcro strap.
4. Fold the pump and tubing sections of the carrying case together and align the inner Velcro tabs to secure.
5. Secure the bag by closing the lid and fastening the button.

The carrying case strap can be worn over the shoulder or around the waist.

To wear over the shoulder, attach the strap clips to the carrying case. To wear as a belt, thread the strap through the two loops located on the back of the case.

Shoulder strap

Belt strap

Operating instructions

Invia Motion NPWT pump settings

Change pressure level

Press [] and hold for 3 seconds to change the pressure level. The set pressure value symbol on the display will start to blink. You now have 5 seconds to press [] for a new pressure setting, otherwise the pump remains at the previous setting. Each time the button is pressed the next available pressure will appear. When the appropriate pressure level appears on display, the pump will adopt it automatically after 5 seconds. The pressure may only be changed if instructed by a physician or qualified healthcare professional.

Available pressure levels

mmHg	-40	-60	-80	-100	-125	-150	-175

Change therapy mode

The Invia Motion NPWT pump is able to deliver constant or intermittent (5 min. on/2 min. off) suction.

Constant mode

Intermittent mode

Press [] and hold for 3 seconds to change the therapy mode. The therapy mode symbol on the display will start to blink. You now have 5 seconds to press [] to change the mode. Otherwise the pump remains at the previous mode. Each time the button is pressed the next available therapy mode will appear. When the appropriate mode symbol appears, the pump will adopt it automatically after 5 seconds. The therapy mode may only be changed if instructed by a physician or a qualified healthcare professional.

Pump in constant mode

Pump in intermittent mode (on)

Pump in intermittent mode (off)

Mute / Information

Press [] to obtain information about the runtime used or to mute the acoustic signal.

Change Invia Motion canister/tubing set

The Invia Motion canister/tubing set has to be changed on the basis of a visual check or according to the instructions on the display but at a minimum of once a week.

1. Wash hands and prepare clean paper towel.
2. Remove the Invia Motion NPWT pump from the carrying case. Switch the Invia Motion NPWT pump off by pressing [] and holding for 3 seconds. Press only once.
3. Position the clamp next to the quick-connector and close the clamp on the pump tubing and close the clamp on the dressing tubing.


Clamp open

Clamp closed

4. Disconnect the pump tubing from the dressing tubing by pressing on the sides of the quick-connector. Have a paper towel or tissue available to cover the tip of the used tubing to protect against wound fluid leakage.

5. Detach the full canister/tubing set by rotating the canister to the left.

6. Wash hands.
7. Open the sterile packaging of a new Invia Motion canister/tubing set.

CAUTION

Do not use an Invia Motion canister/tubing set if the sterile packaging is damaged or has expired.

Make sure the correct canister is used: **REF** 087.4016.

8. Remove the canister/tubing set from the packaging and attach it to Invia Motion NPWT pump.


For dressing applications, please refer to instructions for use provided with Invia dressing kits.

9. Connect the pump tubing to the dressing tubing. Push together until you hear a click.

Make sure both parts are aligned correctly as shown in the picture below.

10. Open the clamp on the dressing tubing. Also make sure the clamp on the pump tubing is open.
11. Switch the Invia Motion NPWT pump on by pressing [] and slide it back into the carrying case.
12. Discard the used canister/tubing set according to local regulations.

Power supply and battery management

The Invia Motion NPWT pump is operated by an internal rechargeable NiMH battery and while attached to a power source. Make sure the battery compartment on the back side of the pump is closed. Do not remove battery from device as long as therapy is ongoing. If the rechargeable battery is removed more than 30 minutes the Invia Motion NPWT pump will suspend to operate.

For disposal of rechargeable battery after therapy, refer to chapter Disposal.

While in use and connected to a power source, the battery is recharged. The available battery time from a single charge depends upon the runtime of the motor. This is influenced by the size of the wound, air leaks in the dressing, the set pressure and the mode of therapy. The pump only runs if the pressure drops below the set pressure. If the Invia Motion NPWT pump is running continuously, a maximum of 4 hours of battery time is expected before it is necessary to be recharged. The expected battery time generally exceeds 10 hours. If the Invia Motion NPWT pump is not in use, the battery must be charged approximately once every 6 months to ensure optimal functioning.

Before attaching the Invia Motion NPWT pump to a power source, check to ensure your local power supply is the same as the voltage indicated on the specification plate located on the back of the pump.


Battery low

Approximately 30 minutes before the battery is fully discharged, an acoustic signal will sound and the battery symbol [] will start to blink. The acoustic signal can be muted, but the battery symbol continues to blink until the battery is fully discharged. The set pressure is maintained during the 30 minutes, but the battery should be recharged as soon as possible.

Battery fully discharged

The battery symbol [] blinks on the display, an acoustic signal sounds and the negative pressure is not maintained. Connect the pump to the power supply. When the pump is not connected to the power supply within 10 minutes, it switches off.

Connecting power source

Plug the Invia Motion power supply into a wall outlet and open the small protective cover to access the DC socket. Plug the coaxial power supply into the DC socket.

Moving bars on the battery display will indicate charging is occurring [].

1. Battery is fully discharged

2. Battery is being charged (moving bars)

3. Battery is fully charged

Disconnecting power source

When the battery is fully charged (after approximately 4 hours) the following non-blinking symbol will appear on the display []. Charging of the battery can be slowed when exposed to high temperatures. To disconnect the Invia Motion NPWT pump from the electrical source, hold the end of the plug and remove it from the pump.

Close the DC socket by pressing the small protective cover over the opening.

Notifications and troubleshooting


WARNING

Negative Pressure Wound Therapy should run 24 hours a day without interruption. If therapy is discontinued for more than 2 hours, the dressing should be replaced and therapy restarted by a healthcare professional.

If the Invia Motion NPWT pump detects an irregularity in the operation of the pump, an acoustic signal sounds and a notification will appear on the display.

1. An acoustic signal sounds.
2. If the Invia Motion NPWT pump is in the carrying case, open the lid of the carrying case and release the Velcro.


3. Press [] to mute the acoustic signal for 60 seconds (the button can be pressed every time the signal sounds, with the exception of "Battery fully discharged" and after the pump time of 15 or 60 days plus 72 hours of safety time is used up).

Every acoustic signal can be permanently muted by switching Invia Motion NPWT pump off (press [] and hold for 3 seconds). Turn the pump back on to continue treatment.

4. Check the display for one of the following visual symbols.

Canister full


Acoustic signal will sound when the Invia Motion NPWT pump detects a full canister and canister full symbol will appear on the display.

Troubleshooting:

Change the Invia Motion canister/tubing set.

Blockage in tubing

Acoustic signal will sound when the Invia Motion NPWT pump detects blockage in tubing and blockage symbol will appear on the display.

Troubleshooting:

Inspect Invia Motion canister and dressing tubing to ensure tubing is not twisted, kinked or clamped. Verify that the correct canister is used: **REF 087.4016**.

If blockage is identified along canister tubing, switch Invia Motion NPWT pump off and follow instructions in chapter "Operating instructions" to change the Invia Motion canister/tubing set. If the used canister is not the correct one, switch the Invia Motion NPWT pump off, remove the canister. Attach a new canister:

REF 087.4016.

If blockage is identified along dressing tubing, switch Invia Motion NPWT pump off and contact your healthcare professional to replace dressing and restart therapy.

If Invia Motion NPWT pump detects blockage in the intermittent mode, it will continue to run in the on-state until blockage is resolved.

The notification symbol and the acoustic signal may still be active, even if situation is resolved.

Wait a short time for the pump to react.


WARNING

Invia Motion pump does not detect blockage along tubing connected to the drain adapter or the curved arm of the Y-connector.

1. Instruct the patient to frequently check dressing tubing connected to the drain adapter or the curved arm of the Y-connector to make sure tubing is not kinked or blocked.
2. If blockage is identified along tubing connected to the drain adapter or tubing connected to the curved arm of the Y-connector, instruct the patient to contact his healthcare professional.


Leak in system

Acoustic signal will sound when the Invia Motion NPWT pump detects a leak and leakage symbol will appear on the display.

Troubleshooting:

Inspect dressing site for leaks and make sure the Invia Motion canister/tubing set quick-connector is fully engaged. Check that canister is secured to the pump. If a leak is identified at the dressing site, press down on dressing to seal.

The notification symbol and the acoustic signal may still be active, even if situation is resolved.

Wait for some time to let the pump react. If acoustic signal and leakage symbol do not disappear, switch Invia Motion NPWT Pump off and follow instructions in chapter "Operating instructions" to change the Invia Motion canister/tubing set.

☞ If leak is not resolved, switch Invia Motion NPWT pump off and contact your healthcare professional to replace dressing and restart therapy.

Battery low



Troubleshooting:

Charge the battery.

Approximately 30 minutes before the battery is fully discharged, an acoustic signal will sound and the battery symbol [] will start to blink. The acoustic signal can be muted, but the

battery symbol continues to blink until the battery is fully discharged. The set pressure is maintained, but the battery should be recharged as soon as possible.

Battery fully discharged

Troubleshooting:

Charge battery

The battery symbol [] blinks on the display, an acoustic signal sounds and the negative pressure is not maintained. Connect the pump to the power supply. When the pump is not connected to the power supply within 10 minutes, it switches off.

Info used time


The Invia Motion NPWT pumps have a runtime of 15 or 60 days. After one hour of use, the timer is activated and will continue even when the pump is turned off. During the first hour of use, the Invia Motion NPWT pump can be turned on and off without starting the timer. If the timer function is used to readout the therapy time, one hour needs to be added to the displayed runtime. Press [] for information about the used operation time.

The example below displays 6 days and 23 hours of operation time used. The remaining runtime of the Invia Motion – 15 Days version would therefore be 8 days and 1 hour.

Runtime expired

After 15 or 60 days, a blinking symbol on the display alike [] will indicate the expiration of the Invia Motion NPWT pump.

Example for 15 days version (blinking)

 The Invia Motion NPWT pump runtime versions of 15 or 60 days provide an additional safety usage time of 72 hours. During this time the device can be operated as usual. An acoustic signal will occur after 72 hours and the negative pressure will not be maintained.

Overview notifications

	Visual symbol on display	Notification Description	Acoustic signal	Troubleshooting	Remarks/Potential cause
Notification	 <p>Blinking! + acoustic signal!</p>	Leak in system	Can be muted for 60 seconds by pressing  _i	Check the system for leaks.	Dressing could be leaking, check that it is properly sealed.
	 <p>Blinking! + acoustic signal!</p>	Blockage	Can be muted for 60 seconds by pressing  _i	Check the system for blockage. Verify that the correct canister is used: REF 087.4016.	The dressing or canister tubing could be twisted, kinked or clamped. If the Invia Motion pump detects blockage in the intermittent mode, it will continue to run in the on-state until blockage is resolved.
	 <p>Blinking! + acoustic signal!</p>	Canister full	Can be muted for 60 seconds by pressing  _i	Change the Invia Motion canister/ tubing set.	The filter could get wet by improper positioning of the pump (upside down or display side down).
	 <p>Blinking! + acoustic signal!</p>	Battery low	Can be muted permanently by pressing  _i	Connect the Invia Motion pump to the power supply.	The visual symbol blinks all the time until battery is fully discharged. The remaining time of the battery is approximately 30 minutes.
	 <p>Blinking! + acoustic signal!</p>	Battery fully discharged	Can be muted by switching off the pump.	Connect the Invia Motion pump to the power supply.	The remaining time of the battery is approximately 10 minutes.
	 <p>Blinking! (Example for 15 days version)</p>	Pump time expired	The blinking signal can be suppressed permanently by pressing  _i	Dispose of the Invia Motion NPWT pump. Attach another device if necessary.	The Invia Motion NPWT pump provides an additional safety usage time of 72 hours. During this time the Invia Motion pump can be operated as usual. A signal will occur after 72 hours and the negative pressure will not be maintained.

Overview notifications

Notification  Blinking! + acoustic signal!	Battery missing	Can be muted for 60 seconds by pressing  ;	Ensure that battery compartment cover is mounted and closed properly.	Do not remove battery during ongoing therapy.
Internal fault  Acoustic signal!	Internal fault	Can be muted by switching off the pump.	Switch the Invia Motion NPWT pump off, release and refix canister/tubing set. Turn the Invia Motion pump back on. If the fault message repeats, contact Medela customer service.	If the fault message repeats, contact Medela customer service.

Contact your local Medela customer service representative for assistance with product operations.

Cleaning and disinfection

Cleaning

WARNING
 Unplug the pump before cleaning. Close the small protective cover on the DC socket.

The Invia Motion NPWT pump and carrying case are for a single patient's entire therapy. Periodic cleaning is recommended.
 Wipe off with a damp cloth or use a non-abrasive detergent.

Disinfection

Invia Motion can be disinfected with the disinfecting agent group "alcohol".

Sterilization

The Invia Motion NPWT pump and accessories cannot be sterilized. The canister/tubing set is a disposable item and may not be reused or resterilized. Replace after every use!

The Invia Motion canister/tubing set is a single patient use product and is not intended to be reprocessed. Reprocessing may cause the loss of mechanical, chemical and/or biological characteristics.

WARNING
 Reuse on another patient might cause cross contamination.

General notes

Comply with facility or institutional guidelines on hygiene and infection control. Dispose of fluids such as blood and fluids, along with all parts contaminated with them in accordance with facility requirements. Wear suitable gloves for cleaning/disinfection and disposal.

Disposal

Turn the Invia Motion NPWT pump off by pressing [] and holding for 3 seconds.

Disconnect the canister/tubing set from the pump.

Dispose of the canister/tubing set and pump separately and according to local regulations.

When the Invia Motion NPWT pump is ready for disposal, take care to dispose of the Invia Motion NPWT pump and its accessories in accordance with your local regulations.

Disposal of electronic equipment

The symbol above means that the electrical and electronic equipment must not be disposed as normal household refuse. The correct disposal of this device protects and prevents possible damage to the environment or human health. For more information about proper disposal, contact Medela or your local healthcare professional. This symbol is only valid in the European Union. Please dispose according to specified institution protocol or state and local regulations of electrical and electronic equipment.

 Do not remove battery from device as long as therapy is ongoing.

Battery can be disposed separately after completion of therapy:

1. Cut the adhesive label.
2. Open battery cover.
3. Remove rechargeable battery.
4. Dispose rechargeable battery according to local regulations.

Warranty and maintenance

Warranty

The warranty period for the Invia Motion NPWT system expires after 15 or 60 days of usage plus 72 hours of safety time, if used in accordance with the instructions for use. The manufacturer is not liable for any damage or consequential damage caused by incorrect operation, inappropriate usage or unauthorized persons using the device.

Maintenance

The Invia Motion NPWT pump is maintenance free and will not require service.

If an Invia Motion NPWT pump fails within the warranty period due to a manufacturing defect, the pump will be replaced. The original Invia Motion NPWT pump will need to be returned to the supplier.

Safety-related checks

The Invia Motion NPWT system is a device in electrical protection Class II. The safety-related checks are confined to visual inspection of the pump housing and power supply for damage. This check must be performed prior to each use.

Devices of protection Class II do not have a protective earth conductor. There is no need to check the earth leakage current.

The Invia Motion NPWT pump enclosure is made entirely of insulated material. Tests of the enclosure leakage current using common measuring instruments will therefore not reveal measurable values.

The Invia Motion NPWT pump does not have patient circuits or functional earth connections.

Signs and symbols

This symbol indicates safety related tip.

This symbol indicates a CAUTION or WARNING associated with the device (see page 7).

This symbol indicates a class II device.

This symbol indicates to not dispose the device together with household refuse (for EU only).

This symbol indicates the date of manufacture.

This symbol indicates the name and the address of the manufacturer.

This symbol indicates the device is sterilized using ethylene oxide.

This symbol indicates a prescription device. CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician (for US only).

This symbol indicates a type BF applied part.

This symbol indicates manufacturer's catalog number.

This symbol indicates manufacturer's serial number.

This symbol indicates manufacturer's batch code.

This symbol indicates the protection against the ingress of dripping water and solid foreign objects of 12,5mm and greater.

This symbol indicates that the device should not be used after the date shown.

This symbol indicates to follow the instruction for use.

This symbol indicates the compliance with additional USA and Canada safety requirements for medical electrical equipment.

This symbol indicates to not use the device if package is damaged.

This symbol indicates the number of items

This symbol indicates the direct current 9 VDC socket.

This symbol indicates a single use device.
Do not reuse the device.

This symbol indicates the temperature limitation for operation, transport and storage.

This symbol indicates MR unsafe.

This symbol indicates the humidity limitation for operation, transport and storage.

This symbol indicates the atmospheric pressure limitation for operation, transport and storage.

This symbol indicates to keep the device dry.

This symbol indicates to handle the fragile device with care.

This symbol indicates to keep the device away from sunlight.

This symbol indicates that the marked item or its material is part of a recovery or recycling process.

This symbol indicates the correct upright position of the transport package.

0123 This symbol indicates that the device is in conformance with the Medical Device Directive 93/42/EEC.

Technical specifications

medium vacuum
-175 mmHg

low flow

350 g
Without canister

IP22

H x W x D
99 x 98 x 52 mm

ISO 9001
ISO 13485
CE (93/42/EEC), IIa

Operating
conditions

Transport/
Storage
conditions

Invia Motion power supply

Model:

UE08WCP-090088SPA

Input: 100-240VAC,
0.4A, 50-60 Hz

Output: 9VDC, 0.88A

Length of power supply
cable: 172 cm

Invia Motion Pump

[VDC] 9.0

[W] 8.0

Invia Motion battery

Type: NiMH

Invia Motion canister/tubing set

Material: Polypropylene, Silicon, ABS

Sterilization Method: Ethylene Oxide

Electromagnetic compatibility (EMC)

Invia Motion is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 and IEC 60601-1-2:2014 4th Edition according clause 7 and 8.9.

Electromagnetic compatibility (EMC, IEC 60601-1-2:2007 3rd Edition and IEC 60601-1-2:2014 4th Edition, Table 1)

Electromagnetic Emissions		
The Invia Motion NPWT system is intended for use in the electromagnetic environment specified below. The customer or the user of Invia Motion NPWT system should assure that they are used in such environment.		
Emission Tests	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Invia Motion NPWT system uses RF energy only for their internal function. Therefore, its RF emissions are very low and are not likely to cause any interference with nearby electronic equipment.
RF emissions CISPR 11	Class B	The Invia Motion NPWT system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

WARNING

The Invia Motion NPWT system should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Invia Motion NPWT system should be observed to verify normal operation in the configuration in which it will be used.

WARNING

Invia Motion is a medical device that requires special safety precautions and must be installed and placed in operation in accordance with the attached EMC information.

WARNING

Do not use other accessories than those specified or sold by the manufacturer as replacement parts for internal components as it may result in increased emissions or decreased immunity of the Invia Motion pump.

WARNING

EMC

HF (high-frequency) surgical equipment, radio networks or the like can influence the operation of the device and may not be operated in combination with the Invia Motion pump.

WARNING

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Invia Motion. Otherwise, degradation of the performance of this equipment could result.

Electromagnetic compatibility

(EMC, IEC 60601-1-2:2014 4th Edition, Table 2)

Electromagnetic immunity The Invia Motion pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Invia Motion pump should assure that it is used in such environment.			
Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	The Invia Motion NPWT system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines not applicable	The Invia Motion NPWT system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV line-to-earth	± 1 kV differential mode ± 2 kV line-to-earth not applicable	The Invia Motion NPWT system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % U_T for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U_T for 1 cycle 70 % U_T for 25 cycles at 50 Hz Single phase: at 0° 0 % U_T for 250 cycles at 50 Hz	0 % U_T for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U_T for 1 cycle 70 % U_T for 25 cycles at 50 Hz Single phase: at 0° 0 % U_T for 250 cycles at 50 Hz	The Invia Motion NPWT system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. If the user of the Invia Motion pump requires continued operation during power mains interruptions, it is recommended that the Invia Motion pump be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	The Invia Motion NPWT system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
NOTE: U_T is the a.c. mains voltage prior to application of the test level.			

Electromagnetic compatibility
(EMC, IEC 60601-1-2:2014 4th Edition, Table 4)

Electromagnetic immunity The Invia Motion pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Invia Motion pump should assure that it is used in such environment.			
Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	10 Vrms 150 kHz to 80 MHz	10 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the Invia Motion pump including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 0.35 \sqrt{P}$ $d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz $d = 0.7 \sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Invia Motion pump is used exceeds the applicable RF compliance level above, the Invia Motion pump should be observed to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as reorienting or relocating the Invia Motion pump.			
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.			

Electromagnetic compatibility
(EMC, IEC 60601-1-2:2014 4th Edition, Table 6)

Recommended separation distance between portable and mobile RF communications equipment and the Invia Motion pump			
The Invia Motion pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Invia Motion pump can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Invia Motion pump as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter	Separation distance according to frequency of transmitter		
	m		
W	150 kHz to 80 MHz $d = 0.35 \sqrt{P}$	80 MHz to 800 MHz $d = 0.35 \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$
0.01	0.035	0.035	0.07
0.1	0.11	0.11	0.22
1	0.35	0.35	0.7
10	1.11	1.11	2.21
100	3.5	3.5	7.0
For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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Medela AG