DIRECTORY

English	1	Português	28	Eesti	55
Français	4	Česky	31	Lietuvių	58
Deutsch	7	Magyar	34	Latviešu	61
Español	10	Polski	37	Türkçe	64
Italiano	13	Slovensky	40	Русский	67
Nederlands	16	Norsk	43	Srpski	70
Dansk	19	Suomi	46	中文	73
Svenska	22	Български	49	繁體中文 (台灣)	76
Ελληνικά	25	Română	52	한국어	79


English

Flow-Through Housing Model 93505

For single use only

Concept/Description

When used in conjunction with the Injectate Temperature Probe, the Flow-Through Housing (Model 93505) allows direct temperature measurement of the thermal indicator at the proximal lumen infusion site. Measurement of the injectate temperature at the infusion site helps reduce errors frequently encountered in standard injection techniques (for example, injectate warming caused by handling the syringe). The temperature of the injectate is fed directly by the injectate temperature probe into a patient monitor for curve analysis and determination of cardiac output. An appropriate computation constant must be entered into the computer to determine the correct cardiac output. A 3-way stopcock should be placed proximal to the flow-through housing (Figure 3) for standard catheter maintenance and fluid withdrawal.

Instructions For Use

Connect Housing to Catheter and Patient Monitor (Figure 2)

Note: Use Aseptic Technique

1. Verify that the fluid path is completely free of air. Attach the flow-through housing to the catheter injectate lumen stopcock. (Figure 3)
2. Insert the reusable injectate temperature probe into the flow-through housing (Figure 4). Secure the system to minimize movement with respect to the patient.

3. Connect the injectate temperature probe cable to the "Injectate Probe" lead of the computer cable.

To Remove Flow-Through Temperature Probe (Figure 4)

Grasp the flow-through housing, gently press both tabs, and withdraw the temperature probe.

Edwards, Edwards Lifesciences, the stylized E logo, Swan and Swan-Ganz are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.


Figure 2. Flow-Through Housing and Temperature Probe

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1. Select Injectate volume and catheter size on monitor to determine computation constant or enter a computation constant corresponding to injectate volume, temperature range, and catheter model being used.

Note: Refer to the catheter's instructions for use for computation constants.

2. Prepare the monitor for operation according to the directions in the Operations Manual.

To Measure Cardiac Output

1. Attach a prefilled syringe (D_2W cold injectate) to the flow-through housing.
2. Turn the stopcock at the catheter injectate hub to close the IV flush and open the fluid path between the syringe and catheter.
3. Press start on patient monitor.
4. Immediately perform injections at a rate of at least 2.5 ml/sec.

5. Perform cardiac output determinations according to the instructions in the appropriate cardiac output computer operations manual.

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6. At the conclusion of measurements, return the catheter stopcock to its original (flush) position.

Injectate Temperature

Ensure that the injectate temperature falls within 6°C to 16°C for cold injectate and 18°C and 25°C for room temperature injectate.

Lower or higher injectate temperatures will affect cardiac output accuracy.

Precautions

Care should be taken to avoid air entrapment during setup and use. The flow-through housing and syringe should be secured to minimize movement with respect to the patient.

How Supplied/Recommended Use

This product is supplied sterile and nonpyrogenic, and will remain so as long as the package is unopened and undamaged. Do not resterilize. It is recommended for single patient use only. It is recommended that the product be changed (replaced) for each patient according to hospital procedures for changing tubing sets containing dextrose solutions.

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Inside the U.S. and Canada
(24 hours):800.822.9837
Outside the U.S. and Canada
(24 hours):949.250.2222
In the UK:0870 606 2040-Option 4
In Ireland:01 8211012 Option 4

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Warning: This device is designed, intended, and distributed for single use only. **Do not re-sterilize or reuse this device.** There are no data to support the sterility, nonpyrogenicity, and functionality of the device after reprocessing.

Because of continuing product improvement, prices, specifications, and model availability are subject to change without notice.

This product is manufactured and sold under one or more of the following US patent(s): US Patent No. 6,248,077; 6,336,902; and corresponding foreign patents.

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11/16
192306006A
©Copyright 2016, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences Services GmbH
Edisonstrasse 6
85716 Unterschleissheim
Germany

Telephone	949.250.2500 800.424.3278
FAX	949.250.2525

DATA MATRIX
BARCODE FPO

CO-Set+ Cooling Container: 93520

DIRECTORY

English	1	Magyar	5
Français	2	Polski	5
Deutsch	2	Slovensky	5
Español	2	Norsk	5
Italiano	3	Suomi	6
Nederlands	3	Български	6
Dansk	3	Română	6
Svenska	3	Eesti	7
Ελληνικά	4	Lietuvių	7
Português	4	Latviešu	7
Česky	4	Türkçe	7

Cooling Container and Coil ■ Container pour glace et serpentín de refroidissement ■ Kühlbehälter und Spirale ■ Recipiente de Refrigeración y Serpentin ■ Contenitore di raffreddamento e serpentina ■ Koelcontainer en Koelingspoel ■ Kølebeholder og spiral ■ Kylbehållare och slinga ■ Περιέκτης ψύξης και σπειρώμα ■ Recipiente e serpentina de refrigeração ■ Chladicí nádoba a vinutí ■ Hűtőtartály és tekercs ■ Pojemnik chłodzący z węzownią ■ Chładiaca nádoba a cievka ■ Kjølebeholder og spole ■ Jäähdytysaläili ja kierukka ■ Охлаждающий контейнер и намотка ■ Recipient de răcire și bobină ■ Jahutuskonteiner ja -pool ■ Šaldomoji talpyklė ir spirale ■ Dzesēšanas tvertne un spirāle ■ Soğutma Kabi ve Spiral

Lid ■ Couvercle ■ Deckel ■ Tapa ■ Coperchio ■ Deksel ■ Låg ■ Lock
 ■ Κοπάκι ■ Tampa ■ Uzāvēr ■ Fedél ■ Pokrywa ■ Veko ■ Lokk
 ■ Kansi ■ Kanak ■ Capac ■ Kaas ■ Dangtis ■ Vāks ■ Kapak

Mounting Tab ■ Languette de fixation
 ■ Befestigungsvorrichtung ■ Resalte de Montaje
 ■ Linguetta montante ■ Monteerlus ■ Monteringsklap
 ■ Fästansordning ■ Προεξοχή στήριξης
 ■ Lingueta de fixaço ■ Pířchytňý úchyt
 ■ Rögztöfűl ■ Występ montażowy ■ Montážny výstupok
 ■ Feste for lokk ■ Kiinnityskorvake
 ■ Накрайник за монтаж ■ Cap de montare
 ■ Kinnitustipik ■ Surinkimo ašelė
 ■ Montāžas ķepiņa ■ Montaj Tırnağı

Outlet Slot ■ Fente de sortie ■ Auslaßöffnung
 ■ Ranura de Salida ■ Scanalatura della via d'uscita
 ■ Uitgangsoopening ■ afløbshane ■ Utflödeskåra
 ■ Έγκομή εξόδου ■ Ranhura de saída
 ■ Výstupní otvor ■ Kimeneti nyílás ■ Rowek wylotowy
 ■ Výstupná drážka ■ Uitgangsapning
 ■ Ulostuloaukko ■ Изходен отвор
 ■ Fantā de ieşire ■ Väljalaskpesa
 ■ Ištečejimo anga ■ Izeļas atvere ■ Çıkış Yuvası

Thermal Tubing ■ Tubulure thermique
 ■ Thermalschlauch ■ Tubo Térmico ■ Tubo termico
 ■ Thermische slang ■ Varmeslange
 ■ Temperaturslang ■ Θερμικός σωλήνας
 ■ Tubo térmico ■ Tepelné odolná hadička
 ■ Hőscő ■ Orurowanie termiczne ■ Tepelná hadička
 ■ Termisk slang ■ Lämpöletku
 ■ Термална тръбичка ■ Tubulatură termică
 ■ Termovooliikud ■ Terminiai vamzdeliai
 ■ Termocaurule ■ Termal Boru

Spigot ■ Robinet ■ Abblbahn ■ Llave ■ Rubinetto ■ Afloop
 ■ Hane ■ Avtappningskran ■ Kouvola ■ Ligaço de tubo
 ■ Kohout ■ Leesztő csap ■ Kurek czerpalny
 ■ Výpustný ventil ■ Avlop ■ Tulppa ■ Кранче
 ■ Robinet ■ Kraan ■ Čiaupas ■ Krāns ■ Musluk

Model 93520 Cooling Container ■ Modèle 93520 Container pour glace ■ Modell 93520 Kühlbehälter
 ■ Modelo 93520 Recipiente de Refrigeración ■ Modello 93520 Contenitore di raffreddamento ■ Model 93520 Koelcontainer ■ Model 93520 Kølebeholder
 ■ Modell 93520 Kylbehållare ■ Περιέκτης ψύξης μοντέλο 93520
 ■ Recipiente de refrigeração, Modelo 93520 ■ Chladicí nádoba model 93520 ■ 93520-as modellszámú hűtőkonténer
 ■ Pojemnik chłodzący model 93520
 ■ Chładiaca nádoba, model 93520 ■ Kjølebeholder, model 93520
 ■ Mallin 93520 jäähdytysaläili
 ■ Охлаждающий контейнер, модель 93520
 ■ Recipient de răcire model 93520
 ■ Jahutuskonteiner, model 93520
 ■ Modelis 93520 Šaldomoji talpyklė
 ■ Modelis 93520 Dzesēšanas tvertne
 ■ Model 93520 Soğutma Kabi

Cooling Coil ■ Serpentin de refroidissement
 ■ Kühlschleife ■ Serpentin de Refrigeración
 ■ Serpentina di raffreddamento ■ Koelingspoel
 ■ Kølespiral ■ Kylslinga ■ Σπειρώμα ψύξης
 ■ Serpentina de refrigeração ■ Chladicí vinutí
 ■ Hűtőtekercs ■ Węzownica chłodząca
 ■ Chładiaca cievka ■ Kjølespole
 ■ Jäähdytyskierukka ■ Охлаждающая намотка
 ■ Serpentină de răcire ■ Jahutuspool
 ■ Šaldomoji spirale ■ Dzesēšanas spirāle
 ■ Soğutma Spirali

Tubing Connector ■ Connecteur de la tubulure
 ■ Schlauchanschluß ■ Conector de Tubo
 ■ Connettore ai tubi ■ Aansluiting voor de slang
 ■ Slangetiisluutning ■ Slanganslutning
 ■ Σύνδεσμος σωλήνα ■ Conector de tubo
 ■ Plastová spojka
 ■ Csőcsatlakozó ■ Złącze węzownicy
 ■ Pripojenie hadičky ■ Slangekobling
 ■ Letkuiliitin ■ Τρῶθεν конектор
 ■ Conector al tubulaturii
 ■ Voolikute pistmik
 ■ Vamzdelių jungtis
 ■ Caurulu savienotājs
 ■ Tüp Konektörü

Preparing Cooling Container

For Figure 1, see page 8.

1. Pack the container with crushed ice up to the inside ridge in the internal ribs.
2. Seat the cooling coil on the first ridge. For maximum cooling, position the coil with outlet at the bottom as shown.
3. Place the inlet and outlet tubing in the slots provided at the top of the container. **Be sure that the white plastic tubing connector is inside cooling container as shown.**
4. Completely cover the cooling coil with crushed ice.
5. Add iced water until it is visible above the ice level.
6. Place any unused length of tubing into the ice container to provide additional cold solution.
7. Lock the lid on the cooling container.

Note: Use spigot to drain excess water from melting ice before refilling. There is no need to disassemble the system (See page 8, Figure 1).

Note: An optional bracket (Model 93521) is available to mount the container on an I.V. pole with a diameter between 1/2 and 1-1/2 inches (1.27 and 3.81 cm).

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Edwards Lifesciences, the stylized E logo, Edwards and CO-Set+ are trademarks of Edwards Lifesciences Corporation; Edwards Lifesciences and the stylized E logo are registered in the U.S. Patent and Trademark Office.

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PM22SL5x7.2

Edwards

**Edwards Lifesciences Services GmbH**Edisonstr. 6
85716 Unterschleißheim
Germany**Manufacturer** 
Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614-5686 USA
Made in USATelephone 949.250.2500
800.424.3278
FAX 949.250.252512/09
115102004 A
©Copyright 2009, Edwards Lifesciences LLC
All rights reserved.DATA MATRIX
BARCODE FPO →

DIRECTORY

English	1	Polski	8
Français	2	Slovensky	9
Deutsch	2	Norsk	10
Español	3	Suomi	10
Italiano	4	Български	11
Nederlands	4	Română	11
Dansk	5	Eesti	12
Svenska	5	Lietuvių	13
Ελληνικά	6	Latviešu	13
Português	7	Türkçe	14
Česky	7	Português do Brasil	14
Magyar	8		

English

CO-Set+ In-Line Injectate Temperature Probe Model 93522

Description

The model 93522 CO-Set+ In-Line Injectate Temperature Probe is used in conjunction with injectate delivery systems: model 93600 (Cold) CO-Set+ or model 93610 (Room Temperature) CO-Set+ or model 93505 Flow-Through Housing. The CO-Set+ In-Line Injectate Temperature Probe is a three element network consisting of two resistors and one thermistor configured to produce a linear resistance versus temperature curve. The nominal measurement accuracy over the range of 0° to 25 °C is ± 0.3 °C. The CO-Set+ In-Line Injectate Temperature Probe is stationed at the point where the syringe injects the thermal indicator into the catheter's proximal lumen hub. The probe is able to precisely measure the injectate temperature as it is being introduced into the catheter.

Instructions for use

Refer to Instructions For Use provided with the corresponding injectate delivery system (model 93600 CO-Set+ or model 93610 CO-Set+ or model 93505 Flow-Through Housing).

Cleaning

Clean the probe by wiping with a cleaning solution such as isopropyl alcohol.

Edwards Lifesciences, the stylized E logo, Edwards and CO-Set+ are trademarks of Edwards Lifesciences Corporation.

The recommended Use By Date imprinted on the flag label on the cable represents the recommended date the cable should be replaced or returned to Edwards Lifesciences for evaluation.

Please contact Edwards Technical Support or your local Edwards representative for further assistance.

Precaution

Do not immerse the probe's electrical connector in liquids.

How Supplied

The Edwards Lifesciences model 93522 is shipped non-sterile.

Warranty

Edwards warrants that the equipment is fit for the purposes and indications described in the labeling for the period of one (1) year after the date of purchase when used in accordance with the directions for use. Unless the equipment is used in accordance with such instructions, this warranty is void and of no effect. No other express or implied warranty exists, including any warranty of merchantability or fitness for a particular purpose. Edwards' sole obligation and purchaser's exclusive remedy for breach of warranty shall be limited to repair or replacement at Edwards' option. Edwards shall not be liable for proximate, incidental, or consequential damages.

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Inside the U.S. and Canada	
(24 hours):	800.822.9837
Outside the U.S. and Canada	
(24 hours):	949.250.2222
In the UK:	0870 606 2040 - Option 4
In Ireland:	01 8211012 Option 4

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EC REP

Edwards Lifesciences Services GmbH
Edisonstr. 6
85716 Unterschleißheim
Germany

5/14
199098002 A
©Copyright 2014, Edwards Lifesciences LLC
All rights reserved.

Manufacturer 
Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA
Made in Malaysia

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

DATA MATRIX
BARCODE FPO


DIRECTORY

English	1	Português	37	Eesti	73
Français	5	Česky	41	Lietuvių	77
Deutsch	9	Magyar	45	Latviešu	81
Español	13	Polski	49	Türkçe	85
Italiano	17	Slovensky	53	Русский	89
Nederlands	21	Norsk	57	Srpski	93
Dansk	25	Suomi	61	中文	97
Svenska	29	Български	65	繁體中文 (台灣)	100
Ελληνικά	33	Română	69	한국어	104

English

CO-Set+ Closed Injectate Delivery System for Cold Injectate 93600


Figure 2

MS033-35


Figure 3. Injectate, IV Spike, Snap Clamp, and Tubing and Cooling Coils


Figure 4. Syringe and Flow-Through/Check Valve

Caution: This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

For single use only

Description

The CO-Set+ closed injectate delivery system for cold injectate (Model 93600), provides simple and convenient delivery of cold injectate for the measurement of cardiac output by thermodilution. The unit is designed for use with a patient monitor and a Swan-Ganz thermodilution catheter.

Edwards, Edwards Lifesciences, the stylized E logo, CO-Set, CO-Set+, Swan and Swan-Ganz are trademarks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

CO-Set+ closed injectate delivery system for cold injection (shown in Figure 1) is connected between a suitable supply of sterile injectate (D_2W) and a Swan-Ganz thermodilution catheter. After the system is primed and the cooling coil is immersed in an ice bath, the injectate will reach operating temperature within 5 minutes.

The covered ice bath will remain usable for up to 6 hours; for best results, the ice bath should be checked periodically to ensure that ice is still present.

Instructions for Use

Equipment Needed

1. CO-Set+ Closed Injectate Delivery System for Cold Injectate (Model 93600)
2. Swan-Ganz thermodilution catheter
3. Flow-Through Injectate Temperature Probe

4. Cooling Container (Model 93520 or equivalent)
5. Crushed ice
6. Iced water
7. Sterile D_2W in IV container (bag recommended for use as injectate (500 ml minimum)
8. Patient Monitor
9. Connecting cable for patient monitor
10. Empty flush container for collecting discarded injectate
11. Cooling Container IV Pole Mounting Bracket (Model 93521) (optional)
12. Continuous flush device


Figure 5. Cooling Container and Coil


Figure 6. Flow-Through/Check Valve Connected to Catheter Lumen

Connecting Injectate Container to Closed System (See Fig. 3)

Use aseptic technique.

1. Suspend injectate container (IV bag or bottle) at the desired height; do not place it higher than three feet (90 cm) above the catheter insertion site.
2. Remove CO-Set+ closed injectate delivery system for cold injection from its package.
3. Close the snap clamp.
4. Remove the protective cap from the IV spike and insert the nonvented spike into the injectate container.

Note: If a non-vented glass bottle is used, an airway needle or vented extension set is required to vent the bottle.

5. Pull apart the **TUBING COIL (not cooling coil)** to obtain the desired tubing length between the injectate container and cooling container.

WARNING: Do not pull apart the Cooling Coil. Injectate may not cool to appropriate temperature if Cooling Coil is uncoiled. If additional tubing length is needed, pull apart the tubing coil.

Flushing Air from Fluid Path

Use aseptic technique.

1. Connect the 10 ml syringe to the flow-through/check valve (See Fig. 4). **Caution: Do not connect the syringe to the catheter side of the flow-through/check valve.** Damage to the flow-through/check valve can occur if the flow-through/check valve is flushed or pressurized from the wrong side. Disabling back-pressurization of the flow-through/check valve can also occur if the 3-way stopcock between the flow-through/check valve and the catheter is not

opened to the catheter before injecting. If the flow-through/check valve becomes damaged, then it is possible that blood could be aspirated back into the CO-Set+ system syringe. Should this occur, use of the CO-Set+ system (syringe and flow-through/check valve) should be discontinued and the system replaced.

2. Ensure that all the components are securely connected.
3. Open the snap clamp to allow solution to flow from the injectate container. To collect the discarded solution, hold the flow-through/check valve over an empty container.
4. **Slowly** pull the syringe plunger to prime the system; then push the plunger in. **Repeat** 5 or 6 times, or until the system is completely free of air. The system may also be filled by squeezing the IV bag or plastic bottle. Remove the syringe from the flow-through/check valve when filling the system by this method.
5. Return the syringe plunger to the fully depressed position.
6. Close the snap clamp.

Preparing Cooling Container (See Fig. 5)

1. Pack the container with crushed ice up to the inside ridge in the internal ribs.
2. Seat the cooling coil on the first ridge. For maximum cooling, position the coil with the outlet at the bottom as shown.
3. Place the inlet and outlet tubing in the slots provided at the top of the container. Be sure that the white plastic tubing connector is inside the container as shown.
4. Completely cover the cooling coil with crushed ice.

5. Add iced water until it is visible above the ice level.
6. Place any unused length of tubing coil into the ice container to provide additional cold solution.
7. Lock the lid on the cooling container.
8. Let sit for 5 minutes.
9. Flush system once to ensure cold fluid is used for injection.

Note: Use the spigot to drain excess water before refilling with more ice. There is no need to disassemble the system.

Note: An optional bracket (Model 93521) is available to mount the container on an IV pole with a diameter between 1/2 and 1 1/2 inches (1.27 and 3.81 cm).

Connecting System to Catheter and Patient Monitor

Use aseptic technique.

1. Verify the absence of air in the fluid path. Attach the flow-through/check valve to the catheter injectate lumen stopcock (See Fig. 6).
2. Insert the injectate temperature probe into the flow-through/check valve until it snaps into place (See page 4, Fig. 7). Secure the system to minimize movement with respect to the patient.
3. Connect the injectate temperature probe cable to the "Injectate Probe" lead of the monitor cable.

Prepare Patient Monitor for Operation

1. Select injectate volume and catheter size on monitor to determine computation constant or enter the computation constant corresponding to the injectate volume and catheter model

MS033-32


Figure 7. Temperature Probe Connection

MS033-29

Pull (to unsnap)


Figure 8. Removing the Temperature Probe

being used (see package insert accompanying the catheter for computation constant).

2. Prepare the monitor for operation according to the directions in the appropriate monitor operations manual.

CO-Set+ closed injectate delivery system is now ready for use.

Cardiac Output Measurement

1. Open the snap clamp.
2. Turn the stopcock at the catheter injectate hub to close the IV flush and to open the fluid path between the syringe and catheter.
3. Withdraw desired amount (5 to 10 ml) of injectate into syringe.
4. Press start on patient monitor.
5. Immediately perform injections at a rate of at least 2.5 ml/sec.

Note: The syringe plunger and the finger grips on the syringe jacket have been designed to help accommodate a wide range of user hand sizes for efficient and comfortable injection of fluid.

To ensure uniform fluid injections, care should be taken to avoid pushing in the syringe plunger at an angle to the syringe barrel during injections.

6. Perform cardiac output determination according to instructions in the appropriate patient monitor operations manual.

Note: To improve cardiac output reproducibility, space injections approximately one minute apart and increase injectate volume to 10 ml if not already doing so.

7. At the end of the measurements, return the stopcock to its original (flush) position (assuming that the syringe is empty), and close the snap clamp.

Injectate Temperature

Ensure that the injectate temperature falls within the following ranges:

Injectate Volume	Temperature Range
10 ml	6 °C - 12 °C
5 ml	8 °C - 16 °C

Lower or higher temperatures will affect cardiac output accuracy when using the computation constants provided.

Removing the Injectate Temperature Probe (See Fig. 8)

Grasp the flow-through/check valve, and carefully withdraw (unsnap) the temperature probe. Do not squeeze or pinch the area around the probe connector as this will make probe removal from the flow-through/check valve difficult.

Precautions

Care should be taken to avoid air entrapment during setup and use. The flow-through/check valve and syringe should be secured to minimize movement with respect to the patient.

How Supplied/Recommended Use

Contents sterile and nonpyrogenic if package is undamaged or unopened. Do not use if package is opened or damaged. Do not resterilize.

The CO-Set+ closed injectate delivery system is recommended for single patient use only. It is recommended that CO-Set+ closed injectate delivery system be changed (replaced) for each patient according to hospital procedures for changing tubing sets containing dextrose solutions.

Shelf Life

The recommended shelf life is as marked on each package. Resterilization will not extend the shelf life of this product.

Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:
 Inside the U.S. and Canada
 (24 hours):800.822.9837
 Outside the U.S. and Canada
 (24 hours):949.250.2222
 In the UK:0870 606 2040 - Option 4
 In Ireland01 8211012 Option 4

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

WARNING: This device is designed, intended, and distributed for single use only.

Do not resterilize or reuse this device. There are no data to support the sterility, nonpyrogenicity, and functionality of the device after reprocessing.

Prices, specifications, and model availability are subject to change without notice.

Refer to the symbol legend at the end of this document.

Product bearing the symbol:

has been sterilized using Ethylene Oxide.

Alternatively, product bearing the symbol:

has been sterilized using Irradiation.

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COS25LSX7.1

Edwards

Edwards Lifesciences LLC
One Edwards Way
Irvine, CA 92614 USA

Made in Dominican Republic
Edwards Lifesciences AG
Parque Industrial Itabo
Km 18.5 Carr. Sanchez
Haina, San Cristobal, Dominican Republic

11/16
192010008 A
© Copyright 2016, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences Services GmbH
Edisonstrasse 6
85716 Unterschleissheim
Germany

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

DATA MATRIX
BARCODE FPO


DIRECTORY

English	1	Português	33	Eesti	65
Français	5	Česky	36	Lietuvių	68
Deutsch	8	Magyar	40	Latviešu	72
Español	12	Polski	43	Türkçe	75
Italiano	15	Slovensky	47	Русский	79
Nederlands	19	Norsk	50	Srpski	83
Dansk	22	Suomi	54	中文	86
Svenska	26	Български	57	繁體中文 (台灣)	89
Ελληνικά	29	Română	61	한국어	92

English

CO-Set+ Closed Injectate Delivery System for Room Temperature Injectate Model 93610


Figure 2. View of Injectate Container, IV Spike, Snap Clamp, and Tubing


Figure 3. Syringe and Check Valve


Figure 4. Flow-Through Housing Connected to Catheter Lumen

For single use only

Description

The CO-Set+, closed injectate delivery system for room temperature injectate (Model 93610), provides simple and convenient delivery of room temperature injectate for the measurement of cardiac output by thermodilution. The unit is designed for use with a patient monitor and a Swan-Ganz thermodilution catheter.

CO-Set+ closed injectate delivery system (shown in Figure 1) is connected between a suitable supply of sterile injectate (D₃W) and a Swan-Ganz thermodilution catheter.

Instructions for Use

Equipment Needed

1. CO-Set+, closed injectate delivery system for room temperature injectate (Model 93610)
2. Swan-Ganz thermodilution catheter

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3. Flow-Through Injectate Temperature Probe
4. Sterile D₃W in IV container (bag recommended) for use as injectate (500 ml minimum)
5. Patient Monitor
6. Connecting cable for patient monitor
7. Empty flush container for collecting discarded injectate
8. Continuous flush device

Connecting Injectate Container to Closed System (Figure 2)

Use Aseptic Technique

1. Suspend injectate container (IV bag or bottle) at the desired height; do not place it higher than 3 feet (90 cm) above the catheter insertion site.
2. Remove a CO-Set+ system from its package.
3. Close the snap clamp.
4. Remove the protective cap from the IV spike and insert the non-vented spike into the injectate container.

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Flushing Air from Fluid Path

Use Aseptic Technique

1. Connect the 10 ml syringe to the flow-through/check valve (Figure 3).
Caution: Do not connect the syringe to the catheter side of the flow-through/check valve. Damage to the flow-through/check valve can occur if the flow-through/check valve is flushed or pressurized from the wrong side. Disabling back-pressurization of the flow-through/check valve can also occur if the 3-way stopcock between the flow-through/check valve and the catheter is not opened to the catheter before injecting. If the flow-through/check valve becomes damaged, then it is possible that blood could be aspirated back into the CO-Set+ system syringe. Should this occur, use of the CO-Set system (syringe and flow-through/check valve) should be discontinued and the system replaced.
2. Ensure that all the components are securely connected.


Figure 5. Temperature Probe and Flow-Through Housing


Figure 6. Removing the Temperature Probe from Housing

3. Open the snap clamp to allow solution to flow from the injectate container. To collect the discarded solution, hold the flow-through/check valve over an empty container.
4. **Slowly** pull the syringe plunger to prime the system; then push the plunger in. **Repeat** 2 or 3 times, or until the system is completely free of air. The system may also be filled by squeezing the IV bag. Remove the syringe from the flow-through/check valve when filling the system by this method.
5. Return the syringe plunger to the fully depressed position.
6. Close the snap clamp.

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Use Aseptic Technique

1. Verify the absence of air in the fluid path. Attach the flow-through/check valve to the catheter injectate lumen stopcock (Figure 4).
2. Insert the injectate temperature probe into the flow-through/check valve until it snaps into place (Figure 5). Secure the system to minimize movement with respect to the patient.
3. Connect the injectate temperature probe cable to the "Injectate Probe" lead of the computer cable.

Prepare Patient Monitor for operation

1. Select injectate volume and catheter size on patient monitor to determine computation constant, or enter the computation constant corresponding to the injectate volume and catheter model being used (see Computation Constant Table).
2. Prepare the monitor for operation according to the directions in the appropriate monitor operations manual. The CO-Set+, closed injectate delivery system, is now ready for use.

Cardiac Output Measurement

1. Open the snap clamp.
2. Turn the stopcock at the catheter injectate hub to close the IV flush and to open the fluid path between the syringe and catheter.
3. Withdraw desired amount of injectate (5 to 10 ml) into syringe.
4. Press start on the monitor.
5. Immediately perform injection at a rate of at least 2.5 ml/sec.

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To ensure uniform fluid injections, care should be taken to: a) avoid pushing in the syringe plunger at an angle to the syringe barrel during injection and b) ensure that the contamination shield does not become twisted.

6. Perform cardiac output determination according to the instructions in the patient monitor operations manual.

Note: To improve cardiac output reproducibility, space injections approximately 1 minute apart and increase injectate volume to 10 ml if not already doing so.

7. At the end of the measurements, return the stopcock to its original (flush) position (assuring that the syringe is empty) and close the snap clamp.

Injectate Temperature

Ensure that the injectate temperature falls within the range of 18° to 25 °C for 5 or 10 ml injectate. **Lower or higher injectate temperatures will affect cardiac output accuracy when using the provided computation constants for room temperature injectate.**

Removing the Flow-Through Temperature Probe (Figure 6)

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How Supplied

Contents sterile and fluid path nonpyrogenic if package is undamaged or unopened. Do not use if package is opened or damaged. Do not resterilize.

The CO-Set+ system is for single patient use only. It is recommended that the CO-Set+ system be changed (replaced) for each patient according to hospital procedures for changing tubing sets containing dextrose solutions.

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Technical Assistance

For technical assistance, please call Edwards Technical Support at the following telephone numbers:

Inside the U.S. and Canada
(24 hours):800.822.9837
Outside the U.S. and Canada
(24 hours):949.250.2222
In the UK: 0870 606 2040 - Option 4
In Ireland: 01 8211012 Option 4

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Warning: This device is designed, intended, and distributed for single use only. **Do not resterilize or reuse the device.** There are no data to support the sterility, nonpyrogenicity, and functionality of this device after reprocessing.

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This product is manufactured and sold under one or more of the following U.S. patent(s): U.S. Patent No. 6,248,077; 6,336,902; and corresponding foreign patents. Likewise, additional patents pending.

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Computation Constants Room Temperature Injectate System (5 or 10 ml: 18° to 25 °C)

Injectate Volume (ml)	131F7 & 151F7	C144F7	096F6	All 7.5F Catheters	132F5
10	0.608	0.585	0.607	0.595	N/A
5	0.301	0.287	0.301	0.298	0.307

3 ml injectate volume is not recommended.

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11/16
192011008 A
©Copyright 2016, Edwards Lifesciences LLC
All rights reserved.

Edwards Lifesciences Services GmbH
Edisonstrasse 6
85716 Unterschleissheim
Germany

Telephone 949.250.2500
800.424.3278
FAX 949.250.2525

DATA MATRIX
BARCODE FPO

DIRECTORY

English	1	Português	10	Eesti	19
Français	2	Česky	11	Lietuvių	20
Deutsch	3	Magyar	12	Latviešu	21
Español	4	Polski	13	Türkçe	22
Italiano	5	Slovensky	14	Русский	23
Nederlands	6	Norsk	15	Srpski	24
Dansk	7	Suomi	16	中文	25
Svenska	8	Български	17	繁體中文 (台灣)	26
Ελληνικά	9	Română	18	한국어	27

English

CO-Set+ Model 93650 Syringe Assembly

for cold or room temperature injectate

For single use only

Description

The CO-Set+ Syringe Assembly is a 10 ml syringe for use with the Model 93600 (Cold) or Models 93-610 and 93610 (Room Temperature) closed injectate delivery systems. This individually packaged and sterilized syringe assembly is intended to replace the originally supplied syringe assembly if the need arises.

Additionally, the CO-Set+ Syringe Assembly may be used with the Model 93505 Flow-Through Housing in place of the user supplied 10 ml syringe described in its Instruction for Use.

Instructions for Use

The CO-Set+ Syringe Assembly is designed for use with patient monitor and Swan-Ganz thermodilution catheters. **Please refer to the package insert provided with Models 93600, 93-610 or 93610 Closed Injectate Delivery Systems, or the Model 93505 Flow-Through Housing for detailed instructions for use.**


Figure 1. Model 93650 CO-Set+ Syringe Assembly

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Edwards Lifesciences, het gestyliseerde E-logo, Edwards en CO-Set+ zijn handelsmerken van Edwards Lifesciences Corporation; Edwards Lifesciences en het gestyliseerde E-logo zijn geregistreerd bij het Amerikaans "Patent and Trademark Office".
Edwards Lifesciences, det stiliserede E logo, Edwards og CO-Set+ er varemærker fra Edwards Lifesciences Corporation; Edwards Lifesciences og det stiliserede E logo er registreret i U.S. Patent and Trademark Office.
Edwards Lifesciences, den stiliserade E-logotypen, Edwards och CO-Set+ är varumärken som tillhör Edwards Lifesciences Corporation; Edwards Lifesciences och den stiliserade E-logotypen är registrerade hos U.S. Patent and Trademark Office.
Edwards Lifesciences, o E estilizado do logótipo, Edwards e CO-Set+ são marcas comerciais da Edwards Lifesciences Corporation; Edwards Lifesciences e o E estilizado do logótipo estão registados na Entidade de Marcas Registradas e Patentes dos EUA (U.S. Patent and Trademark Office).