STORZ

KARL STORZ — ENDOSKOPE

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Telescopes and HOPKINS® telescopes

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**Important information for
users of KARL STORZ devices
and instruments**

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This instruction manual is designed to provide help with the proper application, cleaning, and sterilization of the telescopes manufactured by KARL STORZ. All required details and all actions on your part are clearly presented and explained. We therefore ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for easy reference.

The manual also describes the recommended testing and reprocessing methods for telescopes.

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KARL STORZ is continuously working on the further development of all products. Please understand that changes to the scope of supply, design, equipment, and technology are possible for this reason. Therefore, no claims may be derived from the information, figures, and descriptions in this manual.

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






- ① Screw base for KARL STORZ fiberoptic light cables and Olympus/Winter & Ibe
- ② Plug-in base, 9 mm diameter, for Wolf fiberoptic light cables
- ③ Light inlet piece without adaptor
- ④ Color-code ring to indicate the telescope's direction of view:
Green = 0 degrees, straight forward
Black = 12 degrees/45 degrees
Red = 30 degrees, forward oblique
Yellow = 70 degrees, lateral
Blue = 90 degrees, lateral
White = 120 degrees, rear view
Gold = 0° – 120° (ENDOCAMELEON® size 10 mm)
Gold = 15° – 90° (ENDOCAMELEON® size 4 mm)
- ⑤ Eyepiece
- ⑥ Sheath
- ⑦ 'Autoclave' marking ring
- ⑧ Type designation and direction of view. The first letter after the number means:
A = 0/6 degrees, straight forward
B = 30 degrees, forward oblique
C = 70 degrees, lateral
D = 90 degrees, lateral
E = 120 degrees, rear view
F = 12 degrees/45 degrees
If present, the letter W stands for wide angle.
- ⑨ Lens and distal outlet of light cable
- ⑩ Focusing wheel for the eyepiece drive
- ⑪ DCI® connection
- ⑫ Adjustment wheel for direction of view
- ⑬ Rotating wheel for filter replacement (possible with telescopes for fluorescence diagnostics)
- ⑭ Illuminator (possible for VITOM® telescopes)
- ⑮ Irrigation port (LUER-Lock)
- ⑯ Working channel

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1	General information	3
2	Safety instructions.....	4
2.1	Symbols and hazard information	4
2.2	Safety instructions.....	4
2.3	Intended use	7
2.4	User qualification.....	7
2.5	Safety precautions when using the telescopes	7
3	Handling.....	8
3.1	Unpacking	8
3.2	Application	8
3.3	Inspection	9
3.4	Precautions.....	11
3.5	Instructions for use for accessories	12
3.6	Connecting to the camera system and light source	12
3.7	Light cable assignment	13
3.8	Focus and image size	13
3.9	Assembling the sheath and telescope	14
4	Reprocessing	15
4.1	Accessories	16
4.2	Disassembly.....	16
4.3	Preparation for cleaning and disinfection.....	16
4.4	Manual precleaning	17
4.5	Manual cleaning	17
4.6	Manual disinfection	18
4.7	Machine cleaning and disinfection.....	18
4.8	Connecting	19
4.9	Assembly, inspection and care	19
4.9.1	Removal of deposits	19
4.9.2	Checking the image waveguides.....	20
4.9.3	Checking the light cables	20
4.9.4	Drying the optical end faces	20
4.10	Packaging systems	21
4.11	Sterilization.....	21
4.12	Hydrogen peroxide (H ₂ O ₂) sterilization – ASP STERRAD®	22
4.13	Hydrogen peroxide (H ₂ O ₂) sterilization – STERIS® AMSCO® V-PRO™1 (only for telescopes without a channel)	23

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1 General information

HOPKINS® rod lens telescopes are reusable and consist of an eyepiece, a connection for light cables with screwed-on adaptors for light cables from other manufacturers, and a sheath made of rust-proof material that encloses the rod lens system and an integrated fiberoptic light cable.

HOPKINS® telescopes are designed to keep the strain on a patient to a minimum. They are put to use hundreds of thousands of times every year.

However, HOPKINS® telescopes are much more sensitive than the devices that a physician would usually use.

They must therefore be handled with care, and checked for correct functioning and for signs of damage before and after every use on a patient. Improper handling or incorrect usage may lead to tissue damage, infection, or instrument breakages.

Observing the following information will ensure that your HOPKINS® telescopes have a long service life.

i NOTE: The telescopes referred to in this manual usually have separate instruction manuals. Other general regulations (intended use, user qualification, safety rules etc.) must be taken from these separate instruction manuals and observed!

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2 Safety instructions

2.1 Symbols and hazard information

Please read this manual carefully and follow its instructions exactly. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the following text should be carefully reviewed to ensure the safe and effective use of this instrument. To make these words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or personnel may be endangered. Failure to observe a Warning could result in injury to the patient or personnel.

CAUTION: A Caution indicates that particular procedures or precautions must be followed to avoid possible damage to the instrument set.

NOTE: A Note indicates special information regarding operation or clarifies important information.

2.2 Safety instructions

WARNING: Risk of injury and damage to the products: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals carefully and always follow the instructions given precisely. Check the functioning of the products used in combination.

WARNING: Risk of injury: Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all of the accessories used in

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combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

WARNING: Risk of injury and risk of damaging the products: The highly concentrated light energy can cause the distal end, the light connection, adjacent components and tissue in front of the light emission window to heat up. This can cause burns to the patient, user and operating accessories.

Avoid direct tissue contact with the distal end, as well as with the light connection of the telescope and the light cable.

Do not rest the endoscope or the light cable on the patient or where it is in direct contact with operating accessories.

Always select the lowest possible light output of the cold light source which still allows optimal illumination of the operating field.

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WARNING: Risk of injury: Patient leakage currents may accumulate if endoscopes and endotherapy devices supplied with energy are used together. This is particularly important if an endoscope with an applied part of type CF is used. In this case, an endotherapy device with a type CF applied part must also be used in order to minimize the total patient leakage current.

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WARNING: Risk of fatal injury: In the case of patients with cardiac pacemakers, it must be ensured that the proximal end (= furthest from the patient) of the telescope is not positioned less than 2 cm from the pacemaker. The device generates a magnetic field at the setting ring, which can cause cardiac pacemakers or other electronic devices to malfunction.

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WARNING: Risk of injury: Incorrect application of medical instruments poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

NOTE: VITOM® telescopes are not endoscopes and must therefore only be used outside the body. The minimum working distance is 20 cm, and therefore significantly greater than with an endoscope. The minimum distance is set using a spacer bar (209180 20), which is available separately.

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2.3 Intended use

Indication: KARL STORZ HOPKINS® telescopes are used for visualization of the operating field during diagnostic and surgical endoscopic interventions.

VITOM® telescopes are used for the extracorporeal display of the operating field during diagnostic and therapeutic procedures in endoscopic and open surgery.

Contraindication: Use of KARL STORZ telescopes and HOPKINS® telescopes is contraindicated if, in the opinion of the responsible physician, such an application would endanger the patient.

Unauthorized conversions or modifications to the instrument are not permitted for safety reasons.

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The information given in this manual only serves to instruct in the correct handling, cleaning and sterilization of the telescopes. It is not intended as an introduction to the techniques of endoscopy.

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2.5 Safety precautions when using the telescopes

The telescopes must only be used according to the medical rules and procedures of endoscopy recognized for the method.

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3.2 Application

WARNING: The telescope is not sterile when delivered and must therefore be cleaned, disinfected and/or sterilized prior to initial use and each subsequent reuse.

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WARNING: Patient leakage currents may accumulate if endoscopes with powered accessories are used. This is particularly relevant when using CF devices. In such a case, the other devices used should also be of type CF in order to minimize the total patient leakage current.

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3.3 Inspection

Inspect the telescope for damage immediately after receipt and before and after each use.

Image waveguide

Examine the optical end faces (distal end and eyepiece) for scratches and surgical or cleaning residues. The faces should be smooth and shiny.

In order to test the image quality, slowly rotate the telescope while looking through it.

If the image is completely or partially impaired, either the rod lens system in the sheath is broken or the lens in the eyepiece is defective. In such cases, the telescope must be replaced.

A cloudy or speckled image can be caused by moisture or disinfectant residue on the optical end faces.

NOTE: When looking through HOPKINS® PDD telescopes, the image appears yellowish. This effect is normal and is caused by the use of special filters.

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Light cable

Hold the distal end in the direction of a light source and inspect the number of dark dots on the light cable connection. Such dark dots indicate that the optical fibers of the light conducting bundle are broken. Isolated broken optical fibers do not noticeably compromise the quality of the image.

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However, the light output does significantly deteriorate if a certain number are broken (from approx. 20 % to 25 %), which therefore also results in reduced image quality. The telescope must not be used any longer.

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
Sheath

Inspect the entire sheath for dents and scratches. The sheath must not be bent. LUER-Lock stopcocks, telescope locks and other movable parts should be smooth-running and undamaged. Check plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their original condition, they must be replaced. This applies in particular to seals.

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3.4 Precautions

As with all optical instruments, the HOPKINS® telescopes are a very sensitive part of the endoscopic equipment. In order to ensure the long-term quality of the telescope, the following information must be observed.

- In order to avoid breaks, telescopes must always be held by the eyepiece and never by the distal end alone.
- Do not bend the sheath under any circumstances. This can result in breaks or cracks in the rod lens system.
- Place the telescopes down with care. Hard impacts, particularly to the distal end, can cause damage to or cracks in the luting and allow liquid, steam, etc. to penetrate. This results in cloudy and/or blurred image areas.
- Therefore, during cleaning, disinfection, sterilization and storage, it must be ensured that the telescopes cannot be damaged by other instruments. Ideally, you should store HOPKINS® telescopes individually or use containers in which the telescopes can be secured (see accessories).

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Special protective sheaths are available in different lengths for the HOPKINS® telescopes. These sheaths are suitable for storage and transport, as well as for sterilization (see accessories).

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CAUTION: HOPKINS® telescopes that do not have the "Autoclave" identification ring must not be steam sterilized. This may cause irreparable damage.

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3.5 Instructions for use for accessories

Before use, make sure that the telescope has been correctly cleaned, disinfected, and sterilized where necessary.

Check that the telescope and accessories are all functioning correctly before use.

Operate the video and light source equipment as per the specifications in the relevant manual.

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To prevent the telescopes from fogging during use, we recommend the use of a warmer for rigid telescopes, Art. no. 10905 B (see page 35).

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3.6 Connecting to the camera system and light source

The camera heads generally have an integrated instrument coupler.

Squeeze the levers (① and ②) of the instrument coupler together and insert the endoscope's eyepiece into the coupler. The coupler automatically locks into place when the eyepiece is inserted.


The endoscope can be locked by pressing lever ② to lever ③.

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Connect the light cable to the outlet of the light source and to the endoscope (twist the knurled screw on the screw base by a quarter turn).

WARNING: Reduce the intensity of the light source so that optimal illumination of the operating field is still achieved. The higher the intensity of the light source, the greater the heat emission at the distal end of the telescope.

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WARNING: Never place the end of the fiberoptic light cable or a telescope on or underneath a sterile cover cloth. High light intensity can cause burns to the patient and/or the sterile cover cloth.

WARNING: Danger of glare! Never look into the free end of a light cable.

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3.7 Light cable assignment

If there is an identification symbol on the light cable and the endoscope (see table below), assign matching combinations according to the symbols shown*.

Light cable diameter	Endoscope diameter
4.8 – 5.0 mm 	6.6 – 12.0 mm 
3.0 – 3.5 mm 	3.0 – 6.5 mm 
2.0 – 2.5 mm 	0.8 – 2.9 mm 

* Special endoscopes may deviate from this

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
3.8 Focus and image size

Operate the video equipment as per the specifications in the relevant manual. Perform a new white balance as required, or when replacing the telescope.

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3.9 Assembling the sheath and telescope

Sheaths, sheath inserts and telescopes are provided with a 0 mark or a line mark.

Open the telescope lock (1) by turning it counterclockwise (as far as it will go, marks in alignment (2)).

When assembling, make sure that the marks are aligned with one another (line mark and 0 mark (3) are lined up). Insert the telescope into the sheath lock and close the telescope lock (1) by turning it in a clockwise direction.

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4 **Reprocessing**

WARNING: Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out, or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

CAUTION: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: Danger of damage to medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices. Only use chemicals approved by KARL STORZ for reprocessing. A complete list can be found on the internet at www.karlstorz.com.

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4.1 Accessories

Accessories required for carrying out reprocessing:
Brushes:

Surface: 27652

Inside lumina 2.5 – 5 mm: 27650 B

Inside lumina 2 – 2.4 mm: 27650 C

4.2 Disassembly

Prior to cleaning and disinfection, the telescope must be separated into its individual components as far as possible. Both adaptors on the light inlet piece (Fig. A), and any adaptors (Fig. B), seals on the working channel and LUER-Lock stopcocks (Fig. C) must be disassembled.

4.3 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the telescope immediately after use. To this end, pre-clean the telescope surface by wiping down and rinsing, for example, and rinse the working channel if necessary. As a general rule, KARL STORZ recommends manual pre-cleaning under cold running water.

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4.4 Manual precleaning

Brushing the surfaces

Depending on the degree of soiling and the degree to which residues have dried, heavy soiling must be cleaned from the surfaces of instruments with the aid of a brush (Art. no. 27652)/sponge under cold running water until any visible contamination has been completely removed (Fig. A).

Brushing the lumina

In order to remove all visible soiling, the working channel (Fig. B), lumina and hollow spaces (e.g. disassembled LUER-Lock stopcocks, Fig. C) of the telescope must be precleaned under cold running water with the aid of a suitable brush. The brush with Art. no. 27650 B is suitable for the LUER-Lock stopcock.

No pretreatment with ultrasound

For technical reasons, the telescope is **not** suitable for ultrasound treatment.

4.5 Manual cleaning

The telescope and its dismantled individual components must be completely immersed in a cleaning solution. By filling the lumina in a targeted manner, it must be ensured that even surfaces with restricted access are covered and no air bubbles are present. After the necessary exposure time, clean the instrument mechanically with the aid of a brush or sponge. Finally, it must be rinsed with cold water to ensure neutralization.

The dismantled individual components must be completely immersed in a cleaning solution. It must be ensured that even surfaces with restricted access are covered and no air bubbles are present. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a cloth or sponge. Finally, it must be rinsed with cold water to ensure neutralization.

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4.6 Manual disinfection

The telescope and its dismantled individual components must be completely immersed in a disinfectant solution. Through disassembling and opening the instrument, and filling the lumina in a targeted manner, it must be ensured that even surfaces with restricted access are covered and no air bubbles are present. At the end of the necessary exposure time as specified by the manufacturer, the medical device must be rinsed several times in order to remove all chemical residues. To this end, water of the best possible quality must be used, taking into account the national regulations. Finally, all of the surfaces, joints, openings, channels and lumina are dried completely with (preferably sterile) compressed air in accordance with the national regulations (Fig. A). The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

4.7 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A_0 value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

NOTE: If necessary, the instrument must be dried off afterwards by hand.

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4.8 Connecting

In order to guarantee effective machine cleaning and disinfection, the instrument must be connected up to the washer and disinfecter (e.g. via adaptor), or it must be ensured that it is thoroughly rinsed out by means of appropriate rinsing lances (Fig. A).

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Put a little cleaning paste on a clean, moistened cleaning rod and apply to the cover window (Fig. B). Polishing removes stubborn deposits which can be subsequently rinsed off with water.

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4.9.2 Checking the image waveguides

Examine the optical end faces (distal tip and eyepiece) for scratches and surgical or cleaning residues. The faces should be smooth and shiny. In order to test the image quality, slowly rotate the telescope while looking through it. If the image is completely or partially impaired, either the rod lens system in the sheath or the lens in the eyepiece is defective. In such cases, the telescope must be replaced.

A cloudy or speckled image can be caused by moisture or disinfectant/chemical residue on the optical end faces. In these cases, carefully clean the optical end faces with a cotton tip applicator soaked in 70 % alcohol (Fig. A).

4.9.3 Checking the light cables

Hold the distal end in the direction of a light source and inspect the number of dark dots on the light cable connection (Fig. B). Such dark dots indicate that the optical fibers of the light conducting bundle are broken. Isolated broken optical fibers do not noticeably compromise the quality of the image. However, the light output does significantly deteriorate if a certain number are broken (from approx. 20 % to 25 %), which therefore also results in reduced image quality. The telescope must not be used any longer.

WARNING: If damage is apparent, the telescope should no longer be used. Telescopes must be replaced if the image is cloudy or if no image or only parts thereof are visible.

4.9.4 Drying the optical end faces

The fiber areas in the light inlet piece in particular as well as the distal and proximal end faces must be dried carefully with 70 % alcohol (Fig. C and D). Disinfectant and cleaning agent residue, e.g., in the light inlet piece, can become burnt on when the light cable is connected, thereby seriously impairing light transmission.

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- Damaged or corroded medical devices must be withdrawn from use.
- Threads and glide surfaces on the telescope must be lubricated locally with instrument oil (Art. no. 27656 B) and the running surfaces on the cock plug with special grease (Art. no. 27657) (Fig. A).
- Dismantled medical devices must be assembled.
- Afterwards, a functional check must be carried out.

- ❗ **NOTE:** The oil used for this purpose must be suitable for the subsequent sterilization procedure (silicone-free and paraffin- or white oil-based).
- ❗ **NOTE:** During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

4.10 Packaging systems

Only standardized and approved packaging materials or systems may be used (EN 868 Parts 2-10, EN ISO 11607 Parts 1 + 2, DIN 58953, Fig. B).

4.11 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The method must be selected taking into account the respective national requirements and in consultation with the device manufacturers.

Telescopes with a ring marked 'AUTOCLAVE' can be steam sterilized without restrictions (Fig. C).

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

Steam sterilization using the fractionated prevacuum procedure

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) at 132 °C – 137 °C with a minimum exposure time of 3 minutes and a maximum of 18 minutes. This procedure is only suitable for thermostable instruments.

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CAUTION: Do not open sterilization containers immediately after the sterilization process has finished. Never pour cold sterile water over the telescopes. If the cooling process is accelerated, the telescopes may be damaged.

NOTE: Greased LUER-Lock stopcocks must be dismantled for sterilization so that steam penetration can be ensured.

Steam sterilization is not permitted for the following instrument groups. Only those methods with a maximum temperature of 65 °C may be used, such as gas and plasma sterilization:

- Non-autoclavable telescopes (cleaning up to 93 °C)
- Non-autoclavable miniature endoscopes (cleaning up to 65 °C)
- Magnifiers with drive mechanism
- Non-autoclavable prismatic light deflectors

4.12 Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

WARNING: Medical devices must be dismantled for sterilization.

WARNING: Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

CAUTION: For more detailed information, please consult the user handbook of the respective device.

NOTE: Refer to the 'STERRAD® Sterility Guide' to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

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 **NOTE:** Sterilization is not possible on surfaces which have been greased and lubricated.

The following STERRAD® sterilization processes have been validated and approved by KARL STORZ for telescopes **without** a channel:

- STERRAD® 50, 100S, 200
- STERRAD® NX®
- STERRAD® 100NX®
- STERRAD® 100NX® DUO Cycle

The following STERRAD® sterilization processes have been validated and approved by KARL STORZ for telescopes **with** a channel:

- STERRAD® NX®
- STERRAD® 100NX®

4.13 Hydrogen peroxide (H₂O₂) sterilization – STERIS® AMSCO® V-PRO™1 (only for telescopes without a channel)

Detailed information on the selection of the appropriate cycle for the various device generations is available from the manufacturer STERIS®.

 **NOTE:** Sterilization is not possible on surfaces which have been greased and lubricated.

4.14 Ethylene oxide sterilization (EO)

The ethylene oxide procedure is validated with 100 % ethylene oxide at 55 °C and a hold time of 30-45 minutes.

 **WARNING:** Medical devices must be dismantled for sterilization.

 **CAUTION:** For gas sterilization with ethylene oxide, the airing times stipulated by the device manufacturer must be observed for instruments, due to gas absorption by the materials. The airing times depend on the processing technique of the EO sterilizer (concentration, process control).

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i **NOTE:** Sterilization is not possible on surfaces which have been greased and lubricated.

4.15 Chemical low-temperature sterilization with peracetic acid – STERIS® System 1® (only for telescopes without a channel)

More detailed information on the selection of the sterilization parameters and the STERIS® Quick Connect Kit (QCK) required to rinse out the lumina is available from the manufacturer STERIS®.

i **NOTE:** Sterilization is not possible on surfaces which have been greased and lubricated.

4.16 Chemical low-temperature sterilization with peracetic acid – STERIS® System 1E® (only for telescopes without a channel)

More detailed information on the selection of the sterilization parameters and the STERIS® Quick Connect Kit (QCK) required to rinse out the lumina is available from the manufacturer STERIS®.

i **NOTE:** Sterilization is not possible on surfaces which have been greased and lubricated.

The following procedures are approved in respect of material compatibility, but must be validated by the on-site operator with regard to efficacy:

Low-temperature steam and formaldehyde process (LTSF)

i **NOTE:** Sterilization with FO is tested and approved by KARL STORZ in respect of material compatibility. However, due to the large number of processes used in the global market, validation in respect of efficacy was not carried out. Therefore, the user must validate his specific processes individually in collaboration with the device manufacturer. KARL STORZ has had

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exclusively the aforementioned sterilization methods validated with regard to a reduction in microorganisms of the order of 10^{-6} SAL.

NOTE: Sterilization is not possible on surfaces which have been greased and lubricated.

4.17 Reprocessing limits

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

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- This reprocessing step is part of the validated process
 - This reprocessing step may be carried out additionally in order to aid cleaning
 - x On the grounds of potential damage, this reprocessing step must NEVER be performed
- Fields which are not filled are not part of the performed validation.

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5 Servicing and repair

Defective instruments must be serviced and repaired exclusively by persons authorized by us; all repair work must employ genuine parts only.

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6 Limitation of liability

KARL STORZ SE & Co. KG shall be liable for the failure or a deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the conditions that all assembly, operations, system expansions, readjustments, modifications, or repairs to the same have been performed by persons authorized by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

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8 Directive compliance

This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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9 Repair program

KARL STORZ maintains a repair and replacement warehouse which is normally adequate to ensure prompt replacement of damaged telescopes and instruments. Under the replacement plan, you will receive a new replacement telescope and will only be charged the repair costs for the defective telescope.

In Germany you can refer repairs direct to:

KARL STORZ SE & Co. KG
Service department
Logistics center
Take-Off Gewerbepark 83
78579 Neuhausen ob Eck, Germany

In other countries please contact your local KARL STORZ subsidiary or authorized dealer.

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10 Important information

For hygiene reasons and to prevent infection, telescopes, instruments and devices must be cleaned and disinfected/sterilized before they are sent for repair. We reserve the right to return contaminated instruments/devices to the sender. Repairs, modifications, or expansions which are not performed by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

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**Important information
for users of KARL STORZ
instruments**

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This manual is intended to serve as an aid in the proper use, cleaning, disinfection and sterilization of video endoscopes manufactured by KARL STORZ. The video endoscopes are sensitive precision instruments, which require particular care during preparation. To prevent damage during preparation, guidance is provided here for the disassembly, cleaning, sterilization and assembly of KARL STORZ video endoscopes. For this reason, read these instructions carefully and keep them as a ready reference.

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



- ① Deflection up/down
- ② Up/down brake
- ③ Programmable remote control buttons
- ④ Suction valve
- ⑤ Irrigation/insufflation valve
- ⑥ Working channel with rubber lip valve/biopsy
valve
- ⑦ Deflection right/left
- ⑧ Left/right brake
- ⑨ Working channel
- ⑩ Nozzle for lens cleaning/irrigation/insufflation
- ⑪ Fiberoptic light carrier
- ⑫ Lens
- ⑬ Water jet channel (only for 13806 NKS/PKS,
13807 NKS/PKS, 13910 NKS/PKS and
13911 NKS/PKS)
- ⑭ Supply tube
- ⑮ Connection for leakage test
- ⑯ LUER lock (connector for jet nozzle, only for
13806 NKS/PKS, 13807 NKS/PKS,
13910 NKS/PKS and 13911 NKS/PKS)
- ⑰ Connection for water bottle
- ⑱ Light source/GASTRO PACK™ light cable
connection
- ⑲ Port for video processor/electronic connector
- ㉑ Protection cap for video processor
connection port
- ㉑ Suction/aspiration port
- ㉒ Air pin

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**Important information for users of
KARL STORZ instruments II**

Pictures of the instruments.....	III
Controls, connectors and their uses.....	V

Safety instructions 3

Intended use	3
User qualification	4
Safety precautions when using the video endoscopes	4
Warnings and cautions.....	5

Handling 10

1.0 Scope of supply	10
Compatible video processor units/light sources.....	11
2.0 Operating the video endoscope.....	13
2.1 Preparation of the light source.....	13
2.2 Preparation of the processor unit.....	13
2.3 Connecting the light conductor, air/water connections and the suction channel.....	13
2.4 Adjusting the angle of the distal tip	14
2.5 Suction	14
2.6 Insufflation and light.....	14
2.7 Water rinsing the lens	15
3.0 Preparation before use	15
3.1 Testing the various functions	15
3.2 Holding the endoscope correctly	18
3.3 Additional recommendations for use.....	19
3.4 Operating the endoscope buttons	20
3.5 Disassembling and assembling the Albarran unit (only video duodenoscopes).....	21

Preparation and care24

General information	24
General cleaning.....	26
Preparation for cleaning.....	27

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Please read these safety instructions carefully before using the instrument. Make certain you are well acquainted with the function and operation of the instrument before initial use on the patient.

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Intended use

The KARL STORZ video endoscopes (13800 xxx, 13801xxx, 13804xxx, 13806xxx, 13807xxx) have been developed for endoscopic diagnosis and treatment in the upper digestive tract, including the oesophagus, the stomach and the duodenum.

The KARL STORZ video duodenoscopes (13880 xxx – 13882 xxx) were developed for the endoscopic diagnosis and therapy of the biliopankreatic system.

The KARL STORZ video rectoscopes (13912xxx), video sigmoidoscopes (13900/13901 xxx) and video colonoscopes (13904xxx - 13907 xxx and 13910/13911 xxx) have been developed for endoscopic diagnosis and therapy in the lower digestive tract, including the rectum and colon.

In addition, these video endoscopes are used for endoscopic diagnosis and therapy.

Videoscopes are applied components, which should only be used with approved supply units (video processor units, light sources).

Contraindications:

No contraindications are known which directly relate to the product. The video endoscope must not be used if, in the opinion of a qualified physician, such an application would endanger the patient, e.g., due to the patient's general condition, or if the endoscopic method as such is contraindicated.

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User qualification

KARL STORZ video endoscopes may only be used by persons with an appropriate medical qualification and who are familiar with the endoscopic technique. The information given in these instructions only serves to instruct in the correct handling and preparation of the video endoscopes.

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Safety precautions when using the video endoscopes

Video endoscopes must be used according to the medical rules and procedures of endoscopy recognized for this method.

NOTE: The patient leakage current limit of the CF application part may be exceeded when CF application parts (application on heart) are used in combination with other BF application parts.

In this case, the reduction of the patient leakage current should be effected using only CF application parts whenever this is possible.

NOTE: The endoscope does not offer any electric insulation against high frequency (HF) currents. Protection against HF currents should be ensured through the use of correct HF instruments. Please read the instructions for use in the instruction leaflet supplied with the HF application part for more information.

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Warnings and cautions

Please read this manual and follow its instructions carefully, to avoid putting your patients, personnel or yourself at risk. The words, **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the instrument. To make the words **Warning** and **Caution** stand out more clearly, they are accompanied by a pictogram.

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WARNING: Check the video endoscope prior to each procedure to ensure that it functions correctly and has been properly cleaned and disinfected.

WARNING: Verify the sound condition of the video endoscope, as well as of accessories used in combination, before every application. Damaged video endoscopes or damaged accessories must not be used.

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WARNING: Remove the EO/pressure compensation cap 11025 XE from the connection piece:

- during every endoscopic intervention
- during manual/machine preparation.

WARNING: Do not allow any electrically active components or instruments to come into contact with the video endoscope.

WARNING: The video endoscope must not be used during the discharge of a defibrillator. Remove the video endoscope from the operation area before discharge occurs.

WARNING: Danger of explosion if flammable gases are used in the immediate proximity of the video endoscope. The video endoscope must not be used in the presence of flammable or combustible materials.

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WARNING: High energy light radiated through endoscopes may give rise to high temperatures in front of the light outlet, at the end faces of the light cable and at the tip of the endoscope.

Prolonged direct contact of light-emitting surfaces on tissue can lead to burns, this must therefore be avoided.

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WARNING: Before removal from the patient, the distal tip of the video endoscope must be returned to the neutral, nondeflected position (L and U (demarcations on the wheels of the deflection mechanism) on top of each other and positioned between the valves) to prevent damage to the instrument and/or injury to the patient.

WARNING: The video endoscope may only be repaired by personnel authorized by KARL STORZ. Unauthorized repairs may impair functioning and safety, and are therefore strictly forbidden. KARL STORZ gives no warranty for video endoscopes that are repaired by unauthorized personnel.

WARNING: Damage to the video endoscope and injury to the patient or user resulting from incorrect operation are not covered by the manufacturer's warranty.

WARNING: Always unplug the video endoscope before carrying out any maintenance and cleaning procedures.

WARNING: Risk of infection due to shipment! Contaminated video endoscopes can cause serious or fatal diseases during transportation or at the recipient's premises. Please cover the carrying case with the protective sheet 13990 SFN before sending contaminated instruments back to KARL STORZ.

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WARNING: If a malfunction of the endoscope should occur during procedure on a patient, stop the procedure immediately. Place the distal tip of the video endoscope in a neutral, nondeflected position and remove the video endoscope slowly and carefully from the patient.

WARNING: The use of HF probes, HF electrodes or other energy-emitting probes is not permitted if explosive gases (for example, explosive/flammable anesthetics) are present.

WARNING: Patient leakage currents can accumulate if endoscopes and electrically powered endoscope accessories are used together.

WARNING: Before turning on the high voltage/ electrical supply, ensure that there is no contact with metal parts of the endoscope or any other electricity-conducting accessory, and that the suction unit has no contact with conductive liquids.

WARNING: HF probes, HF electrodes, laser probes and other energy-emitting probes should only be used under constant, visual, endoscopic control.

Activation should only occur when the tip of the energy-emitting probe is outside the working channel and within the endoscope's field of vision.

CAUTION: Only use accessories (guide probes, cleaning brushes, etc.) that are specifically suited for the intended purpose for use with video endoscopes and which can be employed in line with generally accepted technical standards and safety guidelines.

CAUTION: Only use accessories with flexible ends. Other types of accessories may cause damage to video endoscopes.

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CAUTION: Place the EO/pressure compensation cap 11025 XE on the connection piece prior to:

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CAUTION: Make certain to read the 'Cleaning, sterilization and care of KARL STORZ instruments' manual. Methods of cleaning, disinfection and gas sterilization are explained in detail there.

CAUTION: Do not hit the distal tip of the video endoscope against hard objects, and never use a sharp object to remove debris.

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NOTE: Any damage to the video endoscope resulting from incorrect operation is not covered by the warranty.

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KARL STORZ light sources Xenon 100 or Xenon 100 SCB/ Gastropack 2.0 Operating the video endoscope

2.1 Preparation of the light source

We recommend the use of our light source (Xenon 100 or Xenon 100 SCB) with a 100 W xenon lamp and our portable package solution the Gastropack with economical, light-intensive 24 W Solarc® lamp. Please read the full instructions in the instruction manual of the respective light source.

- 1 Install the light source on a solid, even surface.
2. Insert the power cable into the socket provided.

2.2 Preparation of the processor unit

Refer to the relevant instruction manual: TELECAM® SL II 202130 20, GASTRO PACK® 25 0430 11/ 25 0431 11 or Image 1 hub 22200020-1

2.3 Connecting the light conductor, air/ water connections and the suction channel

1. Insert the supply connector into the light source until it audibly clicks into place (Fig. 1).
2. Connect the video connecting cable to the endoscope. Take note of the white markings on the endoscope and video connecting plug. (Fig. 2).
3. Attach the connector from the water bottle tube to the endoscope supply plug (irrigation bottle connector, Fig. 3). The water bottle holder automatically positions itself correctly and cannot be adjusted manually.
IMPORTANT: To fill the water bottle with sterile water, remove the bottle from the light source and unscrew the lid. Only fill the water bottle 2/3 full with sterile water (1/3 is required to cushion air pressure).
4. Connect the suction tube (patient tube) from the suction pump or central suction system to the supply plug connector.

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2.5 Suction

Switch on the suction pump. To aspirate, press the suction valve (with the red ring) on the endoscope (Fig. 5).

Suction can also be used when the biopsy forceps are placed in the biopsy channel and the rubber lip valve is mounted. During the procedure, the rubber lip valve/biopsy valve seals the forceps channel.

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2.6 Insufflation and light

Start the light source, switch on the air pump.

For insufflation, cover the opening on the blue valve with a finger (do not depress the valve, since this starts water irrigation) (Fig. 6).

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2.7 Water rinsing the lens

Pressing on the air/water valve starts the water flow and/or lens cleaning. Different valve positions allow various combinations (Fig. 7,8).

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3.0 Preparation before use

WARNING: Before being used for the first time, the device must undergo a complete preparation cycle.

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3.1 Testing the various functions

1. Prepare light sources as described in 2.1.
2. Connect supply connector, water bottle and suction tube according to 2.3.
3. Check the supply cable for damage (e.g., breaks, cracks, twisted or squashed areas...) and the supply connectors for loose connections.
4. Examine the surface of the flexible insertion sheath by eye and by touch for cracks, cuts, bulges, foreign bodies or other damage. Do not press too hard on the deflection section.

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- 5a. Testing deflection – Upward/downward control
The large deflection wheel controls up and down movements. To test up and down movements, release the up and down brake. Slowly rotate the large deflection wheel clockwise, the distal end should move downwards. Rotating the deflection wheel counterclockwise causes upward movement.

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- 5b. Testing deflection – Right/left control
The smaller/outer deflection wheel controls right/left movement. To test right/left movement, release the right/left brake. Slowly rotate the small deflection wheel clockwise, the tip of the endoscope should move to the right. Rotating the deflection wheel counterclockwise causes movement to the left.
- 5c. Testing deflection – brake function
The deflection wheels should turn freely when the brakes are released. Release the brakes and test their action (see 5a/b). Activate the brakes by moving the brake lever counterclockwise – ensure that the angled position of the distal endpiece is fixed in the desired position.
Do not twist or bend the deflection piece by hand.

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WARNING: If the deflection and/or brake mechanisms do not function properly, do not use the instrument. Contact KARL STORZ service department.

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6. Check the air supply. Switch on the air pump and hold the distal end of the instrument under sterile water. Cover the air control opening in the air valve (blue ring) with a finger; bubbles should emerge from the air nozzle at the distal end.

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7. Check the water supply: Depress the air/water valve (blue ring) fully. Water should flow out of the water nozzle (distal) and rinse the lens clean. **IMPORTANT:** When you first actuate the valve, it may take a couple of seconds for water to start flowing.

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8. Test the suction function. With the unit correctly installed and the suction pump switched on, immerse the distal end in a container of sterile water once again. Depress the suction valve (red ring) fully. Water should be suctioned through the suction channel of the instrument via the connected suction pump and into the collecting bottle of the suction equipment.

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CAUTION: Accessories such as biopsy forceps that are kinked or bent will not function properly, are difficult to guide through the channel and may cause severe damage to the channel. They must therefore be replaced and not used further!

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Holding the endoscope correctly

With the control handle of the instrument in the left hand, the thumb, index and middle fingers can easily control the large deflection wheel. In this position, the suction valve can be depressed with either the index or the ring finger. The air/water valve can also be covered with these fingers.

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3.3 Additional recommendations for use

1. Use a lubricant from the KARL STORZ list of approved and recommended lubricants on our website: www.karlstorz.com under 'Hygiene' for the distal end of the insertion tube.
2. Observe the lumen the whole time. If you become disoriented, withdraw the instrument gently until vision is restored.

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3.4 Operating the endoscope buttons

Both buttons have two functions, which are initiated through short and long presses of the key respectively. Pressing and holding down both buttons at the same time for approx. 3s opens a settings menu and allows you to select the button assignment freely (i.e., Freeze/Zoom or Print (ACC1)/contrast). For more detailed information on button assignment please consult the camera control unit manual TELECAM® SL II, GASTRO PACK™ or Image 1 hub.

The freeze mode can be canceled by a short press on the freeze button.

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NOTE: In case of 'Red-out', a blurred red picture on the screen, caused by the distal endoscope tip resting directly on the mucosa, withdraw the instrument gently until vision is restored.

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NOTE: If the angling of the tip can no longer be properly controlled the examination must be ended immediately. Release the brakes and place the tip in the neutral position. Position the deflection wheels so that 'U' and 'L' are between the valves (see page 13/15), then carefully withdraw the instrument back, whilst continuing to monitor the procedure.

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3.5 Disassembly and assembly of the Albarran unit (video duodenoscopes only)

As conventional duodenoscopes are difficult to clean and maintain due to their design, particularly in the region of the Albarran unit, KARL STORZ has developed a modular Albarran system.

The patented Albarran module is easy to disassemble, which facilitates thorough cleaning and disinfection. In addition, the Albarran module can be autoclaved. It is no longer necessary to disassemble the entire endoscope to eliminate technical malfunctions in the Albarran unit. As a rule, in this case, repairs can be carried out by disassembling the Albarran unit.

Disassembling the Albarran unit

1. Remove the cable as shown in Fig. 1 by turning the rotary knob on the proximal endoscopic housing. Unscrew the rotary knob (2) from the proximal housing counterclockwise.
2. Remove the collet chuck (11), which is now visible (Fig. 2).
3. Unlock the Albarran module from the distal endoscope head (Fig. 3).
The quick locking screw (6) of the Albarran module unit is indicated by two red dots. The red dots are turned counterclockwise towards each other using a special assembly/ disassembly key (7) Art. no. 13991 B. The Albarran module is now unlocked (Fig. 4).
4. Removal of the Albarran module with Albarran wire (Fig. 5).
Separate the Albarran module (3) from the endoscope with your thumb and carefully remove the wire (4) from the endoscope.

NOTE: Do not bend or kink the wire!

5. Removal of the Albarran wire sheath (Fig. 6)
Angle the distal tip of the endoscope by turning the small deflection wheel clockwise ('R').

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
3. Attach the Albarran module (Fig. 8) by turning the quick locking screw (6) clockwise with the assembly/disassembly key (7).

The two red dots must now be opposite one another (Fig. 9).

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
4. Testing the module locking:
Locking of the Albarran module must be tested by trying to press it off the head of the endoscope with your thumb.

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
5. Place the shaft of the endoscope in a straight position on a table.
6. Lock the control slide (12) into the endoscope housing as per Fig. 11. At the same time move the Albarran control lever (13) backward as far as it will go.

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
7. Introduce the collet chuck (11) into the proximal housing with a twisting action until it engages and cannot be turned any further (Fig. 11).

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8. Screw on and tighten the rotary knob (2) as seen in Fig. 12.
9. Testing Albarran function: Grip the Albarran lever at the distal end with your thumb and manipulate the proximal Albarran control 13 simultaneously (Fig. 12).

You should be able to feel a significant movement pressure of the Albarran lever with your thumb. If so, the Albarran unit is functioning correctly.

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NOTE: This documentation covers only the handling and care of the Albarran control lever.

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Cleaning and disinfecting the Albarran module, see page 40.

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General information

WARNING: The stated cleaning, disinfection and sterilization methods alone do not guarantee successful sterilization of the endoscopes. This can only be achieved if personnel follow recognized and validated preparation methods.

WARNING: The local occupational safety regulations must be followed and applied for all cleaning and disinfection work on contaminated endoscopes and instruments.

WARNING: Remember that sterilization is only then successful if

- the instruments have been freed of all organic material and chemical residues.
- work proceeds according to recognized and validated methods.
- any deviation from recommended procedures must be validated and tested by the user.
- the recommended sterilization parameters are maintained. These are only valid with cleaning, disinfection and sterilization equipment that is properly maintained and calibrated.

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CAUTION: The manual 'Cleaning, sterilization and care of KARL STORZ Instruments' must be followed. There, the cleaning, disinfection and sterilization procedures are explained in detail. The list of KARL STORZ approved materials for the preparation of video endoscopes can be found on the Internet at 'www.karlstorz.com' or 'www.karlstorz.de' under 'Hygiene' and also in the appendix of this manual.

CAUTION: When preparing and using the solutions, follow the manufacturer's instructions paying close attention to proper dilution and immersion times. Disregarding the instructions can lead to damage. Bear in mind the microbiological range of action of the chemicals used.

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CAUTION: Never clean video endoscopes with metal brushes.

CAUTION: Constantly changing between different preparation methods causes particular wear on the materials and must therefore be avoided. A once-only change in methods, however, e.g., when acquiring new cleaning equipment, should be carried out in consultation with the manufacturers of the endoscope, the chemicals and the washer/ disinfectant.

The method change is then, as a rule, completely harmless.

CAUTION: Do not clean video endoscopes in an ultrasound bath.

CAUTION: Do not steam sterilize (autoclave) video endoscopes.

CAUTION: Never immerse video endoscopes in normal saline solution, since even brief contact can lead to pitting and corrosion.

CAUTION: For manual preparation, video endoscopes may only be rinsed with microbiologically clean/sterile water.

CAUTION: Decontamination in washers and disinfectors is only possible using valid procedures.

CAUTION: Special attention should be given to the following sterilization conditions:

Temperature:	5°C – 65°C
Humidity:	15% – 100%

CAUTION: Even parts that do not come into direct contact with the patient (e.g., video connecting cables) should be thoroughly disinfected following an intervention.

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General cleaning

Cleaning should only be performed by staff who are appropriately trained. All safety regulations should be observed (gloves, protective glasses, mouth and nose protection, clothing, TBA 250...).

Thorough cleaning of the endoscope and all its accessories after every use ensures problem-free working and the best results for disinfection or sterilization procedures.

The basic guidelines and recommendations for the preparation of instruments (for example RKI, ESGE/ESGENA) should be observed and the effectiveness verified.

CAUTION: Before immersion in cleaning solution, check the instrument for leaks. Otherwise fluid penetrating the video endoscope could cause damage.

CAUTION: Do not autoclave the endoscope or accessories (with the exception of accessories that are specifically designated as autoclavable, e.g., forceps or Albarran module).

NOTE: The endoscope must be cleaned immediately after every use. If biological material is allowed to remain and dry in the instrument, complete or partial blockage of the channel may result. Thorough cleaning is a prerequisite for effective disinfection and sterilization. This will avoid clogs/ blockages and the resulting repairs and instrument failures. For reasons of material compatibility, only cleaning and disinfection chemicals approved by KARL STORZ should be used (see current list in the annex, p. 52 or on our website www.karlstorz.de under 'Hygiene').

NOTE: The endoscope must undergo a complete preparation cycle prior to the initial use.

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Preparation for cleaning

1. The preliminary cleaning of the endoscope is performed at the place of operation. The main cleaning and disinfection is carried out in a separate preparation room.
Immediately after use, remove any heavy soiling (secretions, etc.) and organic material from the outer surfaces of the insertion tube with a disposable cloth moistened with active cleaning solution. Then submerge the distal end in a container with active cleaning solution and activate the suction to rinse the working channel.
2. Put the cleaning valve into the air/water housing. (Do not insert this valve during normal use!)
3. Check that the water bottle is filled (max 2/3). Start the air supply (insufflation pump). Water and air are now passed through both the air channel and the water channel. (The process should run for at least 10 seconds).
4. Switch off video system.
5. Pull out the video plug from the video connection cable at the supply connector of the endoscope when the video processor unit is switched off. Put on the sealing cap for the video cable connector and screw it on tightly.
6. Remove the endoscope from the video processor unit and light source. Place the cleaning valve in an active cleaning solution.

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Leakage test

Before every immersion in liquid, test the endoscope for leaks.

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NOTE: Ensure the video cable connector sealing cap is secured tightly (blue rubber inlet seal should be slightly moistened with instrument oil from time to time).

2. The connection tube of the leakage tester and the ventilation connector on the endoscope must be dry. Connect the leakage tester connector to the supply connector and lock by turning clockwise 90°.

NOTE: The rubber seal should also be treated occasionally with instrument oil.

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NOTE: During inflation, make sure that the red pressure relief button is completely closed.

3. Inflate the leakage tester until the maximum test pressure of 160 - 180 mmHg is reached. The deflection rubber on the distal deflection part of the endoscope expands slightly. Move the deflection part in all four directions as far as possible and check the manometer to see whether the pressure remains constant or falls.

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4. Wait five minutes. If the pressure does not fall by more than 10 mmHg, submerge the endoscope in the cleaning solution. Place instrument with closed leakage tester in the cleaning solution and repeat the process of utilizing the deflection wheels to deflect the tip of the scope fully in all four directions. If this test is also negative and no air bubbles rise in the liquid, the leakage tester can be removed. To do this, take the supply connector out of the liquid and press down on the red button of the leakage tester to vent the endoscope.

NOTE: Never release the pressure in the liquid. To this end, remove the plug from the liquid and actuate the red pressure relief lever, before using the leakage tester again reset the red lever to the initial position. The test pressure must fall to 0. Remove the leakage tester.

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If a pressure loss occurs, this means that the video endoscope leaks and may no longer be used or immersed or prepared by machine.

In this event, fill out the accompanying repairs form (see Appendix), roughly clean and remove any organic soiling/residue from the endoscope, mark it as not disinfected and return it to the Service Department for repair or, alternatively, contact KARL STORZ. Line the storage/carrying case with the enclosed transport protection sheet 13990 SFN.

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WARNING: Risk of infection due to shipment. Contaminated instruments can cause aerogenous or contact infections during transportation or at the recipient's premises.

Shipment without the plastic film packaging must be avoided.

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
Kanalsystem eines KS Standard Video-
Endoskopes

Channel system of a KS standard video
endoscope

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
Channel system of a KS therapy video
endoscope with additional irrigation
channel (water jet channel)

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
Channel system of a KS therapy video endoscope with additional irrigation channel (water jet channel) and 2 working channels

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Cleaning/disinfection of the video endoscope: manual cleaning

Please adhere to the general standards, e.g., ESGENA guidelines or DEGEA/RKI recommendations for cleaning, disinfection and sterilization, as well as to the instructions of the chemical manufacturer.

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1. Remove the rubber lip valve/biopsy valve, suction valve and air/water valve and place in an active cleaning solution and clean using a soft disinfected brush (13991 VB).

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2. Wipe the housing, insertion shaft, and supply tube carefully with a lint-free cloth soaked in an enzymatic cleaner (use a suitable cleaning agent or disinfectant in accordance with the DGHM list and the list of approved chemicals, see annex). Avoid shocks and pulling movements so as not to damage the insertion shaft. Rinse the endoscope carefully with water of at least drinking quality (microbiologically sterile water).

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3. Immerse the endoscope in the active cleaning solution and flush the channels using a syringe and the supplied irrigation tubes (13991 SS, 13991 YS). The channels must be filled without any bubbles, to guarantee a complete moistening of the surfaces.

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

CAUTION: Only immerse instrument parts once they have successfully passed a leakage test.

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4. When submerged, use the brush supplied (Art. No. 27651 H/G) (see scope of supply page 10) to brush below the surface of the liquid in one direction starting from the supply plug (control body --> biopsy channel --> distal end).

! CAUTION: When brushing from the suction valve to the biopsy channel the brush comes against the biopsy channel (see channel system picture) and cannot be pushed any further into the biopsy channel!

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5. Now guide the brush through the biopsy channel to the distal end of the shaft.

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Before retracting the brush head that emerges at the distal end, clean it using another brush. Do not brush backward and forward, rather brush in one direction only.

! **CAUTION:** Cleaning with brushes must be carried out with the instrument immersed.

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
6 Thoroughly clean out the entire suction channel (sheath part and supply part) and the biopsy channel a number of times with a brush. Observe the exposure time for the active cleaning solution.

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7. Connect the Y rinsing tube (13991 YS and 13991 SS) and using a syringe, flush fresh cleaning fluid through the channels using a syringe. Afterwards, rinse with microbiologically clean/sterile water.

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8. Suitable brushes

Art. no.	Designation	Used for
27651 F	Cleaning brush dia. 0.5 mm	Working channel, Albarran channel
27651 G	Cleaning brush dia. 3 mm	Working channel
27651 H	Cleaning brush dia. 5 mm	
13991 VB	Valve brush dia. 12 mm	Valve seat

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9. Carefully clean all particles of debris from the distal end. It may be necessary to clean the air/water valve with a brush (13991 VB) with soft bristles, in order to remove hardened contaminants. Never remove debris with a needle or other sharp objects, because these can cause damage.

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10. Remove the endoscope from the cleaning solution and rinse with microbiologically pure/sterile water.

11. Carefully clean all valves with the disinfected cleaning brush (13991 VB) and then rinse. The valves should be cleaned in an ultrasound bath.

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Once the cleaning work has been completed, clean the cleaning instruments, for example the brushes, of residue and disinfect/sterilize them. Never reuse damaged or kinked brushes.

CAUTION: Do not autoclave the endoscope or its accessories (with the exception of accessories that are specifically designated as autoclavable, e.g., forceps, cleaning brushes, albarran module, etc.).

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
Manual disinfection

After careful rinsing, fully immerse the endoscope in the disinfecting solution.

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! CAUTION: When preparing and using the solutions follow the chemical manufacturers' instructions paying close attention to proper concentration and exposure times. Prolonged immersion can lead to damage. Bear in mind the microbiologic range of action of the chemicals used.

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Using the irrigation tubes and a syringe with LUER connector, rinse out the channels a number of times with disinfecting solution (see illustration).

Between rinses, allow the specified time (chemical manufacturers' instructions) for the disinfecting solution to take effect in the channel system.

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Then rinse the channels and surfaces of the video gastroscope with microbiologically pure/clean water.

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12. After the final step in the cleaning process has been completed, dry the entire endoscope. Thoroughly dry off all external surfaces of the endoscope, preferably with sterile compressed air at a max. of 0.5 bar/7.25 psi, thereby removing any residues of cleaning solution.

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Notes on cleaning special endoscopes

The picture opposite shows the disassembled Albarran module. The two upper parts should be cleaned with brushes (13991 VB) and then in an ultrasound bath with cleaning solution.

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Cleaning of the Albarran module

The following components can be cleaned in the ultrasound bath:

- Albarran module with cable
- Albarran cable guide
- Collet chuck
- Rotary knob


Carefully clean the distal end of the Albarran module with a soft brush (13991 VB).

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To enable better rinsing of the Albarran module's wire sheath (Albarran spring), the irrigation tube 13991 PS can be connected to it.

The lower part of the endoscope (in the above picture) is cleaned as described previously.

NOTE: The Albarran module of the video endoscopes can also be steam sterilized at 134 °C ± 3 °C.

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Starting from the proximal housing area, brush the working channel using the cleaning brush (Art. no. 27651 H/G) and cleaning agent.

Brush the disassembled Albarran channel with cleaning brush 27651 F supplied with the instrument.

Before retracting the brush head that emerges at the distal end, clean it using another brush or compressed air.

After the brush cleaning, the channels must be thoroughly rinsed with germfree, completely demineralized or microbiologically pure/sterile water.

Then disinfect the endoscope. Manual disinfection: see p.38

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Preparation of the video endoscope by machine

All flexible KARL STORZ endoscopes can be connected to approved washer/disinfectors for semi- or fully automatic cleaning and disinfection. Leakage testing, manual cleaning and (where necessary) rinsing must be carried out beforehand. The manufacturer of the cleaning machine must be consulted when selecting the preparation programs. Usually, special programs for flexible endoscopes are available which are suited for decontamination and compatible with the materials. For machine preparation, make sure that the flexible endoscopes are not heated above 65°C.

! **CAUTION:** Check the videoscope for leaks before every disinfection in semi- or fully automatic equipment.

! **CAUTION:** Ensure that only KARL STORZ approved agents are used for procedures involving peracetic acid.

! **CAUTION:** The STERRAD process is not suitable for sterilization of video endoscopes. This process destroys video endoscopes.

Lay the endoscope carefully in the tray or basket of the machine. Connect the irrigation tubes supplied with the cleaning adaptor set 13991 CK to the channels. Proceed with the endoscope cleaning program according to the recommendations of the machine manufacturer.

Cleaning adaptor set 13991 CK

The cleaning adaptor set 13991CK is used for validated preparation of KARL STORZ video endoscopes from series 13xxx NKS/PKS (acc. to ISO 15883 part 4). The cleaning adaptor set is not intended for any other purpose as this may damage it.

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- 1 Connection tube 13991 CL for leakage tester port ⑮ of the videoscope
- 2 Connector (13991CL) for vent port
- 3 LUER connector (13991CL) for washers/disinfectors for endoscopes

- 4 Cleaning adaptor 13991CR for air/water (A/W) channel ⑩ of the videoscope
- 5 Connector (13991CL) for A/W channel (air pin ②)
- 6 Connector (13991CL) for water channel (water bottle connection ⑦)
- 7 2x LUER connectors (13991CR) for washers/disinfectors (1x LUER for air channel, 1 LUER for water channel)

- 8 Cleaning adaptor 13991CS for suction channel (aspiration) ⑪ of the videoscope
- 10 Connector (13991CL) for suction channel at the supply connector ⑪
- LUER connector (13991CS) for washers/disinfectors

- 11 Connection adaptor 13991SS for biopsy/working channel ⑨ of the videoscope (with 2x double channel endoscopes)
- 12 Connector (13991CL) for biopsy/working channel ⑨
- 13 LUER connector (13991SS) for washers/disinfectors

- 14 Channel separator 13991CD for suction valve ④ and irrigation/insufflation valve (A/W) ⑤ of the videoscope
- 15 Valve pin (13991CL) for irrigation/insufflation valve (A/W) ⑤
- 16 Closure (13991CL) for suction valve ④

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WARNING: This product is delivered unsterilized and must therefore be cleaned, disinfected or, if necessary, sterilized prior to initial use and each subsequent reuse.

WARNING: Risk of infection by improper preparation: Preparation may only be carried out by experienced and trained staff.

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- Visual checks of the connection tubes (patency) and the O-ring (wear)
- check the connections are correctly and securely in place
- check for damage

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Assembly

Plug the connector (2) of the connection tube 13991 CL into the vent port (15) of the videoscope at the supply connector and secure by turning clockwise.

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Plug the cleaning adaptor 13991 CR (5) as far as it can go into the connector (22) of the air pin of the videoscope attachment.

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Then plug connector (6) as far as it will go into the 'Connection for water bottle' (17) of the videoscope.

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Plug the cleaning adaptor 13991CS with connector (9) as far as it will go into the suction channel (21) at the supply connector of the videoscope.

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If necessary, remove the suction and irrigation/insufflation valves (4, 5) from the videoscope.

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Insert the channel separator 13991 CD (14) with the valve pin (15) into the free opening of the irrigation/insufflation valve (4) as far as it will go so that the lock (16) lies on top of the free opening of the suction valve (4).

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Press the slide mechanism of the channel separator 13991CD (14) downward and slide it forward until the channel separator 13991 CD (14) audibly locks into place. Check that the channel separator is securely in place.

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NOTE: If necessary, remove the cap of the biopsy/working channel (6) of the video-scope.

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
Connect all open LUER connectors of the attached connection adaptors to the respective connections on the washer/disinfector.

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! **CAUTION:** When the program has finished, the conditions in the preparations machine present the videoscopes with a risk of corrosion.

! **CAUTION:** Immediately after the end of the program cycle, the video endoscopes are hot as they are heated up during chemical thermal preparation. Consequently the instrument's material properties are different and, as such, the prepared endoscope should be allowed to cool down for at least 20 minutes at room temperature before being used again.

! **CAUTION:** Ideally, remove the video endoscopes from the machine immediately after the program has finished. Then dry off the optical surfaces, if necessary, with a cotton swab soaked with 70% alcohol.

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Store the endoscopes freely suspended, dust-free and dry. Store valves separately.

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Trocknung

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Recommended procedure for the preparation of video endoscopes

You also will find an illustrated page of the preparation procedures in the middle of this manual. There is also a poster of the preparation process available from your local sales representative or it can be requested per e-mail directly from KARL STORZ (info@karlstorz.de)

Preliminary cleaning directly after the examination (see page 27)

- The outer surfaces should be wiped down from top to bottom with a lint-free disposable cloth that is soaked with active cleaning solution.
- Suction in ample amounts of cleaning solution through the distal end.
- Apply the air-water supply for a few seconds every cycle
- Connect the cleaning adapter
- Remove the flexible endoscope from the light source

CAUTION: Close the video plug with the sealing cap.

Cleaning in the preparation place

- Carry out leakage test
- Immerse the flexible endoscope
- Clean the channels with special cleaning brushes
- Clean the complete air/water system thoroughly (in the case of duodenoscopes, the Albarran unit must also be disassembled, see page 21)
- Rinse the channels

Disinfection

- Immerse the endoscope and rinse channels

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Preparation in the washer/disinfector

- The machine must fulfill all the cleaning, disinfection and material compatibility requirements (validated process).
- Make sure that the environmental conditions with regard to temperature and pressure are not exceeded: +60°C (tolerance +/- 5%) and 1 bar for the channels, 0.3 bar leakage test.
- Recommendation: Before placing the endoscope in the machine, clean all the channels (working channels, water channels, Albarran channels, air/water channels) whether or not they have been used.
- We recommend sticking to one type of cleaning agent.

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Cleaning of accessories

1. Thorough cleaning is required before disinfection or sterilization. Carefully clean the biopsy forceps and brushes under the surface of a cleaning solution using a soft, disinfected brush. Take special care to remove residues such as blood, excrement and mucous from the biopsy forceps and needles – an ultrasound bath is recommended.
2. Rinse in microbiologically pure/sterile water and dry (using sterile compressed air).
3. Disinfect accessories (work procedures such as cleaning).
4. Always test instruments after cleaning/disinfection for proper functioning and any damage.
5. Before sterilization, carefully oil moveable parts. Use only suitable, steam-resistant, physiologically safe, residue-free oil (e.g., 27656 B). After oiling, open and close the forceps several times. For details on preparation, please read the instruction leaflet of the respective accessory.

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Preparation of water bottle 13992 BS

The water bottle may be cleaned/disinfected manually or mechanically.

1. Empty the bottle
2. Carefully clean the container and lid with an active cleaning solution. Immerse the container and lid in the solution.
3. Attach the water bottle lid connector to the irrigation adapter 13991SW and rinse out with a cleaning liquid. Immerse in disinfectant solution (observe the chemical manufacturer's instructions on exposure time and concentration) making sure that all air bubbles have been removed. On removal rinse off thoroughly using microbiologically pure water.
- By machine: Place and secure in the washer and disinfectant. Choose between either chemical disinfection at 65°C or thermal disinfection at 93°C.
4. Sterilize. Steam sterilization is recommended.

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Check preparation procedure

Before assembly it should be checked that all parts and surfaces are undamaged, clean and dry. Replace defective parts or send them in for repair. Before sending them in for repair, all contaminated parts must be cleaned/disinfected and sterilized or clearly labeled as unsterilized.

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Check plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their original condition, they must be replaced. This especially applies to seals and rubber caps.

Test the instruments for proper operation (deflection of the distal end, brake lever). Check the lens for a clear view, the instrument for appropriate deflection and all channels for patency (see page 11).

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CAUTION: All deviations from the recommended sterilization parameters must be validated by the user.

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Gas sterilization

Sterilization methods using ethylene oxide is mild and therefore the recommended method of sterilization.

NOTE: Observe national laws and regulations for the gas sterilization method.

CAUTION: If the instrument is to be sterilized, the entire system must be ventilated by fitting the EO cap on the pressure compensation part to ensure pressure compensation.

WARNING: Before commencing gas sterilization the instrument must be completely dry or else sterile conditions will not be achieved.

- Place the instrument securely on a sterilization tray.

	20:80 EO/ CO	10:90 EO/ HCFC
Temperature:	55°C	55°C
Relative humidity:	50%	50%
Vacuum:	533 mmHg effective	533 mmHg effective
Outlet-pressure:	69 kPa	97 kPa
	0.70 Kg/cm, 10 psi	0.98 Kg/cm, 14 psi
EO-Concentration:	450 mg/l	600 mg/l
Conditioning:	1 hour	1 hour
Gas exposure time:	5 hours	5 hours
Aeration:	12 h at 55°C	12 h at 55°C

An aeration period of 72 hours is required after EO sterilization. An aeration chamber (desorption chamber) can be used to reduce the aeration period to 12 hours. During sterilization the temperature must not exceed 55°C.

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For gas sterilization with ethylene oxide, the following airing times must be observed for instruments due to gas absorption of the materials:

	at 20°C	at 42°C	at 55°C
Telescopes	24 hrs	18 hrs	12 hrs
Rubber items	7 days	4 days	3 days
Video endoscopes			
Plastics	5 days	2 days	30 hrs

(Minimum time in simple transparent packaging and an application time on the patient of more than 30 minutes).

Maximum acceptable levels (per 21 CFR 812.100; US Federal Register, Vol. 43, No. 122, June 23 1978) of residues following EO sterilization are as follows:

Ethylene oxide:	25 ppm*/250 ppm**
Ethylene chlorohydrine:	25 ppm*/250 ppm**
Ethylene glycol:	250 ppm*/5000 ppm**

** for instruments in contact with skin and mucous

membrane,

* for instruments in contact with blood.

WARNING: The efficiency of the sterilization should be checked to see that it complies with the German Medical Devices Act (MPG).

CAUTION: Attention must be paid to the notes on cleaning, sterilization and storage in the 'Cleaning, sterilization and care of KARL STORZ instruments' manual. The cleaning, disinfection and sterilization procedures are explained therein in detail.

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Storage

CAUTION: Do not expose video endoscopes to direct sunlight or store in an environment with excessive heat.

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The videoscopes should be stored hanging vertically, with the distal end and supply connector facing down, in well ventilated cupboards designed for this purpose.
Store dismantled valves separately.

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Maintenance

Preventive maintenance is not essential. Regular maintenance can contribute to identifying potential problems before they become serious, enhancing the instrument's reliability and extending its useful operating life.

Maintenance services can be obtained from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year.

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Important information

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and, if necessary, disinfected/sterilized before they are sent for repair. We reserve the right to return contaminated instruments/equipment to the sender.

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original parts.

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Limitation of liability

As a supplier of this instrument, KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions: All assembly, operations, system expansions, readjustments, modifications, or repairs to same have been performed by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

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Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see 'Subsidiaries' section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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

- ☐ cracked
- ☐ worn out

Distal tip

- ☐ Lens scraped
- ☐ Insulation disturbed
- ☐ Image transmission error

- ☐ Deflection
- ☐ leaking
- ☐ blocked / damaged

REPAIR DISPATCH NOTE

STORZ
KARL STORZ – ENDOSCOPE

From customer enclosed parts:

Customer	Article / Part	KARL STORZ
☐	Valve / air-water, suction valve	☐
☐	Cleaning brush	☐
☐	Protection wrap	☐
☐	Biopsy forceps	☐
☐	Cleaning tubes	☐
☐	Instrument oil	☐
☐	ETO cap	☐
☐	Bite protector	☐
☐	Leakage tester	☐
☐	Seal cap for working channel	☐
☐	Video connection cable	☐

OS KARL STORZ

Signature

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postfach 230, 78503 Tuttlingen/Germany
Phone: +49 (0)7461 708-0, Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com

KARL STORZ Endoskope Berlin GmbH
Scharnhorststr. 3
10115 Berlin/Germany
Phone: +49 (0)30 30 69090, Fax: +49 (0)30 30 19452

KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzcuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 México D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6° Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
12739 Skärholmen, Sweden
14108 Kungens Kurva, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania
Phone: +370 68567000
E-Mail: LT-LV-info@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytøften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopia France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Adria Eos d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Verovškova c. 60A
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Jinjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenheng District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

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This instruction manual is intended to serve as an aid in the proper handling, cleaning, disinfection and sterilization of the Video Urethro-Cystoscope PDD. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for ready reference.

WARNING: KARL STORZ instruments are provided unsterilized, and must be cleaned, disinfected and/or sterilized prior to the initial use and before each subsequent use.

KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology are possible for this reason. Therefore, no claims may be deduced from the information, figures and descriptions in this manual.

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
- ① Function button (video button)
- ② Suction valve 11900 AV
- ③ Valve holder 11900 VA
- ④ Connecting tube 11900 AS
- ⑤ LUER-Lock connectors (forceps unit) 11272 AZ
- ⑥ Supply connector
- ⑦ Connection for leakage tester or pressure compensation cap
- ⑧ Cap
- ⑨ Video connection cable **202130 70**
- ⑩ Deflection lever
- ⑪ Distal tip
- ⑫ Sheath
- ⑬ Handle
- ⑭ O-ring
- ⑮ O-ring
- ⑯ O-ring
- ⑰ Pressure compensation cap 11025 E
- ⑱ Leakage tester 13242 XL
- ⑲ Grasping Forceps 27023 FE
- ⑳ Biopsy forceps 27023 ZE
- ㉑ Cleaning brush 27651 B
- ㉒ Transport protective foil 13990 SFN
- ㉓ Cleaning adaptor 11900 RA

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10.6	Manual cleaning/disinfection of the Video Urethro-Cystoscope PDD	35
10.7	Machine cleaning/disinfection of the Video Urethro-Cystoscope PDD	40
10.8	Inspection after cleaning/ disinfection	42
10.9	Sterilization	43
10.9.1	Gas sterilization	44
10.9.2	Formaldehyde/Hydrogen peroxide sterilization	45
10.9.3	Sterilization with the STERIS® SYSTEM 1	47
10.10	Overview – Preparation	49
10.11	Assembly of the suction unit	50
10.12	Disinfection and sterilization of the biopsy forceps and the cleaning brush	51
11	Storage	53
11.1	Storage/transport conditions:	53
13	Servicing and repair	54
14	Maintenance	54
15	Limitation of liability	55
16	Warranty	55
17	Directive compliance	55
18	Important information	56
19	Disposal	56

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Please read these safety instructions carefully before using the device. Before using the instrument on the patient it is imperative that you be acquainted with how the instrument operates and is controlled.

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1 Intended Use

The flexible videoscopes/fiberscopes (flex. endoscopes) with special filters (PDD, AF, NIR/ICG) are used for endoscopic imaging during diagnosis and therapy in urology, ENT and bronchoscopy. Flexible endoscopes are designed for short-term use in endoscopic interventions performed through natural orifices.

1.1 Indications:

Use is indicated if fluorescence-supported visualization is helpful in the opinion of the attending physician.

The fluorescence systems with corresponding accessories are suitable for all diagnostic and therapeutic interventions in which the fluorescence imaging aids successful performance of the intervention.

These include, for example:

- The detection of urinary bladder cancer by a cystoscopy using Hexvix® (PDD)
- The detection of cancer in the upper respiratory tract via autofluorescence (AF)
- The visualization of certain structures, organs and tissue in abdominal surgery using indocyanine green (NIR/ICG)

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1.2 Contraindications

Use is contraindicated if, in the opinion of the attending physician, the system is not compatible with successful completion of the planned intervention due to its technical design.

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When using the fluorescence imaging, the medical contraindications apply as listed for the use of suitable marker substances (e.g., porphyria when using ALA derivatives), and use is contraindicated if endoscopic examination techniques contradict the current state of academic knowledge.

2 User qualification

The KARL STORZ Video-Urethro-Cystoscopes PDD may only be used by persons with an appropriate medical qualification and who are acquainted with the cystoscopy technique. The information given in these instructions serves to instruct in the correct handling and cleaning of the Video-Urethro-Cystoscope PDD.

It is not intended as an introduction to the Video-Urethro-Cystoscopy technique.

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3 Safety precautions when using the Video Urethro-Cystoscope PDD

The Video Urethro-Cystoscope PDD must be used according to the recognized medical rules and procedures of endoscopy.

NOTE: Patient leakage currents may accumulate. The patient leakage current limit of the CF application part may be exceeded when CF application parts (application on heart) are used in combination with other BF application parts. In this case, the reduction of the patient leakage current should be effected using only CF application parts as far as this is possible.

NOTE: The endoscope does not offer any electric insulation against high frequency (HF) currents. Protection against high frequency currents should be ensured through the use of correct HF instruments. Please read the instructions for use in the instruction leaflet supplied with the HF application part for more information.

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4 Warnings and cautions

Please read this manual carefully and follow the instructions exactly. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To further emphasize the signal words, they are preceded by a pictogram.

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WARNING: Read the following safety instructions carefully before using the Video Urethro-Cystoscope PDD to avoid putting your patients, personnel or yourself at risk.

WARNING: Inspect the Video Urethro-Cystoscope PDD prior to each surgical procedure to ensure that it functions correctly and that it has been correctly cleaned, disinfected and/or sterilized.

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WARNING: Combinations of medical devices are only then assured to be safe if

- they are identified as such in the respective instruction manual or
- the intended purpose and the interface specifications of the devices used in combination permit this.

Pay careful attention to the instructions and interface specifications of medical products used in combination.

WARNING: The video endoscope is delivered unsterilized and must be cleaned and disinfected/sterilized before the initial use as well as before any further application.

WARNING: The **pressure compensation cap** must be **removed** from the vent port for:

- endoscopic interventions
- manual/machine preparation

WARNING: Only use accessories (guide probes, cleaning brushes etc.) that are specifically suited for the intended purpose and can be employed in line with generally accepted technical standards and safety guidelines.

WARNING: Do not allow any live components or instruments to come into contact with the Video Urethro-Cystoscope PDD.

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WARNING: The Video Urethro-Cystoscope PDD must not be used during the discharge of a cardiac defibrillator. Remove the Video Urethro-Cystoscope PDD from the operation area before discharge takes place.

WARNING: Danger of explosion if flammable gases are used in the immediate proximity of the Video Urethro-Cystoscope PDD. The Video Urethro-Cystoscope PDD must not be used in the presence of flammable or combustible media.

WARNING: High energy light radiated through endoscopes may give rise to high temperatures in front of the light outlet, at the end faces of the light guide and at the tip of the endoscope. Prolonged direct contact of light-emitting surfaces on tissue can lead to burns, this must therefore be avoided at all costs.

WARNING: Even after the proximal connector is removed from the light source, surfaces will remain hot for a period of time (Risk of user being burnt).

WARNING: The **pressure compensation cap** must be **in place** on the vent port for:

- gas sterilization / hydrogen peroxide sterilization
- venting
- shipping

WARNING: The instructions 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' must be followed. There the procedures for cleaning, disinfection and sterilization are explained in detail.

WARNING: The Video Urethro-Cystoscope PDD may only be repaired by personnel authorized by KARL STORZ. Unauthorized repairs may impair functioning and safety, and are therefore strictly forbidden. KARL STORZ gives no warranty for Video Urethro-Cystoscopes PDD that are repaired by unauthorized personnel.

WARNING: Damage to the Video Urethro-Cystoscope PDD and injury to the patient or user resulting from incorrect operation are not covered by the manufacturer's warranty.

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WARNING: Always unplug the Video Urethro-Cystoscope PDD before carrying out any maintenance and cleaning work on it.

WARNING: Risk of infection during shipment! Contaminated Video Urethro-Cystoscopes PDD pose a risk of infection. For this reason, always wipe-down disinfect the Video Urethro-Cystoscope PDD and sterilize with the pressure compensation cap in place before shipping!

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CAUTION: Do not expose the Video Urethro-Cystoscope PDD to direct sunlight or excessive heat, to avoid damage to the device.

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CAUTION: The maximum permissible temperature for decontamination, sterilization, transport and storage is 65°C.

CAUTION: Do not hit the distal tip of the Video Urethro-Cystoscope PDD against hard objects, and never use a sharp object to remove dirt.

CAUTION: Never lay heavy objects on the Video Urethro-Cystoscope PDD.

CAUTION: Constant changing between various preparation methods is especially hard on the materials and must therefore be avoided.

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8 Electromagnetic Compatibility

CE marked equipment has been tested and is found to comply with the EMC limits of EN/IEC 60601-1-2: 2001 CISPR 11 Class B.

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9 Handling

NOTE: Damage to the Video Urethro-Cystoscope PDD resulting from incorrect operation is not covered by the manufacturer's warranty.

WARNING: The Video Urethro-Cystoscope is delivered unsterilized, so it has to be disinfected/sterilized prior to initial use and before each subsequent reuse.

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9.1 Operating

The following instructions outline the correct operation of the Video Urethro-Cystoscope PDD and the corresponding accessories. This is not an introduction to the techniques of cystoscopy. Instructions on endoscopic technique can be found in the relevant medical literature. Suitable accessories for your Video Urethro-Cystoscope PDD, such as biopsy forceps, can be found in the KARL STORZ Urology catalog or extracts from this catalog.

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9.2 Holding the Video Urethro-Cystoscope PDD correctly

Grip the handle of the Video Urethro-Cystoscope PDD in your left or right hand, as shown in the picture. With this kind of control grip, the Video Urethro-Cystoscope PDD can be operated with one hand. You can control the deflection lever with your thumb, while your index finger and/or ring finger operates the function button.

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9.3 Handling the sheath

CAUTION: Do not carry the Video Urethro-Cystoscope PDD by the sheath. Do not pull, clamp or twist the sheath of the Video Urethro-Cystoscope PDD.

CAUTION: The sheaths of the Video Urethro-Cystoscope PDD must never be kinked or coiled up too tightly. It is best to suspend the Video Urethro-Cystoscope PDD for storage or to keep it in a large container.

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CAUTION: Only use the oil 27656 B which has been tested by the manufacturer in order to avoid damaging the Video Urethro-Cystoscope PDD.

WARNING: Do not apply lubricant to the lens, to avoid clouding the image and thereby prevent risks for the patient!

Use a suitable water-soluble lubricant for the distal section of the insertion tube. Observe the lumen the whole time during use. If you lose your orientation, pull the sheath back slightly into the neutral position until vision is restored.

If, during use, the image fails, is disturbed, or becomes blurred, or if vision in general is impeded, put the deflection mechanism into the neutral position and slowly pull the sheath of the Video Urethro-Cystoscope PDD out of the patient's body.

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9.4 Handling the deflection mechanism

NOTE: The neutral position of the distal tip is straight, without any deflection whatsoever.

WARNING: Suitable accessories, such as biopsy forceps, must only be introduced into the instrument channel of the Video Urethro-Cystoscope PDD with the distal tip in the neutral, non-deflected position. Once the instrument has been introduced, carefully actuate the lever gently as far as the limit stop. The deflection lever must never be operated with force or abruptly, to avoid injuring the patient or damaging the device irreparably!

WARNING: Before removing the Video Urethro-Cystoscope PDD from the patient, the distal tip of the Video Urethro-Cystoscope PDD must be put into the neutral position, as otherwise injuries to the patient or damage to the instrument may occur.

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WARNING: If the distal tip can no longer be properly controlled (sticking etc.) via the deflection lever, or if the deflection mechanism develops any other fault, the Video Urethro-Cystoscope PDD must not be used.

WARNING: If a functional defect occurs during use on the patient, discontinue use immediately. Put the distal tip of the Video Urethro-Cystoscope PDD into the neutral, non-deflected position and remove the Video Urethro-Cystoscope PDD slowly and carefully from the patient.

WARNING: The lens at the distal end of the Video Urethro-Cystoscope PDD may become detached, e.g. due to mechanical forces. For this reason, it is essential to check that the lens is undamaged and present before and after every cleaning.

WARNING: The resistance that can be felt when deflecting the tip using the deflection lever increases as the angle of deflection increases. Actuate the deflection gently and carefully, without any sudden or jerky movements, as far as the limit stop.

WARNING: Do not use force to turn the deflection lever, to avoid either injuring the patient or damaging the instrument.

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
The distal tip can be deflected 180° upward or 100° downward by actuating the deflection lever slowly and carefully.

Before use, the user must always have full control over the deflection and must develop a sure feel for the deflection mechanism and in particular for the neutral position.

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9.5 Handling the accessories (forceps, cleaning brushes, HF electrodes etc.)

CAUTION: Only use flexible accessories (forceps, cleaning brushes etc.) from KARL STORZ. Accessories from other manufacturers may damage the Video Urethro-Cystoscope PDD.

CAUTION: Only introduce the instruments from the proximal port of the Video Urethro-Cystoscope PDD.

WARNING: When HF electrodes or laser probes are used, it must be ensured that power is not applied until the active parts have left the instrument channel and are in the user's field of vision, and the site of the operation is clearly recognizable. Activation in the instrument channel may lead to serious injury to the patient and damage to the Video Urethro-Cystoscope PDD.

WARNING: When HF electrodes or lasers are used, it must be ensured that power is only applied when there are no flammable gases in the operating field or in the proximity of the Video Urethro-Cystoscope PDD.

CAUTION: When laser probes are used, it must be ensured that the laser beam cannot hit the distal tip of the Video Urethro-Cystoscope PDD.

CAUTION: Do not force accessories, such as biopsy forceps, through the working channel, to avoid damaging the device. It should be possible to move the accessory easily, and without any hindrance, through the entire channel.

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CAUTION: When the distal tip is deflected, do not introduce any instrument into the Video Urethro-Cystoscope PDD. Return the distal tip to the neutral position before introducing the instrument and then move the deflection lever gently and carefully to the limit stop.

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9.6 Compatible camera control units/light sources

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WARNING: The KARL STORZ Video Urethro-Cystoscope PDD should not be used in combination with camera control units and light sources from other manufacturers. The use of non-compatible units can lead to patient injury or unit damage.

It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces, e.g., voltage-carrying elements of other units, at any time.

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CAUTION: Inserting sharp-tipped accessories into the instrument channel may damage the Video Urethro-Cystoscope PDD.

CAUTION: Foreign bodies or parts to be extracted must not be removed through the instrument channel using instruments. Instead, the entire Video Urethro-Cystoscope PDD must be removed from the patient and the foreign body extracted.

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9.7 Connection to the light sources D-LIGHT C and D-LIGHT C/AF

WARNING: The Video Urethro-Cystoscope PDD can be used with the light sources D-LIGHT C 201336 01-1 and D-LIGHT C/AF 201336 01-133 specified by KARL STORZ for this purpose. The instruction manual of the light sources must be followed.

WARNING: Danger of burns and glare! Never look directly into the light emitted from the light source or the Video Urethro-Cystoscope PDD. The high-intensity light emitted through the Video Urethro-Cystoscope PDD may cause glare or high temperatures, and thereby burns, at the distal light outlet.

WARNING: Danger of burns and fire hazard! Never place the end of the fiberoptic cable or the distal tip of the Video Urethro-Cystoscope PDD on or underneath a sterile cover cloth. The intensity of the light can cause burns to the patient and/or the sterile cover cloth.

WARNING: Reduce the intensity of the light source so that optimal illumination of the operating field is still achieved. The higher the intensity of the light source, the greater the heat emitted at the tip of the Video Urethro-Cystoscope PDD.

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NOTE: For the correct settings for the light sources, please see the operating instructions of the devices used.

1. Install the light source on a firm, level surface.
2. Plug the power cord into the specified power socket.
3. Plug the supply connector into the light source.

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9.8 Connection to the camera control unit

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CAUTION: The Video Urethro-Cystoscope PDD can be used with the camera control unit TELECAM SL II FI 20 2130 20-1 specified by KARL STORZ for this purpose. The instruction manual for the camera control units must be followed.

CAUTION: Take care with the connectors on the Video Urethro-Cystoscope PDD! The connectors on the Video Urethro-Cystoscope PDD or video connection cable can be damaged if they are incorrectly plugged into the Video Urethro-Cystoscope PDD. Make sure that the pins are correctly aligned.

CAUTION: Before using the devices, ensure that all connectors (Video Urethro-Cystoscope PDD to camera control unit) are connected to one another. Connecting the connectors to a device in operation may cause malfunctions or damage to the Video Urethro-Cystoscope PDD.

NOTE: The Video-Urethro-Cystoscope PDD can be connected to the TELECAM SL II FI with the video connecting cable 202130 70.

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Install the camera control unit TELECAM SL II FI on a firm, level surface.

1. Plug the power cord into the specified power socket.
2. Unscrew and remove the cap for the video cable connector.
3. Plug the video connection cable into the round port in the supply connector.
4. Connect the video connection cable to the camera control unit.
5. Press the power switches on the devices.

NOTE: Always set the exposure setting of the TELECAM SL II FI camera control unit, in conjunction with the Video Urethro-Cystoscope PDD, to 'Brightness Low'.

NOTE: For the correct settings for the camera control unit, please see the operating instructions of the devices used.

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9.9 User button

The video-urethro-cystoscope PDD has three control buttons which can be assigned specific functions in the camera control unit's device menu. The assignment set is retained after switching off and on again. Instructions on programming the control keys can be found in the instruction manual of the camera control unit.

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9.10 Connection to the suction unit

CAUTION: The Video Urethro-Cystoscope PDD can be used with the ENDOMAT® SELECT suction pump recommended by KARL STORZ. The instruction manual for the suction roller pumps must be followed.

1. Install the suction roller pump on a firm, level surface.
2. Plug the power cord into the specified power socket.
3. Connect the Video Urethro-Cystoscope PDD to the suction devices with the suction tube.

To commence drainage, press the suction button on the Video Urethro-Cystoscope PDD.

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9.11 Preparation prior to use

NOTE: Before the device is used for the first time, it must be subjected to a complete preparation cycle (without preliminary cleaning) including leakage testing.

WARNING: Every time the Video Urethro-Cystoscope PDD is used, you must first check that it, and any accessories used in combination with it, are in perfect condition. Damaged Video Urethro-Cystoscopes PDD or damaged accessories must not be used.

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9.12 Inspecting the surfaces

WARNING: The surface of the flexible sheath must be inspected visually and by carefully feeling it, to detect any cracks, cuts, dents, foreign bodies or other forms of damage.

Every time the Video Urethro-Cystoscope PDD is used, you must first check that its surface, and above all the surface of the sheath and distal tip, is in perfect condition and that there are no sharp edges. Video Urethro-Cystoscopes PDD with sharp edges must not be used. Otherwise serious injuries may be caused to the patient or the user.

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CAUTION: Do not squeeze the distal tip of the Video Urethro-Cystoscope PDD too tightly.

Inspect the supply tube for damage (such as breaks, cracks, twisted or squashed sections...) and check the supply connector for loose connections.

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Examine the surface of the flexible sheath visually and by feeling it, to detect any cracks, cuts, dents, foreign bodies or other forms of damage.

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9.13 Checking the deflection

WARNING: Actuate the deflection lever gently and carefully as far as the limit stop. The deflection lever must never be operated with force or abruptly, to avoid injuring the patient or damaging the device! Move the deflection lever smoothly and slowly!

WARNING: Before use, always check that the deflection lever and the distal tip connected to it function correctly. Defective Video Urethro-Cystoscopes PDD must not be used.

The distal tip is deflected by means of the deflection lever. Turn the deflection lever slowly; the end piece should move downward. Turn the deflection lever in the opposite direction to produce an upward movement.

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9.14 Inspecting the light source

Prepare the light source as described on page 14 and connect to the Video Urethro-Cystoscope PDD. Emission of light at the distal end of the Video Urethro-Cystoscope PDD must be recognizable. The light intensity can be adjusted manually via the light source.

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9.15 Inspecting the camera control unit

Prepare the camera control unit as described on page 13 resp. 15 and connect to the Video Urethro-Cystoscope PDD. If an image is recognizable, carry out a white balance as explained in the instruction manual for the TELECAM SL II FI camera control unit. The function buttons 2 and 3 can be freely configured. When doing this, please observe the instructions for the camera control unit.

The function button can be individually programmed. Programming of the function button is described in the instructions for the TELECAM SL II FI camera control unit. The factory setting permits a still image by pressing the button once briefly. To return to the current transmission, press the button again briefly. Or to print out the still image (only possible if the system configuration includes a printer), press the button longer (5-7 seconds).

WARNING: If the image fails or there is interference during use in the patient, put the distal tip into the neutral position and remove the Video Urethro-Cystoscope PDD from the patient's body gently and carefully. If the Video Urethro-Cystoscope PDD is used in the patient after the image has failed or when there is interference, this may result in serious injury to the patient.

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9.16 Checking the suction function

When the Video Urethro-Cystoscope PDD has been installed correctly and the suction unit has been connected (see page 18), immerse the distal end in a container filled with water.

Press the suction button (red ring) in fully. The water must flow through the device and into the collecting container.

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9.17 Inspecting the instruments and the instrument channel

CAUTION: It must be possible to move the forceps easily, and without any hindrance, through the entire channel.

CAUTION: Accessories, such as forceps, which are bent or twisted may no longer function properly, are difficult to pass through the channel, and may cause serious damage to the channel. They must therefore be replaced!

Check the instrument channel by introducing the appropriate instrument (max. diameter of instruments: 1.7 mm).

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10 General information

WARNING: Video Urethro-Cystoscopes PDD must be rinsed out immediately after an intervention in order to prevent fluids from becoming dried on. Dried-on fluids may jeopardize decontamination and any sterilization required.

Furthermore, it prevents coagulum/blockage and consequent repairs and equipment failure. For the preliminary cleaning, suitable enzymatic cleaners should preferably be used, as they will not attack the surface.

WARNING: Before and after using the Video Urethro-Cystoscope PDD, always ensure that it has been cleaned, disinfected and/or sterilized correctly and also correctly stored.

WARNING: The listed cleaning, disinfection and sterilization methods alone do not guarantee successful sterilization. This can only be achieved if personnel follow recognized and validated preparation methods.

WARNING: The local occupational safety guidelines must be observed and applied when cleaning contaminated Video Urethro-Cystoscopes PDD and instruments. Aerosols should be removed by means of a suction device.

WARNING: Remember that sterilization is only successful if

- the instruments have been freed of all organic materials and cleaning residues. Successful sterilization can only be ensured on visually clean, disinfected surfaces.
- only recognized and validated methods are used. The reparation processes for Video Urethro-Cystoscopes PDD must be validated on-site by the operator.
- any deviation from the recommended procedures is validated by the user.
- the recommended sterilization parameters are maintained. These are only valid with sterilization equipment that is properly maintained and calibrated.

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CAUTION: The instructions 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' must be followed. There the procedures for cleaning, disinfection and sterilization are explained in detail. A list of KARL STORZ approved products for preparing videoscopes can be found on the Internet at 'www.karlstorz.com' and in the appendix to this instruction manual.

CAUTION: When preparing and using the solutions, the solution manufacturer's instructions concerning dilution and immersion time must be strictly adhered to. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiologic range of action of the chemicals used.

CAUTION: Never clean Video Urethro-Cystoscopes PDD with a metal brush.

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CAUTION: Video Urethro-Cystoscopes PDD must not be cleaned in an ultrasound bath.

CAUTION: Video Urethro-Cystoscopes PDD must not be steam sterilized (autoclaved) under any circumstances.

CAUTION: Video Urethro-Cystoscopes PDD must never be immersed in normal saline solution, as even brief contact can result in pitting and corrosion.

CAUTION: Video Urethro-Cystoscopes PDD may only be rinsed using microbiologically clean/sterile water.

CAUTION: Disinfection in an endoscopy rinsing machine (RDGE) is only possible using suitable methods.

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CAUTION: Special attention should be given to the following sterilization conditions:
Temperature: 5°C – 65°C
Humidity: 15% – 100%

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- gas sterilization / hydrogen peroxide sterilization
- shipping
- venting

Preparation must only be carried out by authorized, trained personnel. The guidelines of the Employer's Liability Insurance Association and equivalent organizations must be observed when cleaning and disinfecting contaminated medical devices.

Thorough cleaning of the Video Urethro-Cystoscope PDD and all accessories after every use ensures interruption-free working and optimal disinfection and sterilization results.

The generally applicable guidelines and recommendations for instrument preparation (e.g. RKI, ESGE/ESGENA) must be followed.

NOTE: Disinfection prior to manual cleaning is ineffective, and results in proteins and organic residues becoming fixed on/in the instrument, thus preventing correct disinfection and sterilization.

CAUTION: Before completely immersing the instrument in disinfectant solution, check that it is leakproof, as otherwise penetration of liquid may cause damage.

Note: It is recommended to store non-sterilized Video Urethro-Cystoscopes PDD in a dustproof place following cleaning and disinfection, preferably in a suspended position or in a tray specifically designed for sterilizing and storage (39403AS).

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10.1 Preparation for cleaning

Wipe the housing, sheath and supply tube thoroughly with a lint-free disposable cloth, preferably moistened with an enzymatic cleaner. Do not use alcohol for wipe-down disinfection of the sheath.

Do not use pushing and pulling movements, to avoid damaging the sheath.

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Immerse the distal end in a container filled with cleaning solution and actuate the suction to rinse out the working channel. Suck air and cleaning fluid alternately through the distal end. Repeat the procedure several times.

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Remove the supply connector of the Video Urethro-Cystoscopes PDD from the light source.

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Pull the video connection cable off the Video Urethro-Cystoscopes PDD and the camera control unit.

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Carry out a leakage test, regardless of whether manual or machine preparation takes place.

Leakage test successful: Continue

Leakage test ambiguous: Hold the leakage tester and immerse the Video Urethro-Cystoscope up to the anti-kink protection. Bubble check: The absence of bubble formation indicates that the Video Urethro-Cystoscope is leakproof.

Leakage test unsuccessful: Perform wipe-down disinfection, fill in the card and return the Video Urethro-Cystoscope to the manufacturer.

For every further step, it is essential to ensure that the Video Urethro-Cystoscope is completely vented and that the the manometer has been removed from the connection to the leakage tester.

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10.2 Leakage test

! WARNING: Before and after each use, the Video Urethro-Cystoscope must be tested for leakage. Leaking Video Urethro-Cystoscopes must not be used. They must be sent to KARL STORZ for repair or replacement.

! WARNING: Always test the Video Urethro-Cystoscope for leakage before starting cleaning or disinfection, and in particular before immersing the Video Urethro-Cystoscope. Leaking Video Urethro-Cystoscopes must not be immersed, but must be wipe-down disinfected and then sent to the manufacturer with an accompanying letter for repair or replacement of the Video Urethro-Cystoscope. Additional gas sterilization must be carried out in accordance with the Infection Protection Act.

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The connecting tube of the leakage tester and the tester connection on the Video Urethro-Cystoscope must be dry! Attach the tester connection cap firmly to the leakage tester connection and lock it by turning 25° clockwise.

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Pump up the leakage tester until the maximum test pressure of 160 mmHg is reached. The rubber sleeve on the distal deflection part of the Video Urethro-Cystoscope inflates slightly.

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Wait five minutes. If the test pressure does not drop more than 10 mmHg, the device is watertight and can be immersed. You can now begin the cleaning and disinfection procedure.

If the results of the leakage testing are ambiguous, we recommend immersing the videoscope with the working sheath in the cleaning solution up to the anti-kink protection. The absence of bubble formation indicates that the Video Urethro-Cystoscope is leakproof. Vent and remove the leakage tester, continue preparation procedure.

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If a pressure loss occurs, this means that the Video Urethro-Cystoscope leaks and must not be used or immersed. In this case, send the Video Urethro-Cystoscope for repair or contact KARL STORZ.

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WARNING: Risk of infection due to shipment. Contaminated instruments can cause serious or fatal infectious diseases during transportation or at the recipient's premises. Gas sterilize the Video Urethro-Cystoscope with the pressure compensation cap mounted before shipment (Please attach a sterilization indicator for documentation purposes). Enclose a written note about any special features (such as use on TB patients) and the disinfection measures taken, to protect personnel.

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10.3 Cleaning/disinfection of the suction unit

Dismantling the suction unit

! **CAUTION:** Two O-rings (Art. no. 5925300 and 5925600) on the suction valve and one O-ring (Art. no. 5928000) on the suction cylinder can be lost. Each time the unit is dismantled and cleaned, disinfected and sterilized, inspect the O-rings and, if necessary, install new ones.

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Unscrew and remove the suction valve by turning counterclockwise.

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Remove the LUER-Lock connectors (forceps unit) by pulling them off the Video Urethro-Cystoscope.

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
Pull the valve holder out of the recess. The suction unit has now been removed completely, and comprises:

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- ② 11900 AV Suction valve (with O-Ring 5925300 and 5925600)
- ③ 11900 VA Valve holder
- ④ 11900 ASE Connecting tube
- ⑤ 11272 AZ LUER-Lock connector (forceps unit) (with O-Ring 5928000)

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10.4 Cleaning the LUER-Lock stopcocks

Clean the inlet and outlet nozzles with a brush; use
water and a compressed air gun if necessary.

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To continue cleaning, dismantle the stopcocks.

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Carefully clean the stopcocks with a disinfected
brush.

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Rinse with water and blow dry with sterile
compressed air.

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Assembly

Lightly grease the running surfaces of the
stopcocks (special grease for stopcocks; Art. No.
27657).

NOTE: This thin layer of grease not only
serves to lubricate the metal surfaces
which are in contact with one another, but
also makes the stopcocks leakproof. Oil is
not suitable as a sealant.

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Insert the stopcocks and fix in place using the
spring cap. Tighten slightly.
Check the function of the stopcocks.
Dismantle the stopcocks again for sterilization.

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10.5 Manual cleaning/disinfection of the suction unit

The dismantled suction unit can be first cleaned in an ultrasound bath.

Clean the surface thoroughly with a soft disinfected brush, removing all soiling, tissue remnants and secretions, and rinse all the surfaces with cleaning solution.

Then clean out all the openings, channels and lumens of the dismantled suction unit thoroughly with a soft brush, removing all soiling, tissue remnants and secretions, and rinse out all the openings with cleaning solution.

Immerse the dismantled suction unit in a disinfectant solution. It must be ensured that all the lumens are filled with solution. Remove the dismantled suction unit from the disinfectant solution after the exposure time specified by the disinfectant manufacturer, and rinse thoroughly with microbiologically clean/sterile water for at least a minute.

Sterilization of the suction unit is described in the section 'Sterilization'.

When the cleaning procedure has been completed, remove residues from the cleaning equipment, such as reusable brushes, and disinfect. Do not reuse damaged brushes, as this could cause damage to the working channel.

Then rinse all the openings, channels and lumens of the individual parts of the suction unit thoroughly with microbiologically clean/sterile water until the water is visibly clear.

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10.6 Manual cleaning/disinfection of the Video Urethro-Cystoscope PDD

- !** **CAUTION:** The pressure compensation cap must be removed from the vent port for:
- endoscopic interventions
 - immersing in cleaning/disinfection solution

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- !** **CAUTION:** The pressure compensation cap must be in place on the vent port for:
- gas sterilization/hydrogen peroxide sterilization
 - venting
 - shipping

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Pull off the video connection cable. Put on the cap for the video cable connector and screw tight. Carry out leakage test. For this, see the section 'Leakage test' (p. 27).

- !** **CAUTION:** Never use a needle or other sharp object for removing dirt particles from the device, as this may cause serious damage.

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Carefully clean the distal end to remove all dirt particles.

Wipe the housing, sheath, and supply tube thoroughly with a lint-free disposable cloth, preferably moistened with an enzymatic cleaner. Do not use pushing and pulling movements, to avoid damaging the sheath. Subsequently, rinse the Video Urethro-Cystoscope PDD thoroughly with microbiologically clean/sterile water.

CAUTION: Only pass the cleaning brush 27651 AL through the instrument channel in one direction, and only pull it back once the brush head has emerged at the distal end. The brush must not be moved to and fro in the channel, to avoid damage.

Introduce the brush 27651 AL supplied into the instrument channel and slide it forward to free the instrument channel of soiling, tissue remnants and secretions.

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WARNING: The lens at the distal end of the Video Urethro-Cystoscope PDD may become detached, e.g. due to mechanical forces. For this reason it is essential to check that the lens is present before and after every cleaning.

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When the brush head emerges at the distal end, clean it with another brush or compressed air before pulling the brush back. Then pull the brush back through the instrument channel and repeat the brushing procedure several times.

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To rinse out the channel, fill a syringe with a cleaning solution.

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Carefully press the cleaning solution through the instrument channel, so that it flows out of the distal end of the channel.

NOTE: For easier rinsing of the channel, Cleaning Adaptor 11900 RA can also be connected to the working channel.

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WARNING: Observe the concentration required for effective cleaning, as specified by the chemical manufacturer. Then rinse with microbiologically clean/sterile water until clear water flows out. Repeat the procedure several times.

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After making sure all parts are completely covered, leave to stand. In doing so, observe the manufacturer's specifications concerning exposure time. Then rinse the disinfectant solution away using microbiologically clean/sterile water until the water runs clear. Repeat the procedure several times.

CAUTION: The pressure must not exceed 50 kPa (0.5 bar/7.25 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.

Dry the instrument channel and all the openings, channels and lumens of the Video Urethro-Cystoscope PDD and its surface thoroughly, preferably with sterile compressed air at a pressure not exceeding 0.5 bar/7.25 psi.

CAUTION: Test the components of the device for leakage before immersing them.

NOTE: We recommend using a plastic basin for manual preparation, to avoid scratching. A plastic basin also prevents electrolytic corrosion, which can occur if different metals are present in a solution.

The leakage tester must be connected and the Video Urethro-Cystoscope PDD pumped up to 160 mmHg. Immerse the Video Urethro-Cystoscope PDD in the prepared solution up to the anti-kink protection.

Observe for 5 minutes to see whether bubbles rise regularly, which would indicate a leak. If this is the case, remove the device and contact KARL STORZ.

Remove the Video Urethro-Cystoscope PDD from the solution after the exposure time specified by the chemical manufacturer, and rinse with microbiologically clean/sterile water for at least a minute.

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10.7 Machine cleaning/disinfection of the Video Urethro-Cystoscope PDD

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- endoscopic interventions
- complete immersion in solutions

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All flexible KARL STORZ Video Urethro-Cystoscopes PDD can be connected to devices for semi- or fully automatic cleaning and disinfection.

Before doing so, preliminary cleaning and a manual leakage test must have been carried out as described above. The choice of programs must be agreed with the machine manufacturer. In most cases, special programs for flexible endoscopes are available, which are suitable for cleaning without damaging the material. Remember that for machine cleaning of flexible Video Urethro-Cystoscopes PDD, the temperature must not rise above 65°C.

The suction channel of the Video Urethro-Cystoscope PDD must – following disassembly of the suction unit – be connected to the washer with the aid of Cleaning Adaptor 11900 RA.

The individual parts of the disassembled suction unit can, preferably when connected to the cleaning machine, be chemically disinfected by machine.

WARNING: Before you place Video Urethro-Cystoscopes PDD in the washer, clean all parts of the Video Urethro-Cystoscope PDD following the instructions given above for preliminary cleaning, irrespective of whether they have been used.

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We recommend a manual leakage test irrespective of whether the machine is equipped with an automatic leak tester.

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The machine must fulfill all the demands made of it in terms of process validation, cleaning performance, disinfection and material compatibility. Ensure that the required temperature and pressure (+60°C, with +/- 5% tolerance, and 1 bar for channels, 0.3 bar leakage test) are not exceeded.

It is not recommended to change between manual and automatic preparation.

! **CAUTION:** Before cleaning in automatic or semi-automatic devices, always check the Video Urethro-Cystoscope PDD for leakage. Carefully place the Video Urethro-Cystoscope PDD in the tray of the machine. Connect the rinsing tubes to the channels and, if available, connect the leakage tester of the machine to the Video Urethro-Cystoscope PDD via the adaptor. Use a cleaning program for flexible endoscopes, following the machine manufacturer's instructions.

! **CAUTION:** After the program has run, the Video Urethro-Cystoscopes PDD are at a high risk of corrosion due to the ambient conditions in the preparation machine. Once the program has ended, remove the Video Urethro-Cystoscopes PDD from the machine immediately.

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After machine preparation, clean/dry the distal lens and eyepiece, and the light port on the supply connector with a lens cleaner or a cotton tip applicator soaked in 70% isopropyl alcohol.

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10.8 Inspection after cleaning/ disinfection

Before assembling, visually inspect all the parts and surfaces of the Video Urethro-Cystoscope PDD for damage and soiling and to ensure that they are dry.

It is essential to ensure that there is no moisture residue under the cap or on the connector unit.

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WARNING: Before reusing the device, replace any defective parts or send them for repair.

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If residues are still present after cleaning, prepare the parts again, give them special treatment and separate out any parts that require repair.

Inspect plastic parts for discoloration, porosity and flexibility. If the parts differ significantly from when they were new, they must be replaced. This applies in particular to seals and rubber caps. Subject the instrument to a functional test. Check the telescope for clear view, instruments for steerability and all channels for patency. Lubricate the O-rings of the suction unit with instrument oil 27656 B.

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10.9 Sterilization

WARNING: Sterilization of the Video Urethro-Cystoscope PDD may be necessary for hygiene-related reasons and to prevent infections. Due to the high risk classification, at least a high-level disinfection to eliminate all micro-organisms that are capable of reproduction (with the exception of some spores) is required. This is usually achieved by means of disinfection with products based on Aldehyde or peracetic acid. However, these products should not be used for cleaning. Manual cleaning and disinfection are performed in a separate step.

WARNING: All deviations from the recommended sterilization parameters must be validated by the user.

WARNING: Before sterilization, all parts of the instrument must be thoroughly cleaned, rinsed and dried. All the parts of the KARL STORZ Video Urethro-Cystoscopes PDD can be disinfected chemically with a disinfectant solution. Suitable products are those specially recommended and approved by KARL STORZ for videoscopes.

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- endoscopic interventions
- immersing in liquid/chemical disinfection.

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- chemical sterilization (i.e. hydrogen peroxide / ETO- / FO)
- venting
- shipping

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10.9.1 Gas sterilization

Sterilization methods using formaldehyde or ethylene oxide are gentle and therefore recommended methods of sterilization.

-  **NOTE:** Observe national laws and regulations for the gas sterilization method.
-  **CAUTION:** If the instrument is to be sterilized, the entire system must be opened by fitting the EO cap on the pressure compensation part to ensure pressure compensation.
-  **WARNING:** Before commencing gas sterilization, the instrument must be cleaned, disinfected and completely dry or else successful sterilization will not be achieved.
-  **NOTE:** The suction unit should preferably be disassembled before being gas sterilized.
-  **CAUTION:** Attention must be paid to the information on preparation, sterilization and storage in the 'Cleaning, sterilization and care of KARL STORZ instruments' instructions. The methods for cleaning, disinfection and sterilization are explained there in detail.

- Place the instrument securely on a sterilization tray.
- Ethylene oxide gas sterilization has been validated with the following parameters:
- | | |
|---------------------------------|--|
| Gas mixture | EO: HCFC* |
| Temperature (°F, °C) | 130°F ± 5°F; 54°C ± 2°C |
| Relative humidity | 60 ± 20% |
| Pressure (p _b , kPa) | 58 ... 70
(0.56 ... 0.7 bar/8 ... 10 psi) |
| Exposure time | 120 mins |
| EO concentration | 600 ± 30 mg/L |
- * Gas mixture consisting of 10% ethylene oxide and 90% chlorotetrafluoroethane (HCFC-124; percent by weight).

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10.9.2 Formaldehyde/Hydrogen peroxide sterilization

The instruments can be used immediately after sterilization with formaldehyde or hydrogen peroxide sterilization. Extra airing time is not necessary.

Remember that for hydrogen peroxide sterilization there are restrictions concerning the lumen size and the material with respect to the instruments which can be sterilized with Sterrad® 100S, Sterrad® 200S, Booster or Sterrad® NX.

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10.9.3 Sterilization with the STERIS® SYSTEM 1

i **NOTE:** Video Urethro-Cystoscopes PDD are compatible with the STERIS® SYSTEM 1.

The STERIS® SYSTEM 1 is a sterile processing method for immersible surgical and diagnostic devices. The STERIS® SYSTEM 1 utilizes a liquid chemical process (primarily buffered peracetic acid) for the rapid, low temperature destruction of microorganisms on the surfaces of surgical instruments and devices. Please consult the operator's manual for safe handling instructions for the sterilant.

! **WARNING:** Failure to thoroughly clean Video Urethro-Cystoscopes PDD may result in an ineffective sterile process.

! **WARNING:** Failure to properly position the Video Urethro-Cystoscopes PDD so that all surfaces will be exposed to the sterilant or overloading the processing container may result in an ineffective sterile process and/or damage to the devices.

! **WARNING:** Always verify that the sterilant container is empty after the cycle is complete. If the sterilant container is NOT empty, the load CANNOT be considered sterile.

! **WARNING:** Endoscopes are NOT sterilized and/or adequately rinsed when a sterile processing cycle is CANCELLED.

! **CAUTION:** Remove the red pressure compensation cap (#11025E) from the vent port before using the STERIS® SYSTEM 1.

i **NOTE:** The latest relevant national rules on the use of the STERIS® SYSTEM 1 must always be observed.

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1. Be sure that the Video Urethro-Cystoscope PDD is completely clean using the methods described in the cleaning section of this manual.
2. Visually inspect and leak test the Video Urethro-Cystoscope PDD to ensure that it is working properly before processing in the STERIS® processor.
3. Position the Video Urethro-Cystoscope PDD in the STERIS® tray to ensure that all surfaces will be exposed to the sterilant.
4. Please consult the STERIS® operator's manual for complete operating instructions.
5. The Video Urethro-Cystoscope PDD is ready to use immediately at the completion of the STERIS® cycle.
6. Contact STERIS® Corporation for the most current STERIS® SYSTEM 1 processing options for KARL STORZ flexible endoscopes.

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Statements as to the efficiency of sterilization for the respective type of instrument (depending on the length and diameter of the lumens) can only be made by the manufacturer.

CAUTION: Attention must be paid to the information on preparation, sterilization and storage in the 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' instructions. The methods for cleaning, disinfection and sterilization are explained there in detail.

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
10.11 Assembly of the suction unit

CAUTION: Two O-rings (Art. no. 5925300 and 5925600) on the suction valve ② can be lost. Each time the unit is dismantled and cleaned, disinfected and sterilized, inspect the O-rings and, if necessary, install new ones.

CAUTION: The O-ring (Art. no. 5928000) on the suction cylinder can be lost. Each time the unit is dismantled and cleaned, disinfected and sterilized, inspect the O-ring and, if necessary, install a new one.

- ② 11900 AV Suction valve (with O-Ring 5925300 and 5925600)
- ③ 11900 VA Valve holder
- ④ 11900 AS (11900 ASE) Connecting tube
- ⑤ 11272 AZ LUER-Lock connectors (with O-Ring 5928000)

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
The suction unit must be connected as shown in the pictures above. Only then install it on the Video Urethro-Cystoscope PDD.

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Insert the valve holder into the recess and press it home.

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Fix the suction cylinder by pressing onto the port of the instrument channel.

CAUTION: The connecting tube and the nozzle for the connecting tube of the suction cylinder must be seated in the guide of the Video Urethro-Cystoscope PDD.

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Screw on the suction valve by turning clockwise.

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10.12 Disinfection and sterilization of the biopsy forceps and the cleaning brush

Thorough cleaning is required before disinfection and sterilization. Free the biopsy forceps from any visual soiling and clean in an ultrasound bath. Take special care to remove residues such as blood, excrement and mucus from the biopsy forceps and cleaning brush. Disinfection can be performed either manually or using a machine.

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Rinse in microbiologically clean/sterile water and dry (compressed air can be used). Before sterilizing, carefully oil the moving parts. For this, only use suitable steam-resistant, biocompatible, residue-free oils (e.g. 27656 B).

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After oiling, open and close the forceps several times.

Disinfect accessories and, if necessary, sterilize with EO gas. After each cleaning, check that the instruments function correctly and are not damaged. The biopsy forceps 11001 KL and cleaning brush 27651 AL can be autoclaved.

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11 Storage

CAUTION: Do not expose the Video Urethro-Cystoscope to direct sunlight or excessive heat.

CAUTION: The sheaths of the Video Urethro-Cystoscope must never be kinked or coiled up too tightly.

CAUTION: The distal tip of the Video Urethro-Cystoscope PDD must be in the neutral position during storage.

CAUTION: Do not keep the Video Urethro-Cystoscope PDD in the carrying case.

CAUTION: The sheaths of the Video Urethro-Cystoscope must never be kinked or coiled up too tightly. It is best to suspend the Video Urethro-Cystoscope for storage or to keep it in a large container. Video Urethro-Cystoscopes should be packaged and stored in sterile cloth or suitable containers in accordance with EN 868, and must be absolutely dry before being put into storage. It is most suitable to store the Video Urethro-Cystoscope suspended vertically with the distal end pointing down.

The room where they are stored should be at a normal temperature and humidity.

Store dismantled valves and caps separately.

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15 Limitation of liability

KARL STORZ SE & Co. KG shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the conditions that all assembly, operations, system expansions, readjustments, modifications, or repairs to same have been performed by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

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The medical device must always be sent to your local subsidiary (see 'Subsidiaries' section), even during the warranty period.

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18 Important information

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and sterilized/disinfected before they are sent for repair. We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

KARL STORZ gives no warranty on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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19 Disposal

To dispose of Video Urethro-Cystoscopes PDD, no special measures are necessary. National laws and regulations must be observed.

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KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongguiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

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Flexibles Video-Choledochoskop Serie 11292 VP/VPU

INSTRUCTION MANUAL

Flexible Video Cysto-Urethroscope series 11272 VN/VP/VPU/VNU

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This instruction manual is intended to serve as an aid in the proper handling, cleaning, disinfection and sterilization of the video endoscope. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for ready reference.

 **WARNING:** KARL STORZ instruments are provided unsterilized, and must be cleaned, disinfected and/or sterilized prior to the initial use and before each subsequent use.

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



**Controls, displays, connectors
and their uses
– Video Cysto-Urethroscope –**

- ① User buttons
- ② Light source connector
- ③ Video connector
- ④ Connection for leakage tester or pressure compensation cap
- ⑤ Protective cap for video connector
- ⑥ Distal tip
- ⑦ Sheath
- ⑧ Handpiece
- ⑨ Working channel with LUER-lock connectors
- ⑩ Deflection lever
- ⑪ Biopsy Forceps 27023 ZE
- ⑫ Cleaning Brush 27651 B
- ⑬ Leakage Tester 13242 XL
- ⑭ Pressure Compensation Cap 11025 E
- ⑮ Grasping Forceps, for small fragments
- ⑯ Video Connecting Cable

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- ① User buttons
- ② Light source connector
- ③ Video connector
- ④ Connection for leakage tester or pressure compensation cap
- ⑤ Protective cap for video connector
- ⑥ Distal tip
- ⑦ Sheath
- ⑧ Handpiece
- ⑨ Working channel with LUER-lock connectors
- ⑩ Deflection lever
- ⑪ Biopsy Forceps 11001 SL
- ⑫ Cleaning Brush 27651 B
- ⑬ Leakage Tester 13242 XL
- ⑭ Pressure Compensation Cap 11025 E
- ⑮ Grasping Forceps 11002 SS
- ⑯ Applicator and Guiding Tube 11272 X
- ⑰ Applicator, for use with Guiding Tube 11272 X
- ⑱ Stone Basket 27023 VK

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**Important information for users of
KARL STORZ instruments.....I**

**Controls, displays, connectors
and their usesII**

Video Cysto-Urethroscope II
Video CholedochoscopeIV

Safety information3

Normal use.....3
User qualification4
Warnings and cautions.....4
Safety precautions.....5

Handling9

Unpacking the device.....9
Basic equipment9
Operating conditions10
Electromagnetic Compatibility.....10
Operation11
Holding correctly11
Handling the sheath11
Handling the deflection mechainsim.....12
Handling the accessories.....14
Compatible camera control units/
light sources.....15
Connection to the TELE PACK/TELE PACK X16
Connection to the light source.....18
Connection to the camera control unit.....19
User button21
Preparation prior to use21
Inspecting the surfaces.....21
Checking the deflection22
Inspecting the light source.....23
Inspecting the camera control unit.....23
Inspecting the instruments and
the instrument channel24

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Cleaning, Disinfection, Care and Sterilization.....	25
General information	25
Preparation for cleaning.....	28
Leakage test.....	29
Cleaning the LUER-Lock stopcocks.....	31
Manual cleaning/disinfection of the video endoscope.....	33
Disinfection.....	36
Machine cleaning/disinfection of the video endoscope.....	38
Inspection after cleaning/disinfection	40
Chemical sterilization	41
Gas sterilization	41
Formaldehyde/Hydrogen peroxide sterilization	43
Sterilization with the Steris® System 1	44
Cleaning the forceps and the cleaning brush.....	46
Overview – Preparation.....	47
Storage	48
Maintenance	49
Servicing and repair.....	49
Maintenance	49
Limitation of liability.....	49
Warranty.....	50
Directive compliance	50
Important information	51
Disposal	51
Appendix	52
Troubleshooting.....	52

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Normal Use

Flexible KARL STORZ video endoscopes are optoelectronic instruments for medical, diagnostic and therapeutic use on human patients.

Video endoscopes are applied parts, which should only be used with supply units approved by KARL STORZ (camera control units, light sources).

Video Cysto-Urethrosopes:

Video cysto-urethrosopes are used to aid endoscopic visualization during diagnosis and therapy of the bladder and ureter.

Indication:

Examination, diagnosis and/or therapy in conjunction with endoscopic accessories in urology.

Video Choledochoscopes:

Video choledochoscopes are used in diagnostic and operative choledochoscopy to illuminate the bile duct.

Indication:

Examination, diagnosis and/or therapy in conjunction with endoscopic accessories in choledochoscopy.

Contraindications:

There are no known contraindications relating directly to the product.

The video endoscopes must not be used if, in the opinion of a qualified physician, such an application would endanger the patient, e.g. due to the patient's general condition, or if the endoscopic method as such is contraindicated.

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User qualification

The KARL STORZ video endoscopes may only be used by persons with an appropriate medical qualification and who are acquainted with cystourethroscopic and/or choledochoscopic techniques. The information given in these instructions only serves to instruct in the correct handling and preparation of the video endoscopes. It is not intended as an introduction to the techniques of cystoscopy and/or choledochoscopy.

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Warnings and cautions

Please read this manual carefully and follow the instructions exactly. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To further emphasize the words **Warning** and **Caution**, they are preceded by a pictogram.

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CAUTION: A Caution indicates that particular service procedures or safety precautions must be followed to avoid any damage to the product.

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Safety precautions

The video endoscopes may only be used according to the medical rules and procedures of endoscopy recognized for this method.

i NOTE: The patient leakage current limit of the CF application part may be exceeded when CF application parts (application on heart) are used in combination with other BF application parts. In this case, the reduction of the patient leakage current should be effected using only CF application parts as far as this is possible.

i NOTE: The video endoscope does not offer any electric insulation against high frequency (HF) currents. Protection against high frequency currents should be ensured through the use of correct HF instruments. Please read the instructions for use in the instruction leaflet supplied with the HF application part for more information.

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WARNING: Read the following safety instructions carefully before using the video endoscope to avoid putting your patients, personnel or yourself at risk.

WARNING: Inspect the video endoscope prior to each surgical procedure to ensure that it functions correctly and that it has been correctly cleaned, disinfected and/or sterilized.

WARNING: Every time the video endoscope is used, you must first check that it, and any accessories used in combination with it, is in perfect condition. Damaged video endoscope or damaged accessories must not be used.

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WARNING: The pressure compensation cap must be removed from the vent port for:

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WARNING: Do not allow any live components or instruments to come into contact with the video endoscopes.

WARNING: The video endoscope must not be used during the discharge of a cardiac defibrillator. Remove the video endoscope from the operation area before discharge takes place.

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WARNING: Danger of explosion if flammable gases are used in the immediate proximity of the video endoscope. The video endoscope must not be used in the presence of flammable or combustible media.

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WARNING: The pressure compensation cap must be in place on the vent port for:

- Chemical sterilization
- Transport/shipping

WARNING: The instructions 'Cleaning, Sterilization and Care of KARL STORZ Instruments' must be followed. There the procedures for cleaning, disinfection and sterilization are explained in detail.

WARNING: The video endoscope may only be repaired by personnel authorized by KARL STORZ. Unauthorized repairs may impair functioning and safety, and are therefore strictly forbidden.

KARL STORZ gives no warranty for video endoscopes that are repaired by unauthorized personnel.

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WARNING: Always unplug the video endoscope before carrying out any maintenance and cleaning work on it.

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WARNING: Risk of infection during shipment! Contaminated video endoscopes can cause serious or fatal diseases during transportation or at the recipient's premises. Defective video endoscopes with attached pressure compensation cap must be wiped down with disinfectant and/or gas sterilized before shipment.

WARNING: If a malfunction should occur during application on a patient, stop the application immediately. Place the distal tip of the video endoscope in the neutral, nondeflected position and remove the video endoscope slowly and carefully from the patient.

CAUTION: Only use accessories (guide probes, cleaning brushes etc.) that are specifically suited for the intended purpose and can be employed in line with generally accepted technical standards and safety guidelines.

CAUTION: Do not expose the video endoscope to direct sunlight or excessive heat, to avoid damage to the device.

CAUTION: Video endoscopes must never be exposed to X-rays, to avoid damage.

CAUTION: Do not steam sterilize (autoclave) the video endoscopes.

CAUTION: The maximum permissible temperature for decontamination, sterilization, transport and storage is 65 °C.

CAUTION: Do not hit the distal tip of the video endoscope against hard objects, and never use a sharp object to remove dirt.

CAUTION: Never lay heavy objects on the video endoscope.

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Check for missing items and evidence of shipping damage.
File any complaints immediately with KARL STORZ or the supplier.

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Basic equipment

Video Cysto-Urethroscope

Video Cysto-Urethroscope 11272 VP/VN/VPU/VNU with the following accessories:

- 1 x Carrying Case 27677 VC
- 1 x Instruction manual
- 1 x Biopsy Forceps 27023 ZE
- 1 x Grasping Forceps 27023 FE
- 1 x Video Connecting Cable 20213070
- 1 x Pressure Compensation Cap 11025 E
- 1 x Leakage Tester 13242 XL
- 1 x LUER Adaptor 27014 Y
- 1 x Cleaning Brush 27651 B

Video Choledochoscope

Video Choledochoscope 11292 VP/VN with the following accessories:

- 1 x Carrying Case 27677 VC
- 1 x Instruction manual
- 1 x Biopsy Forceps 11001 SL
- 1 x Grasping Forceps 11002 SS
- 1 x Stone Basket 27023 VK
- 1 x Applicator and Guiding Tube 11272 X
- 1 x Applicator for use with 11272 X
- 1 x Pressure Compensation Cap 11025 E
- 1 x Leakage Tester 13242 XL
- 1 x LUER Adaptor 27014 Y
- 1 x Seal for Working Channel (package of 10) 27550 N
- 1 x Cleaning Adaptor 27001 RA
- 1 x Cleaning Brush 27651 B

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Electromagnetic Compatibility

CE marked equipment has been tested and is found to comply with the EMC limits of EN/IEC 60601-1-2 :2001 CISPR 11 Class B.

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Handling

NOTE: Damage to the video endoscope resulting from incorrect operation is not covered by the manufacturer's warranty.

WARNING: The video endoscope is delivered unsterilized, so it has to be disinfected/sterilized prior to initial use and before each subsequent reuse.

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Operating

The following instructions outline the correct operation of the video endoscope and the corresponding accessories. This is not an introduction to the techniques of video endoscopy. Instructions on endoscopic technique can be found in the relevant medical literature. Appropriate accessories for your video endoscope, e.g. biopsy forceps, can be found in the catalogs for your respective specialty.

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Holding the video endoscope correctly

Grip the handle of the video endoscope in your left or right hand, as shown in the picture. With this kind of control grip, the video endoscope can be operated with one hand.

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Handling the sheath

CAUTION: Do not carry the video endoscope by the sheath. Do not pull, clamp or twist the sheath of the video endoscope.

CAUTION: The sheaths of the video endoscope must never be kinked or coiled up too tightly. It is best to store unsterile video endoscopes in a suspended position.

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Use a suitable water-soluble lubricant for the distal section of the sheath. Observe the anatomy the whole time during use. If you lose your orientation, pull the sheath back slightly into the neutral position until vision is restored.

If, during use, the image fails, is disturbed, or becomes blurred, or if vision in general is impeded, put the deflection mechanism into the neutral position and slowly pull the sheath of the video endoscope out of the patient's body.

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Handling the deflection mechanism

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WARNING: Suitable accessories, such as biopsy forceps, must only be introduced into the instrument channel of the video endoscope with the distal tip in the neutral, non-deflected position. Once the instrument has been introduced, carefully actuate the lever gently as far as the limit stop. The deflection lever must never be operated with force or abruptly, to avoid injuring the patient or damaging the device irreparably!

WARNING: Before removing the video endoscope from the patient, the distal tip of the video endoscope must be put into the neutral position, as otherwise injuries to the patient or damage to the instrument may occur.

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WARNING: If the distal tip can no longer be properly controlled (sticking etc.) via the deflection lever, or if the deflection mechanism develops any other fault, the video endoscope must not be used.

WARNING: If a functional defect occurs during use on the patient, discontinue use immediately. Put the distal tip of the video endoscope into the neutral, non-deflected position and remove the video endoscope slowly and carefully from the patient.

WARNING: The lens at the distal end of the video endoscope may become detached, e.g. due to mechanical forces. For this reason, it is essential to check that the lens is undamaged and present before and after every cleaning.

CAUTION: The resistance that can be felt when deflecting the tip using the deflection lever increases as the angle of deflection increases. Actuate the deflection gently and carefully, without any sudden or jerky movements, as far as the limit stop.

CAUTION: Do not use force to turn the deflection lever, to avoid either injuring the patient or damaging the instrument.

The distal tip can be deflected 210° upward or 140° downward by actuating the deflection lever slowly and carefully. Before use, the user must always have full control over the deflection and must develop a sure feel for the deflection mechanism and in particular for the neutral position.

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Handling the accessories (forceps, cleaning brushes, HF electrodes etc.)

! WARNING: When HF electrodes or laser probes are used, it must be ensured that power is not applied until the active parts have left the instrument channel and are in the user's field of vision, and the site of the operation is clearly recognizable. Activation in the instrument channel may lead to serious injury to the patient and damage to the video endoscope.

! WARNING: When HF electrodes or laser probes are used, it must be ensured that power is only applied when there are no flammable gases in the operating field or in the proximity of the video endoscope.

! CAUTION: Only use flexible accessories (forceps, cleaning brushes etc.) from KARL STORZ. Accessories from other manufacturers may damage the video endoscope

! CAUTION: Only introduce the instruments from the proximal port of the video endoscope.

! CAUTION: When laser probes are used, it must be ensured that the laser beam cannot hit the distal tip of the video endoscope.

! CAUTION: Do not force accessories, such as biopsy forceps, through the working channel, to avoid damaging the device. It should be possible to move the accessory easily, and without any hindrance, through the entire channel.

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- !** **CAUTION:** When the distal tip is deflected, do not introduce any instrument into the video endoscope. Return the distal tip to the neutral position before introducing the instrument and then move the deflection lever gently and carefully to the limit stop.
- !** **CAUTION:** Inserting sharp-tipped accessories into the instrument channel may damage the video endoscope.
- !** **CAUTION:** Foreign bodies or parts to be extracted must not be removed through the instrument channel using instruments. Instead, the entire video endoscope must be removed from the patient and the foreign body extracted.

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Compatible camera control units/light sources

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- i** **NOTE:** The KARL STORZ video endoscope should not be used in combination with camera control units and light sources from other manufacturers. The use of non-compatible units can lead to patient injury or unit damage.
- i** **NOTE:** It must be ensured that the electrically non-insulated parts of the video endoscope do not come into contact with conductive surfaces, e.g., voltage-carrying elements of other units, at any time.

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Connection to the TELE PACK/ TELE PACK X

-  **NOTE:** The video endoscope can be used with the TELE PACK/TELE PACK X. The instruction manual of the TELE PACK/TELE PACK X must be followed.
-  **CAUTION:** The video endoscope must be fully connected to the TELE PACK/TELE PACK X before the TELE PACK is switched on.
-  **WARNING:** Risk of burns and glare! Never look directly into the light emitted from the TELE PACK/TELE PACK X or the video endoscope. The high-intensity light emitted through the video endoscope may cause glare or high temperatures, and thereby burns.
-  **WARNING:** Danger of burns and fire hazard! Never place the end of the fiberoptic cable or the distal tip of the video endoscope on or underneath a sterile cover cloth. The intensity of the light can cause burns to the patient and/or the sterile cover cloth.
-  **WARNING:** Reduce the intensity of the TELE PACK/TELE PACK X so that optimal illumination of the operating field is still achieved. The higher the intensity of the light source, the greater the heat emitted at the tip of the video endoscope.
-  **WARNING:** The keyboard of the TELE PACK/TELE PACK X must not come into contact with the patient or be placed in the proximity of the patient.
-  **WARNING:** It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces, e.g., voltage-carrying elements of other units, at any time.

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Unscrew the cap for the video cable connector and remove it from the supply connector

Insert the video connection cable into the round port of the supply connector and the square port of the TELE PACK/TELE PACK X.

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TELE PACK/TELE PACK X einschalten.

i NOTE: Always put the exposure setting of the TELE PACK/TELE PACK X, in conjunction with the video endoscope, to 'Auto Border'. For this, see the TELE PACK/TELE PACK X instruction manual.

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Connecting the video endoscope

1. Insert the video endoscope's light connector as far as it will go into the light port.

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1. Insert the video endoscope's camera connector as far as it will go into the TELE PACK/TELE PACK X's camera connection.

2. Turn the device on by pressing the  button.

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Connection to the light source

- NOTE:** The video endoscope can be used with KARL STORZ light sources.
- WARNING:** The instruction manual of the specified light source must be followed.
- WARNING:** Danger of burns and glare! Never look directly into the light emitted from the light source or the video endoscope. The high-intensity light emitted through the video endoscope may cause glare or high temperatures, and thereby burns, at the distal light outlet.
- WARNING:** Danger of burns and fire hazard! Never place the distal tip of the video endoscope on or underneath a sterile cover cloth. The intensity of the light can cause burns to the patient and/or the sterile cover cloth.
- WARNING:** Reduce the intensity of the light source so that optimal illumination of the operating field is still achieved. The higher the intensity of the light source, the greater the heat emitted at the tip of the video endoscope.

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Connection to the camera control unit

CAUTION: The video endoscope can be used with the camera control unit TELECAM DX II 20 2330 20/TELECAM SL II 20 2130 20/ TELECAM SL II 20 FI 20213020-1 and KARL STORZ IMAGE 1 HUB™ 22200020-1/KARL STORZ IMAGE 1 HD 22202020-1XX/KARL STORZ IMAGE 1 HUB™ HD 22201020-1XX specified by KARL STORZ for this purpose. The instruction manual of the camera control unit must be followed. The complete operating instructions for the camera control unit can be found in the instruction manual.

CAUTION: Take care with the connectors on the video endoscope! The connectors on the video endoscope or video connection cable can be damaged if they are incorrectly plugged into the video endoscope. Make sure that the pins are correctly aligned.

CAUTION: The video endoscope must be fully connected to the TELE PACK/TELE PACK X before the TELE PACK/TELE PACK X is switched on.

WARNING: The keyboard of the camera control unit must not come into contact with the patient or be placed in the proximity of the patient.

NOTE: The video endoscope can be connected to the IMAGE 1 HUB™ HD/IMAGE 1 HUB™/IMAGE 1 HD with either the connecting cable 222000 77 (PAL) or 222001 77 (NTSC).

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1. Install the camera control unit on a firm, level surface.
2. Plug the power cord into the specified power socket.
3. Unscrew and remove the cap for the video cable connector.
4. Plug the video connection cable into the round port in the supply connector.
5. Connect the video connection cable to the camera control unit.

NOTE: Always put the exposure setting of the camera control unit, in conjunction with the video endoscope, to "Brightness Low".

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User button

The video endoscope has three control buttons which can be assigned specific functions in the camera control unit's device menu. The assignment set is retained after switching off and on again. Instructions on programming the control keys can be found in the instruction manual of the camera control unit or in the TELE PACK/TELE PACK X instruction manual.

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Preparation prior to use

NOTE: Before the device is used for the first time, it must be subjected to a complete preparation cycle, excluding preliminary cleaning.

WARNING: Every time the video endoscope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition. Damaged video endoscopes or damaged accessories must not be used.

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CAUTION: Do not squeeze the distal tip of the video endoscope too tightly.

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Every time the video endoscope is used, you must first check that its surface, and above all the surface of the sheath and distal tip, is in perfect condition and that there are no sharp edges. video endoscopes with sharp edges must not be used. Otherwise serious injuries may be caused to the patient or the user.

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Inspect the umbilical tube for damage (such as breaks, cracks, twisted or squashed sections...) and check the supply connector for loose connections.

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WARNING: Before use, always check that the deflection lever and the distal tip connected to it function correctly. Defective video endoscopes must not be used.

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
Inspecting the light source

Prepare the light source or TELE PACK/ TELE PACK X as described above and connect to the video endoscope. Emission of light at the distal end of the video endoscope should be recognizable. The light intensity can be adjusted manually via the light source.

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Inspecting the camera control unit

WARNING: If the image fails or there is interference during use in the patient, put the distal tip into the neutral position and remove the video endoscope from the patient's body gently and carefully. If the video endoscope is used in the patient after the image has failed or when there is interference, this may result in serious injury to the patient.

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Prepare the camera control unit or TELE PACK/ TELE PACK X as described above and connect to the video endoscope.

CAUTION: Ensure the CCU of the TELE PACK/TELE PACK X is off before connecting the video endoscope!

If an image is recognizable, carry out a white balance as explained in the instruction manual for the TELE PACK/TELE PACK X or the camera control unit.

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The user buttons can be individually programmed. Programming of the user button is described in the respective instructions for the camera control unit or the TELE PACK.

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Inspecting the instruments and the instrument channel

CAUTION: It should be possible to move the biopsy forceps easily, and without any hindrance, through the entire channel.

CAUTION: Accessories, such as biopsy forceps, which are bent or twisted may no longer function properly, are difficult to pass through the channel, and may cause serious damage to the channel. They must therefore be replaced!

Check the instrument channel by introducing the appropriate instrument (max. diameter of instruments: 5 Fr. or 1.7 mm).

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General information

WARNING: The video endoscope must be rinsed out immediately after an intervention in order to prevent fluids from becoming dried on. Dried-on fluids may impede preparation and possibly jeopardize any required sterilization. Furthermore, it prevents coagulum/blockage and consequent repairs and equipment failure. For the preliminary cleaning, suitable enzymatic cleaners should preferably be used (observing the required concentration), as they will not attack the surface.

WARNING: Before and after using the video endoscope, always ensure that it has been cleaned, disinfected and/or sterilized correctly and also correctly stored.

WARNING: The cleaning, disinfection and sterilization procedures described do not guarantee a reliable and complete preparation cycle when implemented on their own. This can only be achieved if personnel follow recognized and validated preparation methods.

WARNING: The local occupational safety guidelines must be observed and applied when cleaning contaminated video endoscopes and instruments. Aerosols should be removed by means of a suction device.

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- the instruments have been freed of all organic materials and chemical residues. Sterilization is only effective on clean, disinfected surfaces!
- only recognized and validated methods are used.
- any deviation from the recommended procedures is validated by the user.
- the recommended sterilization parameters are maintained. These are only valid with sterilization equipment that is properly maintained and calibrated.

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CAUTION: The instructions 'Cleaning, Sterilization and Care of KARL STORZ Instruments' must be followed. There the procedures for cleaning, disinfection and sterilization are explained in detail. A list of KARL STORZ approved products for preparing video endoscopes can be found on the Internet at "www.karlstorz.com".

CAUTION: When preparing and using the cleaning and disinfectant solutions, follow the manufacturer's instructions paying close attention to proper concentration and exposure times. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiologic range of action of the chemicals used.

CAUTION: Never clean video endoscopes with a metal brush.

CAUTION: Constant changing between various preparation methods is especially hard on the materials and must therefore be avoided. However, a one-off change of method, due to procurement of new cleaning equipment, for example, is possible.

CAUTION: Video endoscopes must not be cleaned in an ultrasound bath.

CAUTION: Video endoscopes must not be steam sterilized (autoclaved) under any circumstances.

CAUTION: Video endoscopes must never be immersed in normal saline solution, as even brief contact can result in pitting and corrosion.

CAUTION: Video endoscopes may only be rinsed using microbiologically clean/sterile water.

CAUTION: Decontamination in instrument washers is only possible if suitable procedures are followed.

CAUTION: The maximum permitted temperature is 65 °C for decontamination, sterilization, transportation and storage.

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- Chemical sterilization
(except for Steris® System 1)
- transport/shipping

Preparation must only be carried out by authorized, trained personnel. All occupational safety regulations must be observed (gloves, safety goggles masks, clothing...). Thorough cleaning of the video endoscope and all accessories after every use ensures interruption-free working and optimal results of disinfection and sterilization.

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CAUTION: Before immersing the instrument in solution, check that it is leakproof, as otherwise penetration of liquid may cause damage.

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NOTE: It is recommended to store unsterile video endoscopes in a dustproof place after preparation, preferably in a suspended position.

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Remove the supply connector of the video endoscope from the light source.

Pull the video connection cable off the video endoscope and the camera control unit.

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Put on the cap for the video cable connector and screw tight.

Ensure that the cap is screwed on properly and tightly.

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Leakage test successful: Continue.

Leakage test ambiguous: Hold the leakage tester and immerse the video endoscope up to the anti-kink protection. Bubble check: The absence of bubble formation indicates that the video endoscope is leakproof.

Leakage test unsuccessful: Perform wipe-down disinfection, fill in the card and return the video endoscope to the manufacturer.

For every further step, it is essential to ensure that the video endoscope is completely vented (setting >0<) and that the manometer has been removed from the connection to the leakage tester.

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Leakage test


! WARNING: Before and after each use, the video endoscope must be tested for leakage. Leaking video endoscopes must not be used. They must be sent to KARL STORZ for repair or replacement.

! WARNING: Always test the video endoscope for leakage before starting cleaning or disinfection, and in particular before immersing the video endoscope. Leaking video endoscopes must not be immersed, but must be wipe-down disinfected and then sent to the manufacturer with an accompanying letter for repair or replacement of the video endoscope. Additional gas sterilization must be carried out in accordance with the Infection Protection Act.

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i NOTE: Before inflating, check that the cap is attached properly and tightly.

The connecting tube of the leakage tester and the tester connection on the video endoscope must be dry! Attach the tester connection cap firmly to the leakage tester connection and lock it by turning 25° clockwise.

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Pump up the leakage tester until the maximum test pressure of 160 mmHg is reached. The rubber sleeve on the distal deflection part of the video endoscope inflates slightly.

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Wait five minutes. If the test pressure does not drop more than 10 mmHg, the device is water-tight and can be immersed. After venting and removing the leakage tester, cleaning and disinfection can begin.

Venting is then complete when the manometer needle points to > 0 <.

If the results of the leakage testing are ambiguous, we recommend immersing the video endoscope with the working sheath in the cleaning solution up to the anti-kink protection. The absence of bubble formation indicates that the video endoscope is leakproof. Vent and remove the leakage tester, continue preparation procedure.

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If a pressure loss occurs, this means that the video endoscope leaks and must not be used or immersed. In this case, send the video endoscope for repair or contact KARL STORZ.

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WARNING: Risk of infection due to shipment. Contaminated instruments can cause serious or fatal infectious diseases during transportation or at the recipient's premises. Gas sterilize the video endoscope with the pressure compensation cap mounted before shipment (Please attach a sterilization indicator for documentation purposes). Enclose a written note about any special features (such as use on TB patients) and the disinfection measures taken, to protect personnel.

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
Cleaning the LUER-Lock stopcocks

Clean the inlet and outlet nozzles with a brush (art. no. 27651 B); use water and a compressed air gun or cleaning gun if necessary.

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
Carefully clean the cock plug with a brush.

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
Clean the inlet and outlet nozzles with a brush.

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
Rinse with water and blow dry with sterile compressed air.

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Assembly

Lightly grease the running surfaces of the cock plug. To this end, spread a pea-sized amount evenly over the cock plug (special grease for stopcocks, Art. No. 27657).


NOTE: This thin layer of grease not only serves to lubricate the metal surfaces which are in contact with one another, but also makes the stopcocks leakproof. Oil is not suitable as a sealant.

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Insert the cock plug and fix in place using the spring cap. Tighten slightly. Check the function of the cock plug. Dismantle the cock plug again for sterilization.

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**Manual cleaning/disinfection of the
video endoscope**

CAUTION: The pressure compensation cap must be removed from the vent port for:

- endoscopic interventions
- immersing in cleaning/disinfection solution

CAUTION: The pressure compensation cap must be in place on the vent port for:

- chemical sterilization (except for Steris® System 1)
- shipping/transport

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Pull off the video connection cable. Put on the cap for the video cable connector and screw tight. Carry out leakage test. For this, see the section 'Leakage test' (p. 29).

CAUTION: Never use a needle or other sharp object for removing dirt particles from the device, as this may cause serious damage.

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Carefully clean the distal end to remove all dirt particles.
All subsequent preparation steps are carried out using the immersion method.
Intermediate rinses and the final rinse must be carried out using microbiologically clean/sterile water. Wipe the housing, sheath, and supply tube thoroughly with a lint-free disposable cloth, preferably moistened with an enzymatic cleaner. Do not use pushing and pulling movements, to avoid damaging the sheath. Subsequently, rinse the video endoscope thoroughly with microbiologically clean/sterile water.

CAUTION: Only pass the cleaning brush 27651 B through the instrument channel in one direction, and only pull it back once the brush head has emerged at the distal end. The brush must not be moved to and fro in the channel, to avoid damage.

Introduce the brush 27651 B supplied into the instrument channel and slide it forward to free the instrument channel of soiling, tissue remnants and secretions.

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WARNING: The lens at the distal end of the video endoscope may become detached, e.g. due to mechanical forces. For this reason it is essential to check that the lens is present before and after every cleaning.

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Before retracting the brush head that emerges at the distal end, clean it using another brush or with a sterile compress. Then pull the brush back through the instrument channel and repeat the brushing procedure several times.

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Carefully press the cleaning solution through the instrument channel, so that it flows out of the distal end of the channel. Repeat the procedure several times if necessary.

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Clean the distal lens, eyepiece and the light port on the supply connector using a lens cleaner or a cotton tip applicator soaked in 70 % isopropyl alcohol.

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CAUTION: Do not apply lens cleaner/ isopropanol to the tube connectors. This may cause damage to the video endoscope due to loss of tube elasticity.

When the cleaning procedure has been completed, remove residues from the cleaning equipment, such as brushes, and disinfect. Do not reuse damaged brushes.

CAUTION: Do not autoclave either the video endoscope or the accessories (except for accessories that are expressly identified as autoclavable, such as forceps etc.).

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**Machine cleaning/disinfection of the
video endoscope**

CAUTION: The **pressure compensation cap** must be **removed** from the vent port for:

- endoscopic interventions
- complete **immersion in solutions**

CAUTION: The **pressure compensation cap** must be **in place** on the vent port for:

- gas sterilization/hydrogen peroxide sterilization
- venting

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All flexible KARL STORZ video endoscopes can be connected to devices for fully automatic cleaning and disinfection. Before doing so, manual cleaning must be carried out as described above. The choice of rinsing programs must be agreed with the machine manufacturer. In most cases, special programs for flexible endoscopes are available, which are suitable for cleaning without damaging the material.

WARNING: Remember that for machine cleaning of flexible video endoscopes, the temperature must not exceed 65°C.

WARNING: Before you place video endoscopes in the washer, clean all parts of the video endoscopes following the instructions given above for manual cleaning, irrespective of whether they have been used.

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The machine must fulfill all the demands made of it in terms of process validation, cleaning performance, disinfection and material compatibility. Ensure that the required temperature and pressure (+65 °C, and 1 bar for channels, 0.3 bar leakage test) are not exceeded. It is not recommended to change between manual and automatic preparation.

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CAUTION: Before cleaning in fully automatic devices, always check the video endoscope for leakage. Carefully place the video endoscope in the tray of the machine. Connect the rinsing tubes to the channels and, if available, connect the leakage tester of the machine to the video endoscope via the adaptor. Use a cleaning program for flexible endoscopes, following the machine manufacturer's instructions.

CAUTION: After the program has run, the video endoscopes are at a high risk of corrosion due to the ambient conditions in the preparation machine. Once the program has ended, remove the video endoscopes from the machine immediately.

NOTE: Cleaning/disinfection should be followed by a final rinse with completely demineralized water.

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
After machine preparation, clean the distal lens, and the light port on the supply connector with a lens cleaner or a cotton tip applicator soaked in 70 % isopropyl alcohol.

When cleaning the distal lens, always wipe with a circular movement.

CAUTION: Do not apply lens cleaner/isopropanol to the tube connectors. This may cause damage to the video endoscope due to loss of tube elasticity.

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**Inspection after cleaning/
disinfection**

Inspect all the parts and surfaces of the video endoscope for damage and to ensure that they are completely clean and dry.

WARNING: Before reusing the device, replace any defective parts or send the video endoscope for repair.

If residues are still present after preparation, prepare the parts again, give them special treatment and separate out any parts that require repair.

Inspect plastic parts for discoloration, porosity and flexibility. If the parts differ significantly from when they were new, they must be replaced. This applies in particular to seals and rubber caps. Subject the instrument to a functional test. Check instruments for steerability and all channels for patency.

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Chemical sterilization

Gas sterilization

Sterilization methods using ethylene oxide are gentle and therefore recommended methods of sterilization.

! WARNING: Before commencing gas sterilization, the instrument must be cleaned, disinfected and completely dry or else successful sterilization will not be achieved.

! WARNING: For gas sterilization with ethylene oxide, the airing times stipulated by the device manufacturer must be observed for instruments, due to gas absorption of the materials. The airing times depend crucially on the processing procedure of the EO sterilizer (concentration, process control).

i NOTE: Observe national laws and regulations for the gas sterilization method.

i NOTE: The suction unit should preferably be disassembled before being gas sterilized.

! CAUTION: If the instrument is to be sterilized, the entire system must be opened by fitting the EO cap on the pressure compensation part to ensure pressure compensation.

! CAUTION: Attention must be paid to the information on preparation, sterilization and storage in the 'Cleaning, sterilization and care of KARL STORZ instruments' instructions. The methods for cleaning, disinfection and sterilization are explained there in detail.

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Formaldehyde/ Hydrogen peroxide sterilization

The instruments can be used immediately after sterilization with formaldehyde or hydrogen peroxide sterilization. Extra airing time is not necessary. Please observe the recommendations concerning validation given in the table on page 47. Remember that for hydrogen peroxide sterilization there are restrictions concerning the lumen size and the material with respect to the instruments which can be sterilized with Sterrad® 100S, Sterrad® 200S with Booster.

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 **WARNING:** Sterility is only assured if the manufacturer's instructions for the Sterrad® procedure are followed.

 **WARNING:** All deviations from the sterilization parameters recommended for the Sterrad® system must be validated by the user.

 **CAUTION:** Attention must be paid to the information on preparation, sterilization and storage in the 'Cleaning, sterilization and care of KARL STORZ instruments' instructions. The methods for cleaning, disinfection and sterilization are explained there in detail.

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Sterilization with the Steris® System 1

 **NOTE:** Video endoscopes are compatible with the Steris® System 1.

The Steris® System 1 is a sterile processing method for immersible surgical and diagnostic devices. The Steris® System 1 utilizes a liquid chemical process (primarily buffered peracetic acid) for the rapid, low temperature destruction of microorganisms on the surfaces of surgical instruments and devices. Please consult the operator's manual for safe handling instructions for the sterilant.

 **WARNING:** Failure to thoroughly clean video endoscopes may result in an ineffective sterile process.

 **WARNING:** Failure to properly position the video endoscopes so that all surfaces will be exposed to the sterilant or overloading the processing container may result in an ineffective sterile process and/or damage to the devices.

 **WARNING:** Always verify that the sterilant container is empty after the cycle is complete. If the sterilant container is NOT empty, the load CANNOT be considered sterile.

 **WARNING:** Endoscopes are NOT sterilized and/or adequately rinsed when a sterile processing cycle is CANCELLED.

 **CAUTION:** Remove the red pressure compensation cap (#11025 E) from the vent port before using the Steris® System 1.

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1. Be sure that video endoscopes are completely clean using the methods described in the cleaning section of this manual.
2. Visually inspect and leak test the video endoscopes to ensure that it is working properly before processing in the Steris® processor.
3. Position the video endoscope in the Steris® tray to ensure that all surface will be exposed to the sterilant. All light cable adaptors should be removed from the light post before processing.
4. Please consult the Steris® operator's manual for complete operating instructions.
5. The Video-Rhino-Laryngoscope is ready to use immediately at the completion of the Steris® cycle.
6. Contact Steris® Corporation for the most current Steris® System 1 processing options for KARL STORZ flexible endoscopes.

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**Cleaning the forceps and
the cleaning brush**

Thorough cleaning is required before disinfection and sterilization. Clean the forceps and brushes thoroughly with (preferably enzymatic) cleaners and a soft brush. Take special care to remove residues such as blood, excrement and slime from the biopsy forceps and cleaning brush. Forceps must be cleaned in the ultrasound bath.

When preparing and using the solutions for manual cleaning follow the manufacturer's instructions paying close attention to proper concentration and exposure time. In the case of flexible forceps, machine decontamination is generally preferable to manual.

Then rinse the channels and surfaces with microbiologically clean/sterile water and dry with compressed air. Before sterilizing, carefully oil the moving parts. For this, only use suitable steam-resistant, biocompatible, residue-free oils (e.g. 27656 B).

After oiling, open and close the forceps several times. Disinfect accessories and, if necessary, sterilize with EO gas. After each cleaning, check that the instruments function correctly and are not damaged. The forceps and cleaning brush 27651 B can be autoclaved.

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Storage

-  **CAUTION:** Do not expose the video endoscope to direct sunlight or excessive heat.
-  **CAUTION:** The distal tip of the video endoscope must be in the neutral position during storage.
-  **CAUTION:** Do not keep the video endoscope in the carrying case.
-  **CAUTION:** The sheaths of the video endoscope must never be kinked or coiled up too tightly. It is best to suspend the video endoscope for storage or to keep it in a large container. Video endoscopes should be packaged and stored in sterile cloths or suitable containers in accordance with EN 868, and must be absolutely dry before being put into storage. It is most suitable to store the video endoscopes suspended vertically with the distal end pointing down. The room where they are stored should be at a normal temperature and humidity. Store dismantled valves and caps separately.

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Maintenance

Preventive maintenance is not essential. Regular maintenance can contribute to identifying potential problems before they become serious, enhancing the instrument's reliability and extending its useful operating life. Maintenance services can be obtained from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year.

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Limitation of liability

As the supplier of this instrument, KARL STORZ shall be liable for the safety, operational reliability and performance of this instrument subject to the following conditions: assembly, system expansion, readjustment, modification, or repair have been performed by persons authorized by KARL STORZ and that the instrument is used in accordance with the instruction manual.

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Disposal

To dispose of video endoscopes, no special measures are necessary. National laws and regulations must be observed.

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KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongquiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

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KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: Fax: +63 2 31745-00
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp
www.karlstorz.com

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14. Cleaning, disinfection, sterilization	31
14.1 Reprocessing Overview 1, German.....	34
14.2 Reprocessing Overview 1, English	35
14.3 Reprocessing Overview 1, Spanish.....	36
14.4 Reprocessing Overview 2, German.....	37
14.5 Reprocessing Overview 2, English	38
14.6 Reprocessing Overview 2, Spanish.....	39
14.7 Reprocessing Overview 3, German.....	40
14.8 Reprocessing Overview 3, English	41
14.9 Reprocessing Overview 3, Spanish.....	42
15. Preliminary cleaning	43
15.1 Disassembling suction valve insert	45
16. Manual reprocessing.....	46
16.1 Manual cleaning.....	46
16.2 Manual disinfection	48
17. Machine reprocessing	50
18. Inspection after reprocessing	58
18.1 Assembling suction valve insert	59
19. Gas sterilization/plasma sterilization	60
20. Chemical sterilization.....	63
21. Storage	64
22. Accessories.....	65
23. Overview: Cleaning brushes.....	68
24. Normatives	70
24.1 Servicing and repair	70
24.2 Limitation of liability.....	70
24.3 Warranty	70
24.4 Directive compliance	70
25. Repair and replacement program... ..	71
25.1 Important information.....	71
Subsidiaries.....	72

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Important note for the users of Broncho-Fiberscopes and Intubation Fiberscopes

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It has been determined that Broncho-Fiberscopes and Intubation fiberscopes (see table next page) incorporate a lumen material that is only approved for use as an external surface material and not included in the STERRAD® sterilization systems FDA-cleared claims for use in flexible endoscope lumens. This could result in lower than expected levels of hydrogen peroxide within the lumen during the exposure phase of the sterilization cycle and compromise sterilization efficacy.

Based on the risk factors discussed above, the probability of a patient infection occurring from a cleaned and disinfected bronchoscope or intubation scope processed in the STERRAD® sterilization systems is defined as "remote". Since bronchoscopes and intubation scopes are semi-critical (B) devices and need not necessarily sterilization, KARL STORZ has validated many sterilization methods.

Available validated sterilization methods include:

- STERRAD® 100S with Booster,
- STERIS® System 1,
- STERIS® System 1E,
- STERIS® AMSCO® V-PRO® maX and
- ethylene oxide (EO) gas.

The instructions „Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments“, Version 1.0.0, No. 96216003DER, describe in detail the workflows and process-relevant parameters of the individual validated procedures. The procedure must be chosen taking into account applicable national requirements and in consultation with the device manufacturers.

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The methodologies for manual and automated cleaning and disinfection are already described in details in this instruction manual. When selecting the disinfectant the range of action should be taken into account with regard to microorganisms. Use only chemicals approved by KARL STORZ for processing. For a complete list, please visit www.karlstorz.com. When preparing and using solution, strictly observe the chemical manufacturer's instructions on concentration, temperature, exposure time, and shelf life.

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
Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, these products are the result of years of experience and great care in manufacture. You and your organization have decided in favor of modern, high-quality items of equipment from KARL STORZ.

This manual is intended to serve as an aid in the proper use, cleaning, disinfection and sterilization of fiberscopes manufactured by KARL STORZ. The fiberscopes are sensitive precision instruments which require particular care during reprocessing. To prevent damage during reprocessing, support is provided here for disassembly, cleaning, sterilization and assembly of the KARL STORZ fiberscopes. For this reason, read these instructions carefully and keep them as a ready reference.

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DCI®-Fiberskope (Direct-Coupled-Interface)

This instruction manual describes three fiberscope versions, distinguished from one another by the imaging system. If a section refers exclusively to one of these versions, it is thus appropriately labeled.

Telescope

Telescope with special filters for PDD/AF

DCI® fiberscopes (Direct Coupled Interface)

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


KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope for supply in form, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, figures and descriptions in this manual.

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
2. Controls

Connectors, accessories and their uses

495 F	Light cable adaptor for Wolf
495 G	Light cable adaptor for KARL STORZ, Olympus, Winter & Ibe
11025 E	Pressure compensation cap (red)
11301 CA	Lip valve, disposable version, pack of 20
13242 XA	Adaptor for leakage test in cleaning and disinfection machines (green)
29100	Plug for instrument channel, pack of 10, reusable
11301 CB	Suction valve insert, reusable version
5917800	Seal ring (suction valve insert)
6885391	Cap (suction valve insert)
5923110	Spring (suction valve insert)
6927791	Housing (suction valve insert)
8567290	Plunger (suction valve insert)
11301 CE	Suction valve insert, disposable version, pack of 20
11301 D1	LED battery-powered light source with fine screw thread
11301 D3	LED battery-powered source with quick-release thread
6927691	Three-way double stopcock adaptor with LUER-Lock
6011590	LUER cap
600007	Barbed tube, thin, with LUER-Lock

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13272 Mouth piece

11301 CF LIPP tube holder

11301 BM Adaptor for leakage test in Belimed cleaning machines

11301 CD Adaptor for machine cleaning the working channel (green)

11301 HH Adaptor for leakage test in BHT cleaning machines

11301 FF2 Adaptor for leakage test in ETD2/ETD3 cleaning machines with nonreturn valve

11301 GG2 Adaptor for cleaning the working channel in ETD2/ETD3 machines

11301 KK2 Connecting adaptor for cleaning the working channel in ETD3 cleaning machines

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
- 20260030 DCI® adaptor, for connecting endoscopes with DIN eyepiece to DCI® camera heads
- 20260031 DCI® eyepiece adaptor for connecting DCI® telescopes to camera heads with mounts for DIN eyepieces

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- ① Focusing ring (if available)
- ② Suction valve insert (if available)
Press = suction activated
Release = suction deactivated
- ③ Valve nozzle (if available)
- ④ Up/down control
- ⑤ Valve for pressure compensation
and leakage tester connection
- ⑥ Suction outlet
- ⑦ Light inlet piece
- ⑧ Instrument channel
- ⑨ Objective lens
- ⑩ Light outlet
- ⑪ Instrument in instrument channel

- ⑫ PDD filter changer:
ALA for special tumor marker
based on aminolevulinic acid
W for white light
AF autofluorescence mode

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3. Intended Use

Fiberscopes are used in diagnostic and surgical endoscopy. They must only be used by appropriately trained and qualified technical staff.

Indications: Examination, diagnosis and/or therapy in conjunction with endoscopic accessories in various medical fields.

Contraindications: No contraindications are known which directly relate to the product. The fibroscope must not be used if, in the opinion of a qualified physician, such an application would endanger the patient, e.g., due to the patient's general condition, or if the endoscopic method as such is contraindicated.

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3.1 Safety precautions when using the fiberscopes

Fiberscopes must be used according to the medical rules and procedures of endoscopy recognized for this method.

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4. User qualification

KARL STORZ fiberscopes may only be used by persons with an appropriate medical qualification and who are acquainted with the endoscopic technique.

The information given in these instructions only serves to instruct in the correct handling, cleaning, disinfection and sterilization of the fiberscopes. It is not intended as an introduction to the endoscopic technique.

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5. Warnings and cautions

Please read this manual and follow its instructions carefully. The words, **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the text should be carefully reviewed to ensure the safe and effective operation of this product. To make the words **Warning** and **Caution** stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or safety precautions must be followed to avoid possible damage to the product.

NOTE: A Note indicates special information about operating the product, or clarifies important information.

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WARNING: Before using the fiberscope, read the following safety instructions carefully, to avoid putting your patients, personnel or yourself at risk.

WARNING: Inspect the fiberscope prior to each surgical procedure to ensure that it is functioning correctly and that it has been correctly cleaned, disinfected and (where applicable) sterilized.

WARNING: This equipment may not be used during the discharge of a cardiac defibrillator. Remove all endoscopic equipment from patient contact and from the surgical site.

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WARNING: Verify the sound condition of the fibroscope as well as of accessories used in combination before every application. Damaged fiberscopes or damaged accessories must not be used.

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WARNING: Lengthy direct contact between the endoscope tip and surrounding tissue may cause tissue damage due to the potentially high temperatures on the tip of the endoscope and must thus be avoided.

WARNING: The fibroscope is delivered unsterilized, and must therefore be cleaned, disinfected and (where applicable) sterilized prior to initial use and each subsequent reuse.

WARNING: Remove the red pressure compensation cap from the vent port

- for every endoscopic intervention
- before immersion in liquid/chemical sterilization.

WARNING: Before the fibroscope is removed from the patient, the locking mechanism ⑤ must be released and the distal tip of the fibroscope put in the neutral, non-angled position. Otherwise the instrument may be damaged or the patient injured.

WARNING: If a malfunction should occur during application on a patient, stop the application immediately. Place the distal tip of the fibroscope in the neutral, non-angled position and remove the fibroscope slowly and carefully from the patient.

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CAUTION: Do not steam sterilize (autoclave) the fiberscopes.

CAUTION: The maximum permissible temperature for cleaning, disinfection, sterilization, transport and storage is 65 °C.

CAUTION: Never immerse fiberscopes in normal saline solution, since even only brief contact can lead to pitting and corrosion.

CAUTION: Material damage: if the outer casing of the fiberscope is cleaned with alcohol, this will damage the material.

CAUTION: Fiberscopes must not be cleaned in an ultrasound bath under any circumstances.

CAUTION: Place the red pressure compensation cap on the vent port before

- gas sterilization/plasma sterilization
- shipping

CAUTION: Only use KARL STORZ accessories (guide probes, cleaning brushes etc.). Accessories from other manufacturers may damage the fiberscopes.

CAUTION: Only use guide probes with flexible ends. Other types of guide probe may cause considerable damage to fiberscopes.

CAUTION: Make certain to read the 'Cleaning, sterilization and care of KARL STORZ instruments' instructions. Methods of cleaning, disinfection and sterilization are explained in detail there.

NOTE: Any damage to the fiberscope resulting from incorrect operation is not covered by the guarantee.

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NOTE: If the instrument is in the working channel, its additional rigidity may limit angling.

NOTE: Separate instruction manuals are available for flexible videoscopes; they are therefore not described in this brochure.

NOTE: Federal law restricts this device to sale by or on the order of a physician (USA only).

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6. Unpacking

Check for missing items and evidence of shipping damage. File any complaints to KARL STORZ or the supplier immediately.

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6.1 Basic equipment

1 Fiberscope	Art. no.
1 Suction valve insert, disposable version, pack of 20 inserts	11301 CE
1 Leakage tester	13242 XL
1 Pressure compensation cap	11025 E
1 Plug for instrument channel, pack of 10	29100
2 Adaptor for cleaning the working channel in cleaning machines (green)	11301 CD
1 Cleaning brush for 2.8 x 50 and for 2.8 x 70	27651 A*)
for 3.7 x 65	27651 AL*)
for 5.0 x 65	27651 B
for 5.6 x 65	27651 C
for 6.4 x 65	27651 C
1 Biopsy forceps (only bronchofiberscopes)	11001 KL
1 Grasping forceps (only bronchofiberscopes)	11002 KS
1 Three-way double stopcock adaptor with LUER-Lock for O ₂ connection (only intubation fiberscopes)	6927691
1 LIPP tube holder (only intubation fiberscopes)	11301 CF
1 Mouth piece (only gastro-, intubation and bronchofiberscopes)	13272
*) For the fiberscopes 11278 A/AI/AU/AIU use the cleaning brush 27651 A (27651 AL is not suitable for these fiberscopes).	

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6.1.1 Accessories

- | | |
|--|----------------------|
| 1 Suction valve insert,
reusable version | 11301 CB |
| 1 Battery light source | 11301 D |
| 1 Adaptor for leakage test in
cleaning and disinfection
machines (green) | 13242 XA |
| 4 Cleaning adaptors for
cleaning machines | 11301 GG2/FF2/HH/KK2 |

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7. Initial operation

WARNING: The fiberscope is delivered unsterilized, and must therefore be cleaned, disinfected and (where applicable) sterilized prior to initial use and each subsequent reuse.

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7.1 Fitting the suction valve

Place the suction valve insert ② on the valve nozzle ③ and turn until the projections on the side of the insert lie over the nozzle recesses.

Push gently until the suction valve insert clicks into the valve nozzle.

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7.2 Fitting the tube holder (only intubation fiberscopes)

Push the tube holder onto the fiberscope from distal and secure it by turning the two parts in opposite directions.

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7.3 Fitting the lip valve (only bronchofiberscopes)

Place the disposable lip valve onto the instrument channel ⑧ and tighten gently by turning it in a clockwise direction (plastic part).

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8. Handling the fiberscopes

The optical section of KARL STORZ fiberscopes consists of a fiberoptic bundle in a flexible and watertight casing, a light guide, as well as objective and eyepiece lenses. Like all optical systems, these parts are very delicate. For this reason, careful handling of the fiberscopes is essential for a long service life and a consistent quality. Therefore be aware of the following points:

- Never kink or tightly coil the fiberscopes and light guides, so that no optical fibers break off. It is therefore best to suspend the fiberscopes for storage or keep in large containers.
- Do not jerk the up/down control ④ or move forcefully.
- Do not insert the instruments or cleaning brushes forcefully through the channels. This may cause damage to the casing, especially if the tip of the instrument is deflected.
- Carefully protect the objective lens ⑨ at the distal end from damage. Do not hit. Never remove soiling with sharp-edged objects.
- Never store or disinfect the fiberscope together with sharp-edged objects to prevent damage to the exterior.
- Never place heavy objects on top of the fiberscope.

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NOTE: Every fiberscope supplied by KARL STORZ includes a photo showing the condition of the fiber bundle upon delivery. A copy of this photo is kept on file in our quality control office. The wedge-shaped marking (at 12 o'clock) serves as a guide for the user.

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9. Using the fiberscopes

The following instructions are recommended for the proper operation of the fiberscope and accessories. This is not a reference for endoscopic techniques.

Consult the appropriate medical literature for information on endoscopic techniques.

You can find accessories for your fiberscope (catheters, flexible scissors and forceps, guide wires, stone baskets, lithotripter probes, laser fibers, electrodes etc.) in the appropriate KARL STORZ catalog.

WARNING: Always adjust the light source to the minimum illumination intensity necessary to achieve optimum visualization of the operating field. The higher the intensity setting of the light source, the greater the heat energy that will be generated at the tip of the fiberscope.

WARNING: Never lay the end of the fiber optic cable or the distal tip of the fiberscope on or under a sterile cover cloth. The light intensity may lead to burns on the patient and/or sterile cover cloth.

- When connecting and operating light sources and suction/irrigation equipment, be certain to read the appropriate instruction manuals for the equipment. The irrigation equipment recommended by KARL STORZ is set to a maximum pressure of 54 kpa (400 mmHg/0,54 bar/7.8 psi).
- Endoscopic examinations must be carried out under aseptic conditions.
- When introducing the fiberscope into the patient, the distal tip of the fiberscope must be in a straight, nondeflected position (control lock-in can be felt and heard). The fiberscope must be handled according to the recognized medical rules and procedures for endoscopic methods.

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- If no instrument is inserted, close the instrument channel proximally with a plug (Art. no. 29100).
- If an instrument is to be inserted, position the disposable lip valve (Art. no. 11301 CA) to seal off the instrument channel proximally.
- Instruments especially designed for this may be inserted through the instrument channel of the fiberscopes. Make certain the accessories have a suitable length and diameter. Unsuitable accessories may damage the fiberscopes. For questions on accessories, please contact KARL STORZ directly or an authorized dealer.

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NOTE: When the lip valve is first used, the instrument is subjected to increased resistance. This resistance will decrease the more often the lip valve is used.

- When inserting an instrument into the instrument channel, the fiberscope must be in a straight, nondeflected position.

CAUTION: When inserting instruments, increasing resistance is felt with increasing deflection. Forced introduction of instruments can lead to damage to the fiberscope or the instrument channel.

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CAUTION: If undue resistance is encountered when inserting an instrument, do not attempt to force the instrument. Withdraw the instrument and apply a water soluble lubricant to the tip of the instrument and attempt to insert again. Prior to removal of the instrument from the patient, the distal tip must be in a straight position.

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- When using forceps, scissors or other accessories that can be opened or closed, ensure that the instrument is completely closed prior to insertion into the instrument channel. Do not attempt to open the jaws until the distal tip is completely free of the channel. Opening the jaws inside the instrument channel can lead to damage to the fiberscope.

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WARNING: When using HF electrodes or laser probes, make sure that the active parts are clear of the instrument channel and in full view of the operating surgeon before power is applied. Activation within the instrument channel may damage the fiberscope and/or injure the patient.

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WARNING: Make certain to read the instructions for the laser unit and/or HF generator.

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CAUTION: If accessories with sharp tips, such as EHL probes or laser probes, are inserted into the instrument channel of a deflected fiberscope, this can damage the fiberscope.

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CAUTION: When using a stone basket, the stone must not be retracted through the instrument channel but instead the entire fiberscope must be removed. Once the stone is removed, the stone basket may then be taken out of the instrument channel.

CAUTION: The fiberscopes must be cleaned and then reprocessed immediately after the intervention in order to prevent incrustations. Dried-on incrustations may jeopardize successful decontamination and disinfection.

- A leakage test must be carried out after every intervention, to prevent damage during reprocessing.

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CAUTION: Use only KARL STORZ accessories (guide probes, cleaning brushes etc.). Accessories from other manufacturers may damage the fiberscopes.

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NOTE: If the instrument is in the working channel, its additional rigidity may limit maneuverability.

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
10. Preparation of the fiberscope

Before using the fiberscope on a patient, exterior damage and operational problems must be eliminated.

Check the outer sheath of the fiberscope. Be especially aware of visible damage to the casing, particularly in the area of the fiberscope tip.

Look through the eyepiece and make certain that you have a clear image. The grid (Moiré effect) should be easily recognizable. If the image is not clear but cloudy, clean the eyepiece and lens with a cotton applicator soaked in alcohol, maintaining hygienic conditions at all times.

DCI® fiberscope image quality can be controlled on screen. If necessary, clean image waveguides, light cables and eyepieces as described above.

If it is still cloudy, it is possible that the objective lens is scratched or that there is moisture inside the telescope. The fiberscope must not be used in this case, but instead sent in for repairs.

! CAUTION: Before every application on the patient, check the instrument for leaks.

Check the fiberscope for leaks before every application on a patient.

Make certain to read the specifications on page 43 of this instruction manual.

10.1 Connect these supply lines

- If necessary, remove the adaptor 495 G or 495 F from the light inlet piece and connect the fiberoptic light cable or battery-powered light source to the light inlet piece ⑦.

i NOTE: Not necessary for battery-powered light sources with coarse threads (11301 D3).

! WARNING: The instruction manuals of the supply units used must be observed.

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- suction to suction outlet ⑥
- or O₂ supply (under anesthesia) to instrument channel ⑧ via the appropriate adaptor

CAUTION: The instruction manuals of the supply units used must be observed.

Switch on the respective supply units.

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Insert the DCI® fiberscope into the DCI® camera head until it locks into place.

Insert the light cable (Art. no. 495 DV) into the DCI® camera head as far as it will go.

WARNING: The instruction manuals of the supply units used must be observed.

Connect the light cable and the video signal cable up to the appropriate supply units.

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10.2 Controlling the distal tip

The up/down movement of the distal tip on the fiberscope can be adjusted at the control ④: move smoothly and slowly.

Some fiberscopes have:

- contrapositive deflection direction (urology);
- a locking mechanism ⑮ for securing the fiberscope tip (gynecology).

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Models with locking mechanism ⑮:

To secure the deflection position of the distal tip, turn the locking mechanism ⑮ slowly to distal until resistance can be felt.

CAUTION: Forcing the locking mechanism ⑮ further in the same direction may damage the locking mechanism ⑮.

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! WARNING: Before the fiberscope is removed from the patient, the locking mechanism 15 must be released and the distal tip of the fiberscope put in the neutral, non-deflected position. Otherwise the instrument may be damaged or the patient injured.

To release the locking mechanism 15 from the secure position, turn it slowly to proximal.

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10.3 Suction

Press the suction valve 2 to start suction (when the suction pump has been connected and switched on).

! CAUTION: There should be no instrument in the suction channel during suction. The instrument channel inlet should be closed so that no secondary air is drawn in.

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10.4 Insufflation

Under anesthesia, connect the O₂ hose to the instrument channel:

- using a barbed tube (Art. no. 600007) on the instrument channel or
- if you are using the three-way double stopcock adapter with LUER-Lock (Art. no. 6927691), attach the hose to the tube holder on the angled part of the adapter. This ensures that the straight-line passage of the adapter remains free for instrument insertion.

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10.5 Focusing

The image is focused with the focusing ring 1.

i NOTE: When connecting photo or video systems recommended by KARL STORZ, the fiberscope focusing ring must be at the zero position (mark).

DCI® fiberscopes are focused using the adjustment wheel of the DCI® camera (see page 28).

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10.6 Instrument channel

Flexible instruments can be inserted through the instrument channel ⑧ as required.

If the instrument channel is not being used: replace the sealing cap.

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10.7 PDD version

Fiberscopes intended for use in Photodynamic Diagnosis (PDD) or autofluorescence (AF) are illuminated with a special diagnosis light (D-Light C). Adjust the PDD filter changer ⑫ for special filters designed for this light and possibly for marker substances:

ALA	for a special tumor marker based on aminolevulinic acid
W	for white light
AF	autofluorescence mode

NOTE: If the filter changer is set to 'ALA' or 'AF', a yellowish image is produced when viewing through the integrated special filter.

NOTE: No filter changer is available on PDD fiberscope art. no. 11272 CI. The filter (ALA) is permanently installed.

For optimal illumination, use the fluid light cable 495 FP (2.5 m) or 495 FO (1.8 m) between the light source and PDD fiberscope.

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11. DCI® fiberscope with camera

11.1 Assembly

Slide the DCI® into the snap-in DCI® coupling mechanism of the DCI® camera until the coupling snaps in automatically.

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Connecting the special light cable 495 DC or 459 DV to the DCI® camera head:

Insert the light cable plug into the DCI® camera head until it locks into place.

While doing this, hold the light cable by the handle piece. Never pull on the cable!

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NOTE: The entrance and exit faces of the light cable and endoscope must be kept clean if optimal light delivery is to be maintained. Clean the surfaces with a cotton applicator soaked in cleaning solution, then disinfect it with a cotton applicator soaked in 70 % isopropyl alcohol and dry.

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DCI® fiberscopes are focused using the adjustment wheel on the DCI® camera.

The camera control unit is controlled using the two programmable control keys (silver, blue) on the DCI® camera head.

NOTE: Further controls are described in the instruction manual of the camera control unit used.

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11.2 Disassembly

Hold the DCI® camera and the DCI® fiberscope still. Push the release button on the DCI® camera: the 'snap-in' DCI® coupling mechanism is released.

Pull the light cable out of the DCI® camera head.

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12. DIN eyepiece on DCI® fiberscope

The DCI® eyepiece adaptor 20260031 ensures the connection of DCI® telescopes to camera heads with DIN eyepiece mounts.

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12.1 Assembly

Slide the DCI® fiberscope into the 'snap-in' DCI® coupling mechanism of the DCI® eyepiece adaptor until the coupling snaps in automatically.

Connect a light cable or a battery-powered light source.

WARNING: The instruction manuals of the supply units used must be observed.

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DCI®

12.2 Disassembly

Remove the light cable or battery-powered light source.

Hold the DCI® eyepiece adaptor and the DCI® fiberscope. Press the release button on the DCI® eyepiece adaptor. The 'snap-in' DCI® coupling mechanism is released.

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13. DCI® camera on DIN eyepiece

13.1 Assembly

Put the focusing ring (if available) of the fiberscope in the zero position (mark).

Slide the DCI® adaptor 20260030 into the 'snap-in' DCI® coupling mechanism of the DCI® camera until the coupling snaps in automatically.

Connect the standard light cable at the light inlet piece ⑦ of the fiberscope.

NOTE: The special light cable 495 DC or 495 DV cannot be used and may have to be removed from the DCI® camera.

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Opening the instrument coupling: turn the outer ring against the spring force, insert the fiberscope DIN eyepiece and let go of the outer ring. The instrument coupling closes automatically as a result of the spring force.

CAUTION: No further tightening of the outer ring to position the eyepiece is required. This may, among other things, cause damage to the instrument coupler!

13.2 Disassembly

Hold the DCI® camera and the fiberscope still. Open the instrument coupler for the DCI® adaptor and remove the fiberscope.

Hold the DCI® camera and the DCI® adaptor still. Push the release button on the DCI® camera: the 'snap-in' DCI® coupling mechanism is released. Remove the DCI® adaptor.

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14. Cleaning, disinfection, sterilization

Cleaning, disinfection, care and sterilization

! **CAUTION:** The manual 'Cleaning, sterilization and care of KARL STORZ instruments' must be followed. There, the cleaning, disinfection and sterilization procedures are explained in detail.

! **CAUTION:** For cleaning and disinfection, use a product from the KARL STORZ approved chemicals list. The current list of approved products can be found on our website (www.karlstorz.com).

! **CAUTION:** When preparing and using the solutions, follow the chemical manufacturers' instructions paying close attention to proper concentration and exposure times. Prolonged immersion and incorrect concentrations may result in damage. Bear in mind the microbiological range of action of the chemicals used.

! **WARNING:** The stated cleaning, disinfection and, where necessary, sterilization methods alone do not guarantee proper reprocessing or sterility. This can only be achieved if the staff follows recognized and validated reprocessing methods.

! **WARNING:** For all cleaning and disinfection of contaminated endoscopes and instruments, observe the guidelines of the Employers' Liability Insurance Association and equivalent organizations.

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WARNING: Remember that sterilization can only be performed successfully if

- all organic and chemical residue on the instruments has been eliminated. Sterilization is only effective on clean, disinfected surfaces!
- work proceeds according to recognized and validated methods.
- any deviation from recommended procedures must be validated by the user.
- the recommended sterilization parameters are maintained. These are only valid with sterilization equipment that is properly maintained and calibrated.

CAUTION: Never clean fiberscopes with a metal brush.

CAUTION: Constantly changing between different reprocessing methods causes particular wear on the materials and must therefore be avoided. A once-only change in methods however, e.g., when acquiring new cleaning equipment, is perfectly safe.

CAUTION: Fiberscopes must not be cleaned in an ultrasonic bath.

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CAUTION: Never immerse fiberscopes in normal saline solution, since even brief contact can lead to pitting and corrosion.

CAUTION: Final rinsing during manual reprocessing of fiberscopes should be carried out using microbiologically pure/sterile water.

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Reprocessing Overview 1										
Product	Cleaning			Disinfection	Sterilization					Remarks
	manual	Cleaning machine	Cleaning machine		Steam	Gas		Plasma	Chemical	
		60 °C	93 °C		134 °C, fractionated	E0	F0	STERRAD	STERIS	
110xx ..., 111xx ..., 112xx ..., 113xx ... Fiberscopes	✓	✓	-	✓	-	✓	✓	✓	✓	
495 F/G; Light cable adaptor	✓	✓	✓	✓	✓	✓	✓	✓	✓	
11025 E; Pressure compensation cap	-	-	-	-	-	✓	✓	✓	-	
11301 CD; Adaptor for cleaning working channel in cleaning machines	✓	✓	✓	✓	-	✓	✓	✓	✓	
11301 CF; LIPP tube holder	¹ ✓	¹ ✓	¹ ✓	✓	✓	✓	✓	✓	✓	¹) dismantle for cleaning
11301 D1, 11301 D3; Battery light source (accessory)	⁵ ✓	-	-	⁵ ✓	-	✓	✓	✓	✓	⁵) Battery compartment with closed cover; liquid <15 cm deep; no pressure impact: rinse only with unpressurized liquids; do not dry with compressed air
11301 FF2; Adaptor for leakage test in ETD2/ETD3 cleaning machines	✓	✓	✓	✓	✓	✓	✓	✓	✓	

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Reprocessing Overview 2										
Product	Cleaning			Disinfection	Sterilization					Remarks
	manual	Cleaning machine	Cleaning machine		Steam	Gas		Plasma	Chemical	
		60 °C	93 °C			E0	F0	STERRAD	STERIS	
11301 GG2 ; Adaptor for cleaning irrigation and working channel in ETD2/ETD3 cleaning machines	✓	✓	✓	✓	✓	✓	✓	✓	✓	
11301 HH ; Adaptor for leakage test in BHT cleaning machines	✓	✓	✓	✓	✓	✓	✓	✓	✓	
11301 KK2 ; Adaptor for cleaning working channel in ETD3 cleaning machines	✓	✓	✓	✓	✓	✓	✓	✓	✓	
13242 XA ; Adaptor for leakage test in cleaning und disinfection machines	² ✓	² ✓	² ✓	✓	✓	✓	✓	³ ✓	✓	²) only PH neutral cleaners ³) change in color possible
13242 XL ; Leakage tester	-	-	-	-	-	-	-	-	-	
13272 ; Mouth piece	✓	✓	✓	✓	✓	✓	✓	✓	✓	
29100 ; Plug for instrument channel	✓	✓	✓	✓	✓	✓	✓	✓	-	
600007 ; Barbed tube, thin, with LUER-Lock	✓	✓	✓	✓	✓	✓	✓	✓	✓	
6011590 ; LUER cap	✓	✓	✓	✓	✓	✓	✓	✓	✓	
6927691 ; Three-way double stopcock adaptor with LUER-Lock	✓	✓	✓	✓	✓	✓	✓	✓	✓	Disassemble stopcocks before cleaning and sterilizing

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14.6 Reprocessing Overview 2, Russian
Preliminary cleaning

- Procedure not appropriate!
- ✓ Procedure appropriate!

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Reprocessing Overview 3										
Product	Cleaning			Disinfection	Sterilization					Remarks
	manual	Cleaning machine	Cleaning machine		Steam	Gas		Plasma	Chemical	
		60 °C	93 °C			EO	FO	STERRAD	STERIS	
11301 CB; Suction valve insert, reusable version	4✓)	4✓)	4✓)	✓	✓	✓	✓	✓	✓	4) dismantle for cleaning
5917800; Seal ring (suction valve insert)	✓	✓	✓	✓	✓	✓	✓	✓	✓	
5923110; Spring (suction valve insert)	✓	✓	✓	✓	✓	✓	✓	✓	✓	
6885391; Cap (suction valve insert)	✓	✓	✓	✓	✓	✓	✓	✓	✓	
6927791; Housing (suction valve insert)	✓	✓	✓	✓	✓	✓	✓	✓	✓	
8567290; Plunger (suction valve insert)	✓	✓	✓	✓	✓	✓	✓	✓	✓	
20260030; DCI® adaptor	✓	✓	-	✓	-	✓	✓	✓	✓	
20260031; DCI® eyepiece adaptor	✓	✓	-	✓	-	✓	✓	✓	✓	

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15. Preliminary cleaning

Remove any heavy soiling (secretions, etc.) from the outer casing of the fiberscope immediately after use. Use a disposable cloth moistened with cleaning solution (no solutions predominantly containing alcohol). Submerge the distal end in a container with cleaning solution and activate the suction to clean the working channel.

CAUTION: Before immersing the instrument in cleaning and disinfectant solution, check that it is leakproof, as otherwise penetration of liquid may cause damage.

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Connect the leakage tester (art. no. 13242 XL) to the valve for pressure compensation ⑤ and build up pressure until the pointer is in the blue marked range. Watch the pointer for approx. 30 seconds: the pressure must not drop during this time. Press the red button on the pressure valve and release the pressure. Remove the test unit from the fiberscope.

NOTE: During inflation, make certain that the red pressure relief button is completely closed.

If there is a loss in pressure, the fiberscope has a leak and must neither be used nor immersed in liquid. In this event, send the fiberscope in for repairs.

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WARNING: Risk of infection: In order to prevent infection, it is strictly forbidden to ship contaminated medical devices! All medical devices must be decontaminated on-site, to avoid contact and aerogenous infections among staff. To eliminate the risk of infection, before the fiberscope is dispatched it should be wiped-down disinfected and sterilized in the NT procedure with the pressure compensation cap in place. We reserve the right to return contaminated instruments/equipment to the sender.

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- Remove all instruments from the instrument channel ⑧.
- Remove adapter and other removable parts from the light inlet piece ⑦ and instrument channel ⑧.
- Remove suction valve insert ② from the valve nozzle ③: press the projections on the side together and pull out the suction valve insert.
- Dispose of disposable parts removed. If necessary, disassemble other parts further and reprocess separately.

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15.1 Disassembling suction valve insert

ⓘ NOTE: Only the multi-use suction valve insert (silver-coloured cap, metal plunger) is suitable for reprocessing.

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
16. Manual reprocessing

! **NOTE:** Before cleaning, remove both adaptors from the light inlet piece ⑦ to permit free access to the fiber surface.

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16.1 Manual cleaning

Submerge the instrument in a cleaning solution approved by KARL STORZ and clean it using a cleaning agent (brush or cleaning gun, Art. no. 27660). Avoid splashing with the cleaning agent.

! **CAUTION:** Pass the brush through the channel in only one direction and pull back only after the brush head has completely left the distal end. The brush must not be pulled back and forth within the channel, since this can lead to damage.

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Remove the brush from the channel carefully. Then thoroughly rinse the channels repeatedly using cleaning and disinfection agents. If required, add a cleaning agent booster.

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Detach the adaptor, clean the stopcocks and grease (see the 'Cleaning, sterilization and care of KARL STORZ instruments' instructions).

The adaptor 11301 CD can be used for manual rinsing of the suction/irrigation channel.

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Rinse to remove the cleaning solution.

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Clean the eyepiece and objective lens with a soft cloth, sponge or cotton swab, soaked with 70 % alcohol.

CAUTION: Cleaning with the special cleaning paste should only be carried out if the image as viewed through the endoscope is visibly cloudy or if there are visible deposits on the optical surfaces, it should certainly not be used routinely during reprocessing.

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Stubborn soiling and deposits may be removed with cleaning paste (art. no. 27661). To do this, a small amount of cleaning paste is placed on a clean, moistened cleaning swab and distributed over the glass surface. Through polishing, even stubborn deposits can be removed.

NOTE: The instrument must be reprocessed again following treatment with the cleaning paste!

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Finally, dry the surfaces of the fibers and the optical end faces with a cotton-tipped applicator soaked in 70 % alcohol.

NOTE: The fiber surface in the light inlet piece ⑦ must be dried off especially thoroughly with alcohol. Residue from disinfectants and detergents in the light inlet piece could bake on if the fiber optic light cable is connected and severely impair light transmission.

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16.2 Manual disinfection

All parts of the KARL STORZ fiberscope can be disinfected. For this purpose, materials approved by KARL STORZ and recommended for use with endoscopes should be used. For a complete list, please visit www.karlstorz.com.

For disinfection, we recommend a plastic basin (Art. no. 27647 V), in order to avoid scratches on the optical end faces in particular. A plastic basin also prevents electrolytic corrosion, which can occur if different metals are present in a solution.

For disinfection, submerge the fiberscope completely. Ensure that the lumina of the channels and sheaths are completely filled with the solution. There must be no more air (air bubbles) in the channels. Make sure that all the stopcocks, the valve insert and the valve cap are open or have been disassembled. Immerse the valve insert, adaptors etc. separately.

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Rinse out the instrument channel ⑧ with a syringe filled with disinfectant solution. Follow the chemical manufacturers' instructions for concentration and exposure times.

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For each rinse, use new microbiologically pure/sterile water. It is absolutely necessary to rinse the fiberscope with microbiologically pure/sterile water in order to prevent potential toxic impact from any remaining chemical residue, which could cause the instrument to become conductive. This is also required to optimize the results of sterilization and avoid damage to the surface.

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Dry the outside of the instruments with a lint-free, cloth and/or sterile compressed air.

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Remove any water in the instrument and suction channels using sterile compressed air.

CAUTION: The pressure must not exceed 50 kPa (0.5 bar/7.25 psi). The cleaning gun might therefore have to be operated with a pressure reducer.

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NOTE: If, according to the high risk classification and/or the Infection Protection Act, sterilization of flexible fiberoptic scopes is not required, the instruments should be stored safely in a hanging position free from dust (distal tip pointing down), so that any remaining moisture in the channels can drain out.

NOTE: For flexible fiberoptic scopes without channels, a hanging storage position is not essential.

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17. Machine reprocessing

After appropriate preliminary cleaning

After the appropriate preliminary cleaning (see p. 43), all KARL STORZ fiberscopes may be connected to devices for semi-automatic or automatic reprocessing. The choice of reprocessing programs must be agreed with the machine manufacturer. There are usually special programs available for fiberscopes which are suitable for cleaning and which are gentle on materials. For machine reprocessing, make sure that the fiberscopes are not heated above 65 °C.

CAUTION: Check the fiberscope for leaks before every reprocessing in semi- or fully automatic equipment.

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Leakage tester connection on washer/disinfector:

- ETD3/ETD2 washers:
Adaptor for leakage test (Art. no. 11301 FF2, see pp. 52/55)
- BHT washer/disinfector:
Adaptor for leakage test (Art. no. 11301 HH)
- Belimed washer/disinfector:
Adaptor for leakage test (Art. no. 11301 BM)
- On the fiberscope using the green adaptor included in delivery (Art. no 13 242 XA);
On the machine: consult the manufacturer of the washer/disinfector.

If the reprocessing machine cannot be connected to the leakage testing unit for any reason, check for leakage manually as described on page 43 or check the adaptor connection using the 13242 XA.

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For machine reprocessing, attach the green cleaning adapter (Art. no. 11301 CD) for the reprocessing of suction and irrigation channels in cleaning machines to a reprocessing tube and position on the valve connections. Use the cleaning adapter (Art. no. 11301 GG2, see page 8) for reprocessing in ETD2/ETD3 machines. In addition, use the special connecting adapter (Art. no. 11301 KK2, see page 8) for cleaning in ETD3 machines. For fiberscopes with a working channel without a suction valve, connect the tube of the washer and disinfectant up to the working channel ⑧. To avoid pressure losses in the washer and disinfectant, lock unused connections with the LUER-Lock (Art. no. 6011590).

WARNING: Please observe the instruction manual of the washer and disinfectant.

Lay the fiberscope carefully in the tray of the machine. Connect the rinsing tubes to the channels. Proceed with the endoscope cleaning and disinfection program according to the regulations of the machine manufacturer.

CAUTION: After the program has finished, the fiberscopes (including all adapters except 11301 CD) are highly susceptible to corrosion due to the ambient conditions in the reprocessing machine.

Remove the fiberscopes from the machine as soon as possible after the program has finished. Then dry the optical surfaces with a cotton-tipped applicator soaked in 70 % alcohol.

NOTE: If, according to the high risk classification and/or the Infection Protection Act, sterilization of flexible fiberscopes is not required, the instruments should be stored safely in a hanging position free from dust (distal tip pointing down), so that any remaining moisture in the channels can drain out.

NOTE: For flexible fiberscopes without channels, a hanging storage position is not essential.

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
Connection for machine reprocessing in ETD3 machines

Dismantling the adaptor 11301 GG2: Press the button corresponding to each tube and pull the tube off simultaneously.

1. Slide the free tube end ① of the adaptor (Art. no. 11301 GG2) onto the olive connector of the green adaptor (Art. no. 11301 CD).

ⓘ NOTE: For fiberoptic scopes with a working channel without a suction valve, pull the tube not required, e.g. ① out of the metal block.

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2. For machine reprocessing, attach the green adaptor (Art. no. 11301 CD) to the valve connections of the fiberoptic scope.

3. Attach the tube with LUER-lock to the fiberoptic instrument channel ⑧ and secure it.

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4. Insert the other tube ends in the metal part of the adaptor (Art. no. 11301 GG2).

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5. Connect the adaptor (Art. no. 11301 FF2) to the 'valve for equalizing pressure and connecting leakage tester'.

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6. Connect the leakage tester to the ETD3 using the adaptor 11301 FF2.

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7. Connect the cleaning adaptor for cleaning the instrument channel in EDT3 cleaning machines (Art. no. 11301 GG2) to the connector of the supply tube. Then connect up the special connecting adaptor for ETD-3 machines (Art.-Nr. 11301 KK2).

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8. Lay the instrument carefully in the tray of the ETD3 machine.

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9. Connect the rinsing tubes to the ETD3 machine.

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WARNING: To avoid larger pressure losses in the EDT3 machine, connections which are not in use should be locked with LUER locks.

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
10. Screw the adaptor 11301 KK2 onto the designated connections of the ETD3 machine, according to the figure.

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11. Use the cleaning and disinfection program for fiberoptic scopes according to the ETD3 manufacturer's instructions.

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After cleaning has been completed, remove the green adaptor (11301 CD) for cleaning the working channel by pressing together the projections at the side.

! **NOTE:** The green adaptor for cleaning the working channel may be left on the tube for later use.

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Remove the connection from adapter 11301 FF2.

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
! **CAUTION:** The fiberoscope could be damaged by liquid which may enter the device via an adaptor for leakage testing (11301 FF2) left on the fiberoscope after machine reprocessing. Therefore, you should always dismantle the adaptor (11301 FF2) immediately after cleaning is completed!

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To disassemble adaptor 11301 GG2: Press respective button and pull off tube simultaneously.

! **NOTE:** If, according to the high risk classification and/or the Infection Protection Act, sterilization of flexible fiberoscopes is not required, the instruments should be stored safely in a hanging position free from dust (distal tip pointing down), so that any remaining moisture in the channels can drain out.

! **NOTE:** For flexible fiberoscopes without channels, a hanging storage position is not essential.

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
Connection for machine reprocessing in ETD2 machines

Push free end of the adaptor tube (Art. no. 11301 GG2) onto the barbed tube (Art. no. 11301 CD).

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
Place the adaptor 11301 CD on the valve nozzle ③ of the fiberscope.

Insert the other end of the tube into the metal part of adaptor 11301 GG2.

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Place the LUER-Lock on the working channel ⑧ of the fiberscope.

Insert the other end of the tube into the metal part of adaptor 11301 GG2.

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Insert the supply tube with the fitting for ETD2 cleaning machines into the metal part of adaptor 11301 GG2.

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Attach adaptor 11301 FF2 to the valve ⑤ for pressure compensation and the leakage tester connection of the fiberscope.

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Connect leakage tester connection of the ETD2 cleaning machine.

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Lay the endoscope carefully in the machine basket. Connect the irrigation tubing to the channels. Proceed with the endoscope cleaning program according to the regulations of the machine manufacturer.

CAUTION: After the program run, the fiberscopes are greatly at risk of corrosion in the reprocessing machine.

The fiberscopes (including all adapters except 11301 CD) must be removed from the machine as soon as possible once the program has finished.

Then if necessary dry off the optical surfaces with a cotton swab soaked with 70 % alcohol.

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
After cleaning has been completed, remove the green adaptor (11301 CD) for cleaning the working channel by pressing together the projections at the side.

NOTE: The green adaptor for cleaning the working channel may be left on the tube for later use.

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Remove the leakage tester connection from the adaptor 11301 FF2.

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CAUTION: The fiberscope could be damaged by liquid which may enter the device via an adaptor for leakage testing (11301 FF2) left on the fiberscope after machine reprocessing. Therefore, you should always dismantle the adaptor (11301 FF2) immediately after cleaning is completed!

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18. Inspection after reprocessing

Before assembly it should be checked that all parts and surfaces are undamaged, free from any organic or chemical residue as well as completely clean and dry.

Replace defective parts or send in for repair. Before being sent in for repair, all parts must be cleaned, wipe-down disinfected and sterilized in the NT procedure if necessary.

If residue is still present after reprocessing, reprocess the parts again, subject them to special treatment, discard or repair. Remove stubborn deposits on glass surfaces with a special cleaning paste (Art. no. 27661).

The part must then be submitted to a complete reprocessing cycle.

Check plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their new condition, they must be replaced. This applies in particular to seals and rubber caps.

Test the instruments for proper operation: check the telescope for a clear view, the instrument for steerability and all canals for patency.

WARNING: Damage to the tube: The use of ULTRASTOP as an anti-fog solution on the optically distal end faces of the flexible endoscopes can damage the tube connector.

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19. Gas sterilization/plasma sterilization

Sterilization using formaldehyde or ethylene oxide and plasma sterilization are gentle and therefore recommended methods of sterilization.

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NOTE: Fiberscopes which are to be used for interventions in areas of the body which lack microbial colonies, must be sterilized in sterile packaging and stored in a closed cupboard free from contamination and dust.

CAUTION: If the instrument is to be sterilized, the entire system must be opened by fitting the red cap on the pressure compensation part to ensure pressure compensation.

WARNING: Before commencing gas sterilization, the instrument must be completely dry or else sterile conditions will not be achieved.

- Place the instrument securely on a sterilization sieve.

Ethylene oxide gas sterilization has been validated under the following parameters:

Gas mixture	EO: HCFC*
Temperature (°F, °C)	130 °F ± 5 °F; 54 °C ± 2 °C
Relative air humidity	60 ± 20%
Pressure (p _a , kPa)	58 ... 70 (0.56 ... 0.7 bar/8 ... 10 psi)
Exposure time	120 mins
EO concentration	600 ±30 mg/L

* Gas mixture consisting of 10 % ethylene oxide and 90 % chlorotetrafluoroethane (HCFC-124; percent by weight)

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	at 20 °C	at 42 °C	at 55 °C
Telescopes	24 hrs	18 hrs	12 hrs
Rubber items	7 days	4 days	3 days
Fiberscopes/ Plastics	5 days	2 days	30 hrs

(minimum time in simple transparent packaging and an application time on the patient of more than 30 minutes.)

Maximum acceptable levels (per 21 CFR 812.100; Federal Register, Vol. 43, No. 122, June 23 1978) of residues following EO sterilization are as follows:

Ethylene oxide:	25 ppm* / 250 ppm**
Ethylene chlorohydrine:	25 ppm* / 250 ppm**
Ethylene glycol:	250 ppm* / 5000 ppm**

** for instruments in contact with skin and mucous membrane,
* for instruments in contact with blood.

The instruments can be used immediately after sterilization with formaldehyde or plasma sterilization. Extra airing time is not necessary.
Remember if you use plasma sterilization that there are limits regarding the size of the lumen and material in regard to instruments which may be sterilized using STERRAD®.

Sterilization using STERRAD NX is also possible.

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Statements as to the efficiency of sterilization for the respective type of instrument (depending on the length and diameter of the lumina) may only be made by the manufacturer.

WARNING: Sterility is only assured if the manufacturer's instructions for the STERRAD® procedure are followed.

WARNING: All deviations from the sterilization parameters recommended for the STERRAD® system must be validated by the user.

CAUTION: Attention must be paid to the information on reprocessing, sterilization and storage in the 'Cleaning, sterilization and care of KARL STORZ instruments' instructions. The methods for cleaning, disinfection and sterilization are explained there in detail.

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20. Chemical sterilization

Chemical sterilization (STERIS®) is a special method which is particularly suited for thermolabile goods such as fiberscopes. Except for telescopes with a HAMOU® eyepiece drive mechanism, all KARL STORZ telescopes are suitable for STERIS® sterilization. Color changes to the surface of aluminum alloys could occur during sterilization, which do not however affect the function of the instrument. It must be pointed out in this context that anodized aluminum parts are not always recognizable as such.

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WARNING: Sterility is only guaranteed if the manufacturer's instructions for the STERIS System have been followed. This also applies to use of the sterilization containers.

WARNING: Chemical sterilization is not possible on greased or oiled surfaces.

WARNING: Parts with a lumen must be reprocessed together with the STERIS Quick Connect kit. All parts must be thoroughly cleaned before sterilization using the STERIS PROCESS™, as also applies to the other sterilization methods.

CAUTION: Attention must be paid to the information on cleaning, sterilization and storage in the 'Cleaning, sterilization and care of KARL STORZ instruments' instructions. The methods for cleaning, disinfection and sterilization are explained there in detail.

CAUTION: For miniature endoscopes with a removable eyepiece, the eyepiece must either be sealed with a plug or firmly attached to the endoscope before sterilization.

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21. Storage

CAUTION: Do not expose fiberscopes to direct sunlight or excessive heat.

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Sterilized fiberscopes should be packaged and stored in sterile cloths or suitable containers in accordance with EN 868, and must be absolutely dry before being put into storage.

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Storage must take place in facilities with normal temperature and air humidity levels.

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You can find special holders for this in the 'Cleaning, sterilization and care of KARL STORZ instruments' instructions.

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22. Accessories

Flexible forceps

Carefully wipe the outer surfaces of the flexible forceps with a disposable cloth soaked in disinfectant. Clean the forceps jaws with a brush.

Next, clean in an ultrasound bath. Afterwards, rinse with microbiologically pure/sterile water.

Manual/machine cleaning/disinfection

Both processes can run parallel to decontaminate flexible forceps.

Exception: If the flexible forceps are thermostable, they can be cleaned and disinfected in the 93 °C program.

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For handles which can be dismantled:

Dismantle handle as follows:

- Undo the thumbscrew on the handle and remove the rear part of the handle.
- Undo the thumbscrew on the connecting piece and remove the front part of the handle to the rear.

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Also oil the joints of the jaws by holding the instrument so that the oil can flow into the instrument channel and wet the pulling wire. Then reassemble in reverse order.

The sterilization of flexible forceps is compulsory!

NOTE: Only use special oils whose lubricating effect remains even in autoclaves (art. no. 27656).

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Papillotome/Polypectomy snarn

Dismantle the device completely for reprocessing.

- Unscrew thumb handle.
- Undo contact socket.
- Unscrew grip piece.
- Remove inner guide sleeve.
- Pull papillotomy or polypectomy loop out of the plastic probe at the front.

Then reassemble the polypectomy snare in reverse order.

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Both processes can run parallel to decontaminate papillotomy/polypectomy snares.

Exception: If the papillotomy/polypectomy snares are thermostable, they can be cleaned and disinfected in the 93 °C program.

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Wipe the exterior of the instruments with a cloth soaked in disinfectant. Clean the forceps jaws with a disinfected brush.

Then clean in an ultrasound bath.

Afterwards, rinse with microbiologically pure/sterile water.

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
Assembly of papillotomes:

- Screw inner guide sleeve into the grip piece
- Fit contact socket.
- Insert loop into plastic probe as follows: insert long end of loop into the rear opening of the probe and push to the rear, keeping the probe stretched.
- Insert shorter end of loop into the front opening of the probe and slide both parts towards the rear keeping them parallel.

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In polypectomy snares, the loop end is already shaped. Insert the proximal loop end into the plastic probe from the distal end and push towards rear.

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Pull the end of the loop out of the plastic probe at the rear as far as possible, preferably arranging the probe in a circle. The free end of the loop must be long enough to extend below the locking screw.

Insert the end of the loop into the handle and lock it in place with the locking screw.

Lock the probe connecting piece to the lower end of the handle.

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
Flexible injection cannula

For reprocessing, pull the needle out of the flexible casing.

Rinse the cladding, then immerse it in solution and clean in the ultrasound bath.

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Both processes can run parallel to decontaminate papillotome/polypectomy snares.

Exception: If the papillotome/polypectomy snares are thermostable, they can be cleaned and disinfected in the 93 °C program.

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Rinse out the cannula using a syringe with disinfectant solution.

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24. Normatives

24.1 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original parts.

24.2 Limitation of liability

As the supplier of this instrument, KARL STORZ shall be liable for the safety, operational reliability and performance of this instrument subject to the following conditions:

assembly, system expansion, readjustment, modification, or repair have been performed by persons authorized by KARL STORZ and that the instrument is used in accordance with the instruction manual.

24.3 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see “Subsidiaries” section), even during the warranty period.

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24.4 Directive compliance

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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25. Repair and replacement program

KARL STORZ maintains a repair and replacement warehouse which is normally adequate to ensure prompt replacement of damaged telescopes and instruments. Under the repair and replacement plan, you receive an identical as-new instrument and are only charged the repair costs for the defective instrument. For telescopes, you receive a warranty of 1 year, and for instruments 6 months.

For fiberscopes and equipment, individual repair is necessary. To bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Straße 34
D-78532 Tuttlingen/Germany

In other countries please contact your local KARL STORZ branch or authorized dealer.

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25.1 Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights.

KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

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**1 Important information
for users of KARL STORZ
medical devices**

**It is recommended to check the suitability of
the products for the intended procedure prior
to use.**

Thank you for your expression of confidence in
the KARL STORZ brand name. Like all of our
other products, this product is the result of years
of experience and great care in manufacture.
You and your organization have decided in favor
of a modern, high-quality medical device from
KARL STORZ.

This manual is designed to help with the correct
assembly and operation of the BETTOCCHI[®]
B.I.O.H.[®] compact hysteroscope.

Please read these instructions carefully. Keep this
manual in a convenient and conspicuous location
and in its protective casing close to the device.

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**3 Controls, displays,
connectors, and their uses**

- ① B.I.O.H.® outer sheath 4 mm (26252 BO)
- ② Working channel
- ③ Compact hysteroscope 26252 BC
(without outer sheath and valves)
- ④ Suction and irrigation valve (26252 BW)
- ⑤ Valve connector for connecting to tubing set
031317-10 (mtp) or monobloc (26252 BM)
- ⑥ Light port
- ⑦a Monobloc adapter (26252 BM) for connecting
the B.I.O.H.® compact hysteroscope and a
tubing set
- ⑦b Tubing set 031317-10 (mtp)
- ⑧ Assembly Aid for B.I.O.H.® Valves
(26252 BWM)

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4 Symbols employed

4.1 Safety instructions

Symbol	Explanation
	WARNING: A Warning indicates that the personal safety of the patient or physician may be involved.
	CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.
	NOTE: A Note indicates special information about operating the device or clarifies important information.

4.2 Pictograms

Symbol	Explanation
	Consult instructions for use
	Manufacturer
	Do not reuse
	Catalogue number
	Serial number
	Number of products in the product packaging

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	Keep dry
	Fragile, handle with care
	Unsteril
	Federal (USA) law restricts this device to sale by or on the order of a physician.
	CE mark. The CE mark is a marking required for certain marketable industrial products in accordance with EU law. It is composed of the CE logo and, where applicable, the four-figure identification number of the responsible testing center (when the product belongs to a certain risk category).

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5.1 Intended use

Telescope:

The telescope is used for endoscopic viewing and for creating a working or irrigation channel during gynecology. Telescopes are designed for short-term use in endoscopic interventions or surgical invasive interventions.

Sheath:

The sheath is used to create a working or irrigation channel during gynecology. Sheaths are designed for short-term use in endoscopic interventions or surgical invasive interventions.

5.2 Indications

The use of hysteroscopes is indicated if, in the view of the attending physician, hysteroscopy is indicated.

This includes, amongst other things, the diagnosis and treatment of:

- AUB (abnormal uterine bleeding)
- Metrorrhagia/Menorrhagia
- Sterility
- Endometrial hyperplasia
- Polyps
- Myomas
- Intrauterine adhesions (e.g., Asherman's syndrome)
- Uterine malformations (e.g., septum)

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5.3 Contraindications

The use of hysteroscopes is contraindicated if, in the opinion of the attending physician, the surgical method as such is contraindicated, or if the patient is not able to undergo surgery or anesthesia due to his or her general condition.

Hysteroscopes must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

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- Pregnancies
- Acute inflammatory situations (cervicitis, colpitis)
- For operative hysteroscopy: unclear histology, exclusion of malignant tumor
- In the case of endometrial ablation: myoma, polyp, desire to have children
- For synechiae: freedom from symptoms without the desire for future fertility

Relative contraindications:

- Severe heart failure due to potential volume overload
- Second half of the menstrual cycle due to not being able to reliably exclude the possibility of pregnancy and the thick endometrium

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6 Safety

6.1 Explanation of warnings and cautions

Please read this manual carefully and follow its instructions precisely. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the system.

WARNING: Risk of injury and damage to products: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals carefully and always follow the instructions given precisely. Check the functioning of the products used in combination.

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NOTE: The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

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WARNING: Risk of injury and risk of damaging the products: If applied parts are used outside of the field of vision there is a risk that tissue and accessories could be damaged unintentionally. Always hold the applied parts of the active electrode, laser and other instruments which transmit energy pointed at the target in the field of vision during use.

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- Avoid direct tissue contact with the distal end, as well as with the light port of the telescope and the light cable.
- Do not rest the endoscope or the light cable on the patient or where it is in direct contact with operating accessories.
- Always select the lowest possible light output of the cold light source which still allows optimal illumination of the operating field.

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WARNING: Risk of injury: Patient leakage currents may accumulate if endoscopes and endotherapy devices supplied with energy are used together. This is particularly important if an endoscope with an applied part of type CF is used. In this case, an endotherapy device with a type CF applied part must be used in order to minimize the total patient leakage current.

WARNING: Before use, check the expiration date and ensure that the packaging is undamaged. If the expiration date has passed or the packaging is damaged, the device must not be used.

WARNING: If CJD (Creutzfeldt Jakob disease) is suspected or indicated, KARL STORZ recommends disposing of the device in accordance with national regulations and not using it again.

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7 Initial operation

7.1 Unpacking

Carefully remove the BETTOCCHI[®] B.I.O.H.[®] compact hysteroscope and the accessories from the packaging.

Check for missing items and for evidence of shipping damage.

If necessary: notify the manufacturer or the supplier.

7.2 Scope of supply

	26252 BC	BETTOCCHI [®] B.I.O.H. [®] compact hysteroscope
	26252 BO	Outer sheath
	2x 26252 BW	Suction and irrigation valve
	1x 26252 BWM	Assembly Aid for B.I.O.H. [®] Valves
	26252 BM	Monobloc adapter
	100010-10	<i>Endoscopic Seal</i> , for working channels, pack of 10

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8 Disassembly

1. If necessary, remove instruments from the working channel ②.
 2. Unscrew the light cable and light cable adapter.
 3. Unscrew the valves out of the handle for cleaning (in non-engaged position).
- ⓘ NOTE:** When doing so, note the 2 different O-ring sizes for valve 26252 BW. There is one O-ring for the sheath.
4. Replace damaged sealing rings.
 5. Remove the outer sheath ① by holding it at the swivel nut and pulling it out forward in a straight line in a controlled manner.
 6. Check the O-ring at the sheath mount and replace if necessary.
 7. If a single-use endoscopic seal (article no. 100010-10) is included in the delivery, remove it from the working channel ② and dispose of it correctly.

The following replacement seals can be found in the spare parts box, article no. 26252 DS:

10 x 7966595 O-ring (for valve 26252 BW)
10 x 5918900 O-ring (for valve 26252 BW)
10 x 5922000 O-ring (for telescope sheath mount)

8. Press the locking mechanism on the handle and remove the monobloc adapter ⑦a or the disposable tubing set ⑦b from the valve connection ⑤.
9. Check the O-rings on the monobloc adapter ⑦a and replace if necessary.

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9 Reprocessing

9.1 Overview table

Article number			Precleaning			Cleaning/disinfection			Sterilization													
			manual						machine			STERRAD®										
Disassembly	Manual wipe-down disinfection	Immersion in cold water	Dissolving incrustations	Brushing the surfaces	Brushing the lumina	Rinsing out with a pressure water gun	Ultrasound treatment	Manual cleaning	Manual disinfection	Chemical disinfection	Thermal disinfection	Connecting	Steam sterilization	STERRAD® 50, 100S, 200	STERRAD® NX®	STERRAD® 100NX®	STERRAD® 100NX® DUO Cycle	AMSCO® V-PRO™1	Ethylene oxide (ETO)	Low-temperature steam and formaldehyde process (LTSF)	STERIS® System 1®	STERIS® System 1E®
26252 BW			•			•					•		•									
Compact hysteroscope 26252 BC	•		•		◦	•					•	•	•									
26252 BO											•	•	•									
26252 BM											•	•	•									
26252 BWM					◦						•		•									

- This reprocessing step is part of the validated process.
- This reprocessing step may be carried out additionally in order to aid cleaning.
The sterilization method is approved in respect of material compatibility, however must be validated by the operator on site in respect of effectiveness.
- x On the grounds of potential damage, this reprocessing step must not be performed.

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CAUTION: When reprocessing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

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NOTE: The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

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- Brush: 27652
- Adapter: 26252 BWA, 26252 BM
- Tray: 39501 XC

9.4 Disassembly

Prior to cleaning and disinfection, the medical device must be separated into its individual components as far as possible (see Section 8 'Disassembly' on page 13).

9.5 Preparation for cleaning and disinfection

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9.6 Manual precleaning

9.6.1 Immersion in cold water

To dissolve dried organic residues, the medical device must be immersed in cold water (10-20°C) for at least 5 minutes.

9.6.2 Brushing the surfaces

Depending on the degree of soiling and the degree to which residues have dried, heavy soiling must be cleaned from the surfaces of instruments with the aid of a brush/sponge under cold running water until any visible contamination has been completely removed.

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9.6.3 Rinsing out with a pressure water gun

The working channel, lumina and hollow spaces must be rinsed using a pressure water gun and a suitable rinsing attachment. The valve 26252 BW must be rinsed through the cleaning hole (see arrow). For rinsing, use pulses of at least 3 x 5 seconds from a jet of water with a pressure of at least 2.5 bar. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

9.7 Machine cleaning and disinfection

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ⓘ NOTE: If necessary, the instrument must be dried off afterwards by hand.

9.8 Connecting

In order to ensure effective machine cleaning and disinfection, the full tray (art. no. 39501 XC) must be connected up to the washer/disinfector.

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9.9 Filling the cleaning tray

1. Screw the valve seat adapter (art. no. 26252 BWA) into the valve holes on the handle.
2. Attach the tubes on the wire tray onto the cleaned monobloc (art. no. 26252 BM) and the LUER connector on the working channel.
3. Insert the BETTOCCHI[®] B.I.O.H.[®] compact hysteroscope into the tray.
4. Undo the outer sheath and place it on the compact hysteroscope in the correct position in the tray.
5. Position the suction and irrigation valves into the designated holders in the tray.
6. Place the equipped tray into the washer and disinfectant and connect the LUER connectors on the front side of the tray.

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- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Threads and glide surfaces must be lubricated locally with instrument oil.
- Dismantled medical devices must be assembled (see section 10 'Assembly' on page 23).
- Afterwards, a functional check must be carried out.

i NOTE: The oil used for this purpose must be suitable for the subsequent sterilization procedure (silicone-free and paraffin- or white oil-based).

i NOTE: Use items from the catalog 'Care, Sterilization and Storage Techniques' during care procedures.

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9.12 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'.

The method must be selected taking into account the respective applicable national requirements and in consultation with the device manufacturers. The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

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The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) at 132 °C – 137 °C with a minimum exposure time of 4 minutes and a maximum of 18 minutes. This procedure is only suitable for thermostable instruments.

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10 Assembly

1. Screw the light cable adapter onto the light inlet piece ⑥.
2. Click on the outer sheath ①.
3. Screw the valves into the handle.
4. If necessary, attach the single-use endoscopic seal onto the proximal end of the working channel ②.
5. Check that the instrument is functioning correctly.

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11 Servicing

11.1 Disposal

For disposal, no special measures are necessary. Disposal must be carried out in accordance with the respective applicable local and national laws and regulations.

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11.3 Application and storage conditions

	Temperature	Rel. humidity
Storage	-20°C ... +60°C	10 % ... 90 %
Application	+10°C ... +40°C	30 % ... 70 %
Air pressure	700 hPa – 1060 hPa	

11.4 Directive compliance

This medical device bears a CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

A code number after the CE mark indicates the responsible notified body.

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com
KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za
TOO KARL STORZ Endoscopy Kasachstan
6, Saryarka str. BG „Arman”, office 314
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz
KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai – United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com
KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in
KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, MPlaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com
KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au
KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk
KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongquiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn
KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn
KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn
KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn
KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg
KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg
KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjail-Ro, Bundang-Gu
Seongnam, 13556, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com
KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw
KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 31745-00
E-Mail: philippines@karlstorz.com
KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp
www.karlstorz.com

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**1 Important information
for users of
KARL STORZ devices**

It is recommended to check the suitability of the products for the intended procedure prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to aid you with proper connection and operation of the VITOM® 3D in conjunction with the instruction manual for the IMAGE1 S. All the necessary details and actions are clearly explained. We thus ask that you read these manuals carefully before proceeding to work with the equipment. Insert in the protective wallet and keep available for ready reference in a convenient and conspicuous location near the equipment.

Please refer to the IMAGE1 S instruction manual for the individual functions of the CCU that can be used with the VITOM® 3D exoscope.

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
**2.1.1 Controls, displays,
connectors and their uses**

- ① Rotary knob for horizontal image adjustment
- ② Sheath support for connection to holding system
- ③ Rocker switch – Scroll up/down in menu/programmed functions
- ④ Control button: Open menu, selection/programmed functions
- ⑤ Handle
- ⑥ Connection for fiberoptic light cable
- ⑦ Stereo lens
- ⑧ Light outlet – integrated illumination unit
- ⑨ Video connection cable

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2.2 IMAGE1 PILOT

2.2.1 Controls, displays, connectors and their uses

- ① Control buttons 1 – 4 for programmed functions
- ② LED for operating status display
- ③ 3D wheel for control of zoom, focus and ROI as well as navigation through the menu
- ④ Right control button: Cancel, jump back from a menu and Esc
- ⑤ Left control button: Open the camera menu, Setup menu and selection
- ⑥ KSLOCK quick-release coupling (male) for connection to a KARL STORZ holding system

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

3 Symbols employed

3.1 Symbols on the product

	Follow instructions for use
	Type CF device/ Defib. resistant
	CE mark

3.2 Symbols on label and packaging


	Consult instructions for use
	Catalogue number
	Serial number
	Number of products in the product packaging
	Fragile, handle with care
	Keep dry
	Temperature limit
	Federal (USA) law restricts this device to sale by or on the order of a physician.

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1	Important information for users of KARL STORZ devices	III
2	Images of the device	IV
2.1	VITOM® 3D exoscope	IV
2.1.1	Controls, displays, connectors and their uses	V
2.2	IMAGE1 PILOT.....	VI
2.2.1	Controls, displays, connectors and their uses	VI
3	Symbols employed.....	VII
3.1	Symbols on the product.....	VII
3.2	Symbols on label and packaging	VII
4	General information	4
4.1	Description of the device.....	4
5	Safety instructions	5
5.1	Warnings and cautions.....	5
5.2	Applications	9
5.2.1	Intended use.....	9
5.2.2	Indication	9
5.2.3	Contraindications.....	9
5.3	Patient profile	9
5.4	User profile.....	9
5.5	Intended conditions of use	10
5.5.1	Use.....	10
5.5.2	Other intended conditions.....	10
5.6	Safety precautions when operating the instrument	11
6	Installation and operating instructions	12
6.1	Unpacking the equipment	12
6.2	Components of the VITOM® 3D system	12
6.2.1	Necessary basic equipment.....	12
6.2.2	Optional accessories	12
6.3	Compatible camera control units (CCU).....	13
6.4	Before each use.....	13

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8.1.1	Care and handling	33
8.2	General warnings	35
8.3	Manual wipe-down disinfection	36
8.4	Assembly, inspection and care.....	37
8.5	Reprocessing limits	37
9	Reprocessing for US customers only.....	38
9.1	General care instructions.....	38
9.1.1	Care and handling	38
9.2	General warnings	41
9.3	Manual wipe-down disinfection	42
9.4	Inspection after cleaning.....	42
9.5	Reprocessing limits	42
10	Servicing	43
10.1	Maintenance and safety check.....	43
10.1.1	Maintenance	43
10.1.2	Safety check.....	43
10.2	Servicing and repair.....	44
10.3	Disposal	44
10.4	Repair program	45
10.5	Limitation of liability	46
10.6	Warranty.....	46
11	Technical description.....	47
11.1	Technical data VITOM® 3D.....	47
11.2	IMAGE1 PILOT technical data	48
11.2.1	Standard compliance	49
11.2.2	Directive compliance (for VITOM® 3D)	49
12	Information on electromagnetic compatibility (EMC).....	50
13	Subsidiaries	63

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4 General information

4.1 Description of the device

The VITOM® 3D exoscope is a high-resolution stereoscopic visualization system with integrated zoom, focusing and image alignment mechanism for visualization of microsurgical interventions. The applications are comparable to those of a surgical microscope. With the VITOM® 3D exoscope, the displayed site is enlarged on a 3D screen.

The VITOM® 3D exoscope is compatible with the IMAGE1 S via the D3-LINK™.

The VITOM® 3D exoscope is positioned for use with the aid of a holding arm in such a way that the user can visualize the site.

The VITOM® 3D exoscope is controlled using the control unit supplied (IMAGE1 PILOT). The control unit is secured to the operating table using a second holding arm in such a way that the user can reach and operate it with ease.

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5 Safety instructions

5.1 Warnings and cautions

Please read this manual carefully and follow its instructions precisely. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in the instruction manual, please read the text that follows carefully, in order to ensure safe and efficient operation of the device.

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CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

NOTE: Notes contain additional information on the safe operation of the device.

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 **WARNING:** Risk of injury and damage to the products: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals carefully and always follow the instructions given precisely. Check the functioning of the products used in combination.

 **WARNING:** Risk of injury: Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all of the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

 **WARNING:** Risk of injury: Incorrect application of medical instruments poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

 **WARNING:** Risk of injury and risk of damaging the products: If applied parts are used outside of the field of vision there is a risk that tissue and accessories could be damaged unintentionally. Always hold the applied parts of the active electrode, laser and other instruments which transmit energy pointed at the target in the field of vision during use.

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WARNING: Risk of injury and risk of damaging the products: The highly concentrated light energy can cause the distal end, the light port, adjacent components and tissue in front of the light emission window to heat up. This can cause burns to the patient, user and operating accessories.

- Avoid direct tissue contact with the distal end, as well as with the light port of the telescope and the light cable.
- Never allow the exoscope and the fiberoptic light cable to rest on the patient or come in direct contact with surgical accessories.
- Always select the lowest possible light output of the cold light source which still allows optimal illumination of the operating field.

WARNING: Risk of injury: Optical radiation poses a risk of injury to eyes. Never look into the free end of either a connected fiberoptic light cable or an exoscope. Always wear standardized protective equipment when using lasers.

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WARNING: Risk of injury: Patient leakage currents may add if powered endotherapy devices are used simultaneously. This is particularly important if an endoscope with a Type CF applied part is used. In that case, an endotherapy device with a Type CF applied part must be used to minimize the total patient leakage current.

WARNING: The video system may stop working during use. An adequate replacement system must be kept available.

WARNING: Risk of electric shock: Do not open the device! Have servicing carried out by authorized personnel only. The removal of covers by unauthorized staff voids the warranty.

WARNING: Risk of fatal injury: In the case of patients with cardiac pacemakers it must be ensured that the proximal end (= furthest from the patient) of the telescope is not positioned less than 2 cm from the pacemaker. The device generates a magnetic field at the setting ring which can cause cardiac pacemakers or other electrical devices to malfunction.

WARNING: Before use, check the expiration date and ensure that the packaging is undamaged. If the expiration date has passed or the packaging is damaged, the product must not be used.

CAUTION: Protect the device from water splashes and other liquids.

CAUTION: Do not store liquids on or above the device.

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5.2 Applications

5.2.1 Intended use

The VITOM® 3D is used for the short-term, non-invasive visualization of the operating field during microsurgical, open and endoscopic interventions.

5.2.2 Indication

Extracorporeal visualization of the operating field during microsurgical, open and endoscopic interventions.

5.2.3 Contraindications

Use of the VITOM® 3D exoscope is contraindicated, if, in the opinion of the attending physician, the VITOM® 3D system is not suitable for visualization of the operating field, or the general condition of the patient does not permit the use of the VITOM® 3D system.

As KARL STORZ light sources are not intended for use in ophthalmology, the use of the VITOM® 3D system is also contraindicated in this field.

5.3 Patient profile

The use of the VITOM® 3D is not limited to a certain patient profile (sex, age, weight etc.). The attending physician must assess a patient's health for the relevant application and declare it suitable.

5.4 User profile

- Adequate and recognized skills in the relevant application
- Adequate powers of comprehension to rationally assess the surgical situation in hand
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- Have been trained thoroughly at least once in the use of the device

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- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic) and that cannot be compensated using auxiliary means.

5.5 Intended conditions of use

5.5.1 Use

The device is intended for use in hospitals and doctors' offices. The technical data and ambient conditions are described in the instruction manual. It is the user's responsibility to make sure the device is safe and operating properly before using the device. During treatment with the VITOM® 3D exoscope the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital signs and the anesthetic.

i NOTE: Any treatment may only be performed if visual observation of the action of the device is ensured.

5.5.2 Other intended conditions

Frequency of use:	one or more times a day
Length of use:	from a few minutes to several hours a day
Place of installation:	Positioning on a suitable holding system
Mobility:	can be positioned and moved in combination with a holding arm
Combination:	can be used on the patient at the same time as other equipment required for the operation
Control:	can be controlled via the IMAGE1 PILOT

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**5.6 Safety precautions when
operating the instrument**

It is the user's responsibility to make sure the instrument is safe and operates properly before using it.

During treatment in conjunction with the VITOM® 3D exoscope, the patient must be treated and kept under observation using the usual medical care. This includes monitoring sterile application conditions where required by the type of intervention.

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6.2.2 Optionales Zubehör

- 1 VITOM® 3D-Illuminator (TH 002)

6 Installation and operating instructions

6.1 Unpacking the equipment

Carefully remove the VITOM® 3D exoscope from its packaging. Check for missing items and evidence of damage during shipment. Please file any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these can come in handy if the device has to be transported.

6.2 Components of the VITOM® 3D system

6.2.1 Necessary basic equipment

Camera system and monitor:

- 1 IMAGE1 S CONNECT module (TC 200 EN)
- 1 IMAGE1 S D3-LINK™ module (TC 302)
- 1 32" 3D monitor (TM 323)
- 1 3D polarization glasses (9800 GF)

VITOM® 3D exoscope:

- 1 VITOM® 3D exoscope (TH 200)
- 1 IMAGE1 PILOT (TC 014)
- 1 Sterile cover for VITOM® 3D exoscope (TH 001)
- 1 Sterile cover for IMAGE1 PILOT (041150-20)

Recommended light source:

- 1 Power LED 300 (TL 300)
- 1 Light cable with a diameter of 4.8 mm (e.g., 495 TIP or 495 VIT)

Holding arms:

- 1 Suitable KARL STORZ holding system for VITOM® 3D exoscope (e.g., 28272 HS with 28272 HP)
- 1 VITOM® 3D clamping jaw with ball joint for KARL STORZ holding systems (28272 VTK)
- 1 Suitable KARL STORZ holding system for the IMAGE1 PILOT (e.g., 28272 HA with 28172 HR)

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6.3 Compatible camera control units (CCU)

CCU

IMAGE1 S CONNECT (TC 200)
from SW version **2.6** in combination
with IMAGE1 S D3-LINK™ (TC 302)
from SW version **2.6**

Please check the software version as follows:

1. Open the Setup menu (Setup key on the dashboard or Shift+Enter key on the keyboard).
2. Select "Module Information".
3. Check the software version.

If the above-mentioned or a higher software version (TC 200 – 2.6, TC 302 – 2.6) is not installed on your IMAGE1 S camera control unit, you can install the software version using the USB stick optionally provided or contact your KARL STORZ service partner.
If it is necessary to upgrade the hardware of the compatible Link module (TC 302), the system will show you this. For more information on this, please contact your KARL STORZ service partner.

6.4 Before each use

Before each use, the instrument must be inspected for damage, especially around the distal end (cover glass, light outlet, ...).
Before each use, the video connecting cable and light cable must be inspected for damage, especially breaks and cuts.
The use of instruments with damaged cable must be discontinued. Cables must be replaced (light cable to be replaced by the user or the Service Dept., and the video connection cable by the Service Dept.).

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6.5 General information on commissioning of the VITOM® 3D exoscope

The VITOM® 3D exoscope is always held over the site by a holding system during a practical application. The holding arm should be positioned in such a way that the operating theater procedures are not interrupted and the user has an optimal and ergonomic view of the 3D monitor. If possible, the monitor should therefore be positioned directly opposite the user at eye level.

- i NOTE:** The VITOM® 3D exoscope including the equipment connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations and standards.
- i NOTE:** Always observe the instructions for the holding system used and its suitability for use with the VITOM® 3D exoscope.

6.6 Attaching the VITOM® 3D exoscope to the holding arm

For operation of the VITOM® 3D exoscope, a KARL STORZ holding arm (e.g., endoscopic holding system) which can carry a load of at least 2 kg is required.

Depending on the holding system, the VITOM® 3D exoscope is attached to the KARL STORZ holding system using a clamping jaw (28272 VTK or 28272 VTP).

Attach the clamping jaw to the holding system correctly and ensure that the clamping jaw is firmly connected to the holding system (see Fig. A: Clamping jaw correctly attached to the holding system and Fig. B: Clamping jaw incorrectly attached to the holding system).

- i NOTE:** Observe the relevant instruction manuals for use of the clamping jaw 28272 VTK or 28272 VTP on the holding system (Art. no.: 97000154 or 97000155).

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Then open the clamping jaw 28272 VTK or 28272 VTP by simultaneously pressing the unlock buttons. Place the VITOM® in the clamping jaw and pay particular attention to correct orientation (see Fig. C).

NOTE: Due to the circumferential groove, the VITOM® 3D can only be placed in the clamping jaw in one direction.

Hold the VITOM® 3D in position and close the clamping jaw with your other hand until it clicks into place audibly (see Fig. D). Check whether the VITOM® 3D is safely held in the clamping jaw before you let it go. The clamping jaw ensures that the VITOM® 3D can rotate along its axis. This permits a comfortable position during the intervention (see Fig. E).

Depending on the holding system and clamping jaw, the VITOM® 3D can be swiveled in all directions either via the ball joint on the clamping jaw (Art. no.: 28272 VTK) or via the integrated ball joint in the POINT SETTER.

NOTE: Before every intervention, check whether the VITOM® 3D exoscope can be positioned using the holding system selected in such a way that optimal visualization of the site is ensured.

6.7 Connecting the light cable to the VITOM® 3D exoscope

The light cable is mounted by fitting to the light port provided and turning clockwise. Check that it is firmly connected (see Fig. F). Light cables with a fiber diameter of 4.8 – 5.0 mm are recommended for the VITOM® 3D.

CAUTION: Please note that the circular symbol on the VITOM® 3D and the circular symbol on the light cable must correspond (see following table).

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Overview: Recommended combinations*

Light cable diameter	
4.8 – 5.0 mm	⊙
3.0 – 3.5 mm	⊙
2.0 – 2.5 mm	⊙
Endoscope/exoscope diameter	
6.6 – 12.0 mm	⊙
3.0 – 6.5 mm	⊙
0.8 – 2.9 mm	⊙

*Special endoscopes/exoscopes may deviate slightly

**6.8 Fitting the VITOM® 3D and
holding arm with a sterile cover**

The VITOM® 3D exoscope cannot be sterilized and must therefore be fitted with a sterile cover for the duration of the intervention. The cover TH 001 has been specially designed for the VITOM® 3D exoscope and features a plastic frame with a transparent screen in the distal area.

Remove the sterile cover from its packaging (see Fig. G). Insert the plastic frame over the cover glass and light outlet on the distal end of the VITOM® 3D exoscope (see Fig. H). Then pull the cover over the VITOM® 3D exoscope and holding arm (see Fig. I).

NOTE: The cover TH 001 has been specially designed for the VITOM® 3D. The use of other sterile covers is prohibited.

NOTE: Observe the instruction manual with article number 97000169 for handling instructions for the sterile cover.

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6.9 Attaching the IMAGE1 PILOT

The VITOM® 3D exoscope system is always operated with the IMAGE1 PILOT (TC 014). It is attached to the operating table using a KARL STORZ articulated stand (e.g., base 28172 HR and articulated stand 28272 HA) and positioned in such a way that the user can reach the control unit with ease.

The holding arm is attached to the operating table and fixed in position as described in the corresponding instruction manual with Art. no. 96116335DF.

The IMAGE1 PILOT is equipped with a KSLOCK quick-release coupling (male) on its underside and can thus be easily attached to KARL STORZ articulated stands (see Fig. J).

To this end, open the KSLOCK quick-release coupling (female) on the articulated stand by turning in the 'OPEN' direction and insert the IMAGE1 PILOT so that it clicks into place audibly. Make sure that the IMAGE1 PILOT is correctly locked in (see Fig. L: fully locked in/ Fig. K: not fully locked in).

Then lock the KSLOCK quick-release coupling (female) on the articulated stand by turning in the 'CLOSE' direction (see Fig. M).

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Then connect the USB cable of the control unit to a USB socket on the IMAGE1 S CONNECT module (see Fig. N). The message 'IMAGE1 PILOT connected' appears on the screen. The LED for display of the function status must be lit up constantly.

i NOTE: The IMAGE1 PILOT calibrates the 3D Wheel position on power up. Therefore, do not press any keys or touch the 3D Wheel while attaching the USB connection, or powering up the IMAGE1 S system.

6.10 Fitting the IMAGE1 PILOT with a sterile cover

The IMAGE1 PILOT cannot be sterilized and is therefore fitted with a sterile cover during the intervention (041150-01).

Remove the sterile cover from its packaging and use it to cover the control unit and articulated stand (see Fig. O).

Make sure that the elastic insert on the tip of the sterile cover is placed directly over the IMAGE1 PILOT. The elastic insert increases the grip when operating the IMAGE1 PILOT (see Fig. P).

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6.11 Connecting the VITOM® 3D to a light source

Plug the light cable into the connecting socket on the light source as far as it will go. To disconnect the light cable pull on the connector (see Fig. Q), never on the cable itself.

6.12 Connecting the VITOM® 3D to the D3-LINK™

Insert the video connecting cable connector into the connection socket on the D3-LINK™.

NOTE: Do not use moist cables.

NOTE: The VITOM® 3D exoscope is first initialized after being connected. It takes approx. 20 seconds until an image is displayed.

6.13 Switching on the CCU, light source and monitor

NOTE: Heed the instruction manuals for the IMAGE1 S, the light source, the light cable and the monitor.

6.14 Optimal light efficiency

NOTE: The entrance and exit faces of the light cable and VITOM® 3D must be kept clean if optimal light transmission is to be maintained. The light cable must first be unscrewed (see Fig. R). Then clean the entrance and exit faces with a cotton swab soaked in cleaning solution and then wipe with a cotton swab moistened with 70% isopropyl alcohol. Please note that both the glass surface inside the light port of the VITOM® 3D (see arrow, Fig. S) and on the light cable itself (see Fig. T) must be cleaned. Please observe the information on cleaning the light cable connection on page 32.

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6.15 Disassembly of the VITOM® 3D system after use

Switch off the equipment after using the VITOM® 3D exoscope. Remove and dispose of the sterile covers of the VITOM® 3D exoscope and IMAGE1 PILOT.

Remove the VITOM® 3D exoscope and IMAGE1 PILOT from the holding system, perform a wipe-down disinfection (also see chapter 8 'Reprocessing') and keep both safe.

NOTE: The handle on the VITOM® 3D exoscope can become very hot during operation, especially on the sides and on the light cable connection (over 50°C under certain circumstances). This is not a malfunction. It is recommended to allow the system to briefly cool down after use prior to disassembly of the holding system. To this end, please remove the camera head cable from the D3-LINK™ and switch the light source to standby mode.

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
7 Operation of the VITOM® 3D exoscope system

7.1 Setting the horizon

There is a rotary knob on the distal end of the VITOM® 3D with which the stereo lens can be rotated through +/- 90° (see Fig. U). This enables flexible positioning of the 3D exoscope and correction of the horizontal image alignment during the application. To set the horizon, turn the black rotary knob in the desired direction (see Fig. V). If the image has to be rotated further, the VITOM® 3D image can be rotated digitally through 180° using the 'Image Flip' function. This makes it possible to rotate an image through 360°. If necessary, the 'Image Flip' function can be assigned to one of the 4 programmable buttons on the IMAGE1 PILOT or to the camera head buttons. Also see chapter 7.4 'Programmable function buttons on the IMAGE1 PILOT'.

NOTE: To avoid damage to the VITOM® 3D exoscope, do not turn the knob beyond the limit stops.

7.2 Controlling the optical functions of the VITOM® 3D via the IMAGE1 PILOT

Most of the functions of the VITOM® 3D exoscope can be remote controlled easily via the IMAGE1 PILOT control unit. To do so, you do not have to hold the VITOM® 3D exoscope.

The IMAGE1 PILOT can be most easily operated using your finger tips (see Fig. W).

The most important functions are explained below.

NOTE: For reasons of simplicity, a sterile cover is not displayed in this image. Always use a sterile cover for practical application in a sterile environment.

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7.2.1 Adjusting the magnification

The VITOM® 3D has an infinitely variable zoom which can be conveniently adjusted via the IMAGE1 PILOT.

To magnify the image, press the 3D Wheel of the IMAGE1 PILOT straight down (see Fig. X).

To reduce the image, pull on the 3D Wheel of the IMAGE1 PILOT.

The more you press the 3D Wheel, the faster you will be able to zoom in or out. This ensures that the desired magnification can be set precisely.

7.2.2 Adjusting the focus

The working distance of the VITOM® 3D can be variably adjusted from 20 – 50 cm. The 3D Wheel enables focusing of the image at any time.

To bring a lower level (greater working distance) into focus, turn the 3D Wheel clockwise.

To bring a higher level (shorter working distance) into focus, turn the 3D Wheel counterclockwise (see Fig. Y).

The more you turn the 3D Wheel, the faster you will be able to focus. This ensures that the required focus point can be set with great precision.

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7.2.3 Selectable region of interest during magnification (ROI)

Whenever a magnification greater than 1x is selected, the selected region of interest can be moved upward, downward, to the left or to the right. This facilitates centering of the magnified image without having to move the holding arm.

To move the region of interest, push the 3D Wheel in the required direction (see Fig. Z and a).

The more you press the 3D Wheel, the faster you will be able to move the region of interest.

i NOTE: The region of interest can only be moved in the image level. This does not change the angle of view of the VITOM® 3D relative to the object. The region of interest can only be moved with a magnification setting of >1x. The movement range is limited. The image cannot be moved beyond the margins of the fully zoomed out image.

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
By pressing the 3D Wheel on the IMAGE1 PILOT, an overview screen (minimap) is displayed in the top right-hand corner which provides the user with visual feedback on the zoom area and the position of the selectable region of interest (ROI) (see Fig. b).

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The large rectangle shows the available overview area. As soon as you use the zoom, a small filled rectangle which represents the selectable region of interest appears. The size of the rectangle varies with the selected zoom level. If the region of interest is moved, the rectangle also moves and shows the position within the available overview area. Under the overview area, a bar shows the relative zoom level. The bar is filled as the magnification factor increases.

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7.3 Control of the IMAGE1 S using the IMAGE1 PILOT

The IMAGE1 PILOT permits direct access to the menu of the IMAGE1 S so that the surgeon can operate all the functions from the sterile field as necessary (see Fig. c).

As soon as the menu of the IMAGE1 S is called up using the 'ENTER' button, you can navigate through the menu, adjust settings or call up functions by moving the 3D Wheel. You can exit the menu by pressing the ESC button.

The following table offers an overview of the functions.

i NOTE: As long as the IMAGE1 S menu is active, the VITOM® 3D functions (e.g., zoom and focus) cannot be controlled.

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7.4 Programmable function buttons on the IMAGE1 PILOT

The 4 buttons on the IMAGE1 PILOT can be programmed in a similar way as the camera head buttons and can be assigned functions (see Fig. d and e). In doing so, certain functions can be performed quickly or rapid access to menus can be provided.

The programmable buttons can be assigned the following functions:

- White Balance
- Light source intensity
- Capture
- Video
- S technologies
- Brightness
- Nothing
- Source swap
- AEL Activate
- Image Flip

7.5 Advanced Exposure Levelling (AEL) function

The function **Advanced Exposure Levelling (AEL)** stabilizes the Exposure control in situations with extreme illumination conditions. and prevents the VITOM® 3D from under- or overexposing the image, for example when a shiny instrument is inserted.

The effect is shown schematically in the following images:

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Well exposed image through automatic exposure control.

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By inserting the shiny instrument, the automatic exposure control dims and the operating field is displayed too dark.

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The exposure is stabilized when the AEL function is activated. The reflection on the instrument may become larger, however the operating field is displayed well exposed.

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The AEL function can be assigned to one of the programmable buttons on the IMAGE1 PILOT. As soon as the button is pressed, the AEL function is active. The AEL logo appears on the screen. The exposure does not change even when inserting instruments.

The AEL function remains active until one of the following events takes place:

- AEL button is pressed again
- The region of interest (ROI) is considerably changed using the IMAGE1 PILOT
- The zoom is considerably changed using the IMAGE1 PILOT

After deactivating the AEL function, the VITOM® 3D returns to the normal exposure mode.

7.6 Extended PIP-Swap function

The PIP-Swap function can be assigned to one of the programmable buttons on the IMAGE1 PILOT. As soon as more than one Camera LINK module is connected to the IMAGE1 Platform, the PIP-Swap function permits switching between the cameras defined as Primary and Secondary Source. For example, this may be beneficial during endoscopically assisted microsurgery (EAM) to quickly switch from the VITOM® 3D to an endoscopic camera.

Pressing the function button to which PIP-Swap is assigned switches the image between the primary source and secondary source.

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7.7 Operation of the VITOM® 3D exoscope without IMAGE1 PILOT

The VITOM® 3D exoscope can also be operated without the IMAGE1 PILOT. This is necessary if no IMAGE1 PILOT is connected or if the IMAGE1 PILOT fails.

If one of these events occurs, the VITOM® 3D exoscope can alternatively be operated via the touchpad.

NOTE: For safety reasons, always connect the keyboard to the IMAGE1 S.

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
	Zoom
	By moving your finger upward on the touchpad, the image view is magnified.
	By performing a downward movement on the touchpad, the image view is reduced.
	Focus
	Move to the left on the touchpad to focus in the distance.

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Move to the right on the touchpad to focus at close range.

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
Moving the ROI

By holding the left mouse button depressed and moving your finger on the touchpad, the magnified region of interest can be moved in the x and y direction.

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8 Reprocessing not validated for US customers

! WARNING: The light cable must be removed for reprocessing of the VITOM® 3D. To do so, unscrew the light cable in a counterclockwise direction. The video connecting cable is firmly connected to the VITOM® 3D and must not be removed for reprocessing.

i NOTE: Instructions for US customers can be found in Chapter 9.

8.1 General care instructions

i NOTE: Always disconnect the power cord from the Camera Control Unit before cleaning

We recommend storing the VITOM® 3D in its original packaging when not in use. Do not coil the video connecting cable of the VITOM® 3D into tight loops or cause it to kink while storing. Do not expose the VITOM® 3D unprotected to direct sunlight or excessive heat.

When returning the VITOM® 3D for service, the user must ensure the equipment being returned is clean and (wipe-down) disinfected. The entrance and exit faces of the light cable and exoscope must always be kept clean (see corresponding instruction manual) to ensure optimal reproduction of the image.

! CAUTION: Chemical residue on the card-edge connector can cause failure of the camera electronics.

! CAUTION: Heat from the light source can cause chemical residue on the light cable and on the light entrance of the exoscope to burn. This can damage the light fibers and have a detrimental effect on the lighting power and the picture quality.

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8.1.1 Care and handling

To maximize the life of the VITOM® 3D, please follow the guidelines given below for proper handling and care.

DO NOT PLUG THE VIDEO CONNECTING CABLE INTO THE PROCESSOR RECEPTACLE IF IT IS WET. The connector must always be completely dry and free of residue. Where necessary, remove any moisture with a low lint disposable cloth before plugging in the device. Nor should the fiberoptic light cable be inserted into the receptacle of the light source if it is wet. The connector must always be free of residue and any remaining moisture. If necessary, remove remaining moisture with a low lint disposable cloth. The connectors of the video connecting cable and of the light cable must also be completely dry before connecting!

NEVER PICK UP THE VITOM® 3D BY ITS CABLE! Always pick up the set by the instrument handle first, followed by the cable and connector. The cable contains tiny electronic conductors and should be treated gently and carefully.

PROTECT THE VITOM® 3D FROM DIRECT SUNLIGHT OR EXCESSIVE HEAT.

HANDLE THE SYSTEM WITH CARE, avoiding extreme impacts to the device. When transporting, use the original packaging to prevent damage to the device.

DO NOT STORE THE VITOM® 3D WITH CABLE COILED TIGHTLY OR GATHERED INTO FOLDS. Store the cable in loops at least 15 cm in size.

DO NOT RAPIDLY PULL THE CABLE WHEN UNCOILING. Harmful 'kinks' can result. Always allow the cable to uncoil gently and loosely.

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- This reprocessing step is part of the validated process.
- This reprocessing step may be carried out additionally in order to aid cleaning. The sterilization method is approved in respect of material compatibility, however must be validated by the operator on site in respect of effectiveness.
- ✗ On the grounds of potential damage, this reprocessing step must not be performed.

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8.2 General warnings

 **WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

 **NOTE:** Risk of infection: Always use the medical device with the sterile cover (TH 001), which is available to order separately.

 **WARNING:** Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

 **WARNING:** When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

 **CAUTION:** When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

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Reprocessing (worldwide, does not apply to USA)

CAUTION: Danger of damage to medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices. Only use chemicals approved by KARL STORZ for reprocessing. A complete list can be found on the Internet at www.karlstorz.com.

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8.3 Manual wipe-down disinfection

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications concerning material compatibility and necessary exposure time for the disinfectant must be observed. At the end of the necessary exposure time, wipe the surface with a dry low lint cloth. Clean and wipe down the cover slips of the VITOM® 3D exoscope for the telescope and lighting, making sure that there are no stripes. Ensure that the lenses are not scratched.

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9 Reprocessing for US customers only

WARNING: The light cable must be removed for reprocessing of the VITOM® 3D. To do so, unscrew the light cable in a counterclockwise direction. The video connecting cable is firmly connected to the VITOM® 3D and must not be removed for reprocessing.

9.1 General care instructions

NOTE: Always disconnect the power cord from the Camera Control Unit before cleaning

We recommend storing the VITOM® 3D in its original packaging when not in use. Do not coil the video connecting cable of the VITOM® 3D exoscope into tight loops or cause it to kink.

Do not expose the VITOM® 3D unprotected to direct sunlight or excessive heat.

Prior to returning for service, the VITOM® 3D must be thoroughly cleaned per validated procedures.

The entrance and exit faces of the light cable and exoscope must always be kept clean (see corresponding instruction manual) to ensure optimal reproduction of the image.

CAUTION: Chemical residue on the card-edge connector can cause failure of the camera electronics.

CAUTION: Heat from the light source can cause chemical residue on the light cable and on the light entrance of the exoscope to burn. This can damage the light fibers and have a detrimental effect on the lighting power and the picture quality.

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NOTE: Risk of infection: Always use the medical device with the sterile cover (TH 001), which is available to order separately.

WARNING: Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

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CAUTION: National laws and regulations must be observed.

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9.3 Manual wipe-down disinfection

After each use, the device must be wiped down with a low level disinfection wipe with 0.1 – 1% Sodium Hypochlorite Wipe. In addition, clean the distal lens with a cotton tip applicator soaked in lens cleaner or 70% isopropyl alcohol.

9.4 Inspection after cleaning

The cleaned medical device must be visually inspected to verify that it is clean, complete, free of damage, and dry.

- If residues or contamination are still present, the medical device must be wipe-down disinfected again.

9.5 The device must be free from any visible signs of damage (cuts, dents, or scratches, etc.) that may trap residual debris. If the device shows any damage, contact KARL STORZ for repairs.

Reprocessing limits

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

CAUTION: Inspect the video connecting cable and the fiberoptic light cable for breaks or cracks before and after disinfection, as well as prior to subsequent use. Breaks or cracks will allow fluids to enter the cables and cause damage. Where a VITOM® 3D has a damaged fiberoptic light cable or damaged video connecting cable, the relevant cable must be replaced.

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10.1 Maintenance and safety check

10.1.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

10.1.2 Safety check

Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, the following checks must be carried out on this device at least once a year by persons who, due to their training, knowledge and practical experience, are able to properly carry out such metrological checks and who are not subject to any instructions as far as such checking activities are concerned.

Visual inspection:

- Check the device and accessories for any mechanical damage which may impair functionality.
- Check that the labels relevant to safety are legible.

Electric measurements:

- Protective grounding measured according to IEC 62353 (current edition), and meets the limits specified in the standard.
- Touch Current, Earth Leakage Current, and Patient Leakage currents, measured according to IEC 62353 (current edition), and meets the limits specified in the standard.

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Test for proper functioning:

- A functional check must be carried out in accordance with the instruction manual.
- The metrological check must be entered in the unit log-book and the results of the check must be documented.
- If the device does not function reliably and/or safely, it must be repaired.

10.2 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

CAUTION: Do not open the instrument or the lens. Opening these sealed assemblies affects the camera's tightness and voids all warranties. The covers on the camera processor must only be removed by authorized personnel.

10.3 Disposal

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the device as electronic scrap.

Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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10.4 Repair program

Usually to bridge the repair period, you will receive a device on loan which you then return to KARL STORZ as soon as you receive the repaired device.

In the case of repairs, please contact your local KARL STORZ branch or authorized dealer.

In Germany, you can refer repairs direct to:
KARL STORZ SE & Co. KG
Abt. Reparaturservice
Take-off-Gewerbepark 83
78579 Neuhausen

Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among personnel). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, modifications, or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights. KARL STORZ gives no warranty on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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10.5 Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the device is connected and operated meet the applicable laws and standards; and
- the device has been used in accordance with its operating instructions at all times.

10.6 Warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see 'Subsidiaries' section), even during the warranty period.

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11 Technical description

11.1 Technical data VITOM® 3D

Dimensions (LxWxH):	approx. 351 x 45 x 80 mm
Weight:	approx. 950 g
Working distance:	20 – 50 cm
Buttons/Input options	
• 3 camera head buttons	
• Rotary knob for mechanical image alignment:	+/- 90°
Interfaces	
• Video interface for the IMAGE1 S D3-LINK™ module	
• Interface for fiberoptic light cable	
• Interface for attaching a holding system using the clamping jaw 28272 VTK or 28272 VTP	
Operating temperature:	5°C ... 30°C
Storage	
Ambient temperature:	-10°C ... +60°C
Humidity:	5 ... 95% (non-condensing)
Air pressure:	500 ... 1060 hPa

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11.2 IMAGE1 PILOT technical data

Dimensions (LxWxH):	approx. 139 x 93 x 77 mm
Weight:	approx. 354.6 g
Buttons/Input options	
<ul style="list-style-type: none"> • 1 3D Wheel with 6 DoF (Degrees of Freedom) • 4 programmable buttons • 2 buttons for menu navigation 	
Interfaces	
<ul style="list-style-type: none"> • USB interface for operation with the IMAGE1 S CONNECT • KSLOCK for fast and secure fixation on KARL STORZ holding systems 	
Operating temperature:	5°C ... 30°C
Storage	
Ambient temperature:	-10°C ... +60°C
Humidity:	5 ... 85% (non-condensing)
Air pressure:	500 ... 1060 hPa

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11.2.1 Standard compliance

for VITOM® 3D exoscope (TH 200) in connection with the KARL STORZ IMAGE1 S CCU

According to IEC 60601-1, UL 60601-1, CAN/CSA C22.2 No. 601.1-M90:

- Degree of protection against electric shock: Defibrillation-protected Type CF 

The Type CF defib protection is realised in the IMAGE1 S™ D3-Link module (TH 302). The VITOM® 3D is always operated with this module.

11.2.2 Directive compliance (for VITOM® 3D)

In accordance with the Medical Device Directive (MDD):

This medical device belongs to Class I

This medical device bears the CE mark in accordance with MDD 93/42/EEC.

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12 Information on electromagnetic compatibility (EMC)

WARNING: Medical electrical equipment needs special precautions regarding electromagnetic compatibility (EMC). Observe the EMC instructions in this chapter during installation and operation.

The VITOM® 3D exoscope model TH 200 corresponds to IEC 60601-1-2 [CISPR 11 Class A] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC. The employed limit values provide a basic measure of safety from typical electromagnetic interference that is to be expected in a medical environment. The VITOM® 3D exoscope model TH 200 is a Group 1 device (according to CISPR 11). Group 1 contains all of the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this chapter provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If electromagnetic interference arises while using the device, the user can eliminate it by taking the following actions:

- reorient or relocate the receiving device;
- increase the separation between the equipment;
- connect the devices to different electrical circuits.

In case of additional questions, please contact your area representative or our Service Department.

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Information on electromagnetic compatibility (EMC)

WARNING: The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

WARNING: The use of accessories and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the device. The accessories and cables listed below have been shown to comply with the requirements of IEC 60601-1-2. When using accessories and cables other than those specified here, it is the responsibility of the user to ensure that they comply with IEC 60601-1-2.

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Table 1 Guidance and manufacturer's declaration – electromagnetic emissions		
The VITOM® 3D exoscope model TH 200 is intended for use in the environment specified below. The customer or user of the VITOM® 3D exoscope model TH 200 should make sure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions acc. to CISPR 11	Group 1	The VITOM® 3D exoscope model TH 200 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions acc. to CISPR 11	Class A	The VITOM® 3D exoscope model TH 200 is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions acc. to IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions acc. to IEC 61000-3-3	Complies	

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Table 2 Guidance and manufacturer's declaration – electromagnetic immunity			
The VITOM® 3D exoscope model TH200 is intended for use in the electromagnetic environment specified below. The customer or user of the VITOM® 3D exoscope model TH200 should make sure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst acc. to IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines acc. to IEC 61000-4-11	$<5\% U_T^*$ (>95% dip in U_T) for 0.5 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $<5\% U_T$ (>95% dip in U_T) for 5 sec	$<5\% U_T^*$ (>95% dip in U_T) for 0.5 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $<5\% U_T$ (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field acc. to IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
* Note: U_T is the a.c. mains voltage prior to application of the test level.			

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Table 4 Guidance and manufacturer's declaration – electromagnetic immunity – for equipment and systems that are not life-supporting			
The VITOM® 3D exoscope model TH200 is intended for use in the electromagnetic environment specified below. The customer or user of the device should make sure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF disturbance acc. to IEC 61000-4-6	3 V _{eff} 150 kHz to 80 MHz	3 V _{eff}	Portable and mobile RF communications equipment should be used no closer to any part of the VITOM® 3D exoscope model TH200, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5 / 3] \sqrt{P}$
Radiated RF disturbance acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = [3.5 / 3] \sqrt{P}$ 80 MHz to 800 MHz $d = [7 / 3] \sqrt{P}$ 800 MHz to 2.5 GHz
			Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.			
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

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Table 6 Recommended separation distances between portable and mobile RF communications equipment and the VITOM® 3D exoscope model TH 200			
<p>The VITOM® 3D exoscope model TH 200 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.</p>			
Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters [m] can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer.</p> <p>Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
6, Saryarka str. BC „Arman“, office 314
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, M Plaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongqiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 610041, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjail-Ro, Bundang-Gu
Seongnam, 13556, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 31745-00
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp
www.karlstorz.com

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1	Important information1
2	Pictures of the instruments2
3	Controls, displays, connectors and their uses.....4
4	Symbols employed5
5	Safety instructions.....6
5.1	Intended use6
5.2	Indications6
5.3	Contraindications.....7
6	Patient population7
6.1	User qualification7
7	Warnings and cautions8
7.1	Safety precautions when using flexible intubation videoscopes8
8	Handling14
8.1	Unpacking14
8.2	Basic equipment for endoscopes with a working channel14
8.3	Accessories (not included in the set).....15
8.4	Basic equipment for endoscopes without a working channel15
8.4.1	Accessories15
8.5	Operating conditions16
8.6	Electromagnetic Compatibility16
8.7	Handling16
9	Operation16
9.1	Holding the device correctly17
9.2	Fitting the suction valve (only necessary with intubation videoscopes with a working channel)17
9.3	Fitting the tube holder.....17
9.4	Handling the sheath.....17
9.5	Handling the deflection mechanism18
9.6	Handling the working channel (does not apply to 11301 ABX)20

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11.8.1 Hydrogen peroxide (H ₂ O ₂) sterilization - ASP STERRAD®	44
11.8.2 Hydrogen peroxide (H ₂ O ₂) sterilization - STERIS® V-PRO™	44
11.8.3 Ethylene oxide sterilization (EO)	45
11.8.4 Chemical low-temperature sterilization with peracetic acid – STERIS System 1E®	45
11.9 Limits of reprocessing.....	46
11.10 Storage.....	46
11.11 Reprocessing of spare parts.....	47
11.11.1 General warnings	48
11.11.2 Accessories.....	49
11.11.3 Manual wipe-down disinfection	49
11.11.4 Manual precleaning	50
11.11.5 Disassembly	50
11.11.6 Machine cleaning and disinfection	50
11.11.7 Connecting.....	51
11.11.8 Assembly, inspection and care	51
11.11.9 Packaging systems	51
11.11.10 Sterilization	51
11.11.11 Limits of reprocessing.....	53
12 Technical data.....	54
12.1 Flexible intubation videoscope 11301 BN/11303 BN:	54
12.2 Flexible intubation videoscope 11302 BD:	55
12.3 Flexible intubation videoscope 11301 AB:	56
13 Servicing and repair	57
14 Maintenance.....	57
15 Limitation of liability	57
16 Warranty	57
17 Directive compliance.....	58
18 Important information	58
19 Disposal.....	58

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1 Important information

It is recommended that the suitability of the products for the intended procedure be checked prior to use.

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This instruction manual is intended to serve as an aid in the proper use, cleaning, disinfection and, if necessary, sterilization of the intubation video-scope. All the necessary details and actions are clearly explained. We therefore ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for easy reference.

WARNING:

Risk of infection:

KARL STORZ instruments are not sterile when delivered, and must be cleaned, disinfected and, if necessary, sterilized prior to initial use and each subsequent use.

KARL STORZ is continuously working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, figures and descriptions in this manual.

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**3 Controls, displays,
connectors and their uses**

- ① Connecting cable
- ② Connection for suction valve/irrigation adaptor
11301 CB1/11301 CE1/11301 CD1 (only for
endoscopes with working channel)
- ③ Connection for leakage tester or pressure com-
pensation cap (pressure compensation valve)
- ④ Working channel with lip valve 29100 (only for
endoscopes with working channel)
- ⑤ Sheath
- ⑥ Handle
- ⑦ Deflection lever
- ⑧ Distal tip
- ⑨ Leakage tester 13242 XL
- ⑩ Pressure compensation cap 11025 E
- ⑪ Monitor 8403 ZX
- ⑫ Second connection for videoscope
- ⑬ Primary connection for videoscope
- ⑭ Mains connection
- ⑮ Reset button
- ⑯ Holder for SD memory card
- ⑰ Status bar
- ⑱ Pilot lamp
- ⑲ On/Off switch button
- ⑳ Charging pilot lamp (orange)
- ㉑ Softkeys
- ㉒ Function symbols
- ㉓ USB port
- ㉔ HDMI port
- ㉕ Video recording Start/Stop
- ㉖ White balance
- ㉗ Single image capture
- ㉘ Tube holder
- ㉙ Suction valve

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Control of pollution caused by electronic information products (China RoHS)

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UL approval mark, valid in the USA and Canada

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Degree of protection in accordance with International Protection Codes: Protection against contact with tools or similar and small foreign bodies, protection against vertically dripping water

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Rated ambient temperature/transducer temperature class

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5 Safety instructions

Please read these safety instructions carefully before using the instrument. Before using the instrument on the patient for the first time, it is imperative that you be acquainted with how the instrument operates and is controlled.

5.1 Intended use

The flexible intubation videoscopes 11301 ABX, 11302 BDX, 11301 BNX and 11303 BNX are used in conjunction with the C-MAC® monitor 8403 ZX for endotracheal/nasal intubation. Intubation endoscopes are designed for transient use in invasive procedures through a body orifice. Intubation endoscopes are used to inspect the upper and lower airways, to check the tube position with double lumen tubes and for monitoring during PCT.

NOTE: The flexible intubation videoscope 11301 BNX can also be used with the C-MAC® monitor 8402 ZX.

Flexible intubation videoscopes are applied parts, which must only be used with KARL STORZ supply units which are intended for this purpose.

5.2 Indications

The use of intubation endoscopes is indicated if, in the opinion of the attending physician, a fiberoptic endotracheal intubation or inspection of the oropharyngeal, supra- or subglottic region is indicated or PCT is performed.

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- Endotracheal intubation with a difficult respiratory tract
- Endoscopic inspection and monitoring of the mouth, nose and throat region.
- Monitoring of percutaneous tracheostomy
- Anesthesiological bronchoscopy

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5.3 Contraindications

The use of the intubation endoscope is contraindicated if, in the opinion of the attending physician, the surgical method as such is contraindicated or if the patient is not able to undergo surgery or anesthesia due to his or her general condition.

Intubation endoscopes must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

No other contraindications are known.

6 Patient population

Men and women as well as children and adolescents. Use is not limited to a certain patient profile (sex, age, weight, etc.).

NOTE: The attending physician decides which patients to use the product on.

6.1 User qualification

KARL STORZ flexible intubation videoscopes may only be used by persons with an appropriate medical qualification and who are acquainted with the flexible laryngoscopy technique. The information given in these instructions only serves to instruct in the correct handling and reprocessing of the flexible intubation videoscope. It is not intended as an introduction to the techniques of flexible laryngoscopy.

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7 Warnings and cautions

Please read this manual carefully and follow the instructions exactly. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this instrument. To further emphasize the words **Warning**, **Caution** and **Note**, they are preceded by a pictogram.

WARNING:

'Warning' indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

ATTENTION:

'Attention' indicates that particular service procedures or safety precautions must be followed to avoid any damage to the product.

NOTE:

'Notes' indicate special information about operating the instrument, or clarify important information.

7.1 Safety precautions when using flexible intubation videoscopes

Flexible intubation videoscopes must only be used according to the medical rules and procedures of endoscopy recognized for this method.

NOTE: The flexible intubation videoscopes do not offer any electric insulation against high frequency (HF) currents. The use of HF instrument sets can result in image faults and damage to the endoscope.

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⚠ ATTENTION:

Damage to product: Whenever the flexible intubation videoscope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition.

- Damaged flexible intubation videoscopes or damaged accessories must not be used.

⚠ WARNING:

Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended use and interface specifications of the devices used in combination permit this.

The instruction manuals and interface specifications of the medical devices used in combination must be observed precisely.

⚠ WARNING:

Risk of infection:

- The flexible intubation videoscopes are not sterile when delivered, and must be cleaned and disinfected/sterilized prior to initial use and each subsequent reuse.

⚠ WARNING:

The pressure compensation cap ⑩ must be removed from the pressure compensation valve ③ during:

- endoscopic interventions
- immersion in liquid, e.g., disinfectant solution and for machine reprocessing

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Safety instructions

⚠ ATTENTION:

Damage to product: The pressure compensation cap ⑩ must be placed on the pressure compensation valve ③ during:

- gas sterilization (EO, FO, hydrogen peroxide)
- transport and shipping

⚠ WARNING:

Risk of burns and glare: The LED light emitted through the flexible intubation videoscope may cause glare or high temperatures, and thereby burns, at the distal light outlet.

- Never look directly into the light emitted by the flexible intubation videoscope.

⚠ WARNING:

Risk of explosion: If flammable anesthetics are used in the immediate proximity of the flexible intubation videoscope.

- The flexible intubation videoscope must not be used in the presence of flammable or combustible media.

⚠ ATTENTION:

Only use accessories that are specifically suited for use with videoscopes in accordance with the intended use and which can be employed in line with generally accepted technical standards and safety guidelines.

⚠ ATTENTION:

Damage to product: The flexible intubation videoscope may only be repaired by personnel authorized by KARL STORZ. Unauthorized repairs may impair functioning and safety, and are therefore strictly forbidden. KARL STORZ gives no warranty for flexible intubation videoscopes that are repaired by unauthorized personnel.

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NOTE: Follow the 'Cleaning, Sterilization and Care of KARL STORZ Instruments' instructions. The cleaning, disinfection and sterilization procedures are explained therein in detail.

NOTE: Damage to the flexible intubation videoscope and injury to the patient or user resulting from incorrect operation are not covered by the manufacturer's warranty.

⚠️ WARNING:

Always unplug the flexible intubation videoscope before carrying out any maintenance and cleaning procedures on it.

⚠️ WARNING:

Risk of injury: The flexible intubation videoscope must not be used during the discharge of a cardiac defibrillator.

- Remove the flexible intubation videoscope from the operation area before discharge takes place.

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- Switch off the flexible intubation videoscope before putting it down, and do not place it in the vicinity of the patient (as per IEC 60601-1-1 or chapter 16 of the 3rd edition of IEC 60601-1).

⚠️ WARNING:

Risk of infection: In order to prevent infection and protect personnel, it is strictly forbidden to ship contaminated medical devices.

- To avoid contact and aerogen infections, medical devices must be decontaminated on-site and wrapped in protective film.

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⚠️ WARNING:

Danger of electric shock:

- Have servicing on these devices carried out by authorized personnel only. Any opening of the device by unqualified persons will result in the forfeiting of all warranty claims. Exceptions are those that are described in the instruction manual.

⚠️ WARNING:

Only KARL STORZ devices may be used in combination with the flexible intubation videoscope. Disregarding this warning can endanger the patient and user.

⚠️ WARNING:

Risk of injury: If the flexible intubation videoscope malfunctions during use on the patient, discontinue use immediately.

- Place the distal tip of the flexible intubation videoscope in the neutral, non-deflected position and remove the flexible intubation videoscope slowly and carefully from the patient.

⚠️ ATTENTION:

Damage to product:

- Only use accessories (guide probes, etc.) from KARL STORZ. Accessories from other manufacturers can damage the flexible intubation videoscope.

⚠️ ATTENTION:

Damage to product:

- Flexible intubation videoscopes must never be steam sterilized (autoclaved).

⚠️ ATTENTION:

Damage to product:

- The maximum permissible temperature for decontamination, sterilization, transport and storage is 65 °C.

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⚠ ATTENTION:

Damage to product:

- Do not allow any live components or instruments to come into contact with the flexible intubation videoscope.

⚠ ATTENTION:

Damage to product:

- Also do not place the flexible intubation videoscope on other devices that become very hot.
- Do not expose the device to direct sunlight or excessive heat, to avoid damage to the device.

⚠ ATTENTION:

Damage to product:

- Flexible intubation videoscopes must never be exposed to X-rays, to avoid damage.

⚠ ATTENTION:

Damage to product:

- Do not hit the distal tip of the flexible intubation videoscope against hard objects.
- Do not remove soiling with sharp-edged objects under any circumstances.

⚠ ATTENTION:

Damage to product:

- Never lay heavy objects on the flexible intubation videoscope.

⚠ ATTENTION:

Damage to product:

- Never clamp or fold the sheath of the flexible intubation videoscope. This can result in damage to the videoscope.

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WARNING:

Risk of infection:

- The flexible intubation videoscope is not sterile when delivered, and must therefore be properly cleaned, disinfected and (where applicable) sterilized prior to initial use and each subsequent reuse.

Check for missing items and evidence of shipping damage.

File any complaints immediately with KARL STORZ or the supplier.

8.2 Basic equipment for endoscopes with a working channel

The 11302 BDXX, 11301 BNXX and 11303 BNXX sets consist of:

- 1x Flexible intubation videoscope
- 11302 BDXX 4.0 x 65 cm or
- 11301 BNXX 5.5 x 65 cm or
- 11303 BNXX 5.5 x 65 cm
- 1x Case 27677 SM
- 1x Instruction manual
- 1x Leakage tester 13242 XL
- 1x Pressure compensation cap 11025 E
- 1x Suction valve insert for single use, pack of 20 11301 CE1
- 1x Lip valve for instrument channel, pack of 10 29100
- 2x Adaptors for cleaning the working channel in cleaning machines (black) 11301 CD1
- 1x Cleaning brush
- 1x LIPP tube holder (intubation fiberscope only) 11301 CFX
- 1x Bronchoscope insertion tube, size 4 10309
- 1x Bronchoscope insertion tube, size 2 10310
- 1x Protective film 13990 SFN
- 1x Protective cap 8401 YZ

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8.3 Accessories (not included in the set)

- 1x Suction valve insert, reusable 11301 CB1
- 1x Adaptor for leakage test in washers and disinfectors (green) 13242 XA
- 1x 11008 C 'Mainz Adaptor' mask attachment, for endoscopic inspection with simultaneous respiration, sterile packaging, pack of 5

8.4 Basic equipment for endoscopes without a working channel

The 11301 ABXK set consists of:

- 1x 11301 ABXK flexible intubation videoscope 3.0 x 51.5 cm
- 1x Case 27677 SL
- 1x Instruction manual
- 1x Leakage tester 13242 XL
- 1x Pressure compensation cap 11025 E
- 1x LIPP tube holder (intubation fiberscope only) 11301 CFX
- 1x Bronchoscope insertion tube, size 4 10309
- 1x Bronchoscope insertion tube, size 2 10310
- 1x Protective film 13990 SFN
- 1x Protective cap 8401 YZ

8.4.1 Accessories

- 1x Adaptor for leakage test in washers and disinfectors (green) 13242 XA

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Ambient temperature: 15 °C – 40 °C
Humidity: 15 % – 80 %

8.6 Electromagnetic Compatibility

This device has been tested and is found to comply with the EMC limit values as per EN/IEC 60601-1-2 CISPR 11 Class B.

8.7 Handling

The flexible intubation videoscope must be used according to the medical rules and procedures recognized for laryngoscopy.

NOTE: Any damage to the flexible intubation videoscope resulting from incorrect operation is not covered by the manufacturer's warranty.

WARNING:

Risk of infection: The flexible intubation videoscope is not sterile when delivered, and must therefore be properly reprocessed prior to initial use and each subsequent reuse.

ATTENTION:

Damage to product:

- In order to avoid damage to the exterior, do not store or disinfect the flexible intubation videoscope together with sharp objects.

9 Operation

The following instructions outline the correct operation of the flexible intubation videoscope and the corresponding accessories. This is not an introduction to the techniques of laryngoscopy. Instructions on endoscopic technique can be found in the relevant medical literature.

The flexible intubation videoscope is used for endotracheal intubation and inspection of the oropharynx.

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9.1 Holding the device correctly

Take the handle of the flexible intubation videoendoscope in your right or left hand and wrap your fingers around it to secure, as shown in the Figure. The thumb is located on the deflection lever ⑦. With this kind of control grip, the flexible intubation videoendoscope can be operated with one hand. (Fig. 1 Intubation videoendoscope with working channel; Fig. 2 Intubation videoendoscope without working channel).

9.2 Fitting the suction valve (only necessary with intubation videoendoscopes with a working channel)

Place the suction valve ② on the connection and rotate until the side nubs on the insert are positioned above the recesses on the connection.

Engage the suction valve ② in the connection by applying gentle pressure.

9.3 Fitting the tube holder

Push the tube holder ② 11301 CFX from distal onto the cone of the flexible intubation videoendoscope until the tube holder is securely in place.

9.4 Handling the sheath

⚠ ATTENTION:

Damage to product:

- The flexible intubation videoendoscope must not be carried by the sheath. Do not pull, clamp or twist the sheath of the flexible intubation videoendoscope.

⚠ ATTENTION:

Damage to product:

- The sheath of the flexible intubation videoendoscope must never be clamped or coiled up too tightly. It is best to store the flexible intubation videoendoscope in a suspended position or to keep it sterilized in a tray specifically designed for sterilizing and storage (39406AS).

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9.5 Handling the deflection mechanism

The deflection lever ⑦ can be used to angle the distal tip of the flexible intubation videoendoscope. To do so, the deflection lever needs to be moved upwards (tip angled downwards) or downwards (tip angled upwards) using the thumb. When the distal tip is moved back into the neutral position, there will be a click. This ensures that the user is aware that the distal tip is back in the neutral position.

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⚠️ WARNING:

Risk of injury: When inserting the flexible intubation videoendoscope into the patient, the distal tip of the flexible intubation videoendoscope must be in the straight, non-deflected position.

⚠️ WARNING:

Risk of injury:

- The deflection lever must never be operated with force or abruptly to avoid injuring the patient or damaging the device!

⚠️ WARNING:

Risk of injury:

- Before removing the flexible intubation videoendoscope from the patient, the distal tip must be put into the neutral position, as otherwise injuries to the patient or damage to the instrument may occur.

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⚠️ WARNUNG:

Risk of injury:

- If the distal tip can no longer be properly controlled (sticking etc.) via the deflection lever, or if the deflection mechanism develops any other fault, the flexible intubation videoscope must not be used.

⚠️ WARNUNG:

Risk of injury:

- If the flexible intubation videoscope malfunctions during use on the patient, discontinue use immediately.
- Place the distal tip of the flexible intubation videoscope in the neutral, non-deflected position and remove the flexible intubation videoscope slowly and carefully from the patient.

⚠️ ATTENTION:

Damage to product: The resistance that can be felt when deflecting the distal tip using the deflection lever increases as the angle of deflection increases.

- Actuate the deflection gently and carefully, without any sudden or jerky movements, as far as the limit stop.

⚠️ ATTENTION:

Damage to product:

- The lens at the distal tip of the flexible intubation videoscope may become detached, e.g., due to mechanical forces. For this reason, it is essential to check that the lens is undamaged and present before and after every decontamination.

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9.6 Handling the working channel (does not apply to 11301 ABX)

Do not force instruments or cleaning brushes through the channels. This can lead to damage to the sheath, especially if the tip of the sheath is deflected.

NOTE: Check instruments, cleaning brushes and other accessories for compatibility before inserting these in the working channel.

ATTENTION:

Damage to product: Unsuitable or damaged instruments may damage the working channel or cause the endoscope to malfunction.

- Check instruments, cleaning brushes and other accessories for sharp edges and corners before inserting them in the working channel.

When connecting and operating suction/irrigation devices, pay attention to the instruction manuals for the particular devices. The irrigation devices recommended by KARL STORZ are set to a maximum pressure of 54 kPa (400 mmHg/0.54 bar/7.8 psi).

- If no instrument is inserted, seal the proximal end of the working channel with a plug (Art. no. 29100).
- If an instrument is to be inserted, attach the lip valve (Art. no. 11301 CA) for single use to seal the proximal end of the instrument channel.
- The working channel of the flexible intubation videoscope means that specially developed instruments can be introduced. Make sure that the accessories are of the appropriate length and diameter. Unsuitable accessories can damage fiberscopes and flexible intubation videoscopes. If you have any questions concerning accessories, please contact KARL STORZ directly or an authorized specialist dealer.

NOTE: When using the lip valve, the instrument is subjected to increased resistance.

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- The flexible intubation videoscope must be in a straight, undeflected position when inserting accessories into the working channel.

⚠ ATTENTION:

Damage to product: The resistance that can be felt when inserting the instruments increases as the angle of deflection increases.

- *If the instruments are inserted with force, this can damage the flexible intubation videoscope or the instrument channel.*

Forceps, scissors and other accessories whose jaw can be opened or closed must always be inserted into the instrument channel with closed jaws. The jaws may only be opened when the distal tip of the instrument is completely clear of the channel. Opening the jaws in the instrument channel can cause damage to the flexible endoscope.

⚠ ATTENTION:

Damage to product: Accessories from other manufacturers may damage the flexible endoscope.

- *Only use accessories (guide probes, etc.) from KARL STORZ.*

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10 Initial operation

The flexible intubation videoscopes (FIVE) can be operated with the following monitors:

- C-MAC® monitor 8403 ZX (11301BNX also with C-MAC® monitor 8402 ZX)
- also with additional monitors via the camera control unit C-HUB® II 20 2903 01

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10.1 Connection to the C-MAC® monitor

The flexible intubation videoscope can be used together with the KARL STORZ C-MAC® monitor 8403 ZX. The instruction manual for the monitor must be observed (see manual 96076008DER).

Older monitor models (8402 ZX/8401 ZX) are not compatible.

WARNING:

Risk of injury: It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces or voltage-carrying elements of other devices at any time.

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⚠️ WARNING:

Risk of injury: Optical radiation poses a risk of injury to eyes.

- *Never look directly into the light output at the distal tip.*

Plug the connecting cable ① of the flexible intubation videoscope into one of the connection sockets ⑫ or ⑬ on the monitor.

The connecting cable ① can be connected or disconnected during operation. If necessary, connect the cord of the power supply up to the connection socket ⑭ on the C-MAC® monitor and plug it into the power supply. The power cord can be fixed to the serial port using the union nut.

The battery is charged automatically as soon as the C-MAC® monitor is connected to the power supply. This is indicated by the lighting up of the orange LED ⑳.

NOTE: The flexible intubation videoscopes 11301 ABX, 11302 BDX and 11303 BNX can only be connected to the C-MAC® monitor 8403 ZX.

NOTE: The flexible intubation videoscope 11301 BNX can be connected to both the C-MAC® monitor 8402 ZX and the C-MAC® monitor 8403 ZX.

10.2 Switching the flexible intubation videoscope On/Off

Switching On/Off: Pressing the on/off button ⑱ on the C-MAC® monitor switches on the flexible intubation videoscope with the integrated LED light source. Pressing the button a second time switches the unit off.

⚠️ WARNING:

Danger of glare:

- *Do not look directly into the light beam.*

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 WARNING:

Risk of injury: When in operation, the integrated LED light source heats up the flexible intubation videoscope.

- *Switch off the flexible intubation videoscope before putting it down, and do not place it in the vicinity of the patient (as per IEC 60601-1-1).*

NOTE: The C-MAC® monitor can be switched on or off by pressing the On/Off button ⑨. The flexible intubation video endoscope itself does not have an On/Off button.

Pressing and holding the On/Off button switches off the monitor together with the connected endoscope. If the On/Off button is actuated only briefly during operation, the monitor will switch to standby mode. This is indicated by the flashing LED on the monitor.

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10.3 Functions of the C-MAC® monitor 8403 ZX

A more detailed description of the operation and functions of the monitor 8403 ZX can be found in the corresponding instruction manual supplied with every monitor (see manual 96076008DER).

If you lose the instruction manual, you can also download it from our homepage at www.karlstorz.com.

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10.4 Connection to the C-HUB® II

The C-HUB® II (20290320) is a plug & play device for processing and transmitting image and control signals to external devices in medical video endoscopy. The flexible intubation videoscope can be used in combination with the C-HUB® II and an external monitor.

To do so, insert the connector ① of the flexible intubation videoscope into the serial port of the C-HUB® II.

A more detailed description of the operation and functions of the C-HUB® II can be found in the corresponding manual (manual 96206529DER).

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NOTE: The flexible intubation videoscopes 11301 ABX, 11302 BDX and 11303 BNX can only be connected to the C-HUB® II 20290301. Older models are not compatible. The flexible videoscope 11301 BNX can be used with both the C-HUB® (20290101) and the C-HUB® II (20290320).

NOTE: Observe the instruction manual for the C-HUB® II (Art. no. 96206234DER).

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1. Connect the power supply unit to the C-HUB® II.
2. Connect the video connecting cable from the flexible intubation videoscope to the video input on the front of the C-HUB® II.
3. Connect the monitor to the S-Video or HDMI output on the rear of the C-HUB® II.

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WARNING:

Risk of injury: The flexible intubation videoscopes from KARL STORZ are not intended for use with light sources.

- *The use of non-compatible units can lead to patient injury or device damage.*

Focus

The image focus is adjusted automatically. Manual adjustment is not necessary.

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
10.5 Function buttons of the endoscope

The endoscope is equipped with three buttons. These are used to activate the functions displayed as symbols on the endoscope.

The following functions are available:

- Video recording ②⑤:
Starts or stops the recording of a video stream on the SD memory card inserted in the receptacle ①⑥.

The video recording symbol appears in the top left-hand corner of the live image during video recording.

- White balance ②⑥:
Enables the user to carry out a manual white balance. To do this, point the tip of the videoscope at a white surface and press the button for 1 second.

NOTE: With the flexible intubation videoscopes, the performance of the white balance procedure is visualized by a progress bar in the menu bar. The white balance can also be performed via the monitor menu.

- Single image capture ②⑦:
Saves the live image on the screen to the SD memory card.

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10.6 Preparation before use

NOTE: Before the device is used for the first time, it must be subjected to a complete reprocessing cycle, excluding precleaning.

⚠ ATTENTION:

Damage to product: Damaged flexible intubation videoscopes or damaged accessories must not be used.

- Whenever the flexible intubation videoscope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition.

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10.7 Checking the suction (does not apply to 11301 ABX)

Suction can be triggered by pressing the suction valve ② (when the suction pump is connected and switched on).

NOTE: There should be no instruments in the suction channel whilst suctioning is being performed. The instrument channel inlet should be sealed (Art. no. 29100) to prevent additional air from being drawn in.

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10.8 Checking the insufflation (does not apply to 11301 ABX)

Connect the O₂ tube to the instrument channel:

- Use either the olive connector (Art. no. 600007) on the instrument channel or,
- if you are using the three-way double stopcock with LUER lock (Art. no. 6927691), connect the tube to the tube holder on the angled part of the adaptor. This keeps the straight opening of the adaptor free for instrumentation.

WARNING:

Observe the instruction manuals of the supply units used.

- *Switch on the respective supply units.*

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
⚠ ATTENTION:

Damage to product:

- The surface of the flexible sheath must be inspected visually and by carefully feeling it, to detect any cracks, cuts, dents, foreign bodies or other forms of damage.

⚠ ATTENTION:

Damage to product:

- Do not squeeze the distal tip of the flexible intubation videoscope too tightly.

Every time the flexible intubation videoscope is used, you must first check that its surface, and above all the surface of the sheath and distal tip, is in perfect condition and that there are no sharp edges. Intubation videoscopes with sharp edges must not be used.

They may cause serious injury to the patient or user.

Also check the connecting cable for damage (such as fractures, cracks, twisted or crushed areas, etc.) and the connecting plug for loose connections.

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10.10 Inspecting the deflection

WARNING:

Risk of injury: The deflection lever must never be operated abruptly or with force, to avoid injuring the patient or damaging the device!

- Actuate the deflection lever gently and carefully as far as the limit stop.
- Move the deflection lever smoothly and slowly!

ATTENTION:

Damage to product: Defective flexible intubation videoscopes must not be used.

- Before each use, check the functionality of the deflection lever and the connected distal tip.

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The distal tip will deflect in accordance with the direction of actuation (upwards or downwards).

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10.11 Inspecting the image quality

⚠️ WARNING:

Risk of injury: If the image fails or there is interference during use, this can result in injury to the patient.


- Put the distal tip into the neutral position and remove it from the patient's body gently and carefully.
- If the intubation flexible videoscope continues to be used in the patient after the image has failed or when there is interference, this may result in serious injury to the patient.

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Switching off: Pressing the On/Off switch (19) for an extended period (approx. 3 seconds) turns the C-MAC® monitor off. The green indicator light (18) goes out.

NOTE: When the monitor is switched off, the integrated light source in the flexible intubation videoscope also switches off.

NOTE: Further information on the operation and functions of the C-MAC® monitor 8403 ZX can be found in the corresponding instruction manual supplied with every C-MAC® monitor (Art. no. 96076008D).

If you lose the instruction manual, you can also download it from our homepage (www.karlstorz.com).

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11 Reprocessing

NOTE: Remove the connector from the monitor or the C-HUB® II. Remove the tube holder, suction valves, lip valves and other accessories from the endoscope, and reprocess them separately.

- This reprocessing step is part of the validated process.
- This reprocessing step may be carried out additionally in order to aid cleaning. The sterilization method is approved in respect of material compatibility; however, it must be validated by the operator on site in respect of effectiveness.
- x On the grounds of potential damage, this reprocessing step must not be performed.

NOTE: Empty fields are not part of the validated process.

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WARNING:

Risk of infection:

- Only reprocess the medical device as per the documents provided with the device or the reprocessing instructions issued by KARL STORZ. These can be downloaded from www.karlstorz.com.

WARNING:

Risk of infection from non-sterile delivery.

- Reprocess the device before each use. Check the medical device for visible soiling before use.

ATTENTION:

Damage to product:

- Before reusing the device, replace any defective parts or send the flexible endoscope for repair.

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ATTENTION:

Damage to product:

- Never clean flexible endoscopes with metal brushes.

ATTENTION:

Damage to product: Constant changing between various reprocessing methods is hard on the materials and should be avoided. However, a one-off change of method, due to procurement of new cleaning equipment, for example, is possible.

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⚠ ATTENTION:

Damage to product:

- Flexible endoscopes must never be steam sterilized (autoclaved).

⚠ ATTENTION:

Damage to product:

- To prevent damage to the flexible endoscope and accessories, it must be ensured that all connections are completely dry prior to use.

NOTE: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

NOTE: Observe the chemical manufacturer's specifications concerning concentration and exposure times. Bear in mind the microbiological range of action of the chemicals.

NOTE: Damage to product. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be found at www.karlstorz.com.

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NOTE: We recommend using a plastic basin for disinfection to avoid damaging the endoscope. Furthermore, use of a plastic basin prevents electrolytic corrosion which can occur if there are different metals in a solution.

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Accessories required for carrying out reprocessing: pressure compensation cap, manual leakage tester

Place of application	Art. number brush/adaptor
Clean surface	27652
Working channel	27651 H (11303 BNX, 11301 BNX), 110910-01 (11302 BDX)
Suction channel	11301 CD1, 27651 H (11303 BNX, 11301 BNX), 110910-01 (11302 BDX)
Pressure compensation valve	13242 XA

11.2 Precleaning directly after use

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the endoscope immediately after use.

Wipe-clean the exterior surfaces of the endoscope with a damp cloth.

Rinse out the working channel in order to prevent drying.

NOTE: To avoid damaging the sheath, do not use pushing and pulling movements.

NOTE: Do not use disinfectants for precleaning; these cause fixation of proteins and organic residues.

11.3 Leakage test

The endoscope should be tested to ensure it is leakproof after use and before reprocessing.

11.3.1 Performing the leakage test – dry

The connection tube of the leakage tester and the vent port on the endoscope must be dry. Connect the leakage tester connector securely to the supply connector vent port and lock by turning the bayonet catch clockwise 25°.

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NOTE: The rubber seal in the connection adaptor of the leakage tester should be treated regularly with instrument oil.

Inflate the leakage tester until the maximum test pressure of 160–180 mmHg (blue area) is reached. The deflection rubber on the distal end may expand slightly. Using the deflection mechanism on the control body, move the deflection part in all possible directions as far as possible and observe the pointer on the manometer. The test pressure can drop by up to 10 mmHg until the endoscope is completely filled with air.

After stabilizing the pressure on the manometer, wait at least 2 minutes to see whether the pressure remains constant or drops.

11.3.2 Performing the leakage test – in liquid

! CAUTION

Damage to product: To avoid damage from the penetration of liquid, a test in the dry state must be performed in advance.

NOTE: Any lumina must be covered with no air bubbles by filling in a targeted manner.

Immerse the device in solution so that it is completely covered. It must be ensured that the leakage tester does not become wet. Perform the leakage test as described in ‘Leakage test – dry’ and note whether bubbles emerge along the device. If, after an observation period of at least 3 minutes, no bubbles can be seen emerging (continuously) from the device and, in addition, the manometer does not show a fall in pressure, the device is leakproof.

Remove the endoscope from the liquid and press the vent button on the manometer. The test pressure must drop to 0 before the leakage tester can be removed.

NOTE: Never release the pressure in the liquid.

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⚠ ATTENTION:

Damage to product: If a pressure loss occurs, this means that the endoscope leaks and may no longer be used, immersed in liquid or reprocessed by machine.

- In this case, fill out the accompanying repairs form and send the endoscope, which has been cleaned by wiping to remove organic soiling/residue, for repair or, alternatively, contact KARL STORZ. The endoscope must not be labeled as disinfected. Line the storage/carrying case with the protective film 13990 SFN.

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⚠ ATTENTION:

Damage to product:

- Never remove dirt with a needle or other sharp objects, since these can cause severe damage.

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11.4.2 Manual disinfection

The endoscope must be completely immersed in a disinfectant solution. The working channel must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time (as specified by the manufacturer), the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

11.4.3 Manual High Level Disinfection (HLD)

The endoscope must be completely immersed in a high level disinfectant solution. The working and suction channel must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time (as specified by the manufacturer), the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

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11.4.4 Drying

Finally, the surfaces and channels/lumina are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

ATTENTION:

Damage to product:

- The pressure must not exceed 100 kPa (1 bar/14.5 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.

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11.5 Machine cleaning and disinfection

The following method for machine decontamination has been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

11.5.1 Manual preparation

The endoscope must be completely immersed in a cleaning solution. The working channel must be covered with no air bubbles by filling in a targeted manner. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a suitable brush. The channel is brushed until the brush head is visibly clean. Never brush back and forth within the channel.

Finally, it must be rinsed with cold water to ensure neutralization.

11.5.2 Connecting

In order to ensure effective machine cleaning and disinfection, the instrument must be connected to the suction and working channel of the endoscope in the washer and disinfecter.

ATTENTION:

Damage to product:

- During cleaning and disinfection the pressure compensation cap must not be in place on the pressure compensation valve.

11.5.3 Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65°C and pressure of 100 kPa (1 bar / 14.5 psi) must not be exceeded. The relevant national requirements and the chemical manufacturer's instructions must be taken into account when using this method.

The selection of a suitable endoscope holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

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- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.

NOTE: Clean the distal lens with a cotton tip applicator soaked in lens cleaner or 70% isopropyl alcohol. When cleaning the distal lens, always wipe with a circular movement.

- Damaged or corroded medical devices must be withdrawn from use.

NOTE: Inspect plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their original condition, they must be replaced. This especially applies to seals and rubber caps.

- If there is still residual moisture on or in the endoscope, it must be dried again manually. The surfaces (e.g., connection) and channels/lumina must be completely dried using (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

ATTENTION:

Damage to product:

- The pressure must not exceed 100 kPa (1 bar/14.5 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.
- Afterwards, a functional check must be carried out.

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NOTE: Check for a clear view and steerability.

11.7 Packaging

Only standardized and approved packaging materials and systems may be used (EN 868 Parts 2-10, EN ISO 11607 Parts 1 + 2, DIN 58953).

11.8 Sterilization

WARNING:

Risk of infection:

- Before sterilization, all parts of the endoscope must be thoroughly cleaned, disinfected, rinsed and then dried as otherwise sterility cannot be achieved.

WARNING:

Risk of infection:

- Sterilize the medical device in a disassembled state.

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ATTENTION:

Damage to product:

- The pressure compensation cap must be in place on the pressure compensation valve for sterilization and transport to ensure pressure compensation.

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11.8.1 Hydrogen peroxide (H₂O₂) sterilization - ASP STERRAD®

⚠ ATTENTION:

Damage to product: Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

NOTE: For detailed information, please consult the user handbook of the respective device.

NOTE: Refer to the ‘STERRAD® Sterility Guide’ to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

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The following STERRAD® sterilization processes have been validated and approved by KARL STORZ for this medical device:

- STERRAD® 100S Long Cycle (without channel)
- STERRAD® 100S Long Cycle with Boosters (with channel)
- STERRAD® NX® Advanced Cycle
- STERRAD® 100NX® Flex Cycle
- STERRAD® 100NX® DUO Cycle

11.8.2 Hydrogen peroxide (H₂O₂) sterilization - STERIS® V-PRO™

NOTE: For detailed information, please consult the user handbook of the respective device.

NOTE: The ‘STERIS Low Temperature Reprocessing Device Compatibility Matrix’ can be used to verify whether the medical device can be sterilized using the various STERIS V-PRO® devices.

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The following STERIS V-PRO® sterilization methods have been validated and approved for this medical device by KARL STORZ:

- V-PRO® maX Flexible Cycle
- V-PRO® 60 Flexible Cycle
- V-PRO® 1 Plus Lumen Cycle
- V-PRO® 1

11.8.3 Ethylene oxide sterilization (EO)

The ethylene oxide procedure is validated with 100% ethylene oxide at 55°C and a hold time of 30-45 minutes.

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11.8.4 Chemical low-temperature sterilization with peracetic acid – STERIS System 1E®

Detailed information on the selection of the sterilization parameters and the STERIS® Quick Connect Kit (QCK) required to rinse out the lumina is available from the manufacturer STERIS®.

ATTENTION:

Damage to product:

- The pressure compensation cap must not be in place on the pressure compensation valve.

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WARNING:

Risk of infection:

- Improperly reprocessed medical devices may pose a risk of infection for patients, users and third parties and may cause malfunction of the medical device. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

WARNING:

Risk of infection:

- These medical devices are not supplied in sterile condition. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

WARNING:

Risk of infection:

- When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

ATTENTION:

Damage to product:

- When preparing and using solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

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 ATTENTION:

Damage to product:

- Danger of damage to medical devices: The use of chemicals that have not been approved by KARL STORZ or chemical manufacturers may cause damage to the medical devices. Only use chemicals approved for the medical device for reprocessing. A helpful guide is available on the Internet at www.karlstorz.com.

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NOTE: The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003) can be downloaded or requested by visiting www.karlstorz.com.

11.11.2 Accessories

Accessories required for carrying out reprocessing:

Brushes: 27652, 27650B

11.11.3 Manual wipe-down disinfection

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications must be observed. At the end of the necessary exposure time, wipe the surface with a dry low lint cloth.

Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, pre-clean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual pre-cleaning under cold running water.

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11.11.4 Manual precleaning

11.11.5 Disassembly

Prior to cleaning and disinfection, the medical device must be separated into its individual components as far as possible, and/or must be opened.

Brushing the surfaces

Any visible contamination or heavy soiling must be removed by cleaning the surfaces under cold running water with the aid of a brush or sponge.

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11.11.6 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual, 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003):

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Thermal disinfection is preferred. The relevant national requirements and the A₀ value must be taken into account when using this method.

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11.11.7 Connecting

In order to ensure effective machine cleaning, disinfection and rinsing, if necessary, the instrument must be connected up to the washer and disinfectant.

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- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Disassembled medical devices must be assembled.
- Afterwards, a functional check must be carried out.

NOTE: During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

11.11.9 Packaging systems

Only standardized and approved packaging materials and systems may be used (EN 868 Parts 2-10, EN ISO 11607 Parts 1 + 2, DIN 58953).

11.11.10 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003). The method must be selected taking into account the respective applicable national requirements and in consultation with the device and product manufacturers.

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Steam sterilization in the fractionated prevacuum procedure

NOTE: Sterilize greased components in a disassembled state.

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) with the following parameters:

- 132°C - 137°C with a minimum exposure time of 4 minutes and maximum of 18 minutes.

NOTE: This validation is applicable for reusable tube sets which are open at both ends and with a maximum length of 200 cm.

Or:

- 134°C - 137°C with a minimum exposure time of 3 minutes and maximum of 18 minutes.

Hydrogen peroxide (H₂O₂) sterilization - ASP STERRAD®

 WARNING:

Risk of infection:

- Medical devices must be sterilized in disassembled condition.

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The following STERRAD® sterilization methods have been validated and approved for this medical device by KARL STORZ:

- STERRAD® 50, 100S, 200 without booster
- STERRAD® NX® Advanced Cycle
- STERRAD® 100NX® Flex Cycle
- STERRAD® 100NX® DUO Cycle

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12 Technical data

12.1 Flexible intubation videoscope 11301 BNX/11303 BNX:

Direction of view:	0°
Angle of view:	100°
Deflection:	up 140°, down 140°
Working length:	65 cm
Outer diameter:	5.5 mm
Working channel:	2.1 mm
Illumination:	Integrated LED light source
Electrical power consumption:	1 W
Sensor:	CMOS technology
Weight	400 g
Drip water protection:	IPX8
Classification:	Ila

Applied part: Provides Type BF protection when
used in conjunction with a KARL STORZ camera
control unit with a BF symbol on the connection
port.

Operating/storage conditions

Operating conditions:	15°C to +40°C 15% to 80% rel. humidity
Storage conditions:	-0°C to +55°C 15% to 90% rel. humidity
Atmospheric pressure:	700 to 1080 hPa

ATTENTION:

The specifications in relation to atmospheric pres-
sure can only be guaranteed if the pressure com-
pensation cap is fitted to the pressure compensa-
tion valve during transportation.

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12.2 Flexible intubation videoscope 11302 BDX

Direction of view:	0°
Angle of view:	100°
Deflection:	up 140°, down 140°
Working length:	65 cm
Outer diameter:	4.0 mm
Working channel:	1.5 mm
Illumination:	Integrated LED light source
Electrical power consumption:	1 W
Sensor:	CMOS technology
Weight:	385 g
Drip water protection:	IPX8
Classification:	Ila

Applied part: Provides Type BF protection when used in conjunction with a KARL STORZ camera control unit with a BF symbol on the connection port.

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ATTENTION:

The specifications in relation to atmospheric pressure can only be guaranteed if the pressure compensation cap is fitted to the pressure compensation valve during transportation.

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17 Directive compliance

In accordance with the Medical Devices Directive
(MDD): 11302 BDX, 11301 BNX and 11303 BNX
are Class II a medical devices.

11301 ABX: Class I medical device.

This medical device bears a CE mark in accor-
dance with the Medical Device Directive (MDD)
93/42/EEC. A code number after the CE mark
indicates the responsible notified body.

18 Important information

To avoid the spread of infectious illnesses among
hospital staff or among KARL STORZ employees,
telescopes, instruments and equipment must
be cleaned and sterilized/disinfected before they
are sent for repair. We reserve the right to return
contaminated instruments/devices to the sender.

Repairs, changes or expansions which are not
carried out by KARL STORZ or by experts autho-
rized by KARL STORZ will invalidate all warranty
rights.

KARL STORZ gives no warranty on the correct
functioning of equipment or instruments which
have been repaired by unauthorized third parties.

19 Disposal

This device has been marked in accordance with
the European Directive on waste electrical and
electronic equipment (WEEE). At the end of its
useful operating life, dispose of the device as elec-
tronic scrap. Please ask KARL STORZ SE & Co.
KG, a KARL STORZ subsidiary or your authorized
dealer about the collection point in your area.
Within the scope of application of this Directive,
KARL STORZ SE & Co. KG is responsible for the
proper disposal of this device.

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29100 Sealing cap,
black sealing cap, for LUER lock irrigation port,
black, pack of 10 autoclavable

11301 CD1 Irrigation adaptor,
for machine cleaning, reusable, for use with
flexible intubation videoscopes 11301 BNХ and
11302 BDХ

11301CE1 Suction valve,
for single use, pack of 20, for use with flex-
ible intubation videoscopes 11301 BNХ and
11302 BDХ

10309/10310 Bronchoscope insertion tube,
size 4, with integrated mouthpiece, for single
use, sterile, insertion length 85 mm, made from
EVA, pack of 10, for use with flexible intubation
endoscopes

11301 CFX Tube holder,
for use with flexible intubation videoscopes
11301 ABX, 11302 BDХ, 11301 BNХ,
11303 BNХ and 11304 BCХ

11025 E Pressure compensation cap,
for venting during gas and plasma sterilization

13242 XL Leakage tester,
with blower and pressure gauge

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22 Electromagnetic Compatibility (EMC) information

⚠️ WARNING:

Medical electrical devices are subject to special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this appendix during installation and operation.

The flexible intubation videoscopes 11301ABX, 11301BNX, 11302 BDX and 11303 BNX in combination with the C-MAC® monitor 8403 ZX, comply with EN/IEC 60601-1-2: [CISPR 11 Class B] and therefore meet the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

The employed limit values provide a basic measure of safety from typical electromagnetic interference that is to be expected in a medical environment. The flexible intubation videoscopes 11301ABX, 11301 BNX, 11302 BDX and 11303 BNX in combination with the C-MAC® monitor 8403 ZX, are Group 1 devices (as per CISPR 11).

Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this Appendix provide information to the customer or user that is essential in determining the suitability of the device or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the device or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If electromagnetic interference arises while using the device, the user can eliminate it by taking the following actions:

- reorient or relocate the receiving device
- increase the separation between the individual devices
- connect the devices to different electrical circuits.

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Electromagnetic Compatibility

⚠ **WARNING:**

The use of portable and mobile RF communications equipment may affect this or other medical electrical equipment.

⚠ **WARNING:**

The use of accessories other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased interference immunity of the flexible intubation videoscopes 11301ABX, 11301BNX, 11302 BDX and 11303 BNX in combination with the monitor 8403ZX. The accessories listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using accessories other than those specified here, it is the responsibility of the operator to ensure that they comply with EN/IEC 60601-1-2.

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Accessories and cables which have been shown to comply with EN/IEC 60601-1-2:

Type	Shielded	Length [m]	Ferrite	Use
Charger ET27-30-0004515	No	2.0	No	Charging the internal battery of the C-MAC® monitor 8403 ZX.

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Table 2 Guidelines and manufacturer's declaration – Electromagnetic Immunity			
The flexible intubation videoscopes 11301ABX, 11301BNX, 11302 BDX and 11303 BNX in combination with the monitor 8403 ZX are designed for operation in an electromagnetic environment as defined below. The user should make sure that it is used in such an environment.			
Interference immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – Guidelines
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	Complies ± 6 kV contact discharge ± 8 kV air discharge	Floors should be made of wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity must be at least 30%.
Electrical fast transient/burst acc. to IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment. If image distortion occurs, it may be necessary to further filter the power quality by using an AC line filter.
Surges acc. to IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	$<5\% U_T^*$ (>95% dip in U_T) for 0.5 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $<5\% U_T$ (>95% dip in U_T) for 5 seconds	Complies $<5\% U_T^*$ (>95% dip of the U_T) for 0.5 cycle Complies $40\% U_T$ (60% dip of the U_T) for 5 cycles Complies $70\% U_T$ (30% dip of the U_T) for 25 cycles Complies $<5\% U_T$ (>95% dip of the U_T) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power supply interruptions, it is recommended that the device be powered from an uninterruptible power supply.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	3 A/m	Complies 3 A/m	The magnetic fields at the mains frequency should be at the values typical of business or hospital environments.
* Note: U_T is the a.c. mains voltage prior to application of the test level.			

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Table 4
**Guidance and manufacturer's declaration – electromagnetic immunity –
for equipment and systems that are not life-supporting**

The flexible intubation videoscopes *11301ABX*, *11301BNX*, *11302 BDX* and *11303 BNX* in combination with the C-MAC® monitor 8403 ZX are designed for operation in an electromagnetic environment as defined below. The user of the device should make sure that it is used in such an environment.

Interference immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – Guidelines
			Portable and mobile RF communications equipment should be used no closer to any part of the flexible intubation videoscopes <i>11301ABX</i> , <i>11301BNX</i> , <i>11302 BDX</i> and <i>11303 BNX</i> in combination with the C-MAC® monitor 8403 ZX including cables than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	$d = [3,5/3]\sqrt{P}$ 150 kHz to 80 MHz
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = [3,5/3]\sqrt{P}$ 80 MHz to 800 MHz $d = [7/3]\sqrt{P}$ 800 MHz to 2.5 GHz
			Where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^{a)} , should be less than the compliance level in each frequency range ^{b)} . Interference may occur in the vicinity of equipment marked with the following symbol: 

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^{a)} Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to stationary transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the device is used exceeds the above compliance levels, the device should be monitored to ensure proper function. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

^{b)} Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

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Table 6

Recommended separation distances between portable and mobile RF telecommunication devices and the flexible KARL STORZ intubation videoscope 11302 BDX in combination with the C-MAC® monitor 8403 ZX

The flexible intubation videoscopes 11301ABX, 11301BNX, 11302 BDX and 11303 BNX in combination with the C-MAC® monitor 8403 ZX are designed for operation in an electromagnetic environment in which the RF disturbances are controlled. The user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication devices (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according the transmitter frequency		
	150 kHz to 80 MHz $d = \left[\frac{3,5}{3} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3,5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.1	0.1	0.2
0.1	0.4	0.4	0.7
1	1.2	1.2	2.3
10	3.7	3.7	7.4

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters [m] can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com
KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za
T00 KARL STORZ Endoscopy Kasachstan
6, Saryarka str. BC „Arman“, office 314
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz
KARL STORZ ENDOSKOPE
East Mediterranean and Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com
KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in
KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, M Plaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com
KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au
KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk
KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongguiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn
KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn
KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn
KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn
KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg
KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg
KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjail-Ro, Bundang-Gu
Seongnam, 13556, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com
KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw
KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com
KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp
www.karlstorz.com

KARL STORZ SE & Co. KG

Dr.-Karl-Storz-Straße 34
78532 Tuttlingen

Postfach 230
78503 Tuttlingen

Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com

STORZ

KARL STORZ — ENDOSKOPE

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SILVER SCOPE® Video-Gastroskope,
Modelle 1382xxx

INSTRUCTION MANUAL

SILVER SCOPE® Video gastroscopes,
Models 1382xxx

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Important information for users of KARL STORZ instruments

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, these products are the result of years of experience and great care in manufacturing. You and your organization have decided in favor of modern, high-quality items of equipment from KARL STORZ.

This manual is intended to serve as an aid in the proper use, cleaning, disinfection and sterilization of videoscopes manufactured by KARL STORZ. The videoscopes are sensitive precision instruments, which require particular care during preparation. To prevent damage during preparation, guidance is provided here for the disassembly, cleaning, sterilization and assembly of KARL STORZ videoscopes. For this reason, read these instructions carefully and keep them as a ready reference.

KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope of supply, design, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, illustrations and descriptions in this manual.

It is recommended to check the suitability of the product for the intended procedure prior to use.

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1	Symbols employed	4
2	Images Gastroscope	7
3	Images Accessories	9
4	Safety instructions	11
5	Intended use	11
6	User qualification	12
7	Warnings and cautions	12
8	Safety precautions when using the videoscopes	13
9	Handling	18
9.1	Compatible video processor units/light sources	18
9.2	Connecting to a light source, for example the KARL STORZ XENON 100 SCB	19
9.3	Operating the videoscope	20
9.3.1	Preparation of the light source	20
9.3.2	Preparation of the processor unit	20
9.3.3	Connecting the light conductor, air/water connections and the suction channel	20
9.3.4	Deflection of the distal tip	21
9.3.5	Suction	21
9.3.6	Insufflation and light	22
9.3.7	Water irrigation for lens cleaning	22
9.4	Preparation before use	22
9.4.1	Testing the various functions	23
9.4.2	Handling the endoscope correctly	26
9.4.3	Additional recommendations for use	27
9.4.4	Operating the endoscope buttons	28
10	Reprocessing	29
10.1	General information	29
10.2	General cleaning	31
10.2.1	Preparation for cleaning	32
10.2.2	Leakage test	33
10.3	Channel system of a KARL STORZ SILVER SCOPE® video gastroscope	36
10.4	Channel system of a KARL STORZ SILVER SCOPE® Interventional video gastroscope 13826 PKS/NKS	37

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10.5	Cleaning/disinfection of the videoscope: manual cleaning.....	38
10.6	Manual disinfection	43
10.7	Machine preparation of the videoscope	45
10.8	Preliminary cleaning directly after the examination (see page 31).....	51
10.9	Cleaning in the preparation place.....	51
10.10	Preparation in the washer/disinfector (AER)	52
10.11	Cleaning of accessories.....	52
10.12	Preparation of water bottle 13992 BS	53
10.13	Checking preparation procedure.....	53
10.14	Sterilization	54
10.14.1	Ethylene oxide sterilization (EO)	55
10.15	Limits of reprocessing.....	55
10.16	Storage.....	55
10.17	Storage.....	56
11	Storage/transport conditions:	57
12	Disposal.....	57
13	Maintenance	58
14	Important information	58
15	Servicing and repair	59
16	Limitation of liability	59
17	Warranty	59
18	Directive compliance.....	59
19	Appendix: Troubleshooting	60
20	Subsidiaries.....	64

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1 Symbols employed

Symbol	Meaning
	Consult instructions for use
	Serial number
	Manufacturer
	Fragile, handle with care
	Keep dry
	Temperature limit
	Keep away from sunlight
	Do not reuse
	Applied part type BF, defibrillation proof, IEC 60601-1
	Potential equalization connector
	DANGER: Risk of explosion if used in the presence of flammable anesthetics.
	Protective earth (ground)
	Alternating current

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Symbol	Meaning
	Device is subject to the requirements of the WEEE Directive, 2012/19/EU.
	CAUTION: To reduce the risk of electrical shock, do not remove cover. Refer servicing to qualified service personnel. Do not store liquids on or above the unit. Keep out of reach of patients.
	On
○	Off
	Video input
	Video output
	Color video camera
	Color video monitor
	White balance
	Output

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
Symbol	Meaning
	Input
	Setup Wizard
	Signifies that this feature/option is available in the Options Menu
	Signifies that this feature/option can be controlled via the User Menu
	Signifies that this feature/option can be controlled via the buttons on the videoscope
	Federal (USA) law restricts this device to sale by or on the order of a physician.

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


Images
SILVER SCOPE® Series
Gastrosopes

- ① Up/down deflection wheels
- ② Up/down deflection brake
- ③ Programmable remote control buttons
- ④ Suction valve
- ⑤ Air/water valve
- ⑥ Rubber lip valve/biopsy valve
- ⑦ Right/left deflection wheels
- ⑧ Right/left deflection brake
- ⑨ Nozzle for lens cleaning/irrigation/insufflation
- ⑩ Light cable
- ⑪ Lens
- ⑫ Supply tube
- ⑬ Connection for leakage test
- ⑭ Connection for water bottle
- ⑮ Light cable connection
- ⑯ Port for video processor/electronic connector
- ⑰ Cap with connection for leakage tester
- ⑱ Suction/aspiration port
- ⑲ Air pin
- ㉑ Lanyard
- ㉑ Working channel
- ㉒ Additional irrigation channel with LUER lock connection (only 13826 PKS/NKS)
- ㉓ Distal water-jet channel (only 13826 PKS/NKS)

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5 Intended use

The KARL STORZ video gastroscopes are optical instruments for medical, diagnostic and endoscopic treatment in the upper gastrointestinal tract of humans.

They are inserted orally through the mouth or transnasally through the nose, via the oesophagus into the patient's stomach. In the process, they can be advanced to within the duodenum.

Gastroscopy: Examination of the stomach.

Videoscopes are applied parts, which should only be used with approved supply units (video processor units, light sources).

Contraindications:

No contraindications are known which directly relate to the product. The videoscope must not be used if, in the opinion of the responsible physician, such an application would endanger the patient, e.g., due to the patient's general condition, or if the endoscopic method as such is contraindicated.

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6 User qualification

KARL STORZ videoscopes may only be used by persons with an appropriate medical qualification and who are familiar with the endoscopic technique. The information given in these instructions only serves to instruct in the correct handling and preparation of the videoscopes.

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7 Warnings and cautions

Please read this manual and follow its instructions carefully, to avoid putting your patients, personnel or yourself at risk. The words, **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the instrument. To make the words stand out more clearly, they are accompanied by a pictogram.

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CAUTION: A Caution indicates that particular service procedures or safety precautions must be followed to avoid possible damage to the instrument.

NOTE: A Note indicates special information about operating the instrument, or clarifies important information.

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8 Safety precautions when using the videoscopes

Videoscopes must be used according to the medical rules and procedures of endoscopy recognized for this method.

WARNING: Check the videoscope prior to each procedure to ensure that it functions correctly and has been properly cleaned and disinfected.

WARNING: Verify the sound condition of the videoscope, as well as of accessories used in combination, before every application. Damaged videoscopes or damaged accessories must not be used.

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WARNING: Remove the EO/pressure compensation cap 11025 XE from the connection piece during:

- endoscopic interventions
- manual/machine preparation.

WARNING: Do not allow any electrically active components or instruments to come into contact with the videoscope.

WARNING: The videoscope must not be used during the discharge of a defibrillator. Remove the videoscope from the operation area before discharge occurs.

WARNING: Danger of explosion if flammable gases are used in the immediate proximity of the videoscope. The videoscope must not be used in the presence of flammable or combustible materials.

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WARNING: High-energy light radiated through endoscopes may give rise to high temperatures in front of the light outlet, at the end faces of the fiberoptic light cable and at the distal tip of the endoscope.

Prolonged direct contact of light-emitting surfaces on tissue can lead to burns; this must therefore be avoided.

Temperatures remain high even after the proximal connector is removed from the light source (risk of burns to the user from the light cable).

WARNING: Before removal from the patient, the distal tip of the videoscope must be returned to the neutral, nondeflected position (L and U on top of each other and positioned between the valves) to prevent damage to the instrument and/or injury to the patient.

WARNING: The videoscope may only be repaired by personnel authorized by KARL STORZ. Unauthorized repairs may impair functioning and safety, and are therefore strictly forbidden. KARL STORZ gives no warranty for videoscopes that are repaired by unauthorized personnel.

WARNING: Damage to the videoscope and injury to the patient or user resulting from incorrect operation are not covered by the manufacturer's warranty.

WARNING: Always unplug the videoscope before carrying out any maintenance and cleaning procedures.

WARNING: Risk of infection due to shipment! Contaminated videoscopes can cause serious or fatal diseases during transportation or at the recipient's premises. Please cover the carrying case with the protective sheet 13990 SFN before sending contaminated instruments back to KARL STORZ.

WARNING: If a malfunction of the videoscope should occur during an intervention on the patient, stop the intervention immediately. Place the distal tip of the videoscope in a neutral, nondeflected position and remove the videoscope slowly and carefully from the patient.

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WARNING: Before turning on the high voltage/electrical supply, ensure that there is no contact with metal parts of the endoscope or any other electricity-conducting accessory, and that the suction unit has no contact with conductive liquids.

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CAUTION: Only use accessories (guide probes, cleaning brushes, etc.) that are specifically suited for the intended use with videoscopes and which can be employed in line with generally accepted technical standards and safety guidelines.

CAUTION: Only use accessories with flexible ends. Other types of accessories may cause damage to videoscopes.

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CAUTION: Place the EO/pressure compensation cap 11025XE on the connection piece prior to:

- gas sterilization
- shipping

CAUTION: Do not expose the videoscope to direct sunlight or excessive heat, to avoid damage to the device.

CAUTION: Do not steam sterilize (autoclave) the videoscopes.

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CAUTION: The 'Preparation and Care' manual must be followed. The procedures for cleaning, disinfection and gas sterilization are explained therein in detail.

CAUTION: Do not hit the distal tip of the videoscope against hard objects, and never use a sharp object to remove dirt.

CAUTION: Never lay heavy objects on the videoscope.

CAUTION: Do not use any accessories with sharp edges or points in or on the endoscope.

NOTE: The patient leakage current limit of the CF application part may be exceeded when CF application parts (application on heart) are used in combination with other BF application parts. In this case, the reduction of the patient leakage current should be effected using only CF application parts whenever this is possible.

NOTE: The endoscope does not offer any electric insulation against high frequency (HF) currents. Protection against HF currents should be ensured through the use of the correct HF instruments. Please read the instructions for use in the instruction leaflet supplied with the HF application part for more information.

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NOTE: Any damage to the videoscope resulting from incorrect operation is not covered by the warranty.

NOTE: Constant changing between various preparation methods is especially hard on the materials and must therefore be avoided.

NOTE: Follow local governing ordinances and recycling plans regarding disposal or recycling of device components.

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9 Handling

9.1 Compatible video processor units/light sources

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NOTE: The KARL STORZ videoscope should not be used in combination with video processor units and light sources from other manufacturers. The use of non-compatible units can lead to patient injury or unit damage. It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces or voltage-carrying elements of other units at any time.

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
**9.2 Connecting to a light source,
for example the KARL STORZ
XENON 100 SCB**

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CAUTION: The correct videoscope adaptor must be used for the video processor.

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Connection of supply connector with light source and water bottle

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9.3 Operating the videoscope

9.3.1 Preparation of the light source

Please read the full instructions in the instruction manual of the respective light source.

1. Install the light source on a solid, even surface.
2. Insert the power cable into the socket provided.

9.3.2 Preparation of the processor unit

Refer to the relevant instruction manual.

9.3.3 Connecting the light conductor, air/water connections and the suction channel

1. Insert the supply connector into the light source until it audibly clicks into place (Fig. 1).
2. Connect the videoscope adaptor cable (VEA) to the endoscope. When doing so, take note of the white markings on the endoscope and the VEA cable (Fig. 2).
3. Attach the connector from the water bottle tube to the endoscope supply plug (irrigation bottle connector) (Fig. 3). The water bottle holder automatically positions itself correctly and cannot be adjusted manually.
IMPORTANT: To fill the water bottle with sterile water, remove the bottle from the light source and unscrew the lid. Only fill the water bottle 2/3 full with sterile water (1/3 is required to cushion air pressure).
4. Connect the suction tube (patient tube) from the suction pump or central suction system to the supply plug connector (Fig. 4).

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9.3.4 Deflection of the distal tip

Rotating the large deflection wheel, labeled 'U/D', (up/down) deflects the distal tip. Rotating clockwise (or forwards) lowers the tip. Rotating counterclockwise (or backwards) raises the tip.

Right and left movements are controlled using the smaller deflection wheel (labeled 'L/R', left/right). Rotating clockwise (or forwards) = right / rotating counterclockwise (or backwards) = left (Fig. 5).

WARNING: Do not rotate the wheels forcibly, this could damage the instrument or cause injury to the patient! Before attempting to introduce or remove the instrument, release the deflection wheel brakes!

9.3.5 Suction

Switch on the suction pump. To perform suction, press the suction valve (marked red) on the endoscope (Fig. 6a).

Suction can also be used when the biopsy forceps are placed in the biopsy channel and the rubber lip valve is mounted. During the procedure, the rubber lip valve/biopsy valve seals off the forceps or the accessories.

With the interventional gastroscope (13826 PKS/NKS), the suction channel can also be determined after the suction valve with the two lugs has been positioned in the direction of the distal end.

Should the pressed-down valve be aligned such that the mark points towards the distal end, suction will be performed via both working channels by pressing the suction valve (Fig. 6b).

Should the pressed-down valve be turned anticlockwise, suction will be performed via the left working channel (dia. 2.8 mm) by pressing the suction valve (Fig. 6c).

Should the pressed-down valve be turned clockwise, suction will be performed via the right working channel (dia. 3.8 mm) by pressing the suction valve (Fig. 6d).

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9.3.6 Insufflation and light

Start the light source, switch on the insufflation pump.
For insufflation, cover the opening on the blue valve with a finger (do not depress the valve, since this starts water irrigation) (Fig. 7).

NOTE: To avoid patient discomfort, use insufflation sparingly. Excessive insufflation may cause cramps, spasms and injuries!

9.3.7 Water irrigation for lens cleaning

Pressing on the air/water valve starts the water flow and/or lens cleaning (Fig. 8).

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9.4.1 Testing the various functions

1. Prepare light source as described in 9.3.1.
2. Connect supply connector, water bottle and suction tube according to 9.3.3.
3. Check the supply cable for damage (e.g., breaks, cracks, twisted or squashed areas...) and the supply connectors for loose connections.
4. Examine the surface of the flexible insertion shaft visually and by touch for cracks, cuts, bulges, foreign bodies or other damage.
5. Testing deflection – upward/downward control

The large deflection wheel controls up and down movement. To test up and down movement,

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release the up and down brake. Slowly rotate the large deflection wheel clockwise (or forwards), the distal end should move downwards. Rotating the deflection wheel counterclockwise (or backwards) causes upward movement.

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6. Testing deflection – right/left control
7. The smaller/outer deflection wheel controls right/left movement. To test right/left movement, release the right/left brake. Slowly rotate the small deflection wheel clockwise or forwards, the distal tip should move to the right. Rotating the deflection wheel counterclockwise (or backwards) causes movement to the left.
8. Testing deflection – brake function
The deflection wheels should turn freely when the brakes are released.
9. Release the brakes and test their action (see 5a/b). Activate the brakes by moving the brake lever counterclockwise – ensure that the deflected position of the distal endpiece is fixed in the desired position.

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Do not twist or bend the deflection piece by hand.

WARNING: If the deflection and/or brake mechanisms do not function properly, do not use the device. Contact the KARL STORZ service department.

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10. Check the air supply. Switch on the insufflation pump and hold the distal end of the instrument under sterile water. Cover the air control opening in the air valve (blue ring) with a finger; bubbles should emerge from the air nozzle at the distal end.

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11. Check the function of the water supply:
Depress the air/water valve (blue ring) fully.
Water should flow out of the water nozzle (distal) and rinse the lens clean.
IMPORTANT: When you first actuate the valve, it may take a couple of seconds for water to start flowing.

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12. Test the suction function. With the unit correctly installed and the suction pump or central suction system switched on, immerse the distal end in a container of sterile water once again. Depress the suction valve (red ring) fully. Water should be suctioned through the suction channel of the instrument via the connected suction pump and into the collecting bottle of the suction equipment.

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13. Test the biopsy channel by introducing biopsy forceps.

CAUTION: It should be possible to guide the biopsy forceps smoothly and without problems through the entire channel.

CAUTION: Full (>180°) deflection can lead to limitation in movement, depending on the accessories. If necessary, adjust the deflection slightly.

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Handling the endoscope correctly

Hold the control handle of the endoscope in the left hand in such a way that the thumb, index and middle fingers can easily control the large deflection wheel. In this position, the suction valve can be depressed with either the index or the ring finger. The air/water valve can also be covered with these fingers. The remote control button (e.g., Stills) is also within reach.

The small, outer deflection wheel can be controlled with the fingers of the left hand, or alternatively with the thumb and fingers of the right hand.

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9.4.3 Additional recommendations for use

1. Use a lubricant from KARL STORZ (see list on our website www.karlstorz.com under 'Hygiene') for the distal end of the insertion tube.
2. Observe the lumen the whole time. If you become disoriented, withdraw the instrument gently until vision is restored.

NOTE: Do not put lubricant on the lens to avoid clouding!

NOTE: Do not use oil or petroleum, silicone, Vaseline or alcohol based media, to avoid damage to the endoscope material.

NOTE: If debris causes the lens to cloud during the procedure, this can normally be rectified by a combination of water and insufflation. Excessive water or air can be removed by applying simultaneous pressure on the suction valve.

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9.4.4 Operating the endoscope buttons

Depending on the video processor used, the buttons can have two functions, which are initiated by means of short and long presses of the button. For more detailed information on button assignment, please consult the manual for the camera control unit used. A list of compatible camera control units can be found in section 9.1.

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NOTE: In case of 'red-out', a blurred red picture on the screen, caused by the distal endoscope tip resting directly on the mucosa, withdraw the instrument gently until vision is restored.

WARNING: Pushing the instrument forwards despite a 'red-out' effect may result in serious injury to the patient.

NOTE: If the angling of the distal tip can no longer be properly controlled the examination must be ended immediately. Release the brakes and place the tip in the neutral position. Position the deflection wheels so that 'U' and 'L' are between the valves (see page 13/15), then carefully withdraw the instrument, whilst continuing to monitor the procedure.

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WARNING: Remember that decontamination is only then successful if

- the instruments have been freed of all organic material and chemical residues.
- work proceeds according to recognized and validated methods.
- any deviation from recommended procedures are validated and tested by the user.
- the recommended sterilization parameters are maintained. These are only valid with cleaning, disinfection and sterilization equipment that is properly maintained and calibrated.

CAUTION: The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' must be followed. The cleaning, disinfection and sterilization procedures are explained therein in detail. The list of KARL STORZ approved materials for the preparation of videoendoscopes can be found on the Internet at 'www.karlstorz.com' or 'www.karlstorz.de' under 'Hygiene' and also in the appendix of this manual.

CAUTION: When preparing and using the cleaning and disinfection solutions, follow the manufacturer's instructions paying close attention to proper concentration and exposure times. Disregarding the instructions can lead to damage. Bear in mind the microbiological range of action of the chemicals used.

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CAUTION: Never clean videoscopes with metal brushes.

CAUTION: Constantly changing between different preparation methods causes particular wear on the materials and must therefore be avoided. However, a once-only change in methods, e.g., when acquiring new cleaning equipment, should be carried out in consultation with the manufacturers of the endoscope, the chemicals and the washer/disinfector.

This change in method is then usually completely harmless.

CAUTION: Do not clean videoscopes in an ultrasound bath.

CAUTION: Do not steam sterilize (autoclave) videoscopes.

CAUTION: Never immerse videoscopes in normal saline solution, since even brief contact can lead to pitting and corrosion.

CAUTION: During manual preparation, videoscopes may only be rinsed with microbiologically clean/sterile water.

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10.2 General cleaning

Cleaning should only be performed by authorized, trained personnel. All safety regulations should be observed (gloves, protective glasses, mouth and nose protection, clothing, German Technical Rules for Biological Agents (TRBA 250).

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The national guidelines and recommendations for the preparation of instruments (for example RKI, ESGE/ESGENA, ASGE, LANGS, SAGES etc.) should be observed and the effectiveness verified.

CAUTION: Before immersion in cleaning solution, check the instrument for leaks. Otherwise fluid penetrating the videoscope could cause damage.

CAUTION: Do not autoclave the endoscope or accessories (with the exception of accessories that are specifically designated as autoclavable, e.g., forceps).

NOTE: The endoscope must be cleaned immediately after every use. If organic material is allowed to remain and dry in the instrument, complete or partial blockage of the channel may result. Thorough cleaning is a prerequisite for effective preparation. This will avoid coagulation residue / blockages and the resulting repairs and instrument failures. For reasons of material compatibility, only chemicals approved by KARL STORZ should be used for cleaning and disinfection (see on our website www.karlstorz.com under "Hygiene").

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

10.2.1 Preparation for cleaning

1. The preliminary cleaning of the endoscope is performed at the place of operation. The main cleaning and disinfection is carried out in a separate preparation room.
2. Immediately after use, remove any heavy soiling (secretions, etc.) and organic material from the outer surfaces of the insertion tube with a disposable cloth moistened with cleaning solution. Then submerge the distal end in a container with cleaning solution and activate the suction to rinse the working channel.
3. Actuate the suction valve and suction at least 100 ml of cleaning solution through the suction channel. (Fig. 14). Suction up at least 200 ml with the interventional gastroscope (13826 PKS/NKS). The suction valve must then be operated in all three positions as described in section 9.3.5 (Left --> left working channel dia. 2.8 mm, right --> right working channel dia. 3.8 mm and middle --> through both channels). Make sure that all the channels are free from blockages.
4. Use the cleaning valve 13991 SRV (Fig. 13) in place of the air/water valve. (Do not use this valve during examination!).
5. Check that the water bottle is filled (max 2/3). Start the air supply (insufflation pump). Water and air are now passed through both the air channel and the water channel, depending on the position of the cleaning valve. (The process should run for at least 10 seconds).
Air channel = Position 1
Water channel = Position 2
6. It is important that both valve settings ('Actuated'/'Not actuated') are used multiple times for several seconds. This is the only way to ensure that both channels are rinsed through.
7. Turn off video system.

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8. Pull out the video connector from the videoscope adaptor (VEA) at the supply connector of the endoscope when the video processor unit is switched off. Put on the sealing cap for the video cable connector and lock it with a 1/4 turn.
9. Disconnect the endoscope from the video processor unit and light source. Place the cleaning valve in a cleaning solution.
10. Then transport in a protected case to the site of reprocessing.

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10.2.2 Leakage test

Preparing for the leakage test – dry

NOTE: Place the cap on the videoscope adaptor (VEA) connector and screw it on tightly. Ensure that the cap is secured properly (Fig. 16).

1. The connection tube of the leakage tester and the vent port on the endoscope must be dry. Connect the leakage tester connector to the supply connector and lock by turning clockwise 90° (Fig. 17).

NOTE: The rubber seal in the connection adaptor of the leakage tester 13242 LX should be treated occasionally with instrument oil.

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Before every immersion in liquid, test the endoscope for leaks. Remove all valves.

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1. Inflate the leakage tester until the maximum test pressure of 160–180 mmHg (blue area) is reached. The deflection rubber on the distal deflection part of the endoscope expands slightly. Move the deflection part with the deflection wheels on the control piece in all four directions as far as possible and check the manometer to see whether the pressure remains constant or falls (dry test).
2. Wait at least 1 minute.

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Performing the leakage test – in liquid

1. Immerse the device in water so that it is completely covered. It must be ensured that the leakage tester is not immersed but rather kept dry. Perform the leakage test as described above, note whether bubbles emerge from the device. If after an observation period of at least 3 minutes no more bubbles can be seen emerging from the device and, in addition, the manometer does not show a fall in pressure, the device is leakproof and can be used.

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If this test is also negative and no bubbles can be seen rising in the liquid, the leakage tester can be removed.

To do this, remove the supply connector from the liquid and press the vent button on the manometer to vent the endoscope.

NOTE: Never release the pressure in the liquid. For this, remove the connector from the liquid and press the vent button. The test pressure must fall to 0. Then remove the leakage tester.

If a pressure loss occurs, this means that the videoscope leaks and may no longer be used, immersed or reprocessed by machine.

CAUTION: In this case, fill out the accompanying repairs form (see Appendix), roughly clean and remove any organic soiling/ residue from the endoscope, mark it as not disinfected and return it to the Service Department for repair or, alternatively, contact KARL STORZ. Line the storage/carrying case with the enclosed transport protection sheet 13990 SFN.

WARNING: Risk of infection during shipment! Contaminated instruments can cause aerogenous or contact infections during transportation or at the recipient's premises. Shipping without the plastic film packaging must be avoided.

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**10.3 Kanalsystem eines
KARL STORZ SILVER SCOPE®
Video-Gastroskops**


(13820 PKS/NKS, 13821 PKS/NKS,
13823 PKS/NKS)

**10.3 Channel system of a
KARL STORZ SILVER SCOPE®
video gastroscope**

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
10.4 Channel system of a
KARL STORZ SILVER
SCOPE® Interventional video
gastroscope 13826 PKS/NKS

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
10.5 Cleaning/disinfection of the videoscope: manual cleaning

Please adhere to the national guidelines, e.g., ESGENA guidelines, DEGEA/RKI recommendations, ASGE, LANGS, SAGES etc. for cleaning, disinfection and sterilization, as well as to the instructions of the chemical manufacturer concerning concentration and exposure times.

1. Remove the rubber lip valve/biopsy valve, suction valve and air/water valve and place in a cleaning solution and clean using a soft disinfected brush (13991 VB) or disposable brush.

2. Wipe the housing, insertion shaft, and supply tube carefully with a lint-free cloth soaked in a cleaner (for suitable cleaning agents or disinfectants approved by KARL STORZ, see www.karlstorz.com). Avoid shocks and pulling movements so as not to damage the insertion shaft. Then rinse with microbiologically pure/sterile water.

3. Immerse the endoscope in the cleaning solution and rinse the channels using a syringe and the supplied irrigation tubes -13820 PKS/NKS, 13821 PKS/NKS, 13823 PKS/NKS (Fig. 23a, 23b): 13991 SS, 13991 YS; -13826 PKS/NKS (Fig. 24a, 24b): 13991 DKSA, 13991 SWJ, 13991 JASA).
4. The channels must be filled without any bubbles, to guarantee a complete covering of the surfaces.

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
- When submerged, use a suitable brush (preferably disposable) (see accessories, page 9/10) to brush below the surface of the liquid in one direction starting from the supply plug (supply plug --> control body, control body --> distal tip, biopsy channel --> distal end) (Fig. 25a: 13820 PKS/NKS, 13821 PKS/NKS, 13823 PKS/NKS; Fig. 25b and 25c: 13826 PKS/NKS)). Repeat these procedures at least 3x.

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- Clean the brush head (if you are using reusable brushes) at the distal end before retracting it. Do not brush backward and forward, rather brush in one direction only.

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

CAUTION: Cleaning with brushes must be carried out with the instrument immersed!

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7. Thoroughly clean out the entire suction channel (shaft part and supply part) and the biopsy channel a number of times with a brush. Observe the exposure time for the cleaning solution.

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
8. Carefully clean all particles of dirt from the distal tip. It may be necessary to clean the air/water valve with a brush (13991 VB) with soft bristles, in order to remove hardened contaminants. Never remove debris with a needle or other sharp objects, because these can cause damage.

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
9. Connect the irrigation tubes:
- 13820 PKS/NKS, 13821 PKS/NKS, 13823 PKS/NKS: 13991 SS, 13991 YS;
- 13826 PKS/NKS: 13991 DKSA, 13991 SWJ, 13991 JASA) and pump fresh cleaning fluid through the channels using a syringe. Afterwards, rinse with microbiologically clean/sterile water.

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

Suitable brushes

Art. no.	Designation	Used for
110915-01	Cleaning brush dia. 3 mm	Working channel 2 mm
110910-01	Cleaning brush dia. 5 mm	Working channel from 2.8 to 4.2 mm
13991 VB	Valve brush	Valve seat
110917-01	Cleaning brush dia. 1.5 mm	Water-jet channel 13826 PKS/NKS

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- Remove the device from the cleaning solution and rinse with microbiologically pure/sterile water.
 - Carefully clean all valves with the disinfected cleaning brush (13991 VB) and then rinse.
 - The valves should also be cleaned in an ultrasound bath.
 - Flow control for the interventional gastroscope 13826 PKS/NKS:
- Remove the irrigation tubes from the device
 - Close the working channel inlets on the handle using a biopsy valve/rubber lip valve ((Fig. 30a) 13991 DA).
 - Connect the irrigation adaptor 13991 CS to the device's supply sheath (Fig. 30b)
 - Connect the suction cleaning valve ((Fig. 30c) 13991 ARV) in the valve seat of the suction valve on the handle in two subsequent positions so that it engages perceptibly in place each time.

CAUTION: The valve 13991ARV must be set in the direction of the right working channel. The marker on the valve then points towards the right working channel (Fig. 33d). Then valve must then be set in the direction of the left working channel. The marker on the valve points towards the left working channel (Fig. 33e). At least 50 ml of microbiologically clean / sterile water must be pumped through the irrigation tube in both valve settings.

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It must be checked that the working channels are not blocked and that the amount of liquid pumped through emerges again at the distal tip.

Once the cleaning work has been completed, clean the cleaning instruments, for example the brushes and irrigation tubes and disinfect/sterilize them. Never reuse damaged or kinked brushes.


CAUTION: Do not autoclave the endoscope or its accessories (with the exception of accessories that are specifically designated as autoclavable, e.g., forceps).

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
10.6 Manual disinfection

1. After careful rinsing, fully immerse the endoscope in the disinfecting solution.

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! **CAUTION:** When preparing and using the solutions, follow the chemical manufacturers' instructions paying close attention to proper concentration and exposure times. Prolonged immersion can lead to damage. Bear in mind the microbiologic range of action of the chemicals used.

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2. Using the irrigation tubes and a syringe with LUER connector, rinse out the channels a number of times with disinfecting solution.

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3. Then rinse the channels and surfaces of the video gastroscope with microbiologically pure/sterile water.

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4. After the final step in the cleaning process has been completed, dry the entire endoscope. Thoroughly dry off all external surfaces of the endoscope, preferably with sterile compressed air at a max. of 0.5 bar/7.25 psi, thereby removing any residues of cleaning and disinfection solution.

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If at all possible, store the videoscope in a dry cupboard or hang it up in a storage case for safe storage.

NOTE: The valves of the videoscope must not be inserted; they must be stored separately.

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10.7 Machine preparation of the videoscope

All flexible KARL STORZ endoscopes can be connected to approved washer/disinfectors for semi- or fully automatic cleaning and disinfection. Observe the national laws and regulations for cleaning and disinfection procedures. Leakage testing, manual cleaning and (where necessary) rinsing with microbiologically pure water must be carried out beforehand. The manufacturer of the cleaning machine must be consulted when selecting the preparation programs. Usually, special programs for flexible endoscopes are available which are suited for decontamination and compatible with the materials. For machine preparation, make sure that the flexible endoscopes are not heated above 65°C.

CAUTION: Check the videoscope for leaks before every disinfection in semi- or fully automatic equipment.

CAUTION: Ensure that only KARL STORZ approved agents are used for procedures involving peracetic acid.

Lay the endoscope carefully in the tray or basket of the machine. Connect the irrigation tubes supplied with the cleaning adaptor set 13991 CK or the irrigation tubes of the manufacturer of the washer/disinfectant to the endoscope.

Cleaning adaptor set 13991 CK

The cleaning adaptor set 13991 CK is used for validated preparation of KARL STORZ videoscopes from series 13xxx NKS/PKS (acc. to ISO 15883 part 4). The cleaning adaptor set is not intended for any other purpose as this may damage it.

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

- 1** {
- 1 Connection tube 13991 CL for leakage tester connection ⑤ of the videoscope
 - 2 Connector (13991 CL) for vent port
 - 3 LUER connector (13991 CL) for washers/disinfectors for endoscopes

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- 4** {
- 4 Cleaning adaptor 13991 CR for air/water (A/W) channel ⑦ of the videoscope
 - 5 Connector (13991 CL) for A/W channel (air pin ②)
 - 6 Connector (13991 CL) for water channel (water bottle connection ⑦)
 - 7 2x LUER connectors (13991 CR) for washers/disinfectors (1x LUER for air channel, 1 LUER for water channel)

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- 8** {
- 8 Cleaning adaptor 13991 CS for suction channel (aspiration) ② of the videoscope
 - 9 Connector (13991 CL) for suction channel on the supply connector ②
 - 10 LUER connector (13991 CS) for washers/disinfectors

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- 11** {
- 11 Connection adaptor 13991 SS for biopsy/working channel ⑨ of the videoscope (with 2x double channel endoscopes)
 - 12 Connector (13991 CL) for biopsy/working channel ⑨
 - 13 LUER connector (13991 SS) for washers/disinfectors

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- 14** {
- 14 Channel separator 13991 CD for suction valve ④ and irrigation/insufflation valve (A/W) ⑤ of the videoscope
 - 15 Valve pin for irrigation/insufflation valve (A/W) ⑤
 - 16 Closure for suction valve ④

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WARNING: This product is not sterile when delivered and must therefore be cleaned, disinfected and if necessary, sterilized prior to initial use and each subsequent reuse.

WARNING: Risk of infection by improper preparation: Preparation may only be carried out by experienced and trained personnel.

CAUTION: Only use chemicals which have been tested for material compatibility by KARL STORZ. If you use other chemicals, get written confirmation from the chemical manufacturers that the material is compatible with KARL STORZ endoscopes.

CAUTION: If the preparation chemicals are changed, it is essential that you consult the manufacturer of the chemicals and of the washer/disinfector.

CAUTION: Please refer to our operating instructions for details on the cleaning preparation under the section: Cleaning, disinfection and sterilization.

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The following must be checked before use:


- visual checks of the connection tubes (patency) and the O-ring (wear)
- check the connections are correctly and securely in place
- check for damage

In the event of noticeable problems or complaints, please contact our service department or send in the cleaning adaptor set for technical examination.

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Assembly cleaning adaptor set

Plug the connector ② of the connection tube 13991 CL into the vent port ⑮ of the videoscope at the supply connector and secure by turning clockwise.

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Plug the cleaning adaptor 13991 CR ⑤ as far as it can go into the connector ② of the air pin of the videoscope.

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Then plug connector ⑥ as far as it will go into the 'Connection for water bottle' ⑰ of the videoscope.

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Plug connector ⑨ of the cleaning adaptor 13991 CS as far as it will go into the suction channel ⑳ at the supply connector of the videoscope.

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
Press the slide mechanism of the channel separator 13991 CD (14) downward and slide it forward until the channel separator 13991 CD (14) audibly locks into place. Check that the channel separator is securely in place.

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! NOTE: If necessary, remove the cap of the biopsy/working channel (6) of the videoscope.


Connect the connection adaptor 13991 SS with the connector (12) to the biopsy/working channel (6) of the videoscope.

! WARNING: Danger of infection due to inadequate rinsing caused by kinked tubes (decreases cleaning performance): avoid kinking of the irrigation tubes.

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
! CAUTION: When the program has finished, the videoscopes are at a risk of corrosion due to the conditions in the preparation machine.

! CAUTION: Immediately after the end of the program cycle, the videoscopes are still hot, as they have been heated up during chemical thermal preparation. Consequently the instrument's material properties are different and, as such, the prepared endoscope should be allowed to cool down for at least 20 minutes at room temperature before being used again.

! CAUTION: Ideally, remove the videoscopes from the machine immediately after the program has finished. Then dry off the optical surfaces, if necessary, with a cotton swab soaked with 70% alcohol.

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Store the endoscopes freely suspended, dust-free and dry. Store valves separately

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10.8 Preliminary cleaning directly after the examination (see page 31)

- The outer surfaces should be wiped down from top to bottom with a lint-free disposable cloth that is soaked with cleaning solution.
- Ample amounts of cleaning solution should be suctioned through the distal end.
- Apply the air-water supply for a few seconds every cycle
- Connect the cleaning adaptor
- Disconnect the flexible endoscope from the light source

CAUTION: Close the video connector with the cap.

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Trocknung

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10.9 Cleaning in the preparation place

- Carry out leakage test
- Immerse the flexible endoscope
- Clean the channels with special cleaning brushes
- Clean the complete air/water system thoroughly
- Rinse the channels

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10.10 Preparation in the washer/ disinfectant (AER)

- The machine must fulfill all the cleaning, disinfection and material compatibility requirements (validated process).
- Make sure that the environmental conditions with regard to temperature and pressure are not exceeded (+60 °C (tolerance +/- 5%) and 1 bar for the channels, 0.3 bar leakage test).
- Recommendation: Before placing the endoscope in the machine, clean all of the channels (working channels, water channels, air/water channels) regardless of whether or not they have been used.
- We recommend sticking to one type of cleaning agent.

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10.11 Cleaning of accessories

1. Thorough cleaning is required before disinfection or sterilization. Carefully clean the biopsy forceps and brushes whilst they are immersed in a cleaning solution using a soft, disinfected brush. Take special care to remove residues such as blood, excrement and mucous from the biopsy forceps and needles – an ultrasound bath is recommended.
2. Rinse in microbiologically pure/sterile water and dry (using sterile compressed air).
3. Disinfect accessories (work procedures as with cleaning).
4. Always test instruments after cleaning/disinfection for proper functioning and possible damage.
5. Before sterilization, carefully oil movable parts. Use only suitable, steam-resistant, physiologically safe, residue-free oil (e.g., 27656 B). After oiling, open and close the forceps several times. For details on preparation, please read the instruction leaflet of the respective accessory.

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10.13 Checking preparation procedure

Before assembly it should be checked that all parts and surfaces are undamaged, clean and dry.

Replace defective parts or send them in for repair. Before sending them in for repair, all contaminated parts must be cleaned/disinfected and sterilized or clearly labeled as unsterilized.

If residues are still present after cleaning, prepare the parts again, subject them to special treatment, discard or repair.

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Test the instruments for proper operation (deflection of the distal end, brake lever). Check the lens for a clear view, the instrument for appropriate deflection and all channels for patency (see page 9).

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10.14 Sterilization

Sterilization of instruments to be used on patients is not necessary when coming into contact with the skin and mucous membrane, as well as when introduced through natural, unsterile orifices.

When used on patients via sterile tissue, a transmucosal access or piercing of the mucosa, sterilization of the instruments is necessary.

WARNING: Before sterilization, all parts of the endoscope must be thoroughly cleaned, disinfected, rinsed and then dried as otherwise sterility cannot be achieved.

WARNING: Risk of infection: Sterilize the medical device in a disassembled state.

CAUTION: Flexible endoscopes must never be steam sterilized (autoclaved).

CAUTION: The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.

NOTE: All deviations from the recommended sterilization parameters must be validated by the user.

NOTE: When using reprocessing equipment, please observe the equipment manufacturers' instructions.

NOTE: Sterilization is not possible on surfaces which have been lubricated and oiled.

CAUTION: The STERRAD process is not suitable for sterilization of videoscopes. This process damages the videoscope!

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The procedure must be selected taking into account the respective applicable national requirements and in consultation with the device manufacturers.

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The following sterilization procedures have been validated and approved by KARL STORZ for this medical device:

10.14.1 Ethylene oxide sterilization (EO)

The ethylene oxide procedure is validated with 100% ethylene oxide at 55°C and a hold time of 30 to 45 minutes.

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10.16 Storage

The storage method must be selected taking into account the respective applicable national requirements.

KARL STORZ only recommends dust-protected, suspended storage for decontaminated endoscopes.

Horizontal storage is possible for sterilized endoscopes. KARL STORZ recommends the use of a specially designed tray for sterilization and storage.

CAUTION: Do not expose flexible endoscopes to direct sunlight or excessive heat.

CAUTION: The distal tip of the flexible endoscope must be in the neutral position during storage.

CAUTION: Do not store the flexible endoscope in the case.

CAUTION: The sheath of the flexible endoscope must never be kinked or coiled up too tightly.

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14 Important information

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and, if necessary, disinfected/sterilized before they are sent for repair. We reserve the right to return contaminated instruments, units or endoscopes to the sender.

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KARL STORZ gives no warranty on the correct functioning of units, endoscopes or instruments which have been repaired by unauthorized third parties.

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15 Servicing and repair

Defective endoscopes must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original parts.

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16 Limitation of liability

As a supplier of this endoscope, KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:
All assembly, operations, system expansions, readjustments, modifications, or repairs to same have been performed by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

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17 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see “Subsidiaries” section), even during the warranty period.
Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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18 Directive compliance

According to Medical Device Directive:
This medical product belongs to Class IIa
This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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1 Preface

Hygiene must be treated as an integral responsibility and a fundamental component of the whole process of medical treatment.

This principle begins with the design of instruments. Alongside the technical solution, effective cleaning, disinfection and sterilization represent the central factors in the development of all instruments.

Particularly in the case of instruments that are used in the field of neurosurgery or in direct contact with the central nervous system (CNS), particular attention and care must be paid to reprocessing with respect to the potential risk of transmission of Creutzfeldt-Jakob disease (CJD or vCJD).

The objective is to create an environmentally-friendly and long-lasting reusable instrument, the materials, design and quality of which allow reprocessing using a suitable process that takes all applicable standards and laws into consideration, while maintaining the full functionality of the instrument.

As the supplier of these medical devices, KARL STORZ shall be liable for the safety, reliability and performance of the devices only where installation, maintenance and repairs have been carried out exclusively by personnel authorized by KARL STORZ and the device has been used as directed in the instruction manual.

WARNING: KARL STORZ telescopes are not sterile when delivered, and must be cleaned, disinfected and/or sterilized prior to initial use and each subsequent reuse.

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

- ① Screw base, for KARL STORZ
fiberoptic light cables and Olympus/
Winter & Ibe
- ② Plug-in base, 9 mm in diameter, for Wolf
fiberoptic light cables
- ③ Light inlet piece, without adaptor
- ④ Color-coded ring to indicate the telescope's
direction of view:
Green = 0°, straight forward
- ⑤ Eyepiece cup
- ⑥ Model designation
- ⑦ Coupling cone for adaptation of
KARL STORZ accessories
- ⑧ Lens and distal outlet of light cable

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2 Terms used

AKI	Instrument Reprocessing Working Group
BfArM	German Federal Institute for Drugs and Medical Devices
DGKH	German Society for Hospital Hygiene
DGSV	German Society for Sterile Supply
DIN	German Institute for Standardization
EN	European Committee for Standardization
H ₂ O ₂	Hydrogen peroxide
ISO	International Organization for Standardization
MPBetreibV	German Medical Device Operator Ordinance
OP	Surgical department
PPE	Personal protective equipment
WD	Washer and disinfectant
RT	Room temperature (20 – 25 °C)
US	Ultrasound
CD water	Completely demineralized water
WHO	World Health Organization
CNS	Central nervous system

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3 Intended use, indications and contraindications

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Intended use

The telescope is used for endoscopic imaging during neurosurgery in direct contact with the central nervous system. The telescope is designed for short-term use in surgically invasive interventions.

Indications

Endoscopic intracranial neurosurgery

- Obstructive hydrocephalus (ventriculostomy)
- Septostomy, foraminoplasty, aqueductoplasty
- Positioning and removal of ventricular catheters and stents
- Resection and fenestration of intracranial cysts
- Resection and biopsy of intraventricular tumors
- Removal of blood clots
- Treatment of aneurysms

Endoscopic skull base and pituitary surgery

- Resection and treatment of tumors
- Resection and treatment of cysts
- Vascular diseases

Endoscopy-assisted microsurgery (EAM)

- Resection and treatment of tumors
- Resection and treatment of cysts
- Cranial nerve compression syndromes
- Vascular diseases

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- Overly small foramen of Monro
- Overly narrow ventricular system
- Overly thick and non-transparent floor of the third ventricle
- Vascularized tumors

Moreover, use of the instrument is contraindicated if, in the opinion of the responsible physician, the health of the patient is endangered through such use, due to the general condition of the patient, or if the endoscopic method itself is contraindicated.

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Use of the instrument is contraindicated if, in the opinion of the responsible physician, the health of the patient is endangered through such use, due to the general condition of the patient or if the endoscopic method itself is contraindicated.

Endoscopy-assisted microneurosurgery (EAM)

Use of the instrument is contraindicated if, in the opinion of the responsible physician, the health of the patient is endangered through such use, due to the general condition of the patient or if the endoscopic method itself is contraindicated.

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* corresponds to the maximum outer diameter here

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Accessories

see catalog NEUROSURGERY
(Art. no. 96020000DER)

Spare parts

see catalog NEUROSURGERY
(Art. no. 96020000DER)

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4 Safety instructions

4.1 Symbols and hazard information

Please read this manual carefully and follow its instructions exactly. The words Warning, Caution and Note convey special meanings. Wherever they are used in this manual, the following text should be carefully reviewed to ensure the safe and effective use of this instrument. To make these words stand out more clearly, they are accompanied by a pictogram.

WARNING: A warning indicates that the personal safety of the patient or personnel may be endangered. Failure to observe a warning could result in injury to the patient or personnel.

CAUTION: A caution indicates that particular procedures or precautions must be followed to avoid possible damage to the instrument set.

NOTE: A note indicates special information regarding operation or clarifies important information.

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WARNING: Risk of injury: Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all of the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

WARNING: Risk of injury and risk of damaging the products: The highly concentrated light energy can cause the distal end, the light connection, adjacent components and tissue in front of the light emission window to heat up. This can cause burns to the patient, user and surgical accessories. Avoid direct tissue contact with the distal end, as well as with the light connection of the telescope and the light cable. Do not rest the endoscope or the light cable on the patient or where they are in direct contact with surgical accessories. Always select the lowest possible light output of the cold light source which still allows optimum illumination of the surgical field.

WARNING: Risk of injury: Patient leakage currents may accumulate if endoscopes and endotherapy devices supplied with energy are used together. This is particularly important if an endoscope with a Type CF applied part is used. In this case, an endotherapy unit with a Type CF applied part must be used in order to minimize the total patient leakage current.

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WARNING: Risk of injury: Overloading the instrument by exerting too much force may cause the medical device to break, bend, and malfunction, and consequently injure the patient or user. Do not overload the instruments. Do not bend bent instruments back to their original position.

WARNING: Risk of injury: Incorrect use of medical instruments poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

WARNING: Risk of injury and risk of damaging the products: If applied parts are used outside of the field of vision, there is a risk that tissue and accessories could be damaged unintentionally. Always hold the applied parts of the active electrode, laser and other instruments that transmit energy pointed at the target in the field of vision during use.

NOTE: The telescopes referred to in this manual usually have separate instruction manuals. Other general regulations (Intended Use, User qualification, Safety precautions, etc.) must be taken from these separate instruction manuals and observed!

NOTE: For the meanings of the symbols printed on the packaging, please refer to the 'Packaging symbols' accompanying instruction leaflet, mat. no. 96216316DF. It can be downloaded from www.karlstorz.com or requested directly from info@karlstorz.com.

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4.3 User qualification

KARL STORZ telescopes and HOPKINS® telescopes may only be used by persons with an appropriate medical qualification and who are acquainted with endoscopic techniques.

The information given in this manual only serves to instruct in the correct handling, cleaning and sterilization of the telescopes. It is not intended as an introduction to the endoscopic technique.

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5 Handling

5.1 Unpacking

Check for missing items and evidence of shipping damage. File any complaints immediately with KARL STORZ or the supplier.

5.2 Application

WARNING: Risk of infection: These instruments are not sterile when delivered. The use of non-sterile instruments poses a risk of infection for patients, users and third parties. Inspect instruments for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the instruments before initial use and before and after every subsequent use using validated procedures.

WARNING: The instruction manuals and interface specifications for the medical devices used in combination must be observed precisely.

WARNING: Risk of injury: Patient leakage currents may accumulate if endoscopes and endotherapy devices supplied with energy are used together. This is particularly important if an endoscope with a Type CF applied part is used. In this case, an endotherapy unit with a Type CF applied part must be used in order to minimize the total patient leakage current.

NOTE: Check that the telescope, the system components used in combination and the overall system are all working correctly prior to use (see Section 5.3).

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In order to examine the image quality, slowly rotate the telescope while looking through it.

If the image is completely or partially impaired, either the rod lens system in the sheath or the lens in the eyepiece is defective. In such cases, the telescope must be replaced.

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A cloudy or speckled image can be caused by moisture or disinfectant/chemical residue on the optical end faces.

In these cases, carefully clean the optical end faces with a cotton tip applicator soaked in 70% alcohol.

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Fiberoptic cable

Hold the distal end in the direction of a light source and inspect the number of dark dots on the light cable connection. These dark dots indicate that the optical fibers of the light conducting bundle are broken. Isolated broken optical fibers do not noticeably compromise the quality of the image.

As of a certain number of broken optical fibers, the light intensity deteriorates considerably. Depending on the application and the combination of light source, light cable and telescope used, the light intensity may then no longer be sufficient. Check that the available light intensity still ensures safe working prior to use. Otherwise, the telescope should no longer be used.

NOTE: If damage is apparent, the telescope should no longer be used. Telescopes must be replaced if the image is cloudy or if no image or only parts thereof are visible.

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Sheath

Inspect the entire sheath for dents and scratches. The sheath must not be bent. LUER-LOCK stopcocks, telescope locks and other movable parts should be smooth-running and undamaged. Check plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their original condition, they must be replaced. This applies in particular to seals.

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5.4 Precautions

As with all optical instruments, the telescopes are a very sensitive part of the endoscopic equipment. In order to ensure the long-term quality of the telescope, the following information must be observed.

- In order to avoid breaks, telescopes must always be held by the eyepiece and never by the distal end alone.
- Do not bend the sheath under any circumstances. This can result in breaks or cracks in the rod lens system.
- Place the telescopes down with care. Hard impacts, particularly to the distal end, can cause damage to or cracks in the luting and allow liquid, steam, etc., to penetrate. This results in cloudy and/or blurred image areas.
- Therefore, during cleaning, disinfection, sterilization and storage, it must be ensured that the telescope cannot be damaged by other instruments. Telescopes are therefore best stored individually, or use containers in which the telescopes can be secured (see catalog entitled HYGIENE – Care, Sterilization and Storage Techniques, Art. no. 96211004DER).

Special storage options that are suitable both for storage and transport as well as for cleaning and sterilization are available for telescopes (see catalog entitled NEUROSURGERY, Art. no. 96020000 DER).

NOTE: If damage is apparent, the telescope should no longer be used. Telescopes must be replaced if the image is cloudy or if no image or only parts thereof are visible.

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5.5 Instructions for use

Before use, make sure that the telescope has been correctly reprocessed (see Section 7).

Check that the telescope, the system components used in combination and the overall system are all working correctly prior to use (see Section 5.3).

Operate the video and light source equipment as per the specifications in the relevant instruction manual.

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5.6 Connecting to the camera system and light source


The KARL STORZ camera heads have an integrated instrument coupler. Squeeze the levers of the instrument coupler together and insert the endoscope's eyepiece into the coupler. The coupler automatically locks into place when the eyepiece is inserted.

NOTE: The eyepiece cup is designed according to ISO/TS 18339. All camera heads with an interface that is compatible with these standards can generally be connected. However, proper functioning must still be checked.

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Connect the light cable to the outlet of the light source and to the endoscope.

NOTE: The light connection is designed according to ISO/TS 18339. All light cables with an interface that is compatible with these standards can generally be connected using appropriate adaptors if necessary (Art. no.: 495 F, 495 G). However, proper functioning must still be checked.

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WARNING: Before each use and after any configuration change, the operator must ensure that the telescope is providing an image in real time and with the correct color reproduction, and that the image is correctly aligned.

WARNING: Risk of injury and risk of damaging the products: The highly concentrated light energy can cause the distal end, the light connection, adjacent components and tissue in front of the light emission window to heat up. This can cause burns to the patient, user and surgical accessories.

- Avoid direct tissue contact with the distal end, as well as with the light connection of the telescope and the light cable.
- Do not rest the endoscope or the light cable on the patient or where they are in direct contact with surgical accessories.
- Always select the lowest possible light output of the cold light source which still allows optimum illumination of the surgical field.

WARNING: Risk of injury: Optical radiation poses a risk of injury to eyes. Never look into the free end of either a connected light cable or endoscope. Always wear standardized protective equipment when using lasers.

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
5.7 Light cable assignment

If there is an identification symbol on the light cable and the endoscope (see table below), assign matching combinations according to the symbols shown.

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5.8 Focus and image size

Operate the video equipment as per the specifications in the relevant instruction manual.

ⓘ NOTE: KARL STORZ HX camera heads do not have an optical zoom.

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
5.9 White balance

Operate the video equipment as per the specifications in the relevant instruction manual. Perform a new white balance as required, or after replacing the telescope.

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
5.10 Assembling the sheath and telescope

Sheaths, sheath inserts and telescopes are provided with a 0 mark or a line mark. Open the telescope lock (1) by turning it counterclockwise (as far as it will go, line marks in alignment (2)). When assembling, ensure that the marks are aligned with one another (line mark (2) and 0 mark (3) are lined up). Insert the telescope into the sheath lock; close the telescope lock (1) by turning it in a clockwise direction.

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6 Reprocessing

**6.1 Overview of reprocessing
procedures**

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- This reprocessing step is part of the validated procedure
- This reprocessing step may be carried out additionally in order to aid cleaning
- ✕ On the grounds of potential damage, this reprocessing step must NEVER be performed
- The process step is not part of the performed validation

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6.2 Overview of important points

WARNING: If CJD/vCJD (Creutzfeldt-Jakob disease) is suspected or indicated, isolate and if necessary dispose of the product in accordance with national regulations and do not use it again.

WARNING: Risk of infection: For instruments belonging to this risk class, the reprocessing requirements are raised. Only the processes mentioned here have been validated. Incorrectly reprocessed instruments expose patients, users and third parties to a risk of infection as well as to the risk of an instrument malfunction.

WARNING: Please note that special process chemicals with prion-inactivating properties must be used for instruments used for neurosurgery or for contact with the central nervous system.

CAUTION: When preparing and using solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiologic range of action of the chemicals used.

CAUTION: Danger of damage to products: Instruments may be damaged by using chemicals that have not been approved by KARL STORZ. Only use chemicals approved by KARL STORZ for reprocessing. A complete list can be found on the Internet at www.karlstorz.com.

CAUTION: Do not use fixing agents when rinsing off (e.g., no alcohol, aldehyde, peracetic acid) and ensure that this is not carried out at temperatures above room temperature.

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CAUTION: Where cleaning agents or disinfectants which have not been approved are used, KARL STORZ cannot provide a guarantee in case of damage.

NOTE: The biological and functional safety of the KARL STORZ instruments used for neurosurgery and for contact with the CNS has been proven for a limited number of cycles by means of reprocessing in accordance with the processes described in this manual (see Section 7). Detailed information on this matter can be found in the relevant instrument-specific manual.

- Observe all Warnings and Notes for your own protection as well as for the protection of the personnel around you.
- Failure to observe this manual poses the risk of damage to the instrument, increased instrument wear and invalidation of the warranty.
- Before brand-new and repaired instruments are sterilized for the first time, they must be subjected to at least one complete reprocessing cycle.
- Instruments which have been used on patients must be taken for cleaning/disinfection without delay in order to avoid the unnecessary drying of contaminants on the instrument, as this may have an extremely negative impact on cleaning in hard-to-reach areas.
- Only use validated procedures/processes for machine cleaning/disinfection.
- The operator is responsible for the on-site validation of the reprocessing processes.
- Only process chemicals (with demonstrable prion-inactivating properties) that have been approved by KARL STORZ with respect to material compatibility and are suitable with regard to the range of action may be used.
- Only use recommended and suitable cleaning accessories (e.g., brushes).

CAUTION: Never use metal brushes or metal sponges for cleaning or repeated cleaning.

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- To avoid unrinsed areas, do not overload washer and disinfectant units (WDUs).
- Microbiologically perfect (germ-free) CD water must be used for the final rinse following cleaning and chemical disinfection. The instruments must be residue-free upon submission for packaging and sterilization.
- Remove instruments immediately after machine cleaning and disinfection in the washer and disinfectant (risk of corrosion).
- Prior to the functional check and specific care measures, each instrument must be visually inspected for purity/cleanliness and, if necessary, subjected to repeated cleaning.
- Only clean, disinfected and dry instruments which have been cared for and are in full working order may be packed and sterilized.
- Drying: The level of dryness of the instruments to be sterilized can have a significant influence on the effectiveness of the sterilization (depending on the method). It must be ensured that the instruments are absolutely dry when submitted for sterilization. This must be carried out whilst taking the various influential factors, such as the load, the geometry of the instruments and management of the process on site, into account.
- Each instrument must undergo a functional check before packaging and sterilization. Disassembled instruments must be assembled for checking and remain in this state for packaging and sterilization. Separate information is provided for instruments that must be sterilized in a disassembled state.
- Defective instruments must be removed from circulation and sent to the manufacturer for repair.

i NOTE: Repairs carried out by third parties may lead to invalidation of the warranty.

- Before being returned to the manufacturer, any instruments which require a repair must – as far as possible – be subjected to a complete reprocessing cycle. The corresponding documentary proof must be enclosed with the delivery.

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NOTE: Disinfection is not assured without cleaning first! Sterilization is not assured without disinfecting first!

- Only the validated procedure may be used for sterilization.

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6.3 Basic information about the reprocessing procedures

As a fundamental principle, always ensure that cleaning is carried out first, and disinfection afterwards in a separate procedure! Rationale for the separation of the cleaning and disinfection processes: Disinfectants may contain fixing components and therefore impede cleaning, which must be considered very critically with regard to a potential prion risk. Thorough cleaning reduces the number of microorganisms and potentially high-risk tissue residues/liquids, and facilitates the disinfection process that follows and makes it more effective.

6.4 Cleaning and disinfection procedures

When selecting the instrument cleaner, the range of action must be considered in relation to its prion-inactivating properties. Highly diverse contaminants such as blood, mucus and bone constituents make specific demands of the process chemicals. The chemical manufacturer's specifications (e.g., concentration, exposure time, temperature) for machine cleaning must be followed (ISO 15883).

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6.4.1 Cleaning, disinfection and sterilization of thermostable instruments and accessories

WARNING: The stated cleaning, disinfection and sterilization procedures alone are not enough to guarantee the required level of hygienic safety for the instruments. This can only be achieved if personnel follow recognized and validated reprocessing methods.

The instruments and accessories for endoscopic surgery are, as a rule, sensitive precision instruments which require particular care during reprocessing, particularly when they are used in contact with the CNS or for neurosurgery to minimize the risk of iatrogenic prion transmission. Only trained personnel (according to §4 of the German Medical Device Operator Ordinance) are in a position to carry out disassembly, cleaning, disinfection, inspection/care, assembly and sterilization procedures in such a way that proper functioning of the instruments is retained.

Some instruments become contaminated on surfaces which cannot be seen and are hard to reach. Therefore, where technically feasible or compatible with the application, these instruments have an irrigation channel or are developed for easy disassembly.

WARNING: When carrying out any work on contaminated instruments, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety (protective equipment to minimize the risk of infection) must be observed. Aerosols must be removed by suitable means, such as a suction device.

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6.5 Sterilization methods

- Steam sterilization in the fractionated prevacuum procedure (ISO 17665)
- The applicability of other procedures based on chemical reagents is subject to country-specific regulations.

6.6 Water and steam quality

For reprocessing in relation to endotoxins, KARL STORZ recommends using a water quality based on Ph. Eur. <0169>.

6.6.1 Water quality for cleaning, disinfection and rinsing

Generally speaking, CD water of an appropriate microbiological quality is used for the final rinse, in accordance with the 'Guideline Compiled by the German Society for Hospital Hygiene, the German Society for Sterile Supply and the Instrument Reprocessing Working Group for the Validation and Routine Monitoring of Automated Cleaning and Thermal Disinfection Processes for Medical Devices as well as Advice on Selecting Washer-Disinfectors' (ISO 15883-1).

Conductivity:	≤ 15 µS
pH value:	5 – 7
Total hardness:	≤ 0.02 mmol CaO/l
Salt content:	≤ 10 mg/l
Phosphate: (as P ₂ O ₅)	≤ 0.5 mg/l
Silica: (as SiO ₂)	≤ 1 mg/l
Chloride:	≤ 2 mg/l

Tab. 1: In accordance with the Guideline, the following values are recommended

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The types of water listed below are free from microorganisms and are suitable for the final rinse:

- Sterile-filtered water
- Sterile water

These types of water are not free from microorganisms and are therefore not suitable for the final rinse:

- Distilled water
- Deionized water
- Utility water
- Aqua purificata

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Tab. 2: EN 285:2009-08 Anhang B

6.6.2 Steam quality for sterilization

For sterilization, water quality must comply with EN 285:2009-08 Appendix B for feed water and condensate. This compliance should be checked regularly according to recognized analytical methods.

Feed water:	Conductivity at 25°C ≤ 5 µS/cm pH value 5 to 7.5
Condensate:	Conductivity at 25°C ≤ 3 µS/cm pH value 5 to 7

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6.7 Handling instruments if prion contamination is suspected/confirmed

The reprocessing of medical devices, e.g., when (variant) Creutzfeldt-Jakob disease ((v)CJD) is suspected, must be based on the respective applicable national regulations. In clinically verified cases, the medical devices used must generally be destroyed.

In Germany, at least, a regulation stipulates that, in suspected cases, the medical devices to be reprocessed must be labeled and stored in isolation until the diagnosis is clarified (according to the current recommendation set down by the Commission for Hospital Hygiene and Infectious Disease Prevention/the German Federal Institute for Drugs and Medical Devices*). Where the suspicion is unfounded, conventional reprocessing as described in Section 6 is performed; where the suspicion is confirmed, the instrument set used must be incinerated.

*: Hygiene requirements during the reprocessing of medical devices (2012)

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6.8 Carrying out reprocessing

The cleaning, disinfection and sterilization procedures that are based on the guidelines set down by the WHO as regards minimizing the risk of iatrogenic prion transmission through surgical instruments used for neurosurgery or contact with the CNS are mentioned here.

NOTE: Machine cleaning/thermal disinfection must be carried out, followed by fractionated steam sterilization at 134-137°C with a hold time of at least 18 mins.

Instruments used for contact with the CNS or for neurosurgery must be reprocessed separately and should not be mixed with instruments from other specialist departments/disciplines.

Proper and effective cleaning is essential to the success of subsequent disinfection and sterilization. It must therefore be carried out with due care, particularly with regard to minimizing the risk of iatrogenic prion transmission through surgical instruments.

Efficient cleaning is the most effective way of minimizing the risk of infection.

Personnel must have the appropriate training and specialist knowledge required for reprocessing. Personnel must have in-depth knowledge of the instruments to enable them to carry out the necessary disassembly and care procedures as well as the functional check. The instrument-specific instruction manuals must be considered in detail on a case-by-case basis.

Machine processes must be validated and used.

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Personnel wear nationally defined protective equipment (PPE). Other protective measures, such as structural requirements, must be observed on a country-specific basis.

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Instrument-specific work steps, processes and parameters relating to reprocessing are included in the relevant accompanying instruction leaflet and must be observed separately.

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6.9 Preparation for cleaning/ disinfection at the site of use

The following basic measures help to ensure that the value and functionality of the instruments are maintained. They make reprocessing significantly easier.

- Heavy soiling must be removed immediately after use whilst still in the operating/procedure room, e.g., by wiping off external contaminants and rinsing out working channels before placing (not 'throwing') into the disposal container.
- Instruments which are found to be defective during use must be visibly labeled and subjected to reprocessing as well before sending them to the manufacturer for repair.

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* in Anlehnung an die Infection Control Guidelines for Transmissible Spongiform Encephalopathies der World Health Organization (WHO)

6.10 Transportation to the reprocessing site

- i** **NOTE:** The instruments must be kept moist for transportation to the reprocessing site in order to prevent or reduce the chance of potentially prion-loaded, organic residues becoming dried on as far as possible.*
- i** **NOTE:** If such residues do dry on, the effectiveness of the validated processes cannot be guaranteed.
- i** **NOTE:** The instruments should not be mixed with instruments from other applications/disciplines.
- i** **NOTE:** In order to ensure safe transportation and storage, specific holding systems are available for delicate instruments (see catalog HYGIENE – Care, Sterilization and Storage Techniques, Art. no. 96211004DER).

6.11 Disassembly and preparation for cleaning/disinfection at the reprocessing site

After delivery of the instruments, unload the transport system immediately.

- !** **WARNING:** Use splash protection for eyes, mouth and nose as well as the body (PPE). Aerosols must be removed by suitable means, such as a suction device.

6.11.1 Disassembly

Prior to cleaning and disinfection, the telescope must be separated into its individual components as far as possible. Both adaptors on the light inlet piece and any available adaptors, seals on the working channel and LUER-LOCK stopcocks must be disassembled.

* Based on the Infection Control Guidelines for Transmissible Spongiform Encephalopathies set down by the World Health Organization (WHO)

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6.11.2 Manual pre-cleaning Brushing the surfaces

Any visible contamination or heavy soiling of the surfaces must be removed by cleaning them under cold running water with the aid of a brush (Art. no. 27652).

No pretreatment with ultrasound

For technical reasons, the medical device is not suitable for ultrasound treatment.

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6.12 Machine cleaning and disinfection

Only washers and disinfectors that have been qualified and validated in accordance with EN ISO 15883 may be used. As a rule, the most important processes (standardized programs) are preinstalled in the washer and disinfectant by the unit manufacturer. The advantage of machine cleaning/disinfection lies in the reduction of risk to personnel, the standardization of the processes and the creation of a basis for a validated reprocessing process.

NOTE: The temperature of the water flowing into the unit should not exceed 45°C; higher temperatures lead to the coagulation of protein and may thereby impair the cleaning performance.

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We recommend using completely demineralized water of an appropriate microbiological quality in all stages of the program.

i NOTE: For reprocessing in relation to endotoxins, KARL STORZ recommends using a water quality based on Ph. Eur. <0169>.

The selection of the process chemicals and their ratios must also be observed precisely. The use of care and rinse agents in the machine is discouraged, as deposits may form on the instruments or plastics may be damaged by the agents.

The A₀ concept (see ISO 15883) must be taken into consideration in the on-site process validation carried out by the user.

During the cleaning process, it must be ensured that cleaning is only carried out using a demonstrably prion-activating (prion-destabilizing/prion-inactivating/prion-decontaminating) cleaning agent that is approved for this purpose. It is preferable to carry out this cleaning work in an alkaline environment with a pH value > 10 (in accordance with the current recommendation set down by the Commission for Hospital Hygiene and Infectious Disease Prevention/the Federal Institute for Drugs and Medical Devices in Germany).

Thermal disinfection is carried out at temperatures of at least 93°C and with an exposure time of 5 minutes which is appropriate for reaching the required A₀ 3000 value and for achieving a level of bacteria reduction of at least 4 log levels.

6.13 Connecting telescopes to the channel

In order to guarantee effective machine cleaning and disinfection, telescopes must be connected to the channel in the washer and disinfectant, or it must be ensured that it is thoroughly rinsed out by means of appropriate rinsing lances (see the section entitled "Recommended storage and reprocessing accessories").

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6.14 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Dismantled medical devices must be assembled.
- Afterwards, a functional check must be carried out.

Removing deposits

Due to foreign particles (inert gases) in the steam, over time stubborn deposits may form on the cover slips during autoclaving. These deposits can be easily removed using the cleaning paste (Art. no. 27661).

Put a little cleaning paste on a clean, moistened cleaning rod and apply to the cover window. Polishing removes stubborn deposits which can be subsequently rinsed off with water.

In order to achieve a smooth surface, polishing wool (Art. no. 27658) can be used if necessary.

WARNING: Following cleaning, a complete reprocessing cycle must be carried out.

CAUTION: The optionally available special cleaning paste should only be used for cleaning when the endoscopic image is visibly cloudy; under no circumstances should the paste be used on a regular basis.

Inspecting the optical system

Examine the optical end faces (distal tip and eyepiece) for scratches and surgical or cleaning residues. The faces should be smooth and shiny. In order to examine the image quality, slowly rotate the telescope while looking through it. If the image is completely or partially impaired, either the rod lens system in the sheath or the lens in the eyepiece is defective. In such cases, the telescope must be replaced.

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A cloudy or speckled image can be caused by moisture or disinfectant/chemical residue on the optical end faces. In these cases, carefully clean the optical end faces with a cotton tip applicator soaked in 70% alcohol (Fig. A).

Inspecting the fiberoptic cables

Hold the distal end in the direction of a light source and inspect the number of dark dots on the light cable connection. These dark dots indicate that the optical fibers of the light conducting bundle are broken. Isolated broken optical fibers do not noticeably compromise the quality of the image. However, the light intensity does significantly deteriorate if a certain number are broken (from approx. 20% to 25%), which therefore also results in reduced image quality. The telescope should then no longer be used.

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Manual (post-)drying

Dry surfaces with medical-grade compressed air for min. 5 seconds and at 2.5 bar until no more residual moisture appears or is visible.

The fiber areas in the light inlet piece in particular as well as the distal and proximal end faces must be post-dried carefully with 70% alcohol. Disinfectant and cleaning agent residue, e.g., in the light inlet piece, can become burnt on when the light cable is connected, thereby seriously impairing light transmission.

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NOTE: During care procedures, use items from the catalog 'HYGIENE – Care, Sterilization and Storage Techniques' (Art. no. 96211004DER).

6.15 Packaging and labeling

Only standardized packaging materials and systems may be used. Containers, for example, are considered rigid packaging and fleece sheets or film/paper, for example, are considered soft packaging*.

General, more detailed information about the 'Packaging and labeling' section can be found in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no.: 96216003DER).

*As part of sterilization validation, the following packaging material (clear paper/film pouch – ISO 11607, CE with STEAM-, EO and FORM indicator, stericlean® VP Group, Art. no. 4FKFSxxxxxx) was used double-packed.

6.16 Steam sterilization

Steam sterilization with an 18-minute exposure time at 134-137°C is a partially suitable risk minimization method according to current scientific knowledge. It must be used in conjunction with other procedures that are rated as more than just partially effective, such as the cleaning/disinfection process step.

NOTE: National laws and regulations must be observed.

6.16.1 Steam sterilization in the prefractionated prevacuum

The following standards must be observed to meet the requirements placed on the steam sterilizer with fractionated prevacuum and to correctly carry out sterilization using saturated steam:

- EN 285 Large steam sterilizers (from 1 StE upwards) – requirements and test procedures
- EN 13060 Small steam sterilizers
- ISO 17665-1 for routine monitoring

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For the purpose of prion prevention, these telescopes must be sterilized using the fractionated prevacuum procedure (ISO 17665-1) at 134°C-137°C with a minimum exposure time of 18 minutes.

NOTE: Greased components must be disassembled for sterilization, so that steam penetration is ensured.

Chronological pressure and temperature profile of a fractionated prevacuum procedure

7 Validated number of reprocessing cycles

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application. The biological and functional safety of the HOPKINS® telescopes has been proven for 50 cycles by means of reprocessing in accordance with the processes described in this manual.

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8 Operation and storage conditions

	Temperature	Rel. humidity
Storage	-20°C ... +60°C	10% ... 90%
Application	+10°C ... +40°C	30% ... 75%
Air pressure	700 hPa – 1060 hPa	

9 Servicing and repair

Defective telescopes, instruments and devices must be serviced and repaired exclusively by persons authorized by us; all work must employ original parts.

10 Limitation of liability

KARL STORZ SE & Co. KG shall be liable for the failure or deterioration in the safe operation, operational reliability, and performance of the this telescope, instrument or device only subject to the conditions that all assembly, operations, system expansions, readjustments, modifications or repairs to the same have been performed by persons authorized by KARL STORZ, and the instrument has been used in accordance with its operating instructions at all times.

11 Warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see the section 'Subsidiaries'), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall invalidate any warranty.

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12 Directive compliance

This Medical Device bears a CE mark in accordance with the Medical Devices Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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In Germany, you can refer repairs direct to:

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14 Important information

For hygiene reasons and to prevent infection, telescopes, instruments and devices must be cleaned and sterilized/disinfected before they are sent for repair. We reserve the right to return contaminated instruments/devices to the sender. Repairs, modifications, or expansions which are not performed by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

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KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 56381188, Fax: +86 10 56381199
E-Mail: info@karlstorz.com.cn

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KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Jinjiang Road, Wuhou District,
Chengdu, 610041, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenyhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

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1	Important information	1
2	Symbols employed	3
3	Intended use	6
3.1	Indications	6
3.2	Contraindications.....	7
4	User qualification.....	8
5	Warnings and cautions	8
6	Safety instructions.....	9
6.1	Unpacking	18
6.2	Basic equipment.....	19
6.3	Operating conditions	20
7	External system components.....	20
8	Preparation before use	21
9	Handling	21
9.1	Operation	21
9.2	Holding the device correctly	21
9.3	Operating the suction valve	22
9.4	Handling the sheath.....	22
9.5	Handling the deflection mechanism	23
9.6	Handling the working channel.....	24
9.7	Handling the accessories (forceps, cleaning brushes, HF electrodes etc.)	26
9.8	Inspecting the surfaces.....	28
9.9	Inspecting the deflection.....	29
9.10	Connecting to IMAGE1 S™	29
9.11	Focus	29
9.12	Function buttons of the endoscope	30
9.13	Switching the flexible cystoscope HD VIEW ON/OFF	31
10	Reprocessing	32
10.1	General warnings.....	32
10.2	Accessories	33
10.3	Precleaning directly after use	34
10.4	Leakage test.....	35
10.5	Disassembly	37

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10.6	Manual disinfection	40
10.6.1	Manual High Level Disinfection (HLD)	40
10.6.2	Drying	41
10.7	Machine cleaning and disinfection	41
10.7.1	Manual preparation	41
10.7.2	Connecting	42
10.7.3	Machine cleaning/chemical-thermal disinfection	42
10.8	Assembly, inspection, care and packaging	43
10.8.1	Packaging	44
10.9	Sterilization	44
10.10	Limits of reprocessing	46
10.11	Storage	46
11	Reprocessing of LUER port with stopcocks	48
11.1	Accessories	49
11.2	Preparation for cleaning and disinfection	49
11.3	Manual precleaning	49
11.4	Machine cleaning and disinfection	50
11.5	Connecting	51
11.6	Assembly, inspection and care	51
11.7	Packaging systems	51
11.8	Sterilization	52
11.9	Limits of reprocessing	52
12	Technical data	53
13	Servicing and repair	55
14	Maintenance	55
15	Limitation of liability	55
16	Warranty	55
17	Directive compliance	56
18	Important information	56
19	Disposal	56
20	Electromagnetic Compatibility user information	57

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This instruction manual is intended to serve as an aid in the proper handling, cleaning, disinfection and, if appropriate, sterilization of the flexible cystoscope HD VIEW. All the necessary details and actions are clearly explained. We therefore ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for easy reference.

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2 Symbols employed

	Manufacturer
	Consult instructions for use
	Catalogue number
	Number of products in the product packaging
	Serial number
	Non-sterile
	Fragile, handle with care
	Keep dry
	Temperature limit
	Federal (USA) law restricts this device to sale by or on the order of a physician.

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- ① Connector
- ② Suction valve
- ③ Connection for leakage tester or pressure compensation cap (vent port)
- ④ Handle
- ⑤ Irrigation and working channel
- ⑥ Deflection lever
- ⑦ Distal tip
- ⑧ Working channel connector seal 11014 DS
- ⑨ Leakage tester 13242 XL
- ⑩ Suction valve, sterile, for single use 091011-20
- ⑪ Suction valve, reusable 11301 CB1
- ⑫ Pressure compensation cap 11025 E
- ⑬ LUER port with stopcocks 11014 L
- ⑭ Double LUER port 11014 LU
- ⑮ Cleaning adaptor for working channel 11014 RA
- ⑯ Cleaning adaptor 11301 CD1
- ⑰ Right navigation button
- ⑱ Left navigation button
- ⑲ Menu button

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3 Intended use

Flexible cysto-urethrosopes (fiber endoscopes and videoscopes) are used for endoscopic imaging of the urethra and the bladder for the purposes of diagnosis and treatment. Flexible cysto-urethrosopes with an integrated working channel allow for the use of flexible auxiliary instruments such as forceps, stone baskets etc., coagulation electrodes and laser fibers for lithotripsy or for the ablation of benign/malignant tissue. With flexible cysto-urethrosopes, the urethra and the entirety of the bladder can be visualized/treated. The special feature of the flexible cysto-urethroscope lies in the fact that it can be used percutaneously in the kidneys or suprapubically in the bladder via an operating sheath.

Flexible cysto-urethrosopes (fiber endoscopes and videoscopes) are designed for temporary or short-term use in endoscopic interventions performed through natural orifices or surgical invasive interventions.

3.1 Indications

Description of all known indications:

The use of cysto-urethrosopes is indicated if urethro-cystoscopy is indicated in the opinion of the attending physician.

This includes, amongst other things:

- Clarification of micro and macro hematuria
- Clarification/biopsy in the event of suspected bladder tumors
- Control urethro-cystoscopy as per TUR-B or P
- Clarification in particular in the event of suspected tumor infiltration due to adjacent organ cancer (e.g., colon carcinoma, cervical carcinoma)
- Clarification or treatment in particular in the event of suspected vesicoenteric or obstetric fistulas
- Urethral or vesical foreign body removal
- Clarification of recurring urinary tract infections
- Clarification of LUTS/BES, subvesical obstructions, BOO, neurogenic micturition disorder

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Intended use

- Visual assessment and possibly treatment (e.g., injection) of the orifices of the ureters (localization, condition, etc.)
- Visual clarification of the bladder sphincter in the event of incontinence
- Visual assessment and possibly treatment of the urethra (tumor, stricture/stenosis, fistula)
- Visual assistance for indwelling catheter system, e.g., with urethral trauma in the event of a pelvic fracture
- Retrograde pyelography, ureteral splinting

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3.2 Contraindications

The use of the cysto-urethrosopes is contraindicated if, in the opinion of the attending physician, the surgical or examination method as such is contraindicated or if the patient is not able to undergo surgery or anesthesia due to his or her general condition.

Cysto-urethrosopes must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

The following contraindications also apply:

- Acute urethritis
- Acute epididymitis
- Acute prostatitis
- Acute/florid cystitis/urinary tract infection
- Possible coagulation disorders
- Restless, agitated patients

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4 User qualification

The flexible cystoscope HD VIEW from KARL STORZ may only be used by persons with an appropriate medical qualification and who are acquainted with the flexible cystoscopy technique. The information given in these instructions only serves to instruct in the correct handling and reprocessing of the flexible cystoscope HD VIEW. It is not intended as an introduction to the techniques of flexible cystoscopy.

5 Warnings and cautions

HAZARD SYMBOL: This is a hazard symbol. It warns of the risk of injury. To prevent injuries, all the measures marked with the hazard symbol must be followed.

Signal words:

WARNING! designates a possible imminent risk. If this is not avoided, it **could** lead to death or serious injuries.

ATTENTION! designates a possibly harmful situation. If it is not avoided, the medical devices could be damaged.

NOTE! indicates special information on the operation of the device or important information.

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6 Safety instructions

The flexible cystoscope HD VIEW must only be used according to the medical rules and procedures of endoscopy recognized for this method.

WARNING:

Risk of injury and damage to the product:

- It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces or voltage-carrying elements of other devices at any time.

WARNING:

Risk of burns and glare: The LED light emitted through the flexible videoscope may cause glare or high temperatures, and thereby burns, at the distal light outlet.

- Never look directly into the light emitted by the flexible videoscope.

WARNING:

Risk of explosion: If flammable anesthetics are used in the immediate proximity of the flexible videoscope.

- The flexible videoscope must not be used in the presence of flammable or combustible media.

WARNING:

Risk of injury: When in operation, the integrated LED light source heats up the flexible videoscope.

- Switch off the flexible videoscope before putting it down, and do not place it in the vicinity of the patient (as per chapter 16 of the 3rd edition of IEC 60601-1).

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⚠ WARNING:

Risk of injury: High energy light radiated through endoscopes may give rise to high temperatures in front of the light outlet, at the end faces of the fiberoptic light cable and at the tip of the endoscope.

- *Prolonged direct contact of light-emitting surfaces on tissue can lead to burns; this must therefore be avoided at all costs.*

⚠ WARNING:

Risk of infection: The flexible videoscopes are not sterile when delivered, and must be cleaned and disinfected/sterilized prior to initial use and each subsequent reuse.

⚠ WARNING:

Risk of injection: To avoid infection and maintain functionality, the chapter on reprocessing in the instruction manual must be observed. The cleaning, disinfection and sterilization procedures are explained therein in detail.

⚠ WARNING:

Risk of infection: Before the device is used for the first time, it must be subjected to a complete reprocessing cycle, excluding precleaning.

⚠ WARNING:

Risk from electromagnetic incompatibility: Medical electrical devices are subject to special precautions regarding Electromagnetic Compatibility (EMC). Observe the stated EMC instructions during installation and operation.

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⚠️ **WARNING:**

The use of accessories other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased interference immunity of the flexible videoscope. The accessories listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using accessories other than those specified here, it is the responsibility of the operator to ensure that they comply with EN/IEC 60601-1-2.

⚠️ **WARNING:**

Risk of injury: The videoscope does not offer any electric insulation against high frequency (HF) voltages. Protection against HF voltage must be achieved using the appropriate HF instrument set.

➤ Please read the instructions for use in the instruction leaflet supplied with the HF application part for more information.

⚠️ **WARNING:**

Combinations of medical devices are only assured to be safe if they are listed as such in the respective instruction manuals or the intended use and interface specifications of the products used in the combination permit this.

➤ The instruction manuals and interface specifications of the medical devices used in combination must be observed precisely.

⚠️ **WARNING:**

The pressure compensation cap must be removed from the vent port for:

- endoscopic interventions,
- immersing in liquid, e.g., disinfectant solution and for machine reprocessing

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Safety instructions

⚠️ **WARNING:**

The pressure compensation cap must be in place on the vent port for:

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⚠️ **WARNING:**

Risk of injury: The flexible videoscope must not be used during the discharge of a cardiac defibrillator.

- Remove the flexible videoscope from the operation area before discharge takes place.

⚠️ **WARNING:**

Risk of injury: Only KARL STORZ devices may be used in combination with the flexible videoscope. Disregarding this warning can endanger the patient and user.

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ATTENTION: Do not allow any live components or instruments to come into contact with the flexible videoscope.

NOTE: Do not hit the distal tip of the flexible videoscope against hard objects, and never use a sharp object to remove dirt.

NOTE: Never place heavy objects on the flexible videoscope.

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Safety instructions

NOTE: Never clamp or fold the sheath of the flexible videoscope. This can result in damage to the videoscope.

NOTE: Constant changing between various reprocessing methods is especially hard on the materials and must therefore be avoided.

NOTE: The flexible videoscope must not be carried by the sheath. Do not pull, clamp or twist the sheath of the videoscope.

NOTE: The sheath of the flexible videoscope must never be clamped or coiled up too tightly. It is best to store the video videoscope in a suspended position or to keep it sterilized in a tray specifically designed for sterilizing and storage (39406AS).

WARNING:

Risk of injury from clouding of the lens.

- Do not apply lubricant to the lens, to avoid clouding the image and thereby prevent risks for the patient!

WARNING:

When inserting the flexible videoscope into the patient, the distal tip of the flexible videoscope must be in the straight, non-deflected position.

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 WARNING:

Risk of injury: Before inserting or removing the flexible videoscope, the distal tip must be put into the neutral position, as otherwise injuries to the patient or damage to the instrument may occur.

 WARNING:

Risk of injury: If the distal tip can no longer be properly controlled (sticking etc.) via the deflection lever, or if the deflection mechanism develops any other fault, the flexible videoscope must not be used.

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- Place the distal tip of the flexible videoscope in the neutral, non-deflected position and remove the flexible videoscope slowly and carefully from the patient.

NOTE: The resistance that can be felt when deflecting the distal tip using the deflection lever increases as the angle of deflection increases. Actuate the deflection gently and carefully, without any sudden or jerky movements, as far as the limit stop.

 WARNING:

Risk of injury: The flexible videoscopes from KARL STORZ are not intended for use with light sources. The use of non-compatible devices can lead to patient injury or device damage.

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WARNING:

Risk of injury: Whenever the flexible videoscope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition. Damaged flexible videoscopes or damaged accessories must not be used.

NOTE: There should be no instruments in the suction channel whilst suctioning is being performed. The instrument channel inlet should be sealed to prevent additional air from being drawn in.

NOTE: Do not squeeze the distal tip of the flexible videoscope too tightly.

WARNING:

Risk of injury: Actuate the deflection lever gently and carefully as far as the limit stop. The deflection lever must never be operated abruptly or with force, to avoid injuring the patient or damaging the device.

- Move the deflection lever smoothly and slowly!

WARNING:

Risk of injury: Before each use, the functionality of the deflection lever and the connected distal tip must be checked.

- Defective flexible videoscopes must not be used.

WARNING:

Risk of injury: If the image fails or there is interference during use in the patient, put the distal tip into the neutral position and remove the device from the patient's body gently and carefully. If the flexible videoscope continues to be used in the patient after the image has failed or when there is interference, this may result in serious injury to the patient.

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Safety instructions

WARNING:

Risk of injury: Only insert suitable accessories, e.g., biopsy forceps, into the instrument channel of the videoscope when the distal tip is in a neutral, non-deflected position.

- > Once the instrument has been introduced, carefully actuate the lever gently as far as the limit stop.
- > The deflection lever must never be operated abruptly or with force, to avoid injuring the patient or damaging the device.

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NOTE: Do not apply force to insert accessories, e.g., biopsy forceps, through the working channel, to prevent damaging the device. It should be possible to move the accessory easily, and without any hindrance, through the entire channel.

NOTE: When the distal tip is deflected, do not insert instruments into the videoscope. Only insert the instrument when the distal tip is in the neutral position and deflect the deflection lever gently and carefully up to the limit stop.

NOTE: Inserting sharp-tipped accessories into the instrument channel may damage the videoscope.

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⚠️ WARNING:

Risk of injury: Foreign bodies and items to be extracted must not be removed via the instrument channel with the aid of instruments. Instead, the entire videoscope must be taken out of the body and the foreign body must then be removed.

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NOTE: Accessories such as biopsy forceps that are kinked or bent will not function properly, are difficult to guide through the channel and may cause severe damage to it. They must therefore be replaced.

⚠️ WARNING:

Ensure that the working channel is completely dry prior to sterilization to avoid damage to the working channel.

⚠️ WARNING:

Always unplug the flexible videoscope before carrying out any maintenance and cleaning procedures on it.

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6.2 Basic equipment

Flexible cystoscope HD VIEW 11272 VHK/VHU with following accessories:

Set 11272VHK

11272 VH	Flexible cystoscope HD VIEW
11301 P	Seal for suction valve seat
11014 L	LUER port with stopcocks
11014 DS	Working channel connector seal
11014 RA	Cleaning adaptor for working channel
11301 CD1	Irrigation adaptor

Set 11272VHUK

11272 VHU	Flexible cystoscope HD VIEW
11301 P	Seal for suction valve seat
11014 LU	Double LUER port
11014 DS	Working channel connector seal
11014 RA	Cleaning adaptor for working channel
11301 CD1	Irrigation adaptor

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Optional accessories

1 x grasping forceps for small fragments 27023 FE
1 x biopsy forceps 27023 ZE
1 x stone basket 27023 KF
1 x fiber fixation kit 11014 YK
1 x coagulating electrode 4 Fr. 27723 T
1 x seal for instrument ports 27550C-10
1 x LUER adaptor with seal 27014 Y
1 x endoscopic seal 100010-10
1 x cleaning adaptor for instrument ports 27001 RA
1x single-use cleaning brushes 110940-50
1 x cleaning brush 27651 B

6.3 Operating conditions

Flexible cystoscope HD VIEW:

Ambient temperature: 15°C ... 30°C
Humidity: 15% – 80%

7 External system components

Further information on external system components can be found in the IMAGE1 S™ (96206286) instruction manual.

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
8 Preparation before use

1. Assemble LUER attachment. To do this, insert the lug of the LUER attachment into the recess on the working channel. Push and turn clockwise until it audibly clicks into place.

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2. Assemble the suction valve/valve seal.

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9 Handling

NOTE: Any damage to the flexible cystoscope HD VIEW resulting from incorrect operation is not covered by the manufacturer's warranty.

9.1 Operation

The following instructions outline the correct operation of the flexible cystoscope HD VIEW and the corresponding accessories. This is not an introduction to the techniques of cystoscopy. Instructions on endoscopic technique can be found in the relevant medical literature.

9.2 Holding the device correctly

Grip the handle of the flexible cystoscope HD VIEW in your left or right hand, as shown in the picture. With this kind of control grip, the flexible cystoscope HD VIEW can be operated with one hand.

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9.3 Operating the suction valve

When the suction valve is assembled, a suction pump must be connected to the suction valve via a tube. Suction is activated upon pressing the valve button.

9.4 Handling the sheath

ATTENTION: Risk of damage to the sheath:

- Do not carry the flexible cystoscope HD VIEW by the sheath.
- Do not pull, clamp or twist the sheath of the flexible cystoscope HD VIEW.
- The sheath of the flexible cystoscope HD VIEW must not be kinked or coiled up too tightly.
- It is best to store the flexible cystoscope HD VIEW in a suspended position or to keep it sterilized in a tray specifically designed for sterilizing and storage (39403 AS).

WARNING:

Risk of injury: Do not apply lubricant to the lens, to avoid clouding the image and thereby prevent risks for the patient!

If, during use, the image fails, is disturbed, or becomes blurred, or if vision in general is impeded, put the deflection mechanism into the neutral position and slowly pull the sheath of the flexible cystoscope HD VIEW out of the patient's body.

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Deflection functions as follows with type 11272 VH endoscopes:

If the lever is pushed upwards, the tip is lowered.

If the lever is pushed downwards, the tip is deflected upwards.

This mechanism is inverted with model 11272 VHU.

ATTENTION: When the flexible videoscope is inserted into the patient, the distal tip must be in a straight, non-deflected position.

⚠️ WARNUNG:

The deflection lever must never be operated with force or abruptly to avoid injuring the patient or damaging the device!

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⚠️ WARNUNG:

Risk of injury and damage from faulty deflection mechanisms: If the distal tip can no longer be properly controlled (sticking etc.) via the deflection lever, or if the deflection mechanism develops any other fault, the flexible videoscope must not be used.

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⚠️ WARNING:

Risk of injury: If the flexible videoscope malfunctions during use on the patient, discontinue use immediately.

- *Place the distal tip of the flexible videoscope in the neutral, non-deflected position and remove slowly and carefully from the patient.*

⚠️ WARNING:

Risk of damage: The lens at the distal end of the flexible videoscope may become detached, e.g., due to mechanical forces.

- *For this reason, it is essential to check that the lens is undamaged and present before and after every decontamination.*

9.6 Handling the working channel

Do not force instruments or cleaning brushes through the channels. This can lead to damage to the working channel, especially if the tip of the instrument is deflected.

When connecting and operating suction/irrigation devices, pay attention to the instruction manuals for the particular devices. The irrigation devices recommended by KARL STORZ are set to a maximum pressure of 54 kpa (400 mmHg/0.54 bar/7.8 psi).

- The LUER adaptor 27014 Y or the seal 27550 C can be used to seal the instrument channel when inserting instruments.
- Instruments can be introduced through the instrument channel of the flexible cystoscope HD VIEW. Make sure that the accessories are of the appropriate length and diameter. Unsuitable accessories can damage the flexible cystoscopes HD VIEW. If you have any questions concerning accessories, please contact KARL STORZ directly or an authorized specialist dealer.

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- The flexible cystoscope HD VIEW must be in the straight, non-deflected position when inserting an instrument into the working channel.

Forceps, scissors and other accessories whose jaw can be opened or closed must always be inserted into the instrument channel with closed jaws. The jaws may only be opened when the distal tip of the instrument is completely clear of the channel. Opening the jaws in the instrument channel can cause damage to the videoscope.

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ATTENTION: Risk of damage: Only use accessories (guide probes, cleaning brushes, etc.) from KARL STORZ. Accessories from other manufacturers may damage the videoscope.

9.7 Handling the accessories (forceps, cleaning brushes, HF electrodes etc.)

WARNING:

Risk of injury: When HF electrodes or laser probes are used, it must be ensured that power is not applied until the active parts have left the instrument channel and are in the user's field of vision, and the site of the operation is clearly recognizable. Activation in the instrument channel may lead to serious injury to the patient and damage to the videoscope.

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ATTENTION: Risk of damage: Only use
accessories (guide probes, cleaning brushes,
etc.) from KARL STORZ. Accessories from other
manufacturers may damage the videoscope.

NOTE: Do not apply force to insert accessories,
e.g., biopsy forceps, through the working channel,
to prevent damaging the device. It should be
possible to move the accessory easily, and without
any hindrance, through the entire channel.

NOTE: When the distal tip is deflected, do not
insert instruments into the videoscope. Only insert
the instrument when the distal tip is in the neutral
position and deflect the deflection lever gently and
carefully up to the limit stop.

NOTE: Inserting sharp-tipped accessories into the
instrument channel may damage the videoscope.

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9.8 Inspecting the surfaces

⚠ WARNING:

Risk of injury: The surface of the flexible sheath must be inspected visually and by carefully feeling it, to detect any cracks, cuts, dents, foreign bodies or other forms of damage.

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Before each use, verify that the device is in good condition by checking the surface of the flexible cystoscope HD VIEW, especially the sheath area and the distal tip, for sharp edges. Cystoscopes HD VIEW with sharp edges must not be used.

Otherwise they may cause serious injury to the patient or user.

Check the connecting cable for damage (such as fractures, cracks, twisted or crushed areas, etc.) and the connecting plug for corrosion.

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⚠️ WARNING:

Risk from damage: Faulty flexible videoscopes must not be used.

- Before each use, check the functionality of the deflection lever and the distal tip.

The distal tip is deflected by actuating the deflection lever slowly.

The distal end will deflect in accordance with the direction of actuation (forward or backward).

9.10 Connecting to IMAGE1 S™

The flexible cystoscope HD VIEW is connected to the camera control unit for operation.

For operation, an IMAGE1 S™ system with an IMAGE1 S Connect® and an IMAGE1 X-Link are required. Both camera modules must be connected to each other.

The flexible cystoscope HD VIEW is then connected to the IMAGE 1X-Link.

9.11 Focus

The image focus is adjusted automatically. Manual focusing is not possible/necessary.

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10.1 General warnings

WARNING:

Risk of infection: Only reprocess the medical device as per the documents provided with the device or the reprocessing instructions issued by KARL STORZ. These can be downloaded from www.karlstorz.com.

WARNING:

Risk of infection from non-sterile delivery. Reprocess the device before each use. Check the medical device for visible soiling before use.

WARNING:

Before reusing the device, replace any defective parts or send the flexible endoscope for repair.

ATTENTION: Risk from damage to the device: The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.

ATTENTION: Risk from damage to the device: Never clean flexible endoscopes with metal brushes.

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ATTENTION: Risk from damage to the device: Flexible endoscopes must never be steam sterilized (autoclaved).

ATTENTION: Risk from damage to the device: To prevent damage to the flexible endoscope and accessories, it must be ensured that all connections are completely dry prior to use.

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NOTE: Damage to the device. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be found at www.karlstorz.com

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NOTE: To avoid damage from the penetration of liquid, a test in the dry state must be performed in advance.

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Immerse the device in solution so that it is completely covered. It must be ensured that the leakage tester does not become wet. Perform the leakage test as described in 'Leakage test – dry' and note whether bubbles emerge along the device. If, after an observation period of at least 3 minutes, no bubbles can be seen emerging (continuously) from the device and, in addition, the manometer does not show a fall in pressure, the device is leakproof.

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Cleaning can then begin.

If a pressure loss occurs, this means that the endoscope leaks and may no longer be used, immersed in liquid or reprocessed by machine.

ATTENTION: In this case, fill out the accompanying repairs form (see Appendix or Art. no. 96096004DER), mark the endoscope, which has been cleaned by wiping to remove organic soiling/residue, as not disinfected and return it to the Service Department for repair or, alternatively, contact KARL STORZ. Line the storage/carrying case with the transport protection sheet 13990 SFN.

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10.5 Disassembly

The flexible cystoscope HD VIEW must be completely dismantled for cleaning/disinfection.

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1. Dismantle the suction valve by pressing and turning.

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2. Dismantle the LUER attachment by pressing and turning

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3. Dismantle the stopcock (if present) on the LUER attachment.

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5. Clean the channels and lumina of the LUER attachment with a suitable brush.

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6. Dismantle the reusable suction valve 11301 CB1 by unscrewing the metal cap.

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8. Brush the working channel connector.

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9. Brush the working channel and suction channel.

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10. Connect the cleaning adaptor 11301 CD1 to the valve seat.

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10.6 Manual disinfection

The endoscope must be completely immersed in a cleaning solution. The suction and working channel must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

10.6.1 Manual High Level Disinfection (HLD)

The endoscope must be completely immersed in a high level disinfectant solution. The working and suction channel must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

The endoscope can be disinfected with the following high level disinfection solutions which have been tested by way of example:

- 2.4% glutaraldehyde solution (e.g., CIDEX® Activated Dialdehyde Solution)
- 0.55% ortho-phthalaldehyde solution (e.g., CIDEX® OPA Solution)
- 2.0% accelerated hydrogen peroxide solution (e.g., Revital-Ox® Resert™ or Revital-Ox Resert XL HLD)

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ATTENTION: The pressure must not exceed 100 kPa (1 bar/14.5 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.

10.7 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003DER):

10.7.1 Manual preparation

The endoscope must be completely immersed in a cleaning solution. The working and suction channel must be covered with no air bubbles by filling in a targeted manner. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a suitable brush. The channels are brushed until the brush head emerges and is visibly clean. Before retracting the brush head that emerges at the distal end, clean it. Never brush back and forth within the channel.

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NOTE: KARL STORZ recommends using the adaptors 11301CD1 and 11014RA for simpler filling of the lumina.

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10.7.2 Connecting

In order to ensure effective machine cleaning and disinfection, the endoscope's working and suction channel must be connected to the washer and disinfectant using the recommended adaptors 11301 CD1 and 11014 RA.

ATTENTION: The pressure compensation cap must not be in place on the vent port for cleaning and disinfection.

10.7.3 Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65°C must not be exceeded. The relevant national requirements and the chemical manufacturer's instructions must be taken into account when using this method.

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**Hydrogen peroxide (H₂O₂) sterilization -
STERIS® AMSCO® V-PRO™1**

NOTE: The 'STERIS Low Temperature Reprocessing Device Compatibility Matrix' can be used to verify whether the medical device can be sterilized using the various STERIS V-PRO® devices.

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ATTENTION: Do not expose flexible endoscopes to direct sunlight or excessive heat.

ATTENTION: The distal tip of the flexible endoscope must be in the neutral position during storage.

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Reprocessing

ATTENTION: Do not store the flexible endoscope in the case.

ATTENTION: The sheath of the flexible endoscope must never be kinked or coiled up too tightly.

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11 Reprocessing of LUER port

⚠️ WARNING:

Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments (art. no. 96216003DER)' and the accompanying documentation.

⚠️ WARNING:

Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

⚠️ WARNING:

When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

ATTENTION: When preparing and using solutions, follow the chemical manufacturer's instructions, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

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ATTENTION: Danger of damage to medical devices: The use of chemicals that have not been approved by KARL STORZ or chemical manufacturers may cause damage to the medical devices. Only use chemicals approved for the medical device for reprocessing. A helpful guide is available on the Internet at www.karlstorz.com.

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NOTE: The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' (Art. no. 96216003) can be downloaded or requested by visiting www.karlstorz.com.

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Brushes: 27652, 27650B

Adaptor: 39107ALK (only for 11014L)

11.2 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

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Brushing the lumina

The working channels, lumina and hollow spaces must be precleaned under cold running water with the aid of a suitable brush.

Rinsing out with a pressure water gun

The working channel, lumina and hollow spaces must be rinsed using a pressure water gun and a suitable rinsing attachment. For rinsing, pulse at least 3 times for 5 seconds from a jet of water with a pressure of at least 2.5 bar. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

11.4 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003):

Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65°C must not be exceeded. The relevant national regulations and the chemical manufacturer's instructions must be taken into account when using this method.

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NOTE: If necessary, the instrument must be dried off manually afterwards.

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11.5 Connecting

In order to ensure effective machine cleaning, disinfection and rinsing, the instrument must be connected in the washer and disinfectant.

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- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Joints, threads and glide surfaces must be treated locally with instrument oil and/or special lubricant. (only valid for 11014L)
- Dismantled medical devices must be assembled.
- Afterwards, a functional check must be carried out.

NOTE: The oil used for this purpose must be suitable for the subsequent sterilization procedure (silicone-free and paraffin- or white oil-based).

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12 Technical data

Flexible cystoscope HD VIEW 11272 VH/VHU

Direction of view:	0°
Angle of view:	100°
Deflection:	up 210°, down 140°
Working length:	37 cm
Outer diameter:	16 Fr.
Working channel:	7 Fr.
Illumination:	Integrated LED light source
Electrical power consumption:	1 W
Sensor:	CMOS technology
Weight:	573 g
Classification:	Ila

Operating/storage conditions

Operating conditions:	15°C to +30°C
Operating conditions:	15% to 80% rel. humidity
Storage conditions:	0°C to +65°C
Storage conditions:	15% to 90% rel. humidity

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13 Servicing and repair

Defective instruments must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original parts only.

14 Maintenance

Preventive maintenance is not essential. However, regular maintenance can contribute to identifying potential problems before they become serious, enhancing the instrument's reliability and extending its useful operating life.

Maintenance services can be requested from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year.

15 Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this instrument only subject to the conditions that all assembly, system expansions, readjustments, modifications, or repairs to same have been performed by persons authorized by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

16 Warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see the section 'Subsidiaries'), even during the warranty period.

Opening of the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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11272 VH/VHU:

This medical device belongs to Class II a

This medical device bears a CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

18 Important information

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and sterilized/disinfected before they are sent for repair. We reserve the right to return contaminated instruments/devices to the sender.

Repairs, modifications or expansions which are not performed by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

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19 Disposal

This device has been marked in accordance with the European Directive 2012/19/EC on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the device as electronic scrap. Please ask KARL STORZ SE & Co. KG, a KARL STORZ subsidiary, or your authorized dealer about the collection point in your area. Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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20 Electromagnetic Compatibility user information

The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is suitable for use within close proximity to an active HF electro-surgical device in professional healthcare facility environment.

Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

WARNING:

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING:

The use of an accessory or cable with the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ other than those specified in this manual may result in increased emissions or decreased immunity of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™. When using an accessory or cable with the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ other than those specified in this manual, it becomes the responsibility of the user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ to determine compliance with IEC 60601-1-2 when using this item.

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 WARNING:

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

NOTE: The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If it is used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio communication service. The user might need to take mitigation measures, such as relocating or re-orienting equipment.

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Table 1 Guidelines and manufacturer's declaration – electromagnetic immunity			
The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is intended for use in an electromagnetic environment specified below. The user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ should make sure that it is used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1kV ± 1 kV for input and output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage Dip:</u> Dip to 0% for 1 cycles @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100% for 250/300 cycles	<u>Voltage Dip:</u> Dip to 0% for 1 cycles @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100% for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ requires continued operation during power supply interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50Hz/60Hz	30 A/m at 50Hz/60Hz	If image distortion occurs, it may be necessary to position the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	

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Table 3

**Guidance and manufacturer's declaration – electromagnetic immunity –
for equipment and systems that are not life-supporting**

The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is designed for operation in an electromagnetic environment as defined below. The user of the device should make sure that it is used in such an environment.

Interference immunity tests	EN/IEC 60601 test level	Conformity level	Electromagnetic environment – guidelines
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the $d = 1.2 \sqrt{P}$ 80 MHz to < 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

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Note 2: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects and persons.

^a) Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to stationary transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is used exceeds the above compliance levels, the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ should be monitored to ensure proper function. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

^b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4 Emission class and group		
Guidelines and manufacturer's declaration – electromagnetic emissions		
The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is intended for use in the electromagnetic environment specified below. The customer or user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ should ensure that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ uses RF energy only for its internal function. The customer or user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ should make sure that it is used in such an environment. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class-A	The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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Table 5 Recommended separation distances between portable and mobile RF telecommunication devices and the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™			
<p>The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ as recommended below, according to the maximum output power of the communications equipment.</p>			
Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>NOTE: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.</p> <p>NOTE: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.</p>			

HF surgical instrument compatibility

This equipment has been evaluated for compatibility with high-frequency surgical equipment.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@kseas.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Carneros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 México D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
127 39 Skärholmen, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kęstucio st. 59 / Lenktoji st. 27
08124 Vilnius, Lithuania
Phone: +370 5 272 0448
Mobile: +370 685 67 000
E-Mail: info-lt-lv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
84124 Høje, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopia Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopia France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Adria Eos d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Cesta v Gorice 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: + 386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongquiao Road, Xuhui District,
200030, Shanghai, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

www.karlstorz.com

KARL STORZ SE & Co. KG

Dr.-Karl-Storz-Straße 34
78532 Tuttlingen

Postfach 230
78503 Tuttlingen

Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com