VARI
LINDAL GROUP

Coordinata da Lindal Group Holding GMBH – Socio Unico

INSTRUCTION MANUAL

The KEMP HFA pMDI valve is designed to be used with the open end of the valve stem facing in a generally downward direction.

The KEMP HFA pMDI valve is suitable for use with a range of inhalation actuators, suitable for pulmonary delivery of drugs. For the latest recommendations on actuators contact V.A.R.I S.p.A. In standard operation the valve stem is depressed relative to the valve ferrule, until the full travel position is reached. In order to achieve the full travel state of the valve a force equal to or greater than 30 newtons is typically required.

The operating sequence of the The KEMP HFA pMDI valve, can be described as a number of progressive stages.

Stage 1

Prior to any depression of the valve stem the valve is at rest. This position allows the filling of the metered chamber with the liquid phase.

Stage 2

Initial depression of the valve stem, brings the internal cylindrical portion of the stem in to a sealing contact with the flexible flange section of the thermoplastic elastomer metering chamber. At this stage the metered volume to be delivered is isolated from the bulk contents of the canister.

Stage 3

Continued depression of the valve stem relative to the ferrule moves the exit hole in the valve stem in to the isolated volume of the metering chamber. When the exit hole enters the metering chamber the contents are released into the atmosphere. It is recommended that the valve is held in the fully depressed position for a time greater than 2.0 seconds.

Stage 4

When the depressing force on the valve stem is removed the stainless steel coil spring urges the valve stem back towards the rest position. The exit orifice in the valve stem is first covered by the wall of the inner gasket, sealing off the flow-path to atmosphere. The contents of the canister is forced in to the metering chamber, via the filling paths in the inner section of the valve stem. Typical time for the dose chamber to refill may be as much as 5 seconds.

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