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GB. Dental carbide burs

Instruction manual

Dental carbide burs are designed for treatment of tooth fundamental substance, fillings and metal structures of prosthetic dentures when attaching to a straight, angled or turbine dental handpiece. Scope of application - dentistry. Dental carbide burs are designed for operation and storage at a temperature from plus 10 to plus 35°C, air relative humidity should not exceed 80% at a temperature plus 25°C. Indoor air should be free of impurities that cause corrosion. Dental carbide burs can be used for at least 300 minutes of machine running time. They are supplied non-sterile.

SAFETY MEASURES

- Dental carbide burs can only be used by qualified medical personnel characterized with sufficient experience and good manual skills;
 - do not use corroded and mechanically damaged tools;
 - do not use a tool with a nominal diameter smaller than the maxilla alveolar canal diameter;
 - touching the tool working part with a hand is not allowed when working with hard-alloy dental drills to exclude cuts against sharp edges;
 - tool rotation speed should not exceed 300,000 rpm for the turbine dental handpiece. Should not exceed 30,000 - 160,000 rpm for straight and angled dental handpiece.
- Use the tool pursuant to its purpose. Dental carbide burs are made of non-toxic material. Used tool disposal is executed in the manner usual for this kind of tools.

PREPARATION FOR WORK

Working part is made of hard alloy, shank end - of stainless steel.

Dental carbide burs must be consistently disinfected, cleaned and sterilized before use. For this purpose you should use solutions with neutral pH.

Disinfection: Place the tool in the disinfectant solution before initial and repeated use.

Observe the curing period and recommendations specified by the disinfectant manufacturer. Then the tools must be washed with running water.

Cleaning: Dental carbide burs must be subjected to ultrasonic cleaning within 10 minutes, then be rinsed in distilled water within 0.5 min. Dry the tool at a temperature of 85°C until the moisture complete disappearance.

Sterilization: Before sterilization, put tools on sterilization stands or trays. Sterilize tools with steam in the autoclave: the curing period should make 5 minutes at a temperature of 134°C and a pressure of 0.21 MPa, or the curing period should make 20 minutes at a temperature of 131°C and a pressure of 0.11 MPa.

When sterilizing a large number of tools in one cycle, make sure that the autoclave maximum loading is not exceeded. After sterilization and cooling, place tools in sealed containers that ensure the tool sterility preservation and eliminate accidental mechanical damage, tool removal should be executed with sterile tweezers fixed at the shank end.

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