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GB. DENTAL DIAMOND BURS

Instruction for Use

DENTAL DIAMOND BURS are used with either straight, counterangular or turbine handpieces.

DENTAL DIAMOND BURS are used for operating at temperature from +10 to +25°C, relative humidity is no more than 80% at +25°C.

Diamond heads should have being kept in closed area at temperature from +5 to +35°C and at relative humidity no more than 98% at +25°C. The air in a room must not consist of any pollution with corrosive effect.

SAFETY INSTRUCTIONS

- only qualified medical personnel should use carbide cutters;
- don't use bent, damaged, overworn, dull, nonconcentric rotative instruments, instruments with overworn or crushed diamond layer;
- don't bend the instrument and don't use it as a lever during any work;
- for cooling the working area use diamond heads with the handpieces with enough amount of cooling liquid;
- check the fixation of diamond head in handpiece and the working efficiency of the handpiece;
- use the form of working part and the type of grain according to specific purpose;
- to start using the instrument it must achieve the working speed;
- recommended pressure: 0,3 – 2 N.
- utilization must be in regular way for such kind of products.

PREPARATION

The heads are made of stainless steel.

Diamond dental heads are supplied in non-sterile condition. Before use they must be disinfected, cleaned and sterilized. Only solutions with neutral pH value must be used for this purpose. Put on gloves and protective eyeglasses when treating dirty instruments.

Disinfection: After the use place the heads immediately in the disinfecting solution. Follow the instructions concerning the duration and other recommendations given by the manufacturer. Then wash the heads with running water.

Cleaning: To remove the hard debris of the dental enamel wash the heads in an ultrasonic washer during 10 minutes or immerse them in a special solution that dissolves the organic debris of dental tissues. Then wash the heads in the distilled water during 0.5 minute. If necessary repeat the cleaning procedure. Dry the heads at the temperature plus 85°C till the moisture disappears completely.

Sterilization: Before the sterilization place the heads on sterilization supports or trays. Sterilize the heads with steam in an autoclave:

- at the temperature plus 134°C and the pressure 0.21 MPa during 5 minutes, OR

- at the temperature plus 121°C and the pressure 0.11 MPa during 20 minutes.

When sterilizing many instruments during one cycle be sure that the maximum loading of the autoclave has not been exceeded.

After the sterilization and cooling place the heads in containers ensuring the safety and sterility of the cutters. The extraction of needed diamond head should be done with sterile forceps on the shank.

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