AMS Ambicor™

Penile Prosthesis

Instructions for Use

English	AMS Ambicor™ <i>Penile Prosthesis</i> Instructions For Use 3
Français	AMS Ambicor™ <i>Prothèse pénienne</i> Mode d'emploi 21
Deutsch	AMS Ambicor™ <i>Schwellkörperimplantat</i> Gebrauchsanweisung 39
Italiano	AMS Ambicor™ <i>Protesi peniena</i> Istruzioni per l'uso 57
Español	AMS Ambicor™ <i>Prótesis de pene</i> Instrucciones de uso 75
Português	AMS Ambicor™ <i>Prótese Peniana</i> Instruções de Utilização 93

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Penile Prosthesis

INSTRUCTIONS FOR USE

NOTE: Refer to the Operating Room Manual for further information on the AMS Ambicor Device and its implantation.

Brief Device Description

The AMS Ambicor Penile Prosthesis is a closed fluid-filled system consisting of two implanted components: a pair of cylinders implanted in the corpora cavernosa, and a pump implanted in the scrotum. Kink-resistant tubing connects the components. The device is delivered prefilled with normal saline and preconnected. The cylinders are inflated as fluid is pumped from the reservoirs, which are located at the proximal end of the cylinders, to the main portion of the cylinder, creating an erection. They are deflated as fluid is transferred back to the reservoirs, making the penis flaccid once again. This device contains solid silicone elastomer. These devices are for men who, after appropriate patient history and diagnostic evaluations as well as discussions with the urologist about other alternative treatment methods, are determined to be suitable candidates for implantation surgery.

Indications for Use

The AMS Ambicor Penile Prosthesis is intended for use in the treatment of chronic, organic, male erectile dysfunction (impotence).

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Warnings

1. **Implantation of the device will make latent natural or spontaneous erections as well as other interventional treatment options, impossible.**
2. Men with diabetes, spinal cord injuries, or open sores may have an increased risk of infection associated with a prosthesis.
3. Failure to evaluate and promptly treat erosion may result in a substantial worsening of the condition leading to infection and loss of tissue.
4. Implantation of a penile prosthesis may result in penile shortening, curvature or scarring.
5. This device contains solid silicone elastomer. The risks and benefits of implanting this device in patients with documented sensitivity to silicone should be carefully considered.
6. Pre-existing abdominal or penile scarring or contracture may make surgical implantation more complicated or impractical.

Precautions

Surgery Related

1. Spontaneous unintended inflation or deflation of the cylinders may occur and may result in unintended partial or full erections.
2. Migration of the device components can occur if the cylinders are improperly sized, or if the pump is not positioned properly.
3. Removal of an implanted prosthesis without timely reimplantation of a new prosthesis may complicate subsequent reimplantation or may make it impossible.
4. Improper measurement technique, positioning or sizing may reduce cylinder life.
5. Unsuccessful outcomes have been reported due to improper surgical technique, anatomical misplacement of components, improper sizing of components, or tubing kinks.

Device Related

1. Some of the materials used in the construction of this device have been shown to cause minor irritation when implanted in animals. Therefore, implantation of this device may cause minor irritation or discomfort in some patients.
2. Do not use product that has damaged or open packaging, as sterility may be compromised.

Patient Related

1. A thorough preoperative consultation should include a discussion between patient and physician of all available treatment options and their risks and benefits.
2. Adequate patient manual dexterity and strength are required for proper device inflation and deflation.
3. Mental or psychological conditions, such as senile dementia, may inhibit the patient's successful operation of the prosthesis.
4. Trauma to the pelvic area, such as impact injuries associated with sports (e.g. bicycle riding), can result in damage of the implanted device and/or surrounding tissues. This damage may result in the malfunction of the device and may necessitate surgical correction, including replacement of the device.
5. The implantation of this device should only be considered in patients whom the physician determines are adequate surgical candidates.

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A clinical trial was conducted to determine the safety and effectiveness of the AMS Ambicor inflatable penile prostheses. A total of 179 patients were implanted with follow-up out to 18 months or longer for 141 patients. The Adverse Device Effects, detailed in the table below, were noted during the duration of this clinical trial for all enrolled patients.

AMS Ambicor Clinical Trial Device Related Adverse Events*

ADE	# Patients (%)	# Events	Mean Onset time in Days (Range in Days)
Urogenital Pain (Typically Associated with Healing Process)	38 (21.2%)	41	83 (0 – 548)
Patient Dissatisfaction (With Length, Ability to Use and Nonspecific Reasons)	22 (12.3%)	25	269 (0 – 1105)
Mechanical Malfunction	11 (6.1%)	11	394 (0 – 1508)
Penile Sensation Change	11 (6.1%)	11	266 (10 – 573)
Urogenital Hematoma	9 (5.0%)	9	9 (2 – 34)
Urogenital Edema	8 (4.5%)	9	9 (0 – 29)
Infection/Device Infection	6 (3.4%)	8	207 (32 – 606)
Auto-inflation	4 (2.2%)	4	110 (27 – 224)
Penile Curvature	4 (2.2%)	4	323 (0 – 999)

ADE		# Patients (%*)	# Events	Mean Onset time in Days (Range in Days)
Abnormal Ejaculation		3 (1.7%)	3	256 (40 – 487)
Application Site Reaction (Wound Separation, Delay in Cutaneous Closure)		3 (1.7%)	3	26 (0 – 46)
Phimosis/Paraphimosis		3 (1.7%)	3	21 (3 – 38)
Post-Operative Bleeding		3 (1.7%)	3	4 (2 – 6)
Seroma		3 (1.7%)	3	242 (34 – 596)
Urinary retention		3 (1.7%)	3	121 (0 – 361)
Urogenital Ecchymosis		3 (1.7%)	3	123 (6 – 351)
Adhesion		2 (1.1%)	2	42 (34 – 50)
Auto-deflation		2 (1.1%)	2	486 (212 – 760)
Nocturia		2 (1.1%)	2	245 (105 – 385)
Penile Rash		2 (1.1%)	2	272 (43 – 501)

ADE	# Patients (%*)	# Events	Mean Onset time in Days (Range in Days)
Renal Failure	2 (1.1%)	2	476 (396 – 555)
Urethral Laceration	2 (1.1%)	2	1 (0 – 2)
Urinary Tract Infection	2 (1.1%)	2	283 (4 – 561)
Urination Impaired	2 (1.1%)	2	225 (64 – 385)
Urogenital Erythema	2 (1.1%)	2	193 (13 – 373)
Other urogenital	23 (12.8%)	24	97 (0 – 727)
Other non-urogenital	10 (5.6%)	14	176 (0 – 556)

* Percentages based on total number of patients implanted (n=179).

The following “other” urogenital adverse events each occurred in one patient: abnormal sexual function, allergic reaction, BPH, burning urination, crossover from left to right corpora, device migration, dribbling, dysuria, epididymitis, extrusion, incontinence, induration in penile area, lower urinary tract symptoms, malposition, meatal excoriation, micro-hematuria, penile scar tissue, prerenal azotemia, scab on corpora, soft glans, tunica perforation, urosepsis, wound dehiscence and wound infection.

Resolution of Device Related Adverse Events

ADE	# Patients (%)	# Events	Method of Resolution				Continuing
			Surgical	Medication	Other**	No Medical Intervention	
Urogenital Pain	38 (21.2 %)	41	3	9	1	29	8
Patient Dissatisfaction	22 (12.3 %)	25	3	3	3	16	14
Mechanical Malfunction	11 (6.1 %)	11	7	0	1	3	1
Penile Sensation Change	11 (6.1 %)	11	0	2	0	10	4
Urogenital Hematoma	9 (5.0 %)	9	0	1	3	6	0
Urogenital Edema	8 (4.5 %)	9	0	2	2	5	0
Infection/Device Infection	6 (3.4 %)	8	6	1	0	1	1
Auto-inflation	4 (2.2 %)	4	0	0	0	4	2

ADE	МИ в Росздравнадзоре ru info@nevacert.ru	Method of Resolution						
		# Patients (%)	# Events	Surgical	Medication	Other**	No Medical Intervention	Continuing
Penile Curvature		4 (2.2 %)	4	2	0	0	2	2
Abnormal Ejaculation		3 (1.7 %)	3	0	1	0	2	1
Application Site Reaction		3 (1.7 %)	3	0	1	0	2	0
Phimosis/Paraphimosis		3 (1.7 %)	3	1	2	0	0	0
Post-operative Bleeding		3 (1.7 %)	3	0	0	0	3	0
Seroma		3 (1.7 %)	3	0	0	1	2	0
Urinary Retention		3 (1.7 %)	3	1	1	2	0	1
Urogenital Ecchymosis		3 (1.7 %)	3	0	0	0	3	0
Adhesion		2 (1.1 %)	2	1	0	0	1	0

ADE	# Patients (%)	# Events	Method of Resolution				Other**	No Medical Intervention	Continuing
			Surgical	Medication					
Auto-deflation	2 (1.1 %)	2	0	0			1	1	0
Nocturia	2 (1.1 %)	2	0	1			0	1	2
Penile Rash	2 (1.1 %)	2	0	1			1	0	1
Renal Failure	2 (1.1 %)	2	0	0			1	1	1
Urethral Laceration	2 (1.1 %)	2	0	0			1	1	0
Urinary Tract Infection	2 (1.1 %)	2	0	2			0	0	0
Urination Impaired	2 (1.1 %)	2	0	1			1	0	2
Urogenital Erythema	2 (1.1 %)	2	0	2			0	0	0

** Other methods of resolution consist of standard medical treatments.

The following are adverse events associated with this product: Abscess, Bleeding, Erosion, Excessive Fibrous Capsular Growth, Exposure to Biohazardous Material, Genital Changes, Granulomas, Hemorrhage, Immune-Related Connective Tissue Disorders, Improper Size, Inguinal Hernia, Ischemia, Necrosis, Pain (which may be prolonged or severe), Perforation or Injury of Bladder, Corpus Cavernosum, Nerve, Tunica, or Urethra, Prolonged Procedure, Super Sonic Transporter (SST) (as a result of improper sizing), Ulceration, Unretrieved Device Fragment, Vascular Compromise, Ventral Chordee.

The following events occurred during this clinical study but were classified by the Investigators as not related to the device or the implantation surgery: urethral stricture (9 cases), prostate cancer (4 cases), BPH (3 cases), incontinence (3 cases), urinary retention (3 cases), urinary tract infection (3 cases), abnormal prostate exam/PSA (2 cases), abnormal sexual function (2 cases), balanitis (2 cases), hematuria (2 cases), pain (2 cases), renal calculus (2 cases), abnormal ejaculation (1 case), bladder instability (1 case), bladder spasms (1 case), blood in urine (1 case), dysuria (1 case), edema (1 case), epididymal cyst (1 case), eroded artificial sphincter (1 case), erythema (1 case), infection (1 cases), inguinal hernia (1 case), nocturia (1 case), prostate nodule (1 case), prostatitis (1 case), pyuria (1 case), renal colic (1 case), renal insufficiency (1 case), testicular cyst (1 case), testicular hydrocele (1 case), trichomonas (1 case), urethral mucosal atrophy (1 case), urination impaired (1 case), urgency (1 case), and urosepsis (1 case).

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There were four patient deaths during the course of the trial. No deaths that occurred during the duration of the clinical study were attributed to the device implantation or use.

A total of 14 patients underwent revision surgeries in the study. Information on device revisions is described in the “Clinical Studies” section.

Clinical Studies

A clinical trial was undertaken to demonstrate that the AMS Ambicor Inflatable Penile Prosthesis provides an erection that is suitable for intercourse and has acceptable rates of surgical revision and of significant clinical events associated with the implantation and use of these devices. This trial was also designed to demonstrate the implantation of these devices does not negatively impact the sexual satisfaction, psychological well being, self-esteem or quality of life of patients who receive these devices. It was a prospective, multi-center cohort trial in which the patients served as their own control.

One hundred and seventy-nine (179) male patients, over 21 years of age, were enrolled and implanted in this study. All patients with diagnosed organic erectile dysfunction were eligible for enrollment, if they did not present with a history of allergy/sensitivity to silicone, preexisting autoimmune or connective tissue diseases or active urogenital infection.

All safety-related data, diagnoses and health status evaluations were captured on detailed case report forms. The Investigators' professional evaluation of the erections provided by the IPPs after implantation and their suitability for intercourse was the primary efficacy endpoint. The number of surgical revisions performed and reported by the Investigators was the primary safety endpoint. Patient self-evaluations on four validated outcome instruments were the secondary efficacy endpoints (concerning quality of life, self-esteem, and sexual satisfaction and functioning).

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This clinical trial provided the following results through the 18-month evaluation for the first 141 patients to reach this post-surgical follow-up. Some patients were studied for longer than 18 months.

Physician Assessment of Device Function:

The function of 141 implanted devices was evaluated on or after the 18-month follow-up exam. Of these 141 devices, 140 could be inflated (99.3%). The device that could not be inflated was surgically replaced. All of the 140 devices that were successfully

inflated were determined to provide an erection suitable for intercourse. This assessment of device function does not include the following information regarding device malfunctions: Of the 140 devices found to be functioning properly on or after 18 months, 2 required revision surgery prior to this evaluation to correct a mechanical malfunction. These and other cases of revision surgery are discussed further in the next section.

Surgical Revisions:

The incidence of revisions was evaluated in the 141 patients with follow-up out to at least 18-months, and an additional 6 patients who were explanted prior to 18 months and subsequently discontinued from the study. (A revision is considered any urogenital surgical intervention that is related to the function, placement or site reaction to the implanted device.) Of these 147 patients, 14 (9.5%) experienced a total of 15 revision surgeries, and 133 (90.5%) have not been revised.

The average time to the first revision surgery was 12 months (ranging from 0 to 50 months). Of the fifteen (15) revision surgeries, there were six (6) revisions due to “Mechanical Malfunction”; five (5) due to “Infection”; one (1) due to “Migration/ Extrusion” (pump migration and cylinder aneurysm/bulging/extrusion); one (1) due to “Malposition”; one (1) due to “Phimosis”; and one (1) due to “Patient Dissatisfaction”.

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No decline in overall patient quality of life, self-esteem, psychological well being, sexual functioning or sexual satisfaction was demonstrated during this clinical trial.

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The psychological well being of these patients, as judged by the Brief Symptom Inventory (BSI) was determined to be equivalent to the pre-implant state. The self-esteem of the patients, as evaluated by the Rosenberg Self-Esteem Scale, was also equivalent to the pre-implant state. The sexual functioning and sexual satisfaction did improve from the pre-implant state, as measured by the Sexual History Form (SHF).

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Implantation of a penile prosthesis may result in penile shortening, curvature or scarring. The prosthetic erection may differ from the patient's original, natural erection in that it may be shorter, less firm, have less girth, and reduced sensations. Realistic cosmetic expectations should be communicated to the patient and should include the potential for skin scarring, scrotal deformity, pump bulge in the scrotum, lack of concealability and other possible adverse events. Patients should also be aware that penile prostheses are not considered to be lifetime implants.

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Patients who undergo revision surgery may notice a change in the character of their erection compared to their previous implant, which may include differences in sensation, length, girth, rigidity, and/or flaccidity. It is also important that the physician discusses with the patient the possibility of an allergic reaction to the materials in the device (See Silicone Information).

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Scientific literature has included reports of adverse events and other observations in patients with implantable silicone devices. As reported, these events/observations indicate “allergic-like” symptoms and in other cases a symptom complex associated with immunological disorders. No casual relationship has been established between these events and silicone elastomer or fluorosilicone lubricant.

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Extensive testing has been conducted on all materials that comprise the prostheses in the AMS Ambicor. This testing has indicated no toxicological response attributable to the materials. However, some of the materials caused minor irritation when implanted in animals.

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Non-clinical testing has demonstrated the

Ambicor product line is MR Conditional.
The device can be scanned safely under the following conditions:

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MRI-Related Heating

Non-clinical testing has demonstrated the Ambicor product line produced the temperature rises during MRI performed for 15 minutes of scanning in the respective MR systems which would not pose a hazard to the human subject.

Static Magnetic Field	1.5 Tesla ^a	3.0 Tesla ^b
Highest Temperature Change	≤ + 0.4° C	≤ + 1.9 ° C
(a) 1.5T - 64 MHz MR System (General Electric Healthcare, Milwaukee, WI) (b) 3.0T MR Excite, General Electric Healthcare, software version 14X.M5		

Artifact Information

Non-clinical testing has demonstrated that the Ambicor product line may compromise the MR image quality if the area of interest is relatively close to the position of the implant. The maximum image artifact produced by a MR gradient echo pulse sequence was a “moderate” localized signal void in size and shape of the implant. Optimization of MR imaging parameters to compensate for the presence of the device may be necessary.

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For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

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In the United States

Before returning any components, whether explanted or unused (sterile or nonsterile), customers must fill out the Return Goods Form located on the last page of the Patient Information Form. Follow all of the instructions on the form carefully, and be sure that the components have been thoroughly cleaned before returning them to AMS.

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This page intentionally left blank.

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TABLE OF CONTENTS

GENERAL INFORMATION

Detailed Device Description	4
Device Characteristics	4
Two-piece Device	4
Easy to Inflate and Deflate	4
Device Function	4
Inflation	4
Deflation	4
Device Components	5
Cylinder Sizes	5
AMS Ambicor Packaging	5
Sterilization	5
Sterilization of Components	5
AMS Tools	6

OPERATING ROOM INSTRUCTIONS

Operating Room Instructions	7
Preoperative Set-up Materials	7
Documentation	7
Patient Preparation	7
Unpacking the Device	7
Intraoperative Procedures	8
Sizing	10
Method	10
Inflate/Deflate Test	10

POSTOPERATIVE PROCEDURES

Postoperative Procedures	11
Immediately Postoperative	11
After the Patient is Released from the Hospital	11
Evaluating Long-Term Function and Placement	11

PROCEDURE ILLUSTRATIONS

Procedure Illustrations	12
-------------------------------	----

BRIEF SUMMARY

Brief Summary	14
---------------------	----

GENERAL INFORMATION

DETAILED DEVICE DESCRIPTION

The AMS Ambicor™ Penile Prosthesis is a closed, fluid-filled system consisting of a pair of cylinders, which are implanted in the corpora cavernosa, and a pump which is implanted in the scrotum. All components are connected by kink-resistant tubing. The device is delivered prefilled with normal saline and preconnected. The cylinders are inflated as fluid is pumped from the reservoirs, which are located in the proximal ends of the cylinders, into the main cylinder body, creating an erection. They are deflated as fluid is transferred back to the reservoirs, making the penis flaccid once again.

DEVICE CHARACTERISTICS

Two-piece Device

This prefilled and preconnected two-piece (pump and cylinders) device eliminates the need to fill and connect the device during surgery.

Easy to Inflate and Deflate

The inflation pump in the scrotum is easy for the patient to locate and requires few squeezes to inflate the device, because each squeeze displaces a large amount of fluid. To deflate the device, the patient simply bends the penis towards the scrotum, holds the position for 6-12 seconds, and then releases.

DEVICE FUNCTION


Inflation

The inflation pump is placed in the scrotum. To inflate the cylinders, the patient should stabilize the pump with one hand and use the other hand to squeeze and release the pump bulb several times to make the cylinders stiff. When the cylinders are fully inflated, the pump bulb will be hard and can no longer be compressed.


Deflation

To deflate the device, the patient places his thumb under the shaft of the penis to act as a fulcrum and his fingers on top of the shaft of the penis. Using the fingers of the same or opposite hand, the patient should bend his penis down over his thumb (fulcrum) towards his scrotum to a 55-65 degree angle, making sure that both cylinders are bent. The cylinders should be held in this position for approximately 6-12 seconds and then released. This will open the valves which allow the fluid to flow back into the reservoirs and pump. The patient may also deflate the device by bending the cylinders upwards; this can be done by placing his thumb on top of the shaft of his penis and following the same instructions as above.

Inflation

Deflation

DEVICE COMPONENTS

Cylinder Sizes

The AMS Ambicor penile prosthesis is provided in the following cylinder sizes:

Sizes (diameter by length)

12.5 mm x 14 cm

12.5 mm x 16 cm

12.5 mm x 18 cm

14 mm x 16 cm

14 mm x 18 cm

14 mm x 20 cm

15.5 mm x 18 cm

15.5 mm x 20 cm

15.5 mm x 22 cm

AMS Ambicor Packaging

The prosthesis is packaged in a fluid-filled inner pouch. The prosthesis is preconnected, prefilled, and sterilized in a pouch within a pouch system which is placed in a dustcover box. The proximal cylinder ends are secured in a protective cushion to minimize movement within the pouch.

Rear tip extenders (RTEs) are packaged in the fluid-filled pouch with the Ambicor device. These RTEs include one pair each of 0.5 cm, 1 cm, 2 cm, and 3 cm RTEs.

Ambicor RTEs are also available in a separate RTE kit that includes one pair each of 0.5 cm, 1 cm, 2 cm, and 3 cm RTEs, which is provided sterile.

Each Ambicor prosthesis package also contains:

- 2 Keith needles
- Patient Information Form (PIF)
- One mailing envelope to return completed PIF to AMS
- Patient Medical Identification Card
- Ambicor Instructions for Use

NOTE: For warnings, precautions and contraindications please refer to the AMS Ambicor Instructions for Use provided in the product packaging.

STERILIZATION

Sterilization of Components

American Medical Systems sterilizes the Ambicor Penile Prosthesis and RTEs. Under normal storage conditions, the devices will remain sterile until the expiration date as long as the sterile barriers of the packaging remain intact. To protect the integrity of the packaging and the function of the prosthesis, store the sterilized devices on a protected shelf or in a cabinet. The environment should be clean, dry, and near room temperature. For maximum protection during storage, leave the Ambicor prosthesis and RTEs in their packaging. Inspect the packaging for damage and check the expiration date before use.

CAUTION: Do not reuse, reprocess, or resterilize the AMS Ambicor Penile Prosthesis or AMS Ambicor RTEs. They are intended for single use only.

CAUTION: Do not use product that is past its expiration date.

GENERAL INFORMATION (CONTINUED)

AMS Tools

American Medical Systems has surgical instruments that can be used during the surgery to help facilitate the surgeon's implantation of the penile prosthesis. For sterilization information, refer to the instructions provided with the tools.

The following tools can be ordered from AMS.

Non-sterile.

These tools may be resterilized and reused.

- AMS Sizer
- Furlow Insertion Tool
- AMS Closing Tool

Sterile

- Proximal Tool (dual use for cylinder insertion and closing)
- Disposable Dilators
- SKW Retractor Kit

CAUTION: Do not resterilize or reuse the sterile tools. They are intended for single use only.

CAUTION: Do not use product that is past its expiration date.

OPERATING ROOM INSTRUCTIONS

OPERATING ROOM INSTRUCTIONS

The following instructions are intended as a guide for the surgeon.

CAUTION: This device is to be used only by physicians who are knowledgeable and have received appropriate training regarding the use of inflatable penile prostheses. This Manual is not intended to be a complete reference.

Preoperative Set-up Materials

The hospital should provide those instruments normally required for a urological surgical procedure. In addition to the AMS Ambicor penile prosthesis, you will need the following sterile setup:

Instruments:

- 2 Keith needles provided in Ambicor packaging for passing cylinder traction sutures through glans
- Corporal length measuring tool.
- Furlow Insertion Tool for passing traction sutures through glans.
- Corporal dilators (7mm-17mm)
- 1 sterile surgical stand (Mayo stand) or a stainless steel tray
- 1 basin of sterile normal saline
- Antibiotic solution for irrigation (optional)
- AMS Closing Tool (optional)

Documentation

- One AMS Ambicor Penile Prosthesis Instructions for Use brochure
- One Patient Information Form (PIF)
- One mailing envelope for returning the completed PIF to AMS
- One AMS Ambicor Patient Medical Identification Card

Patient Preparation

Before the surgery, doctors should take adequate steps to limit the risk of post-operative infection.

Once the patient is in the operating room, the abdominal and genital area should be shaved. Following the shave, the area should be scrubbed with povidone-iodine soap for ten minutes or in accordance with the approved hospital preoperative scrub procedure.

Establish the sterile field, drape, and prepare the patient according to the physician's instructions. Throughout the procedure, the surgical site should be flushed with copious amounts of broad-spectrum antibiotics.

The patient should be positioned so that a penoscrotal incision can be made.

Unpacking the Device

NOTE: Do not unpack the prosthesis until after the surgeon dilates and measures both corpora cavernosa intraoperatively.

NOTE: The adhesive label at one end of the dustcover box and the small, peelable labels on the foil pouches show the part and serial/lot numbers and the size of the components. Record this information on and affix this label to the Patient Information Form (PIF).

Keep the sterile pouches in the dustcover boxes until the components are in the operating room. Remove the pouches from the dustcover box in the operating room by opening the unlabeled end of the dustcover box.

OPERATING ROOM INSTRUCTIONS (CONTINUED)

To remove the Ambicor device from its packaging:

1. Remove the inner fluid-filled foil pouch from the outer pouch.

CAUTION: To avoid the possibility of cutting the Ambicor prosthesis or RTEs, do not cut the package open with scissors.

2. Place the sterile, fluid-filled foil pouch onto a sterile, lint free Mayo stand.

CAUTION: Cloth towels placed on the Mayo stand may transfer lint to the AMS components.

3. Hold the foil pouch over a kidney basin containing sterile normal saline. Carefully tear the foil pouch at the side notch. Empty the Ambicor device, its RTEs, and the fluid into the kidney basin. Add sterile normal saline as needed to cover the components in the basin.

NOTE: Do NOT leave the prosthesis exposed to air. Submerge the prosthesis in the fluid from the foil pouch or in sterile normal saline immediately after removing it from the foil pouch. Keep the prosthesis submerged until implantation to prevent air from entering the device.

4. Remove the tips of the Ambicor cylinders from the protective end cushion.
5. Discard the end cushion.
6. Leave the device and RTEs in the basin until ready to implant.

Intraoperative Procedures

Surgical Procedures

Assemble the necessary equipment and position the patient for a penoscrotal surgical approach. The penoscrotal approach leaves the incision well hidden and provides convenient access to the corpora cavernosa.

Establish the sterile field, drape and prepare the patient according to the physician's instruction. Throughout the procedure the surgical site may be flushed with copious amounts of broad-spectrum antibiotic.

The following description is an overview of the penoscrotal surgical approach.

Penoscrotal Approach

1. To begin, place a Foley catheter to facilitate identification of the urethra.
2. Some physicians use the SKW Retractor System to place the penis on "stretch" to provide exposure of the corpora.
3. Make a 2 cm to 3 cm incision through the median raphe of the scrotum at the penoscrotal angle. Some physicians may prefer a high scrotal incision for better proximal corporal access. Then, laterally retract the corpus spongiosum to avoid damaging the urethra.
4. Dissect through to Buck's fascia to expose the tunica albuginea. Place stay sutures to use as a reference point when measuring the corpora. Then, make an incision into one of the corpora cavernosa.
5. Dilate the proximal corpus (crus) and the distal corpus to create a space for inserting a penile cylinder. Dilate the corpus cavernosum to 13-16 mm.

Note: the available Ambicor cylinders diameters are 12.5 mm, 14 mm, and 15 mm.

OPERATING ROOM INSTRUCTIONS (CONTINUED)

After dilating one corpus cavernosum, incise and dilate the adjacent corpus cavernosum following the same procedure.

6. To select the cylinders and rear tip extenders that will fit the patient's anatomy, measure each corpus proximally and distally using the Furlow Insertion Tool or appropriate measuring device (See Sizing section for choosing cylinder size.). As a general rule, the corporotomy is best placed when two-thirds of the total corporal measurement is distal to the incision and one-third is proximal. This facilitates the placement of the cylinders, and may avoid the need to extend the corporotomy during the procedure.
7. Select the appropriate size cylinders, apply rear tip extenders (if necessary), and use the insertion tool(s) to introduce the cylinders into the corpora cavernosa, as follows:
 - a. When inserting the cylinder distally, use the Furlow Insertion Tool to put the penis on mild stretch. The tool should be palpable beneath the glans.
 - b. First thread the pulling suture at the front of the cylinder through a Keith needle, then place the Keith needle in the insertion tool.
 - c. Insert the Furlow insertion tool distally into the corpora cavernosa and pass the needle through the glans.
 - d. Inflate the cylinders fully, then insert the proximal end of the cylinder.
 - e. Deflate the cylinders by bending both evenly to 55-65 degrees and holding them for 6-12 seconds.
 - f. Position the distal ends of the cylinders by pulling the traction sutures.
 - g. Use the Proximal Tool to place the proximal cylinder ends into the corpora.

During insertion, be sure that the Furlow Insertion Tool is in the ipsilateral corpora at the distal penis. Because the intracavernosal septum may be inconsistent distally, it is easy to cross over to the contralateral side.

If you think you may have crossed over, place a dilator into the other side. If the cylinder does cross over, no repair is necessary. Simply remove the cylinder and place it correctly.

NOTE: The corporal incision may need to be extended to ensure the input tube exits directly from the corporotomy.

Once both cylinders are implanted, close the tunica albuginea. Place the AMS Closing Tool or the closing end of the Proximal Tool (or other suitable instruments) over the cylinder to protect it from inadvertent needle injury, so that the corporal body can be closed. This should be done with meticulous attention to hemostasis.

8. Use blunt dissection to form a pocket in the most lateral and dependent portion of the scrotum.
9. Insert the pump into the scrotal pocket. The tubing between the pump and cylinders should not be palpable to the patient.
10. After both the cylinders and pump are implanted, check the function of the prosthesis by inflating and deflating the device.
11. Some physicians close the dartos in two layers with a running 2-0 absorbable suture. Close the skin.
12. Apply a wound dressing. Leave the device partially inflated. Some choose to tape the penis to the abdomen overnight. In addition, a drain may be placed for 12 to 24 hours.

For postoperative instructions, see the section entitled "Postoperative Procedures".

OPERATING ROOM INSTRUCTIONS (CONTINUED)

Sizing

Below is the recommended method of selecting cylinder sizes for the AMS Ambicor penile prostheses.

Method

This method allows the tubing to exit directly from the corporotomy. Follow the formula described below to select the appropriate cylinder length and number of rear tip extenders. If necessary, extend the length of the corporotomy.

1. Calculate the Total Corporal Length (Distal + Proximal)

Example

Distal Corporal Length 11 cm

Proximal Corporal Length +8 cm

Total Corporal Length 19 cm

2. Subtract 2 cm from the Total Corporal Length to obtain an Adjusted Measurement.

Example

Total Corporal Length 19 cm

 -2 cm

Adjusted Measurement 17 cm

3. Select the closest cylinder size that is shorter than or equal to the Adjusted Measurement.

Example

Adjusted Measurement 17 cm

Selected Cylinder Length 16 cm

4. Subtract the Selected Cylinder Length from the Total Corporal Length to determine the length of rear tip extenders required to fit the patient.

Example

Total Corporal Length 19 cm

Selected Cylinder Length -16 cm

Rear Tip Extender Length 3 cm

Inflate/Deflate Test

After all components are implanted, inflate the device to check the quality of the erection and deflate to evaluate flaccidity. The penis should lie down close to the body when deflated. There may be some swelling that precludes a good flaccid result.

POSTOPERATIVE PROCEDURES

POSTOPERATIVE PROCEDURES

Immediately Postoperative

After the surgery, some physicians partially inflate the cylinders for the first 24 hours. This will aid hemostasis. The physician may place a closed system drain in the abdomen to drain excess fluid from the incision site.

After 24 hours, remove the dressing and completely deflate the cylinders. Support the penis on the abdomen for four to six weeks to obtain a straight erection. Taping the penis to the abdomen is the usual process.

After the Patient is Released from the Hospital

After the patient has returned home and the swelling from the surgery has subsided, the physician may ask the patient to pull down on the pump located in the scrotum to properly position it. Positioning the pump makes it easier for the patient to locate the pump.

The frequency of positioning the pump is the physician's decision. Some physicians have their patients position the pump several times daily.

To position the pump in the scrotum, a patient should:

1. Locate the pump in the scrotum.
2. Grasp the pump firmly and carefully pull the pump down in the scrotum. The patient should gently pull the pump into a position close to the outer scrotal wall.

After three to six weeks, the physician may instruct the patient to begin cycling the device for the first time. To cycle the device, the patient inflates and deflates the prosthesis several times. The cylinders must be fully inflated before they can be deflated.

It may be painful for the patient the first few times that he inflates and deflates the device.

However, after the postoperative healing period, the pain should subside. Instruct the patient to inflate and deflate the prosthesis several times daily.

Four to six weeks postoperatively, instruct the patient that it is possible to begin using the prosthesis to have gentle intercourse. To determine if the patient is ready to use the device:

1. Check the incision site to be sure that it has healed properly. There should be no redness, swelling, or drainage. Any of these things may indicate that an infection is present and the infection should be treated promptly with antibiotics.
2. Ask the patient about pain when cycling the device and observe the patient inflating and deflating the device.

After determining that the patient knows how to operate the device and that the device is functioning correctly, inform the patient that it is possible to begin intercourse.

Evaluating Long-Term Function and Placement


After the postoperative healing period, the physician should continue to have contact with the patient at least on an annual basis to evaluate the function of the device. During the annual evaluation, ask the patient how the device is functioning and if he has noticed any changes in the function, for example, the cylinders losing rigidity. Also check the patient for signs of infection or erosion.

If the patient is having mechanical difficulty with the device, or if there is infection or erosion present, revision surgery may be necessary.


PROCEDURE ILLUSTRATIONS

PROCEDURE ILLUSTRATIONS


1


2


3


4


5


6


7

8


9

PROCEDURE ILLUSTRATIONS (CONTINUED)


10


11


12


13

14

15

BRIEF SUMMARY

BRIEF SUMMARY

The AMS Ambicor Penile Prosthesis is intended for use in the treatment of chronic, organic, male erectile dysfunction (impotence). These devices are contraindicated in patients who have active urogenital infections or active skin infections in the region of surgery. Implantation will make latent natural or spontaneous erections, as well as other interventional treatment options, impossible. Men with diabetes, spinal cord injuries, or open sores may have an increased risk of infection. Failure to evaluate and treat device erosion may result in infection and loss of tissue. Implantation may result in penile shortening, curvature, or scarring. Possible adverse events include, but are not limited to, urogenital pain (usually associated with healing), patient dissatisfaction, mechanical malfunction, auto-inflation, penile curvature or sensation change, urogenital hematoma, urogenital edema, and infection.

Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, and potential adverse events. Rx Only.

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0086

AMS™

©2016 Boston Scientific Corporation or its affiliates. All rights reserved. All trademarks are the property of the respective owners.

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Penile Prosthesis Product Line

Instructions For Use

NOTE: Refer to the Operating Room Manual for further information on the AMS 700™ product line and their implantation.

Device Description

The AMS 700 Series Inflatable Penile Prosthesis product line includes the AMS 700 LGX™ Preconnect, AMS 700 CX Preconnect, AMS 700 LGX, AMS 700 CX, and the AMS 700 CXR Preconnect Penile Prosthesis. These configurations are available with and without InhibiZone™, an antibiotic impregnation of rifampin (rifampicin) and minocycline hydrochloride (HCl).^{*} Use of an InhibiZone treated penile prosthesis results in a significant reduction in the rate of revisions due to infection in patients receiving a first time or a revision AMS 700 device compared to an untreated AMS 700 device (this benefit is also seen among high risk diabetic patients receiving a first time AMS 700 implant). This clinical data demonstrates that the use of InhibiZone reduces device-related infections. These prostheses are closed fluid-filled systems consisting of a pair of cylinders, optional rear tip extenders (RTEs), a Momentary Squeeze (MS) pump, and a fluid reservoir. All components are connected by kink-resistant tubing (KRT). The cylinders are inflated as fluid is pumped from the reservoir, creating an erection. They are deflated as fluid is transferred back to the reservoir, making the penis flaccid once again. This device contains solid silicone elastomer. These devices are for men who, after appropriate patient history and diagnostic evaluations as well as discussions with the urologist about other alternative treatment methods, are determined to be suitable candidates for implantation surgery.

^{*} not available in all markets.

Indications for Use

The AMS 700 Series Inflatable Penile Prosthesis product line is intended for use in the treatment of chronic, organic, male erectile dysfunction (impotence).

Contraindications

The implantation of this device is contraindicated in patients who have active urogenital infections or active skin infections in the region of surgery.

The implantation of the InhibiZone version of this device is contraindicated in patients with known allergy or sensitivity to rifampin (rifampicin), or to minocycline HCl, or other tetracyclines.

The implantation of products with InhibiZone is contraindicated in patients with systemic lupus erythematosus because minocycline HCl has been reported to aggravate this condition.

Warnings

1. Implantation of the device will make latent natural or spontaneous erections, as well as other interventional treatment options, impossible.
2. Men with diabetes, spinal cord injuries, or open sores may have an increased risk of infection associated with a prosthesis.
3. Failure to evaluate and promptly treat erosion may result in a substantial worsening of the condition leading to infection and loss of tissue.
4. Implantation of a penile prosthesis may result in penile shortening, curvature or scarring.
5. This device contains solid silicone elastomer. The risks and benefits of implanting this device in patients with documented sensitivity to silicone should be carefully considered.
6. Pre-existing abdominal or penile scarring or contracture may make surgical implantation more complicated or impractical.
7. If a hypersensitivity reaction develops to a device coated with InhibiZone, the penile prosthesis should be removed and the patient treated appropriately.

Precautions

Surgery Related

1. Improper reservoir placement or filling technique can result in spontaneous unintended inflation or deflation of the cylinders that may result in unintended partial or full erections.
2. Migration of the device components can occur if the cylinders are improperly sized, if the pump or the reservoir is not positioned properly, or if the tubing lengths are incorrect.
3. Removal of an implanted prosthesis without timely reimplantation of a new prosthesis may complicate subsequent reimplantation or may make it impossible.
4. Improper measurement technique, positioning or sizing may reduce cylinder life.
5. Unsuccessful outcomes have been reported due to improper surgical technique, anatomical misplacement of components, improper sizing and filling of components, or tubing kinks.
6. Implantation of AMS 700 LGX cylinders in patients with Peyronie's disease may not provide a satisfactory result.

Device Related

1. AMS Quick Connect Sutureless Window Connectors, provided in the AMS 700 Accessory Kit, should not be used in revision procedures involving previously implanted component tubing. The Quick Connect Sutureless Window Connectors may be less effective.
2. The materials used in the construction of this device have been shown to cause minor irritation when implanted in animals. Therefore, implantation of this device may cause minor irritation or discomfort in some patients.
3. Devices in the AMS 700 product line should be filled with sterile, normal saline. Some patients may have a hypersensitivity to contrast media.
4. Do not use product that has damaged or open packaging, as sterility may be compromised.
5. Devices with InhibiZone should not come into contact with ethyl alcohol, isopropyl alcohol or other alcohols, acetone or other nonpolar solvents. These solvents may remove the antibiotics from the device.
6. InhibiZone components should not be soaked in saline or other solutions prior to implantation. The components may be briefly rinsed or dipped in a sterile solution immediately prior to implant, if desired.
7. CXR RTEs are not compatible with CX or LGX cylinders.
8. Verify proper attachment of RTEs by spinning them once seated. Properly attached RTEs should spin freely without accidental disengagement or material bulging.

9. Do not stack the CX, LGX or CXR Snapcone RTEs with the exception of the 1.5 cm. The locking ring will not engage with the smooth outer surface of the RTE, which may result in the RTE disconnecting.

Patient Related

1. A thorough preoperative consultation should include a discussion between patient and physician of all available treatment options and their risks and benefits.
2. Adequate patient manual dexterity and strength are required for proper device inflation and deflation.
3. Mental or psychological conditions, such as senile dementia, may inhibit the patient's successful operation of the prosthesis.
4. Trauma to the pelvic or abdominal areas, such as impact injuries associated with sports (eg, bicycle riding), can result in damage of the implanted device and/or surrounding tissues. This damage may result in the malfunction of the device and may necessitate surgical correction, including replacement of the device.
5. The contour, elasticity, and dimension of the tunica albuginea may limit the length and/or diameter expansion of the AMS 700 cylinders.
6. The implantation of this device should only be considered in patients whom the physician determines are adequate surgical candidates.
7. Use of injection therapy concurrently with the penile prosthesis can damage the prosthesis. Patients should not use injection therapy after they receive their implant.

InhibiZone Related

1. InhibiZone does not replace your normal antibiotic protocols. Continue using any prophylactic protocols normally used when implanting an inflatable penile prosthesis.
2. Because the products with InhibiZone are impregnated with a combination of rifampin (rifampicin) (a derivative of rifamycin B) and minocycline HCl (a derivative of tetracycline), the contraindications, warnings, and precautions regarding the use of these antimicrobial agents apply and should be adhered to for the use of this device, although systemic levels of minocycline HCl and rifampin (rifampicin) in patients receiving this device are unlikely to be detected.
3. Use of products with InhibiZone should be carefully considered in patients with hepatic or renal disease, as use of rifampin (rifampicin) and minocycline HCl can cause additional stress on the hepatic and renal systems.
4. Patients who receive a device with InhibiZone and are also taking methoxyflourane should be carefully monitored for signs of renal toxicity.
5. Patients who receive a device with InhibiZone and are also taking warfarin should have their prothrombin time monitored, because tetracyclines have been reported to slow coagulation.
6. Use of products with InhibiZone should be carefully considered in patients using thionamides, isoniazid, and halothane, due to potential hepatic side effects that have been reported in patients using these drugs and higher doses of rifampin (rifampicin).

Adverse Events

A clinical trial was conducted to determine the safety and effectiveness of the AMS 700 Series of inflatable penile prostheses. This trial involved only devices without InhibiZone. A total of 300 patients were enrolled with follow-up out to 5 years for 126 patients. The Adverse Device Effects (ADE), detailed in the following table, were noted during the duration of this clinical trial for all enrolled patients.

Table 1: AMS Clinical Trial Adverse Device Effects

ADE	# Patient (%**)	Mean Onset Time in Days (Range in Days)
Urogenital Pain (Typically Associated with Healing Process)	160 (53.3%)	21 (0 – 876)
Urogenital Edema	106 (35.3%)	8 (0 – 722)
Urogenital Ecchymosis	30 (10.0%)	4 (0 – 150)
Reservoir Encapsulation (persistent in 11/19 cases)	19 (6.3%)	275 (38 – 1731)
Patient Dissatisfaction (With Length, Ability to Use and Nonspecific Reasons)	18 (6.0%)	384 (0 – 1830)
Auto-Inflation	17 (5.7%)	141 (0 – 608)
Mechanical Malfunction (Leaks, Incomplete Inflation/Deflation, Kinking)	13 (4.3%)	905 (0 – 1915)
Urination Impaired (Slow Stream, Split Stream, Voiding Difficulties or Obstructive Symptoms)	11 (3.7%)	239 (0 – 930)
Urogenital Erythema	10 (3.3%)	36 (0 – 320)
Joint Pain, Swelling, or Stiffness	9 (3.0%)	609 (1 – 1592)
Decreased Penile Sensation	7 (2.3%)	124 (0 – 214)
Urogenital Hematoma	7 (2.3%)	4 (0 – 25)
Abnormal Ejaculation (Delayed, Burning, or General Nonspecific Problems)	6 (2.0%)	409 (40 – 1797)
Infection	6 (2.0%)	216 (9 – 716)
Dysuria	5 (1.7%)	231 (2 – 684)
Penile Curvature	5 (1.7%)	144 (0 – 257)
Application Site Reaction (Wound Separation, Delay in Cutaneous Closure)	4 (1.3%)	14 (0 – 30)
Erosion/Extrusion (Pump/Cylinder)	4 (1.3%)	425 (72 – 1066)
Urogenital Inflammation	4 (1.3%)	490 (0 – 1897)
Adhesion of the Pump/Scrotum	3 (1.0%)	13 (10 – 19)
Device Malposition	3 (1.0%)	278 (43 – 574)
Device Migration (Pump/Cylinder)	3 (1.0%)	210 (40 – 548)
Transient Urinary Retention	3 (1.0%)	85 (3 – 248)
Urinary Frequency	3 (1.0%)	277 (99 – 409)
Weakness	3 (1.0%)	1072 (519 – 1592)
Abnormal Sexual Function	2 (0.7%)	239 (128 – 349)
Device Cylinder Aneurysm/Bulge	2 (0.7%)	945 (110 – 1780)
Dizziness	2 (0.7%)	929 (7 – 1850)
Dry Mouth	2 (0.7%)	1721 (1592 – 1850)
Hematuria	2 (0.7%)	902 (13 – 1791)

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ADE	# Patient (%**)	Mean Onset Time in Days (Range in Days)
Low Grade Fever	2 (0.7%)	13 (7 – 18)
Memory Difficulties	2 (0.7%)	1318 (1107 – 1592)
Pelvic Pain	2 (0.7%)	270 (42 – 498)
Rheumatoid Arthritis	2 (0.7%)	281 (189 – 372)
Other	22 (7.0%)	N/A

** Percentages based on total number of patients implanted (300).

Table 2: Resolution of Adverse Device Effects

ADE	Method of Resolution					
	# Patient (%)	Surgical	Medication	Other ¹	No Medical Intervention	Continuing
Urogenital Pain	160 (53.3%)	1%	31%	0%	68%	11%
Urogenital Edema	106 (35.3%)	0%	3%	11%	86%	2%
Urogenital Ecchymosis	30 (10.0%)	0%	0%	0%	100%	0%
Reservoir Encapsulation	19 (6.3%)	5%	0%	0%	95%	58%
Patient Dissatisfaction	18 (6.0%)	0%	0%	0%	100%	56%
Auto-Inflation	17 (5.7%)	0%	0%	0%	100%	35%
Mechanical Malfunction	13 (4.3%)	46%	0%	8%	46%	62%
Urination Impaired	11 (3.7%)	0%	64%	9%	27%	0%
Urogenital Erythema	10 (3.3%)	10%	30%	0%	60%	0%
Joint Pain, Swelling, or Stiffness	9 (3.0%)	0%	11%	11%	78%	67%
Decreased Penile Sensation	7 (2.3%)	0%	0%	0%	100%	72%
Urogenital Hematoma	7 (2.3%)	0%	0%	0%	100%	0%
Abnormal Ejaculation	6 (2.0%)	0%	17%	0%	83%	17%
Infection	6 (2.0%)	67%	33%	0%	0%	17%
Dysuria	5 (1.7%)	0%	60%	0%	40%	0%
Penile Curvature	5 (1.7%)	0%	0%	0%	100%	60%
Application Site Reaction	4 (1.3%)	0%	25%	0%	75%	25%
Erosion/Extrusion	4 (1.3%)	100%	0%	0%	0%	0%
Paresthesia	4 (1.3%)	0%	0%	0%	100%	50%
Urogenital Inflammation	4 (1.3%)	0%	50%	0%	50%	0%
Adhesion of Pump/Scrotum	3 (1.0%)	0%	0%	0%	100%	33%
Device Malposition	3 (1.0%)	67%	0%	0%	33%	0%
Device Migration	3 (1.0%)	100%	0%	0%	0%	0%
Transient Urinary Retention	3 (1.0%)	0%	0%	100%	0%	0%
Urinary Frequency	3 (1.0%)	0%	33%	0%	67%	67%
Weakness	3 (1.0%)	0%	0%	67%	33%	67%
Abnormal Sexual Function	2 (0.7%)	0%	0%	0%	100%	100%
Device Cylinder Aneurysm/Bulge	2 (0.7%)	50%	0%	0%	50%	50%
Dizziness	2 (0.7%)	0%	0%	0%	100%	50%
Dry Mouth	2 (0.7%)	0%	0%	0%	100%	100%
Hematuria	2 (0.7%)	0%	50%	0%	50%	50%
Low Grade Fever	2 (0.7%)	50%	50%	0%	0%	0%
Memory Difficulties	2 (0.7%)	0%	0%	0%	100%	0%
Pelvic Pain	2 (0.7%)	0%	0%	0%	100%	50%
Rheumatoid Arthritis	2 (0.7%)	0%	0%	0%	100%	100%
Other	22 (7.0%)	N/A	N/A	N/A	N/A	N/A

1 Other treatments included back brace, physical therapy, urine culture, ice packs, elevation, hot soaks, hot sitz bath, manual manipulation, patient education, filliforms and followers, foley catheter, ultrasound/CT scan, and cystoscopy.

The following “Other” adverse device effects (in alphabetical order) each occurred in less than 0.5% of the patients: Alopecia, Back Pain, Cellulitis, Depression, Diabetes Mellitus, Epigastric Pain, Eye Disorder, Eye Pain, Fecal Incontinence, Fibrosis, Glans Hypermobile Dorsally, Kidney Calculus, Libido Decreased, Migraine, Necrosis, Phimosi, Photosensitivity Reaction, Pump Fixation, Thickening of the Skin, Urinary Tract Infection, Urinary Urgency, and Vertigo.

The following adverse events (in alphabetical order) have been associated with the use of this product: Bleeding, Cavernosal Crossover, Exposure to Biohazardous Material, Granuloma Formation, Ischemia, Non-Rheumatoid Arthritis Immune-Related Tissue Disorders, Pain (which may be prolonged or severe), Perforation or Injury of Bladder, Corpus Cavernosum, Nerve, Tunica, or Urethra, Prolonged Procedure, Seroma, Thrombosis, Ulceration, Unretrieved Device Fragment, Vascular Compromise, and Penile Trauma.

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There were 18 patient deaths during the course of the trial. No deaths that occurred during the duration of the clinical study were attributed to the device implantation or use.

A total of 22 patients underwent revision surgeries in the 5 year study period. Information on device revisions is described in the Clinical Studies section.

Clinical Studies

A clinical trial was undertaken to demonstrate that the AMS 700 product line provides an erection that is suitable for intercourse and has acceptable rates of surgical revision and of significant clinical events associated with the implantation and use of these devices. This trial included only devices **without** InhibiZone. This trial was also designed to demonstrate the implantation of these devices does not negatively impact the sexual satisfaction, psychological well-being, self-esteem or quality of life of patients who receive these devices. It was a prospective, multicenter cohort trial in which the patients served as their own control. The choice of device model implanted (ie, 700 CX, 700 CX Preconnect, 700 CXM, 700 Ultrex™, 700 Ultrex Plus) was at the discretion of the patient and implanting physician.

NOTE: The AMS 700 MS Pump was not available at the time the clinical study was conducted. However, based on the similarities between the AMS 700 MS Pump and the AMS 700 Inflate/Deflate Pump, the clinical results also apply to this new model.

NOTE: The AMS 700 LGX Preconnect was not available at the time the clinical study was conducted. However, since providing the AMS 700 LGX in preconnected form is not expected to affect the safety and effectiveness of the prosthesis, these clinical results also apply to the new model.

NOTE: The AMS 700 CXR was not available at the time the clinical study was conducted. However, based on the similarities between the AMS 700 CXR and the AMS 700 CXM models, the clinical results also apply to this new model.

NOTE: The AMS Conceal™ Low Profile Reservoir* was not available at the time the clinical study was conducted. However, based on similarities between the AMS Conceal Low Profile Reservoir and the spherical reservoir, the clinical results also apply to this new model.

* not available in all markets.

Three hundred male patients, over 21 years of age, were enrolled in this study. All patients with diagnosed organic erectile dysfunction were eligible for enrollment, if they did not present with a history of allergy/sensitivity to silicone, preexisting autoimmune or connective tissue diseases, or active urogenital infection.

All safety-related data, diagnoses, and health status evaluations were captured on detailed case report forms. The investigators' professional evaluation of the erections provided by the IPPs (Inflatable Penile Prostheses) after implantation and their suitability for intercourse was the primary efficacy endpoint. The number of surgical revisions performed and reported by the investigators was the primary safety endpoint. Patient self-evaluations on 4 validated outcome instruments were the secondary efficacy endpoints (concerning quality of life, self-esteem, and sexual satisfaction and functioning).

This clinical trial provided the following results through the 5 year evaluation for the first 126 patients to reach this postsurgical follow-up.

Physician Assessment of Device Function

One hundred twenty-six devices were evaluated at the 5 year follow-up, of which 123 (97.6%) could be inflated. Of these 123 devices, all (100%) were determined to provide an erection suitable for intercourse. However, it should be noted that this analysis does not include the following information regarding device malfunctions: (i) 3 of the 123 devices found to be functioning properly at the 5 year exam were surgically revised prior to this exam to correct a mechanical malfunction and (ii) 3 additional devices not evaluated at the 5 year follow-up exam were also surgically revised due to mechanical malfunction. These cases of device revision are discussed further in the next section.

Surgical Revisions

The incidence of revisions was evaluated in the 126 patients with follow-up out to 5 years, as well as 16 additional patients who experienced one or more revision surgeries and did not reach the 5-year follow-up exam. (A revision is considered any urogenital surgical intervention that is related to the function, placement, or site reaction to the implanted device.) Of these 142 patients, 22 (15.5%; 95% confidence interval = 21.5%) experienced a total of 26 revision surgeries, and 120 (84.5%) were not revised.

The average time to the first revision surgery was 15 months (ranging from 0.9 months to 60.1 months). Of the 26 revision surgeries, there were 5 revisions due to "Infection," 2 for "Infection/ Erosion," 2 for "Migration/Malposition," 2 for "Erosion," 2 for "Malposition," 7 for "Mechanical Malfunction," 2 for "Fibrous Capsular Complication," 2 for "Reimplantation Following Previous Revision," and 2 due to reasons listed as "Other." The "Other" reasons for revision surgeries included 1 for Cylinder kink/auto-inflation and 1 for Corporal body aneurysm. In 5 of these revision surgeries, no device components were explanted or replaced. The components were manipulated/repositioned but were not removed.

Patient Evaluation of Quality of Life, Self-esteem, Psychological and Sexual Well-Being

In accordance with the study protocol, overall health-related quality of life (using the Medical Outcomes Study Health Survey, MOS-20), self-esteem (using the Rosenberg Self-esteem Scale), psychological well-being (using the Brief Symptom Inventory), and sexual functioning and satisfaction (using the Sexual History Form) were evaluated in patients through 2 years post-implantation. Throughout the 2-year follow-up period, patient quality of life, self-esteem, and psychological well-being were determined to be equivalent to the preimplant state. Sexual functioning and sexual satisfaction, on the other hand, was significantly improved over the preimplant state.

Supplementary Clinical Information

Although it is not feasible to predict exactly how long an implanted penile prosthesis will function in a particular patient, American Medical Systems (AMS), Inc., has assembled a set of data on device removals and revisions to help gain insight into the product performance over time.

The following 2 tables provide an estimate of the long-term rates of device removals and revisions for Ultrex and CX models. The first data set comes from Patient Information Forms (PIFs) submitted to AMS by physicians for surgical procedures requiring parts replacement under AMS warranty (Table 3). All forms reporting devices implanted between January 1993 and December 2000 were included in a life table analysis that was used to calculate the revision rates for each category. Revision surgeries may not be reported to AMS. Therefore, the incidence of surgeries after original implant would likely be underreported, if one were to rely solely on the PIF data. AMS also assembled a second set of data directly from a retrospective review of physicians' medical records (Table 4). These medical records capture each surgery performed by that physician after original implant for any reason.

NOTE: These analyses did not include the AMS 700 with MS Pump, snapcone cylinders, or parylene coated components.

NOTE: These analyses did not include the AMS Conceal Low Profile Reservoir.

PIF Study

Table 3: Revision rates based on PIF data:*

<i>Reason for Removal or Replacement Surgery</i>	<i>700 CX Revision Rate (5 YRS) n=12,080</i>	<i>700 Ultrex Revision Rate (5 YRS) n=20,438</i>
Medical Device Revision	5.6%	4.8%
Регистрация МИ в Росздравнадзоре www.nevacert.ru info@nevacert.ru	2.2%	2.0%
Removal for Erosion	1.4%	1.0%
Migration/Malposition of Component	0.5%	0.5%
Cylinder Aneurysm	0.5%	0.5%
Other Reasons	3.6%	2.9%
ALL REASONS	11.6%	9.8%

* Interpretations of the PIF data may be limited by a number of factors:

- Statistics are based only on surgery data voluntarily reported to AMS by hospitals and physicians in the U.S. as part of the AMS Product Replacement. Because surgeries may not be reported to AMS, the number of patients implanted and the incidence of removal/replacement surgery may actually be higher.
- These statistics pertain only to the incidence of removal/replacement surgery and not to the current functioning of devices that have not been removed.

NOTE: The amount for "ALL REASONS" is lower than the total of the individual percentages due to the incidence of multiple reasons for removal/replacement surgery.

Medical Records Study

Table 4: Revision rates based on medical records data*

Reason for Removal or Replacement Surgery	700 CX Revision Rate (5 YRS) n=512	700 Ultrex Revision Rate (3 YRS) n=155
Mechanical Revision	9.4%	3.2%
Removal for Infection	2.9%	2.0%
Removal for Erosion	2.4%	1.6%
Migration/Malposition of Component	2.8%	0.7%
Cylinder Aneurysm	0.6%	0.0%
Other Reasons	2.1%	4.0%
ALL REASONS	18.1%	10.4%

- * Interpretations of the medical record data may be limited by several factors:
- These percentages reflect the known revisions performed by the original implant physicians.
 - These statistics pertain only to the incidence of removal/replacement surgery and not to the current functioning of devices which have not been removed.

NOTE: The amount for "ALL REASONS" is lower than the total of the individual percentages due to the incidence of multiple reasons for removal/replacement surgery.

Antibiotic Information

The antibiotics present in InhibiZone, minocycline HCl and rifampin (rifampicin), are well characterized and have been in use for years. The dosage present on the penile prostheses is intended to act on organisms that attempt to colonize the device. The AMS 700 components are treated with very low levels of antibiotics. AMS provides numerous completed configurations of the AMS 700 to individualize treatment; however, a complete device (reservoir, pump, and 2 cylinders), regardless of configuration, contains ≤33 mg rifampin and ≤12 mg minocycline HCl which represents less than 2% of oral dose exposure for a complete course of rifampin or minocycline HCl with the maximum dose calculated for the most common device configuration's average concentration plus (1) standard deviation.

The following in vitro data describe the Zones of Inhibition provided by samples from InhibiZone treated devices.

Table 5: In vitro zones of inhibition for device samples* with InhibiZone treatment

Organism	Mean (mm)	S.D. (mm)	Number of Isolates
<i>Staphylococcus epidermidis</i>	22.6	2.9	21
<i>Staphylococcus aureus</i>	17.5	5.0	25
<i>Escherichia coli</i> **	6.5	2.6	24
<i>Enterococcus faecalis</i> **	4.8	6.7	21
<i>Candida albicans</i> **	0.1	0.4	21
<i>Proteus mirabilis</i> **	0.6	1.0	17

- * Obtained using standardized KRT test samples containing 12 µg minocycline and 26 µg rifampin (rifampicin)
- ** The isolates tested were not susceptible to rifampin (rifampicin) and/or minocycline control disks

An animal infection study was conducted using 11 rabbits. Five rabbits were implanted subcutaneously with 6 test samples each and 5 rabbits were implanted subcutaneously with 6 control samples each. One rabbit received 3 test samples and 3 control samples. The test samples were portions of an InhibiZone treated AMS 700 Pump and the control samples were portions of a standard AMS 700 Pump without InhibiZone. All samples were soaked in a 103-104 colony-forming units solution of *Staphylococcus aureus*, Sheretz strain for 8 hours. Samples were then allowed to dry for 30 minutes prior to surgical placement in the rabbit. After 2 days, all samples were removed and observed for growth on the samples. The number of coated samples that were infected was statistically significantly lower than the number of control samples that were infected.

Antibiotic Human Clinical Experience

A post marketing study was conducted to evaluate the effect of the antibiotic surface treatment on reducing the incidence of penile implant infections, which consisted of a comprehensive review of more than 43,000 patients in the AMS PIF database. This study shows a significant improvement in the rate of revisions due to infection for patients with original AMS 700 InhibiZone (IZ) implants or revision AMS 700 IZ implants compared to those who received AMS 700 devices without the IZ treatment. Further, patients with diabetes who received original AMS 700 IZ implants also demonstrate significant improvements in the rate of revisions due to infections, compared to those who received AMS 700 devices without the IZ treatment. The data set included PIF data for patients who underwent an AMS 700 implant surgery in the United States from May 1, 2001, through September 1, 2007. The most common reason for implantation of an AMS 700 device was recorded as "organic" on the PIFs (Table 6). Post prostatectomy, vascular disease and diabetes were also commonly reported as the etiology leading to the need for the AMS 700 implants.

Table 6: Etiology for AMS 700 Implant Patients

Description	IZ		Non-IZ		Mixed	
	#	%	#	%	#	%
Original AMS 700 Implant						
Diabetes Mellitus*	4,977	17.7%	627	17.7%	177	18.1%
Peyronie's Disease	2,553	9.1%	328	9.2%	102	10.4%
Post Prostatectomy	7,215	25.7%	819	23.1%	242	24.7%
Vascular Disease	4,037	14.4%	576	16.2%	120	12.3%
Other	1,107	3.9%	151	4.3%	50	5.1%
Total	8,185	29.2%	1,046	29.5%	288	29.4%
	28,074	100.0%	3,547	100.0%	979	100.0%
Revision or Replacement AMS 700 Implant						
Diabetes Mellitus	1,979	16.2%	349	15.4%	63	15.2%
Peyronie's Disease	1,243	10.2%	208	9.2%	40	9.6%
Post Prostatectomy	2,016	16.5%	330	14.5%	68	16.4%
Vascular Disease	2,566	21.0%	599	26.4%	101	24.3%
Other	1,034	8.5%	184	8.1%	34	8.2%
Organic	3,378	27.7%	602	26.5%	109	26.3%
Total	12,216	100.0%	2,272	100.0%	415	100.0%
TOTAL	40,290		5,819		1,394	

NOTE: None or more than one etiology could be selected for each patient; therefore, the numbers of patients within each subgroup by etiology above may not be identical to the number of patients in each life table analysis.

- * The diabetes subgroup with an original implant was identified and analyzed as a specific high-risk subgroup.

Three patient populations were analyzed in this PIF study: the first group included all patients with an AMS 700 original implant between May 1, 2001, and September 1, 2007. A subgroup of these AMS 700 original implant patients included all patients with diabetes. The third group included all patients with an AMS 700 revision implant (removing or replacing all or part of an earlier implanted device) between May 1, 2001, and September 1, 2007. An event was recorded for any revision surgery with a surgery reason recorded as infection, mechanical malfunction, fluid loss, erosion, patient dissatisfaction, or other. Revision surgeries were tracked through December 1, 2007, when the reporting period was closed. When multiple reasons were listed for 1 revision surgery (eg, infection, device malfunction, erosion, etc), the individual revision surgery was counted and included in the data for each of the indicated revision reasons. Thus, all revisions that report some presence of infection will be captured in the infection category presented in this study (even if multiple other reasons for the revision surgery are listed).

Life tables were constructed and the logrank test was used to compare the device survival distributions between the Treatment (InhibiZone) and Control (non-InhibiZone) groups. All significance tests were performed at an alpha level of 0.05. For original implants, device survival distributions were defined as the time between the implantation of the original device in the study period and the first replacement, revision, or removal of any component of the device due to infection, malfunction, fluid loss, erosion, patient dissatisfaction, or other reason (the Original Implant Group). For revision surgery patients, device survival distributions were created for the time between the implantation of the revision device and the next revision device due to infection, malfunction, fluid loss, erosion, patient dissatisfaction, or other reason (the Revision Surgery Group). The study also included a subgroup of patients who received a penile implant after spinal cord injury. However, this aspect of the study continues to be unknown, since there were too few patients in this group to be able to determine valid conclusions.

Original Implant Patient Group

Among the original implant patients, a total of 31,328 (98.7%) patient devices did not have any events reported due to infection and were censored during the timeframe of this report. A total of 338 (1.2%) device events were due to infection in the IZ treated group while 84 (2.5%) were due to infection in the non-IZ treated group. Life table analyses compared the rates of revision surgery due to infection for original IZ or non-IZ treated implants. In this device survival analysis with infection as the reason for the revision surgery, InhibiZone treated AMS 700 devices had a significantly better survival than non-InhibiZone treated devices among original implant patients ($p < 0.0001$) (Figure 1).


Figure 1: Device survival—freedom from revision due to infection among original implant patients

Diabetic Original Implant Patient Group

Among the subgroup of original implant patients who had diabetes, a total of 5604 diabetic patients included 4977 (88.8%) with IZ treated devices and 627 (11.2%) with non-IZ treated devices. A total of 5508 (98.3%) diabetic patient devices did not have any events due to infection during the timeframe of this report. Only 71 (1.4%) infection device events were in the IZ treated group while 25 (4.0%) events were in the non-IZ treated group. Life table analyses compared rates of revision surgery due to infection for diabetic patients. InhibiZone treated devices had significantly better survivals than non-InhibiZone treated devices among diabetic patients ($p < 0.0001$) (Figure 2).


Figure 2: Device survival—freedom from revision due to infection among diabetic patients

Revision Surgery Implant Group

Among the revision surgery implant patients, a total of 10,769 (97.3%) patient devices did not have any events due to infection and were censored during the timeframe of this report. Only 229 (2.5%) infection device events occurred in the IZ treated group while 66 (3.7%) infection events occurred in the non-IZ treated group. Life table analysis compared the rates of revision surgery due to infection for IZ or non-IZ treated implants among revision surgery implant patients. In this device survival analysis with infection as the reason for the revision surgery, InhibiZone treated AMS 700 devices had a significantly better survival than non-InhibiZone treated devices among revision surgery implant patients ($p=0.0252$) (Figure 3).


Figure 3: Device survival—freedom from revision due to infection among revision surgery implant patients

In conclusion, this human clinical evidence from the InhibiZone PIF study indicates the revisions due to infection were significantly lower among patients who received an AMS 700 InhibiZone treated device when compared to a nontreated AMS 700 device. This reduced rate of revisions subsequent to infection for the AMS 700 InhibiZone treated device when compared to a nontreated AMS 700 device was reported for original as well as revision implants and included the high-risk subgroup of original implant patients who had diabetes. The study supports the following conclusions:

- No significant differences were found between InhibiZone treated and untreated AMS 700 penile prostheses for overall rates of revision due to mechanical malfunction, fluid loss, erosion, patient dissatisfaction, or all other reasons.
- Use of InhibiZone results in a significant reduction in the rate of revision due to device-related infection among penile prosthesis implants.
- Use of InhibiZone results in a significant reduction in the rate of revision due to device-related infection among penile prosthesis implants in diabetic patients.
- Use of InhibiZone results in a significant reduction in the rate of revision due to device-related infection in original and revision penile prosthesis implants.

Patient Counseling Information

Patients should be counseled in order to have a realistic expectation of the physical, psychological, and functional outcome of the implantation. The risks, benefits, and potential adverse events of all available treatment options should be discussed with the patient and considered by the physician and patient when choosing a treatment option.

An appropriate patient history, including history of personality disorders and a diagnostic work-up, should be a part of the patient decision-making process.

Some patients may become dissatisfied by the presence of the prosthetic device in their body. This issue should be discussed with the patient prior to the surgery. Patient dissatisfaction may lead to device removal.

Implantation of a penile prosthesis may result in penile shortening, curvature, or scarring. The prosthetic erection may differ from the patient's original, natural erection in that it may be shorter, less firm, have less girth, and reduced sensations. Realistic cosmetic expectations should be communicated to the patient and should include the potential for skin scarring, scrotal deformity, pump bulge in the scrotum, lack of concealability, and other possible adverse events. Patients should also be aware that penile prostheses are not considered to be lifetime implants.

Improper implantation of a penile prosthesis may not provide rigidity to the glans, which may result in a floppy glans and may result in a lack of rigidity of the corpus spongiosum. Penile flaccidity may be less than it was before implantation.

Patients who undergo revision surgery may notice a change in the character of their erection compared to their previous implant, which may include differences in sensation, length, girth, rigidity, and/or flaccidity.

It is also important that the physician discusses with the patient the possibility of an allergic reaction to the device (see Silicone Information section).

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Silicone Information

This device is composed of a number of materials, including solid silicone elastomers and a fluorosilicone lubricant. Silicone gel is not a component in the materials of this device.

Silicone elastomers have been commonly used in a variety of biomedical devices for over 40 years and are used as a biocompatibility reference against which new materials are tested. Silicone fluids have an extensive history of use in medical devices.

Scientific literature has included reports of adverse events and other observations in patients with implantable silicone devices. As reported, these events/observations indicate "allergic-like" symptoms and in other cases a symptom complex associated with immunological disorders. No causal relationship has been established between these events and silicone elastomer or fluorosilicone lubricant.

There are reports of malignant tumor formation in laboratory animals only, not humans, associated with implants of relatively large size. Many different materials are associated with this effect in animals, silicone elastomers among them. No such effect has been described in humans.

Extensive testing has been conducted on all materials that comprise the prostheses in the AMS 700. This testing has indicated no toxicological response attributable to the materials. However, some of the materials caused minor irritation when implanted in animals.

Silicone elastomer particulate shedding and particulate migrations to regional lymph nodes have been reported in the literature on penile implants. There are no known clinical sequelae to this phenomenon.

Magnetic Resonance Imaging (MRI)
Important Safety Information

MR Conditional

The AMS 700 Series Three-Piece Inflatable Penile Implant was determined to be MR Conditional. Nonclinical testing demonstrated that the AMS 700 Series Three-Piece Inflatable Penile Implant is MR Conditional. A patient with this device can be safely scanned immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3 T or less
- Maximum spatial gradient magnetic field of 720 Gauss/cm or less

MRI-Related Heating

In nonclinical testing, the AMS 700 Series Three-Piece Inflatable Penile Implant produced the following temperature rise during MRI performed for 15 minutes of scanning (ie, per pulse sequence) in the 3-Tesla (3 T/128 MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) magnetic resonance (MR) system:

Highest temperature change +1.9°C

Therefore, the MRI-related heating experiments for the AMS 700 Series Three-Piece Inflatable Penile Implant at 3 T using a transmit/receive radio frequency (RF) body coil at an MR system reported whole body averaged specific absorption rate (SAR) of 2.9 W/kg (ie, associated with a calorimetry measured whole body averaged value of 2.7 W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.9°C.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the AMS 700 Series Three-Piece Inflatable Penile Implant. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	6,244 mm ²	1,589 mm ²	10,295 mm ²	2,779 mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

Inventory Returns and Product Replacement Information

In the United States

Before returning any components, whether explanted or unused (sterile or nonsterile), customers must fill out the Return Goods Form located on the last page of the PIF.

Follow **all** of the instructions on the form carefully, and be sure that the components have been thoroughly cleaned before returning them to AMS.

In all cases, obtaining credit or percentage of credit for a returned component is subject to approval under the terms of the AMS Return Goods Policy and the AMS Product Warranty Policy. For complete information regarding these policies, contact the AMS Customer Service Department.

Outside the United States

Customers outside of the United States should contact their local AMS Representative prior to returning any product.

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How Supplied and Storage

WARNING: Contents supplied STERILE. Do not use if sterile barrier is damaged. If damage is found, or if labeling information appears missing, illegible, or compromised, call your AMS representative for clarification, and do not use the AMS 700 device until integrity of the product has been verified.

For single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the AMS 700 device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the AMS 700 device and/or cause patient infection or cross-infection including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the AMS 700 device may lead to injury, illness, or death of the patient. After use, dispose of the AMS 700 device and packaging in accordance with hospital, administrative, and/or local government policy.

Store product in a clean, dry, dark area at room temperature.

Do not use device beyond the indicated expiration date.

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92162915-01 (2017-08)

92162915-01

Information and instructions for patients considering an

AMS 700TM with MS PumpTM Series

Inflatable Penile Protheses

English

AMSTM

Contact list

American Medical Systems, Inc.

10700 Bren Road West
Minnetonka, MN 55343
U.S.A.
U.S. Toll Free: 1 800 328 3881
Tel: +1 952 930 6000
Tel: +31 20 593 8800

American Medical Systems Europe B.V.

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Table of contents

1	General warnings and cautions	1
1.1	How safe are silicone elastomer prostheses?	1
1.2	Will my prosthesis have to be replaced?	1
1.3	Need for manual dexterity	2
1.4	Possibility of malfunction	2
1.5	Possibility of changes in the penis or scrotum	2
1.6	Destruction of any latent natural erection	2
1.7	Different erection	2
1.8	Less softness when deflated	2
1.9	Possibility of infection	2
1.10	Delay intercourse until your doctor gives the OK	2
1.11	Erosion	3
1.12	Pain	3
1.13	Migration	3
1.14	Allergic reaction	3
2	What is impotence?	4
2.1	How does an erection occur?	5
2.2	Impotence is common	7
2.3	Impotence is treatable	7
3	Causes of impotence	8
3.1	Psychological causes	8
3.2	Physical causes	8
4	Am I impotent?	10
4.1	Impotence self-test	10
5	What to expect at your examination	12
5.1	Physical examination	12
5.2	Tests you may be asked to take	13
6	Treatment options for physical causes of impotence	14
6.1	Medication	14
6.2	Vacuum erection devices	14
6.3	Injectations	14
6.4	Vascular surgery	15
6.5	Surgical implants	15
7	Description of the AMS 700 Penile Prosthesis product line	16
8	What to expect during and after implant surgery	18
8.1	The surgical method	18
8.2	After your surgery	19
8.3	Problems that may develop	20
9	Operating instructions	24
9.1	Inflation	24
9.2	Deflation	25
10	Troubleshooting	26
11	Summary	29
12	Index	30

1 General warnings and cautions

1.1 How safe are silicone elastomer prostheses?

The AMS 700™ with MS Pump™ Penile Prostheses product line devices are made of solid silicone elastomers (a type of rubber). Penile prostheses, or solid silicone elastomers, do not contain silicone gels. Silicone elastomers have been commonly used in many different types of biomedical devices for over 40 years.

Solid silicone elastomers are also used to compare against when a new material is being considered for use in a biomedical device. The new material is tested to see if it is as biocompatible (causes as few problems in living tissue) as silicone elastomers.

Scientific literature has included reports of adverse events in some patients with implantable silicone devices. These adverse events indicate allergic-like reactions or autoimmune-like symptoms. (In an autoimmune reaction, the body's own immune cells may attack some or many of the body's own tissues by mistake.) However, even though these reactions or symptoms were seen in some patients, there has been no proof that the silicone elastomer caused them.

Silicone elastomer may sometimes lose tiny particles off its surface after it has been implanted. Sometimes these particles migrate (move) to lymph nodes in other parts of the body where the particles then stay. (Your lymph nodes are a normal part of your body's defense system against infection.) Medical journals, however, have indicated that particle migration has not resulted in any adverse effects to a patient's health.^{1,2}

Fluorosilicone (a silicone fluid) is also used as a lubricant between internal silicone elastomer and fabric layers in the cylinder to reduce wear. Silicone fluids have a large history of use in medical devices, such as lubricating hypodermic syringes.

1.2 Will my prosthesis have to be replaced?

It is not possible to predict how long an implanted penile prosthesis will function in a particular patient. As with any biomedical device, penile prostheses are subject to wear and eventual failure over time, and therefore should not be considered to be lifetime implants. Other causes can lead to malfunctions of the prosthesis which occur more quickly. Discuss any changes you notice in the function of the prosthesis with your doctor.

Product wear or other mechanical problems, such as unintended inflation or deflation, or difficulty or inability to inflate or deflate may lead to additional surgery to remove or replace the device.

Because it is not possible to tell how long an implanted penile prosthesis will function in a particular patient, AMS has collected information on three sets of data on device removals and revisions to help you understand product performance over time.

The first data set comes from a clinical trial conducted for this penile prosthesis. This data shows that approximately 84% of patients receiving an implant did not require additional surgery (to remove or replace the prosthesis) for 5 years after the original surgical implant.

The second set of data comes from Patient Information Forms (PIFs) submitted to AMS by doctors after surgical procedures for penile prostheses requiring replacement of parts. This data indicates that approximately 89% of patients receiving an implant did not require any additional surgery (to remove or replace the prosthesis) for 5 years after the surgical implant.

¹ Barrett DM, O'Sullivan DC, Malizia AA, Reiman HM and Abell-Aleff PC. "Particle Shedding and Migration from Silicone Genitourinary Prosthetic Devices," *J. Urol.* 146: 319-322. (1991)

² Reinberg Y, *J. Urol.* 750: 694-696. (1993)

The third set of data comes directly from physician implant records. These records contain information on each surgery for any reason performed by a doctor after an original implant. This data set indicates that approximately 82% of patients did not require any additional surgery (to remove or replace the prosthesis for any reason) for 5 years after the original implant.

1.3 Need for manual dexterity

The AMS 700 with MS Pump prostheses require some manual dexterity to inflate and deflate.

1.4 Possibility of malfunction

The possibility of leakage, blockage, or device malfunction exists.

1.5 Possibility of changes in the penis or scrotum

Implantation of a penile prosthesis may cause the penis to become shorter, curve, or be scarred. Scrotal deformity (pump bulge in the scrotum) may also occur.

1.6 Destruction of any latent natural erection

Implanting a penile prosthesis is likely to damage or destroy any remaining natural ability to have an erection.

1.7 Different erection

Your erection with the prosthesis may be different from your original, natural erection. Differences may include a shorter penis, less firmness, less width, and reduced sensation in the penis. Also, because the prosthesis will not extend to the tip of your penis (the glans), this part of your penis may be floppy.

1.8 Less softness when deflated

When your prosthesis is deflated, your penis may not be as soft (flaccid) as it was naturally.

1.9 Possibility of infection

Contact your doctor immediately if there is redness, swelling, and/or heat around the incision area or drainage from the incision. This may indicate an infection.

1.10 Delay intercourse until your doctor gives the OK

If you attempt intercourse before your incision has healed completely, you risk the possibility of infection, pain, or surgical complication.

1.11 Erosion

Contact your doctor immediately if there is thinning of skin or tissue over the prosthesis. This may indicate erosion. Failure to treat erosion can make it worse and lead to infection and loss of tissue.

1.12 Pain

Contact your doctor if you have pain that is very severe or if it lasts longer than expected. Such pain may be a symptom of a medical complication or mechanical device malfunction.

1.13 Migration

Contact your doctor if the surface of any part of your device is visible through your skin or if you cannot locate the MS pump in your scrotum. These symptoms may indicate that a part of your device may have moved within your body or may be moving to the outside of your body.

1.14 Allergic reaction

If you have a AMS 700 with MS Pump with InhibiZone™ and exhibit symptoms of an allergic reaction, such as a rash, hives or difficulty breathing, contact your doctor immediately. Because of this slight risk of allergic reaction, you should be sure and inform your doctor about any allergies you may have to antibiotics prior to your surgery.

2 What is impotence?

Impotence (also known as male erectile dysfunction) is the inability to maintain an erection that is firm enough or that lasts long enough to have successful intercourse.

Impotence:


- is frustrating
- has physical or psychological causes
- is common

Impotence:

- is not premature ejaculation
- is not low sex drive
- is not the need for manual or oral stimulation to get an erection
- is not caused by too much sex
- is not caused by masturbation when you were younger
- is not connected to infertility
- is not “normal” at any age


2.1 How does an erection occur?

For an erection to occur, several parts of the body must work together. The brain sends messages to control the nerves, hormone levels, blood flow and muscles that cause an erection. If anything interferes with one or more of these messages, or if any part of the system doesn't function correctly, an erection will not occur.


The flaccid penis

The shaft of the penis contains the corpora cavernosa, two channels which run the full length of the penis and into the pelvis. These channels are rich in special blood vessels which contain relatively little blood when the penis is flaccid (soft, relaxed, not enlarged). Without sexual stimulation, the penis remains flaccid.


The tumescent penis

When the brain is sexually aroused, it stimulates the nervous system to enlarge the blood vessels in the corpora cavernosa, making room for extra blood. The corpora cavernosa absorb arterial blood flow like a sponge. At this point the penis is swollen, but not yet rigid (hard) enough for intercourse.

The erect penis

As the corpora cavernosa continue to absorb blood, the swelling puts pressure on the veins in the penis. This traps the blood in the corpora cavernosa, making the penis rigid (hard).

2.2 Impotence is common

Most men occasionally have difficulty getting an erection. However, about one in ten men are afflicted by chronic (continuing) impotence.

Impotence can be devastating to those affected by it. Unfortunately, it can also be very difficult to talk about. But the first step is to seek professional help. If your doctor regularly treats impotence problems, you can be sure that any question you ask will be one he has heard before from other patients.

2.3 Impotence is treatable

In more than half of all impotence cases the cause is physical-the result of diabetes, blocked arteries, a hormone problem, or other causes. In other cases, the cause could be psychological-the result of stress or depression. And in many cases, physical causes can produce psychological side effects.

Diagnosing the cause is the first step before recommending a treatment. Your doctor will require a thorough history, physical examination, and laboratory tests to determine whether the impotence has a physical or psychological cause.

The good news is that there are treatments available. And there is usually more than one option. Choices may range from sex counseling or marriage counseling, to medical and surgical treatments. Your doctor can determine which treatment alternative is appropriate for you. Most impotence problems can be treated successfully.

3 Causes of impotence

3.1 Psychological causes

Depression can cause a lack of energy and a reduced sex drive. This may result in an occasional inability to get an erection. If this happens, you may become even more depressed. This can lead to impotence. The first step is to treat the depression.

Stress can also result in an occasional inability to get an erection. Stress can be caused by your job, marital, financial, or other situation. Like depression, this inability to have an erection adds to the stress and can lead to impotence.

Performance anxiety (fear that you will fail if you try to have intercourse) happens to most men once in a while. If it keeps happening it can lead to the inability to have an erection and, ultimately, impotence.

Misinformation about sexuality and about how men should or shouldn't be able to "perform" at a certain age can lead to anxiety and stress, which can lead to impotence.

To treat a psychological cause, your doctor may recommend that you seek treatment from a qualified psychologist, psychiatrist, sex therapist, or marriage counselor. Counseling can often resolve the psychological problem causing impotence or be part of the recommended treatment. You and your partner may wish to go through counseling together.

Even if the problem is physical, there may be psychological side effects. Therefore, counseling may also be part of the recommended treatment for a physical problem.

3.2 Physical causes

Diabetes can cause damage to the nerves or blood vessels that control the flow of blood to the penis. In some cases, keeping your diet and blood sugar under control can help. But permanent damage to these nerves and vessels may result in chronic (continuing) impotence.

Cardiovascular problems, such as hardening of the arteries, can decrease the blood flowing into the penis. This makes it difficult for you to get or keep an erection. In other cases, the veins that keep the blood in the penis during an erection are damaged. If this happens, you cannot keep the erection long enough for sexual intercourse. Impotence can also occur if the nerves that control this flow of blood to the penis are damaged by lack of blood.

Trauma (injury) or pelvic surgery, including cancer surgery in the prostate, bladder, colon or rectal area, can cause impotence. In cancer surgery, the surgeon's most important goal is to remove all of the cancer. Nerves and blood vessels that control erections may be near the cancerous tissue. Sometimes these are damaged in an effort to remove the cancer.

Neurological disorders such as spinal cord injuries can cause impotence. The spinal cord is the relay center for nerve impulses, brain messages and blood flow. When the spinal cord is damaged in certain locations, messages can't get through to the nerves of the penis, causing impotence.

Medications may cause impotence by interfering with the nerve impulses or blood flow to the penis. These medications include some prescriptions for high blood pressure, depression and a number of other conditions. Sometimes a change in the medication or the dosage will decrease the risk of impotence.

CAUTION: Medications should never be changed without the doctor's permission.

Alcoholism changes hormone levels and can lead to permanent nerve damage, causing impotence. This type of impotence may be reversible, depending on the severity of the nerve damage.

Hormone problems are rarely the cause of impotence, but certain diseases can change the balance of hormones which control erections. Kidney failure and liver disease are among these conditions.

4 Am I impotent?

If, after reading this, you are still questioning whether or not you are impotent and what the cause may be, the following questions may indicate whether you should consult a urologist, a doctor who specializes in treating problems associated with the urinary tract. In general, if you answer “yes” to any of the first six questions, you should see a urologist who specializes in the treatment of impotence.

If your current doctor does not regularly treat impotence problems, you might ask for a referral to a urologist who specializes in the diagnosis and treatment of impotence.

When you see the urologist, share the answers from the following self-test. This information will be valuable to the doctor as he or she determines a diagnosis.

Most importantly, you should know that continuing impotence is not a problem that you should have to live with, or that will just go away. It's important to make a medical appointment as soon as the problem becomes apparent. It is important so that you can help avoid the psychological problems that may make it worse. Remember that impotence can be treated in almost every case.

4.1 Impotence self-test¹

Current sexual performance

1. Have you had any difficulty recently in achieving erections?
2. Does this problem occur at least three out of every four times that you attempt intercourse?

Sexual performance trends

3. Have you been having difficulty for longer than one month in achieving erections regularly?
4. Are morning and spontaneous erections becoming less common?
5. Does it take much longer to achieve an erection than in the past?
6. Has it become more difficult to have intercourse in certain sexual positions?

Medical history

7. Have you ever been told you have any form of heart disease, especially hardening of the arteries, peripheral arterial disease (PAD), or hypertension?
8. Have you ever had an operation for heart disease or some other cardiovascular problem?
9. Have you ever been told you have an elevated cholesterol level?
10. Do you ever experience serious pain in the legs when walking?
11. Are you taking drugs on a prescription basis for any other problem?
12. Do you have any known glandular disorder, especially diabetes?
13. Do you have any neurological disorder, such as multiple sclerosis or epilepsy?
14. Have you ever had major surgery in the pelvic area, especially surgery involving the prostate gland or colon?
15. Have you ever had an injury involving the pelvic area, back, spinal cord, or head?
16. Have you ever been treated with radiation therapy for a problem in the pelvic area?

¹ Self-test questions are reprinted and edited with permission from Love Again, Live Again by Steven Morganstern, MD, and Allen Abrahams, PhD.

17. Have you ever had an episode of priapism (an erection which lasts for longer than usual and becomes painful)?

Lifestyle

18. Do you now smoke or did you once smoke for a long period of time?

19. Are you a heavy drinker or a diagnosed alcoholic?

20. Have you used illegal drugs, especially cocaine?

21. Are you a frequent user of drugs you can buy at the drugstore without a doctor's prescription?

22. Are you excessively overweight?

If it is determined that you are impotent, and that your impotence is caused by a physical problem, your doctor can offer several options for treatment. These may include hormone therapy, injections, vacuum devices, penile implants, or vascular surgery.

CAUTION: Medications should never be changed without the doctor's permission.

5 What to expect at your examination


5.1 Physical examination

Your doctor will ask you several questions in order to understand when and under what circumstances you experience signs of impotence. Then your doctor will give you a complete physical exam. This exam is to determine if the blood vessels, nerves, and tissues of your penis are working normally.

Your doctor may begin by feeling for the pulse in your penis and surrounding pelvic area. This will give an indication if the blood supply to your penis is adequate.

Your doctor must also perform a rectal examination to check for prostatitis (a swollen prostate gland). Problems with your prostate can decrease blood flow and feeling in the penis. Prostatitis can also make intercourse uncomfortable.

Your doctor will also check for physical abnormalities such as Peyronie's disease (a curved and painful erection caused by scar tissue within the penis). He will also check your history for previous injuries or surgery in the pelvic area which may have caused nerve damage.

5.2 Tests you may be asked to take

To confirm your diagnosis, other tests can detect hormonal abnormalities, determine blood flow problems, and may help to rule out psychological problems.

Blood tests and urine analysis


These tests are used to measure your hormone levels, cholesterol, and triglycerides (to detect hardening of the arteries), and liver and kidney function. To detect diabetes mellitus, a blood glucose test may also be requested.

Penile blood flow studies

Additional tests may be done to see how effectively blood flows into the penis. One test involves an injection with a drug that increases the blood flow to the penis directly, without stimulating the nerves in the penis. If the blood vessels of your penis are healthy, this injection will produce an erection.

Sleep monitoring

Most men experience at least 3 to 4 erections each night when they are dreaming. If you don't have nighttime erections, perhaps the nerve or blood supply to your penis is inadequate for erections. Your doctor may ask you to measure nighttime erections at home with a simple test that he or she will provide.

6 Treatment options for physical causes of impotence

Depending on your diagnosis, your recommended treatment may be medical or surgical. Medical treatments range from simply changing your prescription drugs, to hormone replacement therapy, antidepressant therapy, and devices or self-injection therapy to produce erections. Surgical treatments include vascular surgery or implants. Your doctor will discuss options that may be appropriate in treating your impotence, as well as the risks and benefits of each option.


6.1 Medication

Changing prescription medications or their dosages may change the side effects, which may be causing your impotence. Hormone replacement therapy may be recommended if you have a hormone deficiency. Antidepressant drugs may be the first course of treatment if you've been diagnosed with severe clinical depression. Drugs that increase the flow of blood to the penis to help cause an erection are also available.

6.2 Vacuum erection devices

These devices are placed around the outside of the penis and draw blood into the penis by creating a vacuum. Then you or your partner place a constriction band (rubber band) around the base of the penis to keep the blood inside the penis until sexual intercourse is completed.

6.3 Injections

Injecting medication directly into the penis prior to intercourse can also produce an erection. If you and your doctor choose this option, you will be taught how to administer the injections yourself.

6.4 Vascular surgery

For a few men, vascular surgery may be indicated to improve blood flow into the penis. Leaking veins may also be surgically repaired. In cases where arterial blockage is reducing blood flow to the penis, an arterial bypass around the blockage may be recommended.

6.5 Surgical implants

Penile implants, or prostheses, may be a choice that lasts longer for a number of impotent men, especially those who've tried psychological and other medical treatments without success. Implants have helped over 300,000 return to an active sex life. Many studies show most patients and their partners are highly satisfied with the results.

Implants are concealed entirely within the body. They require manipulation by you or your partner before intercourse to make the penis firm enough for sexual intercourse. Manipulation is also needed afterwards to return the implant to a relaxed state (make it flaccid).

There are several types of implants to choose from. Differences include manner of operation, naturalness of the erection, and the number of components implanted. In choosing a penile prosthesis you should consider the manual dexterity that is needed to operate each type of device. The best choice of penile prosthesis for you will depend upon your medical condition, your lifestyle and, possibly, the cost of each prosthesis.

Before deciding on a penile prosthesis, you should consult with your doctor about the physical, psychological, cosmetic, and functional outcome of the implantation surgery. You should be sure that you understand the risks and benefits of the surgery.

There are certain instances when you and your physician may decide that a surgical implant is not an appropriate choice for you if:

- ▶ the risks associated with surgery are too high because of your medical condition
- ▶ your medical history includes previous surgery which make an implant impossible
- ▶ you are satisfied with one of the less invasive treatment options discussed in section 6
- ▶ you choose not to be implanted with a silicone elastomer device
- ▶ you want to preserve the internal physical structure of your penis in case someday you may again experience a normal erection

In addition to discussions with your doctor, you may want to discuss the various options with your partner.


7 Description of the AMS Penile Prosthesis product line

The AMS 700 with MS Pump prostheses consist of four parts which are connected by tubing: a reservoir, two cylinders, and the MS pump.

The **reservoir** is implanted in the lower abdomen, under the muscle layer, and is filled with sterile saline solution.

The two **cylinders** are inserted side by side into the corpora cavernosa of the penis.

The **pump** is placed in the scrotum. The lower, rounded part of the pump is the **pump bulb**. The upper part of the pump that is rectangular in shape is the **deflation block**. The parts of the pump are marked in the figure below. Squeezing the deflation button starts to deflate the cylinders. Squeeze the deflation button for about 4 seconds to completely deflate the cylinders.

There are three models of inflatable penile prostheses in the AMS 700 product line. Your doctor will choose the right model for you based upon your anatomy and medical history.

AMS 700 LGX™ with MS Pump™ Penile Prosthesis

The 700 LGX cylinders are designed to expand in both length and width (girth). The amount of cylinder lengthening which actually occurs with use depends on how much the tissue in your penis allows them to lengthen and whether you are able to completely inflate the cylinders. This means that even if the cylinders inside your penis get longer when you fill them with fluid, they cannot make your penis any longer than it is now.

CAUTION: Implantation of a penile prosthesis may cause the penis to become shorter, curve, or become scarred.

AMS 700 CX™ with MS Pump™ Penile Prosthesis

The cylinders expand in width (girth) only and, like all AMS cylinders, they have been designed to resist uneven inflation.

AMS 700 CXR™ with MS Pump™ Penile Prosthesis

The 700 CXR cylinders have the same feature as the 700 CX. However, these cylinders are more appropriate for patients whose anatomy requires shorter, narrower cylinders.

8 What to expect during and after implant surgery

Implantation of a penile prosthesis involves a surgical procedure usually lasting from 30 minutes to 2 hours. The length of your hospital stay depends on your physical condition and the type of prosthesis chosen. You will be able to return to work and everyday activities at the discretion or direction of your physician. (Everyday activities include any type of activity you were able to perform before your implant surgery, for example, exercise, working, bathing.)

8.1 The surgical method

Your doctor should be able to give you a thorough explanation of what will happen during the surgery and the rest of your hospital stay. In general, the procedure begins with some preoperative tests, which may include blood tests, urine analysis, and delivery of antibiotics.

Depending on your physical condition and your doctor's preference, you will either be given a local anesthetic (it numbs only the area where the surgery occurs) or a general anesthesia (which puts you to sleep for the surgery). If you are to be given a general anesthetic, you will be asked to abstain from food or drink for 12 hours before surgery.

There are two kinds of incisions. Your doctor will choose the right one for you:

- ▶ **penoscrotal** (between your penis and scrotum)
- ▶ **infrapubic** (in the lower abdomen above the penis)

Penoscrotal

If your doctor chooses the penoscrotal approach, an incision will be made through the skin between your penis and scrotum. Your doctor will usually implant the cylinders, the MS pump and the reservoir through this one incision. Occasionally, the doctor will also make a small incision in the lower abdominal area to place the reservoir.

To implant the cylinders, your doctor will first dilate (widen) the corpora cavernosa (the two channels in the shaft of your penis which fill with blood when you get a natural erection). Then your doctor will measure this area to choose which size cylinder will best fit your anatomy.

After the cylinders are in place, your doctor will make a space in your scrotum in which to place the pump. He will place the pump so that it is easy for you to reach.

Next, your doctor will make a space for the reservoir in your abdomen and connect all the tubing.

Finally, before closing the incision in your skin, your doctor will inflate and deflate the prosthesis to make sure it is working properly.

Infrapubic

An incision will be made in the skin of your lower abdomen. Your doctor will implant the cylinders, pump, and reservoir through this one incision. The same steps used in the penoscrotal approach described above are used here.

Penoscrotal and infrapubic

After surgery is completed, your doctor will usually deflate the prosthesis and tape your penis against your abdomen. Your doctor may insert a catheter (a hollow tube) into your penis to drain your bladder. The catheter will be removed before you leave the hospital.

Your doctor may also insert a tube into your abdomen to drain excess fluid (blood) from the incision site. This will also be removed usually within 12 to 24 hours after surgery (or when drainage stops).

8.2 After your surgery

After you leave the hospital, you will be encouraged to avoid wearing tight-fitting underwear. This helps to prevent any curvature in your penis that could occur during healing. Likewise, after voiding you should carefully retape your penis to the abdomen in the same position (straight) that it was in prior to untaping.

Many doctors recommend that you wait four to six weeks before having intercourse. This time allows your incision site to heal.

CAUTION: If you attempt intercourse before the incision has healed completely, you risk the possibility of infection, pain, or surgical complications.

You will probably have an appointment with your doctor during this time to be sure you are healing properly. Be sure to discuss these possibilities with your doctor and ask how long after your surgery you should wait before having intercourse.

You may experience pain at the operative site(s) during the early time after surgery and when you first use your prosthesis. In most cases the pain goes away within a few weeks of surgery; however, cases of chronic (continuing) pain have been reported.

Recovery times vary from patient to patient. You will be able to return to work and everyday activities at the discretion or direction of your physician. (Everyday activities include any type of activity you were able to perform preoperatively, for example, exercise, working, bathing.) Your doctor will also talk to you about when you will be able to use your device.

After your doctor says you can begin using the device, follow the operating instructions given in Section 9.

CAUTION: Use of injection therapy into the penis (i.e., injection of medications for impotence) can damage the penile prosthesis. Do not use injection therapy after you receive your implant.

You will also have several postoperative visits, and annual or semi-annual follow-up visits with your doctor after the surgery. During this recovery time and after, take care to avoid trauma to the pelvic or abdominal area. Always keep in mind that you have had a surgical implant and choose your activities wisely. Examples of such trauma may include a seat belt jolt from a car accident, being tackled in contact sports, or slipping and falling on ice. Trauma may damage the prosthesis or surrounding tissues.

CAUTION: Take precautions to avoid trauma to your pelvic and abdominal areas such as might occur during contact sports, slipping on ice, etc.

8.3 Problems that may develop

Implant surgery carries the same types of risks that every surgical procedure involves, including that of infection and those associated with anesthesia. In addition, the outcome of your implant surgery may be unsuccessful. For example, the device may fail to function as intended. If this happens, you may need additional surgery to remove or replace the prosthesis. If the prosthesis must be removed, reimplantation of a new prosthesis may be complicated by the amount of time between the two surgeries. Discuss these possibilities with your doctor.

Infection

Infection can happen after any kind of surgery. Your doctor will try to lower your risk by giving you antibiotics before and after your operation and by flushing (washing out) the surgical site with antibiotics during surgery. Some conditions increase the risk of getting an infection:

- ▶ diabetes
- ▶ a spinal cord injury
- ▶ open sores
- ▶ an existing skin infection near the incision site
- ▶ an existing urinary tract infection

Information about the effectiveness of the antibiotic treatment, referred to as InhibiZone™ (IZ), in preventing infections comes from Patient Information Forms (PIFs), which are records submitted to AMS by doctors after surgical procedures to implant penile prostheses or replace parts of previously implanted devices. More than 43,000 surgeries have been recorded in the AMS PIF database during the period since IZ was introduced over 6 years ago. During this period, approximately 1.2% of patients implanted with an IZ-treated prosthesis required a second surgery to replace or remove the prosthesis because of infection. This rate is compared to patients implanted with prostheses that were not treated with IZ (non-IZ), where approximately 2.5% of patients required a second surgery due to infection.

If you get an infection that cannot be treated successfully with antibiotics, your doctor may have to remove the prosthesis. It may not be possible to implant a new one.

Also, if your prosthesis must be removed because of infection, that infection may leave scars inside your penis. This may also make implanting a new prosthesis difficult.

WARNING: Contact your doctor immediately if you notice any redness, swelling, and/or heat around the incision area or drainage from the incision. These symptoms may indicate an infection.

Allergic reaction

Allergic reactions may occur if you have sensitivity to the antibiotics that were used during the penile prosthesis implantation procedure or, if you are implanted with the AMS 700 series with MS pump Penile Prosthesis with InhibiZone™, to the antibiotics impregnated into the penile prosthesis.

WARNING: Contact your doctor immediately if you notice signs of an allergic reaction such as rash, hives or difficulty breathing.

Erosion

Erosion is when the tissue next to any part of the device is “worn away”. Conditions that can cause erosion include:

- ▶ infection
- ▶ pressure on the tissue, cutting off the blood supply
- ▶ improper sizing
- ▶ tissue damage
- ▶ misplacement of the cylinders, reservoir, or MS pump

Erosion involving the cylinders most often involves:

- ▶ the glans (the tip of the penis)
- ▶ the urethra (the tube inside the penis that carries urine out of the body)
- ▶ the skin of the penis

The pump may erode through the skin of the scrotum.

The reservoir can erode into the bladder or bowel.

Symptoms of erosion into the scrotum or penis may include, after having been symptom free, pain, redness of skin, tenderness over the involved part, changes in skin texture, drainage, and/or being able to see the prosthesis through the skin. Erosion into the bladder may result in pain, tenderness in bladder area, a change in your ability to urinate or a change in the color of urine. Less than 1.5% of patients in the clinical trial experienced erosion of any part of the device.

WARNING: Contact your doctor immediately if you notice any pain, tenderness over the involved part, change in skin texture, drainage, or if you can see the prosthesis through your skin. These symptoms may indicate erosion. Failure to treat the erosion can make it worse and lead to infection and loss of tissue.

Your doctor must evaluate any possible erosion. Sometimes the tissue can be repaired and only part of the prosthesis replaced. Other times the entire device must be removed.

Trauma

Trauma (injury) to the hip or stomach area can cause damage to either the device or the surrounding tissue in your penis, scrotum or abdomen. This can cause the device to malfunction and could require surgery to replace it. Some things you can do to decrease possible damage are:

- ▶ avoid contact sports where you might be tackled
- ▶ take extra precautions when walking on ice to prevent slipping and falling

Pain

It is normal to have some pain in your penis and scrotum immediately after surgery and when you are first using the device. Approximately half of the patients in the 5-year clinical trial reported pain to their doctors after surgery.

WARNING: Contact your doctor if you have pain that is very severe or if it lasts longer than expected. Such pain may be a symptom of a medical condition or mechanical device malfunction.

Some patients have had chronic (continuing) pain with no known medical cause. These cases are not common but are unfortunate when they do happen. Sometimes these patients have chosen to have the device removed because the pain would not go away.

Migration

Migration is the movement of one or both cylinders, the pump, or reservoir within the body space where they were originally placed.

If migration occurs, it can cause pain, psychological/medical complications, or device malfunction. Migration may need to be corrected with surgery.

Causes of migration include:

- ▶ improperly-sized cylinders
- ▶ improper positioning of the MS pump or reservoir
- ▶ incorrect tubing length

Only 1% of the patients in the clinical trial had device parts which migrated.

WARNING: Contact your doctor if any part of the device is visible through your skin or if you cannot locate the MS pump in your scrotum.

Mechanical problems

Product wear (the use of the device over a period of time) or other mechanical problems may occur over time. Surgery may be required to correct the problem.

These problems may include unintentional (spontaneous, with no manipulation on your part) inflation or deflation of the device. Also included are difficulty or inability to inflate or deflate the device.

If such problems occur, first check the user instructions to be sure you are operating the device correctly. If you still have the same problem, contact your doctor.

About 10% of patients in the 5-year clinical trial experienced mechanical problems. The AMS 700 series with MS Pump has a new lock-out valve and Parylene coating design to reduce mechanical problems.

Edema

Edema is when the tissue next to any part of the device is swollen. It is normal for some swelling to occur in the healing period after surgery. About 1/3 of the patients in the clinical trial reported having swelling.

Bruising or discoloration of the penis

It is normal to have some bruising after surgery. About 10% of the patients in the clinical trial reported bruising to their doctors. In these patients, the bruising normally went away in about 13 days.

Obvious redness in penis or scrotum

It is normal to have some redness of the skin around the implant. However, extreme redness should be reported to your doctor. In the clinical trial, less than 5% of the patients reported obvious redness to their doctors. The redness went away in an average of 43 days.

Other

Problems which occurred in less than 4% of patients in the clinical trial were: problems with passing urine, painful joints, decreased penile sensation, hematoma, ejaculation problems, curving of the penis, problems with the surgical incision, numbness, tingling, pain in canal in penis which urine flows through, and pain or burning when urinating.

The following "Other" adverse device effects each occurred in 1.0% or fewer of the patients in the 5-year trial: inability to urinate, weakness, device parts "sticking" to the tissue inside the penis or scrotum, fever, dizziness, incorrect positioning of the device, depression (feeling of sadness or hopelessness), cellulitis (swelling and infection far below skin surface), fibrosis (the formation of fiber-like tissue around parts of the device), hematuria (blood in the urine), necrosis (death of body cells, tissue or skin), abnormal sexual function, urinary frequency (having to urinate a lot), infection in the bladder or canal which urine flows through), changes in the tip of the penis, thickening of the skin, stool (bowel movement) leakage, rheumatoid arthritis (swelling of muscles or tendons around joints), dry mouth, memory problems, baldness, diabetes, eye problems, headaches, vertigo, pelvic and back pain, kidney problems, urinary urgency (feeling of severe need to urinate often), reactions to sunlight and epigastric pain ("heartburn" related pain).

Please ask your doctor for explanation on any of the problems that you do not understand.


Note: All complication rates provided in this section are based on a prior version of the AMS inflatable penile prosthesis. Based on design similarities between the AMS 700 with MS Pump and this prior version, it is expected that the complication rates will be comparable.

9 Operating instructions


9.1 Inflation

When you squeeze and release the pump bulb, the fluid moves from the reservoir, through the pump and into the cylinders. When the cylinders fill with fluid, the penis becomes erect.

1. Feel for the pump in your scrotum.
2. Grasp the tubing above the deflation block with one hand to hold the pump in place.
3. Use the thumb and forefinger of your other hand to locate the soft, rounded pump bulb.
4. Initially, squeeze the pump bulb **quickly** and **firmly**, and then release it to activate the pump. The remaining pump bulb squeezes can be slower. Releasing the pump bulb allows it to refill with fluid. Continue to alternate squeezing and releasing until the cylinders fill with fluid and become erect. Full erection may take as few as 10, or as many as 25 squeezes of the pump bulb.

Note: Do not squeeze the deflation button while squeezing the pump bulb.




9.2 Deflation

When you squeeze the deflation button for about four seconds the fluid begins to leave the cylinders and return to the reservoir. After about four seconds release the deflation button and the cylinders will return to their deflated position. The penis becomes soft.

1. Feel for the pump in your scrotum.
2. Grasp the tubing above the deflation block with one hand to hold the pump in place.
3. With your other hand, locate the raised deflation button on the deflation block.
4. To effectively squeeze the deflation button, you must place a thumb and forefinger on opposite sides of the deflation block.
5. Squeeze the deflation button for about four seconds and then release. The cylinders will continue to deflate and the penis will become soft (flaccid).



NOTE: Do not squeeze the deflation button while squeezing the pump bulb.

6. After cylinders have deflated, you may squeeze your penis to make it more flaccid.


10 Troubleshooting

Signs and symptoms that may develop after surgery

Symptom	Problem	What to do*
Inability to inflate.	Deflation button squeezed while squeezing the pump bulb.	<ol style="list-style-type: none"> 1. Grasp the tubing above the deflation block with one hand to hold the pump in place. 2. Do not squeeze the deflation button while squeezing the pump bulb. 3. Inflate normally.
	Pump bulb does not refill and is dimpled or collapsed.	<ol style="list-style-type: none"> 1. Squeeze the deflation button. The pump bulb should refill. 2. First squeeze of the pump bulb should be quick and firm to activate the pump, then release. 3. Do not squeeze the deflation button while squeezing the pump bulb. 4. Continue to squeeze and release normally to inflate cylinders.
 <p>NOTE: Do not squeeze the deflation button while squeezing the pump bulb.</p>	First pump bulb squeeze was not firm enough or long enough to activate pump bulb.	<ol style="list-style-type: none"> 1. Firmly squeeze the deflation button. 2. First squeeze of the pump bulb should be quick and firm to activate the pump, then release. 3. Do not squeeze the deflation button while squeezing the pump bulb. 4. Continue to squeeze and release normally to inflate cylinders.
	Possible mechanical problem.	Contact your doctor.

* For further assistance in operating your device call the AMS patient liaison at +1 952-930-6261, or 1-800-328-3881, extension 6261, or visit our Web site at www.AmericanMedicalSystems.com.

Signs and symptoms that may develop after surgery (continued)

Symptom	Problem	What to do*
Inability to deflate.	Difficulty locating the deflation button.	<ol style="list-style-type: none"> 1. Grasp the tubing above the deflation block with one hand to hold the pump in place. 2. With your other hand, locate the raised deflation button on the deflation block.
	Cylinders stay inflated after deflation button is squeezed and released.	Squeeze the deflation button for a longer time (at least 5 seconds).
	Deflation button and pump bulb were squeezed at the same time.	<ol style="list-style-type: none"> 1. Squeeze the sides of deflation block for about 4 seconds, then release.
		<p>Squeezing the sides of the deflation block</p>  <ol style="list-style-type: none"> 2. Next, squeeze the deflation button for at least 5 seconds, then release. 3. Cylinders should deflate normally.
	Possible mechanical problem.	Contact your doctor.
Pain. Discharge from incision. Redness, swelling.	Infection.	Contact your doctor.
Cylinder surface can be seen through the skin.	Erosion of cylinders (associated with infection).	Contact your doctor.
Stomach Pain.	Erosion of pump (associated with infection).	Contact your doctor.
Pain.	Post-operative surgical pain.	Contact your doctor.
	Pain with intercourse.	Contact your doctor.
Insufficient flaccidity.	Cylinders stay partially inflated after deflation button is pressed and released.	Squeeze the deflation button for a longer period of time.
	Fluid remaining in cylinders.	Squeeze penis to make more flaccid.

* For further assistance in operating your device call the AMS patient liaison at +1 952-930-6261, or 1-800-328-3881, extension 6261, or visit our Web site at www.AmericanMedicalSystems.com.

Symptom	Problem	What to do*
Unintentional inflation.	Fluid has moved into cylinders.	Follow deflation instructions on page 25.
	Possible mechanical problem.	Contact your doctor.
Unintentional deflation or loss of rigidity.	Cylinders have not been fully inflated.	Follow inflation instructions on page 24.
	Deflation button was unintentionally squeezed.	Follow inflation instructions on page 24.
	Possible mechanical problem.	Contact your doctor.
Penile curving or bulging.	Possible mechanical problem.	Contact your doctor.
Inability to locate pump.	Pump has moved.	Contact your doctor.
Any part of the device is visible through your skin.	Possible mechanical problem.	Contact your doctor.
Pain, skin disruption (opening), leaking of body fluids, bruising.	Trauma.	Contact your doctor.
Skin rash or hives.	Possible allergic reaction.	Contact your doctor.
Difficulty breathing.	Possible allergic reaction.	Contact your doctor.

* For further assistance in operating your device call the AMS patient liaison at +1 952-930-6261, or 1-800-328-3881, extension 6261, or visit our Web site at www.AmericanMedicalSystems.com.

11 Summary

This summary is not meant to replace the complete instructions found in this manual. The entire manual should be read before operating your device.

Your inflatable penile prosthesis requires some manual dexterity to inflate and deflate.

The possibility of leakage, blockage, or device malfunction exists. Discuss any changes you notice in the function of your prosthesis with your doctor.

Contact your doctor immediately if there is:

- ▶ redness
- ▶ swelling
- ▶ heat around the incision area or drainage from the incision (symptoms of erosion)
- ▶ if your pain is very severe or lasts longer than expected
- ▶ if any part of your device is visible through your skin or
- ▶ if you cannot locate the MS pump in your scrotum

12 Index

Topic	Section number
Dexterity	[Section 1.3]
Erosion.....	[Sections 1.11, 8.3]
Impotence.....	[Sections 2.2 - 4.1]
Infection	[Sections 1.9, 8.3]
Intercourse	[Sections 1.10, 8.2]
Malfunction.....	[Section 1.4]
Migration.....	[Sections 1.13, 8.3]
Stomach pain	[Sections 1.12, 8.2, 8.3]
Penile change.....	[Section 1.5]
Silicone elastomer safety	[Section 1.1]

AMSTM

American Medical Systems, Inc.

10700 Bren Road West

Minnetonka, MN 55343

U.S.A.

U.S. Toll Free: 1 800 328 3881

Tel: +1 952 930 6000

Tel: +31 20 593 8800

American Medical Systems Europe B.V.

Haarlerbergweg 23 G

1101 CH Amsterdam Zuid-Oost

The Netherlands

© 2016 Boston Scientific Corporation or its affiliates. All rights reserved.

All trademarks are property of the respective owners.

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
The brain controls all sexual functions, from being aroused to starting and controlling the psychological, hormonal, nerve, and blood flow changes that lead to an erection.

Arteries deliver the extra blood to the penis that causes it to stiffen. Veins then drain the blood out of the penis after intercourse.

Nerve impulses relay signals of arousal and sensation to and from the penis.

Hormones, including testosterone, control the male sex drive. Testosterone is secreted by the testicles.

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User's Guide

Table of Contents

The AMS 700 MS Series Penile Prosthesis.....	1
Using the AMS 700 MS Series Penile Prosthesis.....	3
What to expect after the procedure	5
What symptoms may indicate that I should contact my urologist?	6
Where can I get more information about the AMS 700 MS Series Penile Prosthesis?	7

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AMS
Solutions for Life®

AMS 700 MS Series Penile Protheses

The AMS 700 MS Series Penile Protheses is a fluid-filled device that is implanted entirely within the body and is not visible.

The prosthesis consists of three components connected by tubing (Figure 1): a reservoir, two cylinders, and an inflate/deflate pump.

The reservoir is implanted in the lower abdomen, under the muscle layer, and filled with a sterile saline solution. The **reservoir** holds the fluid until the patient decides to inflate the cylinders.

The two **cylinders** are inserted side by side in the shaft of the penis.

The **pump**, consisting of a pump bulb and a deflation button, is placed in the scrotum (Figure 2). The pump is used to inflate and deflate the cylinders by transferring fluid from the reservoir to the cylinders. The pump bulb is the rounded part at the bottom of the pump. The deflation button is above the bulb on the hard, rectangular shape of the pump.


Figure 1


Figure 2

Using the AMS 700 MS Series Penile Prostheses

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1. Feel for the pump in your scrotum.
2. Grasp the tubing above the deflation block with one hand, to hold the pump in place.
3. Use the thumb and forefinger of your other hand to locate the rounded pump bulb.
4. Initially, squeeze the pump bulb **quickly** and **firmly**, and then release it to activate the pump (Figure 3). The remaining pump bulb squeezes can be slower. Releasing allows the inflation bulb to refill with fluid. Continue to alternate squeezing and releasing until the cylinders fill with fluid and become erect. Full erection may take as few as 10, or as many as 25 squeezes of the pump bulb.

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1. Feel for the pump in your scrotum.
2. Grasp the tubing above the deflation block with one hand, to hold the pump in place.
3. With your other hand, locate the raised deflation button on the deflation block.
4. To effectively squeeze the deflation button, you must place a thumb and forefinger on opposite sides of the deflation block (Figure 4).
5. Squeeze the deflation button for about four seconds and then release. The cylinders will continue to deflate and the penis will become soft (flaccid).
Note: Do not squeeze the deflation button while squeezing the pump bulb.
6. After cylinders have deflated, you may squeeze your penis to make it more flaccid.


Figure 4

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You may experience pain or discomfort at the operative site(s) during the early time after surgery and when you first use your prosthesis. In most cases the pain goes away within a few weeks of surgery; however, cases of chronic (continuing) pain have been reported.

Recovery times vary from patient to patient. You will be able to return to work and everyday activities at the discretion or direction of your urologist.

Your doctor will discuss when you will be able to use your device. Most doctors recommend that you wait four to six weeks before having intercourse. This time allows your incision site to heal and help your body adapt to your prosthesis. You will probably have an appointment with your doctor during the recovery period to be sure you are healing properly. At your doctor's discretion, you may have additional visits or follow-up from this procedure.

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After your doctor says you can begin using the device, follow the operating instructions on pages 3 & 4.

During the recovery time and after, take care to avoid trauma to the pelvic or abdominal area. Always keep in mind that you have had a surgical implant and choose your activities wisely. Examples of such trauma may include a seat belt jolt from a car accident, being tackled in contact sports, or slipping and falling on ice. Trauma may damage the prosthesis or surrounding tissues. Use of injection therapy may also damage the prosthesis.

Additionally, your doctor may ask you to take antibiotics before and after any dental or invasive surgical procedure as a precautionary measure to prevent infection.

What symptoms may indicate that I should contact my urologist?

As with any medical procedure, there is the possibility that signs and symptoms may develop after the procedure.

- Unusual or prolonged pain
- Redness or rash near the incision or scrotum
- Swelling near the incision, scrotum or penis
- Pus coming from the incision or head of the penis
- Hives

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- Any problems with urination

You should contact your doctor immediately if any of these situations occur.

Where can I get more information about the AMS 700 MS Series Penile Prostheses?

Your best resource is your doctor. Your doctor is most familiar with your medical situation and can work with you to make sure you get the best possible results from the penile prosthesis.

You can also call the AMS patient liaison at 1-800-328-3881, extension 6261, or visit our Web site at www.AmericanMedicalSystems.com.

Solutions for Life® For more than 35 years, American Medical Systems has provided world-class medical devices used for treating pelvic health issues. Over the past decade, our reputation for quality and medical efficacy has broadened to encompass both devices and therapies that restore pelvic health for men and women. The medical conditions our solutions address include male and female urinary incontinence, erectile dysfunction, prostate disorders (including BPH), urethral strictures, pelvic organ prolapse and fecal incontinence.

AMS

Solutions for Life®

American Medical Systems, Inc. World Headquarters

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Haarlerbergweg 23 G
1101 CH Amsterdam Zuid-Oost
The Netherlands

92127382-15

© 2017 Boston Scientific Corporation or its affiliates. All rights reserved.
All trademarks are the property of the respective owners.

92127382-15 (2017-11)

AMS 800™

Urinary Control System For Female and Pediatric Patients

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AMS 800™

Urinary Control System For Female and Pediatric Patients

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Brief Device Description

The AMS 800 Urinary Control System is an implantable, fluid filled, solid silicone elastomer device used to treat urinary incontinence. It is designed to restore the natural process of urinary control. The device simulates normal sphincter function by opening and closing the urethra at the control of the patient. The AMS 800 consists of three interconnected components: a cuff, a pump, and a pressure-regulating balloon (PRB). The three components are connected with kink-resistant tubing. The AMS 800 can be implanted at the bladder neck.

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Indications for Use

The AMS 800 is used to treat urinary incontinence due to reduced urethral/bladder outlet resistance (intrinsic sphincter deficiency) in males, females and children. The information presented in this document focuses on implantation of the AMS 800 in females and children. For information relating to the implantation of the device in males, see the AMS 800 Urinary Control System Instructions for Use.

Contraindications

1. This device is contraindicated in patients whom the physician determines to be poor candidates for surgical procedures and/or anesthesia due to physical or mental conditions.
2. This device is contraindicated in patients with urinary incontinence due to or complicated by an irreversibly obstructed lower urinary tract.
3. This device is contraindicated in patients with irresolvable detrusor hyperreflexia or bladder instability.
4. The implantation of the InhibiZone version of this device is contraindicated in patients with known allergy or sensitivity to rifampin or to minocycline HCl or other tetracyclines.

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Implantation of products with InhibiZone is contraindicated in patients with systemic lupus erythematosus because minocycline HCl has been reported to aggravate this condition.

Warnings

1. Patients with urinary tract infections, diabetes, spinal cord injuries, open sores, or skin infections in the region of the surgery have an increased risk of infection associated with a prosthesis. Appropriate measures should be taken to reduce the likelihood of infection. Infection that fails to respond to

^{*} InhibiZone™ Antibiotic Surface treatment is not available in all markets.

antibiotic therapy may result in removal of the prosthesis. Infection followed by explantation of the device may result in scarring which may make subsequent reimplantation more difficult.

2. Erosion may be caused by infection, pressure on the tissue, improper cuff sizing, improper balloon selection, tissue damage, and component misplacement. The cuff may erode around the urethra. The control pump may erode through the labia or scrotum. The pressure-regulating balloon may erode into the bladder. Acute urinary tract infection can interfere with proper functioning of the device and may lead to erosion of the urethra in the cuff area. Failure to evaluate and promptly treat the erosion may result in a substantial worsening of the condition leading to infection and/or loss of tissue.
3. Poor bladder compliance or a small fibrotic bladder may require some measure of intervention including, in some cases, augmentation cystoplasty before implanting the prosthesis.
4. Patients with urge incontinence, overflow incontinence, detrusor hyperreflexia or bladder instability should have these conditions treated and controlled (or resolved) prior to implantation of the device.
5. Do not pass a catheter or any other instrument through the urethra without first deflating the cuff and deactivating the device to prevent potential damage to the urethra or the AMS 800.
6. This device contains solid silicone elastomers. This device does not contain silicone gel. The risks and benefits of implanting this device in patients with documented sensitivity to silicone should be carefully considered.
7. Surgical, physical, psychological, or mechanical complications, if they occur, may necessitate revision or removal of the prosthesis. Removal of the device without timely reimplantation of a new device may complicate subsequent reimplantation. The timing of reimplantation should be determined by the treating physician based on the patient's medical condition and history.
8. Product wear, component disconnection or other mechanical problems may lead to surgical intervention. Mechanical complications may include malfunctioning of the components and leakage of fluid. Any mechanical malfunction that does not permit the transfer of fluid from the cuff to the balloon may result in outflow obstruction. Mechanical events should be evaluated carefully by the treating physician and the patient should consider risks and benefits of treatment options, including revision surgery.
9. Previous patient history of adverse reaction(s) to radiopaque solution precludes its use as a filling medium for the prosthesis. Instead, saline should be used to fill the device.
10. Female patients with persistent incontinence should be evaluated to rule out vesicovaginal fistula, which may have been an unrecognized iatrogenic injury.
11. If a hypersensitivity reaction develops to a device coated with InhibiZone, the cuff and pump should be removed and the patient treated appropriately.

Precautions

Patient Related

1. Patient selection requires thorough preoperative consultation and evaluation by the physician.
2. Patients should be counseled in order to have a realistic expectation of the physical, psychological, and functional

outcome of the implantation of an AMS 800. Although the prosthesis is designed to restore urinary control, some patients continue to have a degree of incontinence after this procedure.

3. Patients may experience pain when the device is activated in the postoperative period and during the period of initial use. Cases of chronic pain associated with device have been reported. Pain with a severity or duration beyond what is expected may require medical or surgical intervention. Patients should be counseled on expected postoperative pain including severity and duration.
4. Tissue fibrosis, previous surgery, or previous radiation therapy in the area of the implant may preclude implantation of a cuff at the urethra.
5. Any progressively degenerative disease, e.g. multiple sclerosis, may limit the future usefulness of the implanted prosthesis as a treatment for the patient's urinary incontinence.
6. Adequate manual dexterity, strength, motivation, and mental acuity are required for proper use of the device.
7. Trauma or injury to the pelvic, perineal or abdominal areas, such as impact injuries associated with sports, can result in damage to the implanted device and/or surrounding tissues. This damage may result in the malfunction of the device and may necessitate surgical correction including replacement of the device. The physician should advise patients of these possibilities and warn them to avoid trauma to these areas.
8. Female patients of child-bearing age must be forewarned that pregnancy is acceptable, but cesarean section may be indicated to minimize the risk of damage to the bladder neck and its surrounding cuff. For those patients who become pregnant, deactivation of the device during the third trimester is recommended to reduce the risk of erosion. Those patients contemplating pregnancy should consider delaying implantation.
9. Children receiving device should be evaluated at regular intervals. Lifelong radiological and urodynamic surveillance of the urinary tract is crucial. Before implantation, the patient and their family should be informed of the complication rate and the need for long-term follow-up.
10. Consideration should be given to the diameter of the implanted occlusive cuff relative to catheters or other trans-urethral devices. When fully deflated, the inside diameter of the smallest occlusive cuff (3.5cm) generally exceeds 28F. Additional clearance is required to accommodate the patient's urethral tissue between the trans-urethral device and the occlusive cuff. Urethral tissue thickness is patient specific and requires a physician's assessment to determine its impact on sizing.

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1. Use of products with InhibiZone should be carefully considered in patients with hepatic or renal disease, as use of rifampin and minocycline HCl can cause additional stress on the hepatic and renal systems.
2. Patients who receive a device with InhibiZone and are also taking methoxyflourane should be carefully monitored for signs of renal toxicity.
3. Patients who receive a device with InhibiZone and are also taking warfarin should have their prothrombin time monitored because tetracyclines have been reported to slow

coagulation.

4. Use of products with InhibiZone should be carefully considered in patients using thionamides, isoniazid and halothane, due to potential hepatic side effects that have been reported in patients using these drugs and higher doses of rifampin.
5. Devices with InhibiZone should not come into contact with ethyl alcohol, isopropyl alcohol or other alcohols, acetone or other nonpolar solvents. These solvents may remove the antibiotics from the device.
6. InhibiZone components should not be soaked in saline or other solutions prior to implantation. The components may be briefly rinsed or dipped into a sterile solution, immediately prior to implant, if desired.
7. InhibiZone does not replace your normal antibiotic protocols. Continue using any prophylactic protocols normally used for urological surgical procedures.
8. Because products with InhibiZone are impregnated with a combination of rifampin and minocycline HCl, the contraindications, warnings and precautions regarding the use of these antimicrobial agents apply and should be adhered to for the use of this device, although systemic levels of minocycline HCl and rifampin in patients receiving this device are unlikely to be detected.

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Device Related

1. If the deactivation valve is closed when the cuff is inflated, fluid cannot transfer from the cuff to the balloon and sustained outflow obstruction may arise as a result:
 - a. In the event of large pressures within the bladder, automatic pressure relief that normally occurs with the device would be prevented. Cycling the device can relieve the outflow obstruction.
 - b. Cycling the device may be difficult if deactivation occurs when the pump bulb is deflated. If unable to cycle the device, squeezing the sides adjacent to the deactivation button will allow fluid to fill the pump bulb and then the pump can be cycled normally.
 - c. Release of the deactivation valve may require greater pressure than that used to cycle the device.
2. System pressure changes may occur over time if you fill the balloon with radiopaque solution of incorrect concentration. Follow the instructions in the Operating

Room Manual to prepare the radiopaque solution with the correct concentration.

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For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of the product and packaging in accordance with hospital, administrative and/or local government policy.

Clinical Information – Females

The AMS 800 is indicated for the treatment of urinary incontinence in females and in children. A literature review was performed to evaluate the safety and efficacy of the AMS 800 in female patients. A total of seventeen documents, specifically reporting clinical information for female patients only, were analyzed. Clinical data, including patient demographics, etiology, complications, and success criteria and rates are reported below.

Patients

Clinical efficacy has been reported for 938 female patients, ages ranging from 9-85, implanted with the AMS 800¹⁻¹⁷. The mean reported follow-up on 785 AMS 800 patients was 3.6 years¹⁻¹⁷. Various patient etiologies were reported, including neuropathic bladder dysfunction^{1,2,3, 4,13}, stress incontinence^{1,4,6,7,8,9,12,15,16,17}, pelvic/perineal trauma^{1,4,12,13}, congenital defects (exstrophy)¹, hysterectomy¹, radiotherapy¹, cystectomy¹, nonneurogenic bladder^{2,13}, intrinsic sphincter deficiency^{2,6}, myelomeningocele^{13,14}, spinal injury¹³, and cerebral palsy¹⁴. Prior surgery or radiotherapy were reported in 354 patients^{2,3,5,6,8,10,12,15,16}.

Complications

The following adverse events have been associated with the use of this product: bladder injury, bleeding, bothered by device, contracture, cuff erosion, cutaneous erosion, deep vein thrombosis, device damage – unrelated procedure, device displacement, device fluid leak, device malfunction, dysuria, early post-operative morbidity, edema, enterocutaneous fistula, exposure to biohazardous material, extrusion, fibrosis, foreign body/unretrieved device fragment, hematuria, herniation, herniation of the device, incisional hernia, infection, infection/erosion, labia hematoma, limited urethral coaptation, loosey, overactive bladder, pelvic abscess, perforation, phlebitis, pressure atrophy at cuff site, pulmonary embolus, pump erosion, sphincter-related problem, superficial wound dehiscence, upper urinary tract deterioration, urethral atrophy, urethral injury, urethral stricture, urge de novo, urinary retention, urinary urgency, vaginal injury, wound dehiscence, wound seroma.

Table 1 summarizes the reported complications for 938 patients receiving an AMS 800 implant. The most common reported complication was bladder injury, which occurred in 5.8% of patients (n=54). There may have been more than one reported

complication for each patient.

This trial did not include AMS 800 devices impregnated with InhibiZone.

Table 1. Complications in Patients Studied

Complication	Percentage	Reference
Bladder injury	5.8%	2,8,11,17
Vaginal injury	4.2%	2,4,8,11,13
Cuff Erosion	3.4%	3,6,7,10,13,17
Infection/erosion	3.3%	1
Device Fluid leak	2.8%	2,3,9,12,13,15,16,17
Incisional hernia	1.5%	2
Urinary retention	1.3%	2,16
Urethral injury	1.2%	2,4,11
Urinary urgency	1.1%	2,6
Labia hematoma	1.0%	2,14
Infection	0.6%	1,4,8,10,17
Pelvic abscess	0.4%	16
Early post-operative morbidity	0.4%	1
Pressure atrophy at cuff site	0.3%	17
Cutaneous erosion	0.3%	4
Bothered by device	0.3%	2
Pulmonary embolus	0.1%	1
Enterocutaneous fistula	0.1%	1
Phlebitis	0.1%	2
Upper urinary tract deterioration	0.1%	13
Superficial wound dehiscence	0.1%	16
Device malfunction	0.1%	1,2,13,16,17
Device displacement	0.1%	3
Device damage – unrelated procedure	0.1%	6
Wound seroma	0.1%	9
Pump erosion	0.1%	13
Sphincter-related problem	0.1%	10
Loose cuff	0.1%	15,16

Success Rates

Patient success rates were reported and analyzed. For those defining clinical success^{2,4,5,6,7,12,13,16,17}, rates of success were categorized into the following criteria:

- Fully and socially continent = 0-1 pad use per day
- Occasional incontinence = Consistent use of 1-2 pads per day
- Fair/improved continence = Use of 3 pads per day
- Not improved = more than 3 pads per day

Success rates are summarized in Table 2.

Table 2. Success Rates

Success Criteria	% patients
Fully and socially continent	86.9%
Occasional incontinence	4.4%
Fair/improved	4.3%
Not improved	7.1

Patient Information Form (PIF) Study

A retrospective analysis was performed for female patients implanted (n=637) with the AMS 800 during the five-year period 2005-2010. The study examined Patient Information Form (PIF) data sent to the manufacturer by the implanting physician for original implants and revision surgeries. A revision is a surgical intervention related to the function, placement, or site reaction to the implanted device. PIF data is required to be on file by the manufacturer in order to be eligible for product replacement. Of the 637 female patients implanted during the period studied, 76 patients had undergone a total of 88 revision surgeries (1.16 revisions per patient). There were a total of 121 reported reasons for revisions. Note that more than one reason may be listed for the same revision. The data in Table 3 presents

the reported reasons for revision and the percentage of each reason for revision.

This trial did not include AMS 800 devices impregnated with InhibiZone.

Table 3. PIF Study Revision Data

Reason for revision	% revisions
Incontinence	19.8%
Not specified	19.8%
Infection	13.2%
Fluid loss	11.6%
Erosion	10.7%
Other medical condition	8.3%
Device out of specification/malfunction	5.8%
Improper size of device	3.3%
Other*	3.3%
Urinary retention	1.7%
Pain	1.7%

* Other includes: hematoma, device migration/malposition, patient dissatisfaction.

Clinical Information – Children

The AMS 800 is indicated for the treatment of urinary incontinence in children of both sexes. The safety and efficacy of the AMS 800 in children was evaluated by conducting a literature review. Four documents, reporting clinical information for patients under the age of 18, were analyzed and a summary of the clinical evidence is reported here.

Clinical efficacy has been reported for 322 children who have undergone implantation of the AMS 800. The mean age at the time of implant ranged from 9.9 to 14 years, and mean follow-up time reported was 5 to 15.4 years. Etiologies were similar to those reported for female patients, with the addition of sacral agenesis^{18,19,21}, lipoma^{18,21}, spinal cord abnormalities^{18,20,21}, sacralcoccygeal teratoma¹⁸ or cyst²¹, accidental injection of bicarbonate into the umbilical artery¹⁸, urethral injury²⁰, idiopathic¹⁹, lipomeningocele²¹, posterior urethral valves²¹, interstitial cystitis²¹, congenital absence of bladder neck²¹, and ureter malfunction²¹.

Reported complications were similar to those reported for female patients. Device revisions and bladder augmentation were the most commonly reported reason for additional surgeries^{18,19,21}. Reported success rates were comparable to those reported for female patients, with a reported range of 67% - 86%^{18,20,21}.

This trial did not include AMS 800 devices impregnated with InhibiZone.

Device Survival

Although it is not possible to predict exactly how long an implanted prosthesis will function in a particular patient, American Medical Systems, Inc. gathered data from two sources on device removals and revisions to help gain insight into product performance over time. Although these sources included information gathered from male patients, it is expected that similar data would result from analyzing data exclusively from and children. Figure 1 presents device survival results from a prospective clinical study and a Bayesian analysis that uses data from a prospective clinical study and a PIF Study to estimate device survival at five years.

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Bayesian Analysis – A Bayesian hierarchical model was used to evaluate device safety in the prospective clinical study. The Bayesian model estimated device survival using historical data

(PIF Study n=12,713) on the AMS 800 and prospective clinical study data (n=85) on the AMS 800. A log-normal distribution fit the AMS historical data. Based on the log-normal model, it was estimated that the five-year revision-free rate for the AMS 800 is approximately 73.8% with 95% CI ranging from 67.3% to 79.6%. The results met the primary safety endpoint for the clinical study of a five-year revision free rate at 75% using a 10% delta with two-sided 95% lower bound greater than 65%.

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Patient Counseling Information

Patients should be counseled in order to have a realistic expectation of the physical, psychological and functional outcome of the implantation. The risks, benefits and potential adverse events of all available treatment options should be discussed with the patient and considered by the physician and patient when choosing a treatment option.

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**Table 4. In Vitro Zone of Inhibition for
Device Samples* with InhibiZone Treatment**

Organism	Mean (mm)	S.D. (mm)	Number of Isolates
Staphylococcus epidermidis	22.6	2.9	21
Staphylococcus aureus	17.5	5.0	25
Escherichia coli**	6.5	2.6	24
Enterococcus faecalis**	4.8	6.7	21
Candida albicans**	0.1	0.4	21
Proteus mirabilis**	0.6	1.0	17

* obtained using standardized kink-resistant tubing test samples containing approximately 12µg minocycline HCl and 26 µg rifampin

** the isolates tested were not susceptible to rifampin and/or minocycline HCl control disks

An animal infection study was conducted using 11 rabbits. Five rabbits were implanted subcutaneously with 6 test samples each and five rabbits were implanted subcutaneously with 6 control samples each. One rabbit received three test samples and three control samples. The test samples were portions of an InhibiZone treated AMS 700 pump and the control samples were portions of a standard AMS 700 pump without InhibiZone. The AMS 700 Pumps used in the 700 Series of Inflatable Penile Prostheses is similar to AMS 800 Pumps used in the 800 Series of Artificial Urinary Sphincters in regards to material composition, adhesive, and InhibiZone application process. All samples were soaked in a 10^3 - 10^4 CFU solution of staphylococcus aureus, Sheretz strain for 8 hours. Samples were then allowed to dry for 30 minutes prior to surgical placement in the rabbit. After 2 days, all samples were removed and observed for growth on the samples. The number of coated samples that were infected was statistically significantly lower than the number of control samples that were infected.

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Scientific literature has included reports of adverse events and other observations in patients with implantable silicone devices. As reported, these events/observations indicate “allergic-like” symptoms and in other cases a symptom complex associated with immunological disorders. No causal relationship has been established between these events and silicone elastomer or fluorosilicone lubricant.

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Magnetic Resonance Imaging (MRI) Important Safety Information

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Static Magnetic Field	1.5 Tesla ^a	3.0 Tesla ^b
Highest Temperature Change	≤ + 0.4 °C	≤ + 2.0 °C
(a) 1.5T - 64 MHz MR System (General Electric Healthcare, Milwaukee, WI) (b) 3.0T MR Excite, General Electric Healthcare, software version 14X.M5		

Artifact Information

Non-clinical testing has demonstrated that the AMS 800 product line may compromise the MR image quality if the area of interest is relatively close to the position of the implant. The maximum image artifact produced by a MR gradient echo pulse sequence was a “moderate” localized signal void in size and shape of the implant. Optimization of MR imaging parameters to compensate for the presence of the device may be necessary.

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A Patient Information Form (PIF) must be filled out and filed with AMS at the time of implant to activate the product warranty. Contact your local AMS Representative prior to returning any product. Before returning any components, whether explanted or unused (sterile or nonsterile), customers must fill out the Return Goods Form located on the last page of the Patient Information Form. Follow **all** of the instructions on the form carefully, and be sure that the components have been thoroughly cleaned before returning them to AMS. Request a Return Kit from the AMS Customer Service Department to return any explanted components to AMS.

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References / Références / Literatur / Bibliografia / Referencias / Referências

- 1 Thomas K, Venn SN, Mundy AR, *Outcome of the Artificial Urinary Sphincter in Female Patients*, Journal of Urology, Apr 2002, Vol. 167, pp. 1720-1722
- 2 Costa P, Mottet N, Rabut B, Thuret R, Naoum KB, Wagner L, *The Use of the Artificial Urinary Sphincter in Women With Type III Incontinence and a Negative Marshall Test*, Journal of Urology, Apr 2001, Vol. 165:4, pp. 1172-1176
- 3 Costa P, Mottet N, Lanfrey P, Sahwi A, Dagues F, Louis JF, Navratil H, *Explantations and Revisions with Artificial Urinary Sphincter AMS 800 in Females*, Journal of Urology, May 1996, Vol. 155:5, Supplement
- 4 Richard F, Lefort JM, Bitker MO, Chartier-Kastler E, Chatelain C, *Female Incontinence with Primary Sphincter Deficiency: Results of Artificial Urinary Sphincter (AMS 800) with Long-Term Follow-Up*, Journal of Urology, 1996, Vol. 155:5, supplement
- 5 Costa P, Mottet N, Sandid ME, Sahwi A, Navratil H, *The Artificial Urinary Sphincter AMS 800: Our Experience in Women with or without Uterine Prolapse*, Acta Urol Bel, May 1995, Vol.63(2), pp. 45-46
- 6 Hadley R, Lioses P, Dickinson M, *Long-Term Follow-Up (2-5 Years) of Transvaginally Placed Artificial Urinary Sphincter by an Experienced Surgeon*, Journal of Urology, Apr 1995, Vol. 153, supplement
- 7 Stone KT, Diokno AC, Mitchell BA, *Just How Effective is the AMS 800 Artificial Urinary Sphincter? Results of Long-Term Follow-Up in Females*, Journal of Urology, Apr 1995, Vol. 153, supplement
- 8 Costa P, Mottet N, Pellec LL, Ducatez C, Boukaram M, El Sandid M, Naoum KB, Louis JF, Navratil H, *Artificial Urinary Sphincter AMS 800 in Operated and Unoperated Women With Type III Incontinence*, Journal of Urology, May 1994, Vol. 151:5, supplement
- 9 Karram MM, Rosenzweig, BA, Bhatia, NN, *Artificial Urinary Sphincter for Recurrent/Severe Incontinence in Women: Urogynecologic Perspective*, Journal of Reproductive Medicine, Oct 1993, Vol. 38, pp. 791-794
- 10 Duncan HJ, Nurse DE, Mundy AR, *Role of the Artificial Urinary Sphincter in Women*, British Journal of Urology, Feb 1992, Vol. 69, pp. 141-143
- 11 Salisz JA, Diokno AC, *The Management of Injuries to the Urethra, Bladder or Vagina Encountered During Difficult Placement of the Artificial Urinary Sphincter in the Female Patient*, Journal of Urology, Nov 1992, Vol. 148, pp. 1528-1530
- 12 Webster GD, Perez LM, Khoury JM, Timmons SL, *Management of Type III Stress Urinary Incontinence using Artificial Urinary Sphincter*, Urology, Jun 1992, Vol. XXXIX, Number 6, pp. 499-503
- 13 Parulkar BG, Barrett DM, *Application of the AS-800 Artificial Urinary Sphincter for Intractable Urinary Incontinence in Females*, Gynecology and Obstetrics, Aug 1990, Vol. 171, pp. 131-138
- 14 Abbassian A, *A New Operation for the Insertion of the Artificial Urinary Sphincter*, Journal of Urology, Sept 1988, Vol. 140, pp. 512-513
- 15 Appell RA, *Techniques and Results in the Implantation of the Artificial Urinary Sphincter in Women With Type III Stress Urinary Incontinence by a Vaginal Approach*, Neurourology and Urodynamics, 1998, Vol. 7, pp. 613-619
- 16 Diokno A., Hollander JB, Alderson TP, *Artificial Urinary Sphincter for Recurrent Female Urinary Incontinence: Indications and Results*, Journal of Urology, Oct 1987, Vol. 138, pp. 778-780
- 17 Light JK, Scott FB, *Management of Urinary Incontinence in Women with the Artificial Urinary Sphincter*, Journal of Urology, Sept 1985, Vol. 134, pp. 1528-1530
- 18 Simeoni J, Guys JM, Mollard P, Buzelin JM, Moscovici J, Bondonny JM, Melin Y, Lortat-Jacob S, Aubert D, Costa F, Galifer B, Debeugny P, *Artificial Urinary Sphincter Implantation for Neurogenic Bladder: A Multi-institutional Study in 107 Children*, British Journal of Urology, 1996, Vol. 78, pp. 287-293
- 19 Kryger JV, Barthold JS, Fleming P, Gonzalez R, *The Outcome of Artificial Sphincter Placement After a 15- Year Follow-up in a Paediatric Population*, British Journal of Urology International, 1999, Vol. 83, pp. 1026-1031
- 20 Castera R, Podesta ML, Ruarte A, Herrera M, Medel R, *10-Year Experience with Artificial Urinary Sphincter in Children and Adolescents*, Journal of Urology, June 2001, Vol. 165, pp. 2373-2376

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- 1 Thomas K, Venn SN, Mundy AR, *Outcome of the Artificial Urinary Sphincter in Female Patients*, Journal of Urology, Apr 2002, Vol. 167, pp. 1720-1722
- 2 Costa P, Mottet N, Rabut B, Thuret R, Naoum KB, Wagner L, *The Use of the Artificial Urinary Sphincter in Women With Type III Incontinence and a Negative Marshall Test*, Journal of Urology, Apr 2001, Vol. 165:4, pp. 1172-1176
- 3 Costa P, Mottet N, Lanfrey P, Sahwi A, Dagues F, Louis JF, Navratil H, *Explantations and Revisions with Artificial Urinary Sphincter AMS 800 in Females*, Journal of Urology, May 1996, Vol. 155:5, Supplement
- 4 Richard F, Lefort JM, Bitker MO, Chartier-Kastler E, Chatelain C, *Female Incontinence with Primary Sphincter Deficiency: Results of Artificial Urinary Sphincter (AMS 800) with Long-Term Follow-Up*, Journal of Urology, 1996, Vol. 155:5, supplement
- 5 Costa P, Mottet N, Sandid ME, Sahwi A, Navratil H, *The Artificial Urinary Sphincter AMS 800: Our Experience in Women with or without Uterine Prolapse*, Acta Urol Bel, May 1995, Vol.63(2), pp. 45-46
- 6 Hadley R, Lioses P, Dickinson M, *Long-Term Follow-Up (2-5 Years) of Transvaginally Placed Artificial Urinary Sphincter by an Experienced Surgeon*, Journal of Urology, Apr 1995, Vol. 153, supplement
- 7 Stone KT, Diokn AC, Mitchell BA, *Just How Effective is the AMS 800 Artificial Urinary Sphincter? Results of Long-Term Follow-Up in Females*, Journal of Urology, Apr 1995, Vol. 153, supplement
- 8 Costa P, Mottet N, Pellec LL, Ducatez C, Boukaram M, El Sandid M, Naoum KB, Louis JF, Navratil H, *Artificial Urinary Sphincter AMS 800 in Operated and Unoperated Women With Type III Incontinence*, Journal of Urology, May 1994, Vol. 151:5, supplement
- 9 Karram MM, Rosenzweig, BA, Bhatia, NN, *Artificial Urinary Sphincter for Recurrent/Severe Incontinence in Women: Urogynecologic Perspective*, Journal of Reproductive Medicine, Oct 1993, Vol. 38, pp. 791-794
- 10 Duncan HJ, Nurse DE, Mundy AR, *Role of the Artificial Urinary Sphincter in Women*, British Journal of Urology, Feb 1992, Vol. 69, pp. 141-143
- 11 Salisz JA, Diokno AC, *The Management of Injuries to the Urethra, Bladder or Vagina Encountered During Difficult Placement of the Artificial Urinary Sphincter in the Female Patient*, Journal of Urology, Nov 1992, Vol. 148, pp. 1528-1530
- 12 Webster GD, Perez LM, Khoury JM, Timmons SL, *Management of Type III Stress Urinary Incontinence using Artificial Urinary Sphincter*, Urology, Jun 1992, Vol. XXXIX, Number 6, pp. 499-503
- 13 Parulkar BG, Barrett DM, *Application of the AS-800 Artificial Urinary Sphincter for Intractable Urinary Incontinence in Females*, Gynecology and Obstetrics, Aug 1990, Vol. 171, pp. 131-138
- 14 Abbassian A, *A New Operation for the Insertion of the Artificial Urinary Sphincter*, Journal of Urology, Sept 1988, Vol. 140, pp. 512-513
- 15 Appell RA, *Techniques and Results in the Implantation of the Artificial Urinary Sphincter in Women With Type III Stress Urinary Incontinence by a Vaginal Approach*, Neurourology and Urodynamics, 1998, Vol. 7, pp. 613-619
- 16 Diokno A., Hollander JB, Alderson TP, *Artificial Urinary Sphincter for Recurrent Female Urinary Incontinence: Indications and Results*, Journal of Urology, Oct 1987, Vol. 138, pp. 778-780
- 17 Light JK, Scott FB, *Management of Urinary Incontinence in Women with the Artificial Urinary Sphincter*, Journal of Urology, Sept 1985, Vol. 134, pp. 76-78
- 18 Simeoni J, Guys JM, Mollard P, Buzelin JM, Moscovici J, Bondonny JM, Melin Y, Lortat-Jacob S, Aubert D, Costa F, Galifer B, Debeugny P, *Artificial Urinary Sphincter Implantation for Neurogenic Bladder: A Multi-institutional Study in 107 Children*, British Journal of Urology, 1996, Vol. 78, pp. 287-293
- 19 Kryger JV, Barthold JS, Fleming P, Gonzalez R, *The Outcome of Artificial Sphincter Placement After a 15- Year Follow-up in a Paediatric Population*, British Journal of Urology International, 1999, Vol. 83, pp. 1026-1031
- 20 Castera R, Podesta ML, Ruarte A, Herrera M, Medel R, *10-Year Experience with Artificial Urinary Sphincter in Children and Adolescents*, Journal of Urology, June 2001, Vol. 165, pp. 2373-2376

- 21 Herndon, CDA, Rink RC, Shaw MBK, Simmons GR, Cain MP, Keafer M, Casale AJ, *The Indiana Experience with Artificial Sphincters in Children and Young Adults*, Journal of Urology, Feb 2003, Vol. 169, pp. 650-654
- 22 Sherlock F, *MR Imaging of Metallic Implants and Materials: A Compilation of the Literature*, AJR, October 1988
- 23 Sherlock F, *MR Imaging and Biomedical Implants, Materials and Devices: An Updated Review*, Radiology, 1991, Vol. 180, pp. 541-550
- 24 Sherlock F, *MR Procedures and Biomedical Implants, Materials and Devices: 1993 Update*, Radiology, 1993, Vol. 189, pp. 587-599
- 25 Sherlock F, *MR Procedures and Metallic Objects: Update 1997*. Philadelphia, Lippincott-Raven, 1997, pp. 101, 110

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92116965-15

AMS 800™

Urinary Control System For Male Patients

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For Male Patients

Instructions For Use

NOTE: For implant procedure information, consult the AMS 800 Urinary Control System Operating Room Manual.

Brief Device Description

The AMS 800 Urinary Control System is an implantable, fluid filled, solid silicone elastomer device used to treat stress urinary incontinence. It is designed to restore the natural process of urinary control. The device simulates normal sphincter function by opening and closing the urethra at the control of the patient. The AMS 800 consists of three interconnected components: a cuff, a pump, and a pressure-regulating balloon (PRB). The three components are connected with kink-resistant tubing. The AMS 800 can be implanted at either the bulbous urethra or the bladder neck.

The AMS 800 Urinary Control System Cuff and Pump are available with InhibiZone™, an antibiotic coating of rifampin (rifampicin) and minocycline hydrochloride (minocycline HCl).^{*} For more information please refer to the Antibiotic Information section.

Indications for Use

The AMS 800 is used to treat urinary incontinence due to reduced outlet resistance (intrinsic sphincter deficiency).

Contraindications

1. This device is contraindicated in patients whom the physician determines to be poor candidates for surgical procedures and/or anesthesia due to physical or mental conditions.
2. This device is contraindicated in patients with urinary incontinence due to or complicated by an irreversibly obstructed lower urinary tract.
3. This device is contraindicated in patients with irresolvable detrusor hyperreflexia or bladder instability.
4. The implantation of the InhibiZone version of this device is contraindicated in patients with known allergy or sensitivity to rifampin or to minocycline HCl or other tetracyclines.
5. The implantation of products with InhibiZone is contraindicated in patients with systemic lupus erythematosus because minocycline HCl has been reported to exacerbate this condition.

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Warnings

1. Patients with urinary tract infections, diabetes, spinal cord injuries, open sores, or skin infections in the region of the surgery have an increased risk of infection associated with a prosthesis. Appropriate measures should be taken to reduce the likelihood of infection.

Infection that fails to respond to antibiotic therapy may

^{*} InhibiZone™ Antibiotic Surface Treatment is not available in all markets

result in removal of the prosthesis. Infection followed by explantation of the device may result in scarring which may make subsequent reimplantation more difficult.

2. Erosion may be caused by infection, pressure on the tissue, improper cuff sizing, improper balloon selection, tissue damage, and component misplacement. The cuff may erode around the urethra or bladder neck. The control pump may erode through the scrotum. The pressure-regulating balloon may erode into the bladder. Acute urinary tract infection can interfere with proper functioning of the device and may lead to erosion of the urethra in the cuff area. Failure to evaluate and promptly treat the erosion may result in a substantial worsening of the condition leading to infection and/or loss of tissue.
3. Poor bladder compliance or a small fibrotic bladder may require some measure of intervention including, in some cases, augmentation cystoplasty before implanting the prosthesis.
4. Patients with urge incontinence, overflow incontinence, detrusor hyperreflexia or bladder instability should have these conditions treated and controlled (or resolved) prior to implantation of the device.
5. Do not pass a catheter or any other instrument through the urethra without first deflating the cuff and deactivating the device to prevent potential damage to the urethra or the AMS 800.
6. This device contains solid silicone elastomers. This device does not contain silicone gel. The risks and benefits of implanting this device in patients with documented sensitivity to silicone should be carefully considered.
7. Surgical, physical, psychological, or mechanical complications, if they occur, may necessitate revision or removal of the prosthesis. Removal of the device without timely reimplantation of a new device may complicate subsequent reimplantation. The timing of reimplantation should be determined by the treating physician based on the patient's medical condition and history.
8. Product wear, component disconnection or other mechanical problems may lead to surgical intervention. Mechanical complications may include malfunctioning of the components and leakage of fluid. Any mechanical malfunction that does not permit the transfer of fluid from the cuff to the balloon may result in outflow obstruction. Mechanical events should be evaluated carefully by the treating physician and the patient should consider risks and benefits of treatment options, including revision surgery.
9. Previous patient history of adverse reaction(s) to radiopaque solution precludes its use as a filling medium for the prosthesis. Instead, saline should be used to fill the device.
10. The implanter should check that there is an adequate bulbospongiosus muscle to surround and support a bulbous urethral cuff implant. Thinner spongiosum typically occurs toward the distal end of the bulbous urethra, and implantation of the cuff where the spongiosum is thin increases the chance of erosion and other complications. This warning is especially important for double cuff implants, where the second cuff is placed distal to the first implanted cuff.
11. If a hypersensitivity reaction develops to a device coated with InhibiZone, the cuff and pump should be removed and the patient treated appropriately.

Precautions

Patient Related

1. Patient selection requires thorough preoperative consultation and evaluation by the physician.
2. Patients should be counseled in order to have a realistic expectation of the physical, psychological, and functional outcome of the implantation of an AMS 800. Although the prosthesis is designed to restore urinary control, some patients continue to have a degree of incontinence after this procedure.
3. Patients may experience pain when the device is activated in the postoperative period and during the period of initial use. Cases of chronic pain associated with device have been reported. Pain with a severity or duration beyond what is expected may require medical or surgical intervention. Patients should be counseled on expected postoperative pain including severity and duration.
4. Tissue fibrosis, previous surgery, or previous radiation therapy in the area of the implant may preclude implantation of a cuff at the bulbous urethra or bladder neck.
5. Any progressively degenerative disease, e.g. multiple sclerosis, may limit the future usefulness of the implanted prosthesis as a treatment for the patient's urinary incontinence.
6. Adequate manual dexterity, strength, motivation, and mental acuity are required for proper use of the device.
7. Trauma or injury to the pelvic, perineal or abdominal areas, such as impact injuries associated with sports, can result in damage to the implanted device and/or surrounding tissues. This damage may result in the malfunction of the device and may necessitate surgical correction including replacement of the device. The physician should advise patients of these possibilities and warn them to avoid trauma to these areas.
8. Consideration should be given to the diameter of the implanted occlusive cuff relative to catheters or other trans-urethral devices. When fully deflated, the inside diameter of the smallest occlusive cuff (3.5cm) generally exceeds 28F. Additional clearance is required to accommodate the patient's urethral tissue between the trans-urethral device and the occlusive cuff. Urethral tissue thickness is patient specific and requires a physician's assessment to determine its impact on sizing.

InhibiZone™ Related

1. Use of products with InhibiZone should be carefully considered in patients with hepatic or renal disease, as use of rifampin and minocycline HCl can cause additional stress on the hepatic and renal systems.
2. Patients who receive a device with InhibiZone and are also taking propofol should be carefully monitored for signs of renal toxicity.
3. Patients who receive a device with InhibiZone and are also taking warfarin should have their prothrombin time monitored because tetracyclines have been reported to slow coagulation.
4. Use of products with InhibiZone should be carefully considered in patients using thionamides, isoniazid and halothane, due to potential hepatic side effects that have been reported in patients using these drugs and higher doses of rifampin.

5. Devices with InhibiZone should not come into contact with ethyl alcohol, isopropyl alcohol or other alcohols, acetone or other nonpolar solvents. These solvents may remove the antibiotics from the device.
6. InhibiZone components should not be soaked in saline or other solutions prior to implantation. The components may be briefly rinsed or dipped into a sterile solution, immediately prior to implant, if desired.
7. InhibiZone does not replace your normal antibiotic protocols. Continue using any prophylactic protocols normally used for urological surgical procedures.
8. Because products with InhibiZone are impregnated with a combination of rifampin and minocycline HCl, the contraindications, warnings and precautions regarding the use of these antimicrobial agents apply and should be adhered to for the use of this device, although systemic levels of minocycline HCl and rifampin in patients receiving this device are unlikely to be detected.

Surgery Related

1. Improper cuff sizing, improper balloon selection, or other causes may result in tissue erosion, migration of components, or continued incontinence.
2. Component migration can occur if the cuff is sized improperly, if the pump or balloon is not positioned correctly, or if the tubing lengths are incorrect. Migration can result in pain, complications, device malfunction and surgical revision.
3. Unsuccessful outcomes may result from improper surgical technique, improper sterile technique, anatomical misplacement of components, improper sizing and/or filling of components.
4. Although reinforced tubing has been designed to be more resistant to tubing kinks, tubing kinks may still result from tailoring the connecting tubing to an improper length during the implant procedure.

Device Related

1. If the deactivation valve is closed when the cuff is inflated, fluid cannot transfer from the cuff to the balloon and sustained outflow obstruction may arise as a result:
 - a. In the event of large pressures within the bladder, automatic pressure relief that normally occurs with the device would be prevented. Cycling the device can relieve the outflow obstruction.
 - b. Cycling the device may be difficult if deactivation occurs when the pump bulb is deflated. If unable to cycle the device, squeezing the sides adjacent to the deactivation button will allow fluid to fill the pump bulb and then the pump can be cycled normally.
 - c. Release of the deactivation valve may require greater force than that used to cycle the device.
2. System pressure changes may occur over time if you fill the balloon with radiopaque solution of incorrect concentration. Follow the instructions in the Operating Room Manual to prepare the radiopaque solution with the correct concentration.

Adverse Events

The following adverse events have been associated with the use of this product: bladder spasms, bleeding, contracture, deep vein thrombosis, delayed wound healing, difficult activation, difficult

deactivation, dysuria, edema, exposure to biohazardous material, extrusion, fibrosis, fistula formation, foreign body/unretrieved device fragment, hematoma, hematuria, herniation, herniation of the device, hydrocele, impaired device function, infection, limited urethral coaptation, migration, nerve injury, overactive bladder, pain/discomfort, patient dissatisfaction, perforation, positional incontinence, recurrent incontinence, swelling, tissue erosion, tissue erosion/infection, urethral atrophy, urethral injury, urethral stricture, urge de novo, urinary retention, wound dehiscence, wound infection.

A prospective clinical study was conducted to demonstrate the safety and efficacy of the AMS 800 Urinary Control System. A total of 87 patients were enrolled in the study and 85 patients were implanted with the device. During the study, 26 patients experienced 43 device related adverse events. Table 1 lists the device related adverse events reported during the study.

This trial involved only devices without InhibiZone.

Table 1: AMS 800 Prospective Clinical Study Device Related Adverse Events

Adverse Event Category	Total Events	Patients with AE	Events Resolved	Interventions*		
				None Reported	Medical**	Surgical
Impaired Device Function	7	6	4	2	2	4
Pain/Discomfort	6	5	4	3	3	1
Delayed Wound Healing	5	5	5	2	3	0
Bladder Spasms	2	2	0	0	2	0
Difficult Activation	2	2	2	1	1	0
Migration	3	3	1	2	0	1
Tissue Erosion	2	2	2	0	0	2
Difficult Deactivation	1	1	1	0	1	0
Infection	2	2	2	0	0	2
Recurrent Incontinence	3	3	3	1	0	2
Fistula Formation	1	1	1	0	0	1
Hematoma	1	1	1	0	1	0
Swelling	2	2	2	0	2	1
Hydrocele	1	1	1	0	1	1
Tissue Erosion/Infection	1	1	1	0	0	1
Patient Dissatisfaction	1	1	1	0	0	1
Positional Incontinence	1	1	0	1	0	0
Wound Infection	1	1	1	0	1	0
Urinary Retention	1	1	1	0	1	0

*Events may have been addressed with more than one type of intervention.

**Medical interventions included: medication, education, frequent device deactivation, dressing changes and catheterization.

How Supplied and Storage

Store device in a clean, dry, dark area at room temperature.

WARNING: Contents supplied STERILE. Do not use if sterile barrier is damaged. If damage is found, call your AMS representative.

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or

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including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Clinical Studies

A prospective, multi-center, non-randomized clinical study was undertaken to demonstrate that the AMS 800 can be surgically implanted without serious adverse sequelae, provides an acceptable level of continence and enhances quality of life. Each patient served as their own control. Efficacy data and safety

data related to adverse events, revision surgery, diagnoses and health status evaluations were captured on case report forms. Patient self-evaluations related to health status and non-illness specific quality of life were measured on two validated outcome instruments. Patient and physician assessments of continence were measured on a recognized, standardized non-validated instrument.

Eighty-seven (87) male patients were enrolled in the study of which 85 patients were implanted with the device during the study. Patients available at the follow-up intervals were 6-months (n=67), 12-months (n=60), 18-months (n=55), and 24-months (n=41). Patients diagnosed with intrinsic sphincter deficiency (ISD) resulting from prostate surgery were eligible for enrollment. Patients with a history of allergy/sensitivity to silicone, pre-existing autoimmune or connective tissue disease or active urogenital infection were excluded from the study.

This trial involved only devices without InhibiZone.

Endpoints

The primary effectiveness endpoint evaluated the effect of the prosthesis on patient quality of life using the *Incontinence Impact Questionnaire*, an incontinence-specific quality of life questionnaire. The primary safety endpoint evaluated the five-year revision-free rate using a Bayesian hierarchical model. The safety endpoint was a five-year revision-free rate equivalent to 75% using a 10% delta with a two-sided 95% lower bound greater than 65%.

Incontinence Impact Scores

The primary effectiveness endpoint was a reduction in Incontinence Impact Score from pre- to post-implant status. Incontinence impact was measured pre- and post-implant at 6, 12, 18, and 24 months. Thirty-nine (39) patient answered the *Incontinence Impact Questionnaire* (IIQ) at 24-month follow-up. The IIQ is a 30-item, self-administered questionnaire designed to assess the impact of urinary incontinence on several subscales including physical, emotional, and social. The IIQ used in the study was developed from a validated instrument.¹ The mean pre-implant score was significantly higher ($p<0.0001$) than mean scores at all follow-up visits. Therefore, the impact of incontinence was reduced for patients following AMS 800 implantation and the primary objective was met.

Physician and Patient Assessment of Continence

Physician assessed continence was 63.6% dry and 34.1% required some additional protection at one-year follow-up (n=43). At two-year follow-up (n=30), 73.3% were dry and 23.3% required some additional protection. Patient assessed continence was 61.7% dry and 36.7% required some additional protection at one year follow-up (n=60). At two-year follow-up (n=41), 65.9% were dry and 31.7% required some additional protection. No significant difference existed between physicians' assessment and the patients' assessment of their incontinence.

Patient Evaluation of Health Status, and Self-esteem

General Quality of Life as measured by the *Health Status Questionnaire* and the *Rosenberg Self-esteem Questionnaire* was evaluated at pre- and post-implant at 6, 12, 18, and 24 months. Thirty-eight (38) patients answered the *Health Status* and *Rosenberg Self-esteem Questionnaires* at 24-month follow-up. The self-administered *Health Status Questionnaire*² was used to assess non-illness specific parameters such as physical functioning, social functioning, energy/fatigue, pain, health perception, and

emotional problems. A high score indicates that overall health was perceived to be high. The mean score was 596 at pre-implant and 612 at two-year follow-up. No significant difference in health status scores was observed during the study. The self-administered *Rosenberg Self-esteem Questionnaire*³ was used to assess changes in patient self-esteem. The range of possible scores is 0-6, with a score of 6 indicating high self-esteem. The mean score at implant was 3.5 and at two-year follow-up was 4.1. The increase in mean score indicates a more positive self-esteem following AMS 800 implant. The device did not have an adverse effect on sexual function. Some patients with improved continence following implant also reported increased sexual activity. The positive impact of the device on patient's lives measured in the clinical study is consistent with results obtained by other authors.^{4,5,6}

Surgical Revisions

A revision is a surgical intervention related to the function, placement, or site reaction to the implanted device. For the 85 patients implanted with the device followed under the prospective clinical study, 14 patients (16.5%) experienced a total of 15 revisions up to 24 months following implant. One patient had two revisions. Three (3) revisions were due to mechanical malfunction. Two (2) revisions were due to recurrent incontinence. Two (2) revisions were due to erosion. Two (2) revisions were due to infection. One (1) revision each (total = 6) was due to migration, pain, erosion/infection, persistent incontinence/patient dissatisfaction, recurring incontinence/malfunction, and infection/pain/urethrocutaneous fistula. Multiple reasons were provided for some revisions. Four of the 14 patients who experienced revisions elected to have the device removed and 10 elected to have the device replaced. The probability of remaining revision-free 24 months following implant was 79.5% (95% CI with 95% lower confidence bound 69.8%) based on the prospective clinical study.

Additional data on the number of surgical revisions and their reasons were collected under two retrospective studies. Each of these studies are briefly described below and comparisons of revision data of these two retrospective studies and the prospective study are presented in Tables 2 and 3.

Patient Information Form (PIF) Study – The PIF study was a retrospective analysis of patients implanted (n=12,713) in the U.S during the five-year period 1995-1999. The study examined PIF data voluntarily sent to the manufacturer by the implanting physician for original implants and revisions. PIF data is required to be on file with the manufacturer in order to be eligible for product replacement. Revision data presented in Table 2 and Table 3 below are based on a total of 2,116 revisions reported for 2,014 patients that required one or more revisions during the five-year period of the study.

Retrospective Clinical Study – The retrospective clinical study was an analysis of patients implanted (n=390) in the U.S. between 1995-1999. The study examined pre- and post-implant medical records and follow-up data collected by questionnaires and physician examinations. Post-implant data was available for 356 patients. The study followed patients for up to ten years (mean: 4.1 years; range: 0.03-10.3 years). The revision data presented in Table 2 and Table 3 below are based on a total of 317 revisions reported on 164 patients that required one or more revisions during the ten-year period of the study.

The data in Table 2 presents the percentage of patients revised

during the specified follow-up period, the average number of revisions conducted on patients requiring a revision and the number of revisions expected per 100 patients for these studies in comparison with the data of the prospective clinical study.

Table 2: Comparison of Revision Data from Three Different Clinical Studies

	Prospective Study (85 pts. over 24 months)	PIF Study (12,713 pts. over 5 years)	Retrospective Study (356 pts. over 9 years)
% pts. revised	16.5% (14/85)	15.8% (2014/12713)	46.1% (164/356)
avg.. # of revisions per pts. revised	1.07 (15/14)	1.05 (2116/2014)	1.93 (317/164)
# of revisions per 100 pts.	18 (15/85)	17 (2116/12713)	89 (317/356)

Table 3 shows revision data stratified by each reported reason for revision from three different studies of male patients implanted with the AMS 800. Under the PIF Study and Retrospective Study multiple reasons were sometimes provided for a single revision. Therefore, in order to stratify this revision data by reason, all occurrences were included and presented as “% reason.” The total number of reasons therefore exceeds the total number of revisions reported for these studies.

Table 3: Reasons for Revision in Three Different Studies

Revision Reason ^a	Prospective Study (n=85)		PIF Study (n=12713)		Retrospective Study (n=356)	
	revisions		reasons ^b		reasons	
	%	#	%	#	%	#
Infection	2.4%	(2)	2.3%	(297)	8.1%	(29)
Infection/erosion	1.2%	(1)	----	----	----	----
Erosion	2.4%	(2)	3.6%	(451)	22.5%	(80)
Recurring Incontinence	2.4%	(2)	5.7%	(724)	42.4%	(151)
Fluid Loss	----	----	2.3%	(298)	9.3%	(33)
Fluid Transfer Impaired	----	----	0.3%	(38)	----	----
Pressure too low	----	----	1.1%	(140)	----	----
Mechanical Malfunction	3.5%	(3)	0.7%	(89)	13.8%	(49)
Migration/ Malposition	3.5%	(3)	0.4%	(46)	4.8%	(17)
Iatrogenic Complications	----	----	0.4%	(51)	0.6%	(2)
Reimplantation/ Replacement	----	----	----	----	3.1%	(11)
Pain	1.2%	(1)	0.2%	(22)	1.4%	(5)
Patient Dissatisfaction	1.2%	(1)	0.2%	(27)	1.7%	(6)
Other ^c	----	----	2.4%	(305)	----	----
Not indicated	----	----	1.9%	(242)	----	----

a Note that some adverse events in the table such as fluid loss, pressure too low, fluid transfer impaired and malposition could fall into the category of mechanical malfunction or iatrogenic error. Since information is not available to place them in either category, they are listed separately.

b Numbers of reasons can vary for the same percentage due to rounding.

Регистрация МИ в Росздравнадзоре, www.nevacert.ru | info@nevacert.ru | inflatable cuff, pressure too high, unable to activate, unable to deactivate, atrophy, difficult to operate, urinary retention, air in the system, hematoma.

Device Survival

Although it is not possible to predict exactly how long an implanted prosthesis will function in a particular patient, American Medical Systems, Inc. gathered data from two sources on device removals and revisions to help gain insight into product performance over time. Figure 1 presents device survival results from the prospective clinical study and a Bayesian analysis that uses data from the prospective clinical study and the PIF

Study to estimate device survival at five years.

Prospective Clinical Study – A device survival curve was calculated from data collected during a prospective clinical study (n=85) with two-year follow-up. Using Kaplan-Meier analysis, the two-year revision-free rate for the AMS 800 was 79.5% (95% CI with 95% lower confidence bound 69.8%).

Bayesian Analysis – A Bayesian hierarchical model was used to evaluate device safety in the prospective clinical study. The Bayesian model estimated device survival using historical data (PIF Study n=12,713) on the AMS 800 and prospective clinical study data (n=85) on the AMS 800. A log-normal distribution fit the AMS 800 historical data. Based on the log-normal hierarchical model, it was estimated that the five-year revision-free rate for the AMS 800 is approximately 73.8% with 95% CI ranging from 67.3% to 79.6%. The results met the primary safety endpoint for the clinical study of a five-year revision free rate at 75% using a 10% delta with two-sided 95% lower bound greater than 65%.


Figure 1: Estimated Device Survival Rates for the AMS 800

Device Use

Retrospective Clinical Study – The study included male AMS 800 patients (n=390) implanted between 1987-1990. Data from this study was used to estimate the device use for patients (n=356) with available data through 1997 (range: 0.03-10.3 years). For the retrospective study, “device use” means the span of time from implant to removal, including revisions. Life table analysis indicated that the probability for a 9-year span of device use was 83.9%. Thirty-three (33) of the 356 patients had their device removed. The remaining 323 patients had a functioning device at last contact.

Patient Counseling Information

Patients should be counseled in order to have a realistic expectation of the physical, psychological and functional outcome of the implantation. The risks, benefits and potential of all available treatment options should be discussed with the patient and considered by the physician and patient when choosing a treatment option.

An appropriate patient history, including history of personality disorders, and diagnostic work-up should be a part of the patient decision making process.

Some patients may become dissatisfied by the presence of the prosthetic device in their body. This issue should be discussed with the patient prior to the surgery. Patient dissatisfaction may lead to device removal. Patients should also be aware that the AMS 800 is not considered to be a lifetime implant.

It is also important that the physician discusses with the patient the possibility of an allergic reaction to the materials in the device (See Silicone Information).

Antibiotic Information

The antibiotics present in InhibiZone, minocycline HCl and rifampin, are well characterized and have been in use for years. The dosage present on the artificial urinary sphincter is intended to act on organisms that attempt to colonize on the device.

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The following in vitro data are available, but their clinical significance is unknown. No clinical studies have been performed to evaluate the effect of the antibiotic surface treatment on reducing the incidence of artificial urinary sphincter implant infections.

Table 4: In vitro Zone of Inhibition for Device Samples* with InhibiZone Treatment

Organism	Mean (mm)	S.D. (mm)	Number of Isolates
Staphylococcus epidermidis	22.6	2.9	21
Staphylococcus aureus	17.5	5.0	25
Escherichia coli**	6.5	2.6	24
Enterococcus faecalis**	4.8	6.7	21
Candida albicans**	0.1	0.4	21
Proteus mirabilis**	0.6	1.0	17

*obtained using standardized kink-resistant tubing test samples containing approximately 12µg minocycline HCl and 26 µg rifampin

**the isolates tested were not susceptible to rifampin and/or minocycline HCl control disks

An animal infection study was conducted using 11 rabbits. Five rabbits were implanted subcutaneously with 6 test samples each and five rabbits were implanted subcutaneously with 6 control samples each. One rabbit received three test samples and three control samples. The test samples were portions of an InhibiZone treated AMS 700 pump and the control samples were portions of a standard AMS 700 pump without InhibiZone. The AMS 700 Pumps used in the 700 Series of Inflatable Penile Protheses is similar to AMS 800 Pumps used in the 800 Series of Artifical Urinary Sphincters in regards to material composition, adhesive, and InhibiZone application process. All samples were soaked in a 10³-10⁴ CFU solution of staphylococcus aureus, Sheretz strain for 8 hours. Samples were then allowed to dry for 30 minutes prior to surgical placement in the rabbit. After 2 days, all samples were removed and observed the samples. The number of coated samples that were infected was statistically significantly lower than the number of control samples that were infected.

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Solid silicone elastomers have been commonly used in a variety of biomedical devices for over 40 years. Silicone fluids have an extensive history of use in medical devices.

Scientific literature has included reports of adverse events and other observations in patients with implantable silicone devices. As reported, these events/observations indicate “allergic-like” symptoms and in other cases a symptom complex associated with immunological disorders. No causal relationship has been established between these events and silicone elastomer or fluoro-silicone lubricant.

There are reports of malignant tumor formation in laboratory animals only associated with implants of relatively large size. Many different materials are associated with this effect in animals, silicone elastomers among them. No such effect has been described in humans.

Extensive testing has been conducted on all materials that comprise the AMS 800. This testing has indicated no toxicological response attributable to the materials. However, some of the materials caused minor irritation when implanted in animals.

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Non-clinical testing has demonstrated the AMS 800 product line is MR Conditional. The device can be scanned safely under the following conditions:

Static Magnetic Field	1.5 Tesla ^a	3.0 Tesla ^b
Spatial Gradient Field	450 Gauss/cm or less	720 Gauss/cm or less
Maximum whole body averaged Specific Absorption Rate (SAR)	1.5 W/kg for 15 minutes of scanning as assessed by calorimetry	2.9 W/kg for 15 minutes of scanning as assessed by calorimetry

MRI-Related Heating

Non-clinical testing has demonstrated the AMS 800 product line produced the temperature rises during MRI performed for 15 minutes of scanning in the respective MR systems which would not pose a hazard to the human subject.

Static Magnetic Field	1.5 Tesla ^a	3.0 Tesla ^b
Highest Temperature Change	≤ + 0.4 °C	≤ + 2.0 °C
<div> <div> Регистрация МИ в Росздравнадзоре www.nevacert.ru info@nevacert.ru </div> <div> (a) 1.5T – 64 MHz MR System (General Electric Healthcare, Milwaukee, WI) (b) 3.0T MR Excite, General Electric Healthcare, software version 14X.M5 </div> </div>		

Artifact Information

Non-clinical testing has demonstrated that the AMS 800 product line may compromise the MR image quality if the area of interest is relatively close to the position of the implant. The maximum image artifact produced by a MR gradient echo pulse sequence was a “moderate” localized signal void in size and shape

of the implant. Optimization of MR imaging parameters to compensate for the presence of the device may be necessary.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	5,800 mm ²	1,956 mm ²	6,096 mm ²	2,650 mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

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A Patient Information Form (PIF) must be filled out and filed with AMS at the time of implant to activate the product warranty. Before returning any components, whether explanted or unused (sterile or nonsterile), customers must fill out the Return Goods Form located on the last page of the Patient Information Form.

Follow **all** of the instructions on the form carefully, and be sure that the components have been thoroughly cleaned before returning them to AMS. Request an AMS Product Return Kit from the AMS Customer Service Department to return any explanted components to AMS.

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References

- Shumaker SA, Wyman JF, Ubersax JS, McClish JA, Fantl JA. *Health-related Quality of Life Measures for Women with Urinary Incontinence: the Incontinence Impact Questionnaire and the Urogenital Distress Inventory*. Quality of Life Research, 3:291-306. 1994.
- Health Outcomes Institute, *Health Status Questionnaire 2.0*, 1993 and Rand Health Services Program, *RAND 36-item Health Survey 1.0*, 1986.
- Rosenberg M. *Society and the Adolescent Self-Image*. Princeton, New Jersey: Princeton University Press, 1965.
- Haab F, Trockman BA, Zimmern PE, Leach GE. Continence and Quality of Life after the Artificial Urinary Sphincter: Minimum 3.5 years follow-up. *Journal of Urology* 1997; 158:435-439.
- Litwiller SE, Kim KB, Fone PD, DeVere White RW, Stone AR. Urinary Incontinence and the Artificial Urinary Sphincter: A Long-term Study of Patient Satisfaction and Criteria for Success. *Journal of Urology* 1996; 156:1975-80.
- Fleshner N, Herschorn S. The Artificial Urinary Sphincter for Post-radical Prostatectomy Incontinence: Impact on Urinary Symptoms and Quality of Life. *Journal of Urology* 1996; 155:1260.

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Références

- Shumaker SA, Wyman JF, Ubersax JS, McClish JA, Fantl JA. *Health-related Quality of Life Measures for Women with Urinary Incontinence: the Incontinence Impact Questionnaire and the Urogenital Distress Inventory*. Quality of Life Research, 3:291-306. 1994.
- Health Outcomes Institute, *Health Status Questionnaire 2.0*, 1993 and Rand Health Services Program, *RAND 36-item Health Survey 1.0*, 1986.
- Rosenberg M. Society and the Adolescent Self-Image. Princeton, New Jersey: Princeton University Press, 1965.
- Haab F, Trockman BA, Zimmern PE, Leach GE. Continence and Quality of Life after the Artificial Urinary Sphincter: Minimum Follow-up. *Journal of Urology* 1997; 158:435-439.
- Litwiller SE, Kim KB, Fone PD, DeVere White RW, Stone AR. Post-Prostatectomy Incontinence and the Artificial Urinary Sphincter: A Long-term Study of Patient Satisfaction and Criteria for Success. *Journal of Urology* 1996; 156:1975-80.
- Fleshner N, Herschorn S. The Artificial Urinary Sphincter for Post-radical Prostatectomy Incontinence: Impact on Urinary Symptoms and Quality of Life. *Journal of Urology* 1996; 155:1260.

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Literatur

- 01 Shumaker SA, Wyman JF, Ubersax JS, McClish JA, Fantl JA. *Health-related Quality of Life Measures for Women with Urinary Incontinence: the Incontinence Impact Questionnaire and the Urogenital Distress Inventory*. Quality of Life Research, 3:291-306. 1994.
- 02 Health Outcomes Institute, *Health Status Questionnaire 2.0*, 1993 and Rand Health Services Program, *RAND 36-item Health Survey 1.0*, 1986.
- 03 Rosenberg M. *Society and the Adolescent Self-Image*. Princeton, New Jersey: University Press, 1965.
- 04 Haab F, Trockman BA, Zimmern PE, Leach GE. Continence and Quality of Life after the Artificial Urinary Sphincter: Minimum 3.5 years follow-up. *Journal of Urology* 1997; 158:435-439.
- 05 Litwiller SE, Kim KB, Fone PD, DeVere White RW, Stone AR. Post-Prostatectomy Incontinence and the Artificial Urinary Sphincter: A Long-term Study of Patient Satisfaction and Criteria for Success. *Journal of Urology* 1996; 156:1975-80.
- 06 Fleshner N, Herschorn S. The Artificial Urinary Sphincter for Post-radical Prostatectomy Incontinence: Impact on Urinary Symptoms and Quality of Life. *Journal of Urology* 1996; 155:1260.

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Bibliografia

- 01 Shumaker SA, Wyman JF, Ubersax JS, McClish JA, Fantl JA. *Health-related Quality of Life Measures for Women with Urinary Incontinence: the Incontinence Impact Questionnaire and the Urogenital Distress Inventory*. Quality of Life Research, 3:291-306. 1994.
- 02 Health Outcomes Institute, *Health Status Questionnaire 2.0*, 1993 and Rand Health Services Program, *RAND 36-item Health Survey 1.0*, 1986.
- 03 Rosenberg M. Society and the Adolescent Self-Image. Princeton, New Jersey: Princeton University Press, 1965.
- 04 Haab F, Trockman BA, Zimmern PE, Leach GE. Continence and Quality of Life after the Artificial Urinary Sphincter: Minimum 3.5 years follow-up. *Journal of Urology* 1997; 158:435-439.
- 05 Litwiller SE, Kim KB, Fone PD, DeVere White RW, Stone AR. Post-Prostatectomy Incontinence and the Artificial Urinary Sphincter: A Long-term Study of Patient Satisfaction and Criteria for Success. *Journal of Urology* 1996; 156:1975-80.
- 06 Fleshner N, Herschorn S. The Artificial Urinary Sphincter for Post-radical Prostatectomy Incontinence: Impact on Urinary Symptoms and Quality of Life. *Journal of Urology* 1996; 155:1260.

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Referencias

- Shumaker SA, Wyman JF, Ubersax JS, McClish JA, Fantl JA. *Health-related Quality of Life Measures for Women with Urinary Incontinence: the Incontinence Impact Questionnaire and the Urogenital Distress Inventory*. Quality of Life Research, 3:291-306. 1994.
- Health Outcomes Institute, *Health Status Questionnaire 2.0*, 1993 and Rand Health Services Program, *RAND 36-item Health Survey 1.0*, 1986.
- Rosenberg M. Society and the Adolescent Self-Image. Princeton, New Jersey: Princeton University Press, 1965.
- Haab F, Trockman BA, Zimmern PE, Leach GE. Continence and Quality of Life after the Artificial Urinary Sphincter: Minimum 3.5 years follow-up. *Journal of Urology* 1997; 158:435-439.
- Litwiller SE, Kim KB, Fone PD, DeVere White RW, Stone AR. *Urinary Incontinence and the Artificial Urinary Sphincter: A Long-term Study of Patient Satisfaction and Criteria for Success*. *Journal of Urology* 1996; 156:1975-80.
- Fleshner N, Herschorn S. The Artificial Urinary Sphincter for Post-radical Prostatectomy Incontinence: Impact on Urinary Symptoms and Quality of Life. *Journal of Urology* 1996; 155:1260.

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Referências

- Shumaker SA, Wyman JF, Ubersax JS, McClish JA, Fantl JA. *Health-related Quality of Life Measures for Women with Urinary Incontinence: the Incontinence Impact Questionnaire and the Urogenital Distress Inventory*. Quality of Life Research, 3:291-306. 1994.
- Health Outcomes Institute, *Health Status Questionnaire 2.0*, 1993 and Rand Health Services Program, *RAND 36-item Health Survey 1.0*, 1986.
- Rosenberg M. *Society and the Adolescent Self-Image*. Princeton, New Jersey: Princeton University Press, 1965.
- Haab F, Trockman BA, Zimmern PE, Leach GE. Continence and Quality of Life after the Artificial Urinary Sphincter: Minimum 3.5 years follow-up. *Journal of Urology* 1997; 158:435-439.
- Litwiler SE, Kim KB, Fone PD, DeVere White RW, Stone AR. *Urinary Incontinence and the Artificial Urinary Sphincter: A Long-term Study of Patient Satisfaction and Criteria for Success*. *Journal of Urology* 1996; 156:1975-80.
- Fleshner N, Herschorn S. The Artificial Urinary Sphincter for Post-radical Prostatectomy Incontinence: Impact on Urinary Symptoms and Quality of Life. *Journal of Urology* 1996; 155:1260.

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- Shumaker SA, Wyman JF, Ubersax JS, McClish JA, Fantl JA. *Health-related Quality of Life Measures for Women with Urinary Incontinence: the Incontinence Impact Questionnaire and the Urogenital Distress Inventory*. Quality of Life Research, 3:291-306. 1994.
- Health Outcomes Institute, *Health Status Questionnaire 2.0*, 1993 and Rand Health Services Program, *RAND 36-item Health Survey 1.0*, 1986.
- Registration of MI in the Rosdraznabazore
www.nevacert.ru | info@nevacert.ru Society and the Adolescent Self-Image. Princeton, New Jersey: Princeton University Press, 1965.
- Haab F, Trockman BA, Zimmern PE, Leach GE. Continence and Quality of Life after the Artificial Urinary Sphincter: Minimum 3.5 years follow-up. *Journal of Urology* 1997; 158:435-439.
- Litwiller SE, Kim KB, Fone PD, DeVere White RW, Stone AR. Post-Prostatectomy Incontinence and the Artificial Urinary Sphincter: A Long-term Study of Patient Satisfaction and Criteria for Success. *Journal of Urology* 1996; 156:1975-80.
- Fleshner N, Herschorn S. The Artificial Urinary Sphincter for Post-radical Prostatectomy Incontinence: Impact on Urinary Symptoms and Quality of Life. *Journal of Urology* 1996; 155:1260.

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©2017 Boston Scientific Corporation or its affiliates.
All rights reserved. All trademarks are the property
of the respective owners.


92116966-15 (2017-05)

92116966-15

AMS™

AMS 800™

Urinary Control System

For Male, Female, and Pediatric Patients

Operating Room Manual

English

TABLE OF CONTENTS

CONTENTS

DISCLOSURES

Indications For Use	1
Contraindications	1
Warnings	1
Precautions	2
Patient-Related	2
InhibiZone-Related	3
Surgery-Related	3
Device-Related	4
InhibiZone TM	
Antibiotic Surface Treatment	4

DEVICE DESCRIPTION

System Description	5
System Operation	5
Components	6
Packaging	6
Individual Components	6
Occlusive Cuff	6
Pressure-Regulating Balloon	7
Control Pump	7
Accessory Kit	8
Connectors	9
Quick Connect Assembly Tool (Optional)	9
Tubing Passers (Optional)	9
Deactivation Kit	9

DEVICE STERILIZATION AND STORAGE

Care and Storage of Presterilized Components	10
Sterilization	10
How to Sterilize AMS Tools	10
Component Resterilization	10

SURGICAL PREPARATION

Operating Room Preparation	11
Preoperative Surgical Team Preparation	11
Preoperative Patient Preparation	11
Supplies and Instrument Requirements	12
Open Components	13
Prepare Hemostats	14
Prepare System Components	14
Filling Solutions	14
Prepare Control Pump	16
Prepare Pressure-Regulating Balloon	17
Prepare Cuff	18

BULBOUS URETHRA CUFF PLACEMENT

Surgical Procedures	20
Bulbous Urethra Cuff Placement -	
Perineal Approach	20
Implant Pressure Regulating Balloon	21
Implant Pump	22
Connect Tubing	22
Deactivate	22
Bulbous Urethra Cuff Placement -	

Transverse Scrotal Approach	23
Implant Pressure Regulating Balloon	25
Implant Pump	26
Connect Tubing	26
Deactivate	26
Attaching a Second Cuff to the System	27
Components and Accessories	28
Choose the Correct Size of Cuff	28
Fill the Components	28
Place a Second Cuff	28
Connect Tubing	29
Deactivate	29

BLADDER-NECK CUFF PLACEMENT

Bladder-Neck Cuff Placement	30
Implant Pressure Regulating Balloon	30
Implant Pump	31
Connect Tubing	31
Deactivate	31
Transvaginal Approach	31

CONNECT TUBING

Connect Tubing	32
Using AMS Quick Connect Sutureless Window	
Connectors	32
Using AMS Suture-Tie Connectors	33
Deactivation Kit	35

DEVICE OPERATION

Deactivate Cuff	36
Activate (Reactivate) Cuff:	
Normal Method	36
Activate (Reactivate) Cuff:	
Optional Methods	37
Side Squeeze Method	37
Cotton Swab Method	37
Valve Block Bending Method	37

POST OPERATIVE CARE

Immediately Postoperative	38
After Release From Hospital	38
Evaluating Long-Term Function and Placement	38

DOCUMENTATION

Completing the	
Patient Information Form	39
Inventory Returns and	
Replacement Information	39
Document Information	39

TROUBLESHOOTING

Total Device	40
Occlusive cuff	40
Pressure-Regulating Balloon	40
Control Pump	40

DISCLOSURES

This section lists the following disclosures:

- Indications For Use
- Contraindications
- Warnings
- Precautions

INDICATIONS FOR USE

The AMS 800 is used to treat urinary incontinence due to reduced urethral/bladder outlet resistance (intrinsic sphincter deficiency) in males, females and children.

CONTRAINDICATIONS

1. This device is contraindicated in patients whom the physician determines to be poor candidates for surgical procedures and/or anesthesia due to physical or mental conditions.
2. This device is contraindicated in patients with urinary incontinence due to or complicated by an irreversibly obstructed lower urinary tract.
3. This device is contraindicated in patients with irresolvable detrusor hyperreflexia or bladder instability.
4. The implantation of the InhibiZone™ version of this device is contraindicated in patients with known allergy or sensitivity to rifampin (rifampicin) or to minocycline hydrochloride (minocycline HCl) or other tetracyclines.
5. The implantation of products with InhibiZone is contraindicated in patients with systemic lupus erythematosus because minocycline HCl has been reported to aggravate this condition.

WARNINGS

1. Patients with urinary tract infections, diabetes, spinal cord injuries, open sores, or skin infections in the region of the surgery have an increased risk of infection associated with a prosthesis. Appropriate measures should be taken to reduce the likelihood of infection. Infection that fails to respond to antibiotic therapy may result in removal of the prosthesis. Infection followed by explantation of the device may result in scarring which may make subsequent reimplantation more difficult.
2. Erosion may be caused by infection, pressure on the tissue, improper cuff sizing, improper balloon selection, tissue damage, and component misplacement. The cuff may erode


Figure 1-1.
AMS 800 implanted
in a male patient


Figure 1-2.
AMS 800 implanted
in a female patient

DISCLOSURES (CONTINUED)

around the urethra or bladder neck. The control pump may erode through the scrotum. The pressure-regulating balloon may erode into the bladder. Acute urinary tract infection can interfere with proper functioning of the device and may lead to erosion of the urethra in the cuff area. Failure to evaluate and promptly treat the erosion may result in a substantial worsening of the condition leading to infection and/or loss of tissue.

3. Poor bladder compliance or a small fibrotic bladder may require some measure of intervention including, in some cases, augmentation cystoplasty before implanting the prosthesis.
4. Patients with urge incontinence, overflow incontinence, detrusor hyperreflexia or bladder instability should have these conditions treated and controlled (or resolved) prior to implantation of the device.
5. Do not pass a catheter or any other instrument through the urethra without first deflating the cuff and deactivating the device to prevent potential damage to the urethra or the AMS 800.
6. This device contains solid silicone elastomers. This device does not contain silicone gel. The risks and benefits of implanting this device in patients with documented sensitivity to silicone should be carefully considered.
7. Surgical, physical, psychological, or mechanical complications, if they occur, may necessitate revision or removal of the prosthesis. Removal of the device without timely reimplantation of a new device may complicate subsequent reimplantation. The timing of reimplantation should be determined by the treating physician based on the patient's medical condition and history.
8. Product wear, component disconnection or other mechanical problems may lead to surgical intervention. Mechanical complications may include malfunctioning of the components and leakage of fluid. Any mechanical malfunction that does not permit the transfer of fluid from the cuff to the balloon may result in outflow obstruction. Mechanical events should be evaluated carefully by the treating physician and the patient should consider risks and benefits of treatment options,

including revision surgery.

9. Previous patient history of adverse reaction(s) to radiopaque solution precludes its use as a filling medium for the prosthesis. Instead, saline should be used to fill the device.
10. Female patients with persistent incontinence should be evaluated to rule out vesicovaginal fistula, which may have resulted from an unrecognized iatrogenic injury.
11. If a hypersensitivity reaction develops to a device coated with InhibiZone, the cuff and pump should be removed and the patient treated appropriately.

PRECAUTIONS

Patient-Related

1. Patient selection requires thorough preoperative consultation and evaluation by the physician.
2. Patients should be counseled in order to have a realistic expectation of the physical, psychological, and functional outcome of the implantation of an AMS 800. Although the prosthesis is designed to restore urinary control, some patients continue to have a degree of incontinence after this procedure.
3. Patients may experience pain when the device is activated in the postoperative period and during periods of initial use. Cases of chronic pain associated with device have been reported. Pain with a severity or duration beyond what is expected may require medical or surgical intervention. Patients should be counseled on expected postoperative pain including severity and duration.
4. Tissue fibrosis, previous surgery, or previous radiation therapy in the area of the implant may preclude implantation of a cuff at the bulbous urethra or bladder neck.
5. Any progressively degenerative disease, e.g. multiple sclerosis, may limit the future usefulness of the implanted prosthesis as a treatment for the patient's urinary incontinence.
6. Adequate manual dexterity, strength, motivation, and mental acuity are required for proper use of the device.

DISCLOSURES (CONTINUED)

7. Trauma or injury to the pelvic, perineal or abdominal areas, such as impact injuries associated with sports, can result in damage to the implanted device and/or surrounding tissues. This damage may result in the malfunction of the device and may necessitate surgical correction including replacement of the prosthesis. The physician should advise patients of these possibilities and warn them to avoid trauma to these areas.
 8. Female patients of child-bearing age must be forewarned that pregnancy is acceptable, but cesarean section may be indicated to minimize the risk of damage to the bladder neck and its surrounding cuff. For those patients who become pregnant, deactivation of the device during the third trimester is recommended to reduce the risk of erosion. Those patients contemplating pregnancy should consider delaying implantation.
 9. Children receiving device should be evaluated at regular intervals. Lifelong radiological and urodynamic surveillance of the urinary tract is crucial. Before implantation, the patient and their family should be informed of the complication rate and the need for long-term follow-up.
 10. Consideration should be given to the diameter of the implanted occlusive cuff relative to catheters or other trans-urethral devices. When fully deflated, the inside diameter of the smallest occlusive cuff (3.5cm) generally exceeds 28F. Additional clearance is required to accommodate the patient's urethral tissue between the trans-urethral device and the occlusive cuff. Urethral tissue thickness is patient specific and requires a physician's assessment to determine its impact on sizing.
4. Use of products with InhibiZone should be carefully considered in patients using thionamides, isoniazid and halothane, due to potential hepatic side effects that have been reported in patients using these drugs and higher doses of rifampin.
 5. Devices with InhibiZone should not come into contact with ethyl alcohol, isopropyl alcohol or other alcohols, acetone or other non-polar solvents. These solvents may remove the antibiotics from the device.
 6. InhibiZone components should not be soaked in saline or other solutions prior to implantation. The components may be briefly rinsed or dipped into a sterile solution, immediately prior to implant, if desired.
 7. InhibiZone does not replace your normal antibiotic protocols. Continue using any prophylactic protocols normally used for urological surgical procedures.
 8. Because products with InhibiZone are impregnated with a combination of rifampin and minocycline HCl, the contraindications, warnings and precautions regarding the use of these antimicrobial agents apply and should be adhered to for the use of this device, although systemic levels of minocycline HCl and rifampin in patients receiving this device are unlikely to be detected.

InhibiZone-Related

1. Use of products with InhibiZone should be carefully considered in patients with hepatic or renal disease, as use of rifampin and minocycline HCl can cause additional stress on the hepatic and renal systems.
2. Patients who receive a device with InhibiZone and are also taking methoxyflourane should be carefully monitored for signs of renal toxicity.
3. Patients who receive a device with InhibiZone

Surgery-Related

1. Improper cuff sizing, improper balloon selection, or other causes may result in tissue erosion, migration of components, or continued incontinence.
2. Component migration can occur if the cuff is sized improperly, if the pump or balloon is not positioned correctly, or if the tubing lengths are incorrect. Migration can result in pain, complications, device malfunction and surgical revision.
3. Unsuccessful outcomes may result from improper surgical technique, improper sterile technique, anatomical misplacement of components, improper sizing and/or filling of components.

DISCLOSURES (CONTINUED)

4. Although reinforced tubing has been designed to be more resistant to tubing kinks, tubing kinks may still result from tailoring the connecting tubing to an improper length during the implant procedure.

Device-Related

1. If the deactivation valve is closed when the cuff is inflated, fluid cannot transfer from the cuff to the balloon and sustained outflow obstruction may arise as a result:
 - a. In the event of large pressures within the bladder, automatic pressure relief that normally occurs with the device would be prevented. Cycling the device can relieve the outflow obstruction.
 - b. Cycling the device may be difficult if deactivation occurs when the pump bulb is deflated. If unable to cycle the prosthesis, squeezing the sides adjacent to the deactivation button will allow fluid to fill the pump bulb and then the pump can be cycled normally.
 - c. Release of the deactivation valve may require greater pressure than that used to cycle the device.
2. System pressure changes may occur over time if you fill the balloon with radiopaque solution of incorrect concentration. Follow the instructions in the Operating Room Manual to prepare the radiopaque solution with the correct concentration.

institution.

The AMS patented antibiotic surface treatment process uses a formulation of minocycline HCl and rifampin.

The AMS 800 components are treated with very low levels of antibiotics. AMS provides numerous completed configurations of the AMS 800 to individualize treatment, however, while the AMS 800 PRB is not IZ treated, a complete device (PRB, pump, and one or two cuffs), regardless of configuration, contains ≤ 6.5 mg rifampin and ≤ 8 mg minocycline HCl. This represents less than 2% of oral dose exposure for a complete course of rifampin or minocycline HCl with the maximum dose calculated from the means and 95% tolerance interval.

In vitro studies with the antibiotic treated device material and susceptible strains of *Staphylococcus epidermidis* and *Staphylococcus aureus* show a microbial “zone of inhibition” around the test material. The clinical significance of these in vitro data is unknown. A limited animal model study suggests that this surface treatment may reduce the potential for bacterial colonization of the treated device.

No clinical studies have been performed to evaluate the effect of the antibiotic surface treatment reducing the incidence of sphincter implantation infections.

INHIBIZONETM ANTIBIOTIC SURFACE TREATMENT

American Medical Systems has a proprietary process to impregnate antibiotics into the tissue-contacting surfaces of the urinary control system. This InhibiZone Antibiotic Surface Treatment innovation is intended to elute the antibiotics from the device surface when exposed to a warm, moist environment. In in vitro testing using susceptible organisms, this elution provided antibiotic action both on the surface and in a zone surrounding the treated device.

Existing prophylactic antibiotic protocols should be maintained as determined by the physician and or

DEVICE DESCRIPTION

SYSTEM DESCRIPTION

The AMS 800 prosthesis consists of three components: an occlusive cuff, a control pump, a pressure-regulating balloon, and connectors. The cuff can be implanted at either the bulbous urethra (**Figure 1-1**) or bladder neck (**Figure 1-2**) in men and adolescent males. In females and children, the cuff is placed at the bladder neck. (**Figure 1-2**)


Figure 1-1.
AMS 800 implanted
with cuff at
bulbous urethra


Figure 1-2.
AMS 800 implanted
with cuff at
bladder neck

SYSTEM OPERATION

The AMS 800 simulates normal sphincter function by opening and closing the urethra at the control of the patient. When the cuff is closed, urine stays in the bladder. (**Figure 2-1a**)

When the patient wishes to void, he or she squeezes and releases the pump, implanted in the scrotum or labium, several times. This causes the fluid in the cuff to move from the cuff into the pressure-regulating balloon. The cuff opens and urine passes through the urethra. (**Figure 2-1b**) The balloon then automatically repressurizes the cuff due to the automatic returning of the fluid by the pressure-regulating balloon to the cuff. (**Figure 2-1c**)

The control pump is designed to allow the urologist to deactivate the implanted device without additional surgery (see instructions on how to Deactivate the Cuff).

System Operation Shown in a Male


Figure 2-1a. Urethra closed by cuff


Figure 2-1b. Opening of cuff and
urination by manipulation of pump


Figure 2-1c. Automatic cuff re-closure

DEVICE DESCRIPTION (CONTINUED)

COMPONENTS

This section provides a brief description of the following AMS 800 Urinary Control System components and accessories:

- ✓ Occlusive Cuff
- ✓ Pressure-regulating Balloon
- ✓ Control Pump
- ✓ Accessory Kit (Including Connectors)
- ✓ Quick Connect Assembly Tool (Optional)
- ✓ Tubing Passers (Optional)
- ✓ Deactivation Kit (Optional)

Cuffs and pumps are available with InhibiZone™ Antibiotic Surface Treatment, which is an antibiotic surface treatment of rifampin and minocycline HCl.

Packaging

All components and accessories are sterilized products, except for the non-sterilized designated tools (Quick Connect Tool and tubing passers). Each sterile component is packaged separately within an inner plastic tray which is sterile and sealed with a Tyvek™ lid. The inner tray lies within the sterile environment of a plastic outer tray which is also sealed by a Tyvek lid. This double package is further protected by an outer box for transportation. Adhesive labels on one end of the outer box and on the Tyvek lid of the outer tray identify the components, their sizes, and serial and lot numbers.

INDIVIDUAL COMPONENTS

Occlusive Cuff

Implanted at either the bulbous urethra (male) or the bladder neck (males, females, children), the occlusive cuff occludes the urethra by applying pressure circumferentially. It is made of silicone elastomer and is available in thirteen sizes ranging from 3.5 cm to 11 cm in length (all cuffs are approximately 1.8 cm wide when deflated). The surgeon determines the proper size to be used in the patient by measuring the circumference of the tissue around the urethra intraoperatively. The cuff tubing is joined to the control pump tubing with a connector.


Figure 2-2.
Occlusive
Cuff

DEVICE DESCRIPTION (CONTINUED)

Pressure-Regulating Balloon

The pressure-regulating balloon, implanted in the prevesical space, controls the amount of pressure exerted by the occlusive cuff. It is also made of silicone elastomer and is provided in the following three ranges.

- 51-60 cm H₂O
- 61-70 cm H₂O
- 71-80 cm H₂O

Upon activation the balloon pressure will be within the stated range. The surgeon usually selects the lowest balloon pressure needed to maintain closure of the bladder neck or bulbous urethra. The balloon's tubing is also linked to tubing from the control pump by a connector.

Control Pump

The control pump is implanted in the soft tissue of the scrotum or labium. It is approximately 1.3 cm wide and 3.5 cm long. The upper part of the control pump (the valve block) contains the resistor and valves needed to transfer fluid between the components. The bottom half of the control pump is a bulb that the patient squeezes and releases to transfer fluid in order to void.

The deactivation button is located in the upper portion of the control pump. When the deactivation button is pressed, it stops fluid from being transferred between the components. This feature allows the physician to leave the cuff open:

- During the postoperative healing period
- For transurethral procedures

The AMS 800 has color-coded tubing to help surgeons make the correct connections between the components:

- Clear tubing connects to cuff
- Black tubing connects to balloon

The tubes are connected with connectors. There are two types of connectors: AMS straight or curved Quick Connect Sutureless Window Connectors. They may also be made with straight or curved suture connectors, which are secured with 3-0 polypropylene (non-absorbable) ties.


Figure 2-3.
Pressure-
Regulating
Balloon


Figure 2-4.
Control Pump

DEVICE DESCRIPTION (CONTINUED)

Accessory Kit

The Accessory Kit for the AMS 800 contains accessory materials necessary for one implant procedure. The Accessory Kit includes the following disposable items:

- Cuff Sizer (**Figure 2-5**)
 - One disposable cuff sizer
- Special Purpose Needles (**Figure 2-6**)
 - Two 15-gauge blunt needles (for filling components)
 - Two 22-gauge blunt needles (for flushing the opening of the tubing and inside of the connectors)
- Protection Tubes
 - Two 30 cm lengths of silicone tubing (for covering the tips of the hemostats used to prepare components—to protect the tubing from damage by the hemostat jaws.)
- Tubing Connectors and Collets
 - Three straight suture-tie connectors
 - Two right angle suture-tie connectors
 - One 3-way suture tie connector (Y-connector)
 - Three straight quick-connect window connectors
 - Two right angle quick-connect window connectors
 - One 3-way quick-connect window connector (Y-connector)
 - Eight collets (on a collet holder, used with Quick Connectors).

The Cuff Sizer measures the urethral circumference. When selecting a cuff size, the cuff length refers to the outside diameter of the cuff when it encircles the urethra.

Note: The inside circumference of the cuff is somewhat smaller than the outside circumference of the cuff.

Use the 15-gauge needles to fill the components. They fit snugly into the lumen of the tubing to prevent fluid from leaking during the filling process. The 22-gauge needles are small enough to fit within the tubing to flush away air and blood before a connection is made.

Place the extra silicone tubing found in the Accessory Kit over the tips of the hemostats used for clamping tubing or handling device components (see Preparing Hemostats).

Figure 2-5. Cuff Sizer

Figure 2-6. Blunt-Tip Needles

DEVICE DESCRIPTION (CONTINUED)

Connectors

Two systems for connecting components and tubing are available for the AMS 800.

- Suture-Tie Connectors
- Quick-Connectors

The Accessory Kit contains three straight and two right angle AMS Suture-Tie Connectors and one 3-way “Y” connector. (**Figure 2-7**) These connectors are attached with permanent, non-absorbable sutures (3-0 polypropylene).

The Accessory Kit also contains AMS Quick Connect Sutureless Window Connectors - three straight, two right angle, and one 3-way “Y” connector (**Figure 2-8**), eight collets (or locking rings) (**Figure 2-9**) and one collet holder. To use the Quick Connectors, a Quick Connection Tool (**Figure 2-10**) must be ordered, which is sterilized at the site.

CAUTION: AMS Quick Connect Sutureless Window Connectors should not be used in revision procedures involving previously implanted component tubing. Changes in the tubing over time may reduce the fixation performance of the Quick Connect Sutureless Window Connectors. The Quick Connect system may be used when all previously implanted components are removed and replaced with new components.

Quick Connect Assembly Tool (Optional)

In order to use the quick-connectors in the Accessory Kit, you must order the AMS Quick Connect Assembly Tool. This is a reusable stainless steel instrument used to close the connectors. The Quick Connect Assembly Tool is shipped non-sterile and may be resterilized. Quick-connectors should not be used on revision surgeries, except when all previously implanted components are removed and replaced with new components.

Caution: Sterilize the AMS Quick Connect Assembly Tool properly before use.

Tubing Passers (Optional)

Tubing Passers (**Figure 2-11**) are used to route the tubing of the components from one incision site to another incision site. The machined end of the passer provides a snug fit in the tubing lumen. Tubing passers are shipped in non-sterile condition and must be sterilized before use.

Caution: Sterilize the tubing passers properly before use.

Deactivation Kit

The Deactivation Kit is an optional package not normally required for an initial implant. It contains 3 stainless steel plugs and 1 straight suture-tie connector (**Figure 2-12**). The 3 plugs are often useful during revision surgical procedures to protect the inside of the preserved components and the fluid from contamination.


Figure 2-7. Suture Tie Connectors

Figure 2-8. Quick Connect Sutureless Window Connectors


Figure 2-9. Collets


Figure 2-10. Quick Connect Assembly Tool

Figure 2-11. Tubing Passers


Figure 2-12. Deactivation Kit

DEVICE STERILIZATION AND STORAGE

CARE AND STORAGE OF PRESTERILIZED COMPONENTS

All the components and accessories of this device are sterilized products. Use before the expiration date of the sterilized products.

CAUTION: *Designated tools such as the Quick Connect Tool and tubing passers are shipped without sterilization. Re-sterilize them before each use.*

To protect the integrity of the packaging and function of the AMS 800 prosthesis, store the presterilized components on a protected shelf or in a cabinet. The environment should be clean, dry, and near room temperature. For maximum protection during storage, leave the component trays within their dustcover boxes.

- Store in a clean, dark, and cool place.
- Avoid water, direct sunlight, high temperatures, humidity, and ultraviolet light.
- Avoid abnormal altitudes, temperatures, humidity, aeration, sunlight, dust, salt and ions in the air.
- Do not store in a place where there are chemicals or gas.

The versions of the AMS 800 components with InhibiZone Antibiotic Surface Treatment are light and temperature sensitive. Care should be taken to store the products according to the instructions on the package.

CAUTION: *Do not store products with InhibiZone above 40°C (104°F).*

CAUTION: *Do not use this product past its expiration date.*

STERILIZATION

How to Sterilize AMS Tools

American Medical Systems does not sterilize the AMS Tubing Passers or the AMS Quick Connect Assembly Tool. These instruments are shipped in steam sterilization packages ready for hospital sterilization.

Before using, sterilize the optional Tubing Passers and Quick Connector Assembly Tool by steam autoclave. For reprocessing information, refer to the instructions manual provided with the tools. Thoroughly clean these instruments after each implant procedure so they will be ready for resterilization before the next implant procedure.

Component Resterilization

DO NOT re-sterilize the sterilized components or accessories. Only designated, non-sterilized tools (Tubing Passers and AMS Quick Connect Tool) can be re-sterilized.

Do not resterilize the AMS Suture-Tie Connectors or AMS Quick Connect Sutureless Window Connectors. Only the two stainless steel instruments used in conjunction with the AMS 800 prosthesis may be resterilized: the AMS Tubing Passers and the AMS Quick Connect Assembly Tool.

SURGICAL PREPARATION

OPERATING ROOM PREPARATION

Preoperative Surgical Team Preparation

Before beginning an AMS 800 implant procedure, the surgeon and the OR staff should be familiar with the device and the equipment needed as well as the steps in the procedure.

Both the surgeon and the OR staff should review the information in this manual prior to surgery. In addition, there should be an OR manual in the operating room during the surgery as a quick reference.

In many cases, it is helpful for the surgeon and/or the OR staff to observe an AMS 800 prosthesis implant surgery. This allows them to become familiar with the OR set-up and the procedure before actually being involved with one.

Just prior to surgery, the surgical team should scrub for ten minutes using povidone-iodine soap or the approved hospital scrub procedure.

Preoperative Patient Preparation

Before the surgery, many doctors prescribe prophylactic antibiotics for the patient. This is considered useful to reduce the risk of infection.

It is also important that the surgeon discuss with the patient the possibility of an allergic reaction to the materials in the device. Scientific literature has included reports of adverse events and other observations in patients with implantable silicone devices. As reported, these events/observations indicate “allergic-like” symptoms and, in other cases, a symptom complex associated with immunological disorders. No causal relationship has been established between these events and silicone elastomer.

Once the patient is in the operating room, the abdominal and genital area should be shaved. Following the shave, the area should be scrubbed with povidone-iodine soap for ten minutes or the approved hospital preoperative scrub procedure.

For bulbous urethra cuff placement using the perineal approach, the patient is placed in the lithotomy position, prepped, and draped for both a perineal and an abdominal incision. For bulbous urethra placement using the transverse scrotal approach, place the patient in the supine position with legs gently abducted on arm boards or other spreaders.

For bladder-neck cuff placement, the approach for male patients requires a set-up similar to that used for suprapubic prostatectomy. The patient is placed in the supine position, prepped, and draped for an abdominal incision. Access to the perineum should be maintained.

For bladder-neck cuff placement in female patients, the patient can be placed in the lithotomy or supine position, depending on the surgical approach, prepped, and draped.

SURGICAL PREPARATION (CONTINUED)

SUPPLIES AND INSTRUMENT REQUIREMENTS

A number of conventional surgical tools and supplies are needed for implantation of the AMS 800 prosthesis. Each surgeon may also have a preference about what should be available.

The following equipment, among other standard operating room materials, should be available:

- Sterile stainless steel tray
- 1000cc graduate container
- 500cc graduate container
- Sponge bowl
- Medicine cup
- Emesis basin
- Two 30cc disposable syringes
- One 10cc disposable syringe
- 8 silicone-shod hemostats
- Straight, clean and sharp scissors
- Hegar dilators
- Babcock clamps
- Asepto™ syringe
- Antibiotic solution
- Catheter
- Umbilical tape
- Vaginal pack
- Sterile saline for rinsing gloves & filling components
- Retractor (transverse scrotal approach)
- Extended nasal speculum (optional)
- Rectal tube (optional)
- Centimeter ruler (optional)

Use a plastic-draped Mayo stand or a stainless steel tray as a station for handling and filling the components of the prosthesis. Be certain that the components do not come into contact with paper or cloth drapes.

Submerge the filled components of the prosthesis in a storage basin containing sterile water until they are implanted.

Position splash basins so that the surgeons can conveniently clean their gloves during the surgical implantation procedure, especially before they make the tubing connections.

CAUTION: Be certain that the components do not come into contact with paper or cloth drapes. Fragments of these can cause possible obstruction of fluid flow if they enter the device.

SURGICAL PREPARATION (CONTINUED)

OPEN COMPONENTS

Always keep the sterile trays in their dust cover boxes when moving the components from storage to the operating room. The components should only be opened in the operating room when the surgeon directs the surgical team to do so.

Note: The AMS 800 Accessory Kit must always be opened first in order to prepare the hemostats and syringes used for preparing and filling components.

To open the components in the operating room, follow this procedure:

1. Remove trays from dust covers by opening boxes at the security tab.

Note: The circulating nurse should record the part and serial/lot numbers as well as the size of the components on the AMS Patient Information Form (PIF). Use the labels provided with the device.

Note: The adhesive label at one end of the dust cover box and the small, removable labels on the side of the plastic trays contain the name of the component (English), the part and serial/lot numbers, as well as the size of the components. This information is also listed on the lid of the outer tray (Tyvek sheet).

2. Remove inner tray from outer tray using following method:

- Peel back lid completely from tray in one motion
- Continue to hold outer tray without touching sterile inner tray.
- Have scrub nurse use his or her index finger (not thumb) to carefully lift inner tray up and out of outer tray.
- Scrub nurse should place inner tray on a plastic-draped Mayo stand in the sterile field.

3. Open inner tray just before preparing components using the following method:

- Remove components from inner trays by peeling back lids of sterile inner trays.
- Carefully remove components from tray.
- Place them in their appropriate positions on plastic-draped Mayo stand.

SURGICAL PREPARATION (CONTINUED)

PREPARE HEMOSTATS

To protect component tubing from being damaged by the mosquito hemostat jaws (**Figure 4-1**), cover the jaws with the silicone tubing provided in the Accessory Kit. Complete the following instructions to cover the jaws:

1. Place tubing on both jaws of hemostats up to box lock.
2. Completely cover all teeth on both jaws of six hemostats.
3. Trim tubing at jaw tip with clean, well-cutting scissors.
4. Reserve scissors for use as tubing scissors throughout procedure.

When using the hemostats, clamp jaws together only to the first click to prevent excessive pressure on tubing. (Do not advance more than one click.)


Figure 4-1. Tubing-Shod Hemostat Jaws

PREPARE SYSTEM COMPONENTS

Filling Solutions

The fluid used to fill the prosthesis must be sterile and completely free of particulate matter. The presence of any foreign materials in the fluid can affect the operation of the prosthesis. The solution must also be isotonic to minimize the transfer of fluid across the silicone membrane, which is semi-permeable. **Normal saline is the recommended isotonic solution to use when filling the prosthesis.**

However, if contrast media is preferred, one of the tested solutions in the table below may be used for filling. If you do not use the contrast media in the mixture proportions, you may alter the isotonicity of the mixtures and promote the formation of particulate matter.

Note: The products listed below are some of the radiographic solutions tested by American Medical Systems for use in AMS devices, only sterile water should be used for dilution. For a complete list contact American Medical Systems.

CAUTION: Do not use sterile saline or lactated Ringer's solution to dilute the contrast solutions.

WARNING: Contrast media are contraindicated if the patient has an iodine allergy.

The presence of particulate matter in the filling solution can affect the operation of the prosthesis. The filling solution must be free of blood and debris at all times. The solution (sterile saline) is isotonic with intracellular fluid to minimize the transfer of fluid across the semi-permeable silicone membrane.

SURGICAL PREPARATION (CONTINUED)

Contrast Media		Dilution		Manufacturer	Validated for InhibiZone Use
Conray 43	30 cc Conray 43	+	60 cc sterile H ₂ O	Mallinckrodt	Yes
Cysto Conray II	60 cc Cysto Conray II	+	15 cc sterile H ₂ O	Mallinckrodt	Yes
Hypaque-Cysto	60 cc Hypaque-Cysto	+	58 cc sterile H ₂ O	Nycomed	No
Isovue 200	60 cc Isovue 200	+	23 cc sterile H ₂ O	Bracco	No
Isovue 300	57 cc Isovue 300	+	60 cc sterile H ₂ O	Bracco	No
Isovue 370	38 cc Isovue 370	+	60 cc sterile H ₂ O	Bracco	No
Omnipaque 180	60 cc Omnipaque 180	+	14 cc sterile H ₂ O	Nycomed	No
Omnipaque 240	60 cc Omnipaque 240	+	38 cc sterile H ₂ O	Nycomed	No
Omnipaque 300	57 cc Omnipaque 300	+	60 cc sterile H ₂ O	Nycomed	Yes
Omnipaque 350	48 cc Omnipaque 350	+	60 cc sterile H ₂ O	Nycomed	No
Telebrix 12	53 cc Telebrix 12	+	47 cc sterile H ₂ O	Laboratoire Guerbel	Yes
Use an equivalent ratio dye with sterile water for a larger total volume.					

When filled, each component should be submerged in a storage basin of normal saline to prevent contact between the component and foreign materials. Be certain that the saline used as filling solution is separate from the storage basin where components are submerged after being filled and prepared.

Those components that are labeled as being treated with InhibiZone Antibiotic Surface Treatment should NOT be submerged in sterile normal saline.

CAUTION: Soaking antibiotic impregnated devices in saline will cause the antibiotics to diffuse off the device into the solution. This will cause the solution to turn orange and will reduce the concentration of antibiotics on the device.

CAUTION: It is important to maintain isotonic status of the fluid. Use sterile saline solution for primary cases. In a revision case, when radiopaque solution was used for the primary case and the surgeon preserved some of the primary components, use the same radiopaque solution at the same density. Device performance may deteriorate if different types of fluid solutions or different density solutions were mixed and the isotonicity was lost or particles were formed.

SURGICAL PREPARATION (CONTINUED)

PREPARE CONTROL PUMP

Complete the following instructions to prepare the control pump. A syringe is not necessary to fill the control pump.

1. Place end of each tubing in a basin with the appropriate filling solution. (**Figure 4-2a**)
2. Hold the pump at a 45 degree angle with the black tube on top.
3. Squeeze and release pump bulb repeatedly until all air in the pump and tubing has been displaced with fluid.

Note: It is important to keep tubing submerged during filling procedure.

Note: If air bubbles remain in pump bulb, continue to squeeze and release pump bulb to remove them; they will be evacuated through black color-coded tubing.

4. While keeping tubing submerged, use tubing-shod hemostats to clamp (one notch only) each tube 4 cm or 5 cm from end. (**Figure 4-2b**)

CAUTION: Do not advance hemostat's ratchet more than one notch. Excessive pressure will permanently damage the tubing.

Note: The clear pump tubing will connect to the clear cuff tubing, and black tubing will connected to black balloon tubing.

5. For a non InhibiZone treated pump, submerge the filled pump in a storage basin containing filling solution until implantation.
6. For a pump treated with InhibiZone Antibiotic Surface Treatment, place the pump onto an empty sterile tray or kidney basin and cover with a sterile drape. Prior to implanting, the pump should be inspected for entrapped air.

CAUTION: Soaking antibiotic impregnated devices in saline will cause the antibiotics to diffuse off the device into the solution. This will cause the solution to turn orange and will reduce the concentration of antibiotics in the device.


Figure 4-2a. Place Ends of Both Tubes Into Filling Solution


Figure 4-2b. Clamp Tubing With Ends In Filling Solution


Figure 4-3a. Aspirate Air

SURGICAL PREPARATION (CONTINUED)

PREPARE PRESSURE-REGULATING BALLOON

Complete the following instructions to prepare the Pressure-Regulating Balloon (PRB):

1. Attach 15-gauge blunt tip needle to 30 cc syringe.
2. Fill syringe with approximately 25 cc of filling solution.
3. Hold air-filled balloon with one hand and squeeze it until it is deflated.
4. Insert needle into end of balloon tubing.
5. Aspirate any air left in balloon until you feel resistance on syringe plunger. **(Figure 4-3a)**

Note: Always hold syringe upright: that is, needle pointing down, plunger pointing up.

6. Fill balloon with 20 cc of filling solution. **(Figure 4-3b)**
7. Rotate balloon until all air bubbles are gathered into one bubble. **(Figures 4-3c and 4-3d)**
8. Hold balloon by tubing to draw air bubble into tubing adapter area.
9. Aspirate air first, and then all fluid with syringe held upright until all air has been removed from balloon. Check that no air remains in balloon or tubing.

CAUTION: Do not over-aspirate pressure-regulating balloon because air can be drawn into system through semi permeable silicone membrane.

10. Hold plunger upright to maintain pressure.
11. Clamp tubing (one notch only) with tubing-shod hemostat 3 cm below the needle. **(Figure 4-3e)**

CAUTION: Do not advance hemostat's ratchet more than one notch. Excessive pressure will permanently damage the tubing.

12. Submerge empty balloon in a storage basin containing filling solution until implantation.

CAUTION: Do not place any hemostats on top of balloon. Any instruments on the balloon can damage it.


Figure 4-3b. Fill Balloon


Figure 4-3c. Rotate Balloon


Figure 4-3d. Collect Air Bubbles Into One


Figure 4-3e. Clamp Tubing

SURGICAL PREPARATION (CONTINUED)

PREPARE CUFF

Complete the following instructions to prepare the occlusive cuff:

1. Attach 15-gauge blunt tip needle to 30 cc syringe.
2. Fill syringe with approximately 10 cc of filling solution.
3. Hold air-filled cuff with one hand and squeeze it until it is deflated.
4. Insert needle into end of cuff tubing.
5. Aspirate any air left in cuff until you feel a slight resistance on the plunger. (**Figure 4-4a**)
6. Fill the cuff with 1 cc to 5 cc of recommended filling solution, depending on cuff size. (**Figure 4-4b**)

CAUTION: Do not over fill the cuff. (Figures 4-4c)
Excessive fluid can stretch cuff material.

Note: The amount of fluid that is needed to fill cuff will vary depending on cuff size - larger cuffs require more fluid than smaller cuffs.

7. Place syringe onto Mayo stand and hold both ends of cuff. Collect all small air bubbles into one large air bubble.
8. Roll end of cuff with your thumb to squeeze large air bubble into tubing adapter. (**Figure 4-4d**)


Figure 4-4a. Aspirate Air


Figure 4-4b. Insert Needle Into Tube


Figure 4-4c. Do Not Overfill Cuff

Figure 4-4d. Roll End of Cuff

SURGICAL PREPARATION (CONTINUED)

9. Hold syringe upright again (needle pointing down and plunger pointing up). Maintain pressure on cuff while pulling back on plunger to first aspirate remaining air bubble from cuff and then fluid.
10. If air remains in either cuff or tubing, repeat above steps.
11. With the air and fluid removed and cuff completely empty, use two tubing-shod mosquito hemostats to double-clamp cuff tubing (one notch only) 3 cm below needle and again 3 cm below first hemostat. **(Figure 4-4e)**
12. For a cuff treated with InhibiZone Antibiotic Surface Treatment, place the cuff onto an empty sterile tray or kidney basin and cover with a sterile drape. Prior to implanting, the cuff should be inspected for entrapped air.
13. For a non InhibiZone treated cuff, submerge the cuff into a kidney basin of sterile normal saline until the surgeon is ready to implant the cuff.


Figure 4-4e. Clamp Cuff

Caution: Soaking antibiotic impregnated devices in saline will cause the antibiotics to diffuse off the device into the solution. This will cause the solution to turn orange and will reduce the concentration of antibiotics on the device.

BULBOUS URETHRA CUFF PLACEMENT

SURGICAL PROCEDURES

There are several approaches to implantation of the AMS 800. It is important for the operating room staff to know which approach the surgeon intends to use because the approach will influence the position of the patient, instrumentation, selection of components, and the surgical sequence. The following brief descriptions give a general idea of the approaches.

BULBOUS URETHRA CUFF PLACEMENT - PERINEAL APPROACH

Complete the following steps to place the cuff at the bulbous urethra:

1. Place a Foley catheter or a 20Fr sound into urethra to help identify it during dissection.
2. Make a midline perineal incision and bluntly dissect bulbocavernous muscle from around the bulbous urethra.

(Figure 4-6)

Note: Use blunt dissection.

3. Place cuff sizer (or Penrose drain) around urethra where cuff is to be implanted. It should fit snugly without constricting urethra.

(Figure 4-7)

Note: If catheter or sound is in urethra, remove it before measuring the urethra.

Note: Do not stretch cuff sizer before use.

Note: Surgeon should use his or her judgment in choosing an appropriate cuff size; the measuring tape only provides approximate measurement of bulbous urethra circumference. The inside circumference of cuff is somewhat smaller than the outside circumference of cuff.

4. Select cuff size that corresponds to measured length.

Note: A cuff length of 4.0 cm or 4.5 cm is typically required for bulbous urethra placement.

Note: Cuff length is outside circumference of cuff when it encircles urethra.

5. Prepare cuff for implantation.
6. Position cuff at the implant site with the mesh backing toward the outside and the inflatable side toward the urethra.
7. Pass prepared cuff, tab first, under urethra.

Note: Tubing should be on the pump side.

CAUTION: To avoid damage to the cuff, grasp the cuff tab with a tubing-shod hemostat.


Figure 4-5. Bulbous Urethra Placement


Figure 4-6. Dissect Bulbocavernous Muscle


Figure 4-7. Measure Urethral Circumference

BULBOUS URETHRA CUFF PLACEMENT (CONTINUED)

8. Pass tubing through tab hole as follows:

- Pass the end of the cuff tubing through the hole until the hemostat meets the hole. Clamp a second hemostat onto the cuff tubing on the opposite side of the hole, and then release the first hemostat, so air does not enter the cuff.
- Pull the remainder of the tubing through the hole and close the cuff by pulling the tab over the tubing adapter (button). Ensure that the edges of the hole fit into the slot of the adapter.
- Rotate the cuff so the adapter is lateral to the urethral midline and position the tubing to avoid contact with the cuff.


Figure 4-8. Separate Linea Alba

Implant Pressure Regulating Balloon

9. Select appropriate pressure-regulating balloon.
10. Make a suprapubic incision, divide rectus fascia transversely, and use a spreading motion to separate the linea alba to reach prevesical space. (**Figure 4-8**)
11. Use blunt dissection to create a space for balloon.
12. Position the balloon in prevesical space.
13. Flush balloon tubing end using 22-gauge needle on 10 cc syringe filled with filling solution.
14. Connect 15-gauge needle on 30 cc syringe filled with filling solution to balloon and unclamp tubing.
15. Fill the balloon with 22 cc of the chosen filling solution.
Note: Larger cuff sizes may require more filling solution. Refer to step 17 for Cuff Pressurization Option.
16. Clamp (one notch only) the tubing approximately 3 cm from the end with a tubing-shod hemostat. (**Figure 4-9**)


Figure 4-9. Fill PRB and Clamp Tubing

CAUTION: Do not advance the hemostat's ratchet more than one notch. Excessive pressure will permanently damage tubing.

17. Cuff Pressurization Option: Larger cuff sizes may require more filling solution, which can be provided by pressurizing the cuff as follows:
 - Flush balloon and cuff tubing.
 - Make temporary connection with a suture-tied connector between balloon and cuff.
 - Unclamp; wait one minute; then reclamp (one notch only).
 - Disconnect balloon and cuff.
 - Flush balloon tubing, unclamp and aspirate filling solution.
 - Refill with 20 cc of filling solution and clamp.

BULBOUS URETHRA CUFF PLACEMENT (CONTINUED)

Implant Pump

18. Use blunt dissection to create a dependent subdartos pouch in the scrotum. **(Figure 4-10)**

Note: Control pump should be placed on same side of the patient as the pressure-regulating balloon.

19. Place pump into scrotal pouch with deactivation button facing outward so that it is palpable. **(Figure 4-11)**

20. Route the tubing to abdominal incision.

Note: The pump tubing should be above rectus muscle and fascia in abdominal incision.

Connect Tubing

Follow the instructions given in the section on connecting the tubing. Check the connection by activating the system (closing/opening the cuff).

Deactivate

The device system must be left in the deactivated mode for four to six weeks following implantation. Deactivate the system by following the instructions given in the section on “Deactivating the System”.

Figure 4-10. Create Pouch


Figure 4-11. Place Pump Into Pouch

BULBOUS URETHRA CUFF PLACEMENT (CONTINUED)

BULBOUS URETHRA CUFF PLACEMENT - TRANSVERSE SCROTAL APPROACH

Complete the following steps to place the device through a scrotal incision:

1. Place the patient supine with his legs gently abducted on arm boards or other spreaders. Knees and hips are not flexed. This avoids placing the urethra on a stretch as occurs in the lithotomy position. **(Figure 4-12)** Place a urethral catheter to drain the bladder to avoid injury during placement of the pressure regulating balloon.
2. Make an upper transverse scrotal incision and deepen it through the subcutaneous tissue. Move the incision up onto the penis and stabilize with a retractor and blunt hooks at 1, 3, 5, 7, 9 and 11 o'clock. These stay hooks secure the scrotal incision onto the penis and help prevent unnecessary dissection in the scrotum. The scrotal incision allows excellent access to the proximal bulbar urethra and retropubic and dartos spaces, and leaves the bulbocavernosus muscle intact. **(Figure 4-13)**
3. Expose the tunica albuginea of both corpora cavernosa sharply. Pass the Metzenbaum scissors proximally along the ventral surface of the tunica to the proximal corpora. **(Figure 4-14)**
4. When deep exposure of the proximal corpora is secured, a Deaver retractor is placed on the side of the urethra for caudal retraction. **(Figure 4-15)**

Figure 4-12.


Figure 4-13.


Figure 4-14.

Figure 4-15.

BULBOUS URETHRA CUFF PLACEMENT (CONTINUED)

5. This is repeated on the contralateral side exposing the scrotal septum. **(Figure 4-16)**

Figure 4-16.

6. The scrotal septum is then sharply dissected off the bulbar urethra. **(Figure 4-17)**

Figure 4-17.

7. By sharply dissecting webs of Buck's fascia binding the diverging corpora cavernosum to the corpora spongiosum, mobilization of the urethra can readily be achieved. **(Figure 4-18)**


Figure 4-18.

8. Because the patient is in the supine position, the urethra is mobile and a right angle clamp can be used to conduct the posterior dissection of the urethra almost under direct vision. The right angle clamp is spread to create sufficient space for the placement of the occlusive cuff. **(Figure 4-19)**

Figure 4-19.

BULBOUS URETHRA CUFF PLACEMENT (CONTINUED)

9. Measure the urethra and place the proper size occlusive cuff around the circumference.
(Figure 4-20)

Figure 4-20.

Implant Pressure Regulating Balloon

10. There are two ways to place the pressure-regulating balloon (PRB).
- Empty the bladder, remove the Scott retractor and displace the scrotal incision to the side of the penis. Place the PRB in the retro-pubic space by locating the inguinal ring and piercing the transversalis fascia. After PRB implantation, narrow the opening by placing an absorbable suture.
(Figure 4-21)


Figure 4-21.

- Displace the scrotal incision over the inguinal area and locate the inguinal ring. Use finger dissection to develop a pouch beneath the rectus but anterior to the transversalis fascia (cephalad to the inguinal ring). This avoids the necessity of piercing the fascia in patients with scarred retroperitoneum from radiation or radical surgery. After PRB implantation, narrow the opening by placing an absorbable suture.

BULBOUS URETHRA CUFF PLACEMENT (CONTINUED)

Implant Pump

11. Evaluate the inferior aspect of the scrotal incision and develop a space underneath the scrotal skin and dartos muscle to serve as a pouch for the pump. Begin the development of the tunnel about 2 cm from the skin edge in order to facilitate eventual tubing and connector concealment. Loosely tie a purse string suture around the opening of the tunnel to secure the pump position. (Figure 4-22)

Figure 4-22.

Connect Tubing

Follow the instructions given in the section on connecting the tubing. Check the connection by activating the system (closing/opening the cuff).

Deactivate

The device system must be left in the deactivated mode for four to six weeks following implantation. Deactivate the system by following the instructions given in the section on “Deactivating the System”.

BULBOUS URETHRA CUFF PLACEMENT (CONTINUED)

ATTACHING A SECOND CUFF TO THE SYSTEM

Clinical literature reports that a small percentage of patients with severe stress incontinence may continue to report some degree of incontinence even after the placement of an artificial urinary sphincter. In such cases, a second cuff may be implanted in males around the bulbous urethra.

The second cuff is connected to the device system through the use of a 3-way connector included in the Accessory Kit. Once the device is ready to be cycled, it will take about twice as long for both cuffs to fill as it does for one cuff to fill from a single balloon.

If the patient is experiencing incontinence secondary to urethral atrophy at the site of the first cuff, the first cuff may also require replacement.

Table 1-1. Double Cuff Combinations

Cuff Sizes	3.5	4.0	4.5	5.0	5.5	6.0	6.5	7.0	7.5	8.0	9.0	10	11
3.5													
4.0													
4.5													
5.0													
5.5													
6.0													
6.5													
7.0													
7.5													
8.0													
9.0													
10													
11													

Indicates cuff combinations that may be used.

BULBOUS URETHRA CUFF PLACEMENT (CONTINUED)

Components and Accessories

- 2 cuffs
- 1 pressure regulating balloon
- 1 control pump
- 1 Accessory Kit
- Tubing passers, Quick Connect Tool (when necessary)

Choose the Correct Size of Cuff

Follow the usual instructions as discussed earlier for cuff selection by properly measuring the bulbous urethra with the Cuff Sizer.

Fill the Components

More filling solution must be added to the system to account for the additional cuff. A single cuff system calls for filling the system with 22 cc, aspirating the pressure regulating balloon and refilling it with 20 cc. To account for the second cuff, the system should be filled with 24 cc, the pressure regulating balloon aspirated, and refilled with 20 cc. When pressurizing the cuffs, wait 60 seconds for both cuffs to pressurize.

Place a Second Cuff

First, deactivate the artificial urinary sphincter and insert a Foley catheter to facilitate urethral palpation. Follow the steps discussed earlier in this Manual for urethral measurement and cuff selection. When placing the second cuff either proximal or distal to the first cuff, allow a 1 to 2 cm gap between the two cuffs to prevent them from rubbing against each other and to maintain vascularization.

If Original Cuff is being replaced

With both new cuffs in place and clamped, flush the tubing and secure the 3-way connector to the two cuff ends with 3-0 nonabsorbable polypropylene suture. Flush the 3-way connector. Connect the pressure regulating balloon, which has previously been filled with 24 cc filling solution. Remove the clamps from the cuff and balloon tubing. Allow 60 seconds for the cuffs to pressurize, and then re-clamp each cuff tubing below the 3-way connector.

Disconnect the balloon, aspirate the fluid from the balloon and refill it with 20 cc. When ready to connect the pump tubing, first flush the 3-way connector to remove air and particulate matter. Follow the previously described steps for connecting the tubing and tying with suture.


Figure 4-23.

BULBOUS URETHRA CUFF PLACEMENT (CONTINUED)

If Original Cuff is Preserved

Find the white/clear tubing leading from the existing cuff to the pump. Clamp the tubing on each side of the current connector. Cut the connector. With the new cuff in place, trim the tubing to the desired length, and flush the end. Connect a syringe filled with 10 cc of filling solution and a 15-gauge blunt tip needle attached to the new cuff tubing. Remove the clamp, and add 1-2cc (depending on the size of the cuff) of filling solution to the new cuff. Re-clamp the tubing and remove the syringe. Flush the tubing, attach the 3-way connector and pump tubing, and suture tie the connector to the pump tubing.

Connect Tubing

Follow the usual instructions for connecting the tubing when using the straight and right-angle connectors. The 3-way connector is suture tied to the cuffs and pump tubing. You may secure the 3-way connector to fascia or subcutaneous tissue to keep it stationary.

Deactivate

The device system must be left in the deactivated mode for four to six weeks following implantation. Deactivate the system by following the instructions given in the section on “Deactivating the System”.

BLADDER-NECK CUFF PLACEMENT






BLADDER-NECK CUFF PLACEMENT

Place a catheter into urethra to drain bladder and to help identify during dissection.

1. To begin the procedure, make a suprapubic incision and dissect around the bladder neck. (**Figure 4-24**)
2. To determine the size of cuff that is needed, measure the circumference of the urethra at the bladder neck with the Cuff Sizer. Place the Cuff Sizer around the urethra at the site where the cuff is to be implanted. It has to fit snugly without constricting the urethra. If a catheter has been placed in the urethra, it must be removed before measuring the bladder neck. (**Figure 4-25**)
3. Select the cuff size that corresponds to the measured length.
4. Prepare the selected cuff for implantation (see instructions for Preparing the Occlusive Cuff).
5. To implant the prepared cuff, pass the cuff, tab first, under the bladder neck. To avoid damage to the cuff, grasp it with a silicone-shod hemostat and grab the tab and pull.

Implant Pressure Regulating Balloon

6. Select the appropriate pressure-regulating balloon.
7. Use blunt dissection to create a space for the balloon in the prevesical space. Position the balloon in the prevesical space. (**Figure 4-26**)
8. Fill the balloon with 22cc of the appropriate filling solution and clamp the tubing 3cm from the end with a silicone-shod hemostat (one notch only).
9. Route the balloon tubing and the cuff tubing to the inguinal area. (**Figure 4-27**)
10. To pressurize the cuff, the cuff tubing and the balloon tubing are temporarily connected using a straight suture-tie connector. If a catheter has been placed in the urethra, make sure it is removed before pressurization. Remove the tubing clamps and wait 30 seconds for the cuff to pressurize.
11. Clamp the cuff tubing and the balloon tubing approximately 3cm from the end with silicone-shod hemostats and remove the connector.
12. Insert a syringe with a 15-gauge needle into the balloon tubing, remove the hemostat, and aspirate all of the remaining fluid from the balloon. Refill it with 20cc of filling solution, clamp the tubing with a silicone-shod hemostat, and remove the syringe.

BLADDER-NECK CUFF PLACEMENT (CONTINUED)

Implant Pump

13. To implant the control pump in the scrotum or labium, use blunt dissection to create a dependent subdartos pouch. The control pump should be placed on the same side where the pressure-regulating balloon was placed. **(Figure 4-28)**
14. Place the pump in the pouch, with the deactivation button facing outward so that it is palpable. Route the pump tubing to the inguinal area. **(Figure 4-29)**

Connect Tubing

Follow the instructions given in the section on connecting the tubing. Check the connection by activating the system (closing/opening the cuff).

Deactivate

The device system must be left in the deactivated mode for four to six weeks following implantation. Deactivate the system by following the instructions given in the section on “Deactivating the System”.

Transvaginal Approach

Some physicians prefer to implant the cuff in their female patients through a transvaginal incision. To begin, place the patient in a standard lithotomy position, prep, and drape. Make an inverted “U” incision in the anterior wall of the vagina. Start blunt dissection around the bladder neck. When dissection is complete, measure the circumference of the bladder neck with the Cuff Sizer to determine the size of cuff that is needed. If a catheter has been placed in the urethra, it must be removed before measuring the bladder neck. Select and prepare the appropriate size cuff and implant it around the bladder neck. If the cuff is being implanted trans-vaginally, a second small suprapubic incision must be made to place the pressure regulating balloon and control pump.


Figure 4-28.


Figure 4-29.

CONNECT TUBING

CONNECT TUBING

AMS Suture-Tie Connectors or AMS Quick Connect Sutureless Window Connectors are used to connect the tubing. In most cases, use the straight connectors. Right angle connectors should always be used when the tubing makes a sharp curve at the point of connection.

CAUTION: *AMS Quick Connect Sutureless Window Connectors should not be used in revision procedures involving previously implanted component tubing. Changes in the tubing over time may reduce the fixation performance of the Quick Connect Sutureless Window Connectors. The Quick Connect system may be used when all previously implanted components are removed and replaced with new components.*

Using AMS Quick Connect Sutureless Window Connectors

1. Using clean, well-cutting scissors, trim tubing lengths squarely to fit patient's anatomy. **(Figure 4-30)**
2. Slide collet ring onto tubing, making sure collet ring teeth face toward tubing end. **(Figure 4-31)**
3. Use 22-gauge needle and 10 cc syringe filled with filling solution to flush connector and tubing **(Figure 4-32)** to remove particulate matter and air.
4. Insert tubing end into connector. **(Figure 4-33)**
 - Firmly push one side of tubing to middle wall of connector.
 - Check placement through connector window.
 - Flush both connector and tubing before final connection.
5. Insert other tubing end into connector. **(Figure 4-33)**
 - Firmly push one side of tubing to middle wall of connector.
 - Check that both tubing ends touch middle.
6. Place straight connector ends into assembly tool jaws.


Figure 4-30. Trim Tubing


Figure 4-31. Slide on Collet


Figure 4-32. Flush Tubing and Connector


Figure 4-33. Insert Tubing Into Connector

CONNECT TUBING (CONTINUED)

7. Squeeze tool handles until closure stop touches opposite handle. (Figure 4-34)

CAUTION: Check tubing before you close assembly tool. Do not trap tubing between assembly tool jaw and connector. Tubing must exit straight from ends of connector, through slots in assembly tool.

CAUTION: After using assembly tool, tubing should bulge through connector window. This indicates that tubing is firmly against middle of connector wall.

8. If using right angle connector:
 - Use tool twice, once on each end. (Figure 4-35)

Note: Squeeze tool handles until closure stops touches opposite handle.

9. When using a 3-way connector, the assembly tool must be used 3 times - once on each end of the connector. The assembly tool must be engaged from the side of each connector end. Advance the tubing into each connector end until the tubing touches the inner connector wall as seen through the connector window. (Figure 4-36)

CAUTION: The long branch of the 3-way connector must be connected to the control pump.

10. After all connections have been made, cycle the device to confirm function, and deactivate the device (see instructions).


Figure 4-34. Closing Assembly Tool


Figure 4-35. Using a right angle connector


Figure 4-36. Quick Connect Assembly Tool (Approach from the side of 3-way connector)

Using AMS Suture-Tie Connectors

All connections using AMS Suture-Tie Connectors are tied with 3-0 nonabsorbable polypropylene.

1. Using clean, well-cutting scissors, trim tubing lengths squarely to fit patient's anatomy.
2. Use a 22-gauge blunt tip needle and 10cc syringe to flush tubing ends with filling solution to remove particulate matter and air before connecting.
3. Push connector ends over tubing so that they meet at connector's center hub; use 22-gauge blunt tip needle to flush connector before completing connection.
4. Use a double-throw overhand surgeon's knot followed by two single throws to attach the tubing to the connector.

Note: Be certain that the tubing ends are touching the center hub.

Note: Suture should crimp, but not cut the tubing.

CONNECT TUBING (CONTINUED)

5. Pass the suture to the opposite side of the connector, and use the same tying technique. (**Figure 4-37**)
6. After all connections have been made, cycle the device to confirm function, and deactivate the device (see instructions).

Figure 4-37. Suture-Tie Connector

CONNECT TUBING (CONTINUED)

DEACTIVATION KIT

The Deactivation Kit is useful during revision surgical procedures. The following example describes explanting only the cuff and preserving the other components.

When the cuff needs to be explanted due to erosion, first clamp the tubing with jaw-protected hemostats at both sides of the connector. Then, cut the tubing and explant the cuff. Unless infected, the pressure-regulating balloon and control pump may be left in the body with the tube plugs inserted into the tubing ends while the tissue around the urethra heals after the cuff is explanted.

By inserting a tubing plug into the end of the control pump tubing (with clear, reinforcement filament) the surgeon can protect the filling solution from contamination by blood or other materials during deactivation.

1. Flush the inside of the tubing end with a 22 gauge blunt needle, then insert a tubing plug from the Deactivation Kit.
2. Suture-tie the tube plug with 3-0 non-absorbable suture. When tying, use one double-throw, overhand surgical knot followed by at least two single throw knots to attach the connector to the tubing. The suture should clamp the tubing, but not so tight that it cuts the tubing.
3. The Tubing Plug connected to the control pump tubing is placed in the same superficial location as the original connectors and the incision is closed.
4. After healing takes place, when the new cuff is inserted, clamp the tubing beneath the plug with a tubing-shod clamp. Remove plug. Make sure the ends of the tubing have been cut squarely. Using a 22 gauge needle and a 10 cc syringe, flush end of tube with system fluid. Reconnect to component using a suture-tie connector. No additional filling or priming is necessary as all the original fluid is maintained within the system.

DEVICE OPERATION

DEACTIVATE CUFF

To deactivate the device, complete the following instructions:

1. Squeeze and release the pump bulb several times to remove all fluid from cuff. **(Figure 4-38)**

Note: The cuff will be empty when the pump remains flat.

2. Allow pump bulb to partially refill (approximately 30 seconds to one minute).

Note: It is recommended that time needed to fill pump or number of pumps needed to empty pump be recorded. This information is helpful postoperatively.

3. When slight indentation in pump bulb is felt, press the deactivation button **(Figure 4-39)**.

Note: It is important to leave a slight indentation in the pump bulb to ensure there is enough fluid in the pump to activate the device later. The urethra must be open, but the pump bulb should not be totally flattened.

Note: You will always feel the deactivation button, even when the device is deactivated. (This is not the type of button that bulges into the base).

4. After you press deactivation button, the pump bulb may feel more firm than usual.

Note: The indentation in the pump bulb will remain partially filled until activated. When the device is deactivated, the cuff will not be inflated (urethra not closed) and your patient will be incontinent. No fluid will move into cuff or pump when it's deactivated.

ACTIVATE (REACTIVATE) CUFF: NORMAL METHOD

To activate (reactivate) the device, complete the following instructions:

5. Push deactivation button a few times to loosen the poppet. **(Figure 4-40)**. Then give the pump bulb a quick, forceful squeeze.

*Note: This will move deactivation poppet back to activated position. After device is activated, pump will fill first and then cuff will refill **(Figure 4-41)**. It will take a few minutes for device to refill and for cuff to close off urethra or bladder neck. When system is activated, the pump may become less firm.*

Note: If you have difficulty activating device, there may not be


Figure 4-38. Squeeze and Release Pump Bulb


Figure 4-39. Press Deactivation Button When Slight Indentation Felt


Figure 4-40. Activate (Reactivate) the System (get out of deactivation status)


Figure 4-41. Cuff and Pump Refilled

DEVICE OPERATION (CONTINUED)

enough fluid remaining in pump to push deactivation button to its activated position. Use an optional method described below if this happens.

ACTIVATE (REACTIVATE) CUFF: OPTIONAL METHODS

If the normal activation method does not work, use one of the following optional methods.

Side Squeeze Method

1. Squeeze sides of control pump adjacent to deactivation button to allow fluid to fill the pump bulb. **(Figure 4-42)**

Note: It may take several minutes for the pump to refill. When enough fluid has returned to the pump bulb, give it a quick, forceful squeeze in order to reactivate the system.

Cotton Swab Method

1. Feel control pump to locate deactivation button.
2. Take a cotton tip swab and apply pressure to area directly behind deactivation button. **(Figure 4-43)**

Note: This should unseat the poppet and allow fluid to fill pump and then cuff.

Valve Block Bending Method

1. Feel the control pump, locate the deactivation button, and place your index finger above it (at tubing side). **(Figure 4-44)**
2. Place tip of your thumb below deactivation button on the opposite side.
3. Place index finger of your other hand on firm portion of pump (valve block portion) in front of deactivation button (toward the pump bulb).
4. Firmly bend pump end down to activate, by using thumbs as a fulcrum.
5. Release after bending.
6. Squeeze and release pump bulb several times to transfer fluid.


Figure 4-42. Side Squeeze Method


Figure 4-43. Cotton Swab Method


Figure 4-44. Fulcrum Method

POST OPERATIVE CARE

Some surgeons use a prophylactic antibiotic prior to surgery and intravenous antibiotics immediately postoperatively. Most choose to send their implanted patients home with a five to ten day course of antibiotics. The following paragraphs provide additional details on postoperative care.

IMMEDIATELY POSTOPERATIVE

After the surgery, deactivate the cuff and introduce a catheter into the urethra prior to closing. The length of time the catheter is left in place is at the discretion of the physician.

After 24 hours, the nursing staff may place ice packs in the region of the pump to reduce postoperative edema. At the discretion of the implanting physician antibiotics may be prescribed. The patient should be advised on the use of absorbent pads or condom catheters, until the device is activated four to six weeks after the surgery. The patient should be advised to avoid undue compression of the cuff area.

AFTER RELEASE FROM HOSPITAL

The patient is usually discharged within one to four days after surgery. After leaving the hospital, the patient should take antibiotics as prescribed by the physician.

The patient must return to the physician's office to activate the device prior to use of the AMS 800 Urinary Control System. The device is normally activated four to six weeks postoperatively. At that time, instruct the patient that it is possible to begin using the prosthesis to urinate.

Activating the device may be difficult if deactivation occurred when the pump bulb was deflated. If unable to cycle the prosthesis, squeezing the sides of the control pump adjacent to the deactivation button will allow fluid to fill the pump bulb and then the pump can be cycled normally. Refer to the activation/deactivation instructions.

The patient may experience some mild discomfort the first few times the prosthesis is being used. To determine that the patient is ready to use the device, check the incision site to be sure that it has healed properly.

There should be no redness, swelling, or drainage. Any of these things may indicate that an infection is present and the infection should be treated appropriately. Ask the patient about tenderness and/or discomfort when cycling the device.

The physician may want to observe the patient for up to an hour in the office to determine if sufficient continence is achieved with the device activated.

Provide patient education with regard to the operation. Patients should carry a permanent patient ID card to inform others of the device in case of emergency and **to prevent catheterization without device deactivation (which could damage the urethra or the device).**

EVALUATING LONG-TERM FUNCTION AND PLACEMENT

After the postoperative healing period, the surgeon should continue to have contact with the patient at least on an annual basis to evaluate the function of the device. During the annual evaluation, the surgeon should ask the patient about how the device is functioning and if he/she has noticed any changes in the function.

If the patient is having mechanical difficulty with the device, or there is infection or erosion present, revision or removal surgery may be necessary. In the event of a revision surgery, follow the same preparation and implantation techniques outlined elsewhere in this manual. Suture tie connectors must be used in any revision surgery where the total device is not explanted.

DOCUMENTATION

COMPLETING THE PATIENT INFORMATION FORM

American Medical Systems requires that a Patient Information Form be completed and returned for each implantation procedure. This is necessary to comply with the AMS Limited Warranty Policy and for the patient to receive an ID card.

The top portion of the form concerns the patient and the procedure. The primary etiology should be specified in detail, and the component information—cuff size, balloon pressure, and serial/lot numbers—should be as complete as possible, using self-adhesive labels included within the packaging.

Promptly return the first copy of the Patient Information Form to American Medical Systems. Retain the remaining copies for the files of the hospital, surgeon, or patient (if needed). Complete the Return Goods Form on the last page of the Patient Information Form if components were explanted.

INVENTORY RETURNS AND REPLACEMENT INFORMATION

A Patient Information Form (PIF) must be filled out and filed with American Medical Systems at the time of implant to activate the product warranty. Before returning any components, whether explanted or unused (sterile or nonsterile), customers must fill out the Return Goods Form located on the last page of the Patient Information Form.

Follow all of the instructions on the form carefully, and be sure that the components have been thoroughly cleaned before returning them to American Medical Systems. Request an AMS Product Return Kit from the AMS Customer Service Department to return any explanted components to American Medical Systems.

In all cases, obtaining credit or percentage of credit for a returned component is subject to approval under the terms of the AMS Return Goods Policy and the AMS Limited Warranty Policy. For complete information regarding these policies, contact the AMS Customer Service Department.

Contact your local AMS Representative prior to returning any product.

DOCUMENT INFORMATION

This document is written for professional medical audiences. Contact American Medical Systems for lay publications.

American Medical Systems periodically updates product literature. If you have questions about the currency of this information, contact American Medical Systems.

TROUBLESHOOTING

Problem

What to Do

Total Device

The device fails to cycle

Check connections between components. If they are correct, change the entire device

Leak in any of the components

Check if there is a leak by pump manipulation or ultrasound. If there is a leak, replace all the components (because body fluid has entered into the system.)

Occlusive cuff

Too tight or too loose around the urethra.

Remove improperly-sized cuff. Remeasure with the urethra-cuff sizer and implant proper size.

Punctured or damaged

Remove and replace with new cuff

Pressure-Regulating Balloon

Punctured during filling

Remove and replace with new pressure-regulating balloon.

Control Pump

Difficulty activating (reactivating) the device

Squeeze and release the sides of the control pump adjacent to the deactivation button to allow fluid to fill the pump bulb. When enough fluid has returned to the pump bulb, give it a quick, forceful squeeze (see the activation section for other methods).

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92116967-15

AMS 800™

Urinary Control System User's Guide for Male Patients

IN THIS GUIDE:

- 2 The AMS 800™ Urinary Control System
- 2 Using the AMS 800™ Urinary Control System
- 6 Special precautions to take with your AMS 800™ Urinary Control System before undergoing a medical procedure
- 6 Deactivating the AMS 800™ Urinary Control System
- 7 Activating the AMS 800™ Urinary Control System
- 9 When should I contact my doctor?
- 9 Where can I get more information about the AMS 800™ Urinary Control System?

The AMS 800™ Urinary Control System

The AMS 800™ Urinary Control System is intended for use in the treatment of male urinary incontinence. This device is not for patients who are determined by their doctor to be poor surgical candidates, have permanent blockage of the lower urinary tract or who have uncontrollable contraction of the bladder. Some AMS 800™ Urinary Control devices may contain an antibiotic (InhibiZone™).

The AMS 800™ Urinary Control System consists of three main components: an occlusive cuff, a control pump, and a pressure-regulating balloon (Figure 1). These are connected by tubing and are filled with a sterile solution. The cuff works much like a normal sphincter, gently squeezing the urethra (urinary canal) closed to keep urine in the bladder.


Figure 1

The device may not be suited for patients who are allergic to the antibiotics contained within the device (rifampin, minocycline or other tetracycline's). Patients with urinary tract infections, diabetes, spinal cord injuries, or skin infections may have an increased risk of infection.

Using the AMS 800™ Urinary Control System

While cuff is closed, urine stays inside of the bladder.

To urinate, the cuff is opened as follows:

1. Feel for the control pump in your scrotum (Figure 2).


Figure 2

2. Stabilize the control pump in the proper place by gently grasping the tubing above the control pump (Figure 3).


Figure 3

3. Use the other hand to squeeze and release the lower, softer part of the control pump (Figure 4). Do this several times until the control pump is flat.


Figure 4

4. When the lower part of the control pump is flat, it indicates that the fluid that filled the cuff has moved to the pressure-regulating balloon (PRB), deflating the cuff (Figure 5). The urethra opens and you can urinate.


Figure 5

That is the all the manipulation required to urinate. The fluid automatically returns from the PRB to refill the cuff. Once full, the cuff again squeezes the urethra closed. It usually takes several minutes for the cuff to be filled by the fluid (urethra to be re-closed), but each person is different. After using your AMS 800™ Urinary Control System a few times, you will know how long it takes for your cuff to refill.

If you feel that you did not empty your bladder completely before the cuff refilled, give the soft part of the control pump a few more squeezes to re-open the urethra.

It is a good idea to get into the habit of emptying your bladder on a regular basis, such as every two to three hours.

It is important to let someone else, such as your spouse, know how to operate the AMS 800™ Urinary Control System to assist you if needed.

Special precautions to take with your AMS 800™ Urinary Control System before undergoing a medical procedure

It is important for healthcare providers to know that you have an AMS 800™ Urinary Control System before you undergo any medical procedure. This is especially true of any procedures which require the insertion of a catheter or any other instrument into the urethra. Such procedures can only be safely performed after your AMS 800™ Urinary Control System has been deactivated (see the following section for deactivation instructions). Failure to deactivate the AMS 800™ Urinary Control System before inserting a catheter or other instrument into the urethra can result in damage to the device and/or your urethra. It is possible that such damage could require removal of the device.

Deactivating the AMS 800™ Urinary Control System

The AMS 800™ Urinary Control System has a deactivation button which stops the flow of the fluid from one component to another. Your doctor may deactivate your device at various times, such as while you are healing following your surgery or before inserting a catheter or other instrument into your urethra.

To deactivate the AMS 800™ Urinary Control System:

1. Squeeze and release the lower part of the control pump several times, just as you would do to urinate.

2. After completely emptying the cuff (the cuff is empty when the pump remains flat), allow the pump bulb to partially refill (about 30-60 seconds). When the control pump is partially refilled, firmly press the deactivation button on the top of the control pump. After pressing the deactivation button, you will still feel the button on the upper part of the pump. Pressing the deactivation button locks out the flow of fluid through the control pump. When the device is properly deactivated, a slight indentation will be felt at the center of the soft part of the pump. The cuff will now remain open (Figure 6).


Figure 6

When the AMS 800™ Urinary Control System is deactivated, the cuff around your urethra will not stop the flow of urine. You should use pads while the AMS 800™ Urinary Control System is deactivated.

What happens if the deactivation button is squeezed accidentally? The device will deactivate and two problems could occur.

Fluid may not refill the cuff, so you will keep leaking urine, or fluid may lock in the cuff so you cannot urinate. If either problem occurs, you will notice that the pump feels different from its normal status. You can activate the device to solve these problems.

Activating the AMS 800™ Urinary Control System

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It is important that the AMS 800™ Urinary Control System remain open while you are healing following your surgery. Closely follow your doctor's instructions about the amount of time before activating the device.

To activate the AMS 800™ Urinary Control System:

3. Use one hand to grasp the tubing above the control pump to hold the control pump in an appropriate place (Figure 7).
4. Use the other hand to squeeze the sides of the upper, hard part of the control pump, on the side opposite the deactivation button. Be sure to use steady, firm pressure for several minutes. This will cause the lower part of the control pump to become hard and filled with fluid.


Figure 7

5. Continue stabilizing the control pump by holding the tubing (DO NOT STABILIZE THE PUMP BY

HOLDING THE UPPER, HARDER PART OF THE PUMP DURING THIS STEP. DOING SO WILL MAKE PROPER ACTIVATION EXTREMELY DIFFICULT). Give the lower part of the pump a quick, forceful squeeze (Figure 8). This step may need to be repeated if you're not successful the first time. Because greater pressure is required for activation than for normal operation, you may experience some discomfort while activating the device. The system should now function normally.


Figure 8

When should I contact my doctor?

- ▶ Fever
- ▶ Swelling near the incision
- ▶ Redness near the incision
- ▶ Any other signs of infection
- ▶ If the surgical area does not heal
- ▶ Prolonged pain or discomfort
- ▶ If the device can be seen or felt outside of the skin
- ▶ Return of leakage
- ▶ If the device is not working correctly

Where can I get more information about the AMS 800™ Urinary Control System?

Your doctor is the best source for getting more information on the AMS 800™ Urinary Control System. Your doctor is most familiar with your medical situation and can work with you to make sure you get the best possible results from the AMS 800™ Urinary Control System. Contact your doctor if you have questions or would like to receive additional information.

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You can also call the AMS patient liaison at 1-800-328-3881, extension 6261, or visit our website at amsmenshealth.com.

American Medical Systems, Inc.
10700 Bren Road West
Minnetonka, MN 55343
U.S.A.
US toll-free: 1 800 328 3881
Tel: + 1 952 930 6000
Tel: + 31 20 593 8800

American Medical Systems
Europe B.V.
Haarlerbergweg 23 G
1101 CH Amsterdam Zuid-Oost
The Netherlands

©2016 Boston Scientific Corporation
or its affiliates. All rights reserved.
All trademarks are the property of
the respective owners.
1001691 Rev B (2016-07)

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Sutureless Window Connectors

Instructions for Use

English **AMS™ Quick Connect**
Sutureless Window Connectors
Instructions for Use..... 3

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

Sutureless Window Connectors

Instructions for Use


AMS Quick Connect Sutureless Window Connector Assembly Instructions


1. Trim tubing squarely.


2. Attach collet rings.
Place one ring on end of each tubing.


3. Flush.
Flush ends of both tubing using appropriate solution.


4. Insert first tubing.
Firmly push tubing to middle wall of connector as seen through the connector window.


5. Connect.
Firmly push second tubing to middle wall. Check windows to make certain both tubing ends are still touching the middle wall.


6. Approximate the collets against the connector prior to applying the assembly tool.

7. Do not hold the connection in the assembly tool. This may contribute to the collet misaligning with the connector and result in collet prong(s) missing the connection.


8a. Assure that the connection is inserted completely into the assembly tool. Slowly squeeze the handle, closing the tool. If a prong of the collet does not enter into the connector prior to completing the connection, release the pressure on the connection and allow the collet to realign with the connector. Complete the connection. The closure stop must touch the opposite handle.


8b. When using a right angle connector, the assembly tool must be used twice—once on each end of the connector. Again, make certain the tubing is touching the middle wall on both sides of the connector. Closure stop of assembly tool must touch opposite handle each time.

CAUTION: Before you close the assembly tool, check the tubing. The tubing must not be trapped between the assembly tool jaw and the connector. The tubes must exit straight from the ends of the connector, through the slots in the assembly tool.


8c. When using a 3-way connector, the assembly tool must be used 3 times—once on each end of the connector. The assembly tool must be engaged *from the side* of each connector end. Advance the tubing into each connector end until the tubing touches the inner connector wall as seen through the connector window.

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
9. Check.

After completing the connection, examine both ends of the connection thoroughly. If prong(s) of the collet do not enter the connector and are protruding, do not close. Remake the connection by cutting out the connection and making a completely new connection.

Incorrect

Perforates tubing.

Correct

Tubes exit straight from the ends of the connector. After using the AMS Quick Connect™ Assembly Tool, the tubing should partially bulge through the window of the connector. This indicates that the tubing is still firmly against the middle wall of the connector.

The AMS Quick Connect System cannot be resterilized. Conventional hospital sterilization will damage the connector components. (The assembly tool can be resterilized.)

NOTE: AMS Quick Connect Sutureless Window Connectors should not be used on revision tubing.

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
Spectra™*Concealable Penile Prosthesis***Instructions for Use****INDICATIONS FOR USE**

The Spectra™ Concealable Penile Prosthesis is a sterile, non-pyrogenic, single-use implant that is intended for use in the treatment of chronic, organic, erectile dysfunction (impotence) in men who are determined to be suitable candidates for implantation surgery.

DEVICE DESCRIPTION

Spectra prostheses are provided as sterile, non-pyrogenic, single-use implants. Do not resterilize.

The Spectra device consists of a pair of cylinders which are surgically implanted, one into each corpus cavernosum, to provide penile rigidity. Each device consists of two cylinders (see Figure 1), and may include rear tip extenders (RTEs) for additional length in some cases. All components consist of implantable, biocompatible materials.

Figure 1: Spectra Cylinder

Each cylinder consists of an articulating section of alternating polyethylene and titanium segments. A cable extends through the center of the articulating segments. The proximal and distal ends of the cable are each connected to a stainless steel spring that is encased in a titanium housing. The proximal rigid section of the 9.5 mm diameter cylinders is 5 cm long. For the 12 mm and 14 mm diameter cylinders this section is 6 cm long. The entire outside surface of the cylinder is made of silicone.

The articulating segments, under compression by the cable and spring assemblies, provide sufficient friction and rigidity for the patient to position the device for concealment or for intercourse.

Spectra cylinders are available in 9.5, 12, and 14 mm diameters. Each cylinder diameter is available in three lengths: 12, 16, and 20 cm.

The total length can be adjusted by adding rear tip extenders (RTEs), comprised of silicone with barium sulfate for radiopacity, to the proximal end of each cylinder. A range of RTE lengths are included to accommodate the patient's total intracorporal length (Figure 2a). The addition of RTEs can extend the cylinder length from 0.5 to 7.5 cm (Figure 2b).

SPECTRA DIMENSIONAL SPECIFICATIONS

Cylinder Diameter	x	Cylinder Length	Cylinder Product Number
9.5 mm	x	12 cm	720054-01
9.5 mm	x	16 cm	720054-02
9.5 mm	x	20 cm	720054-03
12 mm	x	12 cm	720074-01
12 mm	x	16 cm	720074-02
12 mm	x	20 cm	720074-03
14 mm	x	12 cm	720056-01
14 mm	x	16 cm	720056-02
14 mm	x	20 cm	720056-03

Cylinder Diameter	RTE Lengths	RTE Product Number
9.5 mm	0.5, 1, 1.5, 2, 3 cm (kit)	720170-01
9.5 mm	4, 5, 6 cm (kit)	720171-01
12 mm & 14 mm	0.5 cm	72404320
12 mm & 14 mm	1 cm	72404321
12 mm & 14 mm	2 cm	72404322
12 mm & 14 mm	3 cm	72404323
12 mm & 14 mm	4 cm	72404324
12 mm & 14 mm	5 cm	72404325
12 mm & 14 mm	6 cm	72404326
12 mm & 14 mm	1.5 cm stackable	72404330



Figure 2b: Available RTE Stack Lengths for Cylinder Attachment

Note: These RTEs are non-stackable, with the exception of the 1.5 cm stackable RTE (shaded in the figure).

CONTRAINDICATIONS

- The implantation of this device is contraindicated in patients who have active urogenital infections or active skin infections in the region of surgery.
- Implantation of this device is contraindicated in patients whose proximal corporal length measurement is less than the proximal rigid section of the Spectra cylinders (5 cm for the 9.5 mm diameter, and 6 cm for the 12 mm and 14 mm diameter).
- Implantation of this device is contraindicated in patients whose total intracorporal length is not within the range of 12 cm to 27.5 cm.
- The implantation of this device is contraindicated in patients whom the physician determines to be poor candidates due to risks associated with open surgical procedures and /or the patient's medical history (physical and mental conditions), or with sensitivity to silicone materials.
- Implantation of the device is contraindicated in patients who require repeated endoscopic procedures.
- Patients who have compromised tissue and as a result cannot withstand constant pressure should not be implanted with a Spectra penile prosthesis.

WARNINGS AND RISKS

The Spectra prosthesis is designed to be implanted as a pair of matched cylinders. A single implanted cylinder may not be adequate to achieve sexual intercourse and may have an adverse effect on the reliability of the device.

Known and potential complications include, but are not limited to infection, erosion, migration/malposition, extrusion, mechanical malfunction, patient dissatisfaction, adverse tissue reaction, allergic reaction, prolonged or intractable pain, neuropraxia, urinary obstruction, urinary retention, silicone particle migration, and other complications (e.g., post-operative bleeding, hematoma, edema, penile necrosis/gangrene, perforation of the corpora or urethra or crus, cold glans, inability to adequately dilate the corpora, incorrect sizing of the implant, urethrocorporeal fistula, SST deformity, and tearing or ripping of the device during or after implantation).

The complications listed above may necessitate surgical revision or removal of the prosthesis.

Product Life

The Spectra Concealable Penile Prosthesis is intended as a prosthetic device which restores to the patient an important physiological function. As with any penile rigidity implant, this device is subject to wear and eventual failure over time.

It is not possible to predict how long an implanted prosthesis will function in a particular patient.

Patients should be advised that the implant is not considered a lifetime prosthesis and additional surgery for replacement and removal may be necessary.

Currently, there is no data on the revision rate for this device. However, in a study of the Dura II penile prosthesis (similar to this device), the revision rate was reported as 5.1% after two years post-implant for 196 patients. There were no mechanical failures reported.¹

Mechanical Malfunction

Product wear, mechanical, or performance malfunctions may occur that could render the prosthesis inadequate for sexual intercourse or concealment. Malfunctions may include fracture of the prosthesis, difficulty or changes in the ability to position the prosthesis, changes in rigidity or column strength of the prosthesis, or breach of the outer silicone layer which may lead to surgical intervention or removal of the prosthesis.

Mechanical events should be evaluated carefully by the treating physician and the benefits and risks of treatment options, including revision surgery, should be considered.

Infection

As with any surgical implantation, the surgery required to implant this device can result in infection. Men with existing medical conditions such as diabetes, spinal cord injuries, open sores or skin infections in the region of the surgery, or urinary tract infections may increase the risk of prosthetic-associated infections.

Appropriate measures should be taken to reduce the likelihood of infection, such as use of sterile techniques and appropriate antibiotic prophylaxis. The patient should be monitored for infection and treated appropriately.

Infection that fails to respond to antibiotic therapy may result in removal of the prosthesis. Implantation of a new device may be contraindicated at the time of removal for infection.

Infection followed by explantation of the device may result in scarring, which may make subsequent re-implantation more difficult.

Erosion

Erosion, which is the breakdown/disruption of tissue planes adjacent to the device, can occur with the device. Erosion may be caused by infection, pressure, improper sizing, tissue

Регистрация МИ в Росздравнадзоре J Peretsman, JO Bolton, RG Holcomb, PK Reddy, PH Bernhard, SM Eppel, www.nevacert.ru, info@nevacert.ru, AA Melman. "Report of a Multicenter Clinical Evaluation of the Dura-II Penile Prosthesis." J. Urol. 155: 1613-1616, 1996.

damage, and cylinder misplacement. The most frequently reported sites associated with cylinder erosion are the glans, the urethra, or the skin. In any case of erosion, the physician must evaluate and decide whether repair and/or removal is necessary.

Failure to timely evaluate and treat the erosion may result in a substantial worsening of the condition leading to infection and loss of tissue.

Migration

Migration is the movement or displacement of components within the space in which they were implanted and can result in surgical revision, pain, psychological/medical complications (such as floppy glans), or device malfunction.

Migration can occur if the cylinders are improperly sized or if there is malpositioning of the cylinders within the corpora cavernosa.

Silicone

This device is composed of a number of materials, including solid silicone elastomers. Silicone gel is not a component in the materials of this device.

Silicone elastomers have been commonly used in a variety of biomedical devices for over 40 years and are used as a biocompatibility reference against which new materials are tested.

Scientific literature has included reports of adverse events and other observations in patients with implantable silicone devices. As reported, these events/observations indicate “allergic-like” symptoms and in other cases a symptom complex associated with immunological disorders. No causal relationship has been established between these events and silicone elastomers.

There are reports of malignant tumor formation in laboratory animals only, not humans, associated with implants of relatively large size. Many different materials are associated with this effect in animals, silicone elastomers are among them. No such effect has been described in humans.

Extensive testing has been conducted on all materials that comprise the AMS penile prostheses product lines. This testing has indicated no toxicological response attributable to the materials. However, some of the materials caused minor irritation when implanted in animals.

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The risk of surgical revision or device removal is common to all devices. Generally, surgical revision or removal of penile prostheses is performed to address other complications. However, the patient may elect to have the device removed due to patient dissatisfaction unrelated to safety or efficacy. Removal of an implanted prosthesis, for any reason, without timely re-implantation of a new prosthesis may substantially complicate re-implantation or make it inappropriate.

The patient must be informed that a penile prosthesis is subject to wear and that eventual failure is expected over time. It is not considered a lifetime prosthesis. The patient should be made fully aware that additional surgery for replacement and removal may be necessary.

Patient Expectations

Patients should be counseled in order to have a realistic expectation of the physical, psychological, and functional outcome of the implantation of a rigid penile prosthesis.

Implantation of a penile prosthesis may result in penile shortening, curvature, or scarring. The prosthetic erection may differ from the patient's original, natural erection in that it may be shorter, less firm, have less girth, and reduced sensations.

Realistic cosmetic expectations should be communicated to the patient and should include the potential for skin scarring and lack of concealability.

The implantation of a penile prosthesis will not provide rigidity to the glans, and may result in a floppy glans and lack of rigidity of the corpus spongiosum. Penile flaccidity will be less than prior to implantation.

Loss of Latent Erectile Capability

Penile prosthesis implantation may result in the loss of any natural or latent spontaneous erectile capacity.

Device Sizing

Proper sizing of the device is critical to a successful outcome. Improper measurement, inappropriate cylinder size selection, as well as malpositioning of the cylinders within the corpora cavernosa may result in migration or buckling of the cylinders or reduce cylinder life.

Pain

Device implantation may result in pain at the operative sites during the post-implantation period and during periods of initial use. Cases have been reported of chronic pain associated with device implantation.

Pain with a severity or duration beyond that which is expected can be symptomatic of medical

complications or mechanical device malfunction, which may lead to medical or surgical intervention. There are also reports of patients without known medical complications who elected to have a functional device removed due to unresolved pain.

Patients should be counseled on expected post-operative course of pain, including severity and duration, to get a sense of the normal healing process.

Surgical Technique

Unsuccessful outcomes have been reported due to improper surgical technique which compromises the integrity of the device (including cuts or abrasions of the device), anatomical misplacement of the device, or improper sizing of cylinders.

PRECAUTIONS

Surgery Related

Direct contact of surgical instruments to the prosthesis may result in damage, rendering it unsuitable for implantation.

During insertion, **do not** over bend cylinders beyond their natural U-shape as it may damage the prosthesis and shorten its product life.

Do not trim the distal or proximal ends of the cylinders, or the rear tip extenders (RTEs). Trimming will damage the device.

Careful intraoperative sizing is required to ensure proper device operation and to minimize the occurrence of sizing-related complications such as migration and/or extrusion.

Device Related

Implantation of a penile prosthesis that has been in previous contact with or contaminated by body tissue or fluid, regardless of intervening, cleaning, or sterilization, is prohibited.

Spectra RTE Kits, 720170-01 and 720171-01, are **only** compatible with the 9.5 mm diameter Spectra cylinders.

Do not use them with the 12 mm or 14 mm Spectra cylinders. (See "Spectra Dimensional Specifications" section for details.)

Spectra RTEs, 72404320 through 72404326 and 72404330, are **only** compatible with the 12 mm and 14 mm diameter Spectra cylinders. **Do not** use these RTEs with the 9.5 mm Spectra cylinders. (See "Spectra Dimensional Specifications" section for details.)

Do not stack the Spectra RTEs on each other, with the exception of the 1.5 cm RTE. The internal locking ring of the 1.5 cm RTE must not engage with the smooth outer surface of

other size RTEs. Improper attachment will result in the RTE disconnecting and/or bulging.

The articulating segments of the Spectra penile implant may produce vibration and/or sound during movement. This is normal and may be detected physically or audibly by some patients.

The total length of RTEs and the proximal rigid section of the cylinders should not extend beyond the penoscrotal junction.

The AMS Sizers are provided non-sterile.

Patient Related

A thorough preoperative consultation between patient and physician should include a discussion of all available treatment options and their risks and benefits.

Adequate patient education is required to ensure safe, effective patient use.

Uncircumcised patients may have an increased risk of post-operative complications with the sub-coronal approach. Surgeons may wish to discuss performing a circumcision to reduce the risks of post-operative complications associated with this approach.

Patients must have sufficient strength and/or dexterity to position the device.

MAGNETIC RESONANCE IMAGING (MRI) IMPORTANT SAFETY INFORMATION

MR Conditional

Non-clinical testing has demonstrated the Spectra product line is MR Conditional. The device can be scanned safely under the following conditions:

- Static magnetic field of 3.0 Tesla (3.0 T) or less
- Maximum spatial gradient magnetic field of 720 Gauss/cm or less

MRI-Related Heating

Non-clinical testing has demonstrated the Spectra product line produced the temperature rises during MRI performed for 15 minutes of scanning in the respective MR systems which would not pose a hazard to the human subject.

Static Magnetic Field	3.0 Tesla ^a
Highest Temperature Change	≤ +1.6 °C
Maximum whole body averaged Specific Absorption Rate (SAR)	2.9 W/kg for 15 minutes of scanning as assessed by calorimetry
(a) 3.0T 128MHz, General Electric Healthcare, Excite software version G3.0-052B	

Artifact Information

Non-clinical testing has demonstrated that the Spectra product line may compromise the MR image quality if the area of interest is relatively close to the position of the implant. The maximum image artifact produced by a MR gradient echo pulse sequence was a “moderate” localized signal void in size and shape of the implant. Optimization of MR imaging parameters to compensate for the presence of the device may be necessary.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	7,193 mm ²	1,553 mm ²	1,160 mm ²	7,030 mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

PACKAGING

The Spectra penile prosthesis is delivered sterile inside a sealed double barrier tray, which is packaged within a dustcover box. **Do not use the Spectra device if packaging integrity is compromised prior to opening or if it is beyond the expiration (use-by) date.**

At the end of the box, an inventory label lists the cylinder diameter and length, lot and serial numbers, product number, expiration date, and manufacturing date. The included self-adhesive labels can be used to record appropriate information on patient charts and forms.

Spectra cylinders are provided sterile, and they are non-pyrogenic. Do not resterilize. For single-use only.

PRE-OPERATIVE CONSIDERATIONS

The implanting surgeon/urologist should provide the implant candidate with the Spectra Patient Manual prior to surgery to allow the patient sufficient time to review the device information, the procedure, and to discuss this information with his physician and sexual partner.

Selection and education are important prior to implantation of an Spectra Concealable Penile Prosthesis.

Proper training of the surgical team on the following pre-operative considerations will minimize anticipated complications for the patient.

CAUTION: This device is to be implanted only by physicians who are knowledgeable and have received appropriate training regarding the use of penile prostheses. This IFU is not intended to be a complete reference.

Considerations before entering the operating room:

- To reduce the likelihood of infection, follow the use of sterile technique and antibiotic prophylaxis to prepare the patient, the operating room, and the implant.
- The patient's urine should be sterile.
- An antimicrobial shower should be given the night before surgery.
- The patient should receive antibacterial prophylaxis.
- The surgeon should perform a 10-minute hand scrub.

In the operating room:

- Parenteral antibiotics should be used liberally (at the physician's discretion).
- The patient's abdominal and genital areas should be shaved.
- The skin should be prepped with a 10-minute scrub of povidone-iodine soap or the approved hospital pre-operative scrub procedure.

Instruments

Assemble the typical surgical instruments required to perform a subcoronal or penoscrotal penile prosthesis implantation. In addition, a penile prosthesis sizer set, such as the AMS Sizers, and also a set of corporal dilators are required. A retractor may be useful with the surgical implantation.

INTRAOPERATIVE PROCEDURES

Surgical Approaches, Dissection, and Corporotomies

Proper training of the surgical team on the following intraoperative procedures will minimize anticipated complications for the patient.

There are two surgical approaches which could be used to implant the Spectra:

1. Penoscrotal: With this approach, corporotomies 5 to 7 cm long must be made to allow for prosthesis insertion. If difficulties are encountered, corporotomies should be extended to facilitate the insertion of the prosthesis.

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With this approach, corporotomies 2 to 3 cm long must be made to allow for the distal tips of the

prosthesis to be inserted under the glans.

For either approach, dissect through Buck's fascia to expose the tunica albuginea and make appropriate corporotomies in each corpus cavernosum.

1. ***Corporal Dilation/Prosthesis Diameter Selection***
- The Spectra cylinders are available in 3 diameter sizes. Using a series of progressively larger dilators, dilate both corpora proximally and distally to accommodate the prosthesis diameter (9.5, 12, or 14 mm). Dilation of the right and left corpora should be equivalent in diameter and length for optimal device fit. It is recommended to dilate approximately 1 mm beyond the diameter of the device to be implanted.

Note: Use care to avoid crossover through the intracavernosal septum during corporal dilation.

After dilating the corpora, select two dilators whose total diameter equals the total diameter of the cylinders to be implanted. Simultaneously insert the dilators side by side into the proximal ends of the corpora cavernosa to evaluate the overall fit. Repeat this step for the distal ends of the corpora.

Example: If two 12 mm diameter cylinders are to be implanted (for a total diameter of 24 mm), insert an 11 mm dilator into one corpus and a 13 mm dilator into the other corpus.

2. ***Total Intracorporal Length Measurement*** - Note the length of the proximal rigid section of the two Spectra cylinders.

Spectra Cylinder Diameter	Proximal Rigid Section
9.5 mm	5 cm
12 mm	6 cm
14 mm	6 cm

Stretch the penis to approximate an erection. Use the AMS Sizers to measure the proximal and distal corporal lengths and add them together to attain the total intracorporal length. Record these measurements on the Patient Information Form (PIF) supplied with the Spectra device.

AMS Sizers Product Information

(The AMS sizers may be used to measure the intracorporal lengths. Both sizers are packaged together in a non-sterile kit.)

- 1 AMS Sizer 9.5 mm (for use with 9.5 mm diameter cylinder)
 - 1 AMS Sizer 12 mm (for use with 12 mm and 14 mm diameter cylinders)
3. **Length Adjustment** - The total intracorporal length should be used to select the appropriate cylinder and/or RTE(s). To lengthen a pair of cylinders, select the appropriate combination of RTEs (up to two per cylinder) and attach them to the proximal end of the cylinders to equal the extra length desired (see Figures 2a and 2b)

When attaching two RTEs per cylinder, one of them must be the 1.5 cm RTE. The 1.5 cm RTE must be attached to the cylinder's proximal end first to allow for stacking of the second Snap-fit RTE. This stacking order ensures that the internal locking mechanism of the RTEs is engaged and their attachment is secure. Rotate the RTE on the cylinder (or on the stackable RTE) to verify proper attachment. Properly attached RTEs will rotate without disengagement or material bulging. **Ensure that the total length of RTEs and the proximal rigid section of the cylinder does not extend beyond the penoscrotal junction.** This is recommended to maximize the Spectra cylinder length and minimize the RTE length (for best concealment).

4. **Insert Prosthesis** - When using the subcoronal approach, first insert the proximal end of the prosthesis cylinder into the proximal corpus. Then insert the distal tip of the prosthesis under the glans by lifting and retracting the corpus distally (a retractor may facilitate this). At the same time, press the distal tip of the prosthesis into the corpus, under the glans. Repeat insertion steps for the other corpus.

Note: If you cross over through the intracavernosal septum to the contralateral side, remove and place the dilator into the contralateral side and reposition the cylinder on the ipsilateral side.

When the penoscrotal approach is used, corporotomies must be of sufficient lengths (5 to 7 cm long) to minimize overbending or disengagement of RTEs during cylinder insertion.

5. **Intraoperative Testing** - To determine device fit and function, perform the following up/down rigidity test after the device has been inserted: bend the penis down for the concealed position and straighten for the erect position. Then, perform the following test to ensure there is no buckling:

With the penis in the erect position, press against the palm of the hand to confirm correct rigidity. To confirm concealment, position a finger at the

dorsal or ventral base of the penis and curve the penis into the desired position. If the device fit or function is not appropriate, re-evaluate cylinder and/or RTE size selections.

6. **Corporotomy Closure** - Close the incision using an acceptable surgical technique.

POST-OPERATIVE CONSIDERATIONS

Proper training of the surgical team on the following post-operative considerations will minimize anticipated complications for the patient.

- Minimize catheter use.
- Administer routine wound care.
- At four to six weeks postoperative, instruct the patient that he may resume sexual activity after performing an evaluation of:
 - anatomical fit/placement
 - proper device functionality
 - healed incision site
- Educate the patient on the following:
 - how to operate his new penile prosthesis
 - the signs of infection
 - the proper way to urinate during the healing period
 - when he can resume sexual activity
 - the pain that may occur during the postoperative period and during the first few months of use

The patient may refer to the Spectra Patient Manual for further information.

- Recommend the appropriate antibiotic prophylaxis for the patient during the post implant period and during any subsequent surgical or dental procedures.
- Complete the patient ID card and give it to the patient after surgery.

After the postoperative healing period, the physician should continue to meet with the patient on an annual basis (at minimum) to evaluate the function of the device. During the annual evaluation, ask the patient about how the device is functioning and if he has noticed any changes in the function. Also, check the patient for signs of infection, extrusion or erosion.

If the patient is having mechanical difficulty with the device, or if there is infection, erosion or extrusion present, revision surgery may be necessary.

IMPLANT REGISTRATION

American Medical Systems, Inc. (AMS) requires implantation of a penile prosthesis to be registered at the time of implantation. In the United States, complete and return the AMS Patient Information Form (PIF) supplied with each pair of Spectra cylinders. Outside the United States, return the appropriate information to the distributor who sold the prosthesis.

The purpose of registration of an implanted prosthesis is to have a record in the event of a complaint, or to identify and protect the patient in the event of a device recall or warranty claim.

RETURNING INVENTORY

In the United States

Before returning any components, whether explanted or unused (sterile or unsterile), customers must complete the Return Goods Form located on the last page of the Patient Information Form. Follow all of the instructions on the form carefully, and be sure that the components have been thoroughly cleaned before returning them to AMS.

In all cases, obtaining credit or percentage of credit for a returned component is subject to approval under the terms of the AMS Return Goods Policy and the AMS Product Replacement Policy. For complete information regarding these policies, contact the AMS Customer Service Department.

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This document is written for professional medical audiences. Contact AMS for lay publications.

AMS periodically updates product literature. If you have questions about the currency of this information, contact AMS.

HOW SUPPLIED AND STORAGE

WARNING: Contents supplied STERILE. Do not use if sterile barrier is damaged. If damage is found, or if labeling information appears missing, illegible, or compromised, call your AMS representative for clarification, and do not use the Spectra device until integrity of the product has been verified.

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury,

illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

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Do not use device beyond the indicated use by date.

There are no special transportation conditions.

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All rights reserved. All trademarks are property of the
©2017 Boston Scientific Corporation or its affiliates. All rights reserved. All trademarks are the property of the respective owners.
All rights reserved. All trademarks are property of
©2017 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific (Australia) Pty Ltd
**American Medical Systems****Europe B.V.**