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CAPIOX[®] NX Oxygenator with Integrated Arterial Filter and UltraPrime[™] Technology

Applicable for Vacuum Assisted Venous Drainage & Post Operative Chest Drainage

Instructions for Use P.3-P.10

ENGLISH

Read all warnings, cautions and instructions carefully prior to use.

DESCRIPTION OF PRODUCTS

The CAPIOX[®] NX is a membrane oxygenator with microporous polypropylene hollow fibers, to be used as an extracorporeal gas exchange device in which the blood flows on the outside of the fibers and with the ventilating gas flowing through the fibers.

The CAPIOX[®] NX consists of a module for gas exchange with an integrated heat exchanger and Arterial Filter. The CAPIOX[®] NX is available with or without a Hardshell Reservoir with an integrated cardiotomy filter to provide a single-unit construction for ease of use.

The NX19 oxygenator is available with a 4000 mL CAPIOX[®] Hardshell reservoir.

The blood contacting surfaces are coated with X_{coating}[™].

X_{coating}[™] is bio-compatible material that is applied on the blood-contacting surfaces of the device to reduce the adhesion of platelets to the device.

INDICATION FOR USE

The CAPIOX[®] NX19 Oxygenator with Integrated Arterial Filter and UltraPrime[™] Technology is intended to be used to exchange gases between blood and a gaseous environment to satisfy the gas exchange needs of adult and small adult patients during cardiopulmonary bypass surgery.

The integrated arterial filter is intended to filter non-biologic particles and emboli and to facilitate air bubble removal from the blood flowing through the cardiopulmonary bypass circuit. The integrated heat exchanger is used to warm or cool blood and/or perfusion fluid as it flows through the device.

The hardshell reservoir is used to store blood during extra-corporeal circulation from the venous line and the cardiotomy line.

The reservoir contains a venous section that is comprised of a filter and defoamer to facilitate air bubble removal.

The cardiotomy section of the reservoir contains a filter to remove particulate matter and a defoamer to facilitate air bubble removal. The 4-liter reservoir may be used for Vacuum Assisted Drainage procedures and Post-Operative Chest Drainage Procedures.

The CAPIOX[®] NX19 Oxygenator with Integrated Arterial Filter and UltraPrime[™] Technology is for use with patients when the required blood flow rate will not exceed 8.0 L/min.

The CAPIOX[®] NX19 Oxygenator with Integrated Arterial Filter and UltraPrime[™] Technology can be used in procedures lasting up to 6 hours.

PRESSURE AND FLOW LIMITATIONS


Oxygenator Module Blood Flow/Pressure	Min. 2.0 L/min (or 0.5 L/min for up to 2 hours)	Max. 8.0 L/min	Max. 1,000 mmHg (133 kPa) (1.36 kgf/cm ²)
Reservoir Blood Flow/Negative Pressure	Min. 2.0 L/min (or 0.5 L/min for up to 2 hours)	Max. 8.0 L/min (5.0 L/min Cardiotomy)	Max. -150 mmHg (-20 kPa)
Water Pressure	Max. 2 kgf/cm ² (1,470 mmHg) (28 PSI) (196 kPa)		
Gas Flow	Min. 0.5 L/min Max. 20 L/min		

DIRECTION OF BLOOD OUTLET PORT

The oxygenator is available with TWO OUTLET PORT CONFIGURATIONS - LEFT side (W) or RIGHT side (E). Refer to the diagram below to determine which configuration is needed for your circuit.

Type	Left	Right
Oxygenator with Hardshell Reservoir	[REF] : 3CX*NX19RW	[REF] : 3CX*NX19RE
Oxygenator without Hardshell Reservoir	[REF] : 3CX*NX19W	[REF] : 3CX*NX19E
Place inlet port facing North as indicated in the illustration. The blood outlet port orientation is (W)est - LEFT side or (E)ast - RIGHT side. Note: this illustration is a top view of the oxygenator.	Blood inlet port Blood outlet port Left (West)	Blood inlet port Blood outlet port Right (East)

STRUCTURE

WARNINGS

Describe serious adverse reactions and potential safety hazards and limitations in use imposed by them.

WARNINGS

- The CAPIOX® NX19 is designed to operate at blood flow rates within the range of 2.0 L/min (or 0.5 L/min for up to 2 hours) to 8.0 L/min. Do not use any blood flow rates outside this range.
- DO NOT OBSTRUCT GAS VENTING AREAS. Avoid build up of excess pressure in the gas phase to prevent gaseous emboli entering the blood phase.
- Pressure in the blood phase should always be higher than that in the gas phase to prevent gaseous emboli entering the blood phase.
- The gas flow rate should not exceed 20 L/min for the NX19. Excessive gas flow rate will bring about pressure increase in the gas phase, allowing gaseous emboli to enter the blood phase.
- During recirculation, do not use pulsatile flow, clamp the arterial-venous circuit, or stop the blood pump suddenly as these actions may cause gaseous emboli to enter the blood phase from the gas phase due to inertia force.
- When the CAPIOX® NX Oxygenator module is used separately from the Hardshell Reservoir, set the module so that the upper end of the fibers is lower than the blood level in the venous reservoir. This prevents gaseous emboli from entering the blood phase from the gas phase.
- To prevent gaseous emboli from entering the blood phase, make sure that the arterial pump flow rate always exceeds the flow rate of the cardioplegia line and/or ancillary pumps. The blood flow rate of the cardioplegia line and/or ancillary pumps should not exceed 1 L/min.
- Ensure the oxygenator is never subjected to a negative pressure. Negative pressure in the oxygenator can cause gaseous emboli to enter the blood phase from the gas phase.
- When PVC tubing is used in the roller pump, water must flow through the heat exchanger prior to blood pump operation to prevent the buildup of static electricity in the heat exchanger.
- Adequate anticoagulation of the blood is required considering patient comorbidities, including hypercoagulable states, and perfusion technique to prevent it from clotting in the system.
- Once blood is introduced to the device, flow through the oxygenator and reservoir must be maintained to prevent stasis of blood (even if patient flow is stopped).
- When using the centrifugal pump, clamp the arterial line distal to the oxygenator (the patient's side) before stopping the pump. Improper clamping may cause uncontrolled back-flow of blood or migration of gaseous emboli into the blood path of the oxygenator.
- Do not use this product for a period in excess of 6 hours. Excessive use for over 6 hours may lead to plasma leak and thrombi formation, which may compromise the gas exchange performance.
- Temperature of water delivered to heat exchanger should not exceed 42°C (108°F), otherwise blood damage may result.
- Do not exceed a 15°C (27°F) temperature difference between blood and water in heat exchanger to prevent gas dissolved in blood from forming bubbles.
- Avoid introduction of hemostatic agents into the Cardiectomy Reservoir/Oxygenator. Hemostatic agents are known to generate thrombi which may compromise performance of the Reservoir/Oxygenator.

Use CAPIOX® NX Oxygenator with Hardshell Reservoir taking note of the additional warnings below.

- The blood flow into the cardiectomy filter should not exceed the rate of 5 L/min. Excessive blood flow rate may increase the pressure in the cardiectomy filter, resulting in back-flow into any solution or blood administration line connected to the reservoir.
- Minimum operating volume in the reservoir is 150 mL. Set appropriate blood storage level, relative to venous flow rate, to prevent gaseous emboli passing to patient. (Refer to page 7 "Safety Margin")
- The yellow vent port cap need not be removed, since this port assures adequate venting with the cap attached. Do not close the vent port, as this may cause positive pressure in the reservoir, resulting in back-flow into the patient, solution or blood administration lines connected to the reservoir.
- Negative pressure below -150 mmHg (-20 kPa) should not be applied to this reservoir as this may damage it.
- When the level of blood in the reservoir is below the minimum operating level, blood flow at a high flow rate into the cardiectomy filter may generate gaseous micro emboli (GME) which could migrate to the patient.

PRECAUTIONS

Include information regarding any special care to be exercised by the practitioner for the safe and effective use of the device.

PRECAUTIONS

- Do not store product at extreme temperatures and humidity. Avoid direct sunlight.
- Do not use solvents such as alcohol, ether, acetone, etc. These solvents may cause damage if used in or on the device.
- Pressure at the blood inlet of the oxygenating module exceeding 1,000 mmHg (133 kPa) may cause leaks or damage to the device.
- Water pressure at the heat exchanger inlet exceeding 2 kgf/cm² (1,470 mmHg) (28 PSI) (196 kPa) may cause leaks or damage to the device.
- Use caution when removing air during priming and perfusion. Excessive shock to the device - especially with hard objects - can cause damage to the device.
- The responsibility for the use of this device in all cases belongs solely to the physician ordering its use.
- The use-by date can be found on the device labeling. Do not use the device past the use-by date.
- The product is only to be used in the operating room environment.
- Read all the instructions of the drugs and the medical devices to be used in combination with this device prior to use.
- Use of phospholipidic drugs in the priming fluid or their infusion in the patient during CPB may facilitate the hydrophilization (plasma breakthrough) of the microporous membrane.
- This device should only be used by properly trained and qualified personnel.
- This device is sterile and non-pyrogenic in the fluid pathway in the unopened, undamaged package. Inspect the device and package carefully. Do not use if the package and / or device is damaged, or if caps are not in place.
- Dispose of device in accordance with hospital and government regulations for compliance with proper medical waste disposal and to avoid risk of infection.
- This device is intended for single use only. DO NOT REUSE. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility, and functional integrity of this product.
- In case of the need to replace the oxygenator, make sure that a spare unit is available.
- Use aseptic technique in all procedures.
- Perform priming using crystalloid solution which contains no blood, plasma and/or blood derived products.
- Place the CAPIOX® NX Oxygenator module downstream from the blood pump.
- A minimum of 0.5 L/min oxygen gas flow is needed when blood is circulated. Less than 0.5 L/min oxygen gas flow may result in inadequate gas exchange. Ensure that adequate gas flow is established when on bypass.
- Do not allow anesthetic liquids to contact the device. Contact with anesthetic liquids may damage the device.
- When using volatile anesthetics a gas scavenging line should be utilized to prevent application of anesthetic gas to the operating room environment.
- Gas scavenging may impact accuracy of capnography sampling.

- Cardiopulmonary bypass circuit temperature and pressure (distal to the oxygenator module) must be monitored. Additional pressure monitoring of the oxygenator module inlet pressure is recommended as an early indicator of increased pressure drop.
- Stop the gas flow when cardiopulmonary circulation is suspended. During recirculation, check blood gas parameters. Excessive gas flow may cause low PaCO₂, alkalosis, or blood damage.
- Upon patient rewarming, adjust O₂ concentration, gas flow rate and blood flow rate by increasing them as needed based on an increase in patient metabolism. Failure to adjust the gas supply and the blood flow rate appropriately may cause insufficient gas exchange.
- Use of a pre-bypass filter is recommended to remove any particulate matter from the circuitry and priming solution.
- The Sampling Line is detachable. If the line is detached and not in use, make sure to cap the end.
- Use of a recirculation line is recommended in the arterial line after the NX oxygenator to facilitate gaseous emboli removal during priming and emergency situations.
- It is recommended to use an air bubble detector between the reservoir outlet and NX oxygenator inlet. In the event gaseous emboli are detected, determine the source and take corrective action to eliminate further introduction of gaseous emboli into the perfusion circuit.

Use CAPIOX® NX Oxygenator with Hardshell Reservoir taking note of the additional precautions below.

- Use the safety devices, including the level sensor and the air detector, in order to control the volume of the blood in the reservoir and to detect gaseous emboli in the arterial lines.
- When rotating the venous blood inlet port on the top of reservoir, make sure that venous line and thermistor probe cable do not come into contact with any solution or administration line connected to the reservoir to prevent it from kinking.
- When separating the reservoir from the oxygenator module, hold the oxygenator firmly and remove connecting ring.
- To remove the L shaped connector, turn the rotating collar away from you when facing the front of the venous return section of the reservoir. Then pull the L shaped female connector off the manifold.
- The caps of unused ports should be left in place. This avoids contamination and prevents leakage of blood.
- Make sure the caps of unused luer ports are tightened firmly to prevent leakage.
- The cardiotomy filter must be wetted with crystalloid solution prior to use in order to achieve the maximum flow capacity.
- Where appropriate, introduce blood, plasma and/or blood derived products through the quick prime port or any of the luer ports leading to the cardiotomy filter after priming and de-airing of the oxygenator.
- Do not exceed the maximum operating volume of 4000mL in the reservoir. This device is designed to be used with the volume of blood in reservoir below 4000mL.
- If the blood entering the Cardiotomy Reservoir filter inlet appears to spill over versus entering the filter this may indicate that the filter is clogged. Discontinue using the Cardiotomy Reservoir filter. Replace the reservoir or flow the Cardiotomy blood through an auxiliary Cardiotomy reservoir and connect it to the auxiliary port.
- Scientific literature contains reports of altered dose response to some medications, such as nitroglycerin, fentanyl, etc., possibly due to varying degrees of absorption by synthetic materials such as those found in cardiopulmonary circuits. Avoid injection of medications into the Cardiotomy Reservoir filter.
- Do not administer drugs through any line or port leading to the cardiotomy reservoir. The medication delivery could be delayed due to inconsistent blood flow through the cardiotomy reservoir and filter.
- Use of a filtered gas source is recommended to prevent potential debris / contaminants from entering gas inlet port of the oxygenator.
- Ensure cleaning chemicals are flushed from the heater cooler according to the manufacturer's instructions for use.

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

METHOD OF OPERATION

Read **WARNINGS** (page 4) and **PRECAUTIONS** (page 4-5) before use.

The following is a basic description of the use of the CAPIOX® NX Oxygenator and Hardshell Reservoir when no other particular indications apply.

Note: The reservoir volumes illustrated by the scale labels are not absolute volumes. Please keep in mind when considering safe operating volumes.

SETUP

1. Remove CAPIOX® NX from package and examine it for defects. Do not open the cardboard box in the operating room. The cardboard box may create particulate.

CAUTION Do not use if the package or device is damaged (e.g. cracked) or any of the port caps are off.

2. Place CAPIOX® NX in its holder.

CAPIOX® NX Oxygenator with Hardshell Reservoir: Set the holder at the appropriate height to prevent the arterial line and other lines from kinking.

Making sure that the connecting ring is in place, set oxygenator in its holder ([REF]: XX*CXH18R, XX*XH032, 801804 or 801139). If the reservoir is not set into the holder correctly, the amount of fluid in the reservoir cannot be correctly assessed.

When rotating the oxygenator module, be sure to hold the support arm. If the oxygenator module is rotated while holding the oxygenator body, it could be damaged.

Note: When using CAPIOX® NX Oxygenator module and Hardshell Reservoir separately,

- a) Hold the module firmly and separate it from Hardshell Reservoir by removing connecting ring.
- b) Set reservoir in the holder ([REF]: XX*CXH18R, XX*XH032, 801804 or 801139) for CAPIOX® NX Oxygenator with Hardshell Reservoir.
- c) Set the oxygenator module in its holder ([REF]: XX*CXH18, XX*CXH25, XX*CXH15).

CAUTION If the product is dropped during set-up, do not use it. Replace with another device.

3. Connect water lines (1/2" (12.7mm) tubing or couplers) to the water ports of CAPIOX® NX (water ports are bi-directional).

Note: The integrity of the heat exchanger water path has been verified with hydrogen peroxide solution (≤ 330 ppm in water).

Chemicals other than hydrogen peroxide, or hydrogen peroxide solutions >330 ppm, have not been verified for safe use with this product.

4. Start water circulation through heat exchanger and circulate for at least 5 minutes. Check for leaks from water path to blood path and water path to the outside of the device.

WARNING When PVC tubing is used in the roller pump, water must flow through the heat exchanger prior to blood pump operation to prevent the buildup of static electricity in the heat exchanger.

5. Attach venous line to the venous blood inlet port of the venous reservoir.
6. Attach arterial line to the blood outlet port of CAPIOX® NX Oxygenator.
7. Remove the cap from the blood cardioplegia port and attach 1/4" (6.4mm) blood cardioplegia line tubing.

CAUTION If the cardioplegia port is not used, push cap in and twist it to seal tightly or attach a short 1/4" (6.4mm) tubing and clamp for security.

8. If using a roller pump, connect the roller pump inlet line to the blood outlet port of Hardshell Reservoir. Connect the roller pump outlet line to the blood inlet port of the oxygenator. If using a centrifugal pump, connect the tubing from the reservoir outlet port to the pump inlet port and the pump outlet tubing to the oxygenator inlet port.
9. Add a recirculation line with 3/16" (4.8mm) or larger diameter to the arterial line after oxygenator blood outlet port.
10. Push all tubing past the second barb and band all connections in the circuit.

11. Connect 1/4" (6.4mm) gas line to gas inlet port. The NX19 oxygenator is compatible with volatile anesthetics (isoflurane and sevoflurane) in concentrations of up to 5% for duration of up to 6 hours. If volatile anesthetics are used, some method of scavenging the gas from the oxygenator should be considered. The protocol, the concentration, and the monitoring of the anesthetic gases administered to the patient, is under the sole responsibility of the physician in charge of the treatment.

CAUTION Do not allow anesthetic liquids to contact the device. Contact with anesthetic liquids may damage the device. When using volatile anesthetics a gas scavenging line should be utilized to prevent release of anesthetic gas to the operating room environment.

12. Aux gas out port may be used for capnography sampling.

CAUTION Gas scavenging may impact accuracy of capnography sampling

WARNING DO NOT OBSTRUCT GAS VENTING AREA.

13. Thermistor probes can be used with Y.S.I. 400 series* temperature monitor by use of the following parts.

Blue Cable (venous line): **REF** : CX*BP021

Red Cable (arterial line): **REF** : CX*BP022 *Y.S.I. 400 (trade name : Measurement Specialties Inc.)

14. Suction ports of the CAPIOX® Hardshell Reservoir have blue caps. Remove the blue caps, connect the suction lines and vent lines to suction ports on reservoir.

15. For use of CAPIOX® NX Oxygenator module without Hardshell Reservoir, connect venous side male luer lock of sampling system to venous line.

16. Prior to priming the cardiopulmonary bypass circuit, confirm all tubing connections and accessory devices are compatible, secure and correct. Also check the NX for water to blood leaks and adequate water supply prior to priming.

NOTE: For the use of the sampling system separated from Hardshell Reservoir, sampling manifold holder (Code No.: XX*XH051) is available.

CAUTION Ensure that all connected parts including the luer caps, the lock adapters and the port caps are securely affixed. Loose connections may cause contamination or a blood leak. Do not apply the connections with excessive force as this may cause damage to the product or accessory.

PRIMING PROCEDURE

CAUTION CAPIOX® NX Oxygenator module should be primed with crystalloid solution containing no blood, plasma, or blood-derived products.

NOTE: CAPIOX® NX can be primed without a CO₂ flush. Incorporating a CO₂ flush with initial priming of the CAPIOX® NX facilitates de-airing the oxygenator module and cardiopulmonary bypass circuit.

1. Position a tubing clamp on the venous reservoir outlet tubing and introduce crystalloid priming solution through the quick prime port or any of the luer ports leading to the cardiotomy filter.
2. Control prime solution flow between reservoir outlet and oxygenator inlet with tubing clamp to minimize fluid-air mixing. Remove all bubbles in this section of tubing (roller pump head/boot) or centrifugal pump.
3. Make sure the purge line is open and the arterial outlet tubing is clamped distal to the recirculation line.
4. Continue priming the oxygenator either by gravity or by initiating forward pump flow. Start the pump at a low speed to begin recirculation to complete priming of the cardiopulmonary bypass circuit.
5. To ensure the system is completely de-aired and no leaks in the system exist, recirculate the priming solution at a rate of 4L/min or higher and establish a system pressure that mimics clinical conditions. This can be accomplished by partially occluding the arterial line distal to the oxygenator and pressure monitoring transducer location. Check oxygenator and tubing for leakage or air.

WARNING • Do not supply gas during priming without adequate arterial pump flow during recirculation in order to prevent the pressure in the gas phase from becoming higher than the fluid phase.

- CAUTIONS** • Do not use an oxygenator and reservoir that leaks. Replace it with another CAPIOX® NX oxygenator and reservoir.
- Do not use a sampling line, a purge line, or tubing with internal diameter less than 3/16" (4.8mm) as a recirculation line. If used, the oxygenator module can be damaged as a result of excessive positive pressure.
 - Recirculate the priming solution at a rate of 4L/min or higher to facilitate air removal. Failure to remove air from the oxygenator may result in serious injury to the patient.
 - Reservoir volume must be 150mL or greater at all times. (Refer to page 7 "Safety Margin")
 - Returning the priming solution back to the cardiotomy filter when there is an insufficient level of solution contained within the reservoir may generate gaseous emboli. Maintain an adequate level of solution in the reservoir. (Refer to page 7 "Safety Margin")

6. After adequate de-airing, introduce blood or blood-derived products through the quick prime port or any of the luer ports leading to the cardiotomy filter.
7. In preparation for cardiopulmonary bypass, close the sampling line, purge line, and recirculation line to prevent arterial to venous shunting and blood flowing back into the reservoir due to the patient's blood pressure and hydrostatic pressure once the patient is cannulated.
8. After priming, if air bubbles continue to appear make necessary corrections to remove the air.

WARNING During recirculation, do not use pulsatile flow, clamp the arterial-venous circuit, or stop the blood pump suddenly as these actions may cause gaseous emboli to enter the blood phase from the gas phase due to inertia force.

INITIATION OF BYPASS

WARNING The blood phase pressure must always be higher than the gas phase pressure in order to avoid air embolisms.

1. Check for adequate anticoagulation level prior to initiating bypass.
2. Ensure that the de-airing process is complete prior to initiating bypass. If air remains, repeat the "Priming Procedure" to dispel air.
3. Confirm that the gas venting mechanisms are not obstructed prior to starting the gas supply.
4. Remove arterial and venous clamps and gradually increase blood flow. Then initiate gas flow with a V/Q = 1 and FiO₂ = 100%

WARNING Obstructing the gas venting mechanisms may lead to pressure build-up in the gas phase, allowing gaseous emboli to enter the blood phase.

DURING PERFUSION

1. Measure blood gases and make necessary adjustments as follows:

- | | |
|---|--|
| a. Control PaO ₂ by changing concentration of oxygen in ventilating gas using gas blender. | b. Control PaCO ₂ by changing the total gas flow. |
| -To decrease PaO ₂ , decrease FiO ₂ . | -To decrease PaCO ₂ , increase total gas flow. |
| -To increase PaO ₂ , increase FiO ₂ . | -To increase PaCO ₂ , decrease total gas flow. |

WARNINGS • Collect blood only while pump is running, or blood-side pressure will decrease and air bubbles may result.

- A phenomenon called wet lung may occur when water condensation occurs inside fibers of microporous membrane oxygenators with blood flowing exterior to the fibers. This may occur when oxygenators are used for a longer period of time. If water condensation and/or a decrease in PaO₂ and/or an increase in PaCO₂ is noted during extended oxygenator use, briefly increasing the gas flow rate may improve the performance. Increase gas flow rate, to 20 L/min for 10 seconds. Flushing more than once will not improve performance further.

- Do not exceed a V/Q = 2 during normal perfusion.

CAUTION • A minimum of 0.5 L/min oxygen gas flow is needed when blood is circulated. Less than 0.5 L/min oxygen gas flow may result in inadequate gas exchange while on cardiopulmonary bypass.

- Adjust blood flow rate, making sure that the blood level in the reservoir is adequate, to prevent gaseous emboli passing to patient.
- For purging, open the purge line.

CAUTION Do not give any shocks to the product during perfusion.

TERMINATION OF BYPASS

Terminate cardiopulmonary bypass per institutional protocol and surgeon direction.

WARNINGS • Before reducing the blood flow rate, confirm that the purge line is closed and that the sampling line is also closed with arterial side stopcock.

- On stopping pump, stop gas flow at once to ensure gas phase pressure does not exceed fluid phase pressure.
- If circulation is to be resumed, recirculation at a low flow rate is recommended. During recirculation, check blood gas parameters. An excessive gas flow may cause low PaCO_2 , alkalosis, or blood damage.

CAUTION • Prior to resuming bypass, set FiO_2 at 100% to ensure adequate oxygenation. Increased PCO_2 and decreased PO_2 in the patient's blood at the beginning of recirculation may not be able to be recovered without appropriate gas supply.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local, state, federal and international laws and regulations for medical waste disposal.

OXYGENATOR REPLACEMENT

Always have a spare CAPIOX® NX oxygenator for replacement if needed.

- Set up and prime spare oxygenator module as described in SET-UP section. (page 5)
- When necessary, maintain a decreased patient temperature as prescribed by physician.
- Disconnect all accessory devices and connections.
- Double-clamp blood inlet and outlet tubing connected to new pre-primed oxygenator, and cut the tubing between the two clamps.
- Stop circulation, double-clamp blood inlet, cardioplegia, and blood outlet lines connected to old oxygenator, and cut the lines between the two clamps.
- Replace old oxygenator with new one by connecting venous, arterial, and cardioplegia lines to the connectors at blood inlet and outlet connectors of the new oxygenator.
- Connect gas line removed from old oxygenator module to the new module.
- Connect the purge line of the new oxygenator to the luer port leading to the cardiotomy section of the reservoir.
- Connect the sampling line to the new oxygenator.

CAUTIONS • After replacement, open recirculation and purge line for de-airing.
• Close recirculation and purge line prior to reinitiating cardiopulmonary bypass.


WARNING • Band and secure all connections in the circuit to prevent blood loss.

- Start circulation at a low flow rate and start gas supply.
- Connect arterial temperature probe.
- Connect water lines removed from old oxygenator module to the new module and start water supply. Check for leaks from water path to blood path and water path to the outside the device.

SUGGESTED ACCESSORIES

- Thermistor Probe (Thermistor probes can be used with Y.S.I. 400 series* temperature monitor by use of the following parts)
Blue Cable (venous line) : [REF] : CX*BP021 & Red Cable (arterial line) : [REF] : CX*BP022 * Y.S.I. 400 (trade name : Measurement Specialties Inc.)
- Oxygenator and Reservoir Holder
[REF] : XX*CXH18R, XX*XH032, 801804 or 801139
- Oxygenator Holder
[REF] : XX*CXH18, XX*CXH25, XX*CXH15
- Pre-bypass Filter
- Reservoir Level Sensor
- Air Detector
- Manometer


SAFETY MARGIN

* Minimum operating blood level in CAPIOX® Hardshell Reservoir.

The above graph indicates the time in which the blood volume in the reservoir decreases to minimum operating blood level of 150 mL at each flow rate. If the blood volume in reservoir is below 150 mL, bubbles may pass to the oxygenator module.

NX19 PERFORMANCE DATA

Data related to blood cell damage is available upon request.

Vacuum Assisted Venous Drainage

INTENDED USE

The CAPIOX® Hardshell Reservoir is also intended for use in vacuum assisted venous drainage procedures. The use of a controlled vacuum on the hardshell improves venous drainage during minimally invasive surgical techniques or regular bypass surgery.

WARNINGS

- For drug administration in the Hardshell Reservoir, be aware that negative pressure in this line could draw more drug than intended. To prevent non intended drug amounts to be drawn, it is recommended to include only the intended amounts in the syringe.
- High negative pressure may lead to cell damage. Use the lowest vacuum level necessary to achieve desired flow rate.
- Discontinue vacuum assisted venous drainage prior to stopping the arterial blood pump.
- When stopping the pump or when the blood flow rate is low during VAVD, all A-V shunt lines (ex. sampling line, purge line, etc) must be closed to prevent drawing air from the gas phase to the blood phase and preventing retrograde and antegrade blood flow from the patient.
- Check for proper occlusion because an under occluded roller pump is more likely to draw air from the gas phase to the blood phase when VAVD is performed.
- Do not open the auxiliary port on the venous reservoir during VAVD to prevent bubble formation inside the venous filter

When VAVD is performed in conjunction with a centrifugal pump take caution as below.


- The line between the oxygenator and the centrifugal pump must be clamped before stopping the pump. Not clamping the arterial line may draw air into the blood side of the oxygenator from the gas phase. An arterial line one-way valve is recommended between the oxygenator and the reservoir.
- Caution must be taken because the relationship between blood flow rate and pump speed changes depending on the degree of negative pressure applied.
- Avoid obstruction or occlusion of the vacuum line to prevent retrograde air to the patient.

PRECAUTIONS

- Always allow for a vented reservoir when vacuum is not being applied. Vent by declamping the vacuum release line.
- Reduce the negative pressure from the vacuum source slowly when terminating VAVD. A sudden change in pressure may create turbulence of the blood inside the reservoir.
- Do not exceed -150 mmHg (-20 kPa) when using vacuum source for assisted venous drainage technique to prevent reservoir damage.
- Do not connect the VAVD circuit to any other ports than the reservoir vent port.
- Do not leave the ports that are not in use open.
- A moisture trap is required to trap condensation generation.
- Use a sterile VAVD circuit and do not reuse or resterilize it.
- A positive pressure-warning device should be equipped on the reservoir.
- A controlled vacuum regulator is required.
- The use of a negative pressure manometer gauge on the Hardshell Reservoir and a negative pressure relief valve [opening range at -150 mmHg (-20 kPa)] are recommended to prevent excessive negative pressure.
- Use of vacuum may cause liquid level to appear higher than actual due to reservoir compliance. A liquid level of 170 mL on volume scale assures actual minimum volume of 150 mL at a negative pressure of -150 mmHg (-20 kPa).

REQUIRED EQUIPMENT

- A CAPIOX® Hardshell Reservoir
- B Sterile moisture trap
- C Clamp
- D Controlled vacuum regulator [adjustable between 0 to -150 mmHg (-20 kPa)]
- E Sterile manometer gauge line
- F Negative manometer gauge
- G Sterile vacuum release line
- H Sterile positive pressure relief valve [integrated in reservoir with opening range from 0 to 8 mmHg (1.1 kPa)]
- I Sterile gas filter
- J Negative pressure relief valve

SET-UP CONFIGURATION**INSTRUCTIONS FOR USE**

Note : A~J refers to the Required Equipment for use with vacuum assisted venous drainage with letter indications.

1. Prepare CPB circuit set-up as per standard procedure.
2. Verify that blue caps from suction inlet ports and auxiliary port of reservoir are pushed all the way to insure adequate seal.
3. Secure in place all yellow luer caps; all caps are non-vented.
4. Install controlled vacuum regulator D on wall suction source.
5. Attach a 3-way stopcock to the non-filtered luer lock port, and connect the sterile manometer gauge line E.
6. Connect the negative manometer gauge F to the sterile manometer line E.
7. Attach the gas filter I to the tubing connected to the moisture trap B and the "Y" connector.
8. Attach sterile moisture trap B tube to the vent port on the hardshell reservoir.

Initiating Bypass

1. Begin with gravity venous drainage; at this time, the vacuum release line G is unclamped.
2. To initiate vacuum assist venous drainage: set the vacuum regulator D then clamp the vacuum release line G.
3. Monitor the negative pressure inside the Advance Hardshell Reservoir with the negative manometer gauge F.
4. Adjust the negative pressure to optimize venous return. Set vacuum regulator D.

Note: Ensure the suction and vent roller pumps are properly occluded.

Coming off Bypass

Unclamp vacuum release line G and turn off the vacuum regulator. Come off bypass as per standard procedure.

Post-Operative Chest Drainage**Contraindications**

Post-operative chest drainage with subsequent autotransfusion procedures are contraindicated in the following events:

- Gross perforations to the chest wall or an air leak in the lungs.
- Infection or malignancy occurring systemically or in the pericardium, mediastinum, or lungs.
- Suspected or evident gross contamination with foreign material, lymphatic failure, or perforated intestine.
- Presence of the following in the aspiration site: topical hemostatic agents, bactericidal wound irrigants or antibiotics not intended for parenteral administration.
- Open chest and vacuum applied.
- Administration of protamine prior to reservoir being removed from bypass circuit.
- Patient returned to surgery for any reason.
- Use of vented chest tubes that do not have vent flow regulation such as stopcocks.

¹Standards for Blood Banks and Transfusion Services, 16th ed. American Association of Blood Banks

Warnings

A qualified person should assess the quality and suitability of returning any blood that has been collected before reinfusion begins. The safe reinfusion of collected fluids is the sole responsibility of the attending physician.

- It is recommended that a transfusion filter be used in all reinfusion procedures.
- Complications such as the following have been associated with chest drainage and subsequent reinfusion: blood trauma, blood coagulation, coagulopathies, and particulate or air embolism.
- Reinfusion of collected blood/fluid should be done on an hourly basis unless fewer than 50 mL are collected hourly.
- A minimum reservoir level of 20 mL of fluid should always be maintained to prevent air embolus passing to patient.
- Blood that has been in the reservoir for longer than 4 hours should not be transfused.¹
- It is not recommended that autotransfusion continue longer than 18 hours after surgery.²
- Attach a shunt line (the bridge) to connect the filtered and unfiltered sections in the reservoir as indicated in instructions for use.
- When using the reservoir, monitor the attached shunt line for blood in the line. If blood is detected in the line this is an indication that the reservoir filter has clotted. Replace reservoir promptly.
- Filter occlusion during high volume chest drainage may cause blood/fluid to pass through the external shunt line bypassing filtration. All fluids that pass through the shunt line must be filtered prior to reinfusion.
- Accepted medical and nursing care routines must be followed during chest drainage.
- When vacuum is used during chest drainage, do not exceed -150 mmHg (-20 kPa) (-204 cmH₂O).

CAPIOX® Hardshell Reservoir Operating Specifications for Post-Operative Chest Drainage

ITEM	SPECIFICATION
Cardiotomy Maximum Blood Flow Range	5 L/min
Blood Storage Capacity	4,000 mL
Minimum level in reservoir	20 mL

List of Required Materials for Use of Reservoir for Post-Operative Chest Drainage

REFERENCE	QUANTITY	DESCRIPTION	TUBING LENGTH
[a]	6 each	1/4" (6.4mm) I.D. dead-end tubing with cap	2" (51mm) each
[b]	1 each	3/8" (9.5mm) I.D. tubing with caps	12" (305mm)
[c]	1 each	1/2" (12.7mm) x 1/2" (12.7mm) connector with 1/2" (12.7mm) I.D. dead-end tubing with cap.	6" (152mm)
[d]	2 each	3/8" (9.5mm) x 3/8" (9.5mm) x 3/8" (9.5mm) Y connector with 3/8" (9.5mm) I.D. 40 Durometer tubing with caps	96" (2438mm)
[e]	1 each	3/8" (9.5mm) I.D. dead-end tubing with cap	4" (102mm)
[f]	1 each	3/8" (9.5mm) perfusion adapter attached to 4" (102mm) of 1/8" (3.2mm) I.D. tubing that includes a clamp with female luer end, connected to male luer with 48" (1219mm) of 1/8" (3.2mm) I.D. tubing with another female luer end connected to a male luer with 2.5" (64mm) of 3/16" (4.8mm) I.D. tubing; the end of the tubing is covered by a cap.	N/A
[g.1.]	1 each	Pressure relief valve: opening vacuum range: 22-44 mmHg (2.9-5.9 kPa) (30-60 cmH ₂ O) with bacterial filter	N/A
[g.2.]	1 each	Pressure relief valve: opening pressure of approximately 8 mmHg (1.1 kPa) (11 cmH ₂ O) [Integrated]	N/A
[h]	1 each	1/4" (6.4mm) I.D. tubing	As needed
[i]	1 each	Water manometer	N/A
[j]	1 each	1/4" (6.4mm) I.D. tubing	18" (457mm)
[k]	1 each	Male non-vented luer lock cap	N/A
[l]	1 each	Four way stopcock capable of opening three ways simultaneously	N/A

²Page.R, et al, Hard-Shell Cardiotomy Reservoir for Reinfusion of Shed Mediastinal Blood. Ann Thorac Surg 1989; 48:514-7.

INSTRUCTIONS FOR USE

Note: When referring to the instructions for Post Operative Chest Drainage, the bracketed letters refer to the List of Required Materials on page 9.

- Carefully remove tubing and/or blue caps from suction inlet ports of reservoir. Plug five of the suction inlet ports with 1/4" (6.4mm) I.D. deadend tubes [a].
- Attach 3/8" (9.5mm) I.D. tubing [b] from the auxiliary port to one of the suction inlet ports; all suction inlet ports are 1/4" (6.4mm) - 3/8" (9.5mm). This will create a bridge between the filtered and unfiltered sections of the reservoir.
- Clamp the 1/2" (12.7mm) circuit tubing attached to the 1/2" (12.7mm) venous blood inlet port of the reservoir, cut it and connect a 1/2" (12.7mm) connector with 1/2" (12.7mm) I.D. dead-end tubing [c].
- Connect the chest drainage tubes to the 3/8" (9.5mm) X 3/8" (9.5mm) Y connector with 96" (2438mm) of 3/8" (9.5mm) 40 Durometer tubing [d]; if only one chest tube is used, plug one of the 3/8" (9.5mm) ports of the of the Y connector with a 3/8" (9.5mm) dead-end tube [e]. Attach the other end of the 3/8" (9.5mm) X 3/8" (9.5mm) Y connector with 96" (2438mm) of tubing to the vertical port to CR filter. In case 3 or more chest tubes are used, attach an additional 3/8" (9.5mm) X 3/8" (9.5mm) Y connector with 96" (2438mm) of 3/8" (9.5mm) 40 Durometer tubing [d] to one of the suction inlet ports that was capped off by a dead-end tube.
- Secure in place all yellow luer caps; all caps are non-vented.
- Disconnect sampling system from the venous blood inlet port and attach a non-vented luer lock cap [k].
- Attach a 1/4" (6.4mm) dead-end tube [a] to the quick prime port on the reservoir.
- Hold the reservoir parallel to the floor with the venous blood outlet port at the top to prevent blood spillage. Clamp the circuit tubing attached to the outlet port, cut it, then connect a 3/8" (9.5mm) to 1/8" (3.2mm) blood infusion adapter line [f]. Clamp the adapter line and leave the end cap in place until autotransfusion is set up. In case of using 1/4" (6.4mm) adapter, take 1/4" (6.4mm) adapter off first.
- Remove luer cap from non-filtered luer port and attach three-way stopcock. Secure to the stopcock the negative pressure relief valve with attached bacterial filter [g.1.].
- Install the reservoir into the autotransfusion holder (product code: XX*XH032).

CAUTION Maintain the hardshell reservoir below the thorax to facilitate chest drainage.

- Using a 1/4" (6.4mm) I.D. tubing [h], connect the vent port on top of the hardshell reservoir to a regulated vacuum source with 11-15 mmHg (1.5 to 2.0 kPa) (15-20 cmH₂O).

Note: The "DO NOT OBSTRUCT" caution label near the vent port is not applicable for post-operative chest drainage and autotransfusion.
OR

- The use of a water seal/water manometer [i] is highly recommended for chest drainage. If using a water manometer, connect one end of a 1/4" (6.4mm) I.D. tubing [j] to the vent port and the other end to the water manometer.

WARNING High negative pressure may lead to cell damage. Use the lowest vacuum level necessary to achieve desired flow rate.

- Record the time that drainage into the reservoir begins.

WARNING It is recommended that autotransfusion continue no longer than 18 hours after surgery.

- Continually monitor drainage and elapsed time. Determine hourly drainage rate.

WARNING Reinfusion of collected blood/fluid should be done on an hourly basis unless fewer than 50 mL are collected hourly.

- Upon collection of 50 mL of blood, remove the end cap of the blood infusion adapter line and attach the blood infusion adapter to an I.V. Administration set/pump apparatus.

Note: It is recommended that the infusion pump be used with an air detection sensor for all reinfusion procedures.

- Unclamp the blood infusion line and proceed with autotransfusion.

WARNINGS • To reduce risk of air embolism, evacuate all air completely from the infusion line prior to autotransfusion.

- The physician is responsible for determining the optimal time, rate and volume to be autotransfused to the patient.
- The clinician should be aware of potential complications such as air emboli, associated with the use of the I.V. Administration set and the infusion pump. Use the I.V. Administration set and the infusion pump in accordance with the instructions provided by the manufacturers.

SPECIFICATIONS

Product Description

Reorder Number	Description
3CX*NX19RW	CAPIOX® NX19 Oxygenator with Integrated Arterial Filter and UltraPrime™ Technology
3CX*NX19RE	CAPIOX® NX19 Oxygenator with Integrated Arterial Filter and UltraPrime™ Technology
3CX*NX19W	CAPIOX® NX19 Oxygenator with Integrated Arterial Filter and UltraPrime™ Technology
3CX*NX19E	CAPIOX® NX19 Oxygenator with Integrated Arterial Filter and UltraPrime™ Technology

(Refer to "Direction of blood outlet port", page 3)

CAPIOX® NX19 Oxygenator Module

COMPONENT	SPECIFICATIONS	
Housing	Material	Polycarbonate
Fibers	Material Surface area	Microporous Polypropylene Approx. 1.9 m ²
Arterial filter	Material	Polyester screen type Pore size 32 µm
Heat Exchanger	Material Surface area	Polyethylene Terephthalate (PET) Approx. 0.43 m ²
Blood Flow Range	Min. 2.0 L/min (or 0.5 L/min for up to 2 hours) Max. 8.0 L/min	
Priming Volume (Static)	185 mL ± 10%	
Blood inlet port (from pump)	3/8" (9.5 mm)	
Blood outlet port	3/8" (9.5 mm)	
Cardioplegia port	1/4" (6.4 mm)	
Gas inlet port	1/4" (6.4 mm)	
Gas outlet port	5/16" (7.9 mm)	
Water ports	1/2" (12.7 mm) Hansen quick connect fitting	
Maximum Pressure	Blood inlet 1,000 mmHg (133 kPa) (1.36 kgf/cm ²) Water inlet 1,470 mmHg (196 kPa) (28 PSI) (2 kgf/cm ²)	

CAPIOX® Hardshell Reservoir

COMPONENT	SPECIFICATIONS	
Housing Material	Polycarbonate	
Blood Flow Range	2.0 L/min (or 0.5 L/min for up to 2 hours)	
Venous Flow:	Min.	8.0 L/min
	Max.	5.0 L/min
Cardiotomy Inlet:	Max.	8.0 L/min
Combined Flow:	Max.	
Blood Storage Capacity	4,000 mL	
Minimum Operating Volume	150 mL	
Venous Filter	Material Pore Size	Polyester screen type 47 µm
Cardiotomy Filter Material	Polyester depth type	
Defoaming Part Material	Polyurethane foam	
Venous blood inlet port	1/2" (12.7 mm) rotatable	
Blood outlet port (to pump)	3/8" (9.5 mm)	
Suction ports	six 1/4" (6.4 mm)	
Vertical port to CR filter Quick	3/8" (9.5 mm)	
Prime port	1/4" (6.4 mm)	
Vent port	1/4" (6.4 mm)	
Auxiliary port	1/4" — 3/8" (6.4 mm — 9.5 mm)	
Luer ports	- three filtered luer locks to cardiotomy filter - a non-filtered luer lock - two luer locks on venous inlet	
Positive Pressure Relief Valve	0 to 8 mmHg (1.1 kPa)	
Maximum Sustainable Negative Pressure in Reservoir	-150 mmHg (-20 kPa)	

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LIMITED WARRANTY

TERUMO Cardiovascular Systems Corporation warrants that the product has been manufactured and packaged with reasonable care, and will be free from defects under normal and proper use. TERUMO Cardiovascular Systems Corporation will not be liable for any direct, indirect, incidental, special, general or consequential loss, damage, or expense from the use of this product. TERUMO Cardiovascular Systems Corporation's sole obligation under this warranty shall be to repair or replace, at its option, any products that it believes were defective at the time of shipment if notice thereof is received before the expiration date indicated on the package.

THIS WARRANTY IS THE SOLE AND ONLY WARRANTY OF TERUMO CARDIOVASCULAR SYSTEMS CORPORATION AND IS EXPRESSLY IN LIEU OF AND EXCLUDES ANY OTHER EXPRESSED OR IMPLIED WARRANTY OF ANY KIND, INCLUDING THE WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE. IN ADDITION, THE REMEDIES SET FORTH HEREIN ARE THE SOLE REMEDIES AVAILABLE TO ANY PERSON OR ENTITY. No agent, employee or representative of TERUMO Cardiovascular Systems Corporation has any authority to change, amend or alter the foregoing or assume or bind TERUMO Cardiovascular Systems Corporation to any additional liability or Responsibility in connection with this product.

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125 BLUE BALL ROAD, ELKTON, MD 21921, U.S.A.
TERUMO AUSTRALIA PTY LTD
Country of Origin: U.S.A