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RX NC SPRINTER®

Instructions For Use

RAPID EXCHANGE BALLOON DILATATION CATHETER

En. Explanation of possible symbols on product labeling
Ja. 製品ラベルに記されている記号の説明
Fr. Explication des symboles éventuellement présents sur l'étiquette
De. Bedeutung der möglicherweise auf der Produktbeschriftung auftretenden Symbole
It. Spiegazione dei possibili simboli presenti sulle etichette del prodotto

Nl. Verklaring van mogelijke symbolen op het etiket van het product
Pt. Explicação de símbolos que podem existir na etiqueta do produto
El. Επεξήγηση των μεθυσών συμβόλων στις ετικέτες του προϊόντος
Es. Explicación de los símbolos de las etiquetas del producto

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Medtronic

Manufacturer:
Medtronic, Inc.
710 Medtronic Parkway NE,
Minneapolis, MN55432, U.S.A.
www.medtronic.com

Manufactured AT/EC Authorised Representative:
Medtronic Ireland
Parkmore Business Park West,
Galway, Ireland
Tel.: +353-91-708000. Fax: +353-91-757524.

Made in: Medtronic Mexico
S. de R.L. de C.V.
Avenida Paseo del Cuapah 10510
Parque Industrial El Lago
Tijuana, B.C. 22570, Mexico
Tel: +52-664-903-1100. Fax: +52-664-903-1175

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Contents:

- One (1) Balloon Dilatation Catheter
- One (1) Flushing Cannula*
- One (1) Re-Fold Tool*
- One (1) Looper™ Device*
- * Components contained in separate pouch

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Warnings:

- For single patient, single procedure use only. Do NOT resterilize and/or reuse it, as this can potentially result in compromised device performance and increased risk of inappropriate resterilization and cross contamination.
- Do not use the catheter if its package has been opened or damaged.
- To reduce the potential for vessel damage the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.
- PTCA in patients who are not acceptable candidates for coronary artery bypass graft surgery requires careful consideration, including possible hemodynamic support during PTCA, as treatment of this patient population carries special risk.
- When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum as this can potentially result in damage to the vessel wall. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
- Balloon pressure should not exceed the rated burst pressure indicated on the package label for each balloon. The rated burst pressure is based on the results of in vitro testing. Use of a pressure monitoring device is recommended to prevent over pressurization.
- PTCA should only be performed at hospitals where emergency coronary artery bypass graft surgery can be quickly performed in the event of a potentially injurious or life-threatening complication.
- Use only the recommended balloon inflation medium. To prevent the possibility of an air embolus, never use air or any gaseous medium to inflate the balloon.
- Use the catheter prior to the "Use By" date (Expiration Date) specified on the package.

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- Death
- Acute myocardial infarction
- Total occlusion of the coronary artery or bypass graft
- Coronary vessel dissection, perforation, rupture or injury
- Restenosis of the dilated vessel
- Hemorrhage or hematomas
- Unstable angina
- Arrhythmias, including ventricular fibrillation
- Drug reactions, allergic reaction to contrast medium
- Hypertension
- Infection
- Coronary artery spasm
- Arteriovenous fistula
- Embolism

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- Suitable guide wire, refer to label claim.
- 20 cc syringe for balloon preparation.
- Suitable guiding catheter, refer to label claim.
- 10 cc or smaller syringe for manual dye injections.
- Appropriate inflation medium (e.g. 50:50 sterile mixture of a contrast medium and saline.)
- Pressure-indicating inflation device.
- Hemostasis valve.

Directions For Use:

- Preparation
 - Select an appropriate balloon catheter for the target vessel.
 - Remove the device from the sterile packaging.

Note: This device allows the guide wire lumen to be flushed in two(2) formats:

- With device contained within its hoop (Ref. In-Hoop Flushing section).
- With device removed from its hoop:
 - Remove the balloon protector (C) from around the balloon. Removing the balloon protector will also remove the balloon stylette.

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- Fill a 20cc syringe with 20cc of a saline solution.
- Remove the flushing cannula (E) from the accessories pouch (D) and, without removing the cover, attach directly to the syringe.
- Remove the cover of the flushing cannula holder leaving the flushing cannula attached to the syringe.
- Insert the flushing cannula into the guidewire entry port (exchange joint).
- Depress the syringe to flush the guide wire lumen until fluid emerges from the tip (distal end of the balloon catheter).
- Remove the flushing cannula from the guide wire entry port (exchange joint).

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2. Balloon Purging

Purge air from the catheter using a 20 cc syringe filled with 2 to 3 cc of the inflation medium with the balloon catheter pointing downward. Attach an inflation device to the balloon inflation port. Ensure that a meniscus of contrast medium is evident in both the catheter luer connector and the inflation device. Apply negative pressure with the inflation device. A flow of bubbles will move from the balloon catheter into the inflation device. Balloon purging is complete when the flow of bubbles stops. Do not attempt Pre-Inflation technique to purge the balloon lumen.

3. Insertion Technique

- Place the guiding catheter, with a hemostasis valve attached, in the orifice of the target coronary artery.
- Advance a guide wire through the guiding catheter to reach and cross the target lesion. Advance the distal tip of the balloon catheter over the proximal end of the guide wire. Ensure that the guide wire exits the balloon catheter through the guide wire exit port (exchange joint).
- The hemostasis valve should be gradually tightened to control back flow. Excessive valve tightening may affect balloon inflation/deflation time as well as movement of the guide wire.
- Track the balloon catheter over the wire to cross the lesion using the radiopaque marker(s) to locate the balloon across the lesion.

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

5. Removing the Catheter

- Apply negative pressure to the inflation device and confirm that the balloon is fully deflated.
- Withdraw the balloon catheter into the guiding catheter while preserving guide wire position and adjusting the hemostasis valve appropriately. Remove the balloon catheter from the system.
- After the deflated balloon catheter is withdrawn, it should be wiped clean with gauze soaked with sterile normal saline.
- Inspect the balloon catheter integrity.
- If reinserting the same balloon catheter, flush the guide wire lumen by using the flushing cannula to flush saline through the guide wire entry port (exchange joint) until fluid emerges from the tip (distal end of the balloon catheter). Prior to reinsertion, the balloon catheter should be wiped clean with gauze soaked with sterile normal saline. The balloon may be refolded using the Re-Fold Tool as described in the Re-Fold Tool Instructions for use.

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In-Hoop Flushing

The in-hoop flushing cannula (E) enables the guide wire to be flushed without removing the catheter from its packaging. Instructions for use:

- With the balloon catheter contained in its packaging, manually detach the blue tubing (B).
- Remove the balloon protector (C) from the catheter. Removing the balloon protector will also remove the balloon stylette.

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- Fill a 20cc syringe with 20cc of a saline solution.
- Remove the flushing cannula (E) from the accessories pouch (D) and, without removing the cover, attach directly to the syringe.
- Remove the cover of the flushing cannula holder leaving the flushing cannula attached to the syringe.
- Insert the flushing cannula into the guide wire entry port (exchange joint).
- Depress the syringe to flush the guide wire lumen until fluid emerges from the tip (distal end of the balloon catheter). Droplets will be visible in the hoop.
- Remove the flushing cannula from the guide wire entry port.

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Re-Fold Tool

The Re-Fold Tool (F) is an accessory component used to re-fold the catheter.


Instructions for use:

- Deflate the balloon by applying negative pressure to the inflation device and maintain under vacuum.
- Visually inspect the balloon to confirm deflation.
- Remove the Re-Fold Tool from the accessories pouch (D).
- Ensure the balloon protector is not on the straight end of the stylette.
- Carefully load the stylette back through the catheter and past the proximal end of the balloon.
- While holding the catheter just proximal to the Re-Fold Tool, pull the catheter into the Re-Fold Tool until the entire balloon is covered.
- Gently remove the Re-Fold Tool and stylette.
- Inspect the balloon for any potential damage. Discard the balloon catheter if there is any visible damage.

Looper™ Device
The LOOPER™ device (G) is an accessory component used to fasten the catheter to the syringe during use.

Instructions for use:

- Remove the LOOPER™ device (G) from the accessories pouch (D).
- Form the catheter into a single or double loop.
- Press down the wings of the LOOPER™ device, clipping section for retention of the hypotube.
- Release the wings of the LOOPER™ device, clipping section to allow the catheter to be released.

Disclaimer of Warranty

NOTE: ALTHOUGH THE CATHETER, HEREINAFTER REFERRED TO AS "PRODUCT," HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC INC. AND ITS AFFILIATES HAS NO CONTROL OVER CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC INC. AND ITS AFFILIATES, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESSED AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC INC. AND ITS AFFILIATES SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC INC. AND ITS AFFILIATES TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

References

The physician should consult recent literature on current medical practice on balloon dilatation, such as published by American College of Cardiology - American Heart Association.

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Patents:

Protected under one or more of the following patents:
U.S. Patents 7,001,358, 6,190,358 and 6,805,057. U.S. and foreign patents pending.

Compliance Chart:

The Compliance Chart is based on actual data from in vitro testing @ 37° C, rounded to two decimal places. See Page 29 for Compliance Chart.

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EC REP

Authorized Representative in the European Community
Medtronic Ireland
Parkmore Business Park West,
Galway, Ireland.
Tel: +353-91-708000 Fax: +353-91-757524

Australia
Medtronic Australasia Pty Ltd
97 Waterloo Road
North Ryde, NSW 2113
Australia

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

RBP

En	Do Not Exceed Rated Burst Pressure
Ja	最大拡張圧を超えてはいけません
Fr	Ne pas dépasser la pression à l'éclatement
De	Nicht überschreiten den maximalen Nennsturzdruck
It	Non superare la pressione nominale di scoppio
Nl	Verscheidings banddruk niet overschrijden
Pl	Nie przekraczać ciśnienia nominalnego
Da	Do ikke overstige tryk
Es	No superar las valores de presión de ruptura
Fi	Älä ylittää nimellistä purkauspainetta
Cs	Nepřekračujte stanovený tlak prasknutí
El	Ανεβάζετε τον πίεση πέρα από το ονομαστικό
Hu	A levegő viszkozitási nyomaték nem szabad túllépni
Lv	Nepārsniedziet nominālo plūsmas spiedienu
Li	Nie przekraczać nominalnego ciśnienia
Pl	Nie przekraczać znamionowego ciśnienia rozrywania
Sk	Nepřekračujte maximální tlak prasknutí
Tr	Anma Patlama Basıncını Aşmayın
Zh	請勿超過額定爆破壓

En	Maximum Guidewire Diameter
Ja	最大導線ワイヤー径
Fr	Diamètre maximum du fil-guide
De	Maximaler Führungsdraht Durchmesser
It	Diametro massimo del filo guida
Nl	Maximale draaddiameter
Pl	Długość maksymalna przewodnika
Da	Maximum diameter of the guide
Es	Máximo diámetro del alfiler
Fi	Ennen käyttöä tarkistettava
Cs	Maximální průměr vodiče drátu
El	Αντικείμενο μέγιστης διαμέτρου
Hu	Maximális vezeték átmérő
Lv	Katēra vadītāja maksimālā diametrs
Li	Maximale draaddiameter
Pl	Maxymalny promień przewodnika
Sk	Maximálny priemer vodiča drátu
Tr	Maximum Katutuz Tel Çapı
Zh	最大導線直徑

USA

En	For U.S. Audiences Only
Ja	米国のクライアントのみに限定
Fr	Uniquement pour les États-Unis
De	Gilt nur für Leser in den USA
It	Esplicitamente per il mercato statunitense
Nl	Alleen van toepassing voor de VS
Pl	Aplikuje się na Stany Zjednoczone
Da	Må kun bruges i USA
Es	Solo para EE. UU.
Fi	Käyttökäyttöön vain Yhdysvalloissa
Cs	Pouze pro uživatele z USA
El	Αντικείμενο μόνο για τους Ηνωμένους Πολιτείες
Hu	Kizárólag az Egyesült Államok területén
Lv	Lietuvai tik naudoti
Li	Tik naudoti JAV
Pl	Tylko dla odbiorców w USA
Sk	Len pre používateľov v USA
Tr	Yalnızca ABD'de Kullanılır
Zh	仅限于美国用户

En	For Post-Dilatation Of Balloon Expandable Stents
Ja	バルーン拡張式ステントの拡張後に使用のこと
Fr	Pour la post-dilatation de stents expansibles par ballonnet
De	Zur nachträglichen Dilatation ballonexpandierbarer Stents
It	Per espansione post-impiego di stent espandibili tramite palloncino
Nl	Voor postdilatatie van ballon-expandeerbare stents
Pl	Para pós-dilatação de stents expansíveis por balão
Da	For post-dilatation af ballonudvidelige stenter
Es	Para post-dilatación de stents expandibles por balón
Fi	Post-dilataatio
Cs	Pro postdilataci stentů
El	Αντικείμενο post-dilatation
Hu	Balloonra történő dilatáció
Lv	Ar balonu izplešmā stenta izplešana pēc ievietošanas
Li	Vlaanderen: Balloon Dilatatie
Pl	Do balonowego rozprężania przedmiotowych stentów, które są przystosowane do takiego rozszerzania
Sk	Pri stenty rozšírení dilatáciou balónika
Tr	Balonla Genişletilebilir Stentlere Dilatasyon Sonrası İyileştirme
Zh	用于球囊扩张式支架的后续扩张

En	Length
Ja	長さ
Fr	Longueur
De	Länge
It	Lunghezza
Nl	Langte
Pl	Comprimento
Da	Length
Es	Longitud
Fi	Pituus
Cs	Délka
El	Μήκος
Hu	Hossz
Lv	Garums
Li	leng
Pl	Długość
Sk	Džūms
Tr	Uzunluk
Zh	长度

En	Do Not Use If Package Is Damaged
Ja	包装が破損している場合は、使用しないでください
Fr	Né pas utiliser si l'emballage est endommagé
De	Nicht verwenden, wenn die Verpackung beschädigt ist
It	Non utilizzare se l'imballaggio non è integro
Nl	Niet gebruiken als de verpakking beschadigd is
Pl	Nie używać, jeśli opakowanie jest uszkodzone
Da	Do not use if the package is damaged
Es	No utilizar si el envase está dañado
Fi	Älä käyttää, jos pakkaus on vaurioitunut
Cs	Nepoužívejte, jest-li obal poškozen
El	Μην χρησιμοποιείτε, εάν η συσκευασία έχει υποστεί ζημία
Hu	Nem szabad használni, ha a csomagolás sérült
Lv	Nie uzturēt, ja iepakojums ir bojāts
Li	Niet gebruiken, als de verpakking beschadigd is
Pl	Nie używać, jeśli opakowanie jest uszkodzone
Sk	Nepoužívať, ak je obal poškodený
Tr	Ambarın hasarlıysa kullanmayın
Zh	如有包装破损，请勿使用

En	Consult Instructions For Use
Ja	取扱説明書を参照してください
Fr	Consulter le mode d'emploi
De	Gebrauchsanweisung/Bedienungsanleitung beachten
It	Consultare le istruzioni per l'uso
Nl	Read the gebruiksaanwijzing
Pl	Consultar instruções de utilização
Da	Tilgå de brugsanvisninger
Es	Consultar instrucciones de uso
Fi	Katso käyttöohjeet
Cs	Viz návody k použití
El	Τίθεται η χρήση του προϊόντος
Hu	Ismerje a használati útmutatót
Lv	Izlasīt lietošanas pamācību
Li	Zie de gebruiksaanwijzing
Pl	Zapoznać się z instrukcją użytkownika
Sk	Přečítajte si pokyny pre používanie
Tr	Kullanma Talimatlarını Okuyun
Zh	參閱“使用說明”

En	Do Not Reuse
Ja	再使用しないでください
Fr	A usage unique
De	Nicht wiederverwenden
It	Non riutilizzare
Nl	Niet opnieuw gebruiken
Pl	Nie reużywać
Da	Do not reuse
Es	No reutilizar
Fi	Älä käyttää uudelleen
Cs	Nepoužívejte opakovaně
El	Οχι επαναχρησιμοποίηση
Hu	Nem újrafelhasználható
Lv	Atkārtotai lietošanai nav paredzēts
Li	Niet opnieuw gebruiken
Pl	Nie używać ponownie
Sk	Nepoužívať opakovaně
Tr	Yeniden kullanılmayın
Zh	一次使用

En	Use By
Ja	使用期限
Fr	Date de péremption
De	Zu verwenden bis einschließlich
It	Utilizzare entro
Nl	Teg gebruiken tot en met
Pl	Nie używać dłużej niż
Da	Benyttes inden
Es	No utilizar después de
Fi	Käyttökäyttöä varten
Cs	Do expirace
El	Χρησιμοποιήστε πριν από την
Hu	Használati határidőig
Lv	Derīgt līdz
Li	Naudoti iki
Pl	Termin ważności
Sk	Deizm nepākoties pabeigšanas
Tr	Son kullanma tarihi
Zh	有效期

LOT

En	Lot Number
Ja	ロット番号
Fr	Numero de lot
De	Losnummer
It	Numero di lotto
Nl	Partijnummer
Pl	Numero di lotto
Da	Apklausos partijos
Es	Numero de lote
Fi	Partinumeri
Cs	Číslo šarže
Hu	Gyártási szám
Li	Partijnummer
Pl	Partijski numer
Sk	Číslo šarže
Tr	Lot Numarası
Zh	批号

REF

En	Catalog Number
Ja	製品コード
Fr	Référence
De	Katalognummer
It	Numero di catalogo
Nl	Catalogusnummer
Pl	Numer katalogowy
Da	Apklausos katalogo
Es	Numero de referencia
Fi	Kataloginumeri
Cs	Katalogové číslo
El	Katalogos nomos
Hu	Katalógusszám
Li	Katalognummer
Pl	Numer katalogowy
Sk	Katalogový číslo
Tr	Katalog Numarası
Zh	产品目录号

En	Manufacturer
Ja	製造元
Fr	Fabricant
De	Hersteller
It	Fabbricante
Nl	Fabrikant
Pl	Fabrykant
Da	Udbydere
Es	Fabricante
Fi	Fabrikant
Cs	Produkční
El	Παραγωγός
Hu	Herkez
Li	Fabrikant
Pl	Fabrykant
Sk	Produkční
Tr	Üretici
Zh	製造商

STERILE R

En	Sterilized Using Irradiation
Ja	放射線滅菌済み
Fr	Sterilisé par irradiation
De	Sterilisiert mittels Bestrahlung
It	Sterilizzato con radiazioni ionizzanti
Nl	Gesteriliseerd met straling
Pl	Esteryzowany metodą promieniowania
El	Αποστειρωμένο με ακτινοβολία
Es	Esterilizado mediante radiación
De	Sterilisiert mit Bestrahlung
Sw	Steriliserad med strålning
Fi	Steriloitu säteilytyksellä
Cs	Sterilizováno ozařením
Et	Sterileeritud kiirgusega
Hu	Besugáztatással sterilizált
Lv	Sterilizēts, izstarojot apstarošanu
Li	Sterilizzata radioactiva strahlung
Pl	Sterylizowany metodą promieniowania
Sk	Sterilizovaný žiarením
Tr	İrinyasyonla Sterilize Edilmiş
Zh	使用射线灭菌

En	Date of Manufacture
Ja	製造年月日
Fr	Date de fabrication
De	Herstellungsdatum
It	Data di fabbricazione
Nl	Productiedatum
Pl	Data produkcji
El	Ημερομηνία κατασκευής
Es	Fecha de fabricación
De	Fabricationsdatum
Sw	Tillveringsdatum
Fi	Välistuotuspäivä
Cs	Datum výroby
Et	Tootmisupäev
Hu	A gyártás dátuma
Lv	Izgatavotības datums
Li	Produktionsdate
Pl	Data produkcji
Sk	Datum výroby
Tr	İmali Tarihi
Zh	生产日期

EC REP

En	Authorized Representative in the European Community
Ja	EC 認定代理店
Fr	Représentant autorisé dans la Communauté européenne
De	Autorisierter Repräsentant in der Europäischen Gemeinschaft
It	Rappresentante autorizzato nella Comunità europea
Nl	Ovetoordenaar/vertrouwsoordnager in de Europese Gemeenschap
Pl	Reprezentantę autoryzowanego na Wspólnocie Europejskiej
El	Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα
Es	Representante autorizado en la Comunidad Europea
De	Autorisiert repräsentant in der EG
Sw	Autoriserad representant inom EG
Fi	Valtuutettu edustaja Euroopan yhteisön alueella
Cs	Ověřovaný zástupce pro Evropské společenství
El	Επίσημος εκπρόσωπος Ευρωπαϊκής Ομοσπονδίας
Hu	Hivatalos képviselő az Európai Közösségekben
Lv	Pārvietotais pārstāvis Eiropas Kopienā
Li	gegründeter abhänger Europäischer Gemeinschaft
Pl	Autorizowany przedstawiciel we Wspólnocie Europejskiej
Sk	Autorizovaný zástupce v Evropské společnosti
Tr	Ayınca Toprakçılık Yönetim Tesisleri
Zh	歐共體的授權代表

En	Manufactured in
Ja	製造場所
Fr	Lieu de fabrication
De	Hergestellt in
It	Fabbricato in
Nl	Vervaardigd in
Pl	Fabryczno w
El	Κατασκευασμένο σε
Es	Fabricado en
De	Fremstillet i
Sw	Tillverkad i
Fi	Välistuotettu
Cs	Vyrobno v
El	Tootamata
Hu	A gyártás helye
Lv	Rādītājs
Li	Produktionsort
Pl	Miejsce produkcji
Sk	Miesto výroby
Tr	İmali Yeri
Zh	製造地

En	Rapid Exchange Balloon Dilatation Catheter
Ja	迅速交換バルーン拡張カテーテル
Fr	Cathéter de dilatation à ballonnet à échange rapide
De	Rapid Exchange-Balloon Dilatationskatheter
It	Catetere di dilatazione a palloncino a scambio rapido
Nl	Rapid Exchange-balloon dilatatorkatheter
Pl	Katheter do dilatacji balonowej typu szybkiego wymiany
El	Καθετήρας μπαλονότυπου διαστολής, γρήγορης αντικατάστασης
Es	Catéter de balón de dilatación de intercambio rápido
De	Balloon Dilatationskatheter st. typ rapid exchange
Sw	Baljongdilataionskatheter av REX-typ (Rapid Exchange)
Fi	Rapid exchange -pallondilataattikathetri

Cs	Balóntíróó dilataciós katéter gyors cserével
El	Καθετήρας μπαλονότυπου διαστολής
Hu	Gyorsan cserélhető dilatációs ballonkatéter
Lv	Ātrīkšanas baloniņa katērs ar ātru nomaiņu
Li	Grafik (pallong) dilatacijs katērs ar ātru nomaiņu
Pl	Cewnik balonowy rozszerzający przystosowany do szybkiej wymiany (Rapid Exchange)
Sk	Rýchle vymeniteľný balónový dilatovaný katéter
Tr	Hızlı Değişim Balon Dilatasyon Kateteri
Zh	快速交換球囊扩张导管

ENGLISH

Medtronic Balloon Catheter Description:

The NC Sprinter Rapid Exchange Balloon Dilatation Catheter is a Percutaneous Transluminal Coronary Angioplasty (PTCA) Rapid Exchange System. The balloon at the distal end of the catheter can be inflated to a defined diameter at a specific pressure (see labeling). The proximal end of the catheter has a female luer for attachment to an inflation device. The catheter provides a lumen which enables the use of a guide wire to position the catheter. Radiopaque balloon marker(s) enable accurate placement. Shaft markers for brachial and femoral techniques are in place.

Contents:

- One (1) Balloon Dilatation Catheter
- One (1) Flushing Cannula*
- One (1) Re-Fold Tool*
- One (1) Looper™ Device*

* Components contained in separate pouch

Indications:

- The balloon dilatation catheter is indicated for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis for the purpose of improving myocardial perfusion.
- The balloon dilatation catheter (balloon models 2.25 - 5.0mm) is also indicated for post-delivery expansion of balloon expandable stents.

Note: Bench testing was conducted with the NC Sprinter Rapid Exchange Balloon Dilatation Catheter and marketed balloon expandable stents. Consideration should be taken when this device is used with different manufacturers' stents due to difference in stent design.

Contra-Indications:

- The catheter is contraindicated for use in:
 - Unprotected left main coronary artery.
 - Coronary artery spasm in the absence of significant stenosis.

Warnings:

- For single patient, single procedure use only. Do NOT resterilize and/or reuse it, as this can potentially result in compromised device performance and increased risk of inappropriate sterilization and cross contamination.
- Do not use the catheter if its package has been opened or damaged.
- To reduce the potential for vessel damage the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.
- PTCA in patients who are not acceptable candidates for coronary artery bypass graft surgery requires careful consideration, including possible hemodynamic support during PTCA, as treatment of this patient population carries special risk.
- When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum as this can potentially result in damage to the vessel wall. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
- Balloon pressure should not exceed the rated burst pressure indicated on the package label for each balloon. The rated burst pressure is based on the results of in vitro testing. Use of a pressure monitoring device is recommended to prevent over pressurization.
- PTCA should only be performed at hospitals where emergency coronary artery bypass graft surgery can be quickly performed in the event of a potentially injurious or life-threatening complication.
- Use only the recommended balloon inflation medium. To prevent the possibility of an air embolus, never use air or any gaseous medium to inflate the balloon.
- Use the catheter prior to the "Use By" date (Expiration Date) specified on the package.

Precautions:

- Prior to angioplasty, the catheter should be examined to verify functionality and ensure that its size and shape are suitable for the procedure for which it is to be used.
- If the device is kinked, it should not be used.
- The catheter system should be used only by physicians trained in the performance of percutaneous transluminal coronary angioplasty.

- Appropriate anticoagulation, antiplatelet and vasodilator therapy should be administered to the patient.
- When using two guide wires, care should be taken when introducing, torquing and removing one or both guide wires to avoid entanglement. It is recommended that one guide wire be completely withdrawn from the patient before removing any additional equipment.
- Care should be taken not to apply excessive force during preparation or use, as this may damage the device.

Adverse Effects:

Possible adverse effects include, but are not limited to the following:

- Death
- Acute myocardial infarction
- Total occlusion of the coronary artery or bypass graft
- Coronary vessel dissection, perforation, rupture or injury
- Restenosis of the dilated vessel
- Hemorrhage or hematoma
- Unstable angina
- Arrhythmias, including ventricular fibrillation
- Drug reactions, allergic reaction to contrast medium
- Hypertension
- Infection
- Coronary artery spasm
- Arteriovenous fistula
- Embolism

Materials To Be Used In Combination With A Balloon Catheter Include:

- Suitable guide wire, refer to label claim.
- 20 cc syringe for balloon preparation.
- Suitable guiding catheter, refer to label claim.
- 10 cc or smaller syringe for manual dye injections.
- Appropriate inflation medium (e.g. 50:50 sterile mixture of a contrast medium and saline).
- Pressure-indicating inflation device
- Hemostasis valve.

Directions For Use:

1. Preparation

- Select an appropriate balloon catheter for the target vessel.
 - Remove the device from the sterile packaging.
- Note: This device allows the guide wire lumen to be flushed in two(2) formats:
- With device contained within its hoop (Ref. In-Hoop Flushing section).
 - With device removed from its hoop.
- Remove the balloon protector (C) from around the balloon. Removing the balloon protector will also remove the balloon stylette.

- Note: Do not discard the balloon stylette until the end of the angioplasty procedure as it may be required for balloon refolding.
- Fill a 20cc syringe with 20cc of a saline solution.
 - Remove the flushing cannula (E) from the accessories pouch (D) and, without removing the cover, attach directly to the syringe.
 - Remove the cover of the flushing cannula holder leaving the flushing cannula attached to the syringe.
 - Insert the flushing cannula into the guidewire entry port (exchange joint).
 - Depress the syringe to flush the guide wire lumen until fluid emerges from the tip (distal end) of the balloon catheter.
 - Remove the flushing cannula from the guide wire entry port (exchange joint).

Note: Do not discard the flushing cannula until the end of the angioplasty procedure as additional flushing may be required.

2. Balloon Purging

Purge air from the catheter using a 20 cc syringe filled with 2 to 3 cc of the inflation medium with the balloon catheter pointing downward. Attach an inflation device to the balloon inflation port. Ensure that a minimum of contrast medium is evident in both the catheter luer connector and the inflation device. Apply negative pressure with the inflation device. A few of bubbles will move from the balloon catheter into the inflation device. Balloon purging is complete when the flow of bubbles stops. Do not attempt Pre-Inflation technique to purge the balloon lumen.

3. Insertion Technique

- Place the guiding catheter, with a hemostasis valve attached, in the orifice of the target coronary artery.
- Advance a guide wire through the guiding catheter to reach and cross the target lesion. Advance the distal tip of the balloon catheter over the proximal end of the guide wire. Ensure that the guide wire exits the balloon catheter through the guide wire exit port (exchange joint).
- The hemostasis valve should be gradually lightened to control back flow. Excessive valve lightening may affect balloon inflation/deflation time as well as movement of the guide wire.
- Track the balloon catheter over the wire to cross the lesion using the radiopaque marker(s) to locate the balloon across the lesion.


4. Balloon Inflation

- Inflate the balloon to dilate the lesion using standard PTCA techniques.
- After each subsequent inflation, the distal blood flow should be assessed.
- If a significant stenosis persists, successive inflations may be required to resolve the stenosis. Do not exceed the rated burst pressure (see labeling).
- Confirm the results with fluoroscopy.

5. Removing the Catheter

- Apply negative pressure to the inflation device and confirm that the balloon is fully deflated.
- Withdraw the balloon catheter into the guiding catheter while preserving guide wire position and adjusting the hemostasis valve appropriately. Remove the balloon catheter from the system.
- After the deflated balloon catheter is withdrawn, it should be wiped clean with gauze soaked with sterile normal saline.
- Inspect the balloon catheter integrity.
- If reinserting the same balloon catheter, flush the guide wire lumen by using the flushing cannula to flush saline through the guide wire entry port (exchange joint) until fluid emerges from the tip (distal end) of the balloon catheter. Prior to reinsertion, the balloon catheter should be wiped clean with gauze soaked with sterile normal saline. The balloon may be refolded using the Re-Fold Tool as described in the Re-Fold Tool Instructions for use.

Note: Do not discard the balloon stylette until the end of the angioplasty procedure as it may be required for balloon refolding.

In-Hoop Flushing

The in-hoop flushing cannula (E) enables the guidewire lumen to be flushed without removing the catheter from its packaging hoop.

Instructions for use:

- With the balloon catheter contained within its hoop (A) manually detach the blue tubing (B).
- Remove the balloon protector (C) from around the balloon. Removing the balloon protector will also remove the balloon stylette.

Note: Do not discard the balloon stylette until the end of the angioplasty procedure as it may be required for balloon refolding.

- Fill a 20cc syringe with 20cc of a saline solution.
- Remove the flushing cannula from the accessories pouch (D) and, without removing the cover, attach directly to the syringe.
- Remove the cover of the flushing cannula holder leaving the flushing cannula attached to the syringe.
- Insert the flushing cannula into the tip (distal end) of the balloon catheter.
- Depress the syringe to flush the guide wire lumen, until fluid emerges from the guide wire entry port (exchange joint). Droplets will be visible in the hoop.
- Remove the flushing cannula from the distal tip of the balloon.

Note: Do not discard the flushing cannula until the end of the angioplasty procedure as additional flushing may be required.

Re-Fold Tool

This Re-Fold Tool (F) is an accessory component that allows the balloon to be re-folded if required.

Instructions for use:


- Deflate the balloon by applying negative pressure to the inflation device and maintain under vacuum.
- Visually inspect the balloon to confirm that it is fully deflated.
- Remove the Re-Fold Tool from the accessories pouch (D).
- Ensure the balloon protector is not on the balloon stylette and load the non-flared end of the Re-Fold Tool onto the straight end of the stylette.
- Carefully load the stylette back through the distal tip of the catheter and past the proximal end of the balloon.
- While holding the catheter just proximal to the balloon, push the re-fold device over the balloon in a gentle twisting motion until the entire balloon is covered.
- Gently remove the Re-Fold Tool and stylette.
- Inspect the balloon for any potential damage. Discard the balloon catheter if there is any visible damage present on the balloon.

Looper™ Device

The LOOPER™ device (G) is an accessory component for use with Medtronic Rapid Exchange balloon catheters. The device allows the catheter to be fastened in a coiled configuration for ease of handling during use.

Instructions for use:

- Remove the LOOPER™ device from the accessories pouch (D).
- Form the catheter into a single or double loop when required.
- Press down the wings of the LOOPER™ device to open the clipping section for retention of the hypotube (H).
- Release the wings of the LOOPER™ device with the hypotube inside the clipping section to secure the catheter in a coiled configuration.

Disclaimer of Warranty

NOTE: ALTHOUGH THE CATHETER, THEREAFTER REFERRED TO AS "PRODUCT," HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC INC. AND ITS AFFILIATES HAS NO CONTROL OVER CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC INC. AND ITS AFFILIATES, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESSED AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC INC. AND ITS AFFILIATES SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC INC. AND ITS AFFILIATES TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

References

The physician should consult recent literature on current medical practice on balloon dilatation, such as published by American College of Cardiology / American Heart Association.

The exclusion and limitations set out above are not intended to and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court or competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected.

Patents:

Protected under one or more of the following patents:
U.S. Patent 7,001,358.

Compliance Chart:

The Compliance Chart is based on actual data from in vitro testing @ 37°C, rounded to two decimal places. See Page 30 for Compliance Chart.

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