DePuy Synthes

MITEK SPORTS MEDICINE

COMPANIES OF *Johnson & Johnson*

Latarjet EXPERIENCE


Instability Shoulder

System Instruments

DESCRIPTION

LATARJET EXPERIENCE INSTRUMENTS

Important: Please read these instructions for use, and the corresponding surgical techniques carefully before use. Ensure that you are familiar with the appropriate surgical technique.

The instruments for the Latarjet EXPERIENCE Instability Shoulder System consist of sterile and non-sterile instruments. These instruments are to be used in two specific procedures: "Latarjet" and "Bristow". The instruments will be supplied in dedicated kits.

Sterile disposable instruments:

Contents are **STERILE** unless packaging is damaged or open.

Cat. #	Description	Quantity
288238	Latarjet Procedure Kit Components:	
	Glenoid K-wire (380mm) - sterile	3
	Coracoid Positioning Double Cannula	1
	Cannula Obturator	2
	Plug	2
288239	Bristow Procedure Kit Components:	
	Glenoid K-wire (380mm) - sterile	2
	Coracoid Positioning Single Cannula	1
	Cannula Obturator	1
	Plug	1

Non-sterile reusable instruments:

Cat. #	Description	Quantity	
		Latarjet	Bristow
288205	Subscapularis Channeler	1	1
288240	Latarjet, Switching Stick	1	1
2882101	Latarjet, Small Switching Stick	1	1
2882102	Latarjet, Coracoid Osteotome	1	1
288209	Coracoid Screw 3.5mm	2	1
288241	Latarjet, Glenoid Drill	1	1
288210	Reference Tube	1	1
288211	Combo Screwdriver	1	1
288216	Extraction Driver	1	1

Cat. #	Description	Quantity	
		Latarjet	Bristow
288200	Inline Coracoid Drill Guide	1	1
288221	Glenoid Guide	1	1
288206	Coracoid Offset Pin	1	1
288202	Coracoid Top Hat Drill	1	1
288201	Coracoid Drill	1	1
288203	Top Hat Tap	1	1
288208	Glenoid Offset Pin	1	1
288215	Subscapularis Spreader	1	1
288213	Latarjet, Hohmann Retractor	2	2
288214	Glenoid Lever, Double Prong	1	1
288212	Latarjet Humeral Head Retractor	1	1
288207	Glenoid K-wire (380mm)	1	1
288242	Latarjet, Small Probe	1	1

INDICATIONS

The Latarjet EXPERIENCE Instability Shoulder Set is intended to provide to the orthopedic surgeon with a means of bone fixation and to assist in the management of reconstructive surgeries.

CONTRAINDICATIONS

Not for use in conditions that would retard healing, (e.g. blood supply limitation, infection), and conditions which would tend to compromise the patient's ability to heal, such as cognitive state.

WARNINGS

- Read these instructions prior to use of the Latarjet EXPERIENCE Instability Shoulder Set.
- Do not use if the label is missing. It is imperative that the surgeon and operating room staff are fully conversant with the appropriate surgical technique for the use of the Latarjet EXPERIENCE Instability Shoulder Set.
- Some instruments are supplied sterile and are for single use only (reference product list under Description above). They have not been designed to be re-used or re-sterilized. Reprocessing may lead to changes in material characteristics such as corrosion and dulled edges, which may impact the strength of the devices and compromise device performance. Reprocessing of single-use devices can also cause cross-contamination leading to patient infection. These risks may potentially affect patient safety.
- Do not use if sterile packaging appears to be damaged.
- Some instruments are supplied non-sterile (reference product list under Description above). Clean and disinfect the instruments after each procedure as described in the Cleaning/Disinfection/ Sterilization section below; sterilize prior to the first use and every subsequent use.
- The K-wires and Drills in the Set have sharp cutting edges that may puncture skin. Handle with care.

PRECAUTIONS

- Inspect the instruments before use for integrity and compatibility in order to ensure proper functionality and safety.
- Replace all dull instruments.
- Highly alkaline conditions can damage products with aluminum parts.
- Instruments are not compatible with hypochlorite solutions, as these will promote corrosion.

CLEANING AND DISINFECTION OF NON-STERILE INSTRUMENTS

Disassembly:

- Disassemble before cleaning and sterilization.

Warnings:

- Wear appropriate protective equipment: eye protection, gloves, etc.
- Cleaning should be performed immediately after each procedure, before blood, saline and debris are dry.
- As disinfection alone is not adequate, instruments should be sterilized before each procedure.
- Prior to sterilization, instruments should be cleaned by either manual or automated means described below:

Manual Cleaning

1. Rinse with water that has been purified (e.g. demineralization, deionization, reverse osmosis, distillation).
2. Prepare an enzymatic cleaning solution according to the manufacturer's instructions.
3. Submerge the devices in the enzymatic solution for 10 minutes, scrub the submerged instruments with a soft sponge and agitate. Use a pipe cleaner for lumens, crevices and slots.

NOTE: Ultrasonic cleaning will provide the most efficient cleaning.

- Rinse with distilled or demineralized water, agitating for 1 minute; repeat rinse twice.
- Dry devices thoroughly before storage.

Automated Cleaning

1. Perform the pre-cleaning step by flushing an enzymatic detergent through the lumens using a 5mL syringe and needle and using a soft bristled brush to scrub the crevices and slots. Rinse the device until no visible detergent is present then clean using the automated process below:
 - Wash – Cold water, temperature 80°F (27°C) for 35 seconds
 - Wash – Hot water, temperature minimum 140°F (60°C) for 2 minutes
 - Wash – Cold water, temperature 80°F (27°C) for 45 seconds
 - Wash – Hot water, temperature minimum 140°F (60°C) for 2 minutes
 - Thermal Rinse – Hot water, temperature 180°F (82°C) for 1 minute
 - Drying – Hot air, temperature minimum 180°F, (82°C) minimum 6 minutes
2. Inspect the instruments after each cleaning cycle to ensure proper functionality function and safety. If you suspect a problem, set aside the Instrument for repair.
3. After cleaning, disinfection may be performed by immersion in a disinfection solution such as Activated Glutaraldehyde followed by rinsing and drying.

STERILIZATION OF NON-STERILE INSTRUMENTS

Sterilization Warnings

Sterilization is not a substitute for cleaning. Devices must be thoroughly cleaned prior to sterilization.

- Prior to surgical use, the instruments should be cleaned as described below.
- Wrap instruments in accordance with the wrap manufacturer's instructions for use. Local procedures, using standard wrapping techniques such as those described in the current revision of ANSI/AAMI ST46 can be also referenced. In the United States (US), use an FDA-cleared Sterilization Wrap.
- Store in a cool, dry place.
- For reusable instruments that are provided in sets, use the Latarjet sterilization tray for containment during sterilization. Place the instrument in the appropriate location within the sterilization tray. If a specific location is not identified for an instrument, it may be placed in the general purpose area ensuring that the sterilant reaches to all surfaces including difficult to reach areas, lumens, etc. Use a validated, properly maintained and calibrated steam sterilizer.
- The total combined sterilization tray system (devices, internal trays, outer sterilization case and lid) should not exceed 10 kg (22 lbs) weight limit.
- Do not stack tray systems for sterilization.

The following methods and parameters have been validated:

1. Prevacuum (Dynamic Air Removal) Steam Sterilization:

Sterilizer type: Prevacuum/Porous Load

Effective steam sterilization can be achieved using the following cycles:

Cycle Type	Temperature	Minimum Exposure Time
Prevacuum	132°C (270°F)	4 min
Prevacuum (outside the U.S.)	134–137°C	3 min
Drying Time: 30 minutes		

2. Gravity -Displacement Steam Sterilization (for instruments only):

Sterilizer type: Gravity

Cycle Type	Temperature	Minimum Exposure Time
Gravity Displacement	132°C (270°F)	10 minutes
Drying Time: 40 minutes		

Instrument Configuration: Individually Wrapped

INSTRUCTIONS FOR USE: LATARJET PROCEDURE

(Bone graft transplantation, fixed with two cortical screws)

Laterjet procedure steps:

1. Joint evaluation and surgical final decision
2. Shoulder access, exposure and preparation
3. Split of *Subscapularis*
4. Bone graft preparation, drill 2 holes
5. Bone graft osteotomy
6. Glenoid drilling
7. Bone graft attachment and transfer
8. Bone graft -Glenoid final fixation

Warning: Care should be exercised; over-tightening to high torque levels can cause damage to the bone and may lead to compromised holding strength.

INSTRUCTIONS FOR USE: BRISTOW PROCEDURE

(Bone graft transplantation, fixed with one cortical screw)

Bristow Procedure steps:

1. Joint evaluation and surgical final decision
2. Shoulder access, exposure and preparation
3. Split of *Subscapularis*
4. Bone graft preparation, drill 1 hole
5. Bone graft osteotomy
6. Glenoid drilling
7. Bone graft attachment and transfer
8. Bone graft -Glenoid final fixation

Warning: Care should be exercised; over-tightening to high torque levels can cause damage to the bone and may lead to compromised holding strength.

HOW SUPPLIED

Instruments are supplied as follows:

- Reusable Instruments are supplied clean and must be sterilized prior to use.
- Disposable Instruments are supplied sterile and need to be properly disposed of, according to the local disposal procedure.

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DePuy Mitek, Inc.