Contour VL™ SureDrive™

Steerable Ureteral Stent Set

Directions for Use	2
Instrucciones de uso	6
Mode d'emploi	10
Gebrauchsanweisung	14
Istruzioni per l'uso	18
Gebruiksaanwijzing	22
Instruções de Utilização	26

50502707-01

2019-10

Contour VL™ SureDrive™

Steerable Ureteral Stent Set

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

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DEVICE DESCRIPTION

The Contour Variable Length Stent is designed to be placed over a guidewire to help drain urine from the kidney into the bladder. The coils on both distal and proximal pigtails allow one stent to be used in adult ureters of various lengths. This device is coated with exclusive HydroPlus Hydrophilic coating which, when activated, significantly increases the lubricous characteristics of the device. Care should be exercised when handling this coated device (See "Prior to Placement" and "Activate Coating").

The stent is MR safe.

Carefully read **all** the instructions prior to using this device.

Contents

The Contour VL Variable Length Steerable Ureteral Stent Set includes:

- (1) Contour VL Stent (Soft Percuflex Material) with HydroPlus™ Coating
- (1) Compatible Steerable Stent Positioner

INTENDED USE/INDICATIONS FOR USE

The Ureteral Stent is intended to facilitate drainage from the kidney to the bladder via placement endoscopically or fluoroscopically by a trained physician.

CONTRAINDICATIONS

The use of ureteral stents should be reconsidered if the following conditions exist:

- Unrepaired ureteral avulsion

WARNINGS

None known.

PRECAUTIONS

1. Retrieval line indwelling time should not exceed fourteen (14) days to avoid possible line encrustation.
2. This device is supplied sterile and recommended for one-time use **only**.
3. Avoid bending or kinking the stent during or prior to placement as to do so could damage the integrity of the stent.
4. If resistance is encountered during advancement or withdrawal of the stent, **STOP**. Do not continue without first determining the cause of resistance and taking remedial action.
5. Periodic upper urinary tract examinations are recommended to ensure the stent is draining adequately. Evaluation for stent encrustation should be performed periodically.

Note: Where long-term use is indicated, it is recommended that indwelling time for stent (with retrieval line removed) not exceed 365 days*. This stent should be evaluated by the physician on or before 90 days post-placement.

*Biocompatibility data on file.

6. This stent is not intended to be a permanent implant device.
7. The recommendations given herein are meant to serve only as a basic guide to the utilization of this device. Use of this device should be restricted to physician specialists with comprehensive knowledge of the indications, techniques and risks of the procedure. The references given below provide a broad overview of the subject of ureteral stent insertion over a preplaced guidewire.
8. Stent placement may complicate the subsequent work-up/evaluation of unexplained hematuria.
9. Patients who are poor surgical risk should be managed appropriately.
10. Formation of knots in multi-length ureteral stents may occur. This may result in injury to the ureter during removal and/or the need for additional surgical intervention. The presence of a knot should be considered if significant resistance is encountered during attempts at removal.

ADVERSE EVENTS

Adverse events associated with retrograde and antegrade positioned indwelling ureteral stents: Reflux-GU (e.g. ureteral reflux); Occlusion/Obstruction (e.g. catheter, stent); Migration (e.g. dislodgement); Hemorrhage; Infection (e.g. sepsis, peritonitis, urinary tract infection); Perforation (e.g. bladder, ureter, kidney, renal pelvis); Extravasation; Encrustation; Loss of renal function; Edema; Urinary symptoms (e.g. frequency, urgency, incontinence, dysuria, nocturia, hematuria); Pain/discomfort; Stent fragmentation; Fistula; Hydronephrosis; Stone formation; Tissue damage; Erosion.

HOW SUPPLIED

The Ureteral Stent and Steerable Positioner are supplied in sealed pouches and intended for single use only.

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Handling and Storage

Store in a cool, dry, dark place. Rotate inventory so that products are used prior to the sterilization expiration date on package label.

OPERATIONAL INSTRUCTIONS

Placement Prerequisites

1. The involved renal collecting system should be visualized via intravenous, retrograde or antegrade pyelography.
2. Fluoroscopy is recommended for more precise control of stent placement, however, standard radiography may be used.

Prior to Placement

Activate Coating

To activate the HydroPlus™ Coating, soak the stent in sterile water or saline for at least 30 seconds. Keep the stent wet during placement. Use a wet gauze pad to handle the stent during placement, if necessary. **DO NOT** wipe stent with dry gauze or any solvents as it may damage coating.

Recommended Technique

Note: Prior to stent placement, grasp retrieval line knot and pull to release retrieval line from under pigtail straightener (Figure 1).


Figure 1

Check the stent system (components)


Figure 2

Note: Pigtail straightener slides easily over ureteral stent (distal end of stent first) straightening out the pigtail and facilitating introduction of guidewire through stent. After assembly, pigtail straightener is removed.

A. Retrograde

1. Insert flexible end of the guidewire into cystoscope and pass up ureter into renal pelvis.
2. Pass the tapered tip of the stent over the guidewire and through cystoscope while assistant maintains guidewire position. Leave the bladder end of the stent exposed.
3. Load steerable positioner onto guidewire and advance until inner stylet of positioner is securely attached to the bladder end of the stent (Figure 3).


Figure 3

4. Advance stent up ureter with positioner until the thick black band on stent is located at the ureteral orifice. Adjust placement of the stent with steerable positioner (push, turn, pull).
5. While maintaining positioner in place, partially withdraw the guidewire allowing the renal coil to form in the renal pelvis.
6. Confirm stent position in bladder. Adjust stent placement with steerable positioner if necessary.
7. Partially withdraw the guidewire to allow a partial formation of the bladder coil.


Figure 4

8. Disengage the stent from the steerable positioner by rotating the thumb knob clockwise followed by pushing the knob forward (Figure 4). Carefully withdraw guidewire and then positioner.

B. Percutaneous Antegrade

1. Establish entry into involved renal pelvis using a percutaneous access set.
2. Pass flexible end of guidewire down ureter into bladder.
3. Pass appropriate open end ureteral catheter over guidewire to confirm entrance into bladder. Placement of a semi-rigid sheath, of appropriate size, into the access site is a useful adjunct at this point.
4. Pass the tapered tip of the stent over the guidewire. Load steerable positioner onto guidewire and securely attach to stent (Figure 3).
5. Advance stent down ureter with steerable positioner. Adjust stent placement with positioner if necessary.
6. Confirm lower coil to be in bladder and upper coil within renal pelvis.
7. Disengage the steerable positioner from the stent by rotating the thumb knob clockwise followed by pushing the knob forward (figure 4).
8. Remove retrieval line by holding the knot, cutting one strand, and while holding the knot, gently pull on the retrieval line. The guidewire is then removed, stabilizing the stent with steerable positioner.

Note: If a temporary percutaneous nephrostomy drainage tube is to be placed post stent insertion, leave the guidewire in place when removing the positioner.

Retrieval Line Management Options

- A. The line may be left in place for subsequent stent removal. To avoid possible retrieval line encrustation the stent should be left in place for a maximum of fourteen (14) days.
- B. The line may be used as an aid to positioning and then be removed once the stent is in the desired position.

After stent placement, but before the steerable positioner is removed, disengage the steerable positioner from the stent by rotating the

thumb knob clockwise followed by pushing the knob forward (figure 4). Hold the knot of the retrieval line and cut one strand. Gently pull on the knot of the retrieval line to remove it, while maintaining the stent position using the steerable positioner and guidewire.

C. The line may be removed before stent placement

The retrieval line may be cut prior to stent placement. Remove the retrieval line by holding the knot, cutting one strand, and while holding the knot, gently pull out the retrieval line.

For all options, stent removal or exchange should be done at the clinically appropriate time.

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Malecot Nephrostomy Catheter Set Re-Entry™ Malecot Nephrostomy Catheter Set

Directions for Use	2
Instrucciones de uso	5
Mode d'emploi	8
Gebrauchsanweisung	11
Istruzioni per l'uso	14
Gebruiksaanwijzing	17
Instruções de Utilização	20

50549798-01

2020-02

Malecot Nephrostomy Catheter Set

Re-Entry™ Malecot Nephrostomy Catheter Set

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DEVICE DESCRIPTION

The Malecot Nephrostomy Catheter Set is used for percutaneous nephrostomy drainage.

Contents

- (1) Re-Entry Catheter (C-Flex™ Elastomer Material) or Malecot Catheter (C-Flex Elastomer Material)
- (2) Flexible Stylets
- (1) Drainage Bag Connecting Tube

INTENDED USE/INDICATIONS FOR USE

The Malecot Nephrostomy Catheter Set is intended to establish percutaneous nephrostomy access and is indicated for the following:

Diagnostic Indications:

- Antegrade pyelography
- Pressure/perfusion study (Whitaker test)

Therapeutic Indications:

- Nephrostomy catheter drainage
- Perfusion chemolysis of renal stones
- Post-percutaneous nephrolithotomy
- Post-percutaneous resection and coagulation of urothelial tumors

CONTRAINDICATIONS

The Malecot Nephrostomy Catheter Set is not indicated for use in patients with blood clotting anomalies due to coagulopathies or pharmacological anticoagulations.

WARNING

None known.

PRECAUTIONS

1. Bending or kinking during, prior to, or after placement (including external fixation) could damage the integrity of the catheter.
2. If resistance is encountered during advancement or withdrawal of the catheter, STOP. Do not continue without first determining the cause of the resistance and taking remedial action.

Note: Where long-term use is indicated, it is recommended that indwelling time for catheter not exceed 90 days*. This catheter should be evaluated by the physician on or before 90 days postplacement.

*Biocompatibility data on file.

3. Periodic examinations are recommended to assess catheter patency and function.
4. Catheters are not intended to be permanent implant devices.
5. These recommendations are meant to serve only as a basic guide to the utilization of this device. The performance of percutaneous nephrostomy should not be undertaken without comprehensive knowledge of the indications, techniques and risks of the procedure.

ADVERSE EVENTS

- Hemorrhage
- Hematoma
- Fistula
- Edema
- Extravasation
- Loss of renal function
- Perforation
- Infection
- Encrustation
- Pneumothorax
- Catheter occlusion
- Catheter dislodgement
- Colon injury
- Tissue Damage
- Discomfort / pain
- Inflammation
- Hydronephrosis

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Handling and Storage

Store at controlled room temperature. Do not expose to organic solvents, ionizing radiation or ultraviolet light. Rotate inventory so that products are used prior to the sterilization expiration date on package label.

OPERATIONAL INSTRUCTIONS

Prior to Use

Carefully read all instructions **prior** to using this device.

RECOMMENDED PROCEDURE

Placement

1. Advance the tip of the flexible stylet squarely inside the tip of the malecot. Straighten the wings of the malecot by securing the male luer-lock on the stylet to the female luer-lock on the catheter.
2. Pass the tip of the catheter over the guidewire and advance down nephrostomy tract using a unidirectional twisting motion.
3. Unlock the stylet to allow malecot wings to form. Prior to removal of stylet, the proper positioning of the catheter should be confirmed. Remove stylet, secure catheter and attach connecting tube provided.

Withdrawal

1. If desired, the malecot wings may be straightened by inserting the additional stylet provided.

Note: In order to ensure stylet can enter the tip of the catheter - pull back slightly on catheter to straighten malecot. Insert stylet and engage tip. Secure stylet and remove catheter.

2. A guidewire may be passed into collecting system after stylet is engaged and prior to catheter removal.

ATTENTION ATTENDING PHYSICIAN: If patient will not be followed up by you, it is recommended that the "instructions for use" or the section on how to remove catheter be attached to the patient's chart.

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C-Flex is a registered trademark of Saint Gobain Performance Plastics.

REFERENCES

- Castaneda-Zuniga, W.R., Amplatz, K.: Percutaneous Nephrostomy with the Stamey Catheter: A New Introducing Technique. Urol. Clin. North America, 9:1 (Feb.) 1982.
- Castaneda-Zuniga, W.R., Smith, A., Tadavarthy, S.M., Kotula, F., Amplatz, A.: Flexible Trocar for Percutaneous Nephrostomy. AJR 136:434, (Feb.) 1981.
- Khasidy, L.R., Smith, A.D.: The Re-Entry Nephrostomy Catheter for Endourological Applications. Journal of Urology Vol. 133, (Feb.) 1985.

Malecot Nephrostomy Catheter Set Re-Entry™ Malecot Nephrostomy Catheter Set

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Advantage™ System
Advantage Fit™ System
Advantage Blue System
Advantage Fit Blue System

Transvaginal Mid-Urethral Sling

Directions for Use	2
Instrucciones de uso	9
Mode d'emploi	16
Gebrauchsanweisung	23
Istruzioni per l'uso	30
Gebruiksaanwijzing	37
Instruções de Utilização	44

50672518-01

2018-03

TABLE OF CONTENTS

WARNING	3
WARNING	3
DEVICE DESCRIPTION.....	3
INDICATIONS FOR USE.....	3
CONTRAINDICATIONS.....	3
GENERAL WARNING.....	4
POST PROCEDURAL WARNING.....	4
PRECAUTIONS.....	4
ADVERSE EVENTS	5
HOW SUPPLIED.....	5
OPERATION INSTRUCTIONS	6
WARRANTY	8

Advantage™ System

Advantage Fit™ System

Advantage Blue System

Advantage Fit Blue System

Transvaginal Mid-Urethral Sling

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician trained in use of surgical mesh for repair of stress urinary incontinence.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process.

Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

WARNING

This product is intended for use only by clinicians with adequate training and experience in the surgical treatment of stress urinary incontinence (SUI). The physician is advised to consult the medical literature regarding techniques, complications and hazards associated with the intended procedures.

DEVICE DESCRIPTION

The Advantage System, Advantage Fit System, Advantage Blue System and Advantage Fit Blue System are sterile, single use systems, consisting of one delivery device and one mesh assembly. Each mesh assembly is comprised of a polypropylene knitted mesh protected by a disposable plastic sleeve. At the distal ends of the mesh assembly are two dilators designed to be placed over the needle end of the delivery device. The disposable delivery device consists of a handle with a curved needle, and a pusher component. The delivery device is designed to facilitate the passage of the mesh assembly through bodily tissues for transvaginal placement.

INDICATIONS FOR USE

The mesh implant is intended as a suburethral sling for the treatment of stress urinary incontinence resulting from urethral hypermobility and/or intrinsic sphincter deficiency.

CONTRAINDICATIONS

A mesh implant is contraindicated in the following patients:

- Pregnant patients, patients with the potential for future growth or patients that are considering future pregnancies.

- Any patients with soft tissue pathology into which the implant is to be placed.
- Patients with any pathology which would compromise implant placement.
- Patients with any pathology, such as blood supply limitations or infections that would compromise healing.

GENERAL WARNING

The risks and benefits of performing a suburethral sling procedure in the following patients should be carefully considered:

- Careful consideration should be given to performing this procedure for patients with untreated coagulopathies or who are being treated with either anticoagulants or antiplatelet agents.
- Patients with hypertonic bladders or vesico ureteral reflux.
- Take special care in cases of bladder prolapse because of anatomical distortion. If the patient requires a cystocele repair, it should be done prior to the suburethral sling procedure.
- Vaginal and urinary tract infection should be treated prior to a suburethral sling implantation procedure.
- User should be familiar with surgical procedures and techniques involving nonabsorbable meshes.
- Good surgical practices should be followed for management of contamination or infected wounds.
- Mesh is considered a permanent implant. Removal of mesh or correction of mesh related complication may involve multiple surgeries.
- Complete removal of mesh may not be possible and additional surgeries may not always fully correct the complications.

POST PROCEDURAL WARNING

- If subsequent infection occurs, follow appropriate medical intervention practices.
- The patient should be advised that future pregnancies may negate the effects of this procedure and the patient may again become incontinent.

PRECAUTIONS

- The use of polypropylene mesh in urogynecologic procedures such as the treatment of stress urinary incontinence, regardless of the route of delivery (transvaginal, suprapubic or transobturator), has been associated with cases of erosion. Erosion has been reported in bladder, vagina, urethra, ureter, and bowel. Treatment of the erosion may require surgical removal.
- As with all surgical procedures, certain risk factors are known to impact patient outcomes in the pelvic floor which include, but are not limited to, impaired vascularity (e.g. diabetes, smoking status, estrogen status, pelvic floor radiation exposure, etc.), age, pelvic floor myalgia, impaired wound healing (e.g. diabetes, steroid usage, etc.), or active infection in or near the surgical site. The above pathophysiologic conditions must be considered when determining whether the patient is an appropriate candidate for mesh implantation, either by transvaginal, suprapubic or transobturator route.
- Standard surgical practices should be followed for the suburethral sling procedure as well as for the management of contaminated or infected wounds.
- Retropubic bleeding can occur. Check carefully before releasing patient from the hospital.
- Cystoscopy must be performed to confirm bladder integrity.
- Do not remove the protective plastic sleeve covering the mesh implant until proper position has been confirmed.
- Ensure the mesh is placed tension free under the mid-urethra.
- Use of this device should be done with the understanding that subsequent infection may require removal of the mesh.
- Physician should determine when it is suitable for each patient to return to normal activities.

- Patients should be counseled when to resume vigorous activities (heavy lifting, exercise), and intercourse after the procedure.
- Should dysuria, bleeding or other problems occur, the patient should be instructed to immediately contact the physician.
- Do not use any mechanical means of contact with the mesh (such as clips, staples, etc.) within the urethral support region of the mesh as mechanical damage to the mesh may occur.
- Avoid excessive tension on the mesh during handling.

ADVERSE EVENTS

The following adverse events have been reported due to suburethral sling placement, but are not limited to:

- As with all implants, local irritation at the wound site and/or a foreign body response may occur.
- Tissue responses to the mesh implant could include:
 - erosion/exposure/extrusion of the mesh through the vaginal or urethral mucosa, bladder wall or other surrounding tissue
 - scarring/scar contracture
 - device migration
 - fistula formation and inflammation

The occurrence of these events may require surgical intervention and possible removal of the entire mesh.

- Like all foreign bodies, the mesh may potentiate an existing infection.
- Excess tension may cause temporary or permanent lower urinary tract obstruction and retention.
- Allergic reaction has been reported.
- Known risks of surgical procedures for the treatment of incontinence include:
 - pain, ongoing pain (pelvic, vaginal, groin/thigh, dyspareunia)
 - infection
 - detrusor instability
 - complete failure of the procedure
 - voiding dysfunction (incontinence, mild to moderate incontinence due to incomplete urethral support or due to overactive bladder)
 - bruising, bleeding (vaginal, hematoma formation)
 - abscess
 - vaginal discharge
 - dehiscence of vaginal incision
 - edema and erythema at the wound site
 - perforation or laceration of vessels, nerves, bladder, urethra or bowel may occur during placement.

The occurrence of these events may require surgical intervention. In some instances the response to these events may persist as a permanent condition after the intervention.

HOW SUPPLIED


Each of the Systems are sterile, single use systems consisting of one (1) delivery device and one (1) mesh assembly.

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Store at a controlled room temperature. Do not expose to organic solvents, ionizing radiation or ultraviolet light. Rotate inventory so that devices are used prior to the expiration date on the package label.

FIGURE 1

OPERATION INSTRUCTIONS

NOTE: Review image above for part description.

Operational Instructions Prior to Use

Each of the Systems are supplied sterile and are intended for single patient use only. Carefully examine the system to verify that neither the contents nor the sterilized package has been damaged in shipment. **DO NOT USE** if the sterile barrier on product is damaged. Immediately return damaged product to Boston Scientific.


Prepare and drape the patient following standard surgical practice.

The design allows the operator a transvaginal route of placement.

WARNING: Make sure that the bladder is empty prior to initiating the use of this product. Ensure that the bladder, urethra and other important landmarks are properly identified.

Prepare the System for Use

FIGURE 2

Preparing the System for Use (See Figure 2)

1. Orient the Delivery Device handle so that the needle end is positioned away from the user and the needle tip is positioned upward.
2. Load the proximal end of the dilator tube by holding the dilator/sleeve joint and placing it over the distal end of the needle.
3. Slide the dilator tube over the needle until the dilator's proximal end abuts the distal end of the dilator pusher.

STEPS FOR USE

1. After preparation of the lower abdominal and vaginal operative sites, create two small transverse abdominal incisions approximately 0.5 cm to 1 cm on each side of the midline just above the symphysis.
2. Make a 1.0 cm to 1.5 cm vertical midline incision on the anterior vaginal wall at the level of the mid-urethra. Dissect bilaterally to the interior portion of the inferior pubic ramus at the 45° angle off the midline creating a pathway for delivery device placement.
3. Resting the tip of the needle on the palmar surface of the non-dominant index finger, gently introduce the delivery device into the vaginal incision and advance it anterolaterally to the lateral edge of the dissected space and then perforate the pubocervical and endopelvic fascia.

WARNING: Make sure the delivery device needle and mesh assembly pass sufficiently lateral to the urethra and bladder so that neither is injured.


4. After perforating the endopelvic fascia, use your fingertip, and guide the distal end of the needle superiorly along the posterior surface of the pubic bone. The curved part of the needle should rest in the operator's non-dominant hand during advancement of the device.
5. Carefully pass the needle through the space of Retzius and perforate the rectus muscle and rectus sheath. Guide the device into the ipsilateral abdominal incision until the needle tip is exposed through the incision.

WARNING: If excessive resistance is encountered during advancement/withdrawal, stop and determine remedial action prior to proceeding.

6. When the needle tip/dilator tube assembly extends extra-abdominally, advance the dilator pusher on the handle forward toward the needle tip end of the Delivery Device. This will cause the dilator tube to advance beyond the tip of the needle.
7. Grasp the dilator by placing a clamp or hemostat on the free end of the dilator end to temporarily secure it extra-abdominally.
8. While holding the dilator in position, withdraw the needle from inside the dilator and out of the vagina. If the dilator should retract back into the abdomen, advance the needle until the proximal end of the dilator abuts the distal end of the cannula and redeploy the Dilator/Mesh Assembly until the needle tip/dilator tube assembly extends extra-abdominally.
9. Repeat steps 1-8 on the contra lateral side.
10. With both dilator tubes in place, cystoscopy must be performed to confirm bladder integrity. If the blue dilator is seen in the bladder, remove the Dilator/Mesh Assembly. Visually inspect the Dilator/Mesh Assembly for integrity. If the Dilator/Mesh Assembly is intact, reload the Delivery Device and redeploy the Dilator/Mesh Assembly. The bladder must be emptied after cystoscopy.
11. Once desired sling placement is achieved, prepare to remove the protective sleeve from the mesh. See Figure 3 in the Tension Mesh/Sleeve Removal Section.

Tension Mesh/Sleeve Removal

FIGURE 3

Tension Mesh/Sleeve Removal

1. Adjust the mesh/sleeve by pulling upwards on the dilators so that the blue centering tab is centered below the urethra.
2. Appropriately tension the mesh/sleeve according to physician preference.
3. Grasp the blue centering tab and cut the tab through the center of the punch hole (see Figure 3) ensuring that both halves of the blue tab are completely removed from the vaginal canal.
4. Pull upwards on the dilators to remove the sleeve out of the body.

5. Verify the tension of the Mesh and adjust mesh as necessary.
6. Once the desired tension has been achieved, gently push downward on the abdomen, cut the distal ends of the mesh and confirm that those ends retract into the incision.
7. Close the incision in the usual manner.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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DEVICE DESCRIPTION

The Ureteral Stent contains a ureteral stent with attached retrieval line (suture), SureDrive Steerable Stent Positioner, and pigtail straightener.

The Ureteral Stent is a sterile, single use, disposable device that allows internal drainage from the kidney to the bladder.

The stent is MR safe.

Contents

- (1) Contour or Percuflex Plus Ureteral Stent with HydroPlus™ Coating
- (1) Compatible Steerable Stent Positioner

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None known.

PRECAUTIONS

Note: To avoid possible retrieval line encrustation, it is recommended that indwelling time not exceed fourteen (14) days if the retrieval line is left attached to the stent.

1. Recommended for one time use only.
2. Bending or kinking during or prior to placement could damage the integrity of the stent.
3. If resistance is encountered during advancement or withdrawal of the stent, **STOP**. Do not continue without first determining the cause of the resistance and taking remedial action.

- Periodic upper urinary tract examinations are recommended to ensure the stent is draining adequately. Evaluation for stent encrustation should be performed periodically.

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*Biocompatibility data on file.

- Stents are not intended to be permanent implant devices.
- The recommendations given are meant to serve only as a basic guide to the utilization of this set. The insertion of a ureteral stent should not be undertaken without comprehensive knowledge of the indications, techniques and risks of the procedure. The references given below provide a broad overview of the subject of ureteral stent insertion over a preplaced guidewire.
- Stent placement may complicate the subsequent work-up/evaluation of unexplained hematuria.
- Patients who are poor surgical risk should be managed appropriately.

ADVERSE EVENTS

Adverse events associated with retrograde or antegrade positioned indwelling ureteral stents include but are not limited to: Reflux-GU (e.g. ureteral reflux); Occlusion/Obstruction (e.g. catheter, stent); Migration (e.g. dislodgement); Hemorrhage; Infection (e.g. sepsis, peritonitis, urinary tract infection); Perforation (e.g. bladder, ureter, kidney, renal pelvis); Extravasation; Encrustation; Loss of renal function; Edema; Urinary symptoms (e.g. frequency, urgency, incontinence, dysuria, nocturia, hematuria); Pain/discomfort; Stent fragmentation; Fistula; Hydronephrosis; Stone formation; Tissue damage; Erosion.

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Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place. Rotate inventory so that products are used prior to the expiration date on package label.

OPERATIONAL INSTRUCTIONS

Placement Prerequisites

- Prior to stent placement, the renal collecting system should be visualized via intravenous, retrograde or antegrade pyelography.
- A stent of proper length should be available. Ideally, the upper coil should lie within the renal pelvis while the lower coil recurves at the ureteral orifice.
- Fluoroscopy is recommended for more precise control of stent placement, however, standard radiography may be used.

Prior to placement: To activate the HydroPlus™ coating, soak the stent in sterile water or saline for at least 30 seconds. Keep the stent wet during placement. Use a wet gauze pad to handle the stent during placement if necessary. Do not wipe stent with dry gauze or any solvents as it may damage coating.

Note: Prior to stent placement, remove and discard tubing from the kidney coil (Figure 1).

Note: The tube should be removed prior to attempting to pull on the retrieval line. Otherwise the line could become caught on the renal coil as it is being pulled through the tube and in some instances can tear or kink the stent.


Figure 1

Check the stent system (components)


Figure 2

Note: Pigtail straightener slides easily over ureteral stent (first proximally, second distally) straightening out the pigtail and facilitating introduction of guidewire. After introduction of guidewire, pigtail straightener is removed.

Recommended Techniques for Placement with Guidewire

Retrograde Endoscopic Placement – Guidewire Technique:

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Percutaneous Antegrade Placement Technique:

1. Establish entry into involved renal pelvis using a percutaneous access set.
2. Pass flexible end of guidewire down ureter into bladder.
3. Pass appropriate open end ureteral catheter over guidewire to confirm entrance into bladder. Placement of a semi-rigid sheath, of appropriate size, into the access site is a useful adjunct at this point.
4. Pass the tapered tip of the stent over the guidewire. Load steerable positioner onto guidewire and securely attach to stent (Figure 3).
5. Advance stent down ureter with steerable positioner. Adjust stent placement with positioner if necessary.
6. Confirm lower coil to be in bladder and upper coil within renal pelvis.

7. Disengage the steerable positioner from the stent by rotating the thumb knob clockwise followed by pushing the knob forward (figure 4).
8. Remove retrieval line by holding the knot, cutting one strand, and while holding the knot, gently pull on the retrieval line. The guidewire is then removed, stabilizing the stent with steerable positioner.

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Retrieval Line Management Options

- A. The line may be left in place for subsequent stent removal. To avoid possible retrieval line encrustation the stent should be left in place for a maximum of fourteen (14) days.

- B. The line may be used as an aid to positioning and then be removed once the stent is in the desired position.

After stent placement, but before the steerable positioner is removed, disengage the steerable positioner from the stent by rotating the thumb knob clockwise followed by pushing the knob forward (figure 4). Hold the knot of the retrieval line and cut one strand. Gently pull on the knot of the retrieval line to remove it, while maintaining the stent position using the steerable positioner and guidewire.

- C. The line may be removed before stent placement

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EC REP EU Authorized Representative

Boston Scientific Limited
Ballybrit Business Park
Galway
IRELAND

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Recyclable Package

C € 2797

© 2019 Boston Scientific Corporation or its affiliates.
All rights reserved.

Obtryx™ II System

CURVED

Transobturator Sling with
PrecisionBlue™ Design

Directions for Use	2
Instrucciones de uso	9
Mode d'emploi	15
Gebrauchsanweisung	21
Istruzioni per l'uso	27
Gebruiksaanwijzing	33
Instruções de Utilização	39

91077751-01

2016-06

TABLE OF CONTENTS

WARNING 3

WARNING 3

DEVICE DESCRIPTION..... 3

INDICATIONS FOR USE..... 3

CONTRAINDICATIONS..... 3

HOW SUPPLIED..... 4

 Handling and Storage..... 4

DIRECTIONS FOR USE..... 4

 Prior to Use..... 4

WARNING 4

WARNING 4

WARNING 5

WARNING 5

TENSION MESH/SLEEVE REMOVAL..... 5

GENERAL WARNING..... 6

POST PROCEDURAL WARNING 6

ADVERSE EVENTS 6

PRECAUTIONS..... 7

WARRANTY 8

Obtryx™ II System

CURVED

Transobturator Sling with PrecisionBlue™ Design

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WARNING

Contents supplied STERILE using ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

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DEVICE DESCRIPTION

The Obtryx II System is a sterile, single use system consisting of two (2) delivery devices and one (1) mesh assembly. The mesh assembly is comprised of a polypropylene knitted mesh with dilator legs and a center tab. At the distal ends of the dilator legs there are association loops designed to be placed in the needle slot of the distal end of the delivery device. The disposable delivery device consists of a handle with a stainless steel needle. The needle is designed to facilitate the passage of the mesh assembly through bodily tissues for placement through the obturator foramen.

INDICATIONS FOR USE

The mesh implant is intended for use as a suburethral sling for the treatment of stress urinary incontinence resulting from urethral hypermobility and/or intrinsic sphincter deficiency.

CONTRAINDICATIONS

The mesh suburethral sling implant is contraindicated in the following patients:

- Pregnant patients, patients with potential for future growth or patients that are considering future pregnancies.
- Any patients with soft tissue pathology into which the implant is to be placed.

- Patients with any pathology which would compromise implant placement.
- Patients with any pathology, such as blood supply limitations or infections that would compromise healing.

HOW SUPPLIED

The device is supplied sterile. Do not use if package is opened or damaged.

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Handling and Storage

Store at a controlled room temperature. Do not expose to organic solvents, ionizing radiation or ultraviolet light. Rotate inventory so that products are used prior to the expiration date on package label.

DIRECTIONS FOR USE

Prior to Use

Carefully examine the system to verify that neither the contents nor the sterilized package has been damaged in shipment. DO NOT USE if sterile barrier on product is damaged. Immediately return damaged product to Boston Scientific.

The design of the Obtryx™ II System allows the operator a percutaneous approach utilizing a transobturator technique. See Figure 1 for parts description.


Figure 1: Parts Description

Prepare and drape the patient using standard surgical practice.

WARNING

Assure that the bladder is empty prior to initiating the use of this product. Ensure that the bladder, urethra and other important landmarks are properly identified.

Steps to Use

1. Prepare the skin lateral to the inferior pubic ramus and vaginal operative sites.
2. Make a 1.0 cm to 1.5 cm vertical midline incision on the anterior vaginal wall at the level of the mid-urethra. Dissect bilaterally to the interior portion of the inferior pubic ramus at the 45° angle off the midline creating a pathway for delivery device placement.
3. Create a vertical skin incision large enough to insert tip of needle just lateral to the edge of the inferior pubic ramus at the junction where the inferior pubic ramus and the adductor longus muscle meet. Repeat on the contralateral side.

WARNING

If excessive force is encountered during advancement/withdrawal, stop and determine remedial action prior to proceeding.

4. Grasp the device handle and insert one (1) needle through one (1) skin incision, piercing through the obturator muscle and obturator membrane. Turn the handle at the 45° angle medial towards the midline. Place the opposite hand's forefinger into the lateral dissection of the vaginal incision, placing the fingertip on the distal end of the needle. Guide the distal end of the needle around the inferior pubic ramus through the vaginal incision, maintaining contact with the finger.

WARNING

Pay careful attention to avoid the adductor longus tendon with the delivery device.

WARNING

Make sure the delivery device and mesh assembly pass sufficiently lateral to the urethra in order to avoid urethral injury.

- Engage one (1) association loop to the distal end of the needle (see Figure 2) protruding through the vagina.


Figure 2: Association Loop Engagement

- Pull the needle out through the skin incision. Be sure that the mesh assembly is not twisted and lies flat under the urethra, with the blue center tab positioned suburethraly, facing outward.
- Remove the association loop from the needle (see Figure 3).


Figure 3: Association Loop Removal

- Repeat Steps 4-7 on the contralateral side with the second needle.
- Cystoscopy may be performed at this time, to be determined at the physician's discretion.
- Next see section "Tension Mesh/Sleeve Removal."

TENSION MESH/SLEEVE REMOVAL

- Adjust the mesh/sleeve by pulling outwards on the dilators so that the blue center tab is centered below the urethra.
- Appropriately tension the mesh/sleeve according to physician preference.
- Once proper tension is achieved, cut the leader loop that is on the outside of the sleeve that is connecting the dilator leg and sleeve to the mesh. Pull outward on the dilator to remove the sleeve leaving the mesh in place. Repeat on the other side. (See Figure 4).


Figure 4: Tension Mesh/Sleeve Removal

- Grasp the blue center tab and cut the center tab lead located on the side of the center tab to release the tab from the mesh. Remove the center tab and center tab lead from the vaginal canal.

5. Gently pushing downward on the skin incisions, cut the distal ends of the mesh and confirm that those ends retract into the skin incisions.
6. Close all incisions per standard practice.

GENERAL WARNING

The risks and benefits of performing a suburethral sling procedure in the following should be carefully considered:

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POST PROCEDURAL WARNING

- If subsequent infection occurs, follow appropriate medical intervention practices.
- The patient should be advised that future pregnancies may negate the effects of this procedure and the patients may again become incontinent.

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The occurrences of these events may require surgical intervention and possible removal of the entire mesh.

- Like all foreign bodies, the mesh may potentiate an existing infection.
- Excess tension may cause temporary or permanent lower urinary tract obstruction and retention.
- Allergic reaction has been reported.
- Known risks of surgical procedures for the treatment of incontinence include:
 - pain, ongoing pain (pelvic, vaginal, groin/thigh, dyspareunia)
 - infection
 - detrusor instability
 - complete failure of the procedure

- voiding dysfunction (incontinence, mild to moderate incontinence due to incomplete urethral support or due to overactive bladder)
- bruising, bleeding (vaginal, hematoma formation)
- abscess
- vaginal discharge
- dehiscence of vaginal incision
- edema and erythema at the wound site
- perforation or laceration of vessels, nerves, bladder, urethra or bowel may occur during placement.

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PRECAUTIONS

- The use of polypropylene mesh in urogynecologic procedures such as the treatment of stress urinary incontinence regardless of the route of delivery (transvaginal, suprapubic or transobturator), has been associated with cases of erosion. Erosion has been reported in bladder, vagina, urethra, ureter, and bowel. Treatment of the erosion may require surgical removal.
- As with all surgical procedures, certain risk factors are known to impact patient outcomes in the pelvic floor which include, but are not limited to, impaired vascularity (e.g. diabetes, smoking status, estrogen status, pelvic floor radiation exposure, etc.), age, pelvic floor myalgia, impaired wound healing (e.g. diabetes, steroid usage, etc.), or active infection in or near the surgical site. The above pathophysiologic conditions must be considered when determining whether the patient is an appropriate candidate for mesh implantation, either by transvaginal, suprapubic or transobturator route.
- Standard surgical practices should be followed for the suburethral sling procedure as well as for the management of contaminated or infected wounds.
- Bleeding can occur. Check carefully before releasing patient from the hospital.
- Do not remove the protective plastic sleeve covering mesh implant until proper position has been confirmed.
- Ensure the mesh is placed without tension under the mid-urethra.
- Use of this device should be done with the understanding that subsequent infection may require removal of the mesh.
- Physician should determine when it is suitable for each patient to return to normal activities.
- Patients should be counseled when to resume vigorous activities (heavy lifting, exercise), and intercourse after the procedure.
- Should dysuria, bleeding or other problems occur, the patient should be instructed to immediately contact the physician.
- Do not use any mechanical means of contact with the mesh (such as clips, staples etc.) within the urethral support region of the mesh as mechanical damage to the mesh may occur.
- Avoid excessive tension on the mesh during handling.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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DEVICE DESCRIPTION

The Obtryx II System is a sterile, single use system consisting of two (2) delivery devices (one patient right and one patient left) and one (1) mesh assembly. The mesh assembly is comprised of a polypropylene knitted mesh with dilator legs and a center tab. At the distal ends of the dilator legs there are association loops designed to be placed in the needle slot of the distal end of the delivery device. The disposable delivery device consists of a handle with a stainless steel needle. The needle is designed to facilitate the passage of the mesh assembly through bodily tissues for placement through the obturator foramen.

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4. Grasp the device handle for the patient's left side with the right hand. Place the left forefinger into the lateral dissection of the vaginal incision. Place the needle tip into the skin incision perpendicular to the skin with the handle at a 45° angle parallel with the thigh.
5. Putting the left thumb on the outside of the needle curve, apply a downward force, piercing through the obturator muscle and membrane.

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6. Rotate the needle medially around the inferior pubic ramus to meet the left hand forefinger. Guide the needle tip through the vaginal incision.
-

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7. Engage one (1) association loop to the distal end of the needle (see Figure 2) protruding through the vagina.

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8. Pull the needle out through the skin incision. Be sure that the mesh assembly is not twisted and lies flat under the urethra, with the blue center tab positioned suburethrally, facing outward.
9. Remove the association loop from the needle (see Figure 3).


Figure 3: Association Loop Removal

10. Repeat Steps 4-9 on the contralateral side with the second needle.
11. Cystoscopy may be performed at this time, to be determined at the physician's discretion.
12. Next see section "Tension Mesh/Sleeve Removal."

TENSION MESH/SLEEVE REMOVAL

1. Adjust the mesh/sleeve by pulling outwards on the dilators so that the blue center tab is centered below the urethra.
2. Appropriately tension the mesh/sleeve according to physician preference.
3. Once proper tension is achieved, cut the leader loop that is on the outside of the sleeve that is connecting the dilator leg and sleeve to the mesh. Pull outward on the dilator to remove the sleeve leaving the mesh in place. Repeat on the other side. (See Figure 4).

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4. Grasp the blue center tab and cut the center tab lead located on the side of the center tab to release the tab from the mesh. Remove the center tab and center tab lead from the vaginal canal.
5. Gently pushing downward on the skin incisions, cut the distal ends of the mesh and confirm that those ends retract into the skin incisions.
6. Close all incisions per standard practice.

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- Like all foreign bodies, the mesh may potentiate an existing infection.
- Excess tension may cause temporary or permanent lower urinary tract obstruction and retention.
- Allergic reaction has been reported.
- Known risks of surgical procedures for the treatment of incontinence include:
 - pain, ongoing pain (pelvic, vaginal, groin/thigh, dyspareunia)
 - infection

- detrusor instability
- complete failure of the procedure
- voiding dysfunction (incontinence, mild to moderate incontinence due to incomplete urethral support or due to overactive bladder)
- bruising, bleeding (vaginal, hematoma formation)
- abscess
- vaginal discharge
- dehiscence of vaginal incision
- edema and erythema at the wound site
- perforation or laceration of vessels, nerves, bladder, urethra or bowel may occur during placement.

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WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC’s control directly affect the instrument and the results obtained from its use. BSC’s obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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