FAST-CATH™ HEMOSTASIS INTRODUCER

English - Instructions for use	2
Deutsch - Gebrauchsanleitung	4
Français - Mode d'emploi	6
Italiano - Istruzioni per l'uso	8
Español - Instrucciones de uso	10
中文 - 使用说明	12
Ελληνικά - Οδηγίες χρήσης	13
Русский - Инструкция по применению	15
Magyar - Használati utasítás	17
Português - Instruções de Utilização	19
Nederlands - Gebruiksaanwijzing	21
Suomi - Käyttöohjeet	23
Svenska - Bruksanvisning	25
Dansk - Brugsanvisning	27
Norsk - Bruksanvisning	29
Türkçe - Kullanım talimatları	30
한국어 - 사용 지침	31
Polski - Instrukcja użytkowania	32
Srpski - Uputstvo za upotrebu	33
SYMBOL	34

FAST-CATH™

HEMOSTASIS INTRODUCER

Instructions for use

See Product Catalog for Reorder Numbers.

THE FRENCH SIZE SPECIFIED REPRESENTS THE INNER DIAMETER OF THE INTRODUCER SHEATH.

DO NOT ATTEMPT TO INSERT A CATHETER HAVING A DISTAL TIP OR BODY SIZE LARGER THAN THE INTRODUCER SIZE INDICATED.

THE ST. JUDE MEDICAL HEMOSTASIS INTRODUCER IS DESIGNED TO INTERLOCK ONLY WITH ST. JUDE MEDICAL DILATORS AND/OR OBTURATORS. MISUSE MAY RESULT IN SERIOUS COMPLICATIONS.

SEE INDIVIDUAL STERILE PACKAGE LABEL FOR CONTENTS.

SINGLE-USE DISPOSABLE MEDICAL DEVICE. CONTENTS ARE STERILE IF PACKAGE IS UNOPENED AND UNDAMAGED.

DO NOT RESTERILIZE.

CAUTION:

DO NOT alter this device.

Federal law (U.S.) restricts this device to sale by or on the order of a physician.

Do not attempt to use a guidewire over maximum diameter specified on package label.

This device should only be used by physicians thoroughly trained in the technique of angiography and/or the use of St. Jude Medical catheter delivery systems.

Individual patient anatomy and physician technique may require procedural variations.

The closing force of the hemostasis valve may alter or impair the function of some catheters.

Damage to guidewire may result if withdrawn through a metal needle cannula. Cannula must be removed first.

Do not attempt to advance or withdraw guidewire if resistance is felt. Use fluoroscopy to determine cause.

Damage to the valve assembly may occur if inner catheter is withdrawn rapidly.

Advance dilator/sheath assembly with a twisting motion to avoid damage to the sheath or vessel.

STORE IN A COOL, DARK, DRY PLACE.

For Single Use Only! Single-use devices are designed and tested for only one patient application. These are disposable devices and are not designed for reprocessing and reuse. Reuse of designated "single-use" devices creates a risk of patient or user infections due to prior patient use and the difficulty in cleaning the narrow structures at material interfaces following direct blood contact. Contamination or reprocessing cleaning agent residues may lead to adverse patient reactions and may damage the device. Use of non-St. Jude Medical packaging may compromise device functionality and sterility due to compromised protection from shipping and handling damage. The absence of labeling after reprocessing, may lead to misuse of the device and impaired traceability. Reprocessing and reuse may result in patient or user injury, permanent impairment or death.

APPLICATION

The introduction of angiographic catheters, closed end catheters, balloon catheters, and electrodes into a vessel where minimizing blood loss is essential.

SUGGESTED DIRECTIONS FOR USE

INSPECT SHEATH INTRODUCER AND ACCESSORIES FOR DAMAGE, DO NOT USE ANY DAMAGED DEVICES.

BEFORE PROCEEDING, VERIFY THE HEMOSTASIS INTRODUCER IS OF PROPER SIZE TO ALLOW PASSAGE OF DESIRED CATHETER OR OTHER MEDICAL DEVICE.

FLUSH THE DILATOR(S), INTRODUCER SHEATH AND INFUSION SIDEPORT AND STOPCOCK WITH HEPARINIZED NORMAL SALINE.

1. Insert dilator into hemostasis valve. Follow normal accepted practice for vessel puncture, guidewire insertion, vessel dilator and hemostasis introducer use. For introducers with shaped distal tips, always advance dilator/sheath assembly over appropriate sized guidewire to the appropriate anatomical site.

CAUTION: DO NOT ATTEMPT ADVANCEMENT WITHOUT GUIDEWIRE. SEVERE VASCULAR DAMAGE AND/OR INJURY MAY OCCUR.

CAUTION: INSERTION INTO ARTERY MAY CAUSE EXCESSIVE BLEEDING AND/OR OTHER COMPLICATIONS.

CAUTION: DO NOT ALLOW GUIDEWIRE TO INADVERTENTLY ADVANCE TOTALLY INTO PATIENT.

2. When a sideport introducer is used, follow normal practice of using a continuous drip of anticoagulant fluid under pressure through the sideport when the hemostasis introducer is in the vessel.

CAUTION: UPON WITHDRAWAL OF DILATOR, CATHETER OR OTHER INSERTED DEVICE, ASPIRATION THROUGH THE STOPCOCK IS RECOMMENDED TO REMOVE ANY FIBRIN ACCUMULATION.

3. A St. Jude Medical repositioning sleeve may be used with the hemostasis introducer to provide an added measure of environmental protection to catheters (see St. Jude Medical product catalog).
4. Follow manufacturer's recommendations for the catheter or device being introduced via the hemostasis introducer.
5. When using the St. Jude Medical luer-lock hemostasis introducer, a St. Jude Medical luer-lock obturator may be used to seal off the hemostasis aperture and aid in maintaining sheath patency and general management of the introduction site.

OVERFLUSHING

The French size shown on each St. Jude Medical introducer represents the tip inside dimension.

To overflush Catheters or St. Jude Medical Luer-Lock Obturators placed into a St. Jude Medical hemostasis introducer, use introducer that is one (1) French size larger than the French size of the catheter or obturator.

CAUTION: UPON REMOVAL OF HEMOSTASIS INTRODUCER, PROPER PRECAUTIONS SHOULD BE TAKEN TO PREVENT BLEEDING.

INTRODUCER WITH SHAPED DISTAL TIPS

Reinsert the dilator fully into the introducer sheath to aid in straightening the top portion. Then remove the dilator and introducer sheath as a unit.

USA ONLY:

LIMITED WARRANTY AND DISCLAIMER

St. Jude Medical ("SJM") hereby warrants that if any SJM product fails to perform within normal tolerances for a patient due to a defect in materials or workmanship, SJM will provide, at no charge, a replacement SJM product for the patient's use. This limited warranty applies only if each of the following conditions are met.

1. The product was packaged and labeled by SJM.
2. The failed product must be returned to SJM and becomes the property of SJM.
3. The product has not been mishandled, reprocessed or altered in any way.
4. The product was used before the "USE BY" date marked on the packaging of the product.

No representation or warranty is made that a SJM product will not fail. SJM disclaims responsibility for any medical complications, including death, resulting from the use of its products. Except as expressly provided by this limited warranty, SJM IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES BASED ON ANY DEFECT, FAILURE, OR MALFUNCTION OF ITS PRODUCTS, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT, OR OTHERWISE. Some states do not allow the exclusion or limitation of incidental or consequential damages however, so the above limitation or exclusion may not apply to you.

Except as expressly provided by the limited warranty, SJM MAKES NO WARRANTY, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE.

STATE OF CALIFORNIA (USA) ONLY:

WARNING: This product and its packaging have been sterilized with ethylene oxide. The packaging may expose you to ethylene oxide, a chemical known to the State of California to cause cancer or birth defects or other reproductive harm.

WARNING: Products may be manufactured with CFC-113 and CFC12 substances which harm public health and the environment by destroying ozone in the upper atmosphere.

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All rights reserved.
are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies.