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Imported and Marketed in India by:
Fresenius Kabi India Pvt. Ltd.
A-3 MIDC, Warehouse-A, GR. & First Floor
Ranjangaon Ganpati, Tal-Shirur,
Dist. Pune - 412220 India
Customer Care No.: 09158898288

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Imported and Distributed By:
Fresenius Kabi Philippines Inc.
18th Flr Units 5-7 Zuellig Bldg.
Makati Ave. cor. Paseo De Roxas
Makati City, 1225 Philippines
Reg. No. MDR-00110

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Document No.: **472311549**
Revision: **A**
Change Order: **CO-030445**
Title: **PRODUCT INSERT AMICUS EXCHANGE KIT FOR TPE
PROCEDURE X6R2339**

Page: 5 of 5
Last Modified: 23 MAY 2017 MKF
Color: Black

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Document No.: **472312915**
Revision: **A**
Change Order: **CO-034005**
Title: **PRODUCT INSERT AMICUS APHERESIS KIT – SINGLE NEEDLE WITH PLATELET ADDITIVE SOLUTION CONNECTOR X6R2301**

Page: 1 of 7
Last Modified: 01 Aug 2018 JML
Color: Black

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6x

FRESENIUS KABI

(en) AMICUS Apheresis Kit – Single Needle with Platelet Additive Solution Connector

For Extended Platelet Storage, Plasma and Red Cell Collection.* For use with the AMICUS separator.

Store at Room Temperature. Protect from freezing. Avoid excessive heat.

Definition of “Room Temperature”:

“The temperature prevailing in a working area.”

Reference: United States Pharmacopeia, General Notices. United States Pharmacopeial Convention, Inc. 12601 Twinbrook Parkway, Rockville, MD.

Description:

This Apheresis Kit consists of:

- (2) 1000 mL plastic platelet storage containers
- (1) 800 mL plastic plasma storage container
- (1) 600 mL plastic red cell container
- (1) 17 gauge needle
- (1) return line filter with pore size of 220 microns
- (3) antimicrobial filters with pore size of 0.2 microns

Indications and Usage:

This kit is designed for use with the AMICUS separator for the collection of leuko-reduced platelet concentrate, plasma and red cells.

Platelet concentrate may be stored in the platelet storage container(s) for up to 5 days at 20°C to 24°C with continuous gentle agitation, and according to the applicable procedures and regulations. Plasma and red blood cells must be handled (freezing, storage conditions, labeling, etc.) according to the applicable procedures and regulations. Red blood cells must further be processed using a cleared AMICUS red blood cell set. Red blood cells should be stored in an appropriate storage solution and container.

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Warnings:

- Procedures involving extracorporeal circulation are associated with the risks of blood loss, hemolysis, air embolism and/or blood clotting. To minimize these risks, do not deviate from the procedures described in the Operator’s Manual.
- To protect system integrity, do not prime the inlet/return line until immediately before venipuncture.
- The separator should be operated only by a qualified healthcare professional who has been trained in the operation of this instrument.
- Failure to follow instructions prior to venipuncture may result in dilution of the whole blood sample.
- The donor should meet appropriate regulatory requirements for the intended procedure and be in good health on the day of the donation.
- The platelet count of the donor should be taken into consideration when performing any apheresis procedure. Because of the increased platelet collection efficiency of the AMICUS separator, it may be necessary to reduce the volume of whole blood processed, compared to other blood cell separators, to avoid unacceptable low post-donation counts.
- Modification of the functionally closed system or replacement of needle under open system conditions negates any claims to a closed system, unless approved sterile connection methods are employed.
If closed system conditions have been maintained, the platelet product may be stored for up to 5 days. If open system collection conditions have occurred, platelets must be stored no longer than 24 hours. Any loss of sterile closure following collection results in a platelet product which must be stored no longer than 4 hours.
- The use of a sterile connecting device to modify the configuration of the apheresis kit should be done using an operating procedure that is consistent with applicable regulatory guidelines and the manufacturer’s directions for use.
- Phthalates are an integral part of the medical device formulation. Blood or blood components collected, prepared or transfused using product may contain low levels of DEHP and present certain risks. In general, the highest risk is associated with medical procedures involving infants, pregnant or nursing mothers, or requiring high volume to body weight blood infusion (such as trauma patients and neonate replacement transfusion). The use for infusion of blood or blood components with limited DEHP content must be weighed against the risks associated with not performing the medical procedure.
- This product is designed for single use. Biological contamination of this product will result from normal use. Re-use of this product can expose patients or donors to this biological contamination. This exposure can result in adverse reactions, including severe illness and possibly death. This product is not designed for re-sterilization. Re-sterilization of this product may not completely eliminate biological contamination and may result in malfunction.

Adverse Effects:

Adverse donor reactions, identical to those seen with routine blood collection procedures may occur. These include “lightheadedness”, fainting, vomiting, hyperventilation and hematoma formation at the site of phlebotomy. Syncopal reactions due to hypovolemia may occur.

Fresenius Kabi AG
61346 Bad Homburg / Germany
Tel.: +49 (0) 61 72 / 686-0
www.fresenius-kabi.com

Fenwal International, Inc.
Carretera Sanchez Km 18.5
Parque Industrial Itabo
Zona Franca Ind. De S.C.
Haina, Dominican Republic
Shelf life: 24 months

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Complications such as blood loss, hemolysis, air embolism and/or blood clotting may be associated with improper operating conditions.

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Tare Weights:

Below are the approximate tare weights of the platelet container and associated kit components. Refer to the appropriate illustration for the tare weights for platelet product calculations.

Approximate Weights (grams)		
Two (2) platelet storage containers with: Two (2) sample pouch assemblies, Two (2) roller clamps and tubings, One (1) Y-Junction, and One (1) 41.9 cm tubing	77	
Two (2) platelet storage containers with: Two (2) sample pouch assemblies, Two (2) roller clamps and tubings, and One (1) Y-Junction	73	
One (1) platelet storage container with: One (1) sample pouch assembly, One (1) roller clamp and tubing, and One (1) Y-Junction	37	
One (1) platelet storage container with: One (1) sample pouch assembly	31	
One (1) platelet storage container	27	
Plastic tubing	0.12 g/cm	
One (1) plasma storage container (800 mL)	39	
One (1) red cell container (600 mL)	32	

*On AMICUS separators configured for red cell collection.

Symbols with Definitions:

Caution, consult instructions for use

Sterilized by irradiation.
Sterile fluid path.

Non-pyrogenic fluid path

Do not reuse

Do not use if the product sterile
barrier system is compromised

Do not vent

Contains or presence of phthalate:
Diethylhexylphthalate (DEHP)

This marking reflects compliance with
the Directive on Medical Devices.

Code

Lot

Expiration date

Manufacturing date

Manufacturer

Manufacturing facility/
Manufactured by

Single needle apheresis
protocol

Fragile

This way up

Recyclable

47-23-12-915 REV: A

Document No.: 472312915
Revision: A
Change Order: CO-034005
Title: PRODUCT INSERT AMICUS APHERESIS KIT – SINGLE NEEDLE WITH PLATELET ADDITIVE SOLUTION CONNECTOR X6R2301

Page: 2 of 7
Last Modified: 01 Aug 2018 JML
Color: Black

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© 2018 Fresenius Kabi AG.
Semua hak dilindungi undang-undang.

Document No.: 472312918
Revision: A
Change Order: CO-034005
Title: PRODUCT INSERT AMICUS APHERESIS KIT -
DOUBLE NEEDLE WITH PLATELET ADDITIVE
SOLUTION CONNECTOR X6R2302

Page: 1 of 6
Last Modified: 06 Sep 2018 JML
Color: Black

PROOFREADING INSPECTION / RELEASED ARTWORK

Proofreading Approval:

Print Name Signature Date

Code
REF X6R2302

6x

(en) AMICUS Apheresis Kit – Double Needle with Platelet Additive Solution Connector

For Extended Platelet Storage and Plasma Collection. For use with the AMICUS separator.

Store at Room Temperature*. Protect from freezing. Avoid excessive heat*.

Definition of “Room Temperature”:
“The temperature prevailing in a working area.”

Reference: United States Pharmacopeia, General Notices. United States Pharmacopeial Convention, Inc.
12601 Twinbrook Parkway, Rockville, MD.

*IP Definition of “Room Temperature”
Indian Pharmacopoeia, General Notices.
Indian Pharmacopoeia Commission, Ghaziabad, India

Description:

This Apheresis Kit consists of:
(2) 1000 mL plastic platelet storage containers
(1) 800 mL plastic plasma storage container
(1) 17 gauge needle (inlet line)
(1) 18 gauge needle (return line)
(1) Return line filter with pore size of 220 microns
(4) Antimicrobial filters with pore size of 0.2 microns

Indications and Usage:

This kit is designed for use with the AMICUS separator for the collection of leukoreduced platelet concentrate and plasma.

Platelet concentrate may be stored in the platelet storage container(s) for up to 5 days at 20°C to 24°C with continuous gentle agitation, and according to the applicable procedures and regulations. Plasma must be handled (freezing, storage conditions, labeling, etc.) according to the applicable procedures and regulations.

Each container can store up to a maximum of 4.7 x 10¹¹ platelets.

Contraindication:

Use of the AMICUS separator is contraindicated in those cases where adequate anticoagulation cannot be achieved.

Refer to the AMICUS Operator’s Manual for a complete listing of contraindications, adverse effects, warnings and cautions.

Warnings:

- Procedures involving extracorporeal circulation are associated with the risks of blood loss, hemolysis, air embolism and/or blood clotting. To minimize these risks, do not deviate from the procedures described in the Operator’s Manual.
- To protect system integrity, do not prime the inlet and return lines until immediately before venipuncture.
- The separator should be operated only by a qualified healthcare professional who has been trained in the operation of this instrument.
- Failure to follow instructions prior to venipuncture may result in dilution of the whole blood sample.
- The donor should meet appropriate regulatory requirements for the intended procedure and be in good health on the day of the donation.
- The platelet count of the donor should be taken into consideration when performing any apheresis procedure. Because of the increased platelet collection efficiency of the AMICUS separator, it may be necessary to reduce the volume of whole blood processed, compared to other blood cell separators, to avoid unacceptable low post-donation counts.
- Modification of the functionally closed system or replacement of the inlet needle under open system conditions negates any claims to a closed system, unless approved sterile connection methods are employed. Replacement of the return line needle does not compromise the closed system.
If closed system conditions have been maintained, the platelet product may be stored for up to 5 days. If open system collection conditions have occurred, platelets must be stored no longer than 24 hours. Any loss of sterile closure following collection results in a platelet product which must be stored no longer than 4 hours.
- The use of a sterile connecting device to modify the configuration of the apheresis kit should be done using an operating procedure that is consistent with applicable regulatory guidelines and the manufacturer’s directions for use.
- Phthalates are an integral part of the medical device formulation. Blood or blood components collected, prepared or transfused using product may contain low levels of DEHP and present certain risks. In general, the highest risk is associated with medical procedures involving infants, pregnant or nursing mothers, or requiring high volume to body weight blood infusion (such as trauma patients and neonate replacement transfusion). The use for infusion of blood or blood components with limited DEHP content must be weighed against the risks associated with not performing the medical procedure.
- This product is designed for single use. Biological contamination of this product will result from normal use. Re-use of this product can expose patients or donors to this biological contamination. This exposure can result in adverse reactions, including severe illness and possibly death. This product is not designed for re-sterilization. Re-sterilization of this product may not completely eliminate biological contamination and may result in malfunction.

Fresenius Kabi AG
61346 Bad Homburg / Germany
Tel.: +49 (0) 61 72 / 686-0
www.fresenius-kabi.com

Fenwal International, Inc.
Carretera Sanchez Km 18.5
Parque Industrial Itabo
Zona Franca Ind. De S.C.
Haina, Dominican Republic
Shelf life: 24 months

Imported and Marketed in India by:
Fresenius Kabi India Pvt. Ltd.
A-3 MIDC, Warehouse-A, GR. & First Floor
Ranjangaon Ganpati, Tal-Shirur,
Dist. Pune - 412220
India
Customer Care No.: 09158898288
Email ID: fkpi@Fresenius-kabi.com
İthalatçı Firma:
GST MEDİKAL SAĞLIK ÜRÜNLERİ SAN. TİC. A.Ş.
No:73 Gazi Mahallesi 06560 Yenimahalle ANKARA
T: +90 312 222 1 333 F: +90 312 222 2 223
info@gst.com.tr

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Adverse Effects:

Adverse donor reactions, identical to those seen with routine blood collection procedures may occur. These include “lightheadedness,” fainting, vomiting, hyperventilation and hematoma formation at the site of phlebotomy. Syncopal reactions due to hypovolemia may occur.

Reactions unique to apheresis collection procedures may also occur. Allergic symptoms including skin redness, itching, hives, etc. may be seen. Cooling of saline or donor blood may be sufficient to induce chills.

Infusion of anticoagulants containing citrate may lead to donor symptoms of moderate hypocalcemia due to chelation of calcium by unmetabolized citrate. Such donor reactions are usually manifested by a “tingling” feeling, often around the mouth or fingers. Other manifestations may include muscle discomfort, muscle twitching or spasms, and/or the presence of an “unusual taste” in the mouth. Return of the donor’s blood may be temporarily slowed by reducing the citrate infusion rate (CIR) or the procedure may be halted if symptoms occur.

Complications such as blood loss, hemolysis, air embolism and/or blood clotting, may be associated with improper operating conditions.

Directions for Use:

Refer to Operator’s Manual for the AMICUS separator for the installation procedure and complete instructions for use and plasma volumes for platelet storage.

Dispose of waste in appropriate biohazard container or according to local regulatory requirements.

Tare Weights:		
Below are the approximate tare weights of the platelet container and associated kit components. Refer to the appropriate illustration for the tare weights for platelet product calculations.		
Approximate Weights (grams)		
Two (2) platelet storage containers with: Two (2) sample pouch assemblies, Two (2) roller clamps and tubings, One (1) Y-Junction, and One (1) 41.9 cm tubing	77	
Two (2) platelet storage containers with: Two (2) sample pouch assemblies, Two (2) roller clamps and tubings, and One (1) Y-Junction	73	
One (1) platelet storage container with: One (1) sample pouch assembly, One (1) roller clamp and tubing, and One (1) Y-Junction	37	
One (1) platelet storage container with: One (1) sample pouch assembly	31	
One (1) platelet storage container	27	
Plastic tubing	0.12 g/cm	
One (1) platelet storage container	39	

Symbols with Definitions:

- Caution, consult instructions for use
- Sterilized by irradiation. Sterile fluid path.
- Non-pyrogenic fluid path
- Do not reuse
- Do not use if the product sterile barrier system is compromised
- Do not vent
- Contains or presence of phthalate: Diethylhexylphthalate (DEHP)
- This marking reflects compliance with the Directive on Medical Devices
- Code

- Lot
- Expiration date
- Manufacturer
- Manufacturing facility/ Manufactured by
- Double needle apheresis protocol
- Fragile
- This way up
- Recyclable

© 2018 Fresenius Kabi AG.
All rights reserved.

Document No.: 472312918
Revision: A
Change Order: CO-034005
Title: PRODUCT INSERT AMICUS APHERESIS KIT - DOUBLE NEEDLE WITH PLATELET ADDITIVE SOLUTION CONNECTOR X6R2302

Page: 2 of 6
Last Modified: 06 Sep 2018 JML
Color: Black

PROOFREADING INSPECTION / RELEASED ARTWORK

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Referencia: United States Pharmacopeia, Avisos generales. United States Pharmacopeial Convention, Inc. 12601 Twinbrook Parkway, Rockville, MD.

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Indian Pharmacopeia, Avisos generales.

Indian Pharmacopeia Commission, Ghaziabad, India

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