Hip Arthroscopy Instruments

It is important to read the Instructions for Use in its entirety prior to use. Failure to properly follow the instructions may lead to serious surgical consequences.

R_x Only

Manufacturer:

T.A.G. Medical Products Corporation LTD.

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English

Table of Contents

3	Description / Intended Use
3	Contraindications
3-4	Warnings
4	Precautions
4	Using the Handle, Hip SwitchStick
4	Using the Hip CordCutter
5	Replacing the Shear Pin
5	Fragment Removal
5-9	Reprocessing
9	Storage
9	Service Life
10	Symbols

Description / Intended Use

- The Hip Arthroscopy Instruments (hereafter hip instruments) are intended to be used in operative arthroscopic procedures in the hip joint.
- The hip instruments consist of reusable instruments:

REF	Catalog Number	Description
	214009	Hip Suture Passer Up
	214013	Hip Suture Passer Straight
	214131	Hip Probe 90
	214132	5.5 Hip Cannulated SwitchStick
	214133	6.0 Hip Cannulated SwitchStick
	214134	4.0 Hip Cannulated SwitchStick
	214135	5.5 Slotted Cannula
	214136	6.0 Slotted Cannula
	214137	Hip Biter 15 Up
	214138	Hip CordCutter
	214139	Hip Biter Straight
	214142	Hip Curette 20 Bent
	214143	Hip Grasper 20 Left
	214144	Hip Grasper 20 Right
	214145	Hip Knot Manipulator
	214146	Hip Suture Manipulator
	214147	Hip Grasper Straight
	214148	Hip TissLib Blade 20 Down
	214151	Handle, Hip SwitchStick
	214153	5.0 Hip Cannulated SwitchStick
	214154	4.5 Hip Cannulated SwitchStick

Contraindications

The hip instruments are NOT intended for use where one or more of the following conditions exist:

- Pathological conditions of the tissue to be treated.
- Physical conditions that would retard healing, such as blood supply limitation and infection.
- Conditions which tend to limit the patient's ability or willingness to follow instructions during the healing period.

⚠ Warnings

- These instructions are not meant to replace proper training in surgical technique or hip arthroscopy. It is imperative that the surgeon and operating theater staff are fully conversant with operative arthroscopic procedures in the hip joint prior to the hip instruments.
- Inspect the instruments prior to each use and after each cleaning / sterilization cycle for integrity and compatibility in order to ensure proper functionality and safety. Do not use damaged, worn, or bent instruments, as these may not function properly and may harm the patient or user.

- Do not attempt to bend, straighten, or repair the hip instruments, as this will compromise instrument integrity and potentially harm the patient or user. For replacement or repair, contact your sales representative.
- Handle hip instruments with care. Some hip instruments have sharp edges that can puncture skin.
- The hip instruments are provided NON-STERILE and must be cleaned, disinfected, and sterilized prior to the first and every subsequent use. For proper cleaning, disinfection, and sterilization, follow the reprocessing instructions provided in this guide.

Precautions

- Caution: Federal law (USA) restricts these instruments to sale by, or on the order of, a physician.
- Do not use the hip instruments for any purpose other than their intended use, as this may result in damage.
- Do not expose the hip instruments to highly alkaline conditions or hypochlorite solutions, as these will promote corrosion and product damage.
- Do not clean delicate instruments with other instruments, as they may sustain damage from coarse handling.
- End of useful instrument life is generally determined by wear or damage from handling or surgical use. Replace all dull instruments.

Using the Handle, Hip SwitchStick

The Handle, Hip SwitchStick provided in the Hip Arthroscopy Set facilitates inserting a switching stick into the surgical site.

1. Slide the entire handle over the switching stick before inserting the switching stick into the surgical portal.
2. Squeeze the handle to grip the switching stick for insertion or to aid extraction.
3. When removing the handle while the switching stick is within the joint, use care to keep the switching stick in position.
4. Remove the switching stick from the handle prior to reprocessing.

Using the Hip CordCutter

The Hip CordCutter is operated in three steps:

1. Load suture into the end of the device with handles fully extended.
2. Capture the suture with handle partly actuated.
3. Fully actuate the handle fully actuated.

Replacing the Shear Pin

Some instruments with actuating loop handles contain a shear pin that connects the handle to the central operating rod. The shear pin is designed to fail prior to a more serious failure mode, such as the jaws breaking apart from excessive force. If the shear pin breaks, contact Mitek for replacement.

Fragment Removal

In the rare event of instrument breakage, fragments can be located visually or by use of radiographic imaging equipment. Fragments can be removed manually through the incision site.

Reprocessing

These reprocessing instructions are provided in accordance with ISO 17664. While they have been validated by the manufacturer of the medical device as being capable of preparing the device for re-use, it remains the responsibility of the processor to ensure that the reprocessing as actually performed, using equipment, materials, and personnel in the reprocessing facility, achieves the desired result. This normally requires validation and routine monitoring of the process.

Limitations on Reprocessing

1. Repeated processing has minimal effect on instrument life and function. End of useful instrument life is generally determined by wear or damage from handling or surgical use.
2. Do not leave the instruments in solutions longer than necessary, as this may accelerate normal product aging.
3. Damage incurred by improper processing will not be covered by the warranty.

Containment and Transportation

1. Reprocess the device as soon as reasonably practical following use.
2. Transport the device in an appropriate manner to prevent damage.

Warnings

1. The following cleaning, disinfection and sterilization methods have been validated. Failure to follow reprocessing instructions as provided may result in incomplete decontamination.
2. Wear appropriate protective equipment (gloves, eye protection, etc.).
3. The Hip CordCutter must be disassembled prior to cleaning, disinfection, and sterilization, or incomplete decontamination will result.

Disassembling the Hip CordCutter

The Hip CordCutter must be disassembled prior to cleaning and sterilization.

1. Depress the button on the top of the handles.
2. Pull the central cutting rod out of the cannulated shaft.
3. Prior to use, reassemble the Hip CordCutter carefully. The proximal end of the central rod must be oriented as such that the notch is facing the handle. Verify that handle and trigger actuation moves the central rod.

Cleaning

Warnings

- Cleaning should be performed as soon as possible⁽¹⁾ after each procedure, before blood, saline and debris are dry.
- Do not use metallic brushes, scrub pads, or other abrasive cleaning aids, as these are likely to damage the instruments.
- Do not use highly alkaline solutions or harsh chemicals such as chlorine, caustic soda, and organic or ammoniated acids or solvents (such as acetone), as these are likely to damage the instruments.
- Follow the instructions and warnings issued by the suppliers of any cleaning and disinfection agents and equipment used.
- Complex devices, such as those with tubes, hinges, retractable features, and matte or textured surface finishes, require special attention during cleaning. Manual pre-cleaning of such device features is required before automated cleaning.

⁽¹⁾ 30 minutes was used for the Cleaning Validation.

Manual cleaning

1. Rinse instruments with distilled or demineralized water.
2. Prepare an enzymatic detergent⁽²⁾ according to the manufacturer's instructions.
3. Submerge the instruments in the enzymatic detergent solution for 20 minutes. Use a 5 ml syringe and needle to flush enzymatic detergent through lumens, and use a soft bristled brush to scrub the lumens, crevices, and slots.
4. Rinse instruments with distilled or demineralized water, agitating for one minute; repeat rinse twice.
5. Dry all instruments thoroughly before storage.
6. Inspect the instruments after each cleaning cycle to ensure proper functionality and safety:
 - 6.1 If you suspect a problem, do not use the instrument.
 - 6.2 Inspect all components for cleanliness: activate all parts to their maximum positions while inspecting for blood or debris. Place a pipe cleaner down any cannulation or crevice, and then inspect the pipe cleaner for blood or debris. If blood or debris is present, repeat the above cleaning procedure.

⁽²⁾ Cleaning efficacy was validated using ENZOL® at one ounce per gallon of water at 35-40°C (95-104°F)

Automated Cleaning

1. Prepare an enzymatic detergent⁽³⁾ according to the manufacturer's instructions.
2. Perform pre-cleaning by flushing an enzymatic detergent through the lumens using a 5 ml syringe and needle and using a soft bristled brush to scrub the lumens, crevices, and slots.
3. Rinse the instruments with distilled or demineralized water until no visible detergent is present.
4. Place the instruments in the washer in a manner to facilitate draining.
5. Clean in a validated automatic washer using the "INSTRUMENTS" cycle and detergents⁽⁴⁾ intended for use with automated cleaning. Effective automated cleaning can be achieved using the following parameters:
 - 5.1 Cold water wash (minimum of 35 seconds at 27°C / 80°F)
 - 5.2 Hot water wash (minimum of 2 minutes at 60°C / 140°F)
 - 5.3 Cold water wash (minimum of 45 seconds at 27°C / 80°F)
 - 5.4 Hot water wash (minimum of 2 minutes at 60°C / 140°F)
 - 5.5 Thermal rinse (minimum of 1 minute at 82.2°C / 180°F)
 - 5.6 Drying (minimum of 6 minutes at 82.2°C / 180°F)
6. Inspect the instruments after each cleaning cycle to ensure proper functionality and safety:
 - 6.1 If you suspect a problem, do not use the instrument.
 - 6.2 Inspect all components for cleanliness: activate all moving parts to their maximum positions while inspecting for blood or debris. Place a pipe cleaner down any cannulation or crevice, and then inspect the pipe cleaner for blood or debris. If blood or debris is present, repeat the above cleaning procedure.

⁽³⁾ Validation was performed using ENZOL®.

⁽⁴⁾ Validation was performed using Renuzyme Plus® and Tec Wash III®.

Disinfection

Warnings

Disinfection alone is not adequate to render the instruments sterile. Sterilization must be performed before each procedure.

1. After cleaning, disinfection should be performed by immersion in 2.4% CIDEX® ⁽⁵⁾ (or equivalent) with a minimum immersion time of 45 minutes at 25°C / 77°F following the manufacturer's instructions.
2. Rinse with distilled or demineralized water, agitating for 30 seconds; repeat rinse twice.
3. Dry all instruments thoroughly before storage.

⁽⁵⁾ Validation was performed using CIDEX®.

Maintenance, Inspection, and Testing

Inspect the instruments for damage, such as deformed shafts or tips. Ensure proper functionality. If a problem is observed or suspected, contact your sales representative for replacement or repair.

Packaging for Sterilization

Use no more than two layers of sterilization wrap.

Sterilization

Warnings

- Sterilization is not a substitute for cleaning. Instruments must be thoroughly cleaned prior to sterilization, or incomplete decontamination will result.
- The sterilization parameters provided for use with a sterilization tray might differ from those provided for the instruments alone. When sterilizing the instruments inside a tray, select the more stringent set of parameters.
- These devices were not validated for ethylene oxide (EO) or gas plasma sterilization methods. Use only steam (autoclave) sterilization as outlined in these instructions, or incomplete sterilization may result.
- Drying time depends on several variables, including altitude, humidity, type of wrap, preconditioning, size of chamber, mass of load, material of load, and placement in chamber. Users must verify that drying time set in their autoclave yields dry surgical equipment.
- The adequacy of any sterilization procedure must be tested. To substantiate the process and its reliability and reproducibility, it is critical that appropriate process parameters be validated for each facility's sterilization equipment and product/load configuration by persons who have training and expertise in sterilization processes.
- The following autoclave processes have been validated by the manufacturer; however, the sterility assurance level for any given sterilization cycle is dependent on the actual conditions selected by the user and the correct calibration of the user's sterilization equipment and instrumentation. Because these factors are not within the manufacturer's control, the user is responsible for the sterility assurance level.

Gravity - Displacement Steam Sterilization:

❖ Sterilizer Type:	Gravity
❖ Minimum Temperature:	132°C (270°F)
❖ Minimum Exposure Time:	10 minutes
❖ Minimum Dry Time:	40 minutes

Pre-vacuum Steam Sterilization:

❖ Sterilizer Type:	Pre-vacuum
❖ Minimum Temperature:	132°C (270°F)
❖ Minimum Exposure Time:	4 minutes
❖ Minimum Dry Time:	15 minutes

Pre-vacuum Steam Sterilization (alternate cycle):

❖ Sterilizer Type:	Pre-vacuum
❖ Minimum Temperature:	134°C (273°F)
❖ Minimum Exposure Time:	3 minutes ⁽⁶⁾
❖ Minimum Dry Time:	15 minutes

⁽⁶⁾ Exposure time of 18 minutes is permitted (as mandatory in some European countries).

Do not exceed 140°C (284°F) during sterilization.

Storage

- Between uses, store instruments in a clean, dry environment.
- To prevent damage, minimize contact between instruments and with other equipment.

 **Warning:** These storage parameters are not intended to preserve instrument sterility. If necessary, reprocess the instruments again prior to use.

Service Life

The instruments should be replaced after one year.

Symbols

See Instructions for Use

Non-sterile

Manufacturer

Made In Israel

Made In Israel

Lot Number

Authorized Representative in the European Community

R_x Only

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician

Number of units

Catalog number

Date of Manufacture

CE Mark (Class I)

ENZOL® and CIDEX® are registered trademarks of Advanced Sterilization Products.

Renuzyme Plus® and Tec Wash III® are registered trademarks of Getinge USA, Inc.

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Renzyme Plus® dan Tec Wash III® merupakan merek dagang terdaftar dari Getinge USA, Inc.