INSTRUCTION GUIDE FOR A.B DENTAL IMPLANTS

DISCLAIMER

A.B Dental Devices products are intended for use only by certified dentists and authorized personnel with specific implant training. A.B Dental Devices implants are used for two-stage and one-piece implantation processes. The implants and abutments are made of titanium alloy. A.B Dental Devices implants are delivered in sterile, sealed containers. They are supplied with the understanding that only AB Dental Devices surgical instruments, which complement each implant, will be used during surgery. If these conditions are not met, the manufacturer will refuse to accept responsibility.

MANDATORY INITIAL INVESTIGATIONS

Patient examination; Patient's medical history; Clinical examination of patient's hygiene, teeth, occlusion, periodontium; Biological observations; Radiographic evaluation: CT scan, intra-oral X-rays, pan-oral, etc. Lack of adequate practitioner training is one of the major factors influencing the success of implant surgery and subsequent long-term patient health.

DESCRIPTION

A.B Dental Implants [Hereinafter: "A.B."] are made of titanium alloy to be placed surgically into the upper or lower jaw in order to replace the patients dental root and to provide a stable base for the prostheses. The abutment, made of zirconium, plastic or titanium, is used in the prosthetic procedure.

THE PRODUCT

A.B implants are designed for use in dental implant surgery. It replaces the natural tooth root. The implant is surgically inserted into the upper or lower alveolar bone. The implant can restore the injured tooth by connecting abutment to the alveolar bone after the implants osseointegration process. The abutment system is intended for use as an aid in prosthetic restoration.

INSTRUCTION GUIDE

Use in accordance with general implant surgical guides. After implant placement, the surgeon should evaluate bone loss, implant percussion analyses and radiological marginal bone loss. Unstable implant and/or more than 50% bone loss should indicate implant's extraction. When the surgeon uses short type implants, it is recommended to: use two stage prosthetic fabrication, connect the short implant to a regular one, use as wide a diameter abutment as possible, avoid immediate loading and provide a longer period of time for osseointegration. Ultra wide implants 5mm or more, should replace molar teeth only. For over denture, it is recommended to inter-connect the implants. Over denture should be placed following appropriate healing time. There are different surgical kits and protocols for the different implants. For further information please see the catalog. For A.B prosthetic components and instructions - see the prosthetic manual.

CONTRAINDICATIONS

Customary observations should be made of the contraindications associated with implant materials used in oral surgery. First, the patient's general health and suitability for oral surgery must be assessed by the general practitioner. Insufficient height and/or bone width, poor bone quality, poor oral hygiene, heavy smoking and tobacco use, systemic and hematologic disorders and uncontrolled diabetes, can lead to implant failure and are contraindications for use with A.B Implants.

WARNINGS

For safe and effective use of A.B Implants, it is strongly suggested that professional training be undertaken since the surgical techniques required to place dental implants are highly specialized and complex

procedures. Improper patient selection and technique can contribute to implant failure and/or loss of supporting bone. A.B implants are intended for use only in the indicated applications. Dental implants must not be altered in any way. The use of electro-surgical instruments or laser around metallic implants and their abutments is not recommended due to the risk of electric shock and/or burns. Implant mobility, bone loss or chronic infection may indicate implant failure. If implant becomes contaminated by patient's body fluids in any way, the implant cannot be used in any other patient.

PRECAUTIONS

The surgical techniques required to place endosseous dental implants require specialized and complex procedures. Formal training for placement of implants is recommended.

Important: Determine local anatomy and suitability of the available bone for implant placement. Thorough screening of prospective implant candidates must be performed. Visual inspection as well as panoramic and peri-apical radiographs is essential to determine location of anatomical structures, occlusal conditions, periodontal status, and adequacy of bone. Lateral cephalometric radiographs, CT scans, and tomography may also be beneficial. Adequate radiographs, direct palpation, and visual inspection of the implant site are necessary prior to treatment, planning and use of A.B implants.

PRECAUTIONS GUIDELINES

Do not expose the packed product to contaminated materials or tissues. Do not use the product if packaging seal is in question or if packaging label is damaged.

CAUTION FOR PATIENT

Keep good oral hygiene. Do not apply too much bite force until final prosthetic placement.

SIDE EFFECTS

Some of the complications (loss of implant anchorage, prosthesis etc) are possible occurrences after surgery. Lack of quantity or quality of remaining bone, infections, poor patient oral hygiene or cooperation, patient discomfort, implant mobility, local soft tissue degeneration, and unfavorable implant placement or alignment are some potential causes for loss of anchorage.

SURGICAL COMPLICATIONS

The surgical procedure is not without complication, including local swelling, temporary sensitivity, hemorrhage, edema or bleeding. There is a risk for loss of sensation in lower lip and chin (after implant placement into the lower jaw) and around the nose (after implant placement into upper jaw). This loss of sensation can be permanent (rare). Injury to the gingival tissue, allergic reaction of tissues or infections, can occur. Usually, these reactions response to local medical treatment.

STERILIZATION

A.B Implant, implant clip carrier and cover screws are shipped cleaned and after sterile under Gama rays. They are delivered packed and ready for use. Abutments and instruments not delivered sterile should be sterilized prior to use using stream sterilization at 132 °C for 15 minutes.

WARNING

Sterility of the product is guaranteed only if package is complete (not opened nor damaged). Product information: Please refer to the label for product name, code, size, manufacturing date, LOT number, and expiry date.

PRODUCT LIMITED USE

Implants and prosthetic components are for single use only! Not for sterilization and reuse!

GUARANTEE

A.B Dental Devices undertakes, for 15 years after manufacturing date, to replace any implant or restoration component bearing the name A.B Dental Devices. This guarantee does not cover the surgical accessories or tools used with the implant. A.B Dental Devices stipulates that to be covered by this guarantee, A.B Dental Devices implants must have been implanted according to the protocols and procedures described in this manual. The practitioner must have undergone formal training in implant surgery. A.B Dental Devices cannot be responsible for any bodily or material damages whatsoever. This guarantee ensures that if necessary, the vendor will replace the damaged implant or cover the expenses of purchasing the replacement. It is important for the implant clinician to take note of the batch number of the implant, should there be any warrant for replacement. Warranty replacement is effected under condition of following protocol and procedures of implantation, described in instruction for use.

NOTE: For general details on implants, prosthetic components and the correct way of use, see catalog, prosthetic manual and other A.B information sources.

AB Dental Devices Ltd.
Ashdod, Israel
www.ab-dent.com

	Symbol for Use until date
	Symbol for Do not reuse
	Symbol for See instructions for use
	Symbol for method of sterilization using gamma radiation
	Product complies with requirements of directive 93/42/EEC for medical devices
	Symbol for Lot number

CEpartner4U B.V., Esdoornlaan 13, 3951 DB MAARN, The Netherlands.
Tel: +31 343 442 524 Fax: +31 343 442 162

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CEPartner4 B.V., Esdoorlaan 13, 3951 DB MAARN, The Netherlands.