IFU, GB: 210mm x 297 mm

IFU, RU: 210mm x 297 mm

BEROCATH®

LATEX FOLEY BALLOON CATHETER

Description of the product:

Latex Foley Balloon Catheters are flexible tubes intended to be inserted into the urinary bladder via the urethra for drainage and/or collection of urine, irrigation of the bladder and instillation of fluids. The balloon retains the catheter inside the bladder once inflated.

Latex Foley Balloon Catheters consist of shaft, drainage funnel, inflation funnel, irrigation funnel (if present), balloon and valve. Their smooth tapered tip facilitates insertion, and their silicone coated surface reduces allergic reactions. They are color-coded for visualization of the size. Size and balloon capacity are marked both on the package and on the funnel of the catheter.

Materials used:

- Medical grade natural latex, silicone, PVC (valve)

Intended use/Indications:

Latex Foley Balloon Catheters are passed through the urethra during urinary catheterization and into the urinary bladder for drainage of urine and/or instillation of fluids into the bladder. They are used in departments of urology, internal medicine, surgery, obstetrics and gynecology for drainage of urine and medication. They are also used for patients suffering from moving with difficulty or being completely bed-ridden.

Contraindications:

- Not to be used in patients with known hypersensitivity to any of the materials used
- Acute urethritis
- Acute prostatitis
- Acute epididymitis

Adverse reactions:

- Bleeding into or around the catheter
- Drainage of only little or no urine despite sufficient fluid intake
- Leakage of large amounts of urine around the catheter
- Urethral swelling around the catheter
- Urinary tract or renal infections
- Pathologic consistence of the urine (viscous, cloudy, strong smell)
- Urethritis
- Fever, chills

Instructions for use:

- Apply aseptic techniques.
- Lubricate tip and shaft of the catheter. Do not use ointments and lubricants based on petroleum.
- Carefully insert the catheter into the bladder via the urethra. Urine flow normally indicates that the catheter tip has reached the bladder.
- Once there is urine flow, continue catheter insertion for further 3 cm to ensure that the balloon is also positioned inside the bladder.
- Inflate the balloon by using a Luer syringe without needle. The recommended inflation volume is marked on the funnel of the catheter and must not be exceeded.
- For deflation of the balloon, cut off the inflation funnel above the valve or push an empty Luer syringe barrel (without needle and without plunger) into the valve.
- The catheter is to be removed /replaced after 14 days at the latest. The accumulated duration of use can exceed 30 days as the product is intended for long-term use.

Warnings:

- Read instructions and warnings before use, improper use can result in serious or fatal illness or injury.
- The product contains natural rubber latex which may cause allergic reactions.
- Use of this product is restricted to trained physicians.
- Visually inspect and check the product and package integrity prior to use and also the expiry date. Improper transport and handling may cause structural and/or functional damage to device or packaging.
- Do not use after the expiry date.
- Do not use if the product is damaged or irregularly shaped.
- The product is sterile and non-toxic if the package is not opened or damaged. Do not use if the individual package is damaged or has been opened.
- Use the product immediately after opening the individual package.
- The product is to be used under aseptic conditions, aseptic techniques are to be applied.
- The balloon is to be tested for leakage and deflation prior to use.
- Do not use ointments or lubricants based on petroleum. Use medical grade silicone oil or water-soluble lubricants.
- Inflate the balloon with sterile distilled water.
- Do not use a needle to inflate the balloon. Use a Luer syringe without needle.
- Do not over-inflate the balloon. Refer to the unit package or funnel of the catheter for the balloon capacity.
- Do not clamp or kink the catheter shaft. It may damage the catheter and impede deflation.
- The catheter must be removed or replaced after 14 days at the latest.
- For single use and for single patient use only. Discard after single use. Do not re-use, recycle, reprocess or re-sterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in injury, illness or death of the patient. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Storage conditions:

Store in a dark, cool and dry place. Avoid excessive heat, protect from direct sunlight and moisture. Do not expose to ultraviolet and fluorescent light. Store in a manner preventing crushing.

Disposal:

Discard the product after single use in accordance with applicable national and institutional regulations for biological hazardous waste.

BEROMED GMBH HOSPITAL PRODUCTS DISCLAIMS ANY RESPONSIBILITY FOR POSSIBLE CONSEQUENCES RESULTING FROM IMPROPER USE.

 Manufacturer	 Use-by date	 Temperature limit	 Do not re-sterilize
 Catalogue number	 Fragile, handle with care	 Contains or presence of natural rubber latex	 Consult instructions for use
 Batch code	 Keep dry	 Do not use if package is damaged	 Caution
 Date of manufacture	 Keep away from sunlight	 Do not re-use	 Sterilized using ethylene oxide
 0123			

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