
Description / Descripción / Descrizione / Beschreibung / 产品描述 / Beschrijving / Leírás / Kuvaus / Descrição / Açıklama / Popis.

A. Echelle de compression - décompression / Compression-decompression scale / Escala de compresión-descompresión / Scala di compressione-decompressione / Kompressions-Dekompressionsskala / 加压/减压刻度 / Schaalverdeling voor compressie-decompressie / Kompressziós-dekompressziós skála / Kompresso- ja dekompressoasteikko / Escala de compressão-descompressão / Kompresyon dekompresyon skalasi / Kompresni/dekompresni stupnice.

B. Molette de compression - décompression / Compression-decompression wheel / Rueda de compresión-descompresión / Ghiera di compressione-decompressione / Kompressions-Dekompressionsrad / 加压/减压轮 / Compressie-Decompressiewiel / Kompressziós-dekompressziós kerék / Kompresso- ja dekompressiopyörä / Roda de compressão-descompressão / Kompresyon-dekompresyon tekerleği / Kompresni/dekompresni kolečko.

C. Bouton de sécurisation de la décompression / Decompression safety button / Botón de seguridad de descompresión / Ghiera di sicurezza di decompressione / Dekompressionssicherheitsknopf / 减压安全按钮 / Veiligheidsknop voor decompressie / Dekompresszió biztonsági gombja / Dekompression suojauspainike / Botão de segurança da descompressão / Dekompresyon emniyet düğmesi / Dekompresni bezpečnostní tlačítko.

D. Indicateur horaire / Time indicator / Indicador de tempo / Indicatore di tempo / Zeitanzeiger / 时间指示器 / Tijdsindicator / Időjelző / Ajastin / Indicador de tempo / Zaman göstergesi / Ukazatel času.

E. Embase transparente / Transparent base / Base transparente / Base trasparente / Transparentes Unterteil / 透明底座 / Transparente basis / Átlátszó alap / Läpinäkyvä alusta / Base transparente / Šeffaf taban / Průhledná základna.

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TVA N° FR 01 317 883 999

Laura JUVEN
[Signature]

Sterile (sterilised with ethylene oxide). Do not re-sterilize. For single use. Do not use if the packaging is open or damaged. Use before the expiry date shown on the packaging. Read the instructions before use. Store away from light, in a cool and dry place.

Indications

The Seal-One® compression device is used for coronary angiography and angioplasty procedures practiced for the radial approach to compress the puncture site.

Warning

Federal (U.S.A) law restricts this device to use or sale by or on the order of a physician. Re-use of the Seal-One® is strictly prohibited. If the product is reused its functionality would be entirely altered. In addition, the iatrogenic risk from inter-patient cross-transmission would be high.

The Seal-One® must be used by trained medical staff.
 Do not leave the patient unsupervised when using the Seal-One®.
 Do not contaminate the Seal-One® compression pad with blood before positioning the device.
 The Seal-One® device should not be left in place for more than 5 hours without physician approval.
 Excessive compression may lead to arterial occlusion.
 Do not apply pressure to the upper part of the compression pad.
 Do not use the device if damaged or if it hasn't been properly adjusted to the wrist.
 After use, dispose of product and packaging in accordance with hospital administrative and/or local government policy.

Patient Information

The patient should be instructed not to move the Seal-One® system nor modify the settings due to the risk of bleeding and bruising. The patient should call the medical staff immediately if he/she experiences pain or bleeding.

The patient must be informed to avoid any impact in order not to damage the device.

The patient must also be informed not to touch the device to avoid changing the compression or the time indicator and to reduce the risk of removing the device unintentionally.

Operation

The Seal-One® may be used on the right or left wrist. The device placement time should be set on the time indicator, which should subsequently be locked. Wrist strap is secured by a strap which is clipped. The wrist strap can be folded back on itself to shorten it.

Compression is applied by turning the compression knob clockwise.

Decompression is done by pressing simultaneously on the safety button and turning the compression wheel anticlockwise.

It is possible to write on the wrist strap with ballpoint pen or indelible felt pen, to indicate the times at which the device was checked and steps during decompression. The dedicated positions for this are shown by the following symbols at the ends of the wrist strap:

Device placement

Before positioning the Seal-One® on the patient wrist, check that the compression wheel is correctly set to 0 (on the compression scale).

Then, set the device placement time by turning the time indicator. Push the indicator backwards towards the transparent base to lock it. (A)

Withdraw the introducer sheath approximately 2-3cm in order not to interfere with positioning the Seal-One®.

Position the Seal-One® ensuring that the central landmark is 0.5cm proximal to the puncture site. (B)

Adjust and fix the wrist strap and then clip the safety strap (the wrist strap snug and the safety clip must be on the practitioner side). (C)

The wrist strap may be folded back on itself if it is too long.

Arterial compression

Turn the compression wheel clockwise to initiate compression to level 3. (D)
 Check that the compression wheel is back-locked in position by turning it anti-clockwise until it locks. If the wheel is not back-locked during decompression, manually lift the safety button to reset the back-lock. If the back-locked position is inoperative, do not use the

device.
 Remove the introducer sheath completely.
 Increase the compression level until bleeding stops.

Final compression level adjustment

Decompress by pressing the safety button (E) and turning the compression wheel simultaneously anticlockwise. Discontinue decompression once a drop of blood is visible. Compress again until bleeding stops.

Note: if a drop of blood is not visible, position the compression level at 3.

After positioning the Seal-One®, check that the radial artery is not occluded by excessive compression from the device, using a Doppler or pulse oximeter. If so, reduce the compression level.

Under usual conditions of use the Seal-One® radial compression device is designed only to compress the radial artery (the ulnar artery is not usually compressed and venous return is not blocked).

Decompression time is patient dependent and varies according to co-morbidities, procedure duration, anticoagulation treatment and the size of the arterial puncture site.

Adjust the level and duration of pressure according to these different parameters.

The Seal-One® device should not be left in place for more than 5 hours without physician approval. Excessive compression may lead to arterial occlusion.

Decompression steps

The artery may be decompressed in steps. To do this press the safety button (E) and turn the compression wheel simultaneously anticlockwise (F) until the desired lower level of compression is reached. If there is any resistance on the safety button during decompression,

gently turn the compression wheel clockwise.

Bleeding may occur when reducing the pressure. Readjust the compression level and wait before decompressing further.

Upon reaching the zero decompression level, wait for 15 minutes minimum and remove. Before removing the device, ensure that no bleeding is present.

Side effects

Depending on the patient state and level of compression applied with the device, side-effects such as occlusion of the radial artery, hypodermal haematoma, haemorrhage, pain, or numbness may occur. Control bleeding and optimally adjust the level of compression as a result.

If the skin becomes itchy or red during compression, stop using the Seal-One® and treat appropriately.

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S.A.S au capital de 1 682 898,50 €