Poches de transfert MacoPharma

FRANCAIS

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MacoPharma transfer bags

ENGLISH

Device intended for the transfer of blood products (whole blood, platelets, plasma, serum, cryoprecipitate).
The device may be centrifuged with the exception of packs with a nominal volume of greater than 1L.

PRECAUTIONS

- Bring to a temperature between 15°C and 30°C prior to use.
- Avoid contact with sharp objects.
- DO NOT USE if the overwrap and/or bag show any visible signs of deterioration.
- The type of sterilisation will be detailed on the bag label.
- Products sterilised by Irradiation may be slightly coloured without affecting the integrity of the bag.
- DO NOT USE 10 days after opening of the packing.
- All or part of this medical device is made of PVC plasticized with DEHP. According to some studies, DEHP could potentially be harmful to the reproductive system of male foetuses. The prescriber is solely responsible for choosing to use this device on women who are either pregnant or breast-feeding, or on young male infants. Nevertheless, DEHP-plasticized PVC is in compliance with the European Pharmacopoeia.
- Risk linked to the reuse: Septic risk.
- All manufacturing processes and all components in contact with the donors, the users and the blood components dedicated to the patients are latex free

INSTRUCTIONS FOR USE

- Place the bag on a clean area, ideally in a laminar flow cabinet.
- Open the transparent overwrap using the peelable corner. (A)
- Label the bag.
- ✓ **Bags with a spike: (B)**
- Clamp the tubing between the spike and the empty bag.
- Remove the spike protector cover by applying gentle pressure at the base of the cover and twisting. (B)
- Place the spike into the outlet port of the bag containing the product to be transferred.
- Twist the spike to secure it in position. (B')
- Check that the spike is inserted correctly.
- Transfer the product by gravity without pulling on the tubing (the transfer bag should be resting flat on the bench).
- Heat seal the tubing.

✓ Bags with a luer: (C)

- Clamp the tubing between the luer and the empty bag.
- Remove the protector covers from the luers which are to be connected.
- Screw together the luers.
- Check that the luers are connected correctly.
- Transfer the product by gravity without pulling on the tubing (the transfer bag should rest flat on the bench).
- Heat seal the tubing.

✓ Bags for connection by Sterile Docking (SCD): (D)

- Clamp the transfer line tubing.
- Connect the tubes following the SCD manufacturer's instructions.
- Check the integrity of the seal.
- Transfer the product by gravity without pulling on the tubing (the transfer bag should rest flat on the bench).
- Heat seal the tubing.
- ✓ **Processing and Storage**
- The following steps will be dependant on the bag configuration used. Follow the established validated procedures.
- For centrifugation, maintain the bags upright in the centrifuge buckets ensuring the bags and tubing are correctly bagged. Use foam pads and weights if necessary.

- Freezing in temperatures below -40°C is not recommended.
- Handle the bags carefully during transportation, processing and storage.
- ✓ **Opening Outlet Port Protectors (E)**
- Using both hands, tear the port protector. (F)
- Hold down the protector pieces with one hand.(G)
- With the other hand, place the spike into the port.
- Twist the spike to secure it in position.

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PLASMAFLEX filter for plasma filtration

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BSV 8000X

ENGLISH

Intended for the filtration of plasma from whole blood or apheresis before freezing.

WARNING

- Bring to a temperature between 15°C to 30°C before use.
- Avoid contact with sharp objects.
- Inspect overwrap and pack.
- DO NOT USE if they present any visible signs of deterioration.
- NEVER FOLD OR PRESS the filter.
- Temperature of filtration : room temperature (20°C + 2°C)
- Volume of plasma: up to 650ml.
- Risks linked to the reuse: septic risk, loss of efficacy.
- All or part of this medical device is made of PVC plasticized with DEHP. According to some studies, DEHP could potentially be harmful to the reproductive system of male foetuses. The prescriber is solely responsible for choosing to use this device on women who are either pregnant or breast-feeding, or on young male infants. Nevertheless, DEHP-plasticized PVC is in compliance with the European Pharmacopeia.
- No MacoPharma product, component or manufacturing process involves the use of latex or any other natural rubber compound.

INSTRUCTION FOR USE

✓ Set up

- Close the clamp
- Sterile connect the plasma bag to the filtration set.

✓ Filtration

- Hang the plasma bag by the hanging point at the base of the bag allowing the transfer bag to hang freely.
- Open the clamp
- Allow the filter to fill
- Filtration starts automatically
- The filtration is finished when the filter becomes flat.

NOTE: If either face of the filter remains white, the plasma unit has not been leucoreduced and must be discarded.

✓ Removal of air

- Squeeze the air back into the filter
- Clamp and seal the tube below the filter
- For the triple pack set (BRV), in order to aliquot the same amount of plasma in each pack they should be suspended at the same level by the eyelet at the base of the pack.
- Store the plasma between -30°C et -40°C

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ENGLISH

INDICATION

- Intended for the collection of human blood and the preservation of blood components.

WARNING

- Bring to a temperature between +15°C and +30°C before use.
- Avoid contact with sharp objects.
- Inspect overwrap and pack.
- DO NOT USE if the overwrap and pack present any visible signs of deterioration.
- DO NOT USE if the solution(s) is (are) not clear.
- Do not VENT
- Risks linked to the reuse : Septic risk, Coagulation risk, Loss of Efficacy.
- All or part of this medical device is made of PVC plasticized with DEHP. According to some studies, DEHP could potentially be harmful to the reproductive system of male fetuses. The prescriber is solely responsible for choosing to use this device on women who are either pregnant or breast-feeding, or on young male infants. Nevertheless, DEHP-plasticized PVC is in compliance with the European Pharmacopoeia.
- All manufacturing processes and all components in contact with the donors, the users and the blood components dedicated to the patients are latex free

INSTRUCTION FOR USE

✓ Venepuncture:

- 1) Open the transparent overwrap using the corner flap.
- 2) Place blood pack at a distance of ~50cm below donor arm, ideally on an automated blood mixer.
- 3) Apply tourniquet or pressure cuff (60 mm Hg).
- 4) Thoroughly disinfect the site of venepuncture.
- 5) To remove the needle cap: Hold the hub and twist the cap. Remove the cap in a straight line to avoid any damage to the needle tip.
- 6) Carry out venepuncture.

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- 1) Allow the BACTIVAM pouch to fill.
- 2) While filling the BACTIVAM, do not let the BACTIVAM pouch hang freely. Fill the BACTIVAM.
- 3) Close the clamp (b) and break the breakaway cannula (c) on the donation line with a side to side motion.
- 4) Immediately on starting the collection, thoroughly mix the blood/anticoagulant and continually mix throughout the collection.
- 5) Remove the cap of the VACUUM device, by either twisting the cap, or by pushing the handle of the cap with a thumb to release. Invert the BACTIVAM pouch to ensure any air is directed to the bottom of the BACTIVAM.
- 6) Collect the blood samples using the VACUUM and the appropriate sampling tubes.

✓ Blood collection:

- 1) Ensure that the anticoagulant and the blood are mixed thoroughly immediately upon the start of collection, and also at least every 50mL during the donation, ideally using an automated blood mixer.
- 2) Collect the quantity of blood stated on the pack label.
- 3) Once the required quantity has been reached, close the clamp of the donation line (a) and seal the donation tubing, using a permanent closure device, between the blood pack and the Y-junction of any attached devices.
- 4) Release the pressure cuff or tourniquet, and immediately mix the full blood pack.
- 5) If using an in-line filter, do not strip the donation tubing. Seal it close to the pack (~2cm), using a permanent closure device.
- 6) If required (quality control, crossmatch, pool of buffycoats...), strip the donation line, and mix the pack. Repeat several times. Seal and cut the donation tubing as needed.

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- 1) After donation and sampling, move the SECUVAM close to the hub of the needle.
- 2) Remove the donation needle, then gently lower the needle into the barrel of the SECUVAM device until a click is heard.
- 3) The needle and SECUVAM must be held upright. Ensure that no pressure is put on the donation line at any time, to prevent possible spillage.
- 4) The SECUVAM may be pushed into the barrel of the VACUUM.
- 5) The needle, SECUVAM and any other devices should be disposed of in a safe manner.

✓ Pack handling whilst processing:

- 1) Processing depends on the pack configuration. Process the pack as required by the blood transfusion centre.
- 2) Ensure the bags and tubing are correctly packed for centrifugation.
- 3) To allow fluid transfer, break the cannula at the seal, by bending the cannula with a side-to-side motion.
- 4) After transfer, seal the transfer tubing, using a permanent closure device, and cut to separate.

✓ Handling of filters:

- 1) Do not fold or mishandle the filter(s).
 - 2) If they are centrifuged, filters should be kept flat between transfer bags (with the arrow on the filter pointing to the side), or between the inner wall of the pot and the whole blood bag.
- NOTE: If either face of the filter remains white, the blood components have not been leucoreduced and must be discarded.

✓ Handling of port protectors:

- 1) Using both hands tear the port protector as indicated.
- 2) Hold down the port protector pieces with one hand.
- 3) Using the free hand, take the spike and engage into the port.
- 4) Twist the spike to secure fully in port.

✓ Precautions for disposable set

Disposal of bags should be carried out in accordance with regulations applicable to potentially infected, biological waste, to potentially infected, biological waste.

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MacoPharma BLOOD BAGS

ENGLISH

Intended for the collection of human blood and the preservation of blood components

WARNING

- Bring to a temperature between 15°C and 30°C before use.
- Avoid contact with sharp objects.
- Inspect overwrap and pack.
- DO NOT USE if they present any visible signs of deterioration.
- DO NOT USE if the solution(s) is (are) not clear.
- DO NOT USE 10 days after opening of the packing.
- The PVC complies with the European Pharmacopoeia. The plasticizer is DEHP. The biological products stored and/or circulating in these devices could contain phthalates.
- Risks linked to the reuse: septic risk, coagulation risk, loss of efficacy.
- All or part of this medical device is made of PVC plasticized with DEHP. According to some studies, DEHP could potentially be harmful to the reproductive system of male foetuses. The prescriber is solely responsible for choosing to use this device on women who are either pregnant or breast-feeding, or on young male infants. Nevertheless, DEHP-plasticized PVC is in compliance with the European Pharmacopoeia.
- All manufacturing processes and all components in contact with the donors, the users and the blood components dedicated to the patients are latex free.

INSTRUCTION FOR USE

✓ Venepuncture:

- 1) Peel open the transparent overwrap using the corner flap.
- 2) Place blood pack at a distance of ~50cm below donor arm, ideally on an automated blood mixer.
- 3) Apply tourniquet or pressure cuff (60 mm Hg).
- 4) Thoroughly disinfect site of venepuncture.
- 5) Clamp donation tubing.
- 6) Open needle by holding needle cap and hub. Twist the cap and remove it in a straight line to avoid damaging the needle tip.
- 7) Carry out venepuncture.

✓ "Pre-donation" sampling (if using BACTIVAM device):

- 1) Ensure that the donation tubing is clamped (a), and that the clamp (b) between the BACTIVAM pouch and the donation line is open.
- 2) Break the cannula (c) on the BACTIVAM line and allow the pouch to fill.
- 3) Close the clamp (b) on the BACTIVAM line.
- 4) Open the clamp (a) on the donation tubing and allow the donation procedure to continue, ensuring that the blood pack is mixed immediately the blood flows into it.
- 5) Remove the cap of the VACUVAM device. This may be accomplished by twisting the cap, or with one hand by pushing a thumb against the vane on the cap-top. Invert the BACTIVAM pouch so that any air is away from the VACUVAM.
- 6) Collect the blood samples using the VACUVAM and appropriate evacuated sampling tubes.

✓ Blood collection:

- 1) Ensure that the anticoagulant and the blood are mixed thoroughly immediately upon the start of collection, and also at least every 50ml during the donation, ideally using an automated blood mixer.
- 2) Collect the quantity of blood stated on the pack label.
- 3) Once the required quantity has been reached, clamp and seal the donation tubing, using a permanent closure device, between the blood pack and the Y-junction of any attached devices.
- 4) Release the pressure cuff or tourniquet, and immediately mix the full blood pack.
- 5) For packs equipped with the VACUVAM device, carry out "post-donation" sampling (see below E-F).
- 6) If using an in-line filter, do not strip the donation tubing. Seal it close to the pack (~2cm), using a permanent closure device.
- 7) If required (quality control, crossmatch, pool of buffycoats...), strip the donation line, and mix the pack. Repeat several times. Seal and cut the donation tubing as needed.

✓ "Post-donation" sampling (if using VACUVAM device):

- 1) After the blood donation is complete, ensure that the donation tubing is sealed.
- 2) Remove the cap of the VACUVAM sampling device. This may be accomplished by twisting the cap, or with one hand by pushing a thumb against the vane on the cap-top.
- 3) Break the cannula (d) by bending the cannula with a side-to-side motion.
- 4) Insert appropriate sampling tubes into the VACUVAM to take blood samples.
- 5) If using a pack without VACUVAM, use the "cut and drip" method for sampling.

✓ Use of SECUVAM needle protection system:

- 1) After donation and sampling, move the SECUVAM close to the hub of the needle.
- 2) Remove the donation needle, then gently lower the needle into the barrel of the SECUVAM device until a click is heard.
- 3) The needle and SECUVAM must be held upright. Ensure that no pressure is put on the donation line at any time, to prevent possible spillage.
- 4) The SECUVAM may be pushed into the barrel of the VACUVAM.
- 5) The needle, SECUVAM and any other devices should be disposed of in a safe manner.

✓ Pack handling whilst processing:

- 1) The processing depends on the pack configuration. Process the pack as required by the blood transfusion centre.
- 2) Ensure a proper packing of bags and tubing for centrifugation.
- 3) To allow fluid transfer, break the cannula at the seal, by bending the cannula with a side-to-side motion.
- 4) After transfer, seal the transfer tubing, using a permanent closure device, and cut to separate.

✓ Handling of soft housing filters

- 1) Do not fold or mishandle the filter(s).
- 2) If they are centrifuged, filters should be kept flat between transfer bags, with the arrow on the filter pointing to the side.
- NOTE: If either face of the filter remains white, the blood components have not been leucocytosed and must be discarded.

✓ Handling of port protectors

- 1) Using both hands tear the port protector as indicated.
- 2) Hold down the port protector piece with one hand.
- 3) Using the free hand, take the spike and engage into the port.
- 4) Twist the spike to secure fully in port.

✓ Precautions for disposable set

Disposal of bags should be carried out in accordance with regulations applicable to potentially infected, biological waste.

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
LEUCOLAB LCG4b

Instruction for use of the LEUCOLAB LCG4b
Leucocyte reduction filter for Red Blood Cells

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Instruction for use of the LEUCOLAB LCG4b Leucocyte reduction filter for Red Blood Cells

EN

INDICATIONS

The MacoPharma LeucoLab LCG4b High Efficiency Leukocyte Removal Filter is a laboratory filter intended for the removal of leukocytes, platelets and micro-aggregates from a single unit of red cells (prepared from approximately 450 to 500 ml of blood), removed from the fridge.
Sterile fluid path sterilized by beta radiation
STERILE and non-pyrogenic fluid pathway.

PRECAUTIONS

Follow instructions carefully.
Follow normal aseptic techniques.
For single use: do not re-sterilise or re-use.
Avoid contact with sharp objects.
Do not vent.
Do not use 10 days after opening the packing.
Do not use if packaging is damaged or end protectors are loose or displaced.
Macopharma recommends filtration of the red cell concentrate within 10 days after collection for optimal leucoreduction.
Risks linked to the reuse: septic risk, loss of efficacy.
All or part of this medical device is made of PVC plasticized with DEHP. According to some studies, DEHP could potentially be harmful to the reproductive system of male foetuses. The prescriber is solely responsible for choosing to use this device on women who are either pregnant or breast-feeding, or on young male infants. Nevertheless, DEHP-plasticized PVC is in compliance with the European Pharmacopoeia.
All manufacturing processes and all components in contact with the donors, the users and the blood components dedicated to the patients are latex free.

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- Close clamp A, B, and C.
- Insert the spike into the blood bag transfusion port by using a twisting motion. Sterile connection of the tubing to the blood bag is preferred to permit a normal storage period. WARNING: If sterile docking is not used, the product must be used within the time frames established for "open" systems.
- Gently and thoroughly mix blood bag and hang on drip stand upside down.
- Hang up the system so that the leucoreduced RBC bag is positioned 70cm below the LCG4b filter.

TO PRIME FILTER

- Priming of the filter with red cells in additive solution may be done by gravity.
- Open clamp A and C.

FILTRATION

Allow red cells to flow by gravity into transfer bag until you see no more blood dripping in the drip chamber. Allow sufficient time for maximum recovery of RBC.

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- Close clamp C.
 - Open clamp B.
 - Hold transfer bag in upright position and squeeze to expel air through the by-pass tubing, till the blood reaches the by-pass tubing.
 - To obtain maximum red cell recovery, open clamp C and allow air to drain the filter.
 - Seal off the segment tubing.
- For quality control the filtered red cells should be sampled from the side tubing (Pig Tail), not the filtration line, after mixing the filtered red cells thoroughly.

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Instruction for use of the LEUCOLAB

pediatric LCG2b / LCG4

Leucocyte removal filter for Red Cells

ENG

INDICATIONS

The MacoPharma LeucoLab LCG2b / LCG4 High Efficiency Leucocyte Removal Filter is a laboratory filter intended for the removal of leucocytes, and micro-aggregates from a single unit of red cells (prepared from approximately 450 to 500 mL of blood), at either 1 to 6 °C or at room temperature (20-24°C).

Sterilized by electron beam (Beta-irradiation).

STERILE and non-pyrogenic fluid pathway.

Risks linked to the reuse: septic risk, loss of efficacy.

All or part of this medical device is made of PVC plasticized with DEHP. According to some studies, DEHP could potentially be harmful to the reproductive system of male fetuses. The prescriber is solely responsible for choosing to use this device on women who are either pregnant or breast-feeding, or on young male infants. Nevertheless, DEHP-plasticized PVC is in compliance with the European Pharmacopoeia.

All manufacturing processes and all components in contact with the donors, the users and the blood components dedicated to the patients are latex free.

PRECAUTIONS

Follow instructions carefully.

Follow normal aseptic techniques.

For single use: do not re-sterilise or re-use.

Do not use if packaging is damaged or end protectors are loose or displaced.

For a better filtration and a better storage of filtered red cells, it is recommended that the blood used is mixed thoroughly and is less than 4 days old when using LCG2b and less than 15 days old when using LCG4.

PREPARATION OF FILTRATION SET AND BLOOD BAG :

- . Close clamp A, B, and C.
- . Insert the spike into the blood bag transfusion port by using a twisting motion.

Sterile connection of the tubing to the blood bag is preferred to permit a normal storage period. WARNING: If sterile docking is not used, the product must be used within the time frames established for "open" systems.

- . Gently and thoroughly mix blood bag and hang on drip stand upside down.
- . Hang up the system so that the leucoreduced RBC bag is positioned 60cm below the LCG2b/LCG4 filter.

TO PRIME FILTER :

- . Priming of the filter with red cells in additive solution may be done by gravity.
- . Open clamp A and C.

Note: For high hematocrit units you may squeeze the blood bag to prime the filter. Continue squeezing the blood bag only until the filter is primed.

FILTRATION :

Allow red cells to flow by gravity into transfer bag until you see no more blood dripping in the drip chamber (Pre-filter). Allow sufficient time for maximum recovery of RBC.

USE OF THE BY-PASS TUBING :
(To expel air from the transfer bag)

- . Close clamp C.
- . Open clamp B.
- . Hold transfer bag in upright position and squeeze to expel air through the by-pass tubing, till the blood reaches the by-pass tubing.
- . To obtain maximum red cell recovery, open clamp C and allow air to drain the filter.
- . Seal off the segment tubing.
- . For quality control of leucoreduced Red Cell Concentrate, samples must be collected from the lateral sampling tubing and not from the main tubing of the device.

Mix the leucoreduced concentrate thoroughly.

PEDIATRIC SEPARATION :

- . Hang the transfer bag upside down to suspend the pediatric bags freely beneath.
- . Break the break-away canula to allow the blood to flow to the pediatric bags.

Note : The original transfer bag may be used as a 5th Pedi bag.

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ENGLISH

The Maco Pharma LeucoLab LCG2b / LCG4 High Efficiency Leucocyte Removal Filter is a laboratory filter intended for the removal of leucocytes, platelets, and microaggregates from a single unit of red cells or whole blood, at either 4°C or at room temperature (20-24°C).

The above set has the capability of washing the leucodepleted red cells, if required.

Sterilised by Beta-irradiation.
Follow normal aseptic techniques.
STERILE and non-pyrogenic fluid pathway.

PRECAUTIONS

Follow instructions carefully.

For single use; do not re-sterilise or re-use.

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For a better filtration and a better storage of filtered red cells or whole blood, it is recommended that the blood filtered with LCG2b is not older than 4 days and 15 days with LCG4.

The PVC complies with the European Pharmacopoeia. The plasticizer is DEHP. The biological products stored and/or circulating in these devices could contain phthalates.

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INSTRUCTION FOR USE

PREPARATION OF FILTRATION SET AND BLOOD BAG

- Close all the clamps.
- Insert spike 1 into the blood bag transfusion port by using a twisting motion. Sterile connection of the tubing to the blood bag is also possible.
- Gently and thoroughly mix blood bag and hang on drip stand upside down.
- Hang the bag up so that the tubings are vertical between the blood bag and the transfer bag, the supernatant bag may be rested on a shelf below.

TO PRIME FILTER

- Priming of the filter with red cells in additive solution or whole blood may be done by gravity.
- Open clamp A and clamp B

N.B. For high hematocrit units squeeze blood bag and to prime the filter. Continue squeezing blood bag until filter is primed.

FILTRATION

- Allow red cells or whole blood to flow by gravity into transfer bag until blood bag has completely drained.

USE OF THE BY-PASS TUBING (To expel air from the transfer bag)

- Close roller clamp B.
- Open clamp C, make sure that clamp D is closed.
- Hold transfer bag in upright position and squeeze to expel air through the by-pass tubing.
- To obtain maximum red cell recovery, open clamp B and allow air to drain the filter.

ADDING NaCl TO THE LEUCODEPLETED RED CELLS:

- Close clamp C, make sure that clamp D is closed.
- Insert spike 2 into the NaCl bag port by using a twisting motion.
- Hang the NaCl bag on a drip stand upside down.
- Open clamp D, make sure that clamp C is closed.
- Transfer the quantity of NaCl you need to the transfer bag, whilst flushing the blood remaining in the by-pass tubing (thus saving more blood).
- Seal off the tubing above the transfer bag.

WASHING THE LEUCODEPLETED RED CELLS:

- Wash using your standard methods, the supernatant may be separated into the supernatant bag attached by breaking the break-away cannula, and expressing the supernatant to the supernatant bag, and sealing it off.
- Make sure that the supernatant is taken care of according to the regulations and laws applicable in your country, for biological waste.

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


MacoPharma Thromboflex TXP Filter - TDV & TRV Systems

Config. 1

Config 2

INDICATIONS

- The MacoPharma Thromboflex TDV and TRV systems are intended for the leucocyte removal from platelet concentrate (PC) prepared from Buffy-Coat (BC).
- The integrated Oxy-PL storage bag is suitable for the storage of Leucoreduced Platelet Concentrate (LR PC) at $22 \pm 2^\circ\text{C}$ up to 5 or 7 days under constant agitation.

WARNING

- Bring to a temperature between $+15^\circ\text{C}$ and $+30^\circ\text{C}$ before use.
- Avoid contact with sharp objects.
- Inspect overwrap and pack.
- The sterility of the device is ensured unless the overwrap is not opened or damaged
- DO NOT USE if the overwrap and pack present any visible signs of deterioration
- Do not VENT
- Do not add any label on the storage bag
- Do not write on bags
- Use only over labels validated
- Processing must be performed by skilled technicians
- Closure piercing device used must comply with the regulation ISO1135-4
- Risks linked to the reuse : Septic risk and loss of efficacy

• All or part of this medical device is made of PVC plasticized with DEHP, except the storage container made of PVC-BTHC. According to some studies, DEHP could potentially be harmful to the reproductive system of male fetuses. The prescriber is solely responsible for choosing to use this device on women who are either pregnant or breast-feeding, or on young male infants. Nevertheless, DEHP-plasticized PVC is in compliance with the European Pharmacopeia.

- All manufacturing processes and all components in contact with the donors, the users and the blood components dedicated to the patients are latex free.

PRECAUTIONS FOR DISPOSABLE SET

- Disposal of bags should be carried out in accordance with regulations applicable to potentially infected, biological waste.

INSTRUCTIONS

✓ Recommendations before filtration

Note: It is recommended to filter the PC within 24 hours from the whole blood collection.

- Select 4 up to 6 ABO compatible BCs

Note: BCs should rest 2 hours before pooling

- Peel open the transparent overwrap using the corner flap.

- Close all clamps on the system.

Note: All tubing connection should be performed using a sterile docking device.

✓ Pooling

- Connect the BCs to the TDV/TRV system using train / cascade or octopus method and connect the appropriate volume of Platelet Additive Solution (PAS) or Plasma (PLS).
- Check the integrity of the seals.
- Pool the BCs together and wash the empty BC bags with PAS/PLS.
- Seal and disconnect the Pooled-BC.

✓ Centrifugation

- Pack the Pooled-BC with the TDV/TRV system according to MacoPharma packing instructions.
- Centrifuge the Pooled-BC under conditions validated for a separation of Pooled-BC into PC and residue.
- When centrifugation is complete, carefully unload the units from the centrifuge.

Note: Rough handling may contaminate the PC with Red Blood Cells (RBCs).

✓ Separation and Filtration

- Gently load the unit on a component separator using the pins.
- Make sure that the TXP filter is positioned vertically: the inlet tubing downward and the outlet tubing upward.
- Separate the PC from the residue and simultaneously filter the PC.

Note: It is recommended to push the residue into the filter to maximize platelet recovery. Care must be taken not to contaminate the LR PC with RBCs. With the config 2, if there is a red cell trapping, in order not to contaminate the LR PC, close the clamp (a) leading to the storage container, open the clamp (b) and break the break-away cannula with a side-to-side motion. Then, launch the filtration and once the red cell trapping has been diverted in the sampling pouch, close the clamp (b), re-open the clamp (a) and keep filtering the PC.

MacoPharma does not recommend a filtration pressure higher than 250 mbar.

- Heat seal the tubing 1 cm below the filter.
- Disconnect the LR PC and dispose of filter and transfer bag.

✓ Air Removal & Sampling

- Remove the air from the LR PC by squeezing the storage bag and transfer the air into the sampling pouch.
 - If necessary transfer into the sampling bag the appropriate volume of LR PC for quality control.
- Note: Make sure that no foam remains into the LR PC.
- Heat seal the tubing between the «Y» junction and storage bag.
 - Disconnect the sampling pouch.
 - Allow the LR PC unit to rest 1 hour without agitation before being placed under constant agitation.

✓ Storage

- Store the LR PC at $22^\circ\text{C} \pm 2^\circ\text{C}$ under constant agitation.

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Anticoagulant solutions
for automated apheresis

ENGLISH
- Sterile and non pyrogenic anticoagulant solutions for blood anticoagulation.
- MacoPharma anticoagulant solutions are strictly reserved for the automated apheresis procedures.

✓ **COMPOSITION OF THE BAG**
- Plastified PVC bag.

✓ **COMPOSITION OF ACD formula A**
Trisodium Citrate, 2H₂O: 22.0g, Citric Acid, 1H₂O: 8.0g, Glucose, 1 H₂O: 24.5g, PPI Water: 1000mL

✓ **COMPOSITION OF CITRATE 4%**
Citric Acid, 1H₂O: 2.8g, Trisodium Citrate, 2H₂O: 39.6g, PPI Water: 1000mL

✓ **PRECAUTIONS FOR USE**
- Do not use for direct infusion.
- Do not use if solution is cloudy, or if there is any visible sign of deterioration.
- Do not re-use the bag. After opening the bag, discard the non used product.
- Do not vent.
- Do not use if the overwrap and/or bag are damaged.
- Check expiry date before use.
- Risks associated with re-use of the product are: Product contaminated by another blood product, Septic risk, risk of misuse.
- All or part of this medical device is made of PVC plasticized with DEHP. According to some studies, DEHP could potentially be harmful to the reproductive system of male foetuses. The prescriber is solely responsible for choosing to use this device on women who are either pregnant or breast-feeding, or on young male infants. Nevertheless, DEHP-plasticized PVC is in compliance with the European Pharmacopoeia.
- Discard according to the local regulations' on the potentially dangerous biological products.
- All manufacturing processes and all components in contact with the donors, the users and the blood components dedicated to the patients are latex free

✓ **STORAGE CONDITIONS**
- 2 years shelf-life in original sealed overwrap.
- Do not store above +35°C in original overwrap.
- Do not freeze.
- Protect from light during storage.

✓ **DOSAGE**
Anticoagulant solutions are added to whole blood during collection in a proportion defined by the person responsible for the procedure. The following rates are commonly advised by apheresis machines manufacturers : 1/9 for ACD-A solution, 1/15 for citrate 4%. This rate, however, can be changed according to the performance of the apheresis machines.

✓ **Instructions for use**
Installation of the apheresis kit on the separator: the connection of the MacoPharma anticoagulant solution bag to the apheresis system is carried out according to the separator's instructions, by 2 means :
① - break and detach the twist-off situated at the end of the anticoagulant solution, and connect the bag using a spike.
② - connect the male luer, located on the anticoagulant solution. If you use the luer connection:
A 1. Connect the Anticoagulant Solution with the luer connector to the apheresis set.
2. Break the breakable connector in the Anticoagulant Solution line by moving the connector back and forth 90° until an audible "snap" sound is heard in both directions. This separates the top piece of the frangible (arrow-shaped design) from the male luer base.
B 3. Bend tubing 180° to force the arrow piece of the frangible upward into the soft tubing closer to the bag opening and away from the base. This prevents this piece from reseating. To ensure proper flow, check to see if a gap exists between the arrow and base of the frangible. If not repeat step 3.
C 4. Hang the Anticoagulant Solution bag on the IV hook.
5. Follow the valid apheresis protocol of your institution.

✓ **CONTRA-INDICATIONS**
Contra-indications are those of blood collection. For cellular apheresis, check haemostasis. Use only for automated plasmapheresis or cytopheresis procedures. Do not perform apheresis in case of: serious haemopathy with haemostatic trouble, cerebral vascular accidents, pericardial or endocardial cardiopathy, post surgical period after intervention on central nervous system and spinal cord.
✓ **UNDESIRABLE EFFECTS**
Paresthesia : stinging sensations on lips and extremities ; Nausea ; Dysgueusia : weird taste sensation ; Cold sensation.
✓ **PREGNANCY AND LACTATION**
Do not carry out apheresis among pregnant women, nor during lactation (except plasmatic exchange if absolutely (indicated) necessary).
✓ **OVERDOSE**
Tetany possibility ; vomiting ; hypotension ; tachycardia. In case of overdose, control clinic and electrographic signs. If necessary inject calcium until (correction of the troubles) symptoms disappear.
✓ **SPECIAL WARNINGS**
Calcium injection is often necessary during cellular apheresis to limit the drop of the ionised calcium rate.
✓ **EFFECTS ON ABILITY TO DRIVE**
Rest time advised after an apheresis collection.

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