

CE 0123 Rx only

For patent information, see www.fphcare.com/ip

Manufacturer Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.co.nz Web: www.fphcare.com Importer/ Distributor **Australia** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 Fax: +61 3 9871 4998 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **Brazil** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraiso, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China** 代理人/售后服务机构:费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 传真: +86 20 32052132 **Denmark** Tel: +45 70 26 37 70 Fax: +46 83 66 310 **Finland** Tel: +358 9 251 66 123 Fax: +46 83 66 310 **France** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 Email: c.s@fphcare.fr **Germany** Fisher & Paykel Healthcare GmbH, Deutschland, Österreich, Schweiz, Wiesenstrasse 49, D 73614 Schorndorf, Germany Tel: +49 7181 98599 0 Fax: +49 7181 98599 66 **Hong Kong** Tel: +852 2116 0032 Fax: +852 2116 0085 **India** Tel: +91 80 2309 6400 **Ireland** Tel: 1800 409 011 Fax: +44 1628 626 146 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea** Tel: +82 2 6205 6900 Fax: +82 2 6309 6901 **Mexico** Tel: + 52 55 9130 1626 **Norway** Tel: +47 21 60 13 53 Fax: +47 22 99 60 10 **Poland** Tel: 00800 49 119 77 Fax: 00800 49 119 78 **Russia** Tel. and Fax: +7 495 782 21 50 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **Taiwan** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625 **Turkey** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 Fax: +44 1628 626 146 **USA/Canada** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

English (en)

ADULT CIRCUIT INSPIRATORY HEATED WITH MR290 AUTOFEED CHAMBER

For the delivery of breathing gases.

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INTERFACE CONNECTIONS	ISO 5356-1 Conical Connectors
COMPRESSIBLE VOLUME	1.00 L
COMPLIANCE @ 60 cmH₂O	1.64 ± 0.08 mL/cmH ₂ O (including 0.02 mL/cmH ₂ O measurement uncertainty)
RESISTANCE TO FLOW @ 30 L/min	0.74 ± 0.06 cmH ₂ O (including 0.03 cmH ₂ O measurement uncertainty)
MAXIMUM OPERATING PRESSURE	8 kPa
BREATHING TUBE LENGTH	1.6 m
MINIMUM TUBE INTERNAL DIAMETER	20 mm
GAS LEAKAGE @ 60 cmH₂O	<40 mL/min

OVERALL PERFORMANCE AT 20 °C TO 26 °C AMBIENT TEMPERATURE	
	NONINVASIVE MODE
Humidification Output	>12 mg/L
Flow Rate	5 – 120 L/min

SYMBOL DEFINITIONS

	Incorrect water level, replace MR290 chamber	Rx only	Prescription only
	Correct water level in the MR290 chamber	CE 0123	CE Marking 93/42/EEC
	Consult operating instructions		Date of manufacture
	Type BF applied part		Manufacturer
	Not made with phthalates		Use-by date
	Not made with natural rubber latex		Caution/Consult instructions for use
	Single use		7 Days maximum use
	Lot number		Transportation and storage temperature limits
	Reference number		European Union authorised representative

WARNINGS, CAUTIONS, AND NOTES

⚠ WARNINGS

- DO NOT reuse this product. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.
- The use of breathing circuits, chambers, accessories or combinations which are not recommended by Fisher & Paykel Healthcare may result in poor humidification system performance, ventilator malfunction, and harm to the patient/ user.
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious harm or death.
- DO NOT touch the heater-plate or chamber base. Surfaces may exceed 85 °C. Failure to comply may result in a skin burn.

- Failure to comply with the following warnings may impair performance of the device or compromise safety (including potentially cause serious harm):**
- When mounting a humidifier adjacent to a patient, ensure that the humidifier is always positioned lower than the patient.
 -  DO NOT use beyond 7 days maximum duration of use.
 - DO NOT cover the circuit with materials such as blankets, towels or bed linen.
 - DO NOT stretch or milk the tubing.
 - DO NOT soak, wash, or sterilize this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
 - DO NOT use the chamber if the water level rises above the maximum water level line.
 - DO NOT use the chamber if the seals are not intact when received, or if it has been dropped.
 - DO NOT operate the chamber at an angle in excess of 10°.
 - DO NOT spike the water source until the blue caps have been removed. Should the primary float fail, splashing into the circuit may occur if the chamber is being operated in excess of 80 L/min.
 - Visually inspect breathing sets for damage (e.g. a crushed tube or cracked connector) before use, and replace if damaged.
 - DO NOT crush tube.
 - The water source must be at least 50 cm higher than the chamber.
 - DO NOT use CPAP valve with air entrainer.
 - DO NOT use heated wire breathing circuits without gas flow. If gas flow is interrupted, turn the humidifier off.

- Ensure there is a water supply connected to the chamber and that water is present within the chamber.
- DO NOT fill the chamber with water in excess of 37 °C.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

Attention

- Ensure appropriate ventilator or flow source alarms are set before connecting breathing set to patient.
- Avoid prolonged contact of heated tube with patient's skin.
- Perform a pressure and leak test on the breathing system and check for occlusions before connecting to a patient.
- Use USP sterile water for inhalation or equivalent for humidification. DO NOT add other substances to the water.
- Regularly monitor and drain condensate build-up in the circuit.
- Check all connections are tight before use.
- Check that the heater wire is evenly distributed along the circuit and not bunched or kinked.

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
REF 185041098 REV T 2019-06 © 2019 Fisher & Paykel Healthcare Limited

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

REF 185041388 REV N 2019-09 © 2019 Fisher & Paykel Healthcare Limited

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
ADULT BREATHING CIRCUIT REUSABLE HEATED WALL 1.5m


900MR810 / 900MR810E (For use with Fisher and Paykel MR810 humidifier only)

Flow: ≥5 LPM to ≤60 LPM

900MR810E


900MR810

Cleaning/Disinfection:

Use only approved disinfection and cleaning methods as shown:


Step 1. Cleaning:


Step 2. High-level Disinfection:

Use one of the disinfection methods shown:


a) Chemical

b) Pasteurization

c) Autoclave (connectors only)

CE 0123 Rx only

Fisher & Paykel
HEALTHCARE

Manufacturer Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.co.nz Web: www.fphcare.com Importer/ Distributor
Australia (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 Fax: +61 3 9871 4998 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554
Brazil Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo – SP, Brazil Tel: +55 11 2548 7002 **China** 代理人/售后服务机构:费雪派克医疗保健（广州）有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋
301号 电话: +86 20 32053486 传真: +86 20 32052132 **Denmark** Tel: +45 70 26 37 70 Fax: +46 83 66 310 **Finland** Tel: +358 94 1590 355 Fax: +46 83 66 310 **France** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât. F5, BP 512, Villebon-sur-Yvette,
91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 Email: c.s@fphcare.fr **Germany** Fisher & Paykel Healthcare GmbH & Co. KG, Deutschland, Österreich, Schweiz, Wiesenstrasse 49, D 73614 Schorndorf, Germany Tel: +49 7181 98599 0
Fax: +49 7181 98599 66 **Hong Kong** Tel: +852 2116 0032 Fax: +852 2116 0085 **India** Tel: +91 80 2309 6400 **Ireland** Tel: 1800 409 011 Fax: +44 1628 626 146 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115
Korea Tel: +82 2 6205 6900 Fax: +82 2 6309 6901 **Norway** Tel: +47 21 60 13 53 Fax: +47 22 99 60 10 **Russia** Tel. and Fax: +7 495 782 21 50 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10
Switzerland Tel: 0800 83 47 63 Fax: 0800 83 47 54 **Taiwan** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625 **Turkey** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374,
Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 Fax: +44 1628 626 146 **USA/Canada** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

English  Refer to MR810 humidifier user instructions for further set up and correct use information.

Intended Use

The Fisher & Paykel Healthcare 900MR810/900MR810E adult breathing circuits are designed for use with the MR810 humidifier in order to provide therapeutic levels of heat and humidity to a patient's inspired respiratory gases, when using a continuous or intermittent ventilator system or a continuous gas flow system.

Addition of heat and humidity to the supply of cold and dry respiratory gases provided through ventilation is beneficial to prevent drying of the patient airways. The MR810 System is designed for use in hospitals, long term care facilities and homes under the prescription of a qualified medical professional.

Symbol:

Cleaning/Disinfection

Use only approved disinfection and cleaning methods as set out in these instructions:

Discard tubes and all circuit components after 20 cleaning cycles, or six months after first use, whichever occurs earlier.
Do not exceed maximum recommended cleaning temperatures as higher cleaning temperatures may damage the circuit.
Do not use automated washer/cleaner as it may damage the circuit.

Cleaning and disinfection instructions:

Step 1 – Cleaning:

Caution: When cleaning, wear appropriate personal protection equipment, such as: disposable gloves, splash visor/safety glasses and protective gown.

- Dis-assemble the circuit
NOTE: Clean circuit as soon as possible after use.
- Rinse circuit components directly after use with warm tap water (<40°C).
- Fill a sink/container with warm water (<40°C) and Sodium Hypochlorite to 0.5% maximum concentration. Then soak tubing and components for 30–40 minutes.
- Clean the inside of the tube removing all visible debris, if needed, with a soft bottle brush. Do not use hard brushes inside the tube.
- Rinse the inside and outside of the tube thoroughly with warm tap water to remove the remaining Sodium Hypochlorite and debris.
- Hang the tubes vertically upright to remove excess water. Tubes can be air dried up to 75°C.

Warning: Do not clean or dry circuit with solvents or alcohol.

- Proceed to high level disinfection, step 2.

Step 2 – High-level Disinfection:

Use one of the disinfection methods shown:

a. Disinfection with 0.3% orthophthalaldehyde (OPA) chemical:

Caution: Before handling OPA
Wear appropriate personal protection materials, such as: disposable gloves, splash visor/safety glasses and protective gown. Ensure disinfection is performed in a fume cupboard or with fume extraction. Read OPA manufacturers instructions for use.

- Fill sink/container with disinfectant OPA 0.3%, such as Johnson and Johnson Cidex OPA (follow manufacturer instructions for correct preparation of OPA).
Submerge tube and components fully in disinfectant.
Soak for a minimum of 30 minutes. Then remove all parts from the OPA.
- Change gloves.
- Rinse tubes and components with sterile water (follow manufacturer instructions for correct rinsing procedure of OPA).
- Hang the tubes vertically to let them air dry. Tubes can be air dried up to 75°C.
Inspect tube and circuit components, check for any cracks or tears or damage.
When all components are dry re-assemble the breathing circuit.
Store circuit in a dry and dust free environment ready for next use.

b. Disinfection by Pasteurization:

- Ensure pasteurizer is set to control between 75 and 78°C maximum.
NOTE: Temperatures higher than 78°C will damage the breathing tube and shorten its life.
- Remove circuit components, hang the tubes to dry vertically. Tubes can be air dried up to 75°C.
Inspect tube and circuit components, check for any cracks or tears or damage.
When all components are dry re-assemble the breathing circuit.
Store circuit in a dry and dust free environment ready for next use.

c. Autoclave:

Warning: Only the circuit connectors can be autoclaved. Autoclaving will destroy the tubes.

- Ensure the autoclave is set at either 121°C or 136°C.
Leave circuit connectors in autoclave for 30 minutes at 121°C or 4 minutes at 136°C.
Inspect circuit connectors, check for any cracks or damage.
When all connectors are dry, re-assemble the circuit connectors to disinfected tubing.
Store circuit in a dry and dust free environment ready for next use.

Technical Specifications:

Type BF

Interface connections: ISO 5356-1 or ISO 5367-1

Ambient temperature and inlet gas operating range: 18°C – 26°C

Maximum delivered temperature: <40°C

	900MR810	900MR810E
Compliance:	4.6 mL/kPa/m	3.4 mL/kPa/m
Compressible volume:	640 mL	1200 mL
Leakage rate:	<25 mL/min @ 6 kPa	<50 mL/min @ 6 kPa
Resistance to flow at rated flow:	0.4 cmH ₂ O @ 60 ± 1 L/min	1.6 cmH ₂ O @ 60 ± 1 L/min
Flow range:	≥5 to ≤60 L/min *	≥5 to ≤60 L/min *

* Tested with MR370 chamber

Warning:

- Clean circuit prior to use and after each patient use, using approved disinfection methods only. Use of unapproved cleaning methods may damage the circuit and reduce its useable life.
- Must be set up and used by a trained medical professional.
- Inspect circuit before re-use, do not use if the circuit shows signs of deterioration, such as: cracks, tears, or damage.
- Use only with Fisher & Paykel MR810, humidification chambers and accessories.
- Ensure that the humidifier is always positioned lower than the patient.
- Do not fill the chamber with water in excess of 37°C.
- Set appropriate ventilator alarms.
- Ensure connector with patient symbol is located at patient end of circuit.
- Use a watertrap in the expiratory limb when applicable.
- Check all tube connections are tight before use.
- Perform a pressure and leak test on the breathing system, and check for occlusions, before connecting to a patient.
- Ensure the complete circuit functions correctly with required ventilator settings, before connecting to a patient.
- Check circuit tubing and watertraps for condensate regularly, and drain as required.
- Do not use the heated breathing circuit without gas flow. If gas flow is interrupted turn the humidifier off.
- Do not cover the circuit with materials such as blankets, towels or bed linen.
- Surface temperature may reach 43°C. Prolonged contact of circuit tube with exposed patients' skin may result in pressure injury or skin damage.
- Do not crush the circuit tube.
- Humidity output of this system can be compromised if used outside recommended ambient temperature or flow range.
- No modification to this equipment is allowed.
- Disconnect tube by handling end connectors only, do not pull or twist tubing, as this may cause damage.
- Check the water level in the humidification chamber regularly and top up as required.
- Do not expose tubes to excessive UV radiation.
- Discard tubes and all circuit components after 20 cleaning cycles, or six months after first use, whichever occurs earlier. At end of product life, disinfect the circuit, then dispose of as uncontaminated medical waste.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

Warranty Statement:

Fisher & Paykel Healthcare warrants that the Evatherm™ circuit, when used in accordance with its instructions for use, shall be free from defects in workmanship and materials and will perform in accordance with Fisher & Paykel Healthcare's official published product specifications for a period of 90 days from the date of purchase by the end user. This warranty is subject to the limitations and exceptions set out in detail in <http://www.fphcare.com/evathermwarranty>

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

CE 0123 Rx only

F&P Evaqua is a trademark of Fisher and Paykel Healthcare.
For patent information, see www.fphcare.com/ip

Manufacturer Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.co.nz Web: www.fphcare.com Importer/ Distributor **Australia** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 Fax: +61 3 9871 4998 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **Brazil** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraiso, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China** 代理人/售后服务机构 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 传真: +86 20 32052132 **Denmark** Tel: +45 70 26 37 70 Fax: +46 83 66 310 **Finland** Tel: +358 9 251 66 123 Fax: +46 83 66 310 **France** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 Email: c.s@fphcare.fr **Germany** Fisher & Paykel Healthcare GmbH, Deutschland, Österreich, Schweiz, Wiesenstrasse 49, D 73614 Schorndorf, Germany Tel: +49 7181 98599 0 Fax: +49 7181 98599 66 **Hong Kong** Tel: +852 2116 0032 Fax: +852 2116 0085 **India** Tel: +91 80 2309 6400 **Ireland** Tel: 1800 409 011 Fax: +44 1628 626 146 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea** Tel: +82 2 6205 6900 Fax: +82 2 6309 6901 **Mexico** Tel: + 52 55 9130 1626 **Norway** Tel: +47 21 60 13 53 Fax: +47 22 99 60 10 **Poland** Tel: 00800 49 119 77 Fax: 00800 49 119 78 **Russia** Tel. and Fax: +7 495 782 21 50 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **Taiwan** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625 **Turkey** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 Fax: +44 1628 626 146 **USA/Canada** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

English (en)

ADULT VENTILATOR CIRCUIT DUAL HEATED WITH EVAQUA TECHNOLOGY

For delivery of breathing gases to adult patients requiring respiratory support. For use in an Intensive Care environment in conjunction with life support equipment.

TECHNICAL SPECIFICATIONS

Compatible with Fisher & Paykel Healthcare MR850 humidifiers. Refer to humidifier user instructions.

INTERFACE CONNECTIONS ISO 5356-1 Conical Connectors

COMPRESSIBLE VOLUME 1.60 L

COMPLIANCE @ 60 cmH₂O

Without Accessories 218 ± 0.11 mL/cmH₂O

With Accessories 2.6 ± 0.6 mL/cmH₂O
(including 0.02 mL/cmH₂O measurement uncertainty)

RESISTANCE TO FLOW @ 30 L/min

Without Accessories

Inspiratory limb 0.77 ± 0.05 cmH₂O

Expiratory limb 0.2 ± 0.04 cmH₂O

(including 0.03 cmH₂O measurement uncertainty)

With Accessories

Expiratory limb 0.99 ± 0.11 cmH₂O

(including 0.03 cmH₂O measurement uncertainty)

MAXIMUM OPERATING PRESSURE 8 kPa

BREATHING TUBE LENGTH 1.60 m

MINIMUM TUBE INTERNAL DIAMETER 20 mm

GAS LEAKAGE @ 60 cmH₂O <40 mL/min

TIDAL VOLUME >120 mL

OVERALL PERFORMANCE AT 20 °C TO 26 °C AMBIENT TEMPERATURE		
	INVASIVE MODE	NONINVASIVE MODE
Humidification Output	>33 mg/L	>12 mg/L
Flow Rate	5 – 60 L/min	5 – 120 L/min

FILTRATION EFFICIENCY

Viral	>99.99%
Organism	ΦX174 Bacteriophage
Bacterial	>99.999%
Organism	Bacillus subtilis
Mean particle size	3 µm
NaCl	98.04%

SYMBOL DEFINITIONS

 	Incorrect water level, replace MR290 chamber	Rx only	Prescription only
 	Correct water level in the MR290 chamber	CE 0123	CE Marking 93/42/EEC
	Consult operating instructions		Date of manufacture
	Type BF applied part		Manufacturer
	Not made with phthalates		Use-by date
	Not made with natural rubber latex		Caution/Consult instructions for use
	Single use		14 Days maximum use
	Lot number		Transportation and storage temperature limits
	Reference number		European Union authorised representative

WARNINGS, CAUTIONS, AND NOTES

⚠️ WARNINGS

- DO NOT reuse this product. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm, or death.
- The use of breathing circuits, chambers, accessories or combinations which are not recommended by Fisher & Paykel Healthcare may result in poor humidification system performance, ventilator malfunction, and harm to the patient/ user.
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious harm, or death.
- DO NOT touch the heater-plate or chamber base. Surfaces may exceed 85 °C. Failure to comply may result in a skin burn.

Failure to comply with the following warnings may impair performance of the device or compromise safety (including potentially cause serious harm):

- When mounting a humidifier adjacent to a patient, ensure that the humidifier is always positioned lower than the patient.
-  DO NOT use beyond 14 days maximum duration of use.
- DO NOT cover the circuit with materials such as blankets, towels, or bed linen.
- DO NOT stretch or milk the tubing.
- DO NOT soak, wash, or sterilize this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- DO NOT use the chamber if the water level rises above the maximum water level line.
- DO NOT use the chamber if the seals are not intact when received, or if it has been dropped.
- DO NOT operate the chamber at an angle in excess of 10°.
- DO NOT spike the water source until the blue caps have been removed. Should the primary float fail, splashing into the circuit may occur if the chamber is being operated in excess of 80 L/min.

- Visually inspect breathing sets for damage (e.g. a crushed tube or cracked connector) before use, and replace if damaged.
- The water source must be at least 50 cm higher than the chamber.
- DO NOT use a CPAP valve with an air entrainer.
- DO NOT use heated wire breathing circuits without gas flow. If gas flow is interrupted, turn the humidifier off. If the expiratory limb is removed from the circuit for any reason, disconnect it electrically from the humidifier.
- Ensure there is a water supply connected to the chamber and that water is present within the chamber.
- Change filter every 24 hours, or sooner if noticeable deterioration/build-up occurs, following standard hospital procedure.
- DO NOT fill the chamber with water in excess of 37 °C.
- DO NOT use medications containing Tyloxapol (such as Tacholiquin) as this may damage the tubing and lead to a loss of ventilation pressure.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

- Attention**
- Ensure appropriate ventilator or flow source alarms are set before connecting breathing set to patient.
- Fit only the supplied Fisher & Paykel Healthcare circuit hanger with care to avoid circuit damage.
- Avoid prolonged contact of heated tubes with patient's skin.
- Perform a pressure and leak test on the breathing system and check for occlusions before connecting to a patient.
- Use USP sterile water for inhalation or equivalent for humidification. DO NOT add other substances to the water.
- Check breathing circuits for condensation every 6 hours and drain if required.
- Check all connections are tight before use.
- Check that the heater wire is evenly distributed along the circuit and not bunched or kinked.
- When nebulized drugs are used, resistance to flow should be monitored, and the filter replaced following standard hospital procedure.

NOTES

- For use under the supervision of trained medical personnel.
- Dispose of product according to hospital protocol. User may be exposed to breathing tract fluids during disposal.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used to connect to the patient and other equipment before use.

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