
IFU Conox Sensor, March 2020 V1.0

Conox® Sensor

Instructions for Use

Manufacturer Quantum Medical S.L.U.

Av. Ernest Lluch 32, Tower 2,
Tecnocampus, 3rd Floor, Office
3.17, 08302 Mataró, Spain
+34 93 702 19 50
www.quantiummedical.com

Distributor Fresenius Kabi AG

61346 Bad Homburg, Germany
+49 (0) 61 72 / 686-0
www.fresenius-kabi.com

Table of Contents

1	General information	3
1.1	About this manual	3
1.2	Contact address	3
1.3	Disclaimer	3-4
1.4	Copyright	4
2	Introduction	5
2.1	Indications for use	5
2.2	Contraindications	5
2.3	Adverse effects	5
2.4	Warnings	6-7
2.5	Storage and Operation	6
2.6	Disposal	7
2.7	Biocompatibility	7
2.8	Safety symbols	7
2.9	Compatibility table	7
2.10	Notice	7
3	Instructions for use	8
3.1	Skin preparation and connections	8
3.2	Electrodes placement	8
3.3	How to use the Conox Sensor	9
4	Sensor Specifications	11

Quantum Medical S.L.U.
 CIF: B6505794
 Av. Ernest Lluch nº 32
 TecnoCampus, TCM2
 08302 Mataró, Barcelona, Spain

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1. About this manual
This operating manual contains all the information required to ensure both safety and correct functioning of the Conox® Sensor.

Before operating the Conox® monitor read this manual with attention to all warnings and cautions.

Note: Although the Conox is built and tested according to exact specifications, it is not intended to replace the supervision of depth of anesthesia monitoring by medical personnel. The user should become thoroughly familiar with the features, operations and accessories of the monitor and sensor, by exercising vigilance in its utilization.

The warnings, precautions and notes have specific meanings in this manual.

Warning: Advise against certain conditions that could cause fatal or serious injury to an operator, patient or any other person.

Caution: Advise against actions or situations that could damage the equipment described in this manual and/or any other equipment or goods, and/or cause environmental pollution.

Note: Gives information with respect to a function, procedure or unusual point to aid the operator.

Follow recommendations, instructions and warnings throughout this manual without exception.

Manufacturer
Quantum Medical S.L.U.
Av. Ernest Lluch 32, Tower 2,
Tecnocampus, 3rd Floor, Office 3.17.
08302 Mataró, Spain
+34 93 702 19 50
www.quantummedical.com

Distributor
Fresenius Kabi AG
61346 Bad Homburg
Germany
+49 6172 686 0
www.fresenius-kabi.com

1.3. Disclaimer

Manufacturer reserves all rights. No part of this document may be reproduced or published, in any format without written consent by the manufacturer.

The instructions in this manual are intended for trained professional medical personnel and/or authorized personnel to service and maintain the device.

All final decisions on patient treatment are the sole responsibility of the healthcare professional in charge and under no circumstances should the

device be used as the sole parameter for adjusting the anesthetic dose.

not be copied, distributed or amended without the written consent of the manufacturer.

1.4. Copyright

This document is the sole copyright of Quantum Medical S.L.U. and must

Introduction

Indications for use

Conox® Sensor is a system of adhesive electrodes applied to the skin of the patient's forehead to allow detection of the EEG signal. The sensor is applied to either side of the patient's forehead. The EEG signal is transferred from the sensor to the Conox® monitor via the Patient Cable. The sensor is intended for single patient use and is disposable after use.

User group: The device is intended to be used by medical personnel such as doctors and nurses and other healthcare professionals trained on the use of the device.

Target group: The device is intended to be used on adult patients undergoing general anesthesia and sedation.

Environment: The device must only be used in medical care facilities such as hospitals and other healthcare facilities.

Clinical Benefits: The Conox® sensor electrodes are designed with semi-liquid gel to allow easy application and good adhesion between electrodes and the patient's skin for the duration surgical procedure requiring general anesthesia and sedation.

Performance characteristics: The Conox® sensor is a low impedance device and will operate only when the skin impedance is <10kΩ.

2.2. Contraindications

DO NOT place the sensor on injured skin or scarred skin nor on patients with a known allergy to sensor.

DO NOT use in pediatric patients.

DO NOT use in patients with a history of psychiatric or neurological disease, drug abuse or medication known to affect the central nervous system.

DO NOT use device during operation of a CT, MRI, or X-RAY machine. Conducted current may cause burns.

DO NOT use in proximity less than 1.5 meters to radiofrequency and microwave emitting devices.

DO NOT use device in presence of ultrasonic devices such as ultrasonic surgery knives.

DO NOT operate device beyond temperature and humidity limits provided and outside of the defined use environment.

2.3. Adverse effects

The Conox® Sensor may cause skin irritation.

2.4. Warnings and

Precautions

Warning: If skin rash or another unusual symptom develops, stop using the sensor and remove it.

Warning: The Conox Sensor must not be used during defibrillator discharge while in conjunction with Conox monitor or Conox Lite QM8000-M as these devices do not possess protective capabilities during defibrillation.

For more information regarding use during defibrillation refer to Conox® and Conox® monitor IFU.

Warning: The usage time of the device is limited to 24 hours. Do not use the device beyond 24 hours as it will lead to degraded performance of the device.

Warning: Multiple use of this device may result in downgraded performance and harm to the patient.

Warning: Do not use the device past the expiry date written on the label. Electrodes may become dry resulting in degraded performance of the sensor.

Warning: Resting patients in prone position may result in increased pressure between skin and the device leading to skin irritation or burns.

Caution: The device is not intended to be sterilized at any point during its lifetime and shelf-life. Sterilization may compromise biocompatibility and material integrity of the device.

Caution: Impedance levels may change during use of the device. Users should use the signal quality index feature on the Conox® monitor to keep track of signal quality. Refer to Conox® monitor Instruction for Use for more information.

Caution: The device is not to be used with non-Conox equipment. The user is responsible for checking the compatibility between the monitor and the sensor; incompatibility between the sensor and the patient cable may result in degradation of device performance.

For further information refer to compatibility table in section 2.9 and to the Conox® monitor IFU.

Caution: Do not use if packaging is not sealed as this may result in degraded performance and/or harm to the patient from infection.

Caution: Do not use the device if it has been dropped outside of its packaging.

Caution: This device is not intended for operation by laypersons.

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The device should be stored according to the conditions indicated below.

Storage and transport conditions	Temperature: 20°C to 45°C. Pressure: 500 hPa to 1060 hPa. Humidity: 0%-80%, non-condensing Hyperbaric pressure: 86kPa-106kPa
Operating conditions:	Temperature: 5°C to 40°C Relative Humidity: ≤80% Air pressure: 86kPa-106kPa

The Device shall be stored with minimum exposure to sunlight.
not store next to combustible or reactive material.

2.6. Disposal

not dispose of in the dumpster.
Recycle the device and packaging according to international standards.

	Manufacturer
	Expiry date
	Single use
	24-hour continuous use
	Temperature limits
	Humidity limits
	Do not use if damaged packaging
	Nonsterile packaging
	Device is a medical device

2.7. Biocompatibility

The materials of the device in contact with the patient's skin have been verified for biocompatibility and cytotoxicity. The device does not contain any phthalates.

2.8. Safety symbols

The following signs and symbols are used on the Device and labelling.

	Consult user manual
	Recycle device according to international standards
	Latex free
	Lot number

2.10 Compatibility table

Monitor	Patient Cable	Sensor
Conox Lite QM8000-M	Z185075	Z029060
Conox Lite QM8001-M	Z185058	Z029061 Pediatric
Conox QM6000-M	Z185075	Z029060
	Z185065	Z029061 Pediatric
Conox QM7000-M	Z185058	Z029061 Pediatric

2.10 Notice

Any serious incident that occurs in relation to this device must be reported to the manufacturer and the competent authority of the state where the user and/or patient is established.

Instructions for use

The sensor can be placed in either the left or right hemisphere and due to its flexible design, it can be adapted to different head sizes.

Warning:

- If the skin presents any rash or other unusual symptoms, remove the electrodes from the patient
- It is important to take special care in patients with skin problems
- DO NOT place the electrodes on wounds.

Warning:

- The skin should always be prepared with the sandpaper provided with the sensor in order to ensure low electrode impedances. It is not recommended to use alcohol as this may leave a fine film on the skin which increases the electrode-skin impedance.

Note:


- Press the borders of the electrodes in order to stick them well. Pressing the electrode itself would spread the gel on the forehead compromising the sensor performance.
- If the adhesion or the contact of any of the electrodes is not enough after pressing around the Sensor, then use adhesive tape to press/stick the electrode to the skin.

3.6. Skin preparation and connections

For a correct measurement of the EEG signal, it is mandatory to prepare the skin in order to reduce the high impedance that exists in normal conditions. It is recommended to use very fine sandpaper designed specifically to clean the surface layer of the skin of the patient. A specific sandpaper is available for this purpose, contained inside of the sensor pouch. Before placing the electrodes, carefully clean the patient's skin using the sandpaper on the indicated area of the patient's forehead.


3.7. Electrodes placement

The sensor is composed of three electrodes: red, yellow and green. Each electrode shall be placed on the patient's forehead as indicated in the figure, with the read electrode as far as possible from the eyeline:


3.8. How to use the Conox Sensor


To use the Conox Sensor, the following steps should be carried out:


Remove the sensor from the pouch.


Clean the electrodes area with physiological serum and dry it. Do NOT use alcohol. Apply the sandpaper to the areas in which the three electrodes will be placed.

Place the Sensor on the forehead according to the figure. Position the Conox Sensor in either the left or the right side of the forehead.


Connect the Sensor to the Patient Cable as indicated in the figure (a,b,c).

Press the borders of the electrodes in order to stick them well. DO NOT press on the colored area of the electrodes to avoid the gel overflowing.

Attach the sensor so that it doesn't form an arch that could pull of the electrodes, as shown in the figure.

Caution: The Conox® Sensor should only be connected to the correct Conox® patient cable.

Caution: The Conox® monitor will have lead-off detection alarm if the sensor is not correctly connected to the patient cable or is disconnected during use.

4. Sensor specifications

Product details		Specifications
Product type		Sensor
Product reference		Z029061
Dimensions		30cm x 6cm
Single use		Yes
Electrode shape		EEG strip
Gel type		Semi-liquid
Material		Silver / Silver Chloride
Shelf Life		2 years (without opening the packaging)
Package Weight		16g
Packaging size		35cm x 9cm
Packaging material		PVC-free
Application		EEG signal acquisition
Compatible devices		Conox
Regulation – CE mark		The device complies with the requirements of the Medical Device Directive, 93-42-EEC

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