



INFANT CIRCUIT KIT

RT131, RT125, RT225 > 4 L/min

RT130, RT126, RT226 0.3 - 4 L/min

Single Use Only **CE 0123 Rx only**

Fisher & Paykel
HEALTHCARE

Manufacturer Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.co.nz Web: www.fphcare.com Importer/ Distributor **Australia** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 Fax: +61 3 9871 4998 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **Brazil** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo – SP, Brazil Tel: +55 11 2548 7002 **China** 代理人/售后服务机构: 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路 31 号 G12 栋 301 号 电话: +86 20 32053486 传真: +86 20 32052132 **Denmark** Tel: +45 70 26 37 70 Fax: +46 83 66 310 **Finland** Tel: +358 94 1590 355 Fax: +46 83 66 310 **France** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât. F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 Email: c.s@fphcare.fr **Germany** Fisher & Paykel Healthcare GmbH & Co. KG, Deutschland, Österreich, Schweiz, Wiesenstrasse 49, D 73614 Schorndorf, Germany Tel: +49 7181 98599 0 Fax: +49 7181 98599 66 **Hong Kong** Tel: +852 2116 0032 Fax: +852 2116 0085 **India** Tel: +91 80 2309 6400 **Ireland** Tel: 1800 409 011 Fax: +44 1628 626 146 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea** Tel: +82 2 6205 6900 Fax: +82 2 6309 6901 **Norway** Tel: +47 21 60 13 53 Fax: +47 22 99 60 10 **Russia** Tel. and Fax: +7 495 782 21 50 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **Taiwan** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625 **Turkey** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 Fax: +44 1628 626 146 **USA/Canada** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

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REF 185042411 Rev W 2019-08 © 2019 Fisher & Paykel Healthcare Limited

English

Technical Specifications

Interface connections

ISO 5356-1 Conical connectors

Single Heated System (including chamber):

Compliance1.03 mL/cm H₂O

Compressible volume 830 mL

Resistance to flow:2 cm H₂O @ 13 L/min

Maximum chamber operating pressure:8 kPa

Maximum circuit gas leakage rate75 mL/min @ 60 cm H₂O

Circuit length:1.5 m

Latex free

RT131, RT125, RT225:	Compatible with MR810, MR850 and 700 series humidifiers. Dräger Babylog: Suitable for flows < 25 L/min.
RT130, RT126, RT226:	Compatible with MR850 humidifier only. <div><div>MR810</div>DO NOT use with MR810 humidifier.</div>

Attention

- Check all connections are tight before use.
- Perform a pressure and leak test on the breathing system and check for occlusions before connecting to a patient.
- Check that the heater wire is evenly distributed along the circuit and not bunched or kinked.
- Avoid prolonged contact with patient's skin.
- Set appropriate ventilator alarms.

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Warning

7

- This product is intended to be used for a maximum of 7 days.
- DO NOT cover the circuit with materials such as blankets, towels or bed linen.
- DO NOT stretch or milk the tubing.
- DO NOT soak, wash, sterilise or re-use this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- DO NOT use heated wire breathing circuits without gas flow.
- If gas flow is interrupted turn the humidifier off.
- When mounting a humidifier adjacent to a patient ensure that the humidifier is always positioned lower than the patient.

- DO NOT use the chamber if the water level rises above the maximum water level line.
- DO NOT use the chamber if the seals are not intact when received, or if it has been dropped.
- DO NOT operate the chamber at an angle in excess of 10°.
- DO NOT fill the chamber with water in excess of 37 °C.
- DO NOT spike the water source until the blue caps have been removed. Should the primary float fail, splashing into the circuit may occur if the chamber is being operated in excess of 80 L/min.
- DO NOT touch the heater plate or chamber base. Surfaces may exceed 85 °C.
- Ensure there is a water supply connected to the chamber and that water is present within the chamber.
- Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

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REF 185042495 REV R 2019-08 © 2019 Fisher & Paykel Healthcare Limited

English (en)

Technical Specifications

Interface connections:
ISO 5356-1 Conical connectors
Single Heated System (including chamber):
Compliance 1.20 mL/cm H₂O
Compressible volume 831 mL
Maximum chamber operating pressure: 8 kPa
Maximum circuit gas leakage rate 75 mL/min @ 60 cmH₂O
Circuit length 1.5 m
Latex free
Compatible with MR850 and 700 series humidifiers

DO NOT use with MR810 humidifier.

Attention

- Check all connections are tight before use.
- RT227 and RT127 with red restrictor for use with SLE2000 and SLE2000 HFO ventilators only.
- RT228 with yellow restrictor for use with SLE4000 and SLE5000 ventilators only.
- Ensure that the red or yellow restrictor is positioned in the inspiratory limb at the swivel before use.
- Perform a pressure and leak test on the breathing system and check for occlusions before connecting to a patient.
- Check that the heater wire is evenly distributed along the circuit and not bunched or kinked.
- Avoid prolonged contact with patient's skin.
- Set appropriate ventilator alarms.
- Not recommended for use above 100mbar ΔP during high frequency ventilation.

MR290 Water Level Symbols

Incorrect water level, replace MR290 chamber

Correct water level in the MR290 chamber

Use USP sterile water for inhalation or equivalent.

Warning

This product is intended to be used for a maximum of 7 days.

- DO NOT cover the circuit with materials such as blankets, towels or bed linen.
- DO NOT stretch or milk the tubing.
- DO NOT soak, wash, sterilise or re-use this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- DO NOT use heated wire breathing circuits without gas flow. If gas flow is interrupted turn the humidifier off.
- When mounting a humidifier adjacent to a patient ensure that the humidifier is always positioned lower than the patient.

- DO NOT use the chamber if the water level rises above the maximum water level line.
- DO NOT use the chamber if the seals are not intact when received, or if it has been dropped.
- DO NOT operate the chamber at an angle in excess of 10°.
- DO NOT fill the chamber with water in excess of 37 °C.
- DO NOT spike the water source until the blue caps have been removed. Should the primary float fail, splashing into the circuit may occur if the chamber is being operated in excess of 80 L/min.
- DO NOT touch the heater plate or chamber base. Surfaces may exceed 85 °C.
- Ensure there is a water supply connected to the chamber and that water is present within the chamber.
- Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.

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



REF 185042559 REV N 2020-01 © 2020 Fisher & Paykel Healthcare Limited

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<div>English</div> <div>en</div> <div><div>Technical Specifications</div><div>Interface connections ISO 5356-1 Conical connectors</div><div>Single Heated System (including chamber): Compliance0.74 mL/ cmH₂O Compressible volume550 mL</div><div>Resistance to flow:2 cmH₂O @ 13 L/min</div><div>Circuit length1.5 m</div><div>Latex free</div><div>Compatible with MR850 and 700 series humidifiers.</div></div> <div>Attention:</div> <div><div><div><div><div></div></div><div></div></div><div><div><div></div></div><div></div></div></div><div><div><div><div></div></div><div></div></div><div><div><div></div></div><div></div></div></div></div> <div><div>Incorrect water level, replace MR290 chamber</div><div>Correct water level in the MR290 chamber</div></div>

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Warning:

DO NOT use the chamber if the water level rises above the maximum water level line.

DO NOT use the chamber if the seals are not intact when received, or if it has been dropped.

DO NOT operate the chamber at an angle in excess of 10°.

DO NOT fill the chamber with water in excess of 37 °C.

DO NOT spike the water source until the blue caps have been removed. Should the primary float fail, splashing into the circuit may occur if the chamber is being operated in excess of 80 L/min.

DO NOT touch the heater plate or chamber base. Surfaces may exceed 85 °C.

Ensure there is a water supply connected to the chamber and that water is present within the chamber.

Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.

California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

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


USER INSTRUCTIONS

- REF
- RT267 – Infant ventilator circuit dual heated with MR290 autofeed chamber for SLE2000 ventilators
- RT268 – Infant ventilator circuit dual heated with MR290 autofeed chamber for SLE4000 & SLE5000 ventilators

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English 

RT267 INFANT VENTILATOR CIRCUIT DUAL HEATED WITH MR290 AUTOFEED CHAMBER FOR SLE2000 VENTILATORS

RT268 INFANT VENTILATOR CIRCUIT DUAL HEATED WITH MR290 AUTOFEED CHAMBER FOR SLE4000 & SLE5000 VENTILATORS

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INTERFACE CONNECTIONS	ISO 5356-1 Conical Connectors
COMPLIANCE (WITH AND WITHOUT ACCESSORIES) @ 60 cmH₂O	1.30 ± 0.15 mL/cmH ₂ O (including 0.02 mL/cmH ₂ O measurement uncertainty)
COMPRESSIBLE VOLUME:	762 mL
MAXIMUM OPERATING PRESSURE	8 kPa
GAS LEAKAGE @ 60 cmH₂O	<30 mL/min
CIRCUIT LENGTH:	1.6 m

OVERALL PERFORMANCE AT 20 °C TO 26 °C AMBIENT TEMPERATURE	
	INVASIVE MODE
Humidification Output:	>33 mg/L
Flow Rate:	4 ~ 15 L/min

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WARNINGS, CAUTIONS, AND NOTES

 WARNINGS

- DO NOT reuse this product. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.
- The use of breathing circuits, chambers, accessories, or combinations which are not approved by Fisher & Paykel Healthcare may result in poor humidification system performance, ventilator malfunction, and harm to the patient/ user.
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious harm or death.
- DO NOT touch the heater-plate or chamber base. Surfaces may exceed 85 °C. Failure to comply may result in a skin burn.

Failure to comply with any of the following warnings may impair performance of the device or compromise safety (including potentially causing serious harm):

- When mounting a humidifier adjacent to a patient, ensure that the humidifier is always positioned lower than the patient.
-  DO NOT use beyond 7 days maximum duration of use.
- DO NOT cover the circuit with materials such as blankets, towels or bed linen.
- DO NOT stretch or milk the tubing.
- DO NOT soak, wash, or sterilize this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- DO NOT use the chamber if the water level rises above the maximum water level line.
- DO NOT use the chamber if the seals are not intact when received, or if it has been dropped.
- DO NOT operate the chamber at an angle in excess of 10°.
- DO NOT spike the water source until the blue caps have been removed. Should the primary float fail, splashing into the circuit may occur if the chamber is being operated in excess of 80 L/min.
- Visually inspect breathing sets for damage (e.g. a crushed tube or cracked connector) before use, and replace if damaged.
- The water source must be at least 50 cm higher than the chamber.

- DO NOT use heated wire breathing circuits without gas flow. If gas flow is interrupted, turn the humidifier off.
- Ensure there is a water supply connected to the chamber and that water is present within the chamber.
- DO NOT fill the chamber with water in excess of 37 °C.
- DO NOT use medications containing Tyloxapol (such as Tacholiquin) as this may damage the tubing and lead to a loss of ventilation pressure.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

Attention

- Ensure appropriate ventilator or flow source alarms are set before connecting breathing set to patient.
- Avoid prolonged contact of heated tubes with patient's skin.
- Perform a pressure and leak test on the breathing system and check for occlusions before connecting to a patient.
- Use USP sterile water or equivalent for humidification. DO NOT add other substances to the water.
- Check breathing circuits for condensation every 6 hours and drain if required.
- Check all connections are tight before use.
- Check that the heater wire is evenly distributed along the circuit and not bunched or kinked.
- DO NOT use with MR810 humidifier.
- RT267 with red restrictor for use with SLE2000 and SLE2000 HFO ventilators only.
- RT268 with yellow restrictor for use with SLE4000 and SLE5000 ventilators only.
- Ensure that the red or yellow restrictor is positioned in the inspiratory limb at the swivel before use.
- Not recommended for use above 100 mbar ΔP during high frequency ventilation.

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

USER INSTRUCTIONS

- REF
- RT265

– Infant Ventilator Circuit (>4 L/min) Dual Heated with MR290 Autofeed Chamber
- RT266

– Infant Ventilator Circuit (0.3 to 4 L/min) Dual Heated with MR290 Autofeed Chamber

CE0123

Rx only

F&P Evaqua is a trademark of Fisher and Paykel Healthcare.
For patent information, see www.fphcare.com/ip

Manufacturer Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.co.nz Web: www.fphcare.com Importer/ Distributor **Australia** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 Fax: +61 3 9871 4998 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **Brazil** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo – SP, Brazil Tel: +55 11 2548 7002 **China** 代理人/售后服务机构:费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 传真: +86 20 32052132 **Denmark** Tel: +45 70 26 37 70 Fax: +46 83 66 310 **Finland** Tel: +358 9 251 66 123 Fax: +46 83 66 310 **France** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 Email: c.s@fphcare.fr **Germany** Fisher & Paykel Healthcare GmbH, Deutschland, Österreich, Schweiz, Wiesenstrasse 49, D 73614 Schorndorf, Germany Tel: +49 7181 98599 0 Fax: +49 7181 98599 66 **Hong Kong** Tel: +852 2116 0032 Fax: +852 2116 0085 **India** Tel: +91 80 2309 6400 **Ireland** Tel: 1800 409 011 Fax: +44 1628 626 146 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea** Tel: +82 2 6205 6900 Fax: +82 2 6309 6901 **Mexico** Tel: + 52 55 9130 1626 **Norway** Tel: +47 21 60 13 53 Fax: +47 22 99 60 10 **Poland** Tel: 00800 49 119 77 Fax: 00800 49 119 78 **Russia** Tel. and Fax: +7 495 782 21 50 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **Taiwan** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625 **Turkey** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 Fax: +44 1628 626 146 **USA/Canada** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

RT265 **INFANT VENTILATOR CIRCUIT (>4 L/MIN) DUAL HEATED WITH MR290 AUTOFEED CHAMBER**

RT266 **INFANT VENTILATOR CIRCUIT (0.3 TO 4 L MIN) DUAL HEATED WITH MR290 AUTOFEED CHAMBER**

For delivery of breathing gases.

TECHNICAL SPECIFICATIONS

Compatible with Fisher & Paykel Healthcare MR850 humidifiers. Refer to humidifier user instructions.

INTERFACE CONNECTIONS	ISO 5356-1 Conical Connectors
COMPRESSIBLE VOLUME	758 mL
COMPLIANCE (WITH AND WITHOUT ACCESSORIES) @ 60 cmH₂O	
RT265	1.30 ± 0.15 mL/cmH ₂ O
RT266	1.12 ± 0.11 mL/cmH ₂ O
(including 0.02 mL/cmH ₂ O measurement uncertainty)	

RESISTANCE TO FLOW (WITH AND WITHOUT ACCESSORIES) @ 2.5 L/min

RT265 Inspiratory limb	0.09 ± 0.04 cmH ₂ O
RT266 Inspiratory limb	0.10 ± 0.05 cmH ₂ O
RT265 Expiratory limb	0.07 ± 0.04 cmH ₂ O
RT266 Expiratory limb	0.07 ± 0.05 cmH ₂ O
(including 0.03 cmH ₂ O measurement uncertainty)	

MAXIMUM OPERATING PRESSURE	8 kPa
BREATHING TUBE LENGTH	1.6 m
MINIMUM TUBE INTERNAL DIAMETER	10 mm
GAS LEAKAGE @ 60 cmH₂O	<30 mL/min
TIDAL VOLUME	<185 mL

OVERALL PERFORMANCE AT 20 °C TO 26 °C AMBIENT TEMPERATURE		
	INVASIVE MODE	
	RT265	RT266
Humidification Output	>33 mg/L	
Flow Rate	4 – 15 L/min	0.8 – 4 L/min

TYPICAL PERFORMANCE AT 23 °C AMBIENT TEMPERATURE		
	INVASIVE MODE	
	RT266	
Humidification Output	>33 mg/L	
Flow Rate	0.3 – 4 L/min	

Dräger Babylog: Suitable for flows <25 L/min

SYMBOL DEFINITIONS

	Incorrect water level, replace MR290 chamber	Rx only	Prescription only
	Correct water level in the MR290 chamber	CE 0123	CE Marking 93/42/EEC
	Consult operating instructions		Date of manufacture
	Type BF applied part		Manufacturer
	Not made with phthalates		Use-by date
	Not made with natural rubber latex		Caution/Consult instructions for use
	Single use		7 Days maximum use
	Lot number		Transportation and storage temperature limits
	Reference number		European Union authorised representative

WARNINGS, CAUTIONS, AND NOTES

⚠ WARNINGS

- DO NOT reuse this product. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm, or death.
- The use of breathing circuits, chambers, accessories or combinations which are not approved by Fisher & Paykel Healthcare may result in poor humidification system performance, ventilator malfunction and harm to the patient/ user.
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious harm, or death.
- DO NOT touch the heater plate or chamber base. Surfaces may exceed 85 °C. Failure to comply may result in a skin burn.
- Covering breathing tubes with a blanket or heating them in an incubator or with an overhead heater can affect the quality of the therapy or injure the patient.

Failure to comply with any of the following warnings may impair performance of the device or compromise safety (including potentially cause serious harm):

- When mounting a humidifier adjacent to a patient, ensure that the humidifier is always positioned lower than the patient.
-  DO NOT use beyond 7 days maximum duration of use.
- DO NOT cover the circuit with materials such as blankets, towels or bed linen.
- DO NOT stretch or milk the tubing.
- DO NOT soak, wash or sterilize this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- DO NOT use the chamber if the water level rises above the maximum water level line.
- DO NOT use the chamber if the seals are not intact when received, or if it has been dropped.
- DO NOT operate the chamber at an angle in excess of 10°.
- DO NOT spike the water source until the blue caps have been removed. Should the primary float fail, splashing into the circuit may occur if the chamber is being operated in excess of 80 L/min.
- Visually inspect breathing sets for damage (e.g. a crushed tube or cracked connector) before use and replace if damaged.
- The water source must be at least 50 cm higher than the chamber.
- DO NOT use heated wire breathing circuits without gas flow. If gas flow is interrupted, turn the humidifier off.
- RT266: Heated tubes should be kept out of direct contact with the patient's skin. At flow rates below 0.8 L/min, the surface temperature can reach 50 °C.
- Ensure there is a water supply connected to the chamber and that water is present within the chamber.
- DO NOT fill the chamber with water in excess of 37 °C.
- DO NOT use medications containing Tyloxapol (such as Tacholiquin) as this may damage the tubing and lead to a loss of ventilation pressure.
- California residents please be advised of the following: pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

Attention

- Ensure appropriate ventilator or flow source alarms are set before connecting breathing set to patient.
- Avoid prolonged contact of heated tubes with patient's skin.
- Perform a pressure and leak test on the breathing system and check for occlusions before connecting to a patient.
- Use USP sterile water for inhalation or equivalent for humidification. DO NOT add other substances to the water.
- Check breathing circuits for condensation every 6 hours and drain if required.
- Check all connections are tight before use.
- Check that the heater wire is evenly distributed along the circuit and not bunched or kinked.
- DO NOT use with MR810 humidifier.

NOTES

- For use under the supervision of trained medical personnel.
- Dispose of product according to hospital protocol. User may be exposed to breathing tract fluids during disposal.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used to connect to the patient and other equipment before use.

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