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USER'S GUIDE

0604-910978B

First Edition: 26 Oct 2012 Third Edition: 07 Mar 2016

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NOTE

- This hematology control is potential biohazardous material. The MEK-3D hematology control is intended solely for in vitro diagnostic use by trained personnel. The hematology control contains human-sourced and/or potentially infectious components. Components from human donors used in preparation of the hematology control were tested by FDA approved methods for the presence of the antibodies to human immunodeficiency virus (HIV-1/HIV-2) and Hepatitis C virus (HCV) as well as for Hepatitis B virus surface antigen and found to be negative.
- Do not use the hematology control after the expiration date. The measurement data may be inaccurate.
- Store the hematology control in a medical refrigerator.
- Store the hematology control at temperatures between 2 and 8°C (36 and 46°F). Do not freeze the hematology control.
- Allow the hematology control to warm to temperatures between 18 and 30°C (64 and 86°F) and then invert and mix the hematology control. Otherwise the contents may leak from the cap.
- If the cap of the hematology control becomes dirty, wipe with a non-woven towel such as Kimwipes.
- Platelet cannot be measured by precipitation or centrifuge.
- This product does not support measuring WBC differential under the microscope.

General

The MEK-3D hematology control is a tri-level control to measure the performance of Nihon Kohden hematology analyzers. Refer to the assay tables for specific instrument models.

NOTE

The MEK-3D hematology control is designed for Nihon Kohden hematology analyzers which classify three types of leukocyte. Do not use with any other instrument.

Manufacturer

NIHON KOHDEN CORPORATION
1-31-4 Nishi-chai, Shinjuku-ku, Tokyo 161-8560, Japan
Phone +81 (3) 5996-8036 Fax +81 (3) 5996-8100

EC REP

European Representative

NIHON KOHDEN EUROPE GmbH
Raiffeisenstrasse 10, D-61191 Rosbach, Germany
Phone +49 6003 827-0 Fax +49 6003 827-599

NIHON KOHDEN SINGAPORE PTE LTD
1 Maritime Square, #10-34 HarbourFront Centre
Singapore 099253
Phone +65 6376-2210 Fax +65 6376-2264

NIHON KOHDEN AMERICA, INC.
15353 Barranca Parkway, Irvine, CA 92618, U.S.A.
Toll-free +1-800-325-0283
Phone +1 949-580-1555 Fax +1 949-580-1550

The MEK-3D hematology control contains human erythrocytes, simulated leukocytes, and mammalian platelets in a plasma-like fluid. The hematology control is for in vitro diagnostic use.

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Instructions for Use

NOTE

- Dark red supernatant fluid, discoloration of the product or unacceptable results may indicate deterioration. Do not use the hematology control if deterioration is suspected.
 - Confirm that the lot number of the product matches the lot number on the assay sheet.
1. Take the hematology control out of the medical refrigerator and allow it to warm to temperatures between 18 and 30°C (64 and 86°F) for 15 minutes before mixing.
 2. Invert and mix the hematology control.
 - Mix the hematology control until the red cells are completely suspended.
 - Immediately before measuring the hematology control, gently turn it upside down 8 to 10 times.

NOTE

- Do not mix the hematology control roughly.
 - Do not mix the hematology control using a mixer.
 - Incomplete mixing of the hematology control affects the measurement result.
3. Measure the hematology control as instructed in the hematology analyzer operator's manual.
 4. Return the hematology control to the medical refrigerator within 30 minutes of use.

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Hematology Control MEK-5DN/MEK-5DH/MEK-5DL

0604-018444E

First Edition: 06 Jan 2003 Sixth Edition: 07 Mar 2016

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General Handling Precautions

WARNING

- Dispose of the hematology control by following your local laws for disposing of infectious medical waste. Otherwise, it may affect the environment.
- Do not swallow the hematology control. It may be infectious.
- When handling the hematology control, wear rubber gloves to protect yourself from infection.

NOTE

- This hematology control is potential biohazardous material. The MEK-5D hematology control is intended solely for in vitro diagnostic use by trained personnel. The hematology control contains human-sourced and/or potentially infectious components. Components from human donors used in preparation of the hematology control were tested by FDA approved methods for the presence of the antibodies to human immunodeficiency virus (HIV-1/HIV-2) and Hepatitis C virus (HCV) as well as for Hepatitis B virus surface antigen and found to be negative.
- Do not use the hematology control after the expiration date. The measurement data may be inaccurate.
- Store the hematology control in a medical refrigerator.
- Store the hematology control at temperatures between 2 and 8°C (36 and 46°F). Do not freeze the hematology control.
- Allow the hematology control to warm to temperatures between 18 and 30°C (64 and 86°F) and then invert and mix the hematology control. Otherwise the contents may leak from the cap.
- If the cap of the hematology control becomes dirty, wipe with a non-woven towel such as Kimwipes.
- Platelet cannot be measured by precipitation or centrifuge.
- This product does not support measuring WBC differential under the microscope.

General

The MEK-5D hematology control is a tri-level control to measure the performance of Nihon Kohden hematology analyzers. Refer to the assay tables for specific instrument models.

NOTE

The MEK-5D hematology control is designed for Nihon Kohden hematology analyzers which use light scatter technique and classify five types of leukocyte. Do not use with any other instrument.

NIHON KOHDEN

Manufacturer

NIHON KOHDEN CORPORATION

1-31-4 Nishiokichai, Shinjuku-ku, Tokyo 161-8560, Japan
Phone +81 (3) 5996 8036 Fax +81 (3) 5996-8100

European Representative

NIHON KOHDEN EUROPE GmbH

Raiffeisenstrasse 10, D-61191 Rosbach, Germany
Phone +49 6003 827-0 Fax +49 6003 827-599

The MEK-5D hematology control contains human erythrocytes, simulated leukocytes, and mammalian platelets in a plasma-like fluid. The hematology control is for in vitro diagnostic use.

To achieve optimum performance and safety, thoroughly read this guide and the hematology analyzer operator's manual before using the hematology control.

Instructions for Use


NOTE

- Dark red supernatant fluid, discoloration of the product or unacceptable results may indicate deterioration. Do not use the hematology control if deterioration is suspected.
- Confirm that the lot number of the product matches the lot number on the assay sheet.

1. Take the hematology control out of the medical refrigerator and allow it to warm to temperatures between 18 and 30°C (64 and 86°F) for 15 minutes before mixing.

2. Invert and mix the hematology control.

- Mix the hematology control until the red cells are completely suspended.
- Immediately before measuring the hematology control, gently turn it upside down 8 to 10 times.

NOTE

- Do not mix the hematology control roughly.
- Do not mix the hematology control using a mixer.
- Incomplete mixing of the hematology control affects the measurement result.

3. Measure the hematology control as instructed in the hematology analyzer operator's manual.

4. Return the hematology control to the medical refrigerator within 30 minutes of use.

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Opened hematology control should be used within 14 days.

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Some specifications of this product may differ depending on the destination country or region. Therefore descriptions in the Japanese manual and English and other language manuals may also differ.

The CE mark is a protected conformity mark of the European Community. Products marked with this symbol comply with the requirements of the In vitro Diagnostic Medical Device Directive 98/79/EC.

General Handling Precautions

WARNING

Dispose of the hematology calibrator by following your local laws for disposing of infectious medical waste. Otherwise, it may affect the environment.

WARNING

Do not swallow the hematology calibrator. It may be infectious.

WARNING

When handling the hematology calibrator, wear rubber gloves to protect yourself from infection.

NOTE

- This hematology calibrator is potential biohazardous material. MEK-CAL hematology calibrator is intended solely for in vitro diagnostic use by trained personnel. This product contains human-sourced and/or potentially infectious components. Components from human donors used in preparation of this product were tested by FDA approved methods for the presence of the antibodies to human immunodeficiency virus (HIV-1/HIV-2) and Hepatitis C virus (HCV) as well as for Hepatitis B virus surface antigen and found to be negative.
- Do not use the hematology calibrator after the expiration date. The measurement data may be inaccurate.
- Store the hematology calibrator in a medical refrigerator.
- Store the hematology calibrator at temperatures between 2 and 8°C (36 and 46°F).
- Use the hematology calibrator at temperature from 18 to 30°C (64 to 86°F). Otherwise, the measurement data may be inaccurate.
- Mix the hematology calibrator thoroughly before use.
- Platelet cannot be measured by precipitation or centrifuge.

General

The MEK-CAL hematology calibrator is for calibration of Nihon Kohden hematology analyzers. For calibration, refer to the provided assay table for specific instrument models and the Operator's Manual of the hematology analyzer.

NOTE

MEK-CAL hematology calibrator is designed for Nihon Kohden hematology analyzers. Do not use with any other instrument.

Manufacturer
NIHON KOHDEN CORPORATION
1-31-4 Nishiochiai, Shinjuku-ku, Tokyo 161-8560, Japan
Phone +81 (3) 5996-8036 Fax +81 (3) 5996-8100

EC REP European Representative
NIHON KOHDEN EUROPE GmbH
Raiffeisenstrasse 10, D-61191 Rosbach, Germany
Phone +49 6003 827-0 Fax +49 6003 827-599

Reagent Production
NIHON KOHDEN FIRENZE S.r.l.
Via Torta 72/74, 50019 Sesto Fiorentino (FI), Italy
Phone +39 055 3045 1 Fax +39 055 308548

SHANGHAI KOHDEN MEDICAL ELECTRONIC INSTRUMENT CORPORATION
7th Floor, Dawning Centre Tower A No.500, Hongbaoshi Road
Changning District, Shanghai 201103, China
Phone +86 (21) 6270-0909 Fax +86 (21) 6270-9700

MEK-CAL hematology calibrator contains human erythrocytes, simulated leukocytes, and mammalian platelets in a plasma-like fluid. MEK-CAL hematology calibrator is for in vitro diagnostic use.


To achieve optimum performance and safety, thoroughly read this instruction and the hematology analyzer operator's manual before using this hematology calibrator.

Instruction for Use

NOTE

- Dark red supernatant fluid, discoloration of the product or unacceptable results may indicate deterioration. Do not use the hematology calibrator if deterioration is suspected.
- Confirm that the lot number of the product matches the lot number on the assay sheet.

- Take the hematology calibrator out of the medical refrigerator and allow it to warm at a temperature of 18 to 30°C for 15 minutes before mixing.
- Mix the hematology calibrator by holding it horizontally between the palms of the hands. Roll it back and forth for 20 to 30 seconds and occasionally turn it upside down.

NOTE

- Do not mix the hematology calibrator roughly.
- Do not mix the hematology analyzer using a mixer.
- Mix the hematology calibrator until the red cells are completely suspended. Mix thoroughly but do not shake it.
- Immediately before measuring, gently turn the hematology calibrator upside down 8 to 10 times.

NOTE

Incomplete mixing of the hematology calibrator affects the measurement result.

- Calibrate the hematology analyzer using the hematology calibrator as instructed in the hematology analyzer Operator's Manual.
- Return the hematology calibrator to the medical refrigerator within 30 minutes of use.

Expected Results

For calibration values, refer to the assay sheet.

Storage

Store at 2 to 8°C.

NOTE

Do not freeze the hematology calibrator.

Expiration Date

The expiration date is shown on the container label.

Expiration Date After Opening the Package

Opened hematology calibrator should be used within 7 days.

Copyright Notice

The entire contents of this manual are copyrighted by Nihon Kohden. All rights are reserved.