Stingray™

Guidewire with Hydrophilic Coating

Directions for Use	2
Instrucciones de uso	4
Mode d'emploi	6
Gebrauchsanweisung	8
Istruzioni per l'uso	10
Gebruiksaanwijzing	12
Instruções de Utilização	14

Stingray™

Guidewire with Hydrophilic Coating

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The Stingray Guidewires are available in a nominal diameter of 0.014 in (0.36 mm) and a nominal length of 185 cm or 300 cm. In general, the Stingray Guidewires have a coiled distal end where the coil is radiopaque to facilitate advancement into and withdrawal from the vasculature using fluoroscopy. The proximal shaft is coated with polytetrafluoroethylene (PTFE) while the distal shaft in the region of the coil is coated with hydrophilic coating. The narrowed distal tip is intended to aid in guidewire control. Refer to the product label for product specifications (e.g. wire length, length of tip radiopacity, tip diameter).

Contents

Quantity	Material
1	Stingray Guidewire with Hydrophilic Coating

INTENDED USE / INDICATIONS FOR USE

The Stingray Guidewires are intended to facilitate placement of balloon dilatation catheters or other intravascular devices during percutaneous transluminal coronary angioplasty (PTCA). Stingray Guidewires are not to be used in cerebral blood vessels.

When used as part of the system consisting of the CrossBoss™ Catheter, Stingray™ LP Catheter, and Stingray Guidewire, the Stingray Guidewire is indicated for use to facilitate the intraluminal placement of conventional guidewires beyond stenotic coronary lesions (including chronic total occlusions [CTOs]) prior to PTCA or stent intervention.

CONTRAINDICATIONS

None known.

WARNINGS

- Only physicians thoroughly trained in interventional procedures should use the Stingray Guidewires.

PRECAUTIONS

- Use prior to the Expiration Date.
- The Stingray Guidewires should only be used in hospitals where emergency coronary bypass surgery can be immediately performed in the event of a potentially injurious or life threatening complication.
- Before insertion of the guidewire, administer appropriate anticoagulant and vasodilator therapy.
- The Stingray Guidewires should be handled with care. Prior to use and during the procedure, inspect the packaging and guidewire for bends, kinks, or other damage. Discontinue use if the guidewire becomes damaged.
- Exercise care during a procedure to reduce the possibility of accidental breakage, kinking, bending or coil separation.
- Do not attempt to straighten a guidewire that has been kinked.
- To reduce the potential for guidewire breakage, do not advance a kinked guidewire into a vessel or endovascular catheter.
- Do not rotate the guidewire if significant resistance is felt.
- The Stingray Guidewires should only be manipulated under fluoroscopic observation.

ADVERSE EVENTS

Potential adverse events include, but are not limited to, the following:

- Acute myocardial infarction
- Artery spasm
- Embolism
- Hemorrhage or hematoma
- Infection
- Neurological deficit
- Stroke
- Vessel trauma requiring surgical repair or intervention

HOW SUPPLIED

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place.

OPERATIONAL INSTRUCTIONS

Supplies

The following supplies are not provided and need to be available and prepped prior to use of the Stingray Guidewire:

- Guiding catheter
- Angiographic imaging supplies (i.e. radiopaque contrast, manifold, tubing, etc.)
- Sterile heparinized normal saline
- Other supplies as needed to complete the established laboratory protocol

Product Preparation

1. Remove the Stingray™ Guidewire from the sterile packaging.
2. After removing the guidewire from the holder tube, inspect to make sure it is not damaged.
3. If desired and according to standard practice, gently shape the flexible tip of the guidewire.
4. After pulling the guidewire out of the body, wipe it with gauze soaked in heparinized saline.
5. Upon completion of the use of the guidewire, discard according to hospital standards and national legal requirements.

Product Operation

Front Load the Guidewire into the Endovascular Catheter

1. Prep the Stingray Guidewire per preparation instructions.
2. Insert a guidewire insertion tool per manufacturer's instructions as to communicate with the guidewire port of the endovascular catheter.
3. Carefully insert distal guidewire tip through the insertion tool and into the endovascular catheter.
4. Remove the insertion tool by withdrawing it over the proximal end of the guidewire.
5. Advance the Stingray Guidewire through the endovascular catheter under fluoroscopy to the intended target vessel. Make sure the guidewire moves smoothly and torque the wire if necessary.
6. A guidewire torque device may be used per manufacturer's instructions.

Backload the Endovascular Catheter onto the Guidewire

1. Access the vasculature with a guiding catheter and hemostatic Y-adapter per established laboratory protocol.
2. Assure the Stingray Guidewire has been prepped per instructions.
3. Advance the guidewire into the vasculature to the target vessel using fluoroscopic guidance. Make sure the wire moves smoothly and torque the wire if necessary.
4. A guidewire torque device may be used per manufacturer's instructions.
5. Once the target location has been accessed, load the distal tip of the endovascular catheter onto the proximal portion of the guidewire.
6. Loosen the knob on the haemostatic Y-adapter.
7. Track the endovascular catheter over the Stingray Guidewire and into the proximal guiding catheter lumen.
8. Tighten the haemostatic Y-adapter to create a seal around the endovascular catheter. Do not over-tighten.
9. Advance the endovascular catheter over the guidewire using fluoroscopic guidance per established laboratory protocol.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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DEVICE DESCRIPTION

The Stingray Extension Wire is a stainless steel wire with a connector at the distal end. The extension wire has a nominal diameter of 0.014 in (0.36 mm) and is 145 cm in length. The shaft is coated with PTFE (polytetrafluoroethylene). It is exclusively compatible with the extendable Stingray Guidewires. The guide pin on the proximal end of the extendable guidewire is inserted into the center lumen of the connector of the extension wire. Once the guide pin is inserted, the extension wire is rotated in a clockwise direction (with respect to the guidewire).

Contents

Quantity	Material
1	Stingray Extension Wire

INTENDED USE / INDICATIONS FOR USE

The attachment of the Stingray Extension Wire to the extendable guidewire creates a guidewire that can be used to exchange a catheter without removing the guidewire from the vessel. When the exchange is complete, the guidewire extension can be detached and the original guidewire can again be used in a conventional manner.

CONTRAINDICATIONS

None known.

WARNINGS

- Only physicians thoroughly trained in interventional procedures should use the extension wire.

PRECAUTIONS

- Use prior to the Expiration Date.
- The extension wire is exclusively compatible with the extendable Stingray guidewires.
- Do not rotate or manipulate the extended guidewire once it is attached. This may cause the extension wire to detach from the guidewire.
- Do not force the proximal end of the guidewire into the connector on the extension wire. Attachment is achieved by rotating and not by pushing.
- The extension wire should be handled with care. Prior to use and during the procedure, inspect the packaging and extension wire for bends, kinks, or other damage. Discontinue use if the extension wire becomes damaged.
- Do not attempt to straighten an extension wire that has been kinked.
- Catheter exchange should only be done under fluoroscopic observation.

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- No additional supplies necessary

Product Preparation

- Remove the extension wire from the sterile packaging.
- After removing the extension wire from the holder tube, inspect to make sure it is not damaged.

Product Operation

- Remove the torque device (if applicable) from the end of the extendable guidewire.
- Maintain position of the extendable guidewire and confirm its position in the vessel by fluoroscopy.
- Insert the guide pin on the proximal end of the extendable guidewire into the distal end of the extension wire. Once the guide pin is inserted, rotate the extension wire in a clockwise direction (with respect to the guidewire) until the two wires are pulled together securely.
- Check for secure attachment by pulling on both wires.
- Withdraw the catheter to be exchanged, over the extended wire, while maintaining wire position in the vessel.
- Inspect the extended wire to make sure it has not been damaged. Confirm the position of the extendable guidewire in the vessel by fluoroscopy. Prepare the catheter to be exchanged according to the manufacturer's *Directions For Use*. Advance the new catheter over the extended wire while maintaining position of the extendable guidewire and guiding catheter.

NOTE: Keep the system as straight as possible to facilitate the exchange.

NOTE: Do not torque or manipulate the guidewire extension. This may cause the connector to unscrew and may lead to detachment of the wires.

7. Detach the guidewire extension by rotating the extension wire in a counterclockwise direction (with respect to the guidewire) until the guidewire is released.
8. Once the extension wire is detached, the extendable guidewire may be torqued and manipulated per standard techniques.
9. The same extension wire may be used for additional exchanges provided that sterility is maintained.
10. Upon completion of the use of the extension wire, discard according to hospital standards and national legal requirements.

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© 2019 Boston Scientific Corporation or its affiliates. All rights reserved.