
Medtronic

EN

Tracer® Pointer (9735317)

Device description

The Tracer® Pointer is a handheld, sterile, single-use disposable device for use during stereotactic surgical navigation. A computer-assisted surgery system tracks the device in an electromagnetic field to localize patient anatomy (that is, to identify the point on the patient images that corresponds to the point touched on the patient anatomy). The device is compatible with AxiEM™-enabled StealthStation® S7® and i7® systems.

Intended use

The Tracer® Pointer is intended for use in completing Tracer® patient registration with a Medtronic computer-assisted surgery system. The Tracer® Pointer is intended to localize patient anatomy during computer-assisted stereotactic surgery. The handheld, sterile, single-use device may come in contact with central nervous system fluids as required by the surgical procedure and at the surgeon's discretion, based on their experience and training.

The Tracer® Pointer is intended for use only with Medtronic computer-assisted surgery systems using electromagnetic navigation and cranial software.

Indications for use

The Tracer® Pointer is indicated for use in registering patient anatomy and in precisely locating anatomical structures during navigated cranial procedures in adult and pediatric patients. It is indicated for use with Medtronic computer-assisted surgery systems using electromagnetic navigation and cranial software.

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
Printed in the USA

**R_x Only**

Medtronic Navigation, Inc.
826 Coal Creek Circle
Louisville, CO 80027 U.S.A.
Phone: +1 (720) 890-3200
Fax: (720) 890-3500

Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
NETHERLANDS
Tel. 31 45 566 80 00

**EN**

Passive Catheter Introducer (PCI)

- ① PCI tip (9731116)
- ② PCI handle (9730092)
- ③ Assembled handle and tip

PCI tip ① is made of stainless steel and is single-use only.

PCI handle ② is made of anodized aluminum and stainless steel.

The single-use, sterile spheres shown in ③ are sold separately.

Intended Use

The Passive Catheter Introducer (PCI) is intended to enable the use of image-guidance for the placement of ventricular catheters and shunts with a StealthStation® system.

Contraindications

There are no known contraindications for this device.

Side Effects

Potential side effects associated with shunt placement procedures include: minor bleeding, hematoma without neurological compromise, infection, overdraining, and shunt migration.

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9660202 packaged as 9731132

Device description

The Skull-Mount Patient Tracker (SMPT) is a sterile single-use disposable electromagnetic field detector with self-tapping screws. The device is compatible with AxiEM™-enabled StealthStation® S7® and i7® systems and Fusion® systems.

The navigation system's mobile emitter generates a low-intensity magnetic field to locate coils mounted within the tracker. Then, the software uses the tracker as a reference point to display the location of instruments in patient image planes and other anatomical renderings.

Intended use

The Skull-Mount Patient Tracker is intended to be attached to the patient's skull with two screws for the purpose of tracking patient location during image guided surgery. The Skull-Mount Patient Tracker is intended for use only with a Medtronic computer-assisted surgery systems using electromagnetic navigation and cranial and ENT software.

Indications for use

The Skull-Mount Patient Tracker is indicated for use as a stereotactic reference frame to track the position of patient anatomy during navigated cranial and ENT procedures. It is indicated for use with Medtronic computer-assisted surgery systems using electromagnetic navigation and cranial and ENT software.

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Contraindications

Contraindicated for use where stereotactic image-guided surgery is contraindicated.

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
Navigation Pointer (9735319)

Device description

The Navigation Pointer is a handheld, sterile, single-use disposable device for use during stereotactic surgical navigation. A computer-assisted surgery system tracks the device in an electromagnetic field to localize patient anatomy (that is, to identify the point on the patient images that corresponds to the point touched on the patient anatomy). The device is compatible with AxiEM™-enabled StealthStation® S7® and i7® systems.

Intended use

The Navigation Pointer is intended to localize patient anatomy during computer-assisted stereotactic surgery. The handheld, sterile, single-use device may come in contact with central nervous system fluids as required by the surgical procedure and at the surgeon's discretion, based on their experience and training.

The Navigation pointer is intended for use only with Medtronic computer-assisted surgery systems using electromagnetic navigation and cranial software.

Indications for use

The Navigation Pointer is indicated for use in precisely locating anatomical structures during navigated cranial procedures. It is indicated for use with Medtronic computer-assisted surgery systems using electromagnetic navigation and cranial software.

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The Navigation Pointer is indicated for use in the following procedures::

- General tumor resections
- Pediatric tumor resections
- Skull base procedures

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
Australia

Medtronic Australasia Pty Ltd
97 Waterloo Road
North Ryde, NSW 2113
Australia

Printed in the USA

**EN****Stylet, 23 cm (9735428)****Device description**

The stylet is a handheld, sterile, single-use disposable device for stereotactic surgical navigation. A computer-assisted surgery system tracks the device in an electromagnetic field to localize patient anatomy (that is, to identify the point on the patient images that corresponds to the point touched on the patient anatomy). The device is compatible with AxiEM™-enabled StealthStation® S7® and i7® systems running Cranial software.

The navigation system's mobile emitter generates a low-intensity magnetic field to locate tracking coils at the tip of the instrument. Then, the software displays the location of the instrument's tip within multiple patient image planes and other anatomical renderings.

Intended use

The stylet is intended to aid in procedures using electromagnetic surgical navigation, localize patient anatomy, and function down the working lumen of a catheter or electrode. The handheld, sterile, single-use device may come in contact with central nervous system fluids as required by the surgical procedure and at the surgeon's discretion, based on their experience and training.

The stylet is intended for use only with Medtronic computer-assisted surgery systems using electromagnetic navigation and cranial software.

Indications for use

The stylet is indicated for use in precisely locating anatomical structures during navigated cranial procedures. It is indicated for use with Medtronic computer-assisted surgery systems using electromagnetic navigation and cranial software.

The Medtronic computer-assisted surgery system and its associated applications are intended as an aid for precisely locating anatomical structures in either open or percutaneous procedures. Their use is indicated for any medical condition in which the use of stereotactic surgery may be appropriate, and where reference to a rigid anatomical structure, such as the skull, can be identified relative to a CT-based or MR-based model, fluoroscopy images, or digitized landmarks of the anatomy.

The stylet is indicated for use in the following procedures:

- General ventricular catheter placement
- Pediatric ventricular catheter placement
- Tumor resections

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All Rights Reserved
©2015 Medtronic Navigation, Inc.