
DIGITAL THERMOMETER PROLIFE PDT100

Thank you for choosing our products! Having assessed the reliability and quality of the Digital thermometer Prolife **PDT100**, you will definitely stay loyal to our products in future.

The device provides safe, accurate and rapid measurement of body temperature. The temperature can be measured in the mouth (oral measurement method), in the anus (rectal measurement method) and in the armpit (axillary measurement method).

If you have any questions regarding body temperature, please consult your doctor.

Before using the Digital thermometer Prolife **PDT100**, please read the Instruction Manual carefully and keep it for future reference.

FEATURES

- Acoustic signal.
- Automatic self-test.
- Ensures accurate readings with every measurement.
- Automatic notification of an internal error affecting the thermometer proper operation (Err).
- Rapid temperature measurement within 60 seconds.
- Last measurement reading auto save.
- Automatic switch off.
- 100% waterproof when immersed at a depth down to 10 cm within 1 minute.

SAFETY PRECAUTIONS

- This thermometer is intended for clinical and domestic use only. Do not use the thermometer for purposes other than intended.
- Use the thermometer only to measure human body temperature.
- The thermometer must be in full contact with the patient during the measurement.
- The thermometer contains small parts (battery and battery compartment cover). To avoid swallowing or suffocation, do not leave the thermometer near the child unattended. Keep the thermometer out of the reach of children.
- Do not expose the thermometer to direct sunlight and high temperatures, as well as mechanical stress.
- If any malfunction is detected, stop using the thermometer and contact a specialized service center. Under no circumstances should you attempt to repair the device yourself.
- Before starting the measurement, the user and the thermometer should stay in the room with constant environmental conditions for at least 30 minutes to ensure measurement accuracy.
- The materials used are intended for patient contact and comply with ISO 10993-5 and ISO 10993-10 standards. Non-toxic, non-allergenic.

Dispose of the batteries used in this product in accordance with the local waste disposal rules.

WARNING

- Seek medical advice in case of fever or continuous fever, particularly with regard to small children.
- Many factors affect temperature readings, including physical activity, drinking hot or cold beverages before taking a measurement, as well as the measurement method.
- Self-diagnosis based on the results obtained is dangerous and may cause deterioration of health. Consult your doctor to correctly interpret the results obtained.
- Do not measure the temperature with a wet thermometer, as readings may be inaccurate.
- When measuring temperature in your mouth, do not clench the teeth with all your might. Doing so may cause damage to the thermometer and/or injury.
- Remove the battery if the thermometer is not going to be used for a long time.
- Do not use the thermometer in places where strong static electricity or electromagnetic fields are generated.
- Do not use the rectal measurement method in case of rectal disorders. This may cause deterioration in condition.

TURNING ON THE THERMOMETER

To turn on the thermometer, push the ON/OFF button.

A short beep will alert you when the internal check is complete. All symbols will be displayed. Then the °C symbol will appear on the display.

The thermometer is now ready for use.

FUNCTIONAL SELF-TEST

Proper functioning of the thermometer is tested automatically each time it is turned on. If «Err» appears on the display after the self-test is complete, please contact your nearest service center or point of purchase.

TEMPERATURE MEASUREMENT

During the measurement, the current temperature will be constantly displayed, and the °C symbol will flash.

When you stop measuring by the thermometer’s acoustic signal (beep), you can only evaluate whether the reported value falls within the range of «normal» temperatures for the selected measurement method.

WARNING

To obtain the actual temperature value, do not stop measuring immediately after the beep. You should wait during the period of time recommended for each measurement method.

MEMORY FUNCTION

To display the last measurement reading, press and hold the ON/OFF button until the temperature value and the «M» symbol appear. Release the button to switch to the temperature measurement mode.

BODY TEMPERATURE MEASUREMENT METHODS

ORAL (IN THE MOUTH)

Place the thermometer in one of the two pockets under the tongue, to the left or right of the root of the tongue. The measuring sensor must be in good contact with the mucosa. Close your mouth and breath evenly through the nose to prevent the measurement from being influenced by inhaled/exhaled air. The «normal» temperature range is 35.5–37.5 °C. **Recommended measuring time: 60 seconds.**

RECTAL (IN THE ANUS)

This is the most reliable measurement method. It is especially appropriate for infants and small children. Carefully insert the tip of the thermometer 2 to 3 cm into the anal orifice. The «normal» temperature range is 36.6–38 °C. Disinfect the thermometer after use. Recommended measuring time: 60 seconds. **Recommended measuring time: 60 seconds.**

AXILLARY (IN THE ARMPIT)

Place the measuring sensor closer to the centre of the armpit and gently press the arm against the body. The «normal» temperature range is 34.7–37.3 °C.

From a medical point of view, this method, as a rule, does not provide high accuracy of readings due to the contact of the armpit skin with the surrounding air. For more reliable results, it is recommended to wait another 2 minutes after the beep.

BATTERY REPLACEMENT

symbol appears on the display to indicate that the battery must be replaced.

To do this, open the battery compartment cover and remove the discharged battery. Insert a new battery, observing the polarity. Please make sure that the new battery is of the same type as the old one.

CLEANING AND DISINFECTING

To clean and disinfect the thermometer at home, it is recommended to wipe it with a cloth (cotton swab) soaked in 70% ethyl or isopropyl alcohol or rinse with water.

To clean and disinfect the thermometer in a clinical setting, it is recommended to wipe it with disinfectants. For these purposes, you can use wipes soaked in disinfectants, 70% ethyl or isopropyl alcohol. After treatment, remove the disinfectant residues from the surface of the thermometer using wet wipes moistened with water, wipe dry and store them in a closed container.

The thermometer is waterproof and for cleaning purposes can be completely (only undamaged and not dismantled) immersed in a disinfectant solution (70% isopropyl alcohol) for up to 1 minute and a depth down to 10 cm. After treatment, wipe the device with a damp cloth moistened with water

WARNING

Never use abrasive cleaners, solvents, or benzine as this may cause damage to the thermometer.

SPECIFICATIONS

ITEMS	DESCRIPTION
Type	Digital thermometer
Measurement Mode	Direct mode
Measurement Range	32.0 °C~42.9 °C Temp <32.0 °C: «L» Temp >42.9 °C: «H»
Resolution	0.1 °C
Accuracy	±0.1 °C for measuring range between 35.5 °C ~ 42.0 °C ±0.2 °C outside the specified range
Measurement Position	oral, axillary or rectal
Power Source	1.5V battery; type LR41, 192, SR41, 392 or equivalent
Power Consumption	Max. 0.15mW in measuring mode
Battery Life	Approx. 2-year (3 measurements per day)
Battery icon	YES; Icon on the display
Beeper	Ready for measurement Notification stating that the temperature rise rate of the measuring sensor is less than 0.1 °C in 8 seconds
Memory	one: the last measuring result
Auto Power-Off	8 min 40sec ± 1 min after measurement completed.
Transport & Storage conditions	-25 °C to +60 °C; 10 to 95% RH
Self-Test Function	Automatic self-test for every turns-on Displays «Err» if the system is failed
Dimensions	129 mm*21 mm*11 mm
Weight	11 g
Life time	5 years

LIST OF SYMBOLS

SYMBOL	MEANING
	The product meets the basic requirements of the Medical Devices Directive 93/42/EEC.
	The product complies with the basic requirements of TR CU 020/2011 «Electromagnetic compatibility of technical means».
	WEEE (Waste Electrical and Electronic Equipment Directive). The symbol on the product or its package means that this product does not fall under the category of domestic waste. To avoid possible damage to the environment and human health, separate such wastes from others and dispose of them in accordance with accepted standards.
	Warning.
	Read the Instruction Manual before use.
	Type BF applied part.
	Manufacturer.
	Authorized manufacturer's representative in the European Union.

WARRANTY

The warranty period is indicated in the warranty certificate. The warranty period is established from the date of sale, providing that all operation conditions specified in this instruction manual are strictly observed. The warranty is valid subject to availability of a warranty card filled out by a duly authorized representative confirming the date of sale and also of the relevant receipt. .

Warranty and free maintenance service is not provided in case of:

- use of the device in violation of the operating instructions;
- damage caused by intentional or erroneous actions of the consumer;
- evidence of mechanical impact, dents, cracks, chips, etc. on the thermometer housing, evidence of opening the device, disassembly, unauthorized repair, ingress of moisture, exposure to aggressive agents or any other unauthorized alterations of the device, and in other cases of violation of storage, cleaning, transportation and operation rules specified in the Instruction Manual.

WARNING

Follow the instructions precisely to ensure reliable and long-term operation of the device.

In case of abnormal operation of the device, please contact your retailer.

For repair and maintenance, please contact a specialized service center.

The manufacturer reserves the right to make alterations in the design of the device.

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
DIGITAL THERMOMETER PROLIFE PDT150

Thank you for choosing our products! Having assessed the reliability and quality of the Digital thermometer Prolife **PDT150**, you will definitely stay loyal to our products in future.

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Before using the Digital thermometer Prolife **PDT150**, please read the Instruction Manual carefully and keep it for future reference.

FEATURES

- Acoustic signal.
- Automatic self-test. Ensures accurate readings with every measurement. Automatic notification of an internal error affecting the thermometer proper operation (Err).
- Rapid temperature measurement within 60 seconds.
- Flexible tip for comfortable measurement.
- Last measurement reading auto save.
- Automatic switch off.

SAFETY PRECAUTIONS

- This thermometer is for domestic use only. Do not use the thermometer for purposes other than intended.
- Use the thermometer only to measure human body temperature.
- The thermometer must be in full contact with the user during the measurement.
- The thermometer contains small parts (battery and battery compartment cover). To avoid swallowing or suffocation, do not leave the thermometer near the child unattended. Keep the thermometer out of the reach of children.
- Do not expose the thermometer to direct sunlight and high temperatures, as well as mechanical stress.
- If any malfunction is detected, stop using the thermometer and contact a specialized service center. Under no circumstances should you attempt to repair the device yourself.
- Before starting the measurement, the user and the thermometer should stay in the room with constant environmental conditions for at least 30 minutes to ensure measurement accuracy.
- The materials used are intended for user contact and comply with ISO 10993-5 and ISO 10993-10 standards. Non-toxic, non-allergenic.

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TURNING ON THE THERMOMETER

To turn on the thermometer, push the ON/OFF button. A short beep will alert you when the internal check is complete. All symbols will be displayed. Then the °C symbol will appear on the display. The thermometer is now ready for use.

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TEMPERATURE MEASUREMENT

During the measurement, the current temperature will be constantly displayed, and the °C symbol will flash. When you stop measuring by the thermometer's acoustic signal (beep), you can only evaluate whether the reported value falls within the range of «normal» temperatures for the selected measurement method.

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BATTERY REPLACEMENT

symbol appears on the display to indicate that the battery must be replaced. To do this, open the battery compartment cover and remove the discharged battery. Insert a new battery, observing the polarity. Please make sure that the new battery is of the same type as the old one.

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SPECIFICATIONS

ITEMS	DESCRIPTION
Type	Digital thermometer
Measurement Mode	Direct mode
Measurement Range	32.0 °C~42.9 °C Temp <32.0 °C: «L» Temp >42.9 °C: «H»
Resolution	0.1 °C
Accuracy	±0.1 °C for measuring range between 35.5 °C ~ 42.0 °C ±0.2 °C outside the specified range
Measurement Position	oral, axillary or rectal
Power Source	1.5V battery; type LR41, 192, SR41, 392 or equivalent
Power Consumption	Max. 0.15mW in measuring mode
Battery Life	Approx. 2-year (3 measurements per day)
Battery icon	YES; Icon on the display
Beeper	Ready for measurement Notification stating that the temperature rise rate of the measuring sensor is less than 0.1 °C in 8 seconds one: the last measuring result
Memory	one: the last measuring result
Auto Power-Off	8 min 40sec ± 1 min after measurement completed.
Operating conditions	10 °C to 40 °C; 10 to 95% RH
Transport & Storage conditions	-25 °C to +60 °C; 10 to 95% RH
Self-Test Function	Automatic self-test for every turns-on Displays «Err» if the system is failed
Dimensions	127 mm*21 mm*11 mm
Weight	13 g
Life time	5 years

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WARRANTY

The warranty period is indicated in the warranty certificate. The warranty period is established from the date of sale, providing that all operation conditions specified in this instruction manual are strictly observed. The warranty is valid subject to availability of a warranty card filled out by a duly authorized representative confirming the date of sale and also of the relevant receipt. Warranty and free maintenance service is not provided in case of:

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WARNING

Follow the instructions precisely to ensure reliable and long-term operation of the device. In case of abnormal operation of the device, please contact your retailer. For repair and maintenance, please contact a specialized service center. The manufacturer reserves the right to make alterations in the design of the device.

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