**REF**

6912, 6913, 6914, CV-6912, CV-6913, CV-6914, 6932, 6933, 6934

CDI™ H/S Cuvette***CDI™ H/S кювета******Kiveta CDI™ H/S******Kyveta CDI™ H/S******CDI™ H/S-cuvette******CDI™ H/S cuvet******CDI™ H/S küvett******CDI™ H/S -kyveti******Cuvette H/S CDI™******CDI™ H/S Küvette******Κυψελίδα H/S CDI™******CDI™ H/S küvetta******Cuvette H/S CDI™******CDI™ H/S kivete******CDI™ H/S kiuveté******CDI™ H/S-kyvette******Kuweta H/S CDI™******Cuvete de H/S CDI™******Cuvă CDI™ H/S******ΚΙΟΒΕΤΑ CDI™ H/S******Kiveta CDI™ H/S******CDI™ H/S kyveta******Kiveta CDI™ H/S******Cubeta de H/S CDI™******CDI™ H/S-kyvett******CDI™ H/S Küveti******Κювета H/S CDI™******CDI™ H/S 試管******CDI™ H/S 光析管******CDI™ H/S キューベット***

This page is intentionally left blank.

CDI™ H/S Cuvette

Table of Contents

English.....	5
Български	9
Hrvatski	13
Česky	17
Dansk	21
Nederlands	25
Eesti keeles	29
Suomi	33
Français.....	37
Deutsch	41
Ελληνικά.....	45
Magyar	49
Italiano.....	53
Latviešu	57
Lietuvių kalba	61
Norsk	65
Polski.....	69
Português	73
Română.....	77
Русский	81
Srpski	85
Slovenčina.....	89
Slovenščina	93
Español	97
Svenska.....	101
Türkçe	105
Українська	109
简体中文	113
繁體中文	117
日本語.....	121
Graphic Symbols and Definitions	125

This page is intentionally left blank.

English

Instructions for Use

CDI™ H/S Cuvette

**Models 6912, 6913, 6914, CV-6912, CV-6913,
CV-6914, 6932, 6933 and 6934**

THE OPERATOR'S MANUAL FOR THE CDI BLOOD PARAMETER MONITORING SYSTEM SHOULD BE READ AND UNDERSTOOD IN ITS ENTIRETY BEFORE ATTEMPTING TO USE THE CDI H/S CUUVETTE WITH THE CDI BLOOD PARAMETER MONITORING SYSTEM.

Product Description

The CDI Hematocrit/Saturation (H/S) Cuvette is a single use bidirectional flow-through device which, when used with the CDI Blood Parameter Monitoring System, provides continuous, in-line monitoring of oxygen saturation, hemoglobin and hematocrit. It is supplied with or without extension tubing for insertion directly into the extracorporeal circuit at the time of use. It may also be pre-connected in extracorporeal circuits by the tubing pack manufacturer.

The CDI H/S Cuvette contains an optical window that provides a means of consistent optical connection between the CDI H/S Probe and CDI H/S Cuvette. A magnet on the CDI H/S Cuvette ensures the correct connection between the CDI H/S Probe and the CDI H/S Cuvette. The CDI H/S Cuvette contains an Ethylene Oxide sterile and non-pyrogenic fluid pathway in an unopened, undamaged package.


Figure 1 – CDI H/S Cuvette

Intended Use

The CDI H/S Cuvette is a disposable conduit positioned in the blood flow pathway of an extracorporeal bypass circuit. It is a component of the CDI blood parameter monitoring system that is used when continuous monitoring of blood hematocrit, hemoglobin and oxygen saturation is desired. The device may be used in cardiopulmonary bypass procedures lasting up to six (6) hours.

Contraindications

None

WARNINGS:

- CDI H/S CUVETTES ARE INTENDED FOR SINGLE USE ONLY. DO NOT REUSE. USED CUVETTES ARE CONTAMINATED AND CANNOT BE RSTERILIZED.
 - DO NOT RSTERILIZE THE CDI H/S CUVETTE.
 - DO NOT REPROCESS CDI H/S CUVETTES. REPROCESSING MAY COMPROMISE THE STERILITY, BIOCOMPATIBILITY AND FUNCTIONAL INTEGRITY OF THE DEVICE.
 - MAINTAIN ADEQUATE LEVELS OF ANTICOAGULATION DURING EXTRACORPOREAL CIRCULATION BY MONITORING ACTIVATED CLOTTING TIME (ACT) OR OTHER APPROPRIATE MEASUREMENT. INADEQUATE LEVELS OF ANTICOAGULANT MAY LEAD TO CLOTTING IN THE DEVICE CAUSING INACCURATE MEASUREMENTS.
 - THE FOLLOWING MAY POTENTIALLY CAUSE INACCURACIES IN DISPLAYED VALUES:
 - USE OF CERTAIN INTRAVASCULAR DYES SUCH AS INDOCYANINE GREEN (CARDIOGREEN) AND METHYLENE BLUE
 - ELEVATED LEVELS OF BLOOD SUBSTANCES INCLUDING IRREGULAR CELL MORPHOLOGIES, PROTEIN LEVELS, PLASMA FREE HEMOGLOBIN AND BILIRUBIN
 - BLOOD CONDITIONS SUCH AS HEMOGLOBINOPATHIES, THALASSEMIA AND A VARIETY OF ANEMIC CONDITIONS (SICKLE CELL, IRON DEFICIENCY, MACROCYTIC)
- INDEPENDENT EXTERNAL LAB ANALYSIS IS REQUIRED FOR ACCURATE DETERMINATION OF ALL MEASURED PARAMETERS NEEDED TO GUIDE THERAPEUTIC DECISIONS.
- USE CAUTION WHEN ADMINISTERING NOVEL PHARMACOLOGICAL AGENTS WHEN THE USER IS UNFAMILIAR WITH THE POTENTIAL EFFECT OF SUCH AGENTS ON THE CDI H/S CUVETTE.
 - PRESSURE IN THE CIRCUIT EXCEEDING 1,000mmHg (133.3kPa) OR BELOW -150mmHg (-20kPa) MAY CAUSE LEAKS OR DAMAGE TO THE DEVICE.
 - DO NOT USE SOLVENTS SUCH AS ALCOHOL, ETHER OR ACETONE. THESE SOLVENTS MAY CAUSE DAMAGE IF USED IN OR ON THE DEVICE.
 - THE PRESENCE OF AIR BUBBLES ON THE OPTICAL WINDOW OF THE CDI H/S CUVETTE CAN RESULT IN REDUCED MEASUREMENT ACCURACY. INTERMITTENT BUBBLES, ONCE MOVED FROM THE OPTICAL WINDOW, WILL NOT AFFECT THE LONG TERM ACCURACY OF THE SYSTEM.
 - THE CDI MONITORING SYSTEM SHOULD ONLY BE USED WHEN THERE IS BLOOD FLOW THROUGH THE EXTRACORPOREAL CIRCUIT. TO PERFORM ACCURATELY, THE H/S CUVETTE REQUIRES BLOOD FLOW RATES SHOWN IN THE FOLLOWING TABLE:

CDI H/S Cuvette Size	Minimum Flow (Liters per Minute)	Maximum Flow (Liters per Minute)
½"	1.0	7.0
⅜"	0.5	4.0
¼"	0.2	1.5

Table 1

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A BLOOD FLOW RATE OUTSIDE THE RANGES LISTED IN TABLE 1 MAY RESULT IN LESS THAN OPTIMAL SATURATION/HEMATOCRIT READINGS. RESTORATION OF THE FLOW RATE ABOVE THE MINIMUM WILL RETURN THE CDI H/S CUVETTE READINGS TO NORMAL.

- CDI H/S CUVETTE MODELS 6932, 6933 AND 6934: PHTHALATE: THIS DEVICE CONTAINS DEHP. BASED ON ANIMAL STUDIES, SIGNIFICANT EXPOSURE TO DEHP MAY INTERFERE WITH THE NORMAL DEVELOPMENT OF THE MALE REPRODUCTIVE TRACT. FOR CHILDREN, PREGNANT AND NURSING WOMEN, ALTERNATIVE DEVICES MAY BE APPROPRIATE.

CAUTIONS:

- THE CONTENTS OF THE CDI H/S CUVETTE PACKAGE ARE STERILE AND NON-PYROGENIC UNLESS THE UNIT PACKAGE IS DAMAGED OR OPENED. INSPECT THE DEVICE AND PACKAGE CAREFULLY. DO NOT USE IF THE PACKAGE OR DEVICE IS DAMAGED OR THE CAPS ARE NOT IN PLACE.
- DO NOT USE THE CDI H/S CUVETTES AFTER THE DATE PRINTED ON THE PACKAGE LABEL. USE BEYOND THE DATE MAY RESULT IN INACCURATE PERFORMANCE.
- AVOID TOUCHING OR SCRATCHING THE OPTICAL SURFACES OF THE CDI H/S CUVETTE OR THE OPTICAL PROBE. FINGERPRINTS OR SCRATCHES ON THESE SURFACES MAY REDUCE THE ACCURACY OF THE SYSTEM. INSTRUCTIONS FOR CLEANING CAN BE FOUND IN THE MAINTENANCE SECTION OF THE CDI MONITORING SYSTEM OPERATOR'S MANUAL.
- THE CDI H/S CUVETTE SHOULD NOT BE USED OUTSIDE THE HOSPITAL ENVIRONMENT OR EXPOSED TO EXTREMES OF TEMPERATURE AND HUMIDITY
- FEDERAL (USA) LAW RESTRICTS THE DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN OR OTHER LICENSED PRACTITIONER.
- FOR USE BY PROFESSIONALLY TRAINED PERSONNEL ONLY
- THE ACCURACY OF RESULTS IS DEPENDENT ON PERIODIC COMPARISON TO A LABORATORY REFERENCE SAMPLE. VERIFY ITS ACCURACY BY INDEPENDENT MEANS BEFORE INITIATING TREATMENT.

INSTRUCTIONS FOR USE

Circuit Preparation

1. Examine the packaging for any damage. Do not use if the package is damaged.
 2. Remove the product from its package and unwrap it using aseptic, sterile technique.
 3. Remove the end caps from the barbed ports and discard.
 4. Ensure that the cuvette is oriented such that the optical window is pointed downward when inserted into the line. This prevents air bubbles from impacting the parameter readings when connected to the CDI H/S probe.
- Note:** The CDI H/S probe is designed for use on the venous side only.

5. Connect the appropriate size tubing to the CDI H/S Cuvette by pushing the tubing over at least two barbs as shown in Figure 2.


Figure 2 – tube insertion

6. Tie-band securely between the first and second barbs. Tighten the connection.
7. The CDI H/S Cuvette should be primed with a crystalloid solution containing no blood, plasma, or blood-derived products.
8. Prior to starting treatment, clip the CDI H/S probe to the CDI H/S Cuvette by inserting the tab on the cable end of the CDI H/S probe into the receptacle on the CDI H/S Cuvette as shown in Figure 3.


Figure 3 – CDI H/S Probe and Cuvette

9. Press the CDI H/S probe and CDI H/S Cuvette together until the release lever latches.
10. Once the CDI H/S probe is secured, verify that the CDI Blood parameter monitoring system recognizes the CDI H/S Cuvette.

Disposal

The product must be disposed of according to hospital policies and procedures in accordance with regulations for proper disposal of medical waste.

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And CDI™ are trademarks of Terumo and may be registered.
© Copyright Terumo Cardiovascular Systems Corporation
Terumo Cardiovascular Systems Corp.,
Terumo Corporation,
Terumo Australia PTY LTD