COULTER DxH Diluent

REF 628017

TABLE OF CONTENTS

English	2
Français (FR)	4
Deutsch (DE)	6
Italiano (IT)	8
Español (ES)	10
Português Portugal (PT-PT)	12
Dansk (DA)	14
Svenska (SV)	16
Norsk (NO)	18
Ελληνικά (EL)	20
日本語 (JP)	22
Lietuviškai (LT)	24
Magyar (HU)	26
Polski (PL)	28
Čeština (CZ)	30
Slovenčina (SK)	32
한국어 (KO)	34
Türkçe (TR)	36
Русский (RU)	38
Hrvatski (HR)	40
Български (BG)	42
中文 (ZH-TW)	44
Română (RO)	46
Србија (SR)	48
Українська (UK)	50
Português Brasil (PT-BR)	52
Nederlands	54
Tiếng Việt (VI)	56

COULTER DxH Diluent

REF 628017

For *In Vitro* Diagnostic Use
Rx Only in the U.S.A.

INTENDED USE

For use as an isotonic buffered diluent in conjunction with a cyanide-free lytic agent for counting and sizing blood cells on all UniCel DxH Coulter Cellular Analysis Systems. Use only with COULTER DxH Cell Lyse. Refer to your instrument product manuals and/or online help, as applicable. For use as a rinsing agent on the UniCel DxH Slidemaker Stainers.

SUMMARY

Blood cell analysis comprises diluting a whole-blood sample with a solution that functions as a diluent. The diluent provides the ability to analyze portions of the diluted blood sample for different blood cell types and the measurement of hemoglobin. COULTER DxH Diluent also acts as a rinsing agent on all UniCel DxH Coulter Cellular Analysis Systems.

PRINCIPLE

The chemical composition of COULTER DxH Diluent provides an isotonic solution with buffering capacity for cell component stabilization and measurement of hemoglobin. Biocides are included for product preservation against microbial growth. Alkaline salts provide a suitable osmolality so as not to adversely affect cell volume (iso-osmotic).

When used on the Unicel DxH Coulter Cellular Analysis Systems in conjunction with COULTER DxH Cell Lyse, COULTER DxH Diff Pack, and COULTER DxH Retic Pack, the diluent is used for erythrocyte, hemoglobin, platelet, NRBC and reticulocyte measurements, as well as, for the measurement and differentiation of leukocytes into five subpopulations – neutrophils, lymphocytes, monocytes, eosinophils, and basophils.

REACTIVE INGREDIENTS

Sodium Sulfate13.73g/L
Sodium Chloride1.04g/L
Tetracaine HCL0.02g/L
Imidazole2.85g/L

WARNING AND PRECAUTIONS

- Do not inhale and/or ingest.
- Avoid eye and skin contact. In case of eye or skin contact, flush affected area with copious amounts of water for at least 15 minutes.
- DO NOT REUSE CONTAINER.

GHS HAZARD CLASSIFICATION		
COULTER DxH Diluent	DANGER	Causes mild skin irritation.
		Obtain special instructions before use.
		Wear protective gloves, protective clothing and eye/face protection.
		IF exposed or concerned: Get medical advice/attention.
		May produce an allergic reaction.
		Imidazole <1%
		Reaction Mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC# 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC# 220-239-6](3:1) <0.05%
		May damage fertility or the unborn child.

STORAGE, STABILITY, AND DISPOSAL

- Store Coulter DxH Diluent at 2-40°C.
- Do not use product past expiration date.
- Use product at temperatures stated in the instrument product manuals and/or online help.
- Opened containers are stable for 60 days when stored as recommended.
- Dispose of waste product, unused product, and contaminated packaging in compliance with federal, state, and local regulations.

PREPARATION

Coulter DxH Diluent is ready to use. Replace the reagent as directed in your instrument product manuals and/or online help.

IMPORTANT: If product has been partially or completely frozen, allow product to warm to room temperature. Invert the container 16 times to ensure complete mixing prior to placement on the instrument. Install and prime the diluent as directed in your instrument product manuals and/or online help. Verify background counts are acceptable before analyzing patient samples.

PRODUCT AVAILABILITY

COULTER DxH Diluent - 1 x 10 L

628017

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

For additional information or if damaged product is received, call Beckman Coulter Customer Service at 800-526-7694 (USA or Canada) or contact your local Beckman Coulter Representative.

ADDITIONAL INFORMATION

Glossary of Symbols is available at techdocs.beckmancoulter.com (PN C07163)

May be covered by one or more pat. - see www.beckmancoulter.com/patents

Safety Data Sheet is available at techdocs.beckmancoulter.com

Revision History

Revision AM, 06/2017

Changes were made to:

- INTENDED USE
- Added new Bulgarian language

Revision AN, 01/2018

Changes were made to:

- SUMMARY
- PRINCIPLE
- Added new Serbian language
- Added new Vietnamese language

Revision AP, 09/2018

Changes were made to:

- Added new Patent Information Statement
- Added new Dutch, Traditional Chinese, and Ukrainian languages

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Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.