KD MEDICAL GMBH HOSPITAL PRODUCTS®

**Postal Address** PF 060107  
10051 Berlin  
Germany  
**Office** Charlottenstrasse 65  
10117 Berlin  
Germany

**Phone** +49-30-20 39 95 91/92  
**Fax** +49-30-20 39 95 99  
**E-mail** info@kdm-berlin.de  
**http** www.kdm-berlin.de

KD Medical GmbH Hospital Products · Charlottenstrasse 65 · 10117 Berlin · Germany

**Please read instruction for use very carefully and follow usage description Instruction for user of Hypodermic needle for single use**

The Hypodermic needle for single use is manufactured according to ISO 7864 standard.

**[Product name]**

General name: KD-Fine Hypodermic needle for single use

**[Product structure, type and model]**

Product contains Protective cap, needle supporter, needle tube.

Models are identified as color of needle supporter in following sheet:

mm

|                     |               |           |            |            |        |       |            |             |           |            |
|---------------------|---------------|-----------|------------|------------|--------|-------|------------|-------------|-----------|------------|
| Gauge               | 34G           | 33G       | 32G        | 31G        | 30G    | 29G   | 28G        | 27G         | 26G       | 25G        |
| Outer diameter (mm) | 0.18          | 0.2       | 0.23       | 0.25       | 0.30   | 0.33  | 0.36       | 0.4         | 0.45      | 0.5        |
| Color               | Orange        | Black     | Deep green | White      | Yellow | Red   | Blue green | Medium grey | Brown     | Orange     |
| Gauge               | 24G           | 23G       | 22G        | 21G        | 20G    | 19G   | 18G        | 16G         | 15G       | 14G        |
| Outer diameter (mm) | 0.55          | 0.6       | 0.7        | 0.8        | 0.9    | 1.1   | 1.2        | 1.4         | 1.6       | 1.8        |
| Color               | Medium purple | Deep blue | Black      | Deep green | Yellow | Cream | Pink       | White       | Blue-grey | Pale green |

**[Product usage and scope]**

Product is used together with disposable syringe or disposable infusion/transfusion set for medicine injection / infusion or blood collection .

**[Usage method]**

Product is sterile by ETO, no pyrogen, non-toxic and sterile, can be used after opening primary package. Open primary package from peer gap, take out product and fix it onto infusion (transfusion) set, syringe or other Lure conical fittings of administration set. Then, according to clinical nursing science vein infusion operating method, the vein carries out a puncture on the human body.

**[Cautions]**

Product is forbidden to be used in case of primary package broken, protective cap removed, needle tip in abnormal status, exceeding expiry date. Take care when you open primary package to prevent product fall down to floor. Product is for single use only, do not use it repeatedly, Please discard it after usage.

**[Warning]**

Re-use of single-use devices creates a potential risk of patient or user. It may lead to infection and/or impairment of functional capability. Infection and/or limited functionality of the device may lead to injury, illness or death of the patient.

**[Storage and maintenance]**

Product should be stored in clean room without corrosive gas and well ventilated. The relative humidity for product storage place is not more than 80%RH.

**[Identification, label and description]**

On the package, "  " means for single use, "  " means product sterilized by ETO.

**[Validity]**

Product validity is 5 years. Production lot number and expiry date are demonstrated on the package.

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