Ortho File Information

Date Copy Control Approved File:	00
Date Pkg. Engineering Approved File:	NA
Designer Who First Created File:	Colleen Miller
Designer Who Last Modified File:	Kim DePaola
Date File Was Last Modified:	07-22-10
Designer review: (layout and requested changes only)	07-17-10 CR
Dimension on Specification:	400 mm x 640 mm
PMS Colors Used In File:	BLACK
NOTE: ALWAYS RELEASE WITH P107 (first 2-sided page of initial languages)	
07-09-10: KD used approved P1106's, rolled letter tailcode to "I" and made revs per markay from Scott Holdaway for Brexit Project. This will still need some translations.	
07-15-10: KD added languages to Manufacture for section after KD review of file.	
07-16-10: KD used translation spreadsheet sent from LM to finish adding all translations to Sum of Revs and Key to Symbols.	
07-17-10: KD added copyright line to the Sum of Revs per conversation with LM and new markup to come.	
07-22-10 CR signed	

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1. Laboratory biosafety manual. 2nd ed. World Health Organization, Avenue Appia 20, 1211 Geneva 27 Switzerland, 1993.
2. Biotechnology – Laboratories for Research, Development & Analysis – Guidelines for Handling, Inactivating and Testing of Waste. BS EN12740, BSI, 389 Cheswick High Road, London, W4 4AL, 1999.
3. Westhoff CM, Sipherd BD, Toalson LD. Red cell antigen stability in K3EDTA. Immunohematology 1993; 9:109-111.
4. Snyder EL, Hezzey A, Joyner R, Davisson W, Buchholz DH. Stability of red cell antigens during prolonged storage in citrate-phosphate-dextrose and a new preservative solution. Transfusion 1983; 23:165-166.
5. Snyder EL, Hedberg SL, Napychank PA, Roberts C, Kagen L, Aster RA, Quinlan A, Strucaly A, Buchholz DH. Stability of red cell antigens and plasma coagulation factors stored in a non-diethylhexyl phthalate-plasticized container. Transfusion 1993; 33:515-519.

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1. Laboratory biosafety manual. 2nd ed. World Health Organization, Avenue Appia 20, 1211 Geneva 27 Switzerland, 1993.
2. Biotechnology – Laboratories for Research, Development & Analysis – Guidelines for Handling, Inactivating and Testing of Waste. BS EN12740, BSI, 389 Cheswick High Road, London, W4 4AL, 1999.

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RTHO® BRC Set
Blood Bank Reagent Control Set

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Ortho File Information

Date Copy Control Approved File:	00
Date Pkg. Engineering Approval File:	NA
Designer Who First Created File:	Coleen Miller
Designer Who Last Modified File:	Kim DePaola
Date File Was Last Modified:	08-09-19
Designer review: (layoff and requested changes only)	09-12-19 CDR
Dimension on Specification:	400 mm x 600 mm
PMS Colors Used in File:	BLACK & PMS 105

07-15-19 KD used approved P1050, rolled letter tailcode to "e" and made revs per markup from Linda Mano for Breast Project and SDS Hazard Symbol Project. This will still need some translations.

07-16-19 KD added translations that are in Sum of Revs and Key to Symbols. Still awaiting translations for the Hazard stuff.

07-25-19 KD made all "IN Vitro" not oblige to match all other P's.

08-06-19 KD removed changes for SDS Hazard Symbol and left the changes for the Breast Project as these changes will occur separately per email and new markup from LM.

08-09-19 KD added translations to Summary of Revisions in all langs from translation spreadsheet.

08-12-19 CD signed

Anti-Kell (KEL1) Monoclonal IgM (human)

BioClone®

For Slide and Tube Tests

Product Code: 713129

INTENDU USE
Qualitative test with monoclonal antibody for recognition of the Kell antigen on human red blood cells FOR IN VITRO DIAGNOSTIC USE

SUMMARY AND EXPLANATION
Anti-Kell (KEL1) BioClone is used to determine whether red blood cells possess the Kell blood group antigens.

PRINCIPLE OF PROCEDURE
The procedure used with this reagent are based on the principle of agglutination. Normal human red cells, possessing K antigens, will agglutinate in the presence of antibody directed toward the specific antigen.

REAGENTS
Anti-Kell (KEL1) BioClone for Slide and Tube Tests is supplied as a slow-wetain reagent by Ortho-Clinical Diagnostics, Inc. This reagent contains the preservative 0.1% (w/v) sodium azide and is prepared from a human/mouse hybridoma cell line secreting human IgM antibodies.

WARNING
Anti-Kell (KEL1) BioClone (anti-Kell monoclonal antibody) (clone M350) **WARNING: Contains sodium azide. Sodium azide may react with lead and copper plumbing to form highly explosive metal azide. No disposal, flush with a large volume of water to prevent azide buildup.**

STORAGE REQUIREMENTS
Store at 2 to 8 °C. May be at room temperature (15 to 30°C) while in use. Replace cap when not in use.

PRECAUTIONS
1. Handle all blood and materials in contact with blood as capable of transmitting infectious agents. It is recommended that blood and materials in contact with blood be handled using established good laboratory practice.

2. Some components may be considered as hazardous or potentially infectious. Dispose of all materials according to applicable guidelines and regulations.
3. Some components may contain hazardous chemicals (i.e., sodium azide). Refer to the Safety Data Sheet (SDS) available from Ortho-Clinical Diagnostics, Inc. for more information.
4. Do not pipette these reagents by mouth since certain materials that are potentially infectious.
5. Dispose of all waste and materials in accordance with applicable regulations.
6. Do not use reagents beyond their labeled expiration dates.
7. To fully and properly interpret test results, serologic testing is necessary to recognize reagent deterioration.

8. Control slide blood specimens and/or supplementary materials used in the procedures described may interfere with test results.

9. Control sera results may be obtained due to improper technique in performing any diagnostic test.

10. Use of a centrifuge with a lid locking device is recommended.

11. During centrifugation, do not touch the centrifuge rotor or minimum, cause the cells to adhere to the bottom of the test tube but vigorous agitation is necessary before they can be resuspended. During such agitation, weak agglutination or flocculation may be observed.

Undercentrifugation, i.e., failure to use proper force necessary to cause the cells to form a "button" and a clear supernatant may result in a false negative reaction.

SPECIMEN COLLECTION AND PREPARATION
No preparation of the patient/donor is required prior to specimen collection. Blood should be collected by standard medical techniques. Samples should be tested as soon as possible following collection. If a delay in testing occurs, samples should be stored at 2 to 8 °C. Control specimens or blood drawn into serum or EDTA should be tested within 14 days. Donor blood may be tested to the end of the expiration.

Donor blood obtained by finger pricks may be tested directly by the slide method but, to avoid clotting, blood collected in this manner should be tested quickly with the reagent.

REAGENT PREPARATION
No preparation of the reagent is required; use directly from the vial.

PROCEDURE
The procedures identified below are for manual testing only. When using automated equipment, follow the procedures that are contained in the reagent manual.

For tests that require approved validation procedures to demonstrate compatibility of this product on automated systems.

Materials Provided
Anti-Kell (KEL1) Monoclonal IgM (human) BioClone for Slide and Tube Tests

(See Reagent section for product description)

Materials Required But Not Provided

Slide Method

1. Glass slides
2. Pipettes designed to deliver approximately 50 µL
3. Applicator sticks

4. Heated water box qualified for use with blood bank reagents

Tube Method

1. Test tubes, 10 x 75 mm or 12 x 75 mm
2. Pipettes designed to deliver approximately 50 µL
3. Centrifuge

4. Isotonic saline, 0.85-0.9% sodium chloride

Test Procedure

Slide Method

1. Use white blood.

2. On a glass slide, place one drop (approximately 50 µL) of Anti-Kell BioClone reagent.

3. Add one drop (approximately 50 µL) of Anti-Kell reagent to the test tube.

4. With an applicator stick, mix the cell/suspension to the end of the capillary to the test tube.

5. Mix well and incubate at room temperature for a minimum of 15 minutes. If desired, the incubation time may be extended to a maximum of 30 minutes.

6. Centrifuge. Suggested centrifugation: approximately 15 seconds at 3400 rpm (900-1000 rcf) or one minute at 1000 rpm (100-125 rcf).

NOTE: The centrifugal force applied to centrifuged mixtures should be the minimum required to produce a "button" of red cells and a clear supernatant. No more speed or time should be used than is recommended to well cover the wide variety of centrifuge available; each laboratory must calibrate its centrifuge.

7. Read results. If agglutination is observed, the test is positive.

8. Completely resuspend the cells by gentle agitation and examine microscopically for agglutination. Record results.

NOTE: All results should be interpreted immediately following centrifugation and resuspension.

Quality Control Procedures
Quality control is necessary to recognize reagent deterioration. Predicted performance is based on demonstrated reactivity of heterozygous (weak) positive and negative controls performed on each day of use.

Positive Control – Use red blood cells known to possess the antigen toward which the reagent antibody is directed.

Negative Control – Use red blood cells known to lack the antigen toward which the reagent antibody is directed.

INTERPRETATION OF RESULTS
Positive Result (+) – Agglutination of the red blood cells is a positive test result indicating the presence of the corresponding antigen.

Negative Result (-) – No agglutination of the red blood cells is a negative test result and indicates the corresponding antigen is not demonstrable.

LIMITATIONS OF THE PROCEDURE
1. The Test Procedure and Interpretation of Results must be followed closely to assure the accuracy of the test results. Each laboratory should have a program that is based on personnel on the proper use of handling of the product.

2. If the sample is known or observed to contain any interfering substances such as free hemoglobin, lipids or high levels of protein, it is recommended that the red blood cells be washed and resuspended in saline.

3. Agglutinated red cells may yield weaker reactions than those obtained with fresh cells.

EXPECTED RESULTS*
In clinical studies, the blood typing results of donor and patient samples obtained with Anti-Kell (KEL1) Monoclonal IgM (human) BioClone were compared to those obtained with the previous formulation of Anti-Kell BioClone. There were 940 samples tested by the slide method and 546 samples tested by the tube method. Anti-Kell BioClone giving results agreed in 100% of the samples tested.

Percent agreement indicates concordance between the two assays and does not indicate which method gave the correct results.

SPECIFIC PERFORMANCE CHARACTERISTICS*
Anti-Kell (KEL1) Monoclonal IgM (human) BioClone has been tested and found to specifically agglutinate human red cells if the corresponding antigen is present, when used according to the recommended directions for use.

*Data on file at Ortho-Clinical Diagnostics, Inc.

Summary of Revisions

Key to Symbols: Added EC Representative address.

Back Page: Added EC Representative address. Updated copyright information.

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Recapit
Anti-Kell (KEL1) BioClone (anti-Kell monoclonal antibody) (clone M350) **WARNING: Contains sodium azide. Sodium azide may react with lead and copper plumbing to form highly explosive metal azide. No disposal, flush with a large volume of water to prevent azide buildup.**

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1. Laboratory biosafety manual. 2nd ed. World Health Organization, Avenue Appia 20, 1211 Geneva 27 Switzerland, 1993.
2. Biotechnology – Laboratories for Research, Development & Analysis – Guidelines for Handling, Inactivating and Testing of Waste. BS EN12740, BSI, 389 Cheswick High Road, London, W4 4AL, 1999.
3. Westhoff CM, Sipherd BD, Toalson LD. Red cell antigen stability in K3EDTA. Immunohematology 1993; 9:109-111.

4. Snyder EL, Hezzey A, Joyner R, Davisson W, Buchholz DH. Stability of red cell antigens during prolonged storage in citrate-phosphate-dextrose and a new preservative solution. Transfusion 1983; 23:165-166.
5. Snyder EL, Hedberg SL, Napychank PA, Roberts C, Kagen L, Aster RA, Quinlan A, Strucaly A, Buchholz DH. Stability of red cell antigens and plasma coagulation factors stored in a non-diethylhexyl phthalate-plasticized container. Transfusion 1993; 33:515-519.

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Millipore (UK) Ltd.