ENGLISH

INSTRUCTIONS FOR USE

Description of 4BONE

4BONE is a Resorbable Bone Graft Substitute. 4BONE is a microporous and macroporous biphasic calcium phosphate ceramic consisting of 60% Hydroxyapatite (HA) and 40% beta-Tricalcium Phosphate (β-TCP). Hydroxyapatite (HA) - Ca10(PO4)6(OH)2 is chemically comparable to biological apatite crystals. Both in vitro and in vivo, chemically pure HA is considered as the least bioactive of calcium phosphate ceramics. Tricalcium phosphate beta (β-TCP) – Ca3(PO4)2 is recognized to have greater bioactivity, greater extent of dissolution than HA in vitro and in vivo. 4BONE is presented in a porous form required for the biological exchanges particularly for bone ingrowth and mineralization. The balance of a stable component (HA) and a more bioactive and resorbable component (β-TCP) allied to the porosity involves the controlled process of ceramic resorption and bone substitution. When packed into a bony site, 4BONE gradually resorbs and is replaced with bone during the healing process. 4BONE is available in granules packaged in vials or syringes.

Indications

: GF= 'c'e]e[\ [f] kC\ Xc\ X Yfe\ ^ix]e^ dX\ix\ Jf Af Xk^d\ej fi i\ZfejikZ] periodontal or oral/maxillofacial defects. These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. 4BONE can be used with autogenous bone grafting materials.

Typical uses include:

- ✓ Periodontal / Intra-bony defects
- ✓ Extraction sites (implant preparation/placement)
- ✓ Cystic cavities
- ✓ Ridge augmentation
- ✓ Sinus lifts

Contra-indications

4BONE has limited initial mechanical properties. Therefore, this product is contraindicated where the device is intended as structural support. Conditions representing contraindications include also:

- ✓ Osteomyelitis
- ✓ Implantation in necrotic surgical sites
- ✓ Degenerative bone disease
- ✓ Intra-articular implantations
- ✓ Opening meninx
- ✓ As with any surgical procedure, care should be exercised in treating individuals with preexisting conditions that may affect the success of the surgical procedure. This includes (but is not limited to) individuals with long- term steroidal therapy or treatment acting on the calcium or phosphorus metabolism.

Adverse-effects

Possible adverse effects include but are not limited to: Wound complications including haematoma, infection, and other complications that are possible with any surgery.

Preparation

Before use, 4BONE must be hydrated without excess with physiological saline to prevent osmotic damage. If desired, 4BONE may be mixed with patient's whole blood following hydration with physiological saline.

Vial

Remove vial from sterile package using sterile technique. Dispense desired quantity of 4BONE granules into a sterile dish. 4BONE granules should be wetted with physiological saline to facilitate delivery to the surgical site using conventional instruments such as curettes or plastic carriers.

Sterilisation and conservation

✓ 4BONE granules are supplied sterile in vials or syringes

✓ 4BONE is sterilized using beta irradiation. This product should not be re-sterilized. The sterilization process of this product using radiations was validated during its conception. We do not guarantee the sterility of the product if a new sterilisation is performed on it. The non-respect of the single use can conduct to septic issues. Additionally, the resorption process of this bone substitute starts immediately after its implantation on the osseous defect so it can not be re-used

✓ In all operations involving 4BONE granules, appropriate aseptic technique must be observed

✓ 4BONE should be stored at ambient temperature

✓ This device is for single patient use only and should never be reused.

DO NOT USE IF PACKAGE IS OPENED OR DAMAGED OR IF EXPIRATION DATE HAS BEEN EXCEEDED.

Warning And Precautions

K_é g]f [kZ] é]f Y\ _Xe[\ [f i' d'g'xe] \ [Yo j]k'e\ [h'k'À\ [g\ i'c'f'c' _X] e^ i\X [_] é instructions for use

✓ 4BONE is intended only for bony voids or gaps that are not intrinsic to the stability of the bony structure

✓ Maximum contact between the product and the recipient bone must be established

✓ <f ef] f\ i\ iAé : [\ \] Zé

✓ A good immobilization of the 4BONE ceramic is necessary in order to prevent micro-movements which could i\X [Y\ i'f'k' éZ'Xg'c'k'X] fe

✓ Do not leave defect open (e.g. mucosa well sutured and/ or use of appropriate membrane)

✓ The device should be secured to prevent motion and migration, use in areas where the graft can be adequately contained

✓ Do not compromise blood supply to the defect area

✓ In order to avoid destruction of the porous structure, do not compact granules during delivery into the bony gap.

Key to codes used:

	Do not use if package is damaged		Consult instructions for use
	Single use		Batch code
	Use by		Sterilized using irradiation
	Manufacture date (including sterilisation)		Catalogue number
	Manufacturer		Caution: Federal law in USA and Canada restricts this device to sale by or on the order of a physician.

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