Intrapur®

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

















GB Instruction for use of Infusion filters

Intrapur Plus/ Intrapur Lipid/ Sterifix
if higher flow rates are needed
Filtration area: 10 cm²
Filling volume of filter housing: 2.4 ml
Max. Operating pressure of filter housing: 3.1 bar

Intrapur Paed/ Intrapur Paed Lipid/ Sterifix Paed
for medium flow rates but lower dead space
Filtration area: 4.5 cm²
Filling volume of filter housing: 0.7 ml
Max. Operating pressure of filter housing: 3.1 bar

Intrapur Neonat/ Intrapur Neonat Lipid/ Sterifix Neonat
If low flow rates are accepted, but low dead space is required
Filtration area: 1.65 cm²
Filling volume of filter housing: 0.4 ml
Max. Operating pressure of filter housing: 5.2 bar

Product: Infusion filter with Supor® (PES) membrane and position-independent air elimination. For use with pressure infusion equipment acc. to ISO 8536-11 (200 kPa/ 2 bar).

Materials used: see table

Indications: This infusion filter may be used with any infusion set. It is indicated for removal of particulate impurities, microbial contaminations and air bubbles from infusion solutions.

Maximum time of use: see table

Contraindications:
This product should only be used acc. to the instructions.
An infusion filter with 0.2 µm pore size must not be used for giving blood, cellular blood constituents, clotting factors, suspensions or emulsions such as fat or drugs which are not completely dissolved in the carrier solution.
An infusion filter with 1.2 µm pore size must not be used for giving blood, cellular blood constituents, clotting factors. It can be used for suspensions and emulsions such as fat (lipids).
The filter should not be used for sterile filtration of solutions with visible contamination by bacteria or pyrogens.

Warnings:

- Do not use if packaging is damaged. Do not resterilise.
- Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

Precautions:

- The filter helps achieve required standards of hygiene, but it is not a replacement for them. Observe aseptic conditions.
- Drugs may be absorbed by or adsorbed to infusion solution containers, infusion sets and membrane filters. Emergency remedies should be administered as near to the patient as possible.
- Filters reach their maximum effectiveness only when installed close to the patient.
- In lipid emulsions electromechanical destabilization leads to enlarged fat globules which might early occlude the filter.
It is recommended to use the lipid filters (pore size 1.2 µm) in connection with infusion pumps in order to keep a continuous flow rate.

Method of use:

1. Remove the protective caps from the filter inlet and couple to the infusion set (picture 1).
2. Hold the filter below the fluid level (pict.2). Open the roller clamp of the infusion set and fill the tubing and filter completely for deaeration.
3. Close the roller clamp.
4. Ensure that no more air bubbles are present on the patient side (ribbed side) of the filter or in the tubing connector. If air bubbles are still present, open the roller clamp slightly and flick the filter and tubing lightly with the finger until all the air bubbles have come out (pict. 3). Close the roller clamp.
5. Remove the protective cap from the male Luer Lock connector and couple to the vascular access point (pict. 4).
6. Mark the date and time acc. to the hospital specific guidelines for change intervals for infusion filters (pict. 5).

Green Housing

0.2 µm + Filter pore size 0.2 µm, Membrane positively charged, endotoxin retentive

Blue Housing

1.2 µm Filter pore size 1.2 µm, Membrane non-charged, low-protein binding

Colorless Transparent Housing

0.2 µm Filter pore size 0.2 µm, Membrane non-charged, low-protein binding

non-pyrogenic

General information on materials used:
 Liquid Filter Membrane: SUPOR® Polyethersulfone (0.2 µm or 1.2 µm) , hydrophilic
 Air Venting Membrane: PTFE 0.02 µm, hydrophobic
 Filter housing: CYROLITE® (modified polyacrylic)

Application of	via 0.2 µm Filter	via 1.2 µm Filter
aqueous solutions blood colloidal solutions lipids suspensions emulsions infusion solutions of high viscosity (like Glucose 70%)	possible not possible not possible* not possible not possible not possible not possible	possible not possible possible possible depends on particle size possible not possible

*Albumin, starches and gelatine (4, 6, 10 %) can be filtered through 0.2µm membranes

REF / Code No.	Product	Filter characteristics								Flow rate (agua dest.) (ml/ min)	Time of use recommended	Minimum Water Bubble point bar (psi)
		Pore size	Positively Charged	Low protein binding	Bacteria retention (100 %)	Endotoxin retention	Fungi: retention	Particle retention	Air elimination			
4099800	Intrapur Plus	0,2 µm	X	Ø	X	X	X	X	X	> 30	96 h	3.1 (45)
4183916	Intrapur Plus	0,2 µm	X	Ø	X	X	X	X	X	> 30	96 h	3.1 (45)
4099753	Intrapur Paed	0,2 µm	X	Ø	X	X	X	X	X	> 10	96 h	3.1 (45)
4099451	Intrapur Neonat	0,2 µm	X	Ø	X	X	X	X	X	> 2	96 h	3.1 (45)
4184386	Intrapur Neonat	0,2 µm	X	Ø	X	X	X	X	X	> 2	96 h	3.1 (45)
4185226	Intrapur Neonat	0,2 µm	X	Ø	X	X	X	X	X	> 2	96 h	3.1 (45)
4099702	Intrapur Lipid	1.2 µm	Ø	X	Ø	Ø	X	X	X	> 100	24 h	0.3 (5)
4093216	Intrapur Paed Lipid	1.2 µm	Ø	X	Ø	Ø	X	X	X	> 90	24 h	0.3 (5)
4099850	Intrapur Paed Lipid	1.2 µm	Ø	X	Ø	Ø	X	X	X	> 90	24 h	0.3 (5)
4099460	Intrapur Neonat Lipid	1.2 µm	Ø	X	Ø	Ø	X	X	X	> 30	24 h	0.3 (5)
4099303	Sterifix	0,2 µm	Ø	X	X	Ø	X	X	X	> 30	24 h	3.1 (45)
4184637	Sterifix	0,2 µm	Ø	X	X	Ø	X	X	X	> 30	24 h	3.1 (45)
4099354	Sterifix Paed	0,2 µm	Ø	X	X	Ø	X	X	X	> 10	24 h	3.1 (45)
4099257	Sterifix Neonat	0,2 µm	Ø	X	X	Ø	X	X	X	> 2	24 h	3.1 (45)

REF / Code No.	Product	Materials used							
		MBS	PTFE	SUPOR® (PES)	PC	PE	PVC (DEHP-free)	MABS	PUR
4099800	Intrapur Plus	X	X	X	X	X	X	X	
4183916	Intrapur Plus	X	X	X	X	X		X	X
4099753	Intrapur Paed	X	X	X		X	X		
4099451	Intrapur Neonat	X	X	X	X	X		X	X
4184386	Intrapur Neonat	X	X	X	X	X		X	X
4185226	Intrapur Neonat	X	X	X	X	X		X	X
4099702	Intrapur Lipid	X	X	X	X	X	X	X	
4093216	Intrapur Paed Lipid	X	X	X		X			
4099850	Intrapur Paed Lipid	X	X	X		X	X		
4099460	Intrapur Neonat Lipid	X	X	X	X	X		X	X
4099303	Sterifix	X	X	X	X	X	X	X	
4184637	Sterifix	X	X	X	X	X		X	X
4099354	Sterifix Paed	X	X	X		X			
4099257	Sterifix Neonat	X	X	X	X	X		X	X

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B. Braun Melsungen AG,