ULTIMUM™ EV HEMOSTASIS INTRODUCER

English – Instructions for Use	2
العربية – إرشادات الاستخدام	4
Български – Инструкции за употреба	6
Česky – Návod k použití	8
Dansk – Brugsanvisning	10
Deutsch – Gebrauchsanleitung	12
Ελληνικά – Οδηγίες χρήσης	14
Español – instrucciones de uso	16
Eesti – kasutusjuhend	18
Suomi – Käyttöohjeet	20
Français – Mode d'emploi	22
Hrvatski – upute za upotrebu	24
Magyar – Használati útmutató	26
Italiano – Istruzioni per l'uso	28
한국어 – 사용 설명서	30
Lietuviškai – Naudojimo instrukcijos	32
Latviski – lietošanas pamācība	34
Македонски – Упатство за употреба	36
Nederlands – Gebruiksaanwijzing	38
Norsk – bruksanvisning	40
Polski – instrukcja obsługi	42
Português – Instruções de utilização	44
Română – Instrucțiuni de utilizare	46
Русский – инструкция по эксплуатации	48
Slovenčina – Pokyny na používanie	50
Srpski – uputstvo za upotrebu	52
Svenska – Bruksanvisning	54
Türkçe – Kullanım Talimatları	56
Інструкція з експлуатації – українською мовою	58
中文 – 使用说明	60
SYMBOL	62

ULTIMUM™ EV

HEMOSTASIS INTRODUCER

English - Instructions for Use

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THE ABBOTT MEDICAL HEMOSTASIS INTRODUCER IS DESIGNED TO INTERLOCK ONLY WITH ABBOTT MEDICAL DILATORS. MISUSE MAY RESULT IN SERIOUS COMPLICATIONS.

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DEVICE DESCRIPTION

The Ultimium™ EV Hemostasis Introducer consists of a polyethylene introducer and dilator with a radiopaque marker band at the distal end. The inner diameter is identified on the package label. The introducer is fitted with a hemostasis valve for use during catheter introduction and/or exchange over a guidewire. The hub is located on the proximal end of the introducer, where a sideport with a three-way stopcock is provided for aspiration, fluid infusion, blood sampling, and pressure monitoring. The dilator is tapered at the distal tip with an internal lumen designed to accept guidewires that have a maximum diameter of 0.035" (0.889 mm).

The Ultimium™ EV Hemostasis Introducer is packaged in a tray and pouch, then sterilized using ethylene oxide. The Instructions for Use are recyclable. Dispose of used product and packaging in accordance with hospital procedures.

INDICATIONS FOR USE

The Ultimium™ EV Hemostasis Introducers are designed for the introduction of angiographic catheters, closed end catheters, balloon catheters, and electrodes into a vessel where minimizing blood loss is essential.

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Damage to the valve assembly may occur if inner catheter is inserted or withdrawn rapidly.

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STORE IN A COOL, DARK, DRY PLACE.

For Single Use Only. Single-use devices are designed and tested for only one patient application. These are disposable devices and are not designed for reprocessing and reuse. Reuse of designated "single-use" devices creates a risk of patient or user infections due to prior patient use and the difficulty in cleaning the narrow structures at material interfaces following direct blood contact. Contamination or reprocessing cleaning agent residues may lead to adverse patient reactions and may damage the device. Use of non Abbott Medical packaging may compromise device functionality and sterility due to compromised protection from shipping and handling damage. The absence of labeling after reprocessing, may lead to misuse of the device and impaired traceability. Reprocessing and reuse may result in patient or user injury, permanent impairment or death.

SUGGESTED DIRECTIONS FOR USE

INSPECT SHEATH INTRODUCER AND ACCESSORIES FOR DAMAGE, DO NOT USE ANY DAMAGED DEVICES.

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1. Insert dilator into hemostasis valve. Follow normal accepted practice for vessel puncture, guidewire insertion, vessel dilator, and hemostasis introducer use.

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3. Follow manufacturer's recommendations for the catheter or device being introduced via the hemostasis introducer.

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™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third-party trademark, which is property of its respective owner.

Pat. <http://www.abbott.com/patents>

© 2020 Abbott. All Rights Reserved.

600108977

CE
2797

Abbott

ULTIMUM™, ULTIMUM™ EV HEMOSTASIS INTRODUCER

English - Instructions for use	3
Deutsch - Gebrauchsanleitung	5
Français - Mode d'emploi	7
Italiano - Istruzioni per l'uso.	9
Español - Instrucciones de uso	11
中文 - 使用说明	13
Ελληνικά - Οδηγίες χρήσης	14
Русский - Инструкция по применению.	16
Magyar - Használati utasítás	18
Português - Instruções de Utilização	20
Nederlands - Gebruiksaanwijzing	22
Suomi - Käyttöohjeet	24
Svenska - Bruksanvisning	25
Dansk - Brugsanvisning	27
Norsk - Bruksanvisning	29
Türkçe - Kullanım talimatları	30
한국어 - 사용 지침	32
Polski - Instrukcja użytkowania	33
Srpski - Uputstvo za upotrebu	35
SYMBOL	37

REF

407637	407653	407687	407830	407840	407850	407860
407638	407654	407688	407831	407841	407851	407861
407639	407655	407689	407832	407842	407852	C407594
407640	407657	407690	407833	407843	407853	C407595
407641	407658	407691	407834	407844	407854	C407596
407647	407659	407692	407835	407845	407855	C407597
407648	407660	407693	407836	407846	407856	C407598
407649	407661	407694	407837	407847	407857	C407599
407650	407662	407695	407838	407848	407858	C407698
407651	407663	407696	407839	407849	407859	C407699
407652	407664	407697				

ULTIMUM™, ULTIMUM™ EV

HEMOSTASIS INTRODUCER

Instructions for use

See Product Catalog for Reorder Numbers.

THE FRENCH SIZE SPECIFIED REPRESENTS THE INNER DIAMETER OF THE INTRODUCER SHEATH.

DO NOT ATTEMPT TO INSERT A CATHETER HAVING A DISTAL TIP OR BODY SIZE LARGER THAN THE INTRODUCER SIZE INDICATED.

THE ABBOTT MEDICAL HEMOSTASIS INTRODUCER IS DESIGNED TO INTERLOCK ONLY WITH ABBOTT MEDICAL DILATORS AND/OR OBTURATORS. MISUSE MAY RESULT IN SERIOUS COMPLICATIONS.

SEE INDIVIDUAL STERILE PACKAGE LABEL FOR CONTENTS.

SINGLE-USE DISPOSABLE MEDICAL DEVICE. CONTENTS ARE STERILE IF PACKAGE IS UNOPENED AND UNDAMAGED. DO NOT RESTERILIZE.

The instructions for use is recyclable. Dispose of used product and packaging following standard solid biohazard waste procedures

CAUTION

DO NOT alter this device.

Federal law (U.S.) restricts this device to sale by or on the order of a physician.

Do not attempt to use a guidewire over maximum diameter specified on package label.

This device should only be used by physicians thoroughly trained in the technique of angiography and/or the use of Abbott Medical catheter delivery systems.

Individual patient anatomy and physician technique may require procedural variations.

The closing force of the hemostasis valve may alter or impair the function of some catheters.

Damage to guidewire may result if withdrawn through a metal needle cannula. Cannula must be removed first.

Do not attempt to advance or withdraw guidewire if resistance is felt. Use fluoroscopy to determine cause.

Damage to the valve assembly may occur if inner catheter is withdrawn rapidly.

Advance dilator/sheath assembly with a twisting motion to avoid damage to the sheath or vessel.

STORE IN A COOL, DARK, DRY PLACE.

For Single Use Only! Single-use devices are designed and tested for only one patient application. These are disposable devices and are not designed for reprocessing and reuse. Reuse of designated "single-use" devices creates a risk of patient or user infections due to prior patient use and the difficulty in cleaning the narrow structures at material interfaces following direct blood contact. Contamination or reprocessing cleaning agent residues may lead to adverse patient reactions and may damage the device. Use of non Abbott Medical packaging may compromise device functionality and sterility due to compromised protection from shipping and handling damage. The absence of labeling after reprocessing, may lead to misuse of the device and impaired traceability. Reprocessing and reuse may result in patient or user injury, permanent impairment or death.

APPLICATION

The introduction of angiographic catheters, closed end catheters, balloon catheters, and electrodes into a vessel where minimizing blood loss is essential.

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SUGGESTED DIRECTIONS FOR USE

INSPECT SHEATH INTRODUCER AND ACCESSORIES FOR DAMAGE. DO NOT USE ANY DAMAGED DEVICES.

BEFORE PROCEEDING, VERIFY THE HEMOSTASIS INTRODUCER IS OF PROPER SIZE TO ALLOW PASSAGE OF DESIRED CATHETER OR OTHER MEDICAL DEVICE.

FLUSH THE DILATOR(S), INTRODUCER SHEATH AND INFUSION SIDEPORT AND STOPCOCK WITH HEPARINIZED NORMAL SALINE.

1. Insert dilator into hemostasis valve. Follow normal accepted practice for vessel puncture, guidewire insertion, vessel dilator, and hemostasis introducer use. For introducers with shaped distal tips, always advance dilator/sheath assembly over appropriate sized guidewire to the appropriate anatomical site.

CAUTION: DO NOT ATTEMPT ADVANCEMENT WITHOUT GUIDEWIRE. SEVERE VASCULAR DAMAGE AND/OR INJURY MAY OCCUR.

CAUTION: INSERTION INTO ARTERY MAY CAUSE EXCESSIVE BLEEDING AND/OR OTHER COMPLICATIONS.

CAUTION: DO NOT ALLOW GUIDEWIRE TO INADVERTENTLY ADVANCE TOTALLY INTO PATIENT.

2. When a sideport introducer is used, follow normal practice of using a continuous drip of anticoagulant fluid under pressure through the sideport when the hemostasis introducer is in the vessel.

CAUTION: UPON WITHDRAWAL OF DILATOR, CATHETER OR OTHER INSERTED DEVICE, ASPIRATION THROUGH THE STOPCOCK IS RECOMMENDED TO REMOVE ANY FIBRIN ACCUMULATION.

3. An Abbott Medical repositioning sleeve may be used with the hemostasis introducer to provide an added measure of environmental protection to catheters (see Abbott Medical product catalog).
4. Follow manufacturer's recommendations for the catheter or device being introduced via the hemostasis introducer.
5. When using the Abbott Medical luer-lock hemostasis introducer, an Abbott Medical luer-lock obturator may be used to seal off the hemostasis aperture and aid in maintaining sheath patency and general management of the introduction site.

OVERFLUSHING

The French size shown on each Abbott Medical introducer represents the tip inside dimension.

To overflush catheters or Abbott Medical luer-lock obturators placed into an Abbott Medical hemostasis introducer, use introducer that is one (1) French size larger than the French size of the catheter or obturator.

CAUTION: UPON REMOVAL OF HEMOSTASIS INTRODUCER, PROPER PRECAUTIONS SHOULD BE TAKEN TO PREVENT BLEEDING.

INTRODUCER WITH SHAPED DISTAL TIPS

Reinsert the dilator fully into the introducer sheath to aid in straightening the top portion. Then remove the dilator and introducer sheath as a unit.

USA ONLY

LIMITED WARRANTY AND DISCLAIMER

Abbott Medical hereby warrants that if any Abbott Medical product fails to perform within normal tolerances for a patient due to a defect in materials or workmanship, Abbott Medical will provide, at no charge, a replacement Abbott Medical product for the patient's use. This limited warranty applies only if each of the following conditions are met.

1. The product was packaged and labeled by Abbott Medical.
2. The failed product must be returned to Abbott Medical and becomes the property of Abbott Medical.
3. The product has not been mishandled, reprocessed or altered in any way.
4. The product was used before the "USE BY" date marked on the packaging of the product.

No representation or warranty is made that an Abbott Medical product will not fail. Abbott Medical disclaims responsibility for any medical complications, including death, resulting from the use of its products. Except as expressly provided by this limited warranty, ABBOTT MEDICAL IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES BASED ON ANY DEFECT, FAILURE, OR MALFUNCTION OF ITS PRODUCTS, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT, OR OTHERWISE. Some states do not allow the exclusion or limitation of incidental or consequential damages however, so the above limitation or exclusion may not apply to you.

Except as expressly provided by the limited warranty, ABBOTT MEDICAL MAKES NO WARRANTY, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE.

STATE OF CALIFORNIA (USA) ONLY

WARNING: This product can expose you to chemicals including ethylene oxide, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

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