AQT90 FLEX

Empty Tube Kit

RADIOMETER 

Content

en	Empty Tube Kit	4
bg	Набор празни епруветки	5
cs	Souprava prázdných zkumavek	6
da	Empty Tube Kit	7
de	Leerröhrchen-Kit	8
el	Πακέτο Κενών Φιαλιδίων	9
es	Kit de Tubos Vacíos para QC	10
et	Tühjade katsutite komplekt	11
fi	Putkipakkaus	12
fr	Kit Empty Tube	13
hr	Komplet praznih epruveta	14
hu	Üres csőkészlet	15
it	Kit provette vuote	16
ja	エンプティチューブキット	17
lt	Tuščių vamzdelių rinkinys	18
lv	Tukšu cauruļu komplekts	19
nl	Kit met lege buisjes	20
no	Empty Tube Kit	21
pl	Empty Tube Kit (Zestaw pustych probówek)	22
pt	Kit de Tubos Vazios	23
ro	Kit de tuburi goale	24
ru	Набор пустых пробирок	25
sk	Súprava prázdných skúmaviek	26
sr	Komplet praznih epruveta	27
sv	Förpackning med tomma rör	28
tr	Boş test kiti	29
zh	质控管	30

en

Empty Tube Kit

Proprietary name and common name: Empty Tube Kit

Catalog no. 944-230

Intended use

For *in vitro* diagnostic use.

To be used for liquid quality controls on the AQT90 FLEX analyzer.

For single-use only.

NOTICE: Do not use Empty Tubes to store samples.

Liquid quality control (LQC) solution

During tests to verify quality control on the analyzer using the Empty Tube as a container for LQC solutions, serum-based Bio-Rad's Liquichek Cardiac Markers Plus Control LT Levels 1, 2 and 3 were used. If other LQC solutions are used new verification tests must be performed.

- Contents:**
- 50 Empty Tubes
- Components:**
- Outer plastic tube
 - Inner glass tube
 - Rubber septum
 - Opal plastic cap
 - Label

Preparing a tube for use in a LQC test procedure

NOTICE: The procedure is illustrated inside the lid of the Empty Tube Kit box.

1. First remove the opal cap from the outer tube and then remove the rubber septum from the inner tube.
2. Use a pipette to transfer the relevant volume of LQC solution into the inner tube - the volume depends on the number of parameters you wish to analyze:
 - One parameter: 400 µL
 - Two parameters: 450 µL
 - Three parameters: 500 µL
 - Four parameters: 550 µL
 - Five parameters: 600 µL
3. Place the rubber septum over the inner tube.
4. Press it into place.

5. Place the opal cap over the outer tube and press it down until it closes with a 'click'.
The tube is now ready for use in a LQC test procedure on the AQT90 FLEX analyzer. See the AQT90 FLEX operator's manual for details.

WARNING – Risk of infection

Disposal of all waste should be in accordance with local regulations. Refer to recommendations of the manufacturer for disposal of liquid quality control (LQC) material.

Performance characteristics

The rubber septum in the opal cap can withstand a maximum of five needle-penetrations.

Explanations of graphical symbols:

	Mark of compliance with the applicable EU Directives
--	--

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Radiometer, the Radiometer logo, ABL, AQT, TCM, RADIANCE, AQUIRE, PICO, CLINITUBES and QUALICHECK are trademarks of or used under license by Radiometer Medical ApS.

© 2018 Radiometer Medical ApS, 2700 Brønshøj, Denmark.
All rights reserved. 994-079. 201802N.

RADIOMETER 

994-078. 201301D

AQT90 FLEX

Cleaning Solution Tubes

en

[REF 905-843**]**

Intended use

For *in vitro* diagnostic use.

For cleaning the sample-handling part of the liquid-transport system of the AQT90 FLEX analyzer.

Preparation: The tubes are ready for use and must not be opened.

Procedure: See AQT90 FLEX documentation.

Chemical composition: Aqueous solution of sodium hypochlorite < 5 % (weight/volume).

pH: approximately 12.

<div></div>	<div>CAUTION — <i>Risk of toxic gas production</i><div>Do not mix the contents of the Cleaning Solution Tubes with any other product. Failure to follow this advice may result in the production of toxic chlorine gas.</div></div>
<div></div>	<div>CAUTION — <i>Risk of corrosive effects</i><div>Avoid contact with skin and eyes. The solution may be corrosive to skin, eyes and mucous membranes.</div></div>

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