SMARTABLATE™ Irrigation Tubing Set

Instructions for Use SMARTABLATE™ Irrigation Tubing Set	English	page	7
Mode d'emploi Jeu de tubulures d'irrigation SMARTABLATE™	Français	page	9
Gebruiksaanweisung SMARTABLATE™ Spülschlauchsatz	Deutsch	seite	11
Istruzioni per l'uso Set di tubi di irrigazione SMARTABLATE™	Italiano	pagina	13
Instrucciones de uso Juego de tubos de irrigación SMARTABLATE™	Español	página	15
Instruções de Utilização Conjunto de Tubagem de Irrigação SMARTABLATE™	Português	página	17
Gebruiksaanwijzing SMARTABLATE™-irrigatieslangenset	Nederlands	pagina	19
Brugsanvisning SMARTABLATE™ slangesæt til skyllning	Dansk	side	21
Käyttöohjeet SMARTABLATE™ -huuhteluletkusto	Suomi	sivu	23
Bruksanvisning SMARTABLATE™ spolningsslangsets	Svensk	sida	25
Bruksanvisning SMARTABLATE™ irrigasjonsslangesett	Norsk	side	27
Οδηγίες χρήσης Σετ σωλήνωσης καταιονισμού SMARTABLATE™	Ελληνικά	σελίδα	29
Kullanım Talimatları SMARTABLATE™ İrigasyon Tüp Seti	Türkçe	sayfa	31
Инструкция по применению Комплект трубок для промывания SMARTABLATE™	Русский	стр.	33
Instrukcja obsługi Zestaw przewodów irygacyjnych SMARTABLATE™	Polski	strona	35
Návod k použití Sada proplachových hadiček SMARTABLATE™	Česky	strana	37
Használati utasítás SMARTABLATE™ öblítő csőkészlet	Magyar	oldal	39
Návod na používanie Súprava preplachových hadičiek SMARTABLATE™	Slovensky	strana	41
Navodila za uporabo Garnitura izpiralnih cevč SMARTABLATE™	Slovensko	stran	43
Инструкции за употреба Комплект тръби за иригация SMARTABLATE™	Български	страница	45

Instrucțiuni de utilizare Set de tubulatură de irigare SMARTABlate™	Română	pagina	47
Kasutusjuhend SMARTABlate™-i irrigatsioonivoolikute komplekt	Eesti	lehekülg	49
Lietošanas norādījumi SMARTABlate™ irigācijas cauruļsistēmas komplekts	Latviski	lappuse	51
Naudojimo instrukcija SMARTABlate™ irigacijos vamzdelių komplektas	Lietuvių k.	psl.	53
使用说明 SMARTABlate™ 灌注管路系统	简体中文	页码	55
사용 설명서 SMARTABlate™ 세척 튜브 세트	한국어	페이지	57
使用說明 SMARTABlate™ 灌注套管組	繁體中文	頁碼	59
Instruksi Penggunaan Perangkat Selang Irigasi SMARTABlate™	Bahasa Indonesia	halaman	61
Upute za uporabu Komplet cijevi za ispiranje SMARTABlate™	Hrvatski	strana	63
Упутство за употребу SMARTABlate™ Комплет цевчица за иригацију	Српски	страна	65
Упатство за употреба Комплет цевки за иригација SMARTABlate™	Македонски јазик	страница	67

Symbol Definitions / Définitions des symboles / Symboldefinitionen / Definizioni dei simboli / Definiciones de símbolos / Definições dos símbolos / Definitie van symbolen / Symboldefinitioner / Symbolien selitykset / Symbolförklaringar / Symboldefinisijoner / Ορισμοί συμβόλων / Sembol Tanımları / Определение символов / Definicje symboli / Definine symbolů / A szimbólumok magyarázata / Definicie symbolov / Definicije symbolov / Определения на символите / Definiții simbolurilor / Sümböller tähendused / Simbolu definicijas / Simbolių apibrėžys / 标识定义 / 기호 - 정의 / 符號定義 / Definisí Simbol / Opis simbola / Дефиниције симбола / Дефиниции за симболите

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www.necavcert.ru | info@necavcert.ru

M-5276-663.02A

WARNING: This is a controlled proprietary and confidential document. Verify revision is current prior to use.

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
SMARTABLATE™ Irrigation Tubing Set

- Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
- STERILE. Sterilized using irradiation.
- For single use only.
- Do not resterilize.
- Do not use if package is opened or damaged.

DESCRIPTION

The SMARTABLATE™ Irrigation Tubing Set consists of a drip chamber (I) with I.V. spike (II) for connection to an I.V. bag, a tubing section (III) designed for use with an irrigation pump, a standard male luer lock (IV), and an attached 4-way stopcock (V) for connection to an irrigated catheter.

Figure 1 - SMARTABLATE™ Irrigation Tubing Set

INTENDED USE

The SMARTABLATE™ Irrigation Tubing Set is designed for use with the SMARTABLATE™ Irrigation Pump or the nGEN™ Irrigation Pump to deliver irrigation solution at specified flow rates to the CELSIUS™ THERMOCOOL™ Catheter, NAVISTAR™ THERMOCOOL™ Catheter, or other irrigated devices for cooling purposes.

CONTRAINDICATIONS

There are no known contraindications for the SMARTABLATE™ Irrigation Tubing Set.

PRECAUTIONS

- Store in a cool, dry place.
- Inspect the sterile packaging and tubing set prior to use. Do not use if opened or damaged.
- The tubing set has been tested for up to 8 hours of use. If a procedure will extend longer than 8 hours, remove and replace the tubing set at or before 8 hours of use.

WARNINGS

- Do not expose the tubing set to organic solvents.
- The tubing set is designed for use with the SMARTABLATE™ Irrigation Pump or the nGEN™ Irrigation Pump. Refer to the User Manual for the pump for additional precautions and warnings.
- Refer to the instructions for use provided with the CELSIUS™ THERMOCOOL™ Catheter, NAVISTAR™ THERMOCOOL™ Catheter, or other irrigated catheter for recommended procedure parameters.
- This device is packaged and sterilized for single use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness, or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

“USE BY” DATE

Use the device prior to the “Use By” date on the package label.

DISPOSAL

Recycle components, or dispose of the product and its residual elements or waste items in accordance with local laws and regulations.

DIRECTIONS FOR USE

1. Remove the sterile tubing set from its package using standard safe hospital practices. Ensure that the 4-way stopcock (V) is securely connected to the male luer lock (IV). Position the 4-way stopcock (V) handle so that the fluid pathway of the tubing set is “Off”.
2. Remove the I.V. spike cap and luer caps from the 4-way stopcock (V). Insert the I.V. spike (II) into an I.V. solution container. Hang the I.V. solution container with the drip chamber (I) in a vertical position near the pump, and fill the drip chamber (I) about 2/3 full.

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3. Prime the tubing set by rotating the handle of the 4-way stopcock (V) so that the fluid pathway of the tubing set is open and completely fill the tubing set with irrigation fluid. Remove all trapped air and then close the 4-way stopcock (V) by rotating the handle so that the fluid pathway of the tubing set is "Off".
4. Install the tubing set in the pump as described in the User Manual for the pump.
5. To prepare for irrigation, open the 4-way stopcock (V) by rotating the handle such that the side port is "Off".
6. Flush the tubing set as described in the User Manual for the pump. Verify that the tubing set is free of leaks or defects and that all bubbles have been flushed from the tubing set. If bubbles are present, continue flushing until they are passed through the tubing set.
7. Securely connect the male luer on the 4-way stopcock (V) or the male luer lock (IV) on the end of the tubing to an irrigated catheter. Use the pump to flush all air from the catheter before inserting the catheter into the patient.

DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE PRODUCT(S) DESCRIBED HEREIN. UNDER NO CIRCUMSTANCES SHALL BIOSENSE WEBSTER, INC., OR ITS AFFILIATED COMPANIES, BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY APPLICABLE LAW.

WITHOUT LIMITING THE FOREGOING, BIOSENSE WEBSTER, INC. OR ITS AFFILIATED COMPANIES, SHALL NOT BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES, ARISING OUT OF THE REUSE OF ANY PRODUCT(S) LABELED FOR SINGLE USE OR WHERE REUSE IS PROHIBITED BY APPLICABLE LAW.

Descriptions and specifications appearing in Biosense Webster, Inc. printed matter, including this publication, are informational only and meant solely to generally describe the product(s) at the time of manufacture and are not made or given as a warranty of

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BIOSSENSE WEBSTER, the Biosense Webster logo, CELSIUS, NAVISTAR, nGEN, SMARTABULATE, and THERMOCOOL are trademarks of Biosense Webster, Inc.
The third-party trademarks used herein are trademarks of their respective owners.