REFSTAR™ with QWIKPATCH™ External Reference Patch REFSTAR™ PLUS with QWIKPATCH™ External Reference Patch

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REFSTAR™ with QWIKPATCH™ External Reference Patch

REFSTAR™ PLUS with QWIKPATCH™ External Reference Patch

- **NON-STERILE.**
- **For single use only.**
- **Do not use if package is opened or damaged.**

DEVICE DESCRIPTION

The Biosense Webster REFSTAR™ with QWIKPATCH™ or REFSTAR™ PLUS with QWIKPATCH™ External Reference Patch is an integral part of a non-fluoroscopic catheter tip location and electrogram capture technology called CARTO™ and NOGA™. When used with the CARTO™ and NOGA™ Systems, the location of the navigation catheter tip is compared to the location of this reference patch. This reference device consists of a sensor embedded in an adhesive patch, which is connected to a handle or connector.

CONSULT THE LOCAL DISTRIBUTOR OR MANUFACTURER FOR THE APPROPRIATE INTERFACE OR ADAPTOR CABLES.

INDICATIONS

The Biosense Webster REFSTAR™ with QWIKPATCH™ or REFSTAR™ PLUS with QWIKPATCH™ External Reference Patch is indicated for use with Biosense Webster navigation catheters and the CARTO™ and NOGA™ Systems to provide a reference point for catheter tip location.

CONTRAINDICATIONS

There are no known contraindications for the use of the external reference patch.

PATIENT TARGET GROUP

This device is targeted at patients who have been diagnosed with cardiac arrhythmias and are undergoing an electrophysiology procedure.

INTENDED USERS

The intended users of this device are appropriately trained personnel in a fully equipped electrophysiological laboratory.

CLINICAL BENEFIT

The external reference patch provides a means for attaching the patch sensor cable to the patient in order to convey catheter location.

WARNINGS

- Do not immerse the proximal handle or cable connector in fluids; electrical performance could be affected.
- Do not expose the external reference patch to organic solvents such as alcohol.
- Do not autoclave the external reference patch.

PRECAUTIONS

- Do not attempt to use the external reference patch prior to completely reading and understanding these Instructions for Use.
- This device is for external use only.

- Refer to the product label for recommended storage conditions. Inspect the packaging and external reference patch prior to use.
- The external reference patch is intended for single patient use only.
- Do not resterilize and reuse.
- Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare practitioner.
- This device is packaged for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

“USE BY” DATE

Use the device prior to the “Use By” date on the package label.

REPORTING INCIDENTS

Per regulation 2017/745/EU on medical devices, if this device is used in the European Union (or in a country with an identical regime) and a serious incident occurs during the use of this device or as a result of use of this device, report the incident to the manufacturer and/or its authorized representative and to your national authority.


DISPOSAL

Recycle components, or dispose of the product and its residual elements or waste items in accordance with local laws and regulations.

SUGGESTED INSTRUCTIONS FOR USE

Please refer to these Instructions for Use: the Biosense Webster ESOPHASTAR™, NAVISTAR™ or NOGASTAR™ Catheter's Instructions for Use, and the CARTO™ or NOGA™ System's User Manual when using the REFSTAR™ with QWIKPATCH™ or REFSTAR™ PLUS with QWIKPATCH™ External Reference Patch in conjunction with the ESOPHASTAR™, NAVISTAR™ or NOGASTAR™ Catheters and the CARTO™ or NOGA™ System.

- Remove the external reference patch from its package. Affix the patch to the patient's back (as shown below). Connect the proximal handle or connector of the patch to the appropriate equipment.
- The external reference patch should be securely positioned on the patient's back. The patch should be placed below the scapula in the area of the seventh (for CARTO™) or the eighth (for NOGA™) thoracic vertebra. The body (shaft) of the reference device should follow the long axis of the patient towards the feet or head. Both the patch and the body of the reference device should be securely positioned and taped in such a manner as not to subject the device to any abnormal stresses during the procedure which may dislodge the

reference device from its location. The handle or connector of the reference device should be connected to the appropriate connection scheme leading to the CARTO™ or NOGA™ System patient interface unit. The CARTO™ and NOGA™ System software program is equipped with a positional reference location feature. This feature is used to ensure that the reference device's tip is within the boundaries of the mapping field.

If there are any questions regarding the use or performance of this product, please consult with the local distributor or the manufacturer.

ADVERSE REACTIONS

A number of adverse reactions have been documented for electrophysiology procedures including:

heart block, pulmonary vein stenosis, esophageal fistula and/or injury, stroke (cerebrovascular accident), other arrhythmias (outside diagnosis), life threatening arrhythmias, myocardial infarction, cardiac perforation, pericardial effusion, cardiac tamponade, thrombosis, embolism, pulmonary embolism, air embolism, valvular damage, phrenic nerve injury, vagal nerve injury, pericarditis, coronary artery stenosis, vessel perforation (peripheral and/or central), soft tissue injury, persistent atrial communication, device related infection, embolization of components, device entrapment, surgical intervention (additional), skin burns, and localized skin reaction.

The following complications associated with cardiac catheterization have also been reported in the literature:

major bleed, hematoma, reaction to medications, allergic reaction, vascular access complication, damage to vasculature, implanted device interactions, renal artery stenosis, pneumothorax, ST segment changes, fluid overload, urinary catheter complications, hypotension, sepsis, wound infection, respiratory failure, heart injury, renal injury, heart failure, cardiac arrest and death.

DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE PRODUCT(S) DESCRIBED HEREIN. UNDER NO CIRCUMSTANCES SHALL BIOSENSE WEBSTER, INC., OR ITS AFFILIATED COMPANIES, BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY APPLICABLE LAW.

WITHOUT LIMITING THE FOREGOING, BIOSENSE WEBSTER, INC. OR ITS AFFILIATED COMPANIES, SHALL NOT BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES, ARISING OUT OF THE REUSE OF ANY PRODUCT(S) LABELED FOR SINGLE USE OR WHERE REUSE IS PROHIBITED BY APPLICABLE LAW.

Descriptions and specifications appearing in Biosense Webster, Inc. printed matter, including this publication, are informational only and meant solely to generally describe the product(s) at the time of manufacture and are not made or given as a warranty of the prescribed product(s) in any way.

ELECTRONIC INSTRUCTIONS FOR USE

This document is available at www.e-ifu.com.

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BIOSENSE WEBSTER, the Biosense Webster logo, CARTO, ESOPHASTAR, NAVISTAR, NOGA, NOGASTAR, QWIKPATCH and REFSTAR are trademarks of Biosense Webster, Inc.
The third-party trademarks used herein are trademarks of their respective owners.