STORZ

KARL STORZ — ENDOSKOPE

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@ksear.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Carneros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 México D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
127 39 Skärholmen, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kęstucio st. 59 / Lenktoji st. 27
08124 Vilnius, Lithuania
Phone: +370 5 272 0448
Mobile: +370 685 67 000
E-Mail: info-lt-lv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopia Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopia France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Adria Eos d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Cesta v Gorice 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: + 386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongquiao Road, Xuhui District,
200030, Shanghai, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

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
- ① Screw base for KARL STORZ fiberoptic light cables and Olympus/Winter & Ibe
- ② Plug-in base, 9 mm diameter, for Wolf fiber optic light cables
- ③ Light inlet piece without adaptor
- ④ Color-coded ring to indicate the telescope's direction of view
- ⑤ Eyepiece
- ⑥ Sheath
- ⑦ Type designation
- ⑧ 'AUTOCLAVE' marking ring
- ⑨ Lens and distal outlet of light cable
- ⑩ ICG marking

The telescopes named are used within the scope of fluorescence imaging with the NIR/ICG system

1 Intended use

The KARL STORZ ICG Imaging System is intended to provide real-time visible (VIS) and near-infrared (NIR) fluorescence imaging.

2 Indications

Upon intravenous administration and use of ICG consistent with its approved label, the KARL STORZ Endoscopic ICG System enables surgeons to perform minimally invasive surgery using standard endoscopic visible light as well as visual assessment of vessels, blood flow and related tissue perfusion, or at least one of the major extra-hepatic bile ducts (cystic duct, common bile duct and common hepatic duct), using near infrared imaging. Fluorescence imaging of biliary ducts with the KARL STORZ Endoscopic ICG System is intended for use with standard of care white light and, when indicated, intraoperative cholangiography. The device is not intended for standalone use for biliary duct visualization.

Upon interstitial administration and use of ICG consistent with its approved label, the KARL STORZ Endoscopic ICG System is used to perform intraoperative fluorescence imaging and visualization of the lymphatic system, including lymphatic vessels and lymph nodes.

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
3 Contraindications

Use of the KARL STORZ ICG Imaging System is contraindicated in patients suspected of being infected with Creutzfeldt-Jakob disease.

ICG Imaging Agent Contraindications:

- ICG contains sodium iodide and should be used with caution in patients who have a history of allergy to iodides or iodinated contrast agents.
- The KARL STORZ ICG Imaging System should not be used with patients who are known to be sensitive to iodides or iodinated contrast agents.

Please refer to the KARL STORZ ICG Imaging System Manual IM-000012 for complete instructions on using the system. Also refer to the CCU and the video camera head and light source instruction manuals for the individual functions of each respective device that is used with the KARL STORZ ICG HOPKINS® Endoscope.

4 Safety instructions

WARNING: Risk of injury and damage to products: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals and always follow the instructions given precisely. Check the functioning of the products used in combination.

WARNING: Risk of injury: Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all of the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

WARNING: Risk of injury: Incorrect application of medical instruments poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

WARNING: Risk of injury and risk of damaging the products: If applied parts are used outside of the field of vision, there is a risk that tissue and accessories could be damaged unintentionally. Always hold the applied parts of the active electrode, laser and other instruments that transmit energy pointed at the target in the field of vision during use.

WARNING: Risk of injury and risk of damaging the products: The highly concentrated light energy can cause the distal end, the light port, adjacent components and tissue in front of the light emission window to heat up (can exceed 105.8°F [41°C]). This can cause burns to the patient, user and surgical accessories.

- Avoid direct tissue contact with the distal end, as well as with the light port of the telescope and the light cable.
- Do not rest the endoscope or the light cable on the patient or where they are in direct contact with surgical accessories.
- Always select the lowest possible light output of the cold light source which still allows optimal illumination of the operating field.

WARNING: Risk of injury: Optical radiation poses a risk of injury to eyes. Never look into the free end of either a connected light cable or endoscope.

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WARNING: Risk of injury: Patient leakage currents may accumulate if endoscopes and endotherapy devices supplied with energy are used together. This is particularly important if an endoscope with a Type CF applied part is used. In this case, an endotherapy device with a Type CF applied part must be used in order to minimize the total patient leakage current.

WARNING: Risk of injury: Overloading the instrument by exerting too much force may cause the medical device to break, bend, or malfunction, and consequently injure the patient or user. Do not overload the instruments. Do not bend bent instruments back to their original position.

WARNING: This device is not indicated for tumor biopsy and resection in newborns (<1 month) and infants (1 month - 2 years).

WARNING: This device is not indicated for aqueduct exploration in newborns (<1 month) and infants (1 month - 2 years).

WARNING: Always adjust the light source to the minimum illumination intensity necessary to achieve optimum illumination of the endoscopic scene, either by direct vision or couple to a video camera. The higher the light intensity setting of the light source is, the greater the heat energy that will be generated at the distal end of the endoscope. Turn the light source to standby mode or initial mode when the endoscope is not in use.

WARNING: Never place the end of a light cable or endoscope connected to a light cable on or under a surgical drape while the light source is activated. The intensity of the light may cause burns to the patient and/or the surgical drape.

WARNING: Avoid excessively long (> 10 mins) continuous direct tissue contact with the shaft.

5 Precautions

In order to avoid damage to the sheath and rod lens system and thus ensure the long-term quality of the telescope, the following information must be observed:

- For cleaning, sterilization and storage, ideally store the telescopes individually or use containers in which the telescopes can be secured.
- Telescopes must always be held by the eyepiece and never by the distal end alone.
- Do not bend the sheath under any circumstances.
- Place the telescopes down with care. Avoid hard impacts, especially at the distal end.

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8 Reprocessing

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a device of infection as well as the risk that the medical device may malfunction.

WARNING: Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out, or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

WARNING: If CJD (Creutzfeldt Jakob disease) is suspected or indicated, KARL STORZ recommends disposing of the product in accordance with national regulations and not using it again.

CAUTION: When preparing and using solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: Risk of damage to medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices. Only use chemicals approved by KARL STORZ for reprocessing.

CAUTION: National laws and regulations must be observed.

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CAUTION: The optionally available special cleaning paste should only be used for cleaning when the endoscopic image is visibly cloudy, under no circumstances should the paste be used on a regular basis.

CAUTION: Do not use Immediate-Use Steam (FLASH) sterilization for KARL STORZ ICG HOPKINS® II endoscopes.

8.1 Accessories

Accessories required for carrying out reprocessing:

Brushes:

Surface: 27652

8.2 Disassembly

Prior to cleaning, the telescope must be separated into its individual components as far as possible. Both adaptors on the light inlet piece are to be disassembled (see figure 1).

8.3 Manual Pre-Cleaning

Recommended Water Quality: Utility Water

1. Disconnect the telescope from the camera head.
2. Wipe the camera head thoroughly with a soft, lint-free disposable cloth moistened with water or a diluted mild/neutral pH enzymatic cleaner (e.g. Enzol) to remove all heavy soiling.
3. Contain and transport the contaminated endoscope to the decontamination area for cleaning.

NOTE: If immediate cleaning is not feasible, place the endoscope in a container and soak for no longer than 30 minutes in a diluted neutral/mild pH enzymatic cleaning solution (e.g. Enzol).

NOTE: When preparing cleaning solutions, follow the manufacturer's instructions for proper concentration and temperature.

8.4 Manual Cleaning

Recommended Water Quality: Utility Water

CAUTION: WEAR PROTECTIVE GLOVES, CLOTHING, AND A FACE MASK FOR CLEANING OF CONTAMINATED INSTRUMENTS.

1. Remove light cable adaptors before cleaning. Do not clean or rinse endoscopes with other instruments to prevent damage to the endoscope.
2. Thoroughly rinse the endoscope and light cable for a minimum of 2 minutes with water to remove all gross debris.
3. Using a low-lint soft cloth dampened with cleaning solution, wrap the cloth around the light cable and wipe the entire cable portion from end to end, back and forth, for at least two passes.
4. Completely immerse the endoscope and light cable in a diluted neutral/mild pH enzymatic cleaning solution (e.g. Enzol). Keep immersed for at least 5 minutes.
5. While immersed, brush all surfaces of the endoscope with a soft bristle brush (P/N: 27652). Do not brush the glass window of the device, but wipe carefully using a lint-free, soft cloth. Rinse the bristles of the brush with water if heavy debris collects on it. Keep brushing until there is no visible debris on the bristle of the brush (at least 2 minutes).
6. Brush the hard surfaces of the light cable with a soft bristle brush (P/N: 27652). Pay close attention to the stainless steel locking mechanism, turning the locking mechanism left and right, brushing for at least 30 seconds. Do not brush the fiber optic ends of the light cable. Instead, use a lint-free, soft cloth

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to wipe the fiber optic ends of the cable and all flexible areas of the cable. Rinse the bristles of the brush with water if heavy debris collects on it. Keep brushing until there is no visible debris on the bristles of the brush (at least 2 minutes).

7. After the immersion period of 5 minutes, remove the endoscope from the cleaning solution and completely immerse it in a large volume of room temperature water. Gently agitate the device while immersed and keep immersed for at least 1 minute. Discard the water and repeat two more times with fresh water for a total of 3 immersion rinses.
8. Remove the light cable from the cleaning solution and rinse with running cold water for at least 1 minute.
9. Do not brush the glass window of the light cable or endoscope. Wipe with a dampened cloth instead.
10. Dry the endoscope and light cable with a soft, lint-free cloth or clean, dry, oil-free compressed air (<5 psi).
11. After cleaning, inspect the endoscope and light cable for cleanliness and damage. If visible debris is observed, repeat the cleaning steps. Do not use the endoscope if visible signs of damage are observed.
12. Inspect the endoscope and light cable for any visible signs of damage (cuts, dents, or scratches) that may trap residual debris. Do not use the endoscope or light cable if it is damaged or is difficult to use. The devices should be returned to KARL STORZ for repairs or replacement.

8.5 Assembly, inspection and care

The cleaned medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- Damaged or corroded medical devices must be withdrawn from use.
- Afterwards, a functional check must be carried out.

8.6 Life span

The following function tests are to be carried out in order to recognize any loss of material through wear:

1. Check the surface of the medical device for mechanical integrity and changes.
2. Check if the image is transmitted.
3. Check if light is transmitted.

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The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

8.8.1 Steam sterilization using the fractionated prevacuum procedure

1. Place the telescope in an FDA-cleared sterilization tray.
2. Double wrap the loaded tray with an FDA-cleared sterilization wrap. If using a sterilization-sealed container, refer to the manufacturer's instructions for use on using the correct disposable filters for the indicated sterilization method.

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- Use the following cycle parameters for steam sterilization:
 Temperature: 270°F (132°C)
 Exposure Time: 4 minutes
 Dry Time: 30 minutes
- Following cycle completion, crack open the sterilizer door with the container still within the sterilization chamber. Allow the load to cool for 15 minutes with the door cracked open to eliminate any condensation.

8.8.2 STERRAD® Systems

WARNING: STERRAD® Sterilization may cause cosmetic changes to the device that do not necessarily impact the functionality of the device.

CAUTION: All telescopes must be thoroughly DRIED before loading into the STERRAD® System chamber. Loads containing moisture may cause a cycle cancellation.

CAUTION: Use only STERRAD® Instrument trays in the sterilization chamber. These trays are specially designed to allow the plasma to surround the items.

CAUTION: If using wrapped trays, use only polypropylene sterilization wrap and polyolefin pouches. Do not use paper pouches or sterilization wraps containing wood pulp or cotton.

NOTE: Instruments that KARL STORZ has determined to be compatible with the STERRAD® System sterilization process have been validated with at least one hundred STERRAD® cycles.

- Place the telescopes in a STERRAD® compatible tray.
- Double wrap the tray with two single layers of FDA-cleared polypropylene wrap, using a sequential double wrapping technique. If using a sterilization sealed container refer to the manufacturer's instructions for use on using the correct disposable filters for the indicated sterilization method.
- Load the tray into the STERRAD® sterilizer.
- Start the system with the following cycle:
 - STERRAD® 100S: Only one cycle available
 - STERRAD® 100 NX: "Standard" cycle
 - STERRAD® NX "Standard" Cycle
- Please consult the STERRAD® System Operators Manual for detailed instructions for use.

8.8.3 V-PRO® Systems

Amsco V-PRO® Sterilization Systems

- Place the telescopes in a V-PRO® compatible tray.
- Double wrap the tray with two single layers of FDA-cleared polypropylene wrap, using a sequential double wrapping technique. If using a sterilization sealed container refer to the manufacturer's instructions for use on using the correct disposable filters for the indicated sterilization method.
- Load the tray into the V-PRO® sterilizer.
- Start the system with the following cycle:
 - V-PRO® 1: Only one cycle available
 - V-PRO® 1 Plus: "Lumen" or "Non-Lumen" cycle
 - V-PRO® max: "Lumen" or "Non-Lumen" cycle
 - V-PRO® 60: "Lumen" or "Non-Lumen" cycle
- Please consult the V-PRO® System Operators Manual for detailed instructions for use.

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9 Overview: Recommended combinations*

If there is an identification symbol on the light cable and the endoscope (see table), assign matching combinations according to the symbols shown.

Light cable diameter	Endoscope diameter
4.8-5.0 mm 	6.6-12.0 mm 
3.0-3.5 mm 	3.0-6.5 mm 
2.0-2.5 mm 	0.8-2.9 mm 

* Special endoscopes may deviate from this

10 Operating/ storage conditions

	Temperature	Rel. humidity
Storage	-4 °F ... +140 °F	10 % ... 90 %
Application	+50 °F ... +104 °F	30 % ... 75 %
Air pressure	700 hPa - 1060 hPa	

11 Complaint shipments

In the interests of our employees' health, only sterilized products are accepted as return shipments. These should be marked as such; otherwise no further processing may be performed.

12 Disposal

To dispose of these medical devices, no special measures are necessary. National laws and regulations must be observed.

13 Directive compliance

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

14 Packaging symbols

For the meaning of the symbols printed on the packaging, please refer to the first page of this manual.

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1	Important information	1
2	Symbols employed	3
3	Intended use	6
3.1	Indications	6
3.2	Contraindications.....	7
4	User qualification.....	8
5	Warnings and cautions	8
6	Safety instructions.....	9
6.1	Unpacking	18
6.2	Basic equipment.....	19
6.3	Operating conditions	20
7	External system components.....	20
8	Preparation before use	21
9	Handling	21
9.1	Operation	21
9.2	Holding the device correctly	21
9.3	Operating the suction valve	22
9.4	Handling the sheath.....	22
9.5	Handling the deflection mechanism	23
9.6	Handling the working channel.....	24
9.7	Handling the accessories (forceps, cleaning brushes, HF electrodes etc.)	26
9.8	Inspecting the surfaces.....	28
9.9	Inspecting the deflection.....	29
9.10	Connecting to IMAGE1 S™	29
9.11	Focus	29
9.12	Function buttons of the endoscope	30
9.13	Switching the flexible cystoscope HD VIEW ON/OFF	31
10	Reprocessing	32
10.1	General warnings.....	32
10.2	Accessories	33
10.3	Precleaning directly after use	34
10.4	Leakage test.....	35
10.5	Disassembly	37

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10.6	Manual disinfection	40
10.6.1	Manual High Level Disinfection (HLD)	40
10.6.2	Drying	41
10.7	Machine cleaning and disinfection	41
10.7.1	Manual preparation	41
10.7.2	Connecting	42
10.7.3	Machine cleaning/chemical-thermal disinfection	42
10.8	Assembly, inspection, care and packaging	43
10.8.1	Packaging	44
10.9	Sterilization	44
10.10	Limits of reprocessing	46
10.11	Storage	46
11	Reprocessing of LUER port with stopcocks	48
11.1	Accessories	49
11.2	Preparation for cleaning and disinfection	49
11.3	Manual precleaning	49
11.4	Machine cleaning and disinfection	50
11.5	Connecting	51
11.6	Assembly, inspection and care	51
11.7	Packaging systems	51
11.8	Sterilization	52
11.9	Limits of reprocessing	52
12	Technical data	53
13	Servicing and repair	55
14	Maintenance	55
15	Limitation of liability	55
16	Warranty	55
17	Directive compliance	56
18	Important information	56
19	Disposal	56
20	Electromagnetic Compatibility user information	57

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1 Important information

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper handling, cleaning, disinfection and, if appropriate, sterilization of the flexible cystoscope HD VIEW. All the necessary details and actions are clearly explained. We therefore ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for easy reference.

It is recommended that the suitability of the products for the intended procedure be checked prior to use.

KARL STORZ is continuously working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, figures and descriptions in this manual.

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2 Symbols employed

	Manufacturer
	Consult instructions for use
	Catalogue number
	Number of products in the product packaging
	Serial number
	Non-sterile
	Fragile, handle with care
	Keep dry
	Temperature limit
	Federal (USA) law restricts this device to sale by or on the order of a physician.

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


- ① Connector
- ② Suction valve
- ③ Connection for leakage tester or pressure compensation cap (vent port)
- ④ Handle
- ⑤ Irrigation and working channel
- ⑥ Deflection lever
- ⑦ Distal tip
- ⑧ Working channel connector seal 11014 DS
- ⑨ Leakage tester 13242 XL
- ⑩ Suction valve, sterile, for single use 091011-20
- ⑪ Suction valve, reusable 11301 CB1
- ⑫ Pressure compensation cap 11025 E
- ⑬ LUER port with stopcocks 11014 L
- ⑭ Double LUER port 11014 LU
- ⑮ Cleaning adaptor for working channel 11014 RA
- ⑯ Cleaning adaptor 11301 CD1
- ⑰ Right navigation button
- ⑱ Left navigation button
- ⑲ Menu button

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3 Intended use

Flexible cysto-urethrosopes (fiber endoscopes and videoscopes) are used for endoscopic imaging of the urethra and the bladder for the purposes of diagnosis and treatment. Flexible cysto-urethrosopes with an integrated working channel allow for the use of flexible auxiliary instruments such as forceps, stone baskets etc., coagulation electrodes and laser fibers for lithotripsy or for the ablation of benign/malignant tissue. With flexible cysto-urethrosopes, the urethra and the entirety of the bladder can be visualized/treated. The special feature of the flexible cysto-urethroscope lies in the fact that it can be used percutaneously in the kidneys or suprapubically in the bladder via an operating sheath.

Flexible cysto-urethrosopes (fiber endoscopes and videoscopes) are designed for temporary or short-term use in endoscopic interventions performed through natural orifices or surgical invasive interventions.

3.1 Indications

Description of all known indications:

The use of cysto-urethrosopes is indicated if urethro-cystoscopy is indicated in the opinion of the attending physician.

This includes, amongst other things:

- Clarification of micro and macro hematuria
- Clarification/biopsy in the event of suspected bladder tumors
- Control urethro-cystoscopy as per TUR-B or P
- Clarification in particular in the event of suspected tumor infiltration due to adjacent organ cancer (e.g., colon carcinoma, cervical carcinoma)
- Clarification or treatment in particular in the event of suspected vesicoenteric or obstetric fistulas
- Urethral or vesical foreign body removal
- Clarification of recurring urinary tract infections
- Clarification of LUTS/BES, subvesical obstructions, BOO, neurogenic micturition disorder

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Intended use

- Visual assessment and possibly treatment (e.g., injection) of the orifices of the ureters (localization, condition, etc.)
- Visual clarification of the bladder sphincter in the event of incontinence
- Visual assessment and possibly treatment of the urethra (tumor, stricture/stenosis, fistula)
- Visual assistance for indwelling catheter system, e.g., with urethral trauma in the event of a pelvic fracture
- Retrograde pyelography, ureteral splinting

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3.2 Contraindications

The use of the cysto-urethrosopes is contraindicated if, in the opinion of the attending physician, the surgical or examination method as such is contraindicated or if the patient is not able to undergo surgery or anesthesia due to his or her general condition.

Cysto-urethrosopes must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

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- Acute urethritis
- Acute epididymitis
- Acute prostatitis
- Acute/florid cystitis/urinary tract infection
- Possible coagulation disorders
- Restless, agitated patients

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4 User qualification

The flexible cystoscope HD VIEW from KARL STORZ may only be used by persons with an appropriate medical qualification and who are acquainted with the flexible cystoscopy technique. The information given in these instructions only serves to instruct in the correct handling and reprocessing of the flexible cystoscope HD VIEW. It is not intended as an introduction to the techniques of flexible cystoscopy.

5 Warnings and cautions

HAZARD SYMBOL: This is a hazard symbol. It warns of the risk of injury. To prevent injuries, all the measures marked with the hazard symbol must be followed.

Signal words:

WARNING! designates a possible imminent risk. If this is not avoided, it **could** lead to death or serious injuries.

ATTENTION! designates a possibly harmful situation. If it is not avoided, the medical devices could be damaged.

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6 Safety instructions

The flexible cystoscope HD VIEW must only be used according to the medical rules and procedures of endoscopy recognized for this method.

WARNING:

Risk of injury and damage to the product:

- It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces or voltage-carrying elements of other devices at any time.

WARNING:

Risk of burns and glare: The LED light emitted through the flexible videoscope may cause glare or high temperatures, and thereby burns, at the distal light outlet.

- Never look directly into the light emitted by the flexible videoscope.

WARNING:

Risk of explosion: If flammable anesthetics are used in the immediate proximity of the flexible videoscope.

- The flexible videoscope must not be used in the presence of flammable or combustible media.

WARNING:

Risk of injury: When in operation, the integrated LED light source heats up the flexible videoscope.

- Switch off the flexible videoscope before putting it down, and do not place it in the vicinity of the patient (as per chapter 16 of the 3rd edition of IEC 60601-1).

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 WARNING:

Risk of injury: High energy light radiated through endoscopes may give rise to high temperatures in front of the light outlet, at the end faces of the fiberoptic light cable and at the tip of the endoscope.

- *Prolonged direct contact of light-emitting surfaces on tissue can lead to burns; this must therefore be avoided at all costs.*

 WARNING:

Risk of infection: The flexible videoscopes are not sterile when delivered, and must be cleaned and disinfected/sterilized prior to initial use and each subsequent reuse.

 WARNING:

Risk of injection: To avoid infection and maintain functionality, the chapter on reprocessing in the instruction manual must be observed. The cleaning, disinfection and sterilization procedures are explained therein in detail.

 WARNING:

Risk of infection: Before the device is used for the first time, it must be subjected to a complete reprocessing cycle, excluding precleaning.

 WARNING:

Risk from electromagnetic incompatibility: Medical electrical devices are subject to special precautions regarding Electromagnetic Compatibility (EMC). Observe the stated EMC instructions during installation and operation.

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WARNING:

The use of accessories other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased interference immunity of the flexible videoscope. The accessories listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using accessories other than those specified here, it is the responsibility of the operator to ensure that they comply with EN/IEC 60601-1-2.

WARNING:

Risk of injury: The videoscope does not offer any electric insulation against high frequency (HF) voltages. Protection against HF voltage must be achieved using the appropriate HF instrument set.

➤ Please read the instructions for use in the instruction leaflet supplied with the HF application part for more information.

WARNING:

Combinations of medical devices are only assured to be safe if they are listed as such in the respective instruction manuals or the intended use and interface specifications of the products used in the combination permit this.

➤ The instruction manuals and interface specifications of the medical devices used in combination must be observed precisely.

WARNING:

The pressure compensation cap must be removed from the vent port for:

- endoscopic interventions,
- immersing in liquid, e.g., disinfectant solution and for machine reprocessing

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Safety instructions

⚠️ WARNING:

The pressure compensation cap must be in place on the vent port for:

- gas sterilization (EO, FO, hydrogen peroxide)
- transport and shipping

⚠️ WARNING:

Only use accessories that are specifically suited for use with videoscopes in accordance with the intended use and which can be employed in line with generally accepted technical standards and safety guidelines.

⚠️ WARNING:

Risk of injury: The flexible videoscope must not be used during the discharge of a cardiac defibrillator.

- Remove the flexible videoscope from the operation area before discharge takes place.

⚠️ WARNING:

Risk of injury: Only KARL STORZ devices may be used in combination with the flexible videoscope. Disregarding this warning can endanger the patient and user.

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ATTENTION: Do not allow any live components or instruments to come into contact with the flexible videoscope.

NOTE: Do not hit the distal tip of the flexible videoscope against hard objects, and never use a sharp object to remove dirt.

NOTE: Never place heavy objects on the flexible videoscope.

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Safety instructions

NOTE: Never clamp or fold the sheath of the flexible videoscope. This can result in damage to the videoscope.

NOTE: Constant changing between various reprocessing methods is especially hard on the materials and must therefore be avoided.

NOTE: The flexible videoscope must not be carried by the sheath. Do not pull, clamp or twist the sheath of the videoscope.

NOTE: The sheath of the flexible videoscope must never be clamped or coiled up too tightly. It is best to store the video videoscope in a suspended position or to keep it sterilized in a tray specifically designed for sterilizing and storage (39406AS).

WARNING:

Risk of injury from clouding of the lens.

- Do not apply lubricant to the lens, to avoid clouding the image and thereby prevent risks for the patient!

WARNING:

When inserting the flexible videoscope into the patient, the distal tip of the flexible videoscope must be in the straight, non-deflected position.

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 WARNING:

Risk of injury: Before inserting or removing the flexible videoscope, the distal tip must be put into the neutral position, as otherwise injuries to the patient or damage to the instrument may occur.

 WARNING:

Risk of injury: If the distal tip can no longer be properly controlled (sticking etc.) via the deflection lever, or if the deflection mechanism develops any other fault, the flexible videoscope must not be used.

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- Place the distal tip of the flexible videoscope in the neutral, non-deflected position and remove the flexible videoscope slowly and carefully from the patient.

NOTE: The resistance that can be felt when deflecting the distal tip using the deflection lever increases as the angle of deflection increases. Actuate the deflection gently and carefully, without any sudden or jerky movements, as far as the limit stop.

 WARNING:

Risk of injury: The flexible videoscopes from KARL STORZ are not intended for use with light sources. The use of non-compatible devices can lead to patient injury or device damage.

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WARNING:

Risk of injury: Whenever the flexible videoscope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition. Damaged flexible videoscopes or damaged accessories must not be used.

NOTE: There should be no instruments in the suction channel whilst suctioning is being performed. The instrument channel inlet should be sealed to prevent additional air from being drawn in.

NOTE: Do not squeeze the distal tip of the flexible videoscope too tightly.

WARNING:

Risk of injury: Actuate the deflection lever gently and carefully as far as the limit stop. The deflection lever must never be operated abruptly or with force, to avoid injuring the patient or damaging the device.

- Move the deflection lever smoothly and slowly!

WARNING:

Risk of injury: Before each use, the functionality of the deflection lever and the connected distal tip must be checked.

- Defective flexible videoscopes must not be used.

WARNING:

Risk of injury: If the image fails or there is interference during use in the patient, put the distal tip into the neutral position and remove the device from the patient's body gently and carefully. If the flexible videoscope continues to be used in the patient after the image has failed or when there is interference, this may result in serious injury to the patient.

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Safety instructions

WARNING:

Risk of injury: Only insert suitable accessories, e.g., biopsy forceps, into the instrument channel of the videoscope when the distal tip is in a neutral, non-deflected position.

- > Once the instrument has been introduced, carefully actuate the lever gently as far as the limit stop.
- > The deflection lever must never be operated abruptly or with force, to avoid injuring the patient or damaging the device.

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NOTE: Do not apply force to insert accessories, e.g., biopsy forceps, through the working channel, to prevent damaging the device. It should be possible to move the accessory easily, and without any hindrance, through the entire channel.

NOTE: When the distal tip is deflected, do not insert instruments into the videoscope. Only insert the instrument when the distal tip is in the neutral position and deflect the deflection lever gently and carefully up to the limit stop.

NOTE: Inserting sharp-tipped accessories into the instrument channel may damage the videoscope.

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 WARNING:

Risk of injury: Foreign bodies and items to be extracted must not be removed via the instrument channel with the aid of instruments. Instead, the entire videoscope must be taken out of the body and the foreign body must then be removed.

NOTE: For the correct settings for the camera control unit, please see the operating instructions of the devices used.

NOTE: Accessories such as biopsy forceps that are kinked or bent will not function properly, are difficult to guide through the channel and may cause severe damage to it. They must therefore be replaced.

 WARNING:

Ensure that the working channel is completely dry prior to sterilization to avoid damage to the working channel.

 WARNING:

Always unplug the flexible videoscope before carrying out any maintenance and cleaning procedures on it.

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6.2 Basic equipment

Flexible cystoscope HD VIEW 11272 VHK/VHU with following accessories:

Set 11272VHK

11272 VH	Flexible cystoscope HD VIEW
11301 P	Seal for suction valve seat
11014 L	LUER port with stopcocks
11014 DS	Working channel connector seal
11014 RA	Cleaning adaptor for working channel
11301 CD1	Irrigation adaptor

Set 11272VHUK

11272 VHU	Flexible cystoscope HD VIEW
11301 P	Seal for suction valve seat
11014 LU	Double LUER port
11014 DS	Working channel connector seal
11014 RA	Cleaning adaptor for working channel
11301 CD1	Irrigation adaptor

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Optional accessories

1 x grasping forceps for small fragments 27023 FE
1 x biopsy forceps 27023 ZE
1 x stone basket 27023 KF
1 x fiber fixation kit 11014 YK
1 x coagulating electrode 4 Fr. 27723 T
1 x seal for instrument ports 27550C-10
1 x LUER adaptor with seal 27014 Y
1 x endoscopic seal 100010-10
1 x cleaning adaptor for instrument ports 27001 RA
1x single-use cleaning brushes 110940-50
1 x cleaning brush 27651 B

6.3 Operating conditions

Flexible cystoscope HD VIEW:

Ambient temperature: 15°C ... 30°C
Humidity: 15% – 80%

7 External system components

Further information on external system components can be found in the IMAGE1 S™ (96206286) instruction manual.

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
8 Preparation before use

1. Assemble LUER attachment. To do this, insert the lug of the LUER attachment into the recess on the working channel. Push and turn clockwise until it audibly clicks into place.

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2. Assemble the suction valve/valve seal.

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9 Handling

NOTE: Any damage to the flexible cystoscope HD VIEW resulting from incorrect operation is not covered by the manufacturer's warranty.

9.1 Operation

The following instructions outline the correct operation of the flexible cystoscope HD VIEW and the corresponding accessories. This is not an introduction to the techniques of cystoscopy. Instructions on endoscopic technique can be found in the relevant medical literature.

9.2 Holding the device correctly

Grip the handle of the flexible cystoscope HD VIEW in your left or right hand, as shown in the picture. With this kind of control grip, the flexible cystoscope HD VIEW can be operated with one hand.

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9.3 Operating the suction valve

When the suction valve is assembled, a suction pump must be connected to the suction valve via a tube. Suction is activated upon pressing the valve button.

9.4 Handling the sheath

ATTENTION: Risk of damage to the sheath:

- Do not carry the flexible cystoscope HD VIEW by the sheath.
- Do not pull, clamp or twist the sheath of the flexible cystoscope HD VIEW.
- The sheath of the flexible cystoscope HD VIEW must not be kinked or coiled up too tightly.
- It is best to store the flexible cystoscope HD VIEW in a suspended position or to keep it sterilized in a tray specifically designed for sterilizing and storage (39403 AS).

WARNING:

Risk of injury: Do not apply lubricant to the lens, to avoid clouding the image and thereby prevent risks for the patient!

If, during use, the image fails, is disturbed, or becomes blurred, or if vision in general is impeded, put the deflection mechanism into the neutral position and slowly pull the sheath of the flexible cystoscope HD VIEW out of the patient's body.

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9.5 Handling the deflection mechanism

NOTE: The neutral position of the distal tip is straight, without any deflection whatsoever.

Deflection functions as follows with type 11272 VH endoscopes:

If the lever is pushed upwards, the tip is lowered.

If the lever is pushed downwards, the tip is deflected upwards.

This mechanism is inverted with model 11272 VHU.

ATTENTION: When the flexible videoscope is inserted into the patient, the distal tip must be in a straight, non-deflected position.

⚠️ WARNUNG:

The deflection lever must never be operated with force or abruptly to avoid injuring the patient or damaging the device!

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⚠️ WARNUNG:

Risk of injury and damage from faulty deflection mechanisms: If the distal tip can no longer be properly controlled (sticking etc.) via the deflection lever, or if the deflection mechanism develops any other fault, the flexible videoscope must not be used.

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⚠️ WARNING:

Risk of injury: If the flexible videoscope malfunctions during use on the patient, discontinue use immediately.

- *Place the distal tip of the flexible videoscope in the neutral, non-deflected position and remove slowly and carefully from the patient.*

⚠️ WARNING:

Risk of damage: The lens at the distal end of the flexible videoscope may become detached, e.g., due to mechanical forces.

- *For this reason, it is essential to check that the lens is undamaged and present before and after every decontamination.*

9.6 Handling the working channel

Do not force instruments or cleaning brushes through the channels. This can lead to damage to the working channel, especially if the tip of the instrument is deflected.

When connecting and operating suction/irrigation devices, pay attention to the instruction manuals for the particular devices. The irrigation devices recommended by KARL STORZ are set to a maximum pressure of 54 kpa (400 mmHg/0.54 bar/7.8 psi).

- The LUER adaptor 27014 Y or the seal 27550 C can be used to seal the instrument channel when inserting instruments.
- Instruments can be introduced through the instrument channel of the flexible cystoscope HD VIEW. Make sure that the accessories are of the appropriate length and diameter. Unsuitable accessories can damage the flexible cystoscopes HD VIEW. If you have any questions concerning accessories, please contact KARL STORZ directly or an authorized specialist dealer.

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Handling

- The flexible cystoscope HD VIEW must be in the straight, non-deflected position when inserting an instrument into the working channel.

Forceps, scissors and other accessories whose jaw can be opened or closed must always be inserted into the instrument channel with closed jaws. The jaws may only be opened when the distal tip of the instrument is completely clear of the channel. Opening the jaws in the instrument channel can cause damage to the videoscope.

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ATTENTION: Risk of damage: Only use accessories (guide probes, cleaning brushes, etc.) from KARL STORZ. Accessories from other manufacturers may damage the videoscope.

9.7 Handling the accessories (forceps, cleaning brushes, HF electrodes etc.)

WARNING:

Risk of injury: When HF electrodes or laser probes are used, it must be ensured that power is not applied until the active parts have left the instrument channel and are in the user's field of vision, and the site of the operation is clearly recognizable. Activation in the instrument channel may lead to serious injury to the patient and damage to the videoscope.

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ATTENTION: Risk of damage: Only use
accessories (guide probes, cleaning brushes,
etc.) from KARL STORZ. Accessories from other
manufacturers may damage the videoscope.

NOTE: Do not apply force to insert accessories,
e.g., biopsy forceps, through the working channel,
to prevent damaging the device. It should be
possible to move the accessory easily, and without
any hindrance, through the entire channel.

NOTE: When the distal tip is deflected, do not
insert instruments into the videoscope. Only insert
the instrument when the distal tip is in the neutral
position and deflect the deflection lever gently and
carefully up to the limit stop.

NOTE: Inserting sharp-tipped accessories into the
instrument channel may damage the videoscope.

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9.8 Inspecting the surfaces

⚠ WARNING:

Risk of injury: The surface of the flexible sheath must be inspected visually and by carefully feeling it, to detect any cracks, cuts, dents, foreign bodies or other forms of damage.

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Before each use, verify that the device is in good condition by checking the surface of the flexible cystoscope HD VIEW, especially the sheath area and the distal tip, for sharp edges. Cystoscopes HD VIEW with sharp edges must not be used.

Otherwise they may cause serious injury to the patient or user.

Check the connecting cable for damage (such as fractures, cracks, twisted or crushed areas, etc.) and the connecting plug for corrosion.

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⚠️ WARNING:

Risk from damage: Faulty flexible videoscopes must not be used.

- Before each use, check the functionality of the deflection lever and the distal tip.

The distal tip is deflected by actuating the deflection lever slowly.

The distal end will deflect in accordance with the direction of actuation (forward or backward).

9.10 Connecting to IMAGE1 S™

The flexible cystoscope HD VIEW is connected to the camera control unit for operation.

For operation, an IMAGE1 S™ system with an IMAGE1 S Connect® and an IMAGE1 X-Link are required. Both camera modules must be connected to each other.

The flexible cystoscope HD VIEW is then connected to the IMAGE 1X-Link.

9.11 Focus

The image focus is adjusted automatically. Manual focusing is not possible/necessary.

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10.1 General warnings

! WARNING:

Risk of infection: Only reprocess the medical device as per the documents provided with the device or the reprocessing instructions issued by KARL STORZ. These can be downloaded from www.karlstorz.com.

! WARNING:

Risk of infection from non-sterile delivery. Reprocess the device before each use. Check the medical device for visible soiling before use.

! WARNING:

Before reusing the device, replace any defective parts or send the flexible endoscope for repair.

ATTENTION: Risk from damage to the device: The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.

ATTENTION: Risk from damage to the device: Never clean flexible endoscopes with metal brushes.

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ATTENTION: Risk from damage to the device: Flexible endoscopes must never be steam sterilized (autoclaved).

ATTENTION: Risk from damage to the device: To prevent damage to the flexible endoscope and accessories, it must be ensured that all connections are completely dry prior to use.

NOTE: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

NOTE: Observe the chemical manufacturer's specifications concerning concentration and exposure times. Bear in mind the microbiological range of action of the chemicals.

NOTE: Damage to the device. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be found at www.karlstorz.com

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10.2 Accessories

Accessories required for carrying out reprocessing: pressure compensation cap, manual leakage tester

Place of application	Art. number brush/ adaptor
Suction channel	27651 H, 11301 CD1
Working channel	27651 H, 11014 RA
Clean surface	27652
Vent port	13242 XA, 11025 E

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10.3 Precleaning directly after use

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the endoscope immediately after use.

Wipe-clean the exterior surfaces of the endoscope with a damp, single-use cloth.

Rinse out the working channel in order to prevent drying.

NOTE: To avoid damaging the sheath, do not use pushing and pulling movements.

NOTE: Do not use disinfectants for precleaning; these cause fixation of proteins and organic residues.

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
10.4 Leakage test

The endoscope should be tested to ensure it is leakproof after use and before reprocessing.

Performing the leakage test – dry

The connection tube of the leakage tester and the vent port on the endoscope must be dry. Connect the leakage tester connector securely to the supply connector vent port and lock by turning the bayonet catch clockwise 25°.

NOTE: The rubber seal in the connection adaptor of the leakage tester should be treated regularly with instrument oil.

Inflate the leakage tester until the maximum test pressure of 160–180 mmHg (blue area) is reached. The deflection rubber on the distal end may expand slightly. Using the deflection mechanism on the control body, move the deflection part in all possible directions as far as possible and observe the pointer on the manometer. The test pressure can drop by up to 10 mmHg until the endoscope is completely filled with air.

After stabilizing the pressure on the manometer, wait at least 2 minutes to see whether the pressure remains constant or drops.

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Performing the leakage test – in liquid

NOTE: To avoid damage from the penetration of liquid, a test in the dry state must be performed in advance.

NOTE: Any lumina must be covered with no air bubbles by filling in a targeted manner.

Immerse the device in solution so that it is completely covered. It must be ensured that the leakage tester does not become wet. Perform the leakage test as described in 'Leakage test – dry' and note whether bubbles emerge along the device. If, after an observation period of at least 3 minutes, no bubbles can be seen emerging (continuously) from the device and, in addition, the manometer does not show a fall in pressure, the device is leakproof.

Remove the endoscope from the liquid and press the vent button on the manometer. The test pressure must drop to 0 before the leakage tester can be removed.

NOTE: Never release the pressure in the liquid.

Cleaning can then begin.

If a pressure loss occurs, this means that the endoscope leaks and may no longer be used, immersed in liquid or reprocessed by machine.

ATTENTION: In this case, fill out the accompanying repairs form (see Appendix or Art. no. 96096004DER), mark the endoscope, which has been cleaned by wiping to remove organic soiling/residue, as not disinfected and return it to the Service Department for repair or, alternatively, contact KARL STORZ. Line the storage/carrying case with the transport protection sheet 13990 SFN.

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10.5 Disassembly

The flexible cystoscope HD VIEW must be completely dismantled for cleaning/disinfection.

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1. Dismantle the suction valve by pressing and turning.

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2. Dismantle the LUER attachment by pressing and turning

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
3. Dismantle the stopcock (if present) on the LUER attachment.

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5. Clean the channels and lumina of the LUER attachment with a suitable brush.

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6. Dismantle the reusable suction valve 11301 CB1 by unscrewing the metal cap.

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8. Brush the working channel connector.

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9. Brush the working channel and suction channel.

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10. Connect the cleaning adaptor 11301 CD1 to the valve seat.

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10.6 Manual disinfection

The endoscope must be completely immersed in a cleaning solution. The suction and working channel must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

10.6.1 Manual High Level Disinfection (HLD)

The endoscope must be completely immersed in a high level disinfectant solution. The working and suction channel must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

The endoscope can be disinfected with the following high level disinfection solutions which have been tested by way of example:

- 2.4% glutaraldehyde solution (e.g., CIDEX® Activated Dialdehyde Solution)
- 0.55% ortho-phthalaldehyde solution (e.g., CIDEX® OPA Solution)
- 2.0% accelerated hydrogen peroxide solution (e.g., Revital-Ox® Resert™ or Revital-Ox Resert XL HLD)

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10.6.2 Drying

Finally, the surfaces and channels/lumina are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

ATTENTION: The pressure must not exceed 100 kPa (1 bar/14.5 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.

10.7 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003DER):

10.7.1 Manual preparation

The endoscope must be completely immersed in a cleaning solution. The working and suction channel must be covered with no air bubbles by filling in a targeted manner. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a suitable brush. The channels are brushed until the brush head emerges and is visibly clean. Before retracting the brush head that emerges at the distal end, clean it. Never brush back and forth within the channel.

Finally, it must be rinsed with cold water to ensure neutralization.

NOTE: KARL STORZ recommends using the adaptors 11301CD1 and 11014RA for simpler filling of the lumina.

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10.7.2 Connecting

In order to ensure effective machine cleaning and disinfection, the endoscope's working and suction channel must be connected to the washer and disinfectant using the recommended adaptors 11301 CD1 and 11014 RA.

ATTENTION: The pressure compensation cap must not be in place on the vent port for cleaning and disinfection.

10.7.3 Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65°C must not be exceeded. The relevant national requirements and the chemical manufacturer's instructions must be taken into account when using this method.

The selection of a suitable endoscope holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

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
10.8 Assembly, inspection, care and packaging

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- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.

NOTE: Clean the distal lens with a cotton tip applicator soaked in lens cleaner or 70% isopropyl alcohol. When cleaning the distal lens, always wipe with a circular movement.

- Damaged or corroded medical devices must be withdrawn from use.

NOTE: Inspect plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their original condition, they must be replaced. This especially applies to seals and rubber caps.

- If there is still residual moisture on or in the endoscope, it must be dried again manually. The surfaces (e.g., connection) and channels/lumina must be completely dried using (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

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- Afterwards, a functional check must be carried out.

NOTE: Check for a clear view and steerability.

NOTE: During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

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WARNING:

WARNING: Before sterilization, all parts of the endoscope must be thoroughly cleaned, disinfected, rinsed and then dried as otherwise sterility cannot be achieved.

WARNING:

Risk of infection: Sterilize the medical device in a disassembled state.

ATTENTION: Flexible endoscopes must never be steam sterilized (autoclaved).

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Reprocessing

ATTENTION: The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.

NOTE: When using reprocessing equipment, please observe the equipment manufacturers' instructions.

NOTE: All deviations from the recommended sterilization parameters must be validated by the user.

NOTE: Sterilization is not possible on surfaces which have been lubricated and oiled.

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Hydrogen peroxide (H₂O₂) sterilization - ASP STERRAD®

ATTENTION: Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

NOTE: Refer to the 'STERRAD® Sterility Guide' to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

The following STERRAD® sterilization methods have been validated and approved for this medical device by KARL STORZ:

- STERRAD® 100S Long Cycle with Boosters
- STERRAD® NX® Advanced Cycle
- STERRAD® 100NX® Flex Cycle
- STERRAD® 100NX® DUO Cycle

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**Hydrogen peroxide (H₂O₂) sterilization -
STERIS® AMSCO® V-PRO™1**

NOTE: The 'STERIS Low Temperature Reprocessing Device Compatibility Matrix' can be used to verify whether the medical device can be sterilized using the various STERIS V-PRO® devices.

The following STERIS V-PRO® sterilization methods have been validated and approved for this medical device by KARL STORZ:

- V-PRO® maX Flexible Cycle
- V-PRO® 60 Flexible Cycle

10.10 Limits of reprocessing

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

10.11 Storage

The storage method must be selected taking into account the respective applicable national requirements.

KARL STORZ only recommends dust-protected, suspended storage for decontaminated endoscopes.

Horizontal storage is possible for sterilized endoscopes. KARL STORZ recommends the use of a specially designed tray for sterilization and storage.

ATTENTION: Do not expose flexible endoscopes to direct sunlight or excessive heat.

ATTENTION: The distal tip of the flexible endoscope must be in the neutral position during storage.

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Reprocessing

ATTENTION: Do not store the flexible endoscope in the case.

ATTENTION: The sheath of the flexible endoscope must never be kinked or coiled up too tightly.

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11 Reprocessing of LUER port

⚠️ WARNING:

Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments (art. no. 96216003DER)' and the accompanying documentation.

⚠️ WARNING:

Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

⚠️ WARNING:

When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

ATTENTION: When preparing and using solutions, follow the chemical manufacturer's instructions, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

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ATTENTION: Danger of damage to medical devices: The use of chemicals that have not been approved by KARL STORZ or chemical manufacturers may cause damage to the medical devices. Only use chemicals approved for the medical device for reprocessing. A helpful guide is available on the Internet at www.karlstorz.com.

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NOTE: The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' (Art. no. 96216003) can be downloaded or requested by visiting www.karlstorz.com.

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Brushes: 27652, 27650B

Adaptor: 39107ALK (only for 11014L)

11.2 Preparation for cleaning and disinfection

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11.3 Manual precleaning

Disassembly (only valid for 11014L)

Prior to cleaning and disinfection, the medical device must be separated into its individual components as far as possible, and/or must be opened.

Brushing the surfaces

Any visible contamination or heavy soiling must be removed by cleaning the surfaces under cold running water with the aid of a brush or sponge.

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Brushing the lumina

The working channels, lumina and hollow spaces must be precleaned under cold running water with the aid of a suitable brush.

Rinsing out with a pressure water gun

The working channel, lumina and hollow spaces must be rinsed using a pressure water gun and a suitable rinsing attachment. For rinsing, pulse at least 3 times for 5 seconds from a jet of water with a pressure of at least 2.5 bar. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

11.4 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003):

Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65°C must not be exceeded. The relevant national regulations and the chemical manufacturer's instructions must be taken into account when using this method.

Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A_0 value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

NOTE: If necessary, the instrument must be dried off manually afterwards.

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11.5 Connecting

In order to ensure effective machine cleaning, disinfection and rinsing, the instrument must be connected in the washer and disinfectant.

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- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Joints, threads and glide surfaces must be treated locally with instrument oil and/or special lubricant. (only valid for 11014L)
- Dismantled medical devices must be assembled.
- Afterwards, a functional check must be carried out.

NOTE: The oil used for this purpose must be suitable for the subsequent sterilization procedure (silicone-free and paraffin- or white oil-based).

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11.8 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003). The method must be selected taking into account the respective applicable national requirements and in consultation with the device and product manufacturers.

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Steam sterilization using the fractionated prevacuum procedure

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) at 132°C–137°C with a minimum exposure time of 3 minutes and a maximum of 18 minutes. This procedure is only suitable for thermostable instruments.

NOTE: Greased components must be disassembled for sterilization, so that steam penetration is ensured.

11.9 Limits of reprocessing

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

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13 Servicing and repair

Defective instruments must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original parts only.

14 Maintenance

Preventive maintenance is not essential. However, regular maintenance can contribute to identifying potential problems before they become serious, enhancing the instrument's reliability and extending its useful operating life.

Maintenance services can be requested from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year.

15 Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this instrument only subject to the conditions that all assembly, system expansions, readjustments, modifications, or repairs to same have been performed by persons authorized by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

16 Warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see the section 'Subsidiaries'), even during the warranty period.

Opening of the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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17 Directive compliance

In accordance with the Medical Device Directive (MDD):

11272 VH/VHU:

This medical device belongs to Class II a

This medical device bears a CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

18 Important information

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and sterilized/disinfected before they are sent for repair. We reserve the right to return contaminated instruments/devices to the sender.

Repairs, modifications or expansions which are not performed by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

KARL STORZ gives no warranty on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

19 Disposal

This device has been marked in accordance with the European Directive 2012/19/EC on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the device as electronic scrap. Please ask KARL STORZ SE & Co. KG, a KARL STORZ subsidiary, or your authorized dealer about the collection point in your area. Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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20 Electromagnetic Compatibility user information

The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is suitable for use within close proximity to an active HF electro-surgical device in professional healthcare facility environment.

Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

WARNING:

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING:

The use of an accessory or cable with the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ other than those specified in this manual may result in increased emissions or decreased immunity of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™. When using an accessory or cable with the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ other than those specified in this manual, it becomes the responsibility of the user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ to determine compliance with IEC 60601-1-2 when using this item.

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 WARNING:

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

NOTE: The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If it is used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio communication service. The user might need to take mitigation measures, such as relocating or re-orienting equipment.

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Table 1 Guidelines and manufacturer's declaration – electromagnetic immunity			
The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is intended for use in an electromagnetic environment specified below. The user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ should make sure that it is used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1kV ± 1 kV for input and output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage Dip:</u> Dip to 0% for 1 cycles @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100% for 250/300 cycles	<u>Voltage Dip:</u> Dip to 0% for 1 cycles @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100% for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ requires continued operation during power supply interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50Hz/60Hz	30 A/m at 50Hz/60Hz	If image distortion occurs, it may be necessary to position the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	

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Table 3

**Guidance and manufacturer's declaration – electromagnetic immunity –
for equipment and systems that are not life-supporting**

The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is designed for operation in an electromagnetic environment as defined below. The user of the device should make sure that it is used in such an environment.

Interference immunity tests	EN/IEC 60601 test level	Conformity level	Electromagnetic environment – guidelines
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the $d = 1.2 \sqrt{P}$ 80 MHz to < 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects and persons.

^a) Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to stationary transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is used exceeds the above compliance levels, the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ should be monitored to ensure proper function. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

^b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4 Emission class and group		
Guidelines and manufacturer's declaration – electromagnetic emissions		
The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is intended for use in the electromagnetic environment specified below. The customer or user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ should ensure that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ uses RF energy only for its internal function. The customer or user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ should make sure that it is used in such an environment. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class-A	The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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Table 5 Recommended separation distances between portable and mobile RF telecommunication devices and the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™			
<p>The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ as recommended below, according to the maximum output power of the communications equipment.</p>			
Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>NOTE: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.</p> <p>NOTE: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.</p>			

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Thank you for your expression of confidence in the KARL STORZ brand. Like all of our other products, this medical device is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality piece of equipment from KARL STORZ.

This instruction manual enables you to properly install, connect, and operate this medical device.

- Failure to observe and follow this instruction manual and the instruction manuals of medical devices used in combination can result in injury to patients, users and third parties as well as damage to the medical device.
- Please read all relevant instruction manuals and always follow the instructions given. Check the functioning of the medical devices used in combination.
- Take note of all the instructions on the device.
- Keep this instruction manual available for ready reference in a clearly visible location near the device.

KARL STORZ reserves the right to make engineering modifications, in particular in the context of ongoing development and improvement of our devices.

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1	Intended use..... 3
1.1	Intended use..... 3
1.2	Indications..... 3
1.3	Contraindications..... 3
1.4	Target groups..... 4
1.5	User qualification..... 4
1.6	Patient population..... 4
1.7	Intended application areas..... 4
1.8	Intended conditions of use..... 4
1.9	User profile of physician and medical support staff..... 4
2	Safety..... 5
3	Description..... 7
3.1	Device overview..... 7
3.2	Symbols employed..... 8
3.2.1	Safety instructions..... 8
3.2.2	Pictograms/packaging symbols..... 9
4	Initial operation..... 10
4.1	Unpacking..... 10
4.2	Scope of supply..... 10
4.3	Required accessories..... 10
4.4	Preparing the device for use..... 11
5	Operation..... 12
5.1	Operating the user control buttons..... 12
5.2	Using accessories..... 13
5.2.1	Mechanical accessories..... 13
5.2.2	HF accessories..... 13
6	Reprocessing..... 15
6.1	Overview table..... 15
6.2	General safety instructions..... 16
6.3	Accessories..... 17
6.4	Preparation for cleaning and disinfection..... 17
6.5	Manual precleaning..... 17
6.6	No use of ultrasound..... 18
6.7	Machine cleaning and disinfection..... 18
6.8	Assembly, inspection and care..... 18

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6.9	Packaging systems.....	19
6.10	Sterilization	19
6.11	Limits of reprocessing.....	19
6.12	Care and handling of the video connecting cable	20
7	Service, repair and disposal	21
7.1	Warranty	21
7.2	Maintenance	21
7.3	Repairing device parts	21
7.4	Repair program.....	22
7.5	Infection prevention.....	22
7.6	Disposal.....	22
8	Technical description	23
8.1	Technical data.....	23
8.2	Operating/storage conditions.....	23
8.3	Directive compliance.....	23
9	Accessories.....	24
10	Electromagnetic compatibility user information	25
11	Subsidiaries.....	43

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1 Intended use

1.1 Intended use

Mediastinoscopes serve to create access for diagnostics and treatment in the middle mediastinum. Mediastinoscopes are designed for transient use in surgically invasive interventions.

1.2 Indications

- Non-small cell lung cancer/metastases
- Small cell lung cancer/metastases
- Staging of mediastinal lymph node metastases
- Interstitial lung disease
- Esophagectomy

1.3 Contraindications

There are no known contraindications relating directly to the device. The video mediastinoscope must not be used if, in the opinion of the responsible physician, the patient's health is endangered through its use, e.g., due to the patient's general condition or if the endoscopic method as such is contraindicated.

The following contraindications also apply:

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1.4 Target groups

This instruction manual is intended for physicians, medical support staff and sterilization assistants.

1.5 User qualification

Only physicians and medical support staff with a relevant specialist qualification and who have received training on the device may use the video mediastinoscope. Only technically suitable persons may provide training.

1.6 Patient population

Gender: irrelevant
Age: irrelevant
Weight: irrelevant
State of health: suitable for endoscopic interventions, in the physician's judgement

1.7 Intended application areas

Application area of mediastinoscopy, see intended use

1.8 Intended conditions of use

The video mediastinoscope may only be used in hospitals in suitable ambient conditions. For ambient conditions see the technical data.

1.9 User profile of physician and medical support staff

- Recognized medical skills in the field of application (specialist physician, qualified medical staff).
- Adequate powers of comprehension to rationally assess the surgical situation in hand.
- Adequate language skills in at least one of the languages used in the instruction manual
- Be thoroughly trained in the operation and use of the device.
- Knowledge of the contents of the instruction manual.
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic).

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- *Inspect reusable instruments for visible contamination. Reprocess contaminated instruments.*
- *Reprocess reusable instruments before and after each use. When doing so, use validated procedures.*

Combination of system components

Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product.

- *Please read all relevant manuals and always follow the instructions given precisely.*
- *Check the functioning of the products used in combination.*

The video mediastinoscope can only be used as a system in combination with other instruments and devices from KARL STORZ.

- *Suitable system components are listed in Section 4.3 'Required accessories'.*
- *Combination with non-suitable instruments and devices can result in uncontrolled conduct and injury to the patient, user or third parties.*

Improper use

Incorrect application of medical instruments poses a risk of injury for patients.

- *Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.*

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Risks due to damaged parts

Incorrectly assembled and damaged instruments can lead to injuries to the patient.

- *Check the following properties of instruments before and after every use to ensure that they:*
 - are free from damage
 - are in full working order,
 - have no unintentional rough surfaces,
 - have no sharp corners, burred edges,
 - have no projecting parts,
 - are complete
- *Do not leave missing or broken-off components inside the patient.*

Risks from incorrect maintenance and installation

Maintenance and installation errors endanger the safety of physicians, medical support staff and patients. The instrument doesn't work reliably or works with only with limited output.

- *Only staff authorized by KARL STORZ may perform the following work:*
 - Assembly
 - Expansion
 - Readjustments
 - Modifications
 - Repairs
- *The electrical installation of the location in which the medical device is connected and operated must conform to the applicable statutes and standards*

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
- ① Video connecting cable, fixed in place
- ② User control button: for scrolling (downwards) through the menu or activating a camera function
- ③ User control button: for accessing the menu/ making a selection
- ④ User control button: for scrolling (upwards) through the menu or activating a camera function

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3.2 Symbols employed

3.2.1 Safety Instructions

HAZARD SYMBOL: This is a hazard symbol. It warns of the risk of injury. To prevent injuries, follow all the measures marked with the hazard symbol.

Signal words:

WARNING! designates a possible imminent risk. If this is not avoided, it could lead to death or serious injuries.

CAUTION! designates a possible imminent risk. If this is not avoided, it could lead to minor injuries.

ATTENTION! designates a possibly harmful situation. If this is not avoided, the medical devices could be damaged.

NOTE! indicates special information on the operation of the device or important information.

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**3.2.2 Pictograms/
packaging symbols**

	Manufacturer
	Consult instructions for use
	Follow instructions for use
	Catalogue number
	Serial Number
	Number of products in the product packaging
	Keep dry
	Applied part type BF
	This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).
	CE mark
	Federal (USA) law restricts this device to sale by or on the order of a physician.

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- Kompatibler HD-Monitor.

4 Initial operation

NOTE: The video mediastinoscope including the equipment connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations and standards.

4.1 Unpacking

Carefully remove the video mediastinoscope and the accessories from the packaging.

Check for missing items and for evidence of shipping damage.

If necessary: Notify the manufacturer or the supplier.

TIP: If possible, retain the original packing materials for later use; these can be used if the device has to be transported.

4.2 Scope of supply

- Video mediastinoscope 10973 HD
- Instruction manual 96086006 DER

4.3 Required accessories

ATTENTION! Non-compatible accessories can lead to malfunctioning or damage to the video mediastinoscope.

- Only use the compatible accessories designated by KARL STORZ.

The video mediastinoscope 10973 HD requires additional accessories for normal use:

- IMAGE1 S CONNECT® TC 200
- IMAGE1 S X-LINK TC 301

NOTE: These devices have their own instruction manual (Art. no. 96206286DER)

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5 Operating instructions

! WARNING!

Risk of injury and risk of damaging the products: The highly concentrated light energy can cause the distal end, the light port, adjacent components and tissue in front of the light emission window to heat up. This can cause burns to the patient, user and operating accessories.

- Avoid direct tissue contact with the distal end.
- Do not rest the video mediastinoscope on the patient or in direct contact with surgical accessories.

! WARNING!

Risk of injury: Optical radiation poses a risk of injury to eyes.

- Never look into the free end of either a connected light cable or endoscope.
- Always wear standardized protective equipment when using lasers.

5.1 Operating the user control buttons

② Scroll down:

Scroll downwards through menus, options, etc. The button also be assigned a programmable function/feature.

③ Menu/selection:

Call up menu. Once the menu is displayed, this button can be used to select the highlighted option (identical to Enter on a computer keyboard). Hold this button down for seven seconds to call up the 'Options' menu.

④ Scroll up:

Scroll upwards through menus, options, etc. This button also be assigned a programmable function/feature.

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Further information on the programing and operation of camera head button can be found in the instruction manual for KARL STORZ IMAGE1 S™ (Art. no.: 96206286DER).

5.2 Using accessories

The spatula with an opening at the side enables the use of accessories such as forceps, scissors and HF coagulation instruments.

5.2.1 Mechanical accessories

WARNING!

Overloading by exerting too much force may cause the medical device to break, bend, and malfunction, and consequently injure the patient or user.

- Do not overload the instruments.
- Do not bend bent instruments back to their original positions.

5.2.2 HF accessories

WARNING!

Risk of injury and risk of damage to the products: If applied parts are used outside of the field of vision, there is a risk that tissue and accessories could be damaged unintentionally.

- Always hold the applied parts of the active electrode, laser and other instruments that transmit energy pointed at the target in the field of vision during use.

WARNING!

Risk of injury: Patient leakage currents may accumulate if endoscopes and endotherapy devices supplied with energy are used together. This is particularly important if an endoscope with a Type CF applied part is used. In this case, an endotherapy device with a Type CF applied part must be used in order to minimize the total patient leakage current.

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WARNING!

Risk of injury: The electromagnetic interference fields generated during the use of HF surgery can impair the functioning of other electronic devices. In the case of patients with cardiac pacemakers or other active implants, there is the risk of disturbances or damage.

- *If in doubt, seek expert advice from KARL STORZ.*

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- This reprocessing step is a key part of the validated process
- This reprocessing step may be carried out additionally in order to aid cleaning
- X On the grounds of potential damage, this reprocessing step must NEVER be performed
No symbol: This reprocessing step has neither been validated in respect of efficacy nor in respect of material compatibility by KARL STORZ.

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6.2 General safety instructions

WARNING!

Incorrectly reprocessed instruments expose patients, users and third parties to a risk of infection as well as the danger of an instrument malfunction. The product-accompanying documents must be observed.

- An overview of the procedures validated by KARL STORZ can be found in the manual 'Cleaning, disinfection, care and sterilization of KARL STORZ instruments' (Art. no. 96216003 DER).
- This manual can be downloaded or requested by visiting www.karlstorz.com.

WARNING!

The instruments are not sterile when delivered. The use of non-sterile instruments poses a risk of infection for patients, users and third parties.

- Inspect reusable instruments for visible contamination. Reprocess contaminated instruments.
- Reprocess reusable instruments before and after each use. When doing so, use validated procedures.

WARNING!

Risk of infection: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety as well as other country-specific laws must be observed.

WARNING!

If CJD (Creutzfeldt Jakob disease) is suspected or indicated, KARL STORZ recommends disposing of the product in accordance with national regulations and not using it again.

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ATTENTION: Using chemicals not approved by KARL STORZ may damage medical devices.

- Only use chemicals approved by KARL STORZ for reprocessing. A list of approved chemicals can be requested from hygiene@karlstorz.com.

ATTENTION: Prolonged immersion and incorrect concentration of the cleaning/disinfection solutions may result in damage.

- When preparing and using solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life.
- The microbiological range of action of the chemicals used must be taken into account.

NOTE: National laws and regulations must be observed.

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- Brush: 27652

6.4 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

6.5 Manual precleaning Brushing the surfaces

Any visible contamination or heavy soiling must be removed by cleaning the surfaces under cold running water with the aid of a brush or sponge

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6.6 No use of ultrasound

For technical reasons, the medical device is not suitable for ultrasound treatment.

6.7 Machine cleaning and disinfection

NOTE: In order to protect the contacts of the video connecting cable connector, machine cleaning and disinfection is only possible if the protective cap is fitted (see Fig.).

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

• Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A_0 value must be taken into account when using this method. The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

NOTE: If necessary, the instrument must be dried off afterwards by hand.

6.8 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Afterwards, a functional check must be carried out.

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NOTE: Use items from the catalog 'Care, Sterilization and Storage Techniques' during care procedures.

6.9 Packaging systems

Only standardized and approved packaging materials and systems may be used (EN 868 Parts 2-10, EN ISO 11607 Parts 1 + 2, DIN 58953).

6.10 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The method must be selected taking into account the respective applicable national requirements and in consultation with the device manufacturers.

The following sterilization procedures have been validated and approved by KARL STORZ for this medical device:

• Steam sterilization using the fractionated pre-vacuum procedure

NOTE: Sterilize greased components in a disassembled state.

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) with the following parameters:

- 132°C with a minimum exposure time of 4 minutes and maximum of 18 minutes.

NOTE: This validation is applicable for reusable tube sets which are open at both ends and with a maximum length of 200 cm.

- 134°C-137°C with a minimum exposure time of 3 minutes and maximum of 18 minutes.

6.11 Limits of reprocessing

The end of the product's service life and correct functioning are determined by mechanical stress and chemical influences within the scope of reprocessing and application.

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6.12 Care and handling of the video connecting cable

To maximize the life of the video mediascope, please follow the guidelines given below for proper handling and care.

DO NOT PLUG THE VIDEO CONNECTING CABLE INTO THE PROCESSOR RECEPTACLE IF IT IS WET. The connector must always be completely dry and free of residue. Where necessary, remove any residual moisture with a low lint disposable cloth before plugging in the unit.

NEVER PICK UP THE VIDEO MEDIASTINOSCOPE BY ITS CABLE! Always pick up the set by the instrument handle first, followed by the cable and connector. The cable contains tiny electronic conductors and should be treated gently and carefully.

PROTECT THE VIDEO MEDIASTINOSCOPE FROM DIRECT SUNLIGHT AND EXCESSIVE HEAT.

HANDLE THE SYSTEM WITH CARE, avoiding extreme impacts to the device. When transporting, use the original packaging to prevent damage to the device.

DO NOT STORE THE VIDEO-MEDIASTINOSCOPE WITH CABLE COILED TIGHTLY OR GATHERED INTO FOLDS. Store the cable in loops at least 15 cm in size.

DO NOT RAPIDLY PULL THE CABLE WHEN UNCOILING. Harmful 'kinks' can result. Always allow the cable to uncoil gently and loosely.

NOTE: We always recommend drying with sterile compressed air after manual reprocessing and, if necessary, also after machine reprocessing.

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7 Service, repair and disposal

7.1 Guarantee

The warranty conditions provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your competent KARL STORZ subsidiary (see the section entitled 'Subsidiaries'), even during the warranty period.

Opening of the equipment or performance of any repairs or modifications to the device by unauthorized persons shall relieve KARL STORZ of any liability for its performance. Any such opening, repair or modification performed during the guarantee period shall void all warranty.

7.2 Maintenance

No further repair measures are required beyond reprocessing.

The legal requirements as per the Medical Devices Operator Ordinance (MPBetreibV) for devices to be used in combination apply (see Section 4.4). The relevant warnings in the instruction manuals must be observed.

7.3 Repairing device parts

WARNING!

Danger of electric shock. Do not open the device!

- Refer servicing to the manufacturer or the manufacturer's authorized service personnel (see Section 3 of the Medical Device Operator Ordinance).
- Removal of covers by unauthorized personnel will void the equipment's warranty.

Only personnel authorized by the manufacturer may repair device parts. The operations described in this manual are exempt from this rule.

The authorized personnel may only repair device parts which the manufacturer designates as repairable.

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Upon request, the manufacturer will provide authorized personnel with information on repairs, e.g., circuit diagrams, spare parts lists, setting instructions.

The fact of having such documents at one's disposal does not correspond to authorization.

7.4 Repair program

For this medical device, individual repair is necessary. To bridge the repair period, you will receive a device on loan, which must then be returned as soon as you receive the repaired medical device.

In the case of repairs, please contact your local KARL STORZ subsidiary (see section 'Subsidiaries') or authorized dealer.

In Germany, you can refer repairs direct to:

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Take-off Gewerbepark 83
78579 Neuhausen

7.5 Infection prevention

To prevent infections, contaminated devices and instruments must not be returned.

Decontaminate instruments and devices before returning. KARL STORZ returns contaminated instruments to the sender.

7.6 Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the device as electronic scrap.

Please ask KARL STORZ, a KARL STORZ subsidiary (see the section 'Subsidiaries'), or your authorized dealer about the collection point in your area.

Within the scope of application of this Directive, KARL STORZ is responsible for the proper disposal of this device.

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8 Technical description

8.1 Technical data

Direction of view	30°
Angular aperture	50°
Spatula working length	198 mm
Distal outer diameter of spatula	16.85 mm
Distal inner diameter of spatula	14.85 mm
Illumination	Integrated LED source
Operating temperature	5 – 30 °C
Weight with cable	590 g
Dimensions max. without cable in mm (W × H × D)	198 × 144 × 25
Cable length	2.90 m

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This medical device bears a CE mark in accordance with MDD 93/42/EEC.

A code number after the CE mark indicates the responsible notified body.

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9 Accessories
MEDIAFIT forceps

MEDIAFIT Sponge and Dissecting Forceps acc. to MAASSEN, fenestrated

MEDIAFIT Biopsy Forceps, oval jaws, size 8 x 16 mm

MEDIAFIT Biopsy Forceps, oval jaws, size 6 x 12 mm

MEDIAFIT Biopsy Forceps, spoon-shaped jaws

MEDIAFIT Grasping and Dissecting Forceps, upturned jaws

MEDIAFIT Grasping and Dissecting Forceps, long and fenestrated jaws

MEDIAFIT clip forceps, jaws curved to right

MEDIAFIT clip forceps, jaws curved to left

Coagulation Suction Tube

MEDIAFIT Outer Sheath with cleaning connection, diameter 5 mm, length 21 cm

HF instrument for mediastinoscopy

HF instrument for mediastinoscopy

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10 Electromagnetic compatibility user information

The video mediastinoscope 10973 HD together with the IMAGE1 S™ is suitable for use in the direct proximity of an active HF electrosurgical unit in professional health institutions.

Professional health institutions include medical practices, dental practices, facilities for restricted care, independent surgical centers, independent birthing centers, several treatment facilities, hospitals (emergency wards, patient rooms, intensive care wards, operating theaters, outside of the HF-shielded room of an ME system for magnetic resonance imaging).

WARNING!

The use of this device next to or with other devices should be avoided as this can result in incorrect operation. If such use is necessary, this and the other devices should be monitored to ensure that they are functioning correctly.

WARNING!

The use of an accessory or cables for the video mediastinoscope 10973 HD with the IMAGE1 S™ other than those specified in this manual may result in increased emissions or decreased immunity of the video mediastinoscope 10973 HD together with the IMAGE1 S™. When using an accessory or cables with the video mediastinoscope 10973 HD together with the IMAGE1 S™ other than those specified in this manual, it becomes the responsibility of the user to determine compliance of the video mediastinoscope 10973 HD and IMAGE1 S™ with IEC 60601-1-2 when using this item.

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 WARNING!

Portable HF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the video mediastinoscope 10973 HD in conjunction with the IMAGE1 S™, including cables specified by the manufacturer. Otherwise, performance may be impaired.

NOTE: Based on the device's emissions properties, it is suitable for use in industrial areas and hospitals (CISPR 11 Class A). If the device is used in a domestic environment (for which CISPR 11 Class B is normally required), the device may not offer sufficient protection for radio transmission operation. It may be necessary for the user to take measures, e.g., choose a different location for the device or realign it.

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Table 1 Guidelines and manufacturer's declaration – electromagnetic immunity			
The video mediastinoscope 10973 HD together with the IMAGE1 S™ is intended for use in an electromagnetic environment specified below. The user of the video mediastinoscope 10973 HD together with the IMAGE1 S™ should make sure that it is used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1kV ± 1 kV for input and output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage Dip:</u> Dip to 0% for 1 cycles @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100% for 250/300 cycles	<u>Voltage Dip:</u> Dip to 0% for 1 cycles @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100% for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the video mediastinoscope 10973 HD together with the IMAGE1 S™ requires continued operation during power supply interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50Hz/60Hz	30 A/m at 50Hz/60Hz	If image distortion occurs, it may be necessary to position the video mediastinoscope 10973 HD together with the IMAGE1 S™ further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	

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Table 3 Guidance and manufacturer's declaration – electromagnetic immunity – for equipment and systems that are not life-supporting			
The Video-Mediastinoscope 10973 HD together with the IMAGE1 S™ is designed for operation in an electromagnetic environment as defined below. The user of the device should make sure that it is used in such an environment.			
Interference immunity tests	EN/IEC 60601 test level	Conformity level	Electromagnetic environment – guidelines
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the Video-Mediastinoscope 10973 HD together with the IMAGE1 S™, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the $d = 1.2 \sqrt{P} \quad 80 \text{ MHz to } < 800 \text{ MHz}$ $d = 2.3 \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects and persons.			
^{a)} Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to stationary transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the video mediastinoscope 10973 HD together with the IMAGE1 S™ is used exceeds the above compliance levels, the video mediastinoscope 10973 HD together with the IMAGE1 S™ should be monitored to ensure proper function. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.			
^{b)} Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

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Table 4 Guidelines and manufacturer's declaration – electromagnetic emissions		
The video mediastinoscope 10973 HD together with the IMAGE1 S™ is intended for use in the electromagnetic environment specified below. The customer or user of the video mediastinoscope 10973 HD together with the IMAGE1 S™ should ensure that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The video mediastinoscope 10973 HD together with the IMAGE1 S™ uses RF energy only for its internal function. The customer or user of the flexible video mediastinoscope 10973 HD together with the IMAGE1 S™ should make sure that it is used in such an environment. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class A	The video mediastinoscope 10973 HD together with the IMAGE1 S™ is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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Table 5
Recommended separation distances between portable and mobile RF telecommunication devices and the Video-Mediastinoscope 10973 HD together with the IMAGE1 S™

The video mediastinoscope 10973 HD together with the IMAGE1 S™ is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the video mediastinoscope 10973 HD together with the IMAGE1 S™ can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the video mediastinoscope 10973 HD together with the IMAGE1 S™ as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.

NOTE: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.

HF surgical instrument compatibility

This equipment has been evaluated for compatibility with high-frequency surgical equipment.

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KARL STORZ Endoscopy Canada Ltd.
1717 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzcuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Av. Ejercito Nacional No. 453 Piso 2,
Colonia Granada, Alcadia Miguel Hidalgo
C.P. 11520 Ciudad de México
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA – Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Stamveien 1
1483 Hagan, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
127 39 Skärholmen, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se
KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kęstucio st. 59 / Lenktoji st. 27
08124 Vilnius, Lithuania
Phone: +370 5 272 0448
Mobile: +370 685 67 000
E-Mail: info-lt-iv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopie France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Croatia d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Cesta v Gorice 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Interchange 21 Tower, Level 33
399 Sukhumvit Road,
North Klongtoey, Wattana,
10110 Bangkok, Thailand
Phone: +84 28 3823 8000
Fax: +84 28 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, MPlaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd.
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch
Room 701A Building 5 & Room 501 Building 7,
No. 3000 Longdong Avenue, Pilot Free Trade Zone,
201203, Shanghai, People's Republic of China
Phone: +86 21 60339888, Fax: +86 21 60339808
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch
Room 803-805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch
Room 02B & 03 & 04A, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 87321281, Fax: +86 20 87321286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

KARL STORZ Endoscopy Australia Pty. Ltd.
68 Waterloo Road, Macquarie Park NSW 2113
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: karlstorz@karlstorz.com.au
www.karlstorz.com

KARL STORZ SE & Co. KG

Dr.-Karl-Storz-Straße 34
78532 Tuttlingen

Postfach 230
78503 Tuttlingen

Germany

Telefon: +49 7461 708-0
Telefax: +49 7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com