**REF****T401160, T401161, T401162, & T401163**

*Hercules™ 3/360 Universal Stabilizer Arms and Herc Flex™ /
Herc Flex™ 360 Stabilizer Arms*

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1025 000

This page is intentionally left blank.

Hercules™ 3/360
Universal Stabilizer Arms
and
Herc Flex™ / Herc Flex™ 360
Stabilizer Arms

Table of Contents

English.....	5
Česky.....	17
Dansk	29
Nederlands	41
Suomi	53
Français (France)	65
Deutsch	77
Magyar.....	89
Italiano	101
Norsk	113
Polski	125
Português (Brasil).....	137
Русский язык	149
简体中文	163
Español (España).....	175
Svenska.....	187
Türkçe.....	199
Slovensky	211
繁體中文	223
Graphical Symbols & Definitions	235

This page is intentionally left blank.

English

Instructions for Use

Hercules™ 3/360 Universal Stabilizer Arms and Herc Flex™ / Herc Flex™ 360 Stabilizer Arms

Catalogue Number T401160, T401161, T401162, and T401163

CAUTION: Carefully read all instructions prior to use. Observe all contraindications, warnings, and precautions noted in these directions. Failure to do so may result in patient complications.

PRODUCT DESCRIPTION

The Hercules™ Universal Stabilizer Arms and the Herc Flex™ Stabilizer Arms are reusable and are supplied non-sterile. They consist of a distal connection mechanism (either standard Terumo® Ball Quick Connect or Terumo® Shaft Quick Connect), an articulating arm with links, a handle for providing cable tension, and a mechanism to mount the device to a sternal retractor. Refer to Table 1 for available versions of the Hercules™ Universal Stabilizer Arms and Herc Flex™ Stabilizer Arms. See Figure 1 and 2 for diagrams of the Hercules™ 3 Universal Stabilizer Arm and the Herc Flex™ Stabilizer Arm.

WARNING Device must be cleaned, lubricated, and sterilized prior to first use and before subsequent uses.


Figure 1 – Hercules™ 3 Universal Stabilizer Arm


Figure 2 – Herc Flex™ Stabilizer Arm

Table 1 – Hercules™ and Herc Flex™ Product Descriptions

Model	Type of distal connection mechanism	# of Links
Hercules™ 360 Universal Stabilizer Arm – Product Reference T401160	Terumo Ball Quick Connect	19
Hercules™ 3 Universal Stabilizer Arm – Product Reference T401161	Terumo Shaft Quick Connect	18
Herc Flex™ 360 Stabilizer Arm – Product Reference T401162	Terumo Ball Quick Connect	24
Herc Flex™ Stabilizer Arm – Product Reference T401163	Terumo Shaft Quick Connect	23

Note: The packaging configuration consists of one Hercules™ Universal Stabilizer Arm or Herc Flex™ Stabilizer Arm unit (1) per box.

A wide variety of single-use attachments are available from Terumo Cardiovascular Group intended for use with the Hercules™ Universal Stabilizer Arms and Herc Flex™ Stabilizer Arms.

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See Table 2 for potential single-use attachments for use with Hercules™ 360 Universal Stabilizer Arm and Herc Flex™ 360 Stabilizer Arm (Product References: T401160 and T401162).

Table 2 - Hercules™ 360 Universal Stabilizer Arm and Herc Flex™ 360 Stabilizer Arm Attachments

Product Reference	Description
T401261U	Atlas™ 360 Positioner Attachment
T401231U	Titan™ 360 Stabilizer Attachment
T401232U	Titan™ Flex 360 Stabilizer Attachment

See Table 3 for potential single-use attachments for use with Hercules™ 3 Universal Stabilizer Arm and Herc Flex™ Stabilizer Arm. (Product References: T401161 and T401163).

Table 3 - Hercules™ 3 Universal Stabilizer Arm and Herc Flex™ Stabilizer Arm Attachments

Product Reference	Description
T401261	Atlas™ Positioner Attachment
T401231	Titan™ Stabilizer Attachment
T401232	Titan™ Flex Stabilizer Attachment
T401231LS	Titan™ Stabilizer Attachment, Long Shaft
T401244	StableSoft™ II Stabilizer
T401245	StableSoft™ II Long Shaft
T401209S	Assistant with StableSoft™ 3L
T401210S	Assistant with StableSoft™ 3S
T401211S	Assistant with StableSoft™ 2L
T401212S	Assistant with StableSoft™ 2S
T401213S	Assistant with StableSoft™ 3L Flex
T401214S	Assistant with StableSoft™ 3S Flex

INTENDED USE

The Hercules™ Universal Stabilizer Arms and the Herc Flex™ Stabilizer Arms are intended for use in cardiothoracic surgery to support stabilization and positioning of various anatomical structures during various procedures.

CONTRAINDICATIONS

The Hercules™ Universal Stabilizer Arms and the Herc Flex™ Stabilizer Arms are contraindicated for use in patients with local or systemic infection.

WARNINGS AND PRECAUTIONS

WARNINGS

- The Hercules™ Universal Stabilizer Arms and the Herc Flex™ Stabilizer Arms should only be operated by experienced physicians properly trained in beating heart surgery, and the assistance of an experienced anesthetist is required.
- Proper surgical procedure and technique is the responsibility of the attending medical professional (physician). Each surgeon must evaluate the appropriate use of this device, case by case, based upon medical training, experience, and the type of surgical procedure employed.
- Excessive torque might cause cable to fray, snap, or break. If the Hercules™ Universal Stabilizer Arms or the Herc Flex™ Stabilizer Arms do not hold position, they should not be used.
- DO NOT reposition the arm when the Hercules™ Universal Stabilizer Arms or the Herc Flex™ Stabilizer Arms are tightened. It will cause the cable in the Hercules™ Universal Stabilizer Arms or the Herc Flex™ Stabilizer Arms to fray and possibly break.
- Heavy or improper use of devices may cause damage requiring replacement of unit.
- Visually inspect cable for signs of fraying. Do not use fingers to check cable for signs of fraying.
- This unit is not to be reprocessed using hydrogen peroxide gas plasma sterilization, such as STERRAD® sterilizer, which can cause damage or breakage to the nylon bearing component. Please contact the manufacturer of the automated washer/disinfector prior to use to determine if the washer/disinfector processing parameters are compatible with the Hercules™ Universal Stabilizer Arms or the Herc Flex™ Stabilizer Arms.
- Strictly follow the manufacturer's instructions for concentration of the detergent solution and proper use to avoid high acidity or alkaline pH balances, which may cause corrosion and result in breakage.

PRECAUTION

- Cable fatigue and wear can be easily detected, and failure prevented, by examining the cable before each use. Prior to beginning each procedure, while the arm is fully loosened, look at the arm to determine if the cable shows signs of fraying (see Figure 4 for fraying). Do not use if signs of fraying are present, as a cable may suddenly fail during a procedure.

- Do not use chlorine based solutions (such as bleach) on the Hercules™ Universal Stabilizer Arms or the Herc Flex™ Stabilizer Arms, as this can cause corrosion.

CAUTION

- United States Federal Law restricts this device to sale by or on the order of a physician or other licensed practitioner.

INSTRUCTIONS FOR USE

I. PRE-USE INSPECTION

1. Prior to use, the unit should be fully loosened and visually inspected. The entire cable should be inspected by examining the spaces between links. No signs of wear or fraying of the cable should be present. If any wear is found, the device should not be used. See Figure 3 and Figure 4 for examples of normal and fraying cables.


Figure 3 – Appearance of a normal cable

Figure 4 – Appearance of a fraying cable – DO NOT USE

2. Verify function of the device by tightening and loosening the cable by turning the main handle.

II. PROCEDURE

WARNING

Device must be cleaned and sterilized prior to first use and before subsequent uses (Refer to chapter on REPROCESSING INSTRUCTIONS).

1. Attach an accessory device to the distal connection mechanism. There are two different types of attachment interfaces that can be used.
 - a. For connection of a standard Terumo SHAFT QUICK CONNECT attachment, slide the sleeve of the Stabilizer Arm forward and insert the shaft of the attachment into the quick connect with the detent on the shaft (of the attachment) facing the flat area with logo on the sleeve of the Stabilizer Arm (see Figure 5). Return the sleeve to its original position to lock the accessory device into the quick connect. Check for secure connection.



Figure 5 – Attachment of accessory device to the Hercules™ 3 Universal Stabilizer Arm or Herc Flex™ Stabilizer Arm (Product References T401161 and T401163)

- b. For connection of a Terumo BALL QUICK CONNECT attachment, grasp the sleeve of the Stabilizer Arm and press the ball of the attachment into the collet of the Stabilizer Arm (see Figure 6). Check for secure connection.

NOTE: Ensure that the sleeve is grasped during attachment insertion. Failure to grasp the sleeve may result in an inability to insert the attachment.


Figure 6 – Attachment of accessory device to the Hercules™ 360 Universal Stabilizer Arm or the Herc Flex™ 360 Stabilizer Arm (Product References T401160 and T401162)

2. Attach the retractor clamp of the Stabilizer Arm to the retractor by sliding the rear jaw backward and placing the retractor clamp onto the retractor arm. Allow the rear jaw to spring closed against the retractor arm before rotating the retractor clamp handle clockwise to clamp the Stabilizer Arm.
3. Position accessory device as required. Turn the main handle clockwise until the gap between the compression washers closes, as depicted in Figure 7. The closing of the gap between the washers indicates that the recommended cable tension is reached, as shown in the following picture.

CAUTION Do not reposition the arm without loosening the main handle.

Compression Washers


Figure 7 – Example of closed compression washers

III. POST-USE INSPECTION

After every use, the arm should be fully loosened. A visual inspection should be conducted between the last link and the clamp assembly for cable fraying. Cable fatigue and fraying typically occurs in the same location, at the back end of the arm nearest the retractor mounting clamp, making it very easy to detect early. There should be no signs of wear or fraying of the cable. See Figure 8 with a cable demonstrating fatigue and fraying. In

addition, check for broken or missing nylon washer depicted in Figure 9 and Figure 10.

If wear is observed, do not continue use, and contact your local Terumo Account Representative or call Customer Service for replacement. Again, the unit should never be used in the conditions shown in Figure 8, Figure 9, and Figure 10.

Figure 8 – Cable demonstrating fatigue and fraying


Figure 9 – Broken nylon washer


Figure 10 – Missing nylon washer

REPROCESSING INSTRUCTIONS

WARNING

Should cleaning and sterilization methods other than those recommended be performed, validation of the effectiveness of these methods must be performed, as well as verifying the effects on the Hercules™ Universal Stabilizer Arms or the Herc Flex™ Stabilizer Arms materials and functions.

CAUTION It is recommended that the device be reprocessed as soon as is reasonably practical following use (recommended within one hour).

I. POINT OF USE

Remove excess soil with disposable cloth/paper wipe.

II. PREPARATION FOR CLEANING

Ensure the main handle is fully loosened prior to cleaning, and that any attachments are disconnected from the distal end of the Hercules™ Universal Stabilizer Arms or the Herc Flex™ Stabilizer Arms. Perform cleaning and disinfection using either the manual or automated methods described below, followed by inspection, lubrication, and sterilization.

III. CLEANING AND DECONTAMINATION

MANUAL METHOD

1. Rinse thoroughly under reverse osmosis/deionized (RO/DI) water for a minimum of 1 minute to remove gross soil and debris, taking into account the gaps between the links.
2. Immerse the device in a 0.8% enzymatic detergent instrument cleaner for a minimum of 5 minutes. Terumo Cardiovascular Group recommends using ENZOL/Cidezyme, or a suitable equivalent. ENZOL/Cidezyme is a cleaner used to break down proteinaceous organic matter.
3. Actuate and tap the device to ensure penetration of the cleaner and release of trapped air bubbles.
4. Scrub each link and cable thoroughly using a soft brush or lint-free cloth under 0.8% enzymatic cleaner (ENZOL/Cidezyme). Remove all traces of blood and debris. Make sure all moving parts are cleaned thoroughly to prevent debris from interfering with movement.
5. Rinse the device thoroughly for a minimum of 1 minute in RO/DI water to remove dislodged surgical debris and the detergent solution, taking into account the gaps between the links.
6. Flush any hard-to-reach areas under water with a spray gun. Thoroughly flush the cable between each link for 5 seconds. Flush the distal connection mechanism for 15 seconds.
7. Sonicate device in 0.8% ENZOL/Cidezyme cleaning solution at 104°F/40°C for 10 minutes. Prior to sonication, actuate device in the cleaning solution to ensure penetration of the cleaner and release of trapped air bubbles.
8. Rinse the device thoroughly for a minimum of 1 minute in RO/DI water taking into account the gaps between the links.
9. Flush any hard-to-reach areas under water with a spray gun. Thoroughly flush the cable between each link for 5 seconds. Flush the distal connection mechanism for 15 seconds.
10. Dry device completely using compressed air.

CAUTION Do not use chlorine based solutions (such as bleach) on the Hercules™ Universal Stabilizer Arms or the Herc Flex™ Stabilizer Arms, as this can cause corrosion.

CAUTION Incomplete rinsing after cleaning can cause enzymatic cleaner to form a residue on the Hercules™ Universal Stabilizer Arms or the Herc Flex™ Stabilizer Arms that may result in binding of the links.

AUTOMATED CLEANING METHOD

MANUAL PRE-CLEAN

1. Rinse thoroughly under reverse osmosis/deionized (RO/DI) water for a minimum of 1 minute to remove gross soil and debris taking into account the gaps between the links.
2. Immerse the device in a 0.8% enzymatic detergent instrument cleaner for a minimum of 5 minutes. Terumo Cardiovascular Group recommends using ENZOL/Cidezyme, or a suitable equivalent. ENZOL/Cidezyme is a cleaner used to break down proteinaceous organic matter.
3. Actuate and tap the device to ensure penetration of the cleaner and release of trapped air bubbles.
4. Scrub each link and cable thoroughly using a soft brush or lint-free cloth under 0.8% enzymatic cleaner (ENZOL/Cidezyme). Remove all traces of blood and debris. Make sure all moving parts are cleaned thoroughly to prevent debris from interfering with movement.
5. Flush any hard-to-reach areas under water with a spray gun. Thoroughly flush the cable between each link for 5 seconds. Flush the distal connection mechanism for 15 seconds.
6. Sonicate device in 0.8% ENZOL/Cidezyme cleaning solution at 104°F/40°C for 10 minutes. Prior to sonication, actuate device in the cleaning solution to ensure penetration of the cleaner and release of trapped air bubbles.
7. Flush any hard-to-reach areas under water with a spray gun. Thoroughly flush the cable between each link for 5 seconds. Flush the distal connection mechanism for 15 seconds.

AUTOMATED PROCESS

1. Two minute pre-clean with cold tap water. Drain.
2. Five minute clean with deionized water at 55°C using 0.5 % of Neodisher Mediclean forte (Dr. Weigert, Hamburg). Drain.
3. Three minute rinsing and neutralization with deionized water. Drain.
4. Three minute thermal disinfection at 93°C with deionized water. Drain.
5. Ten minute drying.

CAUTION Do not use chlorine based solutions (such as bleach) on the Hercules™ Universal Stabilizer Arms or the Herc Flex™ Stabilizer Arms, as this can cause corrosion.

IV. INSPECTION, MAINTENANCE, AND TESTING

Visually inspect the entire Stabilizer Arm for cleanliness to ensure residue has been removed. If any tissue, blood, pus, or soil is still present, repeat cleaning process.

Ensure the arm is fully loosened, and inspect the entire cable for damage. There should be no signs of wear or fraying of the cable. If wear is observed, return to Terumo Cardiovascular Group for replacement. Cable fatigue and wear is a gradual event that occurs over many arm rigidifying cycles. As the cable begins to fatigue, the outer strands will begin to fray and break. The user can recognize onset of fatigue as continued revolutions of the handle do not rigidify the arm and an audible popping is heard. This sound indicates the individual cable wires breaking. Once cable fraying begins, the cable will eventually sever after numerous continued revolutions of the handle.

In the rare and unfortunate event where a cable completely severs and the articulating links disperse into the field, it is imperative to completely account for the distal assembly and all links per Table 1 of this IFU.

Check for smooth movement of retractor clamp handle, rear jaw, and main handle and articulating arm. Inspect the mated stainless steel and nylon thrust bearing washers near the main handle to verify all are intact. The instrument should be returned for replacement if any defects are found.

V. LUBRICATION

For optimal device performance and to maximize useful life, lubrication of the device using Steris Hinge-Free instrument lubricant (or equivalent) per manufacturer’s instructions is recommended prior to each device sterilization cycle. Ensure the handle is fully loosened to expose the screw prior to lubrication.

VI. STERILIZATION

The device should be sterilized by autoclave sterilization method using the following processing parameters shown in Table 4. Ensure the arm is not tightened during sterilization.

Table 4 – Sterilization Process Parameters

Description	Temperature	Exposure Time	Min. Drying Time
Gravity Displacement (wrapped)	132°C to 138°C	15-20 minutes	30 minutes
Prevacuum (wrapped)	132°C to 138°C	4-20 minutes	20 minutes
Prevacuum (wrapped)	134°C to 138°C	3-20 minutes	20 minutes

STORAGE

The Hercules™ Universal Stabilizer Arms and the Herc Flex™ Stabilizer Arms have the following storage conditions:

- A clean, dust-free location out of direct sunlight and free from risk of high temperature, excessive humidity, and water leakage.
- Do not store the device near chemicals.

LIMITATIONS ON REPROCESSING

Repeated reprocessing has no significant influence on the length of life of the device.

END OF LIFE

Inspect the device prior to each use for evidence of unusual wear and/or damage. If such wear is observed, do not use the device. This device has a shelf life of 48 months from the date of manufacture.

DISPOSAL

At end of life, dispose of product and packaging in accordance with hospital, administrative and/or local, state, federal, and international laws and regulations.

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TERUMO CORPORATION
44-1, 2-Chome, Hatagaya, Shibuya-ku
Tokyo 151-0072, Japan

Terumo® is a registered trademark of TERUMO CORPORATION
Hercules™, Herc Flex™, Titan™, StableSoft™, Assistant with StableSoft™,
Atlas™ are trademarks of Terumo Cardiovascular Systems Corporation
U.S. Patent No. 6,860,668

ENZOL® is a registered trademark of Advanced Sterilization Products,
Division of Ethicon US, LLC.

Hinge-Free® is a registered trademark of STERIS Corporation STERRAD®
is a registered trademark of Advanced Sterilization Products. Advanced
Sterilization Products (ASP), a Division of Ethicon, Inc. a Johnson &
Johnson company.

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AU Sponsor **TERUMO AUSTRALIA PTY LTD**
Country of Origin: U.S.A.