**225031 VAPR®
COOLPULSE® Curve
and
225032 VAPR®
COOLPULSE®
Curve XL
Suction Electrodes**

For use only with the VAPR VUE® Generator

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Manufactured for DePuy Synthes by:

Gyrus Medical PN0011069-AB

VAPR® COOLPULSE® CURVE AND VAPR® COOLPULSE® CURVE XL SUCTION ELECTRODES

Please read all information carefully prior to use, giving special attention to the User Manual provided with the VAPR VUE® Radiofrequency (RF) System. VAPR VUE System components are designed for use as a single system. Failure to follow instructions may lead to improper functioning of the device causing electrical or thermal injury.

SYSTEM DESCRIPTION

The VAPR VUE Radiofrequency (RF) System is designed for arthroscopic surgical procedures. The system consists of a high frequency electrosurgical generator and may include a reusable handpiece with connector cable, disposable electrodes, and alternative disposable electrodes with integrated connector cable. Also available are hand controls, a footswitch, and a wireless footswitch alternative. The components are designed and intended to be operated as a single unit. Use only VAPR® electrodes and accessories with the VAPR VUE RF system.

DEVICE DESCRIPTION

The VAPR Electrodes are soft tissue ablation and hemostasis/coagulation devices intended for use with the VAPR VUE RF System. The VAPR® COOLPULSE® Curve and COOLPULSE® Curve XL Electrodes suction capabilities enhance the efficiency of the electrode and extend the utility of the system by assisting in the removal of bubbles and debris created during activation within the operating site. VAPR electrodes are intended to be used by an orthopedic surgeon in a surgical setting.

INDICATIONS

The VAPR Electrodes for use with the VAPR VUE Radiofrequency System are intended for resection, ablation and excision of soft tissue, hemostasis of blood vessels and coagulation of soft tissue in patients requiring arthroscopic surgery.

CONTRAINDICATIONS

The VAPR VUE RF Electrosurgical System is contraindicated in:

- any non-arthroscopic surgical procedure,
- patients with heart pacemakers or other electronic device implants,
- patients for whom arthroscopic procedures are contraindicated for any reason.

WARNINGS

- The operator should be experienced in arthroscopic surgical techniques.
- Please refer to the user manual for step-by-step instructions regarding the assembly and initial system check of the VAPR VUE system.
- As with all electrosurgical devices, do not use in the presence of flammable anesthetics, oxidizing gases or other flammable substances; an electrosurgical device has the potential for providing a source of ignition.
- Do not use in the presence of flammable anesthetics or oxidizing gases (such as nitrous oxide (N₂O) and oxygen), or close to volatile solvents (such as ether or alcohol), as explosion may occur.
- Do not place instruments near, or in contact with flammable materials (such as gauze or surgical drapes). Instruments that are activated or hot from use may cause a fire.
- Carefully insert and withdraw electrodes to avoid possible damage to the device and/or injury to the patient or surgical personnel.
- Ensure that fluid inflow and outflow is adequate and that the electrode is activated only when surrounded by conductive irrigant solution (e.g., saline or Ringer's Lactate). The use of RF energy with inadequate irrigation may overheat the fluid enough to cause skin burns at or near the access site and may result in damage to the electrode tip assembly.
- Do not activate the instrument when not in contact with target tissue, as this may cause injuries due to capacitive coupling with other surgical equipment.
- **Do not activate the probe for unnecessary and prolonged periods, as irrigant overheating and tissue damage may result. Prolonged probe activation at maximum power while adjusting irrigant flow via the flow control valve to a low flow rate may cause the shaft and suction tubing to reach temperatures of up to 49°C. Do not touch the electrode tip when power is being applied.**
- During an electrosurgical procedure, do not allow the patient to come into direct contact with grounded metal objects. Inadvertent contact with the patient when the VAPR system is activated may result in burns.
- Observe extreme caution when using electrosurgery devices in close proximity to, or in direct contact with, any metal objects. Most arthroscopes and arthroscopic instruments are metal. Do not activate the electrode while any portion of the electrode tip is in contact with metal objects or instruments. Doing so may result in unintended injury to the patient or surgical personnel and/or damage to the electrode and/or other equipment.
- Unnecessary and prolonged activation of the electrode during and between tissue applications may result in unintended injury.
- Do not re-use or re-sterilize. The device is labeled as SINGLE USE. These products are for single-use only and are not designed to be re-used/re-sterilized. Reprocessing may lead to changes in material characteristics such as metallic corrosion and dulled edges, ceramic and plastic deformation, or splitting, which may impact the strength of the device and compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection. These risks may potentially affect patient safety.
- Re-sterilization by autoclave may cause a hand control short-circuit resulting in inadvertent activation when the electrode is connected to the generator.
- Avoid touching the distal tip of the instrument (ceramic insulator or metal active component) with fingers or instruments.
- Inadvertent activation or movement of an active VAPR electrode outside the field of vision may result in patient injury and damage to the electrode tip assembly.
- Using arthroscopic guidance, ensure that the electrode is completely immersed in conductive irrigant solution during use. Failure to completely immerse in irrigant solution so may result in electrode tip damage.
- Avoid bubble accumulation in the joint space during use. The accumulation of bubbles around the electrode diminishes performance and may produce overheating sufficient to damage adjacent structures or the electrode tip assembly.
- Do not wrap handpiece, footswitch or generator power cables around metal objects. Wrapping cables around metal objects may induce currents that could lead to electric shock, fire, or injury to patient and/or surgical personnel.
- Operating the electrode at minimal suction rate, or failure to attach the suction tubing to an adequate suction source may cause thermal injury to the physician or patient.

- Electrodes will wear with normal use, dependent on factors such as length of use, high tissue removal rate, prolonged use against bony surfaces, prolonged use in saline, high power settings, and use with minimal suction or fluid management. Periodically assess electrode tip for wear and proper operation. Replace the electrode if excessive wear is noted.
- Attach the suction adapter to standard hospital suction equipment. Ensure that the vacuum is in the range of 300 mmHg to 700 mmHg.
- Damage to surrounding tissue through iatrogenic injury could occur, as an effect of electrosurgery.
- Attempts to bend electrodes can result in electrode fracture or degradation of electrode performance.
- When not using instruments, place them in a clean, dry, highly visible area not in contact with the patient. Inadvertent contact with the patient may result in burns.
- Inspect instruments and cables for damage prior to each use, especially the insulation of electrode. Inspection may be done visually under magnification or with a high voltage insulation testing device. Insulation failures may result in burns or other injuries to the patient or operator.
- The surface of the active electrode may remain hot enough to cause burns after the RF current is deactivated.

PRECAUTIONS

- Prior to initial use, read and understand all package inserts including warnings, cautions, and Instructions for Use. Safe and effective electrosurgery is dependent not only on equipment design, but also on factors under the control of the user (refer to the user manual).
- Use caution when overriding the default power settings. Use the lowest power setting and the minimum tissue contact time necessary to achieve the appropriate surgical effect.
- When not in use, place the active electrode in a clean, dry, nonconductive, and highly visible area not in contact with the patient. Inadvertent activation while in contact with the patient may result in burns.
- Excessive force may damage the electrode tip assembly.
- Failure to maintain the recommended suction may result in device failure.
- The intensity should be set as low as is necessary to achieve the desired effect.
- Keep the active electrodes clean. Build-up of eschar may reduce the instrument’s effectiveness. Do not activate the instrument while cleaning. Injury to operating room personnel may result.

INSTRUCTIONS FOR USE

1. Read all instructions carefully. Failure to properly follow instructions may lead to electrical or thermal injury and cause improper functioning of the device.
2. The operator should be experienced in arthroscopic surgical techniques.
3. Please refer to the user manual for step-by-step instructions regarding the assembly and initial system check of the VAPR VUE RF system.

BEFORE SURGERY

1. Switch on the VAPR VUE generator. Follow the instructions in the user manual. The generator will display: “Attach Instrument”.
2. Connect a new, sterile electrode to the generator. Follow the instructions in the user manual. The generator will display the ABLATE/COAG screen.
3. For optimal performance, safety, and convenience, the VAPR electrodes automatically set the VAPR VUE generator to the default output settings. The default settings are as follows:

CAUTION: If the VAPR VUE generator display does not match the default settings specified below, consult the user manual for further instructions.

ELECTRODE GEOMETRY	COOLPULSE Curve and COOLPULSE Curve XL	
VAPR SYSTEM	VAPR VUE only	
Default mode	COOLPULSE	
	Basic	Advanced (W)
Default ABLATE Power (W)	4	200
Default COAGULATE Power (W)	3	80
Maximum ABLATE Power (W)	5	260
Maximum COAGULATE Power (W)	4	100
Minimum Cannula Size	6.5 mm	

Use caution when overriding the default power settings. Use the lowest power setting and the minimum tissue contact time necessary to achieve the appropriate surgical effect. When overriding power settings to non-default values, the adjusted values are retained when unplugging and re-plugging the same electrode, when connecting a new electrode of the same type, or when connecting a different electrode model with the same default values. If the generator is switched off and then on, a newly attached electrode will always revert to its default values. Always confirm proper, desired power settings before proceeding with surgery.

The generator display can be set up to display output power in Watts (ADVANCED mode) or in relative numbers (BASIC mode). The equivalent values for ADVANCED and BASIC mode vary by electrode type.

Below is a table for the COOLPULSE Curve and COOLPULSE Curve XL Suction Electrodes:

VAPR COOLPULSE Curve and COOLPULSE Curve XL Suction Electrodes			
ABLATE		COAGULATE	
Basic	Advanced (W)	Basic	Advanced (W)
1	10	1	10
2	100	2	40
3	150	3	80
4	200	4	100
5	260		

NOTE: ADVANCED mode provides smaller increments in output power than BASIC mode.

NOTE: When an instrument is attached, switching between ADVANCED and BASIC User Mode after adjusting the power output settings will reset the power output values to the default settings for the specific electrode type. The MODE button on the generator must be pressed to confirm use.

4. Attach the tubing connector to a suction source.

DURING SURGERY

1. Prior inserting the electrode, close electrode suction outflow by moving the Flow Control Valve as far proximal as possible.
2. Insert the electrode into a fluid-filled joint space through a prepared entry point. Ensure that the electrode is immersed in conductive irrigant solution when in use.
3. Ensure that the flow control valve is open (fully distal).

USE OF HAND CONTROLS AND FLOW CONTROL VALVE

4. To activate the ABLATE mode using the hand controls on the electrode handle (ensure suction is turned on), press and hold any of the YELLOW buttons. The ABLATE power setting on the generator will flash and emit a tone. Always keep the electrode in motion and in contact with tissue. Refer to the VAPR VUE User Manual on how to activate using the footswitch (wired or wireless). The hand controls are de-activated when using the footswitch.

5. Refer to the VAPR VUE User Manual for instructions on how to adjust the ABLATE and COAG power settings on the generator panel. Use the lowest power possible to achieve the desired tissue effect.

NOTE: ABLATE power settings can also be adjusted using the footswitch (wired or wireless). Please refer to the VAPR VUE User Manual.

6. To activate the COAG mode using the hand controls on the electrode handle, press and hold any of the BLUE buttons. The COAG power setting on the generator will flash and emit a tone (at a different frequency than during ABLATE). Always keep the electrode in contact with tissue. Refer to the VAPR VUE User Manual for instructions on how to activate using the footswitch (wired or wireless). Use the lowest power possible to achieve the desired tissue effect.

NOTE: COAG power settings are adjusted using only the buttons on the generator. Please refer to the VAPR VUE User Manual.

CAUTION: Based on the multiple button design, there is an increased risk of inadvertent activation.

7. Adjust the flow rate, using the Flow Control Valve, to optimize the fluid removal rate and bubble extraction, and to assist with bleeding blood vessels. Sliding the Flow Control Valve distally will increase flow (up to ~500ml/min); moving the Valve proximally will decrease flow.

8. Reduce suction when device is not in use.

CAUTION: Inspect the active tip of the electrode for any damage after all faults. If the active tip is loose or dislodged, insert a new electrode. If electrode becomes covered in tissue debris, use a new electrode.

Audible Alarm

As a safety feature, and to notify the user of the presence of a short circuit within the VAPR VUE system, a dual frequency audible alarm indicates either contact between the activating electrode tip and a metal object in the joint space, or an electrical fault in the electrode or handpiece.

If the audio signal continues when the electrode is moved away from a metal object, inspect the electrode and/or handpiece. Discard damaged accessories immediately. Remove excessive debris build-up across the distal tip area with an appropriate wipe.

If the alarm sounds continuously, and there is no contact with a metal object, discard the electrode.

AFTER SURGERY

1. Discontinue suction and remove the electrode from the joint space.
2. Disconnect from suction source. Disconnect the electrode cable from the generator by steadying the generator and pulling on the plug body. DO NOT pull on the cable; this can cause damage to the generator.
3. The generator can be left in idle mode between surgical procedures. Turn off the generator at the conclusion of a surgical day.
4. Dispose of the single use components according to hospital procedure.

HOW SUPPLIED

The VAPR COOLPULSE Curve and COOLPULSE Curve XL Suction Electrodes are provided sterile unless the package is opened or damaged. The electrodes are intended to be disposed of after a single use and are not intended to be re-sterilized (see warning section).

STORAGE

Store at controlled temperature up to 25° (77°F). Avoid excessive heat > 40°C (104°F). Bring to ambient temperature before unpacking /use. Do not use after USE BY date.

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DePuy Synthes

MITEK SPORTS MEDICINE

PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

DePuy Mitek VAPR® Non-suction Electrodes

VAPR® End Effect Electrode

3.5 mm x 160 mm

VAPR® Side Effect Electrode

3.5 mm x 160 mm

VAPR® 90 Degree Hook Electrode

3.5 mm x 160 mm

VAPR® Angled End Effect Electrode

3.5 mm x 160 mm

VAPR® Angled Side Effect Electrode

3.5 mm x 160 mm

VAPR® 2.3 Side Effect Electrode

140 mm

VAPR® 2.3 End Effect Electrode

140 mm

VAPR® 2.3 Side Effect Electrode (Short)

85 mm

VAPR® 2.3 Wedge Electrode (Short)

85 mm

VAPR® 2.3 Wedge Electrode

140 mm

VAPR® Flexible Side Effect Electrode

3.5 mm x 160 mm

VAPR® Flexible End Effect Electrode

3.5 mm x 160 mm

VAPR® Side Effect Electrode

3.5 mm Side Effect with Integrated Handpiece

VAPR® Hook Electrode

3.5 mm Hook with Integrated Handpiece

VAPR® Angled Side Effect Electrode

with Integrated Handpiece, 3.5 mm x 160 mm

VAPR® Flexible Side Effect Electrode

with Integrated Handpiece, 3.5 mm x 160 mm

VAPR® 2.3 Wedge Electrode

with Integrated Handpiece, 140 mm

VAPR® 2.3 Wedge Electrode (Short)

with Integrated Handpiece, 85 mm

VAPR® 2.3 Side Effect Electrode (Short)

with Integrated Handpiece, 85 mm

VAPR® 2.3 End Effect Electrode

with Integrated Handpiece, 140 mm

VAPR® 2.3 Wedge Electrode with Hand Controls

2.3 mm with Integrated Handpiece

VAPR® Side Effect Electrode with Hand Controls

3.5 mm Side Effect with Integrated Handpiece

VAPR® Hook Electrode with Hand Controls

3.5 mm Hook Electrode with Integrated Handpiece

CE 2797

Manufactured for DePuy Mitek by:

Gyrus Medical P/N: 117050-JD

P/N: 104285

Rev: N Revised 01/2020

© DePuy Synthes 2004-2020.

All rights reserved.

ENGLISH

DePuy Mitek VAPR® Non-suction Electrodes

INSTRUCTIONS FOR USE

Please read all information carefully, in particular the user manual provided with the DePuy Mitek VAPR generator, prior to use. The components of the DePuy Mitek VAPR systems are designed for use as a single system. Failure to follow instructions may lead to improper functioning of the device and cause electrical or thermal injury.

SYSTEM DESCRIPTION

The DePuy Mitek VAPR Electrosurgical System is designed for arthroscopic surgical procedures. The system consists of a high frequency electrosurgical generator, a reusable handpiece with connector cable, disposable electrodes, alternative disposable electrodes with integrated connector cable, also available with hand controls, a foot switch and a wireless foot switch alternative. The components are designed and intended to be operated as a single unit. Use only DePuy Mitek VAPR electrodes and accessories with DePuy Mitek VAPR systems.

DEVICE DESCRIPTION

The VAPR Non-suction Electrodes are soft tissue ablation and hemostasis / coagulation devices intended for use with the VAPR system.

INDICATIONS

The DePuy Mitek VAPR Electrodes for use with the VAPR Radiofrequency System are intended for resection, ablation and excision of soft tissue, hemostasis of blood vessels, and coagulation of soft tissue in patients requiring arthroscopic surgery of the knee, shoulder, ankle, elbow, and/or wrist.

CONTRAINDICATIONS

The DePuy Mitek VAPR Electrosurgical System is contraindicated in:

- any non-arthroscopic surgical procedure.
- any arthroscopic procedures where saline or Ringer's Lactate is not used as an irrigant.
- patients with heart pacemakers or other electronic device implants.
- patients for whom arthroscopic procedures are contraindicated for any reason.

WARNINGS

- The operator should be experienced in arthroscopic surgical techniques.
- Please refer to the user manual for step-by-step instructions regarding the assembly and initial system check of the DePuy Mitek VAPR system.
- As with all electrosurgical devices, do not use in the presence of flammable anesthetics, oxidizing gases or other flammable substances, as an electrosurgical device has the potential for providing a source of ignition.
- Carefully insert and withdraw electrodes to avoid possible damage to the device and/or injury to the patient or surgical personnel.
- Ensure that fluid inflow and outflow is adequate and that the electrode is activated only when surrounded by conductive irrigant solution (e.g., saline or Ringer's Lactate). The use of RF energy with inadequate irrigation may overheat the fluid enough to cause skin burns at or near the access site and may result in damage to the electrode tip assembly.
- Do not touch the electrode tip when power is being applied.
- During an electrosurgical procedure, the patient should not be allowed to come into direct contact with grounded metal objects. Inadvertent contact with the patient when the VAPR system is activated may result in burns.

- Observe extreme caution when using electrosurgery in close proximity to, or in direct contact with, any metal objects. The majority of arthroscopes and arthroscopic instruments are metal. Do not activate the electrode while any portion of the electrode tip is in contact with metal objects or instruments. Doing so may result in unintended injury to the patient or surgical personnel and/or damage to the electrode and/or other equipment.
- Avoid unnecessary and prolonged activation of the electrode during and between tissue applications. Failure to do so may result in unintended injury.
- Do not reuse any accessories labeled as SINGLE USE. These products are for single-use only and have not been designed to be re-used/re-sterilized. Reprocessing may lead to changes in material characteristics such as metallic corrosion and dulled edges, ceramic and plastic deformation or splitting which may impact the strength of the device and compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection. These risks may potentially affect patient safety.
- Re-sterilization by autoclave may cause hand control shorts resulting in inadvertent activation when the electrode is connected to the generator.
- Avoid touching the distal tip of the instrument (ceramic insulator or metal active component) with fingers or instruments.
- Inadvertent activation or movement of an active VAPR electrode outside the field of vision may result in patient injury and damage to the electrode tip assembly.
- Using arthroscopic guidance, ensure that the electrode is completely surrounded by conductive irrigant solution during use. Failure to do so may result in electrode tip damage.
- Avoid bubble accumulation in the joint space during use. The accumulation of bubbles around the electrode diminishes performance and may produce overheating sufficient to damage adjacent structures or the electrode tip assembly.
- Do not wrap handpiece, footswitch or generator power cables around metal objects. Wrapping cables around metal objects may induce currents that could lead to electric shock, fire or injury to patient and/or surgical personnel.
- Failure to attach the suction tubing to an adequate suction source may cause thermal injury to the physician or patient.
- DO NOT activate the probe for unnecessary and prolonged periods as irrigant overheating and unintended tissue damage may result. Prolonged probe activation while using the suction feature may result in elevated shaft or suction tubing temperatures.
- Electrodes will wear from normal use, dependent on factors such as length of use, high tissue removal rate, prolonged use against bony surfaces, prolonged use in saline, high power settings, and use with minimal suction or fluid management. Periodically assess electrode tip for wear and proper operation. Replace the electrode if excessive wear is noted.
- For electrodes with suction, attach the suction adapter to standard hospital suction equipment. Ensure that the roller clamp is open, and vacuum is in the range of 300 mmHg to 700 mmHg.
- Damage to surrounding tissue through iatrogenic injury could occur, as an effect of electrosurgery.
- Do not withdraw or insert the electrode when power is being applied.
- Due to low joint-space volume in smaller joints, such as the wrist, it may be necessary to manually lower the power settings or increase flow through the joint in order to maintain adequate inter-joint cooling.

PRECAUTIONS

- Prior to initial use, ensure that all package inserts including warnings, cautions, and Instructions for Use are read and understood. Safe and effective electrosurgery is dependent not only on equipment design but also on factors under the control of the user (refer to the user manual).

- When not in use, place the active Electrode in a clean, dry, non-conductive, and highly visible area not in contact with the patient. Inadvertent activation while in contact with the patient may result in burns.
- Excessive force may damage the electrode tip assembly.
- Ensure that the Electrode is fully seated in the Handpiece prior to use. (Note: Not applicable on Electrodes with integrated cable.)
- Avoid fluid contact with the proximal end of the Electrode/Handpiece interface. (Note: Not applicable on Electrodes with integrated cable.)
- Attempts to bend non-FLEX 3.5 mm Electrodes can result in electrode fracture.

INSTRUCTIONS FOR USE

- Read all instructions carefully. Failure to properly follow instructions may lead to electrical or thermal injury and cause improper functioning of the device.
- The operator should be experienced in arthroscopic surgical techniques.
- Please refer to the user manual for step-by-step instructions regarding the assembly and initial system check of the DePuy Mitek VAPR system.

SHAPING THE FLEX ELECTRODE

- If needed, carefully use your fingers to bend the shaft of the FLEX Electrode to the desired shape.
- The Electrode can be bent a maximum of 45° three separate times at the same point.

SHAPING THE 2.3 MM VAPR ELECTRODE

- If needed, carefully use your fingers to bend the shaft of the 2.3 mm Electrode to the desired shape.
- The Electrode can be bent one time to a maximum angle of 30° from the plane of the shaft.

BEFORE SURGERY – VAPR GENERATOR

1. Switch on the DePuy Mitek VAPR Generator. Follow the instructions in the User Manual so that the Generator displays the symbol “CONNECT CABLE.”
- *2. Connect a sterilized Handpiece to the Generator. Follow the instructions in the User Manual so that the Generator displays the symbol “INSERT ELECTRODE.”
3. If you have selected a VAPR FLEX Electrode or a 2.3 mm VAPR Electrode, bend the Electrode as needed.
4. Connect a new DePuy Mitek VAPR Electrode to the Handpiece. Check that the Handpiece and Electrode are properly aligned before pushing together and tightening the screw lock (See Illustration).
***NOTE:** These steps do not apply when using the DePuy Mitek VAPR Electrode with integrated handpiece. In this case, connect the plug at the end of the integrated cable directly to the generator.
5. For safety and convenience, VAPR Electrodes automatically preset the Generator to default output settings and restrict output levels to preserve integrity. If alternate power levels are desired, adjust the VAPR Generator to the desired settings using the BLUE and YELLOW arrows on the front of the VAPR Generator. Use the lowest power possible to achieve the desired end effect.

The default settings for each Electrode are as follows:

REFERENCE ALPHA CODES IN THE CHARTS BELOW FOR ELECTRODE SETTINGS

A – Default Mode

B – Default Vaporization Power (w)

C – Default Desiccation Power (w)

D – Maximum Vaporization Power (w)

E – Maximum Desiccation Power (w)

3.5 MM VAPR ELECTRODES & VAPR FLEX ELECTRODES					
Electrode Type	A	B	C	D	E
Side Effect	V2	120	90	180	120
End Effect	V2	90	60	130	90
Hook	BV2	120	90	180	120

2.3 MM VAPR Electrodes					
Electrode Type	A	B	C	D	E
Side Effect	V2	60	45	90	60
End Effect	V2	90	60	130	90
Wedge	V2	60	45	90	60

CAUTION: Verify the VAPR Generator display matches the VAPR Electrode default setting as specified above. If the Generator display does not match the Electrode default settings refer to User Manual for further instructions.

Using arthroscopic guidance, carefully insert the non-activated Electrode into the joint cavity portal under direct vision. Wherever possible, use an instrument cannula for the access portal.

CHANGING VAPR ELECTRODES DURING SURGERY

- Carefully remove the Electrode from the joint cavity. Do not withdraw the Electrode while power is being applied.
- Unlock the connector assembly and disengage the Electrode connector and Handpiece. Check that the Generator displays the symbol “INSERT ELECTRODE.” (Note: Not applicable on Electrodes with integrated cable. In such cases, unplug the electrode from the generator.)
- Fit a new sterile Electrode as previously described. If the new Electrode has different default settings to the previous Electrode, check that the Generator display matches the Electrode default settings specified on the package insert. If the new Electrode has identical default settings to the previous Electrode, the Generator will retain the settings previously displayed prior to changing the Electrode. If the Generator display does not conform to the above description refer to User Manual for further instructions.

IF USING ELECTRODE WITH HAND CONTROLS AND VAPR VUE GENERATOR

- To activate the ABLATE mode using the hand controls on the electrode handle (ensure suction is turned on), press and hold the YELLOW button. The ABLATE power setting on the generator will flash and a tone will be heard. Always keep the electrode in motion and in contact with tissue. Refer to the VAPR VUE User Manual on how to activate using the footswitch (wired or wireless).
- The ABLATE power setting may be changed by holding down the BLACK MODE button on the electrode handle. The ABLATE power setting may then be increased by pressing the YELLOW button on the electrode handle. The ABLATE power setting may be decreased by pressing the BLUE button on the electrode handle. Pressing the MODE button again will confirm and hold the selected value. Use the lowest power possible to achieve the desired tissue effect.

NOTE: The mode button on the electrode handle can only change ablation power. Coagulation power can be adjusted only using the buttons on the generator.

- To activate the COAG mode using the hand controls on the electrode handle, press and hold the BLUE button. The COAG power setting on the generator will flash and a tone (at a different frequency than during ABLATE) will be heard. Always keep the electrode in motion and in contact with tissue. Refer to the VAPR VUE User Manual on how to activate using the footswitch (wired or wireless). Use the lowest power possible to achieve the desired tissue effect.

BASIC-ADVANCED EQUIVALENCY TABLE IN VAPR VUE GENERATOR:

3.5 MM VAPR SIDE EFFECT, ANGLE SIDE, HOOK, FLEX ELECTRODES			
ABLATE		COAGULATE	
Basic	Advanced	Basic	Advanced
1	10	1	10
2	35	2	45
3	60	3	60
4	80	4	90
5	120	5	120
6	180		

2.3 MM & 3.5 MM VAPR END EFFECT			
ABLATE		COAGULATE	
Basic	Advanced	Basic	Advanced
1	10	1	10
2	35	2	45
3	50	3	60
4	90	4	90
5	130		

2.3 MM VAPR SIDE EFFECT AND WEDGE ELECTRODES			
ABLATE		COAGULATE	
Basic	Advanced	Basic	Advanced
1	10	1	10
2	35	2	45
3	60	3	60
4	90		

NOTE: In ADVANCED Mode, output power can be adjusted in finer increments than in BASIC Mode.

NOTE: When an instrument is attached, switching between ADVANCED and BASIC User Mode after adjusting the power output settings will reset the power output values to the default settings for the specific electrode type. The Mode button must be pressed to confirm use.

BEFORE SURGERY

1. Insert the electrode into a fluid-filled joint space through a prepared entry point. Ensure that the electrode is surrounded by conductive irrigant solution when in use.

CAUTION: Inspect the active tip of the electrode for any damage after all faults. If the active tip is loose or dislodged, please insert a new electrode.

AFTER SURGERY

1. Remove the electrode from the joint space.
2. Disassemble the VAPR Electrode and VAPR Handpiece. (Note: Not applicable on Electrodes with integrated cable.)
3. Disconnect the electrode cable from the generator by steadying the generator and pulling on the plug body. DO NOT pull on the cable. Pulling the cable may cause damage to the generator.
4. The generator can be left in idle mode between surgical procedures. Turn off the generator at the conclusion of a surgical day.

5. Dispose of the single use components per hospital procedure.
6. Prepare reusable accessories for cleaning and sterilization.

HOW SUPPLIED

The Electrodes used in the DePuy Mitek VAPR System are provided sterile and are disposable for single use. The Electrodes are available in a variety of configurations, default settings, and flexibility. Important: The contents are sterile unless the package is opened or damaged. Do not resterilize.

STORAGE

Store below 25°C (77°F) away from moisture and direct heat. Do not use after use by date.

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VAPR® LD Suction Electrode
4.0 mm x 160 mm

VAPR® LP Suction Electrode
4.0 mm x 160 mm

VAPR® S90 Electrode
4.0 mm, 90° Suction with Integrated Handpiece

VAPR® S50 Electrode
50° Suction with Integrated Handpiece

VAPR® Premiere50 Electrode
50° Suction with Integrated Handpiece

VAPR® 90 Degree Suction Electrode
3.5 mm x 140 mm

VAPR® Premiere90 Electrode
90° Suction with Integrated Handpiece

VAPR® S90 Electrode with Hand Controls
4.0 mm, 90° Suction with Integrated Handpiece

VAPR® Premiere50 Electrode with Hand Controls
3.0 mm, 50° Suction with Integrated Handpiece

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Manufactured for DePuy Mitek by:
Gyrus Medical P/N: 118014-MD

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SYSTEM DESCRIPTION

The DePuy Mitek VAPR Electrosurgical System is designed for arthroscopic surgical procedures. The system consists of a high frequency electrosurgical generator, a reusable handpiece with connector cable, disposable electrodes, alternative disposable electrodes with integrated connector cable, a footswitch and a wireless footswitch alternative. The components are designed and intended to be operated as a single unit. Use only DePuy Mitek VAPR electrodes and accessories with DePuy Mitek VAPR systems.

DEVICE DESCRIPTION

The VAPR Suction Electrodes are soft tissue ablation and desiccation devices intended for use with the VAPR system. They extend the utility of the system by removing bubbles created during activation from the operating site.

INDICATIONS

The DePuy Mitek VAPR Electrodes for use with the VAPR System are intended for resection, ablation and excision of soft tissue, hemostasis of blood vessels and coagulation of soft tissue in patients requiring arthroscopic surgery of the knee, shoulder, ankle, elbow, and wrist.

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WARNINGS

- The operator should be experienced in arthroscopic surgical techniques.
- Please refer to the user manual for step-by-step instructions regarding the assembly and initial system check of the DePuy Mitek VAPR system.
- As with all electrosurgical devices, do not use in the presence of flammable anesthetics, oxidizing gases  flammable substances, as an electrosurgical device has the potential for providing a source of ignition.
- Carefully insert and withdraw electrodes to avoid possible damage to the device and/or injury to the patient or surgical personnel.
- Ensure that fluid inflow and outflow is adequate and that the electrode is activated only when surrounded by conductive irrigant solution (e.g., saline or Ringer's Lactate). The use of RF energy with inadequate irrigation may overheat the fluid enough to cause skin burns at or near the access site and may result in damage to the electrode tip assembly.
- Do not touch the electrode tip when power is being applied.
- During an electrosurgical procedure, the patient should not be allowed to come into direct contact with grounded metal objects. Inadvertent contact with the patient when the VAPR system is activated may result in burns.
- Observe extreme caution when using electrosurgery in close proximity to, or in direct contact with, any metal objects. The majority of arthroscopes and arthroscopic instruments are metal. Do not activate the electrode while any portion of the electrode tip is in contact with metal objects or instruments. Doing

so may result in unintended injury to the patient or surgical personnel and/or damage to the electrode and/or other equipment.

- Avoid unnecessary and prolonged activation of the electrode during and between tissue applications. Failure to do so may result in unintended injury.
- Do not reuse any accessories labeled as SINGLE USE. These products are for single-use only and have not been designed to be re-used/re-sterilized. Reprocessing may lead to changes in material characteristics such as metallic corrosion and dulled edges, ceramic and plastic deformation or splitting which may impact the strength of the device and compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection. These risks may potentially affect patient safety.
- Re-sterilization by autoclave may cause hand control shorts resulting in inadvertent activation when the electrode is connected to the generator.
- Avoid touching the distal tip of the instrument (ceramic insulator or metal active component) with fingers or instruments.
- Inadvertent activation or movement of an active VAPR electrode outside the field of vision may result in patient injury and damage to the electrode tip assembly.
- Using arthroscopic guidance, ensure that the electrode is completely surrounded by conductive irrigant solution during use. Failure to do so may result in electrode tip damage.
- Avoid bubble accumulation in the joint space during use. The accumulation of bubbles around the electrode diminishes performance and may produce overheating sufficient to damage adjacent structures or the electrode tip assembly.
- Do not wrap handpiece, footswitch or generator power cables around metal objects. Wrapping cables around metal objects may induce currents that could lead to electric shock, fire or injury to patient and/or surgical personnel.
- Failure to attach the suction tubing to an adequate suction source may cause thermal injury to the physician or patient.
- DO NOT activate the probe for unnecessary and prolonged periods as irrigant overheating and unintended tissue damage may result. Prolonged probe activation while using the suction feature may result in elevated shaft or suction tubing temperatures.
- Electrodes will wear from normal use, dependent on factors such as length of use, high tissue removal rate, prolonged use against bony surfaces, prolonged use in saline, high power settings, and use with minimal suction or fluid management. Periodically assess electrode tip for wear and proper operation. Replace the electrode if excessive wear is noted.
- For electrodes with suction, attach the suction adapter to standard hospital suction equipment. Permanently Marked Electrode the roller clamp is open, and vacuum is in the range of 300 mmHg to 700 mmHg.
- Damage to surrounding tissue through iatrogenic injury could occur, as an effect of electrosurgery.
- Do not withdraw or insert the electrode when power is being applied.

PRECAUTIONS

- Prior to initial use, ensure that all package inserts including warnings, cautions, and Instructions for Use are read and understood. Safe and effective electrosurgery is dependent not only on equipment design but also on factors under the control of the user (refer to the user manual).
- When not in use, place the active Electrode in a clean, dry, nonconductive, and highly visible area not in contact with the patient. Inadvertent activation while in contact with the patient may result in burns.
- Excessive force may damage the electrode tip assembly.
- Attempts to bend electrodes can result in electrode fracture or degradation of electrode performance.
- Failure to maintain the recommended suction may result in device failure.

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BEFORE SURGERY – VAPR Generator

1. Switch on the DePuy Mitek VAPR Generator. Follow the instructions in the User Manual so that the Generator displays the symbol “CONNECT CABLE.”
- *2. Connect a sterilized Handpiece to the Generator. Follow the instructions in the User Manual so that the Generator displays the symbol “INSERT ELECTRODE.”
3. Connect a new DePuy Mitek VAPR Electrode to the Handpiece. Check that the Handpiece and Electrode are properly aligned before pushing together and tightening the screw lock (See Illustration).

*** NOTE:** These steps do not apply when using the DePuy Mitek VAPR Electrode with integrated handpiece. In this case, connect the plug at the end of the integrated cable directly to the generator.

4. For safety and convenience, VAPR Electrodes automatically preset the Generator to default output settings and restrict output levels to preserve electrode integrity. If alternate power levels are desired, adjust the VAPR Generator to the desired settings using the BLUE and YELLOW arrows on the front of the VAPR Generator. Use the lowest power possible to achieve the desired end effect.

IF USING ELECTRODE with HAND CONTROLS and VAPR VUE Generator

1. To activate the ABLATE mode using the hand controls on the electrode handle (ensure suction is turned on), press and hold the YELLOW button. The ABLATE power setting on the generator will flash and a tone will be heard. Always keep the electrode in motion and in contact with tissue. Refer to the VAPR VUE User Manual on how to activate using the footswitch (wired or wireless).
2. The ABLATE power setting may be changed by holding down the BLACK MODE button on the electrode handle. The ABLATE power setting may then be increased by pressing the YELLOW button on the electrode handle. The ABLATE power setting may be decreased by pressing the BLUE button on the electrode handle. Pressing the MODE button again will confirm and hold the selected value. Use the lowest power possible to achieve the desired tissue effect.

Регистрация МН в Росздравнадзоре
www.peracort.ru | info@tevosort.ru
NOTE: The mode button on the electrode handle can only change ablation power. Coagulation power can be adjusted only using the buttons on the generator.

3. To activate the COAG mode using the hand controls on the electrode handle, press and hold the BLUE button. The COAG power setting on the generator will flash and a tone (at a different frequency than during ABLATE) will be heard. Always keep the electrode in motion and in contact with tissue. Refer to the VAPR VUE User Manual on how to activate using the footswitch (wired or wireless). Use the lowest power possible to achieve the desired tissue effect.
4. Caution should be used when overriding the default power settings. Use the lowest power setting and the minimum tissue contact time necessary to achieve the appropriate surgical effect. If power settings have been overridden to non-default values, the adjusted values will be retained when unplugging and replugging the same electrode, when connecting a new electrode of the same type, or when connecting a different electrode model with the same default values. If the generator is switched off and then on, a newly attached electrode will always revert to its

default values. Always confirm proper, desired power settings before proceeding with surgery.

- The generator display can be set up to display output power in Watts (ADVANCED mode) or in relative numbers (BASIC mode). The equivalent values for ADVANCED and BASIC mode vary by electrode type.

NOTE: In ADVANCED mode, output power can be adjusted in finer increments than in BASIC mode.

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- Attach the tubing connector to a suction source. Ensure that suction is off by closing the clamp on the electrode's tubing.

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CAUTION: If the VAPR generator display does not match the default settings specified below, consult the user manual for further instructions.

ELECTRODE GEOMETRY	3.5mm 90° Suction Electrode			
VAPR SYSTEM	VAPR; VAPR II	VAPR 3	VAPR VUE	
Default mode	V2	V3	V3	
			Adv	Basic
Default Vaporization Power (W)	160	160	160	6
Default Desiccation Power (W)	120	120	120	5
Maximum Vaporization Power (W)	200	200	200	7
Maximum Desiccation Power (W)	160	160	160	6
Minimum Cannula Size	5.5mm ID cannula			

ELECTRODE GEOMETRY	LDS Electrode			
VAPR SYSTEM	VAPR; VAPR II	VAPR 3	VAPR VUE	
Default mode	BV2	V3	V3	
			Adv	Basic
Default Vaporization Power (W)	160	240	240	7
Default Desiccation Power (W)	120	120	120	5
Maximum Vaporization Power (W)	200	260	260	8
Maximum Desiccation Power (W)	160	160	160	6
Minimum Cannula Size	6.3mm ID or DePuy Mitek 6mm cannula P/N# 214660, 214661			

ELECTRODE GEOMETRY	LPS Electrode/ VAPR S ⁹⁰ Electrode/Premiere90 Electrode			
VAPR SYSTEM	VAPR; VAPR II	VAPR 3	VAPR VUE	
Default mode	BV2	V3	V3	
			Adv	Basic
Default Vaporization Power (W)	160	240	240	7
Default Desiccation Power (W)	120	120	120	5
Maximum Vaporization Power (W)	200	260	260	8
Maximum Desiccation Power (W)	160	160	160	6
Minimum Cannula Size	5.5mm ID cannula			

ELECTRODE GEOMETRY	S ⁵⁰ Electrode			
VAPR SYSTEM	VAPR; VAPR II	VAPR 3	VAPR VUE	
Default mode	V2	V2	V2	
			Adv	Basic
Default Vaporization Power (W)	60	60	60	3
Default Desiccation Power (W)	45	45	45	2
Maximum Vaporization Power (W)	90	90	90	4
Maximum Desiccation Power (W)	60	60	60	3
Minimum Cannula Size	Not Applicable, Stab incision only			

ELECTRODE GEOMETRY	Premiere50 Electrode			
VAPR SYSTEM	VAPR; VAPR II	VAPR 3	VAPR VUE	
Default mode	BV2	V2	V2	
			Adv	Basic
Default Vaporization Power (W)	60	160	160	6
Default Desiccation Power (W)	45	45	45	2
Maximum Vaporization Power (W)	90	200	200	7
Maximum Desiccation Power (W)	60	60	60	3
Minimum Cannula Size	Not Applicable, Stab incision only			

DURING SURGERY

1. Close suction outflow of the electrode prior to insertion of the electrode.
2. Insert the electrode into a fluid-filled joint space through a prepared entry point. Ensure that the electrode is surrounded by conductive irrigant solution when in use.
3. Adjust the flow rate, by opening suction outflow, to optimize the fluid removal rate and bubble extraction.
4. Reduce suction when device is not in use.

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Do not allow the electrode to become covered in tissue debris. If this occurs, use a new electrode.

5. As a safety feature and to notify the user of the presence of a short circuit within the VAPR system, a dual frequency audible alarm will sound. This will indicate either contact between the activating electrode tip and metal object in the joint space or an electrical fault in the electrode or handpiece. When the electrode is moved away from a metal object, and the audio signal continues, the electrode and/or handpiece should be inspected. Damaged accessories should immediately be discarded. Extensive debris build up across the distal tip area should be removed with an appropriate wipe.

Where there is no contact with a metal object and the alarm persistently sounds, the electrode must be discarded. If the alarm persists, then the handpiece should be discarded.

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HOW SUPPLIED

The VAPR electrode is provided sterile unless the package is opened or damaged. It is intended to be disposed of after a single use and is not intended to be resterilized (see warning section).

STORAGE

Store below 25°C (77°F) and away from moisture and direct heat. Bring to ambient temperature before unpacking/use. Do not use after USE BY date.

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VAPR® TC Electrode

**VAPR® Temperature Control End
Effect Electrode with Integrated
Handpiece**

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Manufactured for DePuy Mitek by:
Gyrus Medical P/N: 124316-FC

P/N: 106261

Rev: K Revised 01/2020

© DePuy Synthes 2004-2020.

All rights reserved.

ENGLISH

VAPR® TC Electrode

Please read all information carefully, in particular the User Manual that is provided with the DePuy Mitek VAPR Generator prior to use. The components of the VAPR System are designed for use together as a System. Failure to follow instructions may lead to improper functioning of the device and cause electrical or thermal injury.

SYSTEM DESCRIPTION

The DePuy Mitek VAPR Electrosurgical System is designed for arthroscopic surgical procedures. Each system consists of a high frequency electrosurgical Generator, a reusable handpiece with connector cable, disposable electrodes, alternative disposable electrodes with integrated connector cable, and a footswitch. The components are designed and intended to be operated as a single unit. Use only DePuy Mitek VAPR Electrodes and accessories with all DePuy Mitek VAPR Systems.

DEVICE DESCRIPTION

The VAPR Temperature Control Electrode is a soft tissue desiccation device intended for use with the VAPR Systems. Utilization with a VAPR System allows the tip temperature of the electrode to be indicated on the generator display.

INDICATIONS

The DePuy Mitek VAPR Electrodes for use with the VAPR System are intended for coagulation of soft tissues in patients requiring arthroscopic surgery of the knee, shoulder, ankle, elbow, and wrist.

CONTRAINDICATIONS

The DePuy Mitek VAPR Electrosurgical System is contraindicated in any non-arthroscopic surgical procedure and in any arthroscopic procedures where saline or Ringer's Lactate is not used as an irrigant. Use of the system is also contraindicated in patients with heart pacemakers or other electronic device implants, or for patients for whom arthroscopic procedures are contraindicated for any reason.

HOW SUPPLIED

The VAPR Temperature Control Electrode is provided sterile unless the package is opened or damaged. It is intended to be disposed of after a single use.

INSTRUCTIONS FOR USE

BEFORE SURGERY

1. Switch on the DePuy Mitek VAPR Generator. Follow the instructions in the User Manual so that the generator displays the symbol "CONNECT CABLE."

- *2. Регистрация НД в Росздравнадзоре
www.roszdrnadzor.ru/nadzor/registratsiya Connect a sterile Handpiece to the Generator. Follow the instructions in the User Manual so that the Generator displays the symbol "INSERT ELECTRODE."

- *3. Connect a new VAPR Temperature Control Electrode to the Handpiece. Check that the Handpiece and Electrode are properly aligned before pushing together and tightening the screw lock.

***NOTE:** Steps 2-3 do not apply when using the DePuy Mitek VAPR Electrode with integrated handpiece. In this case, connect the plug at the end of the integrated cable directly to the generator.

4. For optimal performance, safety and convenience, the VAPR Temperature Control electrode automatically presets the generator to default output. If alternate power levels are desired with the VAPR unit, adjust the VAPR Generator to the desired settings using the BLUE and YELLOW arrows on the front of the VAPR Generator. Use the lowest possible power to achieve the desired end-effect. When the VAPR TC Electrode is used with the VAPR II, VAPR 3 or VAPR VUE System, vaporization is disabled and the display

is replaced by either the set temperature (SET) or the tip temperature (TIP). If an alternate power level is desired with the VAPR II, VAPR 3 or VAPR VUE Generator, adjust it using the BLUE arrows on the front of the generator. Adjust the Set Temperature using the YELLOW arrows on the front of the VAPR Generator. The default settings are as follows:

VAPR TC Electrode					
Electrode Geometry	2.3 mm End-Effect				
VAPR System	VAPR	VAPR II	VAPR 3	VAPR VUE	
Default mode	V1	N/A	N/A	Adv	Basic
Default Vaporization Power (W)	5	N/A	N/A	N/A	N/A
Default Desiccation Power (W)	20	20	40	40	2
Maximum Vaporization Power (W)	10	N/A	N/A	N/A	N/A
Maximum Desiccation Power (W)	50	50	50	50	3
Default Set Temp (°C)	N/A	65C	65C	65C	65C

CAUTION: If the VAPR Generator display does not match the default setting as specified above consult the User Manual for further instructions.

DURING SURGERY

1. Insert the electrode into a fluid filled joint space through a prepared entry point. Ensure that the electrode is surrounded by conductive irrigant solution when in use.
2. Use the lowest power setting to achieve the desired clinical effect.
3. Use the blue foot pedal to activate the electrode. When used with a VAPR II or VAPR 3 System, the Set Temperature (SET) will be replaced by the Tip Temperature (TIP) and the desiccation power value will flash.
4. If the power level is inappropriately low, the Tip Temperature (TIP) may not approach the Set Temperature (SET).
5. As a safety feature, if the Tip Temperature (TIP) exceeds the Set Temperature (SET) by 8°C, then a distinct alarm will sound.

AFTER SURGERY

1. Remove the electrode from the joint space.
2. Disassemble the VAPR TC Electrode and VAPR Handpiece (if used). Disconnect the Handpiece/ electrode cable from the Generator by pulling on the plug body. DO NOT pull on the cable as it may cause damage to the device.
3. The generator can be left in idle mode between surgical procedures. Turn off the generator at the conclusion of a surgical day.
4. Dispose of the single use components of the System, and prepare the reusable accessories for cleaning and sterilization.

WARNINGS

1. The operator should be experienced in arthroscopic surgical techniques.
2. Please refer to the User Manual for step-by-step instructions regarding the assembly and initial system check of the DePuy Mitek VAPR System.
3. Ensure that fluid inflow and outflow is adequate, and the electrode is only activated when surrounded by conductive irrigant solution (e.g. saline or Ringer's Lactate).
4. Do not touch the electrode tip when power is being applied to the tip.
5. Do not withdraw or insert the electrode when power is being applied to the tip.

6. Activation of the electrode outside the field of vision may result in patient injury.
7. Avoid unnecessary and prolonged activation of the electrode between tissue applications as unintended injury may result.
8. Do not reuse any accessories labeled as Single Use. This product is for single use only. It has not been designed to be re-used/re-sterilized. Reprocessing may lead to changes in material characteristics such as metallic corrosion and dulled edges, ceramic and plastic deformation or splitting which may impact the strength of the device and compromise device performance. Reprocessing of single-use devices can also cause cross-contamination leading to patient infection. These risks may potentially affect patient safety.
9. Avoid touching the distal tip of the instrument (insulator and metal active component) with fingers or instruments.

CAUTIONS

1. Prior to initial use, ensure that all package inserts including warnings, cautions and Instructions for Use are read and understood. Safe and effective electrosurgery is dependent not only upon equipment design but also on factors under the control of the user (refer to the User Manual).
2. As with all electrosurgical devices, do not use in the presence of flammable anesthetics, oxidizing gases or other flammable substances, as an electrosurgical device has the potential for providing a source of ignition.
3. Place accessories in a clean, dry, non-conductive and highly visible area (away from the patient) when not in use. The Electrode tip may remain hot sometime after it has ceased to be activated.
4. Inadvertent contact with the patient when the VAPR System is activated may result in burns. The patient should not be allowed to come in contact with grounded metal objects.
5. Do not wrap the Handpiece cable around metal objects, as this may induce currents that could lead to electric shocks, fires or injury to the patient and or surgical personnel.
6. Do not activate the electrode while it is in contact with metal objects or instruments as unintended injury to the patient may occur.
7. Carefully insert and withdraw Electrodes to avoid possible damage to the device and/or injury to the patient or surgical personnel.

IMPORTANT

This package insert is designed to provide instructions for use of VAPR TC Electrode. It is not intended as a reference to electrosurgical technique.

STORAGE

Store below 25°C (77°F) away from moisture and direct heat. Do not use after 'USE BY' date.

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225028 VAPR® TRIPOLAR 90 SUCTION Electrode

**For use only with the VAPR VUE®
Generator**

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© DePuy Synthes 2012-2020.

All rights reserved.

Manufactured for DePuy Mitek Inc. by:

Gyrus Medical P/N: 195021-AE

VAPR® TRIPOLAR 90 SUCTION Electrode

INSTRUCTIONS FOR USE

Please read all information carefully prior to use, with particular attention to The User Manual provided with the VAPR VUE® Radiofrequency (RF) System. The components of the VAPR VUE System are designed for use as a single system. Failure to follow instructions may lead to improper functioning of the device and cause electrical or thermal injury.

SYSTEM DESCRIPTION

The VAPR VUE Radiofrequency (RF) System is designed for arthroscopic surgical procedures. The system consists of a high frequency electrosurgical generator and may include a reusable handpiece with connector cable and disposable electrodes, alternative disposable electrodes with integrated connector cable, also available with hand controls, a footswitch, and a wireless footswitch alternative. The components are designed and intended to be operated as a single unit. Use only VAPR® electrodes and accessories with the VAPR VUE RF system.

DEVICE DESCRIPTION

The VAPR Electrodes are soft tissue ablation and hemostasis / coagulation devices intended for use with the VAPR VUE RF System. The VAPR® TRIPOLAR 90 SUCTION Electrode has suction capabilities which enhance the efficiency of the electrode and extend the utility of the system by assisting in the removal of bubbles and debris created during activation within the operating site.

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VAPR electrodes are intended to be used by an orthopedic surgeon in a surgical setting.

CONTRAINDICATIONS

The VAPR Electrosurgical System is contraindicated in:

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WARNINGS

- The operator should be experienced in arthroscopic surgical techniques.
- Please refer to the user manual for step-by-step instructions regarding the assembly and initial system check of the VAPR VUE system.
- As with all electrosurgical devices, do not use in the presence of flammable anesthetics, oxidizing gases or other flammable substances, as an electrosurgical device has the potential for providing a source of ignition.
- Carefully insert and withdraw electrodes to avoid possible damage to the device and/or injury to the patient or surgical personnel.
- Ensure that fluid inflow and outflow is adequate and that the electrode is activated only when surrounded by conductive irrigant solution (e.g., saline or Ringer's Lactate). The use of RF energy with inadequate irrigation may overheat the fluid enough to cause skin burns at or near the access site and may result in damage to the electrode tip assembly.
- **Do not activate the probe for unnecessary and prolonged periods as irrigant overheating and tissue damage may result. Prolonged probe activation at maximum power while adjusting irrigant flow via the flow control valve to a low flow rate may cause the shaft and suction tubing to reach temperatures of up to 58°C. Do not touch the electrode tip when power is being applied.**
- During an electrosurgical procedure, the patient should not be allowed to come into direct contact with grounded metal objects. Inadvertent contact with the patient when the VAPR system is activated may result in burns.
- Observe extreme caution when using electrosurgery in close proximity to, or in direct contact with, any metal objects. The majority of arthroscopes and arthroscopic instruments are metal. Do not activate the electrode while any portion of the electrode tip is in contact with metal objects or instruments. Doing so may result in unintended injury to the patient or surgical personnel and/or damage to the electrode and/or other equipment.
- Avoid unnecessary and prolonged activation of the electrode during and between tissue applications. Failure to do so may result in unintended injury.
- Do not reuse any accessories labeled as SINGLE USE. These products are for single-use only and have not been designed to be re-used/re-sterilized. Reprocessing may lead to changes in material characteristics such as metallic corrosion and dulled edges, ceramic and plastic deformation or splitting which may impact the strength of the device and compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection. These risks may potentially affect patient safety.
- Re-sterilization by autoclave may cause hand control shorts resulting in inadvertent activation when the electrode is connected to the generator.
- **Ремонтные МН и Потребительские** of the instrument (ceramic insulator or metal active component) with fingers or instruments.
- Inadvertent activation or movement of an active VAPR electrode outside the field of vision may result in patient injury and damage to the electrode tip assembly.
- Using arthroscopic guidance, ensure that the electrode is completely surrounded by conductive irrigant solution during use. Failure to do so may result in electrode tip damage.
- Avoid bubble accumulation in the joint space during use. The accumulation of bubbles around the electrode diminishes performance and may produce overheating sufficient to damage adjacent structures or the electrode tip assembly.
- Do not wrap handpiece, footswitch or generator power cables around metal objects. Wrapping cables around metal objects may induce currents that could lead to electric shock, fire or injury to patient and/or surgical personnel.
- Operating the electrode at minimal suction rate or failure to attach the suction tubing to an adequate suction source may cause thermal injury to the physician or patient.
- Electrodes will wear from normal use, dependent on factors such as length of use, high tissue removal rate, prolonged use against bony surfaces, prolonged use in saline, high power settings, and use with minimal suction or fluid management. Periodically assess electrode tip for wear and proper operation. Replace the electrode if excessive wear is noted.
- Attach the suction adapter to standard hospital suction equipment. Ensure that the vacuum is in the range of 300 mmHg to 700 mmHg.
- Damage to surrounding tissue through iatrogenic injury could occur, as an effect of electrosurgery.
- Attempts to bend electrodes can result in electrode fracture or degradation of electrode performance.

PRECAUTIONS

- Prior to initial use, ensure that all package inserts including warnings, cautions, and Instructions for Use are read and understood. Safe and effective electrosurgery is dependent not only on equipment design but also on factors under the control of the user (refer to the user manual).
- Use caution when overriding the default power settings. Use the lowest power setting and the minimum tissue contact time necessary to achieve the appropriate surgical effect.
- When not in use, place the active electrode in a clean, dry, nonconductive, and highly visible area not in contact with the patient. Inadvertent activation while in contact with the patient may result in burns.
- Excessive force may damage the electrode tip assembly.
- Failure to maintain the recommended suction may result in device failure.

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BEFORE SURGERY

1. Switch on the VAPR VUE generator. Follow the instructions in the user manual. The generator will display: "Attach Instrument"
2. Connect a new, sterile electrode to the generator. Follow the instructions in the user manual. The generator will display the Ablate/Coag screen.
3. For optimal performance, safety, and convenience, the VAPR electrodes automatically preset the VAPR VUE generator to default output settings. The default settings are as follows:

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Caution should be used when overriding the default power settings. Use the lowest power setting and the minimum tissue contact time necessary to achieve the appropriate surgical effect. If power settings have been overridden to non-default values, the adjusted values will be retained when unplugging and re-plugging the same electrode, when connecting a new electrode of the same type, or when connecting a different electrode model with the same default values. If the generator is switched off and then on, a newly attached electrode will always revert to its default values. Always confirm proper, desired power settings before proceeding with surgery.

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BASIC-ADVANCED equivalency table for the VAPR TRIPOLAR 90 SUCTION Electrode:

VAPR TRIPOLAR 90 SUCTION Electrode			
ABLATE		COAGULATE	
Basic	Advanced (W)	Basic	Advanced (W)
1	10	1	10
2	100	2	40
3	190	3	80
4	220	4	100
5	260	5	130
6	280	6	160
7	300		
8	320		
9	340		
10	380		

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NOTE: When an instrument is attached, switching between ADVANCED and BASIC User Mode after adjusting the power output settings will reset the power output values to the default settings for the specific electrode type. The Mode button on the generator must be pressed to confirm use.

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DURING SURGERY

1. Prior to insertion of the electrode, close suction outflow of the electrode by moving the Flow Control Valve as far proximal as possible.
2. Insert the electrode into a fluid-filled joint space through a prepared entry point. Ensure that the electrode is surrounded by conductive irrigant solution when in use.
3. Ensure that the flow control valve is open (fully distal).

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4. To activate the ABLATE mode using the hand controls on the electrode handle (ensure suction is turned on), press and hold any of the YELLOW buttons. The ABLATE power setting on the generator will flash and a tone will be heard. Always keep the electrode in motion and in contact with tissue. Refer to the VAPR VUE User Manual on how to activate using the footswitch (wired or wireless). The hand controls will be de-activated when using the footswitch.
5. Refer to the VAPR VUE User Manual on how to adjust the ABLATE and COAG power settings on the generator panel. Use the lowest power possible to achieve the desired tissue effect.
NOTE: ABLATE power settings can also be adjusted using the footswitch (wired or wireless). Please refer to the VAPR VUE User Manual.
6. To activate the COAG mode using the hand controls on the electrode handle, press and hold any of the BLUE buttons. The COAG power setting on the generator will flash and a tone (at a different frequency than during ABLATE) will be heard. Always keep the electrode in contact with tissue. Refer to the VAPR VUE User Manual on how to activate using the footswitch (wired or wireless). Use the lowest power possible to achieve the desired tissue effect.

NOTE: COAG power settings are adjusted only using the buttons on the generator. Please refer to the VAPR VUE User Manual.

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7. Adjust the flow rate, using the Flow Control Valve, to optimize the fluid removal rate and bubble extraction and to assist with bleeding blood vessels. Sliding the Flow Control Valve distally will increase flow (up to ~500ml/min), moving the Valve proximally will decrease flow.
8. Reduce suction when device is not in use.

CAUTION: Inspect the active tip of the electrode for any damage after all faults. If the active tip is loose or dislodged, insert a new electrode.

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AFTER SURGERY

1. Discontinue suction and remove the electrode from the joint space.
2. Disconnect from suction source. Disconnect the electrode cable from the generator by steadying the generator and pulling on the plug body. DO NOT pull on the cable as this can cause damage to the generator.
3. The generator can be left in idle mode between surgical procedures. Turn off the generator at the conclusion of a surgical day.
4. Dispose of the single use components according to hospital procedure.

HOW SUPPLIED

The VAPR TRIPOLAR 90 SUCTION Electrode is provided sterile unless the package is opened or damaged. It is intended to be disposed of after a single use and is not intended to be re-sterilized (see warning section).

STORAGE

Store below 25°C (77°F) and away from moisture and direct heat. Bring to ambient temperature before unpacking / use. Does not use after USE BY date.

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VAPR® ARCTIC™ SUCTION ELECTRODES

**225029 – VAPR ARCTIC Suction
Electrode with Dome Tip**

**225030 - VAPR ARCTIC Suction
Electrode with Chisel Tip**

**For use only with the VAPR VUE®
Radiofrequency System**

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VAPR® ARCTIC™ Suction Electrodes

Instructions for Use

Please read all information carefully, in particular, the user manual provided with the VAPR VUE® Radiofrequency (RF) System (IFU-110009), prior to use. The components of the VAPR VUE RF system are designed for use as a single system. Failure to follow instructions may lead to improper functioning of the device and cause electrical or thermal injury.

DESCRIPTION

System Description

The VAPR VUE RF System is designed for arthroscopic surgical procedures. The system consists of a high frequency electrosurgical generator, and may include a reusable handpiece with connector cable and disposable electrodes, and alternative disposable electrodes with integrated connector cable. Other available components include hand controls, a footswitch, and a wireless footswitch alternative. The components are designed and intended for operation as a single unit. Use only VAPR electrodes and accessories with VAPR VUE RF system.

Device Description

The VAPR Electrodes are soft tissue ablation and hemostasis/coagulation devices intended for use with the VAPR VUE RF System. The VAPR ARCTIC Suction Electrodes are single use, one-piece bipolar suction articulating electrodes for use in arthroscopic surgery of the hip. The electrodes have a hand-controlled articulating tip to improve access and suction capabilities that will enhance the efficiency of the electrode and extend the utility of the system by assisting in the removal of bubbles and debris created during activation within the operating site.

INDICATION FOR USE

The VAPR ARCTIC Suction Electrode for use with the VAPR VUE RF System is intended for resection, ablation and excision of soft tissue, hemostasis of blood vessels and coagulation of soft tissue in patients requiring arthroscopic surgery of the hip.

INTENDED USE

The VAPR ARCTIC Suction Electrodes are intended for use in a surgical setting by personnel trained in arthroscopy.

CONTRAINDICATIONS

The VAPR ARCTIC Suction Electrodes are contraindicated in:

- Any non-arthroscopic surgical procedure.
- Any arthroscopic procedure that does not use a conductive irrigant solution (saline or Ringer's Lactate).
- Patients with heart pacemakers or other electronic device implants.
- Patients for whom arthroscopic procedures are contraindicated for any reason.

WARNINGS

- The operator should be experienced in arthroscopic surgical techniques.
- Please refer to the user manual for step-by-step instructions regarding the assembly and initial system check of the VAPR VUE RF system.
- As with all electrosurgical devices, do not use in the presence of flammable anesthetics, oxidizing gases or other flammable substances, as an electrosurgical device has the potential for providing a source of ignition.
- Ensure that fluid inflow and outflow is adequate and that the electrode is activated only when surrounded by conductive irrigant solution (e.g., saline or Ringer's Lactate). The use of RF energy with inadequate irrigation may overheat the fluid enough to cause skin burns at or near the access site and may result in damage to the electrode tip assembly.
- Do not touch the electrode tip when applying power.
- During an electrosurgical procedure, do not allow the patient to make direct contact with grounded metal objects. Inadvertent contact between the patient and grounded metal objects when the VAPR system is activated may result in burns.
- Observe extreme caution when using electrosurgery in close proximity to, or in direct contact with, any metal objects. The majority of arthroscopes and arthroscopic instruments are metal. Do not activate the electrode while any portion of the electrode tip is in contact with metal objects or instruments; avoid ablating tissue that is trapped between the arthroscope and the electrode. Регистрация МИ в Росздравнадзоре
www.evacert.ru | info@evacert.ru extended injury to the patient or surgical personnel and/or damage to the electrode and/or other equipment.
- Do not reuse any accessories labeled as SINGLE USE. These products are for single-use only and have not been designed to be re-used/re-sterilized. Reprocessing may lead to changes in material characteristics such as metallic corrosion and dulled edges, ceramic and plastic deformation or splitting which may affect the strength of the device and compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection. These risks may potentially affect patient safety.
- Avoid touching the distal tip of the instrument (ceramic insulator or metal active component) with fingers or instruments.
- Inadvertent activation or movement of an active VAPR electrode outside the field of vision may result in patient injury and damage to the electrode tip assembly.
- Attempts to insert or remove the electrode when it is in the articulated state may result in injury to the patient or surgeon, or cause damage to the device.
- Using arthroscopic guidance, ensure that the electrode is completely surrounded by conductive irrigant solution during use. Failure to do so may result in electrode tip damage.
- Avoid bubble accumulation in the joint space during use. The accumulation of bubbles around the electrode diminishes performance and may produce overheating sufficient to damage adjacent structures or the electrode tip assembly.
- Do not wrap handpiece, footswitch or generator power cables around metal objects. Wrapping cables around metal objects may induce currents that could lead to electric shock, fire or injury to patient and/or surgical personnel.

- Failure to attach the suction tubing to an adequate suction source may cause thermal injury to the physician or patient.
- DO NOT activate the probe for unnecessary and prolonged periods, as irrigant overheating and unintended tissue damage may result. Prolonged probe activation at maximum settings while using the suction feature may cause the shaft and suction tubing to reach temperatures of up to 58°C.
- Electrodes will wear from normal use, dependent on factors such as length of use, high tissue removal rate, prolonged use against bony surfaces, prolonged use in conductive irrigant solution, high power settings, and use with minimal suction or fluid management. Periodically assess electrode tip for wear and proper operation. Replace the electrode if excessive wear is noted.
- Attach the suction adapter to standard hospital suction equipment. Ensure that the roller clamp is open, and vacuum is in the range of 300 to 700 mmHg.
- Damage to surrounding tissue through iatrogenic injury could occur, as an effect of electrosurgery.
- Do not use CP (COOLPULSE®) Mode to shrink or smooth heat-sensitive tissue such as articular cartilage. When using probes that operate in CP Mode, ablation may occur even when ABLATE power is set to the minimum value of 10 W, resulting in possible unwanted ablation and tissue removal. Always use COAG mode when tissue removal is to be avoided.

PRECAUTIONS

- Prior to initial use, read and understand all package inserts including warnings, cautions, and Instructions for Use. Safe and effective electrosurgery is dependent not only on equipment design but also on factors under the control of the user (refer to the user manual).
- Use caution when overriding the default power settings. Use the lowest power setting and the minimum tissue contact time necessary to achieve the appropriate surgical effect.
- When not in use, place the active Electrode in a clean, dry, nonconductive, and highly visible area not in contact with the patient. Inadvertent activation while in contact with the patient may result in burns.
- Excessive force may damage the electrode tip or articulating shaft assembly.
- These devices are intended to be articulated by means of the hand-control. Attempts to bend the electrodes without using the hand-control can result in electrode fracture, degradation of electrode performance, or patient injury.
- Failure to maintain the recommended suction may result in device failure, patient, or user injury.
- Do not apply excessive force to the handle or hand-control during articulation. Damage to the device may result.
- Use a cannula or slotted cannula for insertion. Avoid percutaneous insertion. Attempts to insert the device without a cannula may damage device or patient.
- Carefully insert and withdraw electrodes to avoid possible damage to the device and/or injury to the patient or surgical personnel.
- Ensure that the articulating shaft is straight during insertion and removal from the joint. Inserting or removing the shaft in the articulated position may damage the shaft or tip assembly and may cause injury to the patient.
- These electrodes are designed for ablation and/or coagulation only, and not for mechanical displacement of tissue through applied force, which may result in bent or detached electrodes.
- Monitor suction flow. If suction is blocked, discard the device.
- Use of pre-warmed irrigation solution with the VAPR ARCTIC electrode may result in tissue damage or thermal injury.
- Do not allow suction line to contact patient; this may cause a thermal injury.

INSTRUCTIONS FOR USE

- Read all instructions carefully. Failure to properly follow instructions may lead to electrical or thermal injury and cause improper functioning of the device.
- The operator should be experienced in arthroscopic surgical techniques.
- Please refer to the user manual for step-by-step instructions regarding the assembly and initial system check of the VAPR VUE RF system.

BEFORE SURGERY

1. Switch on the VAPR VUE generator. Follow the instructions in the user manual. The generator displays “Attach Instrument”.
2. Connect a new, sterile electrode to the generator. Follow the instructions in the user manual. The generator will display the ABLATE/COAG screen.
3. For optimal performance, safety, and convenience, the VAPR electrodes automatically preset the VAPR VUE generator to default output settings. The default settings are as follows:

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4. Use caution when overriding the default power settings. Use the lowest power setting and the minimum tissue contact time necessary to achieve the appropriate surgical effect. If power settings were overridden to non-default values, the adjusted values are retained when unplugging and re-plugging the same electrode, when connecting a new electrode of the same type, or when connecting a different electrode model with the same default values. If the generator is switched off and then on, a newly attached electrode will always revert to its default values. Always confirm proper, desired power settings before proceeding with surgery.

5. The generator display can be set up to display output power in Watts (Advanced mode) or in relative numbers (Basic mode). The equivalent values for Advanced and Basic mode vary by electrode type.

Table 2. Basic-Advanced equivalency table for VAPR ARCTIC electrodes

VAPR ARCTIC with Dome Tip				VAPR ARCTIC with Chisel Tip			
ABLATE		COAGULATE		ABLATE		COAGULATE	
Basic	Advanced (W)	Basic	Advanced (W)	Basic	Advanced (W)	Basic	Advanced (W)
1	10	1	10	1	10	1	10
2	100	2	40	2	120	2	40
3	150	3	80	3	160	3	80
4	200	4	100	4	200	4	100
5	260						

NOTE: In Advanced mode, output power can be adjusted in finer increments than in Basic mode.

NOTE: When an instrument is attached, switching between Advanced and Basic User Mode after adjusting the power output settings will reset the power output values to the default settings for the specific electrode type. The Mode button must be pressed to confirm use.

6. Attach the tubing connector to a suction source. Ensure that suction is off by closing the clamp on the electrode tubing.

DURING SURGERY

1. Close suction outflow of the electrode prior to insertion of the electrode.
2. Using either a cannula or slotted cannula, insert the electrode into a fluid-filled joint space through a prepared entry point. Ensure that the articulating shaft is straight during insertion and removal from the joint, and that the electrode is surrounded by conductive irrigant solution when in use.
3. Adjust the flow rate, by opening suction outflow, to optimize the fluid removal rate and bubble extraction.
4. Depress the control on the handpiece to change the angle at the distal end of the probe.
5. Reduce suction when device is not in use.
6. The hand-control can be used to facilitate straightening of the articulated section as required.

CAUTION:

Inspect the active tip of the electrode for any damage prior to and immediately after use. If the active tip is loose or dislodged upon inspection, replace electrode. In the case of device breakage, remove fragments from the operative site. If the electrode becomes covered with extensive tissue debris, use a new electrode.

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If the audio signal continues when the electrode is moved away from a metal object, inspect the electrode and/or handpiece. Discard damaged accessories immediately. Remove debris build-up across the distal tip area with an appropriate wipe, taking care not to apply force to the articulating portion of the device.

If the alarm sounds persistently, and there is no contact with a metal object, discard the electrode.

AFTER SURGERY

1. Discontinue suction and remove the electrode from the joint space.
2. Disconnect from suction source. Disconnect the electrode cable from the generator by steadying the generator and pulling on the plug body. DO NOT pull on the cable as this can cause damage to the generator.
3. The generator can be left in idle mode between surgical procedures. Turn off the generator at the conclusion of a surgical day.
4. Dispose of the single use components per hospital procedure.

HOW SUPPLIED

The VAPR ARCTIC suction electrodes are provided sterile unless the package is opened or damaged. They are intended to be disposed of after a single use and are not intended to be re-sterilized (see warning section).

STORAGE

Store below 25°C (77°F) and away from moisture and direct heat. Bring to ambient temperature before use. Do not use after USE BY date.

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© DePuy Synthes 2020. All rights reserved.
DePuy International Ltd.