REF

T401209S, T401210S, T401211S, T401212S, T401213S, T401214S

Assistant with StableSoft™

Asistent se StableSoft™

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Assistant met StableSoft™

Avustin ja StableSoft™

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zariadenie so StableSoft™

StableSoft™ 辅助装置

This page is intentionally left blank.

Assistant with StableSoft™

Table of Contents

English.....	5
Česky.....	8
Dansk.....	12
Nederlands	15
Suomi.....	19
Français (France)	22
Deutsch.....	25
Magyar	29
Italiano	32
Norsk.....	35
Polski.....	39
Português (Brasil)	43
Русский язык.....	47
简体中文	51
Español (España).....	55
Svenska	59
Türkçe	63
Slovenčina.....	67
繁體中文	71
Graphical Symbols & Definitions	75

This page is intentionally left blank.

English

Instructions for Use

Assistant with StableSoft™

Attachments for Terumo Stabilizer Arms with Quick Connect
Catalogue Numbers T401209S, T401210S, T401211S, T401212S,
T401213S, T401214S

CAUTION: Before use, please read carefully the instructions for use of Terumo Stabilizer Arms with Quick Connect (Product References T401161 and T401163) in addition to these instructions for use.

PRODUCT DESCRIPTION

The Assistant with StableSoft™ attachments with silicone coating are sterile single use accessory devices for use with the Terumo Stabilizer Arms with Quick Connect (Product References T401161 and T401163). These attachments provide cardiothoracic positioning and stabilization. The malleable metal design with textured silicone coating allows the surgeon to bend the attachments to conform to anatomic shape. The attachment is provided in a sterile package and is intended for single use only. Refer to Table 1 for the available versions of the attachments.

Table 1

Product Reference Number	Amount of fingers	Length of fingers
T401209S	3	Long
T401210S	3	Short
T401211S	2	Long
T401212S	2	Short
T401213S	3	Long
T401214S	3	Short


Figure 1 - Assistant with StableSoft™ 3L
(Product Reference: T401209S)


**Figure 2 - Assistant with StableSoft™ 3L Flex
(Product Reference: T401213S)**

INTENDED USE

The Assistant with StableSoft attachments are manual positioning devices used during cardiothoracic surgery intended for positioning and stabilizing of anatomical structures.

CONTRAINDICATIONS

The Assistant with StableSoft attachments are contraindicated for use in patients with local or systemic infection.

Do not position the Assistant with StableSoft on newly infarcted tissue, aneurysmal tissue, fragile tissue, or directly over a coronary artery.

WARNINGS

- This device is sterilized using gamma irradiation and is intended **FOR SINGLE USE ONLY. DO NOT REUSE. DO NOT RESTERILIZE. DO NOT REPROCESS. REPROCESSING MAY COMPROMISE THE STERILITY, BIOCOMPATIBILITY AND FUNCTIONAL INTEGRITY OF THE DEVICE.**
- If the packaging is damaged, do not use, as sterility of the device may have been compromised.
- Inspect the product thoroughly for damage prior to use. Do not use damaged product.
- Proper surgical procedure and technique is the responsibility of the attending medical professional (physician). Each surgeon must evaluate the appropriate use of this device, case by case, based upon medical training, experience, and the type of surgical procedure employed.
- Repeated bending of the malleable zone of the shaft may eventually cause the shaft to break.
- Excessive flexing of the malleable zone may cause breakage.

CAUTION

United States Federal Law restricts this device to sale by or on the order of a physician or other licensed practitioner.

INSTRUCTIONS FOR USE

1. Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru technique, remove the attachment from its packaging and remove the protective tubing from the attachment shaft. Then, pass the Assistant to the surgical field for attachment to the Terumo Stabilizer Arm with Quick Connect (Product References T401161 and T401163).

2. Attach the Assistant with StableSoft attachment to the Terumo Stabilizer Arm with Quick Connect by sliding the sleeve of the Terumo Stabilizer Arm with Quick Connect forward and insert the shaft of the attachment into the quick connect with the detent on the shaft (of the attachment) facing the flat area with logo on the sleeve of the Terumo Stabilizer Arm with Quick Connect (see Figure 3). Return the sleeve to its original position to lock the accessory device into the quick connect. Check for secure connection.


Figure 3 – Attachment of the Assistant with StableSoft™ to the Terumo Stabilizer Arm with Quick Connect

3. All attachments can be bent along the shaft and the fingers to conform to anatomic requirements. Product References T401213S and T401214S can also be bent along the base of the fingers.

WARNINGS

- Repeated bending of the malleable zone of the shaft may eventually cause the shaft to break.
 - Excessive flexing of the malleable zone may cause breakage.
4. Position the Assistant with StableSoft attachment as required. Turn the main handle of the Terumo Stabilizer Arm with Quick Connect clockwise until the gap between the compression washers closes. For further instructions refer to the instructions for use of the Terumo Stabilizer Arm with Quick Connect.

STORAGE

The Assistant with StableSoft attachments storage conditions:

- A clean, dust-free location out of direct sunlight and free from risk of high temperature, excessive humidity, and water leakage.
- Do not store the device near chemicals.

DISPOSAL

After use, dispose of product and packaging in accordance with hospital, administrative and/or local, state, federal and international laws and regulations.

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Terumo® is a registered trademark of TERUMO CORPORATION
Assistant with StableSoft™ is a trademark of Terumo Cardiovascular Systems Corporation
TERUMO AUSTRALIA PTY LTD
Terumo Corp.
Country of Origin: U.S.A.