ULTRAScore™ Focused Force PTA Balloon

ENGLISH

INSTRUCTIONS FOR USE

CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

Device Description

The ULTRAScore™ Focused Force PTA Balloon consists of a flexible, over-the-wire (OTW) catheter shaft with a semi-compliant balloon fixed at the distal end. For all balloon lengths, radiopaque markers delineate the working length of the balloon and aid in balloon placement. For balloon lengths of 100 mm and greater, two radiopaque markers are positioned on the distal portion of the balloon and one radiopaque marker is positioned on the proximal portion of the balloon to differentiate between the distal and proximal ends of the balloon. The catheter includes an atraumatic tip to facilitate advancement of the catheter to and through the stenosis. Two scoring wires, oriented 180° apart, provide focused force upon dilatation. The ULTRAScore™ Focused Force PTA Balloon is compatible with .014" or .035" guidewires, as denoted by the product labeling. The distal portion of the .014" guidewire compatible catheters is hydrophilically coated. The proximal portion of the catheter includes a female luer lock hub connected to the catheter with a guidewire lumen and an inflation lumen. Packaged with every product is a protective sheath that is positioned over the balloon and must be removed prior to use. A stylet is placed into the tip of the catheter. These products are not made with natural rubber latex.

Figure 1

The GeoAlign® Marking System is a non-radiopaque ruler on the catheter shaft measured from the distal tip. The GeoAlign® markings are designated on the catheter shaft by 1cm increment bands with an accuracy within ±1mm. The distance from the distal catheter tip is labeled in 10cm increments. Thicker bands denote the midway point (5cm) between the labeled distances. The GeoAlign® Marking System is designed to be used as a tool to externally measure the intravascular advancement and/or retraction of the catheter. This can provide an intravascular reference regarding the location of the distal tip of the catheter or an approximate intravascular length measurement between two points. The GeoAlign® Marking System may also facilitate geographic alignment of an adjunctive therapy that includes the same GeoAlign® Marking System.

NOTE: The GeoAlign® Marking System provides an approximation that may not be an exact representation of the actual distance traveled intravascularly and should be confirmed under fluoroscopy.


Figure 2

NOTE: The GeoAlign® Marking System includes non-radiopaque white markings and are designed to be utilized outside the sheath.

Indications for Use

The ULTRAScore™ Focused Force PTA Balloon is intended to dilate stenoses in the iliac, femoral, ilio-femoral, popliteal, infra-popliteal and renal arteries and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This device is also recommended for post dilatation of balloon expandable stents, self-expanding stents, and stent grafts in the peripheral vasculature.

Contraindications

The ULTRAScore™ Focused Force PTA Balloon is contraindicated:

- Where there is the inability to cross the target lesion with a guidewire
- For use in the coronary or neuro vasculature

Warnings

1. Contents supplied STERILE using ethylene oxide (EO). Non-Pyrogenic. Do not use if sterile barrier is opened or damaged. Do not reuse, reprocess or re-sterilize. Use the catheter prior to the "Use By" date specified on the package label.
2. This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminate period of time. The residue of biological material can promote the contamination of the device with

pyrogens or microorganisms which may lead to infectious complications.

3. Do not resterilize. After resterilization, the sterility of the product is not guaranteed because of an indeterminate degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.
4. To reduce the potential for vessel damage or difficulty in deflating, the inflated diameter and length of the balloon should approximate the diameter and length of the vessel just proximal and distal to the stenosis.
5. When the catheter is exposed to the vascular system, the location of the balloon should be confirmed while under high quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated. If resistance is met during manipulation, determine the cause of the resistance before proceeding. Applying excessive force to the catheter can result in tip or catheter breakage, catheter kink, or balloon separation.
6. Do not exceed the RBP recommended for this device. Balloon rupture or difficulty in deflation may occur if the RBP rating is exceeded. To prevent over pressurization, use of a pressure monitoring device is recommended.
7. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

Precautions

1. Carefully inspect the catheter prior to use to verify that catheter has not been damaged during shipment and that its size, shape and condition are suitable for the procedure for which it is to be used. Do not use if product damage is evident.
2. The ULTRAScore™ Focused Force PTA Balloon should only be used by physicians experienced in the performance of percutaneous transluminal angioplasty.
3. It is recommended to consider the use of anti-coagulants, anti-platelet agents, and/or vasodilators in conformance with the accepted standard of practice or institutional guidelines surrounding peripheral endovascular procedures.
4. For ULTRAScore™ .014" guidewire sizes only, in order to activate the hydrophilic coating, wet the ULTRAScore™ balloon and catheter with sterile saline or wipe the balloon catheter with sterile saline saturated gauze immediately prior to its insertion in the body. Do not wipe the balloon catheter with dry gauze.
5. The minimal acceptable sheath French size is printed on the package label. Do not attempt to pass the PTA catheter through a smaller size introducer sheath than indicated on the label.
6. Use the recommended balloon inflation medium (25% contrast medium/75% sterile saline solution). Never use air or other gaseous medium to inflate the balloon.
7. The ULTRAScore™ Focused Force PTA Balloon should be used with caution for procedures involving calcified lesions, stents or synthetic vascular grafts due to the abrasive nature of these lesions.
8. Fully evacuate the balloon prior to withdrawing the system. Larger sizes of ULTRAScore™ Focused Force PTA Balloons may exhibit slower deflation times.
9. If resistance is felt during post procedure withdrawal of the catheter through the introducer sheath, determine if contrast medium is trapped in the balloon with fluoroscopy. If contrast is present, push the balloon out of the sheath and then completely evacuate the contrast before proceeding to withdraw the balloon.
10. If resistance is still felt during post procedure withdrawal of the catheter, it is recommended to remove the balloon catheter and guidewire/introducer sheath as a single unit.
11. Do not continue to use the balloon catheter if the shaft has been bent or kinked.
12. For ULTRAScore™ .014" guidewire sizes only, prior to re-insertion through the introducer sheath, re-activate the hydrophilic coating, and clean the balloon catheter by wiping the balloon catheter with sterile saline saturated gauze and rinsing with sterile saline. Do not wipe the balloon catheter with dry gauze.
13. GeoAlign® Marking System is designed to be used as an additional reference tool to accompany the interventionalist standard operation procedure. The use of fluoroscopic imaging is recommended following positioning of the catheter to the target lesion and prior to balloon deployment.

Potential Adverse Reactions

The complications that may result from a peripheral balloon dilatation procedure include:

- Additional intervention
- Allergic reaction to drugs or contrast medium
- Aneurysm or pseudoaneurysm
- Arrhythmias
- Embolization
- Hematoma
- Hemorrhage, including bleeding at the puncture site
- Hypotension/hypertension
- Inflammation
- Occlusion
- Pain or tenderness
- Pneumothorax or hemothorax
- Sepsis/infection
- Shock
- Short term hemodynamic deterioration
- Stroke

Regina Martinez
Senior RA Specialist
24 February 2022

- Thrombosis
- Vessel dissection, perforation, rupture, or spasm

Directions for Use

Handling & Storage

Store in a cool, dry, dark place. Do not store near radiation or ultraviolet light sources. Rotate inventory so that the catheters and other dated products are used prior to the "Use By" date. Do not use if packaging is damaged or opened.

Equipment Required

- Contrast medium
- Sterile saline solution
- Luer lock syringe/inflation device with manometer (10 ml or larger)
- Appropriate introducer sheath and dilator set
- .014" or .035" guidewire
- Optional: Catheter stabilization device

Dilatation Catheter Preparation

1. Remove the balloon catheter from package. Verify the balloon size is suitable for the procedure and the selected accessories accommodate the catheter as labeled.
2. Remove the stylet and protective sheath by grasping the balloon catheter just proximal to the balloon and with the other hand, gently grasp the mandrel and sheath protector and slide distally off of the balloon catheter.

Precaution: While removing the protective sheath, care should be taken to avoid damaging the balloon catheter.

3. Prior to use, the air in the balloon catheter should be removed. To facilitate purging, select a syringe or inflation device with a 10 ml or larger capacity and fill approximately half of it with the appropriate balloon inflation medium (25% contrast medium/75% sterile saline solution). Do not use air or any gaseous medium to inflate the balloon.
4. Connect a stopcock to the balloon inflation female luer hub on the dilatation catheter.
5. Connect the syringe to the stopcock.
6. Hold the syringe with the nozzle pointing downward, open the stopcock and aspirate for approximately 15 seconds. Release the plunger.
7. Repeat step #6 two more times or until bubbles no longer appear during aspiration (negative pressure). Once completed, evacuate all air from the barrel of the syringe/inflation device.
8. In order to activate the coating, wet the balloon catheter with sterile saline or wipe the balloon catheter with saline saturated gauze immediately prior to its insertion in the introducer sheath.
9. Prepare the wire lumen of the catheter by attaching a syringe to the wire lumen hub and flushing with sterile saline solution.

Precaution: Do not wipe the balloon catheter with dry gauze.

Use of the ULTRAScore™ Focused Force PTA Balloon

1. Backload the distal tip of the ULTRAScore™ Focused Force PTA Balloon over the pre-positioned guidewire and advance the tip to the introduction site.
2. Advance the catheter through the introducer sheath/guide catheter and over the wire to the site of inflation. If the stenosis cannot be crossed with the desired dilatation catheter, use a smaller diameter catheter to pre-dilate the lesion to facilitate passage of a more appropriately sized dilatation catheter.

Note: If using a hydrophilic guidewire, ensure that it is kept hydrated with sterile normal saline at all times.

3. Position the balloon relative to the lesion to be dilated, ensure the guidewire is in place and slowly inflate the balloon with an inflation device.

Note: Balloon lengths 100mm and longer have two marker bands located at the distal portion of the working length of the balloon and one marker band at the proximal portion. The working length of the balloon is measured at the proximal edge of the proximal marker band to the distal edge of the most distal marker band.

Precaution: Do not inflate the balloon past its labeled rated burst pressure.

4. Apply negative pressure to fully evacuate fluid from the balloon. Confirm that the balloon is fully deflated under fluoroscopy and that no contrast is left in the balloon.
5. While maintaining negative pressure and the position of the guidewire, grasp the balloon catheter just outside the introducer sheath/guide catheter and withdraw the deflated dilatation catheter over the wire through the introducer sheath/guide catheter.

Optional: To ensure that the sheath does not move during the procedure, use a catheter stabilization device to stabilize the sheath at the access site.

Optional: When used as a location reference tool, utilize the GeoALIGN® Marking System to assist with measuring how far the catheter has advanced. The GeoALIGN® Marking System provides an approximation that may not be an exact representation of the actual distance traveled intravascularly and should be confirmed under fluoroscopy (Reference Figure 3).


Figure 3

Optional: When used for intravascular length measurement, advance or retract the catheter between Points A and B utilizing the GeoALIGN® Marking System to measure the difference. The GeoALIGN® Marking System provides an approximation that may not be an exact representation of the actual distance traveled intravascularly and should be confirmed under fluoroscopy (Reference Figure 4).


Figure 4

Optional: To facilitate repeat catheter alignment with an additional device with the GeoALIGN® Marking System, ensure the GeoALIGN® markings are at the same location outside the sheath as the initial device.

Note: Ensure the sheath does not move during the withdrawal of the initial device or the advancement of subsequent devices.

Balloon Reinsertion

Precaution: Do not continue to use the balloon catheter if the shaft has been bent or kinked.

Precaution: For ULTRAScore™ .014" guidewire sizes only, prior to re-insertion through the introducer sheath, re-activate the hydrophilic coating, and clean the balloon catheter by wiping the balloon catheter with sterile saline saturated gauze and rinsing with sterile saline. Do not wipe the balloon catheter with dry gauze.

1. Load the balloon catheter onto a guidewire.
2. Advance the balloon catheter over the pre-positioned guidewire to the introduction site and through the introducer sheath. If resistance is encountered, replace the previously used balloon catheter with a new balloon catheter.
3. Continue the procedure according to the "Use of the ULTRAScore™ Focused Force PTA Balloon" section herein.

Warning: After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

Warning: After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

Warranty

Bard Peripheral Vascular warrants to the first purchaser of this product, that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited, to repair or replacement of the defective product, at Bard Peripheral Vascular's sole discretion, or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product is not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

Regina Martinez
Senior RA Specialist
24 February 2022

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BAIRD | PERIPHERAL VASCULAR

Printed on: 23 Jun 2021, 01:14:33 pm; Printed by: 10308611.

Regina Martins
Senior RA Specialist
24 February 2022

Do Not Resterilize

Ne pas restériliser
Nicht erneut sterilisieren
Non rsterilizzare
No reesterilizar
Niet opnieuw steriliseren
Näo reesterilizar
Μην επαναστεριλιζετε
Mä ikke resteriliseres
Fär ej omsteriliseras
Ei saa steriloida uudestaan
Mä ikke resteriliseres
Nie należy sterylizować ponownie
Ürasterilizalmi tilos
Nesterilizujte
Tekrar Sterilize Etmeyiniz
Повторная стерилизация запрещена
請勿重複消毒
재방균하지 마십시오

Keep Dry

Conserver au sec
Trocken lagern
Mantenere asciutto
Manténgase seco
Droog bewaren
Manter seco
Δαμπτρίτ το προϊόν στεγνό
Opbevarer tørt
Förvaras tørt
Suojattava kosteudelta
Mä holdes tørt
Chronić przed wilgocią
Tartsa szárazon
Uchovávaťe v suchu
Kuru Yerde Tutunuz
Хранить в сухом месте
保持乾燥
건조한 상태로 보관

Not Made with Natural Rubber Latex

Fabriqué sans latex naturel
Nicht aus Naturgummlatex hergestellt
Non prodotto con lattice di gomma naturale
Este producto no se fabrica con látex de caucho natural
Niet gemaakt met latex van natuurlijk rubber
Näo contém látex de borracha natural
Δεν κατασκευάζεται από φυσικό χυμό από καύσι
Ikke fremstillet med naturlig gummlatex
Ej tillverkad med naturgummlatex
Valmistuksessa ei ole käytetty luonnongummlatexia
Dette produktet er ikke produsert med naturlig gummlatex
Produkt nie jest wywarzany z lateksu kauczuku naturalnego
Nem tartalmaz természetes gummlatexet
Není vyroben z přírodního pryžového latexu
Doğal Kauçuk Lateksden Üretilmemiştir
Не содержит натуральный каучуковый латекс
製造未採用天然膠乳
천연 고무 라텍스를 사용해 제조되지 않음

Non-Pyrogenic

Apyrogène
Nicht pyrogen
Apirogino
Apirigino
Niet-pyrogeen
Apirogénico
Μη πυρρογόνο
Ikke-pyrogen
Pyrogenfri
Pyrogeenlön
Pyrogenfri
Produkt apirogeny
Progenimenes
Nepyrogeeni
Progenit deñjir
Apyrogenico
無熱原
비발열성

Keep Away From Sunlight

Conserver à l'abri de la lumière du soleil
Von Sonneneinstrahlung fern halten
Tenere al riparo dalla luce solare
Mantener alejado de la luz solar
Uit de buurt van zonlicht houden
Manter afastado da luz solar
Φυλάσσειτ το προϊόν μακριά από το ηλιακό φως
Beskyttes mod sollys
Fär ej utsättas för solljus
Suojattava auringon valolta
Beskyttes mot sollys
Chronić przed promieniowaniem słonecznym
Napfénytől védve tartandó
Nevystavujte působení slunečního záření
Güneş Işığından Uzak Tutunuz
Оберечь от воздействия солнечных лучей
避免日光照射
직사광선에 노출시키지 마십시오

Authorized Representative in the European Community

Représentant autorisé au sein de la Communauté européenne
Berechtigtter Vertreter in der Europäischen Gemeinschaft
Rappresentante autorizzato nella Comunità Europea
Representante autorizado en la Comunidad Europea
Gerechtigde vertegenwoordiger binnen de Europese Gemeenschap
Representante autorizado na Comunidade Europeia
Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα
Autoriseret representant i Det Europæiske Fællesskab
Auktoriserad representant i EG
Valtuutettu edustaja Euroopan yhteisössä
Autorisert representant i EU
Autoryzowany przedstawiciel na terenie Unii Europejskiej
Bejegyztett képviselő az Európai Közösségekben
Autorizovaný zástupce pro Evropskou unii
Avrupa Topluluğu Yetkili Temsilcisi
Уполномоченный представитель в Европейском сообществе
歐洲共同體の授權代表
유럽 공동체 내의 인증받은 대리업체

Sterilized Using Ethylene Oxide

Sterilisé à l'oxyde d'éthylène
Mit Ethylenoxid sterilisiert
Sterilizzato con ossido di etilene
Esterilizado mediante óxido de etileno
Gesteriliseerd met behulp van ethyleenoxide
Esterilizado por óxido de etileno
Αποσπορευτικό με αιθυλενοξείδιο
Steriliseret ved brug af ethylenoxid
Steriliserad med etylenoxid
Steriloitu etyleenoksidilla
Steriliseret med etylenoksid
Produkt sterylizowany tlakiem etylenu
Etikín-oxidál sterilizálva
Sterilizováno etylenoxidem
Etikén Oxid Küllantarak Sterilizte Edilmişdir
Стерилизовано этиленоксидом
使用聚氧乙烷消毒
산화에틸렌으로 멸균 처리됨

Manufacturer

Fabricant
Hersteller
Produttore
Fabricante
Fabrikant
Fabricante
Κατασκευαστής
Producent
Tillverkare
Valmistaja
Produsent
Produsent
Gyártó
Výrobce
Üretici
Производитель
製造商
제조사

Do Not Use If Package Is Damaged

Ne pas utiliser si l'emballage est endommagé
Nicht verwenden, wenn die Verpackung beschädigt ist
Non utilizzare se la confezione è danneggiata
No usar si el envase está dañado
Niet gebruiken wanneer de verpakking beschadigd is
Näo utilizar se a embalagem estiver danificada
Μη χρησιμοποιείτε το προϊόν εάν η συσκευασία του έχει υποστεί ζημιά
Mä ikke bruges, hvis emballagen er beskadiget
Använd inte produkten om förpackningen skadats
Ei saa käyttää, jos pakkaus on vaurioitunut
Skal ikke brukes hvis pakningen er skadet
Nie stosować, jeśli opakowanie jest uszkodzone
Ne használja fel, ha a csomagolás sérült
Nepoužívejte, pokud je obal poškozen
Paket Hasarlıysa Kullanmayınız
Не использовать, если упаковка повреждена
如果包裝已受損，請勿使用
포장이 손상되어 있다면 사용하지 마십시오

Copyright © 2018 C. R. Bard, Inc. All rights reserved.

Copyright © 2018 C. R. Bard, Inc. Tous droits réservés.
Copyright © 2018 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2018 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright © 2018 C. R. Bard, Inc. Todos los derechos reservados.
Copyright © 2018 C. R. Bard, Inc. Alle rechten voorbehouden.
Copyright © 2018 C. R. Bard, Inc. Todos os direitos reservados.
Πνευματικό δικαιώματα © 2018 C. R. Bard, Inc. Με την επιφύλαξη παντός δικαιώματος.
Copyright © 2018 C. R. Bard, Inc. Alle rättigheter förbehålls.
Copyright © 2018 C. R. Bard, Inc. Alla rättigheter förbehållna.
Copyright © 2018 C. R. Bard, Inc. Kaikki oikeudet pidätetään.
Copyright © 2018 C. R. Bard, Inc. Med enerett.
Copyright © 2018 C. R. Bard, Inc. Wszystkie prawa zastrzeżone.
Copyright © 2018 C. R. Bard, Inc. Minden jog fenntartva.
Copyright © 2018 C. R. Bard, Inc. Všechna práva vyhrazena.
Tekst Hakkı © 2018 C. R. Bard, Inc. Tüm hakları saklıdır.
Авторские права © 2018 C. R. Bard, Inc. Все права защищены.
Copyright © 2018 C. R. Bard, Inc. 保留所有權利。
Copyright © 2018 C. R. Bard, Inc. 모든 권리는 당사자가 소유합니다.

Bard and UltraScore are trademarks and/or registered trademarks of C.R. Bard, Inc. or an affiliate.

Bard et UltraScore sont des marques commerciales et/ou déposées de C. R. Bard, Inc. ou d'une de ses filiales.
Bard und UltraScore sind Marken und/oder eingetragene Marken von C. R. Bard, Inc. oder einer Tochtergesellschaft.
Bard e UltraScore sono marchi commerciali e/o registrati di C.R. Bard, Inc. o di una sua affiliata.
Bard y UltraScore son marcas comerciales y/o marcas comerciales registradas de C.R. Bard, Inc. o de una filial.
Bard en UltraScore zijn handelsmerken en/of gedeponeerde handelsmerken van C. R. Bard, Inc. of een geaffilieerd bedrijf.
Bard e UltraScore são marcas comerciais e/ou marcas comerciais registradas da C.R. Bard, Inc. ou de uma filial.
Το Bard και UltraScore είναι εμπορικά σήματα και/ή εμπορικά κατασκευασμένα της C.R. Bard, Inc. ή συνδεδεμένης εταιρείας της.
Bard og UltraScore er varemerker og/eller registrerede varemerker tilhørende C.R. Bard, Inc. eller et tilknyttet selskab.
Bard och UltraScore är varumärken och/eller registrerade varumärken som tillhör C.R. Bard, Inc. eller ett affilierat företag.
Bard ja UltraScore ovat C.R. Bard:n, Inc:n tai sen tytäryhtiön tavaramerkkejä jätai niiden rekisteröimä tavaramerkkejä.
Bard og UltraScore er varemerker og/eller registrerte varemerker som tilhører C.R. Bard, Inc. eller et tilknyttet selskap.
Bard i UltraScore są znakami towarowymi i/lub zastrzeżonymi znakami towarowymi firmy C.R. Bard, Inc. lub spółek zależnych.
A Bard és a UltraScore a C.R. Bard, Inc. vagy valamely leányvállalata védjegyei és/vagy bejegyzett védjegyei.
Bard a UltraScore jsou obchodní známky a/nebo registrované obchodní známky společnosti C.R. Bard, Inc. nebo jejich poboček.
Bard ve UltraScore C. R. Bard, Inc. yiketinini veya bağlı kuruluşunun ticari markaları ve/veya tescilli ticari markalarıdır.
Bard e UltraScore — торговые знаки и/или зарегистрированные торговые знаки компании C. R. Bard, Inc. или ее дочерних компаний.
Bard og UltraScore 是 C. R. Bard, Inc. 或其附屬機構的商標和/或註冊商標。
Bard 및 UltraScore는 C. R. Bard, Inc. 또는 제후사의 상표 및/또는 등록 상표입니다.

Regina Moulton
Senior RA Specialist
24 February 2022

 **Manufacturer:**
Bard Peripheral Vascular, Inc.
1625 West 3rd Street
Tempe, AZ 85281
USA
TEL: 1-480-894-9515
1-800-321-4254
FAX: 1-480-966-7062
1-800-440-5376
www.bardpv.com


BARD | **PERIPHERAL VASCULAR**

CE **2797** **EC REP**
Authorized Representative
in the European Community
BD Switzerland Sàrl
Terre Bonne Park – A4
Route de Crassier 17
1262 Eysins – Switzerland

BARD

PK1374600 Rev. 3 11/19

Printed on: 23 Jun 2021, 01:14:33 pm; Printed by: 10308611.

BARD PERIPHERAL VASCULAR

Regina Maltre
Senior RA Specialist
24 February 2022
7 pages

BARD
Peripheral Vascular
1000 West 10th Street
Tulsa, OK 74103
USA
TEL: 1-800-824-6578
FAX: 1-800-824-6579
www.bard.com

BARD
Peripheral Vascular
1000 West 10th Street
Tulsa, OK 74103
USA
TEL: 1-800-824-6578
FAX: 1-800-824-6579
www.bard.com