## Recommended way to use an oral dosing applicator in combination with a bottle adapter

1. Wash your hands carefully.
2. Open the bottle with medicine.
3. Take the oral syringe out of the packaging.

**Warning:** *Do not use a needle with the oral syringe.*

**Note:** *Use the oral syringe only for the medicine with which it is intended to be used and this only for one sales unit.*

**Warning:** *Mix-up of different medicines may cause problems.*

4. If the oral syringe is delivered together with an adapter: Separate the adapter from the oral syringe .

Insert the adapter into the bottle neck.

Ensure it is well fixed.

5. Take the oral syringe and place the tip in the central adapter opening.
6. Turn the bottle slantwise upside down.

Hold the bottle and the oral syringe barrel in one hand.

Use the other hand to fill the oral syringe.

7. Fill the oral syringe with a small amount of solution by pulling the plunger down.

Then push the plunger upward in order to remove any possible bubble. Tiny bubbles may be normal.

Pull the plunger down to the graduation mark of the dosage indicator corresponding to the quantity in millilitres (mL) or other described dosing unit in the leaflet of the medicine.

**Note:** *Check if the graduation of the oral syringe used, corresponds with the measure units described in the leaflet of the medicine.*

**Warning:** *Use of the wrong oral syringe may lead to the administration of the wrong dose of medicine.*

**Note:** *Avoid to get distracted when measuring the correct dose and when administering the medicine.*

**Warning:** *Scattering can lead to the administration of the wrong dose of medicine.*

8. Turn the bottle the right way up.
9. Twist and remove the oral syringe from the adapter.
10. Empty the contents of the oral syringe as described in the leaflet of the medicine: oral administration or in a glass of water or baby's bottle, by pushing the plunger to the bottom of the oral syringe. Product remaining in the tip of the oral syringe is not intended to be administered.

## Recommended way to use an oral dosing applicator in combination with a bottle adapter

11. For oral administration: Place the tip of the oral syringe into the patient's mouth and point the tip of the oral syringe towards either cheek. Push slowly on the plunger, a **small amount at a time**, until all the solution/ mixture in the oral syringe is given.

**Note:** *Be careful not to harm the lips or the inside of the mouth when placing the tip of the oral syringe into the patient's mouth for administration.*

**Warning:** *This could lead to physical harm.*

12. Close the medicine bottle.

13. Wash the oral syringe immediately after use with water. Use lukewarm or warm water. When not sufficient to clean the oral syringe, use a mild detergent and rinse afterwards with plenty of water. For cleaning thoroughly, it is preferable to take the plunger out of the barrel of the oral syringe before rinsing/cleaning both parts and to place the plunger back into the barrel after drying.

**Warning:** *Do not use bleach or an aggressive cleaning agent.*

### Notes:

- Don't use an oral syringe that blocks.

**Warning:** *Blocking during use can lead to inaccurate dosing. Using brute force during administration to try to unblock the oral syringe can cause high flow of the medicine which could lead to choking of the patient.*

- Don't sterilize the oral syringe.

**Warning:** *The use of a microwave or steam can impact the oral syringe negatively.*

- Serious incidents should be reported to the manufacturer and the local authority

### Storage of the oral syringe:

Keep the oral syringe away from little children.

Protect the oral syringe from sunlight and heat.

Store the oral syringe under clean, cool and dry conditions in the original packaging of the medicine.

### Disposal of the oral syringe

The oral syringe is made of plastic. After rinsing the oral syringe with plenty of water, it can be disposed according to the applicable legislations. The oral syringe must be disposed of at the same time as the medicine.

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