Valleylab™

Monopolar Footswitch • Pédale de commande monopolaire • Monopolær fodkontakt •
Monopolar-Fußschalter • Pedal monopolar • Monopolaarinen jalkakytin • Pedaliera
monopolare • Monopolaire voetschakelaar • Monopolar fotbryter • Pedal de acionamento
monopolar • Monopolär fotomkopplare • Педаль управления монополярным
режимом • **单极脚控开关**

 REF E6008, E6008L

Instructions for use • Mode d'emploi • Brugsanvisning • Gebrauchsanweisung •
Instrucciones de uso • Käyttöohjeet • Istruzioni per l'uso • Gebruiksaanwijzing •
Bruksanvisning • Instruções de utilização • Bruksanvisning •

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E6008L Monopolar Footswitch

1. Indications for Use

This product is intended to activate monopolar electrosurgical energy from the generator to an instrument to be used with Covidien's line of electrosurgical generators.

2. Cable Length

The cable length is 15' (4.6m).

3. Compatible Equipment

The E6008/E6008L is compatible with Covidien or other electrosurgical energy platforms/generators with footswitch capabilities.

May be used with the following energy platforms or generators:

- VLFT10GEN SW v1.0 or later (requires adapter 1017577 for Monopolar 2 receptacle)
- ForceTriad™ SW v1.0 or later
- VLFX8GEN SW v1.0 or later (requires adapter 1017577 for Monopolar 2 receptacle)
- Force FX series
- VLSURGEN

4. Handling and Operating Conditions

Table 1. Environmental Conditions

	Operation and storage	Transport
Ambient temperature range	10 to 40 °C (50 to 104 °F)	-30 to 60 °C (-22 to 140 °F)
Relative humidity	15 to 90% non-condensing	15 to 90% non-condensing

5. Clinical Benefits

The electrosurgical footswitches provide the ability for hands-free activation of an energy device through the energy platform it is connected with. This allows the surgeon the flexibility to choose the activation method that best suits their surgical needs.

6. General Warnings and Cautions

Note: Read all instructions, warnings, and cautions before performing electrosurgery using Covidien equipment or accessories.

Warning: Connect accessories to the proper receptacle. Improper connection may result in inadvertent accessory activation or other potentially hazardous conditions.

Warning: Connect only Covidien-approved footswitches to the Covidien™ energy platform/generator. Using footswitches from other manufacturers may cause equipment malfunction.

Warning: Inspect all accessory cords and connections for breaks, cracks, nicks, and other damage before use. Failure to observe this warning may result in injury to the patient or surgical team.

Caution: Position all cords to avoid contact with the footswitch.

7. Before Surgery

7.1. Setup

Attach the connector plug on the two-pedal footswitch to the appropriate receptacle on the energy

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8. During Surgery

8.1. Instructions

Activate the electrosurgical instrument by stepping on the appropriate CUT (yellow) or COAG (blue) pedal.

9. After Surgery

9.1. Instructions

Disconnect the footswitch from the energy platform/generator.

10. Maintenance (optional)

Warning: Always unplug the footswitch from the energy platform/generator before wiping.

If surface maintenance is required by hospital protocol, the exterior of the footswitch may be wiped with a cloth that has been dampened with a mild detergent or mild disinfecting solution.

11. Product Life Expectancy

The product life has been tested to the equivalence of 1,000 surgical procedures.

12. Disposal

Electrical appliances that are incorrectly disposed in dumps or landfills can leach dangerous substances causing contamination of soil and groundwater, and damaging the environment.

Do not dispose of this device in unsorted municipal waste. Use separate collection facilities according to your facility's policy for electrical appliances

Note: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and necessary regulatory authority (if EU, the competent authority of the Member State) in which the user and/or patient is established.

13. Residual Risk Summary

While every attempt has been made to reduce patient and user risks, all surgeries using this accessory carry some residual risk, even when used by trained physicians. The potential adverse events associated with the use of this footswitch includes, but is not limited to, the following risks:

- Burn
- Fall

14. California Proposition 65 Statement

For information regarding California Proposition 65, please refer to www.covidien.com/caprop65.

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O: indicates that the content of the toxic and hazardous substance in all the Homogenous Materials of the part is below the concentration limit requirement as described in GB/T 26572.

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X: indicates that the content of the toxic and hazardous substance in at least one Homogeneous Material of the part exceeds the concentration limit requirement as described in GB/T 26752 (Information provided in this table is solely for compliance with SJ/T 11364 and the Management Measure for the Control of Pollution by Electronic Information Products.)

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG. TM brands are trademarks of their respective owners. Other brands are trademarks of a Covidien company.

May be covered by U.S. Patents: www.covidien.com/patents

Manufactured for Covidien llc
15 Hampshire Street, Mansfield, MA 02048 USA
www.covidien.com

 Medtronic B.V.,
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Made in USA.

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EC REP Tyco Healthcare UK Ltd., Gosport, PO13 0AS, UK
© 2009 Valleylab All rights reserved.
Made in USA Printed in USA

Force Triad™ Three-Pedal Footswitch REF FT6003

BR	Pedal de acionamento triplo
DE	Drei-Pedal-Fußschalter
ES	Conmutador de tres pedales
FR	Pédale de commande à triple pédale
IT	Interruttore a tre pedali
NL	Voetschakelaar met drie pedalen
SE	Fotomkopplare med tre pedaler
RU	Трехпедальный ножной выключатель
CN	三踏板脚控开关

Instructions for Use

Force Triad™

Three-Pedal Footswitch

for use with the ForceTriad Energy Platform

REF **FT6003** Three-pedal footswitch with 4.6 m (15 ft.) cord

Important

The FT6003 requires ForceTriad software version 3.00 or greater to enable the use of the center Valleylab-mode foot pedal.

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Always unplug the footswitch from the generator before cleaning.

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Connect only Valleylab-approved footswitches to the Valleylab electrosurgical generator. Using footswitches from other manufacturers may cause equipment malfunction.

Do not wrap the accessory cords or patient return electrode cords around metal objects. This may induce currents that could lead to shocks, fires, or injury to the patient or surgical team.

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Do not immerse the footswitch.

Caution

Position all cords to avoid contact with the footswitch pedals.

Setup

1. If you are using a two-pedal footswitch, remove the footswitch pigtail from the monopolar footswitch receptacle on the back of the ForceTriad energy platform.
2. Attach the three-pedal footswitch connector to the monopolar footswitch receptacle on the back of the ForceTriad energy platform.

3. Activate the electrosurgical instrument in the desired mode by stepping on the color-coded monopolar footswitch:

- Yellow—Cut mode
- White—Valleylab mode
- Blue—Coag mode

Notice

The Valleylab mode can be activated only if a valid Valleylab mode enabled device is inserted into the monopolar 1 output. Activation of the Valleylab mode pedal when a non-Valleylab mode enabled device is inserted into the monopolar 1 output will result in no RF output and a single audible tone.

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Valleylab™

REM Patient Return Electrode Cord and Clamp

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Part No. 1016717

COVIDIEN, COVIDIEN with logo and the Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG.

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© 2011 Covidien.

 Manufactured for Covidien LLC, 15 Hampshire Street, Mansfield, MA 02048 USA.
EC REP Covidien Ireland Limited, IDA Business & Technology Park, Tullamore.

Made in Taiwan. Printed in USA.

www.covidien.com

REV 02/2011

Rx
ONLY

Consult
instructions
for use

ME20

REF E0560 REM Patient Return Electrode Cord, Reusable
15' (4.6 m)

REF E0560E Cord Clamp, ERBE, Reusable
15' (4.6 m)


For use with E7509 REM Polyhesive™ Adult Cordless Patient Return Electrodes

Caution


Position surgical electrode cables to avoid contact with the patient or other leads.

Before Surgery

1. Inspect the cord prior to each use to make sure that the clamp, plug, and cable are not broken or damaged. Discard damaged cords.
2. Open the clamp by lifting the latch. Do not force the latch beyond 90° when in the lifted position.
3. Remove the release liner from the return electrode tab area, exposing the foil contacts.

4. Insert the tab completely into the clamp with the foil contacts facing up so that the entire foil contact area is covered by the clamp.
5. Lower the latch until it is level with the body of the clamp.

6. Apply the patient return electrode according to the manufacturer's instructions.
7. Connect the cord plug to the generator.

For Covidien generators: prior to inserting the cord plug into the REM™ generator, turn the generator to the ready (on) mode to ensure that the alarm system is operating properly. Refer to the generator instruction manual for additional information.

Connect the cord plug to the generator to correct the alarm condition (audible alarm silenced, alarm light off). Return the generator to the standby (off) mode until ready for use. If the alarm condition is not corrected, refer to "Correcting REM Alarm Conditions" in the REM Polyhesive Patient Return Electrode instructions for use.

For non-Covidien generators: refer to the generator instruction manual.

After Surgery

1. To avoid skin trauma when removing the return electrode, peel off slowly with one hand while supporting the underlying tissue with the other hand.
2. To release the clamp, lift the latch.

Cord Care

If cleaning is necessary, wipe the cord with a mild cleaning solution or blood dissolving detergent. Rinse the cord thoroughly by wiping with a clean, damp cloth. Do not immerse. Wipe the cord completely dry before reuse.

Do not autoclave. This process will damage the cord.

Detailed Cord Inspection Instructions

Follow your institution's procedure for periodic inspection of the cord by the hospital biomedical department.

1. Inspect the cable for nicks or cuts. Bend the cable back and forth at the entrance to the clamp. Look for broken strands or corrosion on the wire. Repeat this inspection at the plug end. If any of these problems are found, discard the cord.
2. Inspect the clamp and plug for cracks. Check to make sure that the pin on the plug is in place and firmly attached. Discard the cord if there are any cracks or a missing or loose pin.
3. Using a milliohm meter, attach one lead to the clamp end of the cord and the other lead to the plug end. Measure the resistance of the cord while flexing the cable where it enters the clamp and the plug. The resistance must be less than 900 milliohms. Discard the cord if the resistance is higher than 900 milliohms.

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**REF LS0300 Tissue Fusion
Footswitch, Purple**

for use with Covidien LigaSure systems:

- LigaSure Vessel Sealing Generator
- ForceTriad™ Energy Platform

Warning

Electric Shock Hazard Do not connect wet accessories to the generator.

The LS0300 instrument is intended for use ONLY with the Covidien LigaSure vessel sealing generator or ForceTriad energy platform. Use of this instrument with other Covidien generators or with generators produced by other manufacturers could result in injury to the patient or surgical team, or cause damage to the instrument.

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Connect accessories to the proper receptacle. Improper connections may result in inadvertent accessory activation or other potentially hazardous conditions.

Connect only Covidien-approved footswitches to the LigaSure vessel sealing generator or ForceTriad energy platform. Using footswitches from other manufacturers may cause equipment malfunction.

Position instrument cords to avoid contact with the patient or other leads. Do not wrap cords around metal objects. This may induce currents that could lead to shocks, fires, or injury to the patient or surgical team.

Inspect the instrument and cords for breaks, cracks, nicks, or other damage before use. If damaged, do not use. Failure to observe this caution may result in injury or electrical shock to the patient or surgical team or cause damage to the instrument.

Do not immerse the footswitch.

Notice


For the LigaSure instrument connected to the vessel sealing handset receptacle on the front panel, the connected single pedal footswitch activates only the vessel sealing output function.

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Instructions for Use

LigaSure Vessel Sealing Generator

1. Connect the purple connector plug on the single pedal footswitch to the purple receptacle on the back of the generator.


① LigaSure Receptacles

2. Activate the LigaSure instrument by stepping on the pedal.

ForceTriad Energy Platform

1. Connect the purple connector plug on the single pedal footswitch to the purple receptacle on the back of the generator

① LigaSure Receptacles

2. Activate the LigaSure instrument by stepping on the pedal.

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© 2011 Covidien.

 Manufactured for Covidien llc, 15
Hampshire Street, Mansfield, MA 02048 USA.
EC REP Covidien Ireland Limited,
IDA Business & Technology Park, Tullamore.

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
E0510

Laparoscopic Monopolar Cord
10' (3 m)

For use with a maximum voltage of 5600 V (peak).

The 10' (3 m) cord is for use with monopolar laparoscopic devices with shrouded 4 mm male electrical connections.

Footswitching

Do not use if package is opened or damaged.

- Open the package and remove the cord using aseptic technique.
- Remove the band from the cord.
- Attach the sterile cord to the sterile accessory. Refer to the instructions provided with the accessory for details on proper connection method.
- Refer to the generator instruction manual for proper procedures to connect the cord to the generator.
- Always disconnect the cord from the active accessory by pulling the connector. Do not pull on the cable. Improper disconnection may damage the cord.

Warning

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

Connect accessories to the proper receptacle. Improper connection of accessories may result in inadvertent accessory activation or other potentially hazardous conditions. Follow the directions provided with electrosurgical accessories for instructions on proper connection and use.

Electric Shock Hazard Always turn the electrosurgical generator off (STANDBY) when connecting/disconnecting accessories. Failure to observe this precaution or handling of these connections while the generator is activated may result in injury or electrical shock to the patient or operating room personnel.

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Adapters

The E0510 requires an adapter (Covidien E0502 series, or E0017) for connection to the generator. Ensure that the cord connector pin is fully inserted and secure in the adapter.

Caution

When using adapters, ensure that adapters and connectors are fully inserted and that no metal is exposed.

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Part No. 1012144

COVIDIEN, COVIDIEN with logo and the Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG. All other brands are trademarks of a Covidien company.

© 2011 Covidien. All Rights Reserved.

 Covidien Inc, 15 Hampshire Street,
Mansfield, MA 02048 USA.**EC REP** Covidien Ireland Limited,
IDA Business & Technology Park, Tullamore.

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Electric Shock Hazard Ensure that adapter and accessory cord are correctly and securely connected and that no metal is exposed.

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Instructions for Use

The E0502 Series Adapters connect the footswitching accessory receptacle of Covidien generators to the distal (cord) end of various footswitching accessories and sterile disposable cords. E0502 adapters may also be used with other electrosurgical generators.


The following is a partial list of accessories designed to connect with these adapters. All items are single-use disposables and are packaged sterile unless otherwise noted.

Catalog Number	Description
E0503	Monopolar Cord
E0510	Laparoscopic Monopolar Cord
E2003	Pencil, Reusable
E2250H	Pencil, Blade Electrode, Holster
E2504, E2504H	Pencil, Holster
E2505-10FR	Suction Coagulator
E2506-10FR	Suction Coagulator
E2520H	Argon Handset

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
Before Surgery


1. Ensure the generator is OFF.
2. Insert the adapter into the footswitching accessory receptacle of the generator.

3. Press the spring-loaded tab on the adapter cap to open the accessory receptacle.

4. Insert the brass terminal on the accessory cord so the spring-loaded tab presses the terminal against the inside of the adapter when the tab is released. Loose connections may result in arcing within the adapter or lack of function.

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
After Surgery


1. Turn OFF the generator before disconnecting the accessory and removing the adapter.
2. Press the spring-loaded tab on the adapter cap to open the accessory receptacle.
3. Remove the accessory cord from the adapter.

4. Remove the adapter from the generator receptacle if applicable.

Cleaning Recommendations

1. Remove all gross matter (blood, mucus, and tissue):
 - a. Dampen a soft cloth with a cleaning agent.
 - b. Thoroughly wipe all surfaces of the adapter.

2. Dry completely.

Notice

To avoid product damage, do not immerse the adapter in cleaning agent.

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Part No. 1012187

COVIDIEN, COVIDIEN with logo and the Covidien logo are U.S. and internationally registered trademarks of Covidien AG.

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
ME20

REF E0005-3C **Active Adapter**
4 mm

For use with all Covidien electrosurgical generators.


Introduction

The E0005-3C adapter allows the monopolar instrument receptacles on all Covidien electrosurgical generators to receive the 4 mm diameter pin connector from various accessories and cords.

The adapter, provided nonsterile, meets IEC601 specifications.

Before Surgery

1. Ensure that the electrosurgical generator is off.
2. Fully insert the adapter into one of the monopolar instrument receptacles on the front panel of the electrosurgical generator.
3. Fully insert the 4 mm diameter pin on the accessory or cord into the adapter.
4. Turn on the generator.

- ① Generator
② Adapter
③ 4 mm diameter connecting pin

Warning

Connect accessories to the proper receptacle. Improper connection of accessories may result in inadvertent accessory activation or other potentially hazardous conditions.

Electric Shock Hazard Connect the adapter to the electrosurgical generator only while the generator is off or in standby mode. Failure to do so may result in injury or electrical shock to the patient or operating room personnel.

Electric Shock Hazard Ensure that adapter and accessory cord are correctly and securely connected and that no metal is exposed. Contact with exposed metal parts may result in shocks or burns to surgical personnel or the patient. Loose connections or incomplete insertion of the pin connector may result in arcing inside the adapter or lack of function.

After Surgery


1. Turn off the electrosurgical generator.
2. Disconnect the accessory cord.
3. Remove the adapter from the generator receptacle.

Notice

Do not pull on the accessory cord to remove the adapter. Damage to the cord or adapter may result.

Cleaning and Decontamination

1. Remove all gross matter (blood, mucus, and tissue):
 - a. Dampen a soft cloth with a cleaning agent.
 - b. Thoroughly wipe all surfaces. Follow the procedures approved by your institution.

2. Clean all surfaces by wiping with a clean and damp cloth.

Notice

To avoid product damage, do not immerse the adapter in liquid.

3. Rinse the adapter thoroughly by wiping with a clean and damp cloth. Dry the adapter thoroughly.

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For use with a max peak voltage of 1000 V.

Intended Use

Cords are to electrically connect bipolar electrosurgical instruments to Covidien electrosurgical generators.

Footswitching

Warning

This product is non-sterile. Clean and sterilize cords before the first and every use.

DO NOT sterilize the cord using low temperature gas plasma, low temperature peracetic acid, dry heat, soaking, or irradiation. Failure to comply may result in patient injury and/or reduced product performance.

Cleaning

The cleaning efficacy has been validated for using either of the following cleaning methods.

Caution

Use only recommended cleaning solutions. Do not use solvents, lubricants, or other chemicals to clean the reusable cord. Some disinfection solutions are corrosive to the connectors.

Manual Cleaning

- Soak cords in an enzymatic cleaner, such as Prolystica® in accordance with the manufacturer's recommendations.
- Remove gross particulate with a soft-bristle brush or cloth.
- Clean the cord's instrument receptacle by inserting a small, soft-bristle brush to remove particulate.
- Rinse and dry with a lint-free cloth.

Mechanical Cleaning

Clean in a low impingement, Steris model Reliance® 444 single chamber washer/disinfector.

1. Pre-rinse in cold water for 30 seconds.
2. Wash with H-Klenz® II (sodium metasilicate, high alkaline), 4 minutes, 140° F (60° C).
3. Rinse with LabKlenz 250® (neutralization cycle, citric acid), 2 minutes, 150° F (66° C).
4. Rinse with water, 15 seconds, 140° F (60° C).
5. Dry, 15 minutes, 180° F (82° C).

6. Cool before autoclave cycle (15 min. - 30 min).

Sterilization

The sterility of the cord has been validated for the following parameters:

1. Wrap each reusable cord separately in a double layer of CSR sterilization wrap.

Notice

Separating reusable cords during sterilization maximizes the life of the cord. Heat from sterilization softens the cord insulation. Contact with other cords or instruments during sterilization can damage the reusable cord.

2. Use the following guidelines for sterilization.

Steam Sterilization - Wrapped

Temperature	Type	Sterilization Time
275° F (135° C)	Prevac	3 min.
273° F (134° C)	Prevac	3 min.
270° F (132° C)	Prevac	3 min.
250° F (121° C)	Prevac	15 min.
275° F (135° C)	Gravity	10 min.
270° F (132° C)	Gravity	10 min.
250° F (121° C)	Gravity	30 min.

Drying cords is not required. Refer to the documentation from the manufacturer of the sterilization system in use and hospital guidelines for specific dry times.

Steam Sterilization - Unwrapped (Flash)

Temperature	Type	Sterilization Time
270° F (132° C)	Prevac	4 min.
270° F (132° C)	Gravity	10 min.

Drying cords is not required. Refer to the documentation from the manufacturer of the sterilization system in use and hospital guidelines for specific dry times.

These instructions have been validated by Covidien as capable of preparing the E0020V, E0021S, and E0022W footswitching bipolar cords for re-use. It is the responsibility of the processor to ensure that sterilization is performed using equipment, materials, and personnel that will achieve the desired results. This requires validation and routine monitoring of the process. Any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences.

Product Life Expectancy

The product life has been validated for a 75-use cycle using the following cleaning and sterilization processes.

- Manual cleaning followed by vacuum-assisted (Prevac) steam sterilization at 280° F (138° C) for 30 minutes, and a 45-minute dry time
- Mechanical cleaning followed by vacuum-assisted (Prevac) steam sterilization at 280° F (138° C) for 30 minutes, and a 45-minute dry time

Assembly

Warning

Malfunction may be caused by a loose or incorrect connection of the cord to the generator or instrument. Worn cords cannot be repaired. Replace cords on a routine basis to maintain optimum reliability in operating room use.

Warning

Connect accessories to the proper receptacle. Improper connection of accessories may result in potentially hazardous conditions.

Inspect the cord insulation before and after each use for cracks, nicks, cuts, dents, or depressions which may decrease the insulation effectiveness.

Do not use the cord if there are cuts, nicks, or depressions in the insulation. Unintended thermal injury could occur.

Caution

Always disconnect the cord from the instrument or generator by pulling on the connector. Do not pull on the cord when disconnecting from the instrument and generator. Improper disconnection may damage the cord insulation.

Prior to assembly, the reusable cord and instrument should be cleaned and sterilized.

1. Inspect the cord for the following end-of-life indicators:
 - Cracks, nicks, cuts, dents or depressions in the cord insulation and connectors
 - Exposed wires
 - Bent or loose plugs
 - Loose receptacles in connectors

If any of these indicators are found, replace the damaged cord with a sterilized cord in good condition.

2. Connect the reusable cord to the instrument. Ensure that the cord receptacle is fully seated against the instrument's connector.
3. Connect the reusable cord to the bipolar output receptacle on the generator.

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Rx
ONLY

Part No. 1045433

COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG.

Polystica, Reliance H-Klenz, LabKlenz 250 are registered trademarks of STERIS Corp. Other brands are trademarks of a Covidien company.

© 2011 Covidien.

 Manufactured for Covidien llc, 15 Hampshire Street, Mansfield, MA 02048 USA.

 Covidien Ireland Limited, IDA Business & Technology Park, Tullamore.

Made in Germany. Printed in Germany.

www.covidien.com

REV 10/2011

ForceTriad™**Bipolar Resection Footswitch**

- en** Instructions for Use
- fr** Mode d'emploi
- pt** Instruções de uso
- de** Gebrauchsanleitung
- es** Instrucciones de uso
- it** Istruzioni per l'uso
- nl** Gebruiksaanwijzing
- sv** Bruksanvisning
- ru** Инструкция по применению
- zh** 使用说明

Part No. 1043111

COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG.

Other brands are trademarks of a Covidien company.

© 2011 Covidien.

 Manufactured for Covidien Ilc,
15 Hampshire Street, Mansfield, MA 02048 USA.
EC REP Covidien Ireland Limited,
IDA Business & Technology Park, Tullamore.

Made in USA. Printed in USA.

www.covidien.com

REV 07/2011

**Rx
ONLY**Consult
instructions
for use

ME20

0086

FT6009

Bipolar Resection Footswitch

with 15' (4.6 m) cord

For use with the ForceTriad energy platform

Important

The FT6009 requires ForceTriad energy platform software version 3.30 or greater and a Covidien-specified bipolar resection-enabled instrument and cord to enable the use of the bipolar resection function. For more information, see the notice on this page and the *Reference Documentation* section.

Warning

Electric shock hazard: Do not connect wet accessories to the generator.

Always disconnect the footswitch from the generator before cleaning.

Connect accessories to the proper receptacle. Improper connection may result in inadvertent accessory activation or other potentially hazardous conditions.

Connect only Covidien-approved footswitches to the ForceTriad energy platform. Using footswitches from other manufacturers may cause equipment malfunction.

Inspect all accessories, cords, and connections for breaks, cracks, nicks, and other damage before each use. Failure to do so may result in injury or electric shock to the patient or surgical team.

Caution

Position all cords to avoid contact with the footswitch.

Intended Use

The Bipolar Resection Footswitch is intended to be used with the ForceTriad Bipolar Resection function for endoscopically controlled removal (resection) or coagulation of tissue using 0.9% NaCl solution (saline) as the irrigation medium.

Setup

1. Attach the Bipolar Resection Footswitch connector to the LigaSure 2/Bipolar Resection footswitch receptacle on the back of the ForceTriad energy platform.
2. Activate a supported electrosurgical, bipolar resection instrument in the desired mode by stepping on the color-coded Bipolar Resection Footswitch:
 - Yellow (left)—Cut mode (resection)
 - Blue (right)—Coag mode (desiccation)

Notice

Bipolar resection function is supported only if a Covidien-specified bipolar resection instrument is connected to the LigaSure 2/ Bipolar Resection receptacle with a Covidien Bipolar Resection cord. See the *Reference Documentation* section on page 3. Attempted activation of the Bipolar Resection Footswitch with an unsupported device or cord results in no radio frequency (RF) output and a single audible tone.

Cleaning

Warning

Electric shock hazard: Before cleaning, turn off the power source and disconnect the footswitch.

Do not immerse the footswitch.

Notice

An enzymatic cleaner may be used in accordance with the manufacturer's recommendations for manual cleaning. Do not immerse the footswitch or clean it with abrasive cleaning or disinfectant compounds, solvents, or other materials that could damage the footswitch.

1. Disconnect the footswitch from the ForceTriad energy platform.
2. Thoroughly wipe all surfaces of the footswitch and cord with a damp cloth and mild cleaning solution or disinfectant.

The footswitch will withstand the effects of cleaning over time without degrading the enclosure.

Reference Documentation

For more information about using the bipolar resection application, refer to these documents. For specific supported cords and bipolar resection instruments, see the Cords documents.

- *ForceTriad Energy Platform Bipolar Resection Supplement*
- *ForceTriad Bipolar Resection Activation Kit Instructions for Use (FTBRKIT)*
- *ForceTriad Bipolar Resection Cord (FT0021S) Instructions for Use*
- *ForceTriad Bipolar Resection Cord (FT0022W) Instructions for Use*

Troubleshooting

For additional questions, contact a technical service representative:

Covidien Service	
Within the U.S.	1-800-255-8522 option 2
Outside of the U.S.	1-303-530-2300
Web Page	www.valleylab.com/valleylab/international/service-world.html
Email	technicalservice@covidien.com

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Assistenza Covidien

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Pagina Web www.valleylab.com/valleylab/international/serivce-world.html

E-mail technicalservice@covidien.com

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG. [™]** brands are trademarks of their respective owners. Other brands are trademarks of a Covidien company.

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Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Product of USA (from US and foreign materials).

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
E0017 Universal Monopolar Adapter

1. Indications for Use

The adapter is intended to be used with compatible accessories to remove tissue and control bleeding by use of high-frequency electrical current.

2. Compatible Equipment

The E0017 adapter connects the footswitching accessory receptacle of Valleylab™ generators to the distal (cord) end of various footswitching accessories and cords. This adapter is designed to be used with an accessory connector up to 6 mm in diameter.

For use with a maximum peak voltage of 5550 V.

Compatible with Valleylab™ and most electrosurgical generators for footswitching activation.

The following is a partial list of accessories designed to connect with these adapters.

Table 1. Compatible Equipment

Product code	Description
E0503	Monopolar Cord 10' (3 m)
E0510	Laparoscopic Monopolar Cord 10' (3m)
E2504	Pencil 10' (3 m)
E2504H	Pencil, Holster 10' (3 m)

Contact your Covidien representative for ordering information and a complete list of items to meet your operating room needs.

3. Handling and Operating Conditions

Table 2. Handling and Operating Conditions

Condition	Operation and storage	Transport
Ambient temperature range	10°C to 40°C (50°F to 104°F)	-30°C to 60°C (-22°F to 140°F)
Relative humidity	15% to 90% non-condensing	15% to 90% non-condensing

4. Clinical Benefits

The electrosurgical devices used in combination with compatible generators and accessories (adapters, footswitches, and cords) will provide monopolar and/or bipolar energy for cutting and coagulating soft tissue. This provides the capability for the surgeon to balance tissue dissection and hemostatic effects as needed during a given surgical case.

5. General Warnings and Cautions

Read all instructions, warnings, and cautions before performing electrosurgery using Covidien equipment or accessories.

5.1. Warnings

- Warning:** Before use, examine the electrosurgical unit and accessories for defects. Do not use cables or accessories with damaged (cracked, burned, or taped) insulation or connectors.
- Warning:** Connect accessories to the proper receptacle. Improper connection of accessories may result in inadvertent accessory activation or other potentially hazardous conditions.
- Warning:** Ensure that the adapter and accessory cords are correctly and securely connected and that no metal is exposed. Contact with exposed metal parts may result in shocks or burns to surgical personnel or the patient. Loose connections or incomplete insertion of the pin connector may result in arcing inside the adapter or lack of function.
- Warning:** Do not use electrosurgery in the presence of flammable anesthetics or oxidizing gases, such as nitrous oxide (N₂O) and oxygen (O₂) or in close proximity to volatile solvents (such as ether or alcohol), as explosion may occur.
- Warning:** Do not use with HF surgical-mode output settings that result in a peak output voltage that exceeds the rated accessory voltage.

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6. Before Surgery

6.1. Setup

1. Insert the adapter into the footswitching accessory receptacle of the generator.
2. Insert the accessory cord into the receptacle in the middle of the adapter. Adjust the tension screw if necessary. Ensure that no metal is exposed. Refer to [Figure 1](#).
3. Adjust the tension screw so the accessory cord is securely connected. Do not over tighten the tension screw; otherwise, accessory cord damage may result.

Note: A nylon cord may be inserted through the hole in the tab and secured to the generator to permanently attach the adapter.

7. After Surgery

Caution: Do not attempt to pull on the accessory cord to remove the adapter. Improper disconnection may result in damage to the cord insulation or adapter.

7.1. Instructions

1. Loosen the tension screw and remove the accessory cord from the adapter.
2. Remove the adapter from the generator receptacle if required.

Note: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and necessary regulatory authority (if European Union, the competent authority of the Member State) in which the user and/or patient is established.

8. Maintenance (Optional)

Caution: Always unplug the adapter from the energy platform/generator before wiping. If maintenance is required by hospital protocol, the exterior of the adapter may be wiped with a cloth that has been dampened with a mild detergent or mild disinfecting solution.

Note: To avoid product damage, do not immerse the adapter in liquid.

9. Product Life Expectancy

The product life has been tested for 10,220 usage cycles.

10. Disposal

Caution: Do not dispose of electrical appliances as unsorted municipal waste. Use separate collection facilities. Electrical appliances incorrectly disposed of in dumps or landfills can leach dangerous substances causing contamination of soil and ground water, and damaging the environment. Contact your local government or point of sale for information regarding the collection of waste electrical appliances.

11. Residual Risk Summary

While every attempt has been made to reduce patient and user risks, all surgeries using this instrument carry some residual risk, even when used by trained physicians. The potential adverse events associated with the use of adapters (accessories) includes, but are not limited to, the following risks:

- Airway obstruction
- Burn
- Electrical shock

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1. Indications for Use

This product is intended to activate bipolar electrosurgical energy from the generator to an instrument to be used with Covidien's line of electrosurgical generators.

2. Cable Length

The cable length is 15' (4.6m).

3. Compatible Equipment

May be used with Covidien™ electrosurgical energy platforms or generators with bipolar capabilities.

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Table 1. Environmental Conditions		
	Operation and storage	Transport
Ambient temperature range	10 to 40 °C (50 to 104 °F)	-30 to 60 °C (-22 to 140 °F)
Relative humidity	15 to 90% non-condensing	15 to 90% non-condensing

5. Clinical Benefits

The electrosurgical footswitches provide the ability for hands-free activation of an energy device through the energy platform it is connected with. This allows the surgeon flexibility to choose the activation method that best suits their surgical needs.

6. General Warnings and Cautions

Note: Read **all** instructions, warnings, and cautions before performing electrosurgery using Covidien equipment or accessories.
Warning: Connect accessories to the proper receptacle. Improper connections may result in inadvertent accessory activation or other potentially hazardous conditions.
Warning: Connect only Covidien-approved footswitches to the Covidien™ energy platform/generator. Using footswitches from other manufacturers may cause equipment malfunction.
Warning: Inspect all accessory cords and connections for breaks, cracks, nicks, and other damage before use. If damaged, do not use. Failure to observe this warning may result in injury or to the patient or surgical team.
Caution: Position all cords to avoid contact with the footswitch.

7. Before Surgery

7.1. Setup

Attach the connector plug on the single pedal footswitch to the appropriate receptacle on the energy platform/ generator.

8. During Surgery

8.1. Instructions

Activate the electrosurgical instrument by stepping on the pedal.

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10. Maintenance (Optional)

Warning: Always unplug the footswitch from the energy platform/generator before wiping.
If surface maintenance is required by hospital protocol, the exterior of the footswitch may be wiped with a cloth that has been dampened with a mild detergent or mild disinfecting solution.

11. Product Life Expectancy

The product life has been tested to the equivalent of 1,000 surgical procedures.

12. Disposal

Electrical appliances that are incorrectly disposed in dumps or landfills can leach dangerous substances causing contamination of soil and groundwater, and damaging the environment.
Do not dispose of this device in unsorted municipal waste. Use separate collection facilities according to your facility's policy for electrical appliances.
Note: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and necessary regulatory authority (if EU, the competent authority of the Member State) in which the user and/or patient is established.

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- Burn
- Fall

14. Compliance Statements

California Proposition 65 Statement
For information regarding California Proposition 65, please refer to www.covidien.com/caprop65.

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REF VLFT10ADP1

Compatible with:


VLFT10GEN software version 2.0 or greater

- en** Instructions for Use
- fr** Mode d'emploi
- pt** Instruções de uso
- da** Brugsanvisning
- de** Gebrauchsanleitung
- es** Instrucciones de uso
- fi** Käyttöohjeet
- it** Istruzioni per l'uso
- nl** Gebruiksaanwijzing
- no** Bruksanvisning
- sv** Bruksanvisning
- ru** Инструкция по применению
- zh** 使用说明

Valleylab™ FT10

[REF] VLFT10ADP1 **Bipolar Adapter**

For use with compatible instruments rated to handle a maximum peak voltage of 531 V.

Important

The Valleylab™ FT10 Bipolar Adapter requires FT series software version 2.0 or greater to enable use of compatible bipolar instruments.

Not made with natural rubber latex

Indications for Use

The Valleylab FT10 Bipolar Adapter is intended to be used with compatible accessories for cutting and coagulating tissue by use of high-frequency electrical current.

Warning

Ensure the Valleylab FT10 Bipolar Adapter and the connected bipolar instrument are fully inserted prior to use. Activation with a partially inserted connector or missing bipolar instrument may result in injury to the surgical team if the exposed connector is touched during use.

Do not connect the Valleylab FT10 Bipolar Adapter and bipolar instrument while depressing the associated footswitch. Inadvertent activation may occur.

Precaution

Do not use instruments with flying leads with the Valleylab FT10 Bipolar Adapter.

Important

The Valleylab FT10 Bipolar Adapter is not compatible with handswitching bipolar instruments.

Cleaning

Disconnect the Valleylab FT10 Bipolar Adapter from the generator and the connected bipolar instrument. Thoroughly wipe surfaces with a damp cloth and mild cleaning solution or disinfectant.

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Consult
instructions
for use

Part No. PT00047615

COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien. Other brands are trademarks of a Covidien company, TM* brands are trademarks of their respective owner.

© 2016 Covidien.

 Covidien llc,
15 Hampshire Street, Mansfield, MA 02048 USA.
 Covidien Ireland Limited,
IDA Business & Technology Park, Tullamore.

Made in USA. Printed in USA.

www.covidien.com

REV 08/2016

Valleylab™

Accessory Holster

- en** Instructions for Use
- fr** Mode d'emploi
- pt** Instruções de uso
- da** Brugervejledning
- de** Gebrauchsanleitung
- es** Instrucciones de uso
- fi** Käyttöohjeet
- it** Istruzioni per l'uso
- nl** Gebruiksaanwijzing
- no** Bruksanvisning
- sv** Bruksanvisning
- ru** Инструкция по применению
- zh** 使用说明

Not made with natural rubber latex


Do not use if package is opened or damaged

Indications for Use



The Valleylab Accessory Holster is intended to provide a safe place for electrosurgical accessories during the surgical procedure when not actively being used.

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1. Remove the accessory holster from the sterile package using aseptic technique.
2. Attach the accessory holster to the surgical drape with tape or by pulling the drape material through the slot as shown.
3. Discard after use. The accessory holster is not designed to withstand resterilization. **SINGLE USE ONLY – DO NOT RESTERILIZE.**

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

PT00115474

Part No. PT00115474

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company, TM* brands are trademarks of their respective owner.

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
Not made with natural rubber latex

Indications for Use

The Sonicision™ Curved Jaw Cordless Ultrasonic Dissection Device is indicated for soft tissue incisions when bleeding control and minimal thermal injury are desired. The Device can be used as an adjunct to or substitute for electrosurgery, lasers, and steel scalpels in general, plastic, pediatric, gynecologic, urologic, exposure to orthopedic structures (such as spine and joint space), and other open and endoscopic procedures. The Sonicision™ Curved Jaw Cordless Ultrasonic Dissection Device can be used to coagulate isolated vessels up to and including 5 mm in diameter.

The Sonicision 13 cm device is also indicated for use in otorhinolaryngologic (ENT) procedures.

Warning

Safe use of the battery pack requires the user to read, understand, and follow the accompanying instructions.

Warning

Refer to the *Sonicision Curved Jaw Cordless Ultrasonic Dissection System User's Guide* provided with the reusable generator for specific indications, contraindications, warnings, precautions, and the operation instructions.

Do not immerse the battery pack.

Do not disassemble or modify the battery pack.

Inspect the battery pack for damage before use. Do not use a damaged battery pack.

The battery pack is non-sterile and cannot be sterilized. The battery pack will be damaged if sterilization is attempted.

Do not rinse the battery pack under running water or submerge the battery pack in water.

Do not use instrument lubricant on the battery pack.

Precaution

Use of a battery pack that is not fully charged can result in a procedural delay.

The battery pack should be cleaned and disinfected after each use and before charging. Clean and disinfect the battery pack as instructed.

Notice

Battery packs are shipped with a low state of charge and must be charged before use.

The battery pack is designed to be used up to 200 times, after which the battery pack is no longer usable and must be replaced.

Important

Instructions provided in this document are limited to the charging, cleaning, and disinfection of the Sonicision Reusable Battery Pack. Refer to the *Sonicision Curved Jaw Cordless Ultrasonic Dissection System User's Guide* for complete information and important warnings and cautions related to the use of the battery pack.

Cleaning and Disinfecting the Battery Pack

Important

Do not expose the battery pack connector to excessive moisture, conductive materials, or contaminants.

Cleaning the Battery Pack

1. Moisten a lint-free cloth with a disinfectant solution or use a disinfectant wipe.

Use one of the following disinfectants:

- Ethyl alcohol or isopropyl alcohol (70%)
- Hydrogen peroxide wipe (e.g. Oxivir™*)

2. Using a clean disinfectant wipe, wipe all external surfaces of the battery pack until visually clean. Ensure the surface remains thoroughly wet by continuously wiping for the contact time recommended per the manufacturer's instructions.
 - If using 70% ethyl alcohol or isopropyl alcohol, ensure a contact time of 1 minute.
 - Pay special attention to the battery pack handle.
 - Avoid contact with the battery pack connector.
 - Use additional wipes if needed.
3. Use a lint-free cloth to dry all exposed surfaces of the battery pack, taking care to avoid damage to the electrical contacts.
4. Visually inspect the battery pack. If it is not clean, repeat cleaning steps above.


Disinfecting the Battery Pack

1. Moisten a lint-free cloth with a disinfectant solution or use a disinfectant wipe.

Use one of the following disinfectants:

- Ethyl alcohol or isopropyl alcohol (70%)

- Hydrogen peroxide wipe (e.g. Oxivir™*)
2. Using a clean disinfectant wipe, wipe all external surfaces of the battery pack. Ensure the surface remains thoroughly wet by continuously wiping for the contact time recommended per the manufacturer's instructions.
 - If using 70% ethyl alcohol or isopropyl alcohol, ensure a contact time of 1 minute.
 - Pay special attention to the battery pack handle.
 - Avoid contact with the battery pack connector.
 - Use additional wipes if needed.
3. Use a lint-free cloth to dry all exposed surfaces of the battery pack, taking care to avoid damage to the electrical contacts.
2. Place the battery pack into a charging bay and ensure the battery pack is fully seated. The charging cycle will start automatically when the battery pack is inserted correctly into the battery charger.

Charging a Battery Pack

1. Refer to the *Sonicision Battery Charger Instructions for Use* or the *Sonicision Curved Jaw Cordless Ultrasonic Dissection System User's Guide* to set up the battery charger.
3. Note the display for the battery pack status. Refer to the *Sonicision Curved Jaw Cordless Ultrasonic Dissection System User's Guide* or the *Battery Charger Instructions for Use* for charging indication.
4. To remove a battery pack from a battery charger bay, pull the battery pack straight up.

Battery Pack Storage

The service life of the battery pack is monitored by the battery charger. Store the battery pack in the battery charger between uses. If the battery pack is stored off of the charger, the battery pack should be returned to the charger before use to verify functionality and state of charge, and recharge as necessary. Refer to the *Sonicision Curved Jaw Cordless Ultrasonic Dissection System User's Guide* for environmental parameters.

Battery Pack Disposal

Discard or recycle the battery pack as required by local, state, and government regulations.

Covidien is committed to a healthy world and recommends that Sonicision reusable lithium-ion batteries be recycled following decontamination.

Within the U.S. and Canada

Covidien has provided funding to support the recycling of Sonicision reusable lithium-ion batteries. Visit www.call2recycle.org for more information and to receive a prepaid collection kit.

Outside the U.S. and Canada

Covidien recommends that Sonicision reusable lithium-ion batteries be recycled in accordance with your local laws, rules and national regulations.

More Information

Refer to the *Sonicision Curved Jaw Cordless Ultrasonic Dissection System User's Guide* for information regarding the following elements:

- Troubleshooting
- Technical Specifications

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Rx
ONLY

Part No. PT00033131

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien. Other brands are trademarks of a Covidien company, ^{TM*} brands are trademarks of their respective owner.

© 2016 Covidien.

 Manufactured for Covidien llc,
15 Hampshire Street, Mansfield, MA 02048 USA.
 Covidien Ireland Limited,
IDA Business & Technology Park, Tullamore.

Made in USA. Printed in Canada.

www.covidien.com

REV 04/2017

Sonicision™

Reusable Battery Insertion Guide

REF SCBIGA

SCBIGA Compatible with:

SCDA13
SCDA26
SCDA39
SCDA48
SCBA

- en** Instructions for Use
- fr** Mode d'emploi
- pt** Instruções de uso
- da** Brugervejledning
- de** Gebrauchsanweisung
- es** Instrucciones de uso
- fi** Käyttöohjeet
- it** Istruzioni per l'uso
- no** Bruksanvisning
- nl** Gebruiksaanwijzing
- sv** Bruksanvisning
- ru** Инструкции по применению
- zh** 使用说明

Rx
ONLY

Part No. PT00045206

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien. Other brands are trademarks of a Covidien company,™™ brands are trademarks of their respective owner.

© 2016 Covidien.

Manufactured for Covidien LLC,
15 Hampshire Street, Mansfield, MA 02048 USA.
EC REP Covidien Ireland Limited,
IDA Business & Technology Park, Tullamore.

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
Not made with natural rubber latex

Indications for Use

The Sonicision Reusable Battery Insertion Guide is intended to facilitate the aseptic transfer of the non-sterile Sonicision Reusable Battery Pack into the Sonicision Curved Jaw Cordless Ultrasonic Dissection Device.

Refer to the *Sonicision Curved Jaw Cordless Ultrasonic Dissection System User's Guide* for the system indications for use. Assembly, cleaning and sterilization, and additional instructions for use can be found in the *Sonicision Curved Jaw Cordless Ultrasonic Dissection System User's Guide*.

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© 2019 Covidien. All rights reserved.
© 2020 Covidien. All rights reserved.
COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company, ^{TM*} brands are trademarks of their respective owner.
COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien. Other brands are trademarks of a Covidien company, TM* brands are trademarks of their respective owner.