"SPIKE"

IT

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INSTRUCTIONS FOR USE

1. PRODUCT DESCRIPTION

The "SPIKE" endocardial leads have been designed for temporary cardiac electro-stimulation and/or for the registration of endocavitary cardiac potentials. These instructions for use refer to the whole range of SPIKE leads 4, 5 and 6 French, number of poles 2, 3, 4, 6 or 10, various inter-electrode distances, and different distal curves.

Depending on leads body structure, three wide families of devices are identified:

- SPIKE LC S. External body in Pebax: the proximal part is stiffer and the distal part is more flexible. The electrodes are in AISI Stainless Steel. Touch-Proof connection.
- SPIKE ULTRA. External body in Pebax and braiding in AISI Stainless Steel; the proximal part is braided and stiffer and the distal part, free from braiding, is more flexible. Platinum-iridium or AISI stainless steel electrodes are provided. Connections may be Redel® or Touch-Proof.

- SPIKE EASY. External body in Pebax and braiding in AISI Stainless Steel, with peculiar flexibility and handling. Platinum-iridium or AISI stainless steel electrodes are provided. Connections may be Redel® or Touch-Proof.

Further details about technical features and models are available in specific data sheets.

The material used, Pebax®, is an high biocompatible polyamide, which can maintain its chemical-physical characteristics stable inside the cardio-circulatory system over a long time.

A safety wire, made of AISI Stainless Steel, is soldered to the tip electrode and fixed on the distal end in order to avoid accidental detachments.

The connection terminals on the proximal end of the lead allow the connection to devices equipped with specific 2mm protected plugs (DIN 42802-2) or to adapter cables with compatible multipolar circular connectors Redel® type.

The Spike leads should be used in conjunction with external cardiac stimulators (class CF) compliant with standards EN 60601-2-31 and EN 60601-1. Before use check on the instructions for use of the equipments any limitations or requirements that may preclude the use of leads.

2. PACKING CONTENTS

SPIKE leads are packaged individually in double pouch (paper-PE internal / tyvek-PE external) and are identified with a label on which are indicated model, diameter (French), curve type, number of poles, distance between electrodes, lot number, sterilisation expiry date, production date, Ohm resistance and number of units contained in the package. The sale packaging is a carton box, with a label bearing the information given above. Furthermore relative instructions for use are enclosed.

ATTENTION: SINGLE USE DEVICE - DO NOT REUSE - DO NOT RE-STERILIZE

Reusing the product may lead to: possible presence of biological residues that may cause cross-infections; alteration of materials; loss of initial functional features of product.

3. STERILISATION

These devices are gas (ETO) sterilised. Please, do open the package only after a careful check of the integrity of the pouches which could lead to a loss of sterility. Sterility is guaranteed until the expiry date indicated on the package if it has been correctly stored and it is not damaged. If, at the moment of receipt of the package, the integrity of package should appear compromised, the product must be returned to FIAB in its original packing for our internal inspections. FIAB is not responsible in any case for consequences of possible re-sterilisation after product has been used.

4. INDICATIONS AND USERS

Temporary intravenous cardiac pacing. Recording of intravenous cardiac electrical signals. The product is designed to be used for a continuous period of not more than 7 days. The leads may be implanted, connected and put into operation only by experienced and trained medical staff in heart surgery or heart pacing.

5. CONTRAINDICATIONS

Do not use any access to the cardiac chambers where phlebitis or phlebothrombosis are present or suspected.

6. INTRODUCTION WAYS


Following a method which is by now consolidated by a wide experience, SPIKE leads can be introduced in a peripheral vessel by means of introduction systems like Seldinger type or puncture set of adequate dimensions and from here pushed into the cardiac cavities. For the use of the introducer or puncture set, follow the instructions for use supplied by their manufacturer.

7. PRECAUTIONS AND WARNINGS

- External cardiac stimulators and intravenous potentials recording devices must comply with regulations in force.
- Before the beginning of the procedure check the compatibility of the lead's connector with the devices it will be connected to. N.B. Shouldn't the sockets of stimulator or polygraph be compatible with touch proof type pins, FIAB can provide suitable adapters (i.e. FIAB adapters PG922/2TPx). Extension cables with Redel® or Touch-Proof connectors are available.

Further details about technical features and models are available in specific data sheets.

- **N.B. use only with CF type devices.**
- Before connecting review the instructions for use, supplied with the device which will be used
- All electrical main powered devices for patient support must be regularly grounded through equipotential node.
- Work in sterile environment.
- Do not wash the leads with a solvent; if necessary, during the procedure, use a sterile physiologic solution.
- Control packing before use. If wet, broken or damaged, sterility can't be guaranteed.
- Make sure a defibrillator is ready for use
- The shelf-life of the product is the same as the sterilisation expiry date indicated on the label.
- The magnetic fields produced during MRI may cause uncontrolled movements of the medical devices positioned inside the heart chambers
- ***In case of use of adapters, the CONNECTION adapter-lead must be the last one (connect the adapter to the device before); on the contrary, the DISCONNECTION adapter-lead must be the first one. N.B. Free terminals far from the patient must always be protected; dangerous arrhythmias could arise due to induced currents transmitted to the cardiac tissue through the bare pin.***

CAUTION: SPIKE leads are also supplied not pre-shaped, leaving to the user the possibility to give the tip the most convenient curve (type Cournand, Josephson, Damato or personalised to the procedure which will be used). In fact one of the exclusive features of this lead consists in the presence of a thin inner steel wire which permits to maintain a steady chosen shape during the time lead remains inside the cardiac chambers.

During handling of the lead pay attention not to cause sharp angles on the lead body in order not to compromise its manoeuvrability.

It is also suggested to avoid both traction and sharp angling, while modelling the curve, in order not to compromise manoeuvrability, but also to avoid breakage of internal electrical contacts.

7.1. PRE-SHAPING OF THE CURVE

To shape the curve it is suggested to proceed as follows:

- a) Hold the lead firmly between thumb and forefinger of one hand, AT ABOUT 8 CM FROM THE TIP SIDE.
- b) Use the thumb of the other hand as a pivot, while forefinger and middle finger exercise an adequate pressure on the part of the lead on the tip side to start forming the desired curve, AVOIDING SHARP ANGLES.
- c) Move the pivoting-thumb slowly towards the tip WITHOUT EXERCISING TRACTION and taking care to respect the main axis of the lead, repeating the pressure exercised by forefinger and middle finger until desired curve is obtained, starting again, if necessary, the described procedure completely or partially.

8. MODE OF USE

8.1. VENTRICLE STIMULATION

- a) Before starting the procedure, patient must be connected to an electrocardiograph to control possible arising of arrhythmia.
- b) Choose the lead with the diameter (French) and length suitable for the age and body surface of the patient.
- c) Remove the lead from the package and place it in a sterile area. If a straight model is being used, shape it following the instructions contained in the paragraph "pre-shaping the curve", using the curvature normally adopted by the operator, and taking care not to make sharp bends.
- d) Insert the lead into the chosen vein using an introducer of a size suitable for the lead diameter.
- e) Guide the lead by making pushing-pulling or torsion movements. Torque control is best obtained by rotating with fingers the lead body at the output from terminal connection system. Rotation of this point will cause the tip of the lead to rotate in the desired direction.
- f) During these manoeuvres it is recommended to constantly follow the lead tip with X-ray procedure.
- g) Position the tip of the lead at the apex of the right ventricle, making sure that it is stable.
- h) Switch on the external cardiac stimulator
- i) To start stimulation, connect the lead directly to cardiac stimulator: the plug corresponding to the electrode "1" (black) to the "negative" (-) socket of the external pacemaker and the plug "2" to the "positive" (+) socket. In case of use of adapters or extension cable, they must to be previously fixed to the stimulator.
- j) Look for the best stimulation threshold, modifying where necessary the position of the tip in the ventricular chamber. For a continuous stimulation, set output of external pacemaker to a value, which is at least twice as high as the recorded threshold.
- k) In case of prolonged stimulation, secure the lead to the skin under the introduction point with a plaster or with not re-absorbable material suture. Check the threshold at least once a day.
- l) ***N.B. At the end of the procedure, disconnect the lead from the device. In case of use of adapters, the disconnection adapter-lead MUST to be the first one.***

8.2. ATRIAL STIMULATION

- a) Before starting the procedure, patient must be connected to an electrocardiograph to control possible arising of arrhythmia.

- b) Choose the lead with the diameter (French) and length suitable for the age and body surface of the patient.
- c) Remove the lead from the package and place it in a sterile area. If a straight model is being used, shape it following the instructions contained in the paragraph “pre-shaping the curve”. Such manoeuvres enable the curve to be modelled in a “J” configuration. Take care not to make sharp bends.
- d) Insert the lead into the chosen vein using an introducer of a size suitable for the lead diameter.
- e) Guide the lead by making pushing-pulling or torsion movements. Torque control is best obtained by rotating with fingers the lead body at the output of terminal connection system. Rotation of this point will cause the tip of the lead to rotate in the desired direction.
- f) During these manoeuvres it is recommended to constantly follow the lead tip with X-ray procedure.
- g) Position the tip of the lead at the apex of the right auricle, making sure that it is stable.
- h) Switch on the external cardiac stimulator
- i) To start stimulation, connect the lead directly to cardiac stimulator: the plug corresponding to the electrode “1” (black) to the “negative” (-) socket of the external pacemaker and the plug “2” to the “positive” (+) socket. In case of use of adapters or extension cable, they must to be previously fixed to the stimulator.
- j) Look for the best stimulation threshold, modifying where necessary the position of the tip in the atrial chamber. For a continuous stimulation, set output of external pacemaker to a value, which is at least twice as high as the recorded threshold.
- k) In case of prolonged stimulation, secure the lead to the skin at the introduction point with a plaster or with not re-absorbable material suture. Check the threshold at least once a day if the stimulation is made for several days.
- l) **N.B. At the end of the procedure, disconnect the lead from the device. In case of use of adapters, the disconnection adapter-lead MUST to be the first one.**

8.3. ELECTROPHYSIOLOGICAL STUDIES – ENDOCAVITARY POTENTIALS RECORDING

- a) Before starting the procedure, patient must be connected to a polygraph both to control possible arising of arrhythmia and to record cardiac intravenous potentials.
- b) Remove the chosen leads from the packages according to the protocol to be applied and place them in a sterile area.
- c) For introduction technique and positioning in cardiac cavities the study of the potentials will be carried out, see paragraphs VENTRICLE STIMULATION and ATRIAL STIMULATION. If necessary, the curve of each lead will be preformed with desired shape before introduction, referring to the procedure described in paragraph PRE-SHAPING OF THE CURVE.
- d) Connect the pins to the polygraph used to record the local cardiac potentials. Once found the desired signals start recording.
- e) In case of use of adapters or extension cable, they must to be previously fixed to the stimulator.
- f) **N.B. At the end of the procedure, disconnect the lead from the device. In case of use of adapters, the disconnection adapter-lead MUST to be the first one.**

9. POSSIBLE COMPLICATIONS

An inserted lead creates a low resistance direct path to the heart. Therefore it is compulsory that all used devices comply with regulations in force. Particular attention must be paid to assure that all devices near the patient are grounded.

A threatening ventricle arrhythmia mechanically induced by the lead tip may arise (this event is more probable in patients suffering from acute myocardial infarct). It might also be electrically induced by a stimulus, happened in a delicate part of the cycle. Therefore it is necessary to always keep a defibrillator in stand-by mode next to the patient.

The perforation of the atrial or ventricle wall is a very rare but possible event. This causes lack of stimulation and sensing. Rarely in these cases, a ventricle fibrillation could arise.

The risks of using temporary pacing catheters include those risks related to heart catheterization, such as but not limited to:

- Infections, arrhythmias, hematoma, pneumothorax, thromboembolic episodes, death.
- Injury of vessels or cardiac structures, valvular injury.
- Bleeding.
- Cardiac perforation, cardiac tamponade.
- Electrode dislocation with loss of capture.
- Unusual increase of pacing threshold or unusual reduction of sensing threshold.
- Problems in the connection of electrodes with temporary cardiac pacemakers or accidental conductor separation, with impossibility to maintain electrical continuity.
- Peripheral nerve stimulation, injury and pain.

Patient should be monitored and pacing / sensing parameters periodically checked.

10. STORAGE

The contents of the packages are protected against damages caused by normal handling. However the package can be damaged or broken if exposed to strong mechanical stress. The product should be stored in its original packaging in rooms with environment conditions of temperature and relative humidity, specified on the label.

11. WASTE DISPOSAL

Waste coming from hospitals must be disposed according to regulations in force.

12. GUARANTEE – LIABILITY

FIAB SpA guarantees that the product complies with Directive 93/42/CEE. The product has been manufactured according to the procedures of FIAB Quality System certified ISO 13485.

SPIKE leads are used for delicate diagnostic or therapeutic procedures and can fail to function because of various reasons including break of inner electrical connections, noises on the signal caused by electro-magnetic interference, artifacts deriving from malfunctioning of recording systems connected to the lead. Although the attention in the designing process, the selection of suitable materials, the manufacturing and the final checks on the device, the product may therefore be damaged before, during or after introduction in the cardiac chambers due to improper manoeuvres or other manoeuvres connected to the chosen procedures.

No liability may be ascribed to the manufacturer who shall not be held liable for medical costs, direct or indirect damage due to lacking function or failures of leads or accessories, when used differently from the instructions for use. It is recommended to report opportunely any malfunction or defect of the product to FIAB Quality Assurance Service.

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“MYOPACE®” EPIMYOCARDIAL ELECTRODES and “MYOCONNECTOR” PLUG

INSTRUCTIONS FOR USE

1. INTRODUCTION

The epimyocardial “MYOPACE®” electrodes are intended for temporary cardiac electrostimulation, via epimyocardium by connection with external cardiac electrostimulators, are sterile and single use only. MYOPACE® are available in unipolar, bipolar and quadripolar models.

2. DESCRIPTION

The MYOPACE® consists of a multithread conductor in stainless steel with an external insulating sheath.

At one end the conductor is connected with a cardiac needle (excluding hook and loop fixing surface, where not expected) and is exposed in order to create a stimulation surface area typically of 9,4 mm² (stimulating surface is calculated on 15 mm of peeled conductor) in the adult versions.

The stimulating surface is formed in advance:

- in a zig-zag in models “TZ”,
- straight in models “T” and “TIN”,
- in a hook shape in models “TV”
- in a loop shape in models “TL”.

Thanks to their different shape, model “TZ” is particularly suitable to be fixed through the conductor wire in the epimyocardic tissue, and model “T”, “TV” and “TL” are suitable for atrial wall implants where the electrode is preferably placed and fixed with simple holding stitches. In the paediatric version, model “P”, the peeled conductor is straight and, as model “T”, it is placed onto the myocardial tissue and not implanted: the paediatric models “TVP” and “TLP” are analogous to the corresponding adult models “TV” and “TL”. The stimulation surface of model “P” is typically 5,0 mm² corresponding to 8 mm of peeled wire.

The “MYOPACE® TIN and TINP” myocardial electrodes are fixed to the epimyocardial tissue in such a way as to offer both effective anchorage and very low risk operation when extracting the electrode after use. In fact, in this electrode model, the mode of fixing is similar to the function of the barbs of an electro-catheter for final stimulation. The barbs come from the same PE coating out of which three strips – 8 mm (mod. “TIN”) and 5 mm long (mod. “TINP”) – are cut out and lifted in all their thickness. Beneath the root of each strip, the PE coating is abraded as much as needed to form an active electrodic surface similar to “T” and “P” models.

At the other end the conductor is crimped into the eye of a needle (called thoracic needle).

In all MYOPACE® models identified by the suffix “M” in the REF code (eg 50060TINM) you must use the Myoconnector adapter (model 50010) to allow connection to the external stimulator

In the unipolar models, without suffix “M” in the REF code (eg. 50060TIN), the thoracic needle has a circular precut. . This furrow allows for manually cutting off the thoracic needle once exteriorized from the thoracic wall, and introducing the stub (connected to the conductor) into the bushing of a threshold analyzer or an external cardiac stimulator.

In bipolar and quadripolar myocardial leads, without suffix “M” in the REF code (eg. 52060TIN or 54060TIN), the connection to the external stimulator is got by preassembled plugs made with metallic tubes. These plugs are crimped on each wire so to guarantee the mechanical and electrical contact in a

permanent way. From the distal part of the plugs, come the corresponding wires, which then are collected and crimped in one single thoracic needle. In this configuration the thoracic needle is not pre-cut and after being extract, it's removed cutting the wires. So the electrode appears as a bipolar / quadripolar cable equipped with plugs ready for the connection to the external stimulator.

The bushes of the external stimulator must be able to tighten pins with a 2 to 0,5 mm Ø. An insulating hose with a suitable diameter, called Myoprotection (Model 50011), available separately as a sterile accessory, protects when needed the pre-cut needle connection under all operating conditions by accidental contacts.

3. WARNINGS

- Single-use Ethylene-Oxide sterilized product. Do not re-use nor sterilize again.
The re-use of the product may lead to: loss of initial functional features of product; alteration of materials; possible presence of biological residues that may cause cross-infections.
- Before opening the package, check carefully to make sure that it has not been damaged in a way that could compromise the sterility of the contents.
- Sterility is guaranteed up to the expiry date printed on the package. If the latter has been correctly stored, it will not be damaged.
- In case of damages or visible defects, do not utilize the product but return it to FIAB.
- The MYOPACE® electrodes should be used in conjunction with external pacemakers (class CF) compliant with standards EN 60601-1 and EN 60601-2-31. Before use check on the instructions for use of the equipments any limitations or requirements that may preclude the use of leads.
- Control the sensing and pacing parameters daily. Unusual increases in the stimulation threshold may occur that could require resorting to temporary intravenous pacing. Unusual reductions in sensing may also occur. Resorting to temporary intravenous stimulation may be necessary also in this case.
- It is recommended to have a defibrillator on standby, both during the control of the sensing and pacing parameters and also while the electrode is being extracted.
- In order to avoid induced currents be conveyed towards the heart, with the risk of bringing on dangerous ventricular arrhythmias, insulate electrically the pins disconnected from the cardiac stimulator. If applicable, use Myoprotection available as sterile accessory.
- Do not subject to MRI the patient with implanted MYOPACE®.
- The device must be implanted, connected and put into operation by a physician specialized in heart-surgery.
- After implantation and connection to pulse generator, the shaft of the lead shall be fixed in order to prevent possible unexpected traction which, if excessive, could cause conductor breakage or detachment of the connection system.

4. INDICATIONS AND USERS

Single and/or dual chamber epicardial temporary stimulation, for adults or pediatric patients, in connection with temporary pulse generators, during and after cardiac open-chest surgery.

MYOPACE® heart wires may be only implanted, connected and put into operation by experienced and trained medical staff in heart surgery.

5. CONTRAINDICATIONS.

- The product must never be used under any circumstances for permanent cardiac electrostimulation.
- Temporary cardiac electrostimulation is not recommended for Patients with implanted pacemakers.
- Do not use the electrodes for cardioversion.

6. DIRECTIONS FOR USE

After surgery, before closing the chest, place the requested electrodes so to guarantee the correct cardiac activity.

For the single chamber pacing (atrial, ventricular), two electrodes must be positioned along the longitudinal direction, same direction of the cardiac current, in the myocardial tissue of the cardiac chamber which has to be paced (according to the personal experiences of the surgeon). The suggested distance between two electrodes is about 1cm, taking care to find a proper zone where the myocardial tissue is not covered by fat or fibrosis.

These two electrical contacts could be realized with two unipolar electrodes or with one bipolar electrode (or with half quadripolar in the dual-chamber pacing).




In case of dual chamber pacing, it's possible to use the different colours in order to select the cardiac chamber to pace.

For this reason, the quadripolar electrodes, specially designed for dual chamber pacing, typically have two wires of one colour (to identify a chamber) and two of another (to identify the other chamber).

Here below it's describe the procedure which has a general validity, except when it's specified, differentiating the use of each model.

6.1 Preliminary operations

After opening the paper/PE pouch remove the MYOPACE® wrapped on the support from the inner plastic bag and perform the following instructions.





MYOPACE® models without "R" in REF code (e.g. 54220TIN) has to be unreeled starting from heart needle.

1 - grab the support with one hand; with the other hand pull the wire holding it at the beginning of the isolation heart side;

2 - partially unreeel the wire and apply the respective endings to the heart surface, following the instructions provided in the next paragraph (operating mode);

3 - unreeel completely the wire by holding the support in the center so as to allow its rotation, taking care not to practice tractions on the part applied to the heart muscle;

4 - release the thoracic needle and perform the necessary steps to bring it outside the thorax described in the next paragraph (operating mode).

MYOPACE® models with "R" in REF code (e.g. 54220RTIN) has to be unreeled starting from the thoracic needle.

1 - grab the support with one hand; with the other hand pull the thoracic needle and perform the necessary steps to bring it outside the thorax as described in the next paragraph (operating mode);

2 - unreeel the wire by holding the support in the center so as to allow its rotation;

3 - release the heart needles;

4 - apply the respective endings to the heart surface, following the instructions provided in the next paragraph (operating mode).

If the MYOPACE® is wrapped in paper support, take the wire from the heart needle side and unreeel until the thoracic needle is released.

6.2 Operating mode

In zig-zag models ("TZ"), by using the heart needle, introduce the electrode in the epimyocardium, avoiding points where the vascular network is particularly evident, until the uncovered part which

precedes the zig-zag shaped conductor comes out, making sure that also the isolating coating penetrates for roughly 1 mm (fig. 1). Cut the conductor at approx. 1 mm from where it comes out of the tissue. It's common practice to apply an appropriate suture point in the position where the zig-zag comes out of the heart tissue, in order to make more solid the contact of the myocardial electrode. In models "P" and "T" the conductor must merely adhere to the epimyocardium through adequate stitching (fig. 2). The positioning of the electrodic surface must be straight. Avoid making crosses or angles which may hinder electrode extraction. Similarly, in the loop models (models "TL" and "TLP"), the conductor must be adhered to the epimyocardium through a suture placed at the beginning of the loop (fig. 3). The size of the knot will determine the force of traction needed for the removal of the electrode. The ligature must not be too tight, in order that after use it is possible to remove the electrode by gently pulling it. **Attention:** do not perform ligatures over the loop, as it could cause lesions during the removal of the electrode. In the hook models ("TV" and "TVP") as well, the conductors must be adhered to the epimyocardium through a suture placed in proximity to the fold (fig. 4). At the time of removal, gently pulling the wire, the tip comes up, allowing the electrode to slide out of the ligature.

In models "TIN" and "TINP", by using the heart needle, introduce the electrode in the epimyocardium, avoiding points where the vascular network is particularly evident, until the PE tines come out. Gently withdraw the electrode: the tines will then open and offer immediate resistance to traction. Stop the traction and cut the electrode at about 5 mm above the tines base (fig. 5 and 6).

After positioning the electrodes, eliminate the corresponding heart needle by cutting the conductor so as to leave an uncovered 15 mm long conductor in adult models "T" and 8 mm long conductor in paediatric models "P".

Place the "In Line" models in such a way that the proximal ring is on and not in the myocardium.

Use each of the thoracic needles to pierce the thoracic wall in order to exteriorize the body of the corresponding electrode, taking care not to keep it in tension inside the chest. Ensure its stability by making one stitch in not absorbable material at the outlet point of the thoracic wall in order to avoid any accidental traction causing the dislocation of the electrode adhering to the epimyocardial tissue. It is recommended exteriorizing the ventricular electrodes on the left-hand side and the atrial ones on the right-hand side.

WARNING: Always check that the thoracic needle and the exposed conductive part connected to it do not come into contact with active devices or parts until the insertion of the adapter has been completed or the external stimulator is connected.

In case of MYOPACE® "M" version, where Myoconnector is required, immediately after completing the exteriorization, remove the thoracic needle by cutting the exposed wires at about 10 mm in relation to the insulator. Insert the exposed conductor part into the adapter hole and slide the conductor until the insulator has reached about 1 mm from the input hole (Fig.14).

Push up the insulating sheath of the adaptor in order to seal off the conductor (Fig.15).

Cut off the wire leaving 5 mm of conductor exposed out the insulating sleeve adapter (fig.16).

Press the cap on for completing the adapter insulation (Fig. 17). Be sure that naked wire is completely cover by the cap, otherwise cut off any surplus. At this point the pin of the connector may be inserted into the bushing for restoring the contact between the electrode and the stimulator.

Check the correct hold by applying a slight traction on the adapter.

WARNING: For the whole period when the pin is not connected to the temporary stimulator, this should be covered by the cap (Fig. 17).

For the other types of MYOPACE®, in case of unipolar, truncate the thoracic needle at the height of the precut (Fig. 8) and insert it into the bushings on the stimulator.

In the case of bipolar / quadripolar, cut the conductors at the end of the pins to remove the chest needle (Fig. 9) and introduce the pins into the corresponding bushings on the pacemaker.

Seeing that threshold measuring devices and external stimulators are not-normally sterile, they must be kept well away from the operating field via the use of a sterile extension cable (such as cable FIAB series F7818 or F7819).

If there are any doubts about the electrical isolation of the pins connected to the stimulator, use Myoprotection, available as a sterile accessory, as follows.

After the extraction, the contact must be protected with the isolation tubes. According to the different models, the following procedures has to be applied:

Unipolar Electrode

The thoracic needle is introduced in the Myoprotection which is moved along the isolated surface of the electrodes (fig. 7). Then the thoracic needle is broken at the level of the pre-cut (fig. 8).

Bipolar/ Quadripolar Electrode

Cut the conductors at the end of the plugs so to remove the thoracic needle (fig. 9) and then introduce the plugs in the Myoprotection which will be moved along the isolated surface of the electrode.

Any insertion of the Myoprotection must be done before cutting off the plug wire lead (fig. 10).

WARNING: The extension cable must be connected to the threshold measuring device or stimulator prior to the connection of the plugs of the implanted electrodes.

Insert the section of cut-off needle into the bushing of the extension cable (Fig. 11) in observance of the polarity (**N.B.** the negative electrode will be the one located more distal in contact with the myocardium, and the positive electrode the one positioned cranially).

Firmly tighten the bushing and perform a slight traction to ensure the stability of the connection system. If Myoprotection is present, slide it until it stops with the bushing teeth (fig. 12).

WARNING: After having moved the patient from the operating table onto the hospital bed with all the required precautions, check that the entire implantation is in perfect order and that the connections are all in order and stable. Check the threshold values 3 – 4 hours after the implantation and following, once a day. Also check the proper working order of the connections.

WARNING: Heart wires are a direct current path, with low resistance to myocardial. To ensure patient safety, all precautions must be taken to eliminate potentially electrostatic discharges from the operator to the metallic parts of the lead. Whenever the electrodes are disconnected from the pacemaker, it is necessary to protect them to prevent accidental contacts. If using Myoprotection, slide it over the contact tip to about 2 mm (Fig. 13).

7. REMOVING ELECTRODES

Before extracting an electrode from the Patient's chest, check which method has been used for implantation.

In order to prevent the tissutal reaction creating any bridles capable of obstructing the removal of the electrode, the removal must be carried out as soon as there is no longer any need to continue with the electrostimulation, and in any case never more than 20 days from the implantation.

WARNING: Although a simple one, this manoeuvre must be carried out by the heart surgeon who performed the surgery and in any case by a physician specialized in heart-surgery.

When using MYOPACE® unipolar "P", "T", "TZ", "TL" or "TV", remove the stitches which block the conductor on the thoracic wall. Firmly holding the conductor between the thumb and forefinger, exert a slow though constant traction until the entire electrode comes out. Do not pull with force.

When using MYOPACE® unipolar "TIN" or "TINP", after having cut suture points on the chest wall, extract the electrode as follows: exert a slow traction until you feel resistance. Then strongly, though briefly, pull the electrode out. The entire electrode is extracted without offering further resistance.

The retention of electrodes by cutting them flush with the skin is not recommended.

WARNING: In bipolar/quadripolar models extract the MYOPACE® up to the bifurcation of the individual wires, then proceed as described for the unipolar models, one wire at a time: never remove the connections together.

8. POTENTIAL COMPLICATIONS

1. Infections, arrhythmias, pneumothorax, hematoma
2. Bleeding;
3. Displacement of the electrode;
4. Unusual increase in the stimulation threshold or unusual reduction in the sensing;
5. Problems of connection of the electrode with external cardiac stimulators equipped with 2 mm sockets or the accidental detachment of the conductor from the needle, thus blocking electric contact between the electrode and cardiac stimulator;
6. Nerve stimulation;
7. Retained epicardial wire's migration;

9. USE OF THE CONNECTOR MODEL 50010 MYOCONNECTOR

In the event of requiring a connection to fixed 2 mm bushings, or if the plug disconnects from the wire, in order to be able to continue the connection of the conductor to the external pacemaker, it will be necessary to use the adapter Model 50010. It is therefore necessary always keeping a 'Myoconnector' on hand. The steps below must be followed:

Remove the Myoprotection, if present, from the implanted electrode. Carefully scrape the insulation with by the cutting edge of a scalpel in order to expose 1 cm of the end portion of the conductor. Check that all the strands of the conductor are correctly braided and if necessary, try to re-braid them. If there are any problems, repeat this operation on another portion of the cable after cutting off the previously peeled section. Insert the exposed part of the conductor into the specific hole in the adapter and push it in until the insulating part comes within 1 mm from the input hole (Fig. 14).

Push up the insulating sheath of the adaptor in order to seal off the conductor (Fig.15).

Cut off the wire leaving 5 mm of conductor exposed out the insulating sleeve adapter (fig.16).

Press the cap on for completing the adapter insulation (Fig. 17). Be sure that naked wire is completely cover by the cap, otherwise cut off any surplus. At this point the pin of the connector may be inserted into the bushing for restoring the contact between the electrode and the stimulator.

Check the correct hold by applying a slight traction on the adapter.

WARNING: whenever the plug is not connected to temporary pulse generator, it must be protected with its special tube (Fig. 17).

10. STORAGE

The contents of each package is protected against damage associated with its normal use. However, the package may deteriorate or break if exposed to severe mechanical stress. The package should be stored in a place where temperatures are between 0°C and 50°C, with a relative humidity between 10% and 90%.

11. DISPOSAL

Refuse coming from hospitals must be disposed of in accordance with existing regulations.

12. GUARANTEE - NOTICE

FIAB SpA guarantees that the product complies with EU Medical Device regulation in force and was manufactured according to the procedures of the FIAB Quality System, ISO 13485 certified. "MYOPACE®" electrodes are utilized for particularly delicate procedures, and may possibly not carry out their function for various reasons, among which are medical complications and mechanical breakage. Notwithstanding the care involved in designing the components, selecting the material, manufacturing, and pre-sale inspection, the product may easily be damaged before, during, or after insertion, due to incorrect manoeuvres or to manoeuvres connected with the surgical operation. The manufacturer will not be responsible, nor will he reimburse medical expenses, nor direct or indirect damages deriving from non-functioning or anomalies in the electrodes or accessories should they be used differently from what is provided for in the present instructions for use. The Quality Assurance Service should be

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INSTRUCTIONS FOR USE

1. PRODUCT DESCRIPTION

The “SPIKE LC N” endocardial leads have been designed for temporary cardiac electrostimulation. These instructions for use refer to the whole range of SPIKE LC N leads 4, 5 and 6 French, with various inter-electrode distances, and different distal curves.

SPIKE LC N are catheters made of a Pebax tube, with a stiffer proximal part, and a more flexible distal section. Electrodes are in AISI stainless steel. The connection system is a 2 mm pin protected with detachable and integral insulating cap to avoid the possibility to be accidentally disconnected or lost.

Further details about technical features and models are available in the specific data sheet.

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The connection terminals on the proximal end of the lead allow the connection to devices equipped with specific 2mm plugs.

The SPIKE LC N leads should be used in conjunction with external cardiac stimulators (class CF), equipped with tightening bushings, compliant with standards EN 60601-2-31 and EN 60601-1. Before use check on the instructions for use of the equipments any limitations or requirements that may preclude the use of leads.

2. PACKING CONTENTS

SPIKE LC N leads are packaged individually in double pouch (paper-PE internal / tyvek-PE external) and are identified with a label on which are indicated model, diameter (French), curve type, number of poles, distance between electrodes, lot number, sterilisation expiry date, production date, Ohm resistance and number of units contained in the package. The sale packaging is a carton box, with a label bearing the information given above. Furthermore relative instructions for use are enclosed.

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4. INDICATIONS AND USERS

Temporary intravenous cardiac pacing. Recording of intravenous cardiac electrical signals. The product is designed to be used for a continuous period of not more than 7 days. The leads may be

implanted, connected and put into operation only by experienced and trained medical staff in heart surgery or heart pacing

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
6. INTRODUCTION WAYS

Following a method which is by now consolidated by a wide experience, SPIKE LC N leads can be introduced in a peripheral vessel by means of introduction systems like Seldinger type or puncture set of adequate dimensions and from here pushed into the cardiac cavities. For the use of the introducer or puncture set, follow the instructions for use supplied by their manufacturer.

7. PRECAUTIONS AND WARNINGS

- External cardiac stimulators and intravenous potentials recording devices must comply with regulations in force.
- Before the beginning of the procedure check the compatibility of the lead's connector with the devices it will be connected to.
- When the pins are not connected to the equipment, they always must be covered by the protective caps.
- The pin protection cap has to be individually removed only before connecting to the instrument, avoiding to touch with fingers the metal part of the pin, or anyway to put it under tension with accidental contacts.
- The protective cap of each connectors has to be removed individually as each pin has to be connected, to avoid to have uncovered and unconnected both pins at the same time.
- Disconnect one pin at a time reinserting the protection and taking care that once the protective cap is reinserted, no conductive surfaces are left uncovered.
- ***N.B. use only with temporary cardiac pulse generators conformed to EN 60601-2-31 and EN 60601-1 standards (CF type) equipped with tightening bushings.***
- Before connecting review the instructions for use, supplied with the device which will be used
- All electrical main powered devices for patient support must be regularly grounded through equipotential node.
- Work in sterile environment.
- Do not wash the leads with a solvent; if necessary, during the procedure, use a sterile physiologic solution.
- Control packing before use. If wet, broken or damaged, sterility can't be guaranteed.
- Make sure a defibrillator is ready for use
- The shelf-life of the product is the same as the sterilisation expiry date indicated on the label.
- The magnetic fields produced during MRI may cause uncontrolled movements of the medical devices positioned inside the heart chambers
- ***In case of use of adapters, the CONNECTION adapter-lead must be the last one (connect the adapter to the device before); on the contrary, the DISCONNECTION adapter-lead must be the first one. N.B. Free terminals far from the patient must always be protected; dangerous arrhythmias could arise due to induced currents transmitted to the cardiac tissue through the bare pin.***

CAUTION: SPIKE LCN leads are also supplied not pre-shaped, leaving to the user the possibility to give the tip the most convenient curve (type Cournand, Josephson, Damato or personalised to the procedure which will be used). In fact one of the exclusive features of this lead consists in the presence of a thin inner steel wire which permits to maintain a steady chosen shape during the time lead remains inside the cardiac chambers.

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- Hold the lead firmly between thumb and forefinger of one hand, AT ABOUT 8 CM FROM THE TIP SIDE.
- Use the thumb of the other hand as a pivot, while forefinger and middle finger exercise an adequate pressure on the part of the lead on the tip side to start forming the desired curve, AVOIDING SHARP ANGLES.
- Move the pivoting-thumb slowly towards the tip WITHOUT EXERCISING TRACTION and taking care to respect the main axis of the lead, repeating the pressure exercised by forefinger and middle finger until desired curve is obtained, starting again, if necessary, the described procedure completely or partially.

8. MODE OF USE

8.1. VENTRICLE / ATRIAL STIMULATION

- Before starting the procedure, patient must be connected to an electrocardiograph to control possible arising of arrhythmia.
- Choose the lead with the diameter (French) and length suitable for the age and body surface of the patient.
- Remove the lead from the package and place it in a sterile area. If a straight model is being used, shape it following the instructions contained in the paragraph "pre-shaping the curve", using the curvature normally adopted by the operator, "J" configuration too. Take care not to make sharp bends.

- d) Insert the lead into the chosen vein using an introducer of a size suitable for the lead diameter.
- e) Guide the lead by making pushing-pulling or torsion movements. Torque control is best obtained by rotating with fingers the lead body at the output from terminal connection system. Rotation of this point will cause the tip of the lead to rotate in the desired direction.
- f) During these manoeuvres it is recommended to constantly follow the lead tip with X-ray procedure.
- g) In case of ventricular stimulation, position the tip of the lead at the apex of the right ventricle, making sure that it is stable.
- h) In case of atrial stimulation, position the tip of the lead at the apex of the right auricle, making sure that it is stable
- i) Switch on the external cardiac stimulator
- j) Remove individually each protective cap and connect each pin to the pulse generator respecting the correct polarity.
- k) Apply a light pull to both pins in order to verify the correct and reliable tightening and retention of the connection bushings.
- l) To start stimulation, connect the lead directly to cardiac stimulator (switched on): the plug corresponding to the electrode "1" (black) to the "negative" (-) socket of the external pacemaker and the plug "2" to the "positive" (+) socket. In case of use of adapters or extension cable, they must to be previously fixed to the stimulator.
- m) Look for the best stimulation threshold, modifying where necessary the position of the tip in the ventricular or atrial chamber. For a continuous stimulation, set output of external pacemaker to a value, which is at least twice as high as the recorded threshold.
- n) In case of prolonged stimulation, secure the lead to the skin under the introduction point with a plaster or with not re-absorbable material suture. Check the threshold at least once a day and the correct connection between lead and pulse generator.
- o) N.B. At the end of the procedure, disconnect the lead from the device (still on) and reinsert the protective caps on the lead pins.. In case of use of adapters, the disconnection adapter-lead MUST to be the first one.**
- p) Switch off the pulse generators.

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Once positioned in the heart muscle, to ensure the safety of the patient it is necessary to take all precautions to eliminate potential electrostatic discharge from the operator towards the metal parts of the connection pin. Therefore the following operations are recommended:

- When the pin are not connected to the equipment, they always must be covered by the protective caps.
- The pin protection cap has to be individually removed only before connecting to the instrument, avoiding to touch with fingers the metal part of the pin, or anyway to put it under tension with accidental contacts.
- The protective cap of each connectors has to be removed individually as each pin has to be connected, to avoid to have uncovered and unconnected both pins at the same time.
- Disconnect one pin at a time by reinserting the protection and taking care that once the protective cap is reinserted, no conductive surface is left uncovered

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SPIKE LC N leads are used for delicate diagnostic or therapeutic procedures and can fail to function because of various reasons including break of inner electrical connections, noises on the signal caused by electro-magnetic interference, artifacts deriving from malfunctioning of recording systems connected to the lead. Although the attention in the designing process, the selection of suitable materials, the manufacturing and the final checks on the device, the product may therefore be damaged before, during or after introduction in the cardiac chambers due to improper manoeuvres or other manoeuvres connected to the chosen procedures.

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