GeneFinder™ COVID-19 Ag Plus Rapid Test

For in vitro diagnostics use only.
For professional use only.

Please follow the instructions carefully when performing the test.


INTENDED USE

The GeneFinder™ COVID-19 Ag Plus Rapid Test is a lateral flow immunoassay cassette intended for the qualitative detection of antigen from SARS-CoV-2 in nasopharyngeal or nasal swab specimens directly collected, or nasopharyngeal swab collected in transport media, from individuals suspected of COVID-19 by their healthcare provider within seven days of symptom onset.

INTRODUCTION & PRINCIPLE

The GeneFinder™ COVID-19 Ag Plus Rapid Test is based on a lateral flow immunoassay. The cassette contains membranes which are pre-coated with anti-SARS-CoV-2 nucleocapsid protein monoclonal antibodies on the test line. Another anti-SARS-CoV-2 nucleocapsid protein monoclonal antibodies are bound to colloidal gold. When the sample is loaded to the sample inlet, SARS-CoV-2 antibodies, that were conjugated with colloidal gold nano-particles, and SARS-CoV-2 antigen complexes are formed and travel up the strip. The complexes will be captured by coated antibodies on the membrane, and then the line will form a visible red line. The presence of SARS-CoV-2 antigen will be indicated by a visible red test line in T-marked position on the side of the result window. The control line (C) line appears in each result window when sample has flowed through the strip. The control line is used as a procedural control. The control line should always appear if the test procedure is performed properly and the reagents are working as intended.

COMPONENTS

- Materials required, but not provided: timer and biohazard waste container
- Please read the **EXTERNAL QUALITY CONTROL** section for instruction and purposes of using the Positive/Negative control swab.
- This product can be used both for nasopharyngeal and nasal swab specimens.
- Reference No. FBE00572 / FBE00578 / FBE00591 / FBE00579 recommend to use for nasopharyngeal test because it includes sampling swab for nasopharyngeal specimen.

Reference No.	FBE00572	FBE00578	FBE00591	FBE00579
Instructions for use	1	1	1	1
Test cassette	25	10	25	10
Extraction buffer tube	25	10	25	10
Tube filter cap	25	10	25	10
Specimen sampling swab (For nasopharyngeal test)	25	10	25	10
Positive control swab	1	1	0	0
Negative control swab	1	1	0	0
Paper stand	1	1	1	1

STORAGE

- The shelf life is up to 24 months at 2-30°C
- Do not freeze the product or store it above 30°C.

SPECIMEN COLLECTION AND PREPARATION

- Specimens should be collected and handled in a standard manner according to health authority guidelines.
- Specimens should be tested as soon as possible after collection.
- Inadequate specimen collection, handling, and/or transport may yield a false negative result; therefore, training in specimen collection is highly recommended for accurate test results.
- If specimens are stored in a transport medium, the available transport mediums are as follows.


Product	Cat.no	Manufacturer
UTM-RT	330C	COPAN Diagnostics
Clinical Virus Transport Medium (CTM)	UTM-0018	Noble Biosciences, Inc
MD Transport Medium	V1000	IVOLAB Co., Ltd.
ASAN Transport Medium III	AM408-10	ASAN Pharm Co., Ltd

- If it is not used immediately after collecting specimen, it can be stored and used as shown below.


Recommended storage condition		
Temperature	Time	
-18°C	Not available	
2-8°C	4 hours	
15-25°C	1 hour	
-18°C	Freezing and thawing once	
2-8°C	12 hours	
15-25°C	8 hours	

TEST PROCEDURE

(For nasopharyngeal swab)

(For specimens in transport media)

INTERPRETATION OF RESULTS


- A colored band will appear at the C line (control line) of the result window to show that the test is working properly.
- Red bands appear at both C and T line (control and test line). The antigen of the COVID-19 virus is present in the sample.
- The presence of any line no matter how faint the result is considered positive.
- Positive results should be considered in conjunction with the clinical history and other data available.
- A strong positive specimen may cause the control line to appear blurred, but the reading of the results is not affected.

EXTERNAL QUALITY CONTROL (Positive/Negative Control swab)

External positive and negative control swabs are supplied in the kit and should be tested using the same procedure as used for patient specimens. The positive control swab contains recombinant protein as SARS-CoV-2 antigen. The control swabs are specifically formulated and manufactured to ensure performance of the GeneFinder™ COVID-19 Ag Plus Rapid Test, and are used to verify the user's ability to properly perform the test and interpret the results. The positive control swab will produce a positive test result, and has been manufactured to produce a visible test line. The negative control swab will produce a negative test result.

If the correct control results are not obtained, do not perform patient test or report patient results. Contact technical support before testing patient specimens.

Procedure for Positive / Negative control swab

WARNING AND PRECAUTIONS

(General Precautions)

- For professional and in vitro diagnostic use only.
- Do not re-use any contents in the product.
- Do not use if any content package is damaged.
- Follow the instructions for use carefully. If there is not followed, the reliability of the test results cannot be guaranteed.
- Do not use the product beyond the labeled expiry date.
- Do not arbitrarily disassemble or alter the product.
- Perform the test in a clean environment free of contaminants.

(Precautions for handling products)

- Store the product at 2-30°C in its sealed pouch.
- If stored at 2-8°C, keep it at 18-30°C for 30 minutes before opening the product before use.
- Open the pouch immediately before inspection and place the cassette on a clean, flat surface.
- Use within 15 minutes of opening the package to prevent the cassette from being exposed to moisture.
- Avoid any foreign substances contaminating the sample inlet and measuring window.
- If there is evidence of microbial contamination or excessive turbidity in the extraction buffer tube, discard it and contact customer service.
- Do not mix any contents from other products into this product box.
- Do not move the cassette during the test.

(Precautions for handling specimens)

- This product is only used for testing human nasopharyngeal swab specimens.
- Incorrect specimen collection, storage, and transport may yield inaccurate test results.
- Wear personal protective equipment and disposable gloves when working with samples and reagents in order to prevent any infection or contamination.
- After the inspection, clean your hands and disinfect them with disinfectant.
- Decontaminate and dispose of all samples, product, and potentially contaminated materials as if they were infectious waste, in a biohazard waste container.

(Precautions for result interpretation)

- It is not possible to diagnose SARS-CoV-2 infection only with the result of this product, final decisions must be interpreted together with other clinical information and results of RT-PCR by physician.
- Positive test results do not rule out co-infections with other pathogens.
- A false-negative test result may occur if the level of viral antigen in a sample is below the detection limit of the test or if the sample was collected or transported improperly; therefore, a negative test result does not eliminate the possibility of SARS-CoV-2 infection.
- This product determines the presence of SARS-CoV-2 antigens, and the strength of the test line and the concentration of SARS-CoV-2 antigens are irrelevant.
- Monoclonal antibodies may fail to detect or detect with less sensitivity, SARS-CoV-2 viruses that have undergone minor amino acid changes in the target epitope region.
- This product cannot distinguish between SARS-CoV and SARS-CoV-2 infections.
- Compared with the RT-PCR results, this product may have a lower sensitivity in samples taken after 8 days from the symptom onset.
- Samples taken eight days after symptom onset appear are likely to be negative because the amount of antigens decreases depending on the duration of the disease.
- Reading test results earlier than 15 minutes or later than 20 minutes from adding the sample may yield erroneous results. Check the results from 15 to 20 minutes.
- To obtain the most sensitive results, directly test patient specimens without transport media.
- Use of transport media other than those described in the manual may result in incorrect results.
- Samples stored in a nucleic acid-only transport medium (NAPT) cannot be used.
- Clinical performance was evaluated with frozen samples, and test performance may be different with fresh samples.
- Negative test results are not intended to rule out other non-SARS-CoV-2 viral or bacterial infections.
- Positive and negative predictive values are highly dependent on prevalence rates. Positive test results are more likely to represent false-positive results during periods of little/no SARS-CoV-2 activity when disease prevalence is low. False-negative test results are more likely when prevalence of disease caused by SARS-CoV-2 is high.
- The performance of this test has not been evaluated for use in patients without signs and symptoms of respiratory infection and performance may differ in asymptomatic individuals.

PERFORMANCE CHARACTERISTIC

Clinical performance

The performance of the GeneFinder™ COVID-19 Ag Plus Rapid Test was established with 603 nasopharyngeal swabs prospectively collected, and enrolled from individual patients who tested with a RT-PCR method for SARS-CoV-2 infection. All specimens within a pre-specified date range were selected and then sequentially tested in a blinded fashion. The performance of the GeneFinder™ COVID-19 Ag Plus Rapid Test compared to RT-PCR for nasopharyngeal swabs.

	RT-PCR		Total
	Positive	Negative	
GeneFinder™ COVID-19 Ag Plus Rapid Test	145	1	146
	6	451	457
	151	452	603

Sensitivity: 96.03% (95% CI 91.55 - 98.53) Specificity: 99.78% (95% CI: 98.77 - 99.96)

Limit of detection (LoD)

The LoD is 0.75 x 10³ TCID₅₀/mL. The study used "2019-nCoV/Hong Kong/VM200101/2020" strain. The titer of cultured virus was confirmed by TCID₅₀ assay.

Strain	2019-nCoV/Hong Kong/VM200101/2020
Stock titer	3.16 x 10 ⁶ TCID ₅₀ /mL
Conc. (TCID ₅₀ /mL)	3.00 x 10 ⁵ 1.00 x 10 ⁴ 0.75 x 10 ³ 0.38 x 10 ² 0.19 x 10 ¹
Accuracy	100% 100% 100% 91.4% 0%
LoD	0.75 x 10 ³ TCID ₅₀ /mL

Cross-reactivity

Cross-reactivity of the GeneFinder™ COVID-19 Ag Plus Rapid Test was evaluated by testing a panel of high prevalence respiratory pathogens that could potentially cross-react with the GeneFinder™ COVID-19 Ag Plus Rapid Test. Each organism and virus was tested in triplicate. The final concentration of each organism is documented in the following table.

Potential cross-reactant	Concentration	Cross-reactivity
Human coronavirus 229E	1.00 x 10 ⁶ TCID ₅₀ /mL	Negative
Human coronavirus OC43	1.26 x 10 ⁶ TCID ₅₀ /mL	Negative
Human coronavirus NL63	1.00 x 10 ⁶ TCID ₅₀ /mL	Negative
Adenovirus type1	2.57 x 10 ⁶ TCID ₅₀ /mL	Negative
Adenovirus type2	1.15 x 10 ⁶ TCID ₅₀ /mL	Negative
Adenovirus type3	3.80 x 10 ⁶ TCID ₅₀ /mL	Negative
Adenovirus type4	3.39 x 10 ⁶ TCID ₅₀ /mL	Negative
Adenovirus type26	3.80 x 10 ⁶ TCID ₅₀ /mL	Negative
Parainfluenza virus 1	9.12 x 10 ⁶ TCID ₅₀ /mL	Negative
Parainfluenza virus 2	1.00 x 10 ⁶ TCID ₅₀ /mL	Negative
Parainfluenza virus 3	2.57 x 10 ⁶ TCID ₅₀ /mL	Negative
Parainfluenza virus 4a	3.80 x 10 ⁶ TCID ₅₀ /mL	Negative
Parainfluenza virus 4b	1.00 x 10 ⁶ TCID ₅₀ /mL	Negative
Influenza A	4800 ng/mL	Negative
Influenza B	4096 HAunits/mL	Negative
Enterovirus	1.00 x 10 ⁶ TCID ₅₀ /mL	Negative
Respiratory syncytial virus A	1.00 x 10 ⁶ TCID ₅₀ /mL	Negative
Respiratory syncytial virus B	1.00 x 10 ⁶ TCID ₅₀ /mL	Negative
Rhinovirus	1.00 x 10 ⁶ TCID ₅₀ /mL	Negative
MERS-coronavirus	1.00 x 10 ⁶ TCID ₅₀ /mL	Negative
Streptococcus pneumoniae 262	1.00 x 10 ⁶ CFU/mL	Negative
Legionella pneumophila	1.00 x 10 ⁶ CFU/mL	Negative
Candida albicans	1.00 x 10 ⁶ CFU/mL	Negative
Streptococcus pyogenes	1.00 x 10 ⁶ CFU/mL	Negative
Mycoplasma pneumoniae	1.00 x 10 ⁶ CFU/mL	Negative
Haemophilus influenzae type b	1.00 x 10 ⁶ CFU/mL	Negative
Human coronavirus HKU1 (recombinant full-length nucleocapsid protein)	1 µg/mL	Negative

Endogenous/Exogenous interference substances

There was no interference for potential interfering substances listed below. However, abnormally high concentration in sample may cause inaccurately low or high results.

Interference substances	Interference substances	Concentration	Result	
Endogenous	Mucin	0.5% (w/v)	Negative	
	Whole blood (Anticoagulant: K2-EDTA)	5% (v/v)	Negative	
	Whole blood (Anticoagulant: K2-EDTA)	5% (v/v)	Negative	
	Whole blood (Anticoagulant: Li-Heparin)	5% (v/v)	Negative	
	Whole blood (Anticoagulant: Na-Citrate)	5% (v/v)	Negative	
	Whole blood (Anticoagulant: Na-Heparin)	5% (v/v)	Negative	
	Human serum albumin	40mg/mL	Negative	
OTC Nasal (Drop, Gel, Spray)	Phenylephrine hydrochloride	15% (w/v)	Negative	
	Sodium chloride (i.e., NaCl)	5% (w/v)	Negative	
	Oxymetazoline	15% (w/v)	Negative	
	Cromolyn	15% (w/v)	Negative	
	Fluticasone Propionate	5% (w/v)	Negative	
	Budesonide	3.15 mg/mL	Negative	
	Dexamethasone	0.25 mg/mL	Negative	
OTC Homeopathic Nasal Spray	Zincum gluconum (i.e., Zinc)	5% (w/v)	Negative	
	Alkalol	10% (w/v)	Negative	
	Artificially Nasal Ointment	Mupirocin	5 mg/mL	Negative
	Anti-bacterial, systemic	Tobramycin	10 µg/mL	Negative
	Anti-inflammatory medication	Acetylsalicylic acid	5.62 mmol/L	Negative
	Anti-viral drugs	Zanamivir	5 mg/mL	Negative
	Tamiflu (Ezetimibe Phosphate)	2.5 mg/mL	Negative	
Throat Lozenge	Menthol/Berzocaine	1.5mg/mL	Negative	
Sore Throat Phenol Spray	Phenol	15% (w/v)	Negative	

High-dose hook effect

The GeneFinder™ COVID-19 Ag Plus Rapid Test was tested up to 0.75 x 10⁶ TCID₅₀/mL of heat-inactivated SARS-CoV-2 strain, which is a concentration of 100 times the limit of detection (LoD). No high-dose hook effect was observed.

REFERENCE

- Centers for Disease Control and Prevention(2020) <https://www.cdc.gov/coronavirus/2019-ncov/index.html>
- YANG, D.-S. Myeonghaksandanso-ro, Yeondong-myeon, Sejong-si, 30608, South Korea info@company.co.kr
- Nature, 2020 Mar;579(7788): 270-273

SYMBOL

Symbol	Description	Symbol	Description
CE	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical device.	EC REP	Authorized representative in the European community.
ⓘ	Consult instructions for use.	⚠	Do not re-use.
⏰	Use by date.	📄	Caution, consult accompanying documents.
🏠	In vitro diagnostic medical device.	🏭	Manufacturer.
📦	Batch code.	🌡	Temperature limitation.
📋	Catalogue number.	🧪	Contains sufficient for in-vitro tests.

Manufacture

OSANG Healthcare Co., Ltd. 132, Anyangheondong-ro, Dongan-gu, Anyang-si, Gyeonggi-do, Korea (14040) www.osanghc.com
 (EU) Obelis S.A. Bd. General Weis 93, 1030 Brussels, Belgium Tel: +32(2) 732.69.84 Fax: +32(2) 732.60.93 E-Mail: info@obelis.net
 Manufacturer and Authorized representative information for sterile swabs is as below.

Specimen sampling swab

1. Noble Biosciences, Inc. 13-50, Sinbaek-gil, Jeongnam-myeon, Hwaseong-si, Gyeonggi-do, 18521 Korea Tel: +82-31-291-0044
 2. MT Promed Consulting GmbH Altenhofstr. 80, 66386 St. Ingbert, GERMANY Tel: +49-6894 581020
 3. FAIRING, 10-5 Myeonghaksandanso-ro, Yeondong-myeon, Sejong-si, 30608, South Korea info@company.co.kr
 Tel: +82-44-862-9134 E-mail: info@company.co.kr
 4. MT Promed Consulting GmbH Altenhofstrasse 80, 66386 St. Ingbert, Germany www.mt-procons.com
 ISO27262 Rev. 2022-02-08

GeneFinder™ COVID-19 Ag Plus Rapid Test For nasal test

For in vitro diagnostics use only.
For professional use only.

Please follow the instructions carefully when performing the test.


INTENDED USE

The GeneFinder™ COVID-19 Ag Plus Rapid Test is a lateral flow immunoassay intended for the qualitative detection of antigen from SARS-CoV-2 in nasopharyngeal or nasal swab specimens directly collected, or nasopharyngeal swab collected in transport media, from individuals suspected of COVID-19 by their healthcare provider within eight days of symptom onset.

INTRODUCTION & PRINCIPLE

The GeneFinder™ COVID-19 Ag Plus Rapid Test is based on a lateral flow immunoassay. The cassette contains membranes which are pre-coated with anti-SARS-CoV-2 nucleocapsid protein monoclonal antibodies on the test line. Another anti-SARS-CoV-2 nucleocapsid protein monoclonal antibodies are bound to colloidal gold. When the sample is loaded to the sample inlet, SARS-CoV-2 antibodies, that were conjugated with colloidal gold nano-particles, and SARS-CoV-2 antigen complexes are formed and travel up the strip. The complexes will be captured by coated antibodies on the membrane, and then the line will form a visible red line. The presence of SARS-CoV-2 antigen will be indicated by a visible red test line in T-marked position on side of the result window. The control (C) line appears in each result window when sample has flowed through the strip. The control line is used as a procedural control. The control line should always appear if the test procedure is performed properly and the reagents are working as intended.

COMPONENTS

- Materials required, but not provided, timer and biohazard waste container.
- Please read the «EXTERNAL QUALITY CONTROL» section for instruction and purposes of using the Positive/Negative control swab.
- This product can be used both for nasopharyngeal and nasal swab specimens.
- Reference No. FBE00577 / FBE00592 recommend to use for nasal test because it includes sampling swab for nasal specimen.

Reference No.	FBE00577	FBE00592
Instructions for use	1	1
Test cassette	25	25
Extraction buffer tube	25	25
Tube filter cap	25	25
Specimen sampling swab (For nasal test)	25	25
Positive control swab	1	0
Negative control swab	1	0
Paper stand	1	1

STORAGE


- The shelf life is up to 24 months at 2-30°C.
- Do not freeze the product or store it above 30°C.

SPECIMEN COLLECTION AND PREPARATION


- Specimens should be collected and handled in a standard manner according to health authority guidelines.
- Specimens should be tested as soon as possible after collection.
- Inadequate specimen collection, handling, and/or transport may yield a false negative result; therefore, training in specimen collection is highly recommended for accurate test results.
- If it is not used immediately after collecting specimen, it can be stored and used as shown below.

	Recommended storage condition	
	Temperature	Time
Nasal swab	-70°C	Not available
	2 - 8°C	4 hours
	15 - 25°C	1 hour

TEST PROCEDURE

INTERPRETATION OF RESULTS


- A colored band will appear at the C line (Control line), of the result window to show that the test is working properly.
- Red bands appear at both C and T line (Control and Test line). The antigen of the COVID-19 virus is present in the sample.
- The presence of any line no matter how faint the result is considered positive.
- Positive results should be considered in conjunction with the clinical history and other data available.
- A strong positive specimen may cause the control line to appear blurred, but the reading of the results is not affected.

EXTERNAL QUALITY CONTROL (Positive/Negative Control swab)

External positive and negative control swabs are supplied in the kit and should be tested using the same procedure as used for patient specimens. The positive control swab contains recombinant protein as SARS-CoV-2 antigen. The control swabs are specifically formulated and manufactured to ensure performance of the GeneFinder™ COVID-19 Ag Plus Rapid Test, and are used to verify the user's ability to properly perform the test and interpret the results. The positive control swab will produce a positive test result, and has been manufactured to produce a visible test line. The negative control swab will produce a negative test result.

If the correct control results are not obtained, do not perform patient test or report patient results. Contact technical support before testing patient specimens.

Procedure for Positive / Negative control swab

WARNING AND PRECAUTIONS

[General Precautions]

- For professional and in vitro diagnostic use only.
- Do not re-use any contents in the product.
- Do not use if any content package is damaged.
- Follow the instructions for use carefully. If there is not followed, the reliability of the test results cannot be guaranteed.
- Do not use the product beyond the labeled expiry date.
- Do not arbitrarily disassemble or alter the product.
- Perform the test in a clean environment free of contaminants.

[Precautions for handling products]

- Store the product at 2-8°C in its sealed pouch.
- If stored at 2-8°C, keep it at 18-30°C for 30 minutes before opening the product before use.
- Open the pouch immediately before inspection and place the cassette on a clean, flat surface.
- Use within 15 minutes of opening the package to prevent the cassette from being exposed to moisture.
- Avoid any foreign substances contaminating the sample inlet and measuring window.
- If there is evidence of microbial contamination or excessive turbidity in the extraction buffer tube, discard it and contact customer service.
- Do not mix any of the contents.
- Do not move the

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[Precautions for handling specimens]

- This product is only used for testing human nasal swab specimens.
- Incorrect specimen collection, storage, and transport may yield inaccurate test results.
- Wear personal protective equipment and disposable gloves when working with samples and reagents in order to prevent any infection or contamination.
- After the inspection, clean your hands and disinfect them with disinfectant.
- Decontaminate and dispose of all samples, product, and potentially contaminated materials as if they were infectious waste, in a biohazard waste container.
- This product is available for transport medium from which nasopharyngeal swab is extracted.

[Precautions for result interpretation]

- It is not possible to diagnose SARS-CoV-2 infection only with the result of this product, final decisions must be interpreted together with other clinical information and results of RT-PCR by physician.
- Positive test results do not rule out co-infections with other pathogens.
- A false-negative test result may occur if the level of viral antigen in a sample is below the detection limit of the test or if the sample was collected or transported improperly; therefore, a negative test result does not eliminate the possibility of SARS-CoV-2 infection.
- This product determines the presence of SARS-CoV-2 antigens, and the strength of the test line and the concentration of SARS-CoV-2 antigens are irrelevant.
- Monoclonal antibodies may fail to detect or detect with less sensitivity, SARS-CoV-2 viruses that have undergone minor amino acid changes in the target epitope region.
- This product cannot distinguish between SARS-CoV and SARS-CoV-2 infections.
- Reading test results earlier than 15 minutes or later than 20 minutes from adding the sample may yield erroneous results. Check the results from 15 minutes to 20 minutes.
- Samples stored in a transport medium cannot be used.
- Clinical performance was evaluated with frozen samples, and test performance may be different with fresh samples.
- Negative test results are not intended to rule out other non-SARS-CoV-2 viral or bacterial infections.
- Positive and negative predictive values are highly dependent on prevalence rates. Positive test results are more likely to represent false-positive results during periods of little/no SARS-CoV-2 activity when disease prevalence is low. False-negative test results are more likely when prevalence of disease caused by SARS-CoV-2 is high.
- The performance of this test has not been evaluated for use in patients without signs and symptoms of respiratory infection and performance may differ in asymptomatic individuals.

PERFORMANCE CHARACTERISTIC

Clinical performance

The performance of the GeneFinder™ COVID-19 Ag Plus Rapid Test was established with 170 direct nasal swabs prospectively collected and enrolled from individuals who tested with a RT-PCR method for SARS-CoV-2 infection. The performance of the GeneFinder™ COVID-19 Ag Plus Rapid Test compared to RT-PCR for nasopharyngeal swabs.

		RT-PCR		
		Positive	Negative	Total
GeneFinder™ COVID-19 Ag Plus Rapid Test	Positive	65	0	65
	Negative	5	100	105
	Total	70	100	170

Positive Percent Agreement (PPA): 92.86% (95% CI: 84.34% - 96.91%)

Negative Percent Agreement (NPA): 100.00% (95% CI: 96.30% - 100.00%)

Limit of detection (LoD)

The LoD is 0.75 x 10³ TCID₅₀/mL. The study used '2019-nCoV/Hong Kong/VM20001061/2020' strain. The titer of cultured virus was confirmed by TCID₅₀ assay.

Strain	2019-nCoV/Hong Kong/VM20001061/2020				
Stock titer	3.16 x 10 ⁷ TCID ₅₀ /mL				
Conc. (TCID ₅₀ /mL)	3.00 x 10 ⁷	1.50 x 10 ⁷	7.50 x 10 ⁶	3.75 x 10 ⁶	1.88 x 10 ⁶
Accuracy	100%	100%	100%	44.4%	0%
LoD	7.50 x 10 ³ TCID ₅₀ /mL				

Cross-reactivity

Cross-reactivity of the GeneFinder™ COVID-19 Ag Plus Rapid Test was evaluated by testing a panel of high prevalence respiratory pathogens that could potentially cross-react with the GeneFinder™ COVID-19 Ag Plus Rapid Test. Each organism and virus were tested in triplicate. The final concentration of each organism is documented in the following table.

Potential cross-reactant	Concentration	Cross-reactivity
Human coronavirus 229E	1.00 x 10 ⁷ TCID ₅₀ /mL	Negative
Human coronavirus OC43	1.28 x 10 ⁷ TCID ₅₀ /mL	Negative
Human coronavirus NL63	1.00 x 10 ⁷ TCID ₅₀ /mL	Negative
Adenovirus type1	2.87 x 10 ⁷ TCID ₅₀ /mL	Negative
Adenovirus type2	1.15 x 10 ⁷ TCID ₅₀ /mL	Negative
Adenovirus type3	3.80 x 10 ⁷ TCID ₅₀ /mL	Negative
Adenovirus type8	3.39 x 10 ⁷ TCID ₅₀ /mL	Negative
hMPV-27 Type A2	3.80 x 10 ⁷ TCID ₅₀ /mL	Negative
Parainfluenza virus 1	9.12 x 10 ⁷ TCID ₅₀ /mL	Negative
Parainfluenza virus 2	1.00 x 10 ⁷ TCID ₅₀ /mL	Negative
Parainfluenza virus 3	2.57 x 10 ⁷ TCID ₅₀ /mL	Negative
Parainfluenza virus 4a	3.80 x 10 ⁷ TCID ₅₀ /mL	Negative
Parainfluenza virus 4b	1.00 x 10 ⁷ TCID ₅₀ /mL	Negative
Influenza A	4800 ng/mL	Negative
Influenza B	4096 ng/mL	Negative
Enterovirus	1.00 x 10 ⁷ TCID ₅₀ /mL	Negative
Respiratory syncytial virus A	1.00 x 10 ⁷ TCID ₅₀ /mL	Negative
Respiratory syncytial virus B	1.00 x 10 ⁷ TCID ₅₀ /mL	Negative
Rhinovirus	1.00 x 10 ⁷ TCID ₅₀ /mL	Negative
MERS-coronavirus	1.00 x 10 ⁷ TCID ₅₀ /mL	Negative
Streptococcus pneumoniae 242	1.00 x 10 ⁷ CFU/mL	Negative
Legionella pneumophila	1.00 x 10 ⁷ CFU/mL	Negative
Candida albicans	1.00 x 10 ⁷ CFU/mL	Negative
Streptococcus pyogenes	1.00 x 10 ⁷ CFU/mL	Negative
Mycoplasma pneumoniae	1.00 x 10 ⁷ CFU/mL	Negative
Haemophilus influenzae type b	1.00 x 10 ⁷ CFU/mL	Negative

Endogenous/Exogenous interference substances

There was no interference for potential interfering substances listed below. However, abnormally high concentration in sample may cause inaccurately low or high results.

Interference substances	Concentration	Result
Endogenous	Nuclein	1 mg/mL
	Whole blood	5%
	Phenylephrine	7.5%
Nasal Spray or drops	Sodium chloride	5%
	Dexamethasone	10 µg/mL
	Fluticasone Propionate	0.5%
Homeopathic Nasal Spray	Bucladine	3.15 mg/mL
	Dexamethasone	0.25 mg/mL
Antibiotic, Nasal ointment	Mupirocin	0.6 mg/mL
Antibiotic, systemic	Tetracycline	10 µg/mL
Anti-inflammatory medication	Acetylsalicylic acid	0.5 mg/mL
Anti-viral drugs	Zanamivir	5 mg/mL
	Osetamivir Phosphate	1 mg/mL

High-dose hook effect

The GeneFinder™ COVID-19 Ag Plus Rapid Test was tested up to 0.75 x 10⁷ TCID₅₀/mL of heat-inactivated SARS-CoV-2 strain, which is a concentration of 100 times the limit of detection (LoD). No high-dose hook effect was observed.

REFERENCE

- Centers for Disease Control and Prevention (2020) <https://www.cdc.gov/coronavirus/2019-ncov/index.html>
- Peng Zhou, Xing-Lou Yang, Xian-Guang Wang, et al., A pneumonia outbreak associated with a new coronavirus of probable bat origin. Nature. 2020 Mar;579(7798): 270-273

SYMBOL

Symbol	Description	Symbol	Description
CE	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical device.	EC REP	Authorized representative in the European community.
IB	Consult instructions for use.	⊗	Do not re-use.
Use by date.		⚠	Caution, consult accompanying documents.
IVD	In vitro diagnostic medical device.	MA	Manufacturer.
LOT	Batch code.	⌚	Temperature limitation.
REF	Catalogue number.	⚡	Contains sufficient for <= tests.

Manufacturer

OSANG Healthcare Co., Ltd. 132, Anyangcheondong-ro, Dangan-gu, Anyang-si, Gyeonggi-do, Korea (14040) www.osanghc.com
OSANG Healthcare Co., Ltd. 132, Anyangcheondong-ro, Dangan-gu, Anyang-si, Gyeonggi-do, Korea (14040) www.osanghc.com
PENG ZHOU, XING-LOU YANG, XIAN-GUANG WANG, et al., A pneumonia outbreak associated with a new coronavirus of probable bat origin. Nature. 2020 Mar;579(7798): 270-273

Specimen sampling swab

Zhejiang Gongdeng Medical Technology Co., Ltd. No.10 Beiyuan Ave., Huangyan, 318020 Taizhou, Zhejiang, People's Republic of China
Tel: +86-571-84115698 E-mail: info@china-gongdeng.com

OSANG Healthcare Co., Ltd. 132, Anyangcheondong-ro, Dangan-gu, Anyang-si, Gyeonggi-do, Korea (14040) www.osanghc.com
PENG ZHOU, XING-LOU YANG, XIAN-GUANG WANG, et al., A pneumonia outbreak associated with a new coronavirus of probable bat origin. Nature. 2020 Mar;579(7798): 270-273

FA IHC, 10-3 Myeonghaksan-dong, Yeongdeungmyeong, Seoul, 04624, South Korea www.facompany.co.kr

Tel: +82-44-862-9134 E-mail: fa@facompany.co.kr

MT Promed Consulting GmbH Altenhofstrasse 80, 46384 St. Ingbert, Germany www.mt-procons.com

OSANG HEALTHCARE Co., Ltd.

132, Anyangcheondong-ro, Dangan-gu, Anyang-si, Gyeonggi-do 14040

Anyang-si, Gyeonggi-do 14040

CE

CE

CE

CE


CE

Original


COVID-19 Ag Plus Test Procedure

ISC02783 Rev. 2021-09-03


For nasopharyngeal swab


Insert a specimen sampling swab into nasopharynx, and collect the specimens.


1. Stir the specimen in the extraction buffer tube more than 10 times.


2. Squeeze out the swab by squeezing the extraction buffer tube.

3. Place the tube filter cap with tube.


4. Apply 4 drops vertically of the solution into the sample inlet.


5. Read in 15-20 minutes.


For specimens in transport media


If you stored a transport media specimen in the refrigerator, it must be at room temperature for 30 minutes before testing.


1. Take 350µL of the specimen from the transport media with a micropipette.


2. Mix the specimen with the extraction buffer.

3. Place the tube filter cap with tube.


4. Apply 4 drops vertically of the solution into the sample inlet.

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



 Obelis S.A. Bd. General Wahis 53, 1030 Brussels, Belgium

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 Zhejiang Gongdong Medical Technology Co., Ltd. No.10 Beiyuan Ave., Huangyan, 318020 Taizhou, Zhejiang, People's Republic of China

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GeneFinder™ COVID-19 Ag Plus Rapid Test (For nasal test)

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
INTENDED USE
The GeneFinder™ COVID-19 Ag Plus Rapid Test is a lateral flow immunochromatographic assay intended for the qualitative detection of antigen from SARS-CoV-2 in nasopharyngeal or nasal swab specimens directly collected, or nasopharyngeal swab collected in transport media, from individuals suspected of COVID-19 by their healthcare provider within eight days of symptom onset.

INTRODUCTION & PRINCIPLE

The GeneFinder™ COVID-19 Ag Plus Rapid Test is based on a lateral flow immunoassay. The cassette contains membranes which are pre-coated with anti-SARS-CoV-2 nucleocapsid protein monoclonal antibodies on the test line. Another anti-SARS-CoV-2 nucleocapsid protein monoclonal antibodies are bound to colloidal gold. When the sample is loaded to the sample inlet, SARS-CoV-2 antibodies, that were conjugated with colloidal gold nano-particles, and SARS-CoV-2 antigen complexes are formed and travel up the strip. The complexes will be captured by coated antibodies on the membrane, and then the line will form a visible red line. The presence of SARS-CoV-2 antigen will be indicated by a visible red line in T-marked position on side of the result window. The control (C) line appears in each result window when sample has flowed through the strip. The control line is used as a procedural control. The control line should always appear if the test procedure is performed properly and the reagents are working as intended.

COMPONENTS

- Materials required, but not provided: timer and biohazard waste container.
- Please read the «EXTERNAL QUALITY CONTROL» section for instruction and purposes of using the Positive/Negative control swab.
- This product can be used both for nasopharyngeal and nasal swab specimens.
- Reference No. FBE00577 / FBE00592 recommend to use for nasal test because it includes sampling swab for nasal specimen.

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
	Recommended storage condition	
	Temperature	Time
Nasal swab	> 70°C	Not available
	2-8°C	4 hours
	15-25°C	1 hour

TEST PROCEDURE

INTERPRETATION OF RESULTS


- A colored band will appear at the C line (control line), of the result window to show that the test is working properly.
- Red bands appear at both C and T line (Control and test line). The antigen of the COVID-19 virus is present in the sample.
- The presence of any line no matter how faint the result is considered positive.
- Positive results should be considered in conjunction with the clinical history and other data available.
- A strong positive specimen may cause the control line to appear blurred, but the reading of the results is not affected.

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External positive and negative control swabs are supplied in the kit and should be tested using the same procedure as used for patient specimens. The positive control swab contains recombinant protein as SARS-CoV-2 antigen. The control swabs are specifically formulated and manufactured to ensure performance of the GeneFinder™ COVID-19 Ag Plus Rapid Test, and are used to verify the product before testing patient specimens. The positive control swab will produce a positive test result, and has been manufactured to produce a visible test line. The negative control swab will produce a negative test result.

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WARNING AND PRECAUTIONS

- (General Precautions)
- For professional and in vitro diagnostic use only.
- Do not re-use any contents in the product.
- Do not use if any content package is damaged.
- Follow the instructions for use carefully. If there is not followed, the reliability of the test results cannot be guaranteed.
- Do not use the product beyond the labeled expiry date.
- Do not arbitrarily disassemble or alter the product.
- Perform the test in a clean environment free of contaminants.

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- Store the product at 2-30°C in its sealed pouch.
- If stored at 2-8°C, keep it at 18-30°C for 30 minutes before opening the product before use.
- Open the pouch immediately before inspection and place the cassette on a clean, flat surface.
- Use within 15 minutes of opening the package to prevent the cassette from being exposed to moisture.
- Avoid any foreign substance from entering the cassette.
- If there is evidence of mold, do not use the product.
- Do not mix any contents from other products into the product.
- Do not move the cassette.

(Precautions for handling specimens)

- This product is only used for testing human nasal swab specimens.
- Incorporate specimen collection, storage, and transport may yield inaccurate test results.
- Wear personal protective equipment and disposable gloves when working with samples and reagents in order to prevent any infection or contamination.
- After the inspection, clean your hands and disinfect them with disinfectant.
- Decontaminate and dispose of all samples, product, and potentially contaminated materials as if they were infectious waste, in a biohazard waste container.
- This product is available for transport medium from only which nasopharyngeal swab is extracted.

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- It is not possible to diagnose SARS-CoV-2 infection only with the result of this product, final decisions must be interpreted together with other clinical information and results of RT-PCR by physician.
- Positive test results do not rule out co-infections with other pathogens.
- A false-negative test result may occur if the level of viral antigen in a sample is below the detection limit of the test or if the sample was collected or transported improperly. Therefore, a negative test result does not eliminate the possibility of SARS-CoV-2 infection.
- This product determines the presence of SARS-CoV-2 antigens, and the strength of the test line and the concentration of SARS-CoV-2 antigens are irrelevant.
- Monoclonal antibodies may fail to detect or detect with less sensitivity, SARS-CoV-2 viruses that have undergone minor amino acid changes in the target epitope region.
- This product cannot distinguish between SARS-CoV and SARS-CoV-2 infections.
- Reading test results earlier than 15 minutes or later than 20 minutes from adding the sample may yield erroneous results. Check the results from 15 minutes to 20 minutes.
- Samples stored in a transport medium cannot be used.
- Clinical performance was evaluated with frozen samples, and test performance may be different with fresh samples.
- A false-negative test result is not intended to rule out other non-SARS-CoV-2 viral or bacterial infections.
- Positive and negative predictive values are highly dependent on prevalence rates. Positive test results are more likely to represent false-positive results during periods of little/no SARS-CoV-2 activity when disease prevalence is low. False-negative test results are more likely when prevalence of disease caused by SARS-CoV-2 is high.
- The performance of this test has not been evaluated for use in patients without signs and symptoms of respiratory infection and performance may differ in asymptomatic individuals.

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Positive Percent Agreement (PPA): 92.85% (95% CI: 84.34% - 96.91%)
Negative Percent Agreement (NPA): 100.00% (95% CI: 96.30% - 100.00%)

Limit of detection (LoD)

The LoD is 0.75 x 10³ TCID₅₀/mL. The study used '2019-nCoV/Hong Kong/VM2001041/2020' strain. The titer of cultured virus was confirmed by TCID₅₀ assay.

Strain	2019-nCoV/Hong Kong/VM2001041/2020				
Stock titer	3.16 x 10 ⁷ TCID ₅₀ /mL				
Conc. (TCID ₅₀ /mL)	3.00 x 10 ⁷	1.50 x 10 ⁷	7.50 x 10 ⁶	3.75 x 10 ⁶	1.88 x 10 ⁶
Accuracy	100%	100%	100%	44.4%	0%
LoD	7.50 x 10 ³ TCID ₅₀ /mL				

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REFERENCE

- Centers for Disease Control and Prevention (2020) <https://www.cdc.gov/coronavirus/2019-ncov/index.html>
- Peng Zhou, Xing-Lu Yang, Xian-Guang Wang, et al. A pneumonia outbreak associated with a new coronavirus of probable bat origin. Nature. 2020 Mar;579(7778): 270-273

SYMBOL

Symbol	Description	Symbol	Description
CE	This product fulfills the requirements of Directive 98/79/EC on in vitro Diagnostic medical device.	EC REP	Authorized representative in the European community.
IFU	Consult instructions for use.	Do not re-use.	
Use by date.		Caution, consult accompanying documents.	
In vitro diagnostic medical device.		Manufacturer.	
Batch code.		Temperature limitation.	
Catalogue number.		Contains sufficient for <no> tests.	

Manufacturer:
OSANG Healthcare Co., Ltd. 132, Anyangcheondong-ro, Dangan-gu, Anyang-si, Gyeonggi-do, Korea (14040) www.osanghc.com
General Sales: General Sales 53, 1030 Brussels, Belgium Tel: +32 (0) 2 732 59 54 Fax: +32 (0) 2 732 60 03 E-Mail: mail@osanghc.com
Distributor and Authorized representative information for sterile swabs is as below:
Specimen sampling swab
Zhongshan Gongsong Medical Technology Co., Ltd. No. 10 Beiyuan Ave., Huangyan, 318020 Taizhou, Zhejiang, People's Republic of China
E-Mail: info@zhongshang.com
Zhongshan Gongsong Medical Technology Co., Ltd. No. 10 Beiyuan Ave., Huangyan, 318020 Taizhou, Zhejiang, People's Republic of China
E-Mail: info@zhongshang.com
FA INC. 10-5 Myeonghag-dong-ro, Yeongdong-myeon, Sejong-si, 30068, South Korea Tel: +82-44-842-9131 E-Mail: fa@falconcompany.co.kr
MT Promet Consulting GmbH Altenhofstrasse 80, 46386 St. Ingbert, Germany www.mt-procon.com

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OSANG Healthcare Co., Ltd. 132, Anyangcheondong-ro, Dongan-gu, Anyang-si, Gyeonggi-do, Korea (14040) www.osanghc.com
MT Promedt Consulting GmbH Altenhofstrasse 80, 66386 St. Ingbert, Germany