EnSite™ X EP System

Surface Electrode Kit

Instructions for Use

ENSITE-SEK-5-I-01

EN: English	1
AR: العربية	6
BG: Български	11
CS: Čeština	17
DA: Dansk	22
DE: Deutsch	27
EL: Ελληνικά	33
ES: Español	39
ET: Eesti	45
FI: Suomi	50
FR: Français	56
HR: Hrvatski	62
HU: Magyar	68
IT: Italiano	74
KO: 한국의	80
LT: Lietuvių	86
LV: Latviski	92
MK: Македонски	97
NL: Nederlands	103
NO: Norsk	109
PL: Polski	114
PT: Português	120
RO: Română	126
RU: Русский	132
SK: Slovenčina	138
SL: Slovenščina	143
SR: Српски	148
SV: Svenska	154
TR: Türkçe	159
UK: Український	164
VI: Tiếng Việt	170
ZF: 繁體中文	176
ZH: 简体中文	182

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This warranty is given by Abbott Medical. To make inquiries regarding this warranty, use the contact information for either the manufacturer or the Australian Sponsor.

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

The benefits to you given by this warranty are in addition to your other rights and remedies under the Australian Consumer Law.

Abbott Medical warrants to the buyer that the goods, for a period equal to the validated shelf life of the goods (the "Warranty Period"), shall meet the product specifications established by the manufacturer when used in accordance with the manufacturer's instructions for use and shall be free from defects in materials and workmanship. Unless otherwise obligated by law or mandatory consumer guarantees, under this warranty Abbott Medical will, at its option, replace or repair the goods at its factory, if the buyer, at the buyer's expense, returns the goods within the Warranty Period to Abbott Medical to the address below and after Abbott Medical confirms that the goods are defective.

Abbott Medical
Attn: Postmarket Surveillance
One St. Jude Medical Drive
St. Paul, MN 55117 USA
+1 651 756 2000

To the maximum extent permitted by law, if a mandatory term is implied by law or a mandatory consumer guarantee applies to the goods, and the goods are not of a kind ordinarily acquired for personal, domestic or household use or consumption, Abbott Medical's liability for the breach of the term or guarantee is limited to, at Abbott Medical's option, either replacement or repair of the goods or payment of the costs of replacing or repairing the goods.

EXCEPT FOR MANDATORY TERMS IMPLIED BY LAW, MANDATORY CONSUMER GUARANTEES, OR AS EXPRESSLY PROVIDED IN THIS WARRANTY, ABBOTT MEDICAL DISCLAIMS ANY REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED.

See the Terms and Conditions of Sale for further information.

™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third-party trademark, which is property of its respective owner.

Pat. <http://www.abbott.com/patents>

© 2020 Abbott. All Rights Reserved.

EN: English

EnSite™ X EP System Surface Electrode Kit

Instructions for Use

Description

The EnSite™ X EP System Surface Electrode Kit is used to transmit information from the patient to the EnSite™ X EP System for signal conditioning and data inference. The kit contains the following single-use, non-sterile components for use in conjunction with the EnSite™ X EP System:

Kit component	Quantity	Description
EnSite™ X EP System, System Reference Electrode	1	A self-adhesive, non-sterile, single use disposable electrode. This electrode serves as a system electrical reference for impedance and ECG measurements.
EnSite™ X EP System RL ECG Electrode	1	A self-adhesive, non-sterile, single use disposable electrode that serves as the right leg ECG electrode during an EnSite™ X EP System procedure.
EnSite™ X EP System Locating Electrodes	1	The EnSite™ X EP System Locating Electrodes consist of the following: <ul style="list-style-type: none">- EnSite™ X EP System Surface Electrodes (6) – Self-adhesive, non-sterile, single use disposable electrodes that are color-coded with a placement diagram printed on the back of each electrode. The electrodes serve as diagnostic tools that impart low-intensity alternating currents generated by the EnSite™ X EP System to the patient's thorax and sense the resultant electrical potentials which supports the acquisition of impedance signals used for locating the position of electrodes on conventional EP catheters.- Patient Reference Sensor (PRS) Patches (4) - Foam adhesive patches that are intended to hold PRS sensors in a stable position on the patient during an EnSite™ X EP System procedure.

Indications for Use

The EnSite™ X EP Surface Electrode Kit is indicated for use with the EnSite™ X EP System in accordance with the EnSite™ X EP System indications for use.

Clinical Benefit

As accessories these products contribute to the realization of the intended clinical benefit of the primary, compatible EnSite™ X EP Systems for which they support.

Contraindications

There are no known contraindications.

Warnings

WARNING: A Warning indicates a risk of a situation which, if not avoided, could result in death or serious injury to the patient or user.

Precautions

CAUTION: A Caution indicates a risk of a situation which, if not avoided, could result in minor or moderate injury to the patient or user or damage to the equipment.

- This device is intended for one time use only. Any attempt to reuse may compromise the integrity of this system. Adverse effects of reusing components may include, but are not limited to, inaccurate functionality.
- Store in a cool, dry, dark place. Inspect the packaging and surface electrode(s) prior to use.
- Do not attempt to use this product prior to completely reading and understanding these Instructions for Use.
- Do not alter device(s) in any way.
- This device must be operated by, or under the supervision of, an electrophysiologist trained in the operation of this device and supported by other qualified personnel trained in the field of cardiac EP.
- Do not warm surface electrodes.
- Do not use this product if the package seal is broken, the conductive adhesive is dry, the "use by" date has passed, or the gel column is dislodged from the electrode.
- Before applying the surface electrode(s), ensure that the application sites are hair-free, clean, non-irritated, and dry, to prevent receiving low quality signals from body surface electrodes. Dryness is particularly important if the preparation agent is flammable or a potential skin irritant.
- It is recommended to avoid removing the liner until the surface electrode is ready to be placed on the patient. Once the liner has been removed, do not attempt to stick the liner back on the adhesive side of a surface electrode. Doing so may result in the liner becoming stuck to the surface electrode.
- Do not overlap kit contents with RF return electrodes.
- To reduce likelihood of skin damage when removing, slowly pull back surface electrode(s) at a low angle with one hand while holding the skin taut with the other hand.
- Do not attempt to relocate the device(s) after initial application.
- Discontinue use of this product if skin sensitivity is observed.
- When removing the EnSite™ X EP System, System Reference Electrode from the patient's abdomen, remove the electrode from the patient before disconnecting from the EnSite™ X EP System hardware.
- To prevent dryout, surface electrodes and patches should be used immediately after opening the pouch.

Limitations

- Do not place the surface electrode(s) over devices such as pacemakers or ICDs.
- Device(s) should only be connected to the EnSite™ X EP System.

Adverse Events

Skin and subcutaneous tissue disorders including irritation, damage, sensitivity, burns, ecchymosis, inflammation, and allergic skin reaction.

Directions for Use

Skin Preparation

Good skin preparation is essential to all electrophysiology studies. To help minimize noise related to surface electrode placement, perform these steps at application sites before placing surface electrodes and patches on the patient's skin:

1. Trim excess hair.
2. Gently abrade the skin using a gauze pad or like material or gel prep (such as Omni Prep) that contains small particles of gritty material. Make sure all the gel material is removed from the skin.
3. Clean skin surface with soap and water. Avoid cleaning with isopropyl alcohol unless the skin is excessively oily or covered with lotion.
4. Ensure the skin is completely dry and that all alcohol has evaporated.

Placement Guidelines

When placing the EnSite™ X EP System Surface Electrodes on the patient, it is important to place them in the correct position, not only in relationship to each other, but also in relationship to other surface electrodes.

It is recommended to avoid opening the package and removing the liner until the surface electrode is ready to be placed on the patient. Once the liner has been removed, do not attempt to stick the liner back on the adhesive side of the surface electrode. Doing so may result in the liner becoming stuck to the surface electrode. Dispose of liner after removing from surface electrode.


Front (Anterior) and Back (Posterior) Electrodes

The precordial lead V2 ECG electrode placement at the 4th intercostal space and the left sternal border guides the placement of the anterior EnSite™ X EP System Surface Electrode. If both are placed correctly, the anterior EnSite™ X EP System Surface Electrode will be offset slightly to the left of the patient's sternum. It is important that the posterior EnSite™ X EP System Surface Electrode be placed directly opposite of the anterior EnSite™ X EP System Surface Electrode. As a result, it should be offset slightly to the left of the patient's spine. This placement will help maintain appropriate anteroposterior (AP) orientation. Deviation from the recommended location will skew the AP orientation in the map display.

Right and Left Electrodes

The right and left side EnSite™ X EP System Surface Electrodes must also be directly across from each other and orthogonal to the AP surface electrode arrangement. Any other placement may skew side-to-side orientation.

Figure 1. Caudal view of thorax - Correct (left) and incorrect (right) placement

Neck and Left Leg Electrodes


The neck and left leg EnSite™ X EP System Surface Electrodes should be vertically aligned to maintain an orthogonal relationship to the other four EnSite™ X EP System Surface Electrodes. The neck EnSite™ X EP System Surface Electrode should be placed at the midline. Placement of the left leg EnSite™ X EP System Surface Electrode on the inner thigh is recommended to achieve as near of a midline location as possible. Any other placement may skew the vertical orientation.

NOTES:

- None of the EnSite™ X EP System Surface Electrodes should overlap each other or overlap any other surface electrodes.
- Edge-to-edge contact with ECG electrodes should not create interference.
- Avoid direct contact of any EnSite™ X EP System Surface Electrodes or the EnSite™ X EP System, System Reference Electrode with other patches including an RF (radiofrequency) return patch or a defibrillator patch.

Surface Electrode Placement

Figure 2. Surface electrode placement (left) and lead connection to EnSite™ X EP System SurfaceLink Module (right)

Place surface electrodes on patient (left graphic) and connect surface electrode leads to EnSite™ X EP System SurfaceLink Module (right graphic).

1. EnSite™ X EP System, System Reference Electrode
2. EnSite™ X EP System Surface Electrode - Front (white)
3. EnSite™ X EP System Surface Electrode - Back (brown)
4. EnSite™ X EP System Surface Electrode - Right (orange)
5. EnSite™ X EP System Surface Electrode - Left (red)
6. EnSite™ X EP System Surface Electrode - Left Leg (blue)
7. EnSite™ X EP System Surface Electrode - Neck (green)
8. EnSite™ X EP System RL ECG Electrode (placed as right leg ECG electrode during EnSite™ X EP System procedure)
9. ECG electrodes

NOTE: Not shown on left graphic. Place ECG electrodes in standard 12-lead configuration for an EnSite™ X EP System procedure.

NOTE: The surface electrode placement instructions include the following devices that are not provided with this kit.

- ECG electrodes (typically used during an EnSite™ X EP System procedure) – Customer is required to procure nine ECG electrodes compliant to ANSI/AAMI EC12. The EnSite™ X EP System RL ECG Electrode is used as the tenth ECG electrode.
- EnSite™ X EP System ECG Cable (required for EnSite™ X EP System procedure) – model ENSITE-ECGRL-I-01 (radiolucent cable) or ENSITE-ECGR0-I-01 (radiopaque cable)
- EnSite™ X EP System SurfaceLink Module (required for EnSite™ X EP System procedure) – model ENSITE-SRFLINK-I-01

1. Place the EnSite™ X EP System, System Reference Electrode on the patient's abdomen as shown in the figure above.
2. Connect the EnSite™ X EP System, System Reference Electrode lead to the EnSite™ X EP System SurfaceLink Module as shown in the diagram above.

NOTE: The EnSite™ X EP System, System Reference Electrode must be the first patient electrode connected to the EnSite™ X EP System Amplifier at the beginning of a study and the last patient electrode to be disconnected at the end of the study.

3. Place the V2 ECG electrode at the 4th intercostal space, left sternal edge.
4. Place each EnSite™ X EP System Surface Electrode on the patient per the illustration printed on the foam of the EnSite™ X EP System Surface Electrode. See the diagram above for proper placement. The orientation of the EnSite™ X EP System Surface Electrode leads is not important.

NOTES:

- Grasp an area away from the lead wire when removing the liner from EnSite™ X EP System Surface Electrodes.
- It is recommended to place the front EnSite™ X EP System Surface Electrode first.
- Press firmly on the center of the EnSite™ X EP System Surface Electrodes with a circular motion to ensure proper contact.

5. Connect each of the EnSite™ X EP System Surface Electrode leads to the designated ports on the EnSite™ X EP System SurfaceLink Module as shown in the diagram above.

NOTE: The leads are color-coded to the port it must be connected to on the EnSite™ X EP System SurfaceLink Module.

6. Place the EnSite™ X EP System RL ECG Electrode as shown in the diagram above.
7. Connect the EnSite™ X EP System RL ECG Electrode lead to the RL ECG port on the EnSite™ X EP System SurfaceLink Module as shown in the diagram above.
8. Place the remaining ECG electrodes on the patient in a standard 12-lead configuration. Use the EnSite™ X EP System ECG cable to connect the ECG electrodes to ECG ports on the EnSite™ X EP System SurfaceLink Module as shown in the diagram above.

NOTE: Each ECG port on the EnSite™ X EP System SurfaceLink Module is color-coded and labeled with a specific ECG electrode location. Match the color of the ECG cables to the ports and ensure the ECG port label matches the location of the ECG electrode on the patient when connecting each of the ECG leads.

- Refer to the EnSite™ X EP System Amplifier (model ENSITE-AMP-EPM-01) IFU for information on connecting the EnSite™ X EP System SurfaceLink Module to the EnSite™ X EP System Amplifier.

PRS Patch and PRS Sensor Placement

Perform the steps below to place the PRS patches and PRS sensors for an EnSite™ X EP System procedure.

NOTE: The PRS patch placement instructions include the following devices which are not provided with this kit:


- Anterior PRS sensor and cable (model ENSITE-PRSFRT-01)
- Posterior PRS sensors and cables (model ENSITE-PRSBCK-01)

Front PRS Patch and PRS Sensor Placement

One PRS patch is placed on the front of a patient to hold the anterior PRS sensor (PRS-A). PRS-A will be used to detect metal distortion during an EnSite™ X EP System procedure. Placing the PRS patch that will house PRS-A correctly is necessary for proper function.

NOTE: The anterior PRS sensor and cable has one PRS sensor and is specifically intended to be placed on the front of the patient. Do not place on the back of the patient.

- Remove the liner from the bottom of a PRS patch and place on the front of the patient while adhering to the following placement guidelines:
 - Orient the patch so the opening for the PRS cable is 30 degrees from midline and allows the PRS-A cable to be routed over the patient's left clavicle.
 - The PRS patch (and PRS-A) must not be inclined more than 45 degrees relative to the patient table.
 - Place on the patient's chest over the heart.
 - The PRS patch may be placed on top of the front EnSite™ X EP System Surface Electrode if necessary.
 - Removing and repositioning the PRS patches may diminish adhesiveness. Avoid if possible.
- Remove the gold liner from the PRS patch to expose the adhesive.
- Place the black side of the PRS sensor in the PRS patch's adhesive area (the green side of the PRS sensor should be facing up).
- Connect the PRS cable to system hardware as instructed in the EnSite™ X EP System Amplifier IFU.

Back PRS Patch and PRS Sensor Placement

Three PRS patches are placed on the back of a patient to hold the three posterior PRS sensors (PRS-P1, PRS-P2 and PRS-P3). The PRS-Ps are placed in a triangular formation around the heart. The center of this triangle will be the 'zero' of the magnetic coordinate system. Placing the PRS-P sensors correctly is necessary for proper function.

NOTE: The posterior PRS sensors and cables have three sensors that are specifically intended to be placed on the back of the patient. Do not place on the front of the patient.


- Remove the liners from the bottom of the PRS patches and place on the patient's back while adhering to the following placement guidelines:
 - The patient's heart must fall within the boundaries of the triangle formed by the three back PRS patches.
 - It is recommended that PRS patches are placed in an equilateral triangle (all internal angles around 60 degrees).
 - The PRS-P sensors must be placed at least 8 cm from one another. The EnSite™ X EP System software will flag PRS-P sensors as invalid if they are closer than 8 cm from each other.
 - Removing and repositioning the PRS patches may diminish adhesiveness. Avoid if possible.

- Remove the gold liner from each of the PRS patches to expose the adhesive.
- Place one PRS sensor in the adhesive area of each of the three PRS patches.

NOTE: The PRS-P sensors are identical and therefore can be placed in any of the three PRS patches.

NOTE: The PRS-P sensors may be placed with the green or black side up.

- Connect the PRS cable to system hardware as instructed in the EnSite™ X EP System Amplifier IFU.

PRS Sensor and PRS Patch Removal

1. Peel off each PRS patch from the body using the white foam pull tab, leaving the PRS sensor attached.
NOTE: Avoid pulling on the PRS cable to remove the PRS patch from the patient. This can cause damage to the PRS.
2. Release the PRS sensor from the PRS patch by grasping and pulling the white foam pull tab and the colored pull tab in opposite directions.
3. Remove the PRS sensor and dispose of the PRS patch.

NOTE: The PRS sensors are reusable. Store in a safe place after the procedure.

Surface Electrode Removal

Disconnect surface electrode leads from EnSite™ X EP System SurfaceLink Module ports. When removing surface electrodes, slowly pull back at a low angle with one hand while holding the skin taut with the other hand to minimize skin trauma.

NOTE: Soap and water may aid in surface electrode removal.

Technical Specifications

Table 1. Storage and Transport Conditions

Temperature	-25°C to +60°C
Relative Humidity	10% to 90% (non-condensing)
Altitude	0 m to 4267 m

Disposal

- The IFU is recyclable.
- After use, device(s) and packaging should be appropriately classified for disposal, e.g. electrical equipment, non-hazardous waste, etc. and carefully disposed of in compliance with facility procedures and applicable laws and regulations. Discard any unused components after the procedure.

Incident Reporting Notice


If, in the course of use of this device, you have reason to believe that a serious incident occurred, please report it to the manufacturer. For customers in the European Union, report the serious incident to your national authority as well as to the manufacturer.

Limited Warranty

Abbott Medical warrants that its products shall be free from defects in materials and workmanship under normal use. This warranty does not exceed the "Expiration" date stated on any product labeling. The authorized uses and approved methods of use of each of our products is set forth in the related "Instructions for Use" that accompany each product. Abbott Medical disclaims any responsibility and liability for the use of its products in a manner that has not been authorized or approved. Abbott Medical's liability under this warranty is limited to replacing its products. The foregoing warranty excludes and is in lieu of all other warranties whether expressed or implied including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Abbott Medical disclaims any liability for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this product, other than as expressly provided by specific law. Abbott Medical neither assumes nor authorizes any other person to assume for it any other or additional liability for loss, damage, or expense in connection with this product. For more details please review complete Abbott Medical warranty policy available from Abbott Medical or on the back of an Abbott Medical invoice.

Symbols

The symbols below and harmonized symbols may be found on the product or product label. For harmonized symbols, refer to the Universal Symbols Glossary at <https://medical.abbott/manuals>.

Symbol	Description
	Surface Electrode Kit
	System Reference Electrode
	Locating Electrodes
	RL ECG Electrode
	Medical Device
	Unique device identification number
	Importer
	Affixed to this device in accordance with European Council Directive 2012/19/EU. This directive calls for separate collection and disposal of electrical and electronic equipment. Sorting such waste and removing it from other forms of waste lessens the contribution of potentially toxic substances into municipal disposal systems and into the larger ecosystem.
	Conformité Européenne (European Conformity). Affixed in accordance with European Council Directive 93/42/EEC and 2011/65/EU. Hereby, St. Jude Medical declares that this device is in compliance with the essential requirements and other relevant provisions of this directive.

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