Medtronic

FREEZOR® 207F1, 207F3, 207F5

Cardiac CryoAblation Catheter
Cathéter de cryoablation cardiaque
Kardialer Kryoablations-Katheter
Catéter de crioablación cardíaca
Cardiale Cryoablatiekatheter
Catetere per crioablazione cardiaca
Kardialt kryoablasjonskateter
Cardiac CryoAblation-kateter
Καρδιακός καθετήρας κρυοκατάλυσης
Kryokateter
Cateter cardíaco de crioablação
Kardiális krioablációs katéter
Cewnik do krioablacji serca
Srdeční katétr pro kryoablaci
Катетер для внутрисердечной криоабляции
Kardiálny kryoablačný katéter
Kardiyak CryoAblasyon Kateteri

Technical Manual ▪ Manuel technique ▪
Gebrauchsanweisung ▪ Manual técnico ▪
Technische handleiding ▪ Manuale tecnico ▪
Teknisk håndbok ▪ Teknisk manual ▪ Τεχνικό
εγχειρίδιο ▪ Teknisk manual ▪ Manual técnico ▪
Műszaki leírás ▪ Instrukcja techniczna ▪ Technická
příručka ▪ Руководство по эксплуатации ▪
Technická príručka ▪ Teknik El Kitabı

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1 Description

The Freezor Cardiac CryoAblation Catheter is a flexible, steerable catheter used to ablate and map cardiac tissue. It is used together with the CryoConsole and related components. The tip of the Freezor catheter reaches cryoablation temperatures when refrigerant is injected from the CryoConsole to the tip of the catheter. The catheter tip has an integrated type T thermocouple for temperature reading capability. The catheter is introduced into the vasculature by traditional, minimally invasive techniques. The Freezor catheter is available in 3 models, as described in the following table:

Model	Curves available
207F1	47 mm
207F3	53 mm
207F5	58 mm

For details about the CryoConsole and how to use it with the catheter to perform cryoablation procedures, see the *CryoConsole Operator's Manual*.

1.1 Contents of package

The catheter is supplied sterile. The package contains the following items:

- 1 Freezor Cardiac CryoAblation Catheter
- product documentation

2 Indications for use

The Freezor Cardiac CryoAblation Catheter is intended for use in treatment of cardiac arrhythmias.

3 Contraindications

The Freezor Cardiac CryoAblation Catheter is contraindicated in patients with the following conditions:

- active systemic infections
- cryoglobulinemia
- other conditions where the manipulation of the catheter would be unsafe (for example, intracardiac mural thrombus)

4 Warnings and precautions

Cryoadhesion – Do not pull on the catheter, sheath, umbilical cables, or console while the catheter is frozen to the tissue, as this may lead to tissue injury.

Improper connection – Do not connect the cryoablation catheter to a radiofrequency (RF) generator or use it to deliver RF energy. Doing this may cause catheter malfunction or patient harm.

CryoMapping – The impact of CryoMapping with respect to patient outcomes has not fully been characterized.

RF ablation – Before powering up an RF generator or applying RF energy, disconnect the cryoablation catheter from the CryoConsole to avoid an error message and unnecessary catheter replacement.

Pressurized refrigerant – The catheter contains pressurized refrigerant during operation. Release of this gas into the circulatory system due to equipment failure or misuse could result in gas embolism.

Post-ablation period – Closely monitor patients undergoing cardiac ablation procedures during the post-ablation period for clinical adverse events.

Induced arrhythmias – Catheter procedures may mechanically induce arrhythmias.

Fluoroscopy required for catheter placement – The use of fluoroscopy during catheter ablation procedures presents the potential for significant x-ray exposure to both patients and laboratory staff. Extensive exposure can result in acute radiation injury and increased risk for somatic and genetic effects. Only perform catheter ablation after giving adequate attention to the potential radiation exposure associated with the procedure, and taking steps to minimize this exposure. Give careful consideration before using the device in pregnant women.

Catheter positioning around the chordae tendineae – Avoid positioning the catheter around the chordae tendineae, as this increases the likelihood of catheter entrapment within the heart, which may necessitate surgical intervention or repair of injured tissues.

Cryoablation near prosthetic heart valves – Do not pass the catheter through a prosthetic heart valve (mechanical or tissue). The catheter may become trapped in the valve, damaging the valve and causing valvular insufficiency or premature failure of the prosthetic valve.

Transaortic approach – Use adequate fluoroscopic visualization during a transaortic approach to avoid placing the ablation catheter within the coronary vasculature. Catheter placement within the coronary vasculature may cause vascular injury.

Do not resterilize – Do not resterilize this device for purpose of reuse. Resterilization may compromise the structural integrity of the device or create a risk of contamination from the device that could result in patient injury, illness, or death.

For single use only – This device is intended only to be used once for a single patient. Do not reuse, reprocess, or resterilize this device for purpose of reuse. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device that could result in patient injury, illness, or death.

Catheter handling –

- Do not use excessive force to advance or withdraw the catheter, especially if resistance is encountered.
- Use extreme care when manipulating the catheter. Lack of careful attention can result in injury such as perforation or tamponade.
- Do not use the catheter if it is kinked, damaged, or cannot be straightened.
- Straighten the cooling segment before inserting or withdrawing the catheter.
- Do not at any time preshape or bend the catheter shaft or cooling segment. Bending or kinking the catheter shaft may damage internal structures and increase the risk of catheter failure. Prebending of the distal curve can damage the catheter.
- Catheter advancement should be performed under fluoroscopic guidance.

Biohazard disposal – Discard all used catheters and sterile components in accordance with hospital procedures.

Leakage current from connected devices – Use only isolated equipment (IEC 60601-1 Type CF equipment, or equivalent) with the CryoConsole and catheters or patient injury or death may occur.

System compatibility – Use only Medtronic cryoablation catheters, refrigerant tanks, and components with the CryoConsole. The safety and use of other catheters or components has not been tested.

Sterile package inspection – Inspect the sterile packaging and catheter prior to use. If the sterile packaging or catheter is damaged, do not use the catheter. Contact your Medtronic representative.

Catheter integrity – Do not use the catheter if it is kinked or damaged. If the catheter becomes kinked or damaged while in the patient, remove it and use a new catheter. Prior to injecting, the physician should ensure that there is no kink in the catheter.

Required use environment – Cryoablation procedures should be performed only in a fully equipped facility.

Qualified users – This equipment should be used only by or under the supervision of physicians trained in cryoablation procedures.

Fluid incursion – Do not expose the catheter handle or coaxial and electrical connectors to fluids or solvents. If these components get wet, the cryoablation system may not function properly, and connector integrity may be compromised.

Cardioversion/defibrillation during ablation procedure – Disconnect the catheter's electrical connection prior to cardioversion/defibrillation. Failure to do so may trigger system messages indicating a need for catheter exchange.

Anticoagulation therapy – Administer appropriate levels of peri-procedural anticoagulation therapy for patients undergoing left-sided and transseptal cardiac procedure, and for selected patients undergoing right-sided procedures. Administer anticoagulation therapy before and after the procedure according to the hospital standards.

Environmental limits – Perform cryoablation procedures only within the environmental parameters. Operating outside these parameters may prevent the start or completion of a cryoablation procedure.

Embolism risk – Introducing any catheter into the circulatory system entails the risk of air or gas embolism, which can occlude vessels and lead to tissue infarction with serious consequences. Always advance and withdraw components slowly to minimize the vacuum created and therefore minimize the risk of air embolism.

5 Adverse events

Potential adverse events associated with cardiac catheter cryoablation procedures include, but are not limited to, the following conditions:

- Access site complications (hematoma, infection, thrombosis, ecchymosis, AV fistula, bleeding from puncture site, hemorrhage)
- Arrhythmia (atrial fibrillation, atrial flutter, tachycardia)
- Cardiac arrest
- Chest discomfort, pain or pressure
- Coronary artery spasm, dissection, thrombosis
- Death
- Endocarditis
- Heart block, requiring permanent pacemaker
- Hemothorax
- Myocardial infarction
- Perforation of venous cardiac or surrounding tissue
- Pericardial effusion, tamponade
- Pericarditis
- Pleural effusion
- Pneumothorax
- Pseudoaneurysm
- Pulmonary edema
- Pulmonary embolism
- Stroke
- Thrombus
- Transient ischemic attack (TIA)
- Vasovagal reaction

6 Instructions for use

This catheter has two modes, CryoMapping and CryoAblation. During CryoMapping, warmer cryo temperatures and shorter durations are used to transiently halt electrical conduction. This allows verification of the ablation site or determines if the site is not appropriate for ablation. The preset parameters for CryoMapping are -30 °C for a duration of 60 seconds. Because of the complexity of arrhythmia electrophysiology, four (4) different cryomapping techniques may be used to assess catheter tip position.

- Change in the slow pathway ERP
- Change in the maximal AH Interval during pacing
- Termination of the target arrhythmia
- Uninducibility of the target arrhythmia

During cryoablation cardiac tissue is permanently ablated. The catheter has preset refrigerant flow designed to optimize the amount of refrigerant delivered to the tip to remove as much heat as possible. The preset duration is 240 seconds.

There are three ways to use the different modes:

- Start and stop a CryoMapping injection if the desired effect is not seen
- Start a CryoMapping injection and proceed to CryoAblation if the desired effect is seen
- Start CryoAblation without performing CryoMapping

Following therapy delivery, effectiveness can be tested by trying to re-induce the clinical arrhythmia.

Note: Use adequate filtering on the recording system to allow continuous monitoring of the surface electrocardiogram (ECG) during cryo applications.

6.1 Connecting the catheter

For more detailed instructions refer to the CryoConsole Operator's Manual.

1. Connect the connection box to the CryoConsole.
2. Connect the Freezor catheter to a sterile coaxial umbilical cable and a sterile electrical umbilical cable.
3. Connect the coaxial umbilical cable to the CryoConsole and connect the electrical umbilical cable to the connection box.

6.2 CryoMapping and CryoAblation

For more detailed instructions refer to the CryoConsole Operator's Manual.

1. Prior to introducing the Freezor catheter into the patient test the deflection mechanism using the lever on the handle to ensure it is operational.
2. Using aseptic techniques create a vascular access with an 8 Fr (minimum) introducer or sheath and insert the Freezor catheter.
3. Choose the therapy (CryoMapping or CryoAblation). At the CryoConsole screen adjust the treatment time for cryoablation.
4. Under fluoroscopic guidance position the tip of the Freezor at the target site, ensuring good tip contact.

Note: The Freezor catheter tip can be deflected to facilitate positioning by using the thumb knob to vary tip curvature. Pulling the thumb knob back causes the catheter tip to bend; pushing the knob forward causes the tip to straighten.
5. Perform the therapy as appropriate.
6. Wait for the CryoMapping or CryoAblation to finish or press the Stop button on the CryoConsole to terminate the injection before the preset time.
7. Perform additional treatments if required.
8. Remove the catheter from the patient.

7 Specifications

Catheter shaft size	2.5 mm (7 Fr)
Size of recommended introducer	2.7 mm (8 Fr) minimum
Tip length	4 mm
Shaft length	108 cm
Number of electrodes on tip	4
Spacing between electrodes	2 mm, 5 mm, 2 mm
Number of thermocouples	1
Curves available	207F1 - 47 mm 207F3 - 53 mm 207F5 - 58 mm
Environmental parameters	
Recommended transit temperature	-35 °C to 58 °C (-31 °F to 136 °F); up to 85% relative humidity (non-condensing)

Recommended storage temperature	15 °C to 30 °C (59 °F to 86 °F)
Operation	15 °C to 30 °C (59 °F to 86 °F) at altitudes less than 2400 meters (8000 feet) above sea level
Refrigerant flow	
Flow ± 200 sccm	2100 sccm

8 Medtronic limited warranty
For complete warranty information, see the accompanying limited warranty document.

9 Service
Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed on the back cover.

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



Manufacturer
Medtronic CryoCath LP
9000 Autoroute Transcanadienne
Pointe-Claire, Quebec
H9R 5Z8, Canada
+1 514 694 1212
Technical Support +1 877 464 2796

Authorized Representative in the European Community
Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
+31 45 566 8000

Europe/Middle East/Africa
Medtronic International Trading Sàrl
Route du Molliau 31
Case Postale 84
CH-1131 Tolochenaz
Switzerland
+41 21 802 7000

Australia
Medtronic Australasia Pty Ltd
97 Waterloo Road
North Ryde, NSW 2113
Australia

Canada
Medtronic of Canada Ltd
99 Hereford Street
Brampton, Ontario L6Y 0R3
Canada
+1 905 460 3800

Technical manuals
www.medtronic.com/manuals

© Medtronic, Inc. 2019
M994922A001 1A
2019-05-16

* M 9 9 4 9 2 2 A 0 0 1 *

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The following list includes trademarks or registered trademarks of Medtronic in the United States and possibly in other countries. All other trademarks are the property of their respective owners. /