Approved

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en: PLEASE READ CAREFULLY AND KEEP THIS PACKAGE INSERT FOR FUTURE REFERENCE.

BAUSCH + LOMB

ReNu® MPS BAUSCH + LOMB ReNu Sensitive Multi-Purpose Solution Kit (lens case included)

DIRECTIONS: To clean and disinfect, complete these simple steps.
STEP 1: Place 3 drops of ReNu MPS Multi-Purpose Solution on each side of the lens surface and gently rub for 20 seconds.
STEP 2: Thoroughly rinse each side of the lens with ReNu MPS Multi-Purpose Solution.
STEP 3: Place cleaned contact lens in the lens case and fill with ReNu MPS Multi-Purpose Solution. Soak at least 4 hours. Remember to always use fresh solution. Your lenses are now ready to wear. In saline time necessary if any debris remains on contacts, rinse with ReNu MPS Multi-Purpose Solution prior to insertion. Always follow your eye care practitioner's instructions. Based upon your individual tear chemistry and lens-wearing schedule, your eye care practitioner may recommend additional products or procedures.

WEEKLY CLEANING: Protein removal with an enzymatic cleaner may be necessary for some patients. Please consult your eye care practitioner.
Note: Ringer lens care requires removal of protein deposits from your lenses.
STORAGE: You may store your lenses in the unopened lens case until ready to wear, up to a maximum of 30 days.
CONTENTS: A sterile, isotonic solution that contains boric acid, edetate disodium, poloxamine 1%, sodium borate and sodium chloride; preserved with DYMED® (polydimethylsiloxane) 0.0005%.

ACTION:
Cleans, loosens and removes accumulations of film, deposits and debris from soft contact lenses. Destroys harmful micro-organisms on the surface of the lenses. Rinses, stores and rewets lenses before insertion.

INDICATIONS (Uses)
BAUSCH + LOMB ReNu Multi-Purpose Solution is indicated for use in the daily cleaning, rinsing, and storage wear and extended wear soft contact lenses and chemical (not heat) disinfection, as recommended by your eye care practitioner.

CONTRAINDICATIONS (Reasons not to use)
If you are allergic to any ingredient in this product, do not use.
WARNINGS: PROBLEMS WITH CONTACT LENSES AND LENS CARE PRODUCTS COULD RESULT IN CORNEAL INFECTION AND/OR ULCER AND LEAD TO LOSS OF VISION. It is essential that you follow all labeling instructions for proper use of lenses and lens care products, including the lens case.
Eye problems, including contact lens care, can develop rapidly and lead to loss of vision. Daily wear lenses are not indicated for overnight wear and should not be worn while sleeping. Clinical studies have shown the risk of serious adverse reactions is increased when these lenses are worn overnight. Extended wear lenses should be regularly removed for cleaning and disinfection or for disposal and replacement on the schedule prescribed by your eye care practitioner. Clinical studies have shown that there is an increased incidence of serious adverse reactions in extended wear contact lens users as compared to daily wear contact lens users. Studies have also shown that the risk of serious adverse reactions increases the longer extended wear lenses are worn before removal for cleaning and disinfection or for disposal and replacement.

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Studies have also shown that smokers have a higher incidence of adverse reactions. If you experience eye discomfort, excessive tearing, vision changes, redness of the eye, increase in redness or irritation, or any other eye problem, stop using your eye care practitioner. It is recommended that contact lens wearers see their eye care practitioner twice each year, or more frequently if recommended.

You should fill your lens case with fresh solution every time you store your lenses, and use only fresh solution. You should discard your solution immediately after use. Your lenses have been removed from the lens case. You should not bottle, or distill, or use any non-sterile solution. Empty, then clean and rinse your lens case with ReNu MPS Multi-Purpose Solution. Air dry your lens case each time you remove your lenses. In order to permit excess solution to drain, you can flip over your lens case while air drying. Replace your lens case monthly.

Warning: Failure to discard solution from lens case after each use or use of water to care for your lenses may lead to contamination resulting in eye injury and potential loss of vision. See accompanying instructions for additional important safety information.

PRECAUTIONS:
Always follow the product directions for use. Failure to follow product directions may lead to vision loss.

- Always use your eye care practitioner regularly.
- Always wash and dry hands before handling lenses.
- Do not use tap water, bottled water or saliva with lenses or lens case.
- Only use fresh solution to clean and disinfect contact lenses.
- Always discard any remaining solution in your lens case after each disinfection cycle.
- Always replace your lenses and lenses as directed.
- Always replace your lens case monthly.
- Do not use any eye medication in conjunction with this solution unless under medical supervision.
- Replace cap after use.
- Do not use any eye medication in conjunction with this solution unless under medical supervision.
- Keep the bottle tightly closed and heat to not in use.
- Use before expiration date marked on the cap and bottle.
- Discard remaining solution within 90 days after opening.
- Keep out of reach of children.
- Store at room temperature.

IMPORTANT: If irritation or excessive tearing occurs, persists or increases, or if vision is impaired, discontinue use and promptly consult your eye care practitioner.

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REV. 2021-10
154085802W

WebCenter Approval Information

Page: 3 of 4
Project: Asset I54085802W
Document: I54085802W_ReNu MPS_Insert_RU-UA-EN Milan_APPROVED.pdf
Document description:
Version: 9
Version ID: 00002_0000283953
Downloaded: May 30, 2022 at 4:21 PM GMT+03:00
By: Levchenko, Andrey-Regulatory (ANDREY.LEVCHENKO)
Approval state:

Approved

Stage Nr	Approver	Role / Group	Date and Time	Action	Comment
1	Cole, Melissa- Graphics QA (COLEM), Bausch Health Company (Rochester_US)	RA GRAPHICS QC	Dec 9, 2021 at 5:17 PM GMT+03:00	Approved (GQC Review Cycle)	
2	Brambilla, Massimo- Plant Eng_OPs (BRAMBIM), Bausch Health Company (Milan_Italy_EMEA)	PLANT ENGINEER	Dec 9, 2021 at 5:31 PM GMT+03:00	Approved (Stage1)	
2	Ponomareva, Tatiana-Marketing (TATIANA. PONOMAREVA), Bausch Health Company (Moscow_Russia)	MARKETING	Dec 13, 2021 at 12:40 PM GMT+03:00	Approved (Stage1)	
3	Makarova, Elena- Regulatory (ELENA. MAKAROVA), Bausch Health Company (Moscow_Russia)	REGULATORY	Dec 13, 2021 at 2:34 PM GMT+03:00	Approved (Stage2)	
3	Polishchuk, Olga- Regulatory (OLGA. POLISHCHUK), Bausch Health Company (Kyiv_Ukraine_EMEA)	REGULATORY	Dec 13, 2021 at 2:41 PM GMT+03:00	Approved (Stage2)	

WebCenter Approval Information

Page: 4 of 4

Stage Nr	Approver	Role / Group	Date and Time	Action	Comment
3	Brand, Sarah-RA Proofreader (SARAH.BRAND), Bausch Health Company (Bridgewater_US)	RA PROOF READER	Dec 15, 2021 at 10:37 PM GMT+03:00	Approved (Stage2)	
4	Acosta, Rich-RA Project Lead (RICH.ACOSTA), Bausch Health Company (Bridgewater_US)	RA PROJECT LEAD	Dec 16, 2021 at 12:23 AM GMT+03:00	Approved (Stage3)	
5	Chichester, John-Graphics QA (CHICHEJ), Bausch Health Company (Rochester_US)	RA GRAPHICS LEAD ALL	Dec 16, 2021 at 4:06 PM GMT+03:00	Approved (Stage4)	

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SPECIAL INSTRUCTIONS: N/A

PRINT SUPPLIERS: Please refer to Bausch Health's Print Supplier Guidelines

BLACK

WebCenter Approval Information

Page: 3 of 4
Project: Asset I54085800W
Document: I54085800W_ReNu_MultiPlus_Insert_ru-uk-en_MLN_APPROVED.pdf
Document description:
Version: 8
Version ID: 00002_0000285078
Downloaded: May 30, 2022 at 4:27 PM GMT+03:00
By: Levchenko, Andrey-Regulatory (ANDREY.LEVCHENKO)
Approval state:

Approved

Stage Nr	Approver	Role / Group	Date and Time	Action	Comment
1	Cole, Melissa- Graphics QA (COLEM), Bausch Health Company (Rochester_US)	RA GRAPHICS QC	Dec 22, 2021 at 5:30 PM GMT+03:00	Approved (GQC Review Cycle)	
2	Soehner, Christina-RA Proofreader (CHRISTINA. SOEHNER), Bausch Health Company (Bridgewater_US)	RA PROOF READER	Dec 23, 2021 at 12:18 AM GMT+03:00	Approved (Stage2)	
2	Opolonets, Vasylyna- Regulatory (VASYLYNA. OPOLONETS), Bausch Health Company (Kyiv_Ukraine_EMEA)	REGULATORY	Dec 29, 2021 at 1:12 PM GMT+03:00	Approved (Stage2)	
2	Makarova, Elena- Regulatory (ELENA. MAKAROVA), Bausch Health Company (Moscow_Russia)	REGULATORY	Dec 30, 2021 at 5:55 PM GMT+03:00	Approved (Stage2)	
3	Acosta, Rich-RA Project Lead (RICH. ACOSTA), Bausch Health Company (Bridgewater_US)	RA PROJECT LEAD	Dec 30, 2021 at 6:41 PM GMT+03:00	Approved (Stage3)	

WebCenter Approval Information

Page: 4 of 4

Stage Nr	Approver	Role / Group	Date and Time	Action	Comment
4	Chichester, John- Graphics QA (CHICHEJ), Bausch Health Company (Rochester_US)	RA GRAPHICS LEAD ALL	Dec 30, 2021 at 7:50 PM GMT+03:00	Approved (Stage4)	
4	Blankenship, Sharah- Graphics Lead (BLANKES), Bausch Health Company (Rochester_US)	RA GRAPHICS LEAD ALL	Jan 3, 2022 at 5:12 PM GMT+03:00	Approved (Stage4)	

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© 2021 Bausch + Lomb Incorporated or its affiliates
Bausch + Lomb are trademarks of Bausch + Lomb Incorporated or its affiliates.